1 t r r v 1 t r r r r r R C R f r t z 1 t r r f u o f t f r l f u o R

- [illegible]

[illegible][illegible]

- [illegible]

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This IFU is not applicable for Corail AMT stem size 6

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

CORAIL AMT HA (HYDROXYAPATITE) COATED DEVICE DESCRIPTION

Corail AMT HA coated stems are manufactured of titanium alloy.

TOTAL HIP PROSTHESIS

A Total Hip Prosthesis involves the use of an individually packaged metal femoral hip stem, a modular metal or ceramic femoral head, and either a one piece or a metal-backed two piece acetabular component. These two piece acetabular components are comprised of a metal shell which is assembled with either an Ultra High Molecular Weight Polyethylene (UHMWPE), ceramic* or metal** liner. The femoral head and acetabular components are designed to replace the natural articular surface of the hip joint.

PARTIAL HIP PROSTHESIS

Partial hip replacement or hip hemiarthroplasty using a Partial Hip Prosthesis involves:

- either the use of an individually packaged metal femoral hip stem and a modular unipolar metallic femoral head
- or the use of an individually packaged metal femoral hip stem and a modular bipolar head construct.

CAUTION: * Devices used in ceramic-on-ceramic applications have additional full prescribing information that must be reviewed prior to use.

** In the USA and Canada, ceramic heads are not approved for use with metal liners.

INTENDED USE

Corail AMT HA coated stems are intended for use in total and partial hip arthroplasty. The stems are intended only for uncemented use.

Total hip arthroplasty is intended to provide increased patient mobility and reduce pain by replacing the damaged hip joint articulation in patients where there is evidence of sufficient sound bone to seat and support the components.

Partial hip arthroplasty (hip hemiarthroplasty) is intended to provide increased patient mobility and reduce pain by replacing the damaged hip joint articulation in patients where there is evidence of a satisfactory natural acetabulum and sufficient femoral bone to seat and support the femoral stem.

INDICATIONS

HA coated stems of the Corail Hip System are indicated for cementless use only.

TOTAL HIP REPLACEMENT

Total hip replacement is indicated in the following conditions:

1. A severely painful and/or disabled joint from osteoarthritis, traumatic arthritis, rheumatoid arthritis, or congenital hip dysplasia.
2. Avascular necrosis of the femoral head.
3. Acute traumatic fracture of the femoral head or neck.
4. Failed previous hip surgery including joint reconstruction, internal fixation, arthrodesis, hemiarthroplasty, surface replacement arthroplasty, or total hip replacement.
5. Certain cases of ankylosis.

PARTIAL HIP REPLACEMENT (HEMIARTHROPLASTY)

Partial hip replacement is indicated in the following conditions:

1. Acute fracture of the femoral head or neck that cannot be appropriately reduced and treated with internal fixation.

2. Avascular necrosis of the femoral head.
3. Non-union of femoral neck fractures.
4. Certain high subcapital and femoral neck fractures in the elderly.
5. Degenerative arthritis involving only the femoral head in which the acetabulum does not require replacement.
6. Pathology involving only the femoral head/neck and/or proximal femur that can be adequately treated by hip hemiarthroplasty.

CONTRAINDICATIONS

The following conditions are contraindications for total or partial hip replacement (hemiarthroplasty):

1. Active local or systemic infection.
2. Loss of musculature, neuromuscular compromise or vascular deficiency in the affected limb rendering the procedure unjustified.
3. Poor bone quality, such as osteoporosis, where, in the surgeon's opinion, there could be considerable migration of the prosthesis or a significant chance of fracture of the femoral shaft and/or the lack of adequate bone to support the implant(s).
4. Charcot's or Paget's disease.
5. For hip hemiarthroplasty, any pathological condition of the acetabulum, such as distorted acetabuli with irregularities (e.g. from congenital hip dysplasia), protrusio acetabuli (arthrocatadysis), or migrating acetabuli, that would preclude the use of the natural acetabulum as an appropriate articular surface for the partial hip prosthesis.

NOTE: Diabetes, at present, has not been established as a contraindication. However, because of the increased risk for complications such as infection, slow wound healing, etc., the physician should carefully consider the advisability of hip replacement in the severely diabetic patient.

WARNINGS AND PRECAUTIONS

- **HA Coated Implants:** HA Coated implants must not be implanted with cement.
- When changing the head on a femoral stem which is still in place, it is essential to use a metal-metal interface.
- Implants and trial components from different manufacturers or implant systems should never be used together.
- Use only DePuy modular femoral heads with DePuy femoral stems. The taper of the femoral head **MUST** match the taper of the femoral stem. All stems of the Corail Hip System have a 12/14 Articulate Mini Taper and are compatible with all DePuy heads with 12/14 - 5°43' tapers. Please refer to the Surgical Technique Guide for a listing of compatible components. Use only DePuy acetabular components with DePuy heads.
- In the case of Partial Hip Replacement (Hemiarthroplasty), only the following femoral components must be used in association with the Corail stems: Self-Centering Bipolar, Hasting Bipolar and Modular Cathcart Unipolar.

Please refer to the Surgical Technique Guide for a listing of compatible femoral components.

CAUTION: The Hastings Bipolar heads are not cleared for use in the U.S.A.

- The inner diameter of the acetabular component bearing surface must correspond to the femoral head size. Use of a liner with a non-matching hip head size (e.g. 28mm inner diameter liner with a 22mm head) will result in higher stresses, accelerated wear and early failure.
- Use only DePuy acetabular shells with DePuy acetabular liners. Acetabular shell and acetabular liner sizes **MUST** match.
- Hip prosthesis components should never be re-implanted. Even though the implant may appear undamaged. NOTE: the implant may have developed microscopic imperfections which could lead to failure.
- Always use trial prosthesis for trial purposes. Trials must have the same configuration and size as the corresponding components.

- Single use devices manufactured by DePuy shall not undergo or withstand any form of alteration, such as disassembly, cleaning or re-sterilization, after a single patient use. Do not reuse single use devices. Reuse can potentially compromise device performance and patient safety.
- Do not alter or modify implants in any way.
- Patient weight range shall be considered prior to selecting which Corail AMT HA coated stem is appropriate for use. These stems have been labelled accordingly.
- Quality femoral bone stock is critical to ensure fixation of the femoral stem. The Corail AMT HA coated stem should not be used for a revision surgery unless there is sufficient available bone stock. Where there is no sufficient bone stock, it is recommended to use the Corail Revision HA coated stem.
- Consider bone grafting and other reinforcement techniques, to seat and support the femoral stem especially when performing revision surgery.

The following conditions, singularly or concurrently, tend to impose severe loading on the affected extremity, thereby placing the patient at higher risk for failure of the hip replacement:

1. Obesity or excessive patient weight.
2. Manual labor.
3. Active sports participation.
4. High levels of patient activity.
5. Likelihood of falls.
6. Alcohol or drug addiction.
7. Other disabilities, as applicable.

The following conditions singularly or concurrently, tend to adversely affect the fixation of hip replacement implants:

1. Marked osteoporosis or poor bone stock.
2. Metabolic disorders or systemic pharmacological treatments leading to progressive deterioration of solid bone support for the implant (e.g., diabetes mellitus, steroid therapies, immunosuppressive therapies, etc.).
3. History of general or local infections.
4. Severe deformities leading to impaired fixation or improper positioning of the implant.
5. Tumors of the supporting bone structures.
6. Allergic reactions to implant materials (e.g., bone cement, metal, polyethylene).
7. Congenital dysplasia of the hip which may reduce the bone stock available to support the acetabular cup prosthesis in total hip replacement.
8. Tissue reactions to implant corrosion or implant wear debris.
9. Disabilities of other joints (i.e., knees and ankles).

WHEN THE SURGEON DETERMINES THAT HIP REPLACEMENT IS THE BEST MEDICAL OPTION AVAILABLE AND DECIDES TO USE A CORAIL AMT HA COATED PROSTHESIS IN A PATIENT WHO HAS ANY OF THE ABOVE CONDITIONS OR WHO IS SIMPLY YOUNG AND ACTIVE, IT IS IMPERATIVE THAT THE PATIENT BE INSTRUCTED ABOUT THE STRENGTH LIMITATIONS OF THE MATERIALS USED IN THE DEVICE AND FOR IMPLANT FIXATION, AND THE RESULTANT NEED TO SUBSTANTIALLY REDUCE OR ELIMINATE ANY OF THE ABOVE CONDITIONS.

The surgical and postoperative management of the patient must be carried out with due consideration for all existing conditions. Mental attitudes or disorders resulting in a patient's failure to adhere to the surgeon's orders may delay postoperative recovery and/or increase the risk of adverse effects including implant or implant fixation failure.

Excessive physical activity or trauma to the replaced joint may contribute to premature failure of the hip replacement by causing a change in position, fracture, and/or wear of the implants. The patient should be informed that factors such as weight and activity levels may significantly affect wear.

INFORMATION FOR USE

Preoperative

THE SURGEON SHOULD DISCUSS ALL PHYSICAL AND MENTAL LIMITATIONS PARTICULAR TO THE PATIENT AND ALL ASPECTS OF THE SURGERY AND THE PROSTHESES WITH THE PATIENT BEFORE SURGERY. The discussion should include the limitations and possible consequences of joint replacement, and the necessity to follow the surgeon's instructions postoperatively, particularly in regard to patient activity and weight.

Surgeons should not begin the clinical use of any hip prosthesis before they have thoroughly familiarized themselves with its specific implantation technique (preoperative planning and surgical technique). Surgical technique brochures and videos are available from DePuy.

Intraoperative

It is recommended that components at least one size larger and one size smaller than were preoperatively determined be available at surgery to accommodate intraoperative selection of the appropriate size(s).

Protective covers should be left on until the components are ready to be implanted. Do not use any components if they have been dropped, have impacted a hard surface or have received any impact during procedure. Damage to the component may not be visible, but could cause early failure of the prosthesis. Before implanting a femoral head, the male taper on the femoral stem must be free of blood, bone chips or other foreign materials. Foreign material between the head and the femoral stem taper may impede proper seating of the head on the stem. This could affect the performance of the femoral head or the locking mechanism between the femoral head and the femoral stem.

When performing revision surgery, it is imperative that any fibrous membrane encapsulating the failed prosthesis be removed in its entirety. Residual fibrous membranes may interfere with implant fixation. Improper selection, placement, or positioning of implants may result in unusual stress conditions and a subsequent reduction in the functional life of the implant. Note that a femoral stem placed in varus increases the stress on the proximal medial femoral cortex and may lead to loosening of the implant. Increasing the anteversion of an acetabular component in total hip replacement may result in instability and/or dislocation of the joint.

Prior to closure, the surgical site should be thoroughly cleansed of bone chips, ectopic bone, etc. Foreign particles at the metal/plastic or ceramic/plastic interface may cause excessive wear. Range of motion should be thoroughly checked for improper mating, instability, or impingement and corrected as appropriate.

Postoperative

Strict adherence by the patient to the surgeon's instructions and warnings is extremely important. Accepted practices should be followed in postoperative care.

The patient should be released from the hospital with complete written instructions and warnings regarding concerning exercises, care and all restrictions on activity as well as those linked to exposure to magnetic fields.

Partial weight bearing with two crutches and later with one crutch should be continued until muscle function is sufficiently restored so that the operated extremity is no longer overloaded if crutches are discarded; this may take 10 to 12 weeks.

A continuing periodic patient follow-up is recommended. Because of the unknown functional lifetime of the implant, particularly with respect to the maintenance of implant fixation and UHMWPE bearing surfaces, A-P radiographs of the pelvis should be taken at each follow-up and compared with previous radiographs and correlated with the clinical assessment of the patient. If any radiographic changes are observed, such as the occurrence of radiolucencies, bone resorption, or any changes in the position

of an implant, these changes should be closely monitored to determine whether they are static or progressive and the patient treated appropriately.

MRI SAFETY INFORMATION

Corail AMT HA coated stems have not been tested for safety and compatibility in the MR environment. The devices have not been tested for heating or migration in the MR environment. Scanning a patient who has this device may result in patient injury.

The risks associated with other passive implants in an MR environment have been evaluated and are known to include heating, migration, and image artefacts at or near the implant site.

ADVERSE EVENTS AND COMPLICATIONS

The following are generally the most frequently encountered adverse events and complications in hip arthroplasty:

General

1. Change in position of the prosthetic components.
2. Early or late loosening of the prosthetic components.
3. Fatigue fracture of the femoral stem.
4. Wear or fracture of the polyethylene component.
5. Early or late infection.
6. Peripheral neuropathies. Subclinical nerve damage may also occur as a result of surgical trauma.
7. Tissue reactions, osteolysis, and/or implant loosening caused by metallic corrosion, allergic reactions, or the accumulation of polyethylene or metal wear debris or loose cement particles.

Intraoperative

1. Acetabular perforation.
2. Femoral shaft perforation, fissure, or fracture, which may require the use of internal fixation.
3. Trochanteric fracture.
4. Damage to blood vessels (e.g. iliac, obturator and femoral artery).
5. Temporary or permanent nerve damage (e.g. femoral, obturator or isolated peroneal nerve).
6. Subluxation or dislocation of the hip joint due to implant size or configuration selection, positioning of components and/or muscle and fibrous tissue laxity.
7. Lengthening or shortening of the affected extremity.

Early Postoperative

1. Cardiovascular disorders including venous thrombosis, pulmonary embolism and myocardial infarction.
2. Hematoma and/or delayed wound healing.
3. Pneumonia and/or atelectasis.
4. Subluxation or dislocation.

Late Postoperative

1. Trochanteric avulsion from excessive muscular tension, weight-bearing, or inadvertent intraoperative weakening of the trochanter.
2. Aggravation of problems in the ipsilateral or contralateral knee and ankle joints due to leg length discrepancy, femoral medialization and/or muscular deficiencies.
3. Femoral or acetabular fracture due to trauma or excessive loading, particularly in the presence of poor bone stock caused by severe osteoporosis, bone defects from previous surgery, intraoperative reaming procedures, or bone resorption.
4. Bone resorption, which may contribute to the deterioration of fixation and eventual loosening of the implant.
5. Periarticular calcification or ossification, which may lead to a decrease in joint mobility and range of motion.

6. Traumatic arthrosis of the ipsilateral knee secondary to intraoperative positioning, of the extremity during surgery.
7. Subluxation or dislocation.

The incidence and severity of complications in hip replacement are usually greater in revisions than in primary operations. Common problems encountered in revision surgery may include difficulty in placement of the incision, removal of ectopic bone and old bone cement, positioning and fixation of components, and/or obtaining adequate bony support. In general, increased operative time, blood loss, and risk of infection, pulmonary embolus and wound hematoma can be expected with revision procedures.

HOW SUPPLIED

The components are individually packaged and supplied STERILE. All metal components are sterilized using radiation as indicated on the outer package label. Remove from the package using accepted aseptic technique only after the correct size has been determined.

DO NOT RESTERILIZE and DO NOT USE IF THE PACKAGING IS DAMAGED OR BROKEN AND STERILITY MAY BE COMPROMISED.

CAUTION: HA coated stems CANNOT be re-sterilized. Do not use these components if the sterile packaging appears damaged.

Components may not be re-sterilized by the hospital because of the possibility of damaging the articulating and interfacing surfaces of the implant. The package should be opened only after the correct size has been determined, as opened packages may not be returned for credit.

ANY USED OR SOILED COMPONENT HAS TO BE ISOLATED FROM THE ENVIRONMENT. WASTE PROCESSING SHOULD FOLLOW LOCAL PROCEDURE. PREVENTION OF ANY RISK OF CONTAMINATION IS THE RESPONSIBILITY OF THE MEDICAL STAFF.

Initial CE mark: 07/2002 for CORAIL AMT HA Coated stems

Further information is available from your DePuy representative upon request.

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This IFU is applicable only for Corail AMT stem size 6

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TOTAL HIP PROSTHESIS

A Total Hip Prosthesis involves the use of an individually packaged metal femoral hip stem, a modular metal or ceramic femoral head and either a one piece or a metal-backed two piece acetabular component. These two piece acetabular components are comprised of a metal shell which is assembled with either an Ultra High Molecular Weight Polyethylene (UHMWPE), ceramic* or metal** liner. The femoral head and acetabular components are designed to replace the natural articular surface of the hip joint.

CAUTION: * **Devices used in ceramic-on-ceramic applications have additional full prescribing information that must be reviewed prior to use.**

In the USA, ceramic heads are not approved for use with ceramic liners.

**** In the USA and Canada, ceramic heads are not approved for use with metal liners.**

INTENDED USE

Corail AMT HA coated stems are intended for use in total hip arthroplasty. The stems are intended only for uncemented use.

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INDICATIONS

HA-coated stems of the Corail Hip System are indicated for cementless use only.

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1. A severely painful and/or disabled joint from osteoarthritis, traumatic arthritis, rheumatoid arthritis, or congenital hip dysplasia.
2. Avascular necrosis of the femoral head.
3. Acute traumatic fracture of the femoral head or neck.
4. Failed previous hip surgery including joint reconstruction, internal fixation, arthrodesis, surface replacement arthroplasty, or total hip replacement.
5. Certain cases of ankylosis.

CONTRAINDICATIONS

The following conditions are contraindications for total hip replacement:

1. Active local or systemic infection.
2. Loss of musculature, neuromuscular compromise or vascular deficiency in the affected limb rendering the procedure unjustified.
3. Poor bone quality, such as osteoporosis, where, in the surgeon's opinion, there could be considerable migration of the prosthesis or a significant chance of fracture of the femoral shaft and/or the lack of adequate bone to support the implant(s).
4. Charcot's or Paget's disease.

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WARNINGS AND PRECAUTIONS

- **HA Coated Implants:** HA Coated implants must not be implanted with cement.
- When changing the head on a femoral stem which is still in place, it is essential to use a metal-metal interface.
- Implants and trial components from different manufacturers or implant systems should never be used together.
- Use only DePuy modular femoral heads with DePuy femoral stems. The taper of the femoral head **MUST** be matched to the taper of the femoral stem. All stems of the Corail Hip System have a 12/14 Articulate Mini Taper and are compatible with all DePuy heads with 12/14 - 5°43 tapers. Please refer to the Surgical Technique Guide for a listing of compatible components. Use only DePuy acetabular components with DePuy heads.
- The inner diameter of the acetabular component bearing surface must correspond to the femoral head size. Use of a liner with a non-matching hip head size (e.g. 28mm inner diameter liner with a 22mm head) will result in higher stresses, accelerated wear and early failure.
- Use only DePuy acetabular shells with DePuy acetabular liners. Acetabular shell and acetabular liner sizes **MUST** match.

- Hip prosthesis components should never be reimplanted. Even though the implant appears undamaged, the implant may have developed microscopic imperfections which could lead to failure.
- Always use trial prosthesis for trial purposes. Trials must have the same configuration and size as the corresponding components.
- Single use devices manufactured by DePuy shall not undergo or withstand any form of alteration, such as disassembly, cleaning or re-sterilization, after a single patient use. Do not reuse single use devices. Reuse can potentially compromise device performance and patient safety.
- Do not alter or modify implants in any way.
- The Corail AMT Size 6 stem is not indicated for hemi arthroplasty surgery.
- The patient's weight shall be considered prior to selecting which Corail AMT HA coated stem is appropriate for use. Stems with a body weight limit have been labelled accordingly.
- For some specific stems, warning labels with body weight limits are placed on the package in order to limit use to specific maximum patient weights. For the Corail AMT Size 6 stem, the maximum body weight is limited to 60kg or 130lbs.
- For some stems, head offset 'Warning' notices are visible on labels to limit the maximum offset used for the head. For the Corail AMT Size 6 stem, the maximum offset used for the head is limited to 13 mm.
- Quality femoral bone stock is critical to ensure fixation of the femoral stem. The Corail AMT stem should not be used for a revision surgery unless there is sufficient available bone stock. Where there is no sufficient bone stock, it is recommended to use the Corail Revision HA coated stem.
- Consider bone grafting and other reinforcement techniques to seat and support the femoral stem especially when performing revision surgery.

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WHEN THE SURGEON DETERMINES THAT HIP REPLACEMENT IS THE BEST MEDICAL OPTION AVAILABLE AND DECIDES TO USE THIS PROSTHESIS IN A PATIENT WHO HAS ANY OF THE ABOVE CONDITIONS OR WHO IS SIMPLY YOUNG AND ACTIVE, IT IS IMPERATIVE THAT THE PATIENT BE INSTRUCTED ABOUT THE STRENGTH LIMITATIONS OF THE MATERIALS USED IN THE DEVICE AND FOR IMPLANT FIXATION, AND THE RESULTANT NEED TO SUBSTANTIALLY REDUCE OR ELIMINATE ANY OF THE ABOVE CONDITIONS.

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INFORMATION FOR USE

Preoperative

THE SURGEON SHOULD DISCUSS ALL PHYSICAL AND MENTAL LIMITATIONS PARTICULAR TO THE PATIENT AND ALL ASPECTS OF THE SURGERY

AND THE PROSTHESES WITH THE PATIENT BEFORE SURGERY. The discussion should include the limitations and possible consequences of joint replacement, and the necessity to follow the surgeon's instructions postoperatively, particularly in regard to patient activity and weight.

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Protective covers should be left on until the components are ready to be implanted. Do not use any components if they have been dropped, have impacted a hard surface or have received any impact during procedure. Damage to the component may not be visible, but could cause early failure of the prosthesis. Before implanting a femoral head, the male taper on the femoral stem must be free of blood, bone chips or other foreign materials. Foreign material between the head and the femoral stem taper may impede proper seating of the head on the stem. This could affect the performance of the femoral head or the locking mechanism between the femoral head and the femoral stem. When performing revision surgery, it is imperative that any fibrous membrane encapsulating the failed prosthesis be removed in its entirety. Residual fibrous membranes may interfere with implant fixation. Improper selection, placement, or positioning of implants may result in unusual stress conditions and a subsequent reduction in the functional life of the implant. Note that a femoral stem placed in varus increases the stress on the proximal medial femoral cortex and may lead to loosening of the implant. Increasing the anteversion of an acetabular component in total hip replacement may result in instability and/or dislocation of the joint.

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A continuing periodic patient follow-up is recommended. Because of the unknown functional lifetime of the implant, particularly with respect to the maintenance of implant fixation and UHMWPE bearing surfaces, A-P radiographs of the pelvis should be taken at each follow-up and compared with previous radiographs and correlated with the clinical assessment of the patient. If any radiographic changes are observed, such as the occurrence of radiolucencies, bone resorption, or any changes in the position of an implant, these changes should be closely monitored to determine whether they are static or progressive and the patient treated appropriately.

MRI SAFETY INFORMATION

The Corail AMT® Size 6 – HA Coated Stems have not been tested for safety and compatibility in the MR environment. The Corail AMT® Size 6 – HA Coated Stems have not been tested for heating or migration in the MR environment. Scanning a patient who has this device may result in patient injury.

The risks associated with other passive implants in an MR environment have been evaluated and are known to include heating, migration, and image artifacts at or near the implant site.

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Intraoperative

1. Acetabular perforation.
2. Femoral shaft perforation, fissure, or fracture, which may require the use of internal fixation.

3. Trochanteric fracture.
4. Damage to blood vessels (e.g. iliac, obturator and femoral artery).
5. Temporary or permanent nerve damage (e.g. femoral, obturator or isolated peroneal nerve).
6. Subluxation or dislocation of the hip joint due to implant size or configuration selection, positioning of components and/or muscle and fibrous tissue laxity.
7. Lengthening or shortening of the affected extremity.

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Late Postoperative

1. Trochanteric avulsion from excessive muscular tension, weight-bearing, or inadvertent intraoperative weakening of the trochanter.
2. Aggravation of problems in the ipsilateral or contralateral knee and ankle joints due to leg length discrepancy, femoral medialization and/or muscular deficiencies.
3. Femoral or acetabular fracture due to trauma or excessive loading, particularly in the presence of poor bone stock caused by severe osteoporosis, bone defects from previous surgery, intraoperative reaming procedures, or bone resorption.
4. Bone resorption, which may contribute to the deterioration of fixation and eventual loosening of the implant.
5. Periarticular calcification or ossification, which may lead to a decrease in joint mobility and range of motion.
6. Traumatic arthrosis of the ipsilateral knee secondary to intraoperative positioning, of the extremity during surgery.
7. Subluxation or dislocation.

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HOW SUPPLIED

The components are individually packaged and supplied STERILE. All metal components are sterilized using radiation. Remove from the package using accepted aseptic technique only after the correct size has been determined.

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Components may not be resterilized by the hospital because of the possibility of damaging the articulating and interfacing surfaces of the femur. The package should be opened only after the correct size has been determined, as opened packages may not be returned for credit.

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© DePuy Synthes 2017. All rights reserved.
DIST**DePuy Orthopaedics, Inc.**