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

Becton, Dickinson and Company, 1 Becton Drive, Franklin Lakes, NJ 07417 USA (США)

Becton, Dickinson and Company, Bellver Industrial Estate, Bellver Way, Roroborough, Plymouth, PL6 7BP United Kingdom (Великобритания)

Место производства
Becton, Dickinson and Company (BD), 1575 Airport Road, Sumter, SC 29153 USA (США)
Nipro Medical Industries, Ltd, Tatebayashi Plant 2-19-64 Matsubara, Tatebayashi-shi, Gunma, 374-8518, Japan (Япония)
Nipro (Thailand) Corporation Limited, 10/2 Moo 8 Bangnomko, Sena, Phra Nakhon Si Ayutthaya 13110, Thailand (Таиланд)

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

VDP40205RU-03
05/2017

BD Vacutainer®
Safety-Lok™
Single-Use Blood Collection Set
with or without Luer Adapter
for Difficult Vein Access

EN English


For Blood
Collection

For Blood
Collection
Using a
Syringe

For Short
Term IV
Administration
(up to 2 hours)

Intended Use

The BD Vacutainer® Safety-Lok™ Blood Collection Set is a winged blood collection needle and flexible tubing for venipuncture to collect blood specimens from patients. It is also indicated for the intravenous administration of fluids and may be used for any patient population with consideration given to patient size, appropriateness for the solution being infused and duration of therapy. It is also recommended for use in patients with small veins.

Product Description

The BD Vacutainer® Safety-Lok™ Blood Collection Set is a sterile winged blood collection needle with flexible tubing, a female luer adapter, an optional male luer adapter and a protective safety shield for covering the used needle prior to disposal. The male adapter contains a non-patient needle end for puncturing the stopper of an evacuated blood collection tube. Reorder numbers without a male luer adapter are provided with a protective cap on the end of the female luer adapter. The flexible wings are color coded according to needle gauge. The safety shield minimizes the possibility of an accidental needlestick if manually activated following use.

Catalog	Needle Gauge	Needle Length	Tubing Length	Luer Adapter	Color
367246	21G (0.8mm)	0.75" (19mm)	12" (305mm)	None	Green
367247	23G (0.6mm)	0.75" (19mm)	12" (305mm)	None	Light Blue
368383	25G (0.5mm)	0.75" (19mm)	12" (305mm)	None	Dark Blue
367282	21G (0.8mm)	0.75" (19mm)	7" (178mm)	With	Green
367286	21G (0.8mm)	0.75" (19mm)	12" (305mm)	With	Green
367284	23G (0.6mm)	0.75" (19mm)	7" (178mm)	With	Light Blue
367288	23G (0.6mm)	0.75" (19mm)	12" (305mm)	With	Light Blue
367295	25G (0.5mm)	0.75" (19mm)	7" (178mm)	With	Dark Blue

Safety Precaution and Warnings

Precautions:

- Practice Universal Precautions. Use gloves, gowns, eye protection, other personal protective equipment, and engineering controls to protect from blood splatter, blood leakage, and potential exposure to bloodborne pathogens.
- Read directions carefully before use.

Warnings

- Handle all biologic samples and blood collection "sharps" (lancets, needles, luer adapters, and blood collection sets) in accordance with the policies and procedures of your facility. Obtain appropriate medical attention in the event of any exposure to biologic samples (e.g., through a puncture injury), since samples may transmit viral hepatitis, HIV (AIDS), or other infectious diseases. Utilize any safety engineered feature if the blood collection device provides one. Discard all blood collection "sharps" in puncture resistant biohazard containers approved for their disposal.
- Reuse may lead to infection or other illness/injury.

General Instructions for Use

For Blood Collection

- Check to ensure that the Female Luer Adapter (A) is securely attached to the Male Luer Adapter (B) and that the Male Luer Adapter (B) is securely attached on the Holder.
- Remove needle sheath.
- Perform venipuncture by holding wings.
- Withdrawing Safety-Lok™ Blood Collection Sets:** Apply pressure to site using fingers as shown.
 - Grasp the translucent yellow safety shield grip area with the thumb and index finger while at the same time grasping the tubing.
- Safety-Lok™ Shield Activation:**
 - One-Hand Technique:** Hold tubing in hand and advance translucent yellow safety shield with thumb and index finger until a click is heard indicating that the safety shield is locked in place over needle tip.
 - Modified One-Hand Technique Withdrawal Activation:** Withdraw blood collection set by grasping the translucent yellow safety shield grip area with the thumb and index finger. With opposite hand, grasp tubing between thumb and index finger while pushing the yellow safety shield forward until a click is heard indicating that the needle is completely covered and the safety shield is locked in place.
- Dispose of all used materials in an appropriate container and wash hands.

For Blood Collection Using a Syringe

- Remove luer cap and attach syringe. Follow Blood Collection steps 2-6.
Note: Use the BD Vacutainer® Blood Transfer Device to transfer blood from the syringe to the evacuated tube.

For Short Term IV Administration (up to 2 hours)

- Remove luer cap.
- Attach IV line or syringe to the Female Luer Adapter.
Remove needle shield.
- Prime set in accordance with recommended procedure.

Priming Volume		
Gauge	Tubing	Volume (ml)
21G (0.8mm)	12" (305mm)	0.337
23G (0.6mm)	12" (305mm)	0.332
25G (0.5mm)	12" (305mm)	0.331

- Perform venipuncture by holding wings. Administer medication.
 - It is recommended that this device be changed per the institution's procedures, but no longer than 2 hours.
- Safety-Lok™ Shield Activation:** See Section 5 under Blood Collection (above).
- Dispose of all used materials in appropriate container and wash hands.

Expiration Date

Expiration date is indicated on the package. Product expires on the last day of the month and year indicated.

Storage and Transportation Conditions

Keep away from sunlight. This product does not have specific storage temperature requirements.

Technical Services Hotline Number

Russia: 8-800-200-80-40, operating 09:00-18:00 Moscow time, working days.

Becton, Dickinson and Company, 1 Becton Drive, Franklin Lakes, NJ 07417, USA

Becton, Dickinson and Company, Belliver Industrial Estate, Belliver Way, Roborough, Plymouth, PL6 7BP, United Kingdom

Manufacturing Site

Becton, Dickinson and Company (BD), 1575 Airport Road, Sumter, SC 29153, USA
Nipro Medical Industries, Ltd, Tatebayashi Plant 2-19-64 Matsubara, Tatebayashi-shi, Gunma, 374-8518, Japan
Nipro (Thailand) Corporation Limited, 10/2 Moo 8 Bangnomko, Sena, Phra Nakhon Si Ayutthaya 13110, Thailand

Regional Representative Office

LLC Becton Dickinson Vostok
Dvintsev Street, 12, Bldg. 1, Block C, 4th Floor
127018 Moscow, Russia
Tel: +7 (495) 775 8582
Fax: +7 (495) 775 8583

Conformity Mark, Registration Number and Issuance Date

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
Symbol and Mark Key

- EC REP** Authorized Representative
- LOT** Batch Code
- 2** Do Not Reuse
- Do Not Use if Package is Damaged**
- Keep Away from Sunlight**
- Latex Free**
- Manufacturer**
- STERILE EO** Method of Sterilization Using Ethylene Oxide
- Use By**
- Consult Instructions for Use**
- Rx Only (USA)**

© 2017 BD, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.
U.S. Patent No. 5,120,32

VDP40205RU-03
05/2017

BACK

Color Break Information

10.92" x 11" (277.368 x 279.4mm)

3 panels - 2 folds

2 Sided Printing

100% PMS 2755


Ensure Print Setting is set to "Spreads"

Graphical dimensions and locations must adhere to electronic file requirement

BD - Preanalytical Systems • Franklin Lakes, NJ						
REV LEVEL	ALOA #	ECR/O #	DESCRIPTION	DESIGNER	DATE	SAP #
01	00802-01	ECR/O 5-51079	This insert is for the Russian region. Translated in Russian from source document VDP40137-08.	Vasanth	11.18.2015	N/A
02	00869-01	ECR/O 5-62013	Changed RC date of issue and Copyright year. Updated new BD Logo. Added SAP Number. Changed the hot-line working hours to 09:00 – 18:00 in EN and RU.	Vasanth	04.07.2016	500015467
03	00960-01	ECR/O 5-93589	Updated the names in heading. Changed the RC issue date, Trademark year updated, Removed the SAP no., Changed the heading and Sub-heading to FS Albert. Added the BD logo in Russian heading.	Bikshan Roy	16/May/2017	NA
PACKAGING LEVEL			Product Circular (Graphics to update this circular as US/CE/Parent circulars get updates)	DRAWING #	VDP40205RU-03	

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The provided copy is correct and true

Becton, Dickinson and Company
* Becton Drive
Franklin Lakes, New Jersey 07417
USA

Becton, Dickinson and Company (BD), USA
International Regulatory Affairs Manager
Jean-François Corbin

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