STORZ

KARL STORZ — ENDOSKOPE

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**SILVER SCOPE® Video-Duodenoskope,
Modelle 13885 x**

Instruction Manual

**SILVER SCOPE® Video Duodenoscopes,
Models 13885 x**

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Thank you for your confidence in the brand name KARL STORZ. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of modern, high-quality items of equipment from KARL STORZ.

This manual is intended to serve as an aid in the proper use, cleaning, disinfection and, where relevant, sterilization of videoscopes manufactured by KARL STORZ. The videoscopes are sensitive precision instruments which require particularly careful reprocessing. To prevent damage during reprocessing, guidance is provided here for the disassembly, cleaning, disinfection, gas sterilization and assembly of KARL STORZ videoscopes. For this reason, read these instructions carefully and keep them as a ready reference.

KARL STORZ is constantly working on the further development of all products. Please understand that changes to the scope of supply, design, equipment, and technology are possible for this reason. Therefore, no claims may be derived from the information, illustrations, and descriptions in this manual.

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Symbol	Explanation
	Consult instructions for use
	Serial number
	Manufacturer
	Fragile, handle with care
	Keep dry
	Temperature limit
	Keep away from sunlight
	Do not reuse
	Applied part type BF, defibrillation proof, IEC 60601-1
	Potential equalization connector
	DANGER: Risk of explosion if used in the presence of flammable anesthetic gases.
	Ground
	Alternating current

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
Symbol	Explanation
	Input
	Setup Wizard
	Signifies that this feature/option is available in the Options Menu
	Signifies that this feature/option can be controlled via the User Menu
	Signifies that this feature/option can be controlled via the buttons on the videoscope
RxOnly	Federal (USA) law restricts this device to sale by or on the order of a physician.

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- ① Up/down deflection wheels
- ② Up/down control brake
- ③ Programmable remote control buttons
- ④ Suction valve
- ⑤ Air/water valve
- ⑥ Working channel with rubber lip valve/biopsy
valve
- ⑦ Right/left deflection wheels
- ⑧ Right/left deflection brake
- ⑨ Locking screw of Albarran module
- ⑩ Albarran control
- ⑪ Nozzle for lens cleaning/irrigation/insufflation
- ⑫ Light cable
- ⑬ Lens
- ⑭ Albarran lever
- ⑮ Albarran module
- ⑯ Locking device of Albarran module (reusable)
- ⑰ Supply tube
- ⑱ Cap with connection for leakage tester
- ⑲ Water bottle connection
- ㉑ Light cable connection
- ㉑ VEA (videoscope adaptor) connection
- ㉒ Suction/aspiration port
- ㉓ Air pin
- ㉔ Connection for leakage test
- ㉕ Lanyard
- ㉖ Slider clamp

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1 Intended use	3
2 Warnings and cautions	4
3 User qualification	5
4 Safety precautions when using the videoscopes	5
5 Handling	10
5.1 Compatible video processor units/light sources	10
5.2 Connecting to a light source, for example the KARL STORZ XENON 100 SCB	11
5.3 Operating the videoscope	12
5.3.1 Preparation of the light source	12
5.3.2 Preparation of the processor unit	12
5.3.3 Connecting the light conductor, air/water connections and the suction channel	12
5.3.4 Deflection of the distal end	13
5.3.5 Suction	13
5.3.6 Insufflation and light	13
5.3.7 Water irrigation for lens cleaning	14
5.4 Preparation for use	14
5.4.1 Testing the various functions	14
5.4.2 Functional test of the Albarran lever	15
5.4.3 Handling the endoscope correctly	18
5.4.4 Additional recommendations for use	19
5.4.5 Operating the endoscope buttons	20
5.4.6 Disassembly and assembly of the Albarran system	21
5.4.7 Disassembly	21
5.4.8 Assembly	25
5.4.9 Test for proper functioning	27
6 Reprocessing	28
6.1 General	28
6.2 General cleaning procedure	30
6.2.1 Preliminary cleaning	31
6.2.2 Leakproof test	32
6.3 Channel system of a KARL STORZ duodenoscope with a replaceable Albarran module	34
6.4 Procedure for manual processing	35
6.4.1 Rinsing the channels	35
6.4.2 Cleaning the outside surface	37

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6.4.3	Brushing the lumina of the Albarran components.....	38
6.4.4	Brushing other endoscope lumina.....	39
6.4.5	Brushing the reusable Albarran module.....	39
6.4.6	Brushing the endoscope channels.....	40
6.4.7	Rinsing the channels.....	41
6.4.8	Neutralization.....	43
6.4.9	Drying.....	44
6.4.10	Disinfection.....	45
6.5	Reprocessing in the washer/disinfector.....	47
6.5.1	Preliminary cleaning and leakage test.....	48
6.5.2	Manual cleaning prior to machine reprocessing in the washer/disinfector.....	48
6.5.3	Machine reprocessing.....	51
6.6	Cleaning of accessories.....	57
6.7	Reprocessing of water bottle 13992 BS.....	58
6.8	Checking the reprocessing procedure.....	58
6.9	Sterilization.....	59
6.9.1	Gas sterilization.....	59
6.10	Storage.....	62
7	Storage/transport conditions	62
8	Maintenance	63
9	Important information	63
10	Servicing and repair	64
11	Responsibility	65
12	Warranty	65
13	Directive compliance	65

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1 Intended use

The KARL STORZ video duodenoscopes are optical instruments for medical, diagnostic and endoscopic treatment in the upper gastrointestinal tract of humans. Accessories are inserted into the pancreaticobiliary tract via the working channel. The duodenoscopes are inserted orally through the mouth, via the oesophagus and the stomach into the patient's duodenum. ERCP: Examination of the bile ducts, gallbladder and pancreas using a duodenoscope.

Contraindications:

There are no known contraindications relating directly to the product. The videoscope must not be used if, in the opinion of the responsible physician, such an application would endanger the patient, e.g. due to the patient's general condition, or if the endoscopic method as such is contraindicated.

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2 Warnings and cautions

Please read this manual carefully and follow the instructions exactly to avoid putting your patients, personnel or yourself at risk. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, the text should be carefully reviewed to ensure the safe and effective operation of this product. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be endangered. Failure to observe a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or safety precautions must be followed to avoid possible damage to the instrument.

NOTE: A Note indicates special information about operating the instrument, or clarifies important information

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3 User qualification

KARL STORZ videoscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the endoscopic technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the videoscopes.

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WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended use and the interface specifications of the devices used in combination permit this.

The instruction manuals and interface specifications of the medical products used in combination must be observed precisely.

WARNING: The pressure compensation / EO cap 11025 XE must be removed from the connection piece for:

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WARNING: Do not allow any live components or instruments to come into contact with the videoscope.

WARNING: The videoscope must not be used during the discharge of a cardiac defibrillator. Remove the videoscope from the operation area before discharge takes place.

WARNING: Danger of explosion if flammable gases are used in the immediate proximity of the videoscope. The videoscope must not be used in the presence of flammable or combustible media.

WARNING: High energy light radiated through videoscopes may give rise to high temperatures in front of the light outlet, at the end faces of the light guide and at the tip of the endoscope.

Prolonged direct contact of light-emitting surfaces on tissue can lead to burns; this must therefore be avoided at all costs.

Even after the proximal connector is removed from the light source, surfaces will remain hot for a period. (Risk of burns for user on the light cable).

WARNING: Before removal from the patient, the distal tip of the videoscope must be returned to the neutral, non-deflected position (L and U on top of each other and positioned between the valves) to prevent damage to the instrument and/or injury to the patient.

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WARNING: Always unplug the videoscope before carrying out any maintenance and cleaning work on it.

WARNING: Risk of infection during shipment! Contaminated videoscopes may cause serious or life-threatening illnesses during transportation or on delivery to the recipient. Please cover the carrying case with the protective sheet 13990 SFN before sending contaminated instruments back to KARL STORZ.

WARNING: If the videoscope should malfunction during an intervention on the patient, stop the intervention immediately. Place the distal tip in a neutral, non-deflected position and remove the videoscope slowly and carefully from the patient.

WARNING: The use of HF probes, HF electrodes, HF instruments or other energy-emitting probes is not permitted if explosive gases (for example, explosive/flammable anesthetics) are present.

WARNING: Patient leakage currents can accumulate if endoscopes and electrically powered endoscope accessories are used together.

WARNING: Before turning on the high voltage/electrical supply, ensure that there is no contact with metal parts of the endoscope or any other electricity-conducting accessory, and that the suction unit has no contact with conductive liquids.

WARNING: HF probes, HF electrodes, laser probes and other energy-emitting probes should only be used under constant, visual, endoscopic control.

Activation should only occur when the tip of the energy-emitting probe is outside the working channel and within the endoscope's field of vision.

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CAUTION: Only use accessories (guide probes, cleaning brushes, etc.) that are specifically suited for the intended use with videoscopes and can be employed in line with generally accepted technical standards and safety guidelines.

CAUTION: Only use accessories with flexible ends. Other types of accessories may cause damage to videoscopes.

CAUTION: The pressure compensation/EO cap 11025XE must be placed on the connection piece for:

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CAUTION: Do not steam sterilize (autoclave) videoscopes.

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CAUTION: The manual 'Preparation and care' must be followed. The cleaning, disinfection, and gas sterilization procedures are explained there in detail.

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NOTE: The patient leakage current limit of the CF application part may be exceeded when CF application parts (application on heart) are used in combination with other BF application parts. In this case, the reduction of the patient leakage current should be effected using only CF application parts as far as this is possible.

NOTE: The endoscope does not offer any electric insulation against high frequency (HF) voltages. Protection against high frequency voltages should be ensured through the use of correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.

NOTE: Check the videoscope prior to each surgical procedure to ensure that it functions correctly and has been properly cleaned and disinfected.

NOTE: Any damage to the videoscope resulting from incorrect operation is not covered by the guarantee.

NOTE: Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided.

NOTE: Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

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5 Handling

5.1 Compatible video processor units/light sources

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NOTE: The KARL STORZ videoscope should not be used in combination with video processor units and light sources from other manufacturers. The use of non-compatible devices can lead to patient injury or unit damage. It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces, e.g., voltage-carrying elements of other devices, at any time.

WARNING: The use of duodenoscopes 13885 xKS in combination with the GASTRO PACK™ is only recommended to a limited extent, as the high-intensity light emitted at the light outlet, at the end faces of the light guide and at the tip of the endoscope may give rise to high temperatures. Prolonged direct contact of light-emitting surfaces on tissue can lead to burns; this must therefore be avoided at all costs.

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
5.2 Connecting to a light source,
for example the KARL STORZ
XENON 100 SCB

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CAUTION: Depending on the video processor used, a suitable videoscope adaptor must be employed.

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5.3 Operating the videoscope

5.3.1 Preparation of the light source

Please read the full instructions in the instruction manual of the respective light source.

1. Install the light source on a solid, even surface.
2. Plug the power cord into the power socket provided.

5.3.2 Preparation of the processor unit

Refer to the relevant instruction manual for the processor unit concerned.

5.3.3 Connecting the light conductor, air/water connections and the suction channel

1. Insert the supply connector into the light source until it audibly clicks into place (Fig. 1).
2. Connect the videoscope adaptor (VEA) cable to the endoscope. When doing so, take note of the white markings on the endoscope and the VEA cable (Fig. 2).
3. Attach the connector from the water bottle tube to the endoscope supply plug (irrigation bottle connector, Fig. 3). The water bottle holder automatically positions itself correctly and cannot be adjusted manually.
IMPORTANT: To fill the water bottle with sterile water, remove the bottle from the light source and unscrew the lid. Only fill the water bottle 2/3 full with sterile water (1/3 is required to cushion air pressure).
4. Connect the suction tube (patient tube) from the suction pump or central suction system to the supply plug connector (Fig. 4).

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5.3.4 Deflection of the distal end

Rotating the large deflection wheel, labeled 'U/D', (up/down) deflects the distal tip. Rotating clockwise or forwards causes the tip to move downwards. Rotating counterclockwise or backwards causes the tip to move upwards. Right and left movements are controlled using the smaller deflection wheel (labeled 'L/R'). Clockwise (or forwards) = right/counterclockwise (or backwards) = left (Fig. 5).

WARNING: Do not rotate the deflection wheels forcibly; this could damage the instrument or traumatize the patient! Before attempting to introduce or remove the instrument, release the deflection wheel brakes!

5.3.5 Suction

Switch on the suction pump. To perform suction, press the suction valve (with the red ring) on the endoscope (Fig. 6). Suction can also be used when the biopsy forceps are placed in the biopsy channel and the rubber lip valve is mounted. During the procedure, the rubber lip valve/biopsy valve seals off the forceps or the accessories.

5.3.6 Insufflation and light

Start the light source, switch on the insufflation pump. For insufflation, cover the opening on the blue valve with a finger (do not depress the valve, since this starts water irrigation) (Fig. 7).

NOTE: To avoid patient discomfort, use insufflation sparingly. Excessive insufflation may cause cramps, spasms and injuries!

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5.3.7 Water irrigation for lens cleaning

Pressing on the air/water valve starts the water flow and/or lens cleaning. (Fig. 8).

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5.4 Preparation for use

WARNING: Before being used for the first time, the instrument should undergo a complete cleaning cycle, including preliminary cleaning.

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5.4.1 Testing the various functions

1. Prepare light sources as described in 2.1.
2. Connect supply connector, water bottle and suction tube according to 2.3.
3. Inspect the supply tube for damage (such as breaks, cracks, twisted or squashed sections...) and check the supply connector for loose connections.
4. Examine the surface of the flexible insertion sheath visually and by touch for cracks, cuts, bulges, foreign bodies or other damage.

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5. Testing deflection – upward/downward control

The large deflection wheel controls up and down movements. To test up and down movements, release the up and down brake. Slowly rotate the large deflection wheel clockwise or forwards; the distal end should move downwards. Rotating the deflection wheel counterclockwise or backwards causes upward movement.

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6. Testing deflection - right/left control
7. The smaller/outer deflection wheel controls right/left movement. To test right/left movement, release the right/left brake. Slowly rotate the small deflection wheel clockwise or forwards; the **tip of the endoscope** should move to the right. Rotating the deflection wheel counterclockwise or backwards causes upward movement.
8. Testing deflection – brake function
The deflection wheels should turn freely when the brakes are released.
9. Release the brakes and test their action (see 5/6). Activate the brakes by moving the brake lever counterclockwise – ensure that the angled position of the distal endpiece is fixed in the desired position.

Do not twist or bend the deflection piece by hand.

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5.4.2 Functional test of the Albarran lever

If the lever of the Albarran control ⑩ on the handle points toward proximal, the Albarran mechanism is retracted.

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The mechanism must extend when the lever is pushed toward distal.

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10. Check the air supply. Switch on the insufflation pump and hold the distal end of the instrument under sterile water. Cover the air control opening in the air valve (blue ring) with a finger; bubbles should emerge from the air nozzle at the distal end.

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11. Check the water supply: Depress the air/water valve (blue ring) fully. Water should flow out of the water nozzle (distal) and rinse the lens clean. IMPORTANT: When you first actuate the valve, it may take a couple of seconds for water to start flowing.

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12. Test the suction function. With the unit correctly installed and the suction pump or central suction system switched on, immerse the distal tip in a container of sterile water once again. Depress the suction valve (red ring) fully. Water should be suctioned through the suction channel of the instrument via the connected suction pump and into the collecting bottle of the suction equipment.

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13. Test the biopsy channel by introducing biopsy forceps.

CAUTION: It should be possible to guide the biopsy forceps smoothly and without problems through the entire channel.

CAUTION: Full ($>180^\circ$) deflection can lead to limitation in movement, depending on the accessories. If necessary, adjust the deflection slightly.

CAUTION: Accessories such as biopsy forceps that are kinked or bent will not function properly, are difficult to guide through the channel and may cause severe damage to it. They must therefore be replaced and not used further!

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
Handling the endoscope correctly

With the control handle of the endoscope in the left hand, the thumb, index and middle fingers can easily control the large deflection wheel. In this position, the suction valve can be depressed with either the index or the ring finger. The air/water valve can also be covered with these fingers. The remote control button (e.g., Still images) is also within reach.

The smaller projecting deflection wheel can be controlled with the fingers of the left hand, or alternatively with the thumb and fingers of the right hand.

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5.4.4 Additional recommendations for use

1. Use a lubricant approved by KARL STORZ (see list on our website www.karlstorz.com under 'Hygiene') for the distal end of the insertion tube.
2. Observe the lumen the whole time. If you become disoriented, withdraw the instrument gently until vision is restored.

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NOTE: If pollution causes a blurry picture during the procedure, this can normally be rectified by a combination of water and insufflation. Excessive water or air can be removed by applying simultaneous pressure on the suction valve.

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5.4.5 Operating the endoscope buttons

Depending on the video processor unit used, the buttons can have dual functions which are triggered by pressing the button briefly or for an extended period. For more detailed information on button assignment please consult the manual for the camera control unit used, as listed in section 5.1.

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NOTE: If the angling of the tip can no longer be properly controlled, the examination must be ended immediately. Release the brakes and place the tip in the neutral position. Position the deflection wheels so that 'U' and 'L' are between the valves (see page 13/17), then carefully withdraw the instrument back, whilst continuing to monitor the procedure.

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5.4.6 Disassembly and assembly of the Albarran system


As conventional duodenoscopes are difficult to clean and maintain due to their design, particularly around the area of the Albarran lever, KARL STORZ has developed a modular Albarran system.

The patented Albarran system is easy to disassemble, which facilitates thorough cleaning and disinfection. Furthermore, the reusable Albarran module is autoclavable. It is no longer necessary to disassemble the entire endoscope to eliminate technical malfunctions of the Albarran unit. As a rule, repairs can be carried out by disassembling the Albarran unit.

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5.4.7 Disassembly

Undo the locking screw ⑨ on the proximal control element by turning counterclockwise several times...

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
On the proximal end control element, you will find a slider clamp which secures the Albarran slide in place. Release the slider clamp by exerting gentle pressure on the lever with your finger to press in the direction of the arrow.

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Then push the slider clamp upwards and release fully from the Albarran slide.

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Now use your fingers to remove the Albarran slide from the proximal end of the handle.

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The wire of the Albarran module is exposed.

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Reusable Albarran module:

There is a silver-colored locking screw with a red dot (in the 12 o' clock position) at the distal tip of the insertion tube.

Single-use Albarran module:

Follow the instruction leaflet supplied.

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
Undo this screw by turning the Albarran key 180° counterclockwise.

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The red dot must now be located exactly above the second red dot in the 6 o' clock position.

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
The unlocked Albarran module can now be pulled out carefully towards distal. Make sure that the wire at the proximal end is exposed.

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Roll up the thin cable without kinking or bending it.

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With regard to the single-use Albarran module, please refer to the instruction leaflet 97000175.

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
5.4.8 Assembly

Now, on the distal side, push the wire of the Albarran module back into the Albarran channel completely until it emerges at the proximal end of the control element.

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Den beiliegenden Beipackzettel Single-Use-Albarran-Modul 97000175 beachten.

Reusable Albarran module

Using 2 fingers, fix the Albarran cap securely to the insertion tube. There must be no gap (see arrow) between the Albarran cap and the endoscope. Turn the screw 180° clockwise and lock. The two red dots must again be located at 6 and 12 o'clock.

Single-use Albarran module:

Follow the single-use Albarran module 97000175 instruction leaflet supplied.

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NOTE: Make sure that the cap is securely in place and cannot be pushed toward distal.

Place the endoscope's insertion tube flat on the work surface.

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
Push the Albarran slide into the Albarran channel of the proximal end of the handle. The device for securing the clamping hook must face towards the deflection wheels.

To secure the slider clamp, the device must face towards the deflection wheels.

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
Placing the slider clamp on the Albarran slide:
1. Hook in the upper lug by pressing sideways in the direction of the arrow.

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2. Secure the clamping hook in place by pressing on the lever of the hook on the Albarran slide.
3. Secure the underlying lug by pressing laterally on the lever of the hook on the Albarran slide.

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
Now push on the chuck as far as it will go. Make sure that the wire is positioned centrally in the chuck and can still be moved freely.

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
Move the Albarran control ⑩ to the proximal end position and retighten the blue locking screw ⑨ by turning clockwise.

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5.4.9 Test for proper functioning

If the lever of the Albarran control ⑩ on the handle points toward proximal, the Albarran lever ⑭ is open .

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Moving the Albarran control ⑩ towards distal closes the Albarran lever ⑭.
Make sure that the Albarran control runs smoothly and that the end positions of the Albarran lever ⑭ are reached

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6 Reprocessing

6.1 General

WARNING: The stated cleaning, disinfection and sterilization methods alone do not guarantee successful sterilization of the endoscopes. This can only be achieved if personnel follow recognized and validated reprocessing methods.

WARNING: The local occupational safety regulations must be followed and applied for all cleaning and disinfection work on contaminated endoscopes and instruments.

WARNING: Remember that sterilization is only successful if

- the instruments have been freed of all organic material and chemical residues.
- only recognized and validated methods are used.
- any deviation from recommended procedures have been validated by the user.
- the recommended sterilization parameters are maintained. These are only valid with cleaning, disinfection and sterilization equipment that is properly maintained and calibrated.

CAUTION: The 'Cleaning, disinfection, sterilization, and care of KARL STORZ instruments' instructions must be followed. The procedures for reprocessing are explained therein in detail. The list of KARL STORZ approved cleaning/disinfection materials can be found on the Internet at 'www.karlstorz.com' or 'www.karlstorz.de' below Hygiene or in the appendix of this manual.

CAUTION: When preparing and using cleaning and disinfecting solutions, the chemical manufacturers' instructions on concentration and exposure times must be followed precisely. Prolonged immersion and an incorrect concentration can lead to damage. Bear in mind the microbiological range of action of the chemicals used.

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CAUTION: Never clean videoscopes with metal brushes.

CAUTION: Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, should be carried out in consultation with the manufacturers of the endoscope, the chemicals and the washer/ disinfectant. Then this change in method is usually completely harmless.

CAUTION: Videoscopes must not be cleaned in an ultrasonic cleaning bath.

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CAUTION: Never immerse videoscopes in normal saline solution, since even brief contact can lead to pitting and corrosion.

CAUTION: Videoscopes may only be rinsed with microbiologically clean/sterile water during manual reprocessing.

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6.2 General cleaning procedure

Cleaning should only be carried out by authorized trained staff. All safety regulations should be observed (gloves, protective glasses, mouth and nose protection, clothing, German Technical Rules for Biological Agents (TRBA 250)).

Thorough cleaning of the endoscope and all its accessories after every use ensures problem-free working and the best results for disinfection or sterilization procedures.

The national guidelines and recommendations for the reprocessing of instruments (for example RKI, ESGE/ESGENA, ASGE, LANGS, SAGES, etc.) should be observed and the effectiveness verified.

CAUTION: Before immersion in cleaning solution, check the instrument for leaks. Otherwise fluid penetrating the videoscope could cause damage.

CAUTION: Do not use autoclaving for the endoscope or accessories (with the exception of accessories that are specifically designated as autoclavable, e.g., forceps and Albarran module).

NOTE: The endoscope must be cleaned immediately after every use. If organic material is allowed to remain and dry in the instrument, complete or partial blockage of the channel may result. Thorough cleaning is a prerequisite for effective reprocessing. This will avoid clots/blockages and the resulting repairs and instrument failures. For reasons of material compatibility, only chemicals approved by KARL STORZ should be used for cleaning and disinfection (see current list in the Appendix, p. 60 or on our website www.karlstorz.com under 'Hygiene').

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
6.2.1 Preliminary cleaning

1. The preliminary cleaning of the endoscope is performed at the place of operation. The main cleaning and disinfection is carried out in a separate reprocessing room.
2. Immediately after use, remove any heavy soiling (secretions, etc.) and organic material from the outer casing of the insertion tube with a disposable cloth moistened with cleaning solution. Then submerge the distal end in a container with cleaning solution and activate the suction to rinse the working channel. Activate the suction valve and suction at least 100 ml of the cleaning solution through the suction channel.
3. Use the cleaning valve 13991 SRV in place of the air/water valve. Do not use this valve during examination.
4. Check that the water bottle is filled (max 2/3). Start the air supply (insufflation pump). Water and air are now passed through both the air channel and the water channel depending on the position and operation of the cleaning valve. (The process should run for at least 10 seconds).
Air channel = Position 1
Water channel = Position 2
5. It is important that both valve positions (activated/not activated) are used repeatedly for several seconds. This is the only way to ensure that both channels are rinsed through.
6. Switch off video system.
7. Pull out the videoscope adaptor connector from the videoscope adaptor (VEA) at the supply connector of the endoscope when the video processor unit is switched off. Put on the sealing cap for the video cable connector and lock it with a 1/4 turn.

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8. Disconnect the endoscope from the video processor unit and light source. Place the cleaning valve in a cleaning solution.
9. Then protect and transport to the reprocessing station

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6.2.2 Leakproof test

Preparation of the leakage test – dry

i NOTE: Place the cap on the videoscope adaptor (VEA) connector and screw it on tightly. Ensure that the cap is secured properly.

1. The connection tube of the leakage tester and the ventilation connector on the endoscope must be dry! Connect the leakage tester connector to the supply connector and lock by turning clockwise 90°.

i NOTE: The rubber seal in the connector adaptor of leakage tester 13242 LX should be treated occasionally with instrument oil.

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Before every immersion in liquid, test the endoscope for leaks.

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- Inflate the leakage tester until the maximum test pressure of 160 – 180 mmHg (blue area) is reached. The deflection rubber on the distal deflection part of the endoscope expands slightly. Move the deflection wheels on the control piece in all four directions as far as possible and check the manometer to see whether the pressure remains constant or falls (dry test).
- Wait at least 1 minute.

Performing the leakage test – liquid

- Immerse the device in water until it is completely covered. When doing so, ensure that the leakage tester is not immersed under water but remains dry. Perform the leakage test as described under point 3 and check whether bubbles rise up from the device. If no more bubbles can be seen rising from the device after a period of at least 3 minutes and, in addition, there is no loss of pressure on the manometer, the device is leakproof and can be used.

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If a pressure loss occurs, this means that the videoscope leaks and may no longer be used or immersed or reprocessed by machine.

CAUTION: In this event, fill out the accompanying repairs form (see Appendix), roughly clean and remove any organic soiling/residue from the endoscope, mark it as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the enclosed transport protection sheet 13990 SFN.

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
**6.3 Channel system of a
KARL STORZ duodenoscope
with a replaceable Albarran
module**

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6.4 Procedure for manual processing

Please adhere to the national guidelines, e.g., ESGENA guidelines, DEGEA/BfArM KRINKO recommendations, ASGE, LANGS, SAGES etc. for cleaning, disinfection and sterilization, as well as to the instructions of the chemical manufacturer concerning concentration and exposure times.

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1. Fill a basin with cleaning solution.

WARNING: *Risk of material damage/risk of infection: Improper preparation of the cleaning solution may cause damage to the instrument or prevent a complete decontamination of the surface. Refer to the chemical manufacturers' instructions on concentration and exposure times.*

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WARNING: *Risk of infection: Make sure there is sufficient cleaning fluid in the container to cover the endoscope and all disassembled parts completely.*

2. Before the endoscope is immersed in the container, remove the Albarran module and, if relevant, dispose of the single-use Albarran module (13991 AM) (see section 3.5.1)

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
6.4.1 Rinsing the channels

1. Attach the double valve closure 13991 CD to the suction and air/water valves on the handle.

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2. Attach the following adaptors to the corresponding connections:

Adaptor	Connection
13991 SS	Working channel entrance/ biopsy valve ⑥
13991 AKRA	Suction port ②
13991 AKRA	Albarran channel Albarran lever (the Albarran control ⑩ must be closed in order to be able to connect the cleaning adaptor 13991AKRA.)
6011590	LUER Lock cap (3x)

3. Close all adaptors with LUER Lock caps (see Fig. 45-47).

NOTE: The adaptor 13991 AKRA on the Albarran closure must not be closed.

- Completely immerse the endoscope in cleaning solution. Through filling in a targeted manner, it must be ensured that even surfaces with restricted access and any lumina are covered and no air bubbles are present.
- Open one of the LUER Lock caps. Screw a 60 cc syringe onto the adaptor's LUER connection (see Fig. 47).
- During the immersion phase, use the syringe to suck at least 50 ml of cleaning solution from the channel to which the syringe is connected (see Fig. 47).
- Unscrew the syringe and dispose of the suctioned off cleaning solution in a separate container.
- Remove the syringe and screw the LUER Lock closure back on.
- Repeat steps 5 to 10 for all channels.
- Connect the 60 cc syringe to the Albarran channel cleaning adaptor (13991 AKRA). Suction off at least 10 ml of cleaning solution via the channel (see Fig. 48).

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11. Remove all tubes, take off the double valve closure 13991 CD and all cleaning adaptors.

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6.4.2 Cleaning the outside surface.


1. Whilst the endoscope is immersed in the cleaning solution, clean the outside of the entire endoscope with a soft, lint-free cloth. The endoscope must be cleaned completely at least once.
2. The outside surfaces of the following parts must be cleaned at least 3 times with a brush:

Part designation	Brush
Handle including buttons	Cleaning brush 27652
Surfaces between the deflection levers and handle	Cleaning brush 27652
Suction channel entrance	Cleaning brush 27652
Suction valve, air/water valve	Cleaning brush 27652
Albarran control	Cleaning brush 27652
All the disassembled components (i.e., valves, Albarran fastening screw, chuck, double valve closure, cleaning adaptors)	Cleaning brush 27652

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6.4.3 Brushing the lumina of the Albarran components

Brushing the lumina of the Albarran components

Part designation	Brush
Locking screw	Valve brush 13991VB
Cleaning adaptor	Valve brush 13991VB
Collet chuck	Cleaning brush P/N: 110917-50
Albarran slide	Cleaning brush P/N: 110917-50

Using a 60 cc syringe, rinse the lumina of each part of the Albarran components and the cleaning adaptors with at least 10 ml of cleaning solution.

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6.4.4 Brushing other endoscope lumina

The following parts must be brushed at least 3 times. To do this, insert the brushes into each valve and rotate through 360°. Pull out the brushes and remove any visible dirt from the bristles.

Part designation	Brush
Suction valve housing	Valve brush 13991VB
Air/water valve housing (see Fig. 57)	
Working channel entrance/biopsy channel ⑥	
Albarran channel exit (see Fig. 58)	
Endoscope tip (distal)	
All other parts of the distal end which come into contact with the Albarran unit (do not brush the optical lens) (see Fig. 59).	

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6.4.5 Brushing the reusable Albarran module

Carefully clean the distal end of the Albarran module with a soft brush (13991 VB).

The following parts can be cleaned in the ultrasound bath or, alternatively, autoclaved in the steam sterilizer:

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Cleaning the reuseable Albarran module

- Use a soft brush (13991 VB) to carefully brush the Albarran module dismantled at the distal end under cold running water until there is no more visible dirt.

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- Then rinse between the head of the Albarran module and the locking screw, which closes on the duodenoscope when the Albarran module is assembled, three times for at least 5 seconds using a water gun or, alternatively, three times using a disposable syringe filled to the 60-ml mark.

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6.4.6 Brushing the endoscope channels

Insert the brushes far enough into the channel that they protrude out of the latter's distal end. Pull the rest of the cleaning brush through the channel and then rub the brush ends against each other to remove any dirt from the bristles. Repeat this brushing procedure at least three times.

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Valve	Brush
Suction channel	Long cleaning brush: 110910-50
Working channel	
Albarran channel	Long cleaning brush: 110917-50

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6.4.7 Rinsing the channels

1. Place the double valve closure 13991 CD on the suction valve and air/water valve (see Fig. 65).
2. Attach the following adaptors (see section 5.8):

Adaptor	Valve
13991 SS	Working channel entrance/biopsy valve ⑥
13991 AKRA	Suction port ②
13991 CR	Air pin ②③
13991 CR	Water bottle connection ⑨
13991 AKRA	Albarran channel Albarran lever (the Albarran control ⑩ must be closed in order to be able to connect the cleaning adaptor 13991AKRA.)
6011590	LUER LOCK cap (4x)

3. Close all adaptors with LUER Lock caps (see Fig. 64).

ⓘ NOTE: The adaptor 13991 AKRA on the Albarran closure must not be closed.

4. Completely immerse the endoscope in cleaning solution. Through filling in a targeted manner, it must be ensured that even surfaces with restricted access and any lumina are covered and no air bubbles are present.
5. Open one of the LUER Lock caps. Screw a 60 cc syringe onto the adaptor's LUER connection (see Fig. 64-65).
6. During the immersion phase, use the syringe to rinse at least 50 ml of cleaning solution

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through the channel to which the syringe is connected (see Fig. 64-65).

7. Screw the LUER Lock cap back on.
8. Repeat steps 5 to 9 for all channels.
9. Connect the syringe to the Albarran channel cleaning adaptor (13991 AKRA). Rinse at least 30 ml of cleaning solution through the channel (see Fig. 67).
10. Observe the exposure times specified by the cleaning solution manufacturer.
11. After the specified exposure time, remove the endoscope and all the disassembled components from the cleaning solution and place in a container for rinsing. Leave all cleaning adaptors on the endoscope.

NOTE: It is recommended that the cleaning solution be removed from all channels by, for example, blowing through these with air.

12. Once the cleaning work has been completed, clean, disinfect/sterilize the cleaning equipment, e.g., brushes and adaptors. Do not reuse damaged or bent brushes and single-use materials.
13. Bring the empty endoscope for neutralization.

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6.4.8 Neutralization

1. Prepare and disinfect a new container as the neutralization container.
2. Fill the container with at least 10 l of water (room temperature). Make sure there is sufficient rinsing water to cover the endoscope completely.
3. Rinse all the inside and outside surfaces as per the instructions of the cleaning agent manufacturer to ensure complete neutralization.
4. The endoscope must remain immersed during the entire neutralization procedure. Rinse the channels with a cleaning gun.

Alternatively, neutralization can be performed as described below:

1. Attach the following adaptors (see section 6.4.7):

Adaptor	Valve
13991SS	Valve working channel
13991CS	Suction port ②
13991CR	Air pin ③
13991CR	Water bottle connection ①
13991 AKRA	Albarran channel Albarran lever (the Albarran control ⑩ must be closed in order to be able to connect the cleaning adaptor 13991AKRA.)
6011590	LUER Lock cap (4x)

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2. Close all adaptors with LUER Lock caps.

i NOTE: The adaptor 13991 AKRA on the Albarran closure must not be closed.

3. Completely immerse the endoscope in water.
4. Open one of the LUER Lock caps. Screw a 60 cc syringe onto the adaptor's LUER connection.
5. During the immersion phase, rinse at least 50 ml of water through the channel with a 60 cc syringe.
6. Unscrew the syringe.
7. Screw the LUER Lock cap back on.
8. Repeat steps 2 to 7 for all channels.
9. Connect the syringe to the Albarran channel cleaning adaptor (13991 AKRA). Rinse at least 30 ml of water through the channel.

6.4.9 Drying

After successful rinsing, dry the channel system and surfaces with compressed air or a syringe using sterile compressed air.

i NOTE: Perform a visual inspection to ensure that the channel system is not dripping.

Dry the outside of the endoscope as well as all disassembled parts with a new, soft, lint-free cloth.

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6.4.10 Disinfection

1. Prepare and disinfect a new container as the disinfection container.
2. Immerse the endoscope and all disassembled parts in the disinfection solution. The endoscope must remain immersed during the disinfection procedure.
3. Attach the following adaptors (see section 6.4.7):

Adaptor	Valve
13991SS	Working channel entrance/biopsy valve ⑥
13991CS	Suction port ②
13991CR	Air pin ③
13991CR	Water bottle connection ⑱
13991 AKRA	Albarran lever (the Albarran control ⑩ must be closed in order to be able to connect the cleaning adaptor 13991AKRA.)
6011590	LUER Lock cap (4x)

4. Close all adaptors with LUER Lock caps.

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5. Completely immerse the endoscope in disinfectant solution.
6. Open one of the LUER Lock caps. Screw a 60 cc syringe onto the adaptor's LUER connection.

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7. During the immersion phase, rinse at least 50 ml of disinfection solution through the channel with a 60 cc syringe. Repeat the process at least twice.
8. Unscrew the syringe.
9. Screw the LUER Lock cap back on.
10. Repeat steps 2 to 7 for all channels.
11. Between rinses, allow the specified time (chemical manufacturers' instructions) for the disinfection solution to take effect in the channel system.
12. Then rinse the channels and surfaces of the video duodenoscope with microbiologically clean/sterile water.
13. After successful rinsing, dry the channel system and surfaces with compressed air or a syringe. To do so, use sterile compressed.

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6.5 Reprocessing in the washer/ disinfectant

All flexible KARL STORZ endoscopes can be connected to approved washers/disinfectors for semi- or fully automatic cleaning or disinfection. The national laws and regulations for cleaning and disinfection procedures must be observed. Leakage testing, manual cleaning and (where necessary) rinsing with microbiologically pure water must be carried out beforehand. The manufacturer of the cleaning machine must be consulted when selecting the reprocessing programs. Usually, special programs for flexible endoscopes are available which are suitable for decontamination and compatible with the materials in terms of protection. For machine reprocessing, make sure that the flexible endoscopes are not heated above 65°C.

CAUTION: Check the videoscope for leaks before every disinfection in semi- or fully automatic equipment.

CAUTION: Ensure that only KARL STORZ approved chemicals are used for procedures involving peracetic acid.

Lay the endoscope carefully in the tray or basket of the machine. Connect the irrigation tubes supplied with the cleaning adaptor set 13991 CK or the irrigation tubes from the manufacturer of the washer/disinfectant to the endoscope.

- The machine must fulfill all the cleaning, disinfection and material compatibility requirements (validated process).
- Make certain that the environmental conditions with regard to temperature and pressure are not exceeded (+60°C, tolerance +/- 5% and 1 bar for the channels, 0.3 bar leak test).

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- Recommendation: Before placing the endoscope in the machine, clean all the channels (instrument channels, water channels, Albarran channels, air/water channels) whether or not they have been used.
- We recommend sticking to one type of cleaning agent.

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6.5.1 Preliminary cleaning and leakage test

- The outer casing should be wiped down from top to bottom with a lint-free disposable cloth that is soaked in disinfecting agent.
- Suck in ample amounts of cleaning solution through the distal end.
- Apply the air-water supply for a few seconds every cycle
- Connect the cleaning adaptor
- Separate the flexible endoscope from the light source

CAUTION: Close the video plug with the sealing cap.

NOTE: Please note sections 6.2.1 and 6.2.2.

6.5.2 Manual cleaning prior to machine reprocessing in the washer/disinfector

1. Please adhere to the national guidelines, e.g., ESGENA guidelines, DEGEA/RKI recommendations, ASGE, LANGS, SAGES etc. for cleaning, disinfection and sterilization, as well as to the instructions of the chemical manufacturer concerning concentration and exposure times.

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2. Remove the rubber lip valve/biopsy valve, suction valve and air/water valve, place in a cleaning solution and clean using a soft disinfected brush (13991 VB) or disposable brush.

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
3. Wipe the housing, insertion sheath, and supply tube carefully with a lint-free cloth soaked in a cleaner (for suitable cleaning agents or disinfectants approved by KARL STORZ, see www.karlstorz.com). Avoid shocks and pulling movements so as not to damage the insertion shaft. Then rinse the endoscope with microbiologically pure/sterile water.

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4. Immerse the endoscope in the cleaning solution and rinse the channels using a syringe and the supplied irrigation tubes (13991 AKRA, 13991 SS, 13991 YS). The channels must be filled without any bubbles, to guarantee a complete moistening of the surfaces.

CAUTION: Only immerse instrument parts once they have successfully passed a leak-age test!

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

5. When submerged, use a suitable brush (preferably disposable) to brush below the surface of the liquid in one direction starting from the supply connector (supply connector --> valve body, valve body --> distal end, biopsy channel --> distal end). Repeat this three times.

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
6. When using reusable and disposable cleaning brushes, clean the brush head at the distal end before retracting it. Do not move the brushes back and forth; always move them in one direction only.


CAUTION: Cleaning with brushes must be carried out with the instrument immersed!

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
7. Observe the exposure time for the active cleaning solution!

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
8. Carefully clean all particles of dirt from the distal end. It may be necessary to clean the air/water valve with a brush (13991 VB) with soft bristles, in order to remove hardened contaminants. Never remove dirt with a needle or other sharp objects, because these can cause damage.

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9. Connect the Y-irrigation tube (13991 AKRA, 13991 YS and 13991 SS) and pump fresh cleaning fluid through the channels using a syringe. Afterwards, rinse with microbiologically clean/sterile water.

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Cleaning the reusable Albarran module

Use a soft brush (13991 VB) to carefully brush the Albarran module dismantled at the distal end under cold running water until there is no more visible dirt.

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- Albarran module with cable
- Collet chuck
- Button
- Albarran slide

Then rinse between the head of the Albarran module and the locking screw, which closes on the duodenoscope when the Albarran module is assembled, three times for at least 5 seconds using a water gun or, alternatively, three times using a disposable syringe filled to the 60-ml mark.

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6.5.3 Machine reprocessing



6.5.3.4. Cleaning adaptor set 13991 CK

The cleaning adaptor set 13991CK is used for validated reprocessing of KARL STORZ series 13xxx NKS/PKS videoscopes (acc. to ISO 15883 part 4). The cleaning adaptor set is not intended for any other purpose as this may damage it.

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- 1 Connection tube 13991 CL for leakage tester connection ⑮ of the videoscope
- 2 Connector (13991CL) for vent port
- 3 LUER connector (13991CL) for washers/disinfectors for endoscopes

- 4 Cleaning adaptor 13991CR for air/water (A/W) channel ⑩ of the videoscope
- 5 Connector (13991CL) for A/W channel (air pin ⑫)
- 6 Connector (13991CL) for water channel (water bottle connection ⑪)
- 7 2x LUER connectors (13991CR) for washers/disinfectors (1x LUER for air channel, 1 LUER for water channel)

- 8 Cleaning adaptor 13991CS for suction channel (aspiration) ⑪ of the videoscope
- 9 Connector (13991CL) for suction channel on the supply connector ⑪
- 10 LUER connector (13991CS) for washers/disinfectors

- 11 Connection adaptor 13991SS for biopsy/working channel ⑨ of the videoscope (with 2x double channel endoscopes)
- 12 Connector (13991CL) for biopsy/working channel ⑨
- 13 LUER connector (13991SS) for washers/disinfectors

- 14 Channel separator 13991CD for suction valve ④ and irrigation/insufflation valve (A/W) ⑤ of the videoscope
- 15 Valve pin for irrigation/insufflation valve (A/W) ⑤
- 16 Closure for suction valve ④

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WARNING: This product is delivered unsterilized and must therefore be cleaned, disinfected or sterilized prior to initial use and each subsequent reuse.

WARNING: Risk of infection by improper reprocessing: Reprocessing may only be carried out by experienced and trained staff.

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CAUTION: If the reprocessing chemicals are changed, it is essential to consult the manufacturer of the chemicals and of the washer/disinfector.

CAUTION: Please refer to our operating instructions under the section: Cleaning, disinfection and sterilization.

The following checks must be made before use:

- Visual checks of the connection tubes (patency) and the O-ring (wear)
- Check the connections are correctly and securely in place
- Check for damage.

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Assembly of the cleaning set

Plug the connector ② of the connection tube 13991 CL into the vent port ⑮ of the videoscope at the supply connector and secure by turning clockwise.

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Plug the cleaning adaptor 13991 CR ⑤ as far as it can go into the adaptor ② of the air pin of the videoscope.

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Then plug the connector ⑥ as far as it will go into the 'Water bottle connection' ⑰ of the videoscope.

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Plug the connector ⑨ of the cleaning adaptor 13991CS as far as it will go into the suction channel ⑳ at the supply connector of the videoscope.

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If necessary, remove the suction and irrigation/insufflation valves (4, 5) from the videoscope.
Insert the channel separator 13991 CD 14, with the valve pin 15, into the free opening of the irrigation/insufflation valve (4) as far as it will go so that the lock 16 lies on top of the free opening of the suction valve (4).

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
Press the slide mechanism of the channel separator 13991CD 14 downward and slide it forward until the channel separator 13991 CD 14 audibly locks into place. Check that the channel separator is securely in place.

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NOTE: If necessary, remove the cap of the biopsy/work channel 6 of the videoscope.


Connect the connection adaptor 13991 SS with the connector 12 to the biopsy/work channel 6 of the videoscope.

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Connect all open Luer connectors of the attached connection adaptors to the respective connections on the washer/disinfector.

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CAUTION: When the program has finished, the conditions in the reprocessing machine present the videoscopes with a risk of corrosion.

CAUTION: Immediately after the end of the program cycle, the videoscopes are hot as they are heated up during chemical thermal reprocessing. Consequently, the instrument's material properties are different and, as such, the reprocessed endoscope should be allowed to cool down for at least 20 minutes at room temperature before being used again.

CAUTION: Ideally, remove the videoscopes from the machine immediately after the program has finished. Then dry off the optical surfaces, if necessary, with a cotton swab soaked with 70% alcohol.

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Once they have been reprocessed, store the endoscopes freely suspended, dust-free and dry. Store valves separately.

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6.6 Cleaning of accessories

1. Thorough cleaning is required before disinfection or sterilization. Carefully clean the biopsy forceps and brushes whilst they are immersed in a cleaning solution using a soft, disinfected brush. Take special care to remove residues such as blood, excrement and slime from the biopsy forceps and needles – an ultrasound bath is recommended.
2. Rinse in microbiologically clean/sterile water and dry (sterile compressed air can be used).
3. Disinfect accessories (work procedures such as cleaning).
4. Always test instruments after cleaning/disinfection for malfunction and any damage.
5. Before sterilization, carefully oil moveable parts. To do this, only use suitable steam-resistant, physiologically safe, residue-free oils (e.g. 27656 B). After oiling, open and close the forceps several times. For full details on reprocessing, please read the instruction leaflet of the respective accessory.

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6.8 Checking the reprocessing procedure

Before assembly it should be checked that all parts and surfaces are undamaged, clean and dry.

Replace defective parts or send them in for repair. Before sending them in for repair, all parts must be cleaned/disinfected and sterilized or clearly labeled as non-sterile.

If residue is still present after cleaning, reprocess the parts again, subject them to special treatment, discard or repair.

Inspect plastic parts for discoloration, porosity, and flexibility. If the parts differ significantly from when they were new, they must be replaced. This especially applies to seals and rubber caps.

Test the instruments for proper operation (deflection of the distal end, fixing lever). Check the lens for a clear view, the instrument for steerability and all channels for patency (see page 11).

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6.9 Sterilization

Sterilization of instruments to be used on patients is not necessary when coming into contact with the skin and mucous membrane, as well as when introduced through natural, non-sterile orifices.

When used on patients via sterile tissue, a transmucosal access or piercing of the mucosa, sterilization of the instruments is necessary.

CAUTION: All deviations from the recommended sterilization parameters must be validated by the user.

CAUTION: The STERRAD process is not suitable for sterilization of videoscopes. This process damages videoscopes.

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6.9.1 Gas sterilization

Sterilization using ethylene oxide is a gentle and therefore recommended method of sterilization.

NOTE: Observe national laws and regulations for the gas sterilization method.

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CAUTION: If the instrument is to be sterilized, the entire system must be opened by fitting the EO cap on the pressure compensation part to ensure pressure compensation.

WARNING: Before commencing gas sterilization, the instrument must be cleaned, disinfected and completely dried or else successful sterilization will not be achieved.

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- Place the instrument securely on a sterilization tray.

	20:80 EO/CO	10:90 EO/HCFC
Temperature:	55°C	55°C
Relative humidity	50%	50%
Vacuum:	533 mmHg effective	533 mmHg effective
Outlet pressure:	69 kPa 0.70 Kg/cm, l 10 PS	97 kP a0.98 Kg/cm, 14 PSI
EO concentration:	450 mg/l	600 mg/l
Conditioning:	1 hour	1 hour
Gas exposure time	:5 hours	5 hours
Aeration:	12 h at 55°C	12 h at 55°C
An aeration period of 72 hours is required after EO sterilization. An aeration chamber (desorption chamber) can be used to reduce the aeration period to 12 hours. During sterilization, the temperature must not exceed 55°C.		

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For gas sterilization with ethylene oxide, the following airing times must be observed for instruments due to gas absorption of the materials:

	at 20°C	at 42°C	at 55°C
Telescopes	24 hrs	18 hrs	12 hrs
Rubber items	7 days	4 days	3 days
Video- scopes, plastics	5 days	2 days	30 hrs
(Minimum time in simple transparent packaging and an application time on the patient of more than 30 minutes).			
Maximum acceptable levels (as per ISO 10993-7:2008) of residues following EO sterilization are as follows:			
Ethylene oxide:	25 ppm*/250 ppm**		
Ethylene chlorohydrine:	25 ppm*/250 ppm**		
Ethylene glycol:	250 ppm*/5000 ppm**		

** for instruments in contact with skin and mucous membrane,
* for instruments in contact with blood.

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WARNING: The efficiency of the sterilization should be checked to see that it complies with the German Medical Devices Act (MPG).

CAUTION: Attention must be paid to the notes on cleaning, sterilization and storage in the 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' instructions. The cleaning, disinfection, and sterilization procedures are explained there in detail.

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6.10 Storage

CAUTION: Do not expose videoscopes to direct sunlight or in excessive heat.

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The videoscopes should ideally be stored hanging vertically, with the distal end and supply connector facing down, in ventilated cupboards designed for this purpose.
Store the disassembled valves separately.

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8 Maintenance

Preventive maintenance is not essential. However, regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Information regarding maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, we recommend a functional or safety test of the device at least once a year.

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9 Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments, and equipment must be cleaned, disinfected, and, if necessary, sterilized before they are sent for repair. We reserve the right to return contaminated instruments, equipment and endoscopes to the sender.

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Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no guarantee on the correct functioning of equipment, instruments or endoscopes which have been repaired by unauthorized third parties.

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10 Servicing and repair

Defective endoscopes must be repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts only.

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11 Responsibility

As the supplier of this endoscope, KARL STORZ shall be liable for the safety, operational reliability and performance of this instrument subject to the following conditions:

assembly, system expansion, readjustment, modification, or repair have been performed by persons authorized by KARL STORZ and that the instrument is used in accordance with the instruction manual.

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12 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see Chapter 'Subsidiaries'), even during the warranty period.

If the device is opened, repaired or modified by unauthorized persons, this will relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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13 Directive compliance

According to the Medical Device Directive (MDD): Medical product of Class IIa

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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
Insertion sheath
☐ cracked
☐ worn out

Distal tip

☐ Lens scraped
☐ Insufflation disturbed
☐ Image transmission error

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
KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com
KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za
TOO KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BC „Arman“, office314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz
KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai – United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com
KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in
KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com
KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au
KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk
KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 610041, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg
KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg
KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com
KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw
KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: , Fax: +63 2 31745-00
E-Mail: philippines@karlstorz.com
KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

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**Important information
for users of KARL STORZ
instruments**

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, these products are the result of years of experience and great care in manufacturing. You and your organization have decided in favor of modern, high-quality items of equipment from KARL STORZ.

This manual is intended to serve as an aid in the proper use, cleaning, disinfection and sterilization of video endoscopes manufactured by KARL STORZ. The video endoscopes are sensitive precision instruments, which require particular care during preparation. To prevent damage during preparation, guidance is provided here for the disassembly, cleaning, sterilization and assembly of KARL STORZ video endoscopes. For this reason, read these instructions carefully and keep them as a ready reference.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply, design, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, illustrations and descriptions in this manual.

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



- ① Deflection up/down
- ② Up/down brake
- ③ Programmable remote control buttons
- ④ Suction valve
- ⑤ Irrigation/insufflation valve
- ⑥ Working channel with rubber lip valve/biopsy
valve
- ⑦ Deflection right/left
- ⑧ Left/right brake
- ⑨ Working channel
- ⑩ Nozzle for lens cleaning/irrigation/insufflation
- ⑪ Fiberoptic light carrier
- ⑫ Lens
- ⑬ Water jet channel (only for 13806 NKS/PKS,
13807 NKS/PKS, 13910 NKS/PKS and
13911 NKS/PKS)
- ⑭ Supply tube
- ⑮ Connection for leakage test
- ⑯ LUER lock (connector for jet nozzle, only for
13806 NKS/PKS, 13807 NKS/PKS,
13910 NKS/PKS and 13911 NKS/PKS)
- ⑰ Connection for water bottle
- ⑱ Light source/GASTRO PACK™ light cable
connection
- ⑲ Port for video processor/electronic connector
- ㉑ Protection cap for video processor
connection port
- ㉑ Suction/aspiration port
- ㉒ Air pin

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**Important information for users of
KARL STORZ instruments II**

Pictures of the instruments.....	III
Controls, connectors and their uses.....	V

Safety instructions 3

Intended use	3
User qualification	4
Safety precautions when using the video endoscopes	4
Warnings and cautions.....	5

Handling 10

1.0 Scope of supply	10
Compatible video processor units/light sources.....	11
2.0 Operating the video endoscope.....	13
2.1 Preparation of the light source.....	13
2.2 Preparation of the processor unit.....	13
2.3 Connecting the light conductor, air/water connections and the suction channel.....	13
2.4 Adjusting the angle of the distal tip	14
2.5 Suction	14
2.6 Insufflation and light.....	14
2.7 Water rinsing the lens	15
3.0 Preparation before use	15
3.1 Testing the various functions	15
3.2 Holding the endoscope correctly	18
3.3 Additional recommendations for use.....	19
3.4 Operating the endoscope buttons	20
3.5 Disassembling and assembling the Albarran unit (only video duodenoscopes).....	21

Preparation and care24

General information	24
General cleaning.....	26
Preparation for cleaning.....	27

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Please read these safety instructions carefully before using the instrument. Make certain you are well acquainted with the function and operation of the instrument before initial use on the patient.

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Intended use

The KARL STORZ video endoscopes (13800 xxx, 13801xxx, 13804xxx, 13806xxx, 13807xxx) have been developed for endoscopic diagnosis and treatment in the upper digestive tract, including the oesophagus, the stomach and the duodenum.

The KARL STORZ video duodenoscopes (13880 xxx – 13882 xxx) were developed for the endoscopic diagnosis and therapy of the biliopankreatic system.

The KARL STORZ video rectoscopes (13912xxx), video sigmoidoscopes (13900/13901 xxx) and video colonoscopes (13904xxx - 13907 xxx and 13910/13911 xxx) have been developed for endoscopic diagnosis and therapy in the lower digestive tract, including the rectum and colon.

In addition, these video endoscopes are used for endoscopic diagnosis and therapy.

Videoscopes are applied components, which should only be used with approved supply units (video processor units, light sources).

Contraindications:

No contraindications are known which directly relate to the product. The video endoscope must not be used if, in the opinion of a qualified physician, such an application would endanger the patient, e.g., due to the patient's general condition, or if the endoscopic method as such is contraindicated.

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User qualification

KARL STORZ video endoscopes may only be used by persons with an appropriate medical qualification and who are familiar with the endoscopic technique. The information given in these instructions only serves to instruct in the correct handling and preparation of the video endoscopes.

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Safety precautions when using the video endoscopes

Video endoscopes must be used according to the medical rules and procedures of endoscopy recognized for this method.

NOTE: The patient leakage current limit of the CF application part may be exceeded when CF application parts (application on heart) are used in combination with other BF application parts.

In this case, the reduction of the patient leakage current should be effected using only CF application parts whenever this is possible.

NOTE: The endoscope does not offer any electric insulation against high frequency (HF) currents. Protection against HF currents should be ensured through the use of correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.

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Warnings and cautions

Please read this manual and follow its instructions carefully, to avoid putting your patients, personnel or yourself at risk. The words, **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the instrument. To make the words **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

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WARNING: Check the video endoscope prior to each procedure to ensure that it functions correctly and has been properly cleaned and disinfected.

WARNING: Verify the sound condition of the video endoscope, as well as of accessories used in combination, before every application. Damaged video endoscopes or damaged accessories must not be used.

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WARNING: Remove the EO/pressure compensation cap 11025 XE from the connection piece:

- during every endoscopic intervention
- during manual/machine preparation.

WARNING: Do not allow any electrically active components or instruments to come into contact with the video endoscope.

WARNING: The video endoscope must not be used during the discharge of a defibrillator. Remove the video endoscope from the operation area before discharge occurs.

WARNING: Danger of explosion if flammable gases are used in the immediate proximity of the video endoscope. The video endoscope must not be used in the presence of flammable or combustible materials.

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WARNING: High energy light radiated through endoscopes may give rise to high temperatures in front of the light outlet, at the end faces of the light cable and at the tip of the endoscope.

Prolonged direct contact of light-emitting surfaces on tissue can lead to burns, this must therefore be avoided.

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WARNING: Before removal from the patient, the distal tip of the video endoscope must be returned to the neutral, nondeflected position (L and U (demarcations on the wheels of the deflection mechanism) on top of each other and positioned between the valves) to prevent damage to the instrument and/or injury to the patient.

WARNING: The video endoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for video endoscopes that are repaired by unauthorized personnel.

WARNING: Damage to the video endoscope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.

WARNING: Always unplug the video endoscope before carrying out any maintenance and cleaning procedures.

WARNING: Risk of infection due to shipment! Contaminated video endoscopes can cause serious or fatal diseases during transportation or at the recipient's premises. Please cover the carrying case with the protective sheet 13990 SFN before sending contaminated instruments back to KARL STORZ.

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WARNING: If a malfunction of the endoscope should occur during procedure on a patient, stop the procedure immediately. Place the distal tip of the video endoscope in a neutral, nondeflected position and remove the video endoscope slowly and carefully from the patient.

WARNING: The use of HF probes, HF electrodes or other energy-emitting probes is not permitted if explosive gases (for example, explosive/flammable anesthetics) are present.

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WARNING: Before turning on the high voltage/ electrical supply, ensure that there is no contact with metal parts of the endoscope or any other electricity-conducting accessory, and that the suction unit has no contact with conductive liquids.

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CAUTION: Only use accessories (guide probes, cleaning brushes, etc.) that are specifically suited for the intended purpose for use with video endoscopes and which can be employed in line with generally accepted technical standards and safety guidelines.

CAUTION: Only use accessories with flexible ends. Other types of accessories may cause damage to video endoscopes.

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CAUTION: Place the EO/pressure compensation cap 11025 XE on the connection piece prior to:

- gas sterilization
- shipping

CAUTION: Do not expose the video endoscope to direct sunlight or excessive heat, to avoid damage to the device.

CAUTION: Do not steam sterilize (autoclave) the video endoscopes.

CAUTION: The maximum permissible temperature for decontamination, sterilization, transport and storage is 65°C.

CAUTION: Make certain to read the 'Cleaning, sterilization and care of KARL STORZ instruments' manual. Methods of cleaning, disinfection and gas sterilization are explained in detail there.

CAUTION: Do not hit the distal tip of the video endoscope against hard objects, and never use a sharp object to remove debris.

CAUTION: Never lay heavy objects on the video endoscope.

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NOTE: Any damage to the video endoscope resulting from incorrect operation is not covered by the warranty.

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Compatible video processor units/light sources

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NOTE: The KARL STORZ video endoscope should not be used in combination with video processor units and light sources from other manufacturers. The use of non-compatible units can lead to patient injury or unit damage.

It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces, e.g., voltage-carrying elements of other units, at any time.

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**Connecting to a light source illustra-
ted with the KARL STORZ XENON 100**

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KARL STORZ light sources Xenon 100 or Xenon 100 SCB/ Gastropack 2.0 Operating the video endoscope

2.1 Preparation of the light source

We recommend the use of our light source (Xenon 100 or Xenon 100 SCB) with a 100 W xenon lamp and our portable package solution the Gastropack with economical, light-intensive 24 W Solarc® lamp. Please read the full instructions in the instruction manual of the respective light source.

- 1 Install the light source on a solid, even surface.
2. Insert the power cable into the socket provided.

2.2 Preparation of the processor unit

Refer to the relevant instruction manual: TELECAM® SL II 202130 20, GASTRO PACK® 25 0430 11/ 25 0431 11 or Image 1 hub 22200020-1

2.3 Connecting the light conductor, air/ water connections and the suction channel

1. Insert the supply connector into the light source until it audibly clicks into place (Fig. 1).
2. Connect the video connecting cable to the endoscope. Take note of the white markings on the endoscope and video connecting plug. (Fig. 2).
3. Attach the connector from the water bottle tube to the endoscope supply plug (irrigation bottle connector, Fig. 3). The water bottle holder automatically positions itself correctly and cannot be adjusted manually.
IMPORTANT: To fill the water bottle with sterile water, remove the bottle from the light source and unscrew the lid. Only fill the water bottle 2/3 full with sterile water (1/3 is required to cushion air pressure).
4. Connect the suction tube (patient tube) from the suction pump or central suction system to the supply plug connector.

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2.5 Suction

Switch on the suction pump. To aspirate, press the suction valve (with the red ring) on the endoscope (Fig. 5).

Suction can also be used when the biopsy forceps are placed in the biopsy channel and the rubber lip valve is mounted. During the procedure, the rubber lip valve/biopsy valve seals the forceps channel.

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2.6 Insufflation and light

Start the light source, switch on the air pump.

For insufflation, cover the opening on the blue valve with a finger (do not depress the valve, since this starts water irrigation) (Fig. 6).

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2.7 Water rinsing the lens

Pressing on the air/water valve starts the water flow and/or lens cleaning. Different valve positions allow various combinations (Fig. 7,8).

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3.0 Preparation before use

WARNING: Before being used for the first time, the device must undergo a complete preparation cycle.

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3.1 Testing the various functions

1. Prepare light sources as described in 2.1.
2. Connect supply connector, water bottle and suction tube according to 2.3.
3. Check the supply cable for damage (e.g., breaks, cracks, twisted or squashed areas...) and the supply connectors for loose connections.
4. Examine the surface of the flexible insertion sheath by eye and by touch for cracks, cuts, bulges, foreign bodies or other damage. Do not press too hard on the deflection section.

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- 5a. Testing deflection – Upward/downward control
The large deflection wheel controls up and down movements. To test up and down movements, release the up and down brake. Slowly rotate the large deflection wheel clockwise, the distal end should move downwards. Rotating the deflection wheel counterclockwise causes upward movement.

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- 5b. Testing deflection – Right/left control
The smaller/outer deflection wheel controls right/left movement. To test right/left movement, release the right/left brake. Slowly rotate the small deflection wheel clockwise, the tip of the endoscope should move to the right. Rotating the deflection wheel counterclockwise causes movement to the left.
- 5c. Testing deflection – brake function
The deflection wheels should turn freely when the brakes are released. Release the brakes and test their action (see 5a/b). Activate the brakes by moving the brake lever counterclockwise – ensure that the angled position of the distal endpiece is fixed in the desired position.
Do not twist or bend the deflection piece by hand.

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WARNING: If the deflection and/or brake mechanisms do not function properly, do not use the instrument. Contact KARL STORZ service department.

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6. Check the air supply. Switch on the air pump and hold the distal end of the instrument under sterile water. Cover the air control opening in the air valve (blue ring) with a finger; bubbles should emerge from the air nozzle at the distal end.

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8. Test the suction function. With the unit correctly installed and the suction pump switched on, immerse the distal end in a container of sterile water once again. Depress the suction valve (red ring) fully. Water should be suctioned through the suction channel of the instrument via the connected suction pump and into the collecting bottle of the suction equipment.

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CAUTION: Accessories such as biopsy forceps that are kinked or bent will not function properly, are difficult to guide through the channel and may cause severe damage to the channel. They must therefore be replaced and not used further!

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Holding the endoscope correctly

With the control handle of the instrument in the left hand, the thumb, index and middle fingers can easily control the large deflection wheel. In this position, the suction valve can be depressed with either the index or the ring finger. The air/water valve can also be covered with these fingers.

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3.3 Additional recommendations for use

1. Use a lubricant from the KARL STORZ list of approved and recommended lubricants on our website: www.karlstorz.com under 'Hygiene' for the distal end of the insertion tube.
2. Observe the lumen the whole time. If you become disoriented, withdraw the instrument gently until vision is restored.

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3.4 Operating the endoscope buttons

Both buttons have two functions, which are initiated through short and long presses of the key respectively. Pressing and holding down both buttons at the same time for approx. 3s opens a settings menu and allows you to select the button assignment freely (i.e., Freeze/Zoom or Print (ACC1)/contrast). For more detailed information on button assignment please consult the camera control unit manual TELECAM® SL II, GASTRO PACK™ or Image 1 hub.

The freeze mode can be canceled by a short press on the freeze button.

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NOTE: In case of 'Red-out', a blurred red picture on the screen, caused by the distal endoscope tip resting directly on the mucosa, withdraw the instrument gently until vision is restored.

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NOTE: If the angling of the tip can no longer be properly controlled the examination must be ended immediately. Release the brakes and place the tip in the neutral position. Position the deflection wheels so that 'U' and 'L' are between the valves (see page 13/15), then carefully withdraw the instrument back, whilst continuing to monitor the procedure.

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3.5 Disassembly and assembly of the Albarran unit (video duodenoscopes only)

As conventional duodenoscopes are difficult to clean and maintain due to their design, particularly in the region of the Albarran unit, KARL STORZ has developed a modular Albarran system.

The patented Albarran module is easy to disassemble, which facilitates thorough cleaning and disinfection. In addition, the Albarran module can be autoclaved. It is no longer necessary to disassemble the entire endoscope to eliminate technical malfunctions in the Albarran unit. As a rule, in this case, repairs can be carried out by disassembling the Albarran unit.

Disassembling the Albarran unit

1. Remove the cable as shown in Fig. 1 by turning the rotary knob on the proximal endoscopic housing. Unscrew the rotary knob (2) from the proximal housing counterclockwise.
2. Remove the collet chuck (11), which is now visible (Fig. 2).
3. Unlock the Albarran module from the distal endoscope head (Fig. 3).
The quick locking screw (6) of the Albarran module unit is indicated by two red dots. The red dots are turned counterclockwise towards each other using a special assembly/ disassembly key (7) Art. no. 13991 B. The Albarran module is now unlocked (Fig. 4).
4. Removal of the Albarran module with Albarran wire (Fig. 5).
Separate the Albarran module (3) from the endoscope with your thumb and carefully remove the wire (4) from the endoscope.

NOTE: Do not bend or kink the wire!

5. Removal of the Albarran wire sheath (Fig. 6)
Angle the distal tip of the endoscope by turning the small deflection wheel clockwise ('R').

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
3. Attach the Albarran module (Fig. 8) by turning the quick locking screw (6) clockwise with the assembly/disassembly key (7).

The two red dots must now be opposite one another (Fig. 9).

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
4. Testing the module locking:
Locking of the Albarran module must be tested by trying to press it off the head of the endoscope with your thumb.

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
5. Place the shaft of the endoscope in a straight position on a table.
6. Lock the control slide (12) into the endoscope housing as per Fig. 11. At the same time move the Albarran control lever (13) backward as far as it will go.

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
7. Introduce the collet chuck (11) into the proximal housing with a twisting action until it engages and cannot be turned any further (Fig. 11).

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8. Screw on and tighten the rotary knob (2) as seen in Fig. 12.
9. Testing Albarran function: Grip the Albarran lever at the distal end with your thumb and manipulate the proximal Albarran control 13 simultaneously (Fig. 12).

You should be able to feel a significant movement pressure of the Albarran lever with your thumb. If so, the Albarran unit is functioning correctly.

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NOTE: This documentation covers only the handling and care of the Albarran control lever.

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Cleaning and disinfecting the Albarran module, see page 40.

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General information

WARNING: The stated cleaning, disinfection and sterilization methods alone do not guarantee successful sterilization of the endoscopes. This can only be achieved if personnel follow recognized and validated preparation methods.

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WARNING: Remember that sterilization is only then successful if

- the instruments have been freed of all organic material and chemical residues.
- work proceeds according to recognized and validated methods.
- any deviation from recommended procedures must be validated and tested by the user.
- the recommended sterilization parameters are maintained. These are only valid with cleaning, disinfection and sterilization equipment that is properly maintained and calibrated.

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CAUTION: The manual 'Cleaning, sterilization and care of KARL STORZ Instruments' must be followed. There, the cleaning, disinfection and sterilization procedures are explained in detail. The list of KARL STORZ approved materials for the preparation of video endoscopes can be found on the Internet at 'www.karlstorz.com' or 'www.karlstorz.de' under 'Hygiene' and also in the appendix of this manual.

CAUTION: When preparing and using the solutions, follow the manufacturer's instructions paying close attention to proper dilution and immersion times. Disregarding the instructions can lead to damage. Bear in mind the microbiological range of action of the chemicals used.

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CAUTION: Never clean video endoscopes with metal brushes.

CAUTION: Constantly changing between different preparation methods causes particular wear on the materials and must therefore be avoided. A once-only change in methods, however, e.g., when acquiring new cleaning equipment, should be carried out in consultation with the manufacturers of the endoscope, the chemicals and the washer/ disinfectant. The method change is then, as a rule, completely harmless.

CAUTION: Do not clean video endoscopes in an ultrasound bath.

CAUTION: Do not steam sterilize (autoclave) video endoscopes.

CAUTION: Never immerse video endoscopes in normal saline solution, since even brief contact can lead to pitting and corrosion.

CAUTION: For manual preparation, video endoscopes may only be rinsed with microbiologically clean/sterile water.

CAUTION: Decontamination in washers and disinfectors is only possible using valid procedures.

CAUTION: Special attention should be given to the following sterilization conditions:
Temperature: 5°C – 65°C
Humidity: 15% – 100%

CAUTION: Even parts that do not come into direct contact with the patient (e.g., video connecting cables) should be thoroughly disinfected following an intervention.

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General cleaning

Cleaning should only be performed by staff who are appropriately trained. All safety regulations should be observed (gloves, protective glasses, mouth and nose protection, clothing, TBA 250...).

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The basic guidelines and recommendations for the preparation of instruments (for example RKI, ESGE/ESGENA) should be observed and the effectiveness verified.

CAUTION: Before immersion in cleaning solution, check the instrument for leaks. Otherwise fluid penetrating the video endoscope could cause damage.

CAUTION: Do not autoclave the endoscope or accessories (with the exception of accessories that are specifically designated as autoclavable, e.g., forceps or Albarran module).

NOTE: The endoscope must be cleaned immediately after every use. If biological material is allowed to remain and dry in the instrument, complete or partial blockage of the channel may result. Thorough cleaning is a prerequisite for effective disinfection and sterilization. This will avoid clogs/ blockages and the resulting repairs and instrument failures. For reasons of material compatibility, only cleaning and disinfection chemicals approved by KARL STORZ should be used (see current list in the annex, p. 52 or on our website www.karlstorz.de under 'Hygiene').

NOTE: The endoscope must undergo a complete preparation cycle prior to the initial use.

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Preparation for cleaning

1. The preliminary cleaning of the endoscope is performed at the place of operation. The main cleaning and disinfection is carried out in a separate preparation room.
Immediately after use, remove any heavy soiling (secretions, etc.) and organic material from the outer surfaces of the insertion tube with a disposable cloth moistened with active cleaning solution. Then submerge the distal end in a container with active cleaning solution and activate the suction to rinse the working channel.
2. Put the cleaning valve into the air/water housing. (Do not insert this valve during normal use!)
3. Check that the water bottle is filled (max 2/3). Start the air supply (insufflation pump). Water and air are now passed through both the air channel and the water channel. (The process should run for at least 10 seconds).
4. Switch off video system.
5. Pull out the video plug from the video connection cable at the supply connector of the endoscope when the video processor unit is switched off. Put on the sealing cap for the video cable connector and screw it on tightly.
6. Remove the endoscope from the video processor unit and light source. Place the cleaning valve in an active cleaning solution.

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NOTE: Ensure the video cable connector sealing cap is secured tightly (blue rubber inlet seal should be slightly moistened with instrument oil from time to time).

2. The connection tube of the leakage tester and the ventilation connector on the endoscope must be dry. Connect the leakage tester connector to the supply connector and lock by turning clockwise 90°.

NOTE: The rubber seal should also be treated occasionally with instrument oil.

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NOTE: During inflation, make sure that the red pressure relief button is completely closed.

3. Inflate the leakage tester until the maximum test pressure of 160 - 180 mmHg is reached. The deflection rubber on the distal deflection part of the endoscope expands slightly. Move the deflection part in all four directions as far as possible and check the manometer to see whether the pressure remains constant or falls.

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4. Wait five minutes. If the pressure does not fall by more than 10 mmHg, submerge the endoscope in the cleaning solution. Place instrument with closed leakage tester in the cleaning solution and repeat the process of utilizing the deflection wheels to deflect the tip of the scope fully in all four directions. If this test is also negative and no air bubbles rise in the liquid, the leakage tester can be removed. To do this, take the supply connector out of the liquid and press down on the red button of the leakage tester to vent the endoscope.

NOTE: Never release the pressure in the liquid. To this end, remove the plug from the liquid and actuate the red pressure relief lever, before using the leakage tester again reset the red lever to the initial position. The test pressure must fall to 0. Remove the leakage tester.

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If a pressure loss occurs, this means that the video endoscope leaks and may no longer be used or immersed or prepared by machine.

In this event, fill out the accompanying repairs form (see Appendix), roughly clean and remove any organic soiling/residue from the endoscope, mark it as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the enclosed transport protection sheet 13990 SFN.

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WARNING: Risk of infection due to shipment. Contaminated instruments can cause aerogenous or contact infections during transportation or at the recipient's premises.

Shipment without the plastic film packaging must be avoided.

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

Kanalsystem eines KS Standard Video-
Endoskopes

Channel system of a KS standard video
endoscope

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
Channel system of a KS therapy video endoscope with additional irrigation channel (water jet channel) and 2 working channels

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Cleaning/disinfection of the video endoscope: manual cleaning

Please adhere to the general standards, e.g., ESGENA guidelines or DEGEA/RKI recommendations for cleaning, disinfection and sterilization, as well as to the instructions of the chemical manufacturer.

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1. Remove the rubber lip valve/biopsy valve, suction valve and air/water valve and place in an active cleaning solution and clean using a soft disinfected brush (13991 VB).

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2. Wipe the housing, insertion shaft, and supply tube carefully with a lint-free cloth soaked in an enzymatic cleaner (use a suitable cleaning agent or disinfectant in accordance with the DGHM list and the list of approved chemicals, see annex). Avoid shocks and pulling movements so as not to damage the insertion shaft. Rinse the endoscope carefully with water of at least drinking quality (microbiologically sterile water).

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3. Immerse the endoscope in the active cleaning solution and flush the channels using a syringe and the supplied irrigation tubes (13991 SS, 13991 YS). The channels must be filled without any bubbles, to guarantee a complete moistening of the surfaces.

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

CAUTION: Only immerse instrument parts once they have successfully passed a leakage test.

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4. When submerged, use the brush supplied (Art. No. 27651 H/G) (see scope of supply page 10) to brush below the surface of the liquid in one direction starting from the supply plug (control body --> biopsy channel --> distal end).

! CAUTION: When brushing from the suction valve to the biopsy channel the brush comes against the biopsy channel (see channel system picture) and cannot be pushed any further into the biopsy channel!

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5. Now guide the brush through the biopsy channel to the distal end of the shaft.

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Before retracting the brush head that emerges at the distal end, clean it using another brush. Do not brush backward and forward, rather brush in one direction only.

! **CAUTION:** Cleaning with brushes must be carried out with the instrument immersed.

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
6 Thoroughly clean out the entire suction channel (sheath part and supply part) and the biopsy channel a number of times with a brush. Observe the exposure time for the active cleaning solution.

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7. Connect the Y rinsing tube (13991 YS and 13991 SS) and using a syringe, flush fresh cleaning fluid through the channels using a syringe. Afterwards, rinse with microbiologically clean/sterile water.

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8. Suitable brushes

Art. no.	Designation	Used for
27651 F	Cleaning brush dia. 0.5 mm	Working channel, Albarran channel
27651 G	Cleaning brush dia. 3 mm	Working channel
27651 H	Cleaning brush dia. 5 mm	
13991 VB	Valve brush dia. 12 mm	Valve seat

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9. Carefully clean all particles of debris from the distal end. It may be necessary to clean the air/water valve with a brush (13991 VB) with soft bristles, in order to remove hardened contaminants. Never remove debris with a needle or other sharp objects, because these can cause damage.

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10. Remove the endoscope from the cleaning solution and rinse with microbiologically pure/sterile water.

11. Carefully clean all valves with the disinfected cleaning brush (13991 VB) and then rinse. The valves should be cleaned in an ultrasound bath.

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Once the cleaning work has been completed, clean the cleaning instruments, for example the brushes, of residue and disinfect/sterilize them. Never reuse damaged or kinked brushes.

CAUTION: Do not autoclave the endoscope or its accessories (with the exception of accessories that are specifically designated as autoclavable, e.g., forceps, cleaning brushes, albarran module, etc.).

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
Manual disinfection

After careful rinsing, fully immerse the endoscope in the disinfecting solution.

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! CAUTION: When preparing and using the solutions follow the chemical manufacturers' instructions paying close attention to proper concentration and exposure times. Prolonged immersion can lead to damage. Bear in mind the microbiologic range of action of the chemicals used.

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Using the irrigation tubes and a syringe with LUER connector, rinse out the channels a number of times with disinfecting solution (see illustration).

Between rinses, allow the specified time (chemical manufacturers' instructions) for the disinfecting solution to take effect in the channel system.

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Then rinse the channels and surfaces of the video gastroscope with microbiologically pure/clean water.

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12. After the final step in the cleaning process has been completed, dry the entire endoscope. Thoroughly dry off all external surfaces of the endoscope, preferably with sterile compressed air at a max. of 0.5 bar/7.25 psi, thereby removing any residues of cleaning solution.

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Notes on cleaning special endoscopes

The picture opposite shows the disassembled Albarran module. The two upper parts should be cleaned with brushes (13991 VB) and then in an ultrasound bath with cleaning solution.

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Cleaning of the Albarran module

The following components can be cleaned in the ultrasound bath:

- Albarran module with cable
- Albarran cable guide
- Collet chuck
- Rotary knob

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
To enable better rinsing of the Albarran module's wire sheath (Albarran spring), the irrigation tube 13991 PS can be connected to it.

The lower part of the endoscope (in the above picture) is cleaned as described previously.

NOTE: The Albarran module of the video endoscopes can also be steam sterilized at 134 °C ± 3 °C.

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Starting from the proximal housing area, brush the working channel using the cleaning brush (Art. no. 27651 H/G) and cleaning agent.

Brush the disassembled Albarran channel with cleaning brush 27651 F supplied with the instrument.

Before retracting the brush head that emerges at the distal end, clean it using another brush or compressed air.

After the brush cleaning, the channels must be thoroughly rinsed with germfree, completely demineralized or microbiologically pure/sterile water.

Then disinfect the endoscope. Manual disinfection: see p.38

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Then blow the channels dry with compressed air.

CAUTION: The pressure must not exceed 50 kPa (0.5 bar/7.25 psi). The cleaning gun might therefore have to be operated with a pressure reducer.

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Preparation of the video endoscope by machine

All flexible KARL STORZ endoscopes can be connected to approved washer/disinfectors for semi- or fully automatic cleaning and disinfection. Leakage testing, manual cleaning and (where necessary) rinsing must be carried out beforehand. The manufacturer of the cleaning machine must be consulted when selecting the preparation programs. Usually, special programs for flexible endoscopes are available which are suited for decontamination and compatible with the materials. For machine preparation, make sure that the flexible endoscopes are not heated above 65°C.

! **CAUTION:** Check the videoscope for leaks before every disinfection in semi- or fully automatic equipment.

! **CAUTION:** Ensure that only KARL STORZ approved agents are used for procedures involving peracetic acid.

! **CAUTION:** The STERRAD process is not suitable for sterilization of video endoscopes. This process destroys video endoscopes.

Lay the endoscope carefully in the tray or basket of the machine. Connect the irrigation tubes supplied with the cleaning adaptor set 13991 CK to the channels. Proceed with the endoscope cleaning program according to the recommendations of the machine manufacturer.

Cleaning adaptor set 13991 CK

The cleaning adaptor set 13991CK is used for validated preparation of KARL STORZ video endoscopes from series 13xxx NKS/PKS (acc. to ISO 15883 part 4). The cleaning adaptor set is not intended for any other purpose as this may damage it.

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- 1 Connection tube 13991 CL for leakage tester port ⑮ of the videoscope
- 2 Connector (13991CL) for vent port
- 3 LUER connector (13991CL) for washers/disinfectors for endoscopes

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- 4 Cleaning adaptor 13991CR for air/water (A/W) channel ⑩ of the videoscope
- 5 Connector (13991CL) for A/W channel (air pin ②)
- 6 Connector (13991CL) for water channel (water bottle connection ⑦)
- 7 2x LUER connectors (13991CR) for washers/disinfectors (1x LUER for air channel, 1 LUER for water channel)

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- 8 Cleaning adaptor 13991CS for suction channel (aspiration) ⑪ of the videoscope
- 10 Connector (13991CL) for suction channel at the supply connector ⑪
- LUER connector (13991CS) for washers/disinfectors

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- 14 Channel separator 13991CD for suction valve ④ and irrigation/insufflation valve (A/W) ⑤ of the videoscope
- 15 Valve pin (13991CL) for irrigation/insufflation valve (A/W) ⑤
- 16 Closure (13991CL) for suction valve ④

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WARNING: This product is delivered unsterilized and must therefore be cleaned, disinfected or, if necessary, sterilized prior to initial use and each subsequent reuse.

WARNING: Risk of infection by improper preparation: Preparation may only be carried out by experienced and trained staff.

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- Visual checks of the connection tubes (patency) and the O-ring (wear)
- check the connections are correctly and securely in place
- check for damage

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Assembly

Plug the connector (2) of the connection tube 13991 CL into the vent port (15) of the videoscope at the supply connector and secure by turning clockwise.

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Plug the cleaning adaptor 13991 CR (5) as far as it can go into the connector (22) of the air pin of the videoscope attachment.

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Then plug connector (6) as far as it will go into the 'Connection for water bottle' (17) of the videoscope.

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Plug the cleaning adaptor 13991CS with connector (9) as far as it will go into the suction channel (21) at the supply connector of the videoscope.

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If necessary, remove the suction and irrigation/insufflation valves (④, ⑤) from the videoscope.

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Insert the channel separator 13991 CD ⑭ with the valve pin ⑮ into the free opening of the irrigation/insufflation valve ④ as far as it will go so that the lock ⑮ lies on top of the free opening of the suction valve ④.

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Press the slide mechanism of the channel separator 13991CD ⑭ downward and slide it forward until the channel separator 13991 CD ⑭ audibly locks into place. Check that the channel separator is securely in place.

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NOTE: If necessary, remove the cap of the biopsy/working channel ⑥ of the videoscope.


Connect the connection adaptor 13991 SS with the connector ⑫ to the biopsy/working channel ⑥ of the videoscope.

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Connect all open LUER connectors of the attached connection adaptors to the respective connections on the washer/disinfector.

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! **CAUTION:** When the program has finished, the conditions in the preparations machine present the videoscopes with a risk of corrosion.

! **CAUTION:** Immediately after the end of the program cycle, the video endoscopes are hot as they are heated up during chemical thermal preparation. Consequently the instrument's material properties are different and, as such, the prepared endoscope should be allowed to cool down for at least 20 minutes at room temperature before being used again.

! **CAUTION:** Ideally, remove the video endoscopes from the machine immediately after the program has finished. Then dry off the optical surfaces, if necessary, with a cotton swab soaked with 70% alcohol.

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Store the endoscopes freely suspended, dust-free and dry. Store valves separately.

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Trocknung

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Recommended procedure for the preparation of video endoscopes

You also will find an illustrated page of the preparation procedures in the middle of this manual. There is also a poster of the preparation process available from your local sales representative or it can be requested per e-mail directly from KARL STORZ (info@karlstorz.de)

Preliminary cleaning directly after the examination (see page 27)

- The outer surfaces should be wiped down from top to bottom with a lint-free disposable cloth that is soaked with active cleaning solution.
- Suction in ample amounts of cleaning solution through the distal end.
- Apply the air-water supply for a few seconds every cycle
- Connect the cleaning adapter
- Remove the flexible endoscope from the light source

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Cleaning in the preparation place

- Carry out leakage test
- Immerse the flexible endoscope
- Clean the channels with special cleaning brushes
- Clean the complete air/water system thoroughly (in the case of duodenoscopes, the Albarran unit must also be disassembled, see page 21)
- Rinse the channels

Disinfection

- Immerse the endoscope and rinse channels

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Preparation in the washer/disinfector

- The machine must fulfill all the cleaning, disinfection and material compatibility requirements (validated process).
- Make sure that the environmental conditions with regard to temperature and pressure are not exceeded: +60°C (tolerance +/- 5%) and 1 bar for the channels, 0.3 bar leakage test.
- Recommendation: Before placing the endoscope in the machine, clean all the channels (working channels, water channels, Albarran channels, air/water channels) whether or not they have been used.
- We recommend sticking to one type of cleaning agent.

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Cleaning of accessories

1. Thorough cleaning is required before disinfection or sterilization. Carefully clean the biopsy forceps and brushes under the surface of a cleaning solution using a soft, disinfected brush. Take special care to remove residues such as blood, excrement and mucous from the biopsy forceps and needles – an ultrasound bath is recommended.
2. Rinse in microbiologically pure/sterile water and dry (using sterile compressed air).
3. Disinfect accessories (work procedures such as cleaning).
4. Always test instruments after cleaning/disinfection for proper functioning and any damage.
5. Before sterilization, carefully oil moveable parts. Use only suitable, steam-resistant, physiologically safe, residue-free oil (e.g., 27656 B). After oiling, open and close the forceps several times. For details on preparation, please read the instruction leaflet of the respective accessory.

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Preparation of water bottle 13992 BS

The water bottle may be cleaned/disinfected manually or mechanically.

1. Empty the bottle
2. Carefully clean the container and lid with an active cleaning solution. Immerse the container and lid in the solution.
3. Attach the water bottle lid connector to the irrigation adapter 13991SW and rinse out with a cleaning liquid. Immerse in disinfectant solution (observe the chemical manufacturer's instructions on exposure time and concentration) making sure that all air bubbles have been removed. On removal rinse off thoroughly using microbiologically pure water.
- By machine: Place and secure in the washer and disinfectant. Choose between either chemical disinfection at 65°C or thermal disinfection at 93°C.
4. Sterilize. Steam sterilization is recommended.

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Check preparation procedure

Before assembly it should be checked that all parts and surfaces are undamaged, clean and dry. Replace defective parts or send them in for repair. Before sending them in for repair, all contaminated parts must be cleaned/disinfected and sterilized or clearly labeled as unsterilized.

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Check plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.

Test the instruments for proper operation (deflection of the distal end, brake lever). Check the lens for a clear view, the instrument for appropriate deflection and all channels for patency (see page 11).

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Gas sterilization

Sterilization methods using ethylene oxide is mild and therefore the recommended method of sterilization.

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CAUTION: If the instrument is to be sterilized, the entire system must be ventilated by fitting the EO cap on the pressure compensation part to ensure pressure compensation.

WARNING: Before commencing gas sterilization the instrument must be completely dry or else sterile conditions will not be achieved.

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	20:80 EO/ CO	10:90 EO/ HCFC
Temperature:	55°C	55°C
Relative humidity:	50%	50%
Vacuum:	533 mmHg effective	533 mmHg effective
Outlet-pressure:	69 kPa	97 kPa
	0.70 Kg/cm, 10 psi	0.98 Kg/cm, 14 psi
EO-Concentration:	450 mg/l	600 mg/l
Conditioning:	1 hour	1 hour
Gas exposure time:	5 hours	5 hours
Aeration:	12 h at 55°C	12 h at 55°C

An aeration period of 72 hours is required after EO sterilization. An aeration chamber (desorption chamber) can be used to reduce the aeration period to 12 hours. During sterilization the temperature must not exceed 55°C.

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	at 20°C	at 42°C	at 55°C
Telescopes	24 hrs	18 hrs	12 hrs
Rubber items	7 days	4 days	3 days
Video endoscopes			
Plastics	5 days	2 days	30 hrs

(Minimum time in simple transparent packaging and an application time on the patient of more than 30 minutes).

Maximum acceptable levels (per 21 CFR 812.100; US Federal Register, Vol. 43, No. 122, June 23 1978) of residues following EO sterilization are as follows:

Ethylene oxide:	25 ppm*/250 ppm**
Ethylene chlorohydrine:	25 ppm*/250 ppm**
Ethylene glycol:	250 ppm*/5000 ppm**

** for instruments in contact with skin and mucous

membrane,

* for instruments in contact with blood.

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CAUTION: Attention must be paid to the notes on cleaning, sterilization and storage in the 'Cleaning, sterilization and care of KARL STORZ instruments' manual. The cleaning, disinfection and sterilization procedures are explained therein in detail.

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Storage

CAUTION: Do not expose video endoscopes to direct sunlight or store in an environment with excessive heat.

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The videoscopes should be stored hanging vertically, with the distal end and supply connector facing down, in well ventilated cupboards designed for this purpose.
Store dismantled valves separately.

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Maintenance

Preventive maintenance is not essential. Regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

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Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and, if necessary, disinfected/sterilized before they are sent for repair. We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no warranty on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts.

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Limitation of liability

As a supplier of this instrument, KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions: All assembly, operations, system expansions, readjustments, modifications, or repairs to same have been performed by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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

- ☐ cracked
- ☐ worn out

Distal tip

- ☐ Lens scraped
- ☐ Insufflation disturbed
- ☐ Image transmission error

- ☐ Deflection
- ☐ leaking
- ☐ blocked / damaged

- ☐ Supply plug
- ☐ electrical contacts damaged / corrosive
- ☐ leaking
- ☐ Light transmission
- ☐ damaged
- ☐ Illumination power not adequate

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KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

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KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenheng District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

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1. Photographs of the unit.....	6
2. Controls	8
3. Intended Use	11
3.1 Safety precautions when using the fiberscopes	11
4. User qualification.....	12
5. Warnings and cautions	13
6. Unpacking	16
6.1 Basic equipment	16
6.1.1 Accessories	17
7. Initial operation	18
7.1 Fitting the suction valve.....	18
7.2 Fitting the tube holder (only intubation fiberscopes).....	18
7.3 Fitting the lip valve (only bronchofiberscopes)	18
8. Handling the fiberscopes.....	19
9. Using the fiberscopes	20
10. Preparation of the fiberscope	24
10.1 Connect these supply lines	24
10.2 Controlling the distal tip.....	25
10.3 Suction	26
10.4 Insufflation	26
10.5 Focusing	26
10.6 Instrument channel	27
10.7 PDD version.....	27
11. DCI® fiberscope with camera	28
11.1 Assembly	28
11.2 Disassembly	29
12. DIN eyepiece on DCI® fiberscope ...	29
12.1 Assembly	29
12.2 Disassembly	29
13. DCI® camera on DIN eyepiece.....	30
13.1 Assembly	30
13.2 Disassembly.....	30
14. Reprocessing	31
14.1 Accessories	33
14.2 Precleaning	33

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14.3	Leakage test	34
14.3.1	Phase 1: Performing the leakage test – dry	34
14.3.2	Phase 2: Performing the leakage test – in liquid	35
14.4	Disassembly	36
14.5	Manual reprocessing	36
14.5.1	Manual cleaning	36
14.5.2	Manual disinfection	37
14.5.3	Manual High-Level Disinfection (HLD)	38
14.5.4	Drying	38
14.6	Machine reprocessing	39
14.6.1	Manual preparation	39
14.6.2	Connecting	40
14.6.3	Machine cleaning/ chemical-thermal disinfection	40
14.7	Assembly, inspection and care	41
14.8	Packaging	42
14.9	Sterilization	43
14.9.1	Hydrogen peroxide (H ₂ O ₂) sterilization – ASP STERRAD®	43
14.9.2	Hydrogen peroxide (H ₂ O ₂) sterilization – STERIS® V-PRO®	44
14.9.3	Ethylene oxide sterilization (EO)	45
14.9.4	Chemical low-temperature sterilization with peracetic acid – STERIS System 1E®	45
14.10	Limits of reprocessing	46
14.11	Storage	46
14.12	Reprocessing of accessories	47
14.13	Overview: Cleaning brushes	48
15	Standard compliance	49
15.1	Servicing and repair	49
15.2	Responsibility	49
15.3	Warranty	49
15.4	Directive compliance	50
16	Repair and replacement program	50
16.1	Important information	51

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, these products are the result of years of experience and great care in manufacture. You and your organization have decided in favor of modern, high-quality items of equipment from KARL STORZ.

This manual is intended to serve as an aid in the proper use, cleaning, disinfection and sterilization of fiberscopes manufactured by KARL STORZ. The fiberscopes are sensitive precision instruments


which require particular care during reprocessing. To prevent damage during reprocessing, support is provided here for disassembly, cleaning, sterilization and assembly of the KARL STORZ fiberscopes. For this reason, read these instructions carefully and keep them as a ready reference.

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DCI®-Fiberskope (Direct-Coupled-Interface)

This instruction manual describes three fiberoscope versions, distinguished from one another by the imaging system. If a section refers exclusively to one of these versions, it is thus appropriately labeled.

Telescope

Telescope with special filters for PDD/AF

DCI® fiberscopes (Direct Coupled Interface)

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


KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope for supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

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2. Controls

Connectors, accessories and their uses

495 F	Light cable adaptor for Wolf
495 G	Light cable adaptor for KARL STORZ, Olympus, Winter & Ibe
11025 E	Pressure compensation cap (red)
11301 CA	Lip valve, disposable version, pack of 20
13242 XA	Adaptor for leakage test in cleaning and disinfection machines (green)
29100	Plug for instrument channel, pack of 10, reusable
11301 CB	Suction valve insert, reusable version
5917800	Seal ring (suction valve insert)
6885391	Cap (suction valve insert)
5923110	Spring (suction valve insert)
6927791	Housing (suction valve insert)
8567290	Plunger (suction valve insert)
11301 CE	Suction valve insert, disposable version, pack of 20
11301 D1	LED battery-powered light source with fine screw thread
11301 D3	LED battery-powered source with quick-release thread
6927691	Three-way double stopcock adaptor with LUER-Lock
6011590	LUER cap
600007	Barbed tube, thin, with LUER-Lock

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13272 Mouth piece

11301 CF LIPP tube holder

11301 BM Adaptor for leakage test in Belimed cleaning machines

11301 CD Adaptor for machine cleaning the working channel (green)

11301 HH Adaptor for leakage test in BHT cleaning machines

11301 FF2 Adaptor for leakage test in ETD2/ETD3 cleaning machines with nonreturn valve

11301 GG2 Adaptor for cleaning the working channel in ETD2/ETD3 machines

11301 KK2 Connecting adaptor for cleaning the working channel in ETD3 cleaning machines

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- 20260030 DCI® adaptor, for connecting endoscopes with DIN eyepiece to DCI® camera heads
- 20260031 DCI® eyepiece adaptor for connecting DCI® telescopes to camera heads with mounts for DIN eyepieces

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3. Intended Use

Fiberscopes are used in diagnostic and surgical endoscopy. They must only be used by appropriately trained and qualified technical staff.

Indications: Examination, diagnosis and/or therapy in conjunction with endoscopic accessories in various medical fields.

Contraindications: No contraindications are known which directly relate to the product. The fibroscope must not be used if, in the opinion of a qualified physician, such an application would endanger the patient, e.g., due to the patient's general condition, or if the endoscopic method as such is contraindicated.

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3.1 Safety precautions when using the fiberscopes

Fiberscopes must be used according to the medical rules and procedures of endoscopy recognized for this method.

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4. User qualification

KARL STORZ fiberscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the endoscopic technique.

The information given in these instructions only serves to instruct in the correct handling, cleaning, disinfection and sterilization of the fiberscopes. It is not intended as an introduction to the endoscopic technique.

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5. Warnings and cautions

Please read this manual and follow its instructions carefully. The words, **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the text should be carefully reviewed to ensure the safe and effective operation of this product. To make the words **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

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CAUTION: A Caution indicates that particular service procedures or safety precautions must be followed to avoid possible damage to the product.

NOTE: A Note indicates special information about operating the product, or clarifies important information.

WARNING: Before using the fiberscope, read the following safety instructions carefully, to avoid putting your patients, personnel or yourself at risk.

WARNING: Inspect the fiberscope prior to each surgical procedure to ensure that it is functioning correctly and that it has been correctly cleaned, disinfected and (where applicable) sterilized.

WARNING: This equipment may not be used during the discharge of a cardiac defibrillator. Remove all endoscopic equipment from patient contact and from the surgical site.

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WARNING: Verify the sound condition of the fibroscope as well as of accessories used in combination before every application. Damaged fiberscopes or damaged accessories must not be used.

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Pay careful attention to the instructions and interface specifications of medical products used in combination.

WARNING: Lengthy direct contact between the endoscope tip and surrounding tissue may cause tissue damage due to the potentially high temperatures on the tip of the endoscope and must thus be avoided.

WARNING: The fibroscope is delivered unsterilized, and must therefore be cleaned, disinfected and (where applicable) sterilized prior to initial use and each subsequent reuse.

ATTENTION: The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

WARNING: Before the fibroscope is removed from the patient, the locking mechanism ⑮ must be released and the distal tip of the fibroscope put in the neutral, non-angled position. Otherwise the instrument may be damaged or the patient injured.

WARNING: If a malfunction should occur during application on a patient, stop the application immediately. Place the distal tip of the fibroscope in the neutral, non-angled position and remove the fibroscope slowly and carefully from the patient.

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CAUTION: Do not steam sterilize (autoclave) the fiberscopes.

CAUTION: The maximum permissible temperature for cleaning, disinfection, sterilization, transport and storage is 65 °C.

CAUTION: Never immerse fiberscopes in normal saline solution, since even only brief contact can lead to pitting and corrosion.

CAUTION: Material damage: if the outer casing of the fiberoscope is cleaned with alcohol, this will damage the material.

CAUTION: Fiberscopes must not be cleaned in an ultrasound bath under any circumstances.

CAUTION: Only use KARL STORZ accessories (guide probes, cleaning brushes etc.). Accessories from other manufacturers may damage the fiberscopes.

CAUTION: Only use guide probes with flexible ends. Other types of guide probe may cause considerable damage to fiberscopes.

NOTE: Any damage to the fiberoscope resulting from incorrect operation is not covered by the guarantee.

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NOTE: If the instrument is in the working channel, its additional rigidity may limit angling.

NOTE: Federal law restricts this device to sale by or on the order of a physician (USA only).

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6. Unpacking

Check for missing items and evidence of shipping damage. File any complaints to KARL STORZ or the supplier immediately.

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6.1 Basic equipment

1 Fiberscope	Art. no.
1 Suction valve insert, disposable version, pack of 20 inserts	11301 CE
1 Leakage tester	13242 XL
1 Pressure compensation cap	11025 E
1 Plug for instrument channel, pack of 10	29100
2 Adaptor for cleaning the working channel in cleaning machines (green)	11301 CD
1 Cleaning brush for 2.8 x 50 and for 2.8 x 70	27651 A*)
for 3.7 x 65	27651 AL*)
for 5.0 x 65	27651 B
for 5.6 x 65	27651 C
for 6.4 x 65	27651 C
1 Biopsy forceps (only bronchofiberscopes)	11001 KL
1 Grasping forceps (only bronchofiberscopes)	11002 KS
1 Three-way double stopcock adaptor with LUER-Lock for O ₂ connection (only intubation fiberscopes)	6927691
1 LIPP tube holder (only intubation fiberscopes)	11301 CF
1 Mouth piece (only gastro-, intubation and bronchofiberscopes)	13272
*) For the fiberscopes 11278 A/AI/AU/AIU use the cleaning brush 27651 A (27651 AL is not suitable for these fiberscopes).	

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6.1.1 Accessories

- | | |
|--|----------------------|
| 1 Suction valve insert,
reusable version | 11301 CB |
| 1 Battery light source | 11301 D |
| 1 Adaptor for leakage test in
cleaning and disinfection
machines (green) | 13242 XA |
| 4 Cleaning adaptors for
cleaning machines | 11301 GG2/FF2/HH/KK2 |

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7. Initial operation

WARNING: The fiberscope is delivered unsterilized, and must therefore be cleaned, disinfected and (where applicable) sterilized prior to initial use and each subsequent reuse.

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7.1 Fitting the suction valve

Place the suction valve insert ② on the valve nozzle ③ and turn until the projections on the side of the insert lie over the nozzle recesses.

Push gently until the suction valve insert clicks into the valve nozzle.

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7.2 Fitting the tube holder (only intubation fiberscopes)

Push the tube holder onto the fiberscope from distal and secure it by turning the two parts in opposite directions.

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7.3 Fitting the lip valve (only bronchofiberscopes)

Place the disposable lip valve onto the instrument channel ⑧ and tighten gently by turning it in a clockwise direction (plastic part).

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8. Handling the fiberscopes

The optical section of KARL STORZ fiberscopes consists of a fiberoptic bundle in a flexible and watertight casing, a light guide, as well as objective and eyepiece lenses. Like all optical systems, these parts are very delicate. For this reason, careful handling of the fiberscopes is essential for a long service life and a consistent quality. Therefore be aware of the following points:

- Never kink or tightly coil the fiberscopes and light guides, so that no optical fibers break off. It is therefore best to suspend the fiberscopes for storage or keep in large containers.
- Do not jerk the up/down control ④ or move forcefully.
- Do not insert the instruments or cleaning brushes forcefully through the channels. This may cause damage to the casing, especially if the tip of the instrument is deflected.
- Carefully protect the objective lens ⑨ at the distal end from damage. Do not hit. Never remove soiling with sharp-edged objects.
- Never store or disinfect the fiberscope together with sharp-edged objects to prevent damage to the exterior.
- Never place heavy objects on top of the fiberscope.

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NOTE: Every fiberscope supplied by KARL STORZ includes a photo showing the condition of the fiber bundle upon delivery. A copy of this photo is kept on file in our quality control office. The wedge-shaped marking (at 12 o'clock) serves as a guide for the user.

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9. Using the fiberscopes

The following instructions are recommended for the proper operation of the fiberscope and accessories. This is not a reference for endoscopic techniques.

Consult the appropriate medical literature for information on endoscopic techniques.

You can find accessories for your fiberscope (catheters, flexible scissors and forceps, guide wires, stone baskets, lithotripter probes, laser fibers, electrodes etc.) in the appropriate KARL STORZ catalog.

WARNING: Always adjust the light source to the minimum illumination intensity necessary to achieve optimum visualization of the operating field. The higher the intensity setting of the light source, the greater the heat energy that will be generated at the tip of the fiberscope.

WARNING: Never lay the end of the fiber optic cable or the distal tip of the fiberscope on or under a sterile cover cloth. The light intensity may lead to burns on the patient and/or sterile cover cloth.

- When connecting and operating light sources and suction/irrigation equipment, be certain to read the appropriate instruction manuals for the equipment. The irrigation equipment recommended by KARL STORZ is set to a maximum pressure of 54 kpa (400 mmHg/0,54 bar/7.8 psi).
- Endoscopic examinations must be carried out under aseptic conditions.
- When introducing the fiberscope into the patient, the distal tip of the fiberscope must be in a straight, nondeflected position (control lock-in can be felt and heard). The fiberscope must be handled according to the recognized medical rules and procedures for endoscopic methods.

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- If no instrument is inserted, close the instrument channel proximally with a plug (Art. no. 29100).
- If an instrument is to be inserted, position the disposable lip valve (Art. no. 11301 CA) to seal off the instrument channel proximally.
- Instruments especially designed for this may be inserted through the instrument channel of the fiberscopes. Make certain the accessories have a suitable length and diameter. Unsuitable accessories may damage the fiberscopes. For questions on accessories, please contact KARL STORZ directly or an authorized dealer.

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NOTE: When the lip valve is first used, the instrument is subjected to increased resistance. This resistance will decrease the more often the lip valve is used.

- When inserting an instrument into the instrument channel, the fiberscope must be in a straight, nondeflected position.

CAUTION: When inserting instruments, increasing resistance is felt with increasing deflection. Forced introduction of instruments can lead to damage to the fiberscope or the instrument channel.

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CAUTION: If undue resistance is encountered when inserting an instrument, do not attempt to force the instrument. Withdraw the instrument and apply a water soluble lubricant to the tip of the instrument and attempt to insert again. Prior to removal of the instrument from the patient, the distal tip must be in a straight position.

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- When using forceps, scissors or other accessories that can be opened or closed, ensure that the instrument is completely closed prior to insertion into the instrument channel. Do not attempt to open the jaws until the distal tip is completely free of the channel. Opening the jaws inside the instrument channel can lead to damage to the fiberscope.

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WARNING: When using HF electrodes or laser probes, make sure that the active parts are clear of the instrument channel and in full view of the operating surgeon before power is applied. Activation within the instrument channel may damage the fiberscope and/or injure the patient.

WARNING: Make certain to read the instructions for the laser unit and/or HF generator.

CAUTION: If accessories with sharp tips, such as EHL probes or laser probes, are inserted into the instrument channel of a deflected fiberscope, this can damage the fiberscope.

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CAUTION: When using a stone basket, the stone must not be retracted through the instrument channel but instead the entire fiberscope must be removed. Once the stone is removed, the stone basket may then be taken out of the instrument channel.

CAUTION: The fiberscopes must be cleaned and then reprocessed immediately after the intervention in order to prevent incrustations. Dried-on incrustations may jeopardize successful decontamination and disinfection.

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CAUTION: Use only KARL STORZ accessories (guide probes, cleaning brushes etc.). Accessories from other manufacturers may damage the fiberscopes.

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
NOTE: If the instrument is in the working channel, its additional rigidity may limit maneuverability.

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10. Preparation of the fiberscope

Before using the fiberscope on a patient, exterior damage and operational problems must be eliminated.

Check the outer sheath of the fiberscope. Be especially aware of visible damage to the casing, particularly in the area of the fiberscope tip.

Look through the eyepiece and make certain that you have a clear image. The grid (Moiré effect) should be easily recognizable. If the image is not clear but cloudy, clean the eyepiece and lens with a cotton applicator soaked in alcohol, maintaining hygienic conditions at all times.

DCI® fiberscope image quality can be controlled on screen. If necessary, clean image waveguides, light cables and eyepieces as described above.

If it is still cloudy, it is possible that the objective lens is scratched or that there is moisture inside the telescope. The fiberscope must not be used in this case, but instead sent in for repairs.

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10.1 Connect these supply lines

- If necessary, remove the adaptor 495 G or 495 F from the light inlet piece and connect the fiberoptic light cable or battery-powered light source to the light inlet piece ⑦.

NOTE: Not necessary for battery-powered light sources with coarse threads (11301 D3).

WARNING: The instruction manuals of the supply units used must be observed.

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- suction to suction outlet ⑥
- or O₂ supply (under anesthesia) to instrument channel ⑧ via the appropriate adaptor

CAUTION: The instruction manuals of the supply units used must be observed.

Switch on the respective supply units.

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Insert the DCI® fiberscope into the DCI® camera head until it locks into place.

Insert the light cable (Art. no. 495 DV) into the DCI® camera head as far as it will go.

WARNING: The instruction manuals of the supply units used must be observed.

Connect the light cable and the video signal cable up to the appropriate supply units.

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10.2 Controlling the distal tip

The up/down movement of the distal tip on the fiberscope can be adjusted at the control ④: move smoothly and slowly.

Some fiberscopes have:

- contrapositive deflection direction (urology);
- a locking mechanism ⑮ for securing the fiberscope tip (gynecology).

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Models with locking mechanism ⑮:

To secure the deflection position of the distal tip, turn the locking mechanism ⑮ slowly to distal until resistance can be felt.

CAUTION: Forcing the locking mechanism ⑮ further in the same direction may damage the locking mechanism ⑮.

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WARNING: Before the fiberscope is removed from the patient, the locking mechanism 15 must be released and the distal tip of the fiberscope put in the neutral, non-deflected position. Otherwise the instrument may be damaged or the patient injured.

To release the locking mechanism 15 from the secure position, turn it slowly to proximal.

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10.3 Suction

Press the suction valve 2 to start suction (when the suction pump has been connected and switched on).

CAUTION: There should be no instrument in the suction channel during suction. The instrument channel inlet should be closed so that no secondary air is drawn in.

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10.4 Insufflation

Under anesthesia, connect the O₂ hose to the instrument channel:

- using a barbed tube (Art. no. 600007) on the instrument channel or
- if you are using the three-way double stopcock adapter with LUER-Lock (Art. no. 6927691), attach the hose to the tube holder on the angled part of the adapter. This ensures that the straight-line passage of the adapter remains free for instrument insertion.

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10.5 Focusing

The image is focused with the focusing ring 1.

NOTE: When connecting photo or video systems recommended by KARL STORZ, the fiberscope focusing ring must be at the zero position (mark).

DCI® fiberscopes are focused using the adjustment wheel of the DCI® camera (see page 28).

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10.6 Instrument channel

Flexible instruments can be inserted through the instrument channel ⑧ as required.

If the instrument channel is not being used: replace the sealing cap.

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10.7 PDD version

Fiberscopes intended for use in Photodynamic Diagnosis (PDD) or autofluorescence (AF) are illuminated with a special diagnosis light (D-Light C). Adjust the PDD filter changer ⑫ for special filters designed for this light and possibly for marker substances:

ALA	for a special tumor marker based on aminolevulinic acid
W	for white light
AF	autofluorescence mode

NOTE: If the filter changer is set to 'ALA' or 'AF', a yellowish image is produced when viewing through the integrated special filter.

NOTE: No filter changer is available on PDD fiberscope art. no. 11272 CI. The filter (ALA) is permanently installed.

For optimal illumination, use the fluid light cable 495 FP (2.5 m) or 495 FO (1.8 m) between the light source and PDD fiberscope.

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11. DCI® fiberscope with camera

11.1 Assembly

Slide the DCI® into the snap-in DCI® coupling mechanism of the DCI® camera until the coupling snaps in automatically.

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Connecting the special light cable 495 DC or 459 DV to the DCI® camera head:

Insert the light cable plug into the DCI® camera head until it locks into place.

While doing this, hold the light cable by the handle piece. Never pull on the cable!

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NOTE: The entrance and exit faces of the light cable and endoscope must be kept clean if optimal light delivery is to be maintained. Clean the surfaces with a cotton applicator soaked in cleaning solution, then disinfect it with a cotton applicator soaked in 70 % isopropyl alcohol and dry.

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DCI® fiberscopes are focused using the adjustment wheel on the DCI® camera.

The camera control unit is controlled using the two programmable control keys (silver, blue) on the DCI® camera head.

NOTE: Further controls are described in the instruction manual of the camera control unit used.

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11.2 Disassembly

Hold the DCI® camera and the DCI® fiberscope still. Push the release button on the DCI® camera: the 'snap-in' DCI® coupling mechanism is released.

Pull the light cable out of the DCI® camera head.

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12. DIN eyepiece on DCI® fiberscope

The DCI® eyepiece adaptor 20260031 ensures the connection of DCI® telescopes to camera heads with DIN eyepiece mounts.

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12.1 Assembly

Slide the DCI® fiberscope into the 'snap-in' DCI® coupling mechanism of the DCI® eyepiece adaptor until the coupling snaps in automatically.

Connect a light cable or a battery-powered light source.

WARNING: The instruction manuals of the supply units used must be observed.

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DCI®

12.2 Disassembly

Remove the light cable or battery-powered light source.

Hold the DCI® eyepiece adaptor and the DCI® fiberscope. Press the release button on the DCI® eyepiece adaptor. The 'snap-in' DCI® coupling mechanism is released.

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13. DCI® camera on DIN eyepiece

13.1 Assembly

Put the focusing ring (if available) of the fiberscope in the zero position (mark).

Slide the DCI® adaptor 20260030 into the 'snap-in' DCI® coupling mechanism of the DCI® camera until the coupling snaps in automatically.

Connect the standard light cable at the light inlet piece ⑦ of the fiberscope.

NOTE: The special light cable 495 DC or 495 DV cannot be used and may have to be removed from the DCI® camera.

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Opening the instrument coupling: turn the outer ring against the spring force, insert the fiberscope DIN eyepiece and let go of the outer ring. The instrument coupling closes automatically as a result of the spring force.

CAUTION: No further tightening of the outer ring to position the eyepiece is required. This may, among other things, cause damage to the instrument coupler!

13.2 Disassembly

Hold the DCI® camera and the fiberscope still. Open the instrument coupler for the DCI® adaptor and remove the fiberscope.

Hold the DCI® camera and the DCI® adaptor still. Push the release button on the DCI® camera: the 'snap-in' DCI® coupling mechanism is released. Remove the DCI® adaptor.

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14. Reprocessing

General warnings

-  **WARNING:** Risk of infection: Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.
-  **WARNING:** Risk of infection from non-sterile delivery. Reprocess the medical device before each use. Check the medical device for visible soiling before use.
-  **WARNING:** Before reusing the device, replace any defective parts or send the flexible endoscope for repair.
-  **ATTENTION:** The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.
-  **ATTENTION:** Never clean flexible endoscopes with metal brushes.
-  **ATTENTION:** Constant changing between various reprocessing methods is hard on the materials and should be avoided. However, a one-off change of method, due to procurement of new cleaning devices, for example, is possible.
-  **ATTENTION:** Flexible endoscopes must never be steam sterilized (autoclaved).
-  **ATTENTION:** To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.

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-  **NOTE:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.
-  **NOTE:** Observe the chemical manufacturer's specifications concerning concentration and exposure times. Bear in mind the microbiological range of action of the chemicals.
-  **NOTE:** National laws and regulations must be observed.
-  **NOTE:** We recommend using a plastic basin for disinfection to avoid damaging the endoscope. Furthermore, use of a plastic basin prevents electrolytic corrosion which can occur if there are different metals in a solution.
-  **NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

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14.1 Accessories

Accessories required for carrying out reprocessing: pressure compensation cap, manual leakage tester (different adaptor sizes)

Application site	Art. number for brush/adaptor
Surface	27652, dry, low-lint cloth
Working channel	Use a suitable brush, adaptor and syringe
Suction channel	Use a suitable brush, adaptor and syringe
Vent port	Use a suitable adaptor

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14.2 Precleaning

- Remove the entire instrument set from the instrument channel.
- Remove the suction valve set from the valve nozzle: press the projections on the side together and pull out the suction valve insert.
- Dispose of the removed single-use parts.

Dismantle other parts if necessary and reprocess them separately.

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the endoscope immediately after use.

Wipe-clean the exterior surfaces of the endoscope with a damp cloth.

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Rinse out the working channel in order to prevent drying.

i NOTE: To avoid damaging the sheath, do not use pushing and pulling movements.

i NOTE: Do not use disinfectants for precleaning; these cause fixation of proteins and organic residues.

14.3 Leakage test

i NOTE: The rubber seal in the connection adaptor of the leakage tester should be treated regularly with instrument oil.

14.3.1 Phase 1: Performing the leakage test – dry

The connection tube of the leakage tester and the vent port on the endoscope must be dry.

1. Attach the adaptor of the leakage tester securely to the product's vent port and lock in place using the bayonet mount.
2. Using the leakage tester, pump air into the product until a pressure of 160 – 180 mmHg (blue area) is reached.
- ➡ The angled rubber at the distal end can easily expand while the air is pumped into the product.
3. Using the angling mechanism, move the distal end on the product's control body in all directions as far as possible while observing the pressure display of the leakage tester.
- ➡ The air pressure in the product can drop by up to 10 mmHg when moving the distal end.
4. After stabilizing the air pressure in the product, wait at least 2 minutes and check on the leakage tester whether the air pressure in the product drops.
5. If the air pressure drops, the product is leaking and may no longer be used, immersed in liquid or reprocessed by machine. Complete the repairs form. Clean the product by wiping to remove organic soiling/residue and mark it

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as not disinfected. Contact KARL STORZ and send the product to KARL STORZ for repair. Line the carrying case for the product with the transport protective film 13990 SFN.

- ☛ If the air pressure does not drop, the product is leakproof. The leakage test: Phase 2 can be performed with the product.

i **NOTE:** To increase safety, the leakage test: Phase 1 can also be carried out directly after using the product.

14.3.2 Phase 2: Performing the leakage test – in liquid

1. Lay the product in liquid and ensure that the leakage tester is not moist.
2. Fully cover the product with liquid and wet any lumens present without bubbles.
3. Using the leakage tester, pump air into the product until a pressure of 160 – 180 mmHg (blue area) is reached.
- ☛ The angled rubber at the distal end can easily expand while the air is pumped into the product.
4. Using the angling mechanism, move the distal end on the product's control body in all directions as far as possible while observing the pressure display of the leakage tester.
- ☛ The air pressure in the product can drop by up to 10 mmHg when moving the distal end.
5. After stabilizing the air pressure in the product, wait at least 3 minutes and check on the leakage tester whether the air pressure in the product drops or whether bubbles are continuously escaping along the product.

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- ☛ If the air pressure drops or bubbles emerge from the product, the product is leaking and may no longer be used, immersed in liquid or reprocessed by machine. Complete the repairs form. Clean the product by wiping to remove organic soiling/residue and mark it as not disinfected. Contact KARL STORZ and send the product to KARL STORZ for repair. Line the carrying case for the product with the transport protective film 13990 SFN.
- ☛ If the air pressure does not drop and no bubbles emerge from the product, the product is leakproof. Remove the product from the liquid and then press the vent button on the pressure tester. Wait until the air pressure has dropped to 0 mmHg. The product can now be cleaned.

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14.4 Disassembly

Before cleaning and disinfection, remove the adaptor and other detachable parts from the light inlet piece, the instrument channel and the stopcock in their parts.

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14.5 Manual reprocessing

14.5.1 Manual cleaning

1. If relevant, connect the adaptor.
2. Completely immerse the product in the cleaning solution.
3. If relevant, fill the lumen with cleaning solution using the syringe and ensure bubble-free wetting.
4. If relevant, remove the adaptor.
5. Observe the exposure time as specified by the chemical manufacturer during the cleaning process.
6. Brush the surfaces of the handle at least 3 times until no residue can be seen anymore.
7. The mobile parts must be brushed.
8. If relevant, clean the valve seats by brushing them and turning them 360° 3 times until no residue can be seen anymore.

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9. Clean the surface of the sheath with a low-lint cloth until no residue can be seen anymore.
10. If relevant, brush the lumen at least 3 times until no residue can be seen anymore on the brush head. Clean the brush head between each stroke of the brush.
11. If relevant, connect the adaptor.
12. If relevant, irrigate the lumen with cleaning solution 3 times.
13. After the exposure time, remove the product from the cleaning solution.
14. If relevant, irrigate the lumen with air.
15. If relevant, remove the adaptor.
16. Dry the surface with a soft, low-lint cloth.
17. If relevant, connect the adaptor.
18. Completely immerse the product in water of the best possible quality.
19. If relevant, fill the lumen with water using the syringe and ensure bubble-free wetting.
20. If relevant, remove the adaptor.
21. Remove the product from the water.
22. Dry the product.

ATTENTION: The pressure compensation cap must not be in place on the vent port for cleaning and disinfection.

14.5.2 Manual disinfection

The endoscope must be completely immersed in a disinfectant solution. The working channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

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14.5.3 Manual High-Level Disinfection (HLD)

The endoscope must be completely immersed in a high-level disinfectant solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

The endoscope can be disinfected with the following high-level disinfection solutions which have been tested by way of example:

- 2.0% accelerated hydrogen peroxide solution (e.g., Revital-Ox® Resert™ or Revital-Ox Resert XL HLD)

14.5.4 Drying

Finally, the surfaces and channels/lumens are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

ATTENTION: The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

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14.6 Machine reprocessing

The following procedure for machine decontamination has been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'.

14.6.1 Manual preparation

1. If relevant, connect the adaptor.
2. Completely immerse the product in the cleaning solution.
3. If relevant, fill the lumen with cleaning solution using the syringe and ensure bubble-free wetting.
4. If relevant, remove the adaptor.
5. Observe the exposure time as specified by the chemical manufacturer during the cleaning process.
6. Brush the surfaces of the handle at least 3 times until no residue can be seen anymore.
7. The mobile parts must be brushed.
8. If relevant, clean the valve seats by brushing them and turning them 360° 3 times until no residue can be seen anymore.
9. Clean the surface of the sheath with a low-lint cloth until no residue can be seen anymore.
10. If relevant, brush the lumen at least 3 times until no residue can be seen anymore on the brush head. Clean the brush head between each stroke of the brush.
11. If relevant, connect the adaptor.
12. If relevant, irrigate the lumen with cleaning solution 3 times.
13. After the exposure time, remove the product from the cleaning solution.
14. If relevant, irrigate the lumen with air.
15. If relevant, remove the adaptor.
16. Dry the surface with a soft, low-lint cloth.

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17. If relevant, connect the adaptor.
18. Completely immerse the product in water of the best possible quality.
19. If relevant, fill the lumen with water using the syringe and ensure bubble-free wetting.
20. If relevant, remove the adaptor.
21. Remove the product from the water.

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14.6.3 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 60°C and a maximum pressure of 100 kPa (1 bar/14.5 psi) must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

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14.7 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.

NOTE: Clean the distal lens with a cotton tip applicator soaked in lens cleaner or 70% isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement.

- Damaged or corroded medical devices must be withdrawn from use.

NOTE: Inspect plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.

- If there is still residual moisture on or in the endoscope, it must be dried again manually. The surfaces (e.g., connections) and channels/lumina must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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NOTE: After treating with grease, the stopcocks must be sterilized separately using the fractionated steam sterilization process. Low-temperature sterilization is not possible on surfaces which have been lubricated and oiled.

- Afterwards, a functional check must be carried out.

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- Clean the eyepiece and lens with a soft cloth, sponge or cotton-tipped applicator soaked in 70% alcohol.

CAUTION: Cleaning with the special cleaning paste should only be carried out if the image as viewed through the endoscope is visibly cloudy or if there are visible deposits on the optical surfaces, it should certainly not be used routinely during reprocessing.

Stubborn soiling and deposits can be removed with cleaning paste (Art. no. 27661). For this purpose, some cleaning paste is put on a clean, damp cleaning bud and distributed on the glass surfaces. Polishing ensures that even stubborn deposits are removed.

NOTE: The instrument must be reprocessed again following treatment with the cleaning paste!

Finally, dry the surfaces of the fibers and the optical end faces with a cotton-tipped applicator soaked in 70% isopropyl alcohol.

NOTE: The fiber area in the light inlet piece in particular ⑦ must be dried carefully with alcohol. Disinfectant and cleaning agent residue in the light inlet piece can become burnt on when the light cable is connected, thereby seriously impairing light transmission.

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14.9 Sterilization

-  **WARNING:** Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.
-  **WARNING:** Risk of infection: Sterilize the medical device in a disassembled state.
-  **ATTENTION:** Flexible endoscopes must never be steam sterilized (autoclaved).
-  **ATTENTION:** The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.
-  **NOTE:** All deviations from the recommended sterilization parameters must be validated by the user.
-  **NOTE:** When using reprocessing devices, please observe the device manufacturers' instructions.
-  **NOTE:** Sterilization is not possible on surfaces which have been lubricated and oiled.

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

14.9.1 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

-  **ATTENTION:** Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

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NOTE: Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S Long Cycle with/without booster
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

14.9.2 Hydrogen peroxide (H₂O₂) sterilization – STERIS® V-PRO®

NOTE: The 'STERIS Low-Temperature Reprocessing Device Compatibility Matrix' can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.

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14.9.3 Ethylene oxide sterilization (EO)

The ethylene oxide method is validated with 100% ethylene oxide according to the ANSI/AAMI/ISO standard 10993-7.

EO sterilization is validated with the following parameters:

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14.9.4 Chemical low-temperature sterilization with peracetic acid – STERIS System 1E®

Detailed information on the selection of the sterilization parameters and the STERIS® Quick Connect Kit (QCK) required to rinse out the lumina is available from the manufacturer STERIS®.

! **ATTENTION:** The pressure compensation cap must not be in place on the vent port.

- The glide surfaces of the stopcock must be lubricated locally with special grease.

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15. Standard compliance

15.1 Servicing and repair

Defective instruments must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts.

15.2 Responsibility

As the supplier of this instrument, KARL STORZ shall be liable for the safety, operational reliability and performance of this instrument subject to the following conditions:

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15.3 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see chapter 'Subsidiaries'), even during the warranty period.

Opening the device or performance of any repairs or modifications to the device by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair or modification performed during the warranty period shall void all warranty.

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15.4 Directive compliance

This medical device bears the CE mark in accordance with the Medical Devices Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

16. Repair and replacement program

KARL STORZ maintains a repair and replacement warehouse which is normally adequate to ensure prompt replacement of damaged telescopes and instruments. Under the repair and replacement plan, you receive an identical as-new device and are only charged the repair costs for the defective device. For telescopes, you receive a warranty of 1 year, for instruments 1/2 year.

For fiberscopes and devices, individual repair is necessary. To bridge the repair period, you will receive a device on loan which you then return to KARL STORZ as soon as you receive the repaired device.

In Germany, you can refer repairs direct to:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Straße 34
78532 Tuttlingen

In other countries, please contact your local KARL STORZ subsidiary or authorized dealer.

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16.1 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site to avoid contact and aerogenous infections (among personnel).

We reserve the right to return contaminated instruments/devices to the sender.

Repairs, modifications or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no warranty on the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

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KARL STORZ Endoscopy Canada Ltd.
1717 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Av. Ejercito Nacional No. 453 Piso 2,
Colonia Granada, Alcaldia Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA – Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamveien 1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se
KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kęstucio st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lt-iv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.com

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Polska Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Indonesia
Sinarmas MSIG Tower Level 37
Jl. Jend. Sudirman No. Kav. 21
Jakarta Selatan
DKI Jakarta 12920
E-Mail: infoindonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd.
68 Waterloo Road, Macquarie Park NSW 2113
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: karlstorz@karlstorz.com.au
www.karlstorz.com

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This instruction manual is intended to aid you with proper connection and operation of the TIPCAM®1 S 3D in conjunction with the instruction manual for the IMAGE1 S™. All required details and all actions on your part are clearly presented and explained. We thus ask that you read these manuals carefully before proceeding to work with the equipment. Insert in the protective wallet and keep available for ready reference in a convenient and conspicuous location near the equipment.

Please refer to the IMAGE1 S™ instruction manual for the individual functions of the CCU that can be used with the TIPCAM®1 S 3D camera.

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Rx only Federal (USA) law restricts this device to sale by or on the order of a physician.

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
**2 Controls, displays,
connectors, and their uses**

- ① Inserter (sheath) with sensors
- ② Handle
- ③ Control button 1: Scroll up in menu/ programmed functions
- ④ Control button 2: Selection, open menu
- ⑤ Control button 3: Scroll down in menu/ programmed functions
- ⑥ Connection for fiberoptic light cable
- ⑦ Video connecting cable connector to CCU
- ⑧ Fiberoptic light cable (495TIP)
- ⑨ Light cable connector to KARL STORZ light source
- ⑩ Video connecting cable

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Number of products in the product packaging

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CE mark. The CE mark is a marking required for certain marketable industrial products in accordance with EU law. It is composed of the CE logo and, where applicable, the four-figure identification number of the responsible testing center (when the product belongs to a certain risk category).

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Federal (USA) law restricts this device to sale by or on the order of a physician

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1	Images of the unit	IV
2	Controls, displays, connectors, and their uses	V
3	Symbols employed.....	VI
3.1	Symbols for operation	VI
3.2	Label and packaging symbols	VI
4	General information	4
4.1	Device description	4
5	Safety instructions	5
5.1	Warnings and cautions	5
5.2	Application areas (worldwide, does not apply to the USA and Canada).....	12
5.2.1	Intended use.....	12
5.2.2	Indication	12
5.2.3	Contraindications.....	12
5.3	Application areas (for US and Canadian customers only)	13
5.3.1	Indications for use	13
5.3.2	Contraindications.....	13
5.4	User qualification.....	14
5.5	Safety precautions when operating the instrument	14
6	Installation and operating instructions	15
6.1	Unpacking the equipment	15
6.2	Basic equipment	15
6.3	Compatible camera control unit (CCU)	15
6.4	Before every use	15
6.5	Start-up of the TIPCAM®1 S 3D	16
6.6	Connecting the light cable to the TIPCAM®1 S 3D	16
6.7	Connecting the TIPCAM®1 S 3D to the light source	17
6.8	Connecting the TIPCAM®1 S 3D to the D3-LINK.....	17

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6.9	Switching on the CCU, light source and monitor	18
6.10	Optimum light efficiency	18
7	Reprocessing not validated for US and Canadian customers	19
7.1	General care	20
7.1.1	Care and handling	20
7.2	General warnings	23
7.3	Accessories	24
7.4	Disassembly	24
7.5	Preparation for cleaning and disinfection	24
7.6	Manual precleaning	24
7.7	Manual cleaning	25
7.8	Manual disinfection	25
7.9	Machine cleaning and disinfection	25
7.10	Assembly, inspection and care	26
7.11	Packaging systems	27
7.12	Sterilization	27
7.12.1	Accessories	29
7.13	Reprocessing limits	29
8	Reprocessing for US and Canadian customers only	30
8.1	General care	30
8.1.1	Care and handling	30
8.2	General warnings	32
8.3	Disassembly	33
8.4	Pre-cleaning	33
8.5	Manual cleaning	34
8.6	Inspection after cleaning	35
8.7	Sterilization	36
8.8	Limits of reprocessing	36

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9	Maintenance	37
9.1	Maintenance and safety check.....	37
9.1.1	Maintenance.....	37
9.1.2	Safety check.....	37
9.2	Servicing and repair.....	38
9.3	Disposal	38
9.4	Repair program	39
9.5	Limitation of liability	40
9.6	Warranty.....	40
10	Technical description.....	41
10.1	Technical data TIPCAM®1 S 3D	41
10.1.1	Directive compliance (for TIPCAM®1 S 3D)	41
11	Recommended accessories.....	42
12	Information on electromagnetic compatibility (EMC).....	43
13	Subsidiaries	60

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4 General information

4.1 Device description

The TIPCAM®1 S 3D is a high-resolution videoendoscope with two distally integrated video chips.

Since the TIPCAM®1 S 3D does not require an optical imaging system, it is very sturdy.

The TIPCAM®1 S 3D is compatible over D3-LINK with the IMAGE1 S™.

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5 Safety instructions

5.1 Warnings and cautions

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of the device.

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CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the unit.

NOTE: A Note indicates special information about operating the unit.

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WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk. Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

WARNING: Also pay particular attention to the instruction manual for the IMAGE1 S™ CCU no. 96206286DER.

WARNING: Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.

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WARNING: Risk of injury to the patient/user: Damaged or incomplete instruments/accessories may cause malfunctions. Before and after every use, always check the entire system (e.g., the TIPCAM®1 S 3D, all cables, the CCU, monitor and accessories used) to ensure that everything is configured properly, in perfect working order and nothing is missing/damaged. Replace incomplete/damaged instruments/accessories and do not continue to use them.

- Before every use, check that all the instruments required for the operation are functioning correctly.
- Be certain to check the cover glass of the TIPCAM®1 S 3D. In addition, the camera cable and the light cable need to be inspected for damage cracks.
- Do not use wet/moist cables.
- When attaching components to an instrument, the secure position of the components must be checked by gently pulling and/or turning them in both directions before use.

WARNING: Check the instrument and accessories used in combination before every application for unintentional rough surfaces, sharp corners, burred edges, and projecting parts, which can endanger personal safety. Damaged, defective, bent, or incomplete instruments or accessories must not be used further under any circumstances.

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WARNING: Risk of injury: Overloading the instrument by exerting too much force on it may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original positions.

WARNING: Risk of injury and risk of damaging the products: When using applied parts outside of the field of vision there is a risk of unintentionally injuring tissue and damaging accessories. Always keep the applied parts of the active electrode and the laser and other instruments which transmit energy pointed at the target in the field of visions during use.

WARNING: Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

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WARNING: Risk of injury and risk of damaging products: The high level of light energy means that the distal end, the light port, adjacent components and tissue in front of the light outlet window to heat up. This can cause burns to the patient, user and surgical accessories.

- Avoid direct tissue contact with the distal end, as well as with the light port of the telescope and the light cable.
- Avoid excessively long (> 10 min) continuous direct tissue contact with the shaft location.
- Never allow the endoscope to rest on the patient or come into direct contact with surgical accessories.
- Always select the lowest possible light output of the cold light source which still allows optimal illumination of the operating field.

WARNING: Risk of injury: Optical radiation poses a risk of injury to eyes. Never look into the free end of either a connected fiberoptic light cable or an endoscope. Always wear standardized protective equipment when using lasers.

WARNING: The TIPCAM®1 S 3D must be replaced if the image is cloudy or if no image or only parts thereof are visible.

WARNING: Risk of fatal injury: In the case of patients with cardiac pacemakers it must be ensured that the proximal end (= furthest from the patient) of the TIPCAM®1 S 3D is not positioned less than 2 cm from the pacemaker. The device generates a magnetic field at the handle, which can cause cardiac pacemakers or other electronic devices to malfunction.

WARNING: The video system can malfunction during use. An adequate replacement system must be kept available.

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-  **WARNING:** In the event of image failure immediately remove TIPCAM®1 S 3D from the patient.
-  **WARNING:** The use of HF probes, HF electrodes, HF instruments or other energy emitting probes is not permitted if explosive gases (for example, explosive/flammable anesthetics) are present.
-  **WARNING:** Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope applied part type CF is used. In this case an endotherapy device with a type CF applied part must be used in order to reduce the total patient leakage current.
-  **WARNING:** Live parts or instruments may not come into contact with TIPCAM®1 S 3D.
-  **WARNING:** Unintentional grounding of the patient: The metal part of the camera connector is electrically connected to the sheath of the endoscope. During use, make certain that these metal parts do not come into contact with ground potential.
-  **WARNING:** Risk of electric shock: Do not attempt to open the TIPCAM®1 S 3D. Refer servicing to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.

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WARNING: For use only by clinicians trained in laparoscopy!

CAUTION: When in contact with or in close proximity of an active electrosurgical instrument temporary image interference or loss can occur. An upgrade of the TC302 hardware is recommended if the system reports this.

CAUTION: Ultrasonography before laparoscopy may identify clinical conditions that will alter patient management.

CAUTION: Abdominal puncture sites 7 mm or greater (for the introduction of auxiliary instrumentation) may be a source of herniation.

CAUTION: The abdomen can be adequately distended by pressure in the range of 15-20 mmHg. It is seldom necessary to use an abdominal pressure greater than 20 mmHg.

CAUTION: Lasers and electrosurgical probes should not be activated simultaneously. Further, when one is activated, the tip of the other should be completely retracted. This avoids deflection of the energy to the other tip.

CAUTION: During monopolar electrosurgery, inadvertent burns can occur when the appropriate patient return path is obstructed or as a result of capacitive coupling.

CAUTION: Avoid combinations of devices made of both conduction and insulating materials, especially metal trocars with nonconductive trocar anchors, screws, etc.

CAUTION: Periodically check the video system.

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Moisture which has penetrated the system due to damage can lead to faulty insulation paths within the videoendoscope. Rough surfaces or sharp edges can lead to patient injury.

i NOTE: Follow the applicable national ordinances and recycling plans regarding disposal or recycling of device components.

Acquaint yourself with the operation and control of the instrument before initial use of the TIPCAM®1 S 3D.

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5.2 Application areas (worldwide, does not apply to the USA and Canada)

5.2.1 Intended use

The rigid video endoscope is used together with the camera control unit for diagnostic and/or surgical procedures when endoscopic video support is required and is intended for short-term use.

5.2.2 Indication

The rigid video endoscope is intended for use in endoscopy and endoscopic surgery within the peritoneal and thoracic cavity.

5.2.3 Contraindications

The use of the rigid video endoscope is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

The rigid video endoscope must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

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5.3 Application areas (for US and Canadian customers only)

5.3.1 Indications for use

The rigid video endoscope is intended to be used together with the camera control unit during diagnostic and/or surgical procedures when endoscopic video assistance is required. For use in all endoscopy and endoscopic surgery within the peritoneal and thoracic cavity, including gynecological and urological anatomy.

5.3.2 Contraindications

This instrument must not be used for observation or treatment of the heart and must not contact the heart or any area near the heart. In addition, this instrument must not come into contact with any device or therapeutic accessory that contacts the heart or any area near the heart.

Laparoscopy may be contraindicated by pregnancy greater than 16 weeks gestation or abdominal mass of comparable size.

Laparoscopy may also be contraindicated by the following conditions, depending on their severity or extent:

- Class IV cardiac decompensation
- Bowel obstruction
- Ileus
- Cardiac disease
- Intraperitoneal hemorrhage
- Diaphragmatic hernia
- Infection with acute peritonitis
- Previous abdominal surgery
- Obesity
- Thin nulliparous patient
- Chronic obstructive lung disease
- Liver failure with established collateral vessels

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5.4 User qualification

The TIPCAM®1 S 3D may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the TIPCAM®1 S 3D.

5.5 Safety precautions when operating the instrument

It is the user's responsibility to make sure the instrument is safe and operates properly before using the instrument.

During treatment in conjunction with the TIPCAM®1 S 3D, the patient must be treated and observed with the standard medical care. This includes keeping a check on the progress of treatment, as well as monitoring sterile application conditions where required by the type of intervention.

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6 Installation and operating instructions

6.1 Unpacking the equipment

Carefully unpack the TIPCAM®1 S 3D and remove the unit from its packaging. Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

6.2 Basic equipment

- 1 KARL STORZ TIPCAM®1 S 3D
- 1 Instruction manual 96206364DER

6.3 Compatible camera control unit (CCU)

CCU

IMAGE1 S CONNECT® (TC200) in combination with IMAGE1 S D3-LINK® (TC302)

6.4 Before every use

Before every use, the instrument must be checked for damage, especially around the distal end (cover glass, light outlet, ...).

Before each use, the video connecting cable and light cable must be checked for damage, especially breaks and cuts.

Instruments with damaged cables must not be used further under any circumstances. Cables must be replaced (light cable to be replaced by the user or the Service Dept., and the video connecting cable by the Service Dept.).

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6.5 Start-up of the TIPCAM®1 S 3D

! **NOTE:** The TIPCAM®1 S 3D including the equipment connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations and standards.

6.6 Connecting the light cable to the TIPCAM®1 S 3D

The light cable is mounted by fitting to the light port provided and turning clockwise. Check that it is firmly connected (see fig. A).

! **CAUTION:** Please note that the circular symbol on the TIPCAM®1 S 3D and the circular symbol on the light cable must correspond (see following table).

Overview: Recommended combinations*

Diameter light cable	
4.8 – 5.0 mm	⊙
3.0 – 3.5 mm	⊗
2.0 – 2.5 mm	⊙
Diameter endoscope	
6.6 – 12.0 mm	⊙
3.0 – 6.5 mm	⊗
0.8 – 2.9 mm	⊙

*Special endoscopes may deviate slightly

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6.9 Switching on the CCU, light source and monitor

i NOTE: Read the instruction manuals for the IMAGE1 S™, the light source, the light cable and the monitor.

6.10 Optimum light efficiency

i NOTE: The entrance and exit faces of the light cable and TIPCAME1 S 3D must be kept clean if optimal light transmission is to be maintained. The light cable must first be unscrewed (see Fig. C).
Following this, clean the light entrance and exit faces with a cotton swab soaked in cleaning solution and then wipe with a cotton swab moistened with 70 % isopropyl alcohol. Please note that both the glass surface inside the light port of the TIPCAME1 S 3D (see arrow, Fig. D.) and the light cable itself (see Fig. E) must be cleaned in this way. Please note the instructions on cleaning the light cable port on page 26.

! CAUTION: Heat from the light source can cause chemical residue on the light cable and on the light entrance of the endoscope to burn. This can damage the light fibers and have a detrimental effect on the lighting power and the picture quality.

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**Reprocessing
(worldwide, does not apply to
USA and Canada)**

**7 Reprocessing
not validated for US and
Canadian customers**

WARNING: The light cable must be removed for reprocessing the TIPCAM®1 S 3D. To do this, unscrew the light cable in a counterclockwise direction. The video connecting cable is firmly connected to the TIPCAM®1 S 3D and must not be removed for reprocessing.

NOTE: Instructions for US and Canadian customers can be found in chapter 8 on page 30.

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7.1 General care

WARNING: Always pull out power plug before cleaning.

Store the TIPCAM®1 S 3D in the storage/sterilization box available as an optional extra. Do not coil the TIPCAM®1 S 3D video connecting cable into tight loops or cause it to kink while storing. Do not store the TIPCAM®1 S 3D in direct sunlight or excessive heat.

When returning the TIPCAM®1 S 3D for service, the user must ensure the equipment being returned is clean and (wipe-down) disinfected. The light entrances and exit faces of the light cable and the endoscope must always be kept clean (please refer to the corresponding instruction manual) to ensure optimal reproduction of the endoscopic image.

CAUTION: Chemical residue on the cardedge connector can cause camera failure.

CAUTION: Heat from the light source can cause chemical residue on the light cable and on the light entrance of the endoscope to burn. This can damage the light fibers and have a detrimental effect on the lighting power and the picture quality.

7.1.1 Care and handling

To maximize the life of your TIPCAM®1 S 3D please follow the guidelines given below for proper handling and care of the camera head.

DO NOT PLUG THE VIDEO CONNECTING CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS WET. Always ensure the connector is completely dry and free of debris. Where necessary, remove any residual moisture with a low lint disposable cloth before plugging in the unit. Nor should the fiberoptic light cable be inserted into the receptacle of the light source if it is wet. The connector must always be free of residue and any remaining moisture. Where necessary, remove remaining moisture with a low lint, disposable cloth.

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- This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility, however must be validated by the operator on site in respect of effectiveness.
- x On the grounds of potential damage, this reprocessing step must not be performed.

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Reprocessing (worldwide, does not apply to USA and Canada)

7.2 General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

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CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

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CAUTION: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

7.3 Accessories

Accessories required for carrying out reprocessing:

- Brush: 27652

7.4 Disassembly

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible. The light cable must be removed for reprocessing the TIPCAM®1 S 3D. To do this, unscrew the light cable in a counterclockwise direction. The video connecting cable is firmly connected to the TIPCAM®1 S 3D and must not be removed for reprocessing.

7.5 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

7.6 Manual precleaning Brushing the surfaces

Any visible contamination or heavy soiling of the surfaces must be removed by cleaning them under cold running water with the aid of a brush or sponge.

No use of ultrasound

For technical reasons, the medical device is not suitable for ultrasound treatment.

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7.7 Manual cleaning

И **NOTE:** In order to protect the contacts of the video connecting cable connector, manual cleaning and disinfection is only possible if the protective cap is fitted (see Fig. F).

The medical device must be rinsed under cold running water for 2 minutes. Then completely immerse the medical device in a cleaning solution. It must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of the brush for 3 minutes. Finally, it must be rinsed with cold water for 2 minutes to ensure neutralization.

7.8 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. To ensure bubble-free wetting, the lumina in particular must be filled. At the end of the necessary exposure time, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces, joints, openings, channels and lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no.27660) is ideal for this purpose.

7.9 Machine cleaning and disinfection

И **NOTE:** In order to protect the contacts of the video connecting cable connector, machine cleaning and disinfection is only possible if the protective cap is fitted (see Fig. F).

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

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Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A₀ value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off afterwards by hand.

CAUTION: After the program has run, the TIPCAM®1 S 3D is at a high risk of corrosion due to the ambient conditions in the reprocessing machine. Once the program has ended, remove the TIPCAM®1 S 3D from the washer and disinfectant as soon as possible.

7.10 Assembly, inspection and care

CAUTION: After machine reprocessing, clean the distal lens and connector of the video connecting cable and fiberoptic light cable with a lens cleaner or a cotton tip applicator moistened in 70 % isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement. Please pay attention to chapter 6.10 "Optimum light efficiency".

CAUTION: After machine reprocessing, the connection for the fiberoptic light cable must be removed in order to clean the optical fiber area with a cotton tip applicator moistened with 70 % isopropyl alcohol. Then the connector needs to be screwed on again. The connector can be removed by turning it to the left counterclockwise. The connector can be fastened by turning it to the right clockwise (see fig. G – I).

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

i NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

7.11 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

7.12 Sterilization

i NOTE: The protective cap of the video connecting cable connector must be disassembled for sterilization (see Fig. J).

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

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Steam sterilization using the fractionated prevacuum procedure

NOTE: The TIPCAM®1 S 3D is located in its storage/sterilization container.

NOTE: Sterilize greased components in a disassembled state.

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132°C - 137°C with a minimum exposure time of 4 minutes and maximum of 18 minutes.

NOTE: This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.

- or 134°C - 137°C with a minimum exposure time of 3 minutes and maximum of 18 minutes.

CAUTION: Do not immerse the TIPCAM®1 S 3D in cold liquid after removing it from the steam sterilizer. Due to the 'temperature shock', this may result in damage to the camera lens and invalidation of the warranty.

Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

WARNING: Medical devices must be dismantled for sterilization.

WARNING: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

CAUTION: For more detailed information, please consult the user handbook of the respective device.

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i NOTE: Please refer to the 'STERRAD® Sterility Guide' to verify whether the respective medical device can be sterilized using the various STERRAD® devices.

i NOTE: Sterilization is not possible on surfaces which have been greased and lubricated.

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7.12.1 Accessories

Accessories recommended for reprocessing, transport and storage:

- Tray: 39501XTC

7.13 Reprocessing limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

! **CAUTION:** Inspect the video connecting cable and fiberoptic light cable for breaks or cracks before and after disinfection or sterilization, as well as prior to subsequent use. Breaks or cracks will allow fluids to enter the cable and cause damage to it. Where a TIPCAM®1 S 3D has a damaged fiberoptic light cable or damaged video connecting cable, the relevant cable must be replaced.

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Reprocessing for US and Canadian customers only

8 Reprocessing for US and Canadian customers only

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Store the TIPCAM®1 S 3D in the sterilization tray. The tray is an optional accessory that may be purchased from KARL STORZ. Do not coil the TIPCAM®1 S 3D video connecting cable into tight loops or cause it to kink while storing.

Do not store the TIPCAM®1 S 3D in direct sunlight or excessive heat.

Prior to returning for service, the TIPCAM®1 S must be thoroughly cleaned and sterilized per validated procedures. The light entrances and exit faces of the light cable and the endoscope must always be kept clean (please refer to the corresponding instruction manual) to ensure optimal reproduction of the endoscopic image.

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DO NOT PLUG THE VIDEO CONNECTING CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS WET. Always ensure the connector is completely dry and free of debris. Where necessary, remove any residual moisture with a low lint disposable cloth before plugging in the unit. Nor should the fiberoptic light cable be inserted into the receptacle

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**Reprocessing
for US and Canadian
customers only**

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of the light source if it is wet. The connector must always be free of residue and any remaining moisture. Where necessary, remove remaining moisture with a low lint, disposable cloth.

The connectors of the video connecting cable and of the fiberoptic light cable must also be completely dry before connecting!

NEVER PICK UP THE TIPCAM®1 S 3D BY ITS CABLE! Always pick up the set by the camera head first, followed by the cable and connector. The cable contains tiny electronic conductors and should be treated gently and carefully.

PROTECT THE TIPCAM®1 S 3D FROM DIRECT SUNLIGHT OR EXCESSIVE HEAT.

HANDLE THE SYSTEM WITH CARE, avoiding extreme impacts to device. When transporting, use the original shipping box to prevent damage to the device.

DO NOT STORE THE TIPCAM®1 S 3D WITH CABLE COILED TIGHTLY OR GATHERED INTO FOLDS. Always arrange the cable in loose coils of 6 inches (15 cm) or more in diameter.

DO NOT RAPIDLY PULL THE CABLE WHEN UNCOILING. Harmful "kinks" can result. Always allow the cable to uncoil gently and loosely.

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8.2 General warnings

WARNING: Risk of infection: Improperly reprocessed medical devices may pose a risk of infection for patients, users, and third parties and may cause malfunction of the medical device.

WARNING: Risk of infection: These medical devices are not supplied in sterile condition. The use of non-sterile medical devices poses a risk of infection for patients, users, and third parties. Check medical devices for visible soiling. Visible soiling indicates missed or improper cleaning. Reprocess medical devices using validated procedures before initial use and after each subsequent use.

WARNING: During all work on contaminated medical devices, observe the occupational safety guidelines.

CAUTION: When preparing the cleaning solution, follow the manufacturer's instructions for dilution to proper concentration and temperature.

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CAUTION: Country-specific laws and regulations must be observed.

8.3 Disassembly

The light cable must be removed for reprocessing the TIPCAM®1 S 3D. To do this, unscrew the light cable in a counterclockwise direction. The video connecting cable is firmly connected to the TIPCAM®1 S 3D and must not be removed for reprocessing.

8.4 Pre-cleaning

1. Immediately after a procedure, at the point of use (e.g. operating room), wipe the scope with a soft, lint-free, disposable cloth moistened with distilled/demineralized or enzymatic cleaning solution to remove all heavy soiling.
2. Contain and transport the scope to the decontamination area for cleaning.
3. If immediate cleaning is not feasible, place the scope in a container and soak for no longer than 1 hour in a neutral/mild pH enzymatic cleaning solution (e.g. Enzol), diluted with water to proper concentration per manufacturer's instructions.

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8.5 Manual cleaning

1. Thoroughly rinse the scope with distilled/demineralized water at 20 °C (68 °F), using a cleaning pistol. Rinse the scope head for 1.5 minutes and the silicone cable/card edge connector for 30 seconds, for a total rinse of 2 minutes.
2. Completely immerse the scope in a mild/neutral pH enzymatic cleaning solution (e.g. Enzol) at 21 °C (70 °F). Keep immersed for a minimum of 5 minutes.
3. While immersed, using a soft bristle brush (P/N: 27652), brush all surfaces of the scope, including the silicone cable, card edge connector and operating buttons.

i NOTE: KARL STORZ cleaning brushes are provided non-sterile, single-use and must be discarded after each use. Additional brushes may be purchased from KARL STORZ.

4. After immersion, rinse off the detergent with distilled/demineralized water at 20 °C (68 °F), using a cleaning pistol. Rinse the scope head for 1.5 minutes and the silicone cable/card edge connector for 30 seconds, for a total rinse of 2 minutes.
5. Dry the scope with a lint-free soft cloth or clean, dry, compressed air (<14.5 psi).

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8.6 Inspection after cleaning

The cleaned medical device must be visually inspected to verify that it is clean, complete, damage-free, and dry.

- If soiling or residues remain, the medical device must be cleaned manually and undergo another complete cleaning cycle.
- The device must be free from any visible signs of damage (cuts, dents, or scratches, etc.) that may trap residual debris. If the device shows any damage, contact KARL STORZ for repairs.

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8.7 Sterilization

WARNING: Fully clean the endoscope prior to sterilization.

KARL STORZ recommends to sterilize the TIPCAM®1 S 3D via steam sterilization.

1. Place the endoscope into a TIPCAM®1 S 3D SteriTite container. Follow the container's manufacturer's instructions for details on how to properly set up the container.
2. Load the container into the sterilizer.
3. Use the following cycle parameters for steam sterilization:
 - Temperature: 132 °C (270 °F)
 - Exposure Time: 4 minutes
 - Drying Time: 20 minutes
4. Following cycle completion, crack open the sterilizer door with the container still within the sterilization chamber. Allow the load to cool for 15 minutes with the door cracked open to eliminate any condensation.
5. Remove the container from the sterilization chamber and allow to cool on a stainless steel wire rack for an additional 25 minutes.

8.8 Limits of reprocessing

The end of the product's lifespan is decisively determined by wear, reprocessing methods, employed chemicals, and any use-related damage.

Always perform a pre-inspection functional test with the full 3D system prior to using the device. The system should not be used if any of the accessories are damaged or non-functional.

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9.1 Maintenance and safety check

9.1.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be requested from your local representative or from the manufacturer.

9.1.2 Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once a year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

Visual inspection:

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that the inscriptions relevant to safety are legible.

Electric measurements:

- Protective grounding measured according to IEC 62353 (current edition), and meets the limits specified in the standard.
- Touch Current, Earth Leakage Current, and Patient Leakage currents, measured according to IEC 62353 (current edition), and meets the limits specified in the standard.

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Test for proper functioning:

- A functional check must be carried out in accordance with the instruction manual (see Chapter 6.3.12 of the IMAGE1 S™ CCU instruction manual no. 96206286 DER).
- The safety check must be entered in the unit log-book and the results of the check must be documented.
- If the device does not function reliably and/or safely, it must be repaired.

9.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

CAUTION: Do not attempt to open the camera head or lens. Opening these sealed assemblies affects the camera's soakability and voids all warranties. Do not remove the covers on the camera processor. Warranty is voided unless removed by authorized personnel.

9.3 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this unit.

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9.4 Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In the case of repairs, please contact your local KARL STORZ subsidiary or authorized dealer.

In Germany you can refer repairs direct to:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off-Gewerbepark 83
78579 Neuhausen

Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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9.5 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

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The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

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- Type CF defibrillation protection when used in conjunction with a KARL STORZ camera control unit with a CF-defib symbol on the connection port. 

Storage:

- Ambient temperature:
-10 °C (14 °F) ... +60 °C (140 °F)
- Humidity: 5 ... 95 % (non-condensing)
- Atmospheric pressure: 500 ... 1060 hPa

10.1.1 Directive compliance (for TIPCAM®1 S 3D)

According to Medical Device Directive (MDD): This medical product belongs to Class IIa

This medical product bears the CE mark in accordance with MDD 93/42/EEC.

i NOTE: The code number after the CE mark indicates the responsible notified body.

According to IEC 60601-1-2: Please observe the information on electromagnetic compatibility in chapter 12 (pages 43-58).

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11 Recommended accessories

39501XTC*

Tray for cleaning, sterilization and storage of the TIPCAM®1 S 3D and a fiberoptic light cable.

External dimensions (w x d x h):
640 x 220 x 87 mm

495TIP

Fiberoptic light cable, diameter 4,8 mm, length 300 cm, resistant to very high temperatures, for use with TIPCAM®1 S 3D

495NCSC*

Fiber optic light cable, extremely heat-resistant, safety lock, diameter 4,8 mm, length 250 cm, for use with TIPCAM®1 S 3D

495-3D

Light port, spare part for the screw base on the light port, to be used with TIPCAM®1 S 3D (26605AA, 26605BA)

* not available in the US

For recommended accessories not offered in the United States or Canada please contact Technical Support at 800-421-0837

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
correct storage in instrument tray 39501 XTC

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Information on electromagnetic compatibility (EMC)

12 Information on electromagnetic compatibility (EMC)

The TIPCAM®1 S 3D together with the IMAGE1 S™ is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment. Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the TIPCAM®1 S 3D, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this equipment make it suitable for use in industrial areas as well as in hospitals (CISPR 11 Class A) and other professional healthcare environment. If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 – Compliance level for Immunity tests Guidance and manufacturer's declaration – electromagnetic immunity			
The TIPCAM®1 S 3D is intended for use in the electromagnetic environment specified below. The user of TIPCAM®1 S 3D should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0 % for 1 cycles @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0 % for 1 cycles @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the TIPCAM®1 S 3D require continued operation during power mains interruptions, it is recommended that the TIPCAM®1 S 3D be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the TIPCAM®1 S 3D further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	

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Table 3 – Test levels for Radiated and conducted Immunity Tests
Guidance and manufacturer's declaration – electromagnetic immunity

The TIPCAM®1 S 3D is intended for use in the electromagnetic environment specified below. The user of TIPCAM®1 S 3D should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the TIPCAM®1 S 3D, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b .
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the TIPCAM®1 S 3D is used exceeds the applicable RF compliance level above, the TIPCAM®1 S 3D should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the TIPCAM®1 S 3D.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 – Emission class and group		
Guidelines and manufacturer's declaration – electromagnetic emissions		
The TIPCAM®1 S 3D is intended for use in the electromagnetic environment specified below. The customer or user of the TIPCAM®1 S 3D should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The TIPCAM®1 S 3D uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The TIPCAM®1 S 3D is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5

Recommended separation distances between portable and mobile RF communications equipment and the TIPCAM®1 S 3D

The TIPCAM®1 S 3D is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the TIPCAM®1 S 3D can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TIPCAM®1 S 3D as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency surgical equipment. It has been tested according to IEC 60601-2-2 Annex BB.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzucuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Av. Ejercito Nacional No. 453 Piso 2,
Colonia Granada, Alcaldia Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA – Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamveien1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se
KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestučio st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk
KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopia Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.com
KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com
KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info-hrv@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Polska Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
610044, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

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KARL STORZ Endoscopy New Zealand Ltd.
31 Morningside Drive Mt Albert
Auckland, 1025, New Zealand
PO Box 56 511, Dominion Rd
Auckland, 1446, New Zealand
Phone: +64 9 846 6044, Fax: +64 9 846 6808
Toll free: +64 508 84 84 84 (New Zealand only)
E-Mail: sales@karlstorz.co.nz

KARL STORZ Endoscopy Australia Pty. Ltd.
68 Waterloo Road, Macquarie Park NSW 2113
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: karlstorz@karlstorz.com.au
www.karlstorz.com

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 7461 708-0
Telefax: +49 7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com

STORZ

KARL STORZ — ENDOSKOPE

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SILVER SCOPE® Video-Gastroskope,
Modelle 1382xxx

INSTRUCTION MANUAL

SILVER SCOPE® Video gastroscopes,
Models 1382xxx

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This manual is intended to serve as an aid in the proper use, cleaning, disinfection and sterilization of videoscopes manufactured by KARL STORZ. The videoscopes are sensitive precision instruments, which require particular care during preparation. To prevent damage during preparation, guidance is provided here for the disassembly, cleaning, sterilization and assembly of KARL STORZ videoscopes. For this reason, read these instructions carefully and keep them as a ready reference.

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It is recommended to check the suitability of the product for the intended procedure prior to use.

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1	Symbols employed	3
2	Images Gastroscope	6
3	Images Accessories	8
4	Safety instructions	10
5	Intended use	10
6	User qualification	11
7	Warnings and cautions	11
8	Safety precautions when using the videoscopes	12
9	Handling	17
9.1	Compatible video processor units/light sources	17
9.2	Connecting to a light source, for example the KARL STORZ XENON 100 SCB	18
9.3	Operating the videoscope	19
9.3.1	Preparation of the light source	19
9.3.2	Preparation of the processor unit	19
9.3.3	Connecting the light conductor, air/water connections and the suction channel	19
9.3.4	Deflection of the distal tip	20
9.3.5	Suction	20
9.3.6	Insufflation and light	21
9.3.7	Water irrigation for lens cleaning	21
9.4	Preparation before use	21
9.4.1	Testing the various functions	22
9.4.2	Handling the endoscope correctly	25
9.4.3	Additional recommendations for use	26
9.4.4	Operating the endoscope buttons	27
10	Reprocessing	28
10.1	General information	28
10.2	General cleaning	30
10.2.1	Preparation for cleaning	31
10.2.2	Leakage test	32
10.3	Channel system of a KARL STORZ SILVER SCOPE® video gastroscope	35
10.4	Channel system of a KARL STORZ SILVER SCOPE® Interventional video gastroscope 13826 PKS/NKS	36

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10.5	Cleaning/disinfection of the videoscope: manual cleaning.....	37
10.6	Manual disinfection	42
10.7	Machine preparation of the videoscope	44
10.8	Preliminary cleaning directly after the examination (see page 31).....	50
10.9	Cleaning in the preparation place.....	50
10.10	Preparation in the washer/disinfector (AER)	51
10.11	Cleaning of accessories.....	51
10.12	Preparation of water bottle 13992 BS	52
10.13	Checking preparation procedure.....	52
10.14	Sterilization	53
10.14.1	Ethylene oxide sterilization (EO)	54
10.15	Limits of reprocessing.....	54
10.16	Storage.....	54
11	Storage/transport conditions:	55
12	Disposal.....	55
13	Maintenance	56
14	Important information	56
15	Servicing and repair	57
16	Limitation of liability	57
17	Warranty	58
18	Directive compliance.....	58
19	Appendix: Troubleshootin	59

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1 Symbols employed

Symbol	Meaning
	Consult instructions for use
	Serial number
	Manufacturer
	Fragile, handle with care
	Keep dry
	Temperature limit
	Keep away from sunlight
	Do not reuse
	Applied part type BF, defibrillation proof, IEC 60601-1
	Potential equalization connector
	DANGER: Risk of explosion if used in the presence of flammable anesthetics.
	Protective earth (ground)
	Alternating current

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Symbol	Meaning
	Device is subject to the requirements of the WEEE Directive, 2012/19/EU.
	CAUTION: To reduce the risk of electrical shock, do not remove cover. Refer servicing to qualified service personnel. Do not store liquids on or above the unit. Keep out of reach of patients.
	On
○	Off
	Video input
	Video output
	Color video camera
	Color video monitor
	White balance
	Output

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
Symbol	Meaning
	Input
	Setup Wizard
	Signifies that this feature/option is available in the Options Menu
	Signifies that this feature/option can be controlled via the User Menu
	Signifies that this feature/option can be controlled via the buttons on the videoscope
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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


Images
SILVER SCOPE® Series
Gastrosopes

- ① Up/down deflection wheels
- ② Up/down deflection brake
- ③ Programmable remote control buttons
- ④ Suction valve
- ⑤ Air/water valve
- ⑥ Rubber lip valve/biopsy valve
- ⑦ Right/left deflection wheels
- ⑧ Right/left deflection brake
- ⑨ Nozzle for lens cleaning/irrigation/insufflation
- ⑩ Light cable
- ⑪ Lens
- ⑫ Supply tube
- ⑬ Connection for leakage test
- ⑭ Connection for water bottle
- ⑮ Light cable connection
- ⑯ Port for video processor/electronic connector
- ⑰ Cap with connection for leakage tester
- ⑱ Suction/aspiration port
- ⑲ Air pin
- ㉑ Lanyard
- ㉑ Working channel
- ㉒ Additional irrigation channel with LUER lock connection (only 13826 PKS/NKS)
- ㉓ Distal water-jet channel (only 13826 PKS/NKS)

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4 Safety instructions

Please read these safety instructions carefully before using the instrument. Before using the instrument on the patient for the first time, it is imperative that you be acquainted with how the instrument operates and is controlled.

5 Intended use

Gastrosopes are used for medical, diagnostic and endoscopic treatment in the gastrointestinal tract.

Gastrosopes are designed for transient use in endoscopic interventions performed through natural orifices.

Indications

The use of gastrosopes is indicated if, in the opinion of the attending physician, a gastroscopy is indicated.

This includes, amongst other things:

- Inexplicable abdominal pain
- Diagnostics
- Acid reflux disease
- Tumor removal, observation after tumor removal
- Treatment of polyps

Contraindications:

The use of gastrosopes is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated, or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Gastrosopes must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

Videoscopes are applied parts, which should only be used with the designated supply units approved by KARL STORZ (camera control units, light sources).

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6 User qualification

KARL STORZ videoscopes may only be used by persons with an appropriate medical qualification and who are familiar with the endoscopic technique. The information given in these instructions only serves to instruct in the correct handling and preparation of the videoscopes.

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7 Warnings and cautions

Please read this manual and follow its instructions carefully, to avoid putting your patients, personnel or yourself at risk. The words, **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the instrument. To make the words stand out more clearly, they are accompanied by a pictogram.

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NOTE: A Note indicates special information about operating the instrument, or clarifies important information.

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8 Safety precautions when using the videoscopes

Videoscopes must be used according to the medical rules and procedures of endoscopy recognized for this method.

WARNING: Check the videoscope prior to each procedure to ensure that it functions correctly and has been properly cleaned and disinfected.

WARNING: Verify the sound condition of the videoscope, as well as of accessories used in combination, before every application. Damaged videoscopes or damaged accessories must not be used.

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WARNING: Remove the EO/pressure compensation cap 11025 XE from the connection piece during:

- endoscopic interventions
- manual/machine preparation.

WARNING: Do not allow any electrically active components or instruments to come into contact with the videoscope.

WARNING: The videoscope must not be used during the discharge of a defibrillator. Remove the videoscope from the operation area before discharge occurs.

WARNING: Danger of explosion if flammable gases are used in the immediate proximity of the videoscope. The videoscope must not be used in the presence of flammable or combustible materials.

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WARNING: High-energy light radiated through endoscopes may give rise to high temperatures in front of the light outlet, at the end faces of the fiberoptic light cable and at the distal tip of the endoscope.

Prolonged direct contact of light-emitting surfaces on tissue can lead to burns; this must therefore be avoided.

Temperatures remain high even after the proximal connector is removed from the light source (risk of burns to the user from the light cable).

WARNING: Before removal from the patient, the distal tip of the videoscope must be returned to the neutral, nondeflected position (L and U on top of each other and positioned between the valves) to prevent damage to the instrument and/or injury to the patient.

WARNING: The videoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for videoscopes that are repaired by unauthorized personnel.

WARNING: Damage to the videoscope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.

WARNING: Always unplug the videoscope before carrying out any maintenance and cleaning procedures.

WARNING: Risk of infection due to shipment! Contaminated videoscopes can cause serious or fatal diseases during transportation or at the recipient's premises. Please cover the carrying case with the protective sheet 13990 SFN before sending contaminated instruments back to KARL STORZ.

WARNING: If a malfunction of the videoscope should occur during an intervention on the patient, stop the intervention immediately. Place the distal tip of the videoscope in a neutral, nondeflected position and remove the videoscope slowly and carefully from the patient.

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CAUTION: Place the EO/pressure compensation cap 11025XE on the connection piece prior to:

- gas sterilization
- shipping

CAUTION: Do not expose the videoscope to direct sunlight or excessive heat, to avoid damage to the device.

CAUTION: Do not steam sterilize (autoclave) the videoscopes.

CAUTION: The maximum permissible temperature for decontamination, sterilization, transport and storage is 65 °C.

CAUTION: The 'Preparation and Care' manual must be followed. The procedures for cleaning, disinfection and gas sterilization are explained therein in detail.

CAUTION: Do not hit the distal tip of the videoscope against hard objects, and never use a sharp object to remove dirt.

CAUTION: Never lay heavy objects on the videoscope.

CAUTION: Do not use any accessories with sharp edges or points in or on the endoscope.

NOTE: The patient leakage current limit of the CF application part may be exceeded when CF application parts (application on heart) are used in combination with other BF application parts. In this case, the reduction of the patient leakage current should be effected using only CF application parts whenever this is possible.

NOTE: The endoscope does not offer any electric insulation against high frequency (HF) currents. Protection against HF currents should be ensured through the use of the correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.

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NOTE: Any damage to the videoscope resulting from incorrect operation is not covered by the warranty.

NOTE: Constant changing between various preparation methods is especially hard on the materials and must therefore be avoided.

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NOTE: The KARL STORZ videoscope should not be used in combination with video processor units and light sources from other manufacturers. The use of non-compatible units can lead to patient injury or unit damage. It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other units at any time.

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
**9.2 Connecting to a light source,
for example the KARL STORZ
XENON 100 SCB**

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CAUTION: The correct videoscope adaptor must be used for the video processor.

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Connection of supply connector with light source and water bottle

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1. Install the light source on a solid, even surface.
2. Insert the power cable into the socket provided.

9.3.2 Preparation of the processor unit

Refer to the relevant instruction manual.

9.3.3 Connecting the light conductor, air/water connections and the suction channel

1. Insert the supply connector into the light source until it audibly clicks into place (Fig. 1).
2. Connect the videoscope adaptor cable (VEA) to the endoscope. When doing so, take note of the white markings on the endoscope and the VEA cable (Fig. 2).
3. Attach the connector from the water bottle tube to the endoscope supply plug (irrigation bottle connector) (Fig. 3). The water bottle holder automatically positions itself correctly and cannot be adjusted manually.
IMPORTANT: To fill the water bottle with sterile water, remove the bottle from the light source and unscrew the lid. Only fill the water bottle 2/3 full with sterile water (1/3 is required to cushion air pressure).
4. Connect the suction tube (patient tube) from the suction pump or central suction system to the supply plug connector (Fig. 4).

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9.3.4 Deflection of the distal tip

Rotating the large deflection wheel, labeled 'U/D', (up/down) deflects the distal tip. Rotating clockwise (or forwards) lowers the tip. Rotating counterclockwise (or backwards) raises the tip.

Right and left movements are controlled using the smaller deflection wheel (labeled 'L/R', left/right). Rotating clockwise (or forwards) = right / rotating counterclockwise (or backwards) = left (Fig. 5).

WARNING: Do not rotate the wheels forcibly, this could damage the instrument or cause injury to the patient! Before attempting to introduce or remove the instrument, release the deflection wheel brakes!

9.3.5 Suction

Switch on the suction pump. To perform suction, press the suction valve (marked red) on the endoscope (Fig. 6a).

Suction can also be used when the biopsy forceps are placed in the biopsy channel and the rubber lip valve is mounted. During the procedure, the rubber lip valve/biopsy valve seals off the forceps or the accessories.

With the interventional gastroscope (13826 PKS/NKS), the suction channel can also be determined after the suction valve with the two lugs has been positioned in the direction of the distal end.

Should the pressed-down valve be aligned such that the mark points towards the distal end, suction will be performed via both working channels by pressing the suction valve (Fig. 6b).

Should the pressed-down valve be turned anticlockwise, suction will be performed via the left working channel (dia. 2.8 mm) by pressing the suction valve (Fig. 6c).

Should the pressed-down valve be turned clockwise, suction will be performed via the right working channel (dia. 3.8 mm) by pressing the suction valve (Fig. 6d).

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9.4.1 Testing the various functions

1. Prepare light source as described in 9.3.1.
2. Connect supply connector, water bottle and suction tube according to 9.3.3.
3. Check the supply cable for damage (e.g., breaks, cracks, twisted or squashed areas...) and the supply connectors for loose connections.
4. Examine the surface of the flexible insertion shaft visually and by touch for cracks, cuts, bulges, foreign bodies or other damage.

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The large deflection wheel controls up and down movement. To test up and down movement, release the up and down brake. Slowly rotate the large deflection wheel clockwise (or forwards), the distal end should move downwards. Rotating the deflection wheel counterclockwise (or backwards) causes upward movement.

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6. Testing deflection – right/left control
7. The smaller/outer deflection wheel controls right/left movement. To test right/left movement, release the right/left brake. Slowly rotate the small deflection wheel clockwise or forwards, the distal tip should move to the right. Rotating the deflection wheel counterclockwise (or backwards) causes movement to the left.
8. Testing deflection – brake function
The deflection wheels should turn freely when the brakes are released.
9. Release the brakes and test their action (see 5a/b). Activate the brakes by moving the brake lever counterclockwise – ensure that the deflected position of the distal endpiece is fixed in the desired position.

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WARNING: If the deflection and/or brake mechanisms do not function properly, do not use the device. Contact the KARL STORZ service department.

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11. Check the function of the water supply:
Depress the air/water valve (blue ring) fully.
Water should flow out of the water nozzle (distal) and rinse the lens clean.
IMPORTANT: When you first actuate the valve, it may take a couple of seconds for water to start flowing.

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12. Test the suction function. With the unit correctly installed and the suction pump or central suction system switched on, immerse the distal end in a container of sterile water once again. Depress the suction valve (red ring) fully. Water should be suctioned through the suction channel of the instrument via the connected suction pump and into the collecting bottle of the suction equipment.

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CAUTION: Full (>180°) deflection can lead to limitation in movement, depending on the accessories. If necessary, adjust the deflection slightly.

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Handling the endoscope correctly

Hold the control handle of the endoscope in the left hand in such a way that the thumb, index and middle fingers can easily control the large deflection wheel. In this position, the suction valve can be depressed with either the index or the ring finger. The air/water valve can also be covered with these fingers. The remote control button (e.g., Stills) is also within reach.

The small, outer deflection wheel can be controlled with the fingers of the left hand, or alternatively with the thumb and fingers of the right hand.

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9.4.3 Additional recommendations for use

1. Use a lubricant from KARL STORZ (see list on our website www.karlstorz.com under 'Hygiene') for the distal end of the insertion tube.
2. Observe the lumen the whole time. If you become disoriented, withdraw the instrument gently until vision is restored.

NOTE: Do not put lubricant on the lens to avoid clouding!

NOTE: Do not use oil or petroleum, silicone, Vaseline or alcohol based media, to avoid damage to the endoscope material.

NOTE: If debris causes the lens to cloud during the procedure, this can normally be rectified by a combination of water and insufflation. Excessive water or air can be removed by applying simultaneous pressure on the suction valve.

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9.4.4 Operating the endoscope buttons

Depending on the video processor used, the buttons can have two functions, which are initiated by means of short and long presses of the button. For more detailed information on button assignment, please consult the manual for the camera control unit used. A list of compatible camera control units can be found in section 9.1.

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NOTE: In case of 'red-out', a blurred red picture on the screen, caused by the distal endoscope tip resting directly on the mucosa, withdraw the instrument gently until vision is restored.

WARNING: Pushing the instrument forwards despite a 'red-out' effect may result in serious injury to the patient.

NOTE: If the angling of the distal tip can no longer be properly controlled the examination must be ended immediately. Release the brakes and place the tip in the neutral position. Position the deflection wheels so that 'U' and 'L' are between the valves (see page 13/15), then carefully withdraw the instrument, whilst continuing to monitor the procedure.

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WARNING: Remember that decontamination is only then successful if

- the instruments have been freed of all organic material and chemical residues.
- work proceeds according to recognized and validated methods.
- any deviation from recommended procedures are validated and tested by the user.
- the recommended sterilization parameters are maintained. These are only valid with cleaning, disinfection and sterilization equipment that is properly maintained and calibrated.

CAUTION: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' must be followed. The cleaning, disinfection and sterilization procedures are explained therein in detail.

CAUTION: When preparing and using the cleaning and disinfection solutions, follow the manufacturer's instructions paying close attention to proper concentration and exposure times. Disregarding the instructions can lead to damage. Bear in mind the microbiological range of action of the chemicals used.

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CAUTION: Constantly changing between different preparation methods causes particular wear on the materials and must therefore be avoided. However, a once-only change in methods, e.g., when acquiring new cleaning equipment, should be carried out in consultation with the manufacturers of the endoscope, the chemicals and the washer/disinfector.

This change in method is then usually completely harmless.

CAUTION: Do not clean videoscopes in an ultrasound bath.

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CAUTION: During manual preparation, videoscopes may only be rinsed with microbiologically clean/sterile water.

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10.2 General cleaning

Cleaning should only be performed by authorized, trained personnel. All safety regulations should be observed (gloves, protective glasses, mouth and nose protection, clothing, German Technical Rules for Biological Agents (TRBA 250).

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The national guidelines and recommendations for the preparation of instruments (for example RKI, ESGE/ESGENA, ASGE, LANGS, SAGES etc.) should be observed and the effectiveness verified.

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CAUTION: Do not autoclave the endoscope or accessories (with the exception of accessories that are specifically designated as autoclavable, e.g., forceps).

NOTE: The endoscope must be cleaned immediately after every use. If organic material is allowed to remain and dry in the instrument, complete or partial blockage of the channel may result. Thorough cleaning is a prerequisite for effective preparation. This will avoid coagulation residue / blockages and the resulting repairs and instrument failures.

NOTE: Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.

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10.2.1 Preparation for cleaning

1. The preliminary cleaning of the endoscope is performed at the place of operation. The main cleaning and disinfection is carried out in a separate preparation room.
2. Immediately after use, remove any heavy soiling (secretions, etc.) and organic material from the outer surfaces of the insertion tube with a disposable cloth moistened with cleaning solution. Then submerge the distal end in a container with cleaning solution and activate the suction to rinse the working channel.
3. Actuate the suction valve and suction at least 100 ml of cleaning solution through the suction channel. (Fig. 14). Suction up at least 200 ml with the interventional gastroscope (13826 PKS/NKS). The suction valve must then be operated in all three positions as described in section 9.3.5 (Left --> left working channel dia. 2.8 mm, right --> right working channel dia. 3.8 mm and middle --> through both channels). Make sure that all the channels are free from blockages.
4. Use the cleaning valve 13991 SRV (Fig. 13) in place of the air/water valve. (Do not use this valve during examination!).
5. Check that the water bottle is filled (max 2/3). Start the air supply (insufflation pump). Water and air are now passed through both the air channel and the water channel, depending on the position of the cleaning valve. (The process should run for at least 10 seconds).
Air channel = Position 1
Water channel = Position 2
6. It is important that both valve settings ('Actuated'/'Not actuated') are used multiple times for several seconds. This is the only way to ensure that both channels are rinsed through.
7. Turn off video system.

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8. Pull out the video connector from the videoscope adaptor (VEA) at the supply connector of the endoscope when the video processor unit is switched off. Put on the sealing cap for the video cable connector and lock it with a 1/4 turn.
9. Disconnect the endoscope from the video processor unit and light source. Place the cleaning valve in a cleaning solution.
10. Then transport in a protected case to the site of reprocessing.

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10.2.2 Leakage test

Preparing for the leakage test – dry

- i NOTE:** Place the cap on the videoscope adaptor (VEA) connector and screw it on tightly. Ensure that the cap is secured properly (Fig. 16).

1. The connection tube of the leakage tester and the vent port on the endoscope must be dry. Connect the leakage tester connector to the supply connector and lock by turning clockwise 90° (Fig. 17).

- i NOTE:** The rubber seal in the connection adaptor of the leakage tester 13242 LX should be treated occasionally with instrument oil.

Before every immersion in liquid, test the endoscope for leaks. Remove all valves.

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1. Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal deflection part of the endoscope expands slightly. Move the deflection part with the deflection wheels on the control piece in all four directions as far as possible and check the manometer to see whether the pressure remains constant or falls (dry test).
2. Wait at least 1 minute.

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Performing the leakage test – in liquid

1. Immerse the device in water so that it is completely covered. It must be ensured that the leakage tester is not immersed but rather kept dry. Perform the leakage test as described above, note whether bubbles emerge from the device. If after an observation period of at least 3 minutes no more bubbles can be seen emerging from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof and can be used.

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If this test is also negative and no bubbles can be seen rising in the liquid, the leakage tester can be removed.

To do this, remove the supply connector from the liquid and press the vent button on the manometer to vent the endoscope.

NOTE: Never release the pressure in the liquid. For this, remove the connector from the liquid and press the vent button. The test pressure must fall to 0. Then remove the leakage tester.

If a pressure loss occurs, this means that the videoscope leaks and may no longer be used, immersed or reprocessed by machine.

CAUTION: In this case, fill out the accompanying repairs form (see Appendix), roughly clean and remove any organic soiling/residue from the endoscope, mark it as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the enclosed transport protection sheet 13990 SFN.

WARNING: Risk of infection during shipment! Contaminated instruments can cause aerogenous or contact infections during transportation or at the recipient's premises. Shipping without the plastic film packaging must be avoided.

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**10.3 Kanalsystem eines
KARL STORZ SILVER SCOPE®
Video-Gastroskops**


(13820 PKS/NKS, 13821 PKS/NKS,
13823 PKS/NKS)

**10.3 Channel system of a
KARL STORZ SILVER SCOPE®
video gastroscope**

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
10.4 Channel system of a
KARL STORZ SILVER
SCOPE® Interventional video
gastroscope 13826 PKS/NKS

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10.5 Cleaning/disinfection of the videoscope: manual cleaning

Please adhere to the national guidelines, e.g., ESGENA guidelines, DEGEA/RKI recommendations, ASGE, LANGS, SAGES etc. for cleaning, disinfection and sterilization, as well as to the instructions of the chemical manufacturer concerning concentration and exposure times.

1. Remove the rubber lip valve/biopsy valve, suction valve and air/water valve and place in a cleaning solution and clean using a soft disinfected brush (13991 VB) or disposable brush.

2. Wipe the housing, insertion shaft, and supply tube carefully with a lint-free cloth soaked in a cleaner (for suitable cleaning agents or disinfectants approved by KARL STORZ, see www.karlstorz.com). Avoid shocks and pulling movements so as not to damage the insertion shaft. Then rinse with microbiologically pure/sterile water.

3. Immerse the endoscope in the cleaning solution and rinse the channels using a syringe and the supplied irrigation tubes -13820 PKS/NKS, 13821 PKS/NKS, 13823 PKS/NKS (Fig. 23a, 23b): 13991 SS, 13991 YS; -13826 PKS/NKS (Fig. 24a, 24b): 13991 DKSA, 13991 SWJ, 13991 JASA).
4. The channels must be filled without any bubbles, to guarantee a complete covering of the surfaces.

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
- When submerged, use a suitable brush (preferably disposable) (see accessories, page 9/10) to brush below the surface of the liquid in one direction starting from the supply plug (supply plug --> control body, control body --> distal tip, biopsy channel --> distal end) (Fig. 25a: 13820 PKS/NKS, 13821 PKS/NKS, 13823 PKS/NKS; Fig. 25b and 25c: 13826 PKS/NKS)). Repeat these procedures at least 3x.

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- Clean the brush head (if you are using reusable brushes) at the distal end before retracting it. Do not brush backward and forward, rather brush in one direction only.

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

7. Thoroughly clean out the entire suction channel (shaft part and supply part) and the biopsy channel a number of times with a brush. Observe the exposure time for the cleaning solution.

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
8. Carefully clean all particles of dirt from the distal tip. It may be necessary to clean the air/water valve with a brush (13991 VB) with soft bristles, in order to remove hardened contaminants. Never remove debris with a needle or other sharp objects, because these can cause damage.

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

9. Connect the irrigation tubes:
- 13820 PKS/NKS, 13821 PKS/NKS, 13823 PKS/NKS: 13991 SS, 13991 YS;
- 13826 PKS/NKS: 13991 DKSA, 13991 SWJ, 13991 JASA) and pump fresh cleaning fluid through the channels using a syringe. Afterwards, rinse with microbiologically clean/sterile water.

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Suitable brushes

Art. no.	Designation	Used for
110915-01	Cleaning brush dia. 3 mm	Working channel 2 mm
110910-01	Cleaning brush dia. 5 mm	Working channel from 2.8 to 4.2 mm
13991 VB	Valve brush	Valve seat
110917-01	Cleaning brush dia. 1.5 mm	Water-jet channel 13826 PKS/NKS

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10. Remove the device from the cleaning solution and rinse with microbiologically pure/sterile water.

11. Carefully clean all valves with the disinfected cleaning brush (13991 VB) and then rinse.

12. The valves should also be cleaned in an ultrasound bath.

13. Flow control for the interventional gastroscope 13826 PKS/NKS:

- Remove the irrigation tubes from the device
- Close the working channel inlets on the handle using a biopsy valve/rubber lip valve ((Fig. 30a) 13991 DA).
- Connect the irrigation adaptor 13991 CS to the device's supply sheath (Fig. 30b)
- Connect the suction cleaning valve ((Fig. 30c) 13991 ARV) in the valve seat of the suction valve on the handle in two subsequent positions so that it engages perceptibly in place each time.

CAUTION: The valve 13991ARV must be set in the direction of the right working channel. The marker on the valve then points towards the right working channel (Fig. 33d). Then valve must then be set in the direction of the left working channel. The marker on the valve points towards the left working channel (Fig. 33e). At least 50 ml of microbiologically clean / sterile water must be pumped through the irrigation tube in both valve settings.

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It must be checked that the working channels are not blocked and that the amount of liquid pumped through emerges again at the distal tip.

Once the cleaning work has been completed, clean the cleaning instruments, for example the brushes and irrigation tubes and disinfect/sterilize them. Never reuse damaged or kinked brushes.


CAUTION: Do not autoclave the endoscope or its accessories (with the exception of accessories that are specifically designated as autoclavable, e.g., forceps).

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
10.6 Manual disinfection

1. After careful rinsing, fully immerse the endoscope in the disinfecting solution.

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CAUTION: When preparing and using the solutions, follow the chemical manufacturers' instructions paying close attention to proper concentration and exposure times. Prolonged immersion can lead to damage. Bear in mind the microbiologic range of action of the chemicals used.

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2. Using the irrigation tubes and a syringe with LUER connector, rinse out the channels a number of times with disinfecting solution.

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3. Then rinse the channels and surfaces of the video gastroscope with microbiologically pure/sterile water.

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4. After the final step in the cleaning process has been completed, dry the entire endoscope. Thoroughly dry off all external surfaces of the endoscope, preferably with sterile compressed air at a max. of 0.5 bar/7.25 psi, thereby removing any residues of cleaning and disinfection solution.

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If at all possible, store the videoscope in a dry cupboard or hang it up in a storage case for safe storage.

NOTE: The valves of the videoscope must not be inserted; they must be stored separately.

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10.7 Machine preparation of the videoscope

All flexible KARL STORZ endoscopes can be connected to approved washer/disinfectors for semi- or fully automatic cleaning and disinfection. Observe the national laws and regulations for cleaning and disinfection procedures. Leakage testing, manual cleaning and (where necessary) rinsing with microbiologically pure water must be carried out beforehand. The manufacturer of the cleaning machine must be consulted when selecting the preparation programs. Usually, special programs for flexible endoscopes are available which are suited for decontamination and compatible with the materials. For machine preparation, make sure that the flexible endoscopes are not heated above 65°C.

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CAUTION: Ensure that only KARL STORZ approved agents are used for procedures involving peracetic acid.

Lay the endoscope carefully in the tray or basket of the machine. Connect the irrigation tubes supplied with the cleaning adaptor set 13991 CK or the irrigation tubes of the manufacturer of the washer/disinfectant to the endoscope.

Cleaning adaptor set 13991 CK

The cleaning adaptor set 13991 CK is used for validated preparation of KARL STORZ videoscopes from series 13xxx NKS/PKS (acc. to ISO 15883 part 4). The cleaning adaptor set is not intended for any other purpose as this may damage it.

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

- 1** {
- 1 Connection tube 13991 CL for leakage tester connection ⑤ of the videoscope
 - 2 Connector (13991 CL) for vent port
 - 3 LUER connector (13991 CL) for washers/disinfectors for endoscopes

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- 4** {
- 4 Cleaning adaptor 13991 CR for air/water (A/W) channel ⑦ of the videoscope
 - 5 Connector (13991 CL) for A/W channel (air pin ②)
 - 6 Connector (13991 CL) for water channel (water bottle connection ⑦)
 - 7 2x LUER connectors (13991 CR) for washers/disinfectors (1x LUER for air channel, 1 LUER for water channel)

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- 8** {
- 8 Cleaning adaptor 13991 CS for suction channel (aspiration) ② of the videoscope
 - 9 Connector (13991 CL) for suction channel on the supply connector ②
 - 10 LUER connector (13991 CS) for washers/disinfectors

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- 11** {
- 11 Connection adaptor 13991 SS for biopsy/working channel ⑨ of the videoscope (with 2x double channel endoscopes)
 - 12 Connector (13991 CL) for biopsy/working channel ⑨
 - 13 LUER connector (13991 SS) for washers/disinfectors

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- 14** {
- 14 Channel separator 13991 CD for suction valve ④ and irrigation/insufflation valve (A/W) ⑤ of the videoscope
 - 15 Valve pin for irrigation/insufflation valve (A/W) ⑤
 - 16 Closure for suction valve ④

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WARNING: This product is not sterile when delivered and must therefore be cleaned, disinfected and if necessary, sterilized prior to initial use and each subsequent reuse.

WARNING: Risk of infection by improper preparation: Preparation may only be carried out by experienced and trained personnel.

CAUTION: Only use chemicals which have been tested for material compatibility by KARL STORZ. If you use other chemicals, get written confirmation from the chemical manufacturers that the material is compatible with KARL STORZ endoscopes.

CAUTION: If the preparation chemicals are changed, it is essential that you consult the manufacturer of the chemicals and of the washer/disinfector.

CAUTION: Please refer to our operating instructions for details on the cleaning preparation under the section: Cleaning, disinfection and sterilization.

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The following must be checked before use:


- visual checks of the connection tubes (patency) and the O-ring (wear)
- check the connections are correctly and securely in place
- check for damage

In the event of noticeable problems or complaints, please contact our service department or send in the cleaning adaptor set for technical examination.

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Plug the cleaning adaptor 13991 CR ⑤ as far as it can go into the connector ② of the air pin of the videoscope.

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Then plug connector ⑥ as far as it will go into the 'Connection for water bottle' ⑰ of the videoscope.

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Plug connector ⑨ of the cleaning adaptor 13991 CS as far as it will go into the suction channel ㉑ at the supply connector of the videoscope.

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If necessary, remove the suction and irrigation/insufflation valves (4, 5) from the videoscope.

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Insert the channel separator 13991 CD (14) with the valve pin (15) into the free opening of the irrigation/insufflation valve (4) as far as it will go so that the lock (16) lies on top of the free opening of the suction valve (4).

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
Press the slide mechanism of the channel separator 13991 CD (14) downward and slide it forward until the channel separator 13991 CD (14) audibly locks into place. Check that the channel separator is securely in place.

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! NOTE: If necessary, remove the cap of the biopsy/working channel (6) of the videoscope.


Connect the connection adaptor 13991 SS with the connector (12) to the biopsy/working channel (6) of the videoscope.

! WARNING: Danger of infection due to inadequate rinsing caused by kinked tubes (decreases cleaning performance): avoid kinking of the irrigation tubes.

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
! CAUTION: When the program has finished, the videoscopes are at a risk of corrosion due to the conditions in the preparation machine.

! CAUTION: Immediately after the end of the program cycle, the videoscopes are still hot, as they have been heated up during chemical thermal preparation. Consequently the instrument's material properties are different and, as such, the prepared endoscope should be allowed to cool down for at least 20 minutes at room temperature before being used again.

! CAUTION: Ideally, remove the videoscopes from the machine immediately after the program has finished. Then dry off the optical surfaces, if necessary, with a cotton swab soaked with 70% alcohol.

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Store the endoscopes freely suspended, dust-free and dry. Store valves separately

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10.8 Preliminary cleaning directly after the examination (see page 31)

- The outer surfaces should be wiped down from top to bottom with a lint-free disposable cloth that is soaked with cleaning solution.
- Ample amounts of cleaning solution should be suctioned through the distal end.
- Apply the air-water supply for a few seconds every cycle
- Connect the cleaning adaptor
- Disconnect the flexible endoscope from the light source

CAUTION: Close the video connector with the cap.

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Trocknung

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10.9 Cleaning in the preparation place

- Carry out leakage test
- Immerse the flexible endoscope
- Clean the channels with special cleaning brushes
- Clean the complete air/water system thoroughly
- Rinse the channels

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10.10 Preparation in the washer/ disinfecter (AER)

- The machine must fulfill all the cleaning, disinfection and material compatibility requirements (validated process).
- Make sure that the environmental conditions with regard to temperature and pressure are not exceeded (+60 °C (tolerance +/- 5%) and 1 bar for the channels, 0.3 bar leakage test).
- Recommendation: Before placing the endoscope in the machine, clean all of the channels (working channels, water channels, air/water channels) regardless of whether or not they have been used.
- We recommend sticking to one type of cleaning agent.

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10.11 Cleaning of accessories

1. Thorough cleaning is required before disinfection or sterilization. Carefully clean the biopsy forceps and brushes whilst they are immersed in a cleaning solution using a soft, disinfected brush. Take special care to remove residues such as blood, excrement and mucous from the biopsy forceps and needles – an ultrasound bath is recommended.
2. Rinse in microbiologically pure/sterile water and dry (using sterile compressed air).
3. Disinfect accessories (work procedures as with cleaning).
4. Always test instruments after cleaning/disinfection for proper functioning and possible damage.
5. Before sterilization, carefully oil movable parts. Use only suitable, steam-resistant, physiologically safe, residue-free oil (e.g., 27656 B). After oiling, open and close the forceps several times. For details on preparation, please read the instruction leaflet of the respective accessory.

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Replace defective parts or send them in for repair. Before sending them in for repair, all contaminated parts must be cleaned/disinfected and sterilized or clearly labeled as unsterilized.

If residues are still present after cleaning, prepare the parts again, subject them to special treatment, discard or repair.

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Test the instruments for proper operation (deflection of the distal end, brake lever). Check the lens for a clear view, the instrument for appropriate deflection and all channels for patency (see page 9).

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10.14 Sterilization

Sterilization of instruments to be used on patients is not necessary when coming into contact with the skin and mucous membrane, as well as when introduced through natural, unsterile orifices.

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WARNING: Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

WARNING: Risk of infection: Sterilize the medical device in a disassembled state.

CAUTION: Flexible endoscopes must never be steam sterilized (autoclaved).

CAUTION: The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

NOTE: All deviations from the recommended sterilization parameters must be validated by the user.

NOTE: When using reprocessing equipment, please observe the equipment manufacturers' instructions.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

CAUTION: The STERRAD process is not suitable for sterilization of videoscopes. This process damages the videoscope!

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The procedure must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

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10.14.1 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100% ethylene oxide at 55°C and a hold time of 60 minutes.

10.15 Limits of reprocessing

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

10.16 Storage

The storage method must be selected taking into account the respective applicable national requirements.

KARL STORZ only recommends dust-protected, suspended storage for decontaminated endoscopes.

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CAUTION: The sheath of the flexible endoscope must never be kinked or coiled up too tightly.

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14 Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and, if necessary, disinfected/sterilized before they are sent for repair. We reserve the right to return contaminated instruments, units or endoscopes to the sender.

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KARL STORZ gives no warranty on the correct functioning of units, endoscopes or instruments which have been repaired by unauthorized third parties.

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15 Servicing and repair

Defective endoscopes must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts.

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16 Limitation of liability

As a supplier of this endoscope, KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:
All assembly, operations, system expansions, readjustments, modifications, or repairs to same have been performed by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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17 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see “Subsidiaries” section), even during the warranty period.

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18 Directive compliance

According to Medical Device Directive:

This medical product belongs to Class IIa

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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

Notes


Damage report

Hand piece
☐ Function buttons leaking
☐ Chassis leaking

STORZ

Steering wheels
☐ loose
☐ leaking

Supply tube
☐ cracked / damaged

Working channel
☐ blocked
☐ leaking


Insertion sheath
☐ cracked
☐ worn out

Distal tip
☐ Lens scraped
☐ Insufflation disturbed
☐ Image transmission error

Deflection
☐ leaking
☐ blocked / damaged

Supply plug
☐ electrical contacts damaged / corrosive
☐ leaking

Light transmission
☐ damaged
☐ Illumination power not adequate

REPAIR DISPATCH NOTE

From customer enclosed parts:

Customer	Article / Part	KARL STORZ
☐	Valve / air-water, suction valve	☐
☐	Cleaning brush	☐
☐	Protection wrap	☐
☐	Biopsy forceps	☐
☐	Cleaning tubes	☐
☐	Instrument oil	☐
☐	ETO cap	☐
☐	Bite protector	☐
☐	Leakage tester	☐
☐	Seal cap for working channel	☐
☐	Video connection cable	☐

QS KARL STORZ

Signature

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postfach 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0, Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com

KARL STORZ Endoskope Berlin GmbH
Scharnhorststr. 3
10115 Berlin/Germany
Phone: +49 (0)30 30 69090, Fax: +49 (0)30 30 19452

KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kse.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorz.cuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA – Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12739 Skärholmen, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: LT-LV-info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongquiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

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**It is recommended to check the suitability of
the product for the intended procedure prior
to use.**

Thank you for your expression of confidence in
the KARL STORZ brand name. Like all of our
other products, this product is the result of years
of experience and great care in manufacture. You
and your organization have decided in favor of
a modern, high-quality item of equipment from
KARL STORZ.

This instruction manual is intended to serve as an
aid in the proper handling, cleaning, disinfection
and, where required, sterilization of the CMOS
Video-Rhino-Laryngoscope. All essential details
of the equipment and all actions required on your
part are clearly presented and explained. We thus
ask that you read this manual carefully before
proceeding to work with the instrument. Keep this
manual available for ready reference.

WARNING: KARL STORZ instruments
are provided unsterilized, and must be
cleaned, disinfected and/or sterilized
prior to the initial use and before each
subsequent use.

Rx only

*Federal (USA) law restricts this device to
sale by or on the order of a physician.*

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KARL STORZ is constantly working on the further
development of all products. Please appreci-
ate that changes to the scope of supply in form,
equipment and technology are possible for this
reason. Therefore, no claims may be deduced
from the information, figures and descriptions in
this manual.

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1 Important note for users of KARL STORZ devices and instruments	1
2 Pictures of the instruments	4
3 Symbol description	6
4 Safety information	8
5 Intended use	8
5.1 Indications	8
5.2 Contraindications	8
6 User qualification	9
7 Kompatibilität	9
8 Warnings and cautions	10
8.1 Safety precautions when using the CMOS Video-Rhino-Laryngoscopes	10
9 Initial operation	14
9.1 Unpacking the device	14
9.2 Basic equipment	14
9.3 Operating conditions	14
9.4 Electromagnetic compatibility	14
9.5 External system components	15
10 Handling	16
10.1 Operating	16
10.2 Holding correctly	16
10.3 Handling the shaft	17
10.4 Handling the deflection mechanism	18
10.5 Connection to the monitor 8403 ZX	20
10.6 Use of the HDMI output	20
10.7 Functions of the C-MAC® Monitor 8403 ZX	21
10.7.1 Focus	21
10.8 Function buttons of the endoscope	22
10.9 Monitor function keys	22
10.10 Preparation prior to use	23
10.10.1 Inspecting the surfaces	23
10.10.2 Checking the deflection	24
10.11 Switching the CMOS Video-Rhino-Laryngoscope on/off	25

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
10.12	Inspecting the image quality.....	26
10.13	Switching off the C-MAC® Monitor	26
10.14	System reset	27
10.15	Monitor bag 8403 YD.....	27
11	Reprocessing	28
11.1	Overview of reprocessing procedures.....	28
11.2	Accessories.....	30
11.3	Preliminary cleaning directly after use	30
11.4	Leakage test	31
11.4.1	Performing the leakage test – dry.....	31
11.4.2	Performing the leakage test – in liquid	31
11.5	Disassembly.....	32
11.6	Manual reprocessing	32
11.6.1	Manual cleaning	32
11.6.2	Manual disinfection.....	33
11.6.3	Manual High Level Disinfection (HLD)	33
11.6.4	Drying	33
11.7	Automated reprocessing	34
11.7.1	Manual preparation	34
11.7.2	Machine cleaning/chemical-thermal disinfection.....	34
11.8	Assembly, inspection and care.....	34
11.9	Packaging	35
11.10	Sterilization.....	36
11.10.1	Hydrogen peroxide (H ₂ O ₂) sterilization - ASP STERRAD®	36
11.10.2	Hydrogen peroxide (H ₂ O ₂) sterilization – STERIS® AMSCO® V-PRO™	37
11.10.3	Ethylene oxide sterilization (EO).....	37
11.10.4	Chemical low-temperature sterilization with peracetic acid – STERIS® System 1®	37
11.10.5	Low-temperature steam and formaldehyde process (LTSF)	38
11.11	Limits of reprocessing	38

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12 Storage	38
13 System description	40
14 Technical description	41
15 Maintenance	43
15.1 Servicing and repair	43
15.2 Maintenance	43
15.3 Limitation of liability	43
15.4 Warranty	44
15.5 Standard compliance	44
15.6 Directive compliance	45
15.7 Important information	45
15.8 Disposal	45
16 Appendix	46
16.1 Troubleshooting	46
16.2 Accessories	48
16.3 Electromagnetic Compatibility (EMC) Information	49

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

NOTE: The illustrations in this instruction manual show model 11101 CM. However, all descriptions similarly apply to model 11102 CM.

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- ① Connecting cord
- ② Connection for leakage tester or pressure compensation cap (vent port)
- ③ Distal tip
- ④ Shaft
- ⑤ Handle
- ⑥ Deflection lever
- ⑦ Leakage tester 13242 XL
- ⑧ Pressure compensation cap 11025 E
- ⑨ Monitor
- ⑩ Connection 2 for video endoscope
- ⑪ Connection 1 for video endoscope
- ⑫ Mains connection
- ⑬ Reset button
- ⑭ Holder for SD memory card
- ⑮ Status bar
- ⑯ Pilot lamp
- ⑰ On/Off switch
- ⑱ Pilot lamp charging (orange)
- ⑲ Softkeys
- ⑳ Function symbols
- ㉑ USB port
- ㉒ HDMI port
- ㉓ Video recording Start/Stop
- ㉔ White balance
- ㉕ Still image recording
- ㉖ C-HUB®

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Please read these safety instructions carefully before using the device. Before using the instrument on the patient it is imperative that you be acquainted with how the instrument operates and is controlled.

5 Intended use

Video-Rhino-Laryngoscopes are used for endoscopic imaging and image transmission, thus enabling an endoscopic diagnosis of body cavities or hollow organs in ENT medicine.

Video-Rhino-Laryngoscopes are intended for temporary use in both surgical invasive interventions and interventions performed through natural orifices.

5.1 Indications

The use of the Video-Rhino-Laryngoscope is indicated if, in the view of the attending physician, an assessment of the nasal, paranasal sinuses, nasopharynx, laryngeal and esophageal structures is indicated with the flexible endoscope.

5.2 Contraindications

No contraindications relating directly to the product are currently known.

The use of the Video-Rhino-Laryngoscope is contraindicated if, in the opinion of the attending physician, the endoscopic surgical method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

The Video-Rhino-Laryngoscope must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

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6 User qualification

The KARL STORZ CMOS Video-Rhino-Laryngoscope may only be used by persons with an appropriate medical qualification and who are acquainted with the rhino-laryngoscopic technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the CMOS Video-Rhino-Laryngoscope. It is not intended as an introduction to the technique of rhino-laryngoscopy.

7 Compatibility

! **CAUTION!** Non-compatible accessories can lead to malfunctioning or damage to the CMOS video rhino-laryngoscope.

Only use the compatible accessories designated by KARL STORZ.

The CMOS video rhino-laryngoscope is exclusively intended for use with the following KARL STORZ devices.

- **C-MAC® monitor 8403 ZX** (art. no. of the instructions for use: 96076008D)
- **C-HUB® II 20 2903 20** (art. no. of the instructions for use: 96206529D)

The instructions for use of the device used in combination must be followed. The instructions for use can be requested from KARL STORZ or downloaded at www.karlstorz.com.

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8 Warnings and cautions

Please read this manual carefully and follow the instructions exactly. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or safety precautions must be followed to avoid any damage to the product.

NOTE: A note indicates special information about operating the product, or clarifies important information.

8.1 Safety precautions when using the CMOS Video-Rhino-Laryngoscopes

The CMOS Video Rhino-Laryngoscope must be used according to the recognized medical rules and procedures of endoscopy.

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In this case, the reduction of the patient leakage current should be effected using only CF application parts as far as this is possible.

NOTE: The flexible endoscope does not offer any electric insulation against high frequency (HF) currents. Protection against high frequency currents should be ensured through the use of correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.

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Application

WARNING: Every time the CMOS Video-Rhino-Laryngoscope is used, you must first check that it, and any accessories used in combination with it, is in perfect condition. Damaged CMOS Video-Rhino-Laryngoscopes or damaged accessories must not be used.

WARNING: Combinations of medical devices are only then assured to be safe if

- they are identified as such in the respective instruction manual or
- the intended purpose and the interface specifications of the devices used in combination permit this.

WARNING: Danger of explosion if flammable gases are used in the immediate proximity of the CMOS Video-Rhino-Laryngoscope. The CMOS Video-Rhino-Laryngoscope must not be used in the presence of flammable or combustible media.

WARNING: Only the KARL STORZ devices listed in chapter 7 may be used. Accessories from other manufacturers can damage CMOS video rhino laryngoscopes.

WARNING: The CMOS Video-Rhino-Laryngoscope must not be used during the discharge of a cardiac defibrillator. Remove the CMOS Video-Rhino-Laryngoscope from the operation area before discharge takes place.

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WARNING: Danger of glare and burns! Never look directly into the light of the CMOS Video-Rhino-Laryngoscope. The LED light radiated through the CMOS Video-Rhino-Laryngoscope may cause glare or high temperatures at the distal light outlet and therefore lead to burns.

WARNING: If the CMOS Video-Rhino-Laryngoscope should malfunction during application on a patient, stop the application immediately. Place the distal tip of the CMOS Video-Rhino-Laryngoscope in the neutral, nondeflected position and remove

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the CMOS Video-Rhino-Laryngoscope slowly and carefully from the patient.

CAUTION: To avoid damaging the monitor, do not remove the SD card or connecting cord while a video is recording. Only change the SD card when the monitor is switched off.

CAUTION: Do not allow any live components or instruments to come into contact with the CMOS Video-Rhino-Laryngoscope.

WARNING: Damage to the CMOS Video-Rhino-Laryngoscope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.

Reprocessing

WARNING: The CMOS Video-Rhino-Laryngoscope is delivered unsterilized and must be cleaned and disinfected/sterilized before initial and subsequent use.

CAUTION: Do not steam sterilize (autoclave) the CMOS Video-Rhino-Laryngoscopes under any conditions.

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CAUTION: Do not hit the distal tip of the CMOS Video-Rhino-Laryngoscope against hard objects, and never use a sharp object to remove dirt.

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WARNING: The instructions 'Cleaning, Sterilization and Care of KARL STORZ Instruments' must be followed. There the procedures for cleaning, disinfection and sterilization are explained in detail.

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Storage

CAUTION: Do not place the CMOS Video-Rhino-Laryngoscope on devices that become very hot. Do not expose it to direct sunlight or excessive heat, to avoid damage to the device.

CAUTION: Never lay heavy objects on the CMOS Video-Rhino-Laryngoscope.

Maintenance

WARNING: The CMOS Video-Rhino-Laryngoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for CMOS Video-Rhino-Laryngoscopes that are repaired by unauthorized personnel.

WARNING: Danger of electric shock! Refer servicing for this device to authorized personnel only. Any opening of the device by unqualified persons will result in the forfeiting of all warranty claims. Exceptions are those that are described in the instruction manual.

WARNING: In order to prevent infection and protect personnel, it is strictly forbidden to ship contaminated medical devices. To avoid contact and aerogen infections, medical devices must be decontaminated on-site.

WARNING: The **pressure compensation cap** ⑧ must be **removed** from the vent port ② for:

- endoscopic interventions
- immersing in liquid, e.g. disinfectant solution

WARNING: The **pressure compensation cap** ⑧ must be **in place** on the vent port ② for:

- gas and plasma sterilization
- shipping

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9 Initial operation

9.1 Unpacking the device

Check for missing items and evidence of shipping damage.

File any complaints immediately with KARL STORZ or the supplier.

9.2 Basic equipment

CMOS Video-Rhino-Laryngoscope
11101 CM/11102 CM with the following accessories:

- 1x Carrying case 27677 SL
- 1x Instruction manual
- 1x Leakage tester 13242 XL
- 1x Pressure compensation cap 11025 E

Optional accessories:

- Transportation protection sheet 13990 SFN
- Protective cap 8401 YZ
- Adaptor for leakage test in the washer and disinfectant (green) 13242 XA

9.3 Operating conditions

CMOS Video-Rhino-Laryngoscope:

- Ambient temperature: 15 °C – 37 °C
- Humidity: 15 % – 80 %

Monitor 8403 ZX:

- Ambient temperature: ±0 °C – +40 °C
- Humidity: 30 % – 70 %
rel. humidity

9.4 Electromagnetic compatibility

The device 11101 CM has been tested and is found to comply with the EMC limits of EN/IEC 60601-1-2: 2007 CISPR 11 Class B.

The device 11102 CM has been tested and is found to comply with the EMC limits of EN/IEC 60601-1-2: 2014 4th edition CISPR 11 Class B.

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9.5 External system components

Charger

Using the medically approved charger (Art. No. ET27-30-0004515), the Li-Ion battery of the C-MAC® Monitor can be recharged extremely quickly.

Moreover, this charger enables the C-MAC® Monitor to be operated with a discharged Li-Ion battery.

The charger has a green LED display 28, which indicates that the output voltage is correct.

Country specific configuration

The charger features interchangeable primary adapters, which allow the system to be adapted to the country-specific power supply outlets.

- a** Primary adaptor for EUROPE:
ET27-30-0004 516
- b** Primary adaptor for UK:
ET27-30-0004 517
- c** Primary adaptor for US/Japan:
ET27-30-0004 518
- d** Primary adaptor for Australia:
ET27-30-0004 519
- e** Primary adaptor for USA:
ET27-30-0004 520
- f** Adaptor cable-ET27-30-0004 370

Press the button on the power supply unit to change the adaptor (fig. 1). This releases the adaptor which can then be lifted out.

When inserting a new adaptor, this must be placed in the guide (fig. 2) and pushed toward the cable until there is an audible click (fig. 2, see arrow).

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10 Handling

NOTE: Damage to the CMOS Video-Rhino-Laryngoscope resulting from incorrect operation is not covered by the manufacturer's warranty.

WARNING: The CMOS Video-Rhino-Laryngoscope is delivered unsterilized, so it has to be reprocessed prior to initial use and before each subsequent reuse.

CAUTION: Do not store or disinfect the flexible CMOS Video-Rhino-Laryngoscope together with sharp objects in order to avoid causing damage to the outer casing.

NOTE: The C-MAC® monitor has passed a fall and crash test successfully.

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10.1 Operating

The following instructions outline the correct operation of the CMOS Video-Rhino-Laryngoscope and the corresponding accessories. This is not an introduction to the techniques of rhino-/laryngoscopy. Instructions on the endoscopic technique can be found in the relevant medical literature.

The CMOS Video-Rhino-Laryngoscope is intended for visualisation and diagnostic purposes.

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10.2 Holding correctly

Grip the handle of the CMOS Video-Rhino-Laryngoscope in your left or right hand, as shown in the picture. With this kind of control grip, the CMOS Video-Rhino-Laryngoscope can be operated with one hand.

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10.3 Handling the shaft

CAUTION: Do not carry the CMOS Video-Rhino-Laryngoscope by the shaft. Do not pull, clamp or twist the shaft of the CMOS Video-Rhino-Laryngoscope.

CAUTION: The shaft of the Video-Rhino-Laryngoscope must never be kinked or coiled up too tightly. It is best to store the video rhinolaryngoscope in a suspended position or to keep it in a tray specifically designed for sterilizing and storage (39405 AS).

WARNING: Do not apply lubricant to the lens, to avoid clouding the image and thus posing a risk to the patient!

Use a suitable water-soluble lubricant for the distal section of the shaft (For materials recommended and approved insertion by KARL STORZ, please see www.karlstorz.com).

If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the shaft of the CMOS Video-Rhino-Laryngoscope out of the patient's body.

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
10.4 Handling the deflection mechanism

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WARNING: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the CMOS Video-Rhino-Laryngoscope must not be used.

WARNING: The deflection lever must never be operated with force or abruptly to avoid perforation and injury to the patient or damaging the device irreparably.

WARNING: The distal tip of the CMOS Video-Rhino-Laryngoscope must be straight, not angled, when the videoscope is inserted into the patient

WARNING: Before removing the CMOS Video-Rhino-Laryngoscope from the patient, the distal tip must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.

WARNING: If a functional defect occurs with the CMOS Video-Rhino-Laryngoscope during use on the patient, discontinue use immediately. Put the distal tip of the CMOS Video-Rhino-Laryngoscope into the neutral, non-deflected position and remove the CMOS Video-Rhino-Laryngoscope slowly and carefully from the patient.

CAUTION: The resistance that can be felt when deflecting the tip using the deflection lever increases as the angle of deflection increases. Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.

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WARNING: The lens at the distal end of the CMOS Video-Rhino-Laryngoscope may become detached, e.g. due to mechanical forces. For this reason, it is essential to check that the lens is undamaged and present before and after decontamination.

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10.5 Connection to the monitor 8403 ZX

WARNING: It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces, or voltage-carrying elements of other units, at any time.

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Insert the connecting cord ① of the flexible CMOS Video-Rhino-Laryngoscope into either connection socket ⑩ or ⑪ on the monitor.
The connecting cord ① can be connected or disconnected during operation.

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Insert the cable from the power supply unit into the connection socket ⑫ of the C-MAC® Monitor and connect the power supply unit with the mains socket. The power cord can be fixed to the connection socket using the union nut.
The battery is charged automatically as soon as the C-MAC® Monitor is connected to the power supply. This will be indicated by the lighting up of the orange colored LED ⑬.

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10.6 Use of the HDMI output

The device features a port for an external monitor. An HDMI cable can be used to send the video signal to an external monitor, if necessary.

CAUTION: The C-MAC® Monitor must be switched off when an external monitor is connected.

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- NOTE:** The connected monitor must support a resolution of 1024 × 768 or 1280 × 720.
- NOTE:** To connect the C-MAC® Monitor to a KARL STORZ monitor, the HDMI-DVI adaptor with article number **209190 10** is required.

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10.7 Functions of the C-MAC® Monitor 8403 ZX

A more detailed description of the operation and functions of the monitor 8403 ZX can be found in the corresponding instruction manual supplied with every monitor.

In case of loss of the instruction manual, you can also download it from our homepage at www.karlstorz.com.

! WARNING: The KARL STORZ flexible CMOS Video-Rhino-Laryngoscope can not be used in combination with camera control units and light sources. The use of non-compatible units can lead to patient injury or unit damage.

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10.7.1 Focus

The image focus is adjusted automatically. Manual adjustment is not necessary.

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10.8 Function buttons of the endoscope

The endoscope is equipped with 3 keys, which are used to activate the functions displayed as symbols on the endoscope.

The following functions are available:

- Video recording 23: Starts or stops the recording of a video stream on the SD memory card inserted in the receptacle.
- 📷 White balance 24: Enables the user to carry out a manual white balance. Point the tip of the videoscope at a white surface and press the button for 1 sec.
- 📷 Single image capture 25: Saves the live image on the screen on the SD memory card.

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10.9 Monitor function keys

The monitor is equipped with 4 softkeys. These can be used to activate the functions displayed as symbols on the screen.

- ❗ **NOTE:** The **video recording** and **image storage** symbols are only visible if an SD memory card is inserted. The **change the camera unit** symbol is only visible if two video units are connected.

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10.10 Preparation prior to use

NOTE: Before the device is used for the first time, it must be subjected to a complete reprocessing cycle, excluding preliminary cleaning.

WARNING: Every time the CMOS Video-Rhino-Laryngoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged CMOS Video-Rhino-Laryngoscopes or damaged accessories must not be used.

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10.10.1 Inspecting the surfaces

WARNING: The surface of the flexible shaft must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

CAUTION: Do not squeeze the distal tip of the CMOS Video-Rhino-Laryngoscope too tightly.

Every time the CMOS Video-Rhino-Laryngoscope is used, you must first check that its surface, and above all the surface of the shaft and distal tip, are in perfect condition and that there are no sharp edges. CMOS Video-Rhino-Laryngoscopes with sharp edges must not be used. Otherwise serious injuries may be caused to the patient or the user.

Inspect the connecting cord for damage (such as breaks, cracks, twisted or squashed sections...) and check the supply connector for loose connections.

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WARNING: Before use, always check that the deflection lever and the distal tip connected to it function correctly. Defective CMOS Video-Rhino-Laryngoscopes must not be used.

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The distal tip is deflected by actuating the deflection lever slowly.

The distal end will deflect in accordance with the direction of actuation (forward or backward).

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10.11 Switching the CMOS Video-Rhino-Laryngoscope on/off

Switching on/off: The CMOS Video-Rhino-Laryngoscope with the integrated LED light source is switched on by pressing the power button ① on the monitor. Pressing the key a second time switches the unit off.

WARNING: Danger of glare: Do not look directly into the light beam.

WARNING: When in operation, the integrated LED light source heats up the CMOS Video-Rhino-Laryngoscope. Switch off the Video-Rhino-Laryngoscope before putting it down and do not place it in the proximity of the patient (as per IEC 60601-1-1).

NOTE: The function keys ① on the monitor are used for switching the unit on/off. The CMOS Video-Rhino-Laryngoscope itself does not have any function keys.

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Video recording: When the control key 'Video recording' on the C-MAC® Monitor or the flexible CMOS Video-Rhino-Laryngoscope is pressed, the flexible CMOS Video-Rhino-Laryngoscope records a video stream on the SD memory card. Pressing the key a second time stops the recording.

NOTE: The video stream can only be viewed on a PC if an MPEG 4-codec is installed on the PC.

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Single image capture: When the control key 'Single image capture' on the C-MAC® Monitor or the flexible CMOS Video-Rhino-Laryngoscope is pressed, the current image on the monitor is saved on the SD memory card.

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10.12 Inspecting the image quality


WARNING: If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the CMOS Video-Rhino-Laryngoscope from the patient's body gently and carefully. If the CMOS Video-Rhino-Laryngoscope is used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.

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10.13 Switching off the C-MAC® Monitor

Standby: Briefly pressing the power button 17 will set the C-MAC® Monitor to standby operation mode. In the case of brief interruptions, this enables the C-MAC® Monitor to be ready for operation again quickly. The standby mode is indicated by a flashing green pilot lamp 16 above the power button.

After 5 minutes, the system switches itself from standby mode to the off mode.

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Power off: Pressing the power button 17 for more than approx. 3 seconds will switch off the C-MAC® Monitor. The green pilot lamp 16 will go out.

NOTE: When the C-MAC® Monitor is switched off, the integrated light source in the flexible CMOS Video-Rhino-Laryngoscope also switches off.

NOTE: Information on the operation and functions of the C-MAC® Monitor 8403 ZX can be found in the corresponding instruction manual supplied with every C-MAC® Monitor (Art. No. 96076008DER). In case of loss of the instruction manual, you can also download it from our homepage at www.karlstorz.com.

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10.14 System reset

If a system error occurs (unit no longer reacts, although connections and power supply are in order), a system reset has to be implemented. This is the equivalent of rebooting a computer. For this, a thin screwdriver or similar is inserted into the round opening next to the SD card holder (14) and carefully pushed downward in order to activate the reset button (13).

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10.15 Monitor bag 8403 YD

The optionally available protective bag 8403 YD facilitates the mobile use of the monitor.

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- This reprocessing step is part of the validated procedure
- This reprocessing step may be carried out additionally in order to aid the cleaning. The sterilization procedure is approved in respect of material compatibility; however, it must be validated by the operator on site in respect of effectiveness.
- X On the grounds of potential damage, this reprocessing step must NEVER be performed

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Reprocessing

 **WARNING:** Risk of infection: Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.

 **WARNING:** Risk of infection from non-sterile delivery. Reprocess the device before each use. Check the medical device for visible soiling before use.

 **WARNING:** Before reusing the device, replace any defective parts or send the flexible endoscope for repair.

 **ATTENTION:** The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

 **CAUTION:** Never clean flexible endoscopes with metal brushes.

 **CAUTION:** Constant changing between various reprocessing methods is hard on the materials and should be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, is possible.

 **CAUTION:** Flexible endoscopes must never be steam sterilized (autoclaved).

 **CAUTION:** To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.

 **NOTE:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

 **NOTE:** Observe the chemical manufacturer's specifications concerning concentration and exposure times. Bear in mind the microbiological range of action of the chemicals.

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- i NOTE:** Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be found at www.karlstorz.com.
- i NOTE:** National laws and regulations must be observed.
- i NOTE:** We recommend using a plastic basin for disinfection to avoid damaging the endoscope. Furthermore, use of a plastic basin prevents electrolytic corrosion which can occur if there are different metals in a solution.
- i NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

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Accessories: pressure compensation cap, manual leakage tester

Place of application	Art. number brush/ adaptor
Clean surface	27652
Ventilation connector	13242XA

11.3 Preliminary cleaning directly after use

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Wipe-clean the exterior surfaces of the endoscope with a cloth moistened with disinfectant.

- i NOTE:** To avoid damaging the shaft, do not use pushing and pulling movements.
- i NOTE:** Do not use disinfectants for preliminary cleaning, these cause fixation of proteins and organic residues

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11.4 Leakage test

The endoscope should be tested to ensure it is leakproof after use and before reprocessing.

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Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal end may expand slightly. Using the deflection mechanism on the control body, move the deflection part in all possible directions as far as possible and observe the pointer manometer. The test pressure can drop by up to 10 mmHg until the endoscope is completely filled with air.

After stabilizing the pressure on manometer wait at least 2 minutes to see whether the pressure remains constant or drops.

11.4.2 Performing the leakage test – in liquid

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Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in 'Leakage test – dry' and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer

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11.6.2 Manual disinfection

The endoscope must be completely immersed in a cleaning solution. At the end of the exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

11.6.3 Manual High Level Disinfection (HLD)

The endoscope must be completely immersed in a high level disinfectant solution. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

The endoscope can be disinfected with the following high level disinfection solutions which have been tested by way of example:

- 2.4% glutaraldehyde solution (e.g., CIDEX® Activated Dialdehyde Solution)
- 0.55% ortho-phthalaldehyde solution (e.g., CIDEX® OPA Solution)
- 2.0% accelerated hydrogen peroxide solution (e.g., Revital-Ox® Resert™ or Revital-Ox Resert XL HLD)

11.6.4 Drying

Finally, the surfaces are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

CAUTION: The pressure must not exceed 100 kPa (1 bar / 14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

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11.7 Automated reprocessing

The following procedure for machine decontamination has been validated and approved subject to compliance with the process parameters described in the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments":

11.7.1 Manual preparation

The endoscope must be completely immersed in a cleaning solution. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush.

Finally, it must be rinsed with cold water to ensure neutralization.

CAUTION: The pressure compensation cap must not be in place on the vent port for cleaning and disinfection.

11.7.2 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C and a maximum pressure of 100 kPa (1 bar / 14.5 psi) must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

The selection of a suitable endoscope holder, which should ensure that the medical device is thoroughly rinsed out, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off manually afterwards.

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- STERRAD® 100S Short Cycle with Boosters
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

11.10.2 Hydrogen peroxide (H₂O₂) sterilization – STERIS® AMSCO® V-PRO™

NOTE: The “STERIS Low Temperature Reprocessing Device Compatibility Matrix” can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.

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- V-PRO® maX Flexible Cycle
- V-PRO® 60 Flexible Cycle

11.10.3 Ethylene oxide sterilization (EO)

100% ethylene oxide at 55°C and a hold time of 62 minutes acc. to ANSI/AAMI/ISO Norm 10993-7.

The EO sterilization is validated with the following parameters:

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Ventilation:

- 12 hours
- 55 +/- 4 °C

11.10.4 Low-temperature steam and formaldehyde process (LTSF)

KARL STORZ had only the specified sterilization method validated with regard to a reduction in microorganisms on the order of 10⁻⁶ SAL. However, due to the large number of processes used in the global market, validation in respect of efficacy was not carried out. Therefore, the user must validate his specific processes individually in collaboration with the device manufacturer.

WARNING: Risk of infection: If FO sterilization is performed, the method must be validated on site.

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Horizontal storage is possible for sterilized endoscopes. KARL STORZ recommends the use of a specially designed tray for sterilization and storage.

CAUTION: Do not expose flexible endoscopes to direct sunlight or excessive heat.

CAUTION: The distal tip of the flexible endoscope must be in the neutral position during storage.

CAUTION: Do not store the flexible endoscope in the case.

CAUTION: The shaft of the flexible endoscope must never be kinked or coiled up too tightly.

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13 System description

The CMOS Video-Rhino-Laryngoscope from KARL STORZ has a 300 mm pellethane working shaft, 3.70 mm (11101 CM) or 2.9 mm (11102 CM) in diameter. The PEEK (polyetheretherketone) distal tip is 3.70 mm or 2.9 mm in diameter and houses the CMOS video module. The distal tip of the flexible CMOS Video-Rhino-Laryngoscope can be deflected 140° up or 140° down. The handle houses the deflection lever, the function keys and the LED light source.

The CMOS video rhino laryngoscope from KARL STORZ is connected to the recommended monitor 8403 ZX using the integrated 2 m supply cable. Images and videos can be saved using a SD memory card and the control keys on the endoscope or on the monitor, which are assigned, for example, the functions of video recording and image capture.

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C-MAC® Monitor 8403 ZX

Length: 230 mm
Width: 155 mm
Depth: 54 mm
Weight: 1250 g
Drip-water protection: IP54

CAUTION: Drip-water protection can only be ensured if the connection sockets are sealed using the designated rubber seals.

Classification: IEC Class II
Input voltage: 12 V
Current consumption: 1.5 A max.
Battery: Lithium-Ion, 3.6 V 7800 mAh
Resolution: 1280 x 800
Screen width: 150 mm
Screen height: 94 mm
Memory card: SD card
Video out signal: The connected monitor must support the 720p mode.

Operating/storage conditions

Operating conditions: 0 °C to +40 °C
30 % to 70 % rel. humidity
Storage conditions: -10 °C to +60 °C
5 % to 95 % rel. humidity
Atmospheric pressure: 700 to 1080 hPa

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SD card

Throughput: 10 MB/s
Compatibility: Compatible with the Sandisk Ultra II HD card and SDHC cards class 10

Color settings

Basic Setting

Brightness 50 %
Contrast 35 %
Saturation 31 %
Hue 50 %

Endoscope 11101 CM/11102 CM

Direction of view: 0°
Angle of view: 85° (11101 CM)
100° (11102 CM)
Deflection: up 140°, down 140°
Working length: 30 cm
Outside diameter: 3,7 mm (11101 CM)
2,9 mm (11102 CM)
Illumination: integrated LED light source
Electrical power consumption: 1W
Sensor: CMOS technology
Weight 245 g (11101 CM)
325 g (11102 CM)
Drip water protection: IPX8
Classification: IIa
Applied part: Type BF

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15.3 Limitation of liability

KARL STORZ SE & Co. KG shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the conditions that all assembly, system expansions, readjustments, modifications, or repairs to same have been performed by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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15.5 Standard compliance

According to IEC 60601-1, 60601-2-18, UL 60601-1, CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electric shock:
 - In combination with the monitor 8403 ZX: Internally powered during battery operation, Protection class II during line operation
 - In combination with the C-HUB®: Protection class II
- Degree of protection against electric shock: Applied part type BF 
- Type of protection against moisture: Protection against continuous immersion: IPX 8

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15.6 Directive compliance

According to Medical Device Directive (MDD):

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8403 ZX: This medical product belongs to Class I

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15.8 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its service life, dispose of the device as electronic waste. Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information about your local collection point. Within the scope of application of this directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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Stand, height 120 cm, star-foot roller stand with antistatic castors, crossbar 25 cm x diameter 25 mm for positioning the monitor, including storage tray for laryngoscopes, measurements (w x d x h): 30 x 20 x 10 cm

Crossbar, for stands 8401 YA, 50 cm x diameter 25 mm, for positioning the monitor 8403 ZX

Crossbar, for stands 8401 YA, 70 cm x diameter 25 mm, for positioning the monitor 8401 ZX

Multifunctional clamp with VESA 75 standard holder. Suitable for securing rods with a max. diameter of 24 mm and rails with a max. height of 25 mm

Multifunctional cap to protect contacts during reprocessing. The cap is reusable.

Holder for flexible endoscopes, for mounting to standard tubes, incl. installation accessories

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Electromagnetic Compatibility user information

The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: The use of an accessory or cable with the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX other than those specified in this manual may result in increased emissions or decreased immunity of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX. When using an accessory or cable with the CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX other than those specified in this manual, it becomes the responsibility of the user of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX to determine compliance with IEC 60601-1-2 when using this item.

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WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Accessories and cables which have been shown to comply with EN/IEC 60601-1-2:

Cable type	Shielded	Length [m]	Ferrite	Used for
Charger ET27-30-0004515	No	2.0	No	Charging the internal battery

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Table 1 Guidelines and manufacturer's declaration – electromagnetic immunity			
The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is intended for use in an electromagnetic environment specified below. The user of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	If image distortion occurs, it may be necessary to position the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	

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Table 3 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.			
Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
			Portable and mobile RF communications equipment should be used no closer to any part of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	$d = 1.2 \sqrt{P}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the $d = 1.2 \sqrt{P}$ 80 MHz to < 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol:
			
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects and persons.			
^{a)} Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is used exceeds the above compliance levels, the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.			
^{b)} Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Table 4 Guidelines and manufacturer's declaration – electromagnetic emissions		
The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is intended for use in the electromagnetic environment specified below. The customer or user of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX uses RF energy only for its internal function. The customer or user of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class B	The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is suitable for use in all establishments including domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class B	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5 Recommended separation distances between portable and mobile RF telecommunication devices and the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX			
<p>The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.</p> <p>NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.</p>			

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BC „Arman“, office 314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Jinliang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: , Fax: +63 2 31745-00
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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Important note for users of KARL STORZ devices and instruments

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It is recommended to check the suitability of the products for the planned intervention prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper handling, cleaning, disinfection and sterilization of the Video-Uretero-Renoscope FLEX-X^c. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for ready reference.

WARNING: KARL STORZ instruments are delivered unsterilized, and must be cleaned, disinfected and/or sterilized prior to the initial use and before each subsequent use.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

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
Optionales Zubehör:

- 11275 FE
- 11275 ZE
- 110931-50

- ① Control buttons
- ② Supply connector
- ③ Connection for leakage tester or pressure compensation cap
- ④ Distal tip
- ⑤ Sheath
- ⑥ Handpiece
- ⑦ Working channel with LUER-Lock connectors
- ⑧ Deflection lever
- ⑨ Biopsy forceps 27023 ZE
- ⑩ Cleaning brush 27651 AL
- ⑪ Leakage tester 13242 XL
- ⑫ Pressure compensation cap 11025 E
- ⑬ LUER adaptor
- ⑭ Single-use cleaning brush 110931-50

11278 VSK/VSUK comprises:

- 11025 E Pressure compensation cap
- 13242 XL Leakage tester
- 27651 AL Cleaning brush
- 27014 Y LUER adaptor, with seal

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Important note for users of KARL STORZ devices and instruments	I
Pictures of the instruments.....	II
Controls, displays, connectors and their uses.....	III
Safety information	3
Intended Use	3
Indications	4
Contraindications.....	5
User qualification	5
Safety precautions when using the Videoendoscopes	6
Warnings and cautions.....	6
Handling	10
Unpacking the device.....	10
Basic equipment	10
Operating conditions	10
Electromagnetic Compatibility.....	10
Handling	11
Operating	11
Holding the Videoendoscope correctly	11
Handling the sheath	11
Handling the deflection mechanism.....	12
Handling the accessories (forceps, cleaning brushes, HF electrodes etc.) ..	14
Compatible camera control units	15
Connection to the camera control unit	16
Control buttons	17
Preparation prior to use.....	18
Inspecting the surfaces.....	18
Checking the deflection	19
Inspecting the LED function.....	20
Inspecting the camera control unit.....	20
Inspecting the instruments and the instrument channel	21

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Cleaning, disinfection, care and sterilization.....	22
Explanation of Warnings and Cautions.....	22
General warnings	23
Accessories.....	25
Disassembly.....	25
Precleaning directly after use.....	25
Leakage test.....	26
Performing the leakage test – dry	26
Performing the leakage test – in liquid	26
Manual cleaning	27
Manual disinfection.....	28
Machine cleaning and disinfection	28
Manual preparation.....	29
Machine cleaning/chemical-thermal disinfection	29
Assembly, inspection and care	30
Sterilization	31
Maintenance	36
Servicing and repair.....	36
Maintenance.....	36
Limitation of liability.....	36
Warranty.....	37
Directive compliance	37
Important information	38
Disposal	38
Appendix	39
Troubleshooting	39
Subsidiaries.....	42

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Please read these safety instructions carefully before using the instrument. Before using the instrument on the patient it is imperative that you be acquainted with how the instrument operates and is controlled.

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Intended Use **11278 V/VS/VU/VSU**

Flexible uretero-renoscopes (fiber endoscopes and videoscopes) are used for endoscopic imaging of the renal collecting system for the purposes of diagnosis and treatment. Flexible uretero-renoscopes with an integrated working channel allow for the use of flexible auxiliary instruments such as forceps, stone baskets etc., coagulation electrodes and laser fibers for lithotripsy or for the ablation of benign/malignant tissue. The flexible uretero-reno-scope can be used for the imaging/treatment of all renal calices, as well as the pyeloureteral junction or the ureter. The special feature of the flexible uretero-reno-scope lies in the fact that it can also be used percutaneously via an operating sheath, and can be used in antegrade mode to treat all renal calices as well as the pyeloureteral junction or the ureter.

Flexible uretero-renoscopes are designed for short-term use in endoscopic interventions performed through natural orifices or surgical invasive interventions.

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Indications

11278 V/VS/VU/VSU:

The use of ureteroscopes/uretero-renoscopes is indicated if uretero(reno)scopy is indicated in the opinion of the attending physician.

This includes:

- Further diagnostic measures in the event of conspicuous imaging/filling defects in the upper urinary tract shown in the imaging
- Pathological bladder-wash cytology without reference to bladder tumors or from the upper urinary tract
- Unclear haematuria (macro/microhaematuria)
- Diagnosis and treatment of urinary stones in the ureter or renal collecting system.
- Diagnosis and treatment of ureteral strictures or tumors
- Iatrogenically introduced foreign material
- Narrow pyeloureteral junction
- Possibly the percutaneous antegrade treatment of kidney stones, tumors and ureter outflow constrictions.
- UUT (Upper Urinary Tract) tumors or UTUC (Upper Tract Urothelial Cancer)
- Nephrolithiasis of lower calyceal group
- Positive bladder-wash cytology from the upper urinary tract
- Tumors in the upper urinary tract
- Indication retrograde treatment of a transitional cell carcinoma
 - impaired renal function/single kidney
 - bilateral tumor development
 - multiple concomitant diseases/limited life expectancy
 - Tumor stages in which only palliation is possible

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Contraindications**11278 V/VS/VU/VSU:**

The use of uretero(reno)scopes is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Uretero(reno)scopes must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

The following contraindications also apply:

- Urinary tract infection
- Urosepsis
- Untreated/unresolved coagulation disorder
- Contraindication for lithotomy position
- Pregnancy

User qualification

The KARL STORZ Videoendoscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the technique of uretero-renaloscopy. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the Videoendoscopes.

It is not intended as an introduction to the technique of renaloscopy.

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WARNING: Read the following safety instructions carefully before using the Videoendoscope to avoid putting your patients, personnel or yourself at risk.

WARNING: Inspect the videoendoscope prior to each surgical procedure to ensure that it functions correctly and that it has been correctly cleaned, disinfected and/or sterilized.

WARNING: Every time the Videoendoscope is used, you must first check that it, and any accessories used in combination with it, is in perfect condition. Damaged Videoendoscopes or damaged accessories must not be used.

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Pay careful attention to the instructions and interface specifications of medical devices used in combination.

WARNING: The Videoendoscope is delivered unsterilized and must be cleaned and disinfected/sterilized before the initial use as well as before any further application.

WARNING: The **pressure compensation cap** must be **removed** from the vent port for:

- endoscopic interventions
- manual/machine reprocessing

WARNING: Only use accessories (guide probes, cleaning brushes etc.) that are specifically intended for use with Videoendoscopes and can be operated in line with generally accepted technical standards and occupational safety guidelines.

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WARNING: Do not allow any live components or instruments to come into contact with the Videoendoscope.

WARNING: The Videoendoscope must not be used during the discharge of a cardiac defibrillator. Remove the Videoendoscope from the operation area before discharge takes place.

WARNING: Danger of explosion if flammable gases are used in the immediate proximity of the Videoendoscope. The Videoendoscope must not be used in the presence of flammable or combustible media.

WARNING: As a result of the light emitted by the endoscope, the tip of the latter may heat up.

Prolonged direct contact of light-emitting surfaces on tissue can lead to burns, and must therefore be avoided at all costs.

WARNING: The **pressure compensation cap** must be **in place** on the vent port for:

- hydrogen peroxide gas sterilization
- venting
- shipping

WARNING: The instructions 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' must be followed. There the procedures for cleaning, disinfection and sterilization are explained in detail.

WARNING: The Videoendoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for Videoendoscopes that are repaired by unauthorized personnel.

WARNING: Damage to the Videoendoscope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.

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WARNING: Always unplug the Videoendoscope before carrying out any maintenance and cleaning work on it.

WARNING: Risk of infection during shipment! Contaminated Videoendoscopes may cause serious or life-threatening illnesses during transportation or on delivery to the recipient. Defective Videoendoscopes with attached pressure compensation cap must be wiped down with disinfectant and, if necessary, sterilized before shipment. Complete the repairs form.

WARNING: If a malfunction should occur during application on a patient, stop the application immediately. Place the distal tip of the Videoendoscope in the neutral, nondeflected position and remove the video endoscope slowly and carefully from the patient.

CAUTION: Do not expose the Videoendoscopes to direct sunlight or excessive heat, to avoid damage to the device.

CAUTION: Do not steam sterilize (autoclave) the Videoendoscopes.

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CAUTION: Do not hit the distal tip of the Videoendoscope against hard objects, and never use a sharp object to remove dirt.

CAUTION: Never lay heavy objects on the Videoendoscope.

NOTE: Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided. However, if the reprocessing method is changed just once this can be regarded as harmless.

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Electromagnetic Compatibility

This instrument has been tested and complies with the EMC limits as per EN/IEC 60601-1-2: 2001 CISPR 11 Class B.

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Handling

! **NOTE:** Damage to the Videoendoscope resulting from incorrect operation is not covered by the manufacturer's warranty.

! **WARNING:** The Videoendoscope is delivered unsterilized, so it has to be cleaned, disinfected/sterilized prior to initial use and before each subsequent reuse.

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Operating

The following instructions outline the correct operation of the Videoendoscope and the corresponding accessories. This is not an introduction to the techniques of renoscopy. Instructions on endoscopic technique can be found in the relevant medical literature.

Suitable accessories for your Video-uretero-Renoscope FLEX-X^c, such as biopsy forceps, can be found in the UROLOGY catalog or extracts from this catalog.

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Holding the Videoendoscope correctly

Grip the handle of the Videoendoscope in your left or right hand, as shown in the picture. With this kind of control grip, the Videoendoscope can be operated with one hand.

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Handling the sheath

! **CAUTION:** Do not carry the Videoendoscope by the sheath. Do not pull, clamp or twist the sheath.

! **CAUTION:** The sheaths of the Videoendoscope must never be kinked or coiled up too tightly. It is best to suspend the Videoendoscopes for storage.

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WARNING: Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!

Use a suitable water-soluble lubricant for the distal section of the insertion sheath. Observe the lumen the whole time during use. If you lose your orientation, pull the sheath back slightly in the neutral position until vision is restored.

If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the sheath of the Videoendoscope out of the patient's body.

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Handling the deflection mechanism

NOTE: The neutral position of the distal tip is straight, without any deflection whatsoever.

WARNING: Suitable accessories, such as laser fibers and biopsy forceps, must only be introduced into the instrument channel of the Videoendoscope with the distal tip in the neutral, non-deflected position. Once the instrument has been introduced, actuate the deflection lever gently and carefully as far as the limit stop. The deflection lever must never be operated with force or abruptly to avoid injuring the patient or damaging the device!

WARNING: Before removing the Videoendoskop from the patient, the distal tip of the Videoendoscope must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.

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Handling the accessories (forceps, cleaning brushes, HF electrodes etc.)

CAUTION: Only use flexible accessories (forceps, cleaning brushes etc.) from KARL STORZ. Accessories from other manufacturers may damage the Videoendoscope.

CAUTION: Only introduce the instruments from the proximal port of the Videoendoscope.

WARNING: When HF electrodes or laser probes are used, it must be ensured that power is not applied until the active parts have emerged from the instrument channel and are in the user's field of vision, and the site of the operation is clearly recognizable. Activation in the instrument channel may lead to serious injury to the patient and damage to the Videoendoscope.

WARNING: When HF electrodes or laser probes are used, it must be ensured that power is only applied when there are no flammable gases in the operating field or in the proximity of the Videoendoscope.

CAUTION: When laser probes are used, it must be ensured that the laser beam cannot hit the distal tip of the Videoendoscope.

CAUTION: Do not force accessories, such as biopsy forceps, through the working channel, to avoid damaging the device. It should be possible to move the accessory easily, without any hindrance and with closed jaws through the entire channel.

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CAUTION: When the distal tip is deflected, do not introduce any instrument into the Videoendoscope. Return the distal tip to the neutral position before introducing the instrument and then move the deflection lever gently and carefully to the limit stop.

CAUTION: Inserting sharp-tipped accessories into the instrument channel may damage the Videoendoscope.

CAUTION: Foreign bodies or parts to be extracted must not be removed through the instrument channel using instruments. Instead, the entire Videoendoscope must be removed from the patient and the foreign body extracted.

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WARNING: The KARL STORZ Videoendoscope should not be used in combination with camera control units from other manufacturers. The use of non-compatible units may cause injury to the patient or damage to the instrument. It must be ensured that parts of the endoscope that are not electrically insulated do not come into contact with conductive surfaces or live parts of other units, at any time.

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Connection to the camera control unit

CAUTION: The Videoendoscope can be used with the above mentioned camera control units. The instruction manual of the camera control unit must be followed. The complete operating instructions for the camera control unit can be found in the instruction manual.

WARNING: Risk of burns and glare! Never look directly into the light emitted by the Videoendoscope. The high-intensity light emitted by the Videoendoscope at the distal light outlet can cause glare or high temperatures and thereby result in burns.

WARNING: Risk of burns and fire hazard! Never place the distal tip of the Videoendoscope on or under a sterile drape. The intensity of the light can cause burns to the patient and/or the sterile drape.

WARNING: Reduce the light intensity to the extent that optimal illumination of the operation field is still ensured. The higher the intensity of the light, the greater the heat emitted at the tip of the Videoendoscope.

WARNING: The keyboard of the IMAGE1 S must not come into contact with patients or be placed in the vicinity of patients.

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CAUTION: Take care with the Videoendoscope connections! The Videoendoscope connections can be damaged if connected incorrectly. Make sure that the connector is correctly aligned.

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
1. Install the camera control unit on a firm, level surface.
2. Plug the power cord into the specified power socket.
3. Plug the Videoendoscope supply cable into the HD socket of the camera control unit.
4. Activate the image and illumination function (integrated LED) by switching on the camera control unit.

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Control buttons

The Videoendoscope has three control buttons which can be assigned specific functions in the camera control unit's device menu.

- a Scroll up: Scrolls up the options in the menu.
- b Menu/Select: Displays the live menu. As soon as this menu is displayed, this button activates the highlighted option (identical to ENTER on a computer keyboard).
- c Scroll down: Scrolls down the options in the menu.

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The assignment is retained after switching off and on again.
Instructions on programming the control keys can be found in the instruction manual of the camera control unit.

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Preparation prior to use


i NOTE: Before the device is used for the first time, it must be subjected to a complete reprocessing cycle (excluding preliminary cleaning) in accordance with the risk classification.

! WARNING: Every time the Videoendoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged Videoendoscopes or damaged accessories must not be used.

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Lightly grease the running surfaces of the cock plug (special grease for stopcocks; Art. no. 27657).

i NOTE: This thin layer of grease not only serves to lubricate the metal surfaces which are in contact with one another, but also makes the stopcocks leakproof. Oil is not suitable as a sealant.

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Inspecting the surfaces


! WARNING: The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

! CAUTION: Do not squeeze the distal tip of the Videoendoscope too tightly.

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Every time the Videoendoscope is used, you must first check that its surface, and above all the surface of the sheath and distal tip, is in perfect condition and that there are no sharp edges. Videoendoscopes with sharp edges must not be used. Otherwise serious injuries may be caused to the patient or the user.

Inspect the supply cable for damage (such as breaks, cracks, twisted or squashed sections...) and the supply connector for loose connections.

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Checking the deflection

WARNING: Actuate the deflection lever gently and carefully as far as the limit stop. The deflection lever must never be operated with force or abruptly, to avoid injuring the patient or damaging the device! Move the deflection lever smoothly and slowly!

WARNING: Before use, always check that the deflection lever and the distal tip connected to it function correctly. Defective Videoendoscopes must not be used.

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The distal tip is deflected by means of the deflection lever. Turn the deflection lever slowly; the end piece should move downward. Turn the deflection lever in the opposite direction to produce an upward movement.

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Inspecting the LED function

Prepare the IMAGE1 S as described above and connect to the Videoendoscope. Emission of light at the distal end should be recognizable. Light intensity is automatically regulated by the camera control unit.

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Inspecting the camera control unit

WARNING: If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the Videoendoscope from the patient's body gently and carefully. If the Videoendoscope is used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.

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Prepare the camera control unit as described above and connect to the Videoendoscope.

CAUTION: The camera control unit must be switched off before connecting the Videoendoscope!

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The control buttons can be individually programmed.

Programming of the control buttons is described in the instructions for the camera control unit.

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Inspecting the instruments and the instrument channel

CAUTION: It should be possible to move the forceps with closed jaws easily and without any hindrance through the entire channel.

CAUTION: Accessories, such as biopsy forceps, which are bent or twisted may no longer function properly, are difficult to pass through the channel, and may cause serious damage to the channel. They must therefore be replaced!

Check the instrument channel by introducing the appropriate instrument (max. diameter of instruments: 3 Fr.).

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**Explanation of Warnings and
Cautions**

WARNING: A Warning indicates a potential hazard. Failure to observe a Warning may result in injury to the patient, user or third parties.

CAUTION: A Caution indicates that particular procedures or precautions must be followed to avoid possible damage to the medical devices.

NOTE: A Note contains special information regarding reprocessing or important explanations.

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General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: Before reusing the device, replace any defective parts or send the flexible endoscope for repair.

CAUTION: When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

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- !** **CAUTION:** Danger of damage to medical devices: The use of chemicals that have not been approved by KARL STORZ or chemical manufacturers may cause damage to the medical devices. Only use chemicals approved for the medical device for reprocessing. A helpful guide is available on the Internet at www.karlstorz.com.
- !** **CAUTION:** National laws and regulations must be observed.
- !** **CAUTION:** The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.
- !** **CAUTION:** Never clean flexible endoscopes with metal brushes.
- !** **CAUTION:** Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, is possible.
- !** **CAUTION:** Flexible endoscopes must never be steam sterilized (autoclaved).
- !** **CAUTION:** To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.
- i** **NOTE:** The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.
- i** **NOTE:** The valves must be removed in their individual parts prior to cleaning and disinfection.

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If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer does not show a fall in pressure, the device is leak-proof.

Remove the endoscope from the liquid and press the vent button on the manometer. The test pressure must drop to 0 before the leakage tester can be removed.

! **NOTE:** Never release the pressure in the liquid.

Cleaning can then begin.

If a pressure loss occurs, this means that the endoscope leaks and may no longer be used, immersed in liquid or reprocessed by machine.

! **CAUTION:** In this case, fill out the accompanying repairs form, mark the endoscope, which has been cleaned by wiping to remove organic soiling/residue, as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the transport protection sheet 13990 SFN.

Manual cleaning

The endoscope must be completely immersed in a cleaning solution. The working channel must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The channels are brushed until the brush head emerges and is visibly clean. Before retracting the brush head that emerges at the distal end, clean it. Never brush back and forth within the channel. Finally, it must be rinsed with cold water to ensure neutralization.

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CAUTION: Never remove dirt with a needle or other sharp objects, since these can cause severe damage.

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Manual disinfection

The endoscope must be completely immersed in a cleaning solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. At the end of the exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

Drying

Finally, the surfaces and channels are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

CAUTION: The pressure must not exceed 100 kPa (1 bar / 14,5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual, 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003).

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Manual preparation

The endoscope must be completely immersed in a cleaning solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The channels are brushed until the brush head emerges and is visibly clean. Before retracting the brush head that emerges at the distal end, clean it. Never brush back and forth within the channel.

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Connecting

In order to ensure effective machine cleaning and disinfection, the endoscope's working and suction channel must be connected to the washer and disinfectant using the recommended adaptors.

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Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

The selection of a suitable endoscope holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

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- If there is still residual moisture on or in the endoscope, it must be dried again manually. The channels must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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NOTE: Check for a clear view and steerability.

NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

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Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

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The following sterilization procedures have been validated and approved by KARL STORZ for this medical device:

**Hydrogen peroxide (H₂O₂)
sterilization - ASP STERRAD®**

 **WARNING:** Medical devices must be dismantled for sterilization.

 **WARNING:** Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

 **WARNING:** Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

 **CAUTION:** For detailed information, please consult the user handbook of the respective device.

 **CAUTION:** The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

 **NOTE:** Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

 **NOTE:** Sterilization is not possible on surfaces which have been lubricated and oiled.

The following STERRAD® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- STERRAD® 100S Long Cycle with Boosters
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

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**Hydrogen peroxide (H₂O₂)
sterilization - STERIS® AMSCO®V-
PROTM 1**

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CAUTION: For detailed information, please consult the user handbook of the respective device.

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NOTE: The 'STERIS Low Temperature Reprocessing Device Compatibility Matrix' can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.

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The following STERIS V-PRO® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- V-PRO® maX Flexible Cycle
 - V-PRO® 60 Flexible Cycle
- Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100% ethylene oxide at 55°C and a hold time of 30-45 minutes.

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CAUTION: For gas sterilization with ethylene oxide, the airing times stipulated by the device manufacturer must be observed for instruments, due to gas absorption of the materials. The airing times depend on the processing technology of the EO sterilizer (concentration, process control).

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Reprocessing of flexible accessories (Forceps, scissors etc)

Separate reprocessing instructions are available for the preparation of flexible accessories. These can be downloaded from www.karlstorz.com.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts.

Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

Limitation of liability

KARL STORZ SE & Co. KG shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the conditions that all assembly, operations, system expansions, readjustments, modifications, or repairs to same have been performed by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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This instruction manual is intended to serve as an aid in the proper handling, cleaning, disinfection and sterilization of the Video-Uretero-Renoscope FLEX-X[©] and the flexible Video-Choledochoscope. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for ready reference.

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

- ① Control buttons
- ② Supply connector
- ③ Connection for leakage tester or pressure compensation cap
- ④ Distal tip
- ⑤ Sheath
- ⑥ Handpiece
- ⑦ Working channel with LUER-Lock connectors
- ⑧ Deflection lever
- ⑨ Biopsy forceps 27023 ZE
- ⑩ Leakage tester 13242 XL
- ⑪ Pressure compensation cap 11025 E
- ⑫ LUER adaptor
- ⑬ Single-use cleaning brush 110931-50

The set 11292 VSK/VSUK comprises:

- 11025 E Pressure compensation cap
13242 XL Leakage tester
11293 F Applicator and guiding tube
11003 KD Biopsy forceps, flexible, 3 Fr.
11003 KC Grasping forceps, flexible, 3 Fr.
27023 TD Stone basket, sterile, for single use,
2.5 Fr., 120 cm
27014 Y LUER adaptor, with seal
27550 N Seal, for working channel,
pack of 10
27001 RA Cleaning adaptor
110931-50 CleaningBrush, length 125 cm,
brush diameter 1.7 mm,
package of 50

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Important note for users of KARL STORZ devices and instruments	I
Pictures of the instruments.....	II
Controls, displays, connectors and their uses.....	III
Safety information	3
Intended Use	3
Indications	4
Contraindications.....	4
User qualification	5
Safety precautions when using the Videoendoscopes	6
Warnings and cautions.....	6
Handling	10
Unpacking the device.....	10
Basic equipment	10
Operating conditions	10
Electromagnetic Compatibility.....	10
Handling.....	11
Operating	11
Holding the Videoendoscope correctly	11
Handling the sheath	11
Handling the deflection mechanism.....	12
Handling the accessories (forceps, cleaning brushes, HF electrodes etc.) ..	14
Compatible camera control units	15
Connection to the camera control unit	16
Control buttons	17
Preparation prior to use	18
Inspecting the surfaces.....	18
Checking the deflection	19
Inspecting the LED function.....	20
Inspecting the camera control unit.....	20
Inspecting the instruments and the instrument channel	21
Cleaning, disinfection, care and sterilization.....	22
General warnings	22
Accessories	23
Disassembly	24
Precleaning directly after use	24

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Leakage test	24
Manual cleaning	26
Manual disinfection	27
Drying	27
Machine cleaning and disinfection	27
Manual preparation	27
Machine cleaning/chemical-thermal disinfection	28
Assembly, inspection and care	28
Packaging systems	29
Sterilization	29
Limits of reprocessing	33
Storage	33
Maintenance	34
Servicing and repair	34
Maintenance	34
Limitation of liability	34
Warranty	35
Directive compliance	35
Important information	36
Disposal	36
Appendix	37
Troubleshooting	37
Subsidiaries	40

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Intended purpose

Flexible video choledochoscopes serve to illuminate the bile duct and create a working or irrigation channel for diagnostic and surgical choledochoscopy.

Flexible video choledochoscopes are designed for short-term use in endoscopic interventions performed through natural orifices or surgical invasive interventions.

Indications

- Cholecystectomy
- Suspected gall stones
- Cholelithiasis

Contraindications

- Known coagulopathies
- Uremia
- Intrahepatic gallbladder
- Liver abscesses
- Purulent or necrotic cholecystitis
- Gallbladder neck abscess or perforation
- Serious, acute pancreatitis
- Multiple preceding laparotomies
- Morbid obesity
- Liver cirrhosis
- Portal vein thrombosis

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It is not intended as an introduction to the technique of renaloscopy or choledochoscopy.

Safety precautions when using the Videoendoscopes

The Video-Uretero-Renoscopes FLEX-X^c and the flexible Video-Choledochoscopes must be used according to the recognized medical rules and procedures of endoscopy.

NOTE: The patient leakage current limit of the CF applied parts may be exceeded if CF applied parts (application on heart) are used in combination with other BF applied parts.

In this case, the patient leakage current should be reduced by using only CF applied parts as far as this is possible.

NOTE: The video endoscope does not offer any electric insulation against high frequency (HF) currents. Protection against high frequency currents should be ensured through the use of correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF applied part for more information.

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Handling

NOTE: Damage to the Videoendoscope resulting from incorrect operation is not covered by the manufacturer's warranty.

WARNING: The Videoendoscope is delivered unsterilized, so it has to be cleaned, disinfected/sterilized prior to initial use and before each subsequent reuse.

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Operating

The following instructions outline the correct operation of the Videoendoscope and the corresponding accessories. This is not an introduction to the techniques of renoscopy or choledochoscopy. Instructions on endoscopic technique can be found in the relevant medical literature.

Suitable accessories for your flexible Video-Choledochoscope, such as biopsy forceps, can be found in the LAPAROSCOPY catalog or extracts from this catalog.

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Handling the sheath

CAUTION: Do not carry the Videoendoscope by the sheath. Do not pull, clamp or twist the sheath.

CAUTION: The sheaths of the Videoendoscope must never be kinked or coiled up too tightly. It is best to suspend the Videoendoscopes for storage.

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Observe the lumen the whole time during use. If you lose your orientation, pull the sheath back slightly in the neutral position until vision is restored.

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WARNING: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the Videoendoscope must not be used.

WARNING: If a functional defect occurs during use on the patient, discontinue use immediately. Put the distal tip of the Videoendoscope into the neutral, non-deflected position and remove the Videoendoscope slowly and carefully from the patient.

WARNING: The lens at the distal end of the Videoendoscope may become detached, e.g. due to mechanical forces. For this reason, it is essential to check that the lens is undamaged and present before and after every cleaning.

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CAUTION: Do not use force to turn the deflection lever, to avoid either injuring the patient or damaging the instrument.

The distal tip can be deflected 270° upward or downward by continuously actuating the deflection lever slowly and carefully. Before use, the user must always have full control over the deflection and must develop a sure feel for the deflection mechanism and in particular for the neutral position.

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Preparation prior to use


i **NOTE:** Before the device is used for the first time, it must be subjected to a complete reprocessing cycle (excluding preliminary cleaning) in accordance with the risk classification.

! **WARNING:** Every time the Videoendoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged Videoendoscopes or damaged accessories must not be used.

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i **NOTE:** This thin layer of grease not only serves to lubricate the metal surfaces which are in contact with one another, but also makes the stopcocks leakproof. Oil is not suitable as a sealant.

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Inspecting the surfaces

! **WARNING:** The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

! **CAUTION:** Do not squeeze the distal tip of the Videoendoscope too tightly.

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CAUTION: When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological

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- !** range of action of the chemicals used.
- !** **CAUTION:** Danger of damage to medical devices: The use of chemicals that have not been approved by KARL STORZ or chemical manufacturers may cause damage to the medical devices. Only use chemicals approved for the medical device for reprocessing. A helpful guide is available on the Internet at www.karlstorz.com.
- !** **CAUTION:** National laws and regulations must be observed.
- !** **CAUTION:** The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.
- !** **CAUTION:** Never clean flexible endoscopes with metal brushes.
- !** **CAUTION:** Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, is possible.
- !** **CAUTION:** Flexible endoscopes must never be steam sterilized (autoclaved).
- !** **CAUTION:** To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.
- !** **NOTE:** The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.
- !** **NOTE:** The valves must be removed in their individual parts prior to cleaning and disinfection.

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Accessories

Accessories required for carrying out reprocessing: pressure compensation cap, manual leakage tester.

Place of application	Art. number brush/ adaptor
Inlet and outlet nozzle	27650 B, 39107 ALK
Working channel	110931-50
Clean surface	27652
Vent connection	13242 XA, 11025 E

Disassembly

The valves must be removed in their individual parts prior to cleaning and disinfection.

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NOTE: To avoid damaging the sheath, do not use pushing and pulling movements.

NOTE: Do not use disinfectants for precleaning; these cause fixation of proteins and organic residues.

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Leakage test

The endoscope should be tested to ensure it is leak-proof after use and before reprocessing.

Performing the leakage test – dry

The connection tube of the leakage tester and the vent port on the endoscope must be dry. Connect the leakage tester connector securely to the supply connector vent port and lock by turning the bayonet catch 25° clockwise.

NOTE: The rubber seal in the connection adaptor of the leakage tester should be treated regularly with instrument oil.

Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal end may expand slightly. Using the deflection mechanism on the control body, move the deflection part in all possible directions as far as possible and observe the pointer on the manometer. The test pressure can drop by up to 10 mmHg until the endoscope is completely filled with air.

After stabilizing the pressure on manometer, wait at least 2 minutes to see whether the pressure remains constant or drops.

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NOTE: To avoid damage from the penetration of liquid, a test in the dry state must be performed in advance.

NOTE: Any lumina must be covered with no air bubbles by filling in a targeted manner.

Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in 'Leakage test – dry' and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer does not show a fall in pressure, the device is leak-proof.

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! NOTE: Never release the pressure in the liquid.

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! CAUTION: In this case, fill out the accompanying repairs form, mark the endoscope, which has been cleaned by wiping to remove organic soiling/residue, as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the transport protection sheet 13990 SFN.

Manual cleaning

The endoscope must be completely immersed in a cleaning solution. The working channel must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The channels are brushed until the brush head emerges and is visibly clean. Never brush back and forth within the channel.

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Manual preparation

The endoscope must be completely immersed in a cleaning solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The channels are brushed until the brush head emerges and is visibly clean. Never brush back and forth within the channel.

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CAUTION: The pressure compensation cap must be in place on the vent port for cleaning and disinfection.

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Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/equipment to the sender.

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Disposal

To dispose of Videoendoscopes, no special measures are necessary. National laws and regulations must be observed.

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