STORZ

KARL STORZ—ENDOSKOPE

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UM100 PowerUnit / 28200MX MultiDrive-X ARTHRO

INSTRUCTION MANUAL

UM100 PowerUnit / 28200MX MultiDrive-X ARTHRO

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Thank you very much for your expression of trust in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern high-quality piece of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, connection, and operation of the motor system. The 28200MX MultiDrive-X ARTHRO handpiece in combination with the PowerUnit UM100 can be used with the following products:

- Attachment 28200JB
- Attachment 28200QS
- Attachment 28200QB
- Attachment 28200SA

Please read these instructions carefully. Keep this manual in a convenient and conspicuous location and in its protective casing close to the device.

Users of the various application instruments should familiarize themselves with the operating instructions included with the respective devices.

Rx only Federal (USA) law restricts this device to sale by or on the order of a physician.

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




It is recommended to check the suitability of the product for the intended procedure prior to use.

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2 Controls, displays, connectors, and their uses

PowerUnit UM100

- ① ON/OFF button
(function operational/standby)
– Operational: green
– Standby: orange
- ② Connection socket for MultiDrive-X ARTHRO
- ③ Potential equalization
- ④ Power socket
- ⑤ Fuse holder

MultiDrive-X ARTHRO 28200MX

- ⑥ Receptacle for pin driver, chuck/saw
attachment
- ⑦ Release buttons for pin driver, chuck/saw
attachment
- ⑧ Switch lever for mode selection
- ⑨ Speed button
- ⑩ Oscillation button
- ⑪ Detachable connecting cable 28200CC
- ⑫ Connecting cable 28200CC release button
- ⑬ Cleaning adaptor 28200RA1
- ⑭ Pin driver
28200QS / 28200QB
- ⑮ Cleaning adaptor 28200RA2
- ⑯ Jacobs chuck 28200JB
- ⑰ Key for Jacobs chuck 28200JBS
- ⑱ Cleaning adaptor 28200RA3
- ⑲ Oil spray adaptor 28200RA4 for pin driver and
Jacobs chuck
- ㉑ Sagittal saw attachment 28200SA
- ㉒ Oil spray adaptor 28200RA5 for the sagittal
saw attachment

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Connection socket for MultiDrive-X ARTHRO

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This device has been labelled in Accordance with the European Directive on waste electrical and electronic equipment (WEEE).

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Deactivated (lever ⑧ in middle position)

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1	Images of the equipment.....	IV
2	Controls, displays, connectors, and their uses.....	V
3	Symbols employed.....	VI
3.1	Symbols on the device	VI
3.2	Label and packaging symbols.....	VII
4	General information	3
4.1	Description of the system	3
4.2	Special features	4
5	Safety instructions	5
5.1	Explanation of warnings and cautions	5
5.2	Normal use	11
5.3	Intended use – PowerUnit	12
5.4	Intended use, indication/ contraindication – MultiDrive-X ARTHRO and accessories.....	12
5.4.1	Multifunction handpiece/ drill handpiece/drive incl. attachments.....	12
5.4.2	Drill/hollow-core drill/burr.....	13
5.4.3	Drill wire/Kirschner wire	13
5.4.4	Saw attachments.....	13
5.5	User qualifications.....	13
5.6	User profile (specialist physician).....	14
5.7	User profile (assistant).....	14
5.8	Patient profile	14
5.9	Affected body part or tissue type	15
5.10	Use.....	15
5.11	Safety precautions at the site of installation	16
5.12	Safety precautions when operating the system.....	16
6	Operating instructions	17
6.1	Unpacking.....	17
6.2	Basic equipment	17
6.3	Installing and connecting the unit	18
6.4	Connecting the multifunction handpiece.....	19
6.5	MultiDrive-X ARTHRO.....	20
6.6	Initial operation	22
6.6.1	Regulating the speed.....	22
6.6.2	Switching modes	22
6.6.3	Switching the function	23
6.7	Decommissioning	23
7	Maintenance.....	24
7.1	Fuse replacement	24
7.2	Reprocessing of the device	25

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7.2.1	General warnings	25
7.2.2	Manual wipe-down disinfection	26
7.3	Reprocessing of the multifunction handpiece.....	27
7.3.1	General warnings.....	27
7.3.2	Accessories	28
7.3.3	Disassembly.....	29
7.3.4	Preparation for cleaning and disinfection	29
7.3.5	Manual precleaning.....	30
7.3.6	Machine cleaning and disinfection	31
7.3.7	Assembly, inspection and care.....	32
7.3.8	Packaging systems.....	33
7.3.9	Sterilization	33
7.3.10	Reprocessing limits.....	33
7.4	Maintenance and safety check.....	34
7.4.1	Maintenance	34
7.4.2	Safety check	34
7.5	Service and repair.....	36
7.6	Disposal	36
7.7	Repair program.....	37
7.8	Infection prevention.....	37
7.9	Limitation of liability	38
7.10	Manufacturer's warranty	38
8	Technical description	39
8.1	General	39
8.2	Troubleshooting	40
8.3	Technical Data	42
8.4	Performance parameters MultiDrive-X ARTHRO.....	43
8.5	Standard conformity	44
8.6	Directive compliance.....	44
8.7	Technical documentation.....	45
9	Spare parts, recommended accessories	46
9.1	Spare parts	46
9.2	Recommended accessories.....	46
9.3	Accessories for MultiDrive-X ARTHRO 28200MX	46
9.4	Accessories for pin driver.....	47
9.5	Accessories for Jacobs chuck	47
9.6	Accessories for sagittal saw attachment 28200SA	49
10	Electromagnetic Compatibility (EMC) Information	50
11	Subsidiaries.....	67

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4.1 Description of the system

The MultiDrive-X ARTHRO can be used with:

- PowerUnit UM100
- Jacobs chuck 28200JB
- Pin driver 28200QS, for drilling wires dia. 0.6...1.8 mm
- Pin driver 28200QB, for drilling wires dia. 1.8...3.2 mm
- Sagittal saw attachment 28200SA

Das MultiDrive-X ARTHRO system is used in combination with a saw blade, drill and larding wire (see Accessories page 46-49) to perform surgical interventions under endoscopic control as well as open surgery in the field of arthroscopy.

With the corresponding instrument, the system is ideal for the controlled cutting and drilling of bone and tissue during arthroscopic interventions on the knee, shoulder, ankle, elbow, wrist and hip.

The PowerUnit UM100 provides the MultiDrive-X Arthro handpiece with the necessary voltage. The PowerUnit is connected to the handpiece with a connecting cable.

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4.2 Special features

- The multifunction handpiece MultiDrive-X ARTHRO can be operated in 3 modes:
'Clockwise/counterclockwise rotation'
'Oscillation'
'Precision' mode (with reduced speed).
- In Normal mode the rotation speed of the MultiDrive-X ARTHRO is up to 1,500 rpm.
- In Oscillation mode the MultiDrive-X ARTHRO has an oscillation frequency of approx. 3.5 Hz vis-à-vis the speed and generates the oscillation cycles accordingly.
- In 'Precision' mode the rotation speed of the MultiDrive-X ARTHRO is up to 800 rpm. Oscillation mode is identical to Normal mode.

*** NOTE:** KARL STORZ recommends the 'Precision' mode when using the sagittal saw attachment 28200SA.

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5 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

5.1 Explanation of warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words Warning, Caution and Note convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved.

CAUTION: A Caution indicates that particular maintenance procedures or safety precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information.

WARNING: Read this instruction manual thoroughly before operating the system. Read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

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Safety instructions Warnings and cautions

-  **WARNING:** Continuous checks must be performed both before and during use to establish whether the required channel and/or mode is set.
-  **WARNING:** The MultiDrive-X ARTHRO may only be used by persons with an appropriate medical qualification and who are acquainted with the technique.
-  **WARNING:** When replacing fuses, use only fuses of the same rating.
-  **WARNING:** The device may only be operated at the voltage stated on the rating plate.
-  **WARNING:** If the motor does not stop as required, switch the unit off.
-  **WARNING:** Store/keep the instrument at a safe distance from the patient when not in use.
-  **WARNING:** The MultiDrive-X ARTHRO must only be used in medical facilities that conform to applicable national, state and local electrical codes.
-  **WARNING:** Keep out of reach of patients.
-  **WARNING:** The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.
-  **WARNING:** Combinations of medical devices are only assured to be safe if:
 - they are identified as such in the respective instruction manuals or
 - the intended use and interface specifications of the devices used in combination permit this (cf. Clause 16 of IEC 60601-1).

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WARNING: Connect the device to a properly grounded 'Hospital grade' or 'Hospital use' electrical outlet. Check grounds regularly. Routinely inspect the power cord and plug and do not use if the plug or the cord is damaged.

WARNING: To avoid the risk of an electric shock, this equipment must only be connected to a supply mains with protective earth.

WARNING: Test this equipment prior to each surgical procedure to ensure that it functions correctly.

WARNING: Do not use in the presence of flammable anesthetics. There is a risk of explosion.

WARNING: Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.

WARNING: Risk of injury: Incorrectly assembled and damaged instruments can lead to burns or other injuries to the patient or user. Instruments and all of the accessories used in combination with them must be checked and if necessary replaced immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

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 **WARNING:** Risk of infection: These instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination is an indication of no preparation or incorrect preparation. Prepare the instruments before initial use and before and after every subsequent use using validated procedures.

 **WARNING:** Risk of injury: Incorrect application of medical instruments or bending of these instruments by applying excessive side load force poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

 **WARNING:** Risk of injury and risk of damaging the products: When using applied parts outside of the field of vision there is a risk of unintentionally injuring tissue and damaging accessories. Always keep the applied parts of the active electrode and the laser and other instruments which transmit energy pointed at the target in the field of visions during use.

 **WARNING:** The MultiDrive-X ARTHRO in combination with the PowerUnit and the attachments is only designed for intermittent operation and not for continuous operation. (see Technical Data on page 42). For this reason, always turn the motor off again, as it will otherwise become too hot.

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 **WARNING:** Risk of infection: Incorrectly prepared instruments expose patients and users to a risk of infection as well as the danger of an instrument malfunction. Observe the 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

 **NOTE:** The 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' instructions can be downloaded or requested at www.karlstorz.com.

 **WARNING:** When using HF surgical devices, the shaver handle and shaver blade must be put out of reach of the patient.

 **WARNING:** Always unplug the unit before performing any maintenance, including cleaning.

 **WARNING:** To avoid the risk of electrical shock, do not open the device. Refer servicing to the manufacturer or qualified personnel that are authorized by the manufacturer. Removal of the cover by non-authorized personnel will void the warranty.

 **CAUTION:** The MultiDrive-X ARTHRO may only be used with accessories, and disposable items which have been designated for use by KARL STORZ.

 **CAUTION:** Do not store liquids on or above the device. Do not allow fluids to get into the unit.

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CAUTION: If the motor becomes blocked for any length of time, it can lead to overheating and cause damage to the motor. If the blockage cannot be disengaged using oscillating operation (if possible), the motor must be shut off and the blockage must be removed manually.

WARNING: Only use accessories that are suitable for use with MultiDrive-X ARTHRO and can be operated in accordance with the generally accepted rules of technology and occupational safety.

NOTE: A defective sterile glove can result in contamination if it comes into contact with sterile goods and can make replacement necessary.

NOTE: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

NOTE: For questions on material compatibility with chemicals for reprocessing, please contact the Hygiene Customer Service at hygiene@karlstorz.com.

NOTE: If the unit is not functioning as required, please contact the KARL STORZ Service.

NOTE: This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap. Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your distributor for information on your local collection point. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this unit.

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Before using the MultiDrive-X ARTHRO motor system on a patient, it is imperative that the instructions in this manual are understood and are followed carefully.

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5.2 Normal use

The MultiDrive-X ARTHRO motor system is universally applicable in conjunction with the respective, suitable application instrument in the orthopedics field of application. The detailed 'Intended use' of the application instrument can be found in the instructions for use provided for the respective application instrument.

The MultiDrive-X ARTHRO motor system may only be used with accessories, wearing parts and disposable products that are designated by KARL STORZ as suitable for the unit.

The MultiDrive-X ARTHRO motor system is used in conjunction with a shaver blade or burr to perform surgical interventions under endoscopic control in the field of arthroscopy.

With the corresponding blade, the system is designed to provide controlled cutting, shaving, and abrading of bone and tissue during arthroscopic interventions on the knee, shoulder, ankle, elbow, wrist, and hip.

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Use of the MultiDrive-X ARTHRO in fields other than those indicated above is not permitted.

The MultiDrive-X ARTHRO may only be used with accessories, replacement parts, and disposable items which have been designated as suitable for the unit by KARL STORZ.

Unauthorized conversions or modifications to the system are not permitted.

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5.3 Intended use – PowerUnit

The PowerUnit UM100 acts as a power supply unit to provide energy to the MultiDrive-X ARTHRO handpiece 28200MX.

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5.4 Intended use, indication/ contraindication – MultiDrive-X ARTHRO and accessories

5.4.1 Multifunction handpiece/ drill handpiece/drive incl. attachments

5.4.1.1 Intended use

The multifunction handpieces serve to accommodate and drive rotating and oscillating instruments by means of mechanical transmission and conversion of speed and torque, as well as the rotational direction for drill and saw applications. Multifunction handpieces are non-invasive and are designed for transient use in surgical invasive interventions.

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5.4.1.2 Indications

The MultiDrive-X ARTHRO saw and drill handpiece and corresponding accessories are suitable for all therapeutic interventions which require the energy/media supplied by the device for successful operation.

5.4.1.3 Contraindications

- Contact with the central nervous system
- The accessories for the motor system for musculoskeletal applications must not be used if, in the opinion of the responsible physician, the patient's health is endangered through their use, e.g., due to the patient's general condition.
- Further contraindications relating directly to the product group are not known.

5.4.2 Drill / hollow-core drill / burr

i NOTE: Please follow the instruction manual enclosed with the respective application instrument!

5.4.3 Drill wire / Kirschner wire

i NOTE: Please follow the instruction manual enclosed with the respective application instrument!

5.4.4 Saw attachments

i NOTE: Please follow the instruction manual enclosed with the respective application instrument!

5.5 User qualifications

The MultiDrive-X ARTHRO motor system in conjunction with an application instrument may only be used by physicians and/or medical assistants who have been instructed in the use of this device/instrument.

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5.6 User profile (specialist physician)

- Recognized medical skills in the area of arthroscopy (specialist physician)
- Adequate powers of comprehension to assess the situation rationally at all times during the application
- Adequate language skills in one of the languages used in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic)

5.7 User profile (assistant)

- Completed training as medical assistant
- Adequate language skills in one of the languages used in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic)

5.8 Patient profile

The use of the MultiDrive-X ARTHRO is not limited to a certain patient profile (sex, age, weight, etc.). The attending physician must assess a patient's health with regard to the endoscopic intervention and declare it adequate.

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5.9 Affected body part or tissue type

The area of application is limited to certain joints (feet, knees, hips, shoulders, elbows, hands). The joints are accessed via surgically invasive openings.

5.10 Use

The device may only be used in operating rooms in hospitals and medical practices in the ambient conditions given in the technical data (page 42).

- Not on intensive care stations.
- Not in mobile facilities (ambulances, aircraft, ships or trains).
- Position it on a vibration-free, fixed base.
- Visual distance from the physician to the device when operated: max. 2 m.
- The use of the device requires the division of work.
- Frequency of use up to several times a day.
- Length of use from a few minutes to several hours a day.
- Can be moved only by installing it on an equipment trolley.
- Can be swiveled only by installing it on ceiling supply units.
- The patient must not necessarily be under general anesthesia during the intervention (local anesthesia).
- May be used with other devices (light source, shaver, etc.).
- To be operated in ambient light between 100 lux and 1500 lux.
- Ambient temperature in immediate proximity of the device during operation: max. 30°C.
- Sound pressure from background noise < 60 dBA.
- The device is placed outside the sterile area.

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
5.11 Safety precautions at the site of installation

The MultiDrive-X ARTHRO motor system may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

The PowerUnit is equipped with a connector for attaching a ground wire. Connect this in accordance with the state, local and national electrical codes.

The MultiDrive-X ARTHRO motor system is not intended for use in hazardous areas. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

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5.12 Safety precautions when operating the system

WARNING: The system must be tested for operational safety and proper condition prior to each surgical procedure.

During surgery with the MultiDrive-X ARTHRO, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital signs and the anesthetic.

The MultiDrive-X ARTHRO in combination with the PowerUnit device has been tested and is found to comply with the EMC limit values as per IEC 60601-1-2 CISPR 11 Class B.

Please read the Electromagnetic Compatibility information in section 10 (page 50-66).

WARNING: Do not operate the handpiece unless the saw, drill or larding wire is visible.

WARNING: When using HF surgery, the handpiece and the instrument must not come into contact with the patient.

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6 Operating instructions

6.1 Unpacking

Carefully remove the unit and the accessories from the packaging.

Check for missing items and for evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can be used if the unit has to be transported.

6.2 Basic equipment

- 1 PowerUnit UM100
- 1 Power cord 400A
- 1 Handpiece 28200MX
- 1 Connecting cable 28200CC
- 1 Cleaning adaptor 28200RA1

in combination with optional accessories
Attachments according to the application 28200JB/
QS/QB/SA and instrument working inserts
1 Instruction manual 96156021DER

i NOTE: Optional accessories are described in section 'Recommended accessories'.

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6.3 Installing and connecting the unit

NOTE: The MultiDrive-X ARTHRO with PowerUnit including the accessories connected may be used only in medical rooms with electrical installations conforming to the applicable national electrical regulations.

! WARNING: The MultiDrive-X ARTHRO with PowerUnit is not intended for use in potentially explosive atmospheres. Do not operate the system within the hazard zone shown in the diagram while explosive anesthetic gases are in use (Fig. 1).

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Place the PowerUnit on a level surface. The PowerUnit is equipped with a connector for potential equalization ③. The unit's grounding should be installed by a qualified electrician (Fig. 2).

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Before connecting the PowerUnit to the power supply line, verify that the line voltage stated on the unit's identification plate conforms with that of the power supply lines to be used (Fig. 3).

! WARNUNG: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

! WARNUNG: Only use the power cord delivered by KARL STORZ or a similar power cord which has a national inspection seal.

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Connect the power, insert the power cord end into the receptacle ④ (Fig. 4). The ON/OFF button on the PowerUnit lights up orange. The system is in Standby mode.

WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

WARNING: The device must be set up so as to ensure that the power cord can simply be pulled at any time.

The power cord connector must be pulled out to ensure an all-pole disconnection from the mains power supply.

WARNING: The ON/OFF button of PowerUnit must always be in Standby mode during all manipulations of the handpiece other than the intended activation of the handpiece.

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6.4 Connecting the multifunction handpiece

Prepare the multifunction handpiece.

NOTE: Prior to turning the device on, the user must connect the multifunction handpiece to the device.

Connect the connecting cord from the multifunction handpiece to the connection socket ② (multifunction handpiece – Fig. 5).

Align the red point on the connecting cord with the red point on the connection socket to guarantee a proper connection.

NOTE: To unplug the connecting cable, pull on the 'sliding ring' of the plug — do not pull directly on the cable!

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6.5 MultiDrive-X ARTHRO

The MultiDrive-X ARTHRO (Fig. 6) is supplied with energy by the PowerUnit. Chucks, sagittal saws and pin drivers for drilling wires can be attached to the handpiece (Fig. 7-13).

A list of available accessories is shown in a separate section of this manual.

The respective accessory (attachment) is connected to the handpiece by inserting the attachment into the receptacle ⑥ as far as it will go. Make sure that the pins on the attachment are congruent with the recesses on the handpiece ⑦. Click in the attachment: Turn the attachment radially to the handpiece and simultaneously push together.

NOTE: The attachment can be positioned in 8 positions radially on the handpiece ('distance' in each case 45°).

Then check that the lock of the attachment is secure, to this end pull the attachment axially in the distal direction.

Remove the attachment: Press the release buttons ⑦ at the same time and remove the attachment in an axial direction.

Pin driver 28200QS/QB (Fig. 7)

Wire nails with a diameter from 0.6 to 3.2 mm fit in the cannulated drive of the pin driver.

28200QS: dia. 0.6 up to 1.8 mm

28200QB: dia. 1.8 up to 3.2 mm

Prerequisite:

- The correct pin driver is attached to the handpiece.

1. Push the wire in the pin driver.

2. Pull the handle of the pin driver to the handpiece to secure the wire.

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
Jacobs chuck 28200JB (Fig. 8)

Prerequisite:

- The Jacobs chuck is assembled on the handpiece.

1. Open the Jacobs chuck with the respective key (28200JBS).
2. Insert the drill into the chuck as far as it will go (Fig. 9).
3. Tighten with the key (Fig. 10).

Sagittal attachment 28200SA (Fig. 11)

Prerequisite:

- The sagittal attachment is assembled on the handpiece.

1. Insert the saw blade into the attachment receptacle until it clicks into place whilst the button to open the receptacle is held down (Fig. 12-13).
2. To click in release the button.

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6.6 Initial operation


Switch on the PowerUnit using the ON/OFF button ① (Fig. 14). The button lights up green. The system is ready for use.

NOTE: Prior to turning the device on, the user must connect the MultiDrive-X ARTHRO to the PowerUnit.

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6.6.1 Regulating the speed

Regulate the speed: The MultiDrive-X ARTHRO can control the speed by moving the upper button ⑨. The speed (rpm) is proportional to the position of the upper button. Pushing the upper button all the way in toward the handpiece achieves the maximum speed.

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NOTE: No action can be triggered when the lever is in the 'OFF' position.

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6.6.3 Switching the function

6.6.3.1 Functions in the applications (modes) Normal and Precision

- **Clockwise operation**
Actuating the top button ⑨ causes the handpiece to be operated in a clockwise direction (Fig. 19).
The speed can be changed continuously.
- **Counterclockwise operation**
When both buttons (⑨ and ⑩) are actuated at the same time, the handpiece operates in a counterclockwise direction (Fig. 20).
- **Oscillation** 
When the bottom button ⑩ is actuated, the handpiece operates in Oscillation mode (Fig. 21).

i NOTE: Oscillation mode starts with a brief delay (of 0.4 s); this prevents oscillation from starting unintentionally should the user wish to activate counterclockwise operation (and presses the bottom button first).

The user can select the 'Precision' or 'Normal' mode with the side lever ⑧ (see section 6.6.2).

The Precision mode works at reduced speed with clockwise/counterclockwise rotation (Fig. 17).

The 'Precision' mode is recommended when a more sensitive change of speed is required.

The oscillation frequency is identical to the Normal mode.

As a rule, the 'Precision' mode is recommended when using the sagittal saw attachment.

Please note that in this mode the maximum speed is 800 rpm.

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6.7 Decommissioning

Switch off the PowerUnit on the standby button ① (Fig. 22).

Remove connected accessories from the device and, if necessary, 'reprocess'.

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7 Maintenance

7.1 Fuse replacement

Switch the PowerUnit to standby mode and disconnect from the power supply (Fig. 23).

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Remove the power fuse holder ⑤ with a screwdriver or another suitable tool (Fig. 24).

NOTE: 4 A fuses suitable for operation with a 100...240 V supply are inserted in the device at the factory.

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! WARNING: When replacing fuses, use only fuses of the same rating.

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Reconnect the power.
Perform function testing (Fig. 26).

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7.2 Reprocessing of the device

7.2.1 General warnings

WARNING: Risk of fatal injury from electric shock: Unplug from power supply before every cleaning!

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

WARNING: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

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7.2.2 Manual wipe-down disinfection**CAUTION: Damage to the PowerUnit from liquid: Prevent liquid from entering the PowerUnit.**

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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***NOTE:** Please read the instruction manual
provided with the application instrument
carefully!

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- This reprocessing step is part of the validated
procedure.

7.3 Reprocessing of the multifunction handpiece

7.3.1 General warnings

! WARNING: Risk of infection: Incorrectly
reprocessed medical devices expose
patients, users and third parties to a risk
of infection as well as the risk that the
medical device may malfunction. Observe the
'Cleaning, Disinfection, Care, and
Sterilization of KARL STORZ Instruments'
instructions and the accompanying
documentation.

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WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

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CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: National laws and regulations must be observed.

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7.3.2 Accessories

Accessories required for carrying out reprocessing:

- Brush: 27652
- Irrigation adaptor 28200RA1 for multifunction handpiece 28200MX
- Irrigation adaptor 28200RA2 for quick-release attachments 28200QS / QB
- Irrigation adaptor 28200RA3 for the Jacobs chuck 28200JB

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
7.3.3 Disassembly

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible. The attachments and connecting cable must be removed for reprocessing (see Figs. 27-29).

7.3.4 Preparation for cleaning and disinfection

! CAUTION: Risk of burns: Dried residue can cause malfunctions and generate heat. The multifunction handpiece and the reusable accessories should be subjected to a reprocessing process immediately after use in order to avoid the drying out of tissue remnants and fluids.

! WARNING: Reprocessing is only permitted for reprocessable products. If you reprocess disposable products, the product safety cannot be guaranteed as a result of material changes. This can result in injury to the patient. Dispose of disposable products after use. Do not reprocess them.

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the multifunction handpiece by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

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7.3.5 Manual precleaning Brushing the surfaces

Any visible contamination or heavy soiling must be removed from the surfaces of the individual components by cleaning them under cold running water with the aid of a brush or sponge (see Fig. 30 and 31).

Rinsing out with a pressure water gun

The lumina and hollow spaces of the medical devices must be rinsed using a pressure water gun and a suitable rinsing attachment. For rinsing, use pulses of at least 3 x 5 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. no. 27660; see Fig. 32 and 33) is ideal for this purpose.

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7.3.6 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A₀ value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

ⓘ NOTE: If necessary, the instrument must be dried off afterwards by hand.

Connecting

In order to ensure effective machine cleaning and disinfection, the instruments must be connected up to the washer and disinfectant.

- Attach the cleaning adaptor 28200RA1 to the handpiece 28200MX (see Fig. 34).
- Attach the cleaning adaptor 28200RA3 to the Jacobs chuck 28200JB (see Fig. 35).
- Attach the cleaning adaptor 28200RA2 to the pin driver (see Fig. 36).

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40

7.3.7 Assembly, inspection and care

The cleaned and disinfected medical devices must be visually inspected to ensure that they are clean, complete, dry and have not been damaged:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.

Greasing attachments

- Attach the oil spray adaptor 28200RA5 to the sagittal saw attachment (see Fig. 37).
- Place the spray tube of the oil spray on the adaptor 28200RA5. Grease the sagittal saw attachment by means of 3 x 1 second long spray intervals (see Fig. 38).
- Attach the oil spray adaptor 28200RA4 to the Jacobs chuck or the pin driver (see Fig. 39).
- Place the spray tube of the oil spray on the adaptor 28200RA4. Grease the Jacobs chuck or pin driver by means of 3 x 1 second long spray intervals (see Fig. 40).
- Remove excess oil.

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7.3.9 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and following consultation with the device manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

Steam sterilization in the fractionated prevacuum procedure

For sterilization of the medical device in its fully assembled state, a fractionated prevacuum procedure (DIN EN ISO 17665-1) at 134°C–137°C with a minimum exposure time of 3 minutes and a maximum of 18 minutes must be used. This procedure is only suitable for thermostable instruments.

WARNING: *Risk of burns! Parts are hot after steam sterilization. Allow them to cool down!*

WARNING: *KARL STORZ recommends that, in order to assure sterility, the handpiece, cables and attachments should be sterilized in a disassembled state and then reassembled in a sterile field.*

7.3.10 Reprocessing limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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7.4 Maintenance and safety check

7.4.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

7.4.2 Safety check

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/ repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: For limit values, please refer to the current standard.
- Measure earth leakage current in accordance with IEC 62353: For limit values, please refer to the current standard.
- Measure touch current according to IEC 62353: For limit values, please refer to the current standard.
- Measure patient leakage current according to IEC 62353: For limit values, please refer to the current standard.

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Test for proper functioning

The test for proper operation must include the following tests at least:

- Test for proper functioning of the ON/OFF button
- Functional test of handpiece
- Change of direction of rotation

Documentation

The safety check and the results must be recorded in the medical device logbook.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during this safety check/repeat inspection, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety checks, please refer to the latest issue of the relevant Service Manual.

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7.5 Service and repair

Repairing defective devices is only to be performed by personnel that have been authorized by KARL STORZ and using original KARL STORZ parts.

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7.6 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this unit.

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7.7 Repair program

For devices, individual repair is necessary. Usually to bridge the repair period, you will receive a unit on loan, which you will then return to KARL STORZ as soon as you receive the repaired unit. In the case of repairs, please contact your local KARL STORZ subsidiary or authorized dealer. In Germany, you can refer the repairs directly to:
KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

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7.8 Infection prevention

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender. Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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7.9 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

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7.10 Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see "Subsidiaries" page), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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8 Technical description

8.1 General

In the event of the device malfunctioning, the troubleshooting list provided in section 8.2 will assist the user in finding the possible cause.

WARNING: High voltage is present in the unit's interior. The unit must only be opened by an authorized KARL STORZ customer service technician.

WARNING: Maintenance work that is not mentioned in this section must only be performed by technical personnel that are authorized by KARL STORZ.

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8.2 Troubleshooting

Problem Attachment has been released unintentionally Possible cause 'Manual' unlocking Recommended remedy Lock the attachment in place again. Functional test before use.
Problem Instrument totally inoperative. Possible cause Outage in the electrical supply. Recommended remedy Check the electrical supply.
Problem Complete device failure. Recommended remedy Notify customer service.
Problem Power switch lamp does not light up. Possible cause 1 Power supply failure. 2 Defective power fuse. 3 Power switch lamp defective. Recommended remedy 1 Check that there is electricity to the wall outlet. 2 Change fuses as described in the instruction manual. Make sure you fit the correct type of fuse. 3 Notify customer service.

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Problem

- 1 Power switch lamp lights up.
- 2 No indication of rotation.
- 3 No energy output when operating the buttons.
- 4 Speed cannot be regulated.

Possible cause

- 1 Faulty power supply unit.
- 2 Connection to motor interrupted.
- 3a Faulty handpiece.
- 3b Handpiece blocked or jammed.
- 4 Faulty handpiece.

Recommended remedy

- 1, 2, 4 Notify customer service.
- 3a Connect another handpiece.
- 3b Clean the handpiece.

Problem

Increased operating noise at the handpiece

Possible cause

Wear.

Recommended remedy

The intervention can be completed. Send in the handpiece for maintenance.

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8.3 Technical Data

PowerUnit	UM100
Line voltage	100...240 V~
Power frequency	50/60 Hz
Line fuse	2 x T 4.0 AH 250 V, IEC 127 format
Power consumption	150...200 VA
Dimensions	305 x 165 x 280 mm (W x H x D)
Weight	2.1 kg
Operating mode: periodic duty (in connection with application instrument*)	
* Multifunction handpiece	
Duty cycle	ON t = 3 min, OFF t = 2 min (total time 10 min)
Operating conditions	
Temperature	10 °C...30 °C
Humidity	30% ... 70% (non-condensing)
Atmosph. pressure	700 hPa...1080 hPa
Storage/transport conditions (UM100)	
Temperature	-5 °C...60 °C
Relative humidity	5%...70% (non-condensing)
Atmosph. pressure	700 hPa...1080 hPa
Storage/transport conditions (28200MX)	
Temperature	-20 °C...60 °C
Relative humidity	5%...95% (non-condensing)
Atmosph. pressure	500 hPa...1080 hPa
Average lifecycle	
PowerUnit	5 years
MultiDrive-X ARTHRO	2 years/max. 300 reprocessing cycles
IP protection class	IP 21 (UM100)
IP protection class	IP X4 (28200MX)

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Designation	Min. speed [rpm]	Max. speed [rpm]	Accuracy	Operating mode	Adjustable increments
Multifunction handpiece Normal mode	50	1,500	+/- 10%	rotating	~5 rpm
Multifunction handpiece Precision mode	50	800	+/- 10%		
Multifunction handpiece (Normal mode and Precision mode)	3.5 Hz		+/- 10%	oscillating	

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(for UM100/28200MX, 28200CC)

According to IEC 60601-1, UL 60601-1,
CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electrical shock:
Protection class I
- Degree of protection against electric shocks:
Applied part type BF
- Type of moisture protection (UM100):
drip-water protection as per IP 21
- Type of moisture protection
(28200MX/28200CC):
drip-water protection as per IP X4

Conforms with IEC 60601-1-2:

Please read the Electromagnetic Compatibility Information in the section 10 (p. 50-66).

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(for UM100/28200MX, 28200CC)

Medical product of class IIa

This medical product bears the CE mark in
accordance with the Medical Device Directive
(MDD) 93/42/EEC.

NOTE: The code number after the CE mark
indicates the responsible notified body.

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8.7 Technical documentation

On request, the manufacturer will provide available parts lists, descriptions, instructions for adjustments, and other documentation deemed useful for suitably qualified user personnel authorized by the manufacturer when repairing device components that have been designated as repairable by the manufacturer.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

i NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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Spare parts, recommended accessories

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9 Spare parts, recommended accessories

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Part	Order number
Jacobs chuck, up to dia. 6.5 mm	28200JB
Sagittal Saw Attachment, key- less, with pushbutton, for use with MultiDrive-X ARTHRO 28200MX and Saw Blades 28207SLA – SLH	28200SA
Pin Driver, for drilling wires diameter 0.6 – 2.0 mm	28200QS
Pin Driver, for drilling wires diameter 2.0 – 3.4 mm	28200QB

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9.3 Zubehör für MultiDrive-X ARTHRO 28200MX

9.3 Accessories for MultiDrive-X ARTHRO 28200MX

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Part	Order number
Cleaning adaptor, for use with MultiDrive-X ARTHRO handpiece 28200MX	28200RA1
Connecting cable MultiDrive-X ARTHRO, for use with MultiDrive-X ARTHRO and PowerUnit UM100	28200CC

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9.4 Accessories for pin driver

Part	Order number
Drilling Wire , spiral shape, with eyelet, diameter 2.4 mm, length 38 cm, for single use, for use with Bone Drills 28729BA – BH, Collar Burrs 28729BKC – BKF/BLC – BLF and Drills 28729GA – GE	28729D
Larding Wire , pyramidal tip, diameter 2.4 mm, length 32 cm, unsterile, for single use, package of 10, for use with Undercut Drills 28729BA – BH, Collar Burrs BLC – BLF and Drills 28729GA – GE	28729E
Cleaning adaptor , for use with pin driver 28200QS/QB	28200RA2

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**Spare parts,
recommended accessories**

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Collar Burr , tibial, diameter 6 mm, cannulated, collar diameter 7.5 mm x 43 mm, for use with Drilling Wire 28729D and Larding Wire 28729E	28729BLC
Same , diameter 7 mm, collar diameter 8.5 mm	28729BLD
Same , diameter 8 mm, collar diameter 9.5 mm	28729BLE
Same , diameter 9 mm, collar diameter 10.5 mm	28729BLF
Same , diameter 6.5 mm, collar diameter 8 mm	28729BLK
Same , diameter 7.5 mm, collar diameter 9 mm	28729BLL
Same , diameter 8.5 mm, collar diameter 10 mm	28729BLM
Collar Burr , femoral, diameter 6 mm, cannulated, collar diameter 7.5 mm at 23 mm, for use with Drilling Wire 28729D and Larding Wire 28729E	28729BKC
Same , diameter 7 mm, collar diameter 8.5 mm	28729BKD
Same , diameter 8 mm, collar diameter 9.5 mm	28729BKE
Same , diameter 9 mm, collar diameter 10.5 mm	28729BKF
Same , diameter 6.5 mm, collar diameter 8 mm	28729BKK
Same , diameter 7.5 mm, collar diameter 9 mm	28729BKL
Same , diameter 8.5 mm, collar diameter 10 mm	28729BKM
Drill , cannulated, diameter 6 mm for use with Drilling Wire 28729D and Larding Wire 28729E	28729GA
Same , diameter 7 mm	28729GB
Same , diameter 8 mm	28729GC
Same , diameter 9 mm	28729GD
Drill , cannulated, diameter 10 mm, for use with Drilling Wires 28729D/E	28729GE

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Drill , for preparation of the tunnel for BIOPLUG 2870411BP, with laser marking for the depth of the drill hole, length 7 cm	28179KB
Same , length 15 cm	28179TB
Same , for BIOPLUG 2870514BP, length 18 cm	28179TC
Same , for BIOPLUG 2870310BP, length 15 cm	28179TD
Drill , for preparation of the tunnel for MPFL reconstruction, with penetration stop, diameter 4.5 mm, for use with Drill 28728BL	28728BLB
Cleaning adaptor , for use with jacob's chuck 28200JB	28200RA3

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**9.6 Accessories for sagittal saw
attachment 28200SA**

Part	Order number
Saw blade , blade width 7 mm, blade thickness 0.4 mm, cutting thickness 0.6 mm, working length 15 mm, sterile, package of 5	28207SAA
Same , blade width 9.8 mm, working length 18.5 mm	28207SAB
Same , blade width 9.8 mm, working length 25 mm	28207SAC
Same , blade width 7 mm, working length 29 mm	28207SAD
Same , blade width 9.5 mm, working length 31 mm	28207SAE
Same , blade width 5.5 mm, working length 18 mm	28207SAF
Saw blade , blade width 12 mm, blade thickness 0.5 mm, cutting thickness 0.7 mm, working length 34 mm, sterile, package of 5	28207SAG
Same , blade width 5.5 mm, working length 18 mm	28207SAH
Same , blade width 9.5 mm, working length 31 mm	28207SAI

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Electromagnetic Compatibility (EMC) Information

10 Electromagnetic Compatibility (EMC) Information

The PowerUnit UM100 in combination with MultiDrive-X ARTHRO is suitable for use in Professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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Accessories and cables used for EMC Compliance

Cable type	Shielded	Length [m]	Ferrite	Used for
PE (potential equalization)	No	> 3	No	Potential equalization
Power Cord	No	3	No	Power supply (connection of device to AC Mains)
Handpiece MultiDrive-X ARTHRO	No	3.5	No	Energy transfer to motor

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! WARNING: The use of an accessory or cable with the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 other than those specified in this manual may result in increased emissions or decreased immunity of the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100. When an accessory or cable other than those specified in this manual is used with the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100, it is the responsibility of the operator when using the device to determine compliance of the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 with IEC 60601-1-2.

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! WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

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i NOTE: The emission characteristics of this equipment make it suitable for use in professional healthcare environment as well as residential environment (CISPR 11 Class B). This equipment offers adequate protection to radio communication service. In the rare event of interference to the radio communication service, the user might need to take mitigation measures, such as relocating or re-orienting the equipment.

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Table 1 – Compliance level for Immunity tests			
Guidance and manufacturer's declaration – electromagnetic immunity			
The MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 is intended for use in the electromagnetic environment specified below. The user of MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 should make sure that it is used in such an environment			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 require continued operation during power mains interruptions, it is recommended that the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	10 V/m 80 MHz to 6 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	10 V _{rms} on 150 kHz to 80 MHz 1 kHz 80 % AM modulation 10 V _{rms} in ISM bands	

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Table 3 – Test levels for Radiated and conducted Immunity Tests
Guidance and manufacturer's declaration – electromagnetic immunity

The MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 is intended for use in the electromagnetic environment specified below. The user of MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	10 V _{rms} 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b .
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m 80 MHz to 6 GHz	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 is used exceeds the applicable RF compliance level above, the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 – Emission class and group		
Guidelines and manufacturer's declaration – electromagnetic emissions		
The MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 is intended for use in the electromagnetic environment specified below. The customer or user of the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class B	The MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5 Recommended separation distances between portable and mobile RF communications equipment and the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100			
The MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MultiDrive-X ARTHRO 28200MX/PowerUnit UM100 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies. Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.			

HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency surgical equipment. It has been tested according to IEC 60601-2-2 Annex BB.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ Tra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzucuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Av. Ejercito Nacional No. 453 Piso 2,
Colonia Granada, Alcaldia Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamveien1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se
KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestučio st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk
KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B.V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N.V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.com

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Polska Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Interchange 21 Tower, Level 33
399 Sukhumvit Road,
North Klongtoey, Wattana,
10110 Bangkok, Thailand
Phone: +66 2 660 3669
E-Mail: info-th@karlstorz.com

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenheng District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Indonesia
Sinarmas MSIG Tower Level 37
Jl. Jend. Sudirman No. Kav. 21
Jakarta Selatan
DKI Jakarta 12920
E-Mail: infoindonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info-jp@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
68 Waterloo Road, Macquarie Park NSW 2113
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: karlstorz@karlstorz.com.au
www.karlstorz.com

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 7461 708-0
Telefax: +49 7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com