Usage instruction for medical device

Medic Control Peel

Gel for exfoliation BLOCK AGE PEEL GEL for professional use only

A set of BLOCK AGE PEEL GEL for exfoliation from Medic Control Peel consists of Phase № 1 (glass vial) and Phase № 2 (strip-monodose – polypropylene tube).

Ingredients of Phase № 1

Retinoic acid (tretinoin), azelaic acid, hydroxypropylcellulose, isopropyl alcohol.

Ingredients of Phase № 2

Dimethylsulphoxide, hydroxypropylcellulose.

Description

When prepared for the application BLOCK AGE PEEL GEL from Medic Control Peel is a homogeneous transparent mass of yellow color with a characteristic odor without impurities and contains 5 % retinoic acid (tretinoin) and 15 % azelaic acid.

At ambient temperature lower than +18 °C solidification of Phase № 2 is possible. Before mixing and use of phases the set of a material should be kept at the room temperature (not lower than +20 °C) for at least 1 hour. Do not heat using a heat source!

Use

BLOCK AGE PEEL GEL from Medic Control Peel is intended for aesthetic correction of small facial defects: wrinkles, nevus pigmentosus, atonic tissues, post-traumatic and post-operative scars, post-acne scars and other damages of the skin surface. It has a keratolytic, stimulating, lightening and antioxidant effect.

Indications for use

1. Fading skin.
2. Small and large wrinkles ("pouch", "crow's feet", peribuccal and so forth).
3. Post-acne scars.
4. Post-traumatic and post-operative scars.
5. Hyperpigmentation.
6. Solar and senile lentigo.
7. Ephelides (freckles).
8. Actinic hyperkeratosis.
9. Photoaging.
10. Striae (stretch marks).
11. Flabbiness of skin on front abdominal wall and the inner surfaces of the shoulders and hips.

Contraindications

1. Hypersensitivity to the components.
2. Dermatoses in the acute stage, dermatitis of various etiologies, including infectious inflammatory diseases (herpes, pyoderma and so forth).
3. Damaged skin surface.
4. Chronic somatic diseases in acute and decompensated stages.
5. Excessive skin insulation.
6. Pregnancy and lactation.

Method of application

For professional use only!

BLOCK AGE PEEL GEL (Medic Control Peel) is allowed for use by trained and licensed specialists in accordance with current legislation!

Important!

Phase № 2 is a gel with a characteristic odor. It contains dimethylsulphoxide which effectively penetrates into the dermis and through mucous membranes and increases the permeability of all biological membrane to facilitate their penetration by the primary active ingredient: retinoic acid. Therefore, before applying BLOCK AGE PEEL GEL (Medic Control Peel) neither keratolytic products nor other peeling compositions are required to enhance penetration of retinoic acid. Home pre-peel care is not necessary. When mixed phases of BLOCK AGE PEEL GEL (Medic Control Peel) come in contact with water, the smell intensifies. To avoid this use only dry brushes and vessels for the treatment.

Chemical peel

Step 1. Remove make-up including the eye area with MICELLAR CLEANSER for express make-up removal (Mesaltera by Dr. Mikhaylova) or with ETHEREAL CLEANSER soft cleansing emulsion (Mesaltera by Dr. Mikhaylova).

Step 2. Cleansing. Apply CLEANSER MOUSSE (Medic Control Peel) onto the skin, massage lightly and leave for 1-2 minutes. Rinse with water. Dry the skin with a paper towel.

Step 3. Add the content of Phase № 2 polypropylene tube into the Phase № 1 glass vial, then tightly close the lid and shake for 2-3 minutes until you get a

homogeneous transparent yellow-colored mass without foreign impurities.

Step 4. Using a dry brush apply an even layer of BLOCK AGE PEEL GEL (Medic Control Peel) onto the skin.

Step 5. The time of exposure of peeling composition varies from 6 to 10 hours: the more coarse and insensitive the skin and the more pronounced hyperkeratosis and aging signs are, the longer is the time of exposure.

Step 6. After application of BLOCK AGE PEEL GEL (Medic Control Peel) to face or body skin the patient goes home. Warn the patient that during the time of exposure he or she may feel a burning sensation that will subside within 1-2 hours, as well swelling /edema of the skin and growing erythema.

At the end of the exposure period (as determined by the specialist) the patient has to rinse BLOCK AGE PEEL GEL (Medic Control Peel) with water and apply VEGEFARMA cream (Medic Control Peel) or THIRIPHALAN balm (Medic Control Peel) and MEDISCREEN SPF 85 sunscreen cream (Medic Control Peel).

Post-peel care

Apply the VEGEFARMA cream (Medic Control Peel) or THIRIPHALAN balm (Medic Control Peel) daily two times and more to complete the process of exfoliation and protect the skin from ultraviolet radiation by applying MEDISCREEN SPF 85 sunscreen cream (Medic Control Peel) or UV EXPERT MOISTURIZER SPF 50+ (Mesaltera by Dr. Mikhaylova) moisturizing cream serum.

Side effects

Erythema may appear in the treated area (similar to the effects of active tanning) and persist for up to 10 days. Slight edema of the treated area is possible as well as skin dryness and tightening, tenderness to palpation. Perhaps the appearance of cracks in the raised mimic the activity, as well as in the corners of the mouth. The eye area may experience dryness, feeling of sand in the eyes, tearing, scleral injection. Large- or small-lamellar exfoliation (desquamation) may occur within 12-48 hours after peeling depending on the severity of patient's hyperkeratosis and can persist up to 15 days. For 2-3 day after the procedure the patient may experience itching, depending on individual sensitivity to retinoids. Erythema and desquamation may appear on skin areas and mucous membranes where the gel was not applied.

Precautions

The protocol of the procedure is to follow!

Avoid contact with the eyes! In case of contact with the eyes rinse immediately with plenty of water. Do not expose treated areas to intensive heating (sun exposure, artificial tanning, sauna etc.) and cooling.

How supplied

BLOCK AGE PEEL GEL from Medic Control Peel is packaged in a cardboard box with instruction for use, 5 glass vials each containing 2.7 ml of Phase № 1 and 5 strip-monodoses (polypropylene tubes) each containing 2.3 ml of Phase № 2.

Storage conditions

Store in a place protected from light at the temperature between +5 and +25 °C.

Terms of sale

To be sold to specialized beauty salons and clinics.

Not for free sale.

Shelf life

2 years. Do not used after expiry date indicated on the package.

Complaints and claims on quality of medical device are accepted by an authorized representative of the manufacturer.

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Description

BLOCK AGE PEEL CREAM from Medic Control Peel is a homogeneous mass of yellow color with a characteristic odor without impurities.

Use

BLOCK AGE PEEL CREAM from Medic Control Peel is intended for aesthetic correction of small defects: wrinkles, nevus pigmentosis, post-traumatic and post-operative scars, post-acne scars and other damages of the skin surface. It has a keratolytic, stimulating, lightening and antioxidant effect.

Indications for use

1. Fading skin.
2. Small wrinkles ("pouch", "crow's feet", cheek wrinkles).
3. Photoaging.
4. Hyperpigmentation.
5. Post-inflammatory and post-traumatic depigmentation.

Contraindications

1. Hypersensitivity to the components.
2. Dermatoses in the acute stage, dermatitis of various etiologies, including infectious inflammatory diseases (herpes, pyoderma and so forth).
3. Damaged skin surface.
4. Chronic somatic diseases in acute and decompensated stages.
5. Excessive skin insolation.
6. Pregnancy and lactation.

Method of application

For professional use only!

BLOCK AGE PEEL CREAM (Medic Control Peel) is allowed for use by trained and licensed specialists in accordance with current legislation!

Pre-peel preparation

Patient should use CLEANSER MOUSSE (Medic Control Peel) every other day for 7-14 days and then apply PREPEEL ACTIVE cream (Medic Control Peel). Apply MEDISCREEN SPF 85 sunscreen cream (Medic Control Peel) or UV EXPERT MOISTURIZER SPF 50+ (Mesaltera by Dr. Mikhaylova) moisturizing cream serum daily during daylight hours.

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Step 2. Cleansing. Apply CLEANSER MOUSSE (Medic Control Peel) onto the skin, massage lightly and leave for 1-2 minutes. Emulsify. Rinse with water. Dry the skin.

Step 3. Apply LOCAO RICO gel with alpha-beta complex (Mesaltera by Dr. Mikhaylova) onto the skin for 10 minutes. Rinse with water. Dry the skin.

Step 4. Apply BLOCK AGE PEEL CREAM (Medic Control Peel) evenly onto the skin (on the face including the eyelids) and rub the cream with a cotton ball or gloved hands for 5 minutes. The intensity and duration of application depends on the patient's individual characteristics and the skin type: the more coarse and insensitive the skin and the more pronounced hyperkeratosis and aging signs are, the more intense the application.

Step 5. The time of exposure of peeling composition on face varies from 8 to 12 hours: the more coarse and insensitive the skin and the more pronounced hyperkeratosis and aging signs are, the longer is the time of exposure.

Step 6. After application of BLOCK AGE PEEL CREAM (Medic Control Peel) to face or body skin the patient goes home. Warn the patient that during the time of exposure he or she may feel some tingle, slight swelling/

edema of skin and growing erythema.

At the end of the exposure period (as determined by the specialist) the patient rinse BLOCK AGE PEEL CREAM (Medic Control Peel) with water and GENTLE CLEANSER MOUSSE (Mesaltera by Dr. Mikhaylova) and apply VEGEFARMA cream (Medic Control Peel) or THRIPHALAN balm (Medic Control Peel) and MEDISCREEN SPF 85 sunscreen cream (Medic Control Peel) or UV EXPERT MOISTURIZER SPF 50+ (Mesaltera by Dr. Mikhaylova) moisturizing cream serum.

Post-peel care

Apply the VEGEFARMA cream (Medic Control Peel) or THRIPHALAN balm (Medic Control Peel) daily two times and more to complete the process of exfoliation and protect the skin from ultraviolet radiation by applying MEDISCREEN SPF 85 sunscreen cream (Medic Control Peel) or UV EXPERT MOISTURIZER SPF 50+ (Mesaltera by Dr. Mikhaylova) moisturizing cream serum.

Side effects

Erythema may appear in the treated area (similar to the effects of active tanning) and persist for up to 10 days. Slight edema of the treated area is possible as well as skin dryness and tightening, tenderness to palpation. Possible cracks appearance in areas of active mimic activity and in the corners of the mouth, a feeling of dryness and skin tightening, skin tenderness. Perhaps the appearance of cracks in the raised mimic the activity, as well as in the corners of the mouth. The eye area may experience dryness, feeling of sand in the eyes, tearing, scleral injection. Large- or small-lamellar exfoliation (desquamation) may occur within 12-48 hours after peeling depending on the severity of patient's hyperkeratosis and can persist up to 15 days. For 2-3 day after the procedure the patient may experience itching, depending on individual sensitivity to retinoids. Erythema and desquamation may appear on skin areas and mucous membranes where the cream was not applied.

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How supplied

35.0 ml of BLOCK AGE PEEL CREAM from Medic Control Peel is supplied in a plastic vacuum dispenser, packaged in a cardboard box together with instruction for use.

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Terms of sale

To be sold to specialized beauty salons and clinics. Not for free sale.

Shelf life

2 years. Do not used after expiry date indicated on the package

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"Laboratory THOSCANE" LLC, Russia