POLYSITE®

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heavy metals: lead, mercury, cadmium and hexavalent chrome

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3/ Autres exigences/ Other requirements

- Conditionnement : En carton, reliés par multiples de 10 / Packaging : In boxes, joined by multiples of 10
- Conditions de stockage : Conditions standard / Storage conditions: Standard conditions
- Conditions de transport : Conditions standard / Transport conditions : Standard conditions
- Autre(s) / Other(s): N/A

- FR** Kit de Chambre à Cathéter Implantable Haute Pression
EN Pressure Injectable Implantable Port Kit
ES Kit de puerto implantable de alta presión
IT Kit camera impiantabile per infusione ad alta pressione
SV Kit med implanterbar högfödport
NL Kit met inplanteerbare drukinjectiepoort
AR قطع القسطرة القابلة للزرع بالحقن بالضغط
TR Basınçlı Enjekte Edilebilir Vücuda Yerleştirilebilir Port Kiti
CS Sada implantabilního portu pro injekční tlakem

GOIX Julie

28-07-2021

DocuSigned by:
GOIX Julie
 D14286D0C51843C...

DAILLY Emi

28-07-2021

DocuSigned by:
DAILLY Emi
 1F0C39ECDDDA4E2...

NOTICE D'UTILISATION
 INSTRUCTION FOR USE
 INSTRUCCIONES DE USO
 ISTRUZIONI PER L'USO
 BRUKSANVISNING
 GEBRUIKSAANWIJZING
 تعليمات الاستخدام
 KULLANIM TALIMATLARI
 NÁVOD K POUŽITÍ

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INSTRUCTION FOR USE

I/ DESCRIPTION

SEESITE® implantable ports are composed of a radiopaque port and a catheter, which are connected and held together by a connection ring. They come either preassembled or ready for assembly. The port is accessed using a Huber needle, placed percutaneously.

The materials from which SEESITE® implantable ports are made of, can be used with almost all solutions for injection (see chapter VII/ A/) and are compatible with CT-scan and MRI-scan examinations.

II/ INTENDED USE

SEESITE® implantable ports allow repeated and prolonged access to the venous system. SEESITE® implantable ports are used to administer chemotherapy, antibiotics and antiviral drugs. They can also be used for parental nutrition, collection of blood samples and transfusion of blood or blood products. All SEESITE® references can be used for high flow injection of contrast media during CT procedures (compatible with high pressure injectors, see Chapter IX/ B/). The SEESITE® placement kits consist of suitable instruments for inserting implantable ports by ultrasound-guided puncture (refer to section VIII/ A/). All references can be used on adults. Only the references 2005SEE, 2105SEE, 2015SEE and 2115SEE can be used in pediatrics (children > 10 kg).

III/ CONTRAINDICATIONS

Inserting a SEESITE® implantable port is contraindicated in the following situations:

- Neonates and infants < 10 kg;
- Known or suspected history of immunoallergic or hypersensitivity reactions to materials contained in the device (see chapter IV/).
- Presence of local or systemic infection, bacteraemia or septicaemia.
- When the substances used for treatment are incompatible with one of the materials of the device.
- The patient's body size is insufficient to accommodate the size of the implantable device.
- Implantation site irradiated a priori or a posteriori.
- Previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site.

IV/ CONTENT

Each SEESITE® placement kit contains:

- one sachet containing:
 - a SEESITE® port in titanium and POM (polyoxylomethylene) with a silicone septum, silicone filled suture holes and radiopaque marking;
 - a XRO catheter, in silicone or polyurethane;
 - two connection rings (only one for pre-assembled references) with a Phynox thread containing nickel;
 - a straight Huber needle;
 - a PPS CT safety Huber needle with a connection line (compatible with the injection of contrast agent during CT scan examination);
 - a vein pick;
 - a flushing connector, pre-connected to the catheter (except for reference 4019SEE).

- a peelable sheath introducer;
- a metal J guidewire;
- an echogenic introducer needle;
- a tunneler;
- a 10-mL syringe;
- a Raulerson metal guidewire introducer device.

This device is manufactured from latex-free components and materials. It does not contain DEHP.

- a CIV-Flex™ Echonet consisting of:
 - a latex-free protective cover for the ultrasound probe;
 - a sachet containing sterile ultrasound gel;
 - two elastic bands to hold the protective cover;
 - a sterile drape.

V/ POSSIBLE COMPLICATIONS

The non-exhaustive complications most frequently described after insertion of an implantable port include:

1. Early complications:
 - Pneumothorax;
 - Haemothorax and pleural effusion;
 - Air embolism;
 - Arterial injury;
 - Cardiac arrhythmia;
 - Cardiac tamponade;
 - Local pain;
 - Haemorrhage and haematoma;
 - Incorrect position;
 - Intolerance reactions to the device;
 - Brachial plexus injury;
 - Thoracic duct injury.
2. Late complications:
 - Inflammation, skin necrosis over the implantation area;
 - Local or systemic infection: bacteraemia and/or septicaemia;
 - Formation of a fibrin cuff;
 - Superior vena cava erosion and perforation;
 - Cutaneous erosion;
 - Clavicular pain;
 - Catheter leakage and extravasation (drug extravasation);
 - Occlusion of catheter or port;
 - Migration of the implant;
 - Disconnection, rupture and embolisation of the catheter;
 - Rotation of the port;
 - Occlusion and/or rupture of the catheter due to compression between the clavicle and the first rib (pinch off syndrome);
 - Arteriovenous fistula;
 - Thromboembolism/phlebitis, thrombosis.

VII/ MR CONDITIONAL

Pre-clinical testings have demonstrated that SEESITE® Implantable Port are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field from 1.5 to 3.0 Tesla
- Maximum spatial gradient magnetic field of 3,700 Gauss/cm (37 T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode).

Under the scan conditions defined above, the SEESITE® Implantable Port is expected to produce a maximum temperature rise of 2 °C after 15 minutes of continuous scanning.

In pre-clinical testings, the image artifact caused by the device extends approximately 10 mm from the SEESITE® Implantable Port when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

VIII/ WARNINGS AND PRECAUTIONS

The warnings and precautions described below must be rigorously taken into account to avoid any incidents likely to cause device dysfunction and a risk to patient safety.

A/ GENERAL WARNINGS

- Single use device. **DO NOT REUSE.** The SEESITE® implantable ports must never be re-implanted. As it is contaminated by blood a used implantable port represents a biological risk (risk of infection) potentially leading to patient death. In addition, PEROUSE MEDICAL cannot guarantee the good functioning of the device after it has been removed (system patency, septum water tightness, catheter integrity, etc.).
- Sterilized with ethylene oxide. **DO NOT RESTERILIZE.**
- **Do not use the device if the packaging has been damaged or opened or after the expiry date.**
- Store away from light, cold, heat (11 °C < T < 29 °C) and humidity. Under these conditions, the product may be used up to the expiry date indicated on the packaging.
- Strict aseptic measures must be followed when inserting and using the device.
- Do not use medicinal products or solvents which are incompatible with the catheter material. Do not use alcohol to unblock a polyurethane catheter as this may cause deterioration of the catheter in the event of repeated or prolonged contact.
- Perform hand hygiene before and immediately after all clinical procedures, before donning and after removal of gloves.
- Keep hands behind needle at all times during use and disposal.
- Use universal blood and body-fluid precautions in the care of all patients due to the risk of exposure to Human Immunodeficiency Virus (HIV) or other blood borne pathogens.
- Handle and dispose of sharps in sharps container.

B/ INSERTION

- Only qualified doctors can insert implantable access ports.
- Before insertion, purge the device of any air and fill with physiological saline or heparinised solution.
- **PURGE the Raulerson metal guidewire introducer device with physiological saline solution or heparin solution before puncturing the vein, to avoid any risk of injecting air into the vein (risk of air embolism).**
- Take care during insertion of the device to avoid any mechanical damage of the catheter. **Only clamp the catheter with atraumatic clamps or forceps.** Do not use the catheter if it is damaged or if it leaks.

- For ready-to-assemble ports, the catheter and port must always be assembled manually. Follow the connection technique described under VIII/ INSERTION TECHNIQUES closely to avoid damaging the catheter and to ensure correct connection.
- Please read section XI/ before using the CIV-Flex™ Echonet.

Precautions for subclavian catheter insertion

With a subclavian access, catheter may become compressed between the clavicle and the first rib (pinch-off syndrome) with a risk of rupture and embolisation of the catheter. Excessive and/or repeated movements of the upper limb increase this risk.

- A number of signs are suggestive of pinching of the catheter, including the need to mobilize or raise the patient's arm to perform stable infusion.
- To avoid this phenomenon, catheterize the subclavian vein away from the costoclavicular angle.
- Take an x-ray to make sure that the catheter is not pinched.
- Remove the device if signs of pinching occur.

C/ USE AND MAINTENANCE

- Lack of maintenance can lead to occlusion of the device.
- **Use only Huber needles such as POLYPERF® or POLYPERF Safe®, PPS Quick®, PPS Flow+® or PPS CT.** Conventional hypodermic needles should not be used, as they can alter the integrity of the septum.
- Do not reuse Huber needles. Deformation of the bevel of the needle can tear or damage the septum of the port.
- Only use large gauge Huber needles (19 G and 20 G) for blood transfusions, parenteral nutrition and injection of contrast media.
- Select the length of the needle by taking into account the depth of the port, the thickness of the tissues and the thickness of the dressings underneath the curvature of the needle. Needles that are too short or too long may be displaced during patient's movements with a risk of drug extravasation.
- **NEVER TRY TO UNBLOCK THE DEVICE BY PRESSURE AND NEVER USE A SYRINGE OF A VOLUME OF LESS THAN 10 mL to avoid applying excessive pressure within the device.** Excessive pressure can damage the catheter or lead to disconnection of the port.
- If in doubt as to the integrity of the system, do not use the device and opacify it using a contrast medium.
- **Stop injection at the slightest sign of extravasation and immediately take appropriate measures.**
- Find out beforehand about the patient's sensitivity to heparin. Heparin-Induced Thrombocytopenia (HIT) has been reported with use of heparin flush solutions.

D/ CONTRAST MEDIUM INJECTION DURING SCANNER EXAMINATIONS

- When the implantable port is used to inject contrast media, first check the patency of the port (by obtaining blood reflux) and follow the instructions in chapter IX/ B/.
- Ensure patency of catheter prior to each pressure injection to reduce the risk of catheter failure and/or patient complications.
- Obtain a chest x-ray to confirm catheter tip position prior to each pressure injection of contrast medium.
- **Heat the contrast liquid to 37 °C (100 °F) before use. Failure to do so may reduce the flow rate up to 50% and lead to malfunction of the port or injection system.**

- **Always ensure that catheter length does not exceed 25 cm (this type of injection should not be administered by femoral site). The flow rate will be significantly reduced if a longer catheter is used.**
- Always check the maximum flow rate (with 25 cm catheter) for injection of the contrast medium with the selected SEESITE® reference (see table 1) before use with a pressure injector. If the flow rate is higher than stated in the IFU, the pressure injection of contrast media risks damaging the system (catheter rupture, disconnection of port from catheter and cause complications to the patient (extravasation of contrast media, catheter migration).
- **Do not use Huber needles that are incompatible with and which have not been validated for use with a high pressure system.**
- Do not administer the injection with a high pressure injector if the port is at risk of compression or costoclavicular pinching.
- Assess each patient for appropriateness of a pressure injection procedure.
- Pressure injection procedures must be performed by trained personnel well versed in safe technique and potential complications. High pressure injection using the SEESITE® Implantable Infusion Port should be performed by trained clinicians who are knowledgeable about the use of the SEESITE® Implantable Port.
- Do not reuse non-coring needles. Deformation of the bevel of the needle can tear or damage the septum of the port. Only use large gauge, high pressure compatible non-coring needles.
- Never inject contrast medium with a viscosity greater than the one mentioned in table 1.
- Failure to follow this recommendation risks causing failure of the system through excess pressure or obstruction.
- Stop injection and re-assess if the patient complains of pain, or if there is a swelling when the device is flushed or when drugs or contrast media are administered. Evaluate the condition of the device by opacification, correct needle placement, and potential complication such as occlusion, kinking, breakage, pinch-off syndrome, thrombosis or malposition.
- Pressure limit settings on injector equipment may not prevent over pressurizing of occluded or partially occluded catheter.
- Use appropriate administration set tubing between catheter and pressure injector equipment to reduce the risk of catheter failure.
- Follow contrast media manufacturer's specified instructions for use, contraindications, warnings, and precautions.
- Do not exceed the maximum pressure of 325 psi (22.4 bar) power injector equipment to reduce the risk of port, catheter failure and/or catheter tip displacement.
- Do not exceed 10 injections or catheter maximum recommended flow rate (located on product labeling and on the device with radiopaque graduations) to reduce the risk of catheter failure or tip displacement.

E/ DURATION OF IMPLANTATION AND REMOVAL OF THE IMPLANTABLE PORT

- The implantable port must be removed at the end of the treatment.
- The lifespan of the system depends on the duration of the treatment, tailored to each disease and the functional state of the device, particularly the catheter. The most common complication of ageing is fracture.
- Do not leave the device in place when it is no longer necessary.
- Remove the catheter before removing the port. Do not pull on the port while the catheter is still inserted, as this can lead to disconnection of the device and embolisation of the catheter.
- After being removed, the SEESITE® implantable port may be a potential biohazard. Handle and discard in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

VIII/ INSERTION TECHNIQUES

A/ SEESITE® IMPLANTATION

The doctor performing insertion must determine the most appropriate method of introduction for the patient concerned, after having evaluated the risks and advantages of each method. The risk of gas embolism can be avoided by placing the patient in the Trendelenbourg position (body laid flat on the back with feet higher than the head). Please read the warnings and precautions in chapter VII carefully.

Surgical method

1. Purge the port and the catheter with physiological saline or heparinised solution.
2. Perform venotomy.
3. Introduce the catheter into the vein using a vein pick. Suturing the catheter to the vein is useless or even dangerous (risk of venous stripping during removal).
4. Advance the catheter until the junction of the superior vena cava and right atrium.
5. Check catheter tip position by radiography, fluoroscopy or intracavitary ECG (electrocardiography).
6. Prepare a subcutaneous site so that the port is situated about 1 cm under the surface of the skin.
7. Perform trial insertion of the port to check that the dimensions of the site are suitable and that the port is not situated underneath the incision.
8. Tunnel the catheter as far as the pocket.
9. Cut the proximal end of the catheter to the required length, making sure that there is sufficient slack to allow body movements and connection to the port (for ready-to-assemble ports).
10. Thread the connection ring over the catheter (illustration 1).

Illustration 1

11. Mount the catheter onto the port to cover the conical tip (illustration 2).

Illustration 2

12. Push the ring (which advances the catheter) until it locks (illustration 3). The tip of the catheter must not be compressed over the port.


Illustration 3

13. Check that the catheter is connected properly by pulling on it gently.
14. Insert the port into the pocket and check for the absence of kinking of the catheter.

15. Check the patency of the device.
16. Anchor the port in its pocket by using permanent sutures, and then close the incisions.
17. Flush and heparinise the device.

Percutaneous method

1. Purge the port and the catheter with physiological saline or heparinised solution.
2. Puncture the vein with the introducer needle assembled on a syringe.
3. As soon as blood return is achieved, remove the syringe while leaving the needle in place, taking care to close the needle hub to prevent any risk of air embolism.
4. Insert the metal guidewire into the vein. Remove the needle.
5. Enlarge the cutaneous puncture site to allow insertion of the introducer.
6. Thread the introducer onto the guidewire and advance using a rotating motion.
7. Remove the metal guidewire and vein pick from the introducer. Leave the sheath in place.
8. Insert the catheter into the sheath.
9. Advance the catheter until the junction of the superior vena cava and right atrium.
10. Check catheter tip position by radiography, fluoroscopy or intracavitary ECG (electrocardiography).
11. Remove the sheath.
12. Follow points 6 to 17 of the surgical method described above.

Percutaneous method using the ultrasound-guided puncture technique

1. Purge the port and the catheter with physiological saline or heparinised solution.
2. Apply sterile gel to the tip of the probe, then put on the protective cover. The cover is held on the ultrasound probe with the two elastic bands supplied in the kit.
3. Apply a layer of sterile gel to the tip of the probe, on top of the protective cover.
4. Apply the ultrasound probe to the skin's surface above the chosen access, until it is displayed on the ultrasound screen.
5. Connect the Raulerson device to the introducer needle.
6. **BEFORE PUNCTURING THE VEIN**, to avoid any risk of injecting air into the vein (risk of air embolism): **INSERT THE INTRODUCER NEEDLE CONNECTED TO THE RAULERSON DEVICE INTO A VIAL OF PHYSIOLOGICAL SALINE SOLUTION (TOP DOWNWARDS) AND COMPRESS THE RESERVOIR OF THE RAULERSON DEVICE UNTIL ALL AIR HAS BEEN EXPELLED.**
7. Puncture the vein with the introducer needle assembled on the Raulerson device while maintaining compression.
8. As soon as the introducer needle is inserted into the access, the reservoir of the Raulerson device fills with blood (when venous return is achieved). The correct position of the needle tip in the access can be checked on the ultrasound screen.
9. Do not disconnect the Raulerson device from the introducer needle.
10. Insert the metal guidewire into the vein through the Raulerson metal guidewire introducer device and introducer needle. During insertion of the guidewire into the lumen of the Raulerson device, the user may feel some resistance as the guidewire passes through the non-return valve.
11. Remove the introducer needle and Raulerson device assembly.
12. Follow points 1 to 12 of the percutaneous method described above.

Specificities regarding insertion of pre-assembled ports

1. Before inserting the catheter into the vein:
 - prepare the subcutaneous pocket;
 - tunnel the catheter from the compartment;
 - cut the catheter to the required length at its distal extremity.
2. A peelable introducer must be used for insertion via the percutaneous route.

IX/ USE AND MAINTENANCE




Read carefully warnings and precautions listed in Chapter VII/.

A/ GENERAL INFORMATION

1. Before use, inspect and perform asepsis of the injection site.
2. Locate the port by palpation.
3. Insert the needle perpendicularly to the skin and advance as far as the base of the port.
4. Check the patency of the device and correct positioning of the needle by obtaining blood reflux and by injection of physiological saline without excessive pressure or local extravasation.
5. Inject or infuse the drug solution or take blood sample (eliminating the first few mL of blood taken). If several products must be administered, flush with physiological saline after each injection.
6. Rinse the device with physiological saline solution after each use.
7. After taking blood samples rinse with physiological saline solution and if necessary flush with heparin (follow the healthcare establishment's protocol).
8. Remove the needle by maintaining positive pressure to avoid reflux of blood inside the catheter.

B/ INJECTION OF CONTRAST AGENT DURING CT SCAN PROCEDURES

1. Follow stages 1 to 4 above (see chapter IX/ A/).
2. Ensure that the patient is carrying a SEESITE® implantable port. The patient should have his identification card as SEESITE implant carrier and SEESITE® patient diary with him. The SEESITE® implantable port and maximum recommended injection flow rate (with 25 cm catheter) for contrast medium injection can also be identified by x-ray, as shown in the table below.

Maximum recommended flow rate (with 25 cm catheter) for injection of the contrast medium			
			
1 mL/sec (viscosity ≤ 6cP)	3 mL/sec (viscosity ≤ 10cP)	5 mL/sec (viscosity ≤ 10cP)	
One graduation = 1 mL/sec			

3. Use a Huber needle with tubing validated for pressure injection such as PPS CT and refer to the specific instructions for use.
4. Select the appropriate Huber needle gauge to the viscosity of the contrast medium used and to the implantable port reference following the recommendations in table 1.
5. Set the maximum injector pressure to 325 psi or 22.4 bar (2240 kPa) in order to guarantee the reliability of the system.
6. Only use 19G or 20G needles for SEESITE® Standard (4000 series) and Mini (3000 series) and 20G or 22G for SEESITE® Micro (2000 series), in accordance with table 1.

7. Do not inject more than 2mL/s when using a 22G Huber needle.
8. Always ensure that the injection circuit is fully operational (Huber needle and implantable port) by obtaining reflux of blood and injecting 10-20mL of physiological saline solution without difficulties.
9. Rinse the implantable port with 10 to 20mL of 0.9% NaCl before and after use, followed by standard rinsing procedures.

X/ TRACEABILITY

Self-adhesive labels and a patient card are provided in each pack to allow the required follow-up of implanted SEESITE® access ports. One label must be affixed to the patient's medical record, one affixed to the patient diary (also supplied with the SEESITE®), and another to the patient card. The patient should be given the patient diary and patient card following the implantation procedure. Ask the patient to carry these with him/her at all times.

XI/ ADDITIONAL INFORMATION ON THE CIV-FLEX™

A/ APPLICATION:

The cover allows the transducer to be used for CT scans and needle-guided examinations for body surface area, internal cavities and perioperative diagnostic ultrasound, while helping to prevent the transfer of microorganisms, body fluid and particulate matter to the patient and healthcare personnel when re-using the transducer.

B/ WARNINGS:

- According to federal regulations in the United States, the sale and use of this device is restricted to physicians or prescription by a physician.
- Ultrasound training is required before use. Please refer to the device manual for instructions on use on how to use the transducer.
- Single-use components should be disposed of as infectious waste.
- Refer to the system manual for more information on re-treating the transducer between each use.
- Do not re-use, re-treat or re-sterilise the single-use device. Re-use, re-treatment or re-sterilization may cause contamination of the device, patient infection or cross-infection.

Covering the transducer and wrist strap:

Only use water-soluble agents or gels. Petroleum or mineral oil-based substances may damage the cover.

- Disposable items are for single use only. Do not use after the expiry date.
 - Only use disposable components labelled as sterile if the packaging is intact.
 - Always place the cover on the transducer to protect patients and users from cross-contamination.
1. Place a suitable amount of gel inside the cover and/or on the front of the transducer. Not using gel can lead to poor image quality.
 2. Insert the transducer into the cover, using an appropriate sterile technique. Carefully pull the cover over the front of the transducer to eliminate folds and air bubbles, taking care to avoid any tearing.
 3. Attach the cover using the enclosed bands (where necessary).
 4. Inspect the cover to ensure there is no perforation or tearing.

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SEESITE® is a registered trademark of PEROUSE MEDICAL