
IH-1000

User Manual



BIO-RAD

Identification

IH-1000 User Manual

Version 2.0 - 01/2022

Document: H009142

4 001000 / 001000RECOND

3 /



K_lsd_arspcp

DiaMed GmbH

Pra Rond 23

1785 Cressier FR

Switzerland

Tel. +41 (0)26 67 45 111

Fax +41 (0)26 67 45 145

web: www.bio-rad.com



@gm-P_b

3 Boulevard Raymond Poincaré

92430 Marnes-la-Coquette

France

Tel. +33 1 47 95 60 00

Amnwpge fr è

The reproduction, even partially, of this document is forbidden. No part may be copied in any form, and it may not be used, edited nor transmitted by any electronic means (photocopy, photography, magnetic supports or other recording processes), without the written authorization of Bio-Rad. All rights and particularly reproduction, translation, edition, distribution and also industrial property and recording are reserved.

Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

Bmas k c l r Tcpgm l	Qmdru_pc Tcpgm l	B_rc	Kmbgda_rgm l
/,.		07/2010	First edition
/,/	T1,.	03/2011	V3.0 Soft updates
/,0	T2,.	11/2012	V4.0 Soft and troubleshooting updates
/,1		11/2013	Compliance 61010-1 2010, User level and minor updates
/,2	T.2, .2	03/2015	V04.04 Soft updates and minor updates
/,3	T.2, .5	07/2017	New layout Content update according to the software version Revised: <ul style="list-style-type: none"> the Glossary on page 12; the procedures of Managing Liquids and Wastes on page 93; the description of the Options on page 125; the Maintenance Operations on page 148; the list of Examples of Test Execution; Barcode Exception Mask Filter Option on page 130; Added: <ul style="list-style-type: none"> a glossary of abbreviations; a procedure for Removing Sample Racks on page 88
/,4	T.2, .6	03/2020	Content update according to the software version Revised: <ul style="list-style-type: none"> the management of the resources used by the instrument and the tests to be performed; the flowchart describing a normal test performed by the instrument; the criterias of the return of gel cards for 2nd reading; the list of troubleshooting. Added: <ul style="list-style-type: none"> a statement about the limitations to install third-party software in the computers of the instrument; a warning before proceeding to the decontamination procedure.
/,5		01/2021	Revised: <ul style="list-style-type: none"> the section describing the clot detection; the cleaning and the replacement of the needle. Added a statement about overvoltages.
0,.		01/2022	Added complementary information about cleaning solution preparation. User Manual fully revised in accordance to the new European regulation (EU) 2017/746 (IVDR).

N?EC GLRCLRGML?JJW JCDR @J?L I

Table of Contents

[illegible]

Af_nrcp 1	Glqps kclr Mtcptjcu	01
3.1	Introduction	23
3.2	Instrument Overview	24
3.2.1	Front View	24
3.2.2	Pipetting Area	25
3.2.3	Rear View	26
3.2.4	Computer Connecting Panel	27
3.2.5	Computer Protection	28
3.3	Technical Data	29
3.3.1	Performances	29
3.3.2	Commercial Information	29
3.3.3	Dimensions of Instrument	29
3.3.4	Electrical Data	29
3.3.5	Environmental Conditions	29
3.3.6	Emissions	30
3.3.7	Identification Plate	30
3.4	Theory of Operations	31
3.4.1	Reaction Process	31
3.4.2	Test, Interpretation and Result	31
3.5	Process Flow	32
3.6	Features	32
3.6.1	Full Positive Identification	32
3.6.2	Liquid Level Detection	32
3.6.3	Clot Detection	33
3.6.4	Priming and washing the Fluid Circuit	33
Af_nrcp 2	Qmdru_pc Mtcptjcu	13
4.1	Software Overview	35
4.2	IH-1000 Software Main Screen	36
4.2.1	Header strip	37
4.2.1.1	Last Screen or Main Screen Button	37
4.2.1.2	Menu Button	37
4.2.1.3	Data Management	37
4.2.1.4	Miscellaneous Information Area	37
4.2.2	Interactive Instrument Image	38
4.2.3	Search Area	38
4.2.4	Traceability Area	39
4.2.5	Samples and Reagents Loading Area	40
4.2.6	Footer Strip	40
4.2.6.1	Resources Buttons	40
4.2.6.2	Tests Buttons	41
4.2.6.3	Other Buttons	41
4.3	Warning Screen	41
4.4	Authentication Screen	42
4.5	Maintenance Reminder	42
4.6	Errors	42
Af_nrcp 3	Ecrryle Qr_prcb	23
5.1	Starting the Instrument	46
5.1.1	Preliminary Checks	46
5.1.2	Starting Procedure	47
5.2	Exiting the Software	48
5.3	Stop All of the Instrument	49

5.4	Shutdown of the Instrument	49
5.5	Rescue Shutdown	50
5.6	Step by Step Procedure	51

Af_nrcp 4	K_l_ejle Pcqm spacq	35
6.1	Managing Gel Cards	58
6.1.1	Description	58
6.1.2	Gel Cards Storage Drawer Screen	59
6.1.2.1	Access	59
6.1.2.2	Description	59
6.1.3	Opening / Closing the Drawer	60
6.1.4	Loading Gel Cards	62
6.2	Managing Samples and Reagents	64
6.2.1	Sample Racks	64
6.2.1.1	Using Standard Tubes	65
6.2.1.2	Using Low Volume Tubes (for Red Cells Only)	66
6.2.1.3	Using Plunger Tubes	66
6.2.1.4	Using Low Volume Segment Blood	66
6.2.1.5	Labeling	68
6.2.1.6	Emergency Sample Rack	69
6.2.1.7	Using a New Rack	69
6.2.1.8	Specifications for Tubes	69
6.2.1.9	Specifications for Barcodes	70
6.2.1.10	Sample Tubes Dead Volume	70
6.2.1.11	IH-1000 Centrifuged Sample Tubes Minimum Plasma Height	71
6.2.1.12	Red Blood Cells Minimal Height	71
6.2.2	Diluent Rack and Diluent Product	72
6.2.2.1	Diluent 1 Rack Product Instructions	73
6.2.2.2	Inserting a Diluent Rack Product	73
6.2.3	Reagent Racks	74
6.2.3.1	Use of the Cleaning Solution SETUP CLEAN	75
6.2.3.2	Preparation of 15ml of SETUP CLEAN	76
6.2.4	Samples / Reagents Area Screen	77
6.2.4.1	Description	77
6.2.4.2	Rack Status	78
6.2.4.3	Sample Status	78
6.2.4.4	Diluent Status	79
6.2.4.5	Reagent Status	79
6.2.4.6	Details Button	80
6.2.4.7	Sample Barcode Manual Input	83
6.2.4.8	No Cap Confirmation	84
6.2.5	Loading Racks	86
6.2.6	Removing Racks	88
6.2.6.1	Removing Sample Racks	88
6.2.6.2	Removing Reagent or Diluent Racks	90
6.2.6.3	Remove all Racks	90
6.2.6.4	Removing Diluent and Sample Rack manually	91
6.2.6.5	Removing Reagent Racks manually	92

6.3	Managing Liquids and Wastes	93
6.3.1	Description	93
6.3.2	Checking Liquids and Wastes Status	94
6.3.3	Refilling a System Liquid Container	95
6.3.4	Cleaning Liquid Preparation	96
6.3.5	Cleaning of the Containers	97
6.3.6	Replacing a Liquid Waste Container	98
6.3.6.1	Preparation of the Antifoaming Solution	98
6.3.6.2	Replacing the Liquid Waste	98
6.3.7	Priming Cycle	99
6.3.8	Cleaning Cycle	99
6.3.9	Replacing the Gel Cards Waste Bin	100
6.4	Checking Resources Status	101
6.4.1	Missing Resources	101
6.4.2	Resources on Board	102
6.4.2.1	Gel Cards Storage Drawer Area	103
6.4.2.2	Reagent / Diluent Area	103
Af_nrcp 5	K_l_ejle Rcqrq	/ .3
7.1	Performing a Test	106
7.2	Performing an Urgent Test	107
7.3	TESTS Buttons of the Main Screen	107
7.4	Tests in Progress	108
7.4.1	Details of a Test in Progress	110
7.5	Tests ready to start	111
7.5.1	Emergency	112
7.6	Tests without Resources	113
7.7	Tests to Repeat	114
7.8	Samples without Request	117
7.9	Checking Pipetting Area Status	119
7.9.1	Opening the Pipetting Area Window	121
7.9.2	Clot Detection Errors	122
7.9.2.1	Clot Detection	123
Af_nrcp 6	Mnrgm lq	/03
8.1	Options Screen	126
8.2	General Options	128
8.2.1	Barcode Exception Mask Filter Option	130
8.3	Profile Options	132
8.4	Controls Options	134
8.4.1	Return of the Cards in Drawer for Second Reading Option	135
8.4.1.1	Remaining Cards Locations in the Drawer	135
8.5	Disabled Mode Option	137
8.6	Cap Detection	139
8.6.1	Type of Reagent Vials with Detected Cap	139
8.6.2	Type of Standard Tubes with Detected Cap	140
8.6.2.1	Type of Cylindrical Tubes Caps not Detected	141
8.6.3	Low Volume Tube (for Red Cells only) Cap Detection	141
8.6.4	Plunger Cap Detection	141
Af_nrcp 7	K_qlrc l_l ac	/21
9.1	Checking Components Status	144

9.2	Maintenance Screen	146
9.3	User Level Screen	147
9.4	Maintenance Operations	148
9.4.1	Decontamination	148
9.4.2	Periodic Maintenance	149
9.4.2.1	Daily Maintenance	149
9.4.2.2	Weekly Maintenance	149
9.4.2.3	6 and 12 Months Maintenance	158
9.4.2.4	Total Estimated Times and Liquid Volumes for Weekly Maintenance	158
9.4.3	QC Check	159
9.4.4	Removing the ID-Cards Piercing Array Equipped	159
9.5	Needle Cleaning Right/left Pipettor	161
9.6	Needle Change Right/left Pipettor	164
9.7	Rinse Right/Left Pipettor	166
9.8	Reading Station Calibration	166
9.9	General Initialization	167
9.10	Storage	168
9.10.1	Short Term Storage in Laboratory Conditions	168
9.10.2	Other Storage Conditions	168
Af_nrcp / .	Rpms`jcqfmmrql e	/47
10.1	Alarms and Error Management	169
10.2	Software Error Messages	170
10.2.1	All Modules	170
10.2.2	Transport Arm	170
10.2.3	Identification Module	171
10.2.4	Centrifuge	172
10.2.5	Reading Module	173
10.2.6	Pipettors	173
10.2.7	Hydraulic	174
10.2.8	Reagent Rack	175
10.2.9	Gel Cards Piercing Module	175
10.2.10	All Cameras	176
10.2.11	System	176
Af_nrcp //	?nnclbgv	/57
11.1	Equipment and Accessories	179
11.1.1	Equipment	179
11.1.2	Accessories	180
11.2	CE Compliance	180
11.3	Documents / Forms	180
11.3.1	Decontamination Certificate	180
11.3.2	Maintenance Record	182
11.4	Spare Parts Ordering	183
11.5	Device Disposal	183
11.5.1	General Warnings	183
11.5.2	Disposal	183
11.6	List of Documents	184
11.6.1	User Manuals	184
11.6.2	Service Manuals	184
11.6.3	Spare Parts Manual	184
Glbcv		/63

N?EC GLRCLRGML?JJW JCDR @J?L I

1 Generalities

Af_nrcp Mtcptqcu

This chapter contains basic information on the instrument and the document structure.

1.1 Intended Use

The IH-1000 is a fully automated instrument intended for the ID-System used in immunohematology testing for ABO blood grouping (forward and reverse), antigen typing, antibody screening, antibody identification, Direct Antiglobulin Testing (DAT) and compatibility tests in human blood samples.

For in vitro diagnostic use, by trained laboratory personnel. Device for professional user in a laboratory environment only.

Using the IH-1000 is only permitted in conjunction with the corresponding software or in a configuration which is authorized by Bio-Rad.

The use of any material other than the one specified in the User Manual (e.g. non-authorized substances) is forbidden.

0 The instructions contained in the current User Manual must be adhered to with particular reference to the safety instructions.

1.2 Warranty Limitation

Although the instrument and its software have been tested, it is highly recommended to perform a backup of the computer before any installation and use of the software.

Bio-Rad denies any responsibility in case of:

- wrong use of the instrument;
- unauthorized modification (willingly or unwillingly);
- not complying to the instructions contained in the manuals provided with the instruments and software;
- not complying to the safety instructions contained in the manuals;
- damages linked with the use of the instrument, in particular any data loss or any financial loss which could possibly be attached to the use of the software;
- if the instrument is used in a manner not specified by the manufacturer, the protection provided by the instrument may be impaired.

When the instrument is connected to a host, the user takes the entire responsibility for an error-less transmission of the results (hardware, software, firmware, etc...) to this system.

```

i ?lw u_pp_lrw uqjj `c bcc k cb tmgb gd d_sjr qq dms l b rm f_tc `ccl a_sqcb `w k_jrpc_rkc l r*
k qsqc* sl_srfmpgxcb k_g lrc l_lac md qcptgac mp lcejgeclac md pcesj_p k_g lrc l_lac _lb
qcptgac* _aagbc l r_j b_k _ec* g l amppcar qrm p_ec mp sqc md rfc npmb sarq dmp mncp_rgm l q
msrqgbc rfcgp qncagdgcb jg k gr_rgm l q* msrqgbc rfcgp qncagdg a_rgm l q* am l rp_pw rm rfc g l qrp sargm l q
egtcl g l rfqg k _ls_j,

```

Each Bio-Rad's IH-1000 is tested by the manufacturer before shipping.

1.3 Glossary

The following terms, among others, are used in this manual.

ES	Graphical User Interface
B+A _{pb}	Term which designates the gel card
J?L	Local-Area Network
OA	Quality Control
SB	Unique Device Identifier
SNQ	Uninterruptible Power Supply
SQ@	Universal Serial Bus

/,1,/

Ncpqm l q**K_lsd_arspcp**

The manufacturer of the IH-1000 is:

DiaMed GmbH
Pra Rond 23
1785 Cressier FR
Switzerland

Mncp_rmp

The operator is the owner of a IH-1000 when using it as its owner and when transferring it to a third-party.

Ncpqm l l cj

The term “personnel” includes any person who performs any kind of operation with the IH-1000. These individuals should be qualified in accordance with the manufacturer’s requirement.

Rcaf l qa_j ncpqm l l cj

The term «technical personnel» designates the duly trained persons, who are allowed to perform specific tasks on the IH-1000.

For instance, an electrician is designated as technical personnel for the activities linked to wiring the IH-1000 to the electrical network.

Ncpqm l l cj os_jgdga_rqm l q

The personnel qualifications define the minimum requirements which must be met by the authorized personnel.

Kg l mp G l hspw

A Reversible injury that does not require medical treatment.

Kmbcp_rc G l hspw

A Reversible injury that does require medical treatment.

K_hmp G l hspw

The damage that is irreversible, causes handicap or death.

/,1,0

Npmbsar**Glqrps k c l r**

This is the IH-1000 distributed by the manufacturer.

@gm-P_b _srfmpqxcb Ecj A_pb

Gel card based on column agglutination technique including a microtube system consumable where the antigen-antibody binding is fixed in the gel matrix.

■ M l jw ecj a_pbq k_lsd_arspcb `w @gm-P\_b k\_w `c sqcb gl rfc gl qrps k c l r,

Q_k njc

Content of any sample tube.

@gm-P_b _srfmpgxcb Pc_ec l r

A consumable substance or compound that is added to a system in order to bring about a chemical reaction, or is added to see if a reaction occurs. Such a reaction is used to confirm the presence of another substance.

MIjw pc_ec l r q k_l sd_arspcb `w @gm-P\_b k\_w `c sqcb g l r f c g l q r p s k c l r,

/ ,1,1

Q_k njc K_rcpg_j**?@M-P f Rcqrj l e**

Centrifuged blood samples are required for tests which involves red cells EDTA is considered an acceptable anticoagulant for patient and donor samples.

Clotted samples (no anti coagulant in the sample) cannot be used when red cell testing is being performed.

Clotted, grossly hemolyzed, grossly lipemic or grossly icteric samples may result in inaccurate typing or increased “not interpretable” results.

Anti-coagulated blood samples can be tested for up to 48 hours after collection. If the samples are not tested within 24 hours of collection, samples should be stored at 2-8 °C.

The samples must return to room temperature prior to analysis. The use of cold samples may cause dispensing errors or non-specific reactions.

Prior testing, samples must be centrifuged according to the local Good Laboratory Practices, e.g. 10 minutes at 1500g, to obtain a distinct separation between cells and plasma/serum.

? l r j `mbw Qapcc l Rcqrj l e

Plasma or serum can be used for the antibody screening, antibody identification and crossmatch test on the IH-1000. The serum must be in a separate tube without the presence of red cells in the same tube.

1.4 **Typographical Conventions**

The following styles are used in this manual.

/ ,2,/

Bcqapgnrj m l

This style, used in conjunction with illustration numbers, is preceded with the corresponding numbers:

Cv_k njc

- 1** First element.
- 2** Second element.
- 3** Etc...

/,2,0 Am k k_l b

Any software command, button, function key, window, icon, option, tab, checkbox, selection box, article, menu, tool bar, field and section used in this document is represented by a bold font.

Cv_k njc

The **EXIT** command allows to quit the software.

/,2,1 Npmacbspc

Each procedure step to be carried out step-by-step by the user is preceded by a letter.

Cv_k njc

A	Open the drawer.
B	Put the gel card into position as shown.
C	Close the drawer.

/,2,2 Npmacbspc Pcqsjr

A procedure result is shown by the following symbol .

Cv_k njc

A	Click on the Parameters button. <i>Rfc n_p_ k crcp qarcc l q q b qnj_wcb,</i>
----------	---

/,2,3 Aprmq Pcdpc lac

This style helps the user to find complementary information linked to the current subject.

Cv_k njc

See chapter [Cross Reference on page 15](#).

/,2,4 Jgqr md Grc k q

This style is used in order to display a list of elements.

Cv_k njc8

- item 1;
- item 2;
- item 3.

/,2,5

Rpms`jcqfmmrgle

The complete description along with the error code, the description and the corrective actions as follows:

Cppmp ambc	Cppmp bcqapynrgml	Amppcargtc_argmlq

/,2,6

Gk_ecq

Versions on screen shots do not necessarily correspond with the current version of software.

/,2,7

U_plgle

The following styles are used in this manual.

Content of software screens and buttons remain in English printed in bold italic type and are followed by their translation in brackets (e.g. click the **VALIDATE** button to confirm).

/,2,7,/

U_plgle

■ Rm bcqgeI_rc8
_l g k k g l c l r _l b b _l e c p m s q q g r s _r g m l u f g a f * g d l m r _t m g b c b * k _w j c _b r m k _h m p g l h s p w m d r f c
s q c p 9
m p _n m r c l r g _j j w b _l e c p m s q q g r s _r g m l u f g a f * g d l m r _t m g b c b * k _w j c _b r m 8
2 k m b c p _r c g l h s p w m d s q c p m p g l r f c r c q r q ` c g l e d _j q c 9
2 g l h s p w m d s q c p m p g l g k n m p r _l r b c j _w g l b c j g t c p w m d r f c r c q r p c q s j r q ,

/,2,7,0

Lmrc

- 0 Indicates:
- a preferred procedure or recommended use,
 - a general or informative remark.

2 Safety and Handling

Af_nrcp am lrc l r q

This chapter sets out safety instructions to ensure safe and trouble-free operation of the IH-1000 and its associated softwares. It also describes the handling and storage conditions.

2.1 Introduction

0,/,/ Npg l agnjc q

- @cdmpc a_ppwg l e msr _lw mncp_rgm l u f_rqmctcp m l rfc GF-/. . . *gr gq g k ncp_rgtc rm pc_b rfgq af_nrcp _l b dsjjw s l bcpqr _l b gr,

In case of any doubt, consult your Bio-Rad representative.

0,/,0 G k nmpr _l ac md rfc Q_dcrw G l q rpsargm l q

- ?jj rfc q_dcrw g l q rpsargm l q g l rfgq Sqcp K_ls_j _l b m l rfc g l q rps k c l r k s q r `c am k njgcb ugrf g l mpbc p rm npct c l r _aagbc l r q rm ncpqm l q* b_k _ec rm cosgn k c l r mp nmjjsrgm l md rfc cltgpm l k c l r,

In a similar manner, the legal bylaws and the recognized technical rules which apply in the country of use of the IH-1000 must be adhered to.

If during the use of this device or as a result of its use, a serious incident occurs, please report it to Bio-Rad Laboratories and/or its authorized representative and to your national Competent Authority.

0,/,1

Bgqpc_e_pbg_l_e rfc Q_dcrw Psjcq

Disregarding the safety rules, as well as existing legal and technical regulations, may lead to accidents, property damages or environmental pollution.

- Bgqpc_e_pbg_l_e rfc g_lqrpsargmlq dmp sqc egtcl `w rfc k\_lsd\_arspcp k\_w pcbsac rfc jctcj md nprmcargml mddcpb `w rfc g_lqrps kclr,

2.2**Environmental Conditions**

- Rfc GF-/... k_w mljw `c sqcb g\_lbmmpq\_l\_b lctcp g\_l rfc g\_k kcbg\_rc tgaglgrw md n\_rgc\_lrq, Rfc cjcargp\_a\_j q\_dcrw md rfc GF-/... gq es\_p\_lrcb mljw gd rfc cjcargp\_a\_j g\_lqr\_jj\_rgm\_l md rfc `sgjbg_l_e am knjgcq ugrf l_rgm_l_j pcesj_rgm_lq _nnjga_`jc rm `sgjbg_l_e q dmp kcbg_a_j nspnmcq mp j_`mp\_rmpgcq\_l\_b gd rfgq g\_lqr\_jj\_rgm\_l gq dslargml\_g\_l\_e ncpdcarjw, Rfc GF-/... k\_w lmr `c sqcb g_l npc kgqcq qs`hcar rm pgq_i_q md dgpc mp cvnjmqgm_l,
- ?q dmp _lw mrfcp cjcargp_a_j bctgac* rfc GF-/... k_w lmr `c jma_rcb lc_p _u_rcp_r_n, Bm lmr _jjmu _lw jgosgb rm nclcrp_rc g_lqgbc rfc fmsqg_l_e md _lw g_lqrps kclr, G_l_a_qc md k_hmp qnqjj_ecq md jgosgb g_lqgbc _lw g_lqrps kclr* npmaccb g_k kcbg_rcjw _q dmjjmuq8
?, Qugraf mdd rfc GF-/... _r rfc nmucp qugraf,
@, Bgqam_l_lcar rfc njse,
A, Bpw rfc g_lqrps kclr,
B, Ajc_l_l_b bcamlr_kgl_rc rfc g_lqrps kclr,
C, Afcai _jj cjcargp_a_j dslargmlq,
Qcc af_nrcp [Bcamlr_kgl_rgm_l ml n_ec /26](#),

0 For more information, contact your representative.

- Rfc GF-/... ksqr `c icnr \_u\_w dpm k nmrc\_lrq\_j qmspacq md g\_lrcpdcplac, Rfc GF-/... k\_w lmr `c cvnmqcb rm bgpcar qs_ljgef* fc_r* bsqr mp cvacqgqtc fs kgbgrw &sqc mljw g_l _ajc_l_j_`mp_rmpw cltgpm_l kclr',
- Rfc g_lqrps kclr ksqr lmr `c jma_rcb lc_p qmspacq md fgef cjcargp_k_elcrp_a_p_bg_rgm_l ufgaf k_w g_lrcpdcpc ugrf npmnpc mncp_rgm_l, Rfc sqc md rfc g_lqrps kclr g_l _bpw cltgpm_l kclr cqncq_jjw gd qwlrfcrp_a k_rcp_a_jq _pc npcqclr &ajmrfcq* qwlrfcrp_a _pncrq* cra,' k_w a_sqc cjcargpqr_rga bgqaf_pecq_l_b jc_b rm d_jqc pcqsjrj,

0 Assess the electromagnetic environment of the site in which the IH-1000 is located before switching it on.

- Rfqq glqrpskclrf_q`cc lrcqrcb\_lbdmslbrm amknjw ugrf rfc jgkgrq dmp \_Aj\_qq @ bgegr\_j bctgac\* nspqs\_lrrm n\_pr /3 md rfc DAA Psjq, Rfcqc jgkgrq \_pc bcqgelcb rm npmtgbc pc\_qml`jc npmrcargml _e_g lqr f_p kdsj glrcpdcplac gl _pcqgbc lrg_j glqr_jj_rgm l, Rfqq glqrpskclr eclcp_rcq* sqcq_lba_l p_bg_rc p_bgm dpcosclaw clcpew_lb* gd lmr glqr_jjcb _lb sqcb gl _aampb_lac ugrf rfc glqrpsargml* k_w a_sqc f_p kdsj glrcpdcplac rm p_bgm amk kslga_rgm lq, Rfcpc qq lmes_p_lrcr rf_r glrcpdcplac ujjj lmr maasp gl _n_prgasj_p glqr_jj_rgm l, Gd rfqq glqrpskclr bmcq a_sqc f_p kdsj glrcpdcplac rm p_bgm mp rcjctqgml pcacnrgml ufgaf a_l`c bcrckglcb`w rslgle rfc glqrpskclr mdd_lbm l* rfc sqcp qq clamspecb rm rpw rm amppcar glrcpdcplac`w mlc mp kmpc md rfc dmjjmugle kc_qspcq8

?, Pcmpgclr mp pcjma_rc rfc pcacgtgle_lrc l l_,

@, Glapc_qc rfc qcn_p_rgm l`cruccl rfc glqrpskclr_lb pcacgtcp,

A, Aml lcar rfc glqrpskclr glrm_l msrjcr ml agpasgr bgddcplr dpm k rf_r rm ufgaf rfc pcacgtcp qq am l lcarcb,

B, Amlqsjr rfc bc_jcp mp_l cvncpgclac p_bgm-RT rcaflga_l dmp fcjn,

0 This class (B) digital apparatus complies with Canadian ICES-003.

0 The instrument is designed to be safe when temporary overvoltages that may occur between the line conductor and earth in electrical installations.

- Af_lecq mp kmbgdga_rgm lq lmr cvnpcqjw _nnpmtcb`w rfc n\_prw pcqnm lq`jc dmp amknj_lac amsjb tmg b rfc sqcp_lq_srfmpgrw rm mnpc_rc rfc glqrpskclr,

2.3 General Safety Instructions and Signs

EXPLOSION

Rfc GF-/. . . qq lmr cvnjmqgm l pcqqr_lr, Qn_piq k_w a_sqc _dypc mp_l cvnjmqgm l, Bm lmr glqr_jj rfc GF-/. . . gl rfc tgag lgrw md gl dj_k k`jc mp cvnjmqgtc npmb sarq, Bm lmr sqc dj\_k k`jc* cvnjmqgtc mp afc kga_jjw pc_argtc k_rcpq_jq bspgle rfc aclrpgdse_rgm l ugrf qsdgagclr dmpac rm a_sqc _b_lhcp,

INFECTION

C_af ucjj ml c_af ecj a_pb a_l amlr_gl sn rm /03 Öj md`jmmb q\_knjc qmjsrgml, Rfcpc qq \_pgqi md gl dcargml dpm k qigl amlr\_ar ugrf`jmmb, ?ju_wq uc_p npmrcargtc ejmtcq ufcl umpigle* gl _aampb_lac ugrf j`mp_rmpw q_dcrw pcesj_rgm lq,

CONTACT WITH WATER

Rfc GF-/. . . qq lmr npmrcarcb _e_g lqr u_rcp nclcrp_rgm lq, U_rcp nclcrp_rgle rfc glqrpskclr k_w a_sqc _l cjarpgaj bgqaf_pec, Bm lmr glqr_jj rfc GF-/. . . gl rfc tgag lgrw md u_rcp mp gl _lw nj_ac ufcp _pgqi md u_rcp nclcrp_rgm l qq nmqqg`jc,

- P**cn_gp_lb qcptgac mncp_rgm_lq ksqr mljw `c ncpdmp kcb `w_os_jgdgcb qcptgac cleglccp
 _nnmglrcb `w rfc k_lsd_arspcp,
 R_ic glrm_aamslr_jj rfc u_plgleq_lb dmijmu_jj glqrpsargmlq npmtgcb m_l rfc GF-/..._lb
 gl rfc bmaskclr_rgm_l,
 Rfc GF-/... k_w mljw `c amllcarcb rm rfc cjcargpaj nmucp qmspacq qncagdgcb,
 Gr_gq g kncp_rgtc rm sqc rfc jgosgbq qncagdgcb gl rfc k_ls_j rm ajc_l rfc GF-/..., Rm sqc _
 bgddpcplr jgosgb* amllqsjr rfc k_lsd_arspcp_lb _qi dmp fgq _nnpmt_j,
 Rfc sqc md k_rcpg_jq mrfcp rf_l rfmqc bcdglcb gl rfc sqcp%q k_ls_j &sl_srfmpgxcb
 f_x_pbmsq qs`qr\_lacq qsaf\_q gl dj\_k k\_`jc* cvnjmqgic mp afckga_jjw slqr_`jc
 qs`qr\_lacq' gq dmp`gbbcl, ?lw gl dpgle ckclr md rfgq psjc ujjj `c amllqgbcpb `w rfc
 k_lsd_arspcp_q esjrw lcejeciac,
 Lctcp _rrcknr rm sqc qn_pc n_prq mrfcp rf_l rfmqc qsnnjgcb `w rfc k_lsd_arspcp,
- R**fc glqrps kclr k_w mljw `c nj\_acb ml \_dj\_r qspd\_ac a\_n\_`jc md `c_pgile grq ucgefr,
 Bm lmr nj_ac_lwrfgle _`mtc rfc GF-/..._q rfgq k_w f_kncp rfc sqc md rfc glqrps kclr,
 Rfgq _jqm _tmgbq_lw m`hcar djgle ml rfc glqrps kclr,
 Rfc GF-/... ksqr `c jma\_rcb lc\_p`w _k sp_j cjcargpaj qmaicr rm _tmgb sqgle_l cvrc_lqgm_l
 nmucp ampb,
 Rfc GF-/... ksqr `c jma_rcb gl qsaf _k_l lcp rf_r mncp_rgle rfcyp bgqamllcargle bctgacq
 &rfc ML-MDD qugraf_lb rfc qcn_p_`jc njse md rfc nmucp qsnnjw ampb' gq nmqcg`jc_r_jj
 rgkcg, Gl rfc j_rrcp a_qc* rfc qmaicr msrjcr ksqr `c nj_acb lc_p rfc glqrps kclr_lb ksqr
 `c c\_qg\_jw\_aacqcg`jc `w rfc mncp_rmp,
 Clqspc rf_r rfc tclrgj_rgm_l _pmslb rfc GF-/... gq qsddgagclr rm _tmgb_lw cvacqgic fc_rgle
 ufgaf k_w a_sqc k_jdslargmlq_lb-mp_d_gj_spc md rfc glqrps kclr, ?qn_ac md 1..k k_jj
 _pmslb rfc glqrps kclr ksqr `c jcdr ajc\_p md\_lw m`qr_ajc,
 Gd rfc glqrps kclr lccbq rm `c kmtcb* njc_qc amlr_ar wmsp qcptgac cleglccp,
- ?ju_wq** qrmn rfc GF-/... ugrf rfc QFSRBMUL ds_largml, Lctcp sqc rfc eclcp_j nmucp
 qugraf* cvacnr gl a_qc md_l ckceplaw,
- R**fc GF-/... k_w mljw `c mncp_rcb ugrf rfc qmdru_pc qsnnjgcb_lb ugrf rfc _aacqqmpgq
 _nnpmtcb `w rfc k\_lsd\_arspcp &cv\_knjc8 cvrcpl\_j`_pambc pc_bcp',
 Lctcp amllcar_SQ@ bctgac &SQ@ icw* a_kcp_ . ' rm GF-/... ugrfmsr npctgmsqjw qa_l lgle
 gr ugrf _l snb_rcb_lrgtgsq,
 Rfc GF-/... gq bcjgtcpb ugrfmsr _l_lrgtgsq, Rfc sqc md _l_lrgtgsq gq pcam k kclbcb_lb
 qfmsjb pcesj_pjw `c snb\_rcb, Rfgq glqr\_jj\_rgm\_l qfmsjb `c bmlc ml `mrf NAq, ?qi wmsp jma_j
 bgqrpq`srmp,
 K_ic qspc rf_r rfc JcQ gq cosgnncb ugrf rfc lcacq_pw npmrcargml rm clqspc lm_tgpsq_lmp
 _aacqg gq rp_lqkgrcb rm rfc GF-/...,
 K_ic qspc rf_r lm_srm k_rga qmdru_pc snb_rc gq _argt_rcb ml rfc glqrps kclr amknsrpb,

Gr gq g kncp_rgtc rm m`r\_g l rfc k\_lsd\_arspcp\_nnpmt\_j`cdmpc sqg le rfc GF-/... g l
am k`g l_rgm l ugrf mrfcp bctgacq,

0,1,/

Qge lq

The following signs can be found on the IH-1000. Associated dangers are explained in the table.

0 Appropriate safety instructions for specific situations are described along this manual. For each specific procedure, read carefully the complete chapter.
In case of malfunction or any other technical incident for which a remedy has not been described in the User Manual, please contact the manufacturer immediately.

Qge l	Bcqapgnrgm l	Cvnj_l_rgm l
	A_srgm l* am lqsjr _aam k n_lwgle bmas k c l r q	-
	@gmjmega_j f_x_pb	GLDCARML All blood samples should be considered potentially infectious. Contact with skin/mucous membranes must be avoided. Always wear protective gloves when working with gel cards, sample tubes, waste bin and with any instrument, in accordance with laboratory safety regulations. All biological waste must be handled and treated by the user before disposal in accordance with local procedures and directives.
	Fgef tmjr_ec b_lecp	CJCARPMA SRML During maintenance work, use extreme caution at all times when the IH-1000 is powered and uncovered; the instrument must not be left unmonitored. Avoid using non-isolated metallic parts to work inside the instrument such as screwdrivers.
	J_qcp `c_k b_lecp	AJ?QQ 0 J?QCP + NCP K?LCLR CWCQGLHSPW <ul style="list-style-type: none"> • Never look directly into the barcode reader laser beam. • Never modify, open or remove any protections of the laser beam. • Never interrupt the laser beam with a reflecting object. Most metals reflect the laser beam.
	Nglafgle b_lecp	NGLAF NMGLR F?X?PB Keep hands clear during operation.
	@sp lgle b_lecp	FMR QSPD?AC Contact may cause burns. Do not touch.
	Cjcarpmqr_rga qfmai b_lecp	-
	Gppgr_lr npmb sar	GLPPGR?RGLE Avoid contact with skin and eyes.

Qye l	Bcqapgnrgm l	Cvnj_l_rgm l
	Amppmqgtc npmb sar	AMPPMQgTC Causes severe burns. Avoid contact with skin and eyes.
	Bm lmr qnp_w	Do not spray the surface with compressed air or water as this may cause damage.
	K_g l qugraf	Power ON
	K_g l qugraf	Power OFF
	Npmrcargtc c_prfg l e	-

0,1,0

N_ai_eg l e Qye l q

Qw k `mj	Bcqapgnrgm l
	Fragile, handle with care.
	Keep dry.
	Keep upright.
	Maximum and minimum humidity rate.
	Maximum and minimum temperature.
	Do not stock another IH-1000 on top of the first one.

3 Instrument Overview

Автоматизированный анализатор

This chapter is a general presentation of the instrument.

3.1 Introduction

The IH-1000 automated analyzer allows to launch different kinds of tests within seconds. A large number of tests can be processed simultaneously. Multiple tests can be performed without any special interaction from the operator.

Several high-definition cameras allow to get clear pictures and precise tracing of information during the whole process.


The IH-1000 software allows to monitor on-going tests and resources status.

An external Data Management software (IH-Com) is required to create test request and to analyze results.

0 Read carefully the User Manual of the IH-Com Data Management Software.

3.2 Instrument Overview

1,0,/ Dpm I r Tgc u



Dyespc I, IH-1000 - Front view opened

- 1** Right front door: gel cards waste
- 2** Left front door: Liquids containers
- 3** Gel cards storage drawer
- 4** Samples and reagents loading area
- 5** Pipetting area, protected by a window
- 6** Computer power button
- 7** Instrument Power Switch
- 8** Touch-screen
- 9** Keyboard
- 10** Barcode reader
- 11** Samples preparation table
- 12** UPS and batteries

1,0,0

Ngnncrrgle ?pc_




Diagram 0, IH-1000 - Front view - Pipetting area

- 1** Flux compensator
- 2** Left pipettor
- 3** Transport arm
- 4** Right pipettor
- 5** Washing module
- 6** Localization module

1,0,1 Pc_p Tgcu

MIjw _ qcptgac cleglccp gq _jjmucb rm mncl rfc ` _ai bmmppq,




Рис. 1, IH-1000 - Rear view opened

- 1 Gel cards piercing module
- 2 Centrifuges (12 gel cards)
- 3 Automate computer
- 4 Master computer
- 5 Connecting panel
- 6 Power supply
- 7 Reading module
- 8 Incubator at 37°C

1,0,2 Am k n s r c p A m l l c a r g l e N _ l c j

MI jw _ q c p t g a c c l e g l c c p g q _ j j m u c b r m m n c l r f c ` _ a i b m m p q ,

On all other connection (USB, RS232, screen, LAN), connect only safety Extra Low Voltage devices.


? S R M K ? R C A m l l c a r g m l q

- 1 WIRELESS** (USB connection): To connect an external Wi-Fi device on the computer (Automate)
- 2 KEYBOARD** (USB connection): To connect an additional keyboard to manage the computer (Automate)
- 3 MOUSE** (USB connection): To connect an additional mouse to manage the computer (Automate)
- 4 UPS** (USB connection): To connect the UPS
- 5 OPTION** (USB connection): To connect optional USB devices on the computer (Automate)
- 6 SCREEN**: To connect another screen on the computer (Automate)

Do not touch or change the connections on the automated PC.

K ? Q R C P A m l l c a r g m l q

- 7 WIRELESS** (USB connection): To connect an external Wi-Fi device on the computer (Master)
- 8 PRINTER** (USB connection): To connect a printer to the computer (Master)
- 9 LAN** (Ethernet cable): To connect the LAN network cable
- 10 RS232** (Serial port): To connect the RS232 communication cable



Dyescp 2, Computer connecting panel

The IP address used by the IH-1000 to transmit between the R2 Master PC and the D2 Automate are defined as follows:

- R2 Master PC IP address: 192.168.0.1
- D2 Automate IP address: 192.168.0.2

1,0,3 Am k nsrcp Npmrcargm l

An application has been implemented to protect the computer.

It prevents any software installation from an external drive or local drive when logged in as **Standard user**. The installation of software such as Microsoft Office will be blocked.

An **Administrator** user level allows the installation of any software but the software must be copied on the desktop to be able to launch it.

A new component installation may require a patch to update the authorized application deployment.

At the end of the installation, the Field Service Engineer must inform the customer about the cybersecurity risks. He must define with the customer a new access login and password to forbid the free access to the computer.

3.3 Technical Data

1,1,/ Ncpdmp k_lacq

Rfpmsefnsr	Depending on the profile up to 120 gel cards per hour
Jm_b a_n_agrw	<ul style="list-style-type: none"> Up to 180 sample tubes; up to 28 reagent vials; up to 240 gel cards; 2x5L container for system liquid; 2x5L container for liquid waste; 1x5L container for cleaning liquid; 320 gel cards capacity waste bin
Gbc lrgdga_rgm l	Full positive identification (by barcode) of primary sample tubes, reagent vials, diluent products and gel cards including lot number and expiry date control

1,1,0 Am k k cpag_j G l dmp k_rgm l

Mpbcp lsk`cp	001000
--------------	--------

1,1,1 Bg k c l qgm l q md G l q rps k c l r

Ugbrf /	1229mm (without screen holder)
Ugbrf 0	1730mm (with screen holder)
Fcgefr	1691mm
Bcnrf	823mm
Ucgefr	540Kg

1,1,2 Cjcarpga_j B_r_

Nmu cp qsnjw	Connection on a standard system (presence of normal overvoltage)
Tmjr_ec	100/240VAC
Dpcosclaw	50/60Hz
Nmu cp	1350VA
Fc_r c k g qgm l	720W

O The IH-1000 earthing protection is via the power cord therefore the mains supply used must be connected to a protective earth.

▪ Rfc g l q rps k c l r k s q r `c n m u c p c b `w _ l c _ p r f c b k _ g l q n m u c p q s n j w,

▪ Rfc n m u c p a m p b k s q r `c p c n j \_ a c b `w _ n m u c p a m p b u f g a f g q g l _ a a m p b _ l a c u g r f j m a _ j p c e s j _ r g m l q, T m j r _ e c - A s p p c l r q n c a g d g a _ r g m l q _ p c 03 . T - / . ? m p / / . T - / 3 ? ,

1,1,3 Cl t g p m l k c l r _ j A m l b g r g m l q

Mncp_rg l e r c k n c p _ r s p c	+15°C to 27°C
Qrmp_ec \$ r p _ l q n m p r r c k n c p _ r s p c	-20°C to +55°C
Pcj_rgtc f s k g b g r w * l m a m l b c l q g l e	5 to 80% rh (storage, transport and operating)
K_v g k s k p c j _ r g t c f s k g b g r w	80% for temperature up to 31°C, decreasing linearly to 50% at 40°C.
Nmu cp d j s a r s _ r g m l	+/-10%
Rp_l q c l r m t c p t m j r _ e c	Category II
? j r g r s b c k _ v	2000m

1,1,4 C k g q q g m l q

L m g q c	75dB(A) max
N m j j s r g m l j c t c j	2

1,1,5 G b c l r g d g a _ r g m l N j _ r c

Please copy the following information from the identification plate in the fields below.

Type (instrument)

Serial number (SN)

BIO-RAD

IH-1000

REF

001000

100-240 VAC

SN

XXXXXXXXX

50/60 Hz

Max. 1350 VA

DiaMed GmbH
Pra Rond 23, 1785 Cressier FR
Switzerland



FXAA



(01)xxxxxxxxxxxxxxxx
(11)xxxxxx
(21)xxxxxxxxxx

This class (B) digital apparatus complies with Canadian ICES-003

X
XXXXXX

Dgespc 3, Identification plate

Qye l	Bc q a p p n r g m l
4	Manufacturer reference (order number)
5	Serial number
:	Manufacturer
;	Date of manufacture (yyyy-mm)
/	CE Mark See chapter CE Compliance on page 180 .
3	In Vitro Diagnostic medical device
7	See documentation
	WEEE mark See chapter CE Compliance on page 180 .
	Curtis-Straus marking - Canadian & US
	2D barcode (SN + UDI-DI) UDI = Unique Device Identification (GTIN-14)

3.4 Theory of Operations

1,2,/ Pc_argm I Npmacqg

The ID-System is based on a process concerned with the determination of antigen-antibody complexes. The agglomerated particles are separated from non-agglomerated particles by the means of an inner gel. The gel cards used contain 6 microtubes.

The reaction occurs on the gel-filled microtubes. The gel contains specific antibodies or reactives in accordance with the desired reaction.

After adding the sample to test, it is possible to visually observe in the gel column the separation between the agglomerated and non-agglomerated hematies after a centrifugation. Depending on the intensity of the reaction, hematies penetrate the gel to a different extend. It is thus possible to classify the reaction into 5 steps: +++++, +++, ++, +, -.

1,2,0 Rcqr* G lrcpnpcr rcm l l b Pcqsjr

The camera interprets the result of the reaction in the corresponding gel card. The camera image is divided into search windows. The search windows are then analyzed.

Search windows are divided into 5 zones: the lower part for a negative result, the upper part for the positive ++++ results and three intermediate zones for the +, ++, and +++ results.

A computer algorithm then analyses the image and determines the reaction result for each microtube (refer to the IH-Com User Manual for detailed information on reaction results):

Pc_argm pcqsjr	Bcqapnrgm
) *)) *))) *)))	Positive (in 4 steps)
+	Negative
=	No result interpretation possible Check result manually during validation
Jgo	Liquid distribution problem
bn	Double population
uD	Doubtful result / Gel not clear. Please check manually for foreign particles
uP	Nonhomogeneous area above the pellet (red cells or foreign bodies) Please check manually
uN	Check surface of the pellet
uB	No clear cell pellet. Please check manually
Y [	Well not used for the test
J	No liquid distribution detected in the well
U	Well not found for interpretation
C	No reaction available in the well
C	No liquid above the gel

A test result is based on the results of the corresponding microtubes. That result must then be validated by the user on the basis of the digital image provided. Manual corrections are marked with an asterisk “*”.

3.5 Process Flow

Depending on the process, the pipettor will take suitable quantities of sample, diluent and reagent. It will then supply them to the appropriate gel cards. The sampling needle is then automatically sent to the washing module.

- 1** As soon as the drawer is closed, the transport arm detects the kind of gel card loaded in each tray.
- 2** The samples, diluents and reagents are identified by the identification module as soon as the racks are loaded. The transport arm carries them into the pipetting area after identification.
- 3** The transport arm brings the gel cards to the piercing module, then to the pipetting area.
- 4** The pipettor aspirates the sample, the diluent and the reagent, and then supplies it in the appropriate gel card.
- 5** If necessary, the transport arm brings the gel card to the 37°C incubator.
- 6** The transport arm puts the gel card into an available centrifuge. If necessary, counterweighting gel cards are loaded; the number of gel cards in the centrifuge must be even.
- 7** The transport arm takes the centrifuged gel card into the reading module.
- 8** The camera reads the result.
- 9** The gel card is put in the gel cards waste bin.

3.6 Features

1,4,/ Dsjj Nmqqrgt c Gbc l r g d g a _ r g m l

The IH-1000 automatically identifies the resources required by reading the barcodes of the elements inside the instrument and compares them with the ones needed by the selected test program.

It also checks the lot number and the expiry date of the devices used in the system.

1,4,0 J g o s g b J c t c j B c r c a r g m l

The IH-1000 is designed to permanently check the availability of reagent, diluent and solutions in order to avoid any failure.

- 0** Level detection is disabled for any tube with Low volume segment blood (racks of type 4 - Barcode 04 or 14).
See chapter [Using Low Volume Segment Blood on page 66](#).

It also checks and prompts the user when the waste container is full.

1,4,1

Ajmr Bcrcargm I

To prevent any obstruction of the fluid circuit by red cells or serum during the pipetting process, the IH-1000 is able to detect the presence of a blood or serum clot by means of a sensor.

▪ Rfcpc qq Im ajmr bcrcargm I _t_gj_`jc uqrf q k_jj bjsrgm I q : / # ,

0

This detection is also performed in case of using whole blood and not only concentrated red cells.

0

We strongly recommend to only use anticoagulated samples and to not use samples that are older than 48 hours.

See chapter [Clot Detection on page 123](#).

1,4,2

Npg kgle _lb u_qfgle rfc Djsqb Agpasgr

The IH-1000 allows several operations on liquids. The wash function enables to clean internally and externally the dispensing needles and to prime the fluid circuit.

N?EC GLRCLRGML?JJW JC DR @J?L I

4 Software Overview

Аf_nrcp Mtcptgcu

This chapter introduces the IH-1000 control software.

4.1 Software Overview

The IH-1000 software allows to monitor on-going tests and resources status in the instrument.

External Data Management software is required to analyze results.

0 Read carefully the User Manual of the IH-Com Data Management software.

4.2 IH-1000 Software Main Screen

The main screen allows to access most of the instrument functions.




Figure 4. Main screen

- 1 [Header strip on page 37](#)
- 2 [Current operation area](#)
- 3 [Interactive Instrument Image on page 38](#)
- 4 [Footer Strip on page 40](#)
- 5 [Search Area on page 38](#)
- 6 [Traceability Area on page 39](#)
- 7 [Samples and Reagents Loading Area on page 40](#)

The header and footer strips (1) and (4) are always visible and give access to the options available, depending on the current operation.

The current operation area (2) will display the detailed information.

2,0,/ Fc_bcp qrpqn

The header strip has three buttons and displays miscellaneous information. It is identical on all software screens.



Dyespc 5, Header strip

2,0,/,/ J_qr Qapcc l mp K_l Qapcc l @srrm l

On any screen, **MAIN SCREEN** button allows to go back to the main screen.



On the main screen, **LAST SCREEN** button is displayed instead. It allows to go back to the last screen consulted.

2,0,/,0 Kc l s @srrm l

The **MENU** button gives access to the **MENU** screen with the following functions:

- **OPTIONS** (See chapter [Options on page 125.](#))
- **MAINTENANCE** (See chapter [Maintenance Screen on page 146.](#))
- **AUTHORIZED MAINTENANCE** (See chapter [User Level Screen on page 147.](#))
- **EXIT** (See chapter [Exiting the Software on page 48.](#))
- **STOP ALL** (See chapter [Stop All of the Instrument on page 49.](#))
- **SHUTDOWN** (See chapter [Shutdown of the Instrument on page 49.](#))

2,0,/,1 B_r_ K_l _eck c l r

Button to switch to Data Management software (IH-Com).
See IH-Com Data Management User Manual.

2,0,/,2 Kqqacjj_l cmsq G l dmp k _rjm l ?pc _

Status screen showing system date and time, system status, software version and logged in user.

2,0,0 G l r c p _ a r g t c G l q r p s k c l r G k _ e c

The instrument image is divided into 5 areas. Clicking on each area gives access to specific functionalities and information.

Areas become red in case of errors. See chapter [Errors on page 42](#).



Dyspc 6, Interactive instrument image

1 COMPONENTS AREA

Qcc af_nrcp Afcaiyle Amknmlc l r q r _ r s q m l n _ e c / 22,

2 LIQUIDS / WASTE AREA

Qcc af_nrcp K _ l _ eyle Jjosibq _ l b U _ q r c q m l n _ e c 71,

3 PIPETTING AREA

Qcc af_nrcp Afcaiyle Nyncryle ?pc_ Q r _ r s q m l n _ e c / / 7,

4 SAMPLE / REAGENT AREA

Qcc af_nrcp K _ l _ eyle Q _ k n j c q _ l b P c _ e c l r q m l n _ e c 42,

5 GEL CARDS STORAGE DRAWER

Qcc af_nrcp K _ l _ eyle Ecj A _ p b q m l n _ e c 36,

2,0,1 Qc _ p a f ? p c _

Specific test related to a sample can be found by inputting sample barcode.
Sample barcode can be typed with the keyboard or scanned with the barcode reader.

0 Only functional for loaded samples.

2,0,2

Rp_ac_`qjgrw ?pc_

A

Click **INFORMATION** from the main screen to access the **LOG LIST** screen. This screen allows to:

- display the traces in detail;
- filter by Level;
- to print the last 2 pages of traces.

Date / Time	Code	Action	Resources
30/08/2016 14:41	701	Authorized with login : USER , level 2	User="RESTRICTEDUSER"
30/08/2016 14:40	700	Request Display screen OPTIONS	User="RESTRICTEDUSER" Action="OPTI
30/08/2016 14:40	111	Instrument state : Ready	State="Ready"
30/08/2016 14:40	110	User "RESTRICTEDUSER" connected to RUNHCOM EXE	User="RESTRICTEDUSER"
30/08/2016 14:40	104	Loading data	Step=6
30/08/2016 14:40	208	Test 1920 - Phenotype: C,c,E,e,K,ctt (5010) in error: Pipetting abort	Failure="Pipetting abort" Barcode="19
30/08/2016 14:40	208	Test 1920 - Bloodgroup: ABO/Rh (DiaClon) (5001) in error: Pipetting abort	Failure="Pipetting abort" Barcode="19
30/08/2016 14:40	208	Test 1920 - Ab.screening: I,II,III (IAT) (5053) in error: "Needle Left clogged"	Failure="Needle Left clogged" Barco
30/08/2016 14:40	208	Test 1920 - Rev.group: A1,A2,B,O (O on separate Card) (5052/5052) in error: ...	Failure="Needle Left clogged" Barco
30/08/2016 14:40	208	Test 1920 - Direct Antiglobulin Test: (IgG) (5054) in error: "Needle Left clogg...	Failure="Needle Left clogged" Barco
30/08/2016 14:40	103	Test <Controleur.exe>	Step=4
30/08/2016 14:40	102	Loading interface <SystemSwingPlus>	Step=3
30/08/2016 14:40	100	Test <AccessIndirect.exe>	Step=1
30/08/2016 14:40	106	Stop IH-1000 System	
30/08/2016 14:40	106	Stop IH-1000 System	User="RESTRICTEDUSER" Action="EXIT"
30/08/2016 14:40	111	Instrument state : Ready	State="Ready"
30/08/2016 14:40	110	User "RESTRICTEDUSER" connected to RUNHCOM EXE	User="RESTRICTEDUSER"
30/08/2016 14:39	111	Instrument state : Waiting	State="Waiting"

Dyespc 7, Log list screen

B

Click **ARCHIVE** to archive one or several number of trace files within a range of time.

Rfc dmjjm ugle qapcc l qg bqnl wcb,

DATES SELECTION

Please select a beginning and an end date.

Start: 1/1/2011

End: 1/31/2011

OK CANCEL

Dyespc 1., Dates selection screen

C


Select the **Start** and **End** date.

D

Click **OK**.

?jj rp_acq _t _qj _`jc _pc am k npcqqcb rm rfc bcqirmn q l _ lcu dmjbcpl _ kcb ugrf
rfc b_rc md jme dyjcq \$JmeDjycq bb k k wwwwμ,

2,0,3 Q_k njcq _l b Pc_ec l r q Jm_b j l e ?pc_




Dyespc /1/, SAMPLES / REAGENT area

This area shows the current status of each rack loaded into the instrument.
Clicking on this image will open the **SAMPLE REAGENT AREA** screen, which gives more information.

0 The green bar indicates the preferential side (left or right) in the instrument to load additional samples (faster processing).
See chapter [Samples / Reagents Area Screen on page 77](#).

2,0,4 Dmmrcp Qrp gn

The footer strip displays the available options and information according to the current operation.
For the main screen, the following buttons are available:



Dyespc /0/, Footer strip on the main screen

- 1 Resources buttons
- 2 Tests buttons

2,0,4,/ Pcqm sp ac q @srrm l q

The resources buttons give access to the **RESOURCES** management screens.
See chapter [Managing Resources on page 57](#).

0 Resources to be charged and Resources on board are accessible from this menu.
Resources missing for one test in particular and Resources missing for all the tests are accessible only from the **TESTS WITHOUT RESOURCES** screen.

2,0,4,0

Rcqrq @srrm lq

These buttons give access to the **TESTS** management screens.

See chapter [Managing Tests on page 105](#).

The **TESTS** buttons allows to access the screens that list tests according to their state.

0 The tests screens are accessible only if the corresponding button displays a non-nil number of tests (otherwise, the button is disabled).

2,0,4,1

Mrfcp @srrm lq

Depending on the current operation, the following buttons appear on the footer strip.


	<p>@?A l @srrm l</p> <p>The BACK button is available on all the software screens (except main screen). It allows to go back to the previously consulted screen.</p>
	<p>NP&LR @srrm l</p> <p>When displayed, the PRINT button allows to print out the displayed datas.</p>
	<p>L_tge_rgm l icwq</p> <p>On some screens the navigation keys are available. It allows to scroll vertically and/or horizontally in the current screen.</p>

4.3 Warning Screen

This screen has a red triangular sign and is displayed when:

- an error on a test is detected;
- an authentication has failed;
- an acknowledgment or an out of service of a component is required by the user, an acknowledgment of an empty container, of a full liquid waste container or of the full gel cards waste bin is required by the user;
- the gel cards storage drawer is open;
- in some cases a confirmation by inputting the word «yes» is required.

? l _nnpmnpq_rc k cqq_ec q l dmp k q rfc sqcp q l c_af a_qc,



Dgespc / 1, Warning screen

4.4 Authentication Screen

This screen has a yellow padlock and allows special logging:

- to access screens with protected access;
- to carry out specific actions on the automated analyzer.

Dyespc / 2, Authentication required

4.5 Maintenance Reminder

O The weekly maintenance has to be done once a week (after 7 days). Otherwise an orange maintenance reminder (WEEKLY MAINTENANCE NECESSARY) will appear. See chapter [Maintenance on page 143](#).

4.6 Errors

When a fault occurs in the instrument, the concerned area is highlighted in red on the interactive image and a red informative area shows up on the screen (the **SAMPLE / REAGENT** area is always visible, it is colored in red in case of error).

Detailed information on errors are given by clicking either on the red area of the interactive image or on the red message.



Dyespc /3, Main screen - Faults

- 1 Error in the **SAMPLE / REAGENT** area

Qcc af_nrcp [K_I_ejle Q_knjcq_lb Pc_ec lrq ml n_ec 42,](#)

0 The green bar with the logo indicates on which side a sample rack has to be loaded to ensure a faster processing.

- 2 Error in the **COMPONENTS** area

Qcc af_nrcp [Afcaigle Amknmlclrq Qr_rsq ml n_ec /22,](#)

- 3 Error in the **PIPETTING** area

Qcc af_nrcp [Afcaigle Njncrnle ?pc_ Qr_rsq ml n_ec //7,](#)

- 4 Error in the **GEL CARDS STORAGE DRAWER** area

Qcc af_nrcp [K_I_ejle Ecj A_pbq ml n_ec 36,](#)

- 5 Error in the **LIQUIDS / WASTES** area

Qcc af_nrcp [K_I_ejle Jjosibq_lb U_qrcq ml n_ec 71,](#)

N?EC GLRCLRGML?JJW JCDR @J?L I

5 Getting Started

Аф_нрсп Мтсрпгсу

This chapter explains how to start using the instrument, including how to shutdown the instrument.

0 It is necessary to be logged in the Data Management software to be able to perform any operation on the IH-1000.


5.1 Starting the Instrument

3,/,/

Npcjg k g l _pw Afcaiq

Before starting a working day, the following checks and actions must be performed.


- A**Check that the room temperature is between 15°C and 27°C.
- B**Empty the gel cards waste bin (4).
Qcc af_nrcp Pcnj_aq le rfc Ecj A_pbq U_qrc @q l ml n_ec / . . ,
- C**Empty the liquid waste containers (3).
Qcc af_nrcp Pcnj_aq le _ Jqosqb U_qrc Am lr _q lcp ml n_ec 76,
- D**Fill the system liquid containers (2) with system liquid.
Qcc af_nrcp Pcdjjg le _ Qwqrc k Jqosqb Am lr _q lcp ml n_ec 73,
- E**Check that the window (1) of the pipetting area is completely closed.
Qcc af_nrcp Q_ knjq - Pc_ec lrq ?pc_ Qapcc l ml n_ec 55,
Qcc af_nrcp Pckmtc _jj P_aiq ml n_ec 7. ,
- 0**  At start up if racks are detected in the pipetting area, they are listed as being in error (in order to check their validity).
Before starting a working day remove ALL RACKS from the instrument.



Djespc /4, Starting the instrument

3./,0

Qr_prg l e Npmacbspc

AStart the instrument by positioning the **On / Off** switch to “I” (1).

Dgespc /5, Starting the instrument

B

Press and hold the power button (2) for a few seconds.

Rfc am k nsrcp qr_prg,
 ? epcc l jje fr _pms l b rfc nmucp `srrm l &O' q l b q a _rcq rf _r rfc am k nsrcp q q
 nmucp cb m l,

C

Double-click on the icon to start the IH-Com Data Management Software.

Qcc G F-Am k B _r _ K _l _eck c l r Sqcp K _l s _j,
 Rfc G F-Am k B _r _ K _l _eck c l r q l rcp d _ac q q b q q n j _wcb,

DLogin to IH-Com with your **User ID** and **Password**.**E**Start the IH-1000 Software: click on the **IH-1000** button in the IH-Com Data Management interface.

Qcc G F-Am k B _r _ K _l _eck c l r Sqcp K _l s _j,
 Rfc q l q rps k c l r qa _l q _l b q bc l r q d q c q _j j ec j a _pbq nrcq c l r q l rfc bp _ucp bsp q l e
 rfc q l q r q _j q x _r q m l,

0

The IH-1000 Software can also be started by double-clicking the IHM1000 icon on the desktop.

F The IH-1000 software checks if there were errors during the initialization process.

*Gd_ l cppmp qg bcrarc8 rfc ?jcpr qapcc l lcvr qg bqgnj_wcb qfmuq le rfc jqr md
q l qrg_jqx_rgm l cppmpq,
Ajjai MI rm ajmqc rfc GF- / . . . qmdru_pc,*



Dyespc /6, Alert screen

5.2 Exiting the Software

This function shuts down the user interface; the master computer remains started.

A From any screen, click on **MENU**.

B Click on **EXIT**.



O The IH-1000 user interface will be closed immediately but the instrument will keep on working.

C To restart the software, double click the corresponding icon on the desktop or in the Data Management software.

Qcc GF-Amk B_r_ K_l_eckc l r Sqcp K_l s_j,

5.3 Stop All of the Instrument

This button stops as soon as possible all the tests in progress on the instrument.

A	From any screen, click on MENU .
B	<div style="display: flex; align-items: center;">  <div> <p>Click on STOP ALL.</p> <p><i>Rfc_srfclnq_rgm qapccI _nnc_pq,</i></p> </div> </div>
<p>Dyesc 17, Stop All instrument - Authentication</p>	
<p>I @w ajgaig le mI STOP ALL* _jj rcqrq gI npmepcqq _pc jmqr _I b gr gq lmr nmqqg`jc rm am lrg lsc, Rfc mljw qmjsrgmI gq rm qfsr bmu l rfc glqps kcl r ugrf rfc <b>SHUTDOWN</b> `srrmI,</p>	
C	<p>Confirm the STOP ALL procedure:</p> <ul style="list-style-type: none"> Enter your login and password; Click on OK.
D	<p>SHUTDOWN the instrument.</p> <p><i>Qcc af_nrcp Qfsrbmu l md rfc lqpskcl r ml n_ec 27,</i></p>

5.4 Shutdown of the Instrument

This function allows to completely shutdown the instrument.

<p>I ?ju_wq qugraf mdd rfc GF-/... ugrf rfc SHUTDOWN dslargmI, Lctcp sqc rfc eclcp_j nmucp qugraf* cvacnr gl a_qc md _I ck cpeclaw,</p>	
<p>I Gd rfc glqpskcl r gq `jmaicb* qfsrbmu l ugrf rfc Pcqasc Qfsrbmu l `srrmI &pcdcp rm af_nrcp Pcqasc Qfsrbmu l ml n_ec 3.,</p>	
A	From any screen, click on MENU .

B Click on **SHUTDOWN**.



0 The shutdown command is not possible if any test is on-going in the IH-1000. Use **STOP ALL** function first.

See chapter [Stop All of the Instrument on page 49](#).

0 Do not remove diluent racks manually before shutdown. The system will remove automatically all racks. In case of manual remove of the diluent racks before shutdown, there is a risk of losing data.

C Confirm or not the shutting down procedure:

Go WCQ8 rfc q lqpskclr qfsrq bmu l qd lm rcqrq _pc ml+emjle9

Go LM8 rfc npmacbspc qg a_lacjcb,

Diespc 0., Shutdown of the instrument - Confirmation

5.5 Rescue Shutdown

0 Only use the Rescue shutdown button if it is not possible to shutdown the instrument through the IH-1000 software (according explanations in chapter [Shutdown of the Instrument on page 49](#)).

A Double-click on the **Rescue Shutdown** icon on the desktop.



Rfc GF-/. . . Pcqasc Qfsrbmu l qapcc l _nnc_pg,



Diespc 0I, Rescue Shutdown - Confirmation

0 If the icon is not on the desktop, ask your service engineer to install it.

B Confirm the **Rescue Shutdown** procedure:

- Type the password: **Shutdown**
- Click on **OK**.

Rfc ?srm k_rc am k nsrcp qfsrq bmu l,

0 The progress is displayed in the bar.

U_qr s l n j q r f_q q f s r b m u l a m k n j c r c j w l j c q q r f_l 3 k q l s r c q',


C If IH-Com is running, quit the software.

Qcc G F-A m k B_r_ K_l_e c k c l r S q c p K_l s_j,

D Shut down the Master PC.

5.6 Step by Step Procedure

A Start the instrument by positioning the **On / Off** switch to “I” (1).



Dgespc 00, Starting the instrument

B Press and hold the power button (2) for a few seconds.

Rfc am k nsrcp q r_p r q,


? e p c c l j g e f r _ p m s l b r f c n m u c p ` s r m l & O' q l b q a _ r c q r f_r r f c a m k n s r c p q q n m u c p c b m l,

Rfc q l q r p s k c l r q a _ l q _ l b q b c l r q d j c q _ j j e c j a _ p b q n p c q c l r q l r f c b p _ u c p b s p q l e r f c q l q r q _ j q x _ r q m l,

C Double-click on the icon to start the IH-Com Data Management software.

D m p b c r _ q j c b q l d m p k _ r q m l p c d c p r m r f c G F-A m k B_r_ K_l_e c k c l r S q c p K_l s_j,
Rfc G F-A m k B_r_ K_l_e c k c l r q l r c p d _ a c q q b q q n j _ w c b,


D Login to IH-Com with your **User ID** and **Password**.

- E** Start the IH-1000 Software: click on the IH-1000 button in the IH-Com Data Management interface.
- Dmp bcr_gjcb g l dmp k_rgm l pcdcp rm rfc GF-Am k B_r_K_l_eckc l r Sqcp K_Is_j, Rfc g l qrp sk c l r g l g r g_jpxcq _srnk _rja_jjw, Gr qa_lq _lb g bcl r g d c q _jj ecj a_pbq n p c q c l r g l r f c b p _ucp b s p g l e r f c g l g r g_jpx_rgm l,*
- 0** The IH-1000 Software can also be started by double-clicking the IHM1000 icon on the desktop.
- F** Wait until the **System status** displays **Ready** (top right of the main screen)
- Djespc 01, information area - System status*
- G** Make sure that no area on the interactive image is highlighted in red. (red = error in the specific area).
- Gd m l c md r f c _pc_q g q pcb* qcc Cppmpq m l n_ec 20,*
- H** Click on the **DATA MANAGEMENT** button to switch back to the IH-Com Data Management Software.
- I** In the IH-Com, go to the **Work list** screen area.
- G l r f g q q a p c c l _j j r f c q _k n j c q u y r f _l m n c l m p b c p & q c l r `w J g Q m p k _l s _j j w a p c _r c b g l G F - A m k ' _p c b y q n j _w c b,*
Dmp bcr_gjcb g l dmp k_rgm l pcdcp rm rfc GF-Am k B_r_K_l_eckc l r Sqcp K_Is_j,
- 
- Djespc 02, IH-Com - Work List screen area*

J Click on the IH-1000 button to switch back to the IH-1000 Software.

K Place the sample tubes correctly on the foreseen sample rack.


*K_ic qspc rf_r rfc a_nq _pc pc kmtcb,
Dmp kmpc bcr_gjcb g l dmp k _ngm l* qcc Q_knjc P_aiq ml n_ec 42,*



Dyesc 03, Sample rack

L Load the racks.


Qcc af_nrcp Jm_byle P_aiq ml n_ec 64,



Dyesc 04, Loading racks

O A green light (1) indicates that the slot is available.
A blinking red light indicates that a rack is moving: loading or unloading.
A red light indicates that the slot is not available. A rack is already loaded. Sample and diluent racks can be loaded in positions 1 to 18. Reagent racks can be loaded in positions R1 to R4. Engage the rack in the slot and push it slowly until it locks automatically. (2) Shows a well-placed rack.

I *?r_lw rjkc* rfc p_ai ksqrj_w ml rfc `mrrmk md rfc qjmr gl mpbcp rm `c jmaicb _r rfc
amppcar nmqgrml,
Gd rfc p_ai gq nsqfcb rmm d_p* _amjjgqgm l ugrf rfc rp_lqnmpr _p k amsjb maasp,*

M	<p>If the button TESTS WITHOUT RESOURCES is active, one or more resources are missing to start the test.</p> <p><i>Dmp k mpc g / dm* qcc af_nrcp Rcqrq u y r f m s r P c q m s p a c q m l n _ e c / / 1,</i></p>
N	<p>Click on MISSING RESOURCES to see which resources are needed to start the test.</p> <p><i>Dmp k mpc g / dm* qcc af_nrcp A f c a i g l e P c q m s p a c q Q r _ r s q m l n _ e c / . / ,</i></p>
O	<p>Load the required resources.</p> <p><i>Qcc af_nrcp K _ l _ e g l e E c j A _ p b q m l n _ e c 36,</i> <i>Qcc af_nrcp B j j s c l r P _ a i _ l b B j j s c l r N p m b s a r m l n _ e c 50,</i> <i>Qcc af_nrcp P c _ e c l r P _ a i q m l n _ e c 52,</i></p>  <p>Dyespc 05, Resources</p>
P	<p>As soon as all resources needed for a specific test are on board, the test is displayed on the TESTS READY TO START list.</p> <p><i>Dmp k mpc g / dm* qcc af_nrcp Rcqrq p c _ b w r m q r _ p r m l n _ e c / / / ,</i></p>
Q	<p>After a few seconds the tests switches to TESTS IN PROGRESS and the instrument starts the run.</p> <p><i>Dmp k mpc g / dm* qcc af_nrcp Qcc af_nrcp Rcqrq g / N p m e p c q q m l n _ e c / . 6,</i></p>
R	<p>If the test process has finished without any error, the test disappears from the TESTS IN PROGRESS list.</p>


S If the test process has finished with an error, the test is displayed on the **TESTS TO REPEAT** list.

*g l r f _ r a _ q c * q c c a f _ n r c p R c q r q r m P c n c _ r m l n _ e c // 2,*

T If the specific test has been processed without any error, click on the **DATA MANAGEMENT** button to get to the IH-Com Data Management Software and to see the results.

U In the IH-Com, go to the **Results** screen area to be able to validate the result.

D m p b c r _ g j c b g l d m p k _ n y m l p c d c p r m r f c g F - A m k B _ r _ K _ l _ e c k c l r S q c p K _ l s _ j,



Dyescp 06, IH-Com - Results screen area

V When all the tests are processed, shut down the instrument from the IH-1000 software.

*D m p k m p c b c r _ g j c b g l d m * q c c Q f s r b m u l m d r f c g l q r p s k c l r m l n _ e c 27,*

I *? j u _ w q q u g r a f m d d r f c g F - / . . . u g r f r f c S H U T D O W N d s l a r g m l ,*
*L c t c p s q c r f c e c l c p _ j n m u c p q u g r a f * c v a c n r g l a _ q c m d _ l c k c p e c l a w,*

O The shutdown command is not possible if any test is on-going in the IH-1000. Use **STOP ALL** function first.

W

Confirm or not the shutting down procedure:

Go **WCQ8** rfc q lqpskc l r qfsrq bmu l qd l m rcqrq _pc m l +emg l e9

Go **LM8** rfc npmacbspc q q a_ l acjcb,

Dyespc 07, Shutdown of the instrument - Confirmation

6 Managing Resources

Аf_nrcp Mtcptqcu

This chapter presents the resources needed to perform the tests.

6.1 Managing Gel Cards

4, /, / Bcqapgnrgm l




Figure 1. Loading gel cards

- 1** Gel card
- 2** Protection foil
- 3** Gel cards tray
- 4** Gel cards storage drawer

Gel cards are stored in trays, each tray can contain up to 12 Gel cards.

The gel cards trays are stored in the instrument drawer, the gel cards storage drawer can contain up to 20 trays.

0 Take care that the gel cards are protected from dust during the storage.

4./,0 Ecj A_pbq Qrmp_ec Bp_u cp Qapcc l


This screen displays a detailed status of the gel cards loaded into the gel cards storage drawer of the instrument.

4./,0,/ ?aacqq

It can be accessed by the following ways:

- From the main screen, click on the **GEL CARDS STORAGE DRAWER** of the interactive image.
- From any screen where the **GEL CARDS STORAGE DRAWER** button appears.

4./,0,0 Bcqapgnrgm l



Dyespc 1/, Gel cards storage drawer screen

Qw k `mj

Bcqapgnrgm l



Lmp k _j ecj a_pb rp_w

Example:

It contains 12 gel cards out of a possible 12

The type is "Liss/Coombs"

The barcode is 5053

The tray is in location #11 in the drawer

C k nrw rp_w jma_rgm l



Rp_w ugrf _l cppmp

When clicking on the tray with an error, a message appears to explain the error.

Qw k `mj Bcqapnrgm l

Rp_w bq_`jcb



Rmr_j os_lrgw md c_af rwnc md ecj a_pb

Sljmaiq rfc bp_ucp j_raf

GEL CARDS BY TYPE

10	5000 - ABO/Rh
12	5011 - DiaClon Rh-subgroups+K
24	5049 - DiaClon ABO/D + DAT
7	5052 - NaCl Enzyme test and cold a...
15	5053 - LISS/Coombs
24	5096 - DiaClon ABO/Rh for Newborns...

4,/,1

Mnclgle - Ajmqgle rfc Bp_ucp

AClick on **OPEN DRAWER** to request the opening of the drawer.

Rfc q lqrp skclr uqjj dylqgf _jj ml+emgle _arqm lq ugrf ecj a_pbq `cdmpc rfc bp\_ucp a\_l `c mnclcb, Gr r_icq _dcu qcamlbq

Ufcl nmqqj`jc\* rfc bp\_ucp j\_raf yq sljmaicb \_lb rfc bp\_ucp a\_l `c mnclcb k_ls_jjw,

- Ufcl mnclgle rfc bp_ucp* r_ic a_pc rf_r ecj a_pbq _lb rp_wq bm lmr kmtc bspgle rfc mncl_rgm l,

? kccq_ec _nnc_pq ylbqa_ryle rf_r yq pcospcb rm ajmqc rfc ecj a_pbq qrmpec bp_ucp _pc _q qmm l _q nmqqj`jc,

**Despc 10,** Wait for drawer unlocking**O**

To load gel cards, the drawer must be completely pulled out.

B

Load the gel cards.

Qcc Jm_byle Ecj A_pbq ml n_ec 40 dmp bcr_gjcb y l dmp k_rgm l,

C

Close the drawer by pushing it manually until it is completely closed.


Rfc q lqrp skclr qa_lq _lb ybc lrgdcq _jj jm_bcb ecj a_pbq,

- Ufcl ajmqgle rfc bp_ucp* r_ic a_pc rf_r ecj a_pbq _lb rp_wq bm lmr kmtc bspgle rfc mncl_rgm l,

D Click on **OK**.

0 If the gel cards storage drawer is in error, the **OPEN DRAWER** button is in red.

A Warning screen with the error is displayed.



4, /, 2 Jm_bgle Ecj A_pbq

A Prepare trays with gel cards.

1 Afcai rf_r rfc nprmcargm l dmjg &/' md c_af ecj a_pb gq lmr b_k_ecb,
Bm lmr jm_b mnc lcb mp n_prq_jjw sqcb ecj a_pbq, Rfc pgefr dpcc c l b md rfc nprmcargm l dmjg
k sqr `c fmpgxm l r_j mp jmucp* lmr dmjbcb mp snqr_lbgle,
Afcai rf_r rfc `_pambc md c_af ecj a_pb gq ajc_l _l b lmr qap_rafcb, Gd l cacqq_pw* ajc_l
rfc `_pambc ugrf _jg l r dpcc p_e,

0

Make sure that:

- each tray only contains one type of gel cards;
- the barcodes (2) are correctly oriented (the barcode must be visible with the free end of the protection foil (1) showing to the right of instrument);
- the gel of the gel card is in good condition;
- there are no air bubbles present;
- the expiry date of the gel card is valid (checked by the instrument as well).

B Open the drawer.

Qcc af_nrcp *Mnclgle - Ajmqgle rfc Bp_ucp ml n_ec 4.*,

1 Ufcl mnclgle rfc bp_ucp* r_ic a_pc rf_r ecj a_pbq _l b rp_wq bm lmr kmtc bspgle rfc
mncp_rgm l,

If there are unsealed gel cards in the drawer:

)

INFECTION

C_af ucjj ml c_af ecj a_pb a_l am l r_g l sn rm /03 Öj md `jmmb q_k njc qmjsrgm l,
Rfcpc gq _pgqi md g ldcargm l dpm k qig l am l r_ar ugrf `jmmb, ?ju_wq uc_p nprmcargtc ejmtcq
ufcl umpigle* gl _aampb_lac ugrf emmb j_`mp_rmpw np_argacq,

C Put the trays in the drawer.

Rp_wq a_l `c nj\_acb gl \_lw\_t\_gj `jc nj_ac gl rfc bp_ucp,

1 Rm `c \_`jc rm ajmqc rfc bp_ucp* _jj ecj a_pb rp_wq k sqr `c amppcarjw gl qcprcb,

0

If the option **RETURN OF THE CARDS IN DRAWER FOR SECOND READING** is set to **ALWAYS** or **CUSTOMIZED**, the gel cards will be put back to the positions 4, 5, 9, 10, 13, 14, 15, 18, 19, 20.

If **ONLY IF UNSPECIFIED RESULT** is set, the gel cards will be put back to positions 15 and 20.

Make sure that these positions are available.

0

To ensure a proper identification of the empty trays, use the white trays (REF 009888) for the cards returned to the drawer.

D



Close the drawer: manually push it until it is completely closed.

Rfc q l q r p s k c l r u q j j _ s r m k _ r q a _ j j w q a _ l r f c d q p q r e c j a _ p b q m d c _ a f e c j a _ p b r p _ w ,
 ? j j r f c e c j a _ p b q m d m l c r p _ w _ p c _ q q s k c b r m ` c q b c l r q a _ j r m r f c d q p q r m l c ,
 C _ a f e c j a _ p b r w n c q q b m s ` j c a f c a i c b u f c l s q c b , G d r f c r w n c q q b q d d c p c l r r f _ l
 r f c m l c c v n c a r c b * _ u _ p l q l e k c q q _ e c _ n n c _ p q q l r f c **ECJ A?PBQ**
QRMP?EC BP?UCP q a p c c l q r f c p c k _ q l q l e e c j a _ p b q _ p c _ q q s k c b r m ` c m d
 r f c q _ k c r w n c _ q r f c j _ q r e c j a _ p b ,

- U f c l a j m q q l e r f c b p _ u c p * r _ i c a _ p c r f _ r e c j a _ p b q _ l b r p _ w q b m l m r k m t c b s p g l e r f c m n c p _ r g m l ,



Sealed gel cards should not stay longer than 7 days in the IH-1000.

**Diagram 11.** Put the trays in the drawer

6.2 Managing Samples and Reagents

Three types of rack can be loaded in the IH-1000: sample, diluent and reagent racks.

- Racks containing samples or diluents can be loaded in positions 1 to 18 of the rack area.
- Racks containing reagents can be loaded in positions R1 to R4.
- Racks are identified with barcodes.

4,0,/

Q_knjc P_a i q




Figure 12. Sample rack - Description

- 1** Sample tube
- 2** Sample barcode
- 3** Rack identification number and barcode
- 4** Sample tube numbering
- 5** Handle

A sample rack can contain up to 10 sample tubes.

Sample racks can be introduced from position 1 to 18 in the racks area.

The barcode of the tube (1) has to be visible on the left side. If the barcode of a sample can not be read, the instrument shows an alarm of an unreadable or non detected barcode. See chapter [Sample Status on page 78](#).

Racks processing order does neither depend on the position of the racks nor on the racks identification number; the first rack in will be processed first.

The position of the sample in the rack defines which sample will be processed first; the sample in position 1 will be processed before the sample in position 2.

0 Use only centrifuged samples.


0 The reliability of the results depends on correct compliance to the Good Laboratory Practices for reagents and samples.

0 The best performance will be obtained when using anticoagulated samples within 48 hours of collection, however older samples may be used. When testing with donor segments, the samples can be used up to the expiration of the donor unit.


▪ Rfc p_ai q _l b rs`cq ksq r`c ajc_l,
Rfc`\_pam bcq md rfc p\_ai \_l b rs`cq ksq r`c ajc_l _l b lmr qap_rafcb,
Gd l cacqq_pw* ajc_l rfc`_pam bcq ugrf _jg l r dpcc ajmrf,

▪ Lctcp rpw rm jm_b _q_ knjc p_ai g l nmqgrgm l P/ rm P2 md rfc jm_bgle _pc_,

▪ Rm npctcl r os_l rgrw bcrcargm l cppmpq* l m dm_k mp`s`jc
ksq r`c npcqc l r _r rfc qsp d_ac md rfc jg osyb,
Pckmtc q_ knjc rs`cq a\_nq`cdm pc jm_bgle rfc p_ai g l rm
rfc GF-/... ,
Rm npctcl r mtc pdjm u u fcl ngncrrgle* r_ica _pc md rfc
k_vg ksk dgj jgle jctcj md rfc q_ knjc rs`c,



▪ Rfc rs`cq ksq r`c g l qcprcb amppcarjw g l rm rfcgp nj_acq,
Rfc rs`cq ksq r`c tcrpa_j9 afcai rfc amppcar dslargm lgle
md rfc esybg l e qnpg l eq,



0 Unload all the samples at the end of the testing.

4,0,/,/

Sqgle Qr_l b_pb Rs`cq

0 Standard tubes must be used only with racks of type 1 (Barcode 01 or 11) and racks of type 4 (Barcode 04 or 14). See chapter [Labeling on page 68](#).
See chapter [Cap Detection on page 139](#).

4,0,/ ,0 Sqg l e Jmu Tmjs k c Rs `cq & dmp Pcb Acjjq M l jw'


Low volume tubes (for red cells only) must be used for small quantity sample.

Use SARSTEDT products:

- 1 Micro tube 1.5mL ref. 72.703
- 2 Rounded end tube 5mL ref. 55.475.005

1 Rfc k g apm rs `c k sqr `c am k njrcjw g l qcprcb g l rm rfc rs `c,

0 Low volume tubes (for red cells only) must be used only with racks of type 2 (Barcode 02 or 12). See chapter [Labeling on page 68](#).



Dgespc 13, Low volume tubes (for red cells only)

4,0,/ ,1 Sqg l e Njs l ecp Rs `cq

0 Plunger tubes must be used only with racks of type 3 (Barcode 03 or 13). See chapter [Labeling on page 68](#).

1 Gd rfcqc am l b g r g m l q _pc l m r pcqncarc b* rfcpc gq _pgqi rm ap_qf rfc lccb j c,

4,0,/ ,2 Sqg l e Jmu Tmjs k c Qce k c l r @jmb

0 Low volume segment blood must be used with Standard tubes and racks of type 4 (Barcode 04 or 14). See chapter [Labeling on page 68](#).

Segment blood is the blood contained in the pipe linked to the blood bag. This blood is used to check the blood bag.

A Cut segments of this pipe to get the blood inside.



Figure 14, Segment bloods

0 The volume of segment blood in the tube is $200 \mu\text{l} \pm 15\%$ (170 - 230 μl). It corresponds to a length of segment between 7 to 9 centimeters.

▪ Rm_tmgb rm ap_qf rfc lccbjc* _ju_wq sqc rfc q_k c rwn c md rs`cq,

▪ K_i c qspc rm jm_b rfc q_k njcq qfmpjw _drp k_gv l e-fm k mec l_gx l e,

0 The 2nd cell suspension is not allowed for the samples in racks of type 4.

0 When using rack type 4 for segments of blood bags on the IH-1000, an adequate amount of red blood cells is required to prepare a correct red blood cell suspension.

Qncagdyagrcq dmp Rcqrq

- Tests not allowed for segment blood are not processed by IH-1000.
- Only one aspiration or one dilution is allowed per sample tube. If a test requires several aspirations or dilutions, only the first one is processed.
- If resources are missing: Only the tests with enough resources are done, the other tests are deleted.
- Samples centrifuged and/or with plasma required for the test have to be loaded on another type of rack.

Mrfcp_uq_c _l cprmp k cq_q_ec q q b_qnj_wcb q l rfc jme, Ce8 **Rcqr-Npmdjic** lmr _jjmucb dmp qce k c l r
`jmmb q_k njc Yqctcp_j _qngp_rgm l q[.

Qcc af_nrcp Rq_ac_`qjrw ?pc_m l n_ec 17,

If a default profile is configured, it is ignored for this racks type. The samples are displayed in the **SAMPLES WITHOUT REQUEST** list. See chapter [Profile Options on page 132](#).

Level detection is disabled for any sample tube with Low volume segment blood on rack of type 4. See chapter [Liquid Level Detection on page 32](#).

1 K_i c qspc rf_r_jj pcqm spacq _pc ml `m\_pb `cdmpc jm_bg le _p_ai md rwnc 2, ?q mljw ml c _qngp_rgm l a_l `c bml c ncp q\_k njc rs`c* rfc rcqr amsjb `c a_lacjcb mp mljw n_prg_jjw npmacqqcb,

4,0,/,3

J_`cjgle





Diagram 15. Description - Labeling

- 1** The first character defines emergency levels of rack:
 - 0: Not urgent;
 - 1: Urgent.
- 2** The second character defines the type of rack:
 - 1: Standard Spherical bottom tube;
 - 2: Low volume tubes (for red cells only);
 - 3: Plunger tube;
 - 4: Segment blood (in Standard Spherical bottom tube);
 - 5 and 6: Not defined.
- 3** Example:
 - 12 = Urgent, Low volume tube (for red cells only);
 - 03 = Not urgent, Plunger tube.



4,0,/,4 C k cpec law Q_ k njc P_ ai

The emergency sample rack is similar to a standard sample rack except that it has a special barcode that confers a top level priority. If the required resources are available, all the samples of the emergency rack will be processed before the samples of a standard rack. A message on the GUI indicates the expected time to start the test.

The emergency rack is identified with a red label.

Standard samples racks instructions apply for this rack.

See chapter [Sample Racks on page 64](#).

0 To get results faster, try to load all the emergency samples in the same rack (if possible).

4,0,/,5 Sqg le _ Lcu P_ ai

Stickers with barcodes must be put on a new rack so that it can be used in the instrument.




Figure 16, Stickers on samples racks

4,0,/,6 Qncagda_rgm l q dmp Rs`cq


Qf_nc	Cylindrical with spherical, flat or conical bottom
Msrcp bq_k crcp	Between 10 and 17mm
Gllcp bq_k crcp	Between 8 and 15mm
Fcgefr	Between 66 and 100mm
K_rcpu_j	Glass, Plastic or polyethylene

4,0,/,7

Qncagda_rgm l q dmp @_pambcq

Types of barcodes read by the instrument:

- Code 39, 93, 128
- Interleaved 2 of 5
- EAN-8
- CODABAR with control character suppressed
- UCC-EAN 128 with control character suppressed
- ISBT 128 with specific characters
- EAN13 (equal to UPC-A 13)



Dyespc 17, Barcode dimensions

Gr gq pcam k k c l bcb rm sqc qcaspbc ` _pambc rwnccq,
 Dmp l m l + qcaspbc rwnccq & Amb _ ` \_p\* Ambc 17\* Ambc 71\* G lrcpjc \_ tcb 0 md 3\* rfc ` _pambc k sqr
 f _tc _r jc _qr 2 af _p _arcpq g l mpbcb rm npctc l r am l dsqgm l q, Rfc j _ `cj k sqr ` c am k njrcjw
 qrsai _pms l b rfc rs `c,
 Gd rfcqc am l bgrgm l q _pc l mr pcqncarc b* rfcpc gq _pgqi l mr rm bcrcar rfc a _n _l b rm ap _qf rfc
 l ccbjc,

Bcqe l _rgm l	Pcd,	Bcqaq nrgm l	Am k k c l r q
Am jmp		Black with white background	
Ls k `cp md		max. 20	
af _p _arcpq			
Ep _bc		B min. recommended	Print quality
K mbsjc		0.2 to 0.35mm	Width of single bars
P _rgm		from 2.2 to 3	Relation between small and wide bars
@j _l i qn _ac	c	Min. 5x Module size x Ratio For code 128 only: 7mm with a «Module» of 0.35 (0.35x4x5) (Code incl. Up to 4 different bar widths)	e.g.: 4.5mm for a Module of 0.3mm with a Ratio of 3
Fcqe fr md ambc	F	min. 5mm	
Jc l erf md ambc	J	min. 35mm	
Rmr _j jc l erf	j	max. 70mm	
@_pambc nmqrgm l 8 `mrrm k	`	min. 10mm for rounded end tubes min. 13mm for flat end tubes (plunger tubes)	
@_pambc nmqrgm l 8 rmn	r	min. 15mm	

4,0,/,/.

Q_k njc Rs `cq Bc_b Tmjs k c

- Dead volume standard tube 13*75 Terumo: 350µl
- Dead volume plunger tube 11*66 Sarstedt: 500µl
- Dead volume low volume tubes (for red cells only) 1,5ml Sarstedt: 50µl

4,0,/,,/ F- / . . . Ac l r p d s e c b Q _ k n j c R s ` c q K g l g k s k N j _ q k _ F c g e f r

L c a c q _ p w n j _ q k _ f c g e f r & k k '		R s ` c b g _ k c r c p & k k '						
		/ /	/ 0	/ 1	/ 2	/ 3	/ 4	/ 5
B g q r p g ` s r c b t m j s k c & O j'	03	7	6	6	6	6	6	6
	3.	7	7	6	6	6	6	6
	53	7	7	7	6	6	6	6
	/ . .	8	7	7	7	6	6	6
	/ 03	8	8	7	7	7	6	6
	/ 3.	8	8	7	7	7	7	6
	/ 53	9	8	8	7	7	7	7
	0 . .	9	8	8	7	7	7	7
	003	9	9	8	8	7	7	7
	03.	10	9	8	8	7	7	7
	053	10	9	8	8	8	7	7
	1 . .	10	9	9	8	8	7	7
	103	11	10	9	8	8	8	7
	13.	11	10	9	9	8	8	7
	153	11	10	9	9	8	8	8
	2 . .	/ 0	10	10	9	8	8	8
	203	12	11	10	9	9	8	8
	23.	12	11	10	9	9	8	8
	253	12	11	10	9	9	8	8
	3 . .	13	11	10	10	9	9	8
	303	13	12	11	10	9	9	8
	33.	13	12	11	10	9	9	8

0 Example: When using 11 mm diameter tubes and when 400 µl of plasma is needed, the height must be at least 12 mm.


4,0,/,,/ 0 P c b @ j m m b A c j j q K g l g k _ j F c g e f r

B g k c l q g m l & k k '	Q r _ l b _ p b r s ` c		N j s l e c p r s ` c q				j m u t m j s k c r s ` c q & d m p p c b a c j j q m l j w'
B g _ k c r c p	13	16	11	13	15	16	-
P c b @ j m m b a c j j q f c g e f r	6	6	9	10	9	9	5

The values are valid for the tubes of type:

- BD diagnosis;
- Terumo and Greiner for standard tubes;
- Sarstedt for plunger tubes;
- Sarstedt ref: 72.703 for low volume tubes (for red cells only).

4,0,0 Bgjsclr P_ai _lb Bgjsclr Npmbsar



Dyespc 2., Diluent rack - Description

- 1** Diluent rack
- 2** Rack guiding block
- 3** Diluent rack product
- 4** Identification sticker
- 5** Protection foil
- 6** Handle

The diluent rack can hold diluent 1 or diluent 2 products.

The barcode has to be visible in the left side.

The protective foil of the diluent rack product must not be damaged before loading; each unit will be pierced when necessary.

0 It is allowed to read sample barcodes manually, but not diluent barcodes.

! Afcai rf_rfc ` \_pambcq md rfc p\_ai \_lb md rfc bgjsclr p\_ai npmb sar \_pc ajc\_l \_lb lmr qap\_rafcb, Gd l cacqq\_pw\* ajc\_l rfc ` _pambcq ugrf _jg l r dpcc p_e, Bm lmr sqc _ sqcb bgjsclr p_ai &ngcpacb ucjjq' ml _l mrfcp glqrpsk clr,

! Lctcp rpw rm jm_b _ bgjsclr p_ai gl nmqgrgm l P/ rm P2 md rfc jm_bg le _pc_,



When closed, the diluent rack products must not stay more than 7 days in the IH-1000.

4,0,0,/ Bgjsclr / P_ai Npmbsar Glqpsargmlq

U fcl sqgle rfc bgjsclr / p_ai npmbsar* pcqncar a_pcdsjw rfc dmjjmugle glqpsargmlq,

- The diluent 1 rack product must be stored upright, protection foil on top, never up-side down.
- The diluent 1 rack product must be stored between +4°C/8°C when not used.

Rfc bgjsclr ksqr `c qrmpcb \_r pmm k rckncp\_rspc dmp \_r jc\_qr 13 kg lsrcq `cdmpc gr gq sqcb gl rfc glqpskclr,

- Before using the product, apply a little shock against a table to make sure that there is no air bubble trapped in the wells.
- The SETUP CLEAN product (ref. 009902, dilution 1 to 5) is required when using the diluent 1 rack product. It must be loaded in a reagent vial, on reagent rack.

4,0,0,0 Glqcpqgle _ Bgjsclr P_ai Npmbsar

The diluent rack product sticker must be on the same side as the rack guiding block.

The diluent rack product must be completely inserted into the rack.

Afcailrf_rfc bgjsclr p_ai npmbsar qc_rq nrmncpjw gl rfc p_ai* ugrf lmm`qr_ajc,
Gd rfc bgjsclr p_ai npmbsar gq lmr amppcarjw glqpcrcb* rfcpc ugjj `c _ amjjgqgm l ugrf rfc ngncrrmp,





Diagram 21, CAUTION - Height of rack

4,0,1 Pc_eclr P_aiq



Dyesc 20, Reagents rack - Description

- 1 Reagent vial
- 2 Reagent barcode
- 3 Rotating reagent place
- 4 Rack identification barcode
- 5 Reagent vial numbering
- 6 Handle

A reagent rack can contain up to 7 reagent vials.

The reagent vial barcode has to be visible through the aperture of the reagent place.



Reagent racks can be loaded in position R1 to R4.

- 0 It is allowed to read sample barcodes manually, but not reagent barcodes.
If a reagent vial of a reagent set is replaced by a new one, make sure to have the same lot number.


1 Rfc ` \_pambc qrgaicpq md rfc p\_ai \_lb tq\_jq ksqr lmr `c qap_rafc b&c,e, `w rslgle rfc tq\_j gl rfc fmjbcp' gl mpbcp rm `c gblrgdgc b amppcarjw `w rfc glqrps kclr,

2 Dmjjmu rfc glqrpsargm lq md sqc dpm k rfc pc_eclr k_lsd_arspcp,

3 Gl mpbcp rm npctclr os_lgrw bcrcargm l cppmpq* l m dm_k mp `s`jc ksqr `c npcqlr_r rfc qspd_ac md rfc jgosgb,
Rfc acjjq md rfc pc_eclrq ksqr `c gl qsqnc lqgm l `cdmpc rfc tq_j qj m_bcb m l rfc p_ai,
Rfc a_n md rfc pc_eclr tq_jq ksqr `c pckmtcb `cdmpc jm_bg l e rfc p_ai gl rm rfc GF- / ... ,



- Rfc tq_jq ksqr dgr amppcarjw glrm rfcgp nj_acq9 rfcpc ksqr `c lm _lejc nj_w,



- The vials have to be taken out of the fridge at least 30 minutes before they are used on the instrument in order that they reach room temperature. The reagent racks should not be stored in the fridge.



MI+`m\_pb qr\_`ggrw gq8

² 0 b_wq &26 fmspq' dmp pcb `jmmb acjj pc_ec lrq9

² 5 b_wq &/46 fmspq' dmp _lrgqcp _pc_ec lrq &GB+? lrg K* L* Q* q* Dw_* Dw`',

Rfc ml+`m\_pb qr\_`ggrw bmcq lmr glajsbrc rfc r9kc u fcl rfc tq_j gq qrmpcb gl rfc dpgbec, Gr glajsbrcq mljw rfc r9kc u fcl rfc tq_j gq gl rfc glrps kclr,

?drpc rf_r r9kc* rfc pc_ec lrq ksqr `c pc kmtcb _lb bqnmqcb,

4,0,1,/

Sqc md rfc Ajc_lgle QmjsrgmI QCRSN AJC?L

The SETUP CLEAN product (REF 009902, dilution 1 to 5) is required when using Enzymatic reagents and/or antisera reagents classified as contaminant or using the Diluent 1.



Rfc ajc_lgle qmjsrgmI QCRSN AJC?L gq b_l ecpsq,

Pcdcp dmp k mpc gl dmp k_rgmI _`msr rfc npmb sar rm8

² rfc GlrpsargmI q md Sqc &PCD ..77.0'9

² rfc K_rcpg_j Q_dcrw B_r_ Qf ccr,

- Exclusively use the SETUP CLEAN product (REF 009902) for cleaning solution preparation.

Prepared cleaning solution should only be filled into its foreseen 15 ml empty vials (REF 009310). Then, the vials have to be loaded in a reagent rack.

009310	SETUP CLEAN vials 15 ml (x12)
--------	-------------------------------

The SETUP CLEAN solution has to be changed daily. The same vial should not be used longer than 1 month. Therefore REF 009310 contains 12 vials (one vial for every month of the year).

- If a used vial containing SETUP CLEAN solution is outside of the IH-1000, it must be closed with the foreseen cap delivered with the package.

4,0,1,0 Npcn_p_rgm l md /3 k j md QCRSN AJC?L

K_rcpg_j pcosgpcb8

- SETUP CLEAN product with REF 009902;
- SETUP CLEAN empty vials of 15 ml (REF with barcode);
- Demineralized water.

U fcl sqgle bckgl cp_jgxcb u_rcp* k_ic qspc rf_r rfc u_rcp os_jgrw amppcqnm l bq_r jc_qr
rm Ep_b 0 &QM 1474' mp Rwnc &&?QRK B//71+7/',

Bm lmr dgjj k mpc rf_l /3 k j md rfc QCRSN AJC?L qmjsrgm l gl m lc tq_j PCD . .71/ . ,

A Add 3 ml of SETUP CLEAN REF 009902 in a vial.

B Add 12 ml of demineralized water.

The final concentration of this SETUP CLEAN REF 009902 is 20%.

0 This once prepared solution has to be preserved at ambient temperature.
If SETUP CLEAN solution is stored in a dark, light-protected bottle, it must be replaced once a month. If SETUP CLEAN solution is not light-protected (if it is in a transparent bottle), it must be replaced every day.

4.0.2 Q_k njcq - Pc_ec l r q ?pc_ Qapcc l

This screen displays a detailed status of the samples, diluents and reagents loaded in the IH-1000.

?aacqq

A




From the main screen, click on the **SAMPLE / REAGENTS AREA** of the interactive image



From any screen where the **REAGENT / DILUENT AREA** button appears.

4.0.2, / Bcqapgnrgm l



Dyespc 21, Sample / Reagent area screen

- 1 Rack Status (see [page 78](#))
- 2 Sample Status (see [page 78](#))
- 3 Diluent Status (see [page 79](#))
- 4 Reagent Status (see [page 79](#))
- 5 Details Button (see [page 80](#))
- 6 Remove Rack Button
- 7 Remove All Button (see [page 88](#))

4,0,2,0

P_ai Qr_rsq

Qw k `mj	Bcqapnrgm l	U f_r rm bm =
	The rack is surplus.	A rack is already loaded at this place in the pipetting area. Remove the supernumerary rack from the loading area.
	The rack is not identified, or the barcode of the rack is not valid, or one of the tubes is not identified, or the diluent rack is out of date.	Remove the rack and check that the barcode sticker is not damaged. For a diluent rack, check the validity date.
	A sample tube or a reagent vial is in error or an un-loadable rack is detected (not fully pushed in the instrument).	Check the samples tubes or reagents vials errors then remove the rack. Push the rack deeper in the IH-1000.
	The reagent rack has a stirring fault.	Remove and replace the rack.
	The rack is in the restitution position (only for particular cases).	Remove the rack from the instrument.
	The restitution is in progress or the restitution request is memorized.	Wait for process completion.
	Rack background color: gray The rack is in the pipetting position.	Can be selected for further information/actions.
	Rack background color: orange The rack is selected in order to view details or to remove it.	Click DETAIL or REMOVE button
	Rack background color: yellow The rack is out of service.	Remove the rack

4,0,2,1

Q_k njc Qr_rsq

Qw k `mj	Bcqapnrgm l	U f_r rm bm =
	No tube is detected at this position of the rack.	-
	The sample is correctly detected.	-
	Sample correctly detected (donor).	-
	Sample correctly detected (receiver).	-
	A tube is detected with a barcode reading error or duplication.	If manual input option is enabled, remove the rack and check the barcode of the defective tube. In case of duplicate, check both tubes.
	A tube is detected with a cap.	See chapter No Cap Confirmation on page 84 .
	Lack of receiver for sample tube.	Check the loaded crossmatch receivers.
	Lack of donor for sample tube.	Check the loaded crossmatch donors.

Click on the error of a sample rack to display a warning screen describing the nature of the error.

If the sample rack grid is in error, it is possible to display the error details by clicking on the representative error icon.

See chapter [Disabled Mode Option on page 137](#).

0 For particular cases with **Low volume segment blood** (rack type 4), error messages are displayed in the log list.
See chapter [Traceability Area on page 39](#).

4,0,2,2

Bqjsc l r Qr_rsq

Qw k `mj	Bcqapgnrgm l	U f_r rm bm =
	The diluent rack is correctly detected.	-
	The diluent rack is detected with an error.	Remove the rack and replace the diluent product.

Click on the error of a diluent rack to display a warning screen describing the nature of the error.

If the diluent rack grid is in error, it is possible to display the error details by clicking on the representative error icon.

See chapter [Disabled Mode Option on page 137](#).

4,0,2,3

Pc_ec l r Qr_rsq

Qw k `mj	Bcqapgnrgm l	U f_r rm bm =
	The vial is correctly detected.	-
	The vial is detected with an error.	Remove the rack and check that the vial is correctly inserted.
	The vial is detected and is reserved.	-
	The vial is detected with a cap error.	Remove the rack and remove the cap.
	The vial is detected with a barcode reading error.	Remove the rack and check the barcode of the vial: orientation, validity date, barcode.
	The vial is detected with a stirring error.	Remove and replace the rack.

Click on the error of a reagent rack to display a warning screen describing the nature of the error.

Remove or reload the rack to initialize it.

If the rack grid is still in error, it is possible to disable it by clicking on the error icon (with sufficient rights).

See chapter [Disabled Mode Option on page 137](#).

4,0,2,4

Bcr_gjq @srrm l

This button gives access to precise information about the selected rack.


A Click on the rack image to select it (orange background).

B Click on **DETAILS**.

Rfc bcr_gjq qapcc l bgyddcp u fcrfcp _ q_ knjc pc_ec l r mp bjsc l r p_a iq q qjc arc b,*

0 This button is activated only when the rack is in the pipetting position.


Q_k njc - Bcr_gjq Qapcc l



Dyespc 22, Sample details screen

- 1** Position of the rack
- 2** Name (barcode) of the rack
- 3** Remarks/Errors (if any)
- 4** Position of the sample
- 5** Barcode of the sample
- 6** Remarks/Errors of the sample (if any)
- 7** **MANUAL INPUT** button (see [page 83](#))
- 8** **REMOVE RACK** button (see [page 88](#))
- 9** **SWITCH ON / SWITCH OFF** button
- 10** **CAP CORRECTION** button (see [page 84](#))


Pc_ec Ir - Bcr_qlq Qapcc l



Dyespc 23, Reagent details screen

- 1 Position of the rack
- 2 Name (barcode) of the rack
- 3 Remarks/Errors (if any)
- 4 Position of the content in the rack
- 5 Barcode of the content
- 6 Reagent name
- 7 Reagent usable volume (ml)
- 8 Reagent lot number
- 9 Expiry date
- 10 Remarks/Errors of the content (if any)
- 11 **REMOVE RACK** button (see [page 88](#))
- 12 **SWITCH ON / SWITCH OFF** button
- 13 **CAP CORRECTION** button (see [page 84](#))

Bijsc Ir - Bcr_ajq Qapcc l



Djespc 24, Diluent details screen

- 1 Position of the rack
- 2 Name (barcode) of the rack
- 3 Diluent type
- 4 Remarks/Errors (if any)
- 5 Lot number
- 6 Expiry date

0 The contents of diluent racks are not displayed.

- 7 **REMOVE RACK** button (see [page 88](#))

4,0,2,5 Q_k njc @_pambc K_ls_j Glnsr

The **MANUAL INPUT** allows to enter manually a barcode for the samples detected with a barcode error.

0

This button is activated only when a sample is detected with a barcode error.

0

If the barcode is read manually, make sure to put the sample back to the original position.

A



From the main screen, click on the **SAMPLE / REAGENTS AREA**.

B

Select the rack with the sample barcode error.

C

Click on **DETAILS**.

D



Click on **MANUAL INPUT**.
Rfc dmjjm ugle qapcc l _nnc_pq8

SAMPLE / MANUEL INPUT

Emergency

Position	Barcode	
Grid 15	01	
Position	Barcode	Confirm barcode
1	4479	4479
2	4477	4477
3	4476	4476
4	4473	4473
5	...	;
6	4475	4475
7	4472	4472
8	4474	4474
9	4478	4478
10	4480	4480

OK

CANCEL

Dyespc 25, Sample manual input

E

Click on **REMOVE RACK**.

F

Scan or type the sample barcode in the **Barcode** (1) column.

G

Scan or type the sample barcode in the **Confirm Barcode** (2) column.

0

Copy and Paste functionalities are disabled during this operation.


H	Click on OK to validate the input.
I	<p>If the positions of some samples have been switched, an alert screen is displayed.</p> <p><i>Npmacbspc uijj f_tc rm `c pcnc_rcb,</i></p>
J	Place the sample rack back in the instrument (same position). Wait until the instrument has inserted the rack completely.
K	Enter your login and password.
	<p>Dgespc 26, Authentication Required screen</p>
L	Click on the OK button of the AUTHENTICATION REQUIRED screen.

4,0,2,6

Lm A_n Am l dgp k_rgm l

The **NO CAP CONFIRMATION** screen allows to confirm the absence of caps when necessary.

A	 <p>From the SAMPLE or REAGENT / DETAILS screen, click on the CAP CORRECTION button.</p> <p><i>Qcc af_nrcp Q_k njc - Bcr_yjq Qapcc l m l n_ec 6.,</i> <i>Rfc dmjjm uyle qapcc l _nnc_pq,</i></p>
O	This button is activated only when a cap is detected.



Djespc 27, No Cap confirmation screen

- 1 Symbol: Tube is detected with a cap
- 2 Check box activated
- 3 **OK** button
- 4 **CANCEL** button (see [page 83](#))


To confirm the absence of a cap after checking visually:

A Click the **check box** (2) next to the symbol; **Tube is detected with a cap** (1).

0 Only tubes with a cap are active. When a check box is unchecked, the tube returns to its initial state.

B Click **OK** to confirm.

Rfc dmjjmugle qapccI _nnc_pq8



Djespc 3., No Cap confirmation warning screen

C Click **YES** to confirm.


Rfc amppcqnmIbgIe qr_rc md rfc rs`c qg snb_rcb,

4,0,3 Jm_bgle P_aig

The loading procedure is identical for all kinds of racks.

A rack can be loaded at any time in a free position, even when the instrument is performing tests.

A Prepare a rack with samples, diluents or reagents.




Diespc 31, Racks preparation

0 Read carefully the instructions given for each type of rack. See:

- [Sample Racks on page 64](#);
- [Diluent Rack and Diluent Product on page 72](#);
- [Reagent Racks on page 74](#).

B Check where to load the rack:

- A green light (1) indicates that the slot is available.
- A blinking red light indicates that a rack is moving: loading or unloading.
- A red light (3) indicates that the slot is not available. A rack is already loaded.
- Sample and diluent racks can be loaded in positions 1 to 18.
- Reagent racks can be loaded in positions R1 to R4.




Diespc 30, Loading racks

C

Engage the rack in the slot and push it slowly until it locks automatically.

&O' qfmuq _ ucjj nj_acb p_ai,

? a_kcp_ ujj_srmk_rga_jjw gbc lrgdw rfc lcu p_ai _lb rfc rp_lqnmpr_pk ujj
r_ic gr glrm rfc glqrpskc lrdmp nqncrrgle mncp_rgm lq,



Dyesc 31, DANGER - Racks area

! ?r_lw rjkc* rfc p_ai ksqrj_w ml rfc `mrrmk md rfc qjmr gl mpbcp rm `c jmaicb_r rfc
amppcar nmqgrgm l,
Gd rfc p_ai gq nsqfcb rmm d_p* _ amjjgqgm l ugrf rfc rp_lqnmpr_pk amsjb maasp,

! Lctcp nsr_lwrfgle cjqc rf_l _p_ai glqgbc _qjmr,
Rfc a_kcp_ umsjb _srmk_rga_jjw rpw rm gbc lrgdw gr _lb _ amjjgqgm l umsjb maasp,

D

If the identification fails, the loading area becomes red.

**E**

The rack can be unloaded and identified manually using the manual barcode reader.

Qcc af_nrcp Q_knjc @_pambc K_ls_j lnsr ml n_ec 61,

4,0,4 Pckmtgle P_aiq

U_gr slrgj rfc jgefr &/' gq epcc l rm pckmtc rfc p_ai,

Rfc p_aiq am knjrcjw ncpdmp k cb am k c msr md rfc GF-/. . . _srm k_rga_jjw &O', @c _u_pc rm
l mr fgr rfck _aaqbc l r_jjw,

4,0,4,/ Pckmtgle Q_knjc P_aiq

A Select the **SAMPLE / REAGENT AREA** screen or the **SAMPLE / DETAILS** screen.

Qcc8

² [Q_knjcq - Pc_ec l r q ?pc_ Qapcc l m l n_ec 559](#)

² [Q_knjc - Bcr_yj q Qapcc l m l n_ec 6.](#),

B Select the rack or search the sample tube barcode to remove.

Dyesc 32, Search function

Gd l cacqq_pw g l rfc **Q?KNJC - PC?ECLR ?PC?** qapcc l,
? l w dms l b q_knjc rs`c gq fge f j e f rcb g l nspnjc,

C Press **REMOVE RACK**.

Gd l m q_knjc gq n p m a c q q c b * r f c r p _ l q n m p r _ p k n s q f c q r f c p _ a i m s r ,

Gd m l c m p k m p c q_knjcq _pc n p m a c q q c b * _ k c q q _ e c n p m k n r q r f c s q c p r m a m l d y p k
r f c p c k m t _ j ,



Dyesc 33,

D

Confirm the rack removal.

Rfc q lqpskclr dylqfcq rfc njncrryle mncp_rgm l ml rfc asppclr q_knjc, Rfc
mrfc q_knjc _pc lmr njncrrcb,
Rfc q_knjc p_ai q nsqfcb msr,
? kccq_ec bjqnj_wq rfc pck_q lgle rjk c dmp rfc sqcp rm pc+q lqcp rfc q_knjc
p_ai ml rfc q_kc nmqrgm l rm nrmaccb rfc rcqrq &qcc pck_pi q `cjm u',



Dyesc 34,

Pck_pi q8

If the unloaded sample rack is reloaded in the **q_kc** position **ugrfq l** 60 seconds:

- for samples already identified, the pipetting operation restarts;
- for missing samples, the tests which were scheduled are placed in **TESTS TO REPEAT** screen.
See [Tests to Repeat on page 114](#).

If the unloaded sample rack is reloaded in the **q_kc** position **_drp** 60 seconds or in other position, the rack is considered as new rack.

For any missing samples, the status of the scheduled tests is managed as follows:

Rcqr Qr_rsq	Bcqapnrgm l
TESTS IN PROGRESS	No impact.
Q_knjc am k njrcjw njncrrcb &Q_knjc rs`c l m kmpc lccbc b'	The test will be completed.
TESTS IN PROGRESS	The status of the test is set to TEST TO REPEAT .
Q_knjc lmr am k njrcjw njncrrcb &Q_knjc rs`c qrgjj lccbc b'	Reason: Missing resources (samples).
TESTS READY TO START	The test is canceled.
TESTS WITHOUT RESOURCES	The test is canceled.

4,0,4,0 Pckmtgle Pc_eclr mp Bjjsclr P_aiq

0 A reagent rack can be removed only if no reagent vial of the rack is in use.

A Select the **SAMPLE / REAGENT AREA** screen or the **REAGENT / DETAILS** or **DILUENT DETAILS** screen.

Qcc8

² Q_knjcq - Pc_eclrq ?pc_ Qapcc l m l n_ec 559


² Pc_eclr - Bcr_yjq Qapcc l m l n_ec 6/9

² Bjjsclr - Bcr_yjq Qapcc l m l n_ec 60,

B If necessary (in the **SAMPLE / REAGENT AREA** screen), click on the rack image to select it (orange background).

C Press **REMOVE RACK**.

Rfc rp_lqnmpr_p k nsqfcq rfc p_ai msr,



Dyespc 35, Remove racks - light

4,0,4,1 Pckmtc _jj P_aiq

A Click on **REMOVE ALL**.



Rfc rp_lqnmpr_p k nsqfcq _jj p_aiq msr* _jj rfc p_aiq _pc sljm_bcb,

0 It is only possible to remove racks which are not in use.

4,0,4,2 Pckmtgle Bgjsclr_l b Q_knjc P_ai k_ls_jjw

If a rack is blocked inside the loading area in pipetting position, manually remove it according to the following procedure.

SHUTDOWN rfc GF- / . . . `cdmpc pc kmtgle rfc p_aiq k_ls_jjw,

A SHUTDOWN the IH-1000.

Qcc af_nrcp Qfsrbmul md rfc glqpskclr ml n_ec 27,

B Take the rack extractor tool (1).


Gr gq jma_rcb gl rfc pgefrjmucp n_pr md rfc glqpskclr* lcvr rm rfc ecj a_pbq u_qrc
`gl,

C Slide the rack extractor tool under the rack.

D Insert the extremity of the tool in the dedicated hole of the rack (2).

E Lift smoothly the rack and pull it out.


F Put the rack extractor tool back in place.



Dyespc 36, Remove racks - pull out

4,0,4,3 Pckmtgle Pc_eclr P_aiq k_Is_jjw

' @c a_pcdsj lmr rm b_k_ec rfc cjcarpg a_j amlr_arq md rfc pc_eclr p_ai,	
A	Raise the back of the rack about 35 mm.
B	Remove the rack about 15 mm.
C	Rest the back of the rack.
D	Remove completely the rack.




Dyespc 37, Remove racks - Reagent

6.3 Managing Liquids and Wastes

This chapter explains how to manage cleaning liquid, system liquid, liquid waste and gel cards waste.

Аjc_lgle_ib qwqrc k jgosgb amlr_g lcpq pcnj_ack clr ksqr `c bmlc a_pcdsjw_aampbgle rm rfc dmjjmugle npmacbspcq gl mpbcp rm_tmgb nmjjsrgml md rfc fwbp_sjga agpasgr,

4,1,/ Bcqapgnrgm l



Dъesc 4., Liquids and wastes - description

- 1 Cleaning liquid container (black)
- 2 System liquid containers (blue)
- 3 Liquid waste containers (red)
- 4 Containers drawer
- 5 Gel cards waste bin
- 6 Gel cards waste slide
- 7 Gel cards waste bag

- A Push the door to open it.
- B Pull the containers drawer (4) to easily access the liquid containers (1), (2) and (3).
- O The cleaning liquid container (1) is not connected during normal operations; it is used only for maintenance operations.
- C After being processed, used gel cards will fall from the slide (6) into the waste bag (7).

Lctcp nsr f_l bq gl rfc ecj a_pbq u_qrc qjgbc u fgjc rfc gl qrs k clr gq nmucpcb ml,

Gd rfc jgosgb u_qrc gq aml lcarb rm_j_`mp\_rmpw bp\_gl\* tpgdw rfc os\_jgrw md rfc rs`gle_ib k_ic qspc rf_r gr gq amppcarjw aml lcarb,

4,1,0 Afcaig le Jgosgbq _lb U_qrcq Qr_rsq

The **SOLUTIONS / WASTES AREA** screen allows to check the liquids and wastes status.


?aacqq

A

From the main screen, click on the **LIQUIDS / WASTES AREA** of the interactive image.



Bcqapgnrnm l



Dyespc 4/, Solution / wastes area

- 1 Cleaning liquid container area
- 2 System liquid container area
- 3 Liquid waste container area
- 4 Gel card waste bin area
- 5 Maximum gel cards waste capacity of the bin

0 If all items in Solutions / Wastes Area are in order, it is not possible to perform an action.

4,1,1 Pcdgjjg le _ Qwqrck Jgsgb AmIr_g lcp



QmjsrgmI BcamI 7. gq b_l ecpsq,
Pcdcp rm rfc K_rcpg_j Q_dcrw B_r_ Qfccc dmp k mpc g l dmp k _rgmI ,

▪ Sqc ejmtcq dmp _jj mncp _rgmI q,

▪ U fcl sqgle bckglcp_jgxcb u_rcp* k_ic qspc rf_r rfc u_rcp os_jgrw amppcqnmI bq _r jc_qr
rm Ep_b 0 &QM 1474' mp Rwnc GG &?QRK B//71+7I',

▪ R_ic a_pc rm amI lcar rfc amppcar amIr_g lcp rm rfc amppcar fmqc, Pcdcp rm rfc amjmp ambc,


0 The system liquid has to be renewed weekly.

0 Do not mix different solutions of system liquid.

A Check that the second system liquid container is connected to the instrument while it is running.

B Pull the container out of the instrument.

C Remove the front cap (1) of the container.



Despc 40, Refilling system liquid container

D Clean the container with demineralized water.

E Check that the container is empty.

F Shake the Decon 90 bottle to mix the solution.


G Fill the container with system liquid:

3J md bckglcp_jgxcb u_rcp _Ib /3kJ md BcamI 7.,

H Put the container cap back in place.

I Put the container back in place in the instrument.

J Click **CONTINUE** to confirm the message, or **CANCEL**.



Djespc 41, Empty container message

- Гд rfc \$Cknrw amlr_glcpμ k cqq_ec _nnc_pq _e_gl _drp rfc amlr_glcp f_q `ccI pcdgjjcb\* BM LMR \_ai lmujcbec rfc k cqq\_ec \_e\_gl \* rfc fwbp\_sjga agpasgr umsjb `c npg k cb _e_gl _l b rfgq umsjb pcqsjr gl mtcpdjmu gle rfc djsv am knclq_rmp, Amlr_ar _qcptgac cle glccp,

4,1,2 Ajc_lgle Jgosgb Npcn_p_rgm l




Qmjsrgm l Kgapmagbc QO × gq b_l ecprmsq,
Pcdcp rm rfc K_rcpg_j Q_dcrw B_r_ Qf ccr dmp k mpc gl dmp k _rgm l,

- U fcl sqgle bck glcp_jgxcb u_rcp* k_ic qspc rf_r rfc u_rcp os_jgrw amppcqnm l bq _r jc_qr rm Ep_b 0 &QM 1474' mp Rwnc GG &?QRK B//71+71',

- Sqc ejmtcq dmp _jj mncp_rgm l q,

A Remove the front cap of the container.



Djespc 42, Refilling cleaning liquid container

B	Fill the container with cleaning liquid: <ul style="list-style-type: none">• 2l of demineralized water and 8ml of Microcide SQ™. <i>R f q q m j s r q m l q q q r _ ` j c d m p 5 b _ w q ,</i>
----------	---


4,1,3 **Aj c _ l g l e m d r f c A m l r _ g l c p q**

It is necessary to clean a container:

- once a month;
- when a container has not been used for a long time.

I	<i>U f c l s q g l e b c k g l c p _ j g x c b u _ r c p * k _ i c q s p c r f _ r r f c u _ r c p o s _ j g r w a m p p c q n m l b q _ r j c _ q r m E p _ b 0 & G Q M 1474' m p R w n c G G & ? Q R K B / / 71+7 / ' ,</i>
----------	---

A	Wait until the ongoing process is finished.
----------	---

B	Take the container out and remove caps and hose.  D g e s p c 43,
----------	--

C	Empty completely the container.
----------	---------------------------------

D	Rinse it with about 1L of cleaning liquid.
----------	--

E	Shake the container.
----------	----------------------

F	Remove and dispose the cleaning liquid.
----------	---

G	Rinse 4x or more the container with demineralized water.
----------	--

H	Wipe the caps and aspiration hose with Isopropanol, then rinse with demineralized water.
----------	--

4,1,4 Pcnj_age _ Jgosgb U_qrc AmIr_glcp

Jgosgb u_qrc amIr_glcpq amIr_gl amIr_kgl_rcb jgosgbq, F_lbjgle pcosgpc qncaq_j a_pc, Rm_tmgb amIr_kgl_rcb jgosgb qngjjq* jgosgb u_qrc amIr_glcpq ksqrlmr `c bqqamllcarb bspgle sqc,

Sqc ejmtcq dmp _jj mncp_rgmIq,

It is strongly recommended to use the antifoaming product each time the waste container has to be emptied.

4,1,4,/ Npcn_p_rgmI md rfc ?lrgdm_kgle QmjsrgmI

Various antifoaming products can be found worldwide from different manufacturers. Bio-Rad recommends antifoaming product that contains 30% of Silicone.



Rfc qggamlc _lrgdm_kgle _ecIrqq b_lcpmsq, Pcdcp rm rfc K_rcpg_j Q_dcrw B_r_ Qfccc dmp npca_srgmIq rm `cr_icl ufcl umpigle ugrf rfgq rwnc md npmbsar,

For 100 ml of diluted silicone antifoaming agent:

A Add 10 ml of silicone anti-foam agent in a bottle.

B Add 90 ml of demineralized water.

4,1,4,0 Pcnj_age rfc Jgosgb U_qrc

The following procedure explains what to do when a liquid waste container is full.

A Check that the second liquid waste container is connected to the instrument while it is running.

B Take the container.

C Remove the plunger.




Figure 44, Liquid waste container

O

As soon as the aspiration cap is removed, liquid wastes are automatically routed to the second liquid waste container.

D

Dispose the liquid waste in accordance with your local regulations.

E

Add 5 ml of antifoaming solution in the **c k nrw** liquid waste container.

F

Put the container back in place with the caps.

I

R_ic a_pc rm am l l car rfc amppcar am l r_g l cp rm rfc amppcar fmqc, Pcdcp rm rfc amjnp ambc,

G

Click **CONTINUE** to confirm the message, or **CANCEL**.



D

espc 45, Waste full message

4,1,5

Npg k g l e Awajc

The priming cycle is made automatically or with the general initialization function.
See chapter [General Initialization on page 167](#).

O

Make sure to have enough liquid.


4,1,6

Ajc_l g l e Awajc

See chapter [Weekly Maintenance on page 149](#).

4,1,7 Pcnj_agle rfc Ecj A_pbq U_qrc @g l


)	Rfc ecj a_pbq u_qrc `gl amlr_g lq amlr_kgl_rcb npmb sarq, F_lbjgle pcosgpcq qncaq_j a_pc,
!	Sqc ejmtcq dmp _jj mncp_rgm lq,
A	Prepare an empty waste bin with a bag.
B	When the instrument is ready or in idle state: <ul style="list-style-type: none"> • replace whole waste bin; • replace plastic bag in bin.
C	Close the waste bag.
D	Dispose the bag in accordance with your local regulations.



Dyespc 46, Replace waste bin

6.4 Checking Resources Status


From the main screen, two buttons in the footer strip allow to check the resources status.



4,2,/

Kqgqle Pcqm spacq

The **MISSING RESOURCES** button gives access to the following screen, showing all the missing resources.



Dyespc 47, Missing resources

- 1 Gel cards storage drawer
- 2 Reagent / Diluent area

0 When the missing diluent concerns both the right and the left side, the number of the wells missing is displayed for each side.

- 3 Liquids / Wastes area

4,2,0 Pcqm spacq m l @m _pb

The **RESOURCES ON BOARD** button gives access to the following screen, showing all the resources already loaded into the instrument (Gel cards, reagents, diluents) and the status of the **LIQUIDS / WASTES** area.




Fig 5.1. Resources on board screen

- 1 Gel cards storage drawer

Qcc af_nrcp Ecj A_pbq Qmp_ec Bp_ucp ?pc_m l n_ec / .1,

- 2 Reagent / Diluent area

Qcc af_nrcp Pc_ec l r - Bjjscl r ?pc_m l n_ec / .1,

- 3 Liquids / Wastes area

Ajjaig l e q l rf q _pc_ e q tcq _aacq q m rfc JGOSBQ-U?QRCQ ?PC? qapcc l u f q af e q tcq bcr_gjcb q l dmp k _r q m l,


Qcc af_nrcp K_l_e q l e Jjosq bq _l b U_qrcq m l n_ec 71,

4,2,0,/ Ecj A_pbq Qrmp_ec Bp_ucp ?pc_

Clicking on this area gives access to the **GEL CARD STORAGE DRAWER** screen which gives detailed information.

See chapter [Gel Cards Storage Drawer Screen on page 59](#).

- 1 Number of gel cards
- 2 Code of the gel card type
- 3 Gel card type




Djespc 5/, Gel cards information

4,2,0,0 Pc_ecIr - Bjjsclr ?pc_

Clicking on this area gives access to the **SAMPLES / REAGENT AREA** screen which contains more detailed information.

See chapter [Samples / Reagents Area Screen on page 77](#).

- 4 Reagent usable volume (ml)
Sq_`jc tmjs kc ; kc_qspcb tmjs kc + bc_b tmjs kc
- 5 Reagent code
- 6 Reagent type
- 7 Number of available diluent wells
- 8 Diluent code
- 9 Diluent type, with indication of side and quantity (left/right or both)



Djespc 50, Reagent and Diluent information

N?EC GLRCLRGML?JJW JCDR @J?L I

7 Managing Tests


Автоматизация тестирования

This chapter presents how to manage the tests.

7.1 Performing a Test

The following chart shows a standard process flow for performing tests.

0 Before performing a crossmatch test, make sure that there are enough gel cards loaded on the IH-1000.



7.2 Performing an Urgent Test

To perform an urgent test the same process applies. Use urgent racks so that the sample is processed before all other samples.

See chapter [Sample Racks on page 64](#).

7.3 TESTS Buttons of the Main Screen

From the main screen, the following buttons allow to access to several functions about tests.

Buttons are disabled if their function does not apply to any test.

TESTS IN PROGRESS

See chapter [Tests in Progress on page 108](#).

TESTS READY TO START

See chapter [Tests ready to start on page 111](#).

TESTS WITHOUT RESOURCES

See chapter [Tests without Resources on page 113](#).

TESTS TO REPEAT

See chapter [Tests to Repeat on page 114](#).

SAMPLES WITHOUT REQUEST

See chapter [Samples without Request on page 117](#).

7.4 Tests in Progress

This screen displays all the on-going tests in the instrument.

A From the main screen, click on **TESTS IN PROGRESS**.




Figure 51. TESTS IN PROGRESS screen


- 1 Samples special statuses
Emergency sample
Quality control sample
- 2 Sample barcode number
- 3 Test description
- 4 Test start and end time
@j_aid pc_j njkc u fcl rfc rcqr qr_prcb _lb c l bcb,
@jsc8 cqr nj_rcb njkc u fcl rfc rcqr qfmsjb qr_pr _lb d q l q q f,
- 5 Grid
Nmqyrm l md rfc p_ai q l rfc jm_b q l e _pc_
- 6 Position
Nmqyrm l md rfc q_k njc q l rfc p_ai

- 7** State
Nqncrrq le q las`_rql* ac lrpqse_rql* pc_bql e mp am k njrcrb,*
- 8** Sorting strip
Ajjai ml rfc amjsk l l_kc rm qmpr rfc b_r_q
- 9** **DETAILS** button
Rfq` srrml egtcq _aacqq rm rfc bcr_gjq md rfc qcjcarch q_k njc rcqr
Qcc af_nrcp [Bcr_gjq md Rcqr q l Npmepcqq ml n_ec //.](#),

5,2,/ Bcr_gjq md _ Rcqr g l Npmeqcq

A From the **TESTS IN PROGRESS** screen, click on **DETAILS**.

The following screen shows all the ongoing tests for one sample.
Each test uses one or more types of gel cards and can be in a different status. A different color is used for each state.



Dyespc 52, SAMPLE / DETAILS screen


- 1 Samples special statuses
 - Emergency sample
 - Quality control sample
- 2 Test in progress special statuses
 - To repeat
 - Without resources
 - Ready to start
- 3 **BACK** button
- 4 **PREVIOUS** / **NEXT** buttons
- 5 Estimated start time for scheduled test

7.5 Tests ready to start

This screen displays the tests that are ready to start: all the required resources are available in the instrument. They will be processed as soon as possible, depending on the availability of the pipettors, incubator and centrifuges.

If necessary resources are present and if no emergency rack is loaded, the processing order depends on the loading order. The first sample loaded will be processed first.

A From the main screen, click on **TESTS READY TO START**.




D\espс 53, TESTS READY TO START screen

- 1 Information field
Qcc af_nrcp Rcqrq yI Npмepcq mI n_ec /.6,
- 2 **EMERGENCY** button
Qcc af_nrcp Ck cpec law mI n_ec //0,

5,3,/ **Скоро**

By clicking on the **EMERGENCY** button, the status of the sample changes (1).
See chapter [Tests in Progress on page 108](#).

0 The status can only be changed if there are tests ready to start.



The screenshot shows a table with test results. A green callout bubble with two dots points to the 'EMERGENCY' button (represented by a star icon) in the first row of the table.


★	4787000000000	PR16D DC Screening I : IgG-IgA-IgM-C3c-C3d-ctl...
	4788000000000	PR16A DC Screening II : IgG-C3d-ctl (5056)
★	4788000000000	PR16D DC Screening I : IgG-IgA-IgM-C3c-C3d-ctl...
Q	4789000000000	MO33A Bloodgr. + Rev.gr:A-B-AB-D/A1-B+C-c...
	4790000000000	MO33A Bloodgr. + Rev.gr:A-B-AB-D/A1-B+C-c...
	4791000000000	MO33A Bloodgr. + Rev.gr:A-B-AB-D/A1-B+C-c...

Рис. 54. Emergency status

7.6 Tests without Resources

This screen displays the tests that can not be performed because of missing resources.

A From the main screen, click on **TESTS WITHOUT RESOURCES**.



Dyspc 55, TESTS WITHOUT RESOURCES screen

1 Information field

Qcc af_nrcp Rcqrg / Npmpcq m l n_ec / .6,

2 DETAILS button

*Rfjq `srm l ehtcq _aacq rm rfc jqr md kjqql e pcqmspacq dmp rfc qcjcarb rcqr,
Qcc af_nrcp Afcaile Pcqmspacq Qr_rs m l n_ec / . /,*

3 MISSING RESOURCES button

*Rfjq `srm l ehtcq _aacq rm rfc jqr md kjqql e pcqmspacq dmp _jj pcoscqrcb rcqrq,
Qcc af_nrcp Kjqql e Pcqmspacq m l n_ec / . /,*

0 It is also possible to display the **MISSING RESOURCES / DETAILS** screen available for one particular test.


7.7 Tests to Repeat

Tests can encounter errors during process. If so, they are moved to the following **TESTS TO REPEAT** list that informs in which step the instrument was when the error occurred.

0 To be able to repeat failed tests directly from this screen, activate the option **Allow failed tests to be repeated** in the **Option** screen.
Refer also to IH-Com Data Management User Manual.

A From the main screen, click on **TESTS TO REPEAT**.

*Rf;q qapcc l q qnjqr y lrm rum k_l _pc_q8 rfc am k njcrc jqr md rcqrq y l cppmp &/'
_lb rfc pcj_rcb k;qqgle q_ k njc q &O',*



Dyespc 56, TESTS TO REPEAT screen

- 1 TESTS TO REPEAT:** Information field
Qcc af_nrcp Rcqrq y l Npmpcq q m l n_ec / .6,
- 2 MISSING SAMPLES,** for each missing sample, the lower list provides the sample barcode.
- 3 ACTION** button with 2 options: **DELETE** and **REPEAT**.
- 4 SELECT ALL** button
- 5 UNSELECT ALL** button

If the option **Allow failed tests to be repeated** is set, the **REPEAT** function is enabled under the **ACTION** button, allowing the selected profiles to be executed again, after confirmation.

See [General Options on page 128](#).

Otherwise only the **DELETE** button is displayed, to acknowledge the error of all selected items, after confirmation.


If a test is present on the **TESTS TO REPEAT** screen, double click on the test to open a screen displaying explanation of the error.

No test can be repeated if no user is logged in to IH-Com.

A Select one or several test to be repeated.

B  Select **ACTION**.
Cd rfc mnrgm l ?jjmu d_gjcb rcqrq rm `c pcnc\_rcb gq cl\_`jcb,
Pcdcp rm rfc am k kclr_`mtc,

C  Select **REPEAT** to repeat the test.
? u_plgle kccq_ec bqqnj_wq gd q_ knjq_pc kqqqle dmp rfc qcjcrcb rcqrq,




Dyespc 57, Warning message

Gd Im sqcp gq jmeecb g l rm GF-Am k* m l jw rfc `srrm l **A?LACJ** gq \_t\_gj\_`jc,

D

Confirm to restart the test.

? u_p l q / e b q n j _ w q c t c l r s _ j k q q q l e q _ k n j c q ,



Dyespc 6.,

E

Load the required samples.

P c j _ r c b r c q r q q r _ r s q r s p l q r m p c _ b w r m q r _ p r ,


F

Confirm the loading of the samples.

0 Unloaded samples stay in the **TESTS TO REPEAT** area. It is possible to delete them.

7.8 Samples without Request

A From the main screen, click on **SAMPLES WITHOUT REQUEST**.



Dyspc 6 I, SAMPLES WITHOUT REQUEST screen

- 1 Sample list
- 2 Available tests list
- 3 Description of the tests
- 4 **CONFIRM** button
- 5 **SELECT ALL** button
- 6 **UNSELECT ALL** button

Rm_qqmag_rc_rcqr

No test can be associated to the sample if no user is logged in to IH-Com.


A Select a sample in the Samples list (1).

B Select a test in the Tests list (2).

C

Click **CONFIRM**.

? am l dyp k _rjm l k cq q _ec _nnc _pq,



Dyespc 60, Confirm test assignment

Gd l m sqcp yq jmeecb y l rm G F-Am k * m l jw rfc ` srrm l **A?LACJ** yq \_t \_vj ` jc,

D


Click **OK** to confirm the test.

Rfc q _k njc yq pck mtc b dpm k rfc j yqr,

7.9 Checking Pipetting Area Status

This screen allows to monitor the on-going actions in the pipetting area.

A  From the main screen, click on the **PIPETTING AREA** in the interactive image.



Dgespc 61, PIPETTING AREA screen

- 1 37°C Incubator area
- 2 Balancing gel cards area

Rm b9qnj_w ecj a_pbq sqcb dmp `_j lac mp ecj a_pbq n_prj_jjw sqcb,

1 N_prg_jjw sqcb ecj a_pbq k_sqr lmr `c pc kmtcb dpm k rfc g lqrps k c l r_l b pcyl qcprcb,

- ### 3 Information area

Rf_gq_gq_gr_gf_gc_gl_gd_gm_gp_gk_gr_gm_gl_g_p_gc_g_d_gm_gp_gecj_ga_gp_gb_gq_gs_gq_gc_gb_gd_gm_gp_g`_j_gl_ga_gc_gm_gp_gecj_ga_gp_gb_gq_gn_gp_gr_gj_gj_gw_gs_gq_gc_gb_g,
Rf_gc_g`\_j<sub>g</sub>l<sub>g</sub>a<sub>g</sub>c<sub>g</sub>ecj<sub>g</sub>a<sub>g</sub>p<sub>g</sub>b<sub>g</sub>q<sub>g</sub>u<sub>g</sub>u<sub>g</sub>j<sub>g</sub>`c_gl_gn_gm_gq_gy_gr_gm_gl_gq_g/._+/0_gm_gl_gj_gw_g_l_gb_g_b_gc_gd_gs_gj_gr_gr_gw_gn_gc_g*`_p_ga_gm_gb_gc_g_l_gb_gp_gc_gk_g_g_gl_gg_gl_ge_g
r_gg_gk_gc_gu_gu_gj_g`c_gb_gg_gq_gn_g_w_gc_gb_g,
Rf_gc_gg_gl_gq_gr_gp_gs_gk_gc_gl_gr_gu_gu_gj_gq_gr_gm_gp_gc_g_l_gb_gs_gq_gc_gn_gp_gr_gj_gj_gw_gs_gq_gc_gb_gecj_ga_gp_gb_gq_gd_gm_gp_gs_gn_gr_gm_g/3._gk_gg_gl_gs_gr_gc_gq_g_d_gr_gc_gp_gr_gf_gc_ga_gp_gb_gg_gq_g
n_gg_gc_gp_ga_gc_gb_g_l_gb_gr_gf_gc_gl_gb_gg_gq_ga_gp_gb_gr_gf_gc_gk_g,Rf_gc_gr_gw_gn_gc_gm_gd_gecj_ga_gp_gb_g*_t_gg_g_`j_gc_gu_gc_gi_gj_gq_g_l_gb_gr_gg_gk_gc_gp_gc_gk_g_g_gl_gg_gl_ge_gu_gu_gj_g`c_g
b_gg_gq_gn_g_w_gc_gb_g,

- 4** Pipetting area Left and Right
- 5** **SWITCH ON**
- 6** **SWITCH OFF**
- 7** **OPEN** (pipetting area window)
Qcc *MncIqle rfc Nqncrnqle ?pc_ Uq/bmu ml n_ec /0/*
- 8** Status left pipettor
- 9** Status right pipettor
- 10** Door state (displayed in red when open, in main screen the pipetting area is also displayed in red)

5,7,/ Mnc l g l e r f c N g n c r r g l e ? p c _ U g l b m u

- P c _ b r f c d m j j m u g l e g l q r p s a r g m l a m k n j c r c j w ` c d m p c m n c l g l e r f c n g n c r r g l e \_ p c \_ u g l b m u , M l + e m g l e r c q r q u g j j ` c g l t _ j g b ,
G l _ l w a _ q c * g r g q d m p ` g b b c l r m m n c l r f c n g n c r r g l e _ p c _ u g l b m u k _ l s _ j j w ,
- @ c a _ p c d s j u f c l _ a a c q q g l e g l q g b c r f c n g n c r r g l e _ p c _ * c q n c a g _ j j w u g r f l c c b j c q * G B + A _ p b q n g c p a g l e k m b s j c * f c _ r m d g l a s ` _ r m p , ,
- B m l m r f g r r f c l c c b j c ,
- ? j j p _ a i q k s q r ` c p c k m t c b ` c d m p c m n c l g l e r f c n g n c r r g l e _ p c _ u g l b m u * g d r f g q a m l b g r g m l g q l m r p c q n c a r c b r f c p c g q _ p g q i m d a m j j g q g m l u g r f r f c p _ a i q ,
@ g m j m e g a _ j p g q i g l a _ q c m d n p m h c a r g m l ,
- B s p g l e r f c g l g r g _ j g x _ r g m l * g d _ j j p _ a i q f _ t c l m r ` c c l p c k m t c b r f c p c g q _ p g q i r m b _ k _ e c r f c g l q r p s k c l r ,

A Remove all returned racks or racks in waiting position.

Q c c a f _ n r c p P c k m t c _ j j P _ a i q m l n _ e c 7 . ,

B Click on **OPEN** (pipetting area window).

*U f c l a j m a g l e r f c u g l b m u * r f c g l q r p s k c l r n c p d m p k q _ l g l g r g _ j g x _ r g m l _ l b q a _ l q r f c u m p i r _ ` j c _ e _ g l ,*

S q c b e c j a _ p b q _ p c n s r r m r f c u _ q r c ` g l ,


R f c p c o s c q r q r f _ r f _ t c ` c c l q r \_ p r c b ` c d m p c m n c l g l e r f c n g n c r r g l e _ p c _ u g l b m u f _ t c r m ` c q c l r _ e _ g l ,

0 When the pipetting area window is in error, the **OPEN** button is highlighted in red.



All racks removed - Correct!

Dyespc 62,




Racks still loaded - Not correct

5,7,0

Ajmr Bcrgm I Cppmpq


The IH-1000 knows three different clot detection errors:

- Clot detected;




Djespc 63, Clot detected error message

- Defusing during sampling;



Djespc 64, Defusing during sampling error message

- Distribution canceled.



Djespc 65, Distribution canceled error message

0 All the tests with error go back to the **TESTS TO REPEAT** screen. Go to that list and click on the single test to display the error message.
See chapter [Tests to Repeat on page 114](#).

5,7,0,/ Ajmr Bcrcargm l

When a clot detection of a needle occurs, the user must loose the needle from its obstruction.

?l cppmp k cqg_ec gq b qnj_wcb,

The resources present in the pipetting area undergo the following management:

- For the pipettor which the needle is clotted, the tests are stopped and set in “**Tests to repeat**”.
- For the other pipettor, which the needle is not clotted, no new pipetting batch is engaged. The current pipetting continues until the users action.

! @c a_pcdsj u fcl _aacqqg le g l qgbc rfc ngncrrg le _pc_* cqncaq_jjw ugrf lccbjcq* GB+A_pbq ngcpag le kmbsjc* fc_r md g las`_rmp,,

! G l _lw a_qc* gr gq dmp`gbbcl rm mncl rfc ngncrrg le _pc_ ug l bmu k_ls_jjw,

! Sqc ejmtcq dmp _jj mncp_rgm l q,
Bm lmr kmte p_aid* rs`cq* tg_jq* pc_eclr p_aid* ecj a_pbq,
@c a_pcdsj lmr rm kmte rfc ngncrrmp _p k g l rfc fmpgxm l r_j b gpcargm l,
Kmtc rfc lccbjc qsnnmpr _q g l bga_rcb m l rfc nfmrmq* mrfcp ugqc rfcpc gq _pgqi rm a_sqc b_k_ecq m l rfc GF-/... ,

! Sqc _bpw lml umtel ugnc rm ajc_l rfc lccbjc rm_tmgb b_k_ec m l rfc lccbjc am_rg le,

A

Click on the **PIPETTING AREA** of the main screen.

Rfc N gNCRR gLE ?PC? qapcc l _nnc_pq,

Qcc af_nrcp Afcaig le Ngncrrg le ?pc_ Qr_rsq m l n_ec //7,

B

Click **OPEN**.

Rfc g l qrp skcl r mncl q rfc ngncrrg le _pc_ ug l bmu* m lac rfc ngncrrg le _pc_ gq _aacqqg`jc & ngncrrg le gq dg l gqfcb',


Qcc af_nrcp Mnc l g le rfc Ngncrrg le ?pc_ Ug l bmu m l n_ec /0/,

Rfc qapcc l lcvr gq b qnj_wcb,



Dyespc 66, Clot detected - information

C	When the pipetting area window is open, lift slowly the support of the needle.
D	Visually verify the presence of a clot at the tip of the needle.
E	<p>If no clot is observed just leave the needle at the top of the support position.</p> <p><i>Gd _ ajmr bcrcarqm / k cq q _ ec am l n q / sc q rm _ nnc _ p ctc l q d r f c p c q q l m ajmr npcqc l r * a _ j j _ qc p t q ac c l e q / cc p</i></p> <p>If there is a clot, with a non woven wipe clean the tip of the needle and if possible remove manually the clot.</p>
F	<p>Close the pipetting area window.</p> <ul style="list-style-type: none">• All the present gel cards on the pipetting area are detected and thrown to the waste bin.• All the present resources are identified again.• The processing procedure can restart.



Dyespc 67,

8 Options

Af_nrcp Mtcptqcu

This chapter presents the available options to configure the instrument.

0 This configuration is allowed only for specific users; login and password are required.

8.1 Options Screen

The **OPTIONS** screen can be accessed from any screen where the **MENU** button appears.

A

Click on **MENU**.

B



Click on **OPTIONS**.

C

Type your login and password.

D

Click **OK**.

Dyespс 7., Options - Authentication screen



Dgespc 7I, Options screen

1 GENERAL options

Qcc af_nrcp *EcIcp_j MnrqmIq mI n_ec /06,*

2 PROFILE options

Qcc af_nrcp *Npmdjic MnrqmIq mI n_ec /10,*

3 CONTROLS options

Qcc af_nrcp *AmIrpmjq MnrqmIq mI n_ec /12,*

4 DISABLED MODE options

Qcc af_nrcp *Bjq_`jcb Kmbs MnrqmI mI n_ec /15,*

8.2 General Options

Despc 70,

Ajmqc ug l bmu _dr cp

This option defines the time (in seconds) before the software automatically comes back to the main screen. Select **Never** to disable.

0 The user can close and open again the software without any effect.

Qms l b rm nj_w m l c p p mp

This option defines the sound that will be played in case of an error.

Select **No Sound** to disable.

@_pam bc qr _pr d j j r c p

This option allows to specify the amount of digits to suppress at the beginning of the samples barcode.

See chapter [Barcode Exception Mask Filter Option on page 130](#).

@_pam bc c l b d j j r c p

This option allows to specify the amount of digits to suppress at the end of the samples barcode.

See chapter [Barcode Exception Mask Filter Option on page 130](#).

K_i c q s p c r f _r r f c d j j r c p g l e b m c q l m r p c k m t c b g e g r q s q c b `w G F - A m k _l b - m p J G Q d m p _ n p m n c p q _k n j c g b c l r g d g a _r g m l , M r f c p u g q c * G F - A m k k _w b c r c a r _ b g q a p c n _l a w ,

@_pambc CvacnrqmI K_qi&q' Dgjrqp

To define exceptions for digits suppression in **Barcode start filter** and **Barcode end filter** options.
See chapter [Barcode Exception Mask Filter Option on page 130](#).

Qfmu `srrmI B_r_ k_l_eckclr

To enable or disable the Data management button in the main screen (check to enable).

Sqc H_n_lcq @_pambcq

To inform the instrument that resources are with Japanese barcodes (check to activate).

A_n brcarqmI _arjt_rcb

To enable the cap detection function (check to activate).

It is possible to activate the samples cap detection during the racks identification.

O Input for samples is always enabled.

If this option is activated the system will detect caps, see [Cap Detection on page 139](#).

If this option is deactivated, no cap will be detected.

?jjmu d_gjcb rcqrq rm `c pcnc_rcb

This option enables the **REPEAT** button. If all resources are available, failed test can be repeated directly from the **TESTS TO REPEAT** screen (check to enable).

See chapter [Tests to Repeat on page 114](#).

01b acjj qsqnc lqgmI dmp ?@M rcqrq le

This option enforces the pipettors to realize one sample take (and eventually dilution) for each distribution concerning the ABO tests family (check to enable).

O This behavior is legally expected in some countries.
All bloodgroup assays are concerned.

6,0,/

@_pambc Cvacnrgm l K_qi Djrcp Mnrgm l


This option defines exceptions for barcode filtering set with **Barcode start filter** and **Barcode end filter**. See the examples on the next page.

r K_i c qspc rf_r rfc dgjrcpg l e bmcq l mr pc k mtc bgegrq sqcb `w GF-Am k _l b-mp JGQ dmp _ nprmnpc q_k njc gbc l r g d g a _rgm l , Mrfcp u g q c * GF-Am k k_w brcar _ bgqapcn_l aw,

0 This option is available only if **Barcode start filter** and/or **Barcode end filter** options are active (value different of **0**).

A In **Options** screen, select **Customize** to open the **Barcode Exception Mask(s) Filter** option.

Qcc @_pambc Cvacnrgm l K_qi & q' Djrcp m l n_ec /07,
Rfc dmjjm u g l e qapcc l _nnc_pq,



Dyespc 71. Mask exception list for barcode caesura

B Configure the table of exceptions:

- Maximum number of digit: 3
- Maximum number of rules: 6
- Characters allowed:
a b c d e f g h i j k l m n o p q r s t u v w x y z
A B C D E F G H I J K L M N O P Q R S T U V W X Y Z
0 1 2 3 4 5 6 7 8 9
Space ! " # \$ % & ' () * + , - . / : ; < = > ? @ [\] ^ { | } ~

C In **Barcode exception pattern** area, enter required characters and define if they should stay **At the beginning** or **At the end**.

D  Click on the arrow to add the rule to the list.

E Click on **OK** to validate the exception list.

Pck_piq

According to the example of exception table above (Figure 93), the sample barcodes are filtered, but the system does not suppress digits for the samples barcodes:

- starting with **D**, **R** and **LAB**;
- finishing with **QC**.


See the examples below.

Djrcpg le apgrcpq_

- **Barcode start filter** = 2
- **Barcode end filter** = 1
- Exception criteria = see table above

Q_k njc ` _pambcq `cdmpc djrcpg le	Q_k njc ` _pambc _drp djrcpg le
P00001	000
XP0002	000
DP0003	BP0003 &Cvacnrqm I* I mr djrcpcb'
ZP0004	000
RP0005QC	PP0005OA &Cvacnrqm I* I mr djrcpcb'
WP0006	000
LAB0009	J?@0009 &Cvacnrqm I* I mr djrcpcb'

8.3 Profile Options



Dyespc 72,

It is possible to define and to apply a default profile for not allocated samples and not allocated emergency samples.

0 Default profile does not apply to Low volume segment blood (racks of type 4 - Barcode 04 or 14). Samples are displayed in **SAMPLES WITHOUT REQUEST** screen.
See chapter [Using Low Volume Segment Blood on page 66](#).
For details about how to assign a profile, refer to the Data Management software User Manual (Test groups).

?nnjw bcd_sjr npmdgjc dmp lmr _jjma_rcb q_ k njcq

This option allows to choose whether the default profile must apply for not allocated samples or not.
If the default profile is not applied, a specific profile must be assigned for each one.

?nnjw bcd_sjr npmdgjc dmp lmr _jjma_rcb c k cpec law q_k njcq

This option allows to choose whether the default profile must apply for not allocated urgent samples or not.

If the default profile is not applied, a specific profile must be assigned for each one.

For details on how to identify a urgent tube, see [Emergency Sample Rack on page 69](#).

Apply default profile for not allocated samples gq qcr* _q_k njc jm_bcb
ugrf l m qncagdga pcoscqr dmp rcqr ugjj `c _srm k_rga_jjw rcqrcb ugrf rfc bcd_sjr npmdgjc*
mr fcpugqc _qncagdga npmdgjc k sqr `c _qqe lcb dmp c_af q_k njc,

Bcd_sjr npmdgjc dmp lmr _jjma_rcb q_k njcq

This options allows to choose which profile will be assigned to the tubes without test.

This option is active only if the check box **Apply default profile for not allocated samples** is checked.

Bcd_sjr npmdgjc dmp lmr _jjma_rcb c k cpec law q_k njcq


This option allows to choose which profile will be assigned to the not allocated urgent samples.

This option is active only if the check box **Apply default profile for not allocated samples** is checked.

For details on how to identify a urgent tube:

See chapter [Emergency Sample Rack on page 69](#).

8.4 Controls Options



Dyesc 73,

The **Controls** area allows to activate:

Ecj Am l rpmj

0 If this option is activated, the camera checks the gel of the card before the gel card is used (e.g. presence of air bubbles - Gel card unusable, it is sent to the waste bin).

Pcrspl md rfc a_pbq gl bp_u cp dmp qcaml b pc_bg le

Choose between:

Lctcp	No card is returned for second reading.
?ju_wq	All cards are returned for second reading.
MI jw gd slqncagdycb pcqsjr	Gel card with unspecified results. Unspecified results: ?, LIQ, dp, wF, wR, wP, wD, L, W, E See chapter Test, Interpretation and Result on page 31.
Asqrm kpxcb	Only selected results will return the gel cards. Press the button Customize to define the list (see below).


0 If the drawer is open or the gel cards trays are full during card reading, the cards will be thrown away automatically.

ASQRMK6XC

Choosing a customized return of the gel cards enables a customization screen.

- Select available results for given categories or use **All** check box to select or de-select all results.
- Click **OK** to confirm.

*Qcc af_nrcp Pcrsp l md rfc A_pbq g l Bp_ucp
dmp Qcam l b Pc_bg l e Mnrqm l m l n_ec /13,*



Dgespc 74, Customize screen

Bgqrpg`srqm l Am l rpmj &os_jgr_rjtc'

0 The option is permanently activated.

IH-1000 software records all information on tests processed. If a volume control error occurs, interpretation code and image of each well used are transferred to IH-Com Data Management software.

6,2,/ *Pcrsp l md rfc A_pbq g l Bp_ucp dmp Qcam l b Pc_bg l e Mnrqm l*


6,2,/,/ *Pck_g l g l e A_pbq Jma_rqm l q l rfc Bp_ucp*

If the **Return of the cards in drawer for second reading** option is activated, a counter on the main screen displays free card places in the drawer (only if less than 12).



Dgespc 75, Cards locations in the drawer

If only one gel card tray remains empty in the drawer, a warning message is displayed:




Дієспс 76,

8.5 Disabled Mode Option

These options allow to disable some of the instrument components.

A module should be disabled if too many errors occur on it.



Дѣспс 77,

A Click on a component button.

Rfc_srfcIna_rml qapccI _nnc_pq,

Дѣспс / ..., Disable component - Authentication screen

0 A disabled component has a grey button.

B Enter your login and password.

C Click on **OK**.

DEnter <YES> and click **OK** to confirm disabling of component.

D:\espc / . /, Disable component - Confirm screen

- @c a_pcdsj u f_r rm bq_`jc, Gd cqqlrq\_j n\_prq\_pc bq\_`jcb* gr qq lmr nmqq_`jc rm am lrg lsc rm mncp_rc rfc glqps kclr,

8.6 Cap Detection


Rfc GF-/. . . ksqr lmr `c cvnmqcb rm bypcar qsljgefr,
 Gd rfgq amlbgrm l gq lmr pcqncarcb* rfcpc gq _pgqi rm lmr pc_b rfc `_pambc* rm lmr brcar rfc
 a_n _lb rm ap_qf rfc lccbjc,

If the cap detection option is deactivated, no cap will be detected.

If this option is activated, only the valid functions are activated.

6,4,/ RwnC md Pc_ec l r Tq_jq ugrf Brcarc b A_n

All these caps are detected by IH-1000.



Dyespc / .0, Reagent cap detection

- Type A-B-C-D on 5 ml vial
- Type E on 10 ml vial (Other caps can not be introduced into the IH-1000)
- Type E on SETUP CLEAN vial


0 10 ml vials with a dropper can not be loaded.

6,4,0 Rwnс md Qr_ l b_pb Rs`сq ugrf Bсrcarcb A_n

These types of tubes must be charged only on rack with code 01 or 11, and in case of segments samples in rack 04 or 14.

Гд rfcqc am l bgrgm l q _pc l mr pcqncarcb* rfcpc qq _ pgqi rm ap_qf rfc lccbjc,

All these caps are detected by IH-1000.



Dyespc / .1, Cylindrical tubes cap detection

Pcd	@p_l b	K_l sd_arspcp Pcdcpclac &/'	Amjmp	Qgxc &0'
@B /	BD diagnostic	367704	Opaque blue	13*75
@B0	BD diagnostic	367862	Opaque purple	13*75
@B1	BD diagnostic	367614	Opaque red	13*75
Rcp /	Terumo	VF-054SBCS	Opaque bleu	13*75
Rcp0	Terumo	VF-054SDK	Opaque purple	13*75
Rcp1	Terumo	VF-054SAS	Opaque red	13*75
Ep /	Greiner	454322	Opaque light blue	13*75
Ep0	Greiner	454320	Opaque dark blue	13*75
Ep1	Greiner	454036	Opaque dark purple	13*75
Ep2	Greiner	454021	Opaque light purple	13*75
Ep3	Greiner	454027	Opaque red	13*75

- (1) All tubes with the same cap are detectable.
- (2) All tubes higher than 100mm with cap, must not be charged on the IH-1000.

6,4,0,/ Rwnс md Awjg l bpgaj Rs `cq A_nq lmr Bсrcarcb

The translucent colors: red, blue, purple are not detected.

All the other colors are detected.



Dgespc / .2, Cylindrical tubes cap detection

6,4,1 Jm u Tmjs k c Rs `c & dmp Pcb Acjjq m l jw' A_n Bсrcargm l

Not valid function (reliability less than 100%).

Low volume tube (for red cells only) must be used only on rack with code 02 or 12.

▪ Gd rfgq am l bgrgm l gq lmr pcqncarcb* rfcpc gq _pgqi rm ap_qf rfc lccbjc,

6,4,2 Njs l ecp A_n Bсrcargm l

Not valid function (reliability less than 100%).

These types of tubes must be loaded only on rack with code 03 or 13.

▪ Gd rfgq am l bgrgm l gq lmr pcqncarcb* rfcpc gq _pgqi rm ap_qf rfc lccbjc,

N?EC GLRCLRGML?JJW JCDR @J?L I

9 Maintenance


Аф_нрсп Мтсрtꝑcu

This chapter presents the maintenance operation that can be performed on the instrument.

9.1 Checking Components Status

Components status can be checked through the components area screen.

- A**  From the main screen, click on the **COMPONENTS AREA** of the interactive image.



Dyespc / 1.3, Components area screen

- 1 Instrument temperature
Rfc q l q r p s k c l r r c k n c p _ r s p c y q k c _ q s p c b y l r f c n y n c r n y l e _ p c _ ,
- 2 External temperature
Rfc c v r c p l _ j r c k n c p _ r s p c k s q r ` c i c n r ` c r u c c l / 3 B A _ l b 0 5 B A ,
- 3 Back doors status

1 *M l j w q c p t g a c c l e g l c c p q _ p c _ j j m u c b r m m n c l ` _ a i b m m p q ,*

4 Components status:

- **Ерсу ык_ес:** The component is functioning correctly. No user action is required / possible;
- **Рсб ык_ес:** The component is under error status.

Ajjai ml rfc ык_ес rm qfmu rfc _jcpr qapccI _lb _ai lmujbec rfc npm`jck,


5 Components images

Rfc amjmp md rfc ылqrps kclr bcncIb ml rfc amknmlclr qr_rsq cvnj_ылcb _`mtc,

9.2 Maintenance Screen

The maintenance screen allows to perform usual maintenance operations required for the instrument. The **MAINTENANCE** screen can be accessed from any screen where the **MENU** button appears.

- A** Click on **MENU**.
- B** Click on **MAINTENANCE**.
- C** Login with your user name and password to access the maintenance screen.



Djespс / .4, MAINTENANCE screen

- 1** See chapter [Weekly Maintenance](#) on page 149.
- 2** See chapter [Needle Change Right/left Pipettor](#) on page 164.
- 3** See chapter [Rinse Right/Left Pipettor](#) on page 166.
- 4** See chapter [Reading Station Calibration](#) on page 166.
- 5** See chapter [General Initialization](#) on page 167.

9.3 User Level Screen

The User level screen allows to define the minimum user level to access the maintenance actions of the instrument.

0 Only users with a level higher than 1 can access this screen.


The **USER LEVEL** screen can be accessed from any screen where the **MENU** button appears.

A Click on **MENU**.

B Click on **AUTHORIZED MAINTENANCE**.



C Login with your user name and password to access the maintenance screen.



Dyespc / .5, User Level screen

- 1 Weekly Maintenance**
- 2 Needle Change Pipettor (right and left)**
- 3 Rinse Pipettor (right and left)**
- 4 Reading Station Calibration**
- 5 General Initialization**

To enable a **Level 1 User** to perform **Weekly Maintenance**, **General Initialization** or any other maintenance operation, check the corresponding maintenance check box of the **User Level** screen.

If checked, users with level 1 will have access to the maintenance action, otherwise only users with a level higher than 1 will have access to the required maintenance.

0 If maintenance actions require a user level greater than 1, the buttons of the **MAINTENANCE** screen will be grayed for level 1 users.

9.4 Maintenance Operations

Log any maintenance operation in the maintenance record.

See chapter [Maintenance Record on page 182](#).

0 If any variation of performance of the instrument is observed or suspected even after maintenance operations, contact your service engineer.

7,2,/

Bcamlr_kgl_rgmI

0 It is user responsibility to do the appropriate decontamination in case of spillage of hazardous material on or inside of the instrument, with a lint-free cloth soaked in cleaning solution.

1 Qfsrbmulrfcglqrpskclr`cdmpcmnclele rfcngncrgle_pc_uglbmu,

In case of spilling and contamination in the pipetting area or in the drawer:

A	Open access to the contaminated area with the user interface software (pipetting area window or drawer). <i>Qcc Mnclele rfc Ngncrgle ?pc_ Uglbmu ml n_ec /0/ mp Mnclele - Ajmqle rfc Bp_ucp ml n_ec 4.,</i>
B	Shut down the instrument. <i>Qcc af_nrcp Qfsrbmul md rfc glqrpskclr ml n_ec 27,</i>
C	Check that there is no damage to the instrument. <i>cd rfcpc qg_lw b_k_ec* a_jj_ qcptqac cleqlccp_lb qmn sqle rfc glqrpskclr,</i>
D	Take out all possible contaminated parts (racks, gel cards, gel cards trays).
E	Clean the contaminated areas (on or inside of the instrument), samples racks and reagents racks with a lint-free cloth soaked in cleaning solution. <i>Sqc `arcpaqbc* tpaqbc_lb dsleaqbc bqql dcar_lr_q ajc_lgle qmjsrgml,</i>

7,2,0 Ncpqmbga K_glrcl_lac

7,2,0,/ B_gjw K_glrcl_lac

- Check daily the outer surface (is it clean, no blood, no serum...);
- If needed, clean outer surfaces of worktable with any bactericide, viricide and fungicide disinfectant after switching off the system.

7,2,0,0 Uccijw K_glrcl_lac

! Rfc uccijw k_glrcl_lac ksqr `c ncpdmp kcb ctcpw ucci ctc l gd rfc glqpskc l r gq lmr sqcb bspgle _jml e ncpqmb,

0 An orange maintenance reminder (**WEEKLY MAINTENANCE NECESSARY**) appears on the IH-1000 main screen if the user did not perform the weekly maintenance within the last seven days.

0 When a component is disabled, it is not possible to perform the weekly maintenance. Except if only one pipettor is disabled.

Check the expiration date of the solutions.

Prepare a new solution:

- 2l demineralized water, 8 ml Microcide.

Rfqg qmjsrjm l gq qr_`jc dmp 5 b_wq,

Check that there is no leakage around liquid containers.

0 Before performing the weekly maintenance, make sure that the cleaning liquid is filled with min. 2 liters.

0 When starting the instrument: the IH-1000 initializes automatically. It also checks the modules and the worktable automatically.

! @cdmpc uccijw k_glrcl_lac* cknw rfc rum jgsgb u_qrc amlr_glcpq_lb nsr rfck `_ai gl nj_ac,
Sqc ejmtcq dmp _jj mncp_rgm l q,

A From any screen, click on **MENU**.

B Click on **MAINTENANCE**.


Rfc jmeyl qarcc l _nnc_pq,

C Enter your login and password.

D Click on **OK**.

EClick on the button **WEEKLY MAINTENANCE**.

Rfc dmjjm ugle qapcc l mnc lq,



Dyespc 1.6, Weekly maintenance - Confirm screen

0

If **WEEKLY MAINTENANCE** requires a user level greater than 1, this button is grayed for level 1 users.

FClick on **OK** to confirm.

Rfc dmjjm ugle qapcc l mnc lq _lb _l _sbj`jc _j_pk _arjt_rcq8



Dyespc 1.7, Remove ID-Cards piercing array equipped

Gr gq nmqqj`jc rm bqj`jc rfc _sbj`jc \_j\_pk ugrf rfc Qmn \_j\_pk `srrm l,

0

The pipetting area window opens automatically.


G

Remove the ID-Cards piercing array equipped.

Qcc af_nrcp Pckmtyle rfc iB+A_pbq Nycpaile ?pp_w Cosynncb m l n_ec /37,


H

Close the pipetting area window.

Rfc dmjjm ugle qapcc l mnc l q8**Dyespc III.**, Weekly maintenance in progress**I**

Cleaning of ID-Cards piercing array equipped (2):


- The decontamination support (3) must be put in a tank.
- The cleaning liquid must be filled in the tank to the top of this side (1).
- Use the prepared solution.

Qcc af_nrcp Uccijw K_g lrc l lac ml n_ec /27,**Dyespc III.**, ID-Cards piercing array equipped - decontamination support

- *Rfc qmjsrgm l ksq r lmr pc_af rfc snncp qspd_ac &2' md rfc bcamlr_k g l _rgm l qsnmpr,*
Rfc ajc_l g l e jgsgb ksq r lmr ecr g l am l r_ar ugrf rfc snncp kcaf_l ga_j n_pr md rfc ngcpag l e
_pp_w,

- The piercing tool must be put on the support.


?drçp _ dçu kçlsrcq* rfc dmjjmuçle qapcc l mnc lq _ lb _ l _ sbç`jc _ j_p k _ arçt_rcq,



Dçespc / / 0, Replace containers

J

Open the door of the liquid container compartment (3).



Dçespc / / 1, Replace containers

K

Pull out the container drawer (4).

L

Remove the two caps of the cleaning liquid container (1).

M

Place the two hoses from the system liquid (2) into the cleaning liquid container (1).

Qcc af_nrcp K _ l _ eçle Jçosçbq _ lb U_qrcq ml n_ec 71,

NClick on **OK** when containers are changed.

Rfc dmjjmuqle qapccI mncIq_lb_l_sbq`jc_j_pk _arqt_rcq8



Dyespc / /2, Empty the liquid waste containers

O


Empty the two liquid waste containers and put them in place.

Qcc af_nrcp Pcnj_aqle _ Jqosyb U_qrc AmIr_qIcp ml n_ec 76,

PClick **OK** when done,

Rfc dmjjmuqle qapccI mncIq8


Npqkqle md k_qIrcI_lac jqosyb qI npmepcq,



Dyespc / /3, Priming in progress

Q

Гд `mrf jjosyb u_qrc amlr_g/cpq_pc dsjj* rfc dmjjmuyle qapccI mncIq_lb_I
_sbj`jc_j_pk _arqt_rcq8



Dyespc / /4, Both liquid waste containers full

or

Гд mlc jjosyb u_qrc amlr_g/cpq_pc dsjj* rfc dmjjmuyle qapccI mncIq8



Dyespc / /5, One liquid waste container full

R

Empty the full liquid waste container(s) and put it (them) in place.

S

Click **OK**,


?drp_dcu ky/lsrq* rfc kcqq_ec af_lecq_srmk_rja_jjw rm bcamlr_kyl_rymIyl
npмepcq8



Dyespc / /6, Decontamination in progress

Rfc bcamlr_kyl_rymI rlykc_qr_prq rm bcapc_qc,
Bspyle rfc bcamlr_kyl_rymI* ncpdmpk rfc dmjjmuyle _arqmlq R rm Tml_n_ec /33,

T Empty the two system liquid containers (2).



Dgespc /17, Empty the two system liquid containers


U Rinse the containers with demineralized water, then fill them with system liquid.

Qcc af_nrcp Pcdjyile Qwqrck Jqosyb AmIr_qlcp ml n_ec 73,

V Place the containers on the containers drawer.


Qcc af_nrcp K_l_eyle Jqosybq_lb U_qrcq ml n_ec 71,

?drqp rfc 1. kqlsrcq md bcamlr_kql_rgmI rfc dmijmuyle qapccI mnclq _srnk_rja_jjw8*



Dgespc /0., Purge of maintenance liquid in progress

?drqp _dcu kqlsrcq rfc dmijmuyle qapccI mnclq_lb_l_sbq`jc_jpk _arqt_rcq8*




Dgespc /0/, Purge of maintenance liquid in progress

W	Remove the plungers, clean first with Isopropanol, then with demineralized water and put them in the system liquid containers.
X	<p>Press OK to confirm installation of the plungers.</p> <p><i>Rfc dmjjmuyle qapccI mncIq _lb _l _sbj`jc _j_pk _arjt_rcq8</i></p>  <p>Dyespc /00, Priming of system liquid in progress...</p>
Y	Rinse the piercing tool with demineralized water.
I	<i>U fcl sqgle bckglcp_jgxcb u_rcp* k_ic qspc rf_r rfc u_rcp os_jgrw amppcqnmI bq _r jc_qr rm Ep_b 0 &QM 1474' mp RwnC GG &?QRK B//71+7/',</i>
Z	<p>Dry the piercing tool.</p> <p><i>Gd mIc mp `mrf jjosyb u_qrc amlr_gIcpq _pc dsjj*_l _sbj`jc _j_pk _arjt_rcq _lb qapccI mncIq _q bcqapq`cb qI O mI n_ec /32,</i></p>
AA	Empty the full liquid waste container(s) and put them in place.
AB	<p>Click OK,</p> <p><i>Rfc dmjjmuyle qapccI mncIq _lb _l _sbj`jc _j_pk _arjt_rcq8</i> <i>Rfc nqncrryle _pc_ uqIbm u mncIq _srnk_rja_jjw,</i></p>  <p>Dyespc /01, Install the ID-Cards piercing array</p>
AC	<p>Put the ID-Cards piercing array back in position.</p> <p><i>Qcc af_nrcp Pckmtyle rfc iB+A_p bq Njcpaqle ?pp_w Cosynncb mI n_ec /37,</i></p>

AD

Close the pipetting area window.

Гд rfc GB+A_pbg nycpaq le _pp_w qg lmr q lqr_jjcb mp lmr jmaicb amppcarjw* rfc
 dmjjmuq le qapcc l mnc lq _lb _l _sbq`jc _j_p k _anq rcq8




Dyespc /02, ID-Cards piercing array not detected

AE


When the ID-Cards piercing array is correctly installed:

Rfc dmjjmuq le qapcc l mnc lq8 **Afcaig le gl npmepcqg,**



Dyespc /03, Checking in progress

?drqp _dcu kq lsrcq* rfc dmjjmuq le qapcc l mnc lq8 **K_g lrc l _lac qsaaccbcb,**



Dyespc /04, Maintenance succeeded

AFClick **OK**.

0 At the end of the weekly maintenance, clean the gel cards storage drawer with a lint free cloth.

AG

Click on the **DATA MANAGEMENT** button to switch to the IH-Com Data Management software.

AH


Perform a **QC** run.
Pcdcp rm rfc 6F-Am k B_r_ K_l_eck c l r Qmdru_pc dmp bcr_gjcb g l dmp k_rgm l,

I

Gd rfc uccijw k_g lrc l _lac u_q g lrcppsncb*_ u_p lgle k cq q _ec b g qnj_wq rfc cppmp rf_r maasppcb,

AI

Click **OK** to acknowledge the error.
Rfc dmjjm ugle qapcc l mnc l q8



Dyespc /05, Weekly maintenance interrupted

AJ

Click **OK**.

I

? pcb k_g lrc l _lac pc k g l bcp & m l rfc rm n jcd r md rfc qapcc l' g l b g a _rcq rf_r rfc uccijw k_g lrc l _lac ksq r `c bmlc \_e\_g l\* `cdmpc _ps l a_l `c qr\_prcb, Mrfcpugqc gr gq lmr nmqqg `jc rm qr_pr _l w rcqr,



7,2,0,1 4 _l b /0 K m lrfq K_g lrc l _lac

0 The 6 and 12 months maintenance has to be performed by a service engineer.

7,2,0,2 Rmr_j Cqrq k _rcb Rq k cq _l b Jqosgb Tmjs k cq dmp Uccijw K_g lrc l _lac

Cqrq k _rcb r q k c dmp g l q rps k c l r	~80mn
Cqrq k _rcb qwqrc k j qosgb tmjs k c	830ml
Cqrq k _rcb ajc _l g l e j qosgb tmjs k c	710ml

7,2,1 OA Afcai

- 0** A quality control check must be performed after any maintenance operation.
Refer to the Data Management software for detailed information.


7,2,2 Pckmtgle rfc GB+A_pbq Ngcpag le ?pp_w Cosgnncb

- ▮ ?jj p_aiq ksqr `c pc kmtcb `cdmpc mnc lgle rfc ngncrrgle _pc_ ug l bmu, Gd rfgq am l bgrm l gq lmr pcqncarc b* rfcpc gq _pgqi md amjjgqgm l ugrf rfc p_aiq,

- ▮ Sqc ejmtcq dmp _jj mncp_rgm l q,




A Open the pipetting area window (from the software).

Qcc af_nrcp *Mnc lgle rfc Ngncrrgle ?pc_ Ug l bmu ml n_ec /0/*,



Diespc /06, Pipetting area window

- B** To unlock the ID-Cards piercing array equipped, turn the lock (4) 90° towards the back of the instrument.
- 1 R_ic a_srgm l rm lmr bpmn rfc GB+A_pbq ngcpag le _pp_w cosgnncb* rfcpc gq _pgqi rm b_k_ec rfc GF-/... ,
- 1 Gl a_qc md qfmai ml ngcpag le ng lq bspgle rfc uccijw k_g lrc l _lac* tgqs_jjw afcai dmp b_k_ec,
- C** Lower slightly the ID-Cards piercing array equipped.
- D** Pull off to the right.
- E** Remove the ID-Cards piercing array equipped.
- 1 Ufc l pcnj_ag le rfc GB+a_pb ngcpag le _pp_w cosgnncb* r_ic a_pc rm jmai gr amppcarjw, &Gl amlr_ar ugrf rfc qapcu* qcc pcb agpajc',
- F** Perform a **QC** run.
- Pcdcp rm rfc B_r_ K_l_eckclrqmdru_pc dmp bcr_gjcb g l dmp k_rgm l,




Dyespc /07,

9.5 Needle Cleaning Right/left Pipettor

This procedure is valid for both pipettors.

	Dmjmmu rfgq nrmacbspc dmp rfc lccbjc ajc_lgle,
	K_ic qspc rf_r rfc jgefr glqgbc rfc glqps kclr gq qugrafcb mdd bspgle rfc nrmacq, Qcc af_nrcp Mnclgle rfc Nyncrrgle ?pc_ Uglbmu ml n_ec /0/ ,
	?jj p_aiq ksqr `c pckmtcb `cdmpc mnclgle rfc nyncrrgle _pc_ uglbmu, Gd rfgq amlbgrml gq lmr pcqncarcb rfcpc gq _pgqi md amjjgqgm l ugrf rfc p_aiq,
	Sqc ejmtcq dmp _jj mncp_rgm lq,
	?ju_wq nj_ac _rqqsc slbcp rfc lccbjc `cdmpc bqcamllcargle rfc d_p_bqxcb tcgl _q gr bpmnq,
	Sqc _qmdr rqqsc rm _tmgb b_k_ec ml rfc lccbjc am_rgle,
0	The needle cleaning must be performed when the needle is visually soiled.
0	If NEEDLE CHANGE RIGHT PIPETTOR (or left pipettor) requires a user level greater than 1, this button will be grayed for level 1 users.
A	Click on NEEDLE CHANGE RIGHT PIPETTOR . ? kccq_ec _nnc_pq* _qigle rm mncl rfc nyncrrgle _pc_ uglbmu,
B	Go back to the main screen to open the pipetting area window.
C	 From the main screen, click on the pipetting area in the interactive image. Qcc af_nrcp Afcagile Nyncrrgle ?pc_ Qr_rsq ml n_ec //7,
D	Open the pipetting area window (from the software). Dmjmmu rfc nrmacbspc `cjm u rm af_lec rfc q_knjle lccbjc dmp qcpqcg,
E	Place a clean tissue below the needle.

F Disconnect the faradized vein (5) from the needle using the connector (3).



Dyspc / 1.,

G Open the 2 locks (2).

H Slide the needle (1) out of the housings (4).

I Verify the needle is purged.

J Clean the needle with a nonwoven wipe and isopropanol solution.

Uqnc rfc lccbjc dpm k rfc rmn rm rfc rgn md rfc lccbjc,

K Dry the needle with a dry nonwoven wipe:

Uqnc rfc lccbjc dpm k rfc rmn rm rfc rgn md rfc lccbjc,

L	Check for any sign of damage (needle coating, bent needle, etc.). <i>Gd rfc lccbjc q q b_k_ecb* pcnj_ac q, Qcc af_nrcp Lccbjc Af_ lcc Pqefr-jcdr Njncrrmp ml n_ec /42,</i>
M	Locks open (2), needle absent, connector (3) near by.
N	Slide the needle (1) in the housings (4).
O	With the needle in place, screw the faradized vein (5) onto the needle until it clicks.
P	Lock the needle by tilting the two locks (2). <i>Jmai q ajmqcb* lccbjc q q am l lcarcb _ lb jmai cb q l nmqyrqm l ml qrq njncrrmp X_vq qsnmmp,</i>
O	The collar on the needle, when properly inserted, prevents the needle from moving independent of the mechanism.
Q	Clean the needle carefully with a dry nonwoven wipe.
R	Close the pipetting area window. <i>Qcc af_nrcp Mnc l q le rfc Njncrrmp le ?pc_ Uq l bmu ml n_ec /0/,</i>
S	Click OK to validate needle replacement. <i>Rfc lccbjc q q _bhsqrcb _ srm k _ rqa_jjw,</i>
T	Perform a QC run. <i>Pcdcp rm rfc B_r_ K_ l _ ec k c l r qmdru _ pc dmp bcr _ qjcb q l dmp k _ rqm l,</i>
I	<i>Gd rfc lccbjc ajc _ l q le f_q `cc l q lrcppsncrb bsc rm _ nmucp mdd* pc+qr_pr rfc lccbjc ajc _ l q le nrmacbspc,</i>
O	External light can have an influence on the needle adjustment.

9.6 Needle Change Right/left Pipettor

This procedure is valid for both pipettors.

1 ?jj p_aiq ksq r`c pc kmtcb `cdmpc mnc l g l e rfc n g n c r r g l e _pc_ u g l b m u, Gd r f g q a m l b g r g m l g q l m r p c q n c a r c b * r f c p c g q _ p g q i m d a m j j g q g m l u g r f r f c p _ a i q,

2 S q c e j m t c q d m p _ j j m n c p _ r g m l q,

3 ?j u _ w q n j _ a c _ r g q q s c s l b c p r f c l c c b j c `c d m p c b g q a m l l c a r g l e r f c d _ p _ b g x c b t c g l _ q g r b p m n q,

0 If **NEEDLE CHANGE RIGHT PIPETTOR** (or left pipettor) requires a user level greater than 1, this button will be grayed for level 1 users.

A Click on **NEEDLE CHANGE RIGHT PIPETTOR**.

? k c q q _ e c _ n n c _ p q * _ q i g l e r m m n c l r f c n g n c r r g l e _ p c _ u g l b m u,

B Go back to the main screen to open the pipetting area window.

C From the main screen, click on the pipetting area in the interactive image.



Q c c a f _ n r c p A f c a i g l e N g n c r r g l e ? p c _ Q r _ r s q m l n _ e c / / 7,

D Open the pipetting area window (from the software).

Q c c a f _ n r c p M n c l g l e r f c N g n c r r g l e ? p c _ U g l b m u m l n _ e c / 0 /,

E Place a clean tissue below the needle.

F Disconnect the faradized vein (5) from the needle using the connector (3).

G Open the 2 locks (2).

H Slide the needle (1) out of the housings (4).

I Replace the sampling needle (1).

A f c a i d m p _ l w q g e l m d b _ k _ e c & l c c b j c a m _ n g l e * `c l r l c c b j c * c r a , ' ,
G d r f c l c c b j c g q b _ k _ e c b * r _ i c _ l c u m l c ,

J Locks open (2), needle absent, connector near by.


K Slide the needle (1) in the housings (4).

L With the needle in place, screw the faradized vein (5) onto the needle until it clicks.

M Lock the needle by tilting the two locks (2).

J m a i q a j m q c b * l c c b j c g q a m l l c a r c b _ l b j m a i c b g l n m q g r g m l m l g r q n g n c r r m p X _ v g q q s n n m p r ,

0 The collar on the needle, when properly inserted, prevents the needle from moving independent of the mechanism.



Dgespc /1/,

N

Close the pipetting area window.

Qcc af_nrcp MncIq/le rfc Nyncrnq/le ?pc_ Uq/bmu m/ n_ec /0/,

O

Click **OK** to validate needle replacement.

Rfc lccbjc qg_blsqrcb_srm k_rqa_jjw,

P

Proceed the same way with the left pipettor.

Q

Perform a **QC** run.

Pcdcp rm rfc B_r_ K_l_eckclr qmdru_pc dmp bcr_gjcb q/dmp k_rqm l,

! Cd rfc lccbjc pcnj_ackclr f_q `ccl glrcppsncrb `w _dmpacb qfsr bmu l* ncpdmp k rfc lccbjc af_lec nrmacbspc _e_g l &ugrfmsr af_legle rfc lccbjc &Mljw qrcnq ? rm B m/ n_ec /42_ lb L rm N m/ n_ec /43 ksqr `c bmlc'' m/ rfc lcvr pcqr_pr md rfc qwqrc k, Rfgq qg rm k_ic qspc rf_r rfc lccbjc qg amppcarjw jma_jgxcb,

9.7 Rinse Right/Left Pipettor

The button **RINSE RIGHT PIPETTOR** or **RINSE LEFT PIPETTOR** allows to manually launch a wash cycle of the right or left pipettor, additionally to the automatic rinse during the runs.

- 0** If the buttons **RINSE RIGHT PIPETTOR** or **RINSE LEFT PIPETTOR** requires a user level greater than 1, the button is grayed for level 1 users.

9.8 Reading Station Calibration

After 5000 working hours, the user is asked to do a calibration of the reading module.

- 0** A calibration can also be started before the 5000 working hours have passed.


A Put 24 5052 type gel cards (NaCl cards) in the gel cards storage drawer.

- 0** Verify that the gel of the gel cards is in good condition and that there are no air bubbles present. The gel card must be valid, check the expiry date.

B Click on the **READING STATION CALIBRATION** button to start the calibration process.



Rfc dmjjmuyle qapccI q q bjqnj_wcb,




Dyespc / 10, Calibration in progress

- 0** If the **READING STATION CALIBRATION** requires a user level greater than 1, the button is grayed for level 1 users.

- 0** The process can take up to 20 minutes. No action on the IH-1000 software can be performed during the calibration process.

U fcl rfc nрmacq q q d q l q fcb rfc dmjjm u q l e k c q q _ e c q q b q q n j _ w c b ,*



Despc / 11, Calibration finished

C Click on **OK**.

O After calibration, the 24 gel cards go back to the drawer and can be used for further tests.

O If launching conditions are not met, an alert screen with the missing items is displayed.

D Perform a **QC** run.

P c d c p r m r f c B _ r _ K _ l _ e c k c l r q m d r u _ p c d m p b c r _ q j c b q l d m p k _ r y m l ,

9.9 General Initialization

O Do not perform an initialization while the instrument is running. The IH-1000 makes all the initialization automatically.

See [Starting the Instrument on page 46](#) and/or [Errors on page 42](#).

9.10 Storage

When it is planned not to use the instrument during many days, specific operations must be considered.

7, / ., /

Q f m p r R c p k Q r m p _ e c g l J _ ` m p _ r m p w A m l b g r g m l q

If the instrument is not used for few weeks and if it stays in laboratory environment, no wrapping or packing is required.

Uccijw k_g lrc l_lac mncp_rgm lq ksq r `c ncpdmp kcb ctc pw ucci ctcl gd rfc glqrps kclr gq lmr sqcb,

7, / ., 0

M r f c p Q r m p _ e c A m l b g r g m l q

When the instrument is stored outside laboratory conditions or if it is stored for a long period, special procedures must be performed by a service engineer.

0 Refer to a service engineer for handling and storage procedures.

10 Troubleshooting

Af_nrcp mtcptgcu

- Rfc sqcp ksqf f_tc pc_b_lb slbcprmmmb rfgq bmask clr_rgm l`cdmpc a_ppwg le msr_lw
rpms`jcqfmmrg le ml rfc GF-/... ,Gl a_qc md slajc_p gl dmp k_rgm l* njc_qc amlr_ar rfc
bgqrpg`srmp,

10.1 Alarms and Error Management

- 0 When an error message appears on the screen, read it and check whether an explanation is given in the following list. Then click on the OK button and execute the required procedure to solve the problem.

10.2 Software Error Messages

0 For any error messages not listed in the following tables (mainly module and software errors), please contact a service engineer.

/ .,0,/

?jj Kmbsjcq

Cpmp ambc	Cpmp bcqapnrgm l	Amppcargtc _argm lq
"0 . . . / . B	Hardware board communication failure.	Restart or contact a service engineer.
"0 . . . 1. 7	Hardware board failure.	Restart or contact a service engineer.
"0 . . . 1. B	Hardware board communication failure.	Restart or contact a service engineer.
"0 . . . 1. D		
"0 . . . 2. 1	Robotic movement failure.	Try to continue, restart or contact a service engineer.
". 0 . . . 2. ?	System execution failure.	Restart or contact a service engineer.
"0 . . . 2. A	Robotic movement interrupted.	Try to continue, restart or contact a service engineer.
"0 . . . 2. C	Robotic movement failure.	Restart or contact a service engineer.
"0 . . . 2/2		
"0 . . . 6. /@	Action impossible due to previous failure.	Solve previous failure or contact a service engineer.
"0 . . . DD . .	Robotic movement failure.	Restart or contact a service engineer.

/ .,0,0

Rp_l qnmpr ?p k

Cpmp ambc	Cpmp bcqapnrgm l	Amppcargtc _argm lq
"0 . . . 2. 2	Robotic movement failure during initialization.	Restart or contact a service engineer.
"0 . . . 2/ /	Axis drag overshoot.	Try to continue or disable the faulty component and contact a service engineer.
"0 . . . 2/1	Axis positioning error.	Try to continue or disable the faulty component and contact a service engineer.
"0 . . . 6. /7	Card presence sensor not working during check.	Try to continue, restart the system and contact a service engineer.
"0 . . . 60. A	Card transport failure.	Check the presence of the gel card. Try to continue or restart the system.

Cpmp ambc	Cpmp bcqapnrgm l	Ampcargtc _argmlq
"0 . . . 61.2	Centrifuge trapdoor failure.	Try to continue or disable the centrifuge and contact a service engineer.
"0 . . . 61.3	Centrifuge trapdoor failure.	Try to continue or disable the centrifuge and contact a service engineer.

/ . , 0,1

Gbc lrgdga _rgm l Kmbsjc

Cpmp ambc	Cpmp bcqapnrgm l	Ampcargtc _argmlq
"0 . . . 2.2	Robotic movement failure during initialization.	Restart or contact a service engineer.
"0 . . . 2 / /	Axis drag overshoot.	Try to continue or disable the faulty component and contact a service engineer.
"0 . . . 2 / 1	Axis positioning error.	Try to continue or disable the faulty component and contact a service engineer.
"0 . . . 62.0	The rack contains a clogged tube or vial. This error appears if the option "cap detection" is activated.	<ul style="list-style-type: none"> Remove the cap on the sample tube or the vial. See No Cap Confirmation on page 84 or Removing Racks on page 88 If the problem persists, contact a service engineer.
"0 . . . 62.2	The rack barcode is unreadable.	Missing or damaged barcode label: <ul style="list-style-type: none"> Verify if the rack is correctly labeled. If the problem persists, contact a service engineer.
"0 . . . 62.5	Missing or damaged rack barcode label.	Verify if the rack is correctly labeled. If the problem persists, contact a service engineer.
"0 . . . 62.7	Presence of a rack impossible to process at the starting.	The instrument was stopped without using the shutdown function: <ul style="list-style-type: none"> Remove the rack with the software. See chapter Removing Racks on page 88.
"0 . . . 63. /	The reagent vial barcode is unreadable.	Reagent vial incorrectly positioned on the rack: <ul style="list-style-type: none"> Verify if the reagent vial is correctly inserted. Damaged barcode label or symbology out of specification: <ul style="list-style-type: none"> Verify if the reagent vial is correctly labeled. If the problem persists, contact a service engineer.
"0 . . . 63.1	The reagent vial is out of date at the location requested.	<ul style="list-style-type: none"> Remove the rack. See chapter Removing Racks on page 88. Change the reagent vial.
"0 . . . 63.2	The reagent vial is unusable at the location requested.	<ul style="list-style-type: none"> Remove the rack. See chapter Removing Racks on page 88. Change the reagent vial.

Сппмп амбс	Сппмп бсқарғынғым І	Амппсартс_арғым І
"0...63.3	Empty barcode on the reagent vial is refused.	Incorrect positioning of the vial in the rack: See chapter Removing Racks on page 88 . <ul style="list-style-type: none"> • Reposition the reagent vial. • If the problem persists, contact a service engineer.
"0...64./	The diluent wells rack is out of date.	Rack out of date. Remove the rack. See chapter Removing Racks on page 88 .
"0...64.0	The barcode of the diluent rack is duplicated.	Two diluent rack with the same barcode. Remove all racks. See chapter Removing Racks on page 88 . <ul style="list-style-type: none"> • If the problem persists, contact a service engineer.
"0...65./	The sample tube barcode is unreadable at the location requested.	Sample tube incorrectly positioned on the rack: <ul style="list-style-type: none"> • Verify if the sample tube is correctly inserted. Damaged barcode label or symbology out of specification: <ul style="list-style-type: none"> • Verify if the sample tube is correctly labeled. <p>Қд рфс ` _памбс қс_пс_б_к_Іс_ііж* к_іс қспс рм нsr рфс қ_кнјс ` _аі рм рфс мрґег І_ј нмқсғрғм І,</p> <ul style="list-style-type: none"> • If the problem persists, contact a service engineer.
"0...65.1	The sample tube barcode is duplicated.	<ul style="list-style-type: none"> • Take the rack out again and check the barcodes. • If the problem persists, contact a service engineer.
"0...65.2	The sample tube has an embarked clot.	Blood too thick: <ul style="list-style-type: none"> • Check that the blood is not over 48 hours old. • If the problem persists, contact a service engineer.

/.,0,2

Ac І рғдсес

Сппмп амбс	Сппмп бсқарғынғым І	Амппсартс_арғым І
"0...2.2	Robotic movement failure during initialization.	Restart or contact a service engineer.
"0...2//	Axis drag overshoot.	Try to continue or disable the faulty component and contact a service engineer.
"0...2/1	Axis positioning error.	Try to continue or disable the faulty component and contact a service engineer.
"0...61.4	Centrifuge balance failure.	Insert gel cards in the drawer.
"0...DB.2	Speed error.	Try to continue, restart and contact a service engineer.

/ .,0,3 Pc_bg le Kmbsjc

Cppmp ambc	Cppmp bcqapgnrjm l	Amppcargtc _argm lq
"0 . . . 2.2	Robotic movement failure during initialization.	Restart or contact a service engineer.
"0 . . . 2//	Axis drag overshoot.	Try to continue or disable the faulty component and contact a service engineer.
"0 . . . 2/1	Axis positioning error.	Try to continue or disable the faulty component and contact a service engineer.
"0 . . . 67. /	Reading module calibration failure.	Launch a new calibration or contact a service engineer.
"0 . . . 67.0		
"0 . . . 67.1	Expired reading module calibration.	Launch a new calibration or contact a service engineer.
"0 . . . 67.3	Reading module calibration failure.	Launch a new calibration or contact a service engineer.
"0 . . . 67.7		
"0 . . . 67. ?	Reading module failure.	Launch a new calibration or contact a service engineer.
"0 . . . 67. @		
"0 . . . 67. A		
"0 . . . 67. B		
"0 . . . 67//	Gel card detected during lighting control.	Restart system of contact a service engineer.
"0 . . . 67/1	Lost gel card.	Restart system or contact a service engineer.
"0 . . . 6B.3	Reading module failure.	Launch a new calibration or contact a service engineer.
"0 . . . 6B.4		
"0 . . . 6B.5		

/ .,0,4 Ngncrrmpq

Cppmp ambc	Cppmp bcqapgnrjm l	Amppcargtc _argm lq
"0 . . . 2.2	Robotic movement failure during initialization.	Restart or contact a service engineer.
"0 . . . 2//	Axis drag overshoot.	Try to continue or disable the faulty component and contact a service engineer.
"0 . . . 2/1	Axis positioning error.	Try to continue or disable the faulty component and contact a service engineer.
"0 . . . 6. / ?	Shared area conflict.	Restart or contact a service engineer.
"0 . . . 6.5/	Diluent piercing failure.	Verify that the rack is fully inserted in its support.
"0 . . . 6.51		

/ .,0,5

Fwbp_sjga

Cppmp ambc	Cppmp bcqapnrym l	Amppcargtc_argm lq
"0 . . . 2 / /	Axis drag overshoot.	Try to continue or disable the faulty component and contact a service engineer.
"0 . . . 2 / 1	Axis positioning error.	Try to continue or disable the faulty component and contact a service engineer.
"0 . . . 2 / 4	Robotic movement failure.	Restart or contact a service engineer.
"0 . . . 6 . 3	Both pipettors out of order for maintenance.	Contact a service engineer.
"0 . . . 6 . 42	Soiled needle cannot be used.	Check washing means or contact a service engineer.
"0 . . . 6 . 4B	Washing failure.	Restart or contact a service engineer.
"0 . . . 6 @ . 0	Level detection failure.	Try to continue or disable pipettor and contact a service engineer.
"0 . . . 6 @ . 1		
"0 . . . 6 @ . 2	Needle external washing volume error.	Try to continue or disable pipettor and contact a service engineer.
"0 . . . 6 @ . 3		
"0 . . . 6 @ . @	Incorrect state of the hydraulic management.	Check liquid containers and circuits and restart the machine or contact service engineer.
"0 . . . 6 @ / 0	Clogged needle circuit.	Try to continue, follow the needle changing procedure: See chapter Needle Change Right/left Pipettor on page 164 .
"0 . . . 70 . /		or disable the pipettor and contact a service engineer.
"0 . . . 6 @ / 3	Both liquid waste containers are full or disconnected.	Empty containers and check caps are completely screwed.
"0 . . . 6 @ / 5	Pipettor pressure sensor failure.	Try to continue or disable the pipettor and contact a service engineer.
"0 . . . 6 @ / 6	Washing pump failure.	Try to continue or disable the pipettor and contact a service engineer.
"0 . . . 6 C . /	Fluidic calibration failure.	Launch a new calibration or contact a service engineer.
"0 . . . 6 C . 0	Fluidic flow control failure.	Try to restart and contact a service engineer.
"0 . . . 6 C . 1		
"0 . . . 6 C . 2		
"0 . . . 6 C . 3		
"0 . . . 6 C . 4		
"0 . . . 6 C . 5		
"0 . . . 6 C . 6	Fluidic calibration failure.	Launch a new calibration or contact a service engineer.
"0 . . . 6 C . 7		
"0 . . . 6 C . ?	Fluidic software execution failure.	Restart or contact a service engineer.

Cppmp ambc	Cppmp bcqapnrgm l	Amppcargtc _argm lq
"0...6C.A	Fluidic calibration failure.	Launch a new calibration or contact a service engineer.
"0...70.1	Missing decontamination solution.	Load a new cleaning solution SETUP CLEAN vial. If the pipettor status is in default (red), <ul style="list-style-type: none"> • acknowledge the alarm; • load the setup clean; • start a test on the same side as the concerned pipettor.
"0...70//	Clot detected.	Execute needle cleaning process. See chapter Rinse Right/Left Pipettor on page 166 .
"0...CB.5	Clot detected and solved, sampling cancelled.	Try to continue or restart.
"0...CB.6	Permanent clot detected, unable to solve automatically.	Execute appropriate procedure. See chapter Clot Detection on page 123 .
"0...CB.@	Fluidic system failure (pressure sensor).	Try to restart and contact a service engineer.
"0...D4./	Calibration failure I/O.	Launch a new calibration or contact a service engineer.
"0...D4.6		
"0...D4.7		
"0...D4.?	Software execution failure I/O.	Restart or contact a service engineer.
"0...D4.A	Calibration failure I/O.	Launch a new calibration or contact a service engineer.

/.,0,6

Pc_eclr P_ai

Cppmp ambc	Cppmp bcqapnrgm l	Amppcargtc _argm lq
"0...D4.B	Reagent agitating speed out of range.	Use another rack or contact a service engineer.
"0...D4.C	Reagent agitating failure.	Use another rack or contact a service engineer.

/.,0,7

Ecj A_pbq Ngcpag le Kmbsjc

Cppmp ambc	Cppmp bcqapnrgm l	Amppcargtc _argm lq
"0...2.2	Robotic movement failure during initialization.	Restart or contact a service engineer.
"0...2//	Axis drag overshoot.	Try to continue or disable the faulty component and contact a service engineer.
"0...2/1	Axis positioning error.	Try to continue or disable the faulty component and contact a service engineer.

/ .,0,/ . ?jj A_k cp_q

Cpmp ambc	Cpmp bcqapnrym l	Amppcargtc _argm lq
"0 . . . B.0	Reading module calibration failure.	Launch a new calibration or contact a service engineer.
"0 . . . D. .5	Camera communication failure.	Restart system or contact a service engineer.
"0 . . . D. .6	Camera calibration failure.	Launch a new calibration or contact a service engineer.
"0 . . . D. .7	Camera failure.	Launch a new calibration or contact a service engineer.
"0 . . . D. .?		
"0 . . . D. .@		

/ .,0, / / Qwqrc k

Cpmp ambc	Cpmp bcqapnrym l	Amppcargtc _argm lq
"0 . . . / . /	Hardware board communication failure.	Restart or contact a service engineer.
"0 . . . / .0		
"0 . . . / .1		
"0 . . . / .2		
"0 . . . / .3		
"0 . . . / .4		
"0 . . . / .5		
"0 . . . / .6		
"0 . . . / .7		
"0 . . . / .?		
"0 . . . / .@		
"0 . . . / .A		
"0 . . . / .C		
"0 . . . / .D		
"0 . . . / / .		
"0 . . . / / /		
"0 . . . / /0		
"0 . . . / /1		
"0 . . . / /2		
"0 . . . / /3		
"0 . . . 7. /	File access failure.	Restart or contact a service engineer.
"0 . . . 7.0		
"0 . . . 7.1		
"0 . . . 7.2		
"0 . . . 7.3		

Cppmp ambc	Cppmp bcqapnrgm l	Amppcargtc _argm lq
"0...?. / "0...?.0 "0...?.1 "0...?.2 "0...?.3 "0...?.4	Software execution failure.	Restart or contact a service engineer.
"0...C.@ "0...C.A	File access failure.	Restart or contact a service engineer.
"0...60.B	Pipetting door not closed.	Close the pipetting area window.
"0...65.4 "0...65.5	No message - Pictogram on the GUI (software interface). See chapter Sample Status on page 78 .	Check the loaded crossmatch receivers.
"0...7.. / "0...7..0	File access failure.	Restart or contact a service engineer.
"0...7/.@ "0...7//0 "0...7//1 "0...7//2	Inconsistent assay database information.	Check input data or contact a service engineer.
"0...7//3	Missing donor.	Check input data. Enter the missing sample.
"0...D5.0 "0...D5.1 "0...D5.2 "0...D5.3 "0...D5.4	Hardware board communication failure.	Restart or contact a service engineer.

N?EC GLRCLRGML?JJW JC DR @J?L I

11 Appendix

Af_nrcp Mtcptqcu

This chapter covers non-standard operations on the instrument.

11.1 Equipment and Accessories

//,/,/ Cosgn k c l r

Bcqapnrgm l	Os_ l rgrw
IH-1000	1
Reagents rack	4
Diluents rack holder	4
Samples rack	18
Samples rack barcode label set	18
Rack extractor tool	1
Decontamination support	1
Touchscreen	1
Keyboard (language depending on country)	1
Touchscreen arm	1
Power cord	1
Preparation table	1
Containers	5
Gel card waste bin and bag	1
Technical documentation	1
User manual (ref. H009142)	1

//,/,0 ?aacqqmpgcq

Bcqqnrgm l	Os_lrgw
Barcode scanner	1
Hub	1
USB key 8Go	1
UPS	1
UPS battery	1
UPS holder	1

11.2 CE Compliance

This instrument is an in-vitro diagnostic medical device and complies with all applicable European Community Directives and associated harmonized standards, including but not limited to the standards pertaining to the electrical safety as well as the emission and immunity requirements as specified in IEC 61326-2-6.

11.3 Documents / Forms

//,1,/ Bcamlr_kgl_rgm l Acprgdga_rc

0 Please contact your local distributor for the official electronic document.

■ Rfc bcamlr_kgl_rgm l ksq r`c ncpdmp kcb`w_l _srfmpgxcb rcaflgaq_l gl _aampb_lac ugrf rfc glqrpsargm lq amlr_glcb gl rfc amppcqnm lbg l e Qcptgac K_Is_j,


The IH-1000 and its accessories must be cleaned and decontaminated before any of the following intervention:

- performing any maintenance;
- revision task;
- before taking the instrument out of service;
- when preparing for transportation;
- before final elimination.

■ Rfc mpgeg l_j acprgdga_rc ksq r pck_g l gl rfc npcqc l r k_Is_j, Ufcl _acprgdga_rc gq pcosgpcb* k_ic _amnw _lb dgjj gl rfc j_rrcp q_k njc, Dgjj gl _amnw md rfc bcamlr_kgl_rgm l acprgdga_rc _lb nj_ac gr ml rfc glqrps kcl r* gl _amlqngasmsq nj_ac,

The absence of the decontamination certificate may extend the time needed for revision and maintenance.

Bio-Rad will bill any cleaning and decontamination operation.



--

--

--

--

--

--

--

--

[illegible]

//,1,0 K_g|rc|_lac Pcampb

0 The original blank maintenance sheet must remain in the present manual. When a maintenance record is required, make a copy and fill in the latter sample.

Type

Page

Serial number

K_g rc _lac - Bcdcar	?arqm l	Mncp_rmp	B_rc

11.4 Spare Parts Ordering

Please contact your local service engineer.

11.5 Device Disposal

//,3,/ EcIcp_j U_pIgleq

ГІ mрbср rm nрmrcar rfc nсрqmІq Іb rfc cІ tgpml k cІr* Іw g Іqrps k cІr Іb grq _aacqmpgcq k sqr `c bgqnmqcb md g І І \_nnpmnpq\_rc u\_w, Grq gq k Іb\_rmpw rm qrpqarjw \_nnjw j\_uq Іb jma\_j `wj_uq pcj_rqtc rm І І _nnpmnpq_rc bgqnmq_j nрmacbspc,

І ? І g Іqrps k cІr k_w m Іjw `c bgq kms Іrcb Іb bgtgbcб g І rm grq `qga am k nm І c Іrq `w \_os\_jgdqcb rcaflgaq І, N\_prq Іb qs`+n_prq pc k mtcb dpm k І І g Іqrps k cІr k_w Іmr `c pcsqcb sljcqq\_drcp f\_tgle `cc І _nnpmmtcb g І upgrgle `w rfc k Іsd\_arspc, ? Іw dsrspc \_nnjga\_rgm І f\_q rm `c npcagqcjw bcqapq`cb g І rfc pcoscqr,

//,3,0 Bgqnmq_j

The following components can be found in the instrument:

N_aigle k_rcpq_j

N_ncp`m_pb

Nmjwcrfwjcl c dm_k

? Іrq+qr_rqa `e

Aj_qq GG dgpсnрmmdcб nmjwqrwpc І c

G Іqrps k c Іr

Nj_qrqa k_rcpq_jq PVC, POM, nylon, rubber

Cjcarpgaj am k nm І c Іrq Motors, copper wires, plastic insulation

Cjcarpm Іga am k nm І c Іrq Epoxy printed boards

Mnrja_j am k nm І c Іrq Camera lens, camera body

Kcr_jq **Amnnсp:** electrical wires;

?js k j І sk: instrument base, holders, guides for electronics boards, partitions, covers;

Qrcсj: springs, holders, screws, bolts.

11.6 List of Documents

//,4,/ Sqcp K_Is_jq

J_les_ec	GF+... Sqcp K_Is_j	GF+Amk Sqcp K_Is_j
Clejgqf	H009142	H009146
Ecpk_l	H009156	H009147
Dpclaf	H009159	H009151
Gr_jq_l	H009160	H009149
Qn_lgqf	H009157	H009148
Nmprsescqc	H009158	H009150
B_lgqf	H009161	H009153
Qucbgqf	H009162	H009154
Lmpuceq_l	H009163	H009155
Dgllgqf	H009166	H009168
H_n_lcqc	H009167	H009169
Nmjgqf	H009170	H009171
Fsle_pq_l	H009190	H009183
Qcp`q_l	H009191	H009192
Afglcqc	H009195	H009365
Nmprsescqc &p_xgj'	H009204	H009205
Qjmtclq_l	H009216	-
Psqqq_l	H009221	H009222
Axcaf	H009194	H009196

//,4,0 Qcptgac K_Is_jq

J_les_ec	GF+... Qcptgac K_Is_j	GF+Amk Qcptgac K_Is_j
Clejgqf mljw	H009144	H009152

//,4,1 Qn_pc N_prq K_Is_j

J_les_ec	GF+... Qn_pc n_prq K_Is_j
Clejgqf mljw	H009215

Index

- ?**
- Alarm and Error Management
 - Behavior 169
 - Authorized Maintenance 147
- A**
- Components area
 - Checking status 144
 - Computer connecting panel 27
 - Connection panel 27
- B**
- Definition
 - Major Injury 13
 - Minor Injury 13
 - Moderate Injury 13
 - Device disposal 183
 - Diluent
 - Rack removal 90
- D**
- Footer strip
 - Resources buttons
 - Description 101
 - Footer strip 101
 - Resources on board 102
 - Test buttons
 - Footer strip 107
 - In progress
 - Details of a test 110
 - General description 108
 - Ready to start
 - Emergency 112
 - General description 111
 - Tests to repeat (description) 114
 - Without resources (description) 113
 - Tests buttons
 - Samples without request (description) 117
 - Forms (decontamination certificate) 180
- E**
- Gel cards area
 - Interactive instrument image 38
 - Managing ID-Cards 58
 - Replacing the ID-Cards waste bin 100
- G**
- Intended Use 11
- J**
- List of documents 184
- K**
- Menu button
 - Exit (software) 48
 - Maintenance screen 146
 - Daily maintenance 149
 - General initialization 167
 - Needle change (right and left pipettor) 164
 - Needle cleaning (right and left pipettor) 161
 - Operations (decontamination) 148
 - Pipettors (rinse) 166, 167
 - Reading station calibration 166
 - Weekly maintenance 149
 - Options screen 126
 - Controls options 134
 - Disabled Mode 137
 - Cap detection 139
 - General options 128
 - Profile options 132
 - Stop All (instrument) 49
 - User Level screen 147
- L**
- Needle change
 - Procedure (right and left pipettor) 164
 - Needle cleaning
 - Procedure (right and left pipettor) 161
 - Norms
 - Safety 16
- N**
- Pipetting area
 - Clot detection 123
 - Clot detection (errors) 122
 - Interactive instrument image 38
 - Pipetting area window (opening) 121
 - Status (screen) 119
- P**
- Reagent
 - Rack removal 90
 - Rescue shutdown 50
 - Risk
 - Classes definition 16
- Q**
- Safety
 - Environmental conditions 18
 - General instructions 19
 - Importance of the instructions 17
 - Non-compliance with rules 18
 - Safety instructions
 - Warning 16

Sample / Reagent Area

Cap confirmation 84

Diluent

Diluent rack and general description 72

Inserting a diluent rack product 73

Status (diluent) 79

General rack status 78

Interactive instrument image 38

Loading racks 86

Main screen (description) 77

Reagent

Cleaning solution

Preparation (SETUP CLEAN) 76

Use (SETUP CLEAN) 75

Reagent rack and general description 74

Status (reagent) 79

Sample

Barcode specifications 70

Input (Sample barcode manual input) 83

Material (sample) 14

New (sample rack) 69

Status (sample) 78

Type of sample tubes

Labeling 68

Urgent (sample rack) 69

Volumes (sample) 70

Solutions / Wastes area

Containers (cleaning) 97

General description 93

Interactive instrument image 38

Liquid and wastes status 94

Preparation (cleaning liquid) 96

Priming cycle 99

System liquid container (refilling) 95

Waste liquid container (replacement) 98

Step by step procedure 51

Storage 168

R

Troubleshooting

Alarm and Error Management 169

S

Use

Intended 11

Introduction 23

User Level 147

U

Warranty

Limitations 12

[illegible]



**Bio-Rad
Laboratories**

Clinical
Diagnostics Group

Website www.bio-rad.com/diagnostics **Australia** 61-2-9914-2800 **Austria** 43-1-877-8901 **Belgium** +32 (3)710-53-00 **Brazil** +55 (31)3689-6600 **Canada** 1-514-334-4372
China 86-21-61698500 **Czech Republic** 420-241-430-532 **Denmark** +45-4452-1000 **Finland** 358-9-804-22-00 **France** 33-1-47-95-60-00 **Germany** +49 (0)89-318-840
Greece 30-210-7774396 **Hong Kong** 852-2789-3300 **Hungary** +36-1-459-6100 **India** 1800-180-1224 **Israel** 972-3-9636050 **Italy** +39-02-216091 **Japan** 81-3-6361-7070
Korea 82-2-3473-4460 **Mexico** +52 (55)5488-7670 **The Netherlands** +31-318-540666 **New Zealand** 64-9-415-2280 **Norway** +47-23-38-41-30 **Poland** 48-22-3319999
Portugal 351-21-472-7700 **Russia** +7-495-721-1404 **Singapore** 65-6415-3170 **South Africa** 27-11-442-85-08 **Spain** 34-91-590-5200 **Sweden** +46-8-555-127-00
Switzerland +41 (0)26-674-55-05/06 **Taiwan** 886-2-2578-7189 **Thailand** 662-651-8311 **United Kingdom** +44 (0)20-8328-2000

IH-1000

խնդրանշանային փոխարկման սարք



$RH_j dz ds q^w_{ww} d\mathbb{F} \Theta y d^w$

† Լճէն Յէմիլ Յէն Դէն † Գրիգոր Զեօն Եւրոպ

IH-1000 1.0 2.0 - 01/2022

rt54zd3jdzs: H009142

REF 001000 / 001000RECOND

IVD CE



† t6tsqLətsHqIsjdz †

DiaMed GmbH

Pra Rond 23

1785 Cressier FR

Switzerland

uid. +41 (0)26 67 45 111

+41 (0)26 67 45 145

1 800-874-6726 www.bio-rad.com



Bio-Rad

3 Boulevard Raymond Poincaré®

92430 Marnes-la-Coquette

France

Ujdz. +33 1 47 95 60 00

Copyright É

[illegible]

၂။ ရန်ကုန်တိုင်းဒေသကြီး

1.1	1.2	1.3	1.4	1.5	1.6	1.7	2.0
1.0							
1.1	V3.0						
1.2	V4.0						
1.3							
1.4	V04.04						
1.5	V04.07						
1.6	V04.08						
1.7							
2.0							

ԹՃՐՐԷՐՓԼՔԷ

	Քննարկի քննարկ	2
Դժժանություն 1	Ջրի միջավայրի վնասվածություն	11
1.1	Լուծույթի միջավայրի վնասվածություն	11
1.2	Ջրի միջավայրի վնասվածությունից պայմանավորված ռիսկի գնահատում	12
1.3	Մշակույթի վնասվածություն	12
1.3.1	Դժվարություն	13
1.3.2	Լուծույթի միջավայրի վնասվածություն	13
1.3.3	Ջրի միջավայրի վնասվածությունից պայմանավորված ռիսկի գնահատում	14
1.4	Սոցիալական և տնտեսական վնասվածություն	14
1.4.1	Ջրի միջավայրի վնասվածություն	14
1.4.2	Տնտեսական վնասվածություն	15
1.4.3	Լուծույթի միջավայրի վնասվածություն	15
1.4.4	Լուծույթի միջավայրի վնասվածությունից պայմանավորված ռիսկի գնահատում	15
1.4.5	Լուծույթի միջավայրի վնասվածությունից պայմանավորված ռիսկի գնահատում	15
1.4.6	Միջավայրի վնասվածությունից պայմանավորված ռիսկի գնահատում	15
1.4.7	Ջրի միջավայրի վնասվածությունից պայմանավորված ռիսկի գնահատում	15
1.4.8	Քննարկի վնասվածություն	16
1.4.9	Լուծույթի միջավայրի վնասվածություն	16
1.4.9.1	Լուծույթի միջավայրի վնասվածություն	16
1.4.9.2	Լուծույթի միջավայրի վնասվածություն	16
Դժժանություն 2	Ջրի միջավայրի վնասվածությունից պայմանավորված ռիսկի գնահատում	17
2.1	Լուծույթի միջավայրի վնասվածություն	17
2.1.1	Լուծույթի միջավայրի վնասվածություն	17
2.1.2	Լուծույթի միջավայրի վնասվածությունից պայմանավորված ռիսկի գնահատում	17
2.1.3	Լուծույթի միջավայրի վնասվածությունից պայմանավորված ռիսկի գնահատում	18
2.2	Սոցիալական և տնտեսական վնասվածություն	18
2.3	Ջրի միջավայրի վնասվածությունից պայմանավորված ռիսկի գնահատում	19
2.3.1	Լուծույթի միջավայրի վնասվածություն	21
2.3.2	Լուծույթի միջավայրի վնասվածություն	22

[illegible]

6.3	† ԹԵՆՏԵՐ ԻՆ ԿԴԽՆԵՐԿՆԵՐԿ Դ ԵՆՏԵՆՈՒԾԴ	93
6.3.1	Ջ ԴԻՊՈՒԴԵՐ	93
6.3.2	† ԵՆՏԵՐԵՐՈՒՆԵՐԿՆԵՐԿ ԿԴԽՆԵՐԿ Դ ԵՆՏԵՆՈՒՆԵՐ	94
6.3.3	† ՕՐԵՆԴԵՐԿ ԵՆՏԵՐԵՐԿ ԵՆՏԵՐԿ ԿԴԽՆԵՐԿ	95
6.3.4	† ԵՆՏԵՆՏԵՐԵՐՈՒՆԵՐԿ ԿԴԽՆԵՐԿ ԵՆՏԵՐԿ	96
6.3.5	Ջ ԵՆՏԵՐԵՐԿ ԵՆՏԵՐԿ	97
6.3.6	† ՕՐԵՆԴԵՐԿ ԵՆՏԵՐԵՐԿ ԿԴԽՆԵՐԿ ԵՆՏԵՆՈՒՆԵՐԿ	98
6.3.6.1	† ԵՆՏԵՐԵՐԿ ԵՆՏԵՐԿ ԵՆՏԵՐԿ ԵՆՏԵՐԿ	98
6.3.6.2	† ՕՐԵՆԴԵՐԿ ԿԴԽՆԵՐԿ ԵՆՏԵՆՈՒՆԵՐԿ	98
6.3.7	ԿԴԽՆԵՐԿ ԵՆՏԵՐԿ ԿԴԽՆԵՐԿ	99
6.3.8	ԿԴԽՆԵՐԿ ԵՆՏԵՐԿ	99
6.3.9	† ՕՐԵՆԴԵՐԿ ԵՆՏԵՐԵՐԿ ԵՆՏԵՐԿ ԵՆՏԵՐԿ ԵՆՏԵՐԿ	100
6.4	† ԵՆՏԵՐԵՐՈՒՆԵՐԿ ԵՆՏԵՐԿ	101
6.4.1	† ԵՆՏԵՐԵՐՈՒՆԵՐԿ ԵՆՏԵՐԿ	101
6.4.2	† ԵՆՏԵՐԿ ԵՆՏԵՐԿ	102
6.4.2.1	† ԵՆՏԵՐԿ ԵՆՏԵՐԿ ԵՆՏԵՐԿ	103
6.4.2.2	† ԵՆՏԵՐԿ ԵՆՏԵՐԿ / ԵՆՏԵՐԿ	103

4.0.0.7	ԵՆՏԵՐԵՐԿ ԵՆՏԵՐԿ	105
7.1	† ԵՆՏԵՐԿ ԵՆՏԵՐԿ	106
7.2	† ԵՆՏԵՐԿ ԵՆՏԵՐԿ	107
7.3	ԵՆՏԵՐԿ ԵՆՏԵՐԿ ԵՆՏԵՐԿ ԵՆՏԵՐԿ	107
7.4	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	108
7.4.1	† ԵՆՏԵՐԿ ԵՆՏԵՐԿ ԵՆՏԵՐԿ	110
7.5	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	111
7.5.1	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	112
7.6	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	113
7.7	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	114
7.8	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	117
7.9	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	119
7.9.1	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	121
7.9.2	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	122
7.9.2.1	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	123

4.0.0.8	ԵՆՏԵՐԿ	125
8.1	ԵՆՏԵՐԿ	126
8.2	ԵՆՏԵՐԿ	128
8.2.1	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	130
8.3	ԵՆՏԵՐԿ	132
8.4	ԵՆՏԵՐԿ	134
8.4.1	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	135
8.4.1.1	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	135
8.5	ԵՆՏԵՐԿ	137
8.6	ԵՆՏԵՐԿ	139
8.6.1	ԵՆՏԵՐԿ	139
8.6.2	ԵՆՏԵՐԿ	140
8.6.2.1	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	141
8.6.3	ԵՆՏԵՐԿ	141
8.6.4	ԵՆՏԵՐԿ	141

4.0.0.9	ԵՆՏԵՐԿ	143
9.1	ԵՆՏԵՐԿ ԵՆՏԵՐԿ	144

Национальный институт
стандартизации и метрологии

Œ Ÿ Ž Ź R r ı Œ Š š Œ ſ ſ ſ R C

၎င်းတို့သည် အောက်ဖော်ပြပါအတိုင်း ဖြစ်ပေါ်ခဲ့သည်။

၂၀၁၄-၂၀၁၅

[illegible]

Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

1.2

[illegible][illegible]

Տեսչությունը՝ Bio-Rad ձևի վերաբերյալ տեղեկություններ ստանալու համար կոչվում է հ զմ կոչված:

- [illegible]

1. Ինչպիսիք են Կապիտալի Գնահատման Գործընկերության Սահմանափակ ՍՊԸ-ի և Կապիտալի Գնահատման Գործընկերության Սահմանափակ ՍՊԸ-ի (ՕԴԴԳԳՍԸ) միջև կապերը, որոնք կապված են Կապիտալի Գնահատման Գործընկերության Սահմանափակ ՍՊԸ-ի և Կապիտալի Գնահատման Գործընկերության Սահմանափակ ՍՊԸ-ի (ՕԴԴԳԳՍԸ) միջև կապերի հետ։

- [illegible]

[illegible]

1.3

ԿՈՇՏՆԵՐԻ ԼՍԻՏԵՐՈՒՄ

၂။ ဇဝါနိစ္စ^w ၏ ပုံသဏ္ဌာန်တစ်ခုမှာ၊ ဓယုဂ်မညီ တိုက်ရိုက်ပျော်၊ ဝှက်တိုက်ရိုက် လွယ်လွယ် မှတ်မိနိုင်သည်ကို ရှိပါကောင်းသည် :

[illegible]

1.3.1

၂၄၃၀

7 t6tsqL 9tshqls jdz

1 t6tsqL9tsHqfsjdzjdz tft6qBtst60 IH-1000 w9dzw jlsfhw:

DiaMed GmbH
Pra Rond 23
1785 Cressier FR
z 9 j 2 y 0 t e d w

၂၅၂၆၀၁၆၁၆

1. ISM ts'f jte Otses tesd'q f'tshad O Ldz'j o O j's h'w' o d'o H'j dz'u f'te q'eb tes'o IH-1000, q'h f'tesd' u Lz' u' h' q' j'es
f'n O d's f'n s'w's j'd' d's q'd q' f'te iH O' u' h' q' j'es l'e j's' j' f'n l'es tesd' j.

1 jtsdz0dz

[illegible]

ပုဂံသို့ရောက်ရှိသည့်နေ့တွင်

[illegible][illegible]

საბჭოთა კავშირის დაშლის შემდეგ

នៅក្នុងចន្លោះពេលនេះ គ្រឹះស្ថានស្រាវជ្រាវ និង គ្រឹះស្ថានសិក្សា បានប្រើប្រាស់ប្រព័ន្ធគណនេយ្យ ដើម្បីសម្រេចបាននូវលទ្ធផល ដូចខាងក្រោម៖

[illegible][illegible]

ဟဲဇကုတဲာ် နီၤတဲာ် နီၤတဲာ် နီၤတဲာ်

ᄡᆞᆫ ᄢᆞᄱᆞᄳᆞᄴᆞᄰᆞ ᄦᆞᄳᆞᄲᆞᄶᆞᄳᆞ ᄠᆞᄳᆞᄳᆞᄴᆞ^{1w}, ᄤᆞᄳᆞᄳᆞᄴᆞ ᄠᆞᄲᆞᄳᆞᄴᆞ ᄠᆞᄠᆞᄠᆞᄴᆞ^w ᄵᆞᄬᄢᆞᄳᆞᄴᆞ? ᄦᆞᄳᆞᄴᆞ ᄦᆞᄳᆞ.

ሕጋዊነትና ለጥቅም

[illegible]

1.3.2

7 t6tSH zCIs

7 1601B1516

Փնտեսման արդյունքում օգտագործվեցին IH-1000, ԴԵԽԶՕՏՕՂԻՐ 2 ԴԵՏԳԼՈՏԻՎԱՆԵՐՆԵՐ:

၂၂၃၂၀^w နိုတေဝ, ဝေလ်တေဂ်လ်တေဝဇ်ဇ်^w နေဒ်ဖ်ဝဇ်ဂ်၂ Bio-Rad

[illegible]

- Bio-Rad.**

၂၀၁၆

ԿԵՆԴՐՈՒՄԻ Չ ԵՍԶ ԿՐԵՑՈՒՄԻ ՆԵՐՈՒՄԻ ԵՐԱՆՈՒՄԻ:

t j0Gjသ်း, 0၁ststedfL၅၀သ်း r ၃ န်းဒ်တိုဝ်းရို၂ Bio-Rad

[illegible]

Bio-Rad.

1.3.3 $\int \Theta \, ds \, j_{\Theta} \, d\Omega \, dz \, t \, B \, t \, \Theta \, L \, y \, \Theta$

၂. တို့ကဲ့သို့ နေရာအသစ်သို့ ရောက်ရှိလာသည့် အဖွဲ့အစည်းများသည် ABO/Rh-^W ၀ နေရာအသစ်သို့

[illegible][illegible][illegible]

ԱՅԵՈԼՎՐ ԸՆԵՏՈՐ ԻՆՎԵՍՏԻՄԵՆՏԱԿԱՆ ՄԵՏԻՆՎԵՐՋԱԿՈՆ ՆԻՖԻՏԻՆՏՆՈՒԺՔ՝ ԾՆՋԱԺՈՐՋ 48 ՎՈՐԻՏՈ ՄԵՐՈՋ ԻՆՎԵՏԵՈ. ԸՐՈՋՈՐ ՄԵՆ ՈՐԻՆՋԱԺՈՐՋ 24 ՎՈՐԻՏՈ ՄԵՐՈՋ ԻՆՎԵՏԵՈ ԵՄԵՈԼՎԱՆՈ ՆԻՖԻՏԻՆՏՆՈՒԺՈՐՋ ԽՈՒՈՒՈՐ ԸՆՎԵՏԵՈԼՎԱՆՈ ՈՐ ԵՐՈՐՆ ՈՐՄԵՆՈՋԱԿԱՆ, ԵՄԵՈԼՎՐ ԽՏՈՋԱՈՐ ԵՐԼՐ ՄԵՆՈՋ ԻՆՈՐ ՈՐՆՈՋԱԺՈՐՋ ՄԵՐՈՐ ԽԵՄՊԵՐԱՏՈՐ 2-8°C.

[illegible][illegible]

မနုဇာတိသမ္ဘာတော်

[illegible]

1.4 $\mathcal{U}(\mathfrak{g})$ 上的 $\mathfrak{g}$ -模

I dɔfɪnɪst^w h jɔz tɛlɛtɓetɪhɪnɪsə j qɪŋfɪtsɔz j lɛ b sɪŋ^w ŋɔz j hɪz b h q j ətɒdɔdɔz r ɪnsɔdɔz j ətɓɛts ts w tɛdɔz j dɔw:

1.4.1 $\int \frac{1}{x} dx = \ln|x| + C$

[illegible]

၇ တေဝုသုဒ္ဓါး:

- 1** 1 j t e o r 2 1 d j d j d s.
- 2** 1 l s t e t s 2 1 d j d j d s.
- 3** q l s. H.

1.4.7 Դժվարությունների լուծումը

Դժվարությունների լուծումը կապված է ջրի և օդի մաքրման համակարգի, կառուցվածքի և օգտագործվող շինանյութերի չորհրդակցության հետ:

Կառուցվածքի տեսակ	Դժվարությունների լուծում	Գնահատման մեթոդ

1.4.8 Ջրի և օդի մաքրման համակարգի

Կառուցվածքի և օգտագործվող շինանյութերի չորհրդակցության հետ լուծումը կապված է ջրի և օդի մաքրման համակարգի, կառուցվածքի և օգտագործվող շինանյութերի չորհրդակցության հետ:

1.4.9 Դժվարությունների լուծումը

Դժվարությունների լուծումը կապված է ջրի և օդի մաքրման համակարգի, կառուցվածքի և օգտագործվող շինանյութերի չորհրդակցության հետ: Կառուցվածքի և օգտագործվող շինանյութերի չորհրդակցության հետ լուծումը կապված է ջրի և օդի մաքրման համակարգի, կառուցվածքի և օգտագործվող շինանյութերի չորհրդակցության հետ:

1.4.9.1 Դժվարությունների լուծումը

- Դժվարությունների լուծումը կապված է ջրի և օդի մաքրման համակարգի, կառուցվածքի և օգտագործվող շինանյութերի չորհրդակցության հետ:
- Կառուցվածքի և օգտագործվող շինանյութերի չորհրդակցության հետ լուծումը կապված է ջրի և օդի մաքրման համակարգի, կառուցվածքի և օգտագործվող շինանյութերի չորհրդակցության հետ:

1.4.9.2 Դժվարությունների լուծումը

- Դժվարությունների լուծումը կապված է ջրի և օդի մաքրման համակարգի, կառուցվածքի և օգտագործվող շինանյութերի չորհրդակցության հետ:
- Կառուցվածքի և օգտագործվող շինանյութերի չորհրդակցության հետ լուծումը կապված է ջրի և օդի մաքրման համակարգի, կառուցվածքի և օգտագործվող շինանյութերի չորհրդակցության հետ:

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20
 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40

མཚན་ཉིད་ཅེས་པའི་མཚམས་ཀྱི་མཚམས་ཀྱི་

[illegible]

2.1 1945-1949

2.1.1 1 tsqdzqfj r

- 1. 2. 3. 4. 5. 6. 7. 8. 9. 10. 11. 12. 13. 14. 15. 16. 17. 18. 19. 20. 21. 22. 23. 24. 25. 26. 27. 28. 29. 30. 31. 32. 33. 34. 35. 36. 37. 38. 39. 40. 41. 42. 43. 44. 45. 46. 47. 48. 49. 50. 51. 52. 53. 54. 55. 56. 57. 58. 59. 60. 61. 62. 63. 64. 65. 66. 67. 68. 69. 70. 71. 72. 73. 74. 75. 76. 77. 78. 79. 80. 81. 82. 83. 84. 85. 86. 87. 88. 89. 90. 91. 92. 93. 94. 95. 96. 97. 98. 99. 100. 101. 102. 103. 104. 105. 106. 107. 108. 109. 110. 111. 112. 113. 114. 115. 116. 117. 118. 119. 120. 121. 122. 123. 124. 125. 126. 127. 128. 129. 130. 131. 132. 133. 134. 135. 136. 137. 138. 139. 140. 141. 142. 143. 144. 145. 146. 147. 148. 149. 150. 151. 152. 153. 154. 155. 156. 157. 158. 159. 160. 161. 162. 163. 164. 165. 166. 167. 168. 169. 170. 171. 172. 173. 174. 175. 176. 177. 178. 179. 180. 181. 182. 183. 184. 185. 186. 187. 188. 189. 190. 191. 192. 193. 194. 195. 196. 197. 198. 199. 200. 201. 202. 203. 204. 205. 206. 207. 208. 209. 210. 211. 212. 213. 214. 215. 216. 217. 218. 219. 220. 221. 222. 223. 224. 225. 226. 227. 228. 229. 230. 231. 232. 233. 234. 235. 236. 237. 238. 239. 240. 241. 242. 243. 244. 245. 246. 247. 248. 249. 250. 251. 252. 253. 254. 255. 256. 257. 258. 259. 260. 261. 262. 263. 264. 265. 266. 267. 268. 269. 270. 271. 272. 273. 274. 275. 276. 277. 278. 279. 280. 281. 282. 283. 284. 285. 286. 287. 288. 289. 290. 291. 292. 293. 294. 295. 296. 297. 298. 299. 300. 301. 302. 303. 304. 305. 306. 307. 308. 309. 310. 311. 312. 313. 314. 315. 316. 317. 318. 319. 320. 321. 322. 323. 324. 325. 326. 327. 328. 329. 330. 331. 332. 333. 334. 335. 336. 337. 338. 339. 340. 341. 342. 343. 344. 345. 346. 347. 348. 349. 350. 351. 352. 353. 354. 355. 356. 357. 358. 359. 360. 361. 362. 363. 364. 365. 366. 367. 368. 369. 370. 371. 372. 373. 374. 375. 376. 377. 378. 379. 380. 381. 382. 383. 384. 385. 386. 387. 388. 389. 390. 391. 392. 393. 394. 395. 396. 397. 398. 399. 400. 401. 402. 403. 404. 405. 406. 407. 408. 409. 410. 411. 412. 413. 414. 415. 416. 417. 418. 419. 420. 421. 422. 423. 424. 425. 426. 427. 428. 429. 430. 431. 432. 433. 434. 435. 436. 437. 438. 439. 440. 441. 442. 443. 444. 445. 446. 447. 448. 449. 450. 451. 452. 453. 454. 455. 456. 457. 458. 459. 460. 461. 462. 463. 464. 465. 466. 467. 468. 469. 470. 471. 472. 473. 474. 475. 476. 477. 478. 479. 480. 481. 482. 483. 484. 485. 486. 487. 488. 489. 490. 491. 492. 493. 494. 495. 496. 497. 498. 499. 500. 501. 502. 503. 504. 505. 506. 507. 508. 509. 510. 511. 512. 513. 514. 515. 516. 517. 518. 519. 520. 521. 522. 523. 524. 525. 526. 527. 528. 529. 530. 531. 532. 533. 534. 535. 536. 537. 538. 539. 540. 541. 542. 543. 544. 545. 546. 547. 548. 549. 550. 551. 552. 553. 554. 555. 556. 557. 558. 559. 560. 561. 562. 563. 564. 565. 566. 567. 568. 569. 570. 571. 572. 573. 574. 575. 576. 577. 578. 579. 580. 581. 582. 583. 584. 585. 586. 587. 588. 589. 590. 591. 592. 593. 594. 595. 596. 597. 598. 599. 600. 601. 602. 603. 604. 605. 606. 607. 608. 609. 610. 611. 612. 613. 614. 615. 616. 617. 618. 619. 620. 621. 622. 623. 624. 625. 626. 627. 628. 629. 630. 631. 632. 633. 634. 635. 636. 637. 638. 639. 640. 641. 642. 643. 644. 645. 646. 647. 648. 649. 650. 651. 652. 653. 654. 655. 656. 657. 658. 659. 660. 661. 662. 663. 664. 665. 666. 667. 668. 669. 670. 671. 672. 673. 674. 675. 676. 677. 678. 679. 680. 681. 682. 683. 684. 685. 686. 687. 688. 689. 690. 691. 692. 693. 694. 695. 696. 697. 698. 699. 700. 701. 702. 703. 704. 705. 706. 707. 708. 709. 710. 711. 712. 713. 714. 715. 716. 717. 718. 719. 720. 721. 722. 723. 724. 725. 726. 727. 728. 729. 730. 731. 732. 733. 734. 735. 736. 737. 738. 739. 740. 741. 742. 743. 744. 745. 746. 747. 748. 749. 750. 751. 752. 753. 754. 755. 756. 757. 758. 759. 760. 761. 762. 763. 764. 765. 766. 767. 768. 769. 770. 771. 772. 773. 774. 775. 776. 777. 778. 779. 780. 781. 782. 783. 784. 785. 786. 787. 788. 789. 790. 791. 792. 793. 794. 795. 796. 797. 798. 799. 800. 801. 802. 803. 804. 805. 806. 807. 808. 809. 810. 811. 812. 813. 814. 815. 816. 817. 818. 819. 820. 821. 822. 823. 824. 825. 826. 827. 828. 829. 830. 831. 832. 833. 834. 835. 836. 837. 838. 839. 840.

1. 1600-1800-1900-2000-2100-2200-2300-2400-2500-2600-2700-2800-2900-3000-3100-3200-3300-3400-3500-3600-3700-3800-3900-4000-4100-4200-4300-4400-4500-4600-4700-4800-4900-5000-5100-5200-5300-5400-5500-5600-5700-5800-5900-6000-6100-6200-6300-6400-6500-6600-6700-6800-6900-7000-7100-7200-7300-7400-7500-7600-7700-7800-7900-8000-8100-8200-8300-8400-8500-8600-8700-8800-8900-9000-9100-9200-9300-9400-9500-9600-9700-9800-9900-10000-10100-10200-10300-10400-10500-10600-10700-10800-10900-11000-11100-11200-11300-11400-11500-11600-11700-11800-11900-12000-12100-12200-12300-12400-12500-12600-12700-12800-12900-13000-13100-13200-13300-13400-13500-13600-13700-13800-13900-14000-14100-14200-14300-14400-14500-14600-14700-14800-14900-15000-15100-15200-15300-15400-15500-15600-15700-15800-15900-16000-16100-16200-16300-16400-16500-16600-16700-16800-16900-17000-17100-17200-17300-17400-17500-17600-17700-17800-17900-18000-18100-18200-18300-18400-18500-18600-18700-18800-18900-19000-19100-19200-19300-19400-19500-19600-19700-19800-19900-20000-20100-20200-20300-20400-20500-20600-20700-20800-20900-21000-21100-21200-21300-21400-21500-21600-21700-21800-21900-22000-22100-22200-22300-22400-22500-22600-22700-22800-22900-23000-23100-23200-23300-23400-23500-23600-23700-23800-23900-24000-24100-24200-24300-24400-24500-24600-24700-24800-24900-25000-25100-25200-25300-25400-25500-25600-25700-25800-25900-26000-26100-26200-26300-26400-26500-26600-26700-26800-26900-27000-27100-27200-27300-27400-27500-27600-27700-27800-27900-28000-28100-28200-28300-28400-28500-28600-28700-28800-28900-29000-29100-29200-29300-29400-29500-29600-29700-29800-29900-30000-30100-30200-30300-30400-30500-30600-30700-30800-30900-31000-31100-31200-31300-31400-31500-31600-31700-31800-31900-32000-32100-32200-32300-32400-32500-32600-32700-32800-32900-33000-33100-33200-33300-33400-33500-33600-33700-33800-33900-34000-34100-34200-34300-34400-34500-34600-34700-34800-34900-35000-35100-35200-35300-35400-35500-35600-35700-35800-35900-36000-36100-36200-36300-36400-36500-36600-36700-36800-36900-37000-37100-37200-37300-37400-37500-37600-37700-37800-37900-38000-38100-38200-38300-38400-38500-38600-38700-38800-38900-39000-39100-39200-39300-39400-39500-39600-39700-39800-39900-40000-40100-40200-40300-40400-40500-40600-40700-40800-40900-41000-41100-41200-41300-41400-41500-41600-41700-41800-41900-42000-42100-42200-42300-42400-42500-42600-42700-42800-42900-43000-43100-43200-43300-43400-43500-43600-43700-43800-43900-44000-44100-44200-44300-44400-44500-44600-44700-44800-44900-45000-45100-45200-45300-45400-45500-45600-45700-45800-45900-46000-46100-46200-46300-46400-46500-46600-46700-46800-46900-47000-47100-47200-47300-47400-47500-47600-47700-47800-47900-48000-48100-48200-48300-48400-48500-48600-48700-48800-48900-49000-49100-49200-49300-49400-49500-49600-49700-49800-49900-50000-50100-50200-50300-50400-50500-50600-50700-50800-50900-51000-51100-51200-51300-51400-51500-51600-51700-51800-51900-52000-52100-52200-52300-52400-52500-52600-52700-52800-52900-53000-53100-53200-53300-53400-53500-53600-53700-53800-53900-54000-54100-54200-54300-54400-54500-54600-54700-54800-54900-55000-55100-55200-55300-55400-55500-55600-55700-55800-55900-56000-56100-56200-56300-56400-56500-56600-56700-56800-56900-57000-57100-57200-57300-57400-57500-57600-57700-57800-57900-58000-58100-58200-58300-58400-58500-58600-58700-58800-58900-59000-59100-59200-59300-59400-59500-59600-59700-59800-59900-60000-60100-60200-60300-60400-60500-60600-60700-60800-60900-61000-61100-61200-61300-61400-61500-61600-61700-61800-61900-62000-62100-62200-62300-62400-62500-62600-62700-62800-62900-63000-63100-63200-63300-63400-63500-63600-63700-63800-63900-64000-64100-64200-64300-64400-64500-64600-64700-64800-64900-65000-65100-65200-65300-65400-65500-65600-65700-65800-65900-66000-66100-66200-66300-66400-66500-66600-66700-66800-66900-67000-67100-67200-67300-67400-67500-67600-67700-67800-67900-68000-68100-68200-68300-68400-68500-68600-68700-68800-68900-69000-69100-69200-69300-69400-69500-69600-69700-69800-69900-70000-70100-70200-70300-70400-70500-70600-70700-70800-70900-71000-71100-7120

2.1.2 ၂၃၀၂ ခုနှစ် ဝန်ထမ်းစာရင်းနှင့် ဖိတ်ခေါ်ခံရသူများ၏ ဝန်ထမ်းစာရင်း

- [illegible]

[illegible][illegible]

- [illegible]

[illegible][illegible]

- [illegible]

2.3 $\mathbb{Z}_B^h$ ላይ ለፍጥነት ምርመራ የሚያስፈልግ የ $\mathbb{Z}_B$ ላይ ለፍጥነት ምርመራ የሚያስፈልግ የ $\mathbb{Z}_B$ ላይ ለፍጥነት ምርመራ የሚያስፈልግ

- [illegible]

- [illegible]

- [illegible]

2.3.1

[illegible][illegible][illegible][illegible]

| រូបសញ្ញា | និយមន័យ | សេចក្តីពន្យល់ |
|----------|--|--|
| | គ្រោះថ្នាក់ទូទៅ
ឱ្យប្រុងប្រយ័ត្ន
អំពីគ្រោះថ្នាក់
អនាម័យ | - |
| | គ្រោះថ្នាក់ជំងឺ
ឬប្រព័ន្ធភ័យ | ប្រព័ន្ធភ័យ
គ្រោះថ្នាក់ជំងឺ
ឬប្រព័ន្ធភ័យ
គ្រោះថ្នាក់ជំងឺ
ឬប្រព័ន្ធភ័យ |
| | គ្រោះថ្នាក់អគ្គិសនី
ឬគ្រោះថ្នាក់អគ្គិសនី | អគ្គិសនី
គ្រោះថ្នាក់អគ្គិសនី
ឬគ្រោះថ្នាក់អគ្គិសនី |
| | គ្រោះថ្នាក់ភ្លើង
ឬគ្រោះថ្នាក់ភ្លើង | ភ្លើង
គ្រោះថ្នាក់ភ្លើង
ឬគ្រោះថ្នាក់ភ្លើង |
| | គ្រោះថ្នាក់ផ្ទុះ
ឬគ្រោះថ្នាក់ផ្ទុះ | ផ្ទុះ
គ្រោះថ្នាក់ផ្ទុះ
ឬគ្រោះថ្នាក់ផ្ទុះ |
| | គ្រោះថ្នាក់ដង្ហើម
ឬគ្រោះថ្នាក់ដង្ហើម | ដង្ហើម
គ្រោះថ្នាក់ដង្ហើម
ឬគ្រោះថ្នាក់ដង្ហើម |
| | គ្រោះថ្នាក់កំដៅ
ឬគ្រោះថ្នាក់កំដៅ | កំដៅ
គ្រោះថ្នាក់កំដៅ
ឬគ្រោះថ្នាក់កំដៅ |
| | គ្រោះថ្នាក់កាត់
ឬគ្រោះថ្នាក់កាត់ | កាត់
គ្រោះថ្នាក់កាត់
ឬគ្រោះថ្នាក់កាត់ |
| | គ្រោះថ្នាក់ទូទៅ
ឬគ្រោះថ្នាក់ទូទៅ | គ្រោះថ្នាក់ទូទៅ
គ្រោះថ្នាក់ទូទៅ
ឬគ្រោះថ្នាក់ទូទៅ |

| រូបសញ្ញា | ព័ត៌មានបន្ថែម | សេចក្តីពន្យល់ |
|---|---------------------|--|
|  | អនាម័យ មិន អនុញ្ញាត | គ្រប់គ្រង ការប្រើប្រាស់ ឧបករណ៍ ឬ វត្ថុដែល អាច បង្ក ភ្លើង បាន បាន ។
គ្រប់គ្រង ការប្រើប្រាស់ ឧបករណ៍ ឬ វត្ថុដែល អាច បង្ក ភ្លើង បាន បាន ។ |
|  | មិន អនុញ្ញាត បារី | គ្រប់គ្រង ការប្រើប្រាស់ ឧបករណ៍ ឬ វត្ថុដែល អាច បង្ក ភ្លើង បាន បាន ។
គ្រប់គ្រង ការប្រើប្រាស់ ឧបករណ៍ ឬ វត្ថុដែល អាច បង្ក ភ្លើង បាន បាន ។ |
|  | មិន អនុញ្ញាត ភ្លើង | គ្រប់គ្រង ការប្រើប្រាស់ ឧបករណ៍ ឬ វត្ថុដែល អាច បង្ក ភ្លើង បាន បាន ។ |
|  | មិន អនុញ្ញាត ភ្លើង | គ្រប់គ្រង ការប្រើប្រាស់ ឧបករណ៍ ឬ វត្ថុដែល អាច បង្ក ភ្លើង បាន បាន ។ |
|  | មិន អនុញ្ញាត ភ្លើង | - |

2.3.2

၂၀၁၆ ခုနှစ် ဇူလိုင်လ ၁ ရက်နေ့မှ ၃၁ ရက်နေ့အတွင်း

| <p>ՖԵՆԼԱՆՈՎԱԿԱՆԴՅՈՒՆ</p> | <p>ՎԻՃԵԿԻՆԵՐԸ</p> |
|---|---|
|  | <p>Մեռնա՞մեռ? Դեռ չի հասնում, ԵՑԵԹ Ի ՕՆ՝ ԴԻՄ ԻՆ ԵՐԵՎԱՆԵՆԿԱՆԻՆԻՍ՝ Ե՛!</p> |
|  | <p>Վաճաճ՝ ԵՑ ԵՃՕՐԵՎ:</p> |
|  | <p>ՆՃ ԸՕՒՆԵՆՅՈՒՆ՝:</p> |
|  | <p>Դ ՕՏԻՆԴԵՐԵՐԸ՝ ՔՕՄ Դ ԵՐԿԱՐԵՐԵՐԸ՝ ՔՕՄ ԵՄՃԸ՝ ՔՕՄ ԵՃՕՎԱՆԻՆԻՍ՝:</p> |
|  | <p>Դ ՕՏԻՆԴԵՐԵՐԸ՝ ՔՕՄ Դ ԵՐԿԱՐԵՐԸ՝ ՔՕՄ ՆՃԵՐԿՆԵՆԵՆԵՕ:</p> |
|  | <p>ՆՃ ԵՐԵՎԱՆԵՐԵՐԸ՝ ԴԵՐԵՎԵՆԵ ԻԻ-1000 ԴՆՅՈՒՆԻ՝ ԻՆԵՆԵՆԵ ԴԵՐԵՎԵՆԵՕ ԻԻ-1000:</p> |

$\int \nabla^2 r r \nabla t r r \Theta v \Phi J \int r \nabla R r \int$
 $\nabla t R J \int t r$

၂။ ဟိုဝဏ္ဏတို့မှ ငွေဝေခံ၍ ဖြစ်ရသော ရန်ကုန်တော်ဝင် နယ်စပ်တို့တွင် ဖြစ်ရသော။


l o j n j o z o j

ရည်မှန်းချက်အရ လူမှုစီးပွားဝန်ကြီးဌာနမှ နယ်လုပ်ငန်းများကို လုပ်ငန်းခွင်နှင့် လုပ်ငန်းသစ်များ ဖန်တီးပေးရန် အားပေးဆောင်ရွက်ပေးမည် ဟု ဝန်ကြီးက ပြောကြားခဲ့သည်။

IH-1000 | 23

3.2 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ

3.2.1 ԼԵՎ ԻՊԻՏ ԾՈՒՆ



ՏՐԱՆԶԻՐ. 1. ԼԵՎ ԻՊԻՏ ԾՈՒՆ IH-1000 - ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ

- 1 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ
- 2 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ
- 3 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ
- 4 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ
- 5 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ
- 6 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ
- 7 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ
- 8 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ
- 9 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ
- 10 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ
- 11 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ
- 12 ԼԵՎ ԻՊԻՏ ԾՈՒՆ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ ԿՐԵՐԵՆԵՐ

| Icon | Meaning |
|---|--|
|  | The result was transmitted to the LIS. It must still be printed. |
|  | The result was transmitted to the LIS. |
| | The result was validated.
See chapter Additional Validation on page 131 . |
|  | The result was validated. It must still be printed. |
|  | The result was validated and transmitted to the LIS. It must still be printed. |
|  | The result was validated and transmitted to the LIS. |
|  | The result was rejected.
This symbol is only visible if the option Show deleted tests in work list respectively Show rejected in Result list is active. |
|  | The result was temporary rejected. |
|  | The result was excluded from the previous result check. |
|  | The test is processed with the status Priority. |

4.6.5 Display Variants of Results

The number preceding the test description on the right-hand side of Work list and Results windows identifies the number of the instrument on which the test was done.

Specific results are highlighted in four different ways:

| | | | | |
|---|--|---|------------------|-----------------|
| 1 | The displayed result has not been edited.
<i>Display: black, normal.</i> | | | |
| | 2010150401 | 1: Bloodgr. + Rev.gr.:A,B,AB,D/A1,B+C,c,E,e,K,CTL (5041+5011) | 15/04/2010 10:30 | IH-1000 0100011 |
| 2 | The displayed result has been edited, e.g. a reaction intensity was changed.
<i>Display: black, bold.</i> | | | |
| | 201015401 | 1: Ab. screening: I,II,III (IAT) (5053) | 15/04/2010 10:24 | IH-1000 0100011 |
| 3 | The displayed result has not been edited. The result contains a special feature (such as a positive ABS, positive autocontrol, incompatible cross match, etc.)
<i>Display: red, normal.</i> | | | |
| | 2010150419 | 1: Ab.screening: I,II,III,IV,AC (IAT) (5053) | 15/04/2010 12:10 | IH-1000 0100011 |
| 4 | The displayed result has been edited and has special features.
<i>Display: red, bold.</i> | | | |
| | 201015401 | 1: Ab. screening: I,II,III (IAT) (5053) | 15/04/2010 10:24 | IH-1000 0100011 |

4.7 Processing Test Requests

If IH-Com is connected to a LIS, requests are directly sent to the request list in the Work list screen. If patient data are attached, the patients are created in the Patients screen. It is possible to create requests manually at any time.

4.7.1 Manual Request with Patient Data

4.7.1.1 Manual Request for New Patients

A In the Patients window select the button New patient... to create a request for a patient unknown to IH-Com.




Figure 51.

A new request may also be created by selecting in the Patients list area to open the context menu then selecting New patient....

B From the Work list and Results windows select the context menu New request (except QC) in the patient list.




Figure 52.






Figure 53. New request window

| | |
|---|--|
| A | Enter the patient data (1). |
| B | Select the Sender (2) and Cost center (3) from the selection list or enter them manually. |
| C | <p>Enter the Sample number (4) and, if necessary, the case number.</p> <p>If required:</p> <ul style="list-style-type: none"> • Change the blood sample collection time (5) (Default time changes after every hour). • Select the option Priority (7) to define a request as emergency |
| D | If required, print any number of sample labels showing the sample number, first and last names, date of birth and a barcode. |
| E | <p>Select the required number under Sample labels (6).</p> <p><i>Contact Bio-Rad technical support to set a default number of labels to print.</i></p> |
| F | <p>Select the test by placing a check mark to the check box beside the test.</p> <ul style="list-style-type: none"> • If test groups (10) were defined under Configuration / Common / Test profiles, they can also be selected. • If test groups for unknown samples were defined, they are suggested immediately when the Order window is opened. • If cross match tests are requested, enter the required number of blood products. |

- G  Entries are saved, the request is created and a blank Order window for the next patient opens when the Save and next patient button (9) is selected.
- With the button Save (8) the entries will be saved, the request will be created and the input mask will close.

4.7.1.2 Manual Requests for Known Patients

- A  From the Patients window select the button New request for [Last name], [First name] to create a request for a patient known to IH-Com.
- A new request may also be created by selecting New request... in the patient list via context menu.

Č Patient data (1) will already be filled in by the system.







Figure 54. New request window

- B Select the sender (2) from the selection list or input it manually.
- If necessary, enter the case number.

| | | |
|---|---|--|
| C | | Enter the sample number (4). |
| D |  | If required: <ul style="list-style-type: none"> • Change the sample collection time (5). • Select the option Priority (7) to define a request as emergency. |
| E |  | Select the required tests (10) and enter the required number of blood products (11).

If there are defined test profiles (12), they can be used for the request.

If there are defined test groups for known patients, they are automatically suggested. |
| F |  | With the button Save and next patient (9) the request is created and a blank request window for a new patient is opened.

With the button Save (8) the entries will be saved, the request will be created and the input mask will close. |

4.7.2 Manual Request without Patient Data

Requests without patient data can be created in Work list and Results windows.



| | |
|---|--|
| <p>A</p>  | <p>Select the context menu Quick request (QC or without patient data) in the patient list.</p> <p><i>This request can also be used for quality controls.</i></p> |
| <p>B</p>  | <p>Select the required test or the required test group from the Test list.</p> |

Figure 55. New request window: Select Test

- C If a cross match test is selected, enter the number of units to cross match in the Number of crossmatches field (1).
If necessary, select the option Priority (3) to define a request as emergency.





Figure 56. New request window: Selecting cross match test

- D  Select <Enter> to move the cursor to the Sample Barcode (2) field.
- E  Enter samples:
- *With a barcode scanner:* Scan the samples for selected test one by one in the Sample (2) scan field.
Ordered samples are shown in the field on the bottom left.
 - *If entering samples manually:* Confirm each sample number input with the Order now (4) button
- F  When all samples are entered, close the window with the button Close (5).
The last setting for Test is saved and automatically proposed again for the next request.




Figure 57. Warning message for assigning sample number to patient

Č A warning message indicates, if the sample number was assigned to a patient in the last 24 hours.

If the message is confirmed with Yes, the new order is automatically assigned to this patient.

- k The New requests window can also be accessed from the top right of the Patients window. Select the Request without patient... button.

4.7.3 Assigning Blood Products

For manual requests where donor unit numbers are not transmitted from the connected LIS, the donor units must be assigned in IH-Com Client.




Figure 58. Display of not assigned blood products

Non-assigned donor units are identified by the lack of a donor unit number appearing next to the test name in the Work List window.

If donor units need to be assigned, the button Assign donations and the Work List tab flash in the Work List window.

The button changes from  to .

In the Work list: 

- A  Open donation assignment by selecting the Assign donations button.
This window can also be opened using the Assign donations context menu in the patient list. Right click in the Patients list area to access the context menu.




Figure 59. Assignment of donations

A list of all non-assigned cross match requests is displayed.

The one to be assigned is highlighted in purple.

Scan the barcode into the Donation barcode field or input the number manually.

- C  Confirm with <Enter> or the <<Assign button
- D  Check and confirm your entries with the button Save.

- E  Modify incorrect entries by selecting the corresponding line.
Make the assignment again.
Discard all entries with the Cancel button and close the window.
Č The donor unit numbers are displayed next to the test in the Work list.



Figure 60. Display of assigned tests

4.7.4 Quick Request of Crossmatches

The quick request for crossmatches is only possible in the Work List.

The donor unit numbers are assigned immediately.

- A Select the context menu Quick request XM in the patient list or hold down the CTRL key and select Quick request.





Figure 61.

- B Select the required test from the Test list.
Č Only crossmatch tests are displayed in this window. The cursor is positioned in the Sample field when this window opens.

Figure 62. Select test (Example New crossmatch requests)


| | | |
|---|--|--|
| C | ☐ Priority | If necessary, select the option Priority to define the request as an emergency. The check mark must be set before entering the sample number.
<i>Č The next scanned barcode is automatically recognized as Donation(s).</i> |
| D | | Scan the sample in the Sample field.
<i>Č The next scanned barcode is automatically recognized as Donation(s).</i> |
| E | | Scan all donor unit numbers that should be crossmatched for this test individually.
<i>Č Previously scanned donor units are shown in the bottom field.</i> |
| F | | If the sample and donor unit barcodes are entered manually, confirm the input with <Enter>.
<i>Č The cursor moves to the Donation(s) field.</i> |
| G | Donation in list | Confirm your entries with <Enter> or select the Donation in list button. |
| H | Order list of crossmatches | If all donations are assigned, confirm the request with the button Order list of crossmatches.
<i>Č The input window is cleared and the next sample can be entered.</i> |
| I | Close | If no more tests will be requested, exit the window with the button Close.
<i>Č The last setting for Test is saved and automatically suggested for the next request.</i> |
| J |  <p>Figure 63. Warning message for assigning patient to sample</p> <p><i>Č A warning message indicates if the sample number was assigned to a patient in the last 24 hours.</i></p> | |

K

Figure 64. New crossmatch requests window: Assigning patient and sample

If the message is confirmed with Yes, the new request is assigned to this patient.
If the sample is not from this patient, close the window and make a new request.

b Ensure that the correct sample ID is assigned to the patient.

L  Exit this window and make a new request.

4.7.5 Request More

A

From the patient list area, open the context menu and select Request more to request the tests for the currently selected sample again.




Figure 65. Select Request more

Č The corresponding Order window opens. The sample number and patient data, if applicable, are already entered in the window.

B A special type of Request more is Repeat crossmatch.




Figure 66. Select Repeat crossmatch

If Repeat crossmatch is selected in the Patients screen, a new request for sample and donation numbers is created immediately in the request list.

k Repeat crossmatch is mainly intended for ambiguous crossmatch reactions and for re-examination of transfusion reactions.

4.7.6 Repetition of a Request

It is possible to retransmit a request to the device.

A In the Work List window, open the context menu with a right click in the Patients list area and select Retransmit orders for....




Figure 67. Selecting Retransmit orders for...

B Select all requests of a patient, a sample or an individual test.

Č They will be retransmitted to the instrument.

Tests in progress cannot be transmitted again. A new request must be created for these.

4.7.7 Deleting an Order

A To delete an order, right click in the left area of the Work List window to open the context menu.

Select Delete orders for....

Č Depending on the selection, it is possible to delete all requests for a patient, a sample or an individual test.




Figure 68. Selecting Delete orders for...

- B** To delete all open requests, right click in the left area of the Work List window to open the context menu.
Select Delete whole workload.




Figure 69. Selecting Delete whole workload

- C** Individual tests can also be deleted in the Manual work window.
Right click in the Manual Work window to open the context menu and select Remove test.




Figure 70. Selecting Remove test in Manual work window





Figure 71. Selecting Show deleted tests

Č The requests are only hidden. If Show deleted tests is active, they are still visible.

- k** Requests not processed after 48 hours are hidden and after 7 days deleted.
See chapter [Configuration / Common Settings / User Interface on page 38](#).
These times can be changed by your Bio-Rad technical representative.

- D** Requests with already available results (e.g. pipetted samples) cannot be deleted.
These must be rejected as described in section [Rejecting Pipetted Samples on page 130](#)

4.7.8 Use Multiple Selection for Requests

A  To simultaneously retransmit multiple selections to the instrument or to delete several requests at the same time, right click in the Work list window to open the context menu.

B Select Show multi selection.

Č If the option is activated, check boxes are enabled next to the tests in the Work List area.





Figure 72. Example for display of all samples and tests

Example: If the check box in front of a patient is selected, all samples and single tests belonging to this patient are selected.

C  Select the context menu from the Work List area.
From the context menu, select between the options:

- Select all requests via Show multi selection;
- Retransmit selected orders;
- Delete selected orders.

k Reselect Show multi selection to deactivate the option and revert back to the standard context menu items.

PAGE INTENTIONALLY LEFT BLANK

5 Instruments, Manual Work

The chapter gives detailed information about all automated instruments, readers and manual workplace.

5.1 Operating Semi and Fully Automated Instruments

The following instruments can be interfaced with IH-Com:

- IH-500
- IH-1000
- Swing TwinSampler II

The interfacing of instruments with IH-Com (except IH-500) is via a device-specific instrument connector, which operates entirely in the background. The instrument connector is installed on the instrument PC.

The instrument connector transmits requests from IH-Com to the instrument and results from the instrument to IH-Com.

If the instrument provides status updates then these are also transferred via the connector. This allows the status of a request to be automatically changed in the IH-Com Worklist.

See chapter [Transmission Status on page 64](#).

Example: An instrument is sending a status error message that withdraws the status In progress of a request.

- A  When the instrument communicates with IH-Com via instrument connector, the connector is automatically started after PC power up and is present in the system tray (1).
- Č If all the connections have been made between the connector and IH-Com, the icon has a green border.
- Č If one or more connections are missing, the icon has a red border.




Figure 73. IH-Com Master Database Connector

By selecting the icon the instrument connector can be opened.

Č Connection data and log entries are displayed.

If support is required, data can be requested by your Bio-Rad technical representative.




Figure 74. IH-Com Master Database Connector: tab Assays

Č On the Assays tab there is an overview of tests that are installed on this instrument connector.

Only tests for which an appropriate assay has been installed can be processed on this instrument.

Installation is performed by your Bio-Rad technical representative.

5.1.1 IH-500




Figure 75. IH-Com Client: tab Work list for IH-500

IH-500 communicates directly with IH-Com via query mode - no instrument connector is used.

Requests are sent to IH-Com after the sample barcode has been read on the instrument.

When requests are available in IH-Com, the samples can be placed directly on IH-500 as described in the [IH-500 User Manual](#).

Results are automatically transferred to IH-Com when processing has completed.

The instrument software of IH-500 is located on the instrument itself.

5.1.2 IH-1000




Figure 76. IH-Com Client: tab Work list for IH-1000

- A In the Work List window, as soon as the requests have the Transfer to device status icon, samples may be started on the IH-1000.
- B  Switch to the instrument software.
If the instrument software is not yet running, it is started with the IH-1000 button.
- C Place the samples in the instrument and start processing as described in the [IH-1000 User Manual](#).
- D When sample processing has been completed, the results are automatically sent to IH-Com.

5.1.3 Swing TwinSampler II




Figure 77. IH-Com Client: tab Work list for Swing TwinSampler II

- A

In the Work List window, as soon as the requests have the Transfer to device status icon, samples may be started on the Swing TwinSampler II.
- B



Switch to the instrument software.
If the instrument software is not yet running, it is started with the Swing TwinSampler II button.
- C

This button is only displayed if IH-Com has been directly installed on the Swing TwinSampler II.
- D

Place the samples in the instrument and start processing as described in the [Swing TwinSampler II User Manual](#).
- E

When sample processing has completed, results are automatically sent to IH-Com.
- F

Depending up on the test, they can now be centrifuged after the corresponding incubation time, and read by Saxo ID-Reader II or Banjo ID-Reader.


5.2 Operating Reader Devices

The following readers can be interfaced with IH-Com:

- [Saxo ID-Reader II](#)
- [Banjo ID-Reader](#)

The IH-Com Client only displays reader devices that are connected.

Interfacing of a reader device takes place via an instrument-specific instrument connector, which operates entirely in the background. It is usually installed on the instrument PC.

- A  This instrument connector is automatically started after PC power up.
 It is possible to detect it from the icon in the system tray (1).
Č If all the connections have been made between the connector and IH-Com, the icon has a green border.
Č If one or more connections are missing, the icon has a red border.

- B By selecting the icon the instrument connector can be maximized.




Figure 78. IH-Com Reader Centrifuge Connector

Č Connection data and the log entries are displayed.

If support is required, these data can be requested by your Bio-Rad technical representative.

- k Each instrument must be initialized before use.

5.2.1 Saxo ID-Reader II





Figure 79. IH-Com Client tab Saxo2

| | | |
|---|---|---|
| A | | Turn on the Saxo ID-Reader II at the ON/OFF switch. |
| B | | Make sure the instrument lid is closed. |
| C | | In the IH-Com Client select the view Saxo [Serial number] (1). |
| D | | Initialize Saxo by selecting the button Stop/Initialize/Open lid.
<i>Č The respective status is displayed at the top right of the screen (2).</i> |
| E | | Load the instrument and start a run. |
| F |  | According to the required action, select e.g. the button Centrifuge and read.
<i>Č If the warming up time of the Saxo has not yet been completed, IH-Com indicates that with an appropriate message.</i>

<div>Reading ...
Handle all available cards: reading
Warmup process until 09:24 is not completed yet, reading is not possible!</div>
Figure 80.

<i>The cards are centrifuged but not yet not read.</i> |
| G |  | It is possible to read the cards later with Read only if the warming up phase has been completed. |

| | | |
|---|---|--|
| H | | Results are displayed under the tabs Results and Batch Verification.
Verify them as described under Results on page 97 . |
| K | | If an already read gel card is read again by Saxo, reaction strengths will be overwritten, if results are not yet accepted and if the reactions strengths differ from the previous reading result. For each overwritten reaction a comment is created automatically to the well containing the previous reaction strength. This function can be disabled by your Bio-Rad technical representative. |
| I |  | While the instrument is switched on, it is possible to shut down or restart IH-Com.
If the instrument or the instrument PC is to be switched off without IH-Com being shut down, select the button Switch off first. |
| J | | To operate and calibrate Saxo, use the menu Reader devices/Saxo/Calibrate Saxo.
<i>Normal user rights are required.</i> |
| K | | For service settings use the menu Reader devices/Saxo/Saxo settings.
<i>Administrator rights are required.</i> |
| L | | The respective status or action of Saxo is documented in the area below the buttons.
The display can be emptied with the context menu Clear display. |



5.2.2 Banjo ID-Reader




Figure 81. IH-Com Client tab Banjo

- A Turn on the Banjo ID-Reader at the ON/OFF switch.
- B In the IH-Com Client select the view Banjo [Serial number] (1).
- C Initialize the Banjo by selecting any button.
Č The respective state is displayed at the top right of the screen (2).
- D  It is now possible to insert a card in the Banjo ID-Reader.
 To open and close the instrument always use the appropriate button.
 The barcode of the gel card must be on the back for reading.
 Selecting the button Read card starts the reading the gel card.
*Č Reading of the card is displayed in the log screen of the screen view.
 When the card has been read successfully, the instrument opens automatically.*

```
Initialization ended (state = Stopped)
closing cover ...
Reading ...
Analyse id-card
closing cover ...
Reading ...
Analyse id-card
```

Figure 82.

| | |
|---|--|
| E | Results are displayed under the tabs Results and Batch Verification.
Verify them as described under Results on page 97 . |
| k | If an already read gel card is read again by Banjo, reaction strengths will be overwritten, if results are not yet accepted and if the reactions strengths differ from the previous reading result.
For each overwritten reaction a comment is created automatically to the well containing the previous reaction strength. This function can be disabled by your Bio-Rad technical representative. |
| F |  While the instrument is switched on, it is possible to shut down or restart IH-Com.
If the instrument or the instrument PC is to be switched off without IH-Com being shut down, it is essential to first select the button Switch on/off. |
| G | To operate the Banjo ID-Reader normal user rights are required.
The Banjo settings can be accessed with administrator rights using the menu Reader devices/Banjo/Banjo settings. |

5.3 Manual Work

- b** Manual work must be performed with extreme care. Activities should only be performed by professionally qualified staff, who have been trained on working with IH-Com. Apply the four eyes principle.

Open requests are displayed in the Manual work screen. As soon as a request receives the in progress feedback from an instrument, it is hidden.




Figure 83. Manual work window: Display of all open requests

- A** Select the tests that are to be processed manually.
All tests in the Do assign column are selected by default.
- B** Remove the check mark from tests that are not to be processed.
All check marks can be removed at once, with the context menu Remove selection.
If only high priority requests should be processed, select Select Priority Only with the context menu

| Do assign | Sample / Patient | Test | State | Sample time |
|-------------------------------------|---|---|-----------------|------------------|
| ☒ | 2014001 Amberg, Marie Luise *15/09/1945 (F) [ID0815] | Ab screening: I,II,III (IAT) (5053) | Test requested: | 03/11/2014 12:00 |
| ☒ | 2014002 Black, Gary *10/10/1988 [75374] | Ab screening: I,II,III (IAT) (5053) | Test requested: | 03/11/2014 12:00 |
| ☒ | 2014005 Fisher, Fish *09/09/1988 (F) [7834657343467] | Ab screening: I,II,III (IAT) (5053) | Test requested: | 03/11/2014 12:00 |
| ☒ | 2014001 Amberg, Marie Luise *15/09/1945 (F) [ID0815] | Autocontrol: (IAT) (5053) | Test requested: | 03/11/2014 12:00 |
| ☒ | 2014002 Black, Gary *10/10/1988 [75374] | Autocontrol: (IAT) (5053) | Test requested: | 03/11/2014 12:00 |
| ☒ | 2014005 Fisher, Fish *09/09/1988 (F) [7834657343467] | Autocontrol: (IAT) (5053) | Test requested: | 03/11/2014 12:00 |
| ☒ | 2014006 Forster, Elton *07/07/1977 (M) [37843563478346] | Bloodgr.+Rev.gr.:A,B,AB,D/A1,B+C,c.E.e.Kctl (...) | Test requested: | 03/11/2014 12:00 |
| ☒ | 2014001 Amberg, Marie Luise *15/09/1945 (F) [ID0815] | Bloodgroup: A,B,D Confirmation for Patients (Dia... | Test requested: | 03/11/2014 12:00 |
| ☒ | 2014002 Black, Gary *10/10/1988 [75374] | Bloodgroup: A,B,D Confirmation for Patients (Dia... | Test requested: | 03/11/2014 12:00 |

Figure 84. Sorting tests for a better overview

The list can be sorted by column headers.

5.3.1 Assign Samples

5.3.1.1 Assign Samples by Scanning Sample Barcode

It is also possible to select samples to be assigned by scanning the sample barcode.


If all samples have a checkmark in the Do assign column and a sample barcode is scanned, all checkmarks will be removed and only the scanned sample remains selected.

If at least one request is not selected, selection of the scanned sample will be reversed.

Example:

- 1 All samples in the Do assign column are selected.
- 2 The barcode of sample 001 is scanned.
- 3 Now only sample 001 is selected for assignment and all the other checkmarks have been removed from the Do assign column.
- 4 Now the barcode of sample 002 is scanned.
- 5 Consequently, in the Do assign column both sample 001 and sample 002 are selected.
- 6 If the barcode of sample 002 is now scanned again, the selection for sample 002 is removed and only sample 001 is selected for assignment.


5.3.1.2 This workplace only option



If assays are in use at the manual workplace, which are not available for an instrument (for example: ID-PaGIA Heparin/FP4 Antibody Test), select the option This workplace only to display only requests for these tests.

This option can be used if multiple manual workstations are set in the lab.

5.3.2 Assign Gel Cards



Open card assignment using the context menu or the button Assign cards.

Available tests for assignment are displayed in the upper screen.




Figure 85. Assign cards for pipetting window

5.3.2.1 Sorting

It is possible to sort the list by tests or by samples/patients. This sort order only applies within a gel card type.

Example:

For multiple patients, ABS, autocontrol and crossmatch are requested.

These tests are all performed on the same gel card type.

According to the selection it is possible to display the request of that card type in such a way that either:

- all the tests for sample 1 are listed under each other, followed by all the tests for sample 2, etc.
- or first all the ABS tests, followed by all the autocontrols and lastly all the crossmatches.

Sort by test




Figure 86.

5.3.2.2 Documentation of Reagents

- Barcodes must be scanned in the assignment window once a day.
- Only one lot can be scanned per reagent type.
- It is not possible for two lots to be active at the same time.
- For correct lot documentation, in the event of a lot change during the day, the new lot must be scanned as soon as it is going to be used.

Assignment not possible

| State suggested | | | | |
|-----------------|---|--------|---------|--------|
| Reagent | ✓ | Expiry | Checked | QC Lot |
| ID-DiaCell I | | | | |
| ID-DiaCell II | | | | |
| ID-DiaCell III | | | | |

Never scanned reagents are highlighted with a red background.


Assignment possible

| State suggested | | | | |
|-----------------|---|------------|------------------|------------|
| Reagent | ✓ | Expiry | Checked | QC Lot |
| ID-DiaCell I | | 10/11/2014 | 03/11/2014 14... | 06084.82.1 |
| ID-DiaCell II | | 10/11/2014 | 03/11/2014 14... | 06094.82.1 |
| ID-DiaCell III | | 10/11/2014 | 03/11/2014 14... | 06104.82.1 |

Expiry date, date when the reagent was last scanned and whether a valid control exists are displayed for information purposes

If the time period since the last scan has expired, these reagents are highlighted with a yellow background. The period of validity is 24 hours.

If a gel card is scanned before confirming reagents, the following message appears:



All reagents that might be used during the day can be scanned immediately.

If a test has never been used, it may happen that the reagent can only be scanned when that test is to be assigned for the first time.

If a lot changes during the day, only the new lot must be scanned. The other reagents will remain.

5.3.2.3 Gel Card Suggestion

Selected card: LISS/Coombs (Incubation 37°C)

| I | II | III | I | II | III |
|--|----|-----|-----------------------------------|----|-----|
| 2014001 | | | 2014002 | | |
| Amberg, Marie Luise *1509/1945 (F) [00815] | | | Black, Garry *10/10/1988 [75/374] | | |
| | | | Incubation 37°C | | |

Figure 87.

Č According to the selected test in the list of tests to be assigned, the first proposal is displayed in the Selected card screen.

Č Under that all the reagents required for that test are displayed.

k In case of gel card types on which various tests can be performed, a distinction is also made for incubation temperature and incubation time.

Tests with different incubation temperature or incubation time cannot be pipetted onto the same gel card together. Only when the gel card has been read, vacant wells can be used for a different incubation temperature or incubation time.

5.3.2.4 Options of Assigning Gel Cards

Option 1

A Select a test in the list of gel card types for a sample.

Č In the Selected card screen a proposal is displayed for the gel card assignment.

B If the proposal is acceptable, scan a gel card of the required gel card type.

A label with the gel card assignment is printed out.

The status is changed from not assigned to pipetted.




Figure 88. Assigned gel cards

Č In the Assigned cards screen the assigned gel card is displayed with the sample data.

C

| | | | | | |
|---|------|-------------------------------|--|----|-----|
| Selected card: | | LISS/Coombs (Incubation 37°C) | | | |
| | | 5053355011504209909 | | | |
| AC | K006 | K007 | I | II | III |
| 2014005 | | | 2014001 | | |
| Fischer, Fish *09/09/1988 (F) [7634657343467] | | | Amberg, Marie Luise *15/09/1945 (F) [Muller] | | |
| 5053355011504209909 | | | Incubation 37°C | | |

Figure 89. Selected gel card

Č If the well is already pipetted on the gel card, they are highlighted with an orange background and the sample data.

The pipetted status is displayed.


In such a case it is necessary to scan the gel card again in order to achieve an assignment:

- first scan shows the gel card proposal
- second scan accepts it

Option 2

| | |
|---|--|
| A | <p>Irrespective of what is selected in the list of gel card types, scan a gel card of the type that is to be assigned.</p> <p>The next test with the scanned gel card type is automatically selected and a gel card assignment proposed in the Selected cards screen.</p> |
| B | <p>If the proposal is acceptable, scan the gel card again.</p> <p>A label with gel card assignment is printed out and the status is changed from not assigned to pipetted.</p> <p><i>Č In the Assigned cards screen the assigned gel card is displayed with the sample data.</i></p> |
| C | <p>In both options the next test that has to be assigned is automatically selected for the gel card type just scanned.</p> <p><i>Č If there is no other open test, an appropriate message appears.</i></p> |

Further steps

| | |
|----------|---|
| D | <p>Apply the printed label with gel card assignment to the correct side of the gel card (the side opposite that of the barcode) in order to uniquely identify the wells.</p> |
| b | <p>It is essential to adhere to this procedure.</p> <p>A different sequence (first pipette, then assign) can lead to incorrect results.</p> |
| E | <p>If the gel card must be assigned differently than proposed by IH-Com, it is possible by using the context menu.</p> <p>Individual samples can be removed from the gel card by the context menu, which can be opened directly on the gel card proposal.</p> |
| F | <p>With the context menu in the list of tests to be assigned it is possible to assign tests to a certain sample on the gel card with Sample XXX at position Y as long as the test does not occupy all positions on the gel card.</p> <p><i>Example: Rhesus subtype gel cards</i></p>  <p>Figure 90. Selecting Sample XXX at position Y</p> |

- G With New proposal for all tests on this card, it is possible to generate the IH-Com proposal again.




Figure 91. Selecting New proposal for all tests on this card

- k Gel cards that have been processed on one instrument with at least a status: pipetted are taken into account for manual work. The pipetted wells are flagged as used.
- k Gel cards pipetted without the use of IH-Com must not be used. In this case, IH-Com has no information about the wells already used.

- H Label printing can proceed again.
Scan the gel card barcode twice.






Figure 92. Printing label again

5.3.2.5 Remove a Test from a Gel Card

In rare cases it is necessary to delete the assignment for a gel card.

Example:



| Test | State |
|---|----------------------------|
| Ab screening: I,II,III (IAT) (5053) | Test requested: |
| Autocontrol: (IAT) (5053) | Test requested: |
| Bloodgroup: A,B,D Confirmation for Patients (Dia... | Test requested: (pipetted) |
| Crossmatch: IAT(5053) x K004 | XM requested: |

Figure 93. Tests with status Test requested

A gel card has been assigned and not yet pipetted, but in the meantime the test is no longer required.

This gel card can be used for a different sample by opening the gel card assignment, scanning the barcode of the gel card and executing Remove all assignments for this card in the context menu directly on the gel card suggestion.

Subsequently a related message indicates that the gel card assignment is closed.

The status Test requested (1) indicates tests that were removed from the gel card.

These requests can be deleted if necessary. The gel card can now be assigned differently.

| | | |
|---|--------------------|---|
| A | | <p>All assigned gel cards can be pipetted and processed. Saxo ID-Reader or Banjo ID-Reader can be used.</p> <p>See section Saxo ID-Reader II on page 85 or Banjo ID-Reader on page 87</p> <p>In the Manual work screen the manually processed gel cards can be directly read, as in section Reading Results Manually on page 115.</p> |
| B | Assign cards | Select the option Show completed results too. |
| | Assign microplates | |
| C | Select all | Select the required test |
| | Select sample | |
| | Remove selection | |
| D | Read manually | Open the result input window with the context menu Read manually or the button Read manually. |
| | Request more | |

5.3.2.6 Adapting the Gel Card Label to the Label Printer



If in the Manual work screen the Assign cards for pipetting window is opened with administrator rights, a test label can be printed to adapt the printer settings.

A test label is printed out with the currently displayed card proposal and a fictitious card number.

6 Results

The chapter comprehends an explanation about handling and procedures of processing results in IH-Com.

When the first result is transmitted from the instrument to IH-Com (or if it has pipetted status), the sample ID is displayed in the Results window.

Depending upon the settings defined, results may require verification by the user.

See chapter [Configuration / Common Settings / Verification on page 40](#).

Results which have verification missing status are additionally displayed under tab Batch Verification.

- b** Before final test interpretation, results have to be validated by trained operators. This result verification process includes visually reviewing the images and editing (if necessary) the test interpretation. Editing is not limited to indeterminate results; initially positive or negative reactions might also be edited after visual review. It is recommended to activate the 4-eyes-principle. However, this is left to the discretion of the facility.

A

Open Result Verification by:

- using the context menu Edit result;
- selecting the button Edit result;
- double-clicking the sample number or the name in the patient list in the Results window.




Figure 94. Selecting Edit result via context menu

Refer to the sections [Result Verification Window on page 101](#) and [Changing and Confirming Results on page 103](#) for more details on editing results.

Č All tests are displayed in the Search cards window, showing gel card position and sample data.




Figure 95. Search cards window with active button Show result

b If using the search cards functionality, continue to scan the gel cards which have been segregated for second reading in the main window of the IH-Com Client.

B



Select the required test and open Result Verification with button Show result.



If the Show single test immediately option is active, the result validation of the corresponding sample is opened immediately. This will occur only when there is a single result on the scanned gel card.


C The Search cards window can be opened using the Work list / Search card menu.

By opening Result Verification and selecting the button Edit results and when the Automatic result run option is activated, the next patient or sample with open results is displayed immediately after saving the Result Verification. Result verification starts with the selected sample and is performed in sequential order to the end of the sample list.

With Automatic result run the tabs of the Results screen are taken into account, i.e. only the result validation of the samples which corresponds to the current selection is displayed.

k The Automatic result run only applies to results that require verification. For results under the tab Extra check (results to be validated) it is not enabled. See chapter [Additional Validation on page 131](#).

D  When all results have been verified, they are automatically sent to the LIS.

E  It is also possible to send temporary results to the LIS manually:

Select the context menu Send results in the Results window.

Temporary results are marked differently and are identified as such in the LIS.

6.1 Selection Buttons and Filter in the Results Window




Figure 96. Results window: Buttons

| | |
|---------------|---|
| × All | <p>This view is activated by default when the program is started.</p> <p>The following are displayed:</p> <ul style="list-style-type: none">• all results of the current day• results that are not validated regardless of when they were processed• pipetted results for which no reaction intensities were read |
| 1 To Read | <p>All results that have not yet been accepted and with reaction strengths that have not yet been saved are displayed here.</p> |
| 2 To Read 2nd | <p>Results that have been temporarily accepted (e.g. automatically accepted results), but still have irregularities (e.g. discrepancy between blood group and serum crossmatch, missing quality control). Results with reaction strengths that are saved and not accepted.</p> <p>This button is activated only if corresponding results exist.</p> |
| Extra check | <p>Results that have been selected for an extra check in verification. This button is only shown if corresponding results exist.</p> |

The number in brackets on the selection buttons represents the number of samples displayed in the window.




Figure 97.

- 1 If a filter for results is selected under Common Settings, the selection is displayed.
See chapter [Configuration / Common Settings / User Interface on page 38](#).
A click on the selection buttons of the result list displays only results for the selected filter.
- 2 If a Lab location is selected, only results of instruments are displayed that belong to this laboratory location. If an instrument is selected as a filter, only results of this instrument are displayed.

6.2 Result Verification Window

The screenshot shows the 'Result verification - User Admin' window. It contains the following elements:

- 1**: Patient / sample data section at the top left, including patient name, date of birth, and sample information.
- 2**: Audit trail icon (history of changes) next to the patient data.
- 3**: Test time and verification details section.
- 4**: Comment field for the sample and/or the test.
- 5**: Tabs for individual categories of test results: Current results, Crossmatch, Antibody panels, and Other Samples (13 wells).
- 6**: Reaction results section showing gel card images and corresponding results (AC, XM, XM).
- 7**: Result interpretation section with dropdown menus for AB0, Phenotype, Rhesus-D, Kell, ABS, DAT, and Auto ctrl., along with a 'Recalculate Interpretation' button.
- 8**: Magnified view of the selected well with gel card data and comment for the well.

Figure 98. Result verification window

- 1 Patient and sample data
- 2 Audit trail (history of changes)
- 3 Tested by..., time of test, verified by..., validated by...
- 4 Comment for the sample and/or the test
- 5 Tabs for the individual categories of test results
- 6 Reaction results
- 7 Result interpretation
- 8 Magnified view of the selected well with gel card data and comment for the well

k According to the window setting and resolution it is possible that not all wells are shown in this view. If this is the case a scrollbar is displayed.

k In the interpretation section a part of the interpretation can be marked with a superscript 1 (default setting). This means that this part of the interpretation is from a pre-result.

Result Verification is identical for all instruments and for the manual workplace.

Single reaction results are displayed up to a maximum of 5 tabs. Only tabs with reaction results are displayed."

| Tab | Reaction result |
|-----------------|--|
| Current results | All reaction results of the current sample are displayed, if they are not a crossmatch test, a panel test or a titration test |
| Crossmatch | All reaction results from a crossmatch test of the current sample are displayed |
| Antibody panels | All reaction results from a panel test (antibody identification) of the current sample are displayed |
| Titration | All reaction results of a titration test of the current sample are displayed |
| Other samples | All reaction results which do not belong to the current sample but to the same patient are displayed, e.g. pre-results or results detected after the currently marked sample.
The number of existing reaction results from other samples is displayed in brackets.
Images of reaction wells, displayed under this tab, are created in the course of opening the tab.
Related to the number of single reactions loading of the tab can last. In case of a waiting cycle a progress bar is displayed. |

Results under tabs are sorted by test in a sequence that can be modified by your Bio-Rad technical representative. The test sequence is automatically defined at the test import.

If multiple tests are processed for the same sample, it will be sorted by descending testing time.

On tab Other samples sorting takes place by descending sample number, within the same sample number by descending testing time.

Colors representing interpretations

| | |
|--------|--|
| White | Plausible result without edits |
| Green | Result corresponds to the automatic interpretation after editing |
| Orange | Result does not match automatic interpretation after editing |
| Red | Not a plausible result |

k For DAT, Ab Screening and Autocontrol deviations of the representing colors green and orange appear. It depends on whether the test is monospecific or polyspecific.

6.3 Changing and Confirming Results

In the Result verification view results not accepted are marked with the  icon.

A

The reaction strength read by the instrument is displayed under the individual # column.

If your visual assessment of a result differs from that of the reader device, change the reaction strength via the context menu.

Select a column and open the context menu using the right mouse button.




Figure 99. Changing of reaction result via context menu

IH-1000 and IH-500 provide additional well interpretations in case of unclear reactions




Figure 100. Result verification for IH-1000

| | |
|---------|------------------------------------|
| IH-1000 | wF, wR, wP, wD, E, L, W, LIQ or [] |
| IH-500 | wR, wP, E, LIQ |

Other instruments in the event of unclear reactions only issue “?” as reaction strength.

These interpretations are automatically generated and not available for selection.


Č A comment is automatically displayed on the column.

A changed reaction intensity is marked with an asterisk.

B Place the cursor on the corresponding reaction column to display the original instrument result.

A tooltip is opened with following information:

- Test name
- Sample number
- Time of reading
- User
- Rotor position in the instrument (if available)
- Time of first reading
- Current value, etc.



See chapter [Tooltip of a Well](#) on page 108.

C




Figure 101. Audit trail window

All manual changes and original values are recorded in the audit trail.

Open the audit trail window by selecting the  button.

D

To verify results:

Select the context menu on the reaction column.




Figure 102.

E

The following options are available:

- Verify or reject all results of the current view (e.g. tab Current results).#
Rejecting the result prevents validation and transmission to the LIS or;
- Verify each test (e.g. ABD) individually or;
- Accept all results simultaneously by selecting Accept all (All results in the current view will be verified simultaneously.).

✓ A verified result is indicated with .


✗ A rejected result is indicated by a strikethrough.

For crossmatch test results, the preselected treatment of blood sample (transfusion status) can be changed if necessary via the context menu.

Crossmatch samples with negative reactions can be changed from Approve for Transfusion to Block for Transfusion.

Conversely, unreleased samples from Block for Transfusion to Approve for Transfusion can also be changed.

A comment is required to change the transfusion status. An entry to the Daily Journal is automatically added.



b This will release samples for transfusion, even though they have reacted (positive) in the crossmatch test.

F




After all required changes have been made and all results are validated, save your results.

| User rights | Description |
|-------------|--|
| 40 | User can change reaction strengths but cannot verify the results. |
| 50 | User must check and verify results once again and may change reaction strengths again if required. User can change reaction strengths, confirm and release (save) results. |
| 60 | If the appropriate option is enabled, it is possible to perform additional validation, a second verification of results with access rights of at least 60
See chapter Additional Validation on page 131 . |

| | | |
|---|---|---|
| G |  | The individual users are shown in the top portion of the window. |
| H | | Results that have been validated additionally are marked. |
| I |  | The Skip button is available for the automatic result validation.
Skip a sample or cancel result validation by selecting the Cancel button |

6.3.1 Illustrations of Reaction Intensities

IH-Com displays the original images of the readers and automated instruments.



Ideally, the separated column images are exactly the same size as the interpretation fields in Result Verification.

If the barcode is too high or the card is not clean (shown exaggerated in the illustration), the image is reduced.

For this reason it may happen that the reaction columns of a card have different sizes in the Result verification screen.

k IH-Com applies the distance between the last dark pixel of the column edge and the next dark pixel to represent the image of the entire card as individual reaction columns.
Example: At the barcode or at the edge of the next reaction column, the image separation is placed.

k If the original image is not available, due to manually read results or a defective database image, reactions are shown as a graphic.
See chapter [Reading Results Manually on page 115](#).

6.3.2 Color Code for Reactions Illustrations

IH-Com uses the colors of the rectangular area above the individual reaction pictures to identify certain properties.




Figure 103.

| | |
|--------|--|
| Yellow | Results that have not yet been read or verified. |
| Brown | Results: <ul style="list-style-type: none">• with positive reactions in the ABS, self-monitoring or DAT;• with weak reactions (according to assay definition);• whose reaction strength has been modified by a user;• with flags of categories C, D and E from the device (see section Bio-Rad Result Flags on page 22);• with a positive control. |
| Red | Incompatible crossmatch samples. |
| Green | Compatible crossmatch samples. |
| Gray | Accepted, rejected or past results. |

6.3.3 Tooltip of a Well

The tooltip of a well details the steps leading to the result of the single reaction.




Figure 104.

| Section of the tooltip | Description |
|------------------------|--|
| Processing | <p>Displays the time well processing was started. The displayed time is transmitted from the instrument. At the manual workplace this entry corresponds to the time when the gel card has been assigned by scanning the barcode.</p> <p>Displays the user logged in at the instrument where the well was pipetted. If there is no user login at the instrument, instrument and serial number are entered. If pipetting took place manually, the user that has been logged in at the manual workplace at the moment the gel card has been assigned by scanning the barcode, is entered.</p> |
| Well reading | <p>Displays the time of well reading. The displayed time is transmitted from the instrument. At the manual workplace this entry corresponds to the time when the manually entered reaction strength has been saved.</p> <p>Displays the user logged in at the instrument where the well was read. If there is no user login, instrument and serial number are entered.</p> <p>For a reading device the logged in user and the serial number of the reader are entered. For manual reading the user is entered who input the reaction strengths.</p> |
| Machine result | <p>If a reaction strength has been changed manually, the original reaction strength read by the instrument is displayed. If the reaction strength has not been changed no entry is displayed.</p> |
| First reading | <p>The time when the first confirmation of the reaction strength took place.</p> <p>User who confirmed the reaction strength</p> <p>Confirmation of a reaction strength takes place e.g. by batch verification or by selecting Save under Result Verification, without accepting it previously. If confirmation of reaction strengths and accepting takes place within one working step no entry is displayed.</p> |

| Section of the tooltip | Description |
|------------------------|---|
| Last changed | The time when the last modification of the reaction strength took place.
User who initiated the last modification of the reaction strength.
If the reaction strength has not been changed, no entry is displayed. |
| Current value | Current reaction strength |
| Well verification | The time when the well has been verified (accepted).
User who verified (accepted) the well.
If single tests are accepted by batch verification and the overall result of the sample is thereby verified automatically even though the verification differs from the standard value, an entry automatically accepted is displayed. |
| Result verification | Point in time when the sample has been verified (accepted).
User who verified (accepted) the sample. |
| Validation | Point in time when the sample has been validated.
User who validated the sample. |

6.3.4 Confirming Results with Discrepancies

Result verification - User Admin BIO-RAD

Patient / sample data:

Sample: 0802001
 IH-1000 0300999:
 Bloodgr.+Rev.gr.:A,B,AB,D/A1,B+C,c,E,e,K,ctl
 (5041+5011)

Test time: 08/02/2017 13:57
 Processed by: IH-Result Simulator User
 Verified by:
 Validated by:

Anti-D VI-
 50410.00.00/17.02/000003
 Reaction comment:

Current results

| 1 | 2 | 3 | 1 | 2 | 3 | 4 | 5 | 6 |
|---|----|-----|--------|--------|---------|-----------|----|-----|
| I | II | III | Anti-A | Anti-B | Anti-AB | Anti-D VI | A1 | B |
| - | - | - | ++++ | - | ++++ | ++ | - | +++ |

Result:
A Rh D not interpretable CcEE Kell negative ABS negative

AB0: A ABS: ABS negative
 Phenotype: Cc EE DAT:
 Rhesus-D: Rh D not interpretable Auto ctrl.:
 Kell: Kell negative

A discrepancy has been identified

List of discrepancies:
 - Rhesus-D system (Selected value: Rh D not interpretable)

Accept all

Figure 105. Result verification window: A discrepancy has been identified

When IH-Com identifies a discrepancy, the interpretation line in the verification check is highlighted in red and a notification A discrepancy has been identified is displayed.

If the cursor is moved over this message, the discrepant parts of the interpretation are displayed (AB0 system, Rhesus D system, etc.).

If the reason for the discrepancy is not to be found within the current result, but is based on a discrepancy between current and previous result, an appropriate message appears

Example: AB0 system (Preresult)

A preresult discrepancy has a higher priority in the display than the current discrepancy in the same test area.

Example:

If the AB0 system differs from the previous result and at the same time there is a discrepancy in the current test, only the message AB0 system (Preresult) appears for the test area AB0

No double display is shown.

The discrepancy message displayed is automatically refreshed as soon as reaction strengths are changed.

Example:

In case of discrepancy attributed to a doubtful reaction, the notification disappears immediately after it is changed to a valid reaction strength.

A If results are to be verified despite the discrepancy, a comment on the sample must be entered

See chapter [Adding Comments to Results on page 123](#).




Figure 106.

B If a comment was entered, a confirmation request appears when saving the result.

This ensures that a discrepant result cannot be saved accidentally.




Figure 107.

k Results that exhibit a reaction intensity that requires further verification based on the package insert of the gel card are also displayed as being discrepant.

6.3.5 Confirming Results in Case of invalid Quality Controls

If no valid quality control is available at the time of testing, the reaction column is marked with a red QC. The tooltip indicates: There is at least one reagent without valid quality control.




Figure 108.

To accept this result, a comment is mandatory.

Č The Missing quality control(s) window automatically opens.




Figure 109. Missing quality control(s) window

| | |
|---|--|
| A | <p>Enter a comment for each single test.</p> <p><i>The comment Accepted by [user] without valid QC: is automatically included. Enter a comment next to it.</i></p> |
| B | <p>With Use this comment for all tests (not commented yet) (1) comments are assigned to all tests without valid quality control.</p> <p><i>The comment Accepted by [user] without valid QC: is automatically included. Enter a comment next to it.</i></p> |
| C | <p>Save the comment with Save comment.</p> <p><i>It is possible to assign a comment to each test individually and save all the comments with Save comment.</i></p> |
| D | <p>If a sample comment is already saved, it can remain and confirmed.</p> <p><i>The default comment Accepted by [user] without valid QC: appears automatically below.</i></p> |

Comments for invalid quality controls are printed out for each test separately only on the daily journal. Comments are not printed out on the printout of results.

6.3.6 Messages before Reading Results

IH-Com displays messages for certain results when result verification is opened. These messages must be confirmed before the verification window can be viewed.




Figure 110. Example of message

The following examples show a selection of such messages and their background:

| Message | Background |
|---|---|
| D weak or D partial is currently selected. Please verify the Rhesus-D result! | The result for Rhesus D was saved as a result D weak or D partial in the meantime. When assessing the results, the interpretation of the D feature should be carefully considered. |
| The current sample has a rare phenotype. | <p>Rare factors of R, C, D, E, e (phenotype) are mentioned:</p> <ul style="list-style-type: none">• CCddee;• ccddEE;• CCD.EE;• CCddEE;• CCddEe;• CcddEE. <p>NOTE: These characteristic values that generate the message can be adapted. Contact your service representative.</p> |
| A discrepancy has been identified ... | The result has since been saved with a discrepant interpretation. In the assessment of results, the interpretation of the characteristics sought should be carefully considered. |

6.3.7 Reading Results Manually

Results with pipetted status can be read manually.

A

Select the appropriate test and the context menu Read manually.

| Sample | Device / test | Processed at | Processed by |
|------------|---|---------------------------------|--------------|
| 1140869080 | 7: Bloodgroup: A.B.D Confirmation for ... | 23/07/2014 15:40 | saxo |
| 1140869080 | 7: Ab.screening: I,II,III (IAT) (5053) | 23/07/2014 15:40 | saxo |
| 1140869080 | 7: Autocontrol: (IAT) (5053) | 23/07/2014 15:40 | saxo |
| 1140899320 | 6: Identification: 11 test cells (IAT) (50... | 23/07/2014 11:33 | Techno D0108 |
| 1140899320 | 6: DC Screening II: IgG.C3d.ctI (5056) | 23/07/2014 11:33 | Techno D0108 |
| 1140902160 | 7: Ab.scre... | Green, Lima *12/10/1988 [48583] | |

Edit result

Show result details

Read manually

Request more

Figure 111. Selecting Read manually via context menu

It is also possible to select the button Read manually.

B

Reading card 5053111111612111124

| Pos. | Test | Result | Sample / Patient |
|------|------|---------|-----------------------------------|
| 1 | I | + | 1140902160 Green, Alma *12/10/... |
| 2 | II | no res. | 1140902160 Green, Alma *12/10/... |
| 3 | III | no res. | 1140902160 Green, Alma *12/10/... |

1

2

3

4

5

6

+

...

...

++++

+++

++

+

+/-

-

?

DP (Double population)

no res (empty)

Edit reaction comment

Figure 112. Selecting reaction result via context menu

Move the cursor over the appropriate reaction column and select visual interpretation of the reaction grade from the context menu.

C

Save the entry with the button Save.

Enter reaction intensities with the numeric keypad

Reading card 5053111111612111124

| Pos. | Test | Result | Sample / Patient |
|------|------|--------|-----------------------------------|
| 1 | I | + | 1140902160 Green, Alma *12/10/... |
| 2 | II | - | 1140902160 Green, Alma *12/10/... |
| 3 | III | ++++ | 1140902160 Green, Alma *12/10/... |

1

2

3

4

5

6

+

-

++++

| Input | Reaction |
|-------|------------------|
| 0 | 0 |
| 1 | + |
| 2 | ++ |
| 3 | +++ |
| 4 | ++++ |
| 5 | +/- |
| - | 0 |
| + | ++++ |
| d/D | DP (mixed field) |
| ? | ? (questionable) |

Figure 113.

- D

The cursor moves automatically to the next entry field.

Save the entry with the button Save.

If the cursor is placed directly over a well, movement to the next entry field is not automatic.
- E

50530.00.00/19.12/000003

Reaction comment:

In the Result verification view a graphical representation of the reaction strength is displayed.




Figure 114.

Further procedure corresponds with in section [Changing and Confirming Results on page 103](#) described.

6.3.8 Batch Verification

Batch verification allows to verify multiple similar results in one working step.

As soon as verifiable results exist for a sample, they are displayed under the Batch Verification tab. Sorting takes place in groups by gel card type, and then by test and single interpretation.

If a device filter has been set in Common Settings, results only from this device or Lab location are displayed.




Figure 115.

- 1 Start batch verification
- 2 Gel card type
- 3 Name of the processed tests
- 4 Interpretation of single tests
- 5 Number of processed tests and gel cards used

The interpretation (e.g. AB Screening positive) in this case is only an interpretation of a single test for display, grouping and sorting. No **discrepancy** control or alignment with pre-results takes place at this stage.

Immediately after batch verification an overall interpretation of the sample is initiated which can differ from the one displayed.

If crossmatches are processed using different tests (e.g. Coombs and NaCl), the interpretation group all reactions negative may also contain incompatible crossmatches. A modification of the status of transfusion is only possible under Result Verification of the corresponding sample.

If no interpretation text is available (e.g. for antibody panels), groups are displayed in this way: all reactions negative, positive and negative reactions, ? Reactions.

Positive and suspect reaction groups are displayed in multiples, if reaction strengths differ.

- A

Select the result group that is to be verified
- B

Open the Result Verification with the button Open Batch Verification.
Under Result Verification results of the selected group are displayed by gel card number descending.




Figure 116.

- 1

Patient and sample data of the corresponding tests
- 2

User and testing time of the selected test
- 3

Test comment field for the selected test
- 4

Reactions of all tests of the selected result group
- 5

Context menu
- 6

Magnified well of the selected reaction
- C

Accept all displayed reaction within one working step with the button or context menu Accept all displayed wells.
- D

Reject all displayed reactions with the context menu Reject all displayed wells.
- E

Modify reaction strengths by using the context menu

| | | |
|---|---|---|
| F | | By accepting with the button the next group of results are displayed automatically. |
| G | | Eventually enter a comment to the result.
<i>Refer to section Adding Comments to Results on page 123.</i> |
| k | | In batch verification, test comments rather than sample comments should be entered.
Sample comments should be entered in result verification. |
| H |  | If all results of the selected group are accepted, save the procedure with Save.
<i>The Result Verification window will be closed and results are registered as automatically accepted.</i>
<i>Results are then subject to the settings defined in Verification of overall result for further validation or LIS transfer.</i> |

If a result group contains a large amount of samples, only one portion is displayed under Result Verification. After confirming with Save, the Result Verification opens automatically again for the following samples of this result group.

If the Result Verification is opened by the button Batch Verification Run, it will be re-opened automatically for the next result group.

6.3.9 Changing Interpretation

It is possible to change the blood group interpretation determined by IH-Com.

Example:

When no isoagglutinins are demonstrable in a sample for an elderly patient, IH-Com assesses the result as discrepant.

A Modify the specified interpretation in the Result window:




Figure 117. Result verification window: Changing interpretations

k Interpretations can only be changed for current tests.

Example: An interpretation of a rhesus subgroup cannot be changed, if it originates from a previous result.

B Comment on your changes in the Sample comment accordingly.

Č If no comment has been entered, a prompt will be displayed.




Figure 118. Result verification window: Discrepancy tooltip


Č After the change the overall interpretation is highlighted with an orange background.

Č The discrepancy tooltip displays the interpretation indicated by IH-Com and the interpretation selected by the user.


- C If your selection is not saved, it can be undone with Recalculate Interpretation, to restore the interpretation automatically calculated by IH-Com.


Recalculate Interpretation

- D If your selection is to be applied, a sample comment must be entered. See chapter [Adding Comments to Results on page 123](#).



- E If a comment is entered, a confirmation request appears when saving the result. This ensures that a changed interpretation cannot be saved accidentally.





6.3.10 Adding Comments to Results

In IH-Com there are several ways to add comments to results.

Comment on a Sample, a Test or a Crossmatch Test

- A Select any reaction column and open the input mask via context menu Edit sample and test comment.




Figure 120.

Alternatively it is also possible to select the button beside the comment field to open the input screen.

- B Enter a sample comment (1) and/or a test comment (2).




Figure 121. Sample comment window: Inserting predefined comments

If predefined comments are saved, they can be inserted with Insert text block using the context menu.

See how to store predefined comments in chapter [Text Modules \(Pre-defined Comments\)](#) on page 182.

- C Confirm your entry with <OK>.

Č The comments are displayed in a box in the upper right corner of the results check.

If a crossmatch test has been completed, a crossmatch comment can also be included.




Figure 122.

Comment on a Single Reaction

- A Select the reaction column to be edited.
- B Open the input window using the Edit reaction comment context menu.
- C Insert a text.

Predefined comments can be inserted using the context menu.





Figure 123. Reaction comment window: Inserting predefined comments

- D  Confirm your entry with <OK>.
- E Move the cursor over the reaction column.

Č The comment is displayed in the Reaction comment: box.




Figure 124.

- F Confirm your comment input with Save.

k It is possible to input comments at any time.

6.3.11 Sending a Result

Occasionally it is necessary to resend a result to the LIS.




Figure 125.

Retransmit all Results of a Sample

- A Select the context menu Send results in the patient list

Resend the Result of a single Test

- A Select the test in the result list.
- B Select the context menu Send results.

Retransmit Results of Several Patients Simultaneously to the LIS

- A Activate Show multi select from the context menu.
- Č A checkbox displays in the patients list between sample number and name fields.*




Figure 126.

- B Select the patients to be transmitted.
- C Select Send selected results.
- Č The multiple selection in the Results area is automatically closed when the transfer is completed.*
- k** If the results selected for sending are not yet verified, they are sent to the LIS as temporary results.

6.3.12 Discrepancy Check with previous Results

When the result of a test is available in IH-Com, discrepancy checks will automatically be performed:

- Discrepancy check within the current result

E.g. a + or ++ reaction for blood grouping would be displayed in IH-Com as not interpretable and a message will be displayed “A discrepancy has been identified”.

See sections [Test, Interpretation and Result on page 22](#) and [Confirming Results with Discrepancies on page 110](#).

- Discrepancy check with previous results

The current result will be compared with previous results for the same patient or sample.

The comparison with previous results is based upon the following categories:

- blood group and reverse grouping
- RH D including Rh D weak and Rh D partial
- Phenotype including Anti-CDE
- Kell
- each single Antigen incl. Cw and Cellano
- Ab Screening
- DAT
- Autocontrol

It is expected that the current result will be the same as the most recently processed previous result. If not, this will be identified as discrepant in IH-Com.

For example: Previous result A Negative, current result A Positive would be identified as discrepant in IH-Com.

The screenshot shows a web interface for blood test results. At the top, a red banner displays the message "A Rh D not interpretable". Below this, there are several dropdown menus for different blood test parameters: AB0 (set to 'A'), ABS, Phenotype, DAT, Rhesus-D (set to 'Rh D not interpretable'), Auto ctrl., and Kell. A red text message "A discrepancy has been identified" is prominently displayed in the center. To the right of the Rhesus-D dropdown, there is a text label "Antibodies / other antigens:" followed by an empty box. At the top right, there is a button labeled "Recalculate Interpretation". At the bottom right, there is a button labeled "Details ...".

Figure 127.

If a previous result was not interpretable, this will also be applied during the discrepancy checking.

If a test is rejected or a sample marked as a suspect result, it will not be applied during the discrepancy check.

D Weak and D Partial results must be set manually in result verification to be applied as previous results. When a test for the same sample or patient is later processed for Rh D, a warning message appears before result verification is opened:




Figure 128.

If only a part of a result category is defined in a test, the whole category is saved as a current result. Example: ccee is defined as a pre-result, current testing is only for Anti-c. The entire phenotype is applied as the previous result even though Anti-e has not been tested.

When a sample or patient has been tested previously for one of the result categories but the current test does not include this category, the result will be displayed with a superscript 1 in the interpretation.

| Consulted for plausibility check / as pre-result | |
|--|---|
| Within 24 hours | Results without patient data for a sample Later on results for this sample number are determined as a new sample. |
| Within 72 hours | Ab-Screening tests, DAT, Autocontrol
This time period can be changed by your Bio-Rad service representative. |
| Permanently | Results with patient data (except Ab-Screening tests, DAT and autocontrol) for all sample numbers. |

Examples

| Test | Current result | Interpretation |
|---|----------------|--|
| Example 1 | | |
| AB0/D | A D+ | A D+ |
| AB0/D | B D+ | ? D+, manually change to B |
| AB0/D | B D+ | B D+ |
| AB0/D | B D+ | B D+, manually changed to Dweak |
| AB0/D | B D+ | B Dweak |
| AB0/D | B D+ | B Dweak |
| AB0/D | B D- | B Dweak |
| Example 2 | | |
| AB0/D | A D+ | A D+ |
| Rh D | CcEe | A ¹ D ¹ CcEe |
| DAT | DAT+ | A ¹ D ¹ CcEe ¹ DAT+ |
| Anti-c | c+ | A ¹ D ¹ CcEe |
| Anti-c | c? | A ¹ D ¹ RhF n.i. |
| DAT | DAT+ | A ¹ D ¹ RhF n.i. ¹ DAT+ |
| Anti-c | c+ | A ¹ D ¹ RhF n.i. |
| Each line corresponds to a new sample number. | | |

6.3.13 Excluding Results from Pre-Result Check

Sometimes an already verified or validated result should not be no longer be used as a pre-result.

User rights 60 are required for excluding results from the pre-result check.

b Exclusion of pre-results must be performed with careful consideration in order to avoid missed discrepancies.

A Select the test of a sample

Figure 129.

B Select the context menu Suspect result (No longer used as a previous result).

All tests of the sample are set to suspect and are not consulted for a [discrepancy](#) check for not yet verified samples and future tests

User rights of level 50 are sufficient for this procedure.

The interpretation of the as suspect specified sample remains under Result details for traceability

C To deactivate this specification, reselect Suspect result (No longer used as a previous result).

Figure 130.

