
Approved

«21» November 2022

President

Mr. Wang Zhongliang

Lakong Medical Devices Co Ltd., China
No. 13 Jingci Road, Shilong, Dongguan, Guangdong, 523326 China

tel.: +86-(769)-81807881;
fax: +86-(769)-81807111

E-mail: <http://lakong.medincn.com>
Website: <http://www.medincn.com>

INSTRUCTION FOR USE OF MEDICAL DEVICE:

Disposable dental micro-applicator Clean+Safe®

1. Name of the medical device:

Disposable dental micro-applicator Clean+Safe[®], versions:

Disposable dental micro applicator Clean + Safe[®], size Regular (M6500R), color: green/ orange;
Disposable dental micro applicator Clean + Safe[®], size Fine (M6500F), color: yellow/ blue;
Disposable dental micro applicator Clean + Safe[®], size Superfine (M6500SF), color: white/ pink;
Disposable dental micro applicator Clean + Safe[®], size X-Slim (M6500XS), color black.

2. Name and address of the manufacturer:

Lakong Medical Devices Co Ltd, China
Address: No. 13 Jingci Road, Shilong, Dongguan, Guangdong, 523326 China
Made in China
tel.: +86-(769)-81807881; +86 (076)-981807881
fax: +86-(769)-81807111
website: <http://lakong.medincn.com>

Name and address of the manufacturer's authorized representative in the RF

"PROTECO", JSC
192102, Saint Petersburg, Salova street 57, Building 1, Letter "I"
Tel. / fax: +7 (812) 779-30-90
e-mail: info@protecodent.ru

Importer to Republic of Belarus:

PTUE VMK-Dent

Platonova str.1B, premises 41B room 1, 220034 Minsk, Republic of Belarus

3. Purpose and application:

Intended for application of dental substances and materials on teeth of a patient.

Application:

- For professional use in dentistry.

Potential users of the medical device: dentists.

4. Contraindications to use:

- patient's allergy to any component of the micro-applicator;
- impossibility to observe the established procedure of application.

Side effects: none.

The medical device is intended for single use only.


Instruction for use of medical device:

Disposable dental micro-applicator Clean+Safe[®]

5. Main functional elements of the medical device:

The micro-applicator represents a dental hand instrument in the shape of a rod whose bending part ends with a brush-like head of different diameters. Brush fibers do not absorb fluids, solutions, vanishes, which allows application of the whole amount of taken materials. The tip, bendable into different angles, allows application of materials even into hard-to-access areas. Color marking and changeable length of the handle (if necessary, it can be shortened by breaking off the bottom part) exclude cross-infecting. The micro-applicator X-Slim has a finer prolonged head and is intended for endo-canal treatment.

The medical device is non-sterile single-use.

6. Application procedure

- Take a micro-applicator of the needed size. Turn carefully the feeder-cap and take out a micro-applicator. If necessary, bend the micro-applicator into the desired angle (it is not recommended to bend the micro-applicator many times)
- Take the dental material that is necessary to apply and place it (required amount) on the mixing pad (or any other device intended for this)
- Using the micro-applicator, take the necessary amount of the dental material (micro-applicator's fibers retain the dental material but do not absorb it)
- Apply the dental material
- The micro-applicator is intended for single use and cannot be reused.

The sterilization and/ or disinfection of the micro-applicators is not required before the first use.

7. Operating and storage conditions

8.1 Keep in a dry well-ventilated room, avoid exposure to direct sunlight.

Storage conditions

Ambient temperature	-20 - +50°C
Relative humidity	10 to 80%
Air pressure	500 to 1060 hPa

Operation conditions: temperature of use from + 15°C up to +42°C

8. Completeness

No.	Name	Quantity
I	Disposable dental micro-applicator Clean+Safe [®] *	4 tubes of 100 pcs
II	Instruction for use	1 pc.

* – select the needed version within the range of presented variants:

Regular: green, orange

Fine: yellow, blue

Superfine: white, pink

X-Slim: black

9. Shelf life

Shelf life: is not limited, if the storage and transportation conditions are observed

10. Symbols and images used in labeling of the medical device and their description

- Symbols used in labeling:

	Lot number		Temperature range
	Non-sterile		CE marking
	For single use only		- Dispose separately from household waste
	Catalogue No.		Marking of conformity to GOST R

11. Procedures of disposal and destruction of the medical device

The medical device is to be disposed according to the rules of collection, record-keeping and disposal established by the medical institution.

The medical device falls to category of epidemiologically hazardous medical waste: materials and instruments, objects contaminated with blood and/or other biological fluids.

Consumer package, defective device refer to epidemiologically safe waste, similar in composition to household waste.

12. List of international standards

ISO 9001: 2015, ISO 13485:2016, EN 1041:2013, EN ISO 14971:2007, ISO 15223-1:2016, ISO 10993-1:2009, ISO 10993-2:2006, ISO 10993-5:2009, ISO 10993-10:2010, ISO 10993-12:2012.

Instruction for use of medical device:
Disposable dental micro-applicator Clean+Safe[®]

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