
Vision Screener

User Manual

(V100)

Shanghai MediWorks Precision Instruments Co., Ltd

Catalogue

Preface.....	3
■ Performance parameter.....	4
■ Attentions:	4
■ Prohibition:.....	6
■ Main security features of the product.....	6
■ EMC (Electromagnetic Compatibility).....	7
■ Chapter 1 Installation.....	11
1.1 Product list.....	11
1.2 Graphical description of each part of the vision screener.....	11
1.3 Battery installation.....	12
1.4 How to charge.....	13
1.5 Remove and installation of lens cover.....	13
1.6 Power on and off.....	14
1.7 Flip the screen.....	14
1.8 Preparation before screening.....	15
■ Chapter 2 Introduction of the home interface.....	16
2.1 Home interface.....	16
2.2 Add patient.....	17
2.3 Patient list.....	19
2.4 Records detail.....	24
2.5 Report.....	25
2.6 Test.....	27
2.7 Scan the QR code.....	28
2.8 Settings.....	28
2.9 Diagnostic result reference.....	38
■ Chapter 3 Cleaning, protection and maintenance.....	39
3.1 Cleaning method.....	39
3.2 Product cleaning cycle.....	39

3.3 Protection.....	39
3.4 Battery maintenance.....	39
3.5 Product service life.....	40
3.6 Final disposition.....	40
■ Corrective maintenance.....	40
■ Exceptions.....	40

Preface

Thank you for purchasing the vision screener produced by Shanghai MediWorks Precision Instruments Co., Ltd!

Here are some basic information and feature parameters of the vision screener you purchased!

Overview:

✎ This User Manual is a part of the vision screener. The operation instructions and technical instructions are included in this manual.

✎ This User Manual contains instrument operation instruction and technical Instruction, the equipment classification of the vision screener, according to the requirements of IEC 60601-1, is also included in this manual.

✎ According to the IEC 60601-1 standard, the applicable marks have been permanently pasted on the instrument, and instructions are also given in the manual.

✎ The Vision screener is intended to measure the refractive power of the eye by measuring light reflexes from the retina. It also estimates pupil size, pupil distance and eye gaze deviation.

✎ This instrument is a Type I laser product (containing invisible laser radiation).

✎ The laser may injure the human eye. Do not make the infrared light irradiate the human eye directly for more than 2 hours.

✎ The vision screener is a hand-held instrument with a touch screen on the back for displaying and aligning to the patients' eyes. The screen also provides user interface functionality for data input and display. The vision screener provides visual and audible signals to attract the subject's attention to make them stare. The device provides a fast data collect for easy evaluation of the most difficult pediatric subjects of limited compliance. Each evaluation provides simultaneous binocular assessment of the subject's refraction, pupil size and gaze.

■ Product composition

The vision screener consists of the device, lens cover, power adapter, USB cable, bluetooth printer and lithium battery.

■ Scope of application

The vision screener is designed to be used in a medical environment by a person with basic knowledge of ophthalmology to assess the refractive status of the eye for people older than 6 months. Parameters measured by the Vision Screener include the spherical power, cylinder power, angle of cylinder power, pupil size, and pupil distance.

■ Performance parameter

Distance auto-induction	100cm±5cm
Pupil diameter	4.0~9.0mm,±0.1mm
Pupil distance	30~80mm,±1mm
Spherical power:	-9D~+7D,resolution0.25D/0.01D,±0.5D
Cylinder power:	0D~3D,resolution 0.25D/0.01D,±0.5D
Angle of cylinder power:	0~180degree,resolution 1degree,±5 degree
Fixation angle	0-20 degree
Measuring time:	1 Sec.(after focusing)
Data transmission:	WIFI
LCD display screen:	4.3" full LCD touch screen
Power supply:	Two 18650 rechargeable lithium batteries 3.7V/3400mAh
Charging interface:	USB
External power adapters:	Power supply:~100-240V, 50/60HZ,35VA Output:5V,2A
Illumination near-infrared spectral range	800nm~880nm
Attracting way	Sound and flashing light
Applicable group	People older than 6 months

Operating environment	Temperature	+5℃~+40℃
	Relative humidity	≤90%
	Atmospheric pressure	860hPa~1060hPa
Storage environment	Temperature	-40℃~+55℃
	Relative humidity	≤90%
	Atmospheric pressure	860hPa~1060hPa
Transport environment	Temperature	-40℃~+55℃
	Relative humidity	≤90%
	Atmospheric pressure	860hPa~1060hPa

■ Attentions:

- ✎ This product should be used by qualified medical personnel.
- ✎ Before the operation, please read carefully the“patterning, symbols and marks on productssymbolsafety precautions”.

- ✎ This equipment is not used in the oxygen-rich environment.
- ✎ In order to prevent hurting or damaging the device, please do follow the operation instruction.
- ✎ In order to prevent electric shock, do not disassemble and repair this device by yourself. If any problems, please contact Shanghai MediWorks or its authorized distributors.
- ✎ Please do use the appointed charging plug to avoid electric shock and damaging the device.
- ✎ In order to avoid electric shock, please do not insert any metal objects into the holes or slots of the device.
- ✎ If the device fails: smoking or burning smell, in order to prevent fire, immediately turn off the power supply, take out the battery.
- ✎ The device is not waterproof, in order to prevent fire and damaging the device, do not put cups or containers with water beside the device. If the device is wet, do not dry it by baking, microwave, high pressure or ultraviolet.
- ✎ In order to prevent electric shock, fires, short circuit and other danger, please do not place the device in a water or humid environment.
- ✎ In order to prevent loosening of the device due to shaking, when the device is not used, please put the device into the portable box.
- ✎ The device should be placed in a clean environment, and keep far away from the fire source.
- ✎ Please do check if the device is damaged or not before operation. If there are indications that the device has been damaged or any defects, please do contact Shanghai MediWorks Precision Instruments Co., Ltd.
- ✎ The battery is a replaceable consumable. Only the 3.7V 18650 battery that meets the current IEC62336 standard is allowed. The Panasonic NCR18650B battery is recommended. Make sure the battery is installed correctly. Any problems caused by the battery are not related to Shanghai MediWorks Precision Instruments Co., Ltd.
- ✎ When not using the device for a long time, please remove the battery.
- ✎ Because of the strong infrared rays in sunlight, the vision screener cannot be used under strong sunlight and is recommended to use in a darker room. Under the use of strong light, it will not be able to shoot or increase the error due to the oversize change of the examiner's pupil.
- ✎ Shanghai MediWorks Precision Instruments Co.,Ltd will provide the circuit diagram,

components list, description and technical parameters for maintenance service personnel.

 In order to avoid danger, do not use the device while charging.

■ Prohibition:

There are no known contraindications for this product.

■ Safety marks, icons and warning symbols stuck on this instrument.

According to the IEC60601-1 standard and other relevant standards, the following patterning, symbols and marks are used on the products, and their specific meaning is as shown in table 1 below:

No.	mark	Illustration
1		Production date
2		Manufacturer information
3		Refer to operation instruction
4		WEEE mark: please dispose the waste generated by the machine in accordance with relevant laws and regulations
5		CE mark
6		Power switch

■ Main security features of the product

- According to the classification standard of medical electrical equipment of IEC 60601-1 standard, V100 belongs to the ordinary equipment with internal power supply;
- Classified according to the degree of protection against electric shock: the device has no application part;
- According to the degree of protection against harmful liquid: IPX0;
- This instrument is categorized to Class II Type B according to IEC 60601-1-2005 standard, which cannot be used under two circumstances: a flammable anesthetic gas and air mixture, oxygen or nitrous oxide gas and air mixture. This instrument

can be running continuously;

- e) The instrument is not an AP type or APG type equipment.

■ EMC (Electromagnetic Compatibility)

The V100 complies with the International Electrotechnical Commission standards (IEC 60601-1-2) for electromagnetic compatibility as listed in the tables below. Follow the guidance in the tables for use of the V100 in an electromagnetic environment.

Guidance and manufacturer's declaration electromagnetic emissions		
The V100 is intended for use in the electromagnetic environment specified below. The customer or the user of the V100 should assure that it is used in such an environment.		
Emissions test	Compliance level	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The V100 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The V100 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	N/A	

Guidance and manufacturer's declaration electromagnetic immunity			
The V100 is intended for use in the electromagnetic environment specified below. The customer or the user of the V100 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floor should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV	±2 kV for power supply lines ±1 kV	Mains power quality should be that of a typical commercial or hospital environment.

	for input/ output lines	for input/output lines	
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage, dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 sec	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the V100 requires continued operation during power mains interruptions, it is recommended that the V100 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration electromagnetic immunity			
The V100 is intended for use in the electromagnetic environment specified below. The customer or the user of the V100 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms (V1=3)	Portable and mobile RF communications equipment should be used no closer to any part of the V100, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ 150 kHz to 80 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m (E1=3)	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).

			Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the V100 is used exceeds the applicable RF compliance level above, the V100 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the V100. b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the V100 The V100 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the V100 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the V100 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1.2 \sqrt{P}$	80 MHz to 800 MHz $d=1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d=2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.


■ Chapter 1 Installation

1.1 Product list


Before use, please check that the package contains all of the following items. If any missing, please contact with the agencies.

Series No.	Accessory name	Qty
1	Host	1
2	Lens cover	1
3	Power adapter	1
4	Connecting wire	1
5	Portable box	1
6	Printer	1


1.2 Graphical description of each part of the vision screener

Basic composition schematic diagram 1


Basic composition schematic diagram 2

Basic composition schematic diagram 3

1.3 Battery installation


Open the battery cover, load the batteries, and please note that the two batteries are in the same direction. The positive poles of both batteries face the battery cover, as shown below. You will not be able to turn on the machine if the direction is wrong. In order to keep the battery in good condition, please fully charge before first time use.

Battery installation

1.4How to charge

Connect the vision screener to the power adapter with a USB cable. The vision screener will start to charge automatically.

USB socket schematic diagram


Charging mode instruction

Power supply indicator	Status
Shut down and charge	Green light will be always on during charging. The green light will be off when the instrument is fully charged. The green light will be flashing when the battery error occurs.
Turn on and charge	White light is always on during charging. Blue light is on when it is fully charged. The white/blue light is flashing when the battery error occurs.
Running/Standby	Blue light is always on.

The general situation of battery error is that the battery is not installed or the polarity is reversed. Be sure to keep the battery installed correctly.

1.5 Remove and installation of lens cover

As shown below, pinch the buttons on both sides of the lens cover to remove the lens cover. Slightly force the slider to retract, then loose them, the lens cover is installed.


1.6 Power on and off

Power on: Press the power button on the device, the device will be turned on.

Power off:


Method 1. Long press the power button for 3 seconds, the dialog box of shutdown will pop -up on the screen. Click “Shut Down” on the screen to shut down.

Method 2. Long press the power button for about 8 seconds, and the device will be forced to shut down.

1.7 Flip the screen

Under different operating situation, such as the sitting position of the subject and standing position of the tester, the screen can be flipped to a suitable Angle for easier observation and operation of the display screen.


Schematic diagram

1.8 Preparation before screening

Step 1: Take off the lens cover, starting up, entering the shooting interface.

Step 2: The subject should stand or sit comfortably, try to keep his head in a stable position. Make the subject stare at the camera, or you can use light and sound to attract the subject's attention. Move this device to keep 1 meter away from the subject, remind the subject to keep watching, as shown in below picture.

Taking picture schematic diagram


Step 3: Adjust the device slightly until both eyes of the subject are clear on the screen. Keep the device stable until the result comes out.

Note: Harmful radiation exposure may occur with the incorrect steps during operation.

■ Chapter 2 Introduction of the home interface

2.1 Home interface

When the vision screener is powered on, the home interface will appear. The home interface includes status display, quick test, patient list, add patient, QR code scanning and system settings. The status display shows, respectively, the server login prompt (the server login successful state), the WIFI connection prompt (Display when the WIFI has been connected), the printer connection prompt (Display the status that the printer has been connected), the battery indicator, as shown in below picture.

2.1.1 Quick test, patient list, adding patients and QR code scanning

➤ Quick test: In this mode, the user only needs to choose the age bracket of the patient, no need to fill in a specific medical record. The information can be tested. After the test is completed, the user can choose whether to save the test data.

Advantages: Users can quickly test patients, which is more convenient for temporary tests.

➤ Patient list: This mode has functions as below:

1. View, select or search for records that have been created.
2. The records to be tested can be downloaded in batches from the server, and can be uploaded to the server after test.
3. View reports of tested records.
4. Management of medical record information.

Advantage: It can create records in batches. That meets the requirement of creating large numbers of records (such as schools).

➤ Adding patients: In this mode, the user must create medical record at first, and the data will be saved after the test automatically.

Advantage: The case can be created without network, which is more convenient for

individual patients or small scale tests.


➤ QR code

Advantages: Creating new patient case can work in a network-free situation, manual input will not be necessary when creating patient information, which can greatly improve test efficiency.

2.1.2 Operation method

1).Click "Quick Test" button to pop up the age range selection box.

- Select the age range, such as "1-3 years old", enter the test interface(please refer to section 2.6 of this chapter for detail operation). This option will affect the test report, please select the appropriate age range.

2).Click "Patient List" button to enter the records screening interface (refer to section 2.3 of this chapter for detail operation).

3).Click "Add Patients" button to enter the records creation interface (please see section 2.2 of this chapter for detail operation).

4).Click the “QR code” button to enter the QR code scanning interface (for details, please refer to section 2.7 of this chapter).

5).Click “Setting” button to enter the setting interface (please refer to section 2.8 for detail operation).

2.2Add patient

The medical record of the subject to be tested should be created before starting the test, including the basic information of the patient (ID, name, birthday, gender) and the glasses type to be photographed.

Click “Add Patient” in the main interface to enter the records creation interface, as shown in below picture.

The screenshot shows a patient information form with the following fields and options:

- ID:** 000002
- Birthdate:** 06/08/2009 (with a calendar icon)
- First Name:** Test01
- Last Name:** (empty)
- Glasses:** None (selected), Glasses, IOL
- Gender:** M (selected), F
- Buttons:** Save (checkmark icon), Cancel (X icon)

- **ID:** According to the rules of the user's hospital, Letters and figures are acceptable. The maximum length is not more than 10 characters.

- **Birthdate:** Birthdate input can be done either keyboard input or date selection.


2.2.1 Keyboard input: Click birthday input box, the keyboard will appear at the bottom of the interface (by default, numeric input), click on the numeric key to enter the number. The input process will fill the separator automatically. Until you delete the separator manually, then you need to fill it manually.

The screenshot shows a numeric keyboard interface with the following elements:

- Input field:** 2017/
- Language/Theme:** 中文, En, 123
- Keypad:**

1	2	3	←
4	5	6	↶
7	8	9	@,\$...
.	0	/	✓

2.2.2 Date selection: Click "Birthday Selection" button  on the right side of the birthday input box to enter the birth date selection interface, as shown below.

In the date selection box, select year, month, day in turn, then click "OK" button in the interface to complete the date selection.


- Name: Click First Name and Last Name input box, the interface displays the keyboard, click the input key to complete the name input.
- Gender: Click "Gender" radio button to complete gender selection of the patient.
- Glasses: Click any radio button of the glass selection to complete the selection.

After the information input, click "Save" button if it required. It will save to the database, and entering into shooting interface (please refer to section 2.6 of this chapter for detail operation) if the data is correct. The frame of input box will turn red as a prompt, if the input information is missing or the input information format is not correct.

2.3 Patient list

The main functions of the "Patient List" mode include displaying the created records, searching and View records in different conditions, downloading the records from the server, uploading the tested records and deleting the records.

Click "Patient List" button in Home interface to enter the patient records screening interface, as shown in below picture.

2.3.1. Record list: the record list will display the newest 50 records which have been created; user can find the non-display records by the function of “Record Screening”.

2.3.2. Records screening: Search and View records in different conditions, the search condition includes name or ID, test status, as shown below.

- 1) Name or ID: In the search condition box, enter the ID or name of the medical record you want to search for, such as (001 or fn001), the medical record information will be displayed in the medical record list, as shown in below picture.

- 2) Test and Upload status: In the test status, select a test status (such as tested) in the drop-down list, and the records list shows the information of the tested records queried. As shown in below picture.


Two kinds of search conditions above can be used together.

- 3) Lock Search Conditions: Click on "Lock Search Conditions". After switching the interface back, the current search conditions will be retained. If not, the search conditions will be cleared.

2.3.3. To view the details of records: click any item in the records list, and the interface will show the detail information of records (please refer to section 2.4 of this chapter for detail operation).

2.3.4. Tick the medical record: Click "Select" button, a check button appears in the medical record list, as shown in below picture.


➤ Click any item in the records, the item will be in “Selected”. Click item again, it will be reset to “Unselected”.

➤ The user can select all by clicking “Select All” button, all items will be selected after clicking “Select All” button. Click “Select All” button again to reset all items to unselected.

➤ Click “Select” button again, the “Selected” button will be hidden. Then, click any items and the interface will enter the detail information interface.

2.3.5. Records Rank:

Click on the sorting button in the lower left corner of the interface and sort it after selecting the sorting condition in the pop-up sorting condition selection box, as shown in below picture.

2.3.6. Delete the records:


“Delete” button will be available when any records is selected. Click “Delete” button will enter the deletion confirmation interface. Click “OK” to delete selected

records. Click “Cancel” to cancel the deletion.


2.3.7. Download and upload records: Users can download records from the Mediview server or upload the tested records to the server. This operation needs to be implemented under the condition of WIFI connection and successful login to the server firstly (please refer to sections 2.8.6 and 2.8.7).

① Download the records: Under this operation, the medical records to be tested already created on the server will be downloaded. It is convenient for users to create medical records to be tested in batches.

➤ Click “Get” button, if the server is not successfully logged in, a login prompt box will pop up, as shown in below picture. Click “OK” button in the prompt box to enter the login interface (please refer to section 2.8.4), click “Cancel” button to cancel the download.

After login the server successfully, click “Get” button, please be patient and wait for the download to complete. After the download is complete, the "Waiting" prompt disappears and the downloaded medical records are displayed in the medical record list.


Special reminder: If a medical record with the same name and PID is established on both the device side and server side, when the user downloads, the medical record will be overwritten and the test data will be cleared. In order to ensure data security, upload the data to the Mediview server before downloading it.

② Upload records: Click the "Select" button, check the records to be uploaded, and click the "Send" button. The operation is the same as the download. The records can be viewed in the Mediview server after upload finished. If automatic upload is required by the user, check “auto-upload” during login (see section 2.8.4 of this chapter for detail operation). The records which have been tested but not have been uploaded will be uploaded automatically. If the records that have been uploaded to the server need to be uploaded again, they need to be selected manually and uploaded again.

During the download and upload process, the error message will pop out, when the device is disconnected from the server, the user need to try it again or re-login before retry.

2.4 Records detail

The interface shows the detailed information of the records and its relevant operation. The interface is shown as in below picture.

Detail information interface

1. Information display: Including the basic information of the patient, glasses types, test status and the creation time of the records.
2. Edit: Click “Edit” button, the interface enter the records information edit interface (please refer to section 2.2 of this chapter for detail operation). The birth date or glasses type will affect the test results, so the test data will be clean when the birth date or glasses type was changed, you need to test again.


3. Test: Click “Test” button to enter the test interface (please refer to section 2.6 of this chapter for detail operation). This operation will clean the previous data, please upload the data to Mediview server before the operations if user needs to save the previous data.
4. Report: the report button will be available when the records have been tested. Click “Report” button, it will enter report interface (please refer to section 2.5).
5. Back: Click “Back” button, the interface return to record list.

2.5 Report

The report interface shows the report of a existed record or the report of the “Quick Test”.

2.5.1 Report of the “existed medical record”

The report interface contains the name and ID of the patient, data display areas and function buttons, the data display mode is either negative cylinder display or positive cylinder display, as shown below.

4.09 mm ID: 00001 Name: test 4.18 mm

↓ 2° ↓ 1°

← 2° → 1°

OD

DS	DC	Axis
-6.50	-0.50	103
SE	-6.75	

PD 62.90 mm


OS

DS	DC	Axis
-7.00	-0.25	69
SE	-7.25	

Abnormal: Myopia(OD,OS)

Binocular Mono Print Next QR Code Patient List Home

Negative cylinder

4.09 mm ID: 00001 Name: test 4.18 mm

↓ 2° ↓ 1°

← 2° → 1°

OD

DS	DC	Axis
-7.00	+0.50	13
SE	-6.75	

PD 62.90 mm

OS

DS	DC	Axis
-7.50	+0.25	159
SE	-7.25	

Abnormal: Myopia(OD,OS)

Binocular Mono Print Next QR Code Patient List Home

Positive cylinder


1. Negative & positive cylinder selection: Click the cylinder in the setting interface, and return to the interface after the setting is completed (please refer to section 2.8.2 of this chapter for detail operation), the display mode in the interface will be switched to the data display under the corresponding cylinder.
2. Test results: The test results show abnormal items in the test data.
3. Retest: If the user need to retest, click “Binocular” button to Retest (binocular mode as default), click "Mono" button to Retest (mono mode as default).
4. Print: Click “Print” button to print the report.
5. Next: Click on the "Next" button, and the next patient will be selected for testing according to the sorting rules of the record list.
6. QR Code: Click “QR Code” button to scan the QR code.
7. Home: Click “Home” button to enter the main interface.

Attentions:

1. The size of the pupil will be displayed in red when its value is less than 3.5mm.
2. Re-test will clear the original data, if the original data is needed, user can upload to server to save it.

2.5.2 Report of the “Quick test”

In the “Quick Test” report interface the name and ID are empty and the “Save” button is available.


“Quick Test” report interface

1. Save: Click “Save” button to enter the records creation interface. Then you can save the test data (please refer to section 2.2).
2. Other operation: Please refer to section 2.5.1 for other operation.

Attention: the data and the test result in the picture above is only as an example, please do not take it as the reference for diagnose.

2.6 Test


In the shooting interface, the shooting program starts automatically in binocular test mode as default, the interface is shown as below picture.

1. Test:

- Position yourself approximately 1 meter from the subject. Slowly adjust the device to meet both of the subject's eyes. Adjust your distance from the subject until both eyes are clear on the screen.
- Please keep the image clear during the test.
- After the test is completed. It will enter the report interface.

2. Monocular and binocular mode switching: During the shooting process, click "Mono" button to switch to the monocular mode (OD test as default). Then click the left eye to switch to the OS test if required, shown in below picture.

3. Return to homepage: Click “Home” button to enter the home interface.

2.7 Scan the QR code

After clicking the “QR code” button on the interface to enter the scanning interface, then the shooting process starts automatically, the interface is shown as below picture.

- Position yourself approximately 1 meter from the QR code. Adjust your distance from the subject until the QR code is clear on the screen.
- When the scanning is completed, the patient’s scanned information will be added to the local data base automatically. It will enter the records detail interface.
- If the patient exists in the local database (with the same ID and same name), it will automatically read the corresponding patient information and jump to the patient detailed information viewing interface.


2.8 Settings

The “settings” interface contains language setting, report setting, system time setting, server setting and sleeping time setting, WIFI setting, software information (about), music setting, printer setting and other setting function, as shown in below picture.

2.8.1 Language setting

Click “Language” button in the “settings” interface to enter the language setting interface, you can set the language as Chinese or English, as shown below.


- Select language in the language drop-down list.
- Select date display format in the date format drop-down list.
- Click “OK” button to save the setting.

Click “Cancel” button to Cancel and return to the setting interface.

2.8.2 Cylinder setting


Click “Report” button in the “settings” interface to enter the report setting interface, the setting will affect the parameter display of the report interface.(please refer to section 2.5), as shown in below picture.

- Select the cylinder in the cylinder drop-down list, click “OK” button to save setting, and return to the setting interface.
- Click “Cancel” button to cancel and return to the setting interface.

2.8.3 System time setting

Click “System Time” button in the setting interface to enter system time setting interface, as shown in below picture.


➤ In the system time selection box, slide up and down to select the time ("now" refers to the current system time), click "OK" button to save and return to setting interface.

➤ Click “Cancel” button to cancel and return to the setting interface.

2.8.4 Server setting

Click the “Server” button in the “settings” interface to enter the server setting interface, as shown in below picture. An available network (WIFI) is required to login the server.

HTTPS:

➤ Enter correct login information in the message input box.


◆ User name and password: the username and password must be the same as the username and password created on the Mediview server.

◆ IP: The user must know the IP address and port of the server. If the server is in the

LAN, the IP address is the IP address of the computer where the server is located. If the server is in the Internet, the corresponding IP address needs to be added.

- ◆ Port: The default server port is 4430, no user modification is required.
- Automatic upload: Click "Auto" button, the icon turns to green, the tested medical record information and test data will be uploaded to the server automatically.
- HTTPS: Switch transmission mode.
- After the above process is completed, click "Login" button to login. It will return to the last interface after login successfully, and the "logged in" icon will display in the top. The login failure message will display after login failed.

HTTP:

- Input the correct login information in the login information input box.
 - ◆ <http://>: Input the server IP (or domain name) and port (such as 127.0.0.1:8088 or service.com:8088).
 - ◆ Path: Service path.
 - ◆ Timeout: Set the timeout for uploading to service.
- Auto: Click the "Auto" button, the tested result will be automatically uploaded to the server.
- After completing the above process, click the "OK" button to complete the setting.

Attention:

1. In "HTTP" mode, the record interface is displayed in a checkable state. If you need to view the record, cancel the checkable state and click to view it.
2. The format of the information scanned by the QR code is "programcode_idtype_idcode_sex_DOB__name"(such as 20_01_XSWSDEWWEEEEEW1237_1_20130809_firstname, last name).

2.8.5 Sleep time setting

Click “Sleep” button in the “settings” interface to enter the standby time setting interface, as shown in below picture.

- Select the sleep time in the standby time drop-down list.
- Click “OK” button to save setting and return to setting interface.
- Click “Cancel” button to Cancel and return to the setting interface.

2.8.6 WIFI Setting

Click “Wifi” button in the “settings” interface to enter the WIFI interface, the WIFI list displays the searched WIFI, as shown in below picture.

- Click “Search” button, refresh the list and re-display the searched WIFI.
- Click “On” to disable or enable the WIFI function.
- Click “Back” to return to the setting interface.

2.8.7 WIFI Connection

Click the WIFI item to connect in the WIFI list, then it will enter the password input interface, as shown in below picture.


Click password input box then input the password, click “Visible Passwords” button to show the password.


- Click “OK” button to connect WIFI, the "WIFI connecting" icon will display in the top when the connection is successful.
- Click "Back" button to cancel and back.

2.8.8 WIFI Configuration

WIFI configuration is mainly used for the configuration setting of connected WIFI. Click the connected WIFI in the WIFI list, the interface jumps to the WIFI setting interface, as shown in below picture.

- Click "Disconnect" button to disconnect WIFI.
- Select "Static" in the WIFI mode drop-down list, then the edit box will be available, as shown in below picture.


- Enter the correct configuration information in the "WIFI configuration input area", click "Commit" button to finish the configuration.
- Click "Back" button to cancel and return to the WIFI list.

Attention:

1. Click WIFI item to connect the WIFI directly, when no password is set.
2. The input information is only as an example.

2.8.9 Music setting


Click "Music" button in the "settings" interface to enter the music settings interface, as shown in below picture.

- Select the music in the music list and listen to it.
- Click "Back" button to return to the setting interface.

2.8.10 Printer

Click "Printer" button in the "settings" interface to enter the printer setting interface, the interface includes printer connection, printer search and paper size setting function.

- Select the required paper width under paper drop-down list, this setting will affect the display of the printed report. It will be effective once the selection is completed, no need to reconnect the printer.
- Place the printer near the device, click "Search" to obtain the name and address of the printer.
- Click "Connect" to connect the printer and wait for the connection.
- After the successful connection of the printer, the interface enters the setting interface, and the printer icon will display in the top.

Print Examples:

1、Paper-80mm:

Vision Screener (V100)

ID name 11/15/2018

	SE	DS	DC	PS	AXIS	XY
OD	-0.25	0.00	-0.25	3.93	117	3,3
OS	+0.25	+0.25	0.00	3.84	18	0,3
PD	53.69 mm					


2、Paper-58mm:

Vision Screener (V100)

ID name 11/15/2018		
	OD	OS
SE	-0.25	+0.25
DS	0.00	+0.25
DC	-0.25	0.00
PS	3.93	3.84
AXIS	117	18
XY	3, 3	0, 3
PD	53.69 mm	

2.8.11 Disk


Click the "Disk" button in the settings interface to see how much disk is used, as shown in the following figure.

Click on the "Save Experimental Pictures" radio button to save the test pictures during the testing. The saved pictures can be used to continuously improve the accuracy of V100. Click the "Clean up Pictures Button" to clean up the experimental pictures. This functionality only exists on higher version of equipment.

2.8.12 Other settings

Click "Others" button in the "settings" interface to enter the relevant interface, the setting interface includes fixation lamp brightness setting, volume, and shooting mode settings, as shown in below picture.

- Click Fixation Lamp progress bar or "+" "-" to set the brightness of the fixation lamp.
- Click Volume progress bar or "+" "-" to set the volume.
- Click on the "Distance Limit" radio button, and the effective shooting distance is limited to the range of 100 ± 6 mm.
- Click on "Test Number" to select the number of test, After select 3, test result will take the average of 3 tests.
- Click "Mode" radio button, choose fast, high precision mode or high precision 2 for shooting. Hi-Precision2 only exists on higher version equipment.
- Click "OK" button to save the setting.
- Click "Cancel" button to cancel and return to the setting interface.

2.8.13 About

Click "About" button in the "settings" interface to enter about interface, as shown in below picture.

➤ Click "Back" button, it will return to the setting interface.

2.9 Diagnostic result reference

V100 will give a reference result judgment, and the basis for its judgment is as follows:

Name	Age(Month) Range and Comparison				
	[6-12)	[12-36)	[36-72)	[72-240)	240
Anisometropia	1.5	1	1	1	1
Astigmatism	2.25	2	1.75	1.5	1.5
Myopia	2	2	1.25	1	0.75
Hyperopia	3.5	3	2.5	2.5	1.5
PupilDiff	1	1	1	1	1
Gaze (ver.)	8	8	8	8	8
Gaze(nasal)	5	5	5	5	5
Gaze (temp)	8	8	8	8	8
Gaze(asymmetric)	8	8	8	8	8

■ Chapter 3 Cleaning, protection and maintenance

Correct, regular cleaning, protection and maintenance, will help to prolong the life-span of the vision screener. Maintenance should be done every 3 months.

The vision screening device is a precise optical instrument, so be sure to handle it lightly.

3.1 Cleaning method

1. Screen cleaning: If dust sticks to the screen, wipe it off with a soft and clean cloth dipped with pure water, or clean it with a screen cleaning cloth. Do not use organic solvents such as alcohol, acetone or other soluble detergents.

Attention: Do not wipe with hard objects.

2. Clean the plastic parts: To clean plastic parts such as the main unit, use a soft and clean cloth dipped with pure water or other soluble detergents to clean the dirt, then wipe it off with medical alcohol.

3. Clean the dust shield of the lens: Use soft and clean cloth dipped with pure water to wipe gently. Do not use the organic solvent such as alcohol, acetone or other soluble detergent.

Attention: Please do not use any corrosive detergent for wiping, in case the surface is damaged.

3.2 Product cleaning cycle

The vision screener should be used in the relevant clean environment, the main parts which need to be cleaned is specified in 3.1. In order to make sure the vision screener works fine, the regular cleaning should be done, it is suggested to clean every 3 month. After cleaning, wait until the equipment is dry before use.

3.3 Protection

1. When the vision screener is not in use, please put the lens cover on.
2. It should be place in the portable box after use to avoid bumps.

3.4 Battery maintenance

1. After the battery is fully charged, the wire should be pulled out, to prevent that the batteries are overcharged or in danger.
2. When the instrument is not in use for a long time, the battery should be removed and placed in the appropriate position of the portable box.
3. When the battery is not used for a long time, the battery should be charged every 3 months. Generally, the battery needs to be kept at 60%-80%.

3.5 Product service life

The service life of the vision screener is 8 years. Since the date of delivery, Shanghai MediWorks Precision Instruments Co., Ltd. provides 1 year warranty.

3.6 Final disposition

1. Battery disposal

The vision screener has the internal electrical power source. For the sake of environmental protection, the battery should be disposed by the professional department.

2. Final disposition of V100

This mark  is on the label. It means the product cannot be disposed like the domestic garbage. It must be disposed according to the local laws and regulations.

If you have any questions, please connect Shanghai MediWorks or its authorized distributors.

■ Corrective maintenance

If a failure occurs, please check according to the following table for guidance. If the problem is still not resolved, please contact Shanghai MediWorks or its authorized distributors.

Fault	Possible cause	Solution
Black screen after starting up for a while	run out of battery	Charge the battery for using
It can't be detected for a long time	The pupils of human eyes become smaller as a result of the bright light in the environment.	Use in a darker environment
The deviation of detection result is too large.	The environment has strong infrared rays	Use in a darker darkroom

■ Exceptions

- Shanghai MediWorks Precision Instruments Co., Ltd. is not responsible for damage caused by fire, earthquake, third party behavior, other accidents, carelessness of the user, misuse, or use under abnormal conditions.

2. Shanghai MediWorks Precision Instruments Co., Ltd. is not responsible for the deficit, bankrupt, and any loss due to unable to use this instrument.
3. Shanghai MediWorks Precision Instruments Co., Ltd. is not liable for any damage to the operation not described in the instruction.
4. Diagnosis is the responsibility of the doctors, and Shanghai MediWorks Precision Instruments Co., Ltd. is not responsible for the results of those diagnoses.

When you contact with us ,please tell us the following information:

Item No.: V100

Serial No. Refer to the product nameplate.

Purchasing time: Please tell us your purchasing date of the instrument.

Malfunction: Please tell us the malfunction details as much as possible.

Service Life: 8 years

1

Version: 1.3

20220414

Shanghai MediWorks Precision Instruments CO., Ltd.


No.7, MingPu Phase II, No.3279 SanLu Road, MinHang District, Shanghai, China

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