
ENDORET[®]

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www.bti-biotechnologyinstitute.com | bti.implantes@bti-implant.es

PLASMATERM H MODEL

MANUAL OF USE AND MAINTENANCE OF THE PLASMATERM H
FORCED AIR OVEN TO KEEP PLASMA FRACTIONS USED IN ENDORET®
TECHNOLOGY AT BODY TEMPERATURE

en

es

de

fr

it

pt

ro

ru


0123

WARNING:

Before using the oven, please read this manual carefully.

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PLASMATERM H MODEL

MANUAL OF USE AND MAINTENANCE OF THE PLASMATERM H FORCED AIR OVEN TO KEEP PLASMA FRACTIONS USED IN ENDORET® TECHNOLOGY AT BODY TEMPERATURE

INTRODUCTION	5
• Preventive measures	5
• Safety	5
• Symbols	5
CHAPTER I: INSTALLATION	6
1.1. inspection	6
1.2. location	6
1.3. power supply	7
1.4. installation	7
CHAPTER II: DESCRIPTION	7
2.1. applications	7
2.2. functional features	8
2.3. technical specifications*	8
2.4. accessories	8
2.5. spare parts	8
2.6. rs232 port	8
CHAPTER III: OPERATION	8
3.1. front panel	8
3.2. container holding tray	9
3.2.1. load	9
3.2.2. inspection and cleaning	9
3.3. operation of the oven	9
3.4. use of hazardous materials	9
CHAPTER IV: MAINTENANCE	10
4.1. technical support	10
4.2. inspection and cleaning	10
4.3. disposal of used equipment	11
APPENDICES	11
A.1. warranty	11
A.2. manufacturer's details	11

INTRODUCTION

PREVENTIVE MEASURES

In order to guarantee the user's safety and to ensure perfect functioning of the forced air oven, the following guidelines and the different warnings included in the handbook should be carefully observed.

DANGER:

- The **PLASMATERM H** forced air oven has been designed to keep at body temperature Plasma Rich in Growth Factors samples used in the Endoret® technology, so its use involves handling of blood fractions.
- If you need to process Endoret® samples from patients with infectious diseases, you should be aware of the potential contamination and risks associated to their manipulation.
- Please take the necessary preventive measures and use adequate disinfectants to avoid potential contamination.
- For increased safety, potential access to electric parts has been considered in the **PLASMATERM H** oven's design. Make sure that the containers do not surpass maximum plasma capacity to avoid potential spills when entering the oven.
- The **PLASMATERM H** forced air oven includes safety systems which avoid equipment's overheating in the event of failure.
- Please do not lean on the oven nor leave any object over it when running.
- If the oven is not used following BTI's guidelines, protection will not be guaranteed; if the device is used in a way different from the one stated by BTI, equipment protection may be damaged.

CAUTION:

Children, pregnancy and breast-feeding: There are no additional preventive measures other than those standard warnings within the Endoret® technology. However, these particular situations should be considered by the consultant before using the oven.

SAFETY

This device, like any other electrical appliance, is a potential source of electromagnetic radiation. It fulfils with Electromagnetic Compatibility requirements regarding emission and immunity for electrical equipment for

measurement, control and laboratory use. (IEC 61326-1).

In any case, before connecting the equipment please check the electromagnetic environment.

CAUTION:

Do not place the oven when running close to strong electromagnetic fields (e.g. deliberately non-screened RF sources) as they might interfere with its adequate performance.

Check for the oven's mains specifications so that it is not connected to an inadequate power supply unit.

Use just the container holding trays provided by default with the **PLASMATERM H** equipment. In the event of damage or loss, please contact your supplier for their respective replacement or repair. Do not use other components.

For good maintenance of the oven, please follow the indications contained in this manual.

WARNING, HOT SURFACE:

The **PLASMATERM H** forced air oven is an electrical appliance for laboratory use to warm and heat different materials, so surfaces exposed to the user might be hot during use. Standard operation guarantees that temperature does not surpass $37 \pm 1^\circ\text{C}$.

SYMBOLS

The symbols present in the equipment and label are:

SYMBOL	DESCRIPTION
	Warning, risk of danger Caution, read the attached information.
	Biological risk.
	Warning, hot surface.
	Danger, risk of electric shock.
	Protection cable terminal.
	EC label according to Council Directive 2007/47/EC concerning medical devices

SYMBOL	DESCRIPTION
	Selective collection of WEEE (waste from electrical and electronic equipment). The equipment must be disposed of separately from domestic waste.
	Serial number
	Catalogue number
	Batch code
	Manufacturer
	Consult Instructions for Use
Rx only	Federal laws restricts this device to sale by or on the order of a licensed healthcare practitioner
	This product has been approved with the requirements of CAN/CSA-C22.2 No 61010-1, second edition, including Amendment 1, or a later version of the same regulation that includes the same level of testing requirements.

WARNING, DANGER, CAUTION, HOT SURFACE and BIOLOGICAL RISK symbols included in this manual are used to stress critical and major instructions.

 **DANGER:** The symbol followed by DANGER points out that its use can endanger the user's personal safety.

 **CAUTION:** The symbol followed by CAUTION tells the operator about a use that might damage the equipment. Moreover, the presence of these symbols advises the user to check the documents attached to the product to know relevant safety warnings or recommendations.

 **WARNING, HOT SURFACE:** This symbol tells the user that the upper surface of the equipment may be hot.

 **BIOLOGICAL RISK:** This symbol tells the user that if the equipment is used with blood samples from patients with infectious diseases, the user's personal safety can be endangered.

CHAPTER I: INSTALLATION

This chapter includes all the necessary instructions for a correct installation of the forced air oven at workplace.

1.1. INSPECTION

Check packaging conditions; if damaged (wet, torn or broken), please contact your supplier. Once the package is open, take out the power cord and the oven, within which both container holding trays are placed, and verify that the content list matches that of the delivered material. If any spare part is missing, it does not correspond to the referred ones in the content list or if it is faulty, please contact your supplier

 **DANGER:**

Be careful during the oven's carriage, as recipient holding trays are free separate parts not attached to the device. Incorrect handling during transportation might cause physical harm to the equipment.

 **CAUTION:**

When carrying the forced air oven, please take it by the bottom casing. Taking the oven by other different part might damage some equipment's component (See Picture 1).

Picture 1


1.2 LOCATION

The surface on which the oven will be placed has to be flat enough to avoid spills. A minimum safety distance of 90 cm. forward and 30 cm. on either sides and behind must be left. Within this safety distance there should not be other apparatus, walls, panels or any object or obstacle. Thus, user's safety and good air circulation can

be guaranteed. The table or work surface on which the oven will be laid must have an accessible power point at maximum 1 m distance. The electric power supply must comply with the specifications for correct oven operation. Refer to Section "2.3 Technical Specifications".


1.3 POWER SUPPLY

The electrical connector is at the left side of rear panel and it comprises the power switch and the fuse box (See *Picture 2*). Before connecting the electric cable to the oven, check if your power supply values (voltage, frequency, etc.) are appropriate. You can validate the values with those specific for the equipment as shown in the rating plate (*Picture 2*) next to the electrical connector or as described in Section 2.3. Ground power supply is obligatory and the oven must be connected to earth.

① Switch ② Fuse Box

Picture 2

③ Rating plate

1.4 INSTALLATION

To install the forced air oven:

1. Once unpackaged and placed at its work site, make sure that safety distances as indicated in Section 1.2 are observed and that all the oven's rubber supports are correctly laid on the work surface.
2. Check that the rear power switch is in "0" position (off).
3. Connect the heating block making sure that power supply specifications as stated in the rating plate or in Section 2.3 are observed.

The forced air oven is correctly installed and ready for use.

 CAUTION:

To move the oven, see *Picture 1*.

Check for the oven's mains specifications so that it is not connected to an inadequate power supply unit.

CHAPTER II: DESCRIPTION

Technical and operational features of the **PLASMATERM H** forced air oven and target applications for which it has been conceived are dealt with next.

2.1. APPLICATIONS

The PLASMATERM H forced air oven has been designed to keep at body temperature Plasma Rich in Growth Factors samples used in the Endoret® technology, so its use involves handling of blood fractions.

The technology involves plasma extraction and further fractionation to get the so called Plasma Rich in Growth Factors to be administered to patients for bone or tissue regeneration in skeletal muscle damage, oral surgery, skin ulcers and aesthetic mesotherapy.

 BIOLOGICAL RISK  CAUTION

2.2. FUNCTIONAL FEATURES.

The **PLASMATERM H** forced air oven is a light and easy to use product thanks to its external plastic casing coated with a special surface treatment which provides high resistance to external aggression. Container holding trays are in anodised aluminium and are highly resistant and easy to clean.

The front information panel includes two light emitting diodes (LED) that show the process stage at which we are.

The working temperature is reached by means of a fan that takes the air inside the oven, conveys it to an aluminium radiator with an embedded heating coil and returns it to the oven's chamber. A microprocessor allows general control of the oven. Furthermore, it monitors the temperature and the protection systems.

2.3 TECHNICAL SPECIFICATIONS*

Laptop **PLASMATERM H** forced air oven:

- Power supply(**)	100 - 240 V
- Maximum frequency	50/60 Hz
- Maximum intensity	1.5A
- Fuse (T type)	2 x T 2A L (5 x 20mm)
- Maximum capacity	8x60ml
- CISPR 11 Class (***)	1B
- Operating temperature	factory set (37±1°C)
- Dimensions (W/L/H)	302 / 201.2 / 175.5 mm
- Weight	3.8 Kg.
- Warranty	24(Twenty four months)

Recommended conditions for use:

- Temperature range	10°C - 37°C
- Relative humidity	<75%

Storage and transport conditions:

- Moisture	< 95%
- Temperature	-2°C to 50°C

* Tech specifications for the following conditions:

- Room temperature between 10°C and 37°C.
- Power range $\pm 10\%$.
- Relative humidity <75%
- For class I installation.
- Pollution degree I

- Indoor use
- Up to 2000 m height

** Ground power supply is necessary and the oven must be connected to earth.

***Household use is allowed under surveillance of a health-care provider.

2.4 ACCESSORIES

No accessories are delivered with the device.

2.5 SPARE PARTS

To deliver any original spare part, your request must include the code or component description and the serial number of the oven.

CAUTION:

BTI is not responsible for the damages that may cause non original spare parts. Please use original spare parts only.

2.6 RS232 PORT


The **PLASMATERM H** forced air oven has a RS232 interface at the rear **for exclusive use of BTI's technical staff**. Please do not connect any other device to this port.

CHAPTER III: OPERATION

To start using the **PLASMATERM H** forced air oven, first you should get used to it and to its front panel, which will display all the necessary information while it is running.

3.1 FRONT PANEL

Picture 3 shows the front panel. For further information, a table including all the components in the device with a short description of each one is enclosed.

1-Symbol :

Optimum temperature symbol.

2-Temperature indicator :

Indicator of optimum operating temperature: it switches on when optimum temperature is reached

3-State indicator :

Indicator of power connection it switches on when the equipment is in operation

4-Symbol :

Symbol of power connection.

Picture 3

In case of breakdown or malfunction, Indicator No. 2 will flash. Look up Section 4.1 in the manual.

In the event of failure of the processor or power shutdown, all the indicator lights will turn off.

3.2. CONTAINER HOLDING TRAY

3.2.1. LOAD

Check that the tray is free of any object (clips, caps, etc.) before placing plasma containers.

Make sure that the containers are correctly closed before entering the oven to avoid potential spills.

3.2.2. INSPECTION AND CLEANING

It is necessary to periodically check the oven and trays to clean them when necessary so as to remove substances that may alter process reliability.

Cleaning indications are described in Section "4.2 Inspection and Cleaning" of the handbook.

Inspection and placement of the oven and its trays are relevant to prevent potential failures in the equipment.

3.3. OPERATION OF THE OVEN.

Start-up of the oven involves the following stages:

1. Set the power switch placed in the electrical connector in "1" position (on). The power-on indicator (green colour) turns on and the equipment starts running.
2. Optimum operating temperature is reached when the optimum temperature indicator (blue colour) turns on.

It is recommended to switch on the device around 25 minutes before use to guarantee optimum operating temperature.

3. The equipment keeps optimum temperature until the door is opened or the device is switched off.

3.4. USE OF HAZARDOUS MATERIALS.

The **PLASMATERM H** forced air oven has been designed to keep at body temperature Plasma Rich in Growth Factors samples used in the **Endoret®** technology, so its use involves handling of blood fractions.

Consequently, there is risk of pathogen microorganisms in blood fractions, so health service professionals handling the product should wear adequate protection to safeguard any potential contamination.

The **PLASMATERM H** forced air oven is designed for its target use only. Any other use not strictly specified in this manual will be the end user's responsibility.

When the oven or any of its components may need maintenance operations or repairing, they should be previously cleaned and disinfected as described in Section 4.2 of this manual.

If you need to process **Endoret®** samples from patients with infectious diseases, you should be aware of the potential contamination and risks associated to their manipulation.

Please take the necessary preventive measures. Make sure that the containers do not surpass maximum plasma capacity to avoid potential spills when entering the oven and use adequate disinfectants to avoid potential contamination in the event of accidental spill. Reference guidelines on biosafety are provided in the Laboratory Biosafety Manual published by the WHO. Moreover, fulfilment of the legislation in force on exposure to biological agents is compulsory.

DANGER:

The **PLASMATERM H** forced air oven has been designed to keep at body temperature plasma fractions to be used in the **Endoret®** technology. When entering products outside the scope of **Endoret®** technology, please verify the correct temperature for the product before using the device.

CHAPTER IV: MAINTENANCE

This chapter describes routine mechanical and electrical maintenance protocols to be periodically addressed in order to guarantee safe operation of the equipment. The user is responsible for fulfilling with next procedures when required and for keeping the oven in good operating conditions to guarantee accuracy in the results. If your purchased goods need any repair or maintenance operation, do not hesitate to contact our technical service in BTI.

4.1. TECHNICAL SUPPORT.

Technical support will be always provided by skilled staff approved by BTI's technical services. If the oven has been used with pathogenic material and needs technical assistance, the user should proceed with the cleaning and disinfection procedure as described in Section 4.2 to guarantee safety to the operator before any repair or maintenance activity.

- *According to the standard EN 61010-2-010 Part.5.4.5, we advise you to contact BTI's technical service to carry out the necessary verification and maintenance operations in the equipment on an annual recommended basis.*

- *Verification actions will be done by skilled technical staff approved by BTI to provide the Maintenance and Technical Support service.*
- *Keep in mind that correct operation of the PLASMATERM H oven greatly depends on the accomplishment of the periodic maintenance and verification activities established by BTI.*

MALFUNCTION:

- *There is power supply, the power switch is on and the power-on indicator (3) is not lit: Contact BTI's technical service.*
- *There is power supply, the power switch is on, the power indicator (3) is lit and the temperature indicator (2) flashes: Contact BTI's technical service. The equipment might be damaged.*

4.2. INSPECTION AND CLEANING.

Before using any cleaning or disinfection method other than those recommended in this manual, users should confirm with BTI that the methods proposed will not cause any harm to the device.

It is recommended to clean the equipment periodically and in the presence of dirt in order to guarantee correct operation of the oven. The parts described below should be adequately disinfected if an accidental spill or breakage of the containers on the equipment occurs, with the corresponding loss of content. Thorough cleaning and disinfection is recommended at least every three months, or even more often depending on use.

OUTER AND INNER PARTS OF THE OVEN.

Clean the outer shell of the oven with a cloth soaked in water and carefully rub it until complete cleaning. Repeat the operation with a wringed cloth soaked in water. Dry with an absorbent cloth. Do not use abrasives. If it is necessary to disinfect the surface because of accidental spills or container breakage, make sure you previously clean and rinse it with water to remove any dirt and then disinfect it with either a 70% ethanol or 2% glutaraldehyde solution.

DANGER:

When cleaning the different components of the equipment, it should be out of operation and disconnected. Disinfection must be performed when the device is cold.

Observe the prevention measures recommended for the cleaning and disinfection products employed.

CONTAINER HOLDING TRAYS.

Since the container holding tray is in light aluminium alloys, common alkaline detergents may be harmful and cause damage to it. **DO NOT USE ALKALINE DETERGENTS** for cleaning it. The tray can be completely removed and taken out from the device for cleaning. Use a neutral non alkaline detergent or an aluminium-specific acid detergent.


DANGER:

Do not use alkaline detergents to clean container holding trays. Use a neutral non alkaline detergent or an aluminium-specific acid detergent.

If it is necessary to disinfect the surface because of accidental spills or container breakage, make sure you previously clean and rinse it with water to remove any dirt and then disinfect it with either a 70% ethanol or 2% glutaraldehyde solution.

4.3. DISPOSAL OF USED EQUIPMENT.

The presence of this symbol in the product, its spare parts or package indicates that it should not be classified as domestic waste. Thus, dispose the device in the closest collection point for electrical and electronic waste recycling. The European Union and some other European countries have different collection systems for used electrical and electronic appliances. Safe elimination of this waste avoids potential environmental and health risks that might occur if the product is not disposed of correctly. Materials recycling contributes to the conservation of natural resources. Consequently, do not dispose of electrical and electronic appliances jointly with domestic waste. For more detailed information on this product, contact your corresponding municipal office, domestic waste disposal service or the place of purchase.

APPENDICES

A.1 WARRANTY

BTI warrants this equipment against any manufacturing defect for a period of 2 years from the date of purchase in the following cases:

- Coverage of any manufacturing defect, including labour costs for location and replacement of faulty parts by our Technical Service.
- The warranty **DOES NOT COVER** failures that in the opinion of BTI's Technical Service derive from incorrect installation, rough handling, misuse or manipulation by people outside BTI's Technical Service.
- It is deliberately excluded any recognition by BTI of compensation for direct or indirect damage of any nature to any person or thing.
- Any equipment to which the rating plate has been removed or modified is considered out of warranty.

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