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Title	Instruction for use of Lancing Device				
Document no.	QS-0621-241			Version	01
Drafted by	Yang Jingzhi	Reviewed by	Jin Jin	Approved by	Cao Yaping
Draft date	2023-01-10	Review date	2023-01-10	Approval date	2023-01-10

Instruction for use of Lancing Device

1. Purpose

The operation instruction is established for the correct use of the product.

Target market: Belarus

2. Scope of application

It's applicable for the Lancing Device with Ejector of SteriLance.

3 Person responsible

User of this operation instruction.

4 Effective date

It becomes effective since January 12, 2023.


Document Change History

Drafted by	Version	Draft date	Amendment content
Yang Jingzhi	01	01/10/2023	Newly added
CC	Technical Dept, Quality Management Dept, Marketing Dept, International Sales Dept		

Instruction for use of Lancing Device


[Product Name] Lancing Device

[Model] LDE 4/LDE 4-AST


[Intended use] Lancing device is used with disposable blood lancet to collect blood from peripheral parts of human body for blood tests.


[Instructions for Use]


1. Wash hands with soap and clean water, rinse and dry thoroughly. Twist off the End Cap of the Lancing Device.


2. Push Lancet firmly into Holder.


3. Remove Lancet Protective Cover.
Note: Save the Protective Cover for safe disposal of used Lancet.


4. Adjust to the suitable step.


5. Gently pull Priming Barrel until a click is heard.

6. Place the End Cap firmly against finger. Press Trigger Button.

7. Remove End Cap and use the saved Lancet Protective Cover to recap Lancet.

8. Push forward the Ejector to dispose of the used Lancet in appropriate container.

[Important] please read carefully all instructions before testing.


Lancing Depth Guide

LDE 4 lancing device has 9 penetration depths and turn Depth Dial to desired lancing depth.

Lancing Depth Guide: number 1 for soft / thin skin, number 3 for average skin, number 5 for thick / calloused skin, and finer adjustments can be made by setting arrow between numbers.

[Cleaning and disinfecting the Lancing Device and End Cap.]

Clean the surface of the lancing device with a soft cloth dipped with clean water. Please don't immerse the body of lancing device into any liquid to avoid damage. The surface of the lancing device can be wiped and disinfected with 75% ethanol. If there is blood stain inside the cap of lancing device, turn the Depth Dial to arrow " $\leftrightarrow$ " position and pull out the Depth Dial to wipe and clean the blood stain with 75% ethanol (as shown in the figure). After disinfection, align the arrow mark on the Depth Dial with the Depth Indicator arrow and install the Depth Dial back to the End Cap.


[Alternate-site testing (LDE4/AST End Cap) AST End Cap is clear and an optional part and it's only suitable for the blood sugar test strips supporting the alternate-site testing sample.]

Alternate-site testing refers to blood collection from other sites except the fingertips, such as the palm, upper arm, and forearm. Please refer to the instruction of blood glucose meter and consult the doctor before trying the alternate-site testing.

Notes for alternate-site testing:

1. The installation of AST End Cap and use of the Lancet are the same as above-mentioned.

2. Place the End Cap firmly against forearm or other AST sites, press the Trigger Button, apply a firm pressure on the Lancing Device until an adequate blood drop is formed.

3. Avoid collecting blood from skin with blood vessel or nevus to avoid excessive blood flow. Avoid blood collection in areas of tendons and bones, such as joints or anklebones.

[Contraindications] Unknown.

[Caution]

For single patient use, Do not share with another person.

Always use a new, sterile lancet. Lancets are for single use only.

Reuse lancets can dull or bend the tip of the lancet, causing damage to your skin, scarring and greater pain.

Not for clinical or multi-patient use.


[Serious Incident notice]

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

[Symbolic interpretation]

Manufacturer

European
Authorized
Representation

Consult instructions
for use

Date of manufacture

Batch code

Caution

Medical Device

[European Authorized Representation]

Emergo Europe B.V.

Prinsessegracht 20, 2514 AP, The Hague, The Netherlands

[Manufacturer]

SteriLance Medical (Suzhou) Inc .

No.168 PuTuoShan Road, New District, Suzhou, Jiangsu 215153, P.R.China

[Revised date.] January 10, 2023 (Version 01)