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No. 168 PuTuoShan Road, New District,215153 Suzhou,Jiangsu ,P.R.China

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[Intended use] The disposable blood lancet is used with lancing device to collect capillary blood from fingertip or alternate sites for blood glucose testing or other testing utilizing small amounts of blood.

Intended patient population: Adults and children over 2 years old

Intended user: The device is intended to be used by professionals and individuals.

Indications: It's suitable for diabetes or blood glucose detection and capillary blood collection.

[Instructions for Use]

1. Load your lancing device and twist off the protective cap.
2. Consult the Instructions for Use of your lancing device.
3. After use, please discard the used lancet in appropriate container.

Compatibility between Disposable blood lancet and Lancing Device: Compatible with LDE, LDE3, LDE4 series lancing device from SteriLance and most other lancing devices.

[Contraindications] Unknown.

[Caution]

- Do not use if protective cap has been previously removed from lancet.
- Do not store lancet in lancing devices.
- Please keep the product away from the children.
- The lancet is for disposable use and do not reuse the lancet. Lancets become dull and painful when reused.
- Do not use the lancet beyond the use-by date.
- Discard the used lancet in appropriate container.

[Sheft Life] 5 years

[Symbolic interpretation]

Sterilized using
irradiation

Do not re-use

Batch code

Use-by date

Manufacturer

European
Authorized
Representation

Notified Body

Caution

Medical Device

Date of manufacture

[European Authorized Representation]

Emergo Europe B.V.
Prinsessegracht 20, 2514 AP, The Hague, The Netherlands

[Manufacturer]

SteriLance Medical (Suzhou) Inc .
No.168 PuTuoShan Road, New District, 215153 Suzhou, Jiangsu, PEOPLE'S REPUBLIC OF CHINA

[Revised date] September 8, 2022 (Version 03)