SeQure™


Snare System Instructions for use	EN1-7
Schlingenkatheter-System Gebrauchsanweisung	DE8-15
Σύστημα αφαίρεσης με βρόχο Οδηγίες χρήσης	EL16-22
Sistema de lazo Instrucciones de uso	ES23-30
Système d'extraction à l'anse Mode d'emploi	FR31-38
Sistema ad ansa Istruzioni per l'uso	IT39-45
Sistema de Captura (Laço) Instruções de Utilização	PT46-52
Система петлевой ловушки Инструкция по применению	RU53-59
Kısaç Sistemi Kullanım Talimatları	TR60-66

Lifetech Scientific (Shenzhen) Co., Ltd.

	English Deutsch Español Ελληνικά Français Italiano Português Русский Türk dili	Date of manufacture Herstellungsdatum Fecha de fabricación Ημερομηνία κατασκευής Date de fabrication Data di fabbricazione Data de fabrica��o Дата изготовления �retim tarihi
	English Deutsch Español Ελληνικά Français Italiano Português Русский T�rk dili	Manufacturer Hersteller Fabricante Κατασκευαστής Fabricant Produttore Fabricante Производитель � retici
	English Deutsch Español Ελληνικά Français Italiano Português Русский T�rk dili	Use by Verwendbar bis Fecha de caducidad Χρήση έως Date de péremption Utilizzare entro Validade Использовать до Son kullanma
	English Deutsch Español Ελληνικά Français Italiano Português Русский T�rk dili	Serial number Seriennummer ��mero de serie Αριθμός σειράς Num�ro de s�rie Numero di serie Serial N.� Серийный номер Seri numarası

	English Deutsch Español Ελληνικά Français Italiano Português Русский Türk dili	Lot Number Chargennummer Número de lote Αριθμός παρτίδας Numéros de lot Numero di lotto Código de lotes Номер партии Parti numarası
	English Deutsch Español Ελληνικά Français Italiano Português Русский Türk dili	Sterilized using ethylene oxide Mit Ethylenoxid sterilisiert Esterilizado con óxido de etileno Αποστειρωμένο με αιθυλενοξειδίο Stérilisé à l'oxyde d'éthylène Sterilizzato con ossido di etilene Óxido de etileno usado esterilizado Стерилизовано этиленоксидом Etilen oksitle sterilize edilmiştir
	English Deutsch Español Ελληνικά Français Italiano Português Русский Türk dili	Non-pyrogenic Pyrogenfrei Apirógeno Μη πυρετογόνο Non pyrogène Apirogeno Não-piragénico Апирогенно Pirojenik değildir
	English Deutsch Español Ελληνικά Français Italiano Português Русский Türk dili	Do not reuse Nicht wiederverwenden No reutilizer Να μην επαναχρησιμοποιείται Ne pas réutiliser Non riutilizzare Não re-utilize Не использовать повторно İkinci kez kullanmayın

	English Deutsch Español Ελληνικά Français Italiano Português Русский Türk dili	Consult instructions for use Gebrauchsanweisung beachten Consultar las instrucciones de uso Συμβουλευτείτε τις οδηγίες χρήσης Consulter le mode d'emploi Consultare le istruzioni per l'uso Instrução de consulta para utilização Обратитесь к инструкции по применению Kullanım Talimatlarına Başvurun
	English Deutsch Español Ελληνικά Français Italiano Português Русский Türk dili	Keep away from sunlight Von Sonnenlicht fernhalten Mantener alejado de la luz solar Φυλάξτε μακριά από το ηλιακό φως Conserver à l'abri de la lumière Tenere lontano dalla luce del sole Mantenha fora da luz do sol Беречь от воздействия солнечного света Güneş ışığından koruyun
	English Deutsch Español Ελληνικά Français Italiano Português Русский Türk dili	Keep dry Trocken halten Mantener seco Να διατηρείται στεγνό Conserver au sec Conservare in luogo asciutto Mantenha seco Хранить в сухом месте Kuru ortamda saklayın
	English Deutsch Español Ελληνικά Français Italiano Português Русский Türk dili	Catalogue number Katalognummer Número de catálogo Αριθμός καταλόγου Référence Numero di catalogo N.º de catálogo Номер по каталогу Katalog numarası

	English Deutsch Español Ελληνικά Français Italiano Português Русский Türk dili	Do not use if package is damaged Bei beschädigter Verpackung nicht verwenden No utilizar si el paquete está dañado Να μην χρησιμοποιείται αν υπάρχει ζημιά στη συσκευασία Ne pas utiliser si l'emballage est endommagé Non utilizzare se la confezione è danneggiata Não use se a embalagem estiver danificada Не использовать в случае повреждения упаковки Paket hasar görmüşse kullanmayın
	English Deutsch Español Ελληνικά Français Italiano Português Русский Türk dili	Do not re-sterilize Nicht erneut sterilisieren No reesterilizar Να μην επαναποστειρώνεται Ne pas restériliser Non risterilizzare Não re-esterilize Не стерилизовать повторно Tekrar sterilize etmeyin
	English Deutsch Español Ελληνικά Français Italiano Português Русский Türk dili	Authorised representative in the European Community Autorisierte Vertretung in der Europäischen Gemeinschaft Representante autorizado en la comunidad Europea Εξουσιοδοτημένος αν Représentant agréé dans la communauté Européenne Rappresentante autorizzato nella Comunità Europea Representante autorizado na Comunidade Europeia Уполномоченный представитель в Европейском сообществе Avrupa Topluluğu'nda bulunan yetkili temsilci

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1 Device Description

- SeQure™ Snare System, including General Snare System, Petite Snare System and Micro Snare System according to different loop sizes and the length of the snare, consists of a snare, a snare catheter, a snare loader and a torquer. The snare is constructed of a nitinol shaft and a loop, and they are connected together on the tapered section of the nitinol shaft by the shrinking compression force of PTFE heat shrink tubing. The loop is made of braided nitinol wires coiled with a Titanium Nitride coated tungsten wire to enhance radiopaque. The plane of the loop is perpendicular to the nitinol shaft. The pre-formed snare loop can be introduced through catheter without risk of deformation because of shape memory and super-elastic properties of nitinol. The snare catheter made of PE tubing and luer connector has a radiopaque markband at its distal tip to enhance the visibility during the operation. The snare is pre-assembled with the snare loader and the torquer, and then packaged with corresponding catheter. (See Figure 1 to Figure 5)


Figure 1 SeQure™ Snare System

Figure 2 Snare


Figure 3 Snare catheter


Figure 4 Snare loader


Figure 5 Torquer

2 Indications and Usages

- The SeSure™ Snare System is intended for use in the cardiovascular system or hollow viscus to retrieve and/or manipulate foreign objects. Manipulation procedures include retrieval and/or repositioning of intravascular foreign objects such as indwelling central venous access catheters, inferior vena cava filters, coils and guidewire pull-through within the cardiovascular system. The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheath should be employed.

3 Contraindications

- This device is not intended for removing foreign bodies that have become entrapped by tissue growth.
- This device should not be used for fibrin sheath stripping in the presence of atrial or ventricular septal defects or Persistent Foramen Ovale.
- This device is not intended for removal of implanted pacing leads.
- Patients who do not tolerate anticoagulation therapy.

- Patients with an active infection.
- Do not use substances that are incompatible with any of the product's components. (See product components section)

4 Warnings

- Pulling forces applied to catheters during fibrin sheath stripping may damage, stretch, or break indwelling catheters 6 French or smaller in diameter. Do not use excessive pull force when attempting fibrin sheath stripping of catheters 6 French or smaller in diameter.
- Excessive forces applied to retrieve and/or manipulate foreign objects may result in damaging the snare system.
- Caution should be observed when using this device for removal of a large fibrin sleeve in order to minimize risk of pulmonary embolism.
- Contents supplied sterile using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Lifetech representative.
- For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

5 Precautions

- Care should be observed when using this device for removal of a large fibrin sheath in order to minimize risk of pulmonary embolism.

6 Potential Adverse Events

- Adverse Events
 - ✧ Testing has indicated that indwelling catheter damage could occur when attempting fibrin sheath stripping on small bore catheters (See warning).
 - ✧ One incident of pulmonary embolism after fibrin sheath stripping has been published (See warning).
- Potential Complications

Potential complications associated with snare retrieval devices in arterial vasculature include but are not limited to

 - ✧ Embolization
 - ✧ Stroke
 - ✧ Myocardial Infarction (depending on placement)

- ✧ Potential complications associated with snare retrieval devices in venous vasculature include but are not limited to:
- ✧ Pulmonary Embolism
- ✧ Other potential complications associated with snare retrieval devices include but are not limited to:
- ✧ Vessel perforation
- ✧ Device entrapment

7 Product Features

● Product Specifications

The SeSure™ Snare System is available with specifications in the following (Table 1).

Table 1 Product specifications of the SeSure™ Snare System

Spec.	Loop Dia. (D1)/mm	Snare Length (L1)/cm	Catheter Size (D2)	Catheter Length (L2)/cm
Snare System-General				
LT-SG-05	5	120	4 Fr	100
LT-SG-10	10	120	4 Fr	100
LT-SG-15	15	120	6 Fr	100
LT-SG-20	20	120	6 Fr	100
LT-SG-25	25	120	6 Fr	100
LT-SG-30	30	120	6 Fr	100
LT-SG-35	35	120	6 Fr	100
Snare System-Petite				
LT-SP-10	10	65	4 Fr	50
LT-SP-25	25	65	6 Fr	50
Snare System-Micro				
LT-SM-02	2	200	3 Fr	175
LT-SM-04	4	200	3 Fr	175
LT-SM-07	7	200	3 Fr	175

Note: 1Fr=1 French=1/3 mm

● Production Components

The following components may come in contact with skin or blood:

SeSure™ Snare System	
Components name	Material
Snare	Nitinol Tungsten with Titanium Nitride coat PTFE
Snare catheter	HDPE Tantalum Markband Barium Sulfate
Loader	HDPE Barium Sulfate
Torquer	POM Brass

8 Directions for Use

- **Assembling the Snare**
 - ✦ Select the appropriate size SeSure™ Snare System for the location of the foreign body. In general, the snare loop should approximate the size of the vessel in which it will be used.
 - ✦ Remove the snare and snare catheter from their protective holders and inspect for any damage.
 - ✦ Remove the loader and steering handle (Torquer) from the proximal end of the snare shaft.
 - ✦ Load the snare into the snare catheter by inserting the proximal (non-looped) end of the snare into the distal (non-hubbed) end of the snare catheter, until the proximal end of the snare shaft exits the hub and the loop end can be retracted into the distal end of the snare catheter.
 - ✦ Inspect device by extending and retracting the snare loop through the distal end of the catheter 2-3 times, while carefully examining the catheter, radiopaque markband and snare for defects or damage.
 - ✦ When appropriate, the snare/snare catheter can be advanced into the vessel as one unit, assembled as described above.
- **Alternative Procedure**
 - ✦ If the snare catheter is already positioned within the vasculature, the provided loader (located on the proximal end of the snare and just distal to the steering handle) may be used to position the snare loop in the indwelling snare catheter. To do so:
 - ✦ Remove the snare from the protective holder and inspect for any damage.
 - ✦ Move the provided loader (located on the proximal end of the snare and just distal to the steering handle) distally until the loop of the snare is enclosed within the

- tubing portion of the loader.
- ✧ Insert the distal end of the loader into the hub of the indwelling snare catheter until resistance is felt. This will indicate the tip of the loader is properly aligned with the inner lumen of the catheter.
- ✧ Keeping the loader as straight as possible, grasp the shaft of the snare just behind the hub of the loader and advance the snare until it is secure within the lumen of the catheter. The loader can be removed at this time by first removing the steering handle and pulling the loader off the proximal end of the snare's shaft.
- **Snare Assisted Manipulation/Retrieval**
 - ✧ Remove the indwelling balloon or delivery catheter, if present. It may be necessary to change or extend the indwelling guidewire to facilitate balloon removal and to upsize the indwelling guide catheter to accommodate the snare.
 - ✧ If a guidewire is indwelling in the patient at the location of the foreign body, cinch the snare over the proximal end of the wire and advance the snare into the guide catheter or sheath until the distal end of the snare catheter is positioned just proximal to the foreign body.
 - ✧ If no guidewire is present, then pull the snare just into the distal end of the snare catheter and advance in the guide catheter until it is positioned just proximal to the foreign body.
 - ✧ Push the snare shaft forward gently to completely open the loop. The loop is then slowly advanced forward, and around the proximal end of the foreign body.
 - ✧ By advancing the snare catheter, the loop of the snare is closed to grab the foreign body. (Note that attempting to close the loop by pulling the Snare within the Snare catheter will move the loop from its position around the foreign body.)
 - ✧ To manipulate a foreign body, maintain tension on the snare catheter to retain the hold on the foreign body, and move the snare and snare catheter together to manipulate the foreign body to the desired position.
 - ✧ To retrieve a foreign body, maintain tension on the snare catheter and move the snare and snare catheter together (proximally) to or into the guide catheter or sheath. The foreign body is then withdrawn through or together with the guiding catheter or vascular sheath. Withdrawal of large foreign bodies may require the insertion of larger sheaths, or guiding catheters, or a cut down at the peripheral site.
- **Snare Assisted Removal of Fibrin Sheaths from Indwelling Catheters**
 - ✧ From a common femoral vein approach, advance the selected snare to the inferior vena cava or right atrium.
 - ✧ Advance a 035 inch (0.89mm) guidewire through the port (distal or venous port if more than one lumen) of the indwelling catheter and into the right atrium or inferior vena cava.
 - ✧ Encircle guidewire with the snare loop.
 - ✧ Advance the snare over the distal end of the catheter to a position proximal to the

fibrin sheath.

- ✧ Close the snare around the catheter and continue applying light traction while gently pulling the snare down (towards distal end of catheter) over the port(s).
- ✧ Repeat steps 4&5 until the catheter is free of fibrin sheath.
- **Snare Assisted Venous Canalization**
 - ✧ Introduce the snare at a patent venous access site and position in the vasculature at the desired venous access site.
 - ✧ Open snare loop to provide a target to guide a puncture into the desired venous access site.
 - ✧ Introduce guidewire through needle and through snare loop.
 - ✧ Close snare over needle by advancing snare catheter.
 - ✧ Pull the guidewire into inferior vena cava.

9 Expiry Date

- The SeQure™ Snare System is sterilized by ethylene oxide. The shelf life is indicated on the label, and the production date and expiry date are marked on the label. Do not use the overdue product.

10 Package and Label

The SeQure™ Snare System is supplied sterile. All the components of the SeQure™ Snare System are fixed and protected with protective holders, then sealed in two Tyvek1073B dialyzing bags, to which a primary label are stuck. The product is sterilized and put in a box with IFU, patient card, customer feedback form and certificate of compliance. The primary label is stuck on the box.

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