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
1 Device Description

The CeraFlex™ Atrial Septal Defect (ASD) Closure Systems include CeraFlex™ ASD Occluders and the SteerEase™ Introducer. The CeraFlex™ ASD Occluders have two types of occluders, CeraFlex™ ASD Occluder and CeraFlex™ Multi-fenestrated (Cribriform) ASD Occluder as shown in **Figure 1**. The CeraFlex™ ASD Occluders are percutaneous, transcatheter closure devices for the non-surgical closure of atrial septal defects. The CeraFlex™ ASD Occluder is self-expandable double disc device, which is made of a nitinol wire mesh that is shaped into two flat discs and a waist between the two discs. PET membranes sewn into each disc (no PET membrane sewn into the waist for multi-fenestrated ASD occluder) help to seal the hole and provide a foundation for growth of tissue over the occluder after placement. All metallic structures are plated with titanium nitride (TiN) to improve the device's biocompatibility.




Figure 1. The Structure of the CeraFlex™ ASD Occluders

The CeraFlex™ ASD Occluders must be used in combination with the SteerEase™ Introducer. The introducer contains a coil reinforced sheath, a dilator, a loader, a haemostatic valve and a delivery cable. It is used to advance the CeraFlex™ ASD Occluder to the proper position. When the occluder is released from the sheath, the disc expands on each side of the defect. At the same time the waist expands to its pre-determined diameter and closes the hole between the left and right atria. The CeraFlex™ ASD Closure Systems consist of an ASD occluder, a Delivery cable, a Loader, a Haemostatic valve, a Sheath and a Dilator. The system (Refer to **Figure 1** and **Figure 3**) is comprised of:

a) Left atrial disc view of CeraFlex™ ASD Occluder

b) Right atrial disc view of CeraFlex™ ASD Occluder

Figure 2. The Structure of the CeraFlex™ ASD Occluders

- ASD Occluder: The device for closing ASD. The CeraFlex™ ASD Occluder is connected with a delivery cable by a loop connection through holes in the proximal portion of the occluder, adjacent to the proximal hub. The CeraFlex™ ASD Occluder is preassembled with delivery cable.
- Sheath: The sheath is used to advance the device to its desired position within the heart.
- Dilator: The dilator is used to ease penetration of tissues and the vessel wall.
- Hemostatic valve: The haemostatic valve at the proximal end of the sheath minimizes bleeding. The sideport with the flexible extension tube and stopcock is used to flush the system.
- Loader: The loader is used to introduce the occluder with the attached delivery cable into the sheath.
- Delivery cable with Plastic vise/Handle (proximal): The delivery cable is used to advance (push) the occluder through the sheath, holding it in position while the sheath is pulled back to deploy the occluder. The Handle and delivery cable are also used to retrieve and/or reposition the occluder if the size, position or the expansion of the occluder is deemed to be unsatisfactory. The handle is used to facilitate direction control and serve as "release control handle" for disconnecting (releasing) the occluder from the delivery cable.


Figure 3. CeraFlex™ ASD Occluder System and the SteerEase™ Introducer

2 Indications and Usage

The CeraFlex™ ASD Occluder is a percutaneous, transcatheter, closure device intended for the occlusion of atrial septal defects (ASDs) and multi-fenestrated atrial septal defects in secundum position or patients who have undergone a fenestrated Fontan procedure and who now require closure of the fenestration. The CeraFlex™ ASD Occluder applies to closing secundum atrial septal defects, and the CeraFlex™ Multi-Fenestrated (Cribiform) ASD Occluder applies to closing multi-fenestrated atrial septal defects. Patients indicated for ASD closure have echocardiographic evidence of (fenestrated) ostium secundum atrial septal defect and clinical evidence of right ventricular volume overload.

3 Contraindications

- Patients who have extensive congenital cardiac anomaly which can only be adequately repaired by cardiac surgery.
- Patients who have sepsis within one month prior to implantation, or any systemic infection that cannot be successfully treated prior to device placement.
- Patients known to have a bleeding disorder, untreated ulcer or any other contraindications to aspirin therapy unless another anti-platelet agent can be administered for 6 months.
- Patients who have demonstrated intra-cardiac thrombi on echocardiography (especially left atrial or left atrial appendage thrombi).
- Patients whose size (i.e., too small to tolerate TEE probe, catheter size, etc) or condition (active infection, etc.) would cause the patient to be a poor candidate for cardiac catheterization.
- Patients where the margins of the defect are less than 5mm from the coronary sinus, AV valves or right upper lobe pulmonary vein (for CeraFlex™ ASD Occluder).
- Any patients for whom the radius of the device is greater than the distance from the central defect to the aortic root or superior vena cava (for the CeraFlex™ Multi-Fenestrated (Cribiform) ASD Occluder).

4 Warnings

- Patients allergic to nickel may suffer an allergic reaction to the device.
- The CeraFlex™ ASD Closure Systems must only be used in hospitals, which have experience with invasive treatment of congenital heart disease, and used by physicians who are trained in transcatheter defect closure techniques.
- Physicians must be prepared to deal with emergency situations which require removal of embolized devices that result in critical hemodynamic compromise. This includes the availability of an onsite surgeon.
- Embolized devices must be removed. Embolized devices should not be withdrawn through intracardiac structures unless they have been adequately collapsed within a sheath.
- The CeraFlex™ ASD Closure Systems are for single use only. Do not reuse or resterilize. Structural integrity and/or function may be impaired or lost through cleaning, resterilization, or reuse and may cause adverse reactions. Lifetech will not be responsible for any direct or consequential damages or expenses resulting from reuse of any of the components in the ASD closure systems.
- Do not release the CeraFlex™ ASD Occluder from the delivery cable if the device does not conform to its original configuration or if the device position is unstable. Recapture the device and redeploy. If still unsatisfactory, recapture the device and replace with a new device.
- Implantation of this device may not eliminate the need for Coumadin in patients with ASD and paradoxical embolism.
- Implantation requires the use of echocardiographic imaging (TTE, TEE or ICE). Do not inflate the sizing balloon, if used, beyond the cessation of the shunt (i.e., stop-flow) or the visualization of a small waist in the balloon. Do not overinflate. (Balloon Sizing may not be suitable for the CeraFlex™ Multi-fenestrated (Cribiform) ASD Occluder.)
- Do not select a device which size is 1.5 times larger than the ASD diameter measured by the echocardiographic – device prior to balloon sizing (for CeraFlex™ ASD Occluder).
- The Occluder is made of NiTi wires, which is a material sensitive to temperature. Therefore, some occluders might present malformation at low temperature. Keeping the temperature of saline water and operating room above 24°C would help occluders to exhibit the desired shape.

A. Cause for peanut malformation:

- a) Over-sized Occluder is selected.
- b) Retract and release the occluder from loader or sheath repeatedly.
- c) The temperature of saline water or the operating room is too low when an occluder is loaded into the loader. This may result in peanut malformation of the occluder, which cannot recover immediately after the occluder is implanted.

B. Cause for tulip malformation:

- a) When a force is applied on the left disc or when the left disc is being released against the left atrial wall, the top of the left disc might be distorted resulting in tulip malformation.
- b) During retraction of the occluder, if the left disc of the occluder is pushed against the left atrial wall or the wall of saline water container, tulip malformation may occur when the left disc is released afterwards.

C. How to prevent malformation:

- a) Some occluders might present malformation at low temperature. Keeping the temperature in the operating room and the saline water above 24 °C would help occluders to exhibit the desired shape.
- b) Pull the left disc of the occluder axially to extend the occluder, and then retrieve the occluder into the loader.
- c) Do not retrieve the occluder into loader against the wall of saline water container.
- d) Do not retract and release the occluder from loader or sheath repeatedly.
- e) Do not release the occluder against the left atrial wall.
- f) The membrane might be jammed when retrieving the occluder into the loader, which is another possible reason for malformation. Pay attention to the condition of the membrane when retracting the occluder into loader.

D. How to correct the device, when malformation is occurs.

- a) If malformation occurs at the left disc, retrieve the occluder, then push the occluder completely out of the sheath in the left atrium, use the tip of sheath to push occluder until the left disc is rested against the left atrial wall to correct the shape, and then release the occluder in the defect area.
- b) If malformation occurs at the right disc, advance the sheath and gently push the tip of sheath against the right disc to correct the shape.
- c) If a peanut or tulip malformation is observed after the occluder is released, retrieve the occluder and immerse it in warm water (above 24 °C). Extend and retract the occluder repeatedly to relax the NiTi wires until the shape of the occluder is recovers.

5 Precautions

- The use of this device in patients with patent foramen ovale has not been studied.
- The CeraFlex™ ASD Closure Systems are for single use only. Do not reuse or resterilize.
- Patient Selection
Certain patients may be at higher risk for complications such as tissue erosion and device embolization. If higher risk patients have devices implanted, closer follow-up is warranted. Higher risk patients include the following:
 - Patients with deformation of the device at the aortic root
 - Patients with high defects (minimal aortic and superior rims)
 - Patients with IVC rim deficiency (risk of device embolization)
- Pre-procedure, procedure and post-procedure
 - Pre-procedure
Aspirin (3-5 mg/kg/day) is to be started at least 24 hours prior to the procedure. In the rare case of aspirin intolerance, two times 200 mg of Ticlopidin are given. Cephalosporin therapy is optional.

- Procedure

Patient should be fully heparinized throughout the procedure with a minimum active clotting time (ACT) of 250 seconds prior to device insertion.

Transesophageal echocardiography (TEE) or similar imaging equipment is recommended as an aid in placing the ASD Occluders in the proper position. If TEE used, the patient's esophageal anatomy must be adequate for placement and manipulation of the TEE probe.
- Post-procedure

Patients are required to take appropriate endocarditis prophylaxis for 6 months following device implantation. The decision to continue endocarditis prophylaxis beyond 6 months is at the discretion of the physician.

Patients, who had any residual shunt, should be inspected by echocardiographic examination for the residual shunt until the defect closed completely.
- MR conditional

A patient with an implanted ASD Occluders can be scanned safely immediately after placement of the device under the following conditions:

 - Static magnetic field of 3 T or less
 - Spatial gradient magnetic field of 720 G/cm or less
 - Maximum MR system-reported, whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

Note: MRI image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the device. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

6 Potential Adverse Events

Placement of the CeraFlex™ ASD Occluder involves using standard interventional cardiac catheterization techniques. The following adverse events might occur from interventional cardiac catheterization techniques or from device.

- | | |
|--|--------------------------|
| ✧ Allergic reaction | ✧ Arrhythmias |
| ✧ Arterial fistula or left or right atrial fistula | ✧ Thrombus formation |
| ✧ Septal tearing | ✧ Device migration |
| ✧ Injury to the nerve or vessel wall | ✧ Cardiac perforation |
| ✧ Incomplete sealing of the defect | ✧ Infection |
| ✧ Mitral insufficiency | ✧ Cardiac tamponade |
| ✧ Transient ischemic attack | ✧ Valvular regurgitation |
| ✧ Fever | |

7 Product Features

The CeraFlex™ ASD Occluders and the introducers are available with specifications in the following (Table 1).


Figure 4. CeraFlex™ ASD Occluders

Table 1. List of specifications of CeraFlex™ the ASD Closure System

Occluders specification	ΦD2	ΦD1	ΦD	L	Recommended introducers
CeraFlex™ ASD Occluder					
LT-ASDf-06	18	14	6	4	SFA8F-f
LT-ASDf-08	20	16	8	4	SFA8F-f
LT-ASDf-10	22	18	10	4	SFA8F-f
LT-ASDf-12	26	22	12	4	SFA9F-f
LT-ASDf-14	28	24	14	4	SFA10F-f
LT-ASDf-16	30	26	16	4	SFA10F-f
LT-ASDf-18	32	28	18	4	SFA10F-f
LT-ASDf-20	34	30	20	4	SFA12F-f
LT-ASDf-22	36	32	22	4	SFA12F-f
LT-ASDf-24	38	34	24	4	SFA12F-f
LT-ASDf-26	40	36	26	4	SFA12F-f
LT-ASDf-28	42	38	28	4	SFA12F-f
LT-ASDf-30	44	40	30	4	SFA14F-f
LT-ASDf-32	46	42	32	4	SFA14F-f
LT-ASDf-34	50	44	34	4	SFA14F-f

LT-ASDf-36	52	46	36	4	SFA14F-f
LT-ASDf-38	54	48	38	4	SFA14F-f
LT-ASDf-40	56	50	40	4	SFA14F-f
LT-ASDf-42	58	52	42	4	SFA14F-f
CeraFlex™ Multi-Fenestrated ASD Occluder					
LT-ASDf-MF-1818	18	18	-	3	SFA9F-f
LT-ASDf-MF-2518	18	25	-	3	SFA10F-f
LT-ASDf-MF-2525	25	25	-	3	SFA10F-f
LT-ASDf-MF-3030	30	30	-	3	SFA12F-f
LT-ASDf-MF-3525	25	35	-	3	SFA14F-f
LT-ASDf-MF-3535	35	35	-	3	SFA14F-f
LT-ASDf-MF-4040	40	40	-	3	SFA14F-f

8 Individualization of Treatment

● Patients Selection

Device placement should only be attempted in those patients with sufficient rim around the defect to allow stable seating of the device.

- Patients with multiple ASDs
 - Closure of multiple ASDs should only be attempted by physicians who have sufficient experience (more than 10-15 cases) to undertake more technically challenging procedures.
 - When two large ASDs separated by more than 7 mm rim of tissue, then implantation of two devices may be justified.
 - When there are multiple ASDs close to each other, one device may be used to cover all defects by placing the device in the largest defect.
- Device placement and size selection
 - Device placement is required to be done with the assistance of TEE or similar imaging equipment.
 - Select an occluder about 2-4 mm in diameter larger than the ASD defect (refer to **Table 1**).
- Use in specific populations
 - Pregnancy – Care must be taken to minimize the radiation exposure to the fetus and the mother.
 - Nursing mothers – Although appropriate biocompatibility testing has been conducted for this implant device, there has been no quantitative assessment of the presence of leachables in breast milk.

9 Directions for Use

- Perform general or local anesthesia. If TEE used, the patient has to be intubated endotracheally.
- Following percutaneous puncture of the femoral vein, perform a standard right heart catheterization.
- Perform an angiogram in order to demonstrate the atrial communication. Catheterize the left atrium using a 45° LAO position and a 35-45°, cranial angulation; inject contrast medium into the left atrium or into the right superior pulmonary vein.
- Introduce a 0.038-inch exchange “J” tip guidewire into the left atrium. If desired, insert a compliant sizing balloon catheter over the exchange wire into the left atrium and measure the defect size.
- Sizing the defect: measure the diameter of the ASD by using a compliant sizing balloon, which is filled with diluted contrast, placed across the ASD, and inflated until a waist in the balloon is visible on fluoroscopic imaging. The balloon is gradually inflated until shunting is eliminated on the Doppler evaluation.
- For multi-fenestrated ASD, device size and placement (the fenestration in which the device will be placed) is based on the location of the fenestration. The use of transesophageal echocardiography or intracardiac echocardiography is mandatory in positioning the sheath in the most centrally-located fenestration. Once the most centrally-located defect has been crossed, measurements should be made to determine the distance to the outer rim of the farthest fenestration, and the device size selected should be 2:1 of this measurement. The device selected must completely cover all fenestrations.
- Once the diameter of the defect has been determined, select an occlusion device equal to or, if the confirmed size is not available, one size larger than the defect. Select the proper size sheath to accommodate the occluder.
- Remove the balloon catheter leaving the exchange guidewire in place.
- Flush the dilator and sheath with saline solution, and then insert the dilator into the sheath, tighten the clamp nut to ensure the dilator fasten with sheath securely. Advance the sheath with dilator over the exchange guidewire through the septal defect into the left superior pulmonary vein. Remove the exchange guidewire and dilator (Refer to **Figure 5**).
- Clockwise screw the loader into the haemostatic valve until it stops, and then check connection is reliable.
- Infuse the saline solution from the stopcock to flush the haemostatic valve and the loader, while keeping the tip of the loader immersed under the saline solution.
- Carefully check the connection between the occluder and the delivery cable (Refer to **Figure 6**). Immerse the device and loader in saline solution and pull the device into the loader with a quick jerking motion (Refer to **Figure 7**).
- Twist the seal cover of the haemostatic valve to suitable degree (Refer to **Figure 8**), and then repeatedly push and pull to make sure the tension is appropriate.

Warning:

If the seal cover is too tight, then the cable will be hard to push. However, if the seal cover is too loose, there will be a risk of aeroembolism.

- Flush the device via the side arm.

Figure 5. Illustration of delivering sheath into heart


Figure 6. Illustration of CeraFlex™ ASD Occluder System


Figure 7. Illustration of loading the device into a loader


Figure 8. Illustration of securing the seal cover of the haemostatic valve

- Keep the sterile heparinized saline outflowing from the loader to eliminate air in the loader. Simultaneously, insert the loader into the haemostatic valve of the sheath gently until it stops, and then tighten the clamp nut of loader. Advance the occluder into the sheath by pushing (not rotating) the delivery cable. **Do not remove the loader from the sheath when advancing the device.**

Warning:

If there is air in the loader and sheath when advancing the device, it will result in an air embolism. If the loader is removed from the sheath when advancing the device, there will also be risk of an air embolism.

- Under fluoroscopic and echocardiographic guidance, deploy the LA disc and waist and pull the device gently against the atrial septum, which maybe felt and also observed by ultrasound imaging. Then immobilize the delivery cable, pull the sheath back and deploy the right atrial disk. Pull the sheath back by approximately 15 mm. A gentle “forward and backward” motion with the delivery cable verifies a secure position across the atrial septal defect, which can also be observed by TTE or TEE (Refer to **Figure 10**).
- Confirm the device position and evaluate residual shunts. Perform angiography or echocardiography. If the position is unsatisfactory or there is obvious residual shunt, stabilize the delivery cable and re-advance the sheath until the device is completely within the sheath. Reposition the device and deploy it again, or remove the device from the patient. If device malformation has occurred, deploy two discs in the left atrium and it may be necessary to push gently against the atrial wall to recover the original device configuration.

Warning:

When withdrawing the occluder back into the sheath, pull back on the delivery cable, and do not push the sheath. Otherwise, the tip of the occluder maybe stuck outside the sheath, and then the occluder cannot be resheathed.

- Release the device when the position is satisfactory. To release the device, press the blue button and move back, while maintaining the orange button press state.
- Hold the sheath to touch the hub of the occluder, and then pull back the cable gently until it separates from the occluder. Finally remove the delivery cable and sheath from the patient (Refer to **Figure 9**).
- Use TTE or TEE to confirm device placement, assess for residual shunts, obstruction, or regurgitation induced by the device. Perform an angiogram to assess for residual flow through the device.


Figure 9. Illustration of releasing the occluder

1. Press the lock button as the arrow of ① (the color of the button is orange).
2. Move the control button as the arrow of ② (the color of the button is blue).

Figure 10. Illustration of implanting the occluder

10 Post-procedure

- All patients should be kept overnight for observation.
- Anticoagulation treatment with heparin for 24 hours.
- Oral anticoagulation with Aspirin 3-5 mg/(kg.d) for children, and 3 mg/(kg.d) for adults, 6 months; conditionally, plus Plavix 75 mg/d for some adult patients.
- Antibiotic therapy
- Re-examination by transthoracic echocardiography (TTE) at 24 hours prior to discharge.
- Patients with high risk should be followed more closely, including the following:
 - Clinical follow-up with echocardiogram one week following device implantation.
 - Education of patients about the higher risk and the need for echocardiography in case of symptoms (i.e., chest pain or shortness of breath).

11 Expiry Date

The ASD Closure Systems are sterilized by ethylene oxide. The shelf life is indicated on the labels, and the manufacturing date and use by date are marked on the label. **Do not** use an expired product.

12 Package and Label

The CeraFlex™ ASD Closure Systems are supplied sterile. The CeraFlex™ ASD Occluder is pre-connected with delivery cable by loop connection and pre-assembled with its loader and Haemostatic valve. The sheath and dilator are packaged separately. The package content is shown in **Table 2**.

Table 2. The packaging contents of the ASD Closure Systems

Package	CeraFlex™ ASD Closure Systems	
	CeraFlex™ ASD Occluders	SteerEase™ Introducers
Contents	an ASD Occluder, a Delivery Cable, a Loader, and a Haemostatic Valve	a Sheath and a Dilator

The ASD device is fixed and protected on a PETG tray, then sealed into inner and outer Tyvek/filmpouches, where a label and a sterilization indicator are attached. The product is sterilized by ethylene oxide gas and put in a box with the Instructions for Use, a patient card, a customer feedback form, and a certificate of compliance. An outer label is attached to the box.

The sheath and dilator of the SteerEase™ Introducer are fixed and protected on a scale board, then sealed in two dialyzing pouches, where a label and a sterilization indicator are attached. The product is sterilized and put in a box with the Instructions for Use, a customer feedback form, and a certificate of compliance. A label is attached to the box.

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1 Device Description

The CeraFlex™PFO Closure System includes CeraFlex™ PFO Occluder and SteerEase™Introducer. The PFO Occluder is shown in **Figure 1**. The PFO Occluder isa percutaneous, transcatheter closure device for the non-surgical closure of patent foramen ovale (PFO). The PFO Occluder is self-expandable double discs device, which is made of a nitinol wire mesh that is shaped into two flat discs and a waist between the two discs. PET membranes sewn into each disc help to seal the hole and provide a foundation for growth of tissue over the occluder after placement. All metallic structures are plated with titanium nitride (TiN) to improve the biocompatibility.



Figure 1.The Structure of the CeraFlex™ PFOOccluders

The CeraFlex™PFO Occluders must be used in combination with the SteerEase™ Introducer. The introducer contains a coil reinforced sheath, a dilator, a loader, a haemostatic valve and a delivery cable. It is used to advance the PFO Occluder to the proper position. When the occluder is released from the sheath, the disc expands on each side of the defect. The PFO Closure System consists of a PFO occluder, a Delivery cable, a Loader, a Haemostatic valve, a Sheath and a Dilator. The system (Refer to **Figure 1** and **Figure 3**) is comprised of:

- PFO Occluders: The device for closing PFO. The CeraFlex™ PFO Occluder is connected with a delivery cable by a loop connection through holes in the proximal portion of the occluder, adjacent to the proximal hub. The CeraFlex™PFO Occluder is preassembled with delivery cable.

a) Left atrial disc view of CeraFlex™ PFOOccluder

b) Right atrial disc view of CeraFlex™ PFOOccluder

Figure 2.The Structure of the CeraFlex™ PFOOccluder

- Sheath: The sheath is used to advance the device to its desired position within the heart.
- Dilator: The dilator is used to ease penetration of tissues and the vessel wall.
- Haemostatic valve: The haemostatic valve at the proximal end of the sheath minimizes bleeding. The sideport with the flexible extension tube and stopcock are used to flush the system.
- Loader: The loader is used to introduce the occluder with the attached delivery cable into the sheath.
- Delivery cable with Plastic vise/Handle(proximal): The delivery cable is used to advance(push) the occluder through the sheath, holding it in position while the sheath is pulled back to deploy the occluder. The Handle and delivery cable are also used to retrieve and/or reposition the occluder, if the size, position or the expansion of the occluder is deemed to be unsatisfactory. The handle is used to facilitate direction control and serve as "release control handle" for disconnecting(releasing) the occluder from the delivery cable.


Figure 3. CeraFlex™ PFO Occluder System and the SteerEase™ Introducer

2 Indications and Usage

- The CeraFlex™ PFO Occluder is a percutaneous transcatheter occlusion device intended for the non-surgical closure of patent foramen ovale.
 - Indications
 - Confirmed PFO by medical examination;
 - Associated with recurrent migraine/headache or TIA or cryptogenic stroke.

3 Contraindications

- Patients known to have intracardiac thrombi, especially left atrial or left atrial appendage thrombi, demonstrated by echocardiography. Patients known to have a bleeding disorder, untreated ulcer or any other contraindications to aspirin therapy, unless another anti-platelet agent can be administered for 6 months.
- Anatomy in which the CeraFlex™ PFO device of required size would interfere with intra-cardiac structure or intra-vascular structure, such as pulmonary vein, aorta root.
- Patients whose heart or vein size is too small to allow TEE probing or catheterization.
- Patients who are in the condition, e.g. active infection, which would cause them to be poor candidates for cardiac catheterization.
- Patients known to have sepsis within one month prior to implantation, or any systemic infection that cannot be successfully treated prior to device placement.
- Patients whose heart does not have enough tissue to secure the device.
- Patients with hypercoagulation disease.
- Any patient for whom the radius of the device is greater than the distance from the patent foramen ovale to the aortic root or superior vena cava.

4 Warnings

- Patients allergic to nickel may suffer an allergic reaction to the device.
- The CeraFlex™ PFO Occluder should only be used by physicians who are trained in transcatheter defect closure.
- Physicians must be prepared to deal with emergency situations which require removal of embolized devices that result in critical hemodynamic compromise. This includes the availability of an on-site surgeon.
- Any incorrect procedure will cause severe complications which include, but are not limited to septum tearing, femoral vein dissection, nerve and vessel wall injury, cardiac perforation, infection, air bubble in a vessel.
- Do not release the CeraFlex™ PFO Occluder from the delivery cable if the device does not reshape to its original configuration or if the device position is unstable. In such cases pull back the occluder into the sheath, reposition and redeploy. If still unsatisfactory, withdraw the occluder and replace with a new one.
- The device is sterilized with ethylene oxide. Do not use if the sterile barrier has been compromised in any way. If damage is found, call your Lifetech representative.
- The expiry date is indicated on the label. Please store in a cool dry place.
- The PFO Occluder Closure System is sterilized with ethylene oxide. The production date and the expiry date are marked on the labels. Do not use expired products.
- For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

5 Precautions

- The PFO Occluder is for single use only. Do not reuse or resterilize.
- Accurate defect sizing is crucial and mandatory for device selection. Device selection should be appropriate based on defect sizing.
- The PFO Occluder must be used in combination with the SteerEase™ SFA Introducer. Do not use any parts of other deliver set.

- Pre-procedure and in-procedure medication and precautions
 - Aspirin (3-5 mg/kg/day) is to be started at least 24 hours prior to the procedure. In the rare case of aspirin intolerance, twice 200 mg of Ticlopidin can be given. Cephalosporin therapy is optional.
 - Patient should be fully heparanized throughout the procedure with a minimum active clotting time (ACT) of 200 seconds prior to device insertion.
 - Transesophageal echocardiography (TEE) or similar imaging equipment must be used during placement of the device.
- Post-procedure medication and precautions
 - Patients should take appropriate endocarditis prophylaxis for 6 months following device implantation. The decision to continue endocarditis prophylaxis beyond 6 months is at the discretion of the physicians.
 - Patients should be treated with anti-platelet/anticoagulation therapy such as aspirin treatment for 6 months post implantation. The decision to continue anti-platelet/anticoagulation therapy beyond 6 months is at the discretion of the physicians.
- MR conditional
A patient with an implanted CeraFlex™ PFO Occluder can be scanned safely immediately after placement of the device under the following conditions:
 - Static magnetic field of 3 T or less
 - Spatial gradient magnetic field of 720 G/cm or less
 - Maximum MR system-reported, whole-body-averaged specific absorption rate (SAR) at 3 W/kg for 15 minutes scanning

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6 Potential Adverse Events

Adverse events that may occur during the procedure and/or device placements include, but are not limited to:

- ✦ Allergic dye reaction
- ✦ Anesthesia reaction
- ✦ Temporary absence of breathing
- ✦ Heart rhythm irregularities
- ✦ Bruising at the groin
- ✦ Dislodgment of the device
- ✦ An air bubble or clot that blocks blood flow in a vessel
- ✦ Incomplete sealing of the defect
- ✦ Infection
- ✦ Injury to the artery, vein or nerves in the groin
- ✦ Perforation of esophagus, vein or heart

7 Product features

The CeraFlex™ PFO Occluder and the SteerEase™ Introducer are packaged separately. The occluder and the introducers are available with specifications in the following (Table 1).


Figure 4. CeraFlex™ PFO Occluder

Table1. List of specifications of the CeraFlex™ PFO Occluders and recommended introducers

CeraFlex™ PFO Occluders specification	Left disc (diameter/mm) ØD2	Right disc (diameter/mm) ØD1	Waist height/mm L	Recommended introducers
LT-PFOf-1818	18	18	3	SFA9F-f
LT-PFOf-2518	18	25		SFA10F-f
LT-PFOf-2525	25	25		SFA10F-f
LT-PFOf-3025	25	30		SFA12F-f
LT-PFOf-3030	30	30		SFA12F-f
LT-PFOf-3525	25	35		SFA14F-f

8 Directions for Use

- Perform a complete TEE study using the following standard imaging locations; transgastric, frontal 4-chamber and basal short axis views. Longitudinal as well as oblique views will also be obtained from these locations.
- Insert a diagnostic catheter through the puncture of the femoral vein, and then navigate it via the vena cava into the right atrium. Then perform an angiogram to visualize the heart and the PFO.

- Measure the pressure and oxygen content in different chambers of the heart. Sizing the defect: insert a balloon catheter through the right atrium into the PFO, and placed across the defect and inflated with diluted contrast medium until the shunt ceases as observed by echocardiography. The balloon is deflated until flow is seen, and then re-inflated until the shunting ceases. Measurements can then be made using fluoroscopy or echocardiographic imaging. Device size selection: Measure the distance from the defect to the aorta root. Measure the distance from the defect to the SVC rim. Select a device with radius of right disc will not exceed the lesser of these two distances.
- Remove the balloon catheter leaving the exchange guidewire in place.
- Flush the dilator and sheath with saline solution, and then insert the dilator into the sheath, tighten the clamp nut to ensure the dilator fasten with sheath securely. Advance the sheath with dilator over the exchange guidewire through the septal defect into the left superior pulmonary vein. Remove the exchange guidewire and dilator (Refer to Figure 5).
- Clockwise screw the loader into the Haemostatic valve until stops, and then check connection is reliable.
- Infuse the saline solution from the stopcock to flush the haemostatic valve and the loader, while keeping the tip of the loader immersed under the saline solution.
- Carefully check the connection between the occluder and the delivery cable (Refer to Figure 6). Immerse the device and loader in saline solution and pull the device into the loader with a quick jerking motion (Refer to Figure 7).
- Twist the seal cover of the haemostatic valve to suitable degree (Refer to Figure 8), and then repeatedly push and pull to make sure the tension is appropriate.

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
Figure 6. Illustration of CeraFlex™ PFO Occluder System

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
Figure 8. Illustration of clockwise screwing the seal cover of Haemostatic valve

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- Under fluoroscopic and echocardiographic guidance, deploy the LA disc and waist and pull the device gently against the atrial septum, which may be felt and also observed by ultrasound imaging. Then immobilize the delivery cable, pull the sheath back and deploy the right atrial disk. Pull the sheath back by approximately 15 cm. A gentle "forward and backward" motion with the delivery cable verifies a secure position across the atrial septal defect, which can also be observed by TTE or TEE.
- Confirm the device position and evaluate residual shunt. Perform angiography or echocardiographic. If the position is unsatisfactory or there is obvious residual shunt, stabilize the delivery cable and re-advance the sheath until the device is completely within the sheath. Reposition the device and deploy it again, or remove the device from the patient. If device malformation has occurred, advance both discs into the left atrium, and apply gentle pressure against the atrial wall to recover the original disk shape.

Warning:

When withdrawing the occluder back into the sheath, pull back on the delivery cable, and do not push the sheath. Otherwise, the tip of the occluder may be stuck outside the sheath, and then the occluder cannot be resheathed.

- Release the device when the position is satisfactory. To release the device, press the blue button and move back, while maintaining the orange button press state.
- Hold the sheath to touch the hub of the occluder, and then pull back the cable gently until it separates from the occluder. Finally remove the delivery cable and sheath from the patient (Refer to Figure 9).

- Use TTE or TEE to confirm device placement, assess for residual shunt, obstruction. Perform an angiogram to assess for residual flow through the device.

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1. Pressing the lock button as the arrow of ① (the color of the button is orange);
2. Moving the control button as the arrow of ② (the color of the button is blue).

- At the end of the procedure the sheath, diagnostic catheters (if used) and the imaging probe (if used) are removed and the procedure is completed.

9 Post-procedure

- Overnight stay for observation required for all patients
- Antibiotic therapy
- Follow up examinations by echocardiography, ECG and Chest X-ray at 24 hours, 1, 3, 6 and 12 months respectively after the operation.

10 Expiry Date

The PFO Closure Systems are sterilized by ethylene oxide. The shelf life is indicated on the labels, and the manufacturing date and use by date are marked on the label. **Do not** use an expired product.

11 Package and Label

The CeraFlex™ PFO Closure Systems are supplied sterile. The CeraFlex™ PFO Occluder is pre-connected with delivery cable by loop connection and pre-assembled with its loader and Haemostatic valve. The sheath and dilator are packaged separately. The package content is shown in **Table 2**.

Table 2. The packaging contents of the PFO Closure Systems

Package	CeraFlex™ PFO Closure Systems	
	CeraFlex™ PFO Occluders	SteerEase™ Introducers
Contents	a PFO Occluder, a Delivery cable, a Loader, and a Haemostatic valve	a Sheath and a Dilator

The PFO device is fixed and protected on a PETG tray, then sealed into inner and outer Tyvek/film pouches, where a label and a sterilization indicator are attached. The product is sterilized by ethylene oxide gas and put in a box with the Instructions for Use, a implant card, a customer feedback form and a certificate of compliance. An outer label is attached to the box.

The sheath and dilator of the SteerEase™ Introducer are fixed and protected on a scale board, then sealed in two Tyvek pouches, where a label and a sterilization indicator are attached. The product is sterilized and put in a box with the Instructions for Use, a customer feedback form, and a certificate of compliance. A label is attached to the box.

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1 Device Description

The CeraFlex™ Patent Ductus Arteriosus (PDA) Closure Systems including CeraFlex™ PDA Occluder and SteerEase™ Introducer. The CeraFlex™ PDA Occluder is a percutaneous transcatheter occluding device for the non-surgical closure of patent ductus arteriosus defects. The device consists of a knitted nitinol wire mesh, which has been heated and set to expand into a T-shape. The disc on the aorta side is designed to hold the plug at the orifice of the ductus arteriosus, while the waist expands to fill and occlude it. APTFE membrane is sewn inside the disc with nylon thread on the aorta side, and two PTFE membranes are sewn inside the waist to help block the abnormal blood flow from the aorta into the pulmonary artery. All metallic structures are plated with a titanium nitride (TiN) bioceramic coating (applied by vacuum plasma deposition technology), to effectively reduce the release of nickel ions into blood and surrounding tissues in order to mitigate thrombogenicity of the device and to speed up endothelial cell growth on the surface of device.


Figure 1. The CeraFlex™ PDA Occluder

The CeraFlex™ PDA Closure System consists of PDA occluder, delivery cable, loader, Haemostatic valve, sheath and dilator. The system (Refer to **Figure 1** and **Figure 2**) is comprised of:

- **PDA Occluder:** A device for closing PDA. The CeraFlex™ PDA Occluder is connected with delivery cable by a loop connection through the holes. The CeraFlex™ PDA Occluder is preassembled with delivery cable.
- **Sheath:** The sheath is used to advance the device to its appropriate position.
- **Dilator:** The dilator is used to ease penetration of tissues and the vessel wall.
- **Haemostatic valve:** The haemostatic valve at the proximal end of the sheath minimizes bleeding. The sideport with the flexible extension tube and stopcock are used to flush the system.
- **Loader:** The loader is used to introduce the occluder with the attached delivery cable into the sheath.
- **Delivery cable with Plastic vise/Handle (proximal):** The delivery cable is used to advance (push) the occluder through the sheath, hold it in position while the sheath is pulled back to deploy the occluder. The Handle and delivery cable are also used to retrieve and/or reposition the occluder, if the size, position or the expansion of the occluder is deemed to be unsatisfactory. The handle is used to facilitate direction control and serve as "release control handle" for disconnecting (releasing) the occluder from the delivery cable.


Figure 2. CeraFlex™ PDA Occluder and the SteerEase™ Introducer

2 Indications and Usage

The CeraFlex™ PDA Occluder is a percutaneous, transcatheter closure device that is intended for the non-surgical closure of patent ductus arteriosus (PDA).

Indications:

- Patients have PDA or PDA accompanies with other mild cardiac disease;
- Patients' weight more than 6 kgs, and 6 months old and above;
- The narrowest portion of the PDA is 2 mm or larger.

Contraindications:

- Patients have PDA which is needed for patients' survival due to other cardiac anomalies;
- Patients' weight less than 6 kgs, or less than 6 months old;

- Presence of thrombus at the intended site of implant, or documented evidence of venous thrombus in the vessels through which access to the defect is gained;
- Active endocarditic or other infections producing bacterium;
- Patients' vasculature, through which access to the defect, is inadequate to accommodate the appropriate sheath size;
- Patients with pulmonary hypertension higher than 0.4 Rp/Rs or with pulmonary vascular resistance more than 8 woods units.

4 Warnings:

- The device should be removed if more than half of the left pulmonary artery lumen is obstructed by the device, or if more than 3 mm length extends into the pulmonary artery.
- Patients allergic to nickel may suffer an allergic reaction to the device.
- The PDA Closure System must only be used in hospitals, which have experience with invasive treatment of congenital heart disease, and used by physicians who are trained in transcatheter defect closure techniques.
- Physicians must be prepared to deal with emergency cases which require removal of embolized devices that result in a critical hemodynamic compromise. This includes the availability of an onsite surgeon.
- Embolized devices must be withdrawn via the delivery track within a sheath.
- The PDA Closure System is for single use only. Do not reuse or resterilize. Structural integrity and/or function may be impaired or be lost through cleaning, resterilization, or reuse and may cause adverse patient reactions. Lifetech will not be responsible for any direct or consequential damages or expenses resulting from reuse of any of the components in the PDA Closure System.
- The PDA Occluder should not be released from the delivery cable if the device does not match with the defect or if the device position is unstable. Withdraw the device into the introduce sheath and redeploy it. If positioning is still unsatisfactory, withdraw the device completely and replace with a new one.
- Ensure the connection safety before loading the device into the loader.
- Retrieve the device when it is hard to release the device from the delivery cable.

5 Precautions:

- Both the PDA Occluder and the SteerEase™ Introducer are for single use only, and the valid period of sterilization is indicated on the labels. Do not reuse or resterilize.
- Post-procedure
 - Patients should take appropriate endocarditis prophylaxis for 6 months following the device implantation. The decision to continue endocarditis prophylaxis beyond 6 months can be made by physicians.
 - Lung perfusion scan should be done if flow speed is greater than 3m/s, or if the Z-score is -2 for the left pulmonary artery diameter.
- MR conditional
A patient with an implanted PDA Occluders can be scanned safely immediately after placement of the device under the following conditions:
 - Static magnetic field of 3 T or less
 - Spatial gradient magnetic field of 720 G/cm or less
 - Maximum MR system-reported, whole-body-averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning

Note: MRI image quality may be compromised if the scan area is in the same area or relatively close to the position of the device. Therefore, optimization of MR imaging parameters to compensate for the presence of this device is necessary.

- Use in specific populations
- ✧ Pregnancy – Care must be taken to minimize the radiation exposure to the fetus and the mother.
- ✧ Nursing mothers – Although appropriate biocompatibility testing has been conducted for this implant device, there has been no quantitative assessment of the presence of leachables in breast milk.

6 Potential Adverse Events

Placement of the CeraFlex™ PDA Occluders involves using standard interventional cardiac catheterization techniques. The following adverse events might be expected from interventional cardiac catheterization techniques.

- | | |
|--------------------------------------|-----------------------|
| ✧ Allergic reaction | ✧ Arrhythmias |
| ✧ Cardiac tamponade | ✧ Thrombus formation |
| ✧ Septal tearing | ✧ Device migration |
| ✧ Injury to the nerve or vessel wall | ✧ Cardiac perforation |
| ✧ Incomplete sealing of the defect | ✧ Infection |
| ✧ Pulmonary artery dissection | ✧ Fever |

7 Product Features

The CeraFlex™ PDA Closure system is available with specifications in the following table (Table 1).

Table 1.List of the available PDA Occluders and recommended introducers

Spec.	ΦD/mm	ΦD2/mm	ΦD1/mm	L/mm	Recommended introducers
CeraFlex™ PDA Occluders					SteerEase™ Introducers
LT-PDAf-0406	10	6	4	7	SFP6F-f
LT-PDAf-0608	12	8	6	7	SFP7F-f
LT-PDAf-0810	14	10	8	7	SFP7F-f
LT-PDAf-1012	16	12	10	7	SFP8F-f
LT-PDAf-1214	20	14	12	7	SFP9F-f
LT-PDAf-1416	22	16	14	8	SFP9F-f
LT-PDAf-1618	24	18	16	8	SFP10F-f
LT-PDAf-1820	26	20	18	9	SFP12F-f
LT-PDAf-2022	28	22	20	9	SFP12F-f
LT-PDAf-2224	30	24	22	10	SFP14F-f


Figure 3.CeraFlex™ PDA Occluder

8 Directions for Use

- Perform a right heart catheterization in a routine way.
- Use angiography to size the heart defect. There are two options for angiography performance of the patent ductus arteriosus. The first is to introduce an exchange guidewire through the ductus and pass a pigtail catheter with side holes in the communication. Perform a biplane angiogram to opacify the PDA. The second option is to pass a pigtail catheter into the proximal descending aorta via the femoral artery and perform the biplane angiogram to opacify the PDA. Under the angiography, measure ductal size both on the pulmonary and aortic ends of the duct.
- Select the proper size occluder. Select a CeraFlex™ PDA Occluder based on the smallest diameter measured in the PDA. It is recommended to select a device of which the smaller end is at least 2 mm larger than the narrowest portion of the PDA.
- Introduce a 0.038-inch J-tipped exchange guide wire.
- Flush the dilator, sheath with saline solution, and then insert the dilator into the sheath, tighten the clamp nut to ensure the dilator fasten with sheath securely. Advance the introducer sheath with dilator over the exchange wire via the pulmonary artery into the aorta and position the sheath in the descending aorta while withdrawing the dilator.
- Clockwise screw the loader into the Haemostatic valve until it stops, and then check connection is reliable.
- Infuse the saline solution from the stopcock to flush the haemostatic valve and the loader, while keeping the tip of the loader immersed under the saline solution.
- Carefully check the connection between the occluder and the delivery cable (Refer to **Figure 4**). Immerse the device and the loader in saline solution and pull the CeraFlex™ PDA Occluder into the loader with a quick jerking motion (Refer to **Figure 5**).
- Twist the seal cover of the haemostatic valve to suitable degree (Refer to **Figure 6**), and then repeatedly push and pull to make sure the tension is appropriate.

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- Inject sterile heparinized saline from stopcock of the loader to eliminate air from the loader.


Figure 4. Illustration of CeraFlex™ PDA Occluder System


Figure 5. Illustration of loading the device into a loader


Figure 6. Illustration of clockwise screwing the seal cover of haemostatic valve

- Keep the sterile heparinized saline outflowing from the loader to eliminate air of the loader. Simultaneously, insert the loader into the Hemostatic Valve of the sheath gently until stop, and then tighten the clamp nut of loader. Advance the occluder into the descending aorta by delivery cable without any rotation. **Do not remove the loader from the sheath when advancing the device.**

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- Deploy the retention skirt only and pull firmly against the orifice of the PDA. This can be observed under fluoroscopy, or it can be clearly felt as a tugging sensation in synchrony with the aortic pulsation. The position of the device is confirmed with repeated angiograms in the aorta using the pigtail catheter. The device can be adjusted until the retention skirt is well seated at the orifice. Retract the delivery sheath and deploy the cylindrical portion of the device securely in the patent ductus arteriosus while applying slight tension.
- Perform an aortography to verify the correct position of the device. If the position is unsatisfactory, stabilize the sheath and pull the delivery cable until the device is completely within the sheath. It can be repositioned and deployed, or removed completely from the patient.

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- Release the device when the position is satisfactory. To release the device, press the blue button and move back, while maintaining the orange button press state.
- Hold the sheath to touch the hub of the occluder, and then pull back the cable gently until it separates from the occluder. Finally remove the delivery cable and sheath from the patient. (Refer to **Figure 7**)
- Repeat aortography.


Figure 7. Illustration of releasing the occluder

1. Pressing the lock button as the arrow of ① (the color of the button is orange);
2. Moving the control button as the arrow of ② (the color of the button is blue)

Warning:
Remove the device if over 3 mm length extends into the pulmonary artery, or if more than half of the left pulmonary artery lumen is obstructed by the device. In questionable cases, perform TTE before releasing the device with Doppler measurement of left pulmonary artery flow velocity. The device should be removed if the left pulmonary artery flow is greater than 3.0 m/s (or greater than 75% of the LPA velocity before cardiac catheterization).

Note:
The device should not be released from the delivery cable if the device does not match the defect or if device position is unstable. Withdraw the device into the sheath and redeploy. If still unsatisfactory, withdraw the device completely and replace it with a new one.

If the CeraFlex™ occluder becomes detached from the delivery cable during loading, it cannot be used and must be discarded.

9 Post-procedure

- All patients should be kept overnight for observation.
- Antibiotic therapy
- Re-examination by TTE at 24 hours prior to discharge

10 Expiry Date

Both the CeraFlex™PDA Occluder and the SteerEase™ Introducer are sterilized with ethylene oxide. The shelf life is indicated on the label. The manufacturing date and use by date are marked on the label. **Do not** use an expired product.

11 Package and Label

The CeraFlex™ PDA Closure Systems are supplied sterile. The CeraFlex™PDA Occluder is pre-connected with delivery cable by loop connection and pre-assembled with its Loader and Haemostatic valve. The sheath and dilator are packaged separately. The package content is shown in **Table 2**.

Table 2. The packaging contents of the PDA Closure Systems

Package	CeraFlex™ PDA Closure Systems	
	CeraFlex™ PDA Occluders	SteerEase™ Introducers
Contents	a PDA Occluder, a Delivery Cable, a Loader, and a Haemostatic Valve	a Sheath and a Dilator

The PDA device is fixed and protected on a PETG tray, then sealed into inner and outer Tyvek/film pouches, where a label and a sterilization indicator are attached. The product is sterilized by ethylene oxide gas and put in a box with the Instructions for Use, a patient card, a customer feedback form and a certificate of compliance. An outer label is attached to the box.

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Instructions for Use

LAmbre™

LAmbre™ Left Atrial Appendage Closure System

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I. DEVICE DESCRIPTION

The Left Atrial Appendage Closure System is a transcatheter, self-expanding device intended for cardiac structures not involving the septal wall, which require closure or occlusion. The System consists of a Left Atrial Appendage Occluder (**Figure 1**) and a Delivery System (**Figure 2**). The Left Atrial Appendage Occluder consists of an Umbrella and a Cover connecting with each other by a Connecting Sleeve. The Umbrella is made of a nitinol elastic framework with hooks on it to secure the device to the walls of left atrial appendage (LAA). The Cover is made of a nitinol elastic mesh shaped into a flat disc with a recessed hub located at the middle for connecting to the Delivery Cable, which is used to advance the LAA occluder through the Delivery Sheath to the intended position. Both of the Umbrella and the Cover are sewn with PET membranes to prevent entry of blood into LAA after placement. The Delivery System comprises a Delivery Sheath (**Figure 3**), a Dilator (**Figure 4**), a Delivery Cable (**Figure 5**), a Loader (**Figure 6**) and two Haemostatic valves (**Figure 7**).





Figure 1 Left Atrial Appendage Occluder

Work assembly of the Delivery Cable, the Delivery Sheath, the Loader and the Haemostatic Valve


Work assembly of the Delivery Sheath and the Haemostatic Valve


Work assembly of the Delivery Sheath and the Dilator


SL8F/9F/10F45x30-900

Work assembly of the Delivery Cable, the Delivery Sheath, the Loader and the Haemostatic Valve


Work assembly of the Delivery Sheath and the Haemostatic Valve

Work assembly of the Delivery Sheath and the Dilator

SL8F/9F/10F45-900

Figure 2 Delivery System

NOTE: ID=Inner diameter, OD=Outer diameter

Figure 3 Delivery Sheath


Figure 4 Dilator

Figure 5 Delivery Cable


Figure 6 Loader


Figure 7 Haemostatic Valve

Table 1 REF of Left Atrial Appendage Closure System UNIT: mm

REF of Device	D	D1	REF of Delivery System	L2	L3	L4
LT-LAA-1622	16	22	SL8F45×30-900/SL8F45-900(ID=2.7,OD=3.5);	900	920	1240
LT-LAA-1824	18	24	SL9F45×30-900/SL9F45-900(ID=3.0,OD=3.8);			
LT-LAA-2026	20	26	SL10F45×30-900/SL10F45-900(ID=3.4,OD=4.1)			
LT-LAA-2228	22	28	SL9F45×30-900/SL9F45-900(ID=3.0,OD=3.8);			
LT-LAA-2430	24	30	SL10F45×30-900/SL10F45-900(ID=3.4,OD=4.1)			
LT-LAA-1630	16	30	SL10F45×30-900/SL10F45-900(ID=3.4,OD=4.1)			
LT-LAA-2632	26	32				
LT-LAA-2834	28	34				
LT-LAA-3036	30	36				
LT-LAA-3236	32	36				
LT-LAA-3438	34	38				
LT-LAA-3640	36	40				
LT-LAA-1832	18	32				
LT-LAA-2032	20	32				
LT-LAA-2234	22	34				
LT-LAA-2436	24	36				
LT-LAA-2638	26	38				

Note: Generally SL8/9/10F45x30-900 are used in procedure, SL8/9/10F45-900 with one 45°curve could be used in special LAA anatomy as determined by the physician.

II. INDICATIONS

The Left Atrial Appendage Closure System is a percutaneous transcatheter device intended to prevent thrombus embolization from the left atrial appendage (LAA) in patients who have nonvalvular atrial fibrillation.

III. CONTRAINDICATIONS

The Left Atrial Appendage Closure System is contraindicated for the following:

- Patients' LAA anatomy not suitable for the **REF** of the device.
- Patients' with intracardiac thrombus.
- Patients with active endocarditis or other infections causing bacteremia.
- Patients where placement of the device would interfere with any intracardiac or intravascular structures.
- Patients with contraindications to X-Ray and/or trans-esophageal echocardiographic examinations.
- Patients with known hypersensitivity to nickel.
- Other percutaneous, general contraindications for cardiac catheterization. For example, the patient's height(short stature not suitable for TEE probe, catheter size, etc.) or physical condition (active infection, hemorrhagic disease, untreated ulcer, etc.).

IV. WARNINGS

- The Left Atrial Appendage Closure System should only be used by physicians who have been trained in transcatheter techniques. The physician should determine which patients are suitable candidates for the procedure using this device according to the indication. Lifetech will not support and be responsible for any direct or consequential damages or expenses resulting from using this device by untrained personnel.
- Potential for device embolization exists. Physicians should be prepared to deal with urgent situations that may require the removal of the device.
- If resistance occurs during procedure, the physician should check the reason before continuing the procedure.
- During operation, if there is a need for the LAA occluder to be removed or replaced, the device and the Delivery Sheath must both be removed and replaced. Otherwise it may result in the Delivery Sheath and/or device malfunctioning.
- Patients who are allergic to nickel may have an allergic reaction to this device.
- **DO NOT use the device if the packaging sterile barrier is open or damaged.**

- The device is sterilized using ethylene oxide and is for single use only. **DO NOT** reuse or resterilize. Attempts to resterilize the device may result in device malfunction, inadequate sterilization, or patient harm.
- Use the device before the "USE by" date specified on the label. Read all instructions and check carefully of each component before using. Do not use if any damage observed.
- Research has not been conducted for patients under 18 years of age.

V. PRECAUTIONS

- Store the System in a cool, dry place.
- The device should not be sterilized by any autoclaving. The upper limit for relative humidity is 80%, Do not constantly keep exposure to the temperature above 54°C.
- Do not contact with any organic solvent.
- Be careful in handling the System. No squeeze, wet, exposure in the sun during transportation.
- Please contact the sales agent or the manufacturer directly if any questions arised before using the device.
- MRI Compatibility: Through non-clinical testing, a patient with an LAA occluder implanted can be scanned safely immediately after placement of the device under the following conditions:
 - Static magnetic field of 3Tesla or less.
 - Spatial gradient magnetic field of 720G/cm or less.
 - Maximum MR system-reported, whole-body-averaged specific absorption rate (SAR) of 3W/kg for 15 minutes of scanning.

***Note:** MRI image quality may be compromised if the area of interest is in the exact same area of relatively close to the position of the device.*

- The physician should exercise clinical judgment in situations that involve the use of anticoagulants or anti-platelet drugs before, during, and/or after the use of this device.
- The physician should exercise caution if implanting a device in a patient who has an implantable cardioverter defibrillator (ICD) and leads.
- The physician should ensure that the vasculature is adequate for the delivery sheath size being selected.

- The physician should exercise caution if performing arrhythmia ablation at or near the implant site after the device is implanted.
- The physician should use standard interventional cardiovascular catheterization techniques, such as transseptal puncture, during LAA occlusion.
- Adopt a prudent policy for using of the device in the pregnancies and the breastfeeding mothers due to the radiation exposure and medication use during procedure.

VI. COMPLICATIONS

Complications may occur during or after a procedure placing this device include, but are not limited to:

Air embolism, allergic reaction, anesthesia reactions, arrhythmia, bleeding, cardiac arrest, cardiac tamponade, death, device embolization, device migration, embolic event, fever, foreign body embolization, hypertension/hypotension, infection, multi-organ failure, myocardial infarction, cardiac perforation, pericardial effusion, renal failure/dysfunction, epilepsy, significant incomplete LAA occlusion, stroke or transient ischemic attack, thrombus formation on the device, valvular regurgitation/insufficiency, vascular access site injury, vessel trauma/injury, cardiac failure, pseudoaneurysm, pleural effusion.

VII. DIRECTIONS FOR USE

1. Prepare manifold: heparinized saline for flush, pressure monitoring kit, injection contrast.
2. Prepare the patient
 - GA (general anesthesia) with TEE guidance
 - Standard percutaneous puncture at femoral vein. Then insert 0.035" guidewire for the advancement of transseptal sheath.
 - Puncture at infero-posterior site of inter-atrial septum by the standard transseptal puncture kit. After a transseptal puncture is performed, remove the transseptal needle, dilator and guidewire. Maintain a recommended activated clotting time (ACT) of 250 seconds until the procedure is completed by giving intravenous heparin.


Note: *All percutaneous puncture kit and transseptal puncture kit should be sterile before use, and adequately flush by heparinized saline.*

3. Gently advance a stiff guidewire(0.035", length 260cm) into the upper pulmonary vein through the transseptal sheath.
4. Prepare the Delivery System: take out the Delivery System from the sterile packaging with aseptic

technique, flush all the components with heparinized saline prior to use in order to remove air within the system.

5. Remove the transeptal sheath, exchange it with the Delivery Sheath with the Dilator through the stiff guidewire, and then remove the Dilator. Connect the Haemostatic Valve to the end of the Delivery Sheath as shown in **Figure 8**. Allow bleed back from the Haemostatic Valve to remove air in the system before close the Haemostatic Valve.

Note: Remove the air in the Delivery Sheath and the Dilator before use. Carefully and gently advance the Delivery Sheath with the Dilator to prevent trauma to LAA.

- ① Insert the distal of the Haemostatic Valve to the end of the Delivery Sheath
- ② Screw tightly the Haemostatic Valve
- ③ Gently open the Haemostatic Valve for bleed back
- ④ Screw tightly the Haemostatic Valve again afterwards.

Figure 8 Screw the Haemostatic Valve to the end of the Delivery System

6. Open the Haemostatic Valve, insert a 5F, marked pigtail catheter to the left upper pulmonary vein through the stiff guidewire. Remove the stiff guidewire afterwards.
7. Pull back the Delivery Sheath together with the pigtail catheter en-bloc from the pulmonary vein. Use the pigtail catheter to engage into LAA. Advance the Delivery Sheath to the proximal part of LAA through the pigtail catheter. Then screw tightly the Haemostatic Valve.

Note: Position the pigtail catheter at distal part of LAA.

8. LAA angiogram at RAO 20 caudal 20.
9. Use both Transesophageal echocardiography (TEE) and angiography to measure the landing zone for the device.

Note: Generally recommend Transesophageal echocardiography (TEE) and angiography at multi-angle to measure. Choose the maximum of all measurement landing zone sizes for selecting the device.

10. Choose a device that is 3-8mm larger than the landing zone diameter. See table 2 and table 3 to determine if an appropriate device size to occlude the structure is available.

Warning: Do not implant the device if the measurements of the structure to be occluded do not fall within the sizing chart in table 2 and table 3.

Note: The device sizing in table 2 and table 3 is based on angiographic measurements.

Ls: the depth of LAA



Ds: the size of the landing zone

Lo: the distance between the orifice of LAA and the crotch for double lobes

Table 2 The device sizing chart for single lobe LAA UNIT: mm

Measure		Device	
Ds	Ls	D	REF
11.0~12.5	≧ 10	16	LT-LAA-1622
12.5~14.0	≧ 10	18	LT-LAA-1824
14.0~16.0	≧ 10	20	LT-LAA-2026
16.0~18.0	≧ 10	22	LT-LAA-2228
18.0~20.0	≧ 10	24	LT-LAA-2430
20.0~22.0	≧ 10	26	LT-LAA-2632
22.0~24.0	≧ 10	28	LT-LAA-2834
24.0~26.0	≧ 10	30	LT-LAA-3036
24.0~26.0	≧ 10	32	LT-LAA-3236
26.0~28.0	≧ 10	34	LT-LAA-3438
28.0~30.0	≧ 10	36	LT-LAA-3640


Table 3 The device sizing chart for multiple lobes LAA UNIT: mm

				
Measure			Device	
Ds	L0	Ls	D	REF
11.0~12.5	<10	≧ 10	16	LT-LAA-1630
12.5~14.0	<10	≧ 10	18	LT-LAA-1832
14.0~16.0	<10	≧ 10	20	LT-LAA-2032
16.0~18.0	<10	≧ 10	22	LT-LAA-2234
18.0~20.0	<10	≧ 10	24	LT-LAA-2436
20.0~22.0	<10	≧ 10	26	LT-LAA-2638

Note: when $L_0 \geq 10\text{mm}$, the device sizing as table 2. The figure of LAA in this table only shows double lobes. For the multiple lobes the same select method can be used.

11. Prepare Left Atrial Appendage Occluder: take out the LAA occluder from sterile packaging under sterile conditions, completely immerse the occluder in heparinized saline prior to use.
12. As shown in **Figure 9**, screw tightly between the Loader and the Haemostatic Valve. As **Figure 10** insert the Delivery Cable into the Loader until its tip is revealed from the distal the Loader. As **Figure 11** connect the Delivery Cable to the hub of the LAA occluder by screwing in a clockwise manner, then, completely immerse the occluder in heparinized saline and pull the Delivery Cable until the occluder is fully retracted within the Loader.

Note: Exchange for another new LAA Closure System if the hub is loose or not match when connecting the occluder to the Delivery Cable.

- ① Insert the Haemostatic Valve into the proximal of the Loader
- ② Screw tightly the distal nut of the Haemostatic Valve

Figure 9 Screw the Haemostatic Valve onto the Loader



Figure 10 Insert the Delivery Cable into the Loader

- ① Connect the occluder with the Deliver Cable, after tightening, turn back half a turn
- ② Completely immerse the occluder in heparinized saline, then pull the occluder into the Loader completely

Figure 11 Screw the occluder onto the Delivery Cable

13. Remove the air in the system: immerse the device in a bowel of heparinized saline, forcefully flush the device in the Loader via 50cc Luer lock syringe for several times until no air coming out from the system. Then as **Figure 12** screw tightly the adjust nut at the proximal of the Haemostatic Valve.



Figure 12 Screw tightly the adjust nut at the proximal of the Haemostatic Valve after remove the air

***Note:** The product should be fully flushed until no air bubbles exist in the loader. Remember to close the three-way valve and screw tightly the Haemostatic Valve after flush. The following action*

can be repeated to help to remove the air in the system: immerse the distal of the Loader into heparinized saline adequately; push the occluder and pull the occluder again into the Loader.


14. Connect the three-way valve with the pressured heparinized saline. Slight unscrew the Haemostatic Valve. Confirm heparinized saline flow out from the distal of the Loader and the proximal of the Haemostatic Valve.

15. Remove the pigtail catheter. As **Figure 13** withdraw the Haemostatic Valve at the proximal of the Delivery Sheath. As **Figure 14** immediately insert the Loader with the occluder attaching to the Delivery Cable into the Delivery Sheath. Screw tightly the distal nut of the Loader.

- ① Unscrew the distal nut of the Haemostatic Valve
- ② Remove the Haemostatic Valve

Figure 13 Remove the Haemostatic Valve


- ① Insert the Loader with the occluder attaching to the Delivery Cable into the Delivery Sheath
- ② Screw tightly the distal nut of the Loader

Figure 14 Connect the Loader with the Delivery Sheath

Note: During insert the Loader into the Delivery Sheath, keep flushing the loader with continuous pressurized heparinized saline from the manifold and make the heparinized saline flow out of the distal of the Loader. Advance the device under heparinized saline.

16. As **Figure 15** unscrew the adjust nut at the proximal of the Haemostatic Valve. Gently push the Delivery Cable to advance the occluder till it reaches the distal of the Delivery Sheath.

- ① Unscrew the adjust nut at the proximal of the Haemostatic Valve
- ② Gently push the Delivery Cable to advance the occluder till it reaches the distal of the Delivery Sheath.

Figure 15 Advance the occluder

Note: Do not rotate the delivery cable when advancing the occluder. When the LAA occluder nearly reach the distal of the Delivery Sheath, gently pull back about 1cm to determine if the LAA occluder is still attached to the Delivery Cable under fluoroscopic guidance, then advance the occluder again to the tip of the Delivery Sheath.

17. Under DSA position the Delivery Sheath at around 5mm proximal to landing zone, gradually deploy the device by pushing the Delivery Cable forward while hold the Delivery Sheath.

Note: When pushing the Delivery Cable slowly, the Umbrella will open gradually. It is important to note the shouldering of the Umbrella, implying engagement of the stabilizing hooks into the appendage walls. The position of the Umbrella could be verified by contrast injections before and during the deployment of the Umbrella.

18. Withdraw the Delivery Sheath and hold the Delivery Cable at the same time until the Cover of LAA occluder is exposed. Keep a slight pushing tension on the Delivery Cable to make the Cover lean towards the direction of the Umbrella. This maneuver ensures the device would cover the LAA orifice.

Note: Perform echocardiography and angiography to confirm if the device is stable: as **Figure 16**, the Umbrella is compressed and is positioned the distal to the origin of left circumflex artery (by echocardiography); all the parts of the Umbrella are unfolded (by angiography in various projections); and the Cover resumes a concave conformation (by angiography). If not, recapture the device by holding the Delivery Sheath in place while pulling the Delivery Cable to bring the device completely back to the Delivery Sheath. Adjust the position and redeploy. Or exchange for another new LAA occluder with different **REF** may be required.

Warning: Do not advance the Delivery Cable during device recapture. The device can be recaptured

and redeployed for a maximum of 4 times. If the device positioning is still unsatisfactory, remove and replace both the device and the Delivery Sheath with a new system.

19. Perform a tug test to check for device stability.

20. Ensure the orifice of LAA is sealed by the Cover by injecting contrast and TEE.


Figure 16 Status of the Left Atrial Appendage Occluder implanted in LAA

Note: Recapture the device any time if necessary before releasing.

21. When the device placement is optimal it may be released:

- Detach the device from the cable by turning the delivery cable vise in a counter-clockwise manner.
- Remove the delivery cable and the delivery sheath from the patient.

22. Perform echocardiography and angiography to confirm that the device remains in place. Check for any adverse events.

23. Close the right femoral venous wound.

VIII. POST-PROCEDURE INSTRUCTIONS

- Double anti-platelet therapy (Aspirin plus Clopidogrel or the alternate) for initial 3 months after device implantation. Continuation of aspirin indefinitely after 3 months is recommended. Routine standard of care during clinical follow-ups are recommended, including trans-esophageal echocardiography preferably at 3-month follow up, for evaluation of incomplete LAA occlusion and adverse events (including thrombus formation on device surface).

IX. EXPIRY DATE

The Left Atrial Appendage (LAA) Closure System are sterilized by ethylene oxide. The shelf life is indicated on the labels, and the manufacturing date and use by date are marked on the label. **Do not use an expired product.**

X. PACKAGE AND LABEL

- The Left Atrial Appendage (LAA) Closure Systems are supplied sterile. The Occluders and the Delivery Systems are packaged separately.
- The Occluder is packaged and protected in a PETG tray sealed with dialyzing paper, and then the product is sealed with a dialyzing pouch, where a label and a sterile indicator are attached. The product is sterilized and put in a box with IFU, patient card, customer feedback form and quality guarantee card. An outer label is stuck on the box.
- The Delivery System is protected on a support plate, and sealed with double dialyzing pouches, on which a primary label and a sterile indicator are stuck. The product is sterilized and put in a box with IFU, customer feedback form and quality guarantee card. An outer label is stuck on the box.
- Information about product **REF**, sterilization date, expiry date, serial number, batch code, sterilization mark, etc. are indicated on the label.
- The symbols used are as follows.

XI. DISPOSAL

- The carton and IFU are recyclable. Dispose of all packaging materials as appropriate.
- Devices may be returned to Lifetech for disposal.
- Devices may be disposed of per standard solid biohazard waste procedures.

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Lifetech Scientific (Shenzhen) Co., Ltd