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INSTRUCTIONS FOR USE

Tracheal tube holder, blue, one-piece

Applies to holder REF 19-0002

I. FOREWORD

These instructions for use are valid for Tracheal tube holder. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling. Please read the instructions for use carefully before using the product for the first time!

II. INTENDED USE

Tracheal tube holder is a one-piece fixation strap for endotracheal tubes that enables secure fixation in position. Tracheal tube holder can be adapted to different neck sizes. The enclosed adhesive hook-and-loop fastener provides additional fixation.

III. WARNINGS

Tracheal tube holder are single-patient products and are intended for single use only. They must not be cleaned or disinfected.

To prevent the endotracheal tube from accidentally becoming detached, the secure hold of Tracheal tube holder must be checked regularly.

Too tight fixation can cause pressure ulcers at the corners of the mouth or in the mucous membranes of the mouth or nose.

When putting on the Tracheal tube holder, keep the endotracheal tube fixed in place with one hand to avoid the risk of accidental extubation or dislocation, due for instance to an urge to cough induced by irritation of the throat, or to tensile forces.

IV. INSTRUCTIONS

- Check whether the packaging and the product are in a faultless condition.
- Fix the endotracheal tube in place with one hand, and if necessary, remove the Tracheal tube holder that is to be replaced.

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and if necessary, remove the
for the Quality System

Zbigniew Milewski

- If this has not yet been done, affix the adhesive hook-and-loop fastener to the desired position on the endotracheal tube where the Tracheal tube holder is to be fixed in place ❶.
- Place the Tracheal tube holder with the blue side onto the adhesive hook-and-loop fastener on the endotracheal tube and cross it in such a way that the Tracheal tube holder is divided up into a long section and a short section ❷.
- Place the longer end around the patient's neck with the white side of the material facing the patient's skin ❸.
- Connect the ends by means of the enclosed adhesive strip and check the endotracheal tube with the Tracheal tube holder for secure hold ❹.
- If the Tracheal tube holder is too long, it can be shortened to the desired length by cutting it with a pair of scissors ❺.
- A mouth wedge/oropharyngeal tube can additionally be fixed in place by wrapping the Tracheal tube holder around it ❻.

V. SERVICE LIFE

Tracheal tube holder is single-patient products and is intended for single use only. The service life depends on individual handling. The products must be replaced immediately when damaged.

VI. LEGAL NOTICES

The manufacturer SUMI will not accept any liability for damages (especially not for functional deficiencies, injuries, infections, and/or other complications or adverse events) caused by unauthorised product alterations, by repairs not performed by the manufacturer or by improper handling, care (cleaning/disinfection) and/or storage of the products in violation of the provisions of these instructions for use.


This applies – to the extent permitted by law – both to damages to the products themselves caused thereby and to any consequential damages caused thereby.

Product specifications described herein are subject to change without notice.

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Director Representative
for the Quality System

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TRACHEAL TUBE INTRODUCER TRACHEAL TUBE GUIDE INSTRUCTIONS FOR USE

These instructions for use refer to the following products:

09-3617, 09-3627, 09-3610, 09-3620, 09-3717, 09-3727, 09-3817, 09-3827, 09-3810, 09-3820,
09-3117, 09-3127, 09-3110, 09-3120, 09-5617, 09-5627, 09-5610, 09-5620, 09-5717, 09-5727,
09-5817, 09-5827, 09-5810, 09-5820, 09-5117, 09-5127, 09-5110, 09-5120, 09-1610, 09-2010.

Description

1,6 mm for 2,0 - 3,0 intubation tubes 2,0 mm for 2,5 - 4,0 intubation tubes
3,3 mm for 4,0 - 5,5 intubation tubes 5,0 mm for 6,0 - 10,0 intubation tubes
Lengths: 600mm, 700mm, 800 mm and 1000mm
They can be used again after washing and sterilization, maximum 4 times
(it concerns the multi-use guides)

Indications

A difficult intubation

Replacement of the intubation tubes

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Instructions for use

Upon difficult intubations:

- To gently remove the sterile tracheal tube from the packaging and check before use.
- To wash the guide with a sterile solution of water-soluble lubricant right before use
- to connect the patient to oxygen and perform endoscopy upon observance of current medical techniques
- If there is such a necessity, to tentatively shape the guide in order to get an optimal flexure.
- to keep particular caution.
- To insert the guide into the throat and next manipulate so as to insert it into the larynx.
- It is advisable to have confirmed correct arrangement of the guide by correct laryngospasms.
- To leave an endoscope in the patient's throat, and next put a tube on the guide and carefully insert the tube through the larynx into the trachea.
- To pay attention to observance of current approved medical techniques during the intubation.
- The marker of depth of intubation may be used as assistance upon correction of arrangement of the tube.
- Keep the tube and next carefully remove the guide.

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- To gently remove the sterile guide from the packaging and check before use.
- To observe approved medical techniques with regard to correct size of the devices, intubation and extubation.
- To select a point of reference, for instance teeth or nostrils and mark the depth of insertion of intubation tube.
- To insert the guide and to slip it until a moment, when a depth marker of the guide comes level with a depth marker of the tube.
- Keep the guide and remove the intubation tube.
- Take advantage of a relevant marker in order to check arrangement of the tube's guide.
- Place a replaceable intubation tube on the tube's guide and insert deeply until a depth marker of the tube comes level with a depth marker of the guide. Both of them are to be correctly inserted taking into consideration relevant limiters.
- While airing, check once more, whether the new intubation tube is arranged correctly in the trachea.
- Remove the guide from the intubation tube.

Washing, disinfection and sterilization of the multi-use guides

- Insert the used guide into the disinfectant.
- Wash the guide carefully using water solution of soap, removing precisely all visible impurities.
- In order to remove impurities, wash the product in distilled water with a disinfecting agent.
- Dry the guide carefully.
- Check whether the product is not damaged, throw aside a damaged product.
- Pack up the product in new packaging of paper - foil type.
- Close the packaging lightly.
- Carry out sterilization using ethylene oxide with application of relevant sterilization procedures.
- During intervals between procedures the product should be kept in a sterile condition, in packaging of paper - foil type, in dry rooms, upon 5-30° C temperature, in a manner protecting against direct action of solar radiation, vapours of acids and organic solvents as well as other harmful agents.
- The SUMI company does not bear responsibility for the sterilization process.

Precautions

A physician should take proper precautions in order to avoid potential damage of an intubation tube with a guide during intubation or during replacement of a tube.

Warnings

Owing to a danger of respiratory tract infection, mainly by treatment resistant hospital bacterial strains, one should observe current approved medical procedures and techniques. Repeated use of a guide designed for single use may constitute a danger to the patient. Observance of the current approved medical techniques and performing of procedures by skilled medical personnel reduces to a minimum possibility of occurrence of complications e.g.: excoations of epithelium and ulcerations of lips, mouth, gullet, vocal folds, injuries of lips, tongue, gullet, larynx, vocal apparatus of larynx, palate, tonsils and teeth, paresis of hypoglossal or glossal nerves, perforation of gullet, perforation of larynx, bleedings from throat, larynx, stricture of larynx, trachea.

Notes

- Do not use in case of damage of packaging.
- Do not use after expiry date.

Responsibility

The guides of SUMI make undergo physical, biological and chemical investigations. It allows their safe use and top quality. There is recommended making familiar with the instructions for use before using the guides for difficult intubations. The tracheal tubes manufactured by SUMI may be used only by qualified medical personnel. SUMI shall not bear responsibility for consequences caused by inappropriate selection of the size or using the tracheal tubes contrary to their use. In the event of a claim one should keep the guides together with direct packaging.

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Edition no 06 Date of publication: 02.11.2018

AKTUALNY
data 23.08.20

Director Representative
for the Quality System
Zbigniew Milewski

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TRACHEAL TUBE GUIDE WITH VENTILATION LUMEN INSTRUCTIONS FOR USE

Description.

The Assortment of sterile, single-use tracheal tube guides with ventilation lumen be used through mouth, is characterized by the following properties:

1. They are made of medical PVC and PE of special properties providing both stiffness and flexibility upon body temperature of a patient.
2. Tracheal tube guides are provided with 15 mm connector which enables the connection of a set succour breathing and also are equipped with a Luer lock
3. There are placed some markers of depth on the distal end to make it easier replacement and use of the device.
4. They occur in four sizes:
 - 2,6 mm for 3,0 – 6,0 intubation tubes
 - 3,3 mm for 4,0 – 6,0 intubation tubes
 - 3,7 mm for 4,0 – 6,0 intubation tubes
 - 4,0 mm for 5,0 – 11,0 intubation tubes
 - 4,7 mm for 5,0 – 11,0 intubation tubes
 - 5,0 mm for 6,0 – 11,0 intubation tubes
 - 6,0 mm for 7,0 – 11,0 intubation tubes
 - 6,3 mm for 7,0 – 11,0 intubation tubes
 and are of lengths: 600mm and 830 mm
5. They are sterile ones until opening or damaging of their package.
6. They are delivered in a sterile condition.

Set composition:

1. the guide
2. 15 mm connector
3. Luer lock connector
4. instruction for use

Indications

A difficult intubation. Replacement of the intubation tubes

Contraindications

The guide is not recommended for intubation through nose or upon changes of trachea.

Instructions for use.

Upon difficult intubations:

1. To gently remove the sterile tracheal tube from the packaging
2. To wash the guide with a sterile solution of water-soluble lubricant right before use
3. to connect the patient to oxygen and perform endoscopy upon observance of current medical techniques
4. to keep particular caution.
5. To insert the guide into the throat and next manipulate so as to insert it into the larynx.
6. It is advisable to have confirmed correct arrangement of the guide by correct laryngospasms.
7. To leave an endoscope in the patient's throat, and next put a tube on the guide and carefully insert the tube through the larynx into the trachea.
8. To pay attention to observance of current approved medical techniques during the intubation.
9. The marker of depth of intubation may be used as assistance upon correction of arrangement of the tube.
10. Keep the tube and next carefully remove the guide.
11. Remove the endoscope and next connect the tube to the system of respiration.

Upon replacement of an intubation tube:

1. To gently remove the sterile tracheal tube from the packaging take out the wire guide of the introducer and check before use. Push connector in till resistance.
2. To observe approved medical techniques with regard to correct size of the devices, intubation and extubation.
3. To prepare the oxygen set compatible with a connector which is in the connecting set
4. To remove the fixation if the tube is fixed to a patient
5. To select a point of reference, for instance teeth or nostrils and mark the depth of insertion of intubation tube.
6. To insert the guide and to slip it until a moment, when a depth marker of the guide comes level with a depth marker of the tube.
7. To remove the air from the inside of the tube
8. Keep the guide and remove the intubation tube.
9. Take advantage of a relevant marker in order to check arrangement of the tube's guide.
10. In case of hypoxia pass the oxygen through a 15 mm connector placed on the end of the guide or apply the Yet type guides using a Luer lock connector
11. Place a replaceable intubation tube on the tube's guide and insert deeply until a depth marker of the tube comes level with a depth marker of the guide. Both of them are to be correctly inserted taking into consideration relevant limiters.
12. While airing, check once more, whether the new intubation tube is arranged correctly in the trachea
13. Remove the guide from the intubation tube
14. Throw aside the guide. Observe the procedure of particular warnings given in the instructions for use

Warning: The connector has to be disconnected from a guide in case when the tube is placed on a guide

Precautions

1. A physician should take proper precautions in order to avoid potential damage of an intubation tube with a guide during intubation or during replacement of a tube.

Warnings.

1. Owing to a danger of respiratory tract infection, mainly by treatment resistant hospital bacterial strains, one should observe current approved medical procedures and techniques.
2. Repeated use of a guide designed for single use may constitute a danger to the patient.
3. A substance that allows smooth applying can not be in the inner lumen of the guide.
4. The use of the above mentioned substance on the basis of silicone may result in negatively on the connection of the connector with the guide
5. Before the connection of the connector carefully clean the tip and remove the possible remains of the slide-allowing substance

Notes.

1. Do not use in case of damage of packaging.
2. Do not use after expiry date.

Possible complications

Observance of the current approved medical techniques and performing of procedures by skilled medical personnel reduces to a minimum possibility of occurrence of complications. there are possible the following complications (no all of them have been mentioned): excoriations of epithelium and ulcerations of lips, mouth, gullet, vocal folds, injuries of lips, tongue, gullet, larynx, vocal apparatus of larynx, palate, tonsils and teeth, paresis of hypoglossal or glossal nerves, perforation of gullet, perforation of larynx, bleedings from throat, larynx, strictures of larynx, trachea.

Disposal.

The guides are a single product and must be disposed of after use in accordance with the laws in force in the country.

Responsibility.

The guides of SUMI make undergo physical, biological and chemical investigations.

It allows their safe use and top quality.

There is recommended making familiar with the instructions for use before using the guides for difficult intubations.

The tracheal tubes manufactured by SUMI may be used only by qualified medical personnel.

SUMI shall not bear responsibility for consequences caused by inappropriate selection of the size or using the tracheal tubes contrary to their use.

In the event of a claim one should keep the guides together with direct packaging.

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TRACHEAL INTUBATION STYLET INSTRUCTIONS FOR USE

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09-1923	1,9 mm for 2,0 intubation tubes
09-2023	2,0 mm for 2,5 intubation tubes
09-2223	2,2 mm for 3,0 intubation tubes
09-3034	3,0 mm for 3,5 - 5,0 intubation tubes
09-4034	4,0 mm for 5,5 - 7,0 intubation tubes
09-4060	4,0 mm for 5,5 - 7,0 intubation tubes
09-5037	5,0 mm for 6,0 -11,0 intubation tubes
09-5060	5,0 mm for 6,0 -11,0 intubation tubes
09-6037	6,0 mm for 7,0 -11,0 intubation tubes
09-6060	6,0 mm for 7,0 -11,0 intubation tubes

Indications.

Stylets are used to shape tracheal tube for easier intubation

Instructions for use.

1. To gently remove the sterile tracheal intubation stylet from the packaging.
2. Introduce stylet into tracheal tube and shape it according to anatomy of a patient.
3. Intubate patient through mouth or nose paying attention to observance of current approved medical techniques.
4. Hold the tracheal tube in place while taking out the stylet

Precautions

A physician should take proper precautions in order to avoid potential damage to tracheal tube during intubation or during replacement of a tube.

Warnings.

Owing to a danger of respiratory tract infection, mainly by treatment resistant hospital bacterial strains, one should observe current approved medical procedures and techniques. Repeated use of a stylet designed for single use may constitute a danger to the patient. Observance of the current approved medical techniques and performing of procedures by skilled medical personnel reduces to a minimum possibility of occurrence of complications e.g.: excoriations of epithelium and ulcerations of lips, mouth, gullet, vocal folds, injuries of lips, tongue, gullet, larynx, vocal apparatus of larynx, palate, tonsils and teeth, paresis of hypoglossal or glossal nerves, perforation of gullet, perforation of larynx, bleedings from throat, larynx, strictures of larynx, trachea.

Notes.

1. Do not use in case of damage of packaging.
2. Do not use after the expiry date.
3. Do not use in case of damaged stylet

Disposal.

The stylets are a single use product and must be disposed of after use in accordance with the laws in force in the country.

Responsibility.

The stylets made undergo physical, biological and chemical inspections. It allows their safe use and top quality.

It is recommended to read carefully the instructions for use before using the stylets.

The stylets may be used by qualified medical personnel, only.

Producer shall not bear responsibility for consequences caused by inappropriate selection of the size or using the tracheal intubation stylets contrary to their designation.

In event of a claim one should keep the stylets together with direct packaging.

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Cleaning instructions (only for multiple use masks)

Only neutral (pH7) detergents / cleaning solutions should be used for disinfection/sterilization. The mask cannot be sterilized/disinfected by exposure to reaction of formaldehyde, glutaraldehyde or ethylene oxide. Avoid prolonged immersion in chlorhexidine, iodine-based antiseptics and spraying with silicon-based lubricants as it can damage the mask. Cleaning or wiping with the solvent solution, exposure to auxiliary rinsing agents or use of non-neutral detergents (pH7) has a negative effect on the silicone components and on other mask components.

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Necessary equipment:

- Disinfection apparatus for anaesthesia/respiratory equipment. This is usually a cabinet type equipment. Selected disinfection apparatus has to be equipped with a cycle counter.
- A disinfection apparatus should be in accordance with BS2745 Part 1 and 3 and HTM2030.
- Neutral detergent (pH7) preferably with a dosage from a calibrated, automatic closed system.
- Mechanical dryer as a component stage in automatic process or as a separate fan with temperature control in a drying cabinet.

Procedure:

1. Provide fully operative automatic disinfection machine. Protective clothes is recommended (apron, gloves and eye protection, if spraying is possible).
2. Load the mask to the automatic machine so the loading configuration does not interfere with cleaning process.
3. Secure the lock, select a cycle and push the start button.
4. During the cycle, make sure that all stages of the process are complete, remove the content and check the cleanliness of the mask. Pay special attention to the internal surfaces of the ventilation tube and the cuff. There should not be any visible loose particles nor secretion.
5. If the mask is not dry, the water residue should be removed and the mask should be put into a dryer. The mask has to be completely dry before the sterilization.

Scope of actions:

- Cleaning is obtained by continuous spraying, flooding or irrigation of the mask with water and detergents during several stages of the cycle. A typical cycle includes the following stages:
- Cold flushing below 35°C.
 - Hot bath in neutral detergent (pH7) at a temperature app. 55°C.
 - Hot disinfection when the mask surface temperature should amount to minimum 71°C for 3 minutes or up to 90°C, for 1 minute.

Manual cleaning

Necessary equipment

Sink or a vessel that could contain a sufficient volume of water/detergent (pH7) so the mask can be fully immersed. Hot water with a neutral detergent solution (pH7) in proper concentration and temperature. A vessel with water for rinsing. Surface drying - mechanical equipment for drying (a cabinet for drying or industrial hot air dryer). Brushes, including a brush with a soft bristle for external surface and a tube brush for cleaning the internal side of a tube and a spray.

Procedure:

1. Make sure that the cleaning vessel is clean and dry.
2. Protective clothes is recommended - (gloves, apron and eye protection in case of spraying). Fill the vessel with a proper amount of hot water/neutral detergent solution so the mask can be fully immersed.
3. Immerse the mask thoroughly in the solution to remove the remainder of the air. It is important to allow the cleaning solution to access the whole area including the internal surface of the tubes.
4. Brush, wipe, rinse or bath the mask by spraying to remove visible dirt, making sure that this action takes place under the surface of the solution in order to avoid splashing and evaporation. Brush the interior of the ventilation tube, placing the tube brush from the cuff end side carefully to avoid cuff or tube damage.
5. Remove the mask from the cleaning solution and dry it before putting it into a clean rinse vessel or sink.
6. Rinse the mask thoroughly with clean water or a spray, make sure the whole mask is immersed.
7. Take the mask out of the rinse water and dry it. Carefully check the cleanliness of the mask, paying special attention to internal surfaces of the ventilation tube and the cuff. There should be no visible loose particles nor secretion.
8. Put it into a drying cabinet. The mask has to be completely dry before sterilization.
9. Fill in the record card.

Scope of action

This process physically removes dirt and other contaminations, but it does not include the disinfection process. Therefore, there is a risk of infection of the servicing staff during the actions before sterilization.

Sterilization after cleaning

A steam sterilizer is recommended. Remove the air from the cuff completely before sterilization. The cuff has to be emptied completely otherwise the cuff will widen in the autoclave resulting in the cuff and the pilot balloon damage. The mask should be sterilized with steam at the temperature 134°C (+4°C / -0°C). Before putting a product in the machine, make sure the mask is completely dry.

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REF	catalogue number		manufacturer
LOT	batch no.	STERILE EO	sterilized using ethylene oxide
	date of manufacture	LATEX	latex-free
	use-by date		consult instructions for use
CE 1011	the product complies with the requirements of the European Union/Number of the Notified Body participating in the device conformity assessment		caution consult precautions and warnings in instructions for use
			keep dry
	keep away from sunlight		keep away from sunlight
	store at 5-40 °C		store at 5-40 °C

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1. Description

1.1. Kit components.

Laryngeal masks include the following components:

- 15 mm connector
- high diameter ventilation tube
- cuff
- inflating valve – for inflating and deflating a cuff
- pilot balloon – for indication of the cuff inflation level
- inflation line


Figure no. 1

1.2. Necessary equipment.

The following equipment is required to insert the masks:

- anti-bite blockade
- water-based lubricant
- Luer type syringe, size should correspond to the required volume
- a guide facilitating the insertion of reinforced masks

2. Indications.

A laryngeal mask is used for patients in order to enable breathing in the anaesthesia conditions and artificial ventilation in case of a lack of their natural breath.

Compared to the intubation tube, the laryngeal mask is not inserted into a sensitive trachea, reduces the irritation and can be easily and comfortably inserted.

3. Usage instructions.

1. Choose a proper mask size adequate for a patient weight using the following table:

Description	Weight (kg)	Recommended size	Cuff volume (ml)
Newborn infant	up to 5 kg	1	4
Infant	5-10	1 1/2	7
Small child	10-20	2	10
Bigger child	20-30	2 1/2	14
Big child/ small adult	30-50	3	20
Adult	50-70	4	30
Big adult	70>	5	40

3.2. Tests before usage.

Mask tests should be performed before usage, incorrect result of any test indicates that the product is improper for use. It should be destroyed and replaced with another one. A single use laryngeal mask comes in a ready-to-use condition. It only requires taking into account the following 3rd part of the manual.

3.2.1. Visual inspection.

Perform visual inspection of the mask for mask transparency, discoloration and damage.

Do not use a discoloured product of low transparency or a damaged product.

Check the internal tube in order to verify if it is free from any contamination.

Check whether the connector fits the airways and it cannot be disconnected accidentally.

Check the cuff inflation line in order to make sure that there is no leakage or damage.

Check the pilot balloon for any defect.

3.2.2. Inflation and deflation tests.

Step 1: Leakage

Use a Luer type syringe to inflate the cuff so the cuff walls are properly flat.

Check the mask for any wrinkles on the cuff after deflation and preparation for insertion. After cuff deflation to vacuum, inflate the cuff with proper air volume using the following table:

Cuff inflation volume test	
Mask size	Volume* (ml)
1	6
1 1/2	10
2	15
2 1/2	21
3	30
4	45
5	60

* inflate the cuff with these volumes only for the test purposes.

NOTICE: Any cuff deflation tendency should be visible within two minutes and should indicate a leakage.

Step 2: Rupture

When the cuff is inflated, check its shape and symmetry, inflation line and the pilot balloon. Do not use the mask if there is a cuff unevenness or bulge. Do not use if the inflation line has an irregular shape. The pilot balloon should easily lose the pressure and its elliptical shape.

WARNING: Do not use the mask if any of the tests performed before the use indicated any irregularity. Throw away the product and replace it with a new mask. Always have a replacement mask of a proper size for immediate use.

3.3. Insertion technique

Apply general precautions: gloves, eye protection etc.

Deflate the cuff completely using the syringe, pressing the mask bowl to a flat surface.

Apply lubricant only on the back of the cuff. Do not apply excessive lubricant!

Place a patient on his/her back, with a neck pulled away, unless there is a suspicion that a spine can be damaged or in case of a contraindication for changing the patient's position.

In order to insert a mask, hold it as a pen, with an index finger placed outside the edge of the cuff and the tube connection. A black line should be directed towards the patient's nose. Insert the mask to the mouth with the cuff opening directed towards the patient's feet (see fig. 2 and 3)

3.6.1. Mechanical cleaning

Mechanical cleaning using disinfection apparatus for hot washing

Necessary equipment:

- Disinfection apparatus for anaesthesia/ respiratory equipment. This is usually a cabinet type equipment. Selected disinfection apparatus has to be equipped with a cycle counter.
- A disinfection apparatus should be in accordance with BS2745 Part 1 and 3 and HTM2030.
- Neutral detergent (pH7) preferably with a detergent from a calibrated, automatic closed system.
- Mechanical dryer as a component stage in automatic process or as a separate fan stage with temperature control in a drying cabinet.

Procedura:

1. Provide fully operative automatic disinfection machine. Protective clothes is recommended (apron, gloves and eye protection, if spraying is possible).
2. Load the mask to the automatic machine so the loading configuration does not interfere with cleaning process.
3. Secure the lock, select a cycle and push

3.6.2. Manual cleaning

During the cycle, make sure that all stages of the process are complete, remove the mask and check the cleanliness of the mask. Pay special attention to the internal surfaces of the ventilation tube and the cuff. There should not be any visible loose debris or secretion.

If the mask is not dry, the water residue

Figure 2


Figure 3

Press the end of the mask to a tough palate in the mouth and push the mask backward and downwards in one deft move until you feel a resistance. The end of the cuff should be located in the fauces (see Figure 4).

Figure 4

Apply proper size of the syringe and inflate the cuff with proper air volume as shown on the tube.

Successful insertion of the mask into the mouth.

should be removed and the mask should be put into a dryer. The mask has to be completely dry before the sterilization.

Scope of actions:

- Cleaning is obtained by continuous spraying, flooding or irrigation of the mask with water and detergents during several stages of the cycle. A typical cycle includes the following stages:
1. Gold flushing below 35°C.
 2. Hot bath in neutral detergent (pH7) at a temperature approx. 55°C.
 3. Hot disinfection when the mask surface temperature should amount to minimum 71°C for 3 minutes or up to 90°C for 1 minute.

3.6.2. Manual cleaning Necessary equipment

- Sink or a vessel that could contain a sufficient volume of water/detergent (pH7) so the mask can be fully immersed.
- Hot water with a neutral detergent solution (pH7) in proper concentration and temperature.
- A vessel with water for rinsing.
- Surface drying – mechanical equipment for drying (a cabinet for drying or industrial hot air dryer).
- Brushes, including a brush with a soft bristle for external surface and a tube brush for cleaning the internal side of a tube and a spray.

Procedure:

1. Make sure that the cleaning vessel is clean and dry.
2. Protective clothes is recommended – (gloves, apron and eye protection in case of spraying). Fill the vessel with a proper amount of hot water/neutral detergent solution so the mask can be fully immersed.
3. Immerse the mask thoroughly in the solution to remove the remainder of the air. It is important to allow the cleaning solution to access the whole area including the internal surface of the tubes.
4. Brush, wipe, rinse or bath the mask by spraying to remove visible dirt, making sure that this action takes place under the surface of the solution in order to avoid splashing and evaporation. Brush the interior of the ventilation tube, placing the tube brush from the cuff end side carefully to avoid cuff or tube damage.
5. Remove the mask from the cleaning solution and dry it before putting it into a clean rinse vessel or sink.
6. Rinse the mask thoroughly with clean water or a spray, make sure the whole mask is immersed.
7. Take the mask out of the rinse water and dry it. Carefully check the cleanliness of the mask, paying special attention to internal surfaces of the ventilation tube and the cuff. There should be no visible loose particles nor secretion.
8. Put it into a drying cabinet. The mask has to be completely dry before sterilization.
9. Fill in the records card.

Scope of action

This process physically removes dirt and other contaminations, but it does not include the disinfection process. Therefore, there is a risk of infection of the servicing staff during the actions before the sterilization.

4. Sterilization after cleaning

A steam sterilizer is recommended. Remove the air from the cuff completely before sterilization. The cuff has to be emptied completely otherwise the cuff will widen in a pilot balloon damage.

The mask should be sterilized with steam at the temperature 134°C (+4°C/-0°C). Before putting a product in the machine, make sure the mask is completely dry.

5. Additional information.

- A multiple use mask has a colour cuff and a pilot balloon.
- A single use mask has a transparent cuff and a pilot balloon.
- A multiple use mask has a guarantee for 40-time use during a period of 2 years starting from the first sterilization. Excessive use over the limit can negatively affect the use of the mask.
- Each use has to be recorded in the records card.
- SUMI cannot undertake the liability for any product complained as defective if it is returned without the records card properly filled in.

6. Precautions.

1. Read the laryngeal mask user's manual thoroughly. It includes proper warnings describing potential risk during the use and proper precautions to be undertaken in case of difficulties with the product.
2. Do not use the multiple use larynx mask if it was exposed, cleaned or sterilized with disinfecting sterilizing chemicals as glutaraldehyde, ethylene oxide, phenol-based cleaners or cleaners containing iodine. These substances are absorbed by the material of the product and therefore can cause unnecessary risk for a patient or possible deterioration of the product quality.
3. Perform the tests described herein before use. Any test failure indicates that the mask cannot be used and should be destroyed and replaced with another one.
4. In order to avoid injury, force cannot be used at any stage of insertion. Never top over the cuff when the mask is inserted. If proper tightness is not achieved at a recommended cuff inflation volume, then consider other size.
5. A multiple use mask comes with a 40-time use guarantee. Do not use the product after its 40th use. Users who will ignore it will use the product not in accordance with the manufacturer's instructions and they or their employees shall be responsible for such use. A user shall bear the whole liability for any loss, damage or any injuries resulting from the use of the product when the 40-time limit is exceeded.
6. The product manufacturer, distributor or supplier shall not bear any liability for any loss, damage or injuries resulting from the use of the product when the 40-time limit is exceeded.
7. The laryngeal mask cannot be put into the autoclave if the cuff contains the air or moisture as the high temperature and resulting pressure will cause the cuff widening and damage that cannot be repaired. Remove the air from the cuff so the walls are exactly flat before putting the mask into the autoclave.
8. Under no circumstances use the set if the tube or cuff is damaged (e.g. surface crack, leakage).
9. Do not use during the magnetic resonance MRI.
10. Warnings.
1. Due to a risk of airway infection with treatment-resistant hospital strains, the following approved medical procedures and techniques should be obeyed.
2. Resterilization or reuse can create a risk for a patient.
11. Remarks.
1. Store at room temperature.
2. Do not use if the package is damaged.
3. Do not use after expiry date.
12. Possible complications.
1. Throat irritation.
2. Dryness in throat and/or mucous membrane.
3. Side effects due to improper insertion can differ depending on the way of use.
13. Responsibility.
1. The masks manufactured by SUMI undergo physical, biological and chemical tests. Therefore they can be used safely and their highest quality is assured.
2. Read the user's manual thoroughly before using the masks.
3. Products manufactured by SUMI can be used only by a qualified medical staff.

SUMI do not bear the responsibility for consequences resulting from improper size selection, using the product not according to their intended use. In case of complaint, keep the masks with their direct package and the records card (in case of multiple masks).

SUMI cannot undertake the liability for any product complained as defective if it is returned without the records card properly filled in.

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ACTUALNY

2018

Director Representative for the Quality System

Krzysztof Milewski

Description: the HiCapno® is a patented, single-use, high concentration oxygen mask which has specific construction of Exhale/Inhale valve. Two chambers separated by a low resistance membrane together with oxygen bag supplied with O₂, facilitate ventilation of patient with oxygen concentration of up to 98%. Design of exhaust chamber allows for evacuation of the exhaled air to prevent reinhalation of carbon dioxide. At the same time, it provides a possibility of sampling of CO₂ to capnograph on "side stream" via Luer-Lock connector. Mask connected to capnograph provides continuous monitoring of breathed out CO₂ concentration. HiCapno® set consist of oxygen tube, CO₂ sampling line used for delivering sample of exhaled air to capnograph with antibacterial filter.

Indications: Oxygen Therapy

Instructions for use:

1. Using oxygen tubing connect oxygen mask and oxygen source.
2. Set the prescribed oxygen flow.
3. Make sure that the oxygen bag is filled with oxygen (min ½ of bag capacity).
4. Place the mask on the patient's face.
5. Adjust strap tension.
6. Mould nose clip to the patient's face
7. Connect the luer-lock tip in mask with capnograph using sampling line with luer-lock

Precautions: because of the risk of respiratory tract infections, with mostly resistant to treatment of hospital strains, comply with the current accepted medical techniques and procedures. To resterilization or reuse may pose a risk to the patient.

Warning: Product is intended for single use only. Condition for the proper use of HiCapno® is the original fill oxygen bag. Do not block the outlet of the bag.

Comments: Absolutely ensure continuous supply of O₂ by oxygen drain, flow min. 9 to 10 l/min to prevent re-inhalation of CO₂. Proper placement of the mask on a patient's face along with securement on patient head with elastic band, allows maximum tightness. Mask should be stored in a dark place at a temperature of 5-40° C. Do not use in case of damaged packaging. Do not use after the expiry date.

Responsibility: SUMI's products are subjected to physical, biological and chemical research. This allows for their safe use and guarantees the highest quality. We recommend you carefully read the instructions before using the masks. Products manufactured by the SUMI may be used only by qualified medical staff. The manufacturer is not liable for consequences caused by misuse of masks, in accordance with the intended purpose. In case of complaint, please keep the products together with the immediate container

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SUCTION HANDLE WITH FILTER – INSTRUCTION FOR USE

DESCRIPTION.

Suction handle with filter is a device used for effective suctioning of bone chips and debris from the operation field, without the danger of tubing occlusion. Construction of the container and filter, provide easy, undisturbed, flow in even the hardest conditions.

PRODUCT FEATURES.

- Ergonomic, transparent handle provides a firm grip
- Exchangeable filter gathers bone chips, cement and other larger debris
- Rounded distal tip with two side holes ensures effective suctioning
- Additional filter for exchange
- In set with 3 removable tips with lengths: 140mm, 220mm, 260mm, for various cases
- Sterile, in double packing

INDICATIONS.

Suction handle with filter is designated for intra-operative suction of body fluids and debris during the orthopaedic treatments. At the same time, while suctioning, the handle filters the fluids from larger debris, after connecting to the outside source of vacuum.

KIT COMPONENTS.

canister for holding with 2 exchangeable filters and 3 tips with different lengths.

INSTRUCTION FOR USE.

Check if the set is complete

Use in accordance with its designation

To exchange a filter, screw the top back part of the handle

PRECAUTIONS.

Do not use if either package or any element of set, is damaged.

WARNINGS.

1. Due to the risk of infection of airway by some hospital stains, current approved medical procedures and techniques should be observed. Repeated sterilization or use may constitute a risk to the patient.
2. The product is for single use, only.

NOTES.

1. The product should be stored in room conditions
2. Do not use if the package is damaged
3. Do not use after the expiry date
4. Dispose of the used product in accordance with the local standards for utilization

RESPONSIBILITY.

The suction handles manufactured by SUMI Company undergo physical, biological and chemical tests. This allows them to be used safely and guarantees top quality. It is recommended to thoroughly understand the instructions for use prior to use of the suction handles. The suction handles manufactured by SUMI may be used only by qualified medical personnel. SUMI shall not bear responsibility for consequences caused by inappropriate selection of the size or using the suction handles contrary to their use. In the event of a claim one should keep the suction handles together with direct packaging.

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