PL

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GB

NASOPHARYNGEAL TUBES – USER'S MANUAL

These instructions for use refer to the following products:

41-2500, 41-3000, 41-3500, 41-4000, 41-4500, 41-5000, 41-5500, 41-6000, 41-6500, 41-7000, 41-7500, 41-8000, 41-8500, 41-9000

Description.

Range of sterile disposable nasopharyngeal tubes for nasal application made of polyvinyl chloride (PVC) characterized by the following features:

1. Anatomical design with properly shaped curvature.
2. Remain sterile until opening or damaging the package.

Indications.

1. Required mechanical ventilation with positive pressure.
2. Protection against aspiration.
3. Suction of secretion from trachea and bronchi.
4. Application during head, neck and mouth surgical procedures where all auxiliary equipment is to be removed from the operation area – shaped or reinforced tubes.

Symbols:

Instructions for use.

1. Check nostril size and patency.
2. Select proper tube size. Its size should correspond to a distance from the lower jaw angle to a tip of the nose.
3. Carefully remove the sterile tube from the package.
4. Make sure that the tube diameter is not bigger than the nostril diameter.
5. Moisten the tube with water or gel to facilitate insertion. Do not use any substances based on oil e.g. Vaseline.
6. Insert the tube to a right-hand hole (as the tube is designed for a right-hand nostril), and insert towards the back side of the head through the bevel to the nose root or to the septum.
7. In case of resistance withdraw the tube and then gently turn it around and try again.
8. Do not apply excessive force!
9. When the tube restores nostril patency the resistance will decrease. Then insert the tube until its flange stops on the end of the nose.
10. Extubate a patient obeying existing medical techniques
11. Dispose of the tube. Follow the detailed warning procedures given in the instructions for use.

Precautions.

1. Avoid contacts with electrosurgery electrodes or laser beam as PVC will give off toxic smoke or will cause ignition of oxygen enriched environment.
2. If the tube is lubricated before use, it is essential to check if the lubricant does not get into the tube and blocks the tube gap as it can result in the ventilation blockade.

Warnings.

1. Due to a danger of airway infection with hospital bacterial strains resistant to medical treatment, existing and accepted medical procedures and techniques should be obeyed.
2. Resterilization or reuse can be dangerous to a patient.
3. This is a disposable tube for a single use only.

Remarks.

1. Keep the tube packages in room conditions.
2. Do not use if the package is damaged.
3. Do not use after the expiry date.

Possible complications.

Obeying existing accepted medical techniques and performing a treatment by a qualified medical staff minimize danger of complications.

The following complications are

possible (not all of them are listed): epithelium abrasion, infection of sinuses, epiglottitis, trachea, bronchi, lungs, nose injury, nose, throat and larynx bleedings.

Responsibility.

The nasopharyngeal tubes manufactured by SUMI undergo physical, biological and chemical tests. It allows for safe use and guarantees the highest quality.

Read thoroughly the instructions for use before using the nasal-pharyngeal tubes.

The nasopharyngeal tubes manufactured by SUMI can be used by a qualified medical staff only.

SUMI is not responsible for the effects caused by improper selection of size or the use of the nasal-pharyngeal tubes not according to its intended use.

In case of a claim, keep the tubes together with their direct packaging.

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Director Representative
for the Quality System

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