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IFU9509/28122020

1. DEVICE DESCRIPTION:

The Promesa™ BMS – Self-Expanding Nitinol Peripheral Stent System is a GHTF Class D Implantation Medical Device. This device consists of following components –


a) A self-expandable stent made from Nitinol alloy.

b) An over-the-wire (OTW) stent delivery system consisting of PTA catheter that has flexible retractable sheath.

The stent is pre-loaded on the delivery catheter constrained by covering it with the sheath. The stent is deployed by retracting the sheath which allows the stent to expand against the lesion in the arterial lumen to restore its patency.

Promesa™ BMS is permanent implantation device and its schematic diagram is shown in below Fig. 1.

Figure 1- Promesa™ BMS System Schematic Diagram

**1.1. Device Components Description:****1.1.1 Available Stent Lengths & Diameters:**

Available stent lengths & diameters (76 configurations) are shown in table-1 below for Promesa™ BMS;

Table 1 - Promesa™ BMS Size matrix

(OTW 0.035" 80cm)						
Nominal Stent Length (mm)	Unconstrained stent Diameter (mm)					
	7.00	8.00	9.00	10.00		
	Reference Number					
30	PRB0700030A	PRB0800030A	PRB0900030A	PRB1000030A		
40	PRB0700040A	PRB0800040A	PRB0900040A	PRB1000040A		
60	PRB0700060A	PRB0800060A	PRB0900060A	PRB1000060A		
80	PRB0700080A	PRB0800080A	PRB0900080A	PRB1000080A		
100	PRB0700100A	PRB0800100A	PRB0900100A	PRB1000100A		
120	PRB0700120A	PRB0800120A	PRB0900120A	PRB1000120A		
150	PRB0700150A	PRB0800150A	-	-		
180	PRB0700180A	PRB0800180A	-	-		
(OTW 0.035" 120cm)						
Nominal Stent Length (mm)	Unconstrained stent Diameter (mm)					
	5.00	6.00	7.00	8.00	9.00	10.00
	Reference Number					
30	PRB0500030B	PRB0600030B	PRB0700030B	PRB0800030B	PRB0900030B	PRB1000030B
40	PRB0500040B	PRB0600040B	PRB0700040B	PRB0800040B	PRB0900040B	PRB1000040B
60	PRB0500060B	PRB0600060B	PRB0700060B	PRB0800060B	PRB0900060B	PRB1000060B
80	PRB0500080B	PRB0600080B	PRB0700080B	PRB0800080B	PRB0900080B	PRB1000080B
100	PRB0500100B	PRB0600100B	PRB0700100B	PRB0800100B	PRB0900100B	PRB1000100B
120	PRB0500120B	PRB0600120B	PRB0700120B	PRB0800120B	PRB0900120B	PRB1000120B
150	PRB0500150B	PRB0600150B	PRB0700150B	PRB0800150B	-	-
180	PRB0500180B	PRB0600180B	PRB0700180B	PRB0800180B	-	-

Table 3 - Product Description

Stent Material	The stent is made by laser cutting seamless tubing of Nitinol alloy and electro polished. The stent has 06/08 radiopaque markers (03/04 at each end).	
Stent delivery catheter system	The delivery system is an OTW PTA catheter with three coaxial tubes viz. inner lumen, pusher tube and braided sheath. The inner lumen has two platinum-iridium radiopaque markers. The stent is mounted between these markers and hence they define location of the stent. The guide wire lumen within inner lumen covers the entire length of the catheter i.e. from hub to the soft tip. The braided sheath extends from Puller Shaft to the soft tip and covers the entire length of the stent restraining it from expansion.	
Delivery system usable length	The delivery systems are available in two usable lengths viz. 80cm and 120cm.	
Outer braided sheath profile	6F	
Recommended Minimum Sheath/ Introducer	6F*	
Guide wire compatibility	OTW 0.035" (0.89mm)	
Product code format	PRB= Promesa BMS xxxx=nominal stent diameter (mm) yyy=nominal stent length (mm) A=OTW 0.035" 80cm) For Example. PRB0700030A xxxx=0700=Diameter 7.00mm yyy=030=Length 30mm A=80cm catheter length	PRB= Promesa BMS xxxx=nominal stent diameter (mm) yyy=nominal stent length (mm) B=OTW 0.035" 120cm) For Example. PRB0500030B xxxx=0500=Diameter 5.00mm yyy=030=Length 30mm B=120cm catheter length

*see the individual manufacturer specification for (F) equivalent on box label and pouch label.

2. HOW SUPPLIED:

Sterile: Promesa™ BMS is sterilized with Ethylene oxide (ETO) gas and is non-pyrogenic. It is intended for single use only. Do not re-sterilize. Do not use the device if the package is opened or damaged.

Contents: One (1) Promesa™ BMS housed in a protective circular hoop tray which fixed onto a product holding tray, one (1) Instructions for Use, two (2) Stent Implant Cards.

Storage: Store in a cool, dry and dark place in its original packaging below 30° C.

Shelf Life: 3 years.

3. INDICATIONS:

The Promesa™ BMS - Self-Expanding Nitinol Peripheral Stent System is intended for use in the treatment of atherosclerotic diseases of peripheral arteries with reference vessel diameter of 5.00 mm to 7.00 mm for Superficial Femoral Artery (SFA) and 7.00 mm to 10.00 mm for Iliac Artery in patients eligible for Percutaneous Transluminal Angioplasty (PTA) and Stenting procedures.

4. CONTRAINDICATIONS:

The Promesa™ BMS is contraindicated but may not be limited to the following patient types:

- Patients with hypersensitivity or allergies to titanium, cobalt, chromium, nickel, molybdenum, tungsten, or any contrast media.
- Patients in whom anti-platelet and/or anticoagulant therapy are contraindicated.
- Transplant patients

5. WARNINGS:

- This device is intended for single-use only; do not re-sterilize. Do not use if the package is open or damaged. Reuse, reprocessing or re-sterilization may increase the risk of patient infection and risk of compromised device performance.
- DO NOT use after the "Use By" date specified on the label.
- The lesion length should be less than 10mm with respect to the stent length. For example; the lesion length may not be more than 170mm for the stent with 180mm length.
- Persons with known hypersensitivities to nitinol and / or its components (e.g., nickel, titanium) may suffer an allergic reaction to this implant.
- Use of appropriate anticoagulant and / or antiplatelet therapy per standard of care is recommended for use with this stent system.
- Once the stent is partially deployed, it cannot be recaptured using the stent system.
- If resistance is met during advancement of the delivery system, do not force passage. Remove the delivery system and replace with a new device.
- Improper stent size selection may lead to stent migration or stent jumping.
- Patients judged to have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the stent or stent delivery system.
- The product is intended for use of physicians trained and experienced in diagnostic and interventional vascular technique.

6. PRECAUTIONS:**6.1 General Precautions:**

- DO NOT use if the sterile package is opened or damaged or any information provided is obscured.
- DO NOT use if the device is damaged or if the stent is partially deployed.
- Do not try to remove the stent from the introducer system before use.
- Do not expose the delivery system to organic solvents (e.g., alcohol).
- Do not use power injection systems with the delivery system.
- Do not rotate any part of the system during deployment.
- Repositioning of the device after deployment is not possible.
- Re-sheathing / Recapturing of the stent is not possible, once the stent has apposed the vessel wall.
- Do not use the Promesa™ BMS in conjunction with other stents.

6.2 Stent Handling Precautions:

- Do not remove the stent from its delivery system, as removal may damage the stent and / or lead to stent embolization. Stent system is intended to perform as a system.
- Special care must be taken not to handle or in any way disrupt the stent on the delivery system. This is most important during system removal from packaging and device insertion.

6.3 Stent Placement Precautions:

- Advance the delivery system past the lesion and pull back to help remove slack from the system. Removing all slack from the delivery system prior to stent deployment will help ensure accurate stent delivery.
- Do not attempt to pull a partially expanded stent back through the sheath or guiding catheter; dislodgment of the stent from the delivery system may occur.
- Do not expand the stent if it is not properly positioned in the vessel. (See Stent / System Removal – Precautions.)

6.4 Stent / System Removal Precautions :

- Do not attempt to pull a partially-expanded stent back through the introducer sheath or guiding catheter. The stent is not designed for recapturing. The stent is not designed for repositioning once the stent has apposed the vessel. Once the stent is apposed to the vessel, it is not recommended to reposition the stent with the delivery system.
- Should unusual resistance be felt at any time during lesion access or removal of the delivery system post stent implantation, the entire system should be removed together with the introducer sheath or guiding catheter as a single unit. Failure to follow these instructions could result in failure to deploy, difficulties with deployment, partial stent deployment or deployment in an unintended location.
- If any unusual resistance be felt at any time including rotating the handle gear (or thumbwheel), during stent deployment, the entire system should be removed together with the introducer sheath or guiding catheter as a single unit.
- Do not apply excessive force if resistance is felt during delivery system retrieval through the deployed stent. Slightly advance the delivery system and continue withdrawal of the delivery system.

6.5 Post Implant Precautions:

Great care must be exercised when crossing a newly deployed stent with other devices such as another stent delivery system an intravascular Ultrasound (IVUS) catheter, a guide wire or balloon catheter to avoid disrupting the stent geometry.

7. ADVERSE EFFECTS:

Possible Adverse Effects that may occur include, but are not limited to, the following::

- Allergic reaction to anticoagulant and/or antithrombotic therapy or contrast medium
- Allergic reaction to nitinol
- Atheromolization (Blue toe Syndrome)
- Artrial aneurysm
- Artrial rupture
- Artrial thrombosis
- Arteriovenous fistula
- Death
- Embolism
- Hematoma/hemorrhage
- Hypotension/hypertension
- Hypersensitivity reactions
- Infection, local or systemic, including bacteremia or septicemia
- Infection/abscess formation at access site
- Infection and sepsis
- Ischemia requiring intervention (bypass or amputation of toe, foot, or leg)
- Pseudoaneurysm formation
- Renal Failure
- Restenosis of the stented artery
- Stent embolization
- Stent malapposition
- Stent migration
- Stent strut fracture
- Stent thrombosis or occlusion
- Vessel perforation or rupture

- Worsened claudication/rest pain

8. INDIVIDUALIZATION OF TREATMENT:

- The risk and benefits should be considered for each patient before use of Promesa™ BMS. Patient selection factors should include a judgement regarding risk of antiplatelet therapy. Special consideration should be given to those patients with recently active gastritis or peptic ulcer disease.
- Pre-morbid conditions that increase the risk of a poor initial result or the risk of emergency referral for bypass surgery (diabetes mellitus, renal failure and severe obesity), should be reviewed.
- A review of the vessel location, reference vessel size, lesion length, qualitative target lesion characteristics, and implication of acute or subacute thrombosis must also be considered.
- Thrombosis following stent implantation is affected by several baseline angiographic and procedural factors. These include intra-procedural thrombus and dissection following stent implantation. In patients who have undergone peripheral stenting, the persistence of a thrombus or dissection should be considered a marker for subsequent thrombotic occlusion. These patients should be monitored very carefully during the first month after stent implantation.

9. USE IN SPECIAL POPULATIONS:

- The safety and effectiveness of the Promesa™ BMS has not been established in the following patient populations:
- Patient with unresolved vessel thrombus at the lesion site
 - Patients with brachytherapy treatment of the target lesion
 - **Pregnant patient:** There are no adequate and well controlled studies in pregnant women or men intending to father children. Effective contraception should be initiated before implanting Promesa™ BMS and for 12 weeks after implantation. The Promesa™ BMS should be used during pregnancy only if the potential benefit outweighs the potential risk to the embryo or foetus.
 - The safety and effectiveness of using brachytherapy treatment, mechanical arterectomy devices (directional arterectomy catheters, rotational arterectomy catheters) or laser angioplasty catheters to treat in stent stenosis of a Promesa™ BMS has not yet been established.

10. CLINICAL USE INFORMATION:

10.1 Inspection prior to use:

- Carefully inspect the sterile package before opening.
- Do not use if the package has been damaged or opened.
- The product should not be used after the "Use By" Date mentioned on the label.
- If the sterile package appears intact, tear open the sterile pouch and carefully remove the product in the sterile field using aseptic technique.
- Remove the system carefully from tray packing.
- Inspect the stent system for bends, kinks and other damage.
- Do not use if any of the above defects is noted.

10.2 Materials Required:

- Introducer Sheath
- 2-3 Syringes (10-20cc)
- 1000u/500cc, normal heparinised saline (HepNS)
- 0.035" (0.89mm) diameter guide wire
- Contrast diluted 1:1 with normal saline
- Inflation device
- Three-way stopcock
- Torque Device (If applicable)
- Guide Wire Introducer

10.3 Selection of stent size:

- Measure the diameter of the reference vessel (proximal and distal to the lesion or obstruction) and use the largest reference diameter as your basis for choosing the appropriate stent size. Select stent as recommended in table below:

Table-4

Labelled stent diameter (mm)	Reference vessel diameter (mm)
5.00	3.00 - 4.50
6.00	4.00 - 5.50
7.00	5.00 - 6.50
8.00	6.00 - 7.50
9.00	7.00 - 8.50
10.0	8.00 - 9.50

Warning: Improper stent size selection may lead to stent migration or stent jumping.

10.4 Multiple stent Placement:

- If placement of multiple stents is required in a patient, to cover the length of the lesion, the following recommendations should be considered;
- In relation to the lesion site, the distal area of lesion may be stented first, followed by the proximal locations (i.e. second stent should be placed proximally to the previously placed stents).
 - Stent placed in tandem must overlap to allow for complete coverage of the lesion.

10.5 Deployment procedure:

- Prepare vascular access site according to standard practice.

- Predilatation before stent placement is optional and at the discretion of the physician.

A) Delivery system removal

1. After full stent deployment, carefully remove the delivery system under fluoroscopic guidance while observing the radiopaque tip, leaving the guide wire in place.

Note: If resistance is met during the withdrawal of delivery system through the stent, re-advance the braided sheath over the inner lumen to its pre-deployment position. Withdraw the system in this position. If resistance is still encountered while attempting to withdraw the delivery system, the entire system should be removed together with the introducer sheath or guiding catheter as a single unit.

B) Flush the delivery system

Use the 5ml syringe to flush the wire guide lumen with saline by connecting the syringe to the Luer hub (Refer Figure 2 below).

1

Fig. 2 Delivery System Flushing

Note: The time between flushing and stent deployment should be minimized and must be no longer than 5 minutes.

Under fluoroscopy, advance the delivery system over the wire guide through the introducer sheath / guiding catheter until the distal tantalum radiopaque markers on the stent are beyond the target lesion site (Refer Figure below).

2

Fig. 3 Advancing the delivery system across the lesion

C) Alignment of the stent

Before stent alignment, it is important to:

- a. Straighten the proximal part of the delivery system as much as possible (Refer Figure 5 below).
- b. Keep the Handle in a stable position.

1

Fig. 4 Keep the part of the delivery system straight : do not bend

2

The stent alignment must be performed under fluoroscopic control.

Align the radiopaque stent markers to the desired position. The Stent should now be ready to be deployed. (Refer Figure below)

3

Fig. 5 Stent Alignment

D) Deployment of stent

Before stent deployment, it is important to:

- a. Straighten the proximal part of the delivery system as much as possible.
- b. Keep the Handle in a stable position. (Refer Figure below).

1

Fig. 6 Keep the part of the delivery system straight : do not bend

Hold the hub firmly in one hand and with other hand, rotate the thumbwheel in clockwise direction to release the stent from delivery system. Click sound will be audible when stent is released from the delivery system.

2

Fig. 7 Stent Deployment: Rotating the Thumbwheel

Note: Care should be taken to hold the hub stationary and to remove any slack in the stent system to ensure the stent is not stretched or compressed lengthwise during deployment (i.e., so that the stent is deployed to its proper length). Full deployment of the stent length will occur when the distal end of the braided sheath has been retracted past the proximal part of the stent.

As deployment occurs, continue rotating the thumbwheel in a slow, smooth and consistent manner (Refer Figure below) until complete stent is deployed or thumbwheel frees up and puller shaft appears from the proximal end of the handle.

3

Fig. 8 Stent Deployment: Rotating the Thumbwheel till Puller Shaft appears

Note: Once stent deployment has begun, the stent must be fully deployed. Repositioning of the Promesa™ BMS Stent is not possible since the delivery system's pusher coil cannot be re-advanced over the stent once deployment begins. Refer the Multiple Stent Placement section of this Instructions for Use for information on missed lesions

Instructions for deployment of Stent Length above 120mm

Now pull back the blue puller shaft from the proximal end of the handle until it locks with the handle. It will confirm complete stent deployment. (Refer Figure below).


Fig. 9 Stent Deployment: Pull back the puller shaft

E) Delivery system removal

- After full stent deployment, carefully remove the delivery system under fluoroscopic guidance while observing the radiopaque tip, leaving the guide wire in place.
Caution: Do not apply excessive force if resistance is felt during delivery system retrieval through the deployed stent. Slightly advance the delivery system and continue withdrawal of the delivery system.
- Perform the arterial angiogram to verify full deployment of the device. If incomplete expansion exists within the stent at any point along the lesion, post-deployment balloon dilatation (standard PTA) can be performed at the discretion of the physician.
- Check the delivery system after removal from the patient to ensure that the tip and inner catheter have not separated. Maintain wire guide access until this check is completed.
- Remove the guide wire and introducer sheath from patient.
- After use, dispose of the product and packaging in accordance with hospital, administrative and/or local government policy.

10.6 Post procedural:

- If the stent is incomplete expanded along the lesion, post deployment balloon dilatation can be performed. Select a PTA balloon whose inflated diameter matches the reference diameter of the target vessel.
- Obtain angiographic visualization of the stented vessel segment.

11. Disclaimer of Warranty and Limitation of Remedy:

There is no express or implied warranty, including without limitation any implied warranty of merchantability of fitness for a particular purpose, on the Meril Life Sciences Pvt. Ltd. product(s) described in this publication. Under no circumstances shall Meril Life Sciences Pvt. Ltd. be liable for any direct, indirect, incidental or consequential damages. Some jurisdictions do not allow the exclusion of or limitations on an implied warranty. For jurisdictions which do not allow the exclusion or limitations of incidental or consequential damages, some of the above exclusions may not apply. The patient may also have other rights, which vary from jurisdiction to jurisdiction.

Descriptions or specifications in Meril Life Sciences Pvt. Ltd. printed matter, including this publications, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Meril Life Sciences Pvt. Ltd. will not be responsible for any direct, indirect, incidental, or consequential damages resulting from reuse of the product.

Symbols used in labeling

Stent diameter	Stent length	Do not use if package is damaged	Non-pyrogenic
Contains one unit	Keep dry	Reference number	Lot number
Max. guide wire diameter	Use-by date	Consult instructions for use	Manufacturer
Sterilized using ethylene oxide	Date of manufacture	Do not re-use	Do not re-sterilize
Upper limit of temperature	Keep away from sunlight	Caution	
Authorized representative in the European community			

Meril

CE
1783

An ISO 13485 & GMP Certified Company

Manufactured by:
Meril Life Sciences Pvt. Ltd.
Muktanand Marg,
Chala, Vapi 396191,
Gujarat, India.
Web: www.merillife.com

Obelis s.a.
Bd., General Wahis 53,
1030, Brussels, Belgium.
Tel.: +32.2.732.5954
Fax: +32.2.732.6003
E-mail: mail@obelis.net

Customer Care Contact:
Tel.: +91 (260) 240 8000
E-mail: askinfo@merillife.com

IFU9509/28122020