S90, S60, S30, U100, U200, U500

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
The diagram illustrates the electrical architecture of a dental unit, organized into several functional sections:

- Power Section:** Includes a 230V AC input, a fuse (FUSE T10A), a power switch, and a terminal block. A transformer provides multiple secondary voltages (24V, 18V, 6V, 0V, 10V, 12V, 17V).
- Control and Interface Section:** Features the Assistant Touch PCB (U500ZS-03-F5) with a touch screen, and the Optic-control PCB (U500ZL-02-F5) which manages the Optic-linkage PCB (U500ZL-02-F0) and the Dentist Touch PCB (U500ZL-02-F5). It also includes a Sensor PCB and an LBD-01-F5 PCB.
- Actuator and Solenoid Section:** Contains various solenoids (J24, J25, J26, J27, J28, J29, J30, J31, J32, J33, J34, J35, J36, J37, J38, J39, J40, J41, J42, J43, J44, J45, J46, J47, J48, J49, J50, J51, J52, J53, J54, J55, J56, J57, J58, J59, J60, J61, J62, J63, J64, J65, J66, J67, J68, J69, J70, J71, J72, J73, J74, J75, J76, J77, J78, J79, J80, J81, J82, J83, J84, J85, J86, J87, J88, J89, J90, J91, J92, J93, J94, J95, J96, J97, J98, J99, J100, J101, J102, J103, J104, J105, J106, J107, J108, J109, J110, J111, J112, J113, J114, J115, J116, J117, J118, J119, J120, J121, J122, J123, J124, J125, J126, J127, J128, J129, J130, J131, J132, J133, J134, J135, J136, J137, J138, J139, J140, J141, J142, J143, J144, J145, J146, J147, J148, J149, J150, J151, J152, J153, J154, J155, J156, J157, J158, J159, J160, J161, J162, J163, J164, J165, J166, J167, J168, J169, J170, J171, J172, J173, J174, J175, J176, J177, J178, J179, J180, J181, J182, J183, J184, J185, J186, J187, J188, J189, J190, J191, J192, J193, J194, J195, J196, J197, J198, J199, J200, J201, J202, J203, J204, J205, J206, J207, J208, J209, J210, J211, J212, J213, J214, J215, J216, J217, J218, J219, J220, J221, J222, J223, J224, J225, J226, J227, J228, J229, J230, J231, J232, J233, J234, J235, J236, J237, J238, J239, J240, J241, J242, J243, J244, J245, J246, J247, J248, J249, J250, J251, J252, J253, J254, J255, J256, J257, J258, J259, J260, J261, J262, J263, J264, J265, J266, J267, J268, J269, J270, J271, J272, J273, J274, J275, J276, J277, J278, J279, J280, J281, J282, J283, J284, J285, J286, J287, J288, J289, J290, J291, J292, J293, J294, J295, J296, J297, J298, J299, J300, J301, J302, J303, J304, J305, J306, J307, J308, J309, J310, J311, J312, J313, J314, J315, J316, J317, J318, J319, J320, J321, J322, J323, J324, J325, J326, J327, J328, J329, J330, J331, J332, J333, J334, J335, J336, J337, J338, J339, J340, J341, J342, J343, J344, J345, J346, J347, J348, J349, J350, J351, J352, J353, J354, J355, J356, J357, J358, J359, J360, J361, J362, J363, J364, J365, J366, J367, J368, J369, J370, J371, J372, J373, J374, J375, J376, J377, J378, J379, J380, J381, J382, J383, J384, J385, J386, J387, J388, J389, J390, J391, J392, J393, J394, J395, J396, J397, J398, J399, J400, J401, J402, J403, J404, J405, J406, J407, J408, J409, J410, J411, J412, J413, J414, J415, J416, J417, J418, J419, J420, J421, J422, J423, J424, J425, J426, J427, J428, J429, J430, J431, J432, J433, J434, J435, J436, J437, J438, J439, J440, J441, J442, J443, J444, J445, J446, J447, J448, J449, J450, J451, J452, J453, J454, J455, J456, J457, J458, J459, J460, J461, J462, J463, J464, J465, J466, J467, J468, J469, J470, J471, J472, J473, J474, J475, J476, J477, J478, J479, J480, J481, J482, J483, J484, J485, J486, J487, J488, J489, J490, J491, J492, J493, J494, J495, J496, J497, J498, J499, J500, J501, J502, J503, J504, J505, J506, J507, J508, J509, J510, J511, J512, J513, J514, J515, J516, J517, J518, J519, J520, J521, J522, J523, J524, J525, J526, J527, J528, J529, J530, J531, J532, J533, J534, J535, J536, J537, J538, J539, J540, J541, J542, J543, J544, J545, J546, J547, J548, J549, J550, J551, J552, J553, J554, J555, J556, J557, J558, J559, J560, J561, J562, J563, J564, J565, J566, J567, J568, J569, J570, J571, J572, J573, J574, J575, J576, J577, J578, J579, J580, J581, J582, J583, J584, J585, J586, J587, J588, J589, J590, J591, J592, J593, J594, J595, J596, J597, J598, J599, J600, J601, J602, J603, J604, J605, J606, J607, J608, J609, J610, J611, J612, J613, J614, J615, J616, J617, J618, J619, J620, J621, J622, J623, J624, J625, J626, J627, J628, J629, J630, J631, J632, J633, J634, J635, J636, J637, J638, J639, J640, J641, J642, J643, J644, J645, J646, J647, J648, J649, J650, J651, J652, J653, J654, J655, J656, J657, J658, J659, J660, J661, J662, J663, J664, J665, J666, J667, J668, J669, J670, J671, J672, J673, J674, J675, J676, J677, J678, J679, J680, J681, J682, J683, J684, J685, J686, J687, J688, J689, J690, J691, J692, J693, J694, J695, J696, J697, J698, J699, J700, J701, J702, J703, J704, J705, J706, J707, J708, J709, J710, J711, J712, J713, J714, J715, J716, J717, J718, J719, J720, J721, J722, J723, J724, J725, J726, J727, J728, J729, J730, J731, J732, J733, J734, J735, J736, J737, J738, J739, J740, J741, J742, J743, J744, J745, J746, J747, J748, J749, J750, J751, J752, J753, J754, J755, J756, J757, J758, J759, J760, J761, J762, J763, J764, J765, J766, J767, J768, J769, J770, J771, J772, J773, J774, J775, J776, J777, J778, J779, J780, J781, J782, J783, J784, J785, J786, J787, J788, J789, J790, J791, J792, J793, J794, J795, J796, J797, J798, J799, J800, J801, J802,

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S90, S60, S30, U100, U200, U500

**USER MANUAL
FOR
SIGER DENTAL UNIT**

ZHUHAI SIGER MEDICAL EQUIPMENT CO., LTD

Manufacturer: Zhuhai Siger Medical Equipment Co., Ltd.

Manufacture & Registration Address: Building 2, No. 1 Chuangxin Yi Road, Tangjiawan Town, Zhuhai City,
Guangdong Province, P.R. China

Tel: +86-756-3881018

Fax: +86-756-3881028

Website: <http://www.siger.cn>

E-mail: info@siger.cn

Post Code: 519000

Doc No.: SG-S90-RD-57

Version: A1

Date of issue: 11/2021

EU Authorized Representative: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg, Germany

Tel: +49-40-2513175

Fax: +49-40-255726

email: shholding@hotmail.com

Foreword

Before any operation is performed, please read the operation instructions carefully and keep it for future reference. Please strictly follow the operation procedures during operation and perform proper maintenance.

Please read "WARNING" and "NOTE" paragraphs carefully and strictly perform them to avoid any damage to the unit or injury to the operator or patients.

If any malfunction occurs during the use, please contact the local dealer or SIGER without delay.

• Symbols

WARNING

This indicates a situation which may lead to injuries to persons or damages to the unit and/or its parts if the information is ignored.

NOTE

This indicates the important information for the user and/or service personnel.

Guard against damp.

Fragile. Handle with care.

This side up.

Storage temperature limit sign.

Storage humidity limit sign.

Storage Atmospheric Pressure limit sign.

No stowing under other cargo.

See the random file.

It is forbidden to discard the marking of electronic equipment casually, and after the expiry date, the equipment should be handled according to local legal requirements to avoid polluting the environment and injuring the user.

Classification: Type I Class B.

Ground mark.

Opening sign.

Closing sign.

IPX4 Foot switch waterproof rating mark.

• Safety Information

- Please always follow the following safety precautions to avoid injuries to persons or damages to the unit:
- Please follow all the warnings and instructions provided with the unit. If an instruction contradicts with the safety information, please pay attention to safety information, for you may misunderstand the instruction. If you cannot solve the contradiction, please contact our service personnel to seek help.
- Please always shut off the water supply, air supply and the main power before performing any cleaning or maintenance.
- Please do not place the unit on shaky floor, cart and shelf or it may tumble over and badly damaged.
- Please do not install the unit near radiators or heaters.
- This unit is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide, where such gasses may accumulate in concentration (in closed space).
- Please do not place or hang any heavy things on the unit.
- Please do not connect unit with incompatible wires which may affect its performance, or even cause fire or electric shock.
- Please do not place anything in the unit so as to avoid electric shock or fire. Please do not spill any liquid into the inner part of the unit.

- Please do not disassemble the unit so as to avoid electric shock. Contact technical service personnel whenever the unit needs repairing.
- Please do not open or remove the cover plate so as to avoid injuries caused by dangerous voltage or other dangers. An improper assembly may cause electric shock during the following operation.
- Please switch off the power and seek help from technical service personnel if:
 - (1) any break-down or wear-down of the power, plug or cable occurs.
 - (2) unknown liquid spills into the inner unit or water flows into the inner circuits.
 - (3) the unit malfunctions, even if it is operated according to the instructions.
- Please only adjust the control items mentioned in the instructions. Any adjustment of other control items may cause damage to the unit.
- Please avoid operating the unit during thunderstorm or electric shock may occur. If possible, please turn off the power supply during thunderstorm.
- Please do not use damaged or loose plug. An incomplete insertion of plug may cause electric shock or fire.
- Please make sure to use earthed plug and AC outlet. Bad earthing may cause electric shock to the user or damage to the unit.

- ⚠ WARNING!** When treating patients who has pacemakers and/or hearing aids, please always consider the effects the dental unit may have on them.
- ⚠ WARNING!** The maximum load of the dental seat is 150kg. Please never exceed this limit.
- ⚠ WARNING!** The maximum load of the table tray attached to the dentist's instrument board is 3.0kg, which should never be exceeded.
- ⚠ WARNING!** Replace the dental burrs only when the dental instrument has stopped completely, otherwise the cartridge may be damaged and the dental burrs may cause a personal injury.
- ⚠ WARNING!** This unit is used by authorized and fully trained personnel only (doctors or assistants). The manufacturer shall not be responsible for any damage or injury caused by mishandling, neglect or improper use.
- ⚠ WARNING!** Before leaving the surgery, shut off the water supply, air supply and the main power.
- ⚠ WARNING!** Wastewater and waste emissions should comply with national laws and regulations, and the amalgam collection is equipped with an amalgam separator.
- ⚠ WARNING!** The device has a one-way stop valve on the total inlet to prevent backflow, the outlet and other pollutants are greater than 20mm spacing, distilled water system and tap water system completely independent. Before use to ensure its reliability.
- ⚠ WARNING!** The use of gas requirements should comply with national laws and regulations, no gas requirements of laws and regulations of the country, should ensure that the gas source is not oil, water, bacteria and harmful substances.

The equipment is expected to have a service life of five years. Production date see nameplate

Please contact with manufacturer and local competent authority if have any serious incident related with the dental units.

• General Information

Information of optional parts of S/U series is included in this manual which may not suit for your units.

These instructions describe the operation and maintenance procedures for the dental units of S/U series. And all the information provided is available till the publishment of the Manual. All the information, specifications and illustrations herein are subject to change. SIGER reserves the right to make changes without prior notice.

No addition, modification, distribution, republication or sale of this manual, in any form or any manner (electronic, mechanic, or graphic), are allowed unless it is authorized by Siger in a written form in advance.

SIGER reserves the right to make explanation and changes of all the information contained in the instruction.

⚠ NOTE! The guarantee will be disclaimed when the following precautions are not taken:

- (1) Observe the terms specified in the guarantee letter.
- (2) The unit must be used according to the instructions herein.
- (3) During indoor electric wiring, a $3 \times 1 \text{ mm}^2$ power cable should be provided.
- (4) The unit should be installed in a fixed way.
- (5) Repairs, modifications, adjustments and any other operations which involve gaining access to the inside of the unit must be carried out by SIGER technicians or authorized technicians only.
- (6) Replacements for any spare parts and instruments of the unit should be provided by SIGER otherwise the performance and safety may be affected. Regarding the spare parts, please refer to the packing list.

— Installers and maintenance personnel must meet the following requirements:

- (a) Education:
 - At least 18 years old – 9 years intensive reading experience (school).
- (b) Knowledge:
 - Read and understand 'westernized Arabic' numerals when written in Arial font.
 - Read this manual before use.
- (c) Language understanding:
 - Can understand English language in the User Manual and displayed on touchscreen console.
- (d) Permissible impairments:
 - With reading vision impairment or vision corrected to log MAR 0, 2 (6/10 or 20/32)
- (e) Experience
 - Training by SIGER or fully trained for dental unit.

The user must meet the following requirements:

- (a) Person Group and Occupation: Adult-qualified dentist and assistant
- (b) Education:
 - At least 18 years old – 9 years intensive reading experience (school).
- (c) Knowledge:
 - Read and understand 'westernized Arabic' numerals when written in Arial font.
 - Read this manual before use.
 - Understand stomatology.
 - Understand hygiene.
- (d) Language understanding:
 - Can understand English language in the User Manual and displayed on touchscreen console.
- (e) Permissible impairments:
 - With reading vision impairment or vision corrected to log MAR 0, 2 (6/10 or 20/32)
- (f) Experience
 - Fully trained for dental unit.

Operating principle

Dental Units is a combination of dental equipment components used to carry, position and control the devices used in the practice of dentistry.

Contents

Chapter 1 Description of the Unit	1
1.1 Intended Use	1
1.2 Description of the Unit	1
Chapter 2 Control Panels, Display and Buttons of S90/S60/S30/U200 Dental Unit	2
2.1 Operation Panel and Keys Description	3
2.2 Parameter Setup and Display	4
2.3 Status and Functions of Doctor's and Assistant's Control Panel	9
Chapter 3 Control Panels and Buttons of U500/U100 Dental Unit	11
3.1 Description of Control Panels and Buttons	11
3.2 Status and Functions of Doctor's and Assistant's Control Panel	12
Chapter 4 Operation Instruction	13
4.1 Main Switch	13
4.2 Operation Light	13
4.3 Foot pedal	14
4.4 Image	16
4.5 Dental instrument	16
4.6 Air Suction	17
4.7 Saliva Ejector	17
4.8 Height Adjustment of Balance Arm	17
4.9 Instrument Tray	17
4.10 Swivel of Swing Type Instrument Table	17
4.11 Restent	17
4.12 Right Armrest (optional)	17
4.13 Height Adjustment of Assistant's Instrument Board	18
4.14 Distilled Water System	18
4.15 Disinfection of the Pipeline	18
4.16 Auxiliary Devices	18
Chapter 5 Cleaning, Disinfection and Maintenance	19
5.1 Disinfection	19
5.2 Maintenance of the Unit	20
Chapter 6 Dental Unit installation, commissioning and maintenance	20
6.1 Installation Conditions	20
6.2 Requirements for Installation Field	20
6.3 Loading and Unpacking of Equipment	21
6.4 Equipment installation method	22
6.5 Equipment debugging method	22
6.6 Water sampling connection point	23
6.7 Water Air schematic diagram and electrical schematic diagram	24
Chapter 7 Common Faults and Exclusion	35
7.1 Dental Chair	35

12.1 Operating Light.....	35
12.2 Instrument Table.....	35
12.3 Replacement of fuse tube.....	36
Chapter 6 Electromagnetic Compatibility.....	36
Chapter 6 Safety Markings.....	39
Chapter 6 Composition of the Unit.....	40

1. Composition of the Unit

1.1. Composition of the unit

Number of main parts

1. Main unit	1	1. Main unit
2. Control panel	1	2. Control panel
3. Instrument table	1	3. Instrument table
4. Control panel	1	4. Control panel
5. Instrument table	1	5. Instrument table
6. Control panel	1	6. Control panel
7. Instrument table	1	7. Instrument table
8. Control panel	1	8. Control panel
9. Instrument table	1	9. Instrument table
10. Control panel	1	10. Control panel
11. Instrument table	1	11. Instrument table
12. Control panel	1	12. Control panel
13. Instrument table	1	13. Instrument table
14. Control panel	1	14. Control panel
15. Instrument table	1	15. Instrument table
16. Control panel	1	16. Control panel
17. Instrument table	1	17. Instrument table
18. Control panel	1	18. Control panel
19. Instrument table	1	19. Instrument table
20. Control panel	1	20. Control panel
21. Instrument table	1	21. Instrument table
22. Control panel	1	22. Control panel
23. Instrument table	1	23. Instrument table
24. Control panel	1	24. Control panel
25. Instrument table	1	25. Instrument table
26. Control panel	1	26. Control panel
27. Instrument table	1	27. Instrument table
28. Control panel	1	28. Control panel
29. Instrument table	1	29. Instrument table
30. Control panel	1	30. Control panel
31. Instrument table	1	31. Instrument table
32. Control panel	1	32. Control panel
33. Instrument table	1	33. Instrument table
34. Control panel	1	34. Control panel
35. Instrument table	1	35. Instrument table
36. Control panel	1	36. Control panel
37. Instrument table	1	37. Instrument table
38. Control panel	1	38. Control panel
39. Instrument table	1	39. Instrument table
40. Control panel	1	40. Control panel
41. Instrument table	1	41. Instrument table
42. Control panel	1	42. Control panel
43. Instrument table	1	43. Instrument table
44. Control panel	1	44. Control panel
45. Instrument table	1	45. Instrument table
46. Control panel	1	46. Control panel
47. Instrument table	1	47. Instrument table
48. Control panel	1	48. Control panel
49. Instrument table	1	49. Instrument table
50. Control panel	1	50. Control panel
51. Instrument table	1	51. Instrument table
52. Control panel	1	52. Control panel
53. Instrument table	1	53. Instrument table
54. Control panel	1	54. Control panel
55. Instrument table	1	55. Instrument table
56. Control panel	1	56. Control panel
57. Instrument table	1	57. Instrument table
58. Control panel	1	58. Control panel
59. Instrument table	1	59. Instrument table
60. Control panel	1	60. Control panel
61. Instrument table	1	61. Instrument table
62. Control panel	1	62. Control panel
63. Instrument table	1	63. Instrument table
64. Control panel	1	64. Control panel
65. Instrument table	1	65. Instrument table
66. Control panel	1	66. Control panel
67. Instrument table	1	67. Instrument table
68. Control panel	1	68. Control panel
69. Instrument table	1	69. Instrument table
70. Control panel	1	70. Control panel
71. Instrument table	1	71. Instrument table
72. Control panel	1	72. Control panel
73. Instrument table	1	73. Instrument table
74. Control panel	1	74. Control panel
75. Instrument table	1	75. Instrument table
76. Control panel	1	76. Control panel
77. Instrument table	1	77. Instrument table
78. Control panel	1	78. Control panel
79. Instrument table	1	79. Instrument table
80. Control panel	1	80. Control panel
81. Instrument table	1	81. Instrument table
82. Control panel	1	82. Control panel
83. Instrument table	1	83. Instrument table
84. Control panel	1	84. Control panel
85. Instrument table	1	85. Instrument table
86. Control panel	1	86. Control panel
87. Instrument table	1	87. Instrument table
88. Control panel	1	88. Control panel
89. Instrument table	1	89. Instrument table
90. Control panel	1	90. Control panel
91. Instrument table	1	91. Instrument table
92. Control panel	1	92. Control panel
93. Instrument table	1	93. Instrument table
94. Control panel	1	94. Control panel
95. Instrument table	1	95. Instrument table
96. Control panel	1	96. Control panel
97. Instrument table	1	97. Instrument table
98. Control panel	1	98. Control panel
99. Instrument table	1	99. Instrument table
100. Control panel	1	100. Control panel

1.2. Composition of the unit with dental unit

The unit is composed of the main unit, control panel, instrument table, and various accessories. The unit is designed to be used in a dental clinic.

The unit is composed of the main unit, control panel, instrument table, and various accessories. The unit is designed to be used in a dental clinic.

1.3. Composition of the unit with hanging type and hanging type

Chapter 1 Description of the Unit

1.1 Intended Use

The dental units are usually AC-powered, in which a patient sits. The device is intended to properly position a patient to perform dental procedures. A dental operative unit is attached.


Population: adult.

Contraindications: None.

Side-effect: The side-effects of dental unit are the noise, eye fatigue of the dental doctor by the dental light.

1.2 Description of the Unit

1.2.1 Panorama of the unit

Names of main parts

1	HV suction tube	9	Headrest
2	Assistant's control panel	10	Foot pedal
3	Saliva ejector tube	11	Patient's chair
4	Assistant's syringe	12	Left armrest
5	Dental spittoon	13	Doctor's syringe
6	Hydro-group	14	Operating light
7	Instrument table arm	15	Doctor's control panel
8	Instrument table		


Note: syringe means the 3-way syringe

1.2.2 The dental operative unit connected with dental unit

Dental saline handpiece, dental suction handpiece (HV suction and Saliva ejector) and vacuum pump can be connected with dental unit to use.

Note: Dental suction handpiece (HV suction and Saliva ejector) and vacuum pump in the IFU are not individually CE certified by DNV.

1.2.3 Instrument table - swing type and hanging type

The nameplate is posted on the hydro-group. It contains the basic parameters of the unit.

Model: S90

Power Consumption: 900 VA

Safe working load of patient chair: 150 kg

Intermittent operation with intermittent loading

Intermittent operation with intermittent loading

Manufacturer: SINGER Medical Equipment Co., Ltd

Address: Building 2, No. 1 Chuangxin Yi Road, Tangjiawan

City: Dongguan City, Guangdong Province, P.R. China

Representative: SINGER International Holding Corp. GmbH (Europe)

Address: 10, 2057 Hamburg, Germany

CE Marking: MD CE 2460

(01)06974276290011(11)XXXXXX(21)XXXXXXX

Technical parameters

Technical parameters of the unit are as follows:

Type of protection against electric shock	I	Main fuse tube	T10AL250V
Degree of protection against electric shock	Type B	Rated illuminance of operating light	8000-30000Lx
Mode of operation	Intermittent operation	LED operating light	AC12V 8VA
Degree of water-proof	Non-water-tight	Power consumption	900VA
Rated voltage	AC230 V	Frequency	50/60Hz
Working pressure limits	(0.7±0.1) MPa	Air flow rate limit	Greater than 80 NI/min
Humidity limit	Dew point not greater than -20°C at atmospheric pressure	Particulate contamination limit	Not greater than 100 particles per cubic meter for 1 µm to 5 µm particle size
Contamination limit	Max. 0,5 mg/m ³	Air filter net	40 µm
Water pressure limits	200 kPa~600 kPa	Water flow rate limit	Greater than 5 l/min
Water hardness limit	Less than 2,14 mmol/l (<12 °dH)	PH limits	6,5 to 8,5
Water source requirement	Conformance to local drinking water regulations	Water filter net	75µm
Heater	AC24V 120AV	Software version	A0

Working conditions

Normal working		Shipment and Storage	
Temperature	5°C~40°C	Temperature	-40 °C~+55 °C
Relative humidity	≤80%	Relative humidity	≤93%
Atmospheric pressure	70kPa~106kPa	Atmospheric pressure	50kPa~106kPa

Suggested Shelf Life: 10 years from the manufacture date

Working range of the unit (Max length: 2100mm, Max width:1780mm)

Dimensions: Length:190cm, Width: 85cm, Height: 180cm;

Maximum seat height: ≤420mm;

Extension range of the headrest: > 200mm;

Maximum table lifting distance 480 mm;

② maximum seat height: ≥820mm;

④ inclination of the backrest: 105°~185°;

⑥ maximum caster angle of the seat: ≥10°;


⑧ rotation range of the spittoon: ≥180°.

Chapter 2 Control Panels, Display and Buttons of S90/S60/S30/U200 Dental

Unit

2.1 Operation Panel and Keys Description


2.1.1 Illustration of Doctor's and Assistant's Control Panel


S90/S60/S30/U200 Swing tubing
dentist element control panel


90 Swing tubing dentist element
control panel

S90/S60/S30/U200 assistant
control panel

S60/S30/U200 Hanging tubing dentist element control panel

S90 Hanging tubing dentist element control panel

2.1.2 Description of Buttons

Buttons	Description	Buttons	Description
	Increase the current value		Decrease the current value
	Rinsing position		Setup 1
	Emergency position		Setup 2
	Programmed position A		Timer
	Programmed position B		Spittoon rinsing
	Programmed position C		Cup filling (cold water)
	Zero position		Distilled water
	Seat up		Mouthwash heating
	Backrest up		Operating light(ON/OFF)
	Seat down		Optical fiber (ON/OFF)
	Backrest down		X-ray image viewer (ON/OFF)


Parameter Setup and Display

Turning on the machine

Check and make sure the unit is properly installed, then plug in and turn on the power supply. Then the unit conducts self-checking, after which the main menu appears as follows:

Self-check completion screen(for S30/S60/U200)

Self-check completion screen(for S90)

Main menu setup

Main screen setting(for S30/S60/U200)

Long press the set 1  for 4s. After hearing two prompt tones from the buzzer, it enters setting state. The corresponding icon will flash at setting mode. The cup water icon will flash at first. The screen is as below:



By set 1  at the screen, move the cursor to select an icon. The definitions of the icons from up to down, from left to right are as below respectively:

Cup water adding time setting , clock setting , spittoon flushing time setting , chair-spittoon linkage  at mouth rinsing position selection, buzzer on/off  selection, successive circulation.

After the setting, long press the set 1  for 4s. After hearing two prompt tones from the buzzer, it will exit the setting state.

Main screen setting(for S90)

In the main menu, three dentists A, B and C can be selected through button "setup 1" . See as follows:

Or press button "setup 2"  to enter the parameter setup menu as follows:

Press "+" or "-" to select between different icons, and the corresponding parameter of the selected icon will appear on the top right corner.


From left to right, the icon parameters are:

Mouthwash time setup -> Rinsing time setup -> Mouthwash &

Rinsing linkage -> Rinsing & Rinsing Position linkage ->

Clock setup -> Language (Chinese/English) -> Buzzer (ON/OFF)

Factory default: long press the button "setup 2"  for 5 seconds, after an indication of beep, the unit is restored to factory default status.

Mouthwash time setup

Cup water adding time setting (for S30/S60/U200)

Under the setting mode, enter cup water setting screen as below:

Press set 2 , start adding mouth rinsing water. The icon  will flash and start counting simultaneously. The clock display icon will switch


ing automatically.

After the mouth rinsing water reaches the ideal amount, press set 1
 once more to stop adding mouth rinsing water. The icon  will
 appear and the program stop counting simultaneously. The clock
 icon will stop timing automatically.


After the setting, if no setting for other options required, long press set
 for 4s. After hearing two prompt tones from the buzzer, it will exit the
 setting state.

2.2 Cup water adding time setting (for S90)


When the "mouthwash" icon is selected, the menu displays as follows:

Press "setup 1"  to enter the sub-menu of "mouthwash time":

Press "setup 1"  again to fill the cup, and the process starts
 timed.

When the cup filling is completed, press "setup 1"  to stop the water, at the same time, the clock stops
 the mouthwash time displays. Press "setup 2"  to go back to the parameter setup menu.

2.3 Clock setup

2.3.1 Clock setting (for S30/S60/U200)


Press the set 1 , enter clock setting screen as below:

The whole clock icon will flash at the beginning. Press set2 , it
 enter setting state and will shift between hour and minute.
 corresponding icon of hour and minute will flash.


Press the ascending  and descending key  at the operation chair to increase or decrease the value of
 minute that cursor locates. The adjustable range of hour is 24 and 60 for minute. The incremental quantity
 long pressing for rapidly increasing is supported.

2.3.2 Clock setting (for S90)

Select the "clock" icon and the menu displays as follows:

Press button "setup 1"  to enter the following clock setting menu:

Press button "setup 1"  to move the cursor between HH:MM (hour
 minute), and press "+", "-" to set the hour and minute.

When finish the clock setting, press button "setup 2"  to go back to the parameter setup menu.

2.4 Rinsing time setup

2.4.1 Spittoon flushing time setting (for S30/S60/U200)

Press the set 1 , enter spittoon flushing time setting screen as

Press set2 , start spittoon flushing. The icon  will flash and start counting simultaneously. The clock display icon will switch to timing automatically.


After the spittoon flushing reaches the ideal duration, press Set 2  to stop spittoon flushing. The icon  will not appear and the program stop counting simultaneously. The clock display icon will stop timing automatically.


After the setting, if no setting for other options required, long press set1  for 4s. After hearing two prompt tones from the buzzer, it will exit the setting state.

12.5.2 Spittoon flushing time setting (for S90)


Select the "rinsing" icon  and the following menu displays:

Press the "setup 1" button  to enter the rinsing time setup menu:

Press the "setup 1" button  again and the process starts being

When the rinsing is completed, press the "setup 1"  to stop rinsing, at the same time, the rinsing time displays. Press "setup 2"  to go back to the parameter setup menu.

12.5.3 Rinsing and rinsing position linkage

12.5.3.1 Rinsing and rinsing position linkage (for S30/S60/U200)

Press the set 1 , enter chair-spittoon linkage setting screen as


Press set 2 , switch between "On" and "Off". Display the icon  when it is on and only display the icon **SP** when it is off.

After the setting, if no setting for other options required, long press set 1  for 4s. After hearing two prompt tones from the buzzer, it will exit the setting state.


12.5.3.2 Rinsing and rinsing position linkage (for S90)

Select the "rinsing & rinsing position linkage" icon and the following menu displays:

Press button "setup 1"  to switch on or off the linkage:


Press button "setup 2"  to go back to the parameter setup menu.

12.5.4 Buzzer setup

12.5.4.1 Buzzer on/off selection (for S30/S60/U200)

Press set 1 , enter buzzer setting screen as below:

Press set 2 , switch between "On" and "Off". Display the icon  when it is on and display the icon  when it is off.

After the setting, if no setting for other options required, long press set 2  for 4s.



After hearing two prompt tones from the buzzer, it will exit the setting state.

12.1.2 Buzzer on/off selection (for S90)

Select the "buzzer" icon and the menu displays as follows:

Press button "setup 1"  to switch between on and off:

Press button "setup 2"  to go back to the parameter setup menu or it will go back automatically after 30 seconds.

12.1.3 Timer setup

12.1.3.1 Timer function (for S30/S60/U200)

When the main screen is not under setting operation, long press set 2  for 4s. After hearing two prompt tones from the buzzer, it will enter timing program setting state. The clock icon will be switched to timing mode. The icon  will appear.


Press the ascending  and descending key  at the operation chair to increase or decrease the first two numerical values of the clock icon. The adjustable range is 30 min. The incremental quantity is 1 minute. Long pressing for rapidly increasing is supported. Press set 2  after the setting, it will start timing. The first two digits of the clock icon will show minute while the last two digits will show second and start counting down. After reaching the time alarm, the buzzer will make 6 sounds. After the buzzer sound stopped, the icon  will disappear and the timing completes.

During the setting process, long press the set 2  for 4s. After hearing two prompt tones from the buzzer, it will exit the setting state and go back to the main screen.

12.1.3.2 Timer function (for S90)

Press the "timer" button  on the doctor's control panel to enter the following timing menu and press the same button each time to choose between 6 preset time selections - 1min, 5min, 10min, 15min, 20min and 30min, or press "+" or "-" to adjust the time. Press button "setup 1" to start timing. When the time is up, the buzzer buzzes and the timing is completed. Or press button "setup 1" to end timing.

Press button "setup 2"  to go back the parameter setup menu or it will go back automatically after 30 seconds.


12.1.4 Dental chair protection

12.1.4.1 Seat protection function (for S30/S60/U200)

When surgery chair descending, the bottom of the chair will automatically stop as for resistance encountering, and the operating chair will rise slightly. The icon  will appear.

12.1.4.2 Seat protection function (for S90)

When the chair is obstructed while moving down, it will move up automatically for 2 seconds and there will be a "C" flickering on the main menu as shown below:

22.10 Backrest protection

22.10.1 Backrest protection function (for S30/S60/U200)

When surgery chair or backrest descending, the backrest of the chair will automatically stop as for resistance encountering, and the operating chair will rise up. The icon **E02** will appear.

22.10.2 Backrest protection function (for S90)

While the backrest or seat is moving down, if the backrest is obstructed, the chair will stop moving automatically and there will be a "B" flickering on the main menu as shown below:

22.11 Spittoon protection

22.11.1 S30/S60/U200 Spittoon protection function

When the spittoon moved to the front of the operating chair, the movement of the chair was locked and all the movement buttons are invalid. Display icon **E03**.

22.11.2 S90 Spittoon protection function


When the spittoon is at the protection position, the dental chair can not move, and there will be an "S" flickering on the main menu as shown below:

22.12 Communication anomalies (for S30/S60/U200)

At communication anomalies between main operating panel and treatment table control panel, the icon **E04** will appear.

22.13 Mouthwash & rinsing linkage (for S90)

Select the "mouthwash & rinsing linkage" icon, the menu displays as follows:

Press button "setup 1"  to switch on or off the linkage:

Press button "setup 2"  to go back to the parameter setup menu.

22.14 Language setup (for S90)







Select the "language" icon as shown below:

Press button "setup 1"  to choose between Chinese and English:

Then press button "setup 2"  to go back to the parameter setup menu.

22.15 X-brightness setup (for S90)

Press the "X-ray image viewer" button on the doctor's control panel, the viewer is on. At the same time, the following menu appears and the brightness can be adjusted through button "+" or "-". There are 8 levels of brightness to choose.


Press button "setup 2"  to go back to the parameter setup menu or it will go back automatically after 30 seconds.

12.16 Optical fiber brightness setup (for S90)

When an instrument with optical fiber is in the standby mode, press the "optical fiber" button on the doctor's control panel to enter the following brightness setup menu. There are 16 levels of brightness to choose through buttons "+" or "-".

Press button "setup 2"  to go back the parameter setup menu or it will go back automatically after 30 seconds.


12.17 Turbine Setup and Display (for S90)

When the turbine handpiece is in working position as shown below, by using the button "setup 1" , its rotating speed can be roughly selected among 5 levels: 1%, 25%, 50%, 75%, 100%.

The drilling speed can be accurately adjusted 1% by 1% through button "+" and "-". Button "Setup 2"  is used to choose between adjustable rotating speed and locked rotating speed.

The sign  indicates that the rotating speed is adjustable:


The sign  indicates that the rotating speed is locked at the value displayed:

The sign  indicates that the optical fiber is on. Press the button  to turn its brightness.

16 levels of brightness are allowed to select by pressing button "+"

selected power of the scaler. The sign  indicates that the scaler is working at the displayed power, while the sign  indicates that the selected power is the maximum power of the scaler.

12. Buttons and Functions of Doctor's and Assistant's Control Panel

12.1 Seat up/down

- Press the button  on the control panel, the seat would move up.
- Press the button  on the control panel, the seat would move down.

12.2 Backrest up/down

- Press the button  on the control panel, the backrest would move up.
- Press the button  on the control panel, the backrest would move down.

12.3 Programmed positions

Programmed position A: operate the dental chair to the desired position and long press button "A"  for 5 seconds, then the current chair position is automatically saved as a programmed position after an indication of water signal. Short press button "A" , the dental chair would automatically move to the saved position A. At the same time, the operation light turns on.

Programmed position B: the same with programmed position A.

Programmed position C: the same with programmed position A.

12.4 Rinsing position

Setting of rinsing position : operate the dental chair to the needed rinsing position and press button  for 5 seconds, then the current chair is automatically saved as a memory rinsing position after an indication of water signal.

Operating of rinsing position: short press the "rinsing position" button, then the dental chair automatically moves to the saved memory rinsing position at which water flows to the spittoon and fills the cup and the operation light goes out automatically. After that, press the button again to move the chair to the last working

When the operation light turns on automatically.

Emergency position

Setting of emergency position  : operate the chair to the needed emergency position, then long press the "emergency position" button  for 5 seconds. After the indication of an audible signal, the current position is saved as the emergency position.

Operating: short press the "emergency position" button , then the chair moves to the saved emergency position and at the same time, the operation light turns on automatically.

Zero position (exit position)

Setting: operates the chair to the needed zero position, then long press the "zero position" button  for 5 seconds. After an indication of an audible signal, the current position is saved as the zero position.

Operating: short press the "zero position" button , then the chair moves to the saved zero position and operation light goes out.

Spittoon rinsing

Operating: short press the "spittoon rinsing" button  to rinse the spittoon once for the saved time, or long press the button to stop.

Cup filling

Short press the "cup filling" button , the water comes out and flows according to the saved time; then long press the button to stop the water.

Mouthwash heating

Heating key operation (for S30/S60/U200)

Short press mouth rinsing heating key  to start heating. During the operation process, the heating display shows Heating screen icon . When the water temperature reaches the standard value, the icon  of  will be constantly; and when the water temperature not to the standard value, the icon will flash at 1s interval; when not heating, the screen will display the icon .

Heating key operation (for S90)

Short press the "mouthwash heating" button , the indicator light flickers. When the water is heated to the standard temperature, the indicator light is on, which means the heating is accomplished. During heating, the indicator light works in a cyclical way between flickering and lighting.

Dental chair key locking operation (for S30/S60/U200)

During the cleaning and disinfection process by the operators, make the dental chair system keys locked. The operation is as below. When no operation undergoing, long press mouth rinsing heating key  for 4s. After hearing two prompt tones from the buzzer, it will show the icon , and the system will enter the key lock state. In this state, all keys will be invalid when pressed (including the main operating touch key, assistant touch key, electric induction, pull-rod switch, etc.) and the display backlight will be off. Long press mouth rinsing heating key  for 4s again. After hearing two prompt tones from the buzzer, the icon  will disappear. The system will be unlocked and restore to the state before locking.

3-way image viewer(Restore factory settings S90 No icon, no action)

Film viewing illuminator key operation: Short press the film viewing illuminator key . The film viewing illuminator will turn on. It will turn off 3min. later or when it is pressed again.

Long press the film viewing illuminator switch  for 4s. After hearing a buzzer beep, the icon  will flash. The dental chair will ascend to the limit position and then descend to the limit position at the same time. Reset 5 programming memories automatically. After hearing 2 prompt tones from the buzzer, all settings will restore to factory defaults, and the icon  will stop flashing.

12.1.2 Optical fiber

When the instrument is in the standby mode, short press the "optical fiber" button  to turn on/off the optical fiber.

12.1.3 Operation light

Short press the "operation light" button , the light is on. Press the button again to turn it off.

12.1.4 Dry and wet grinding switching of the instrument

When picking up an instrument, pressing the left button on foot control and switch on/off working water. When switching on working water, the setting button signal light is on. (S30/S60/S90/U200 has 'dental instrument with water' icon on the display). When switching off working water, signal light will be off accordingly. (On S30/S60/S90/U200, 'dental instrument with water' icon is disappeared).

12.1.5 System standby mode(for S30/S60/U200)

If the dental chair does not make any action and operation for two minutes, the system will enter standby mode, and the display back light will automatically shut down.

When the dental chair takes any action, it will exit standby mode automatically when the pedal pull rod is being pulled and the instrument being taken up. The display back light will be on.

Note: System standby mode is propitious to protect the display screen, prolong its service life and reduce the power consumption of the control system.

Chapter 3 Control Panels and Buttons of U500/U100 Dental Unit

3.1 Description of Control Panels and Buttons

3.1.1 Illustration of doctor's and assistant's control panel

U100 doctor control panel (swing version)


U100 assistant control panel

U500 main control panel

U500 assistant control panel

U100 doctor control panel (hanging version)

12.1.2 electronic system control panel (touch button), U100/U500 pneumatic system control panel (press button)

3.1.2 Description of buttons

Buttons	Description	Buttons	Description
	Seat up		Rinsing position
	Seat down		Settings

	Backrest up		Spittoon rinsing
	Backrest down		Cup filling
	Programmed position A		Mouthwash heating
	Programmed position B		Operation light (ON/OFF)
	Zero position		Optical fiber (ON/OFF)
	Emergency position		X-ray image viewer (ON/OFF)

Buttons and Functions of Doctor's and Assistant's Control Panel

12.1 Seat up/down

Press the button  on the control panel, the seat would move up.

Press the button  on the control panel, the seat would move down.

12.2 Backrest up/down

Press the button  on the control panel, the backrest would move up.

Press the button  on the control panel, the backrest would move down.

12.3 Programmed positions

Setting of programmed position A  : operate the dental chair to the desired position and press button "A"  for 5 seconds. The current chair position will be automatically saved as a programmed position after an indication of audible signal.

Operating of programmed position A: short press button "A" , the dental chair would automatically move to the saved position A. At the same time, the operation light turns on.

Setting and operating of programmed position B : the same with programmed position A.

12.4 Rinsing position

Setting of rinsing position  : operate the dental chair to the desired rinsing position and press button  for 5 seconds. The current chair will be automatically saved as a memory rinsing position after an indication of audible signal.

Setting of linkage between rinsing position and spittoon rinsing: press the "settings" button  for 5 seconds, the indicator light is on. Press the "rinsing position" button , then the linkage between the chair movement and spittoon rinsing is on after an audible indication and the indicator light goes out.

Operating of rinsing position: short press the "rinsing position" button, then the dental chair will automatically move to the saved memory rinsing position at which water rinses to the spittoon and fills the cup and the operation light goes out automatically. After that, press the button again to return the chair to the last working position and the operation light turns on automatically.

12.5 Emergency position

Setting of emergency position  : operate the chair to the desired emergency position, then long press the "emergency position" button  for 5 seconds. After the indication of an audible signal, the current position will be saved as the emergency position.

Operating: short press the "emergency position" button , then the chair will move to the saved emergency position and at the same time, the operation light turns on automatically.

12.6 Zero position (exit position)

Setting: operates the chair to the desired zero position, then long press the "zero position" button  for 5 seconds. After an indication of an audible signal, the current position is saved as the zero position.

Operating: short press the "zero position" button , the chair will move to the saved zero position and the operation light goes out.

12.7 Spittoon rinsing

Time setting: long press the "settings" button  for 5 seconds, the indicator light will be on. Then short press the "spittoon rinsing" button , when the desired time is reached, press again the "spittoon rinsing" button to stop. The indicator will go out and the setting is accomplished.

Operating: short press the "spittoon rinsing" button  to rinse the spittoon once for the saved time, or press again the button to stop.

12.2.1 Cup filling

Setting: long press the "settings" button  for 5 seconds, the indicator light is on. Then short press the "cup filling" button  to fill up the cup with water. When the desired water volume is reached, press again the button . The indicator light will go out and the setting is accomplished.

Operating: short press the "cup filling" button , the water will come out with desired volume. Press again the button  to stop the water.

12.2.2 Mouthwash heating

Short press the "mouthwash heating" button , the indicator light will flicker. When the water is heated to the preset temperature, the indicator light is on, which means the heating is accomplished. During operating, the indicator light will work in a cyclical way between flickering and lighting.

12.2.3 X-ray image viewer

Operating: short press the "X-ray image viewer" button , the viewer will be on. It will be turned off automatically after 3 minutes or be turned off by pressing again the button.

Factory default: long press the button  for 5 seconds. After an indication of beep, all memory positions are restored to factory default.

12.2.4 Optical fiber

Operating: when the instrument is in the standby mode, short press the "optical fiber" button , the light is on. Fiber optical signal light is on. (not available on U500 pneumatic model and U100). Press the button again to turn off. The Fiber optical signal light is off. (not available on U500 pneumatic model and U100).

Operation light setting function on foot control: long pressing fiber optical light button for five seconds till signal sound "BEE" is heard, the function of lower button on the left side of foot control is changing from last position/programmed B position function to light on/off and last position function. The button functions on main control panel are unchanged. Long pressing again for five seconds, two functions can be exchanged. The users can set the operation light function button freely as their own preference.

12.2.5 Operation light

Operating: short press the "operation light" button , the light will be on. Press the button again to turn it off.

Long press the "operation light" button  for 5 seconds. After an indication of beep, the indicating function of the signal will be off. Another long press can turn it on again.

Chapter 4 Operation Instruction

WARNING! Before any operation, make sure the unit is correctly installed and adjusted according to the Manual of Installation and Adjustment.

4.1 Main Switch

Press the main switch to "I", the green power indicator on the switch will illuminate. By now, the unit is power-on.

4.2 Operation light

4.2.1 LED light


The angle of light handle can be adjusted according to the user's need.

When the unit is turned on, the illuminance of the LED light in default is Level 3.

The sensing range of the LED light sensor is 50mm-120mm. Put your hand close to the sensor and let go after 1.5 second, the light can be turned on/off.

When the light is on, put your hand close to the sensor. The illuminance will go from low to high (level 1 to level 5) (V2 LED light has four levels), and then high to low. Take your hand away when it reaches the desired level.


The reflector type V2 LED operation light has blue light stopping function. When doctor use light sensitive resin for tooth filling, slide switch can turn the light to yellow to prevent curing of sensitive resin. During the normal dental operation, please open up the slide switch.

Assembly and disassembly of the light handle


Assembly: Grab the light with your right hand->grab the left handle with your left hand and press the handle button with your thumb->Put the handle in the slot on the rotating shaft and release the handle button. Same method for right handle.

Disassembly: Grab the light with your right hand->grab the left handle with your left hand and press the handle button with your thumb->Pull out the handle powerfully
Same method for right handle.

Foot pedal

Electronic Multi-functional lever type foot control

Chair movement control: slide the chair movement control, the up and down movements of seat and backrest can be performed.

Left button:

① When no dental instrument is picked up or pulled out and short press the button, the dental chair will move to the rinsing position. Short press it again, the chair will return to the last working position. Long press the button (for 2 seconds) and after an indication of audible signal, the chair will move to programmed position A.

② When a dental instrument is picked up or pulled out, it is an activating/deactivating switch of water spray. When the water spray is activated, the setting button indicator is on (S30/S60/S90/U200 has 'dental instrument with water' icon on the display). When the water spray is deactivated, the setting button indicator is off (S30/S60/S90/U200 has 'dental instrument with water' icon on the display).

Right button:


① When no dental instrument is picked up or pulled out and short press the button, the dental chair will move to the reset position. Long press it (for 2 seconds) and after an indication of audible signal, the chair will move to programmed position B.

② When a dental instrument is picked up or pulled out, it is an air blast button. The air keeps delivering when the button is kept pressing. The air flow will stop when the button is released.

Driving rod: When a dental instrument is taken out or pulled out, the instruments are ready for work. At this time, the foot pedal is available for controlling the speed or power of instruments. From left to right, the instruments will be sped up gradually.

Handle: Lift the handle with your foot to move the foot-control.

Electronic Multi-functional pedal type foot control

Chair movement control: slide the chair movement control, the up and down movements of seat and backrest can be performed.

Left button:

① When no dental instrument is picked up or pulled out and short press the button, the dental chair will move to the rinsing position. Short press it again, the chair will return to the last working position. Long press the button (for 2 seconds) and after an indication of audible signal, the chair will move to programmed position A.

② When a dental instrument is picked up or pulled out, it is an activating/deactivating switch of water spray. When the water spray is activated, the setting button indicator is on (S30/S60/S90/U200 has 'dental instrument with water' icon on the display). When the water spray is deactivated, the setting button indicator is off (S30/S60/S90/U200 has not 'dental instrument with water' icon on the display).

Right button:


① When no dental instrument is picked up or pulled out and short press the button, the dental chair will move to the reset position. Long press it (for 2 seconds) and after an indication of audible signal, the chair will move to programmed position B.

② When a dental instrument is picked up or pulled out, it is an air blast button. The air keeps delivering when the button is kept pressing. The air flow will stop when the button is released.

Foot control pedal: When a dental instrument is taken out or pulled out, the instruments are ready for work. At the time, the foot pedal is available for controlling the speed or power of instruments. From left to right, the instruments will be sped up gradually.

Handle: Lift the handle with your foot to move the foot-control.

Pneumatic Multi-functional lever type foot control

Chair movement control: slide the chair movement control, the up and down movements of seat and backrest can be performed.

Up key operation: Short press the Up key, it will go to mouth rinsing programming position. Press it once more, it will go to mouth rinsing reset. Long press the key (for about 2s), after hearing a "du" from the buzzer, it goes to A programming position.

Down key operation: Short press the key, it will go to chair bit reset function. Long press the key (for about 2s), after hearing a "du" from the buzzer, it goes to B programming position.


Left key operation: Short press the key, it will shift from dry grinding to wet grinding.

Right key operation: When press the key, the dynamic equipment will be as single jet mode. When loosen the key, single jet will stop.

Driving rod: Move from left to right to control the working speed or power increasing of the dynamic equipment.

Lifting yoke operation: It is allowed to hook the lifting yoke by feet to move the pedal.

4.2.4 Pneumatic Multi-functional pedal type foot control

Chair movement control: slide the chair movement control, the up and down movements of seat and backrest can be performed.

Programming position left switch pedal operation: Short press the programming position left switch pedal, and it will go to mouth rinsing position. Press the key once more, it will be mouth rinsing reset. Long press the key (for about 2s), after hearing a "du" from the buzzer, it goes to A programming position.

Programming position right switch pedal operation: Short press the programming position right switch pedal and it will go to chair bit reset function. Long press the key (for about 2s), after hearing a "du" from the buzzer, it goes to B programming position.

Left key operation: Short press the key, it will shift from dry grinding to wet grinding.

Right key operation: When press the key, the dynamic equipment will be as single jet mode. When loosen the key, single jet will stop.

Foot control pedal: Move from up to down to control the working speed or power increasing of the dynamic equipment.

Lifting yoke operation: It is allowed to hook the lifting yoke by feet to move the pedal.

4.2.5 Logic priority of S90 control panel (for chair movement control only)

Logic priority of panels: foot pedal, doctor's control panel, assistant's control panel.

Eg: if you firstly press the "seat down" button on the assistant's panel, then slide the "seat up" button on the foot pedal, the order will be taken from the foot pedal to move the chair upward.

4.2.6 Syringe

Press the water button to spray water. Press the air button to spray air.

Press both of the buttons to spray mist.

Press the outer ring near the spout to release the spout and then pull the spout out to have it disinfected.

4.2.7 Dental instrument

Operating: pick up or pull out a dental instrument and operate the foot pedal, the dental instrument will start to work.

The input air pressure of the dental instrument can be regulated through the air regulating valves under the instrument table. Screw clockwise to turn down the air pressure, or vice versa.

4.2.8 Water volume regulating

Water volume can be regulated through the water regulating valves under the instrument table. Screw clockwise to turn down the water volume, or vice versa. One valve applies to only one instrument.

4.2.9 Spray volume regulating

Spray volume can be regulated through the spray regulating valve under the instrument table. Screw clockwise to turn down the spray volume, or vice versa. The valve applies to all instruments.

NOTE:

Dental instrument is purchased and used at the discretion of professional doctors. It should be used in accordance with the requirements of the CE standard and obtain a CE certificate. The recommend dental instrument see following table:

Name	Model/specification	Manufacturer*
Dental Turbine handpiece	THALIA	TEKNE DENTAL SRL

4.5 HV Suction

The suction system starts working once the HV suction tube is pulled out from the tube holder. When the tube is put back in the holder, the suction system won't stop until 2 seconds later. This delay assures the hygiene of the suction tube. (not available on U500 pneumatic model and U100).

Under static (i.e. no flow) conditions, the maximum vacuum pressure provided to the dental unit at the HV Suction connection port is 25kPa.

4.7 Saliva Ejector


The ejector system starts working once the saliva ejector tube is pulled out from the tube holder. When the tube is put back in the holder, the ejector system will stop.

Under static (i.e. no flow) conditions, the maximum vacuum pressure provided to the dental unit at the Saliva Ejector connection port is 15kPa.

4.8 Height Adjustment of Balance Arm

S90: grab the handle with your right hand, and put your thumb on the air brake sensor; move the handle to a needed position and take off your thumb from the air brake sensor. The arm is thus fixed.

S60, S30, U200, U100: grab the handle with your right hand, and press the air brake with your thumb; move the handle to a needed position and release the air brake. The arm is thus fixed.

U500: grab the handle with your right hand, and press the air brake with your thumb; move the handle to a needed position and release the air brake. The arm is thus fixed.

4.9 Instrument Tray

The maximum load of the instrument tray is 3kg, which should never be exceeded.

The instrument tray of swing type can rotate 180 degrees.

4.10 Drawbars of Swing Type Instrument Table

Take up an instrument and pull it backwards, after you hear a click, the drawbar is locked and the instrument is ready for working.

When the operation is finished, pull the instrument backwards again to unlock the bar, and after you hear another click, the instrument can be put back.

4.11 Headrest

Rotate the clamping knob counterclockwise for one or two circles, then the headrest can be removed at will.

Adjust the headrest to the desired position (to the lower end of the backrest for children treatment), then tighten it clockwise.

4.12 Right Armrest (optional)

During treatment, grab the end of the armrest and lift it up. An audible clicking sound indicates the armrest is ready to be locked. Release the armrest to finally lock it.

After treatment, grab the end of the armrest and lift it to the maximum height, the armrest will be free to fall to minimum height.

4.3.3 Height Adjustment of Assistant's Instrument Board

Hold the assistant's instrument board with both hands and lift it upwards until the desired position is reached. It can be positioned vertically at four levels.

To lower the instrument board, lift it fully upwards until the lock is released and then move downwards to the lowest position. After that, follow step 1 to readjust the height.

4.3.4 Distilled Water System

The distilled water system is composed of a selecting switch and a distilled water tank at the front of the hydro-group. Distilled water in the tank or sanitized tap water is sent to all hydrodynamic instruments on the dentist's and assistant's instrument table, so the tank must be filled with medical distilled water in time.

Water filling:

- (1) Switch the door lever to open the front door.
- (2) Switch the air pressure lever to "Decompress" and wait till the air in the tank is completely released.
- (3) Take out the water tank and screw it off from the lid anticlockwise. Fill up in the tank with distilled water.
- (4) Screw on the tank clockwise and put it back. Switch the air pressure lever to "Plus press".
- (5) Switch the door lever to close the door.

Decompress ⇌ Plus press

Tap water ⇌ Distilled water

NOTE:

The bottle must only be opened when the air in it is completely released.

To screw off the tank, hold the lid with left hand and screw the tank with right hand. Do not screw the lid, otherwise the air and water pipes connected to it will be damaged.

We suggest using distilled water during treatment so as to avoid the negative effect on the unit's pipeline and the service life of the instruments caused by impurities in the tap water.

4.3.5 Disinfection of the Pipeline


Keep the water tank empty and switch the air pressure lever to "Plus press" to let out the waste and water from the pipeline.

Switch the air pressure lever to "Decompress". Remove the white distilled water tank and install the red tank with disinfectant.

Switch the air pressure lever to "Plus press" and operate the dental instruments one by one until disinfectant comes out and lasts for seconds. Keep the unit in non-working status for one or two hours.

Switch the air pressure lever to "Decompress", remove the red disinfectant tank and reinstall the white distilled water tank. Switch the air pressure lever again and clean the internal water line the same way. Thus the disinfection is completed.

It is recommended to use five thousandths of hydrogen peroxide as a disinfectant, disinfected once a day. The pipeline has to be replaced when it was disinfected for 1,900 time. If there is any problem with the pipeline about appearance, please contact Siger to replace on time.

4.3.6 Safety Devices

The bottom of the dental chair is provided with a safety device which stops the chair from moving down in the presence of obstacles. The chair will move up to allow the user to remove the obstacles (for S90, the letter "C" will appear and flicker on the lower part of the display. (S60/S30/U200 has E01 icon on display).

When the backrest is moving down, it will stop automatically in the presence of obstacles (for S90, the letter "B" will appear and flicker on the lower part of the display. (S60/S30/U200 has E02 icon on display).

When the spittoon is in front of the chair, the chair is locked and all movement buttons are unavailable (for S90, the letter "S" will appear and flicker on the lower part of the display. (S60/S30/U200 has E03 icon on display).

Put Syringe samples into heat-sealing sterilization pouches. The Syringe were placed in trays, then transferred to sterilizer cavity.

2) Moist heat sterilization:

Program: Quick B.

Sterilization temperature: 134°C.

Sterilization time: 4 minutes.

Drying time: 15 minutes.

Pressure: 2.1 bar.

The Syringe has to sterilization after each use.

The Maximum sterilization times: 250

5.2 Maintenance of the Unit

NOTE: proper maintenance can extend the service life of the unit!

Check the electric supply every month. Water pressure and air pressure must be in accordance with the unit requirement. Keep the hoses smooth and unblocked.

Check the chair's movement and control panel.

Clean the spittoon after use everyday. Wash the filter, pull out the filter mesh and reinstall them after cleaning to avoid air leakage.

As for the cleaning and lubrication of dental instrument, please refer to the instruction manual of dental instrument.

Water filter: in order to ensure the proper function of the unit, a water filter is installed at the water inlet to filter the impurities in water. After a long-time operation, the filter collector may be filled with contaminant and the water flow may be obstructed. Therefore, the filter collector must be cleaned or replaced timely.

Generally, if any of the following cases happens, the filter collector should be cleaned or replaced:

(1). the filter collector has been used for more than 1 year;

(2). the filter collector is contaminated;

(3). the water coming out becomes turbid.

How to replace the filter collector: open the hydro-group and remove the water filter. Then remove the filter anticlockwise, take off the seal ring and the collector successively. After cleaning or replacing, reinstall all the parts. Please take care of the seal ring to avoid leakage.

Chapter 6 Dental Unit installation, commissioning and maintenance

6.1 Installation Conditions

Caution: The equipment must be installed by our professionals or any personnel authorized by us.

This equipment shall be installed at a clean, dry, cool and well-ventilated location, and sufficient space shall be reserved around the unit to facilitate access during routine maintenance.

The equipment shall not be installed or operated at any location:

A. exposed to humidity, heavy dust, poor ventilation or direct sun; or

B. subject to substantial temperature or humidity change, such as place near air conditioner or heater.


Caution: Improper installation may cause damage to the equipment!

6.2 Requirements for Installation Field


The equipment shall be located based on general layout of the treatment room, lighting condition and operation convenience, and should be installed on flat, horizontal and solid ground at a clean, dry, cool and well-ventilated location.

Standard Layout for S Series Comprehensive Treatment Platform:



The layout is for single-person operation and excludes any auxiliary equipment.. For double-person operation, accommodate additional auxiliary equipment or larger based on actual requirement.

Operating Chair Sizes Diagram

Ground Case Hole Arrangement and Case Sizes Diagram:


The Drain Hose shall be $\Phi 50$ mm PVC hose, and water supply hose and air supply hose are both $\Phi 15$ mm furnished with HP aluminum-plastic complex tube port for connection of ball valve. The power supply shall be 230V and grounded properly, with wire diameter $\geq 1.5 \text{ mm}^2$. The underground hose line shall lie at the same elevation as that of the floor ground as possible, or not above 100 mm higher than the floor. Connection between PC and endoscope (if any) shall be laid through the underground hose line. The VGA cable and AV cable shall be connected to the computer.

Supply voltage 230 VAC
 Air pressure limits $(7 \pm 1) \text{ kPa}$
 Water pressure limits 200~600 kPa
 Maximum wastewater flow rate: 6 L/min

Power frequency 50/60Hz
 Air flow rate limit Greater than 80 NI/min
 water flow rate limit Greater than 5 l/min
 Minimum gradient of the wastewater lines: 5°

Handling and Unpacking of Equipment

The equipment shall be handled to the installation field in accordance with the

Instructions as stated on the label.


The product packaging includes two cases respectively containing operating chair and side case, and the product can be unpacked by removing the locking bolts.

After unpacking, take out the packing list and check the components and accessories against the list. For any question, contact us immediately.

6.4 Equipment installation method

6.4.1 Installation of patient chair and instrument table

Installation of Patient Chair: Remove the front and bottom plate enclosures of the chair, and remove the four bolts locking the chair to the bottom plate, and mount the bottom of the chair onto the position as predetermined, and insert the headrest into the backrest in place.

6.4.2 Installation of Side Case:

- Remove the packaging foam, displace the accessory case and other accessories.
- Remove the external shield at both sides of the side case, and the three bolts connecting the bottom plate and bottom of packing board, and handle the case in balanced status to the installation position as predetermined, joint the two bottom plates of operating chair and side case together in place, mount the connecting plate and fasten it with six socket head screws.
- Fix the bottom plates of side case and operating chair onto the ground with seven 8mm expansion bolts through the fixing holes.
- Connect all plugs of wires on operating chair to the sockets on the circuit board in side case and then close the cover.

6.4.3 Installation of the instrument table

Pull the instrument rotation arm from the holder, keep it balance and move it to hydro-group position

Insert the instrument arm into the corresponding hole on the hydro-group, pay attention to the fixation orientation of the instrument arm. Install the 2 set-bolts and fasten them, then assemble the instrument table. Plug the 2 sets of wires of the instrument table to the corresponding connector on the PCB of the hydro-group, connect the water and air pipe to the corresponding connector in the hydro-group.

6.4.4 Installation of the operation light

Insert lamp rod into the hole in balancing swiveling arm, and thread wire on lamp arm into lamp column, and connect and fix operation lamp head, and then connect lamp arm wire plug to the connector on circuit board in side case, as shown below.

6.4.5 Water, air and electricity installation of ground-group

Put the G1/4" Φ 8 direct connector into the air inlet valve, connect one side of the blue Φ 8PA pipe into the direct connector, and the other side to the general air inlet of the dental unit.

Put the G1/4" Φ 8 direct connector into the air inlet valve, connect one side of the white Φ 8PA pipe into the direct connector, and the other side to the general air inlet of the dental unit.

Insert the drain pipe, suction vent nozzle to the ground hole of Φ 50PVC drain outlet

Connect the power and grounding wire well, plug all the connector on the corresponding socket on the PCB of ground-group, organize all the pipe well.

6.4.6 vacuum pump installation of ground-group


Insert the Φ 25 straight connector into the suction pipe of the vacuum pump, and connect the other side to the negative pressure suction connection pipe of the dental unit. Connect the negative pressure pump cable of the dental unit to the start port of the vacuum pump.

Under static (i.e. no flow) conditions, the maximum vacuum pressure provided by the vacuum pump to the dental unit is 25kPa, the air flowrate not less than 90 l/min.

6.5 Equipment debugging method

Main Air Pressure Adjustment:

As shown below, adjust the inlet pressure regulating valve to 550kPa and 300kPa based on the two pressure meters on operating chair (for S90).

As shown below, adjust the inlet pressure regulating valve to 550kPa and 300kPa based on the two pressure meters on operating chair (for S60/S30).

300kPa

550kPa


As shown below, adjust the inlet pressure regulating valve to 550kPa based on the pressure meters on operating chair (for U500/U200/U100).

550kPa

High-Speed Handheld sprayer air pressure adjustment: Adjust the air pressure to 250 kPa with the knob on the bottom of instrument panel.

Low-Speed Handheld sprayer air pressure adjustment: Adjust low-speed handheld sprayer pressure to 250kPa-300kPa by using the knob on the bottom of instrument panel.

Air adjustment

250-300kPa

250-300kPa


Water sampling connection point

Due to the water supply of dental unit is from domestic portable water, the user can install the input water sampling connection point near the water supply port according to actual needs. The specific method is as follows:

Connect G1/4*Φ8 barrel connector to the water sampling rotary valve of.


Connect the Φ8 three-way connection between the water input check valve and the water filter, and connect one end of the orange Φ8 polyether pipe to the three-way connector, and one end into the water sampling rotary valve.

The other end of the water sampling rotary valve is connected to an orange Φ8 polyether tube of appropriate length.


When sampling, just switch on the water sampling rotary valve, and switch off the water sampling rotary valve after sampling.

Water/Air Supply Schematic Diagram


[illegible]


Water/Air Supply Schematic Diagram


Water-Air schematic diagram of U200, U500 electrical model


Water-Air schematic diagram of U500 pneumatic model


S60/S30/U200 Electric Block Diagram


S90 Electric Block Diagram


U100 Electric Block Diagram

Electrical schematic diagram of U500 pneumatic model

Electrical schematic diagram of U500 hydraulic model

Chapter 7 Common Faults and Exclusion

7.1 Dental Chair

Fault	Cause	Action
No Chair Movement	The fuse of the power supply is damaged.	Replace the same type of fuse.
The up and down function of the seat is abnormal.	The dental instrument is damaged or the wiring is loose	Check the line and contact the engineer.
The function of the backrest is abnormal.	The dental instrument is damaged or the wiring is loose	Check the line and contact the engineer.
Movement Chaos	The chip of the CPU is damp or damaged.	Moisture absorption treatment and contact the engineer.

7.2 Operating Light

Bad condensation	The bulb is not installed on the focal spot	Adjustment
------------------	---	------------

7.3 Instrument Table

Fault	Cause	Action
-------	-------	--------

Not enough air comes from the dental instruments.	Air hoses deformed or filter blocked	Adjust the trachea or clean the filter
Not enough water comes from syringe	Throttle adjustment is too small or too large	Adjust the throttle
The suction failed.	The solenoid valve is damaged	Contact the engineer for replacement.
Hydraulic suction force is too small	O-ring is broken or water pressure is lower than required	Replace the new O-ring or adjust the water pressure
Pneumatic suction force is too small	O-ring is broken or water pressure is lower than required	Replace the new O-ring or adjust the water
None warm water supply	Boiler is damaged.	Contact the engineer for replacement.

7.4 Replacement of fuse tube

Shut down the main power supply, then open the hydro-group box, swing out the cover of the fuse tube with medium size cross screwdriver. Replace a new fuse tube, and operate in the opposite procedure as above.

Fuse Tube	Specification
Main fuse tube	T10AL250V

Chapter 8 Electromagnetic Compatibility

Notice:

- Dental Unit S90, S60, S30, U100, U200, U500 meets the requirement of standard electromagnetic compatibility of IEC60601-1-2.
- The users should install and use the unit according to the information of electromagnetic compatibility of the documents which are shipped with the unit.
- Portable and mobile radio frequency communication equipment may affect the performance of dental unit S90, S60, S30, U100, U200, U500. During the usage of these dental units, please avoid strong electromagnetic interference, such as being close to mobile phones, microwave ovens and so on.
- The detail of the guide and the manufacturer's statement, please refer to the annex.

⚠ Warning

- Dental unit S90, S60, S30, U100, U200, U500 should not be close to or stacked with other equipment during usage. If it must be used close to or stacked with other equipment, you should observe and test the unit whether it could work successfully with its configuration.
- In addition to cables sold as spare parts of the dental unit S90, S60, S30, U100, U200, U500 by the manufacturer, the use of other accessories and cables may cause radio frequency launching increasing and anti-interference reducing of dental unit S90, S60, S30, U100, U200, U500.

Annex:

Guidance and manufacturer's declaration – electromagnetic emissions		
S90, S60, S30, U100, U200, U500 Dental Unit is expected to use in the following electromagnetic environment, buyers or users should ensure that it is in the electromagnetic environments:		
Emission testing	Compliance	Electromagnetic environment – Guide
RF emission CISPR 11	Group 1	S90, S60, S30, U100, U200, U500 Dental Unit uses RF energy only for its internal function. Therefore, its RF emission is very low and is unlikely to interfere with nearby electronic equipment
RF emission CISPR 11	Class B	S90, S60, S30, U100, U200, U500 Dental Unit Suitable for use in all facilities, including direct to home

Harmonic emissions IEC61000-3-2	Compliant	Household facilities and public low-voltage power supply network
Voltage fluctuations / flicker emissions IEC61000-3-3	Compliant	

Guidance and manufacturer's declaration – electromagnetic immunity			
S90, S60, S30, U100, U200, U500 Dental Unit is expected to use in the following electromagnetic environment, buyers or users should ensure that it uses in the electromagnetic environment:			
Immunity test	IEC 60601 Test level	In line with levels	Electromagnetic environment – Guide
Electrostatic discharge (ESD) IEC61000-4-2	±6 kV Contact discharge ±8 kV Air discharge	±6 kV Contact discharge ±8 kV Air discharge	Floors should be wood, concrete or ceramic tiles, if the ground covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC61000-4-4	±2kV On the power supply line ±1kV The input / output lines	±2kV On the power supply line ±1kV The input / output lines	Net power supply should be used with a typical commercial or hospital environment quality
Electrical surge IEC61000-4-5	±1 kV Wire to wire ±2 kV Line to ground	±1 kV Wire to wire ±2 kV Line to ground	Net power supply should be used in a typical commercial or hospital environment quality
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % U_T , 0.5 cycle (U_{of-t} -Shang , >95% SAG) 40 % U_T , 5 cycles (U_{of-t} -Shang , 60% SAG) 70 % U_T , 25 cycle (U_{of-t} -Shang , 30% SAG) <5 % U_T , 5s (U_{of-t} -Shang , >95% SAG)	<5 % U_T , 0.5 cycle (U_{of-t} -Shang , >95% SAG) 40 % U_T , 5 cycles (U_{of-t} -Shang , 60% SAG) 70 % U_T , 25 cycle (U_{of-t} -Shang , 30% SAG) <5 % U_T , 5s (U_{of-t} -Shang , >95% SAG)	Mains should have quality used in a typical commercial or hospital environment. S90, S60, S30, U100, U200, U500 of dental comprehensive treatment machine user requires continuous operation during power outage it is recommended S90, S60, S30, U100, U200, U500 dental comprehensive treatment machine using an uninterruptible power supply or battery powered
Power frequency (50 / 60Hz) magnetic field IEC61000-4-8	3A/m	3A/m, 50/60Hz	Power frequency magnetic fields should have in a typical commercial or hospital environment, level of power frequency magnetic field characteristics of typical sites
Note: U_T test voltage is applied before the AC voltage.			

Guidance and manufacturer's declaration – electromagnetic immunity			
S90, S60, S30, U100, U200, U500 Dental comprehensive treatment machine is expected to use in the following electromagnetic environment, buyers or users shall ensure that its use in the electromagnetic environment:			
Immunity test	IEC 60601 Test level	In line with levels	Electromagnetic environment – Guide
Conducted RF IEC61000-4-6	3 V (Effective value) 150 kHz ~80 MHz	3V (Effective value)	Portable and mobile RF communications equipment should not be closer than the recommended separation distance S90, S60, S30, U100, U200, U500 uses any part of dental unit, including cables. Calculation formula of the distance from the transmitter frequency. Recommended separation distance $d = 1.2\sqrt{P}$

Radiation RF IEC61000-4-3	3V/m 80 MHz ~2.5 GHz	3V/m	$d = 1.2\sqrt{P}$ 80 MHz ~800 MHz $d = 2.3\sqrt{P}$ 800 MHz ~2.5 GHz Type in the: P -According to the transmitter manufacturer's maximum transmitter Rated output power, in Watts (w) as a unit; d -The recommended separation distance in meters (m) ^b . Stationary RF field strength of the transmitter through the investigation of electromagnetic field ^c to determine, in each frequency range ^d should more than meet the low level. Mark the following symbol interference may occur in the vicinity of the device. 
------------------------------	----------------------------	------	---

Note 1: in the 80MHz and 800MHz frequencies, high frequency formula.



Note 2: these guidelines may not be appropriate for all circumstances, electromagnetic propagation affected by absorption and reflection of buildings, objects and human impact.


^a Fixed transmitters, such as radio (cellular / Cordless) telephones and land mobile radio base stations, amateur radio, AM and FM radio and television broadcasting, its strength in theory, cannot be accurately predicted. For evaluating the electromagnetic environment of stationary radio transmitters, electromagnetic site survey should be considered. If measured S90, S60, S30, U100, U200, U500 dental unit places the field strength is higher than the above applicable RF compliance with level, it should be observed S90, S60, S30, U100, U200, U500 dental unit to verify that it works. If observed unusual performance, additional measures may be necessary, such as restructuring S90, S60, S30, U100, U200, U500 the direction or position of dental unit.

^b 150kHz~80MHz the entire frequency range, field strength should be less than 3V/m.

Portable and mobile RF communications equipment and the S90, S60, S30, U100, U200, U500 recommended separation distance between dental comprehensive treatment machine			
S90, S60, S30, U100, U200, U500 Dental Unit is expected to RF radiation disturbance in a controlled electromagnetic environment of use. According to the maximum rated output power of communications equipment, purchasers or users by using the following recommended maintaining portable and mobile RF communications equipment (transmitters) and S90, S60, S30, U100, U200, U500 between dental unit minimum distance to prevent electromagnetic interference.			
Rated maximum output power of transmitter W	Corresponds to the transmitter frequency isolation distance /m		
	150 kHz ~ 80 MHz $d = 1.2\sqrt{P}$	80 MHz ~ 800 MHz $d = 1.2\sqrt{P}$	800 MHz ~ 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
Transmitter maximum nominal output power that is not listed in the previous table, the recommended isolation distance d in meters (m) as a unit, available in the relevant transmitter frequency bar formula to determine, where p is the transmitter maximum nominal output power provided by the transmitter manufacturer, in Watts (w) as the unit.			
Note 1: in 80MHz and 800MHz on the frequency, high frequency formula.			
Note 2: these guidelines may not be appropriate for all situations. Electromagnetic propagation of buildings, objects and the effects of the body's absorption and reflection.			

Chapter 9 Safety Markings

Symboles/ Marks	Description	Symboles/ Marks	Description
	Warnings		Date of manufacture
	Caution		Manufacturer
	Serial number		Unique device identifier
	TYPE B APPLIED PART		Temperature limitations
	Indicates a medical device that needs to be protected from moisture.		Refer to instruction manual/ booklet
	Indicates a medical device that can be broken or damaged if not handled carefully.		Medical device
	Stacking limit by 2		WEEE Directive mark
	Temperature limit		European representative
	Atmospheric pressure limitation		CE marked

	Sterilizable Autoclave		sterilizable in a steam sterilizer (autoclave) at the temperature specified at 134°C
---	---------------------------	--	--

Chapter 10 Composition of the Units

- patient's chair;
- the doctor's stationary instrument unit;
- doctor's mobile Cart instrument unit (if necessary);
- hydroblock with spittoon and assistant tools;
- multifunctional wired pedal;
- multifunctional wireless pedal (if necessary);
- shadowless lamp with mounting to the unit;
- a shadowless lamp with a ceiling mount (if necessary)
- monitor (if necessary)
- monitor holder (if necessary)
- doctor's chair (if necessary)
- assistant chair (if necessary)

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USER MANUAL
FOR
SIGER DENTAL UNIT

Foreword

Before the unit is used, please read the operation instructions carefully, and avoid any wrong operation. Please follow the operation procedures during operation and maintain proper maintenance.

Please read "SAFETY" and "NOTE" paragraphs carefully and strictly perform them to avoid any personal injury or damage to the operator or patients.

If any abnormalities occur during the use, please contact the local dealer or SIGER without delay.

1. Symbols

1.1 WARNING

This indicates a situation which may lead to injuries to persons or damages to the unit, and/or the loss of the information if ignored.

1.2 NOTE

This indicates the important information for the user and/or service personnel.

1.3 General instructions

Please handle with care.

Use side up.

Humidity sensitive.

No stacking under other cargo.

See the random file.

It is forbidden to discard the marking of electronic equipment casually, and after the expiry date, the equipment should be handled according to local laws and requirements to avoid polluting the environment and causing the loss.

1.4 Classification: Type I Class II

1.5 Ground mark

1.6 Warning sign

1.7 Closing sign

1.8 Red switch is compressed during work

2. Safety information

Please always follow the following safety precautions to avoid injuries to persons or damages to the unit.

Please follow all the warnings and instructions provided with the unit. If an instruction contradicts with the safety information, please pay attention to safety information, as you may not understand the instructions. If you cannot solve the contradiction, please contact our service personnel in stock field.

Please always shut off the water supply, air supply and the main power before performing any repairing or maintenance.

Please do not place the unit on shaky floor, cart and shelf or it may tumble over and badly damaged.

Please do not touch the unit near rotation or rollers.

The unit is not suitable for use in the presence of a flammable and/or explosive mixture with air, oxygen or other gases, where such gases may accumulate in concentration (in closed space).

Please do not place or hang any heavy things on the unit.

Please do not connect and use incompatible wires which may affect its performance, or even cause

Manufacturer: Zhuhai Siger Medical Equipment Co., Ltd.

Manufacture & Registration Address: Building 2, No. 1 Chuangxin Yi Road, Tangjiawan Town, Zhuhai City, Guangdong Province, P.R. China

Tel: +86-756-3881018

Fax: +86-756-3881028

Website: <http://www.siger.cn>

E-mail: info@siger.cn

Post Code: 519000

Doc No.: SG-S90-RD-58 Version: A1

Date of issue: 11/2021

EU Authorized Representative: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffeistrasse 80, 20537 Hamburg, Germany

Tel: +49-40-2513175

Fax: +49-40-255726

email: shholding@hotmail.com

Foreword

Before any operation is performed, please read the operation instructions carefully and keep it for future reference. Please strictly follow the operation procedures during operation and perform proper maintenance.

Please read "WARNING" and "NOTE" paragraphs carefully and strictly perform them to avoid any damage to the unit or injury to the operator or patients.

If any malfunction occurs during the use, please contact the local dealer or SIGER without delay.

• Symbols

WARNING

This indicates a situation which may lead to injuries to persons or damages to the unit and/or its parts if the information is ignored.

NOTE

This indicates the important information for the user and/or service personnel.

Guard against damp.

Fragile. Handle with care.

This side up.

Humidity limitation.

No stowing under other cargo.

See the random file

It is forbidden to discard the marking of electronic equipment casually, and after the expiry date, the equipment should be handled according to local legal requirements to avoid polluting the environment and injuring the user.

Classification: Type I Class B

Ground mark

Opening sign

Closing sign

IPX4 Foot switch waterproof rating mark

• Safety Information

- * Please always follow the following safety precautions to avoid injuries to persons or damages to the unit:
- * Please follow all the warnings and instructions provided with the unit. If an instruction contradicts with the safety information, please pay attention to safety information, for you may misunderstand the instruction. If you cannot solve the contradiction, please contact our service personnel to seek help.
- * Please always shut off the water supply, air supply and the main power before performing any cleaning or maintenance.
- * Please do not place the unit on shaky floor, cart and shelf or it may tumble over and badly damaged.
- * Please do not install the unit near radiators or heaters.
- * This unit is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide, where such gasses may accumulate in concentration (in closed space).
- * Please do not place or hang any heavy things on the unit.
- * Please do not connect unit with incompatible wires which may affect its performance, or even cause fire or electric shock.
- * Please do not place any thing in the unit so as to avoid electric shock or fire. Please do not spill any liquid into the inner part of the unit.
- * Please do not disassemble the unit so as to avoid electric shock. Contact technical service personnel whenever the unit needs repairing.
- * Please do not open or remove the cover plate so as to avoid injuries caused by dangerous voltage or other dangers. An improper assembly may cause electric shock during the subsequent operation.
- * Please switch off the power and seek help from technical service personnel if:
 - (1) any break-down or wear-down of the power, plug or cable occurs.
 - (2) unknown liquid spills into the inner unit or water flows into the inner circuits.
 - (3) the unit malfunctions, even if it is operated according to the instructions.
- * Please only adjust the control items mentioned in the instructions. Any adjustment of other control items may cause damage to the unit.
- * Please avoid operating the unit during thunderstorm or electric shock may occur. If possible, please turn off the power supply during thunderstorm.
- * Please do not use damaged or loose plug. An incomplete insertion of plug may cause electric shock or fire.

- * Please make sure to use earthed plug and AC outlet. Bad earthing may cause electric shock to the user or damage to the unit.
- ⚠ **WARNING!** When treating patients who have pacemakers and/or hearing aids, please always consider the effects which the dental unit may have on them.
- ⚠ **WARNING!** The maximum load of the dental seat is 135kg. Please never exceed this limit.
- ⚠ **WARNING!** The maximum load of the table tray attached to the dentist's instrument board is 3.0kg, which should never be exceeded.
- ⚠ **WARNING!** Replace the dental burrs only when the dental instrument has stopped completely, otherwise the cartridge may be damaged and the dental burrs may cause a personal injury.
- ⚠ **WARNING!** This unit is used by authorized and fully trained personnel only (doctors or assistants).
- ⚠ **WARNING!** The manufacturer shall not be responsible for any damage or injury caused by mishandling, neglect or improper use.
- ⚠ **WARNING!** Before leaving the surgery, shut off the water supply, air supply and the main power.
- ⚠ **WARNING!** Wastewater and waste emissions should comply with national laws and regulations, and the amalgam collection is equipped with an amalgam separator.

The equipment is expected to have a service life of five years. Production date see nameplate

Please contact with manufacturer and local competent authority if have any serious incident related with the dental units.

• General Information

Information of optional parts of V series is included in this manual which may not suit for your units.

These instructions describe the operation and maintenance procedures for the dental units of V series. All the information provided is available till the publishment of the Manual. All the information, specifications and illustrations contained herein are subject to change. SIGER reserves the right to make changes without prior notice.

No addition, modification, distribution, republication or sale of this manual, in any form or any manner (electronic, mechanic, or graphic), is allowed unless it is authorized by Siger in a written form in advance.

SIGER reserves the right to make explanation and changes of all the information contained in the Instruction.

⚠ **NOTE!** The guarantee will be disclaimed when the following precautions are not taken:

- (1) Observe the terms specified in the guarantee letter.
- (2) The unit must be used according to the instructions herein.
- (3) During indoor electric wiring, a 3×1mm² power cable should be provided.
- (4) The unit should be installed in a fixed way.
- (5) Repairs, modifications, adjustments and any other operation which involves gaining access to the inside of the unit must be carried out by SIGER technicians or authorized technicians only.
- (6) Replacements for any spare parts and instruments of the unit should be provided by SIGER otherwise the performance and safety may be affected. Regarding the spare parts, please refer to the packing list.

* Installers and maintenance personnel must meet the following requirements:

- a) Education:
 - At least 18 years old – 9 years intensive reading experience (school).
- b) Knowledge:
 - Read and understand 'westernized Arabic' numerals when written in Arial font.
 - Read this manual before use.
- c) Language understanding:
 - Can understand English language in the User Manual and displayed on touchscreen console.
- d) Permissible impairments:
 - Mild reading vision impairment or vision corrected to log MAR 0, 2 (6/10 or 20/32)
- e) Experience
 - Training by SIGER or fully trained for dental unit.

* The user must meet the following requirements:

- a) Person Group and Occupation: Adult-qualified dentist and assistant
- b) Education:
 - At least 18 years old – 9 years intensive reading experience (school).
- c) Knowledge:

- Read and understand 'westernized Arabic' numerals when written in Arial font.
- Read this manual before use.
- Understand stomatology.
- Understand hygiene.
- d) Language understanding:
 - Can understand English language in the User Manual and displayed on touchscreen console.
- e) Permissible impairments:
 - Mild reading vision impairment or vision corrected to log MAR 0, 2 (6/10 or 20/32)
- f) Experience
 - Fully trained for dental unit.

* Operating principle

Dental Units is a combination of dental equipment components used to carry, position and control the devices used in the practice of dentistry.

Contents

Chapter 1 Description of the Unit	5
1.1 Intended Use	5
1.2 Description of the Unit	5
Chapter 2 Operation Instruction	6
2.1 Description of Control Panels and Button	6
2.2 Buttons and Functions of Doctor's and Assistant's Control Panel	7
Chapter 3 Operation Instruction	8
3.1 Main Switch	9
3.2 Operation light	9
3.3 Foot Pedal	9
3.4 Height Adjustment of Balance Arm	11
3.5 Instrument Tray	11
3.6 Pull rods of Swing Type Instrument Table	11
3.7 Headrest	11
3.8 Syringe	11
3.9 Dental instrument	12
3.10 HV Suction	12
3.11 Saliva Ejector	12
3.12 Distilled Water System	12
3.13 Disinfection of the Pipeline	12
3.14 Safety Devices	13
Chapter 4 Cleaning, Disinfection and Maintenance	13
4.1 Disinfection	13
4.2 Cleaning of Suction Tube and Filter	13
4.6 Syringe cleaning	13
4.7 Syringe sterilization	14
4.8 Maintenance of the Unit	14
Chapter 5 Dental unit installation, commissioning and maintenance	14
5.1 Installation Conditions	14
5.2 Requirements for Installation Field	14
5.3 Transportation, package and disassembling	16
5.4 Equipment installation method	16
5.5 Equipment debugging method	17
5.6 Water sampling connection point	17
5.7 Water/Air schematic diagram and electrical schematic diagram	17
Chapter 6 Common faults and exclusion	22
Chapter 7 Electromagnetic Compatibility	22
Chapter 8 Safety Markings	25
Chapter 9 Composition of the Units	26

Chapter 1 Description of the Unit

1.1 Intended Use


The dental units are usually AC-powered, in which a patient sits. The device is intended to properly position a patient to perform dental procedures. A dental operative unit is attached.

Population: adult. **Contraindications:** None.

Side-effect: The side-effects of dental unit are the noise, eye fatigue of the dental doctor by the dental light.

1.2 Description of the Unit

1.2.1 Panorama of the unit

Main parts of dental unit

1	HV suction tube	9	Doctor's control panel
2	Saliva ejector tube	10	Tray
3	Assistant's syringe	11	Headrest
4	Dental spittoon	12	Foot pedal
5	Hydro-group	13	Assistant's control panel
6	Operation light	14	Patient chair
7	Instrument table arm	15	Left armrest
8	Doctor's syringe	16	Instrument table

Note: syringe means the 3-way syringe

1.2.2 The dental operative unit connected with dental unit

Dental turbine handpiece, dental suction handpiece (HV suction and Saliva ejector) and vacuum pump can be connected with dental unit to use.

Note: Dental suction handpiece (HV suction and Saliva ejector) and vacuum pump in the IFU are not individually CE certified by DNV.

1.2.3 Instrument table

V1000 Swing type

V1000 Hanging type

V2000/V5000 Swing type

V2000/V5000 Hanging type

1.2.4 Nameplate

The nameplate is posted on the arm connected to the hydro-group. It contains the basic parameters of the unit.

Dental Unit
 Input: AC230 V 50 Hz
 SN: V1000-XXXXXXX
 Operation Mode: Continuous operation with intermittent loading
 Duty Cycle For Chair: 2 min ON/18 min OFF

MODEL: V1000
 Power Consumption: 900 VA
 Safe working load of patient chair: 150 kg
 Zhuhai Siger Medical Equipment Co., Ltd
 Address: Building 2, No. 1 Chuangxin Yi Road, Tangjiawan
 Town, Zhuhai City, Guangdong Province, P.R. China
 Shanghai International Holding Corp. GmbH (Europe)
 Eifelstrasse 80, 20537 Hamburg, Germany

(01)06974276290073(11)XXXXXX(21)XXXXXX

1.2.5 Technical parameters

The technical parameters of the unit are as follows:

Type of protection against electric shock	I	Main fuse	T10AL250V
Degree of protection against electric shock*	Type B	Rated illuminance of operation light	8000-30000Lx
Mode of operation	Intermittent operation	LED operation light	AC12V 8VA
Degree of water-proof	Non-water-tight	Power consumption	900VA
Input voltage	AC230 V	Frequency	50/60Hz
Air pressure limits	(0.7±0.1) MPa	Air flow rate limit	Greater than 80 NI/min
Humidity limit	Dew point not greater than -20℃ at atmospheric pressure	Particulate contamination limit	Not greater than 100 particles per cubic meter for 1 μm to 5 μm particle size
Oil contamination limit	Max. 0.5 mg/m ³	Air filter net	40 μm
Water pressure limits	200 kPa~600 kPa	Water flow rate limit	Greater than 5 l/min
Water hardness limit	Less than 2,14 mmol/l (<12 °dH)	PH limits	6,5 to 8,5
Water source requirement	Conformance to local drinking water regulations	Water filter net	75μm
Heater	AC24V 120W	Software version	A0

1.2.6 Ambient conditions

Normal working		Shipment and Storage	
Temperature	5 °C~40 °C	Temperature	-40 °C~55 °C
Relative humidity	≤80%	Relative humidity	≤93%
Atmospheric pressure	70kPa~106kPa	Atmospheric pressure	50kPa~106kPa

Suggestion Shelf Life: 10 years from the manufacture date

1.2.7 Movement Range (Max length: 2100mm, Max width:1780mm)

- | | |
|---|--|
| ① Dimensions: Length: 190CM, Width: 85cm, Height: 180cm | ② maximum seat height: ≥820mm |
| ③ minimum seat height: ≤420mm | ④ inclination of the backrest: 105°~185° |
| ⑤ elongation range of the headrest: >200mm | ⑥ maximum caster angle of the seat: ≥10° |
| ⑦ instrument table lifting distance 480 mm | ⑧ rotation range of the spittoon: ≥180° |

Chapter 2 Operation Instruction

2.1 Description of Control Panels and Button

2.1.1 Illustration of Doctor's and Assistant's Control Panel

V2000/V5000 control panel (touch button), V1000 control panel (press button)

V1000 doctor control panel

V1000 assistant control panel

V5000/V2000 main control panel

V5000/V2000 assistant control panel

2.1.2 Description of buttons

Buttons	Description	Buttons	Description
	Seat up		Rinsing position
	Seat down		Settings
	Backrest up		Spittoon rinsing
	Backrest down		Cup filling (cold water)
	Programmed position A		Mouthwash heating
	Programmed position B		Operating light(ON/OFF)
	Zero position		Optical fiber (ON/OFF)
	Emergency position		X-ray image viewer (ON/OFF)

2.2 Buttons and Functions of Doctor's and Assistant's Control Panel

2.2.1 Seat up/down

Press the button on the control panel, the seat would move up.

Press the button on the control panel, the seat would move down.

2.2.2 Backrest up/down

Press the button on the control panel, the backrest would move up.

Press the button on the control panel, the backrest would move down.

2.2.3 Programmed positions

Setting of programmed position A: operate the dental chair to the desired position and press button "A" for 5 seconds. The current chair position will be automatically saved as a programmed position after an indication of audible signal.

Operating of programmed position A: short press button "A" , the dental chair would automatically move to the saved position A. At the same time, the operation light turns on.

Setting and operating of programmed position B : the same with programmed position A.

2.2.4 Rinsing position

Setting of rinsing position: operate the dental chair to the desired rinsing position and press button for 5 seconds. The current chair will be automatically saved as a memory rinsing position after an indication of audible signal.

Setting of linkage between rinsing position and spittoon rinsing: press the "settings" button for 5 seconds, the indicator light is on. Press the "rinsing position" button , then the linkage between the chair movement and spittoon rinsing is on after an audible indication and the indicator light goes out.

Operating of rinsing position: short press the "rinsing position" button , the dental chair will automatically move to the saved memory rinsing position at which water rinses to the spittoon and fills the cup and the operation light goes out automatically. After that, press the button again to return the chair to the last working position and the operation light turns on automatically.

2.2.5 Emergency position

Setting of emergency position : operate the chair to the desired emergency position, then long press

the "emergency position" button  for 5 seconds. After the indication of an audible signal, the current position is saved as the emergency position.

Operating: short press the "emergency position" button , the chair will move to the saved emergency position and at the same time, the operation light turns on automatically.

2.2.6 Zero position (Exit position)

Setting: operate the chair to the desired zero position, then long press the "zero position" button  for 5 seconds. After an indication of an audible signal, the current position is saved as the zero position.

Operating: short press the "zero position" button , then the chair will move to the saved zero position and the operation light goes out.

2.2.7 Spittoon rinsing

Time setting: long press the "settings" button  for 5 seconds, the indicator light will be on. Then short press the "spittoon rinsing" button . When the desired time is reached, press again the "spittoon rinsing" button  to stop. The indicator will go out and the setting is accomplished.

Operating: short press the "spittoon rinsing" button  to rinse the spittoon once for the saved time, or press again the button to stop.

2.2.8 Cup filling

Setting: long press the "settings" button  for 5 seconds, the indicator light is on. Then short press the "cup filling" button  to fill in the cup with water. When the desired water volume is reached, press again the button to stop. The indicator light will go out and the setting is accomplished.

Operating: short press the "cup filling" button , the water will come out with desired volume. Press again the button  to stop the water.

2.2.9 Mouthwash heating

Short press the "mouthwash heating" button , the indicator light will flicker. When the water is heated to the preset temperature, the indicator light is on, which means the heating is accomplished. During using, the indicator light will work in a cyclical way between flickering and lighting.

2.2.10 X-ray image viewer

Operating: short press the "X-ray image viewer" button , the viewer will be on; it will be turned off automatically after 3 minutes or be turned off by pressing again the button.

Factory default: long press the button  for 5 seconds. After an indication of beep, the chair will firstly move up to the highest position then down to the lowest position. Now the 5 programmed positions will be reset to restore factory default.

2.2.11 Optical fiber

Operating: when the instrument is picked up or pulled out, short press the "optical fiber" button , the fiber-optic light will be on. Press the button again to turn it off (this function applies to V2000/V5000).

Function shift of button "O/B" on the foot pedal: long press the "optical fiber" button  for 5 seconds. After an indication of beep, the button "O/B" on the left of the foot pedal will change its function from chair memory position control to the switch of operation light on/off. The function of the button on the control panels remains the same. Long press the same button for 5 seconds again to shift the two functions. The doctor can adjust this function button at his will (this function applies to V1000).

Automatic testing: when the power is off, press the button  and hold it still. Then turn on the power and release the button after an indication of beep. At this time, the chair will move up and down automatically, and other electric equipments such as the operation light and the heater will start working simultaneously to test the condition of the electric parts. Press any button to enter the main program (this function applies to V1000 for after-sales repair).

2.2.12 Operation light

Operating: short press the "operation light" button , the operation light will be on. Press the button again to turn it off.

Long press the "operation light" button  for 5 seconds. The indicating function of audible signal will be off after an indication of beep. Another long press can turn it on again.

Chapter 3 Operation Instruction

⚠ WARNING! Before any operation, make sure the unit is correctly installed and adjusted according to the Manual of Installation and Adjustment.

3.1 Main Switch

Press the main switch to "I", the green power indicator on the switch will illuminate. By now, the unit is power-on.


3.2 Operation light

3.2.1 Halogen Light (Standard configuration for V1000)

Short press the "operation light" button  on the doctor or assistant's control panel, the light can be turned on/off. It should be noted that the automatic control of the operation light in the case of programmed positions is not subject to the button.

There is a manual switch with the halogen light (3 levels):

☉ High intensity; ☉ Low intensity; ☉ Off level

3.2.2 LED light (Standard configuration for V2000/V5000)

The angle of light handle can be adjusted according to the user's need.

When the unit is turned on, the illuminance of the LED light in default is Level 3.

The sensing range of the LED light sensor is 50mm-120mm. Put your hand close to the sensor and let go after 0.5 second, the light can be turned on/off.

When the light is on, put your hand close to the sensor. The illuminance will go from low to high (level 1 to level 6) (V2 LED light has four levels), and then high to low. Take your hand away when it reaches the desired level.


The reflector type V2 LED operation light has blue light stopping function. When doctor use light sensitive resin for tooth filling, slide switch can turn the light to yellow to prevent curing of sensitive resin. During the normal dental operation, please open up the slide switch.

3.2.3 Assembly and disassembly of the light handle

Assembly: Grab the light with your right hand->grab the left handle with your left hand and press the handle button with your thumb->Put the handle in the slot on the rotating shaft and release the handle button. Same method for right handle.


Disassembly: Grab the light with your right hand->grab the left handle with your left hand and press the handle button with your thumb->Pull out the handle powerfully

Same method for right handle.

3.3 Foot Pedal

3.3.1 Electronic Multi-functional lever type foot control

Chair movement control: slide the chair movement control, the movements of seat up and down, backrest up and down can be performed.

Left button:

① When no dental instrument is picked up or pulled out and short press the button, the dental chair will move to the rinsing position. Short press it again, the chair will return to the last working position. Long press

the button (for 2 seconds) and after an indication of audible signal, the chair will move to programmed position A.

② When a dental instrument is picked up or pulled out, it is an activating/deactivating switch of water spray. When the water spray is activated, the setting button indicator is on. When the water spray is deactivated, the setting button indicator is off.

Right button:


① When no dental instrument is picked up or pulled out and short press the button, the dental chair will move to the zero position. Long press it (for 2 seconds) and after an indication of audible signal, the chair will move to programmed position B.

② When a dental instrument is picked up or pulled out, it is an air blast button. The air keeps delivering when the button is kept pressing. The flow of air will stop when the button is released.

Driving rod: When a dental instrument is taken out or pulled out, the instruments are ready for work. At this time, the foot pedal is available for controlling the speed or power of instruments. From left to right, the instruments will be sped up gradually.

Handle: Lift the handle with your foot to move the foot-control.

3.3.2 Electronic Multi-functional pedal type foot control

Chair movement control: slide the chair movement control, the movements of seat up and down, backrest up and down can be performed.

Left button:

① When no dental instrument is picked up or pulled out and short press the button, the dental chair will move to the rinsing position. Short press it again, the chair will return to the last working position. Long press the button (for 2 seconds) and after an indication of audible signal, the chair will move to programmed position A.

② When a dental instrument is picked up or pulled out, it is an activating/deactivating switch of water spray. When the water spray is activated, the setting button indicator is on. When the water spray is deactivated, the setting button indicator is off.

Right button:


① When no dental instrument is picked up or pulled out and short press the button, the dental chair will move to the zero position. Long press it (for 2 seconds) and after an indication of audible signal, the chair will move to programmed position B.

② When a dental instrument is picked up or pulled out, it is an air blast button. The air keeps delivering when the button is kept pressing. The flow of air will stop when the button is released.

Foot control pedal: When a dental instrument is taken out or pulled out, the instruments are ready for work. At this time, the foot pedal is available for controlling the speed or power of instruments. From left to right, the instruments will be sped up gradually.

Handle: Lift the handle with your foot to move the foot-control.

3.3.3 Pneumatic Multi-functional lever type foot control

Chair movement control: slide the chair movement control, the movements of seat up and down, backrest up and down can be performed.

SP/A button: short press the button to move the chair to rinsing position. Short press it again, the chair will return to the last working position. Long press this button for about 2 seconds until an audible signal is heard, the chair will move to programmed position A automatically.

0/B button: short press the button to move the chair to zero position. Long press the button for about 2 seconds until an audible signal is heard, the chair will move to programmed position B automatically.

Left button: when a dental instrument is picked up or pulled out, it's an activating/deactivating switch of

3.9 Dental instrument

Operating: pick up or pull out a dental instrument and operate the foot pedal, the dental instrument will start to work.

The input air pressure of the dental instrument can be regulated through the air regulating valves under the instrument table. Screw clockwise to turn down the air pressure, or vice versa.

3.9.1 Water volume regulating

Water volume can be regulated through the water regulating valves under the instrument table. Screw clockwise to turn down the water volume, or vice versa. One valve applies to only one instrument.

3.9.2 Spray volume regulating

Spray volume can be regulated through the spray regulating valve under the instrument table. Screw clockwise to turn down the spray volume, or vice versa. The valve applies to all instruments.

NOTE:

Dental instrument is purchased and used at the discretion of professional doctors. It should be used in accordance with the requirements of the CE standard and obtain a CE certificate. The recommend dental instrument see following table:

Name	Model/specification	Manufacturer
Dental Turbine handpiece	THALIA	TEKNE DENTAL SRL

3.10 HV Suction

The suction system starts working once the HV suction tube is pulled out from the tube holder. The regulating handspike is used to adjust the suction flow.

For V1000, when the tube is put back in the holder, the suction stops; for V2000/V5000: when the tube is put back in the holder, the suction system won't stop until 3 seconds later. This delay assures the hygiene of the suction tube.

Under static (i.e. no flow) conditions, the maximum vacuum pressure provided to the dental unit at the HV Suction connection port is 25kPa.

3.11 Saliva Ejector

The ejector system starts working once the saliva ejector tube is pulled out from the tube holder.

When the tube is put back in the holder, the ejector system stops.

Under static (i.e. no flow) conditions, the maximum vacuum pressure provided to the dental unit at the Saliva Ejector connection port is 15kPa.

3.12 Distilled Water System

The distilled water system is composed of a selecting switch and a distilled water bottle at the lower back of the hydro-group. Distilled water in the bottle or sanitized tap water is sent to all the dynamic instruments on the dentist's and assistant's instrument table, so the bottle must be filled with medical distilled water in time.

Water filling:

① Switch air pressure lever to "Decompress" position and release the air in the bottle.

② Screw off the bottle anticlockwise and fill up with distilled water.

③ Screw on the bottle clockwise and switch air pressure lever back to "Plus press" position.

④ Adjust the water selecting switch to "distilled water".

NOTE:

The bottle must only be opened when the air in it is completely released.

We suggest using distilled water during treatment so as to avoid the negative effect on the unit's pipeline and the service life of the instruments caused by impurities in the tap water.

3.13 Disinfection of the Pipeline

Keep the water tank empty and switch the air pressure lever to "Plus press" to let out the waste and

water from the pipeline.

Switch the air pressure lever to "Decompress". Remove the white distilled water tank and install the white tank with disinfectant.

Switch the air pressure lever to "Plus press" and operate the dental instruments one by one until disinfectant sprays out and lasts for seconds. Keep the unit in non-working status for one or two hours.

Switch the air pressure lever to "Decompress", remove the red disinfectant tank and reinstall the white distilled water tank. Switch the air pressure lever again and clean the internal water line the same way. Thus the disinfection is completed.

It is recommended to use five thousandths of hydrogen peroxide as a disinfectant, disinfected once a day. The pipeline has to be replaced when it was disinfected for 1,900 time. If there is any problem with the pipeline about appearance, please contact Siger to replace on time.

3.14 Safety Devices

The bottom of the dental chair is provided with a safety device which stops the chair from moving down in the presence of obstacles. The chair will move up to allow users to remove the obstacles.

When the backrest is moving down, it will stop automatically in the presence of obstacles (this applies to V2000/V5000).

Any button can stop the chair from moving.

When the dental instrument stops working, air spray will be continued 2 seconds to activate the anti-suction function. This is effective to prevent cross infection (available on V2000/V5000).

When a dynamic instrument is in use, the dental chair will be under safety protection and no movement is available.

If one dynamic instrument is in working status, none of the others can be activated. This interlocking system enables the user to replace a dental burr on one instrument while another is in use (this applies to V2000/V5000).

Chapter 4 Cleaning, Disinfection and Maintenance

4.1 Disinfection

As for the sterilization of the dental instrument, please refer to the instruction manual for dental instrument.


4.2 Cleaning of Suction Tube and Filter

Take off the suction tube handles from the holder. Pull the tube out from the filter.

Pull out the filter net from the filter and dump the waste, which should be disposed in a centralized way.

Clean the suction tubes and the filter mesh with water and then have them disinfection.

Reinstall the suction tubes and the filter mesh. Once-a-day cleaning is suggested.

Filt

4.3 Cleaning of Dental Spittoon

Take out the lid and the filter net and dump the waste, which should be disposed in a centralized way.

Clean the spittoon, the lid and the filter net and have them disinfection.

Reinstall the lid and the filter mesh. Once-a-day clean is suggested.

4.4 Cleaning of Air Return Tank

Screw off the tank anticlockwise, dump the waste in the tank and clean the noise-absorbing sponge.

Screw on the tank clockwise.

4.5 Other Cleaning

Dental chair does not need special maintenance. Suggest to wipe the upholstery after each treatment.

The external parts which are frequently touched by the doctor, the nurse or the patient must be wiped and disinfected with alcohol or disinfectant once a day.

⚠ NOTE: please do not use abrasives or acid products.

4.6 Syringe cleaning

1) Rinsing

Using daily drinking water rinsing with a brush, removal of organic residues on the 3 way Syringe. then perform disinfectant soaking.

detachable suction system for cleaning

2) Disinfectant soaking

Disinfectant: Medical alcohol.

Put the rinsed devices on medical utensils, Pour into medical alcohol , the 3 way Syringe. should be submerged medical alcohol, soak for 10 minutes, then scrub the surface of the 3 way Syringe with the brush. in the end perform end rinsing.

3) End rinsing

With purified water or distilled water to flush, Clean the residual alcohol. Since then, it can be used to sterilization.

The Syringe has to clean after each use.

The Maximum clean times: 300

4.7 Syringe sterilization

1)Syringe loading:

Put Syringe samples into heat-sealing sterilization pouches. The Syringe were placed in trays , then transferred to sterilizer cavity.

2)Moist heat sterilization:

Program: Quick B.

Sterilization temperature: 134°C.

Sterilization time: 4 minutes.

Drying time: 15 minutes.

Pressure: 2.1 bar.

The Syringe has to sterilization after each use.

The Maximum sterilization times:250

4.8 Maintenance of the Unit

NOTE: proper maintenance can extend the service life of the unit!

Check the electric supply every month. Water pressure and air pressure must be in accordance with the unit requirement. Keep the hoses smooth and unblocked.

Check the chair' s movement and control panel.

Clean the spittoon after use everyday. Wash the filter, pull out the filter net and reinstall them after cleaning to avoid air leakage.

As for the cleaning and lubrication of dental instrument, pleased refer to the instruction manual of dental instrument.

Water filter: in order to ensure the proper function of the unit, a water filter is installed at the water inlet to block the impurities in water. After a long-time operation, the filter collector may be filled with contaminant and the water flow may be obstructed. Therefore, the filter collector must be cleaned or replaced timely.

Generally, if any of the following cases happens, the filter collector should be cleaned or replaced:

①the filter collector has been used for more than 1 year;

②the filter collector is contaminated;

③the water coming out becomes turbid.

How to replace the filter collector: open the hydro-group and remove the water filter. Then remove the filter shell anticlockwise, take off the seal ring and the collector successively. After cleaning or replacing, reinstall all the parts. Please take care of the seal ring to avoid leakage.

Chapter 5 Dental unit installation, commissioning and maintenance

5.1 Installation Conditions

⚠ Caution: The equipment must be installed by our professionals or any personnel authorized by us.

This equipment shall be installed at a clean, dry, cool and well-ventilated location, and sufficient space shall be reserved around the unit to facilitate access during routine maintenance.

The equipment shall not be installed or operated at any location:

A. exposed to humidity, heavy dust, poor ventilation or direct sun; or



B. subject to substantial temperature or humidity change, such as place near air conditioner or heater.

⚠ Caution: Improper installation may cause damage to the equipment!


5.2 Requirements for Installation Field

The equipment shall be located based on general layout of the treatment room, lighting condition and operation convenience, and should be installed on flat, horizontal and solid ground at a clean, dry, cool and well-ventilated location.

Standard Layout for V Series Comprehensive Treatment Platform:

User Manual for V Series Dental Unit

The Drain Hose shall be $\Phi 50$ mm PVC hose, and water supply hose and air supply hose are both $\Phi 15$ mm and furnished with HP aluminum-plastic complex tube port for connection of ball valve. The power supply shall be of 230V and grounded properly, with wire diameter $\geq 1.5 \text{ mm}^2$. The underground hose line shall lie at the same elevation as that of the floor ground as possible, or not above 100 mm higher than the floor. Connection between PC and endoscope (if any) shall be laid through the underground hose line. The VGA cable and AV video cable shall be connected to the computer.

Supply voltage 230 VAC
Air pressure limits $(7 \pm 1) \text{ kPa}$
Water pressure limits 200~600 kPa
Maximum wastewater flow rate: 6 L/min

Power frequency 50/60Hz
Air flow rate limit Greater than 80 NI/min
water flow rate limit Greater than 5 l/min
Minimum gradient of the wastewater lines: 5°

5.3 Transportation, package and disassembling

Transport the equipment to the clinic according to the request on the package.

The equipment can be packed in one case or several cases, the case can be opened by unscrew the bolt. Pick out the packing list, check the parts according to the list and contact the manufacturer if any doubt.

5.4 Equipment installation method


5.4.1 Patient chair installation

Take apart the front cover and bottom cover, unscrew the 4 bolts which fix the chair on the board, then move the chair to the appointed position, please keep the chair balance during the movement.

Assemble the foam part on the chair.

Insert the headrest into the backrest, adjust the tightness of the headrest pulling plate through adjusting the fastening screw.

Using the 3pcs 8mm setscrew to fix the bottom board of the chair on the ground.

5.4.2 Installation of operating light

Insert the light pole into the corresponding hole on the hydro-group crutch, pass the wire which are connected to the light arm through the light pole and connect it to the light, then plug the wire of light arm onto the socket of hydro-group PCB.

5.4.3 Installation of water, air and electricity of ground-group

Put the $G1/4 \times \Phi 8$ direct connector into the air inlet valve, connect one side of the blue $\Phi 8$ PA pipe into the direct connector, and the other side to the general air inlet of the dental unit.

Put the $G1/4 \times \Phi 8$ direct connector into the air inlet valve, connect one side of the white $\Phi 8$ PA pipe into the direct connector, and the other side to the general air inlet of the dental unit.

Insert the drain pipe, suction vent nozzle to the ground hole of $\Phi 50$ PVC drain outlet.

Connect the power and grounding wire well, plug all the connector on the corresponding socket on the PCB of ground-group, organize all the pipe well, then close the hydro-group cover.

Connection method of power cord: Please make sure to turn off the indoor main power before operation. Distinguish the live wire (brown, L), null line (blue, N) and ground wire (yellow-green, G) on the power cord of 2 meters on the device and the power cord, air switch or plug provided indoors by color or sign, and connect firmly according to the corresponding color or position and do insulation treatment.

5.4.4 vacuum pump installation of ground-group

Insert the $\Phi 25$ straight connector into the suction pipe of the vacuum pump, and connect the other side to the negative pressure suction connection pipe of the dental unit. Connect the negative pressure pump cable of the dental unit to the start signal port of the vacuum pump.

Under static (i.e. no flow) conditions, the maximum vacuum pressure provided by the vacuum pump to the dental unit is 25kPa, the air flowrate not less than 90 l/min.

5.5 Equipment debugging method


Main Air Pressure Adjustment:

As shown below, adjust the inlet pressure regulating valve to 550kPa based on the pressure meters on operating chair.

High-Speed Handheld sprayer air pressure adjustment: Adjust the air pressure to 250 kPa with the knob on the bottom of instrument panel.

Low-Speed Handheld sprayer air pressure adjustment: Adjust low-speed handheld sprayer pressure to 250 kPa-300 kPa by using the knob on the bottom of instrument panel.


5.6 Water sampling connection point

Due to the water supply of dental unit is from domestic portable water, the user can install the input water sampling connection point near the water supply port according to actual needs. The specific method is as follows:

Connect G1/4" $\Phi 8$ barrel connector to the water sampling rotary valve of.

Connect the $\Phi 8$ three-way connection between the water input check valve and the water filter, and connect one end of the orange $\Phi 8$ polyether pipe to the three-way connector, and one end into the water sampling rotary valve.


The other end of the water sampling rotary valve is connected to an orange $\Phi 8$ polyether tube of appropriate length.


When sampling, just switch on the water sampling rotary valve, and switch off the water sampling rotary valve after sampling.

5.7 Water/Air schematic diagram and electrical schematic diagram


V1000 Water/Air Supply Schematic Diagram


V2000/V5000 Water/Air Supply Schematic Diagram

V1000 Electric Block Diagram

V2000/V5000 Electric Block Diagram

Chapter 6 Common faults and exclusion

6.1 Dental Chair

Fault	Cause	Action
No Chair Movement	The fuse of the power supply is damaged.	Replace the same type of fuse.
The up and down function of the seat is abnormal.	The dental instrument is damaged or the wiring is loose	Check the line and contact the engineer.
The function of the backrest is abnormal.	The dental instrument is damaged or the wiring is loose	Check the line and contact the engineer.
Movement Chaos	The chip of the CPU is damp or damaged.	Moisture absorption treatment and contact the engineer.

6.2 Operating Light

Bad condensation	The bulb is not installed on the focal spot	Adjustment
------------------	---	------------

6.3 Instrument Table

Fault	Cause	Action
Not enough air comes from the dental instruments.	Air hoses deformed or filter blocked	Adjust the trachea or clean the filter
Not enough water comes from syringe	Throttle adjustment is too small or too large	Adjust the throttle
The suction failed.	The solenoid valve is damaged	Contact the engineer for replacement.
Hydraulic suction force is too small	O-ring is broken or water pressure is lower than required	Replace the new O-ring or adjust the water pressure
Pneumatic suction force is too small	O-ring is broken or water pressure is lower than required	Replace the new O-ring or adjust the water
None warm water supply	Boiler is damaged.	Contact the engineer for replacement.

6.4 Replacement of fuse tube

Shut down the main power supply, then open the hydro-group box, swing out the cover of the fuse tube with medium size cross screwdriver. Replace a new fuse tube, and operate in the opposite procedure as above.

Fuse Tube	Specification
Main fuse tube	T10AL,250V

Chapter 7 Electromagnetic Compatibility

Notice:

- Dental Unit V1000/V2000/V5000/V6000 meets the requirement of standard electromagnetic compatibility of IEC60601-1-2.
- The users should install and use the unit according to the information of electromagnetic compatibility of the documents which are shipped with the unit.
- Portable and mobile radio frequency communication equipment may affect the performance of dental unit V1000/V2000/V5000/V6000. During the usage of these dental units, please avoid strong electromagnetic interference, such as being close to mobile phones, microwave ovens and so on.
- The detail of the guide and the manufacturer's statement, please refer to the annex.

 Warning

- Dental unit V1000/V2000/V5000/V6000 should not be close to or stacked with other equipment during usage. If it must be used close to or stacked with other equipment, you should observe and test the unit whether it could work successfully with its configuration.
- In addition to cables sold as spare parts of the dental unit V1000/V2000/V5000/V6000 by the manufacturer, the use of other accessories and cables may cause radio frequency launching increasing and anti-interference reducing of dental unit V1000/V2000/V5000/V6000.

Annex:

Guidance and manufacturer's declaration – electromagnetic emissions		
V1000/V2000/V5000/V6000 Dental Unit is expected to use in the following electromagnetic environment, buyers or users should ensure that it is in the electromagnetic environments:		
Emission testing	Compliance	Electromagnetic environment – Guide
RF emission CISPR 11	Group 1	V1000/V2000/V5000/V6000 Dental Unit uses RF energy only for its internal function. Therefore, its RF emission is very low and is unlikely to interfere with nearby electronic equipment
RF emission CISPR 11	Class B	V1000/V2000/V5000/V6000 Dental Unit Suitable for use in all facilities, including direct to home House household facilities and public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Compliant	
Voltage fluctuations / flicker emissions IEC61000-3-3	Compliant	

Guidance and manufacturer's declaration – electromagnetic immunity			
V1000/V2000/V5000/V6000 Dental Unit is expected to use in the following electromagnetic environment, buyers or users should ensure that it uses in the electromagnetic environment:			
Immunity test	IEC 60601 Test level	In line with levels	Electromagnetic environment – Guide
Electrostatic discharge (ESD) IEC61000-4-2	±6 kV Contact discharge ±8 kV Air discharge	±6 kV Contact discharge ±8 kV Air discharge	Floors should be wood, concrete or ceramic tiles, if the ground covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC61000-4-4	±2kV On the power supply line ±1kV The input / output lines	±2kV On the power supply line ±1kV The input / output lines	Net power supply should be used with a typical commercial or hospital environment quality
Electrical surge IEC61000-4-5	±1 kV Wire to wire ±2 kV Line to ground	±1 kV Wire to wire ±2 kV Line to ground	Net power supply should be used in a typical commercial or hospital environment quality
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % U_T , 0.5 cycle (U_{of-t} -Shang , >95% SAG) 40 % U_T , 5 cycles (U_{of-t} -Shang , 60% SAG) 70 % U_T , 25 cycle (U_{of-t} -Shang , 30% SAG) <5 % U_T , 5s (U_{of-t} -Shang , >95% SAG)	<5 % U_T , 0.5 cycle (U_{of-t} -Shang , >95% SAG) 40 % U_T , 5 cycles (U_{of-t} -Shang , 60% SAG) 70 % U_T , 25 cycle (U_{of-t} -Shang , 30% SAG) <5 % U_T , 5s (U_{of-t} -Shang , >95% SAG)	Mains should have quality used in a typical commercial or hospital environment. V1000/V2000/V5000/V6000 of dental comprehensive treatment machine user requires continuous operation during power outage it is recommended V1000/V2000/V5000/V6000 dental comprehensive treatment machine using an uninterruptible power supply or battery powered
Power frequency (50 / 60Hz) magnetic field IEC61000-4-8	3A/m	3A/m,50/60Hz	Power frequency magnetic fields should have in a typical commercial or hospital environment, level of power frequency magnetic field characteristics of typical sites
Note: U_T test voltage is applied before the AC voltage.			

Guidance and manufacturer's declaration – electromagnetic immunity			
V1000/V2000/V5000/V6000 Dental comprehensive treatment machine is expected to use in the following electromagnetic environment, buyers or users shall ensure that its use in the electromagnetic environment:			
Immunity test	IEC 60601 Test level	In line with levels	Electromagnetic environment – Guide
Conducted RF IEC61000-4-6	3 V (Effective value) 150 kHz ~80 MHz	3V (Effective value)	Portable and mobile RF communications equipment should not be closer than the recommended separation distance V1000/V2000/V5000/V6000 uses any part of dental unit, including cables. Calculation formula of the distance from the transmitter frequency. Recommended separation distance $d = 1.2\sqrt{P}$
Radiation RF IEC61000-4-3	3V/m 80 MHz ~2.5 GHz	3V/m	$d = 1.2\sqrt{P}$ 80 MHz ~800 MHz $d = 2.3\sqrt{P}$ 800 MHz ~2.5 GHz Type in the: P -According to the transmitter manufacturer's maximum transmitter Rated output power, in Watts (w) as a unit; d -The recommended separation distance in meters (m) ^b . Stationary RF field strength of the transmitter through the investigation of electromagnetic field ^c to determine, in each frequency range ^d should more than meet the low level. Mark the following symbol interference may occur in the vicinity of the device. 
Note 1: in the 80MHz and 800MHz frequencies, high frequency formula.			
Note 2: these guidelines may not be appropriate for all circumstances, electromagnetic propagation affected by absorption and reflection of buildings, objects and human impact.			
^a Fixed transmitters, such as radio (cellular / Cordless) telephones and land mobile radio base stations, amateur radio, AM and FM radio and television broadcasting, its strength in theory, cannot be accurately predicted. For evaluating the electromagnetic environment of stationary radio transmitters, electromagnetic site survey should be considered. If measured V1000/V2000/V5000/V6000 dental unit places the field strength is higher than the above applicable RF compliance with level, it should be observed V1000/V2000/V5000/V6000 dental unit to verify that it works. If observed unusual performance, additional measures may be necessary, such as restructuring V1000/V2000/V5000/V6000 the direction or position of dental unit. ^b 150kHz~80MHz the entire frequency range, field strength should be less than 3V/m.			

Portable and mobile RF communications equipment and the V1000/V2000/V5000/V6000 recommended separation distance between dental comprehensive treatment machine	
V1000/V2000/V5000/V6000 Dental Unit is expected to RF radiation disturbance in a controlled electromagnetic environment of use. According to the maximum rated output power of communications equipment, purchasers or users by using the following recommended maintaining portable and mobile RF communications equipment (transmitters) and V1000/V2000/V5000/V6000 between dental unit minimum distance to prevent electromagnetic interference.	
Rated maximum	Corresponds to the transmitter frequency isolation distance /m

User Manual for V Series Dental Unit

output power of transmitter W	150 kHz ~ 80 MHz $d = 1.2\sqrt{P}$	80 MHz ~ 800 MHz $d = 1.2\sqrt{P}$	800 MHz ~ 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23



Transmitter maximum nominal output power that is not listed in the previous table, the recommended isolation distance d in meters (m) as a unit, available in the relevant transmitter frequency bar formula to determine, where p is the transmitter maximum nominal output power provided by the transmitter manufacturer, in Watts (w) as the unit.

Note 1: in 80MHz and 800MHz on the frequency, high frequency formula.

Note 2: these guidelines may not be appropriate for all situations. Electromagnetic propagation of buildings, objects and the effects of the body's absorption and reflection.

Chapter 8 Safety Markings

Symboles/ Marks	Description	Symboles/ Marks	Description
	Warnings		Date of manufacture
	Caution		Manufacturer
	Serial number		Unique device identifier
	TYPE B APPLIED PART		Temperature limitations
	Indicates a medical device that needs to be protected from moisture.		Refer to instruction manual/ booklet
	Indicates a medical device that can be broken or damaged if not handled carefully.		Medical device
	Stacking limit by 2		WEEE Directive mark

	Temperature limit		European representative
	Atmospheric pressure limitation		CE marked
	Humidity limitation		sterilizable in a steam sterilizer (autoclave) at the temperature specified at 134°C

Chapter 9 Composition of the Units

- patient's chair;
- the doctor's stationary instrument unit;
- doctor's mobile Cart instrument unit (if necessary);
- hydroblock with spittoon and assistant tools;
- multifunctional wired pedal;
- multifunctional wireless pedal (if necessary);
- shadowless lamp with mounting to the unit;
- a shadowless lamp with a ceiling mount (if necessary)
- monitor (if necessary)
- monitor holder (if necessary)
- doctor's chair (if necessary)
- assistant chair (if necessary)

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