Medtronic

Medtronic Burs

Fraises Medtronic / Frese Medtronic /
Knochenfräsen von Medtronic / Fresas
Medtronic / Medtronic boren / Medtronic
bor / Medtronic -poranterät / Medtronic
borr / Brocas Medtronic / Γλύφανα
Medtronic / Wiertła Medtronic / Frézy
Medtronic / Medtronic fúrók / Medtronic
frezler / Borene fra Medtronic

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Product Information & Instructions

Informations et instructions concernant le produit /Informazioni
sul prodotto ed istruzioni per l'uso /Produktinformation und
Gebrauchsanweisung /Instrucciones e información sobre el producto
/ Productinformatie en instructies / Produktoplysninger og vejledning
/ tuotetta koskevat tiedot ja käyttöohjeet / Produktinformation och
instruktioner / Informação do produto e instruções / Πληροφορίες
προϊόντος και οδηγίες / Informacje o produkcie i instrukcje / Údaje o
výrobku a pokyny / Termék információ és használati utasítás / Ürün
Bilgileri ve Talimatlar / Produktinformasjon og instruksjoner

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DEVICE DESCRIPTION

Medtronic burs are interchangeable drill bits and burs intended for use with the Medtronic ENT and Xomed line of otologic drill handpieces. Burs are supplied in a variety of materials for improved cutting performance including: stainless steel with titanium nitride coating, carbide and stainless steel, tool-steel, carbide, and stainless steel with diamond studded heads.

OPERATOR NOTICE — The clinician is strongly recommended to have a thorough understanding of and experience with powered surgical equipment prior to use. Medtronic Xomed provides surgical technique booklets relating to the application of the various cutting blades and accessories for the XPS® 3000 System. The health care professional should seek out and review all other clinical reference materials as dictated by the individual patient's clinical condition.

INDICATIONS FOR USE

Medtronic Xomed burs are indicated for use in incising or removing bone and tissue during general otorhinolaryngology, head and neck, or otoneurological surgery.

CONTRAINDICATIONS

None are known.

WARNINGS

- DO NOT use burs above the speed indicated on the bur label. Exceeding speed limits may cause bur to break.
- Use adequate irrigation. The use of a bur without irrigation may cause an inordinate amount of heat buildup resulting in thermal injury to tissue.
- Always keep the cutting tip away from fingers and loose clothing.
- Operate the blade or bur only after the appropriate anatomical landmarks and the intended surgical site have been confirmed.
- Ensure bur is properly secured in the handpiece prior to surgeon activation.
- DO NOT change burs with handpiece running.
- Test for bur wobble (eccentricity) at desired speed prior to use. Use a bur guard, if needed. Select a new bur if wobble is observed or reduce speed.
- Excessive pressure applied to bur may cause bur fracture which may cause injury.
- DO NOT use dull, damaged, or bent burs. Change burs frequently for optimal performance. Use of dull burs can reduce handpiece cutting effectiveness and cause handpiece temperature to increase.
- DO NOT attempt to resharpen burs. Worn burs should be replaced with new ones frequently to ensure effective cutting and control of drill.
- Carefully inspect burs both prior to and following each use for signs of excessive wear, fragmentation, eccentricities or other defects. Replace any suspicious bur with a new one prior to use.
- Improper priming of Visao®/ Xcalibur® Hi-Speed Drill will result in excessive handpiece temperature.
- DO NOT modify accessories used with any handpiece.
- Bending or prying may break the blade or bur, causing harm to patient or staff.

SURGICAL PRECAUTIONS

- All reusable components of the device must be cleaned and sterilized and all disposable components replaced before using the device system on another patient.
- Use methods at the operative site (airway, throat, head, neck or sinus cavity) to control bleeding that could compromise patient safety during at-risk surgery.
- DO NOT use any parts other than Medtronic Xomed, Inc. components as damage or substandard operation could result.
- Excessive noise from the bur when drilling close to the cochlea or ossicular chain may cause hearing damage.
- The use of a bur without irrigation may cause an inordinate amount of heat build-up resulting in thermal injury to tissue.
- Excessive force may lead to bur bending, breakage or excess vibration. If the bur fractures during use, the resultant sharp surfaces may cause lacerations of the body tissue. A bur which fractures inside a patient may also cause tissue injury and necessitate emergency procedures for extraction.
- Eccentricity of the bur can cause bur vibration and may result in excess tissue and bone destruction and hearing damage.
- When operating or testing the Xcalibur®, Powerforma® and Visao® drill handpiece, ensure bur is properly inserted into the handpiece. Running the handpiece with an unsecured bur can damage the handpiece. Running a handpiece without a bur is not recommended.
- During procedures, operating room technicians and nurses should always check bur function, fit, tightness, and security in the handpiece before handing off to surgeon for use.
- **IMPORTANT:** Always examine operation of each bur in a handpiece before each use. Operating certain burs at high speeds can cause vibration of the bur. The following are suggestions for reducing bur vibration:
 - Reduce handpiece operating speed.
 - Use only burs that are rated for High Speed operation.
 - Use Xcalibur® Hi-Speed Extended 64 mm Bur Guards (REF # 3333645 or 3333655) with 64 mm burs in the Xcalibur® Hi-Speed Drill.
 - Use Xcalibur® Hi-Speed Extended 75 mm Bur Guards (REF # 3333625 or 3333635) with 75 mm burs in the Xcalibur® Hi-Speed Drill.
 - Use Xcalibur® Hi-Torque Extended Bur Guards (REF # 3333620 or 3333630) with 75 mm burs in the Xcalibur® Hi-Torque Drill.
 - Use Visao Hi-Speed Guards (REF # 3334625 or 3334635) with medium and long burs.
 - Operate handpiece at 50% speed when using burs 75 mm or longer.
 - Select a new bur.

INSTRUCTIONS FOR USE

Note: Medtronic Xomed burs are provided sterile for first time convenience. DO NOT use if package is opened or damaged. Broken sterile seal offers no protection against cross-contamination.

Before Surgery:

For Use with the Powerforma® handpieces: BUR ATTACHMENT

1. Using your thumb, pull back on the bur release button.
2. Insert Medtronic Xomed surgical bur into handpiece, giving bur a 1/4 turn as it is seated (or until it seats fully into the handpiece).
3. Allow the bur release button to slide forward.
4. Tug on bur to ensure bur is secured and locked into the handpiece.
5. **IMPORTANT:** Test drill by depressing foot pedal. The system is now ready to use.
6. To remove the bur, pull back on the bur release button, and pull the bur out of the chuck.

For use with the Xcalibur®/ Visao® handpieces: BUR ATTACHMENT

1. Using thumb and index finger, twist handpiece locking mechanism counter clockwise to unlock handpiece.
2. Insert Medtronic Xomed Surgical Bur into handpiece, giving bur a ¼ turn as it is seated (or until it seats fully into the handpiece).
3. Lock and secure bur by twisting handpiece locking mechanism clockwise.
4. Tug on bur to ensure bur is secured and locked into the handpiece.
5. **IMPORTANT:** Test drill by depressing foot pedal prior to use.
6. To remove the bur, twist the handpiece locking mechanism counter clockwise and pull the bur out of the chuck.

TROUBLESHOOTING

IF Burs wobble in Xcalibur®/Visao®/Powerforma® Handpieces...

- Reduce handpiece operating speeds.
- Use burs that are rated for the console speed selected.

- If necessary, use extended bur guards with Xcalibur®/ Visao® drills for 64 mm and 75 mm burs.
- Operate handpiece at 50% of full speed for 64 mm and 75 mm burs.
- Select a new bur.

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Please contact your Medtronic Customer Service Representative for Warranty information for this product.

CLEANING AND STERILIZATION GUIDELINES

Warnings / Precautions	Remove burs from the handpiece before cleaning and sterilizing. Before sterilization, carefully inspect the bur tips, bur flutes, under a microscope or magnifying glass for any irregularities or eccentricities. Discard any burs that show signs of damage or wear. Cold soak in glutaraldehyde, chlorine, or ammonium solutions, or dry heat sterilization is not recommended as damage to the bur may occur.			
Limitations	Discard any burs that show signs of damage or wear.			
Instructions				
Point of Use	Remove burs from the handpiece before cleaning and sterilizing. Promptly and thoroughly rinse instruments with deionized water after each use.			
Containment and Transportation		No particular requirements		
Preparation for decontamination	Promptly and thoroughly rinse instruments with deionized water after each use.			
Cleaning: Automated (Do NOT use ultrasonic washer)	Remove instruments and equipment from any sterilization trays before placing into washer baskets. Orient devices following recommendations of washer/disinfector manufacturers. Use alkaline or neutral pH detergent recommended by washer/disinfector or detergent manufacturers. These products have been validated for effective cleaning using an automatic washer/disinfector cycle consisting of a minimum 44 minutes total time, including a pre-wash, main wash & rinse, and thermal rinse. The thermal rinse shall be at least 10 minutes long at a minimum temperature of 60°C.			
Cleaning: Manual	Soak in lukewarm*, mild* enzymatic detergent, and deionized water for a minimum of two minutes. Then clean ultrasonically in lukewarm* solution of mild* detergent and deionized water for at least 30 seconds. Rinse thoroughly with deionized water and wipe dry. *less than 43°C; pH 7.0 - 8.5 Note: When using an ultrasonic cleaner or a spray washing machine, follow the manufacturer's recommendations, particularly with regard to articulated instruments and positioning of instruments.			
Disinfection	No particular requirements			
Packaging	No particular requirements			
Sterilization (Temperatures are minimum required, times are minimum required)	Cycle:	Gravity	Gravity	
	Temperature:	121-123°C	132°C	
	Time:	30 min	10 min	
	Cycle:	Pre-Vac:	Pre-Vac: (FR/WHO)	Pre-Vac: (UK)
	Temperature:	132°C	134°C	134°C
	Time:	4 min	18 min	3 min
Drying: 8 minutes, or until visibly dry				
STERRAD Sterilization: Not validated				
100% EtO: Not validated				
Maintenance, Inspection and Testing		Discard any burs that show signs of damage or wear.		
Storage	Store in a clean dry place.			
Additional Information		None		

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the desired result. This normally requires validation and routine monitoring of the process.

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EN

Indications

Sinus indications include septoplasty, removal of septal spurs, polypectomy, antrotomy, ethmoidectomy/sphenoidectomy, frontal sinus trephination and irrigation, frontal sinus drill out, endoscopic DCR, trans-sphenoidal procedures, maxillary sinus polypectomy, circumferential maxillary antrotomy, choanal atresia, sphenoidectomy, and medial, lateral, and posterior frontal sinusotomy. Nasopharyngeal/laryngeal indications include adenoidectomy, tracheal procedures, laryngeal polypectomy, laryngeal lesion debulking, tonsillectomy, tonsillotomy for obstructive tonsillar disease, removal of endobronchial lesions and surgical management of recurrent respiratory papillomatosis (RRP).

Head and neck (ENT) indications include soft tissue shaving, rhinoplasty (narrowing of the bony vault and revision of the bony pyramid), removal and shaping of bone during rhinoplasty procedures, removal of adipose tissue (lipo debridement) in the maxillary and mandibular regions of the face, removal of acoustic neuroma, and incision and removal of soft tissue during plastic, reconstructive, and/or aesthetic surgery.

Device description

The StraightShot Microdebrider System is designed to accommodate a number of Otorhinolaryngology procedures, including Functional Endoscopic Sinus Surgery (FESS), Adenoidectomy, removal of laryngeal and vocal cord lesions, Rhinoplasty, dermabrasion, and submental lipectomy. A variety of disposable blades and burs are available for this purpose. These disposables are for use with all Straightshot and all Magnum Microdebrider handpieces.

Contraindications

None are known.

Warnings

- DO NOT change accessory with handpiece running to prevent laceration of user and cross-contamination through compromised glove.
- For Airway Blades: Should the blade become clogged during use, risk assessment has determined 1-5 cc of irrigant could be aspirated by the patient before detection by the user.
- DO NOT modify accessories used with the handpiece. Performance could be diminished with modified accessories.
- DO NOT use burs above the speed indicated on the bur label.
- For high speed drilling (greater than 6,000 rpm) ensure that an XPS High Speed Bur is used.
- Use adequate irrigation from a separate user-provided irrigating source. The use of an accessory without irrigation may cause an inordinate amount of heat buildup resulting in thermal injury to tissue.
- Disposable devices are for single-use only. Do not attempt to sterilize disposable devices. The disposables are packed sterile and not intended for repeat use. To prevent contamination, use only once.
- Insertion of metal objects in blade or bur tip may cause the blade or bur to break leaving fragments in the wound. The fragments may be difficult to remove, causing irritation, inflammation and foreign-body response at surgical site.
- Always keep the cutting tip of the accessory away from fingers and loose clothing. Prevent laceration of user and cross-contamination through compromised glove.
- Bending or prying may break the blade or bur, causing harm to patient or staff.
- Operate the blade or bur only after the appropriate anatomical landmarks and the intended surgical site have been confirmed.
- Blades, burs and irrigation tubing are disposable and intended for single-use only, UNLESS OTHERWISE MARKED.
- Test for wobble at desired speed prior to use. Discontinue use of accessory if tip begins to wobble and replace accessory to prevent unintended tissue removal from patient.
- Blade and bur accessories are available for resection of soft tissue and bone for surgical procedures. Use of accessories depends on the intended application and patient needs. Sharp-cutting powered accessories induce bleeding and removal of significant tissue and bone.
- Use care in application of the moving cutting end to only appropriate anatomical landmarks and the intended surgical site when using XPS accessories. The use of powered reciprocating/rotating instruments may result in vibration-related injury. Use appropriate precautions.
- Employ visualization when using rotating XPS accessories, including use of image guided surgery system if needed. Discontinue powered application in the event of lack of visualization of the surgical site.
- Use methods at the operative site to control bleeding that do not compromise patient safety during at-risk surgery.
- Always inspect the components before and after use for any damage. If damage is observed, do not use damaged part until it is replaced. Damaged parts may deposit metal shavings on

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Do not use any parts other than Medtronic system components as damage or substandard performance could result.

- Remove and discard accessories following local regulations for proper disposal of contaminated materials.
- Do not use accessory if package is opened or damaged. Broken seal offers no protection against cross-contamination.
- Use lock on handpiece to prevent inadvertent rotation of blade or bur during use.
- It is important that the XPS 3000 operator be familiar with the system User's Guide, its precautions, procedures, and safety issues.
- Properly dispose of single-use devices removed from sterile packages. Devices lose sterility upon removal from packaging.
- Excessive pressure applied to bur may cause bur fracture. Should a bur fracture occur during use, extreme care must be exercised to ensure that all fragments of the bur are retrieved and removed from the patient. Unremoved bur fragments may cause tissue damage to the patient.


Surgical precautions

- Straightshot handpieces are intended to operate at speeds greater than 6,000 rpm ONLY when used with the XPS High Speed bur line.
- XPS Blades should be operated in the oscillate mode only. Operating in the forward mode may cause damage to the blade. XPS burs should be operated in the forward mode only.
- To prevent damage to curved blades and burs, disconnect suction tube prior to changing blade or bur during procedure.
- Be sure the blade is fully engaged in the microdebrider and verify the tip is fully engaged with the outer cannula prior to use.

Instructions for use


Blade or bur installation:

1. Using thumb, depress the locking collar (A) on the front of the handpiece.
2. Insert blade or bur with a slight rotating motion until blade or bur is seated.
3. Align blade or bur tip opening to desired position.
4. Release the locking collar.
5. Pull on the blade or bur to ensure engagement and visually check to make sure distal tip of inner blade is in contact with the distal tip of the outer cannula.

Tubing attachment:

1. Adjust the clamp on the irrigation tubing to the off position and spike the bag of irrigant.
2. Push the free end of the irrigation tubing onto the irrigation port on the blade (B) until tubing passes barb on the port and is secure.
3. Attach suction tubing securely to the suction port on the handpiece (C).
4. Secure the suction and irrigation tubing to the handpiece cable with the clip (if provided) and the white clamps included in the blade packaging.
5. During the procedure, it is recommended to periodically submerge the blade in sterile water with suction connected to the handpiece.

After surgery

Blade or bur removal:

1. If the irrigation pump is not used, turn off the irrigation flow using the roller clamp on the irrigation tubing.
2. Remove and discard the tubing following local regulations for proper disposal of contaminated materials.
3. Depress the locking collar of the handpiece and pull the blade or bur out of the handpiece and discard following local regulations for proper disposal of contaminated materials.

Reuse and cleaning

XPS disposable blades and burs:

- Do Not reuse. Remove and discard following local regulations for proper disposal of contaminated materials.

Sterility

XPS disposable blades and burs:

- Disposable blades, burs, and irrigation tubing are for single-use only. Do not attempt to sterilize disposable blades, burs, or irrigation tubing. They are packaged sterile and are not intended for repeated use. To prevent contamination, use only once.

Troubleshooting

Blade or bur malfunctions:

Issue	Solution
Upon inspection, blade or bur appears to be damaged or defective	• Remove and replace defective blade or bur.
Blade or bur is not firmly seated	• Pull back locking collar and reseal the blade or bur.
Регистрация МИ в Росздравнадзоре Вывод на рынок info@mesvacort.ru	• Remove blade from surgical site and submerge the blade tip in sterile water with suction connected to the handpiece to evacuate the obstruction, or use stylet to clear blade. The suction tube may be removed from the handpiece and inserted over the end of the blade and blade activated to evacuate the obstruction from the tip of the blade.
Blade or bur is leaking irrigant	• Check for proper blade/bur insertion by pulling back locking collar, and reseating blade/bur. • Check blade for obstruction. If visible, evacuate obstruction by one of the above methods. • Inspect suction connection on handpiece and at suction cannister, and if suction tubing is disconnected, connect it. • Remove and inspect suction tubing, and if obstructed, remove obstruction.

Issue	Solution
Blade or bur continues to malfunction	- Remove and replace blade or bur. - Call Customer Service.

Customer service information

For further information regarding the use of this product or to report any problems, please contact Medtronic Xomed using the appropriate information provided on the blue and white contact information card packaged with each device; or contact your local distributor.

Limited warranty

A. This LIMITED WARRANTY provides assurance for the customer who purchases a StraighShot Handpiece Accessory (hereinafter the "Product") that should the Product fail to function to Medtronic Xomed's published specifications during the term of this LIMITED WARRANTY (one year from the date of shipment), Medtronic Xomed will either replace, repair, or issue a credit (adjusted to reflect the age of the Product) for the Product or any portion thereof. This LIMITED WARRANTY is extended only to the buyer purchasing the Product directly from Medtronic Xomed or from its affiliate or its authorized distributor or representative.

B. To qualify for this LIMITED WARRANTY, the following conditions must be met:

The Product must be used on or before its "Use By" or "Use Before" date, if applicable.

The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.

Medtronic Xomed must be notified in writing within thirty (30) days following discovery of a defect.

The Product must be returned to Medtronic Xomed within thirty (30) days of Medtronic Xomed receiving notice as provided for in (3) above.

Upon examination of the Product by Medtronic Xomed, Medtronic Xomed shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic Xomed or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.

C. This LIMITED WARRANTY is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Medtronic Xomed be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the Product, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.

D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

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© 2014 Medtronic, Inc.
EN The following are trademarks or registered trademarks of Medtronic, Inc. in the United States and other countries: StraightShot. All other trademarks, service marks, registered trademarks, or registered service marks are the property of their respective owners in the United States and other countries.
Medtronic B.V.