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

INSULIN DELIVERY DEVICE

Model: XLB01A

Please read the instructions carefully prior to use your pen to ensure proper operation.

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Content

Introduction	3
Intended use	3
I. Insert an Insulin Cartridge	3
Attentions	4
Frequently Asked Questions about Inserting the Cartridge	4
II. Pen Priming	4
Attentions	4
Frequently Asked Questions about Priming	4
III. Injection of Required Dose	6
Attentions	6
Frequently Asked Questions about Injection	6
Maintenance and Storage	6
Get to Know Your GanleePen™	7
1. Preparation	8
2. Insert the Cartridge	8
3. Injection Priming (before every injection)	10
4. Dose Injection	13
Description of Graphic Symbols	15
Replacement	16

Introduction

Intended use:

- GanleePen™ is a reusable insulin delivery device designed for single patient use only, used by diabetes patients intended for the self-injection of insulin. GanleePen™ is designed for ease of use. You can give yourself multiple doses from 1 to 60 units of insulin from one Gan&Lee 3mL insulin cartridge (100 units/mL) in each injection. You can dial your dose in one-unit increment at a time. If you dial the wrong dose, you can correct the dose without wasting any insulin.
- GanleePen™ can be used with injection pen needles and Gan&Lee 3mL cartridges or similar cartridges, subject to the availability of relevant studies confirming this possibility.

Important Note:

- Read and follow these instructions completely and carefully before use. Failure to follow the instructions may result in over-dose or under-dose. If any part of your pen is broken or damaged, STOP USING IT. Contact your supplier or your healthcare professional for replacement.
- Blinds or patients with visual problem can use GanleePen™ with the assistance of a trained person with good eyesight, otherwise, do not use it. When children use GanleePen™ for self-injection, please ensure that the child is over 6 years old and operate under the supervision of an adult. Prepare a spare pen for replacement in case this one is damaged or lost.
- If any abnormality is found during use, stop the injection immediately, and check whether the GanleePen™ is normal, whether the needle is clogged, and whether the needle is tightened well.
- Contraindication: Unknown.

Gangan Medical Technology Jiangsu Co., LTD

1st-4th Floor, Building G130, Fifth Phase of Standard Workshop, China Medical
City, Taizhou, 225300 Jiangsu, P.R. China

After-Sale Service Supplier: Gan & Lee Pharmaceuticals

For information contact: info.marketing@ganlee.com

I. Insert an Insulin Cartridge

Attentions

- Your healthcare professional will prescribe and select the most appropriate insulin injection type and injection dosage. Any change of insulin injection on formation, dosage, etc. shall be made only under medical supervision.
- GanleePen™ is used with Gan&Lee 3mL insulin cartridges (100 U/mL) or similar, subject to the availability of relevant studies confirming this possibility.
- Read cartridge label and ensure the pen contains the right insulin cartridge before each injection.

Frequently Asked Questions about Inserting the Cartridge

1. Why doesn't the leadscrew move out when there is no cartridge in the pen?

The leadscrew will not move out when there is not a cartridge in the pen. When the cartridge is inserted in place, the leadscrew can move out when the injection button is pushed down.

2. What should I do if the cartridge holder can't be attached to the pen body?

Check whether the 3mL insulin cartridge is completely inserted into the cartridge holder. Then carefully align the cartridge holder with the pen body, and tighten them together.

II. Pen Priming

Attentions

- GanleePen™ could only be used together with compatible injection pen needle.
- **Do not share GanleePen™ and needle with others.**
- Needle shall be disposed in accordance with local laws and regulations.
- Always use a new needle for each injection to ensure sterility. It will help to avoid leakage of insulin, keep out of air bubbles and reduce needle plugging possibility.

Frequently Asked Questions about Priming

1. Why is it important to prime before every injection?

If you do NOT prime, insulin dose may be over-prescribed or under-prescribed.


Pen priming is beneficial to ensure that pen and needle are working properly. Once the pen is properly primed, insulin will squirt from the needle tip. You may need to prime several times

before you see insulin squirt.


2. Why injection priming shall be repeated several times when a new cartridge is inserted?

Gaps may exist between leadscrew and cartridge plunger. Repeating the priming steps can move leadscrew out to touch the cartridge plunger. When leadscrew pushes the cartridge plunger out, insulin can squirt from the needle tip.

Before priming

After priming

No gap between leadscrew and cartridge plunger

- Priming for several times until insulin will squirt from the needle tip.
- If you are still unable to see insulin squirt, see Question 3.

3. Why should I prime until I see insulin squirt from the needle tip? And why isn't one insulin drop at the tip of the needle considered as properly priming?

Priming moves the leadscrew to contact with the cartridge plunger and gets the air bubbles out of the cartridge. When you prime:

- You may see a drop of insulin on needle when you operate for the first time. This only means the needle is attached and not blocked. You still need to prime the pen.
- You may not find any insulin squirting, it may be because leadscrew is moving forward to eliminate gaps between leadscrew and cartridge plunger.
- If cartridge contains air, insulin may squirt till air are forced out
- Insulin will squirt only when the pen is properly primed.
- If it is hard to press down the injection button, the needle may be clogged. Attach a new needle repeating the priming steps till insulin squirts. If no insulin squirts, **DO NOT USE** this

pen. Contact your supplier or consult your healthcare professional for help.

4. What should I do if there are air bubbles in cartridge?

Priming your pen can remove air. Hold GanleePen™ with needle point upward; then tap the cartridge gently with your finger so any air bubbles will rise up towards the needle.

Repeat priming steps until insulin comes out of the needle tip. A small air bubble may remain in the cartridge after completion of the priming steps. If you have properly primed the pen, the small air bubble will not affect your insulin dose.

III. Injection of Required Dose

Attentions

- Injection button shall be pressed down so that insulin can be delivered into human body.
- Do not try to inject insulin by rotating dose selector.
- Do not change the dose during injecting.

Frequently Asked Questions about Injection

1. Why is it hard to push the injection button when I try to inject?

- The needle may be clogged. Replace the needle.
- It may be caused by a fast push. It may be easier to press down the button more slowly.
- It will be easier to push down the button by using a larger diameter needle.
- Consult your healthcare professional for the suitable needles.
- The injection button may be hard to push down if it gets polluted by insulin, food, drinks or other substances. If the problem is not solved by aforementioned methods, replace the pen.

2. Why does the dose selector fail to go to "0" during injection?


It may happen because the cartridge does not have enough insulin in it for your entire dose. The number displayed in the dose window is the amount you did not receive. Remove the needle and replace the insulin cartridge, and then attach a new needle for pen priming. Finally, inject the dose that you required.

Maintenance and Storage


- Clean the pen cap, pen body and pen case of the GanleePen™ with a damp cloth.
- Keep the pen away from moisture, dust, extreme temperatures and direct sunlight.
- The allowed operating temperature for GanleePen™ is -20-40 °C. No relative humidity requirements.
- The allowed storage temperature for GanleePen™ is -40-70 °C without the insulin cartridge, and the pen should be stored in a no corrosive gas and well ventilated room.
- Do not use alcohol, hydrogen peroxide or bleach to wipe pen body and dose window. Do NOT immerse in or apply any liquid or lubrication such as oil, as this could damage the pen.
- The cartridge can be inserted in GanleePen™ and stored under the storage conditions recommended by the insulin manufacturer.
- The package of GanleePen™ should be placed right side up during transportation. It should be handled with care. Do not expose to direct sunlight or rain. The maximum number of stacked layers is 8 for the same transport package.

Get to Know Your GanleePen™

Injection Pen Needle (supplied separately)


1. Preparation

Wash hands and prepare the following things before each injection:


- GanleePen™
- Gan&Lee 3mL insulin cartridge
- New needle
- Alcohol swab

2. Insert the Cartridge


Pull off the pen cap

2

Twist off the cartridge holder to detach


3

Push back the leadscrew until it stops


Note: If the leadscrew is out, use the finger to push it back. The leadscrew of your new pen might not be out. Only when the pen is attached the cartridge holder, press down the injection button and the leadscrew will move forward.

4

Take a new Gan&Lee 3mL insulin cartridge. Check the name and color of the cartridge label to ensure that it contains the type of insulin you need. **Check the insulin injection** of the formulation, expiration date, appearance **without cracks or damage**.


5

Insert the cartridge

Hold the cartridge holder and insert the cartridge in as illustrated.

6


Install the cartridge holder

Line up the cartridge holder with the pen body straight, and secure together. Then use alcohol swab to wipe rubber seal.


Note: The leadscrew may not move forward if the cartridge holder is not attached properly and you cannot get your required dose.

7

Roll the pen gently 10 times (only for premixed insulin)


8

Swing gently the pen up and down 10 times to mix insulin evenly (only for premixed insulin)


3. Injection Priming (before every injection)

1

Take a new needle, and tear off the paper


2

Attach the needle

Screw the needle straight tightly onto the pen.


3

Remove the external needle guard


Remove the external needle guard and internal needle guard of the needle. Throw the internal needle guard away. Do not discard the external needle guard, so that it can be used to remove the needle after injection.

4

Rotate the dose selector to 2 units


5

Tap the cartridge

Point the needle upwards. Tap the cartridge gently to collect air at the top for removal.

6

Injection priming

Push injection button until you see insulin squirt from the needle tip. The priming is completed when you see insulin squirt.


The injection priming has not been finished if you only see a few drops of insulin.

If you don't see a stream of insulin squirt, repeat 4-6 steps until you see it. Do not try to inject if no insulin appears.

Warning:

Usually, several injection priming may be required for the first use of the cartridge. If you still don't see insulin squirt from the needle tip after several injection priming, the needle may be clogged, change the needle and prime again. If insulin does not appear after priming repeatedly, do not use the pen. Contact your supplier to get a new pen.

4. Dose Injection

Dial the required dose


Make sure the dose window shows "0" before you start.

Rotate dose selector to choose dose you need.

Reverse dose selector to change dose if overdosed.

For example: 18U shown in the figure.

Note: You can select up to 60 units per injection. Always use dose window to see the selected dose. Do not use the insulin scale on the cartridge holder because it only shows the approximately volume left in your cartridge.

Inject dose you need

Insert needle into your skin as per the injection technique suggested by your healthcare professional.

As shown in the figure, place your thumb on the injection button firstly; then slowly press down the injection button until the button cannot be pushed forward.

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Gangan Medical Technology Jiangsu Co., LTD


When the number shown on dose window is "0", you have injected enough medicine.

Hold the injection button for five seconds and then release the push button. Then pull out needle.

Warning: Please do not press the injection button too fast or too hard, it may lead injection too fast to cause pain feeling for the patient.

Notes:


The set dose can exceed the remaining insulin in cartridge. The number on dose window should be "0" after injection. Otherwise, you have not injected enough insulin. Please remember this number. Remove needle and empty cartridge. Insert a new cartridge as the step 1, and make injection priming as the step 2. Then dial dose selector to your required dose amount. Please do not press the injection button too fast during injecting, which may lead injection too fast to cause pain feeling.

Push the external needle guard back

Push the external needle guard carefully without touching the needle. Never try to put the internal needle guard back on the needle (follow your healthcare professional instructions)

4


Remove the needle

Carefully twist the capped needle off and discard the needle, as directed by healthcare professional.

Notes:

To prevent the air from entering the cartridge, do NOT store the pen with needle attached.

5

Cover pen cap

Description of Graphic Symbols

CAUTION***

CE MARKING

BATCH CODE*

DATE OF MANUFACTURE

STACKING LAYER LIMIT**

MANUFACTURER

FRAGILE, HANDLE WITH CARE**

KEEP AWAY FROM SUNLIGHT**

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No copying without permission

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Page 15 of 16

甘甘医疗科技江苏有限公司
Gangan Medical Technology Jiangsu Co., LTD

THIS SIDE UP**

KEEP DRY**

CONSULT INSTRUCTIONS FOR
USE

TEMPERATURE LIMIT**

Note:

*Symbol used only on the carton.

**Symbol used only on the transportation/shipping box.

***Symbol used only on the instructions for use.

Replacement

GanleePen™ can be used for 3 years after the first use. The maximum number of allowable reuses should be no more than 5,000.

Record the first date your GanleePen™ was used here: ____/____/____.

After 3-year use, consult your healthcare professional or visit local pharmacies for a new GanleePen™.

Date of manufacture: Please see package.

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