PMS Steripack®

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PMS Steripack®

Instructions for Use

PRODUCT NAME: Sterilization Reel and Pouch

PRODUCT DESCRIPTION:

PMSSteripack Sterilization Reel & Pouch is constructed from PET/PP copolymer film and Medical Kraft Paper. A process indicator is applied on the surface of the reel & pouch. Indicators show the method of sterilization.

INSTRUCTIONS FOR USE:

- PMSSteripack** Pouches are ready for use. **PMSSteripack** Reels need to be cut to the required length and sealed at one end to form a pouch.
- Do not use if the pack is damaged. Do not reuse.
- Cover sharp ends of instruments with an instrument protector, to prevent puncturing or cutting.
- Place the item that is to be sterilized into the pouch. Do not fill pouches to more than 75% of their capacity.
- Always leave minimum 1,5 cm gap between the material and the each seal edge of the pouch for a smooth airflow.
- Manually remove as much air as possible from the pouch before sealing.
- Seal the pouch with a calibrated and validated sealer according to the manufacturer's recommendations. Strength of the seal must assure a clean peel.
- Avoid heat seal stress during handling, e.g. fill the pouch horizontally, working on a flat surface.
- When loading into the sterilizer place the pouches into the basket in a way that film faces to film and paper faces to paper.
- Sterilization pouches should be put upright in the basket or container and not too tight together, such that a hand can slide in between them.
- Process the load according to the established procedures
- Handle pouches with care when removing from the sterilizer to prevent them from tearing. Store pouches in an area protected from sunlight, water and other liquids.
- Do not use a wet pouch. Repeat the sterilization process again after changing the pouch.
- To open pouch, use a thumb and index finger to hold plastic side of pouch and a thumb and index finger to pull back paper side of pouch. Always open the pouch from the direction of the arrow. For a clean peel open the pouch slowly.

STORAGE CONDITIONS:

Store at max. +35°C and max. 70% RH humidity. Sterilization reel have a shelf life of 5 years when stored at recommended conditions.

The Pouches, either sterilized or not, must be kept away from moisture, excessive airflow, direct sunlight, once the paper is wet, pouch must not be used anymore.

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- Process the load according to the established procedures
- Handle pouches with care when removing from the sterilizer to prevent them from tearing. Store pouches in an area protected from sunlight, water and other liquids.
- Do not use a wet pouch. Repeat the sterilization process again after changing the pouch.
- To open pouch, use a thumb and index finger to hold plastic side of pouch and a thumb and index finger to pull back paper side of pouch. Always open the pouch from the direction of the arrow. For a clean peel open the pouch slowly.

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PMS is an international producer of sterilization consumables and devices. Since 1997, it has been producing efficient, reliable and flexible solutions for healthcare institutions and medical device manufacturers with its wide product range.

Operating in more than 69 countries in five continents and with production plants in international standards in Turkey, India and sales office in Germany, PMS is a reliable partner with its strategic and global collaborations and more competitive and innovative approach.

For further information: www.pmsmedikal.com

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EN 868-5, ISO 11607, ISO 11140
ISO 9001, ISO 13485

© 2014 PMS. All rights reserved. PMSSteripack, PMSSteriTest, PMSSteriSeal and Pouchmate are registered trademarks of PMS Corp. Other products mentioned within this document are trademarks of their respective owners. Disclosing of the contents, reproduction, adaption or translation of this document without prior written confirmation is prohibited. Because of continued product improvements, we reserve the right to make changes without notice.

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