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1. Images of the device	IV
2. Controls, displays, connectors, and their uses	V
3. Symbols employed	VI
3.1 Symbols on the device	VI
3.2 User interface symbols	VII
4. General information	3
4.1 Device description	3
4.2 Property rights	3
5. Safety instructions	4
5.1 Explanation of warnings and cautions	4
5.2 Normal use	8
5.2.1 Indication	8
5.2.2 Contraindication	8
5.3 Intended use	8
5.4 User qualification	8
5.5 Patient profile	9
5.6 User profile	9
5.7 Intended conditions of use	9
5.7.1 Use	9
5.7.2 Other intended conditions	10
5.8 User position	10
5.9 Safety precautions at the site of installation	10
5.10 Safety features	11
5.10.1 Self-test	11
6. Installation and operating instructions	12
6.1 Unpacking the equipment	12
6.2 Basic equipment	12
6.3 Installation	12
6.3.1 Connecting the ground line	13
6.3.2 Connecting the power cord	13
6.3.3 KARL STORZ-SCB	13
6.3.4 Connecting light cable	14
6.3.5 Connecting the light cable to the endoscope or videoscope	14
6.3.6 TL 005 Triple adaptor (optional)	15
6.4 Commissioning the power LED 300 light source	16
6.4.1 Using the automatic light source control	17
6.4.2 Setting the brightness manually without an SCB connection to the camera control unit	18

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WARNING: The electrical installation of the operating room in which the device is connected and operated must comply with the applicable IEC standards.

WARNING: Keep device out of reach of patients.

WARNING: For safety reasons, do not simultaneously touch the device output sockets and the patient.

WARNING: The instruction manuals and interface components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only assured to be safe if

- `sgdx`qd hcdmsdcd`r rrtbg hlm sgd qdrodbslud`
instruction manuals or
- `sgd hmsdmcdc trd`mc hmsdqe`bd rodbheb`shnmr`
of the devices used in combination permit this
"be- Bk`trd 05 ne sqd 2uc Dc- ne iDB 5/5/0.0(-

WARNING: Grounding reliability can only be achieved when the unit is connected to a properly installed 'Hospital Only' or 'Hospital Grade' receptacle (i.e. approved for use in an operating room). If inspection reveals damage.

WARNING: To avoid risk of electric shock, this equipment must only be connected to supply main with protective earth.

WARNING: At the end of its useful operating life, dispose of the unit as electronic scrap.

WARNING: Do not store liquids on or directly above the device. Do not store liquids on or directly above the device.

If liquid penetrates the device, in spite of the evaporation (this also applies for the formation of

WARNING: Dmrtqd rteebidms `hqbqbtk`shnm enq sgd device.

WARNING: It is not permitted to press more than one touch point simultaneously.

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

[illegible]

WARNING: The compatibility of this device with any accessories and/or energized endotherapy devices should be checked in the individual instructions for use.

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WARNING: The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment.

WARNING: To avoid the risk of electrical shock, do not open the unit. Refer servicing to the manufacturer or to personnel authorized by the manufacturer. Removal of covers by unauthorized personnel will void the equipment's warranty.

WARNING: Danger of glare! Never look into the light outlet when the light system is switched on.

WARNING: If the light source is not properly installed, it can cause an internal build-up of heat which leads to a safety shut-off of the light source.

WARNING: Sgd bnqqdbb bnme f tq`shnm ne sgd cdubh must be checked after each software update.

s5YUlls5YUY2zY9 pūōāUōUōKū Iāūōk ū
 I k̄əūōk ōāUōUōū ōĀūIūĀā ĀUōāāūūĀĀUōk āĀūū
 IūōūĀUōUōū ĀūUōūū ōUō ōūōūŌgUōUōū ō āĀūIūĀā
 āūūĀIūĀāIāŌgūū āĀāUōāĀā 3N5 5/5/0-0

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Trd ne sgd cduibd lm ędikr nsgdq sg`m sgnrd
indicated above is not allowed for safety reasons.

3 WARNING: Mn I nchēb`shnm ne sghr equipment is allowed.

The device is intended to create light for endoscopic diagnostic observation and surgical procedures.

The medical device is contraindicated for ophthalmologic interventions. Mandatory protection implemented for the medical device.

For other interventions the responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

Cold light sources are designed to supply light for endoscopic diagnostic and surgical procedures.

Nmrix pt¹ Wedc nbnsnr² mc l dcb³ k⁴ rrlrs⁵ ms rs⁶ eer⁷
are allowed to use the power LED 300 light source
after training. The training has to be conducted
by skilled persons only. The area agency or the
responsible KARL STORZ sales representative is
available for trainings, for advices about further
training alternatives.

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3

- Adequate and recognized skills in the relevant application
- Adequate powers of comprehension to rationally assess the surgical situation in hand
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- Have been thoroughly trained at least once in the use of the device
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic) and that cannot be compensated using auxiliary means

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
Frequency of use: one or more times a day
Length of use: from a few minutes to several hours a day.
Place of installation: Positioning on a level, vibration-free surface
Mobility: can be moved if positioned on a cart
Combination: can be used on the patient at the same time as other equipment required for the operations
Control: can be controlled via the KARL STORZ-SCB.

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
For operation, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of $\pm 45^\circ$ at a distance of approx. 30 – 70 cm to the front panel.

The device may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.
It is not intended for use in hazardous zones.

This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants. The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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6.3.6 TL 005 Triple adaptor (optional)

The TL 005 triple adaptor (see »Recommended accessories«) can be used for the light cables from the following manufacturers:

- KARL STORZ
- Stryker (position II)
- Olympus (position III)

1 NOTE: Richard Wolf light cables can be used after lowering the metal sleeve on the connection piece at the KARL STORZ position.

3 WARNING: Make sure that the light source is switched off when changing the light bulb.

1. Unscrew the light outlet point ③ (standard adaptor) by turning in a counter-clockwise direction (Fig. 8).
2. Place the triple adaptor onto the screw thread and turn to the left until the adaptor can be felt to click into place (Fig. 9).
3. Then screw in the adaptor hand-tight in a clockwise direction (Fig. 10).

1 NOTE: After being screwed in, it may be the case that the adaptor is positioned asymmetrically inside the device (marking not in the 01 position). This asymmetry has no impact on the use of the adaptor, but should be corrected in order to prevent operating errors.

4. To do so, slightly undo the screw at the edge of the adaptor using a Phillips screwdriver (Fig. 11).
5. Rotate the adaptor (Fig. 12) such that the marking on the surface of the adaptor is in the 12 o'clock position (see arrow).
6. Re-tighten the screw (hand-tight) using the Phillips screwdriver (Fig. 12).
7. Select the required working position (KARL STORZ, II or III) by rotating the inner adaptor component (Fig. 13). This working position should then always be in the 12 o'clock position (see note).
8. Insert the light cable into the corresponding socket, and connect (Fig. 14). See also section 6.3.4 and 6.3.5.

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
Press button  to switch to the Settings menu. Further details can be found in the section 'Settings'.

zǎyǔnǚāu úóùèúǎ □+gǎuǐk əŋǎuóǎu ỹ ŷuòq
ùǎǎǎǎùù- ɸùèùùùǎuùùùkũ ǎĩɸùùùù əǎĩĩùðqǎǎq
ĩ òùǎĩù ǎzǎǎǎǎùùùùù-

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Press the Start/Stop button  to switch on the light (the symbol on the Start/Stop button  is orange) (Fig. 20). The light intensity is then controlled automatically via the #L@FD0R + ensuring optimal illumination in every situation.

1 NOTE: The automatic light regulation is deactivated as soon as the MANUAL  button is pressed during operation. The automatic light source control can be switched on again via the camera control unit. Further information can be found in the

[illegible][illegible]

* Depending on the software version of the IMAGE1 S™, the images may differ slightly from the actual appearance.

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1 NOTE: The orange triangles indicate the set setpoint value in relation to the possible setting range.

1 ||ΣΤΖΤΖ zY| 5Aā0ūūPīTκV AĀ0ā0ūū0q0ūū0ū
 ūīTūū0ū0āgāQĀ 0āāĀĀūū000ūū 0āōā00ūū
 00āgū0ūū ūĀ0ūūāāĀ00q0ūū Tūū0ūūūū0ūū0ū
 ḡūūāōā0ūū0ā 0āāĀĀūū0ūū-

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1 NOTE: By pressing the + and – buttons, the setpoint value of the light intensity on the device is set. The device switches the preselected setpoint value at the light outlet only after releasing the button. The intensity

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7.1 Fuse replacement

3

WARNING: The line voltage to the device is not to be unplugged.
power plug has been unplugged.

Switch off the device and remove the power plug from the electrical outlet (Fig. 30).

Remove the fuse holder 7 with a screwdriver or other suitable tool (Fig. 31).

3

WARNING: The fuse must fulfill the requirements for a component of high integrity characteristics.

Insert new fuses (Fig. 32).

	100...240 V~
Line fuse	2 x T 4.0 AH 250 V

Reconnect the power cord.
Test the device for proper operation (Fig. 33).

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3

3

1 NOTE: The manual “Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments” can be downloaded or requested by visiting www.karlstorz.com.

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7.3.2 Safety checks/repeat inspections

Visual inspection

9. Check the device and accessories for any mechanical damage which may impair functionality.
10. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the current standard.
- Patient leakage current measured according to IEC 62353: The limit values can be found in the current standard.

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Test for proper functioning

Carry out the functional check as per the instruction manual (see section 6.4).

Documentation

The safety check/repeat inspection and results must be documented.

- 3 WARNING:** If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been corrected or the device has been serviced.
- 1 NOTE:** For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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7.4 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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7.5 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

- 3 WARNING:** At the end of its useful operating life, dispose of the device as electronic scrap.
- Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.
- Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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8.1 Information signals

8.1.1 Configuration

All of the information signals created, with the exception of those which display the cause of an inoperable device ('Malfunction', see e.g. page 16 Self test), are not self-retaining. This means that the signal is only issued whilst the signal conditions are present.

If the signal conditions are only very brief, there are rules as to the minimum visual/acoustic signal duration, with which the described information signals also comply. This is explained in more detail in the following chapters.

Signals are issued depending on priority. Any low priority signal will be overwritten by a signal with a higher priority, or a low priority information signal is suppressed for as long as there is a signal with a higher priority. If numerous signal conditions with the same priority occur, only the condition most recently detected is signaled.

8.1.2 Visual signals

Information signals are displayed inversely (blue font on a white background) compared to normal `chrol`xr`mc cn mns 0`rg-`

The visual signal is issued for as long as the signal conditions apply. In the case of signal conditions of ``udqx rgnqs c tq`sinm+ ls l`x ad cleqbtis enq sgd t r dq` to recognize which condition triggered the signal. For this reason, the minimum display duration is set at 5s.

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8.1.3 Acoustic signals

The acoustic signals described below satisfy the requirements in terms of distinctiveness.

The acoustic signal is issued for as long as the signal conditions apply. But at least one complete signal sequence is played.

The acoustic signal can be switched on or off in the Settings menu.

But the acoustic signal can only be switched off temporarily by the user. When the line voltage is switched on, the setting is always 'Audio On'.

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


8.1.4 Information signals (general)


The pitch is modulated by ± 2 Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated. In the case of messages which are related to the Safe State, the signal sequence is repeated every 15s, but otherwise it is just played once.

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
On completion of the self-test after switching on, an availability signal sounds. This signal, which is issued following a successful self-test, has the following characteristics:

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- Complete failure of the device.

Possible causes:

- Power supply failure.
- Defective power fuse.
- Power cord connector is not properly connected to device socket.

Remedy:

- Check that there is electricity to the wall outlet.
 - Change fuses as described in the instruction manual. Use only fuses of the correct rating.
- ± O t r g o n v d q b n q c b n m m d b s n q e q l x l m s n c d u i b d socket.

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On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of technical data to personnel duly authorized by the manufacturer for their use in repairing those components of the device that have been designated by their respective manufacturers as repairable.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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9.2 Recommended accessories

Item	Cat. no.
SCB connecting cable , length 100 cm	20 0901 70
Triple adaptor , for use with Power LED 300 cold light fountain in conjunction with KARL STORZ, Olympus, Stryker and Wolf light cables.	TL 005
VE adaptor for halogen light sources	487 VE
All KARL STORZ light cables are suitable for use	(see catalog; corresponding specialities)

9.3 Recommended accessories

1 NOTE: The degree of protection against electric shock at the applied part as type BE `bbnqcmf sn `DB 5/5/0,0 `cdeaqk `snrm, oqnsdbdc(ir ft `q`msddc trlmf sgd KARL STORZ light cables listed below. 

Item	Cat. no.
Fiber optic light cable , length 230 cm, dia. 3.5 mm	495 NA
Same , length 230 cm, dia. 3.5 mm	495 NAS
Same , length 300 cm, dia. 3.5 mm	495 ND

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Table 1 – Compliance level for Immunity tests Guidance and manufacturer's declaration – electromagnetic immunity			
Sgd onvdq KDC 2// hr hmsdmcdc enq trd hm sgd dldbsqn l `fmdsib dmulqnm l dms rodbiedc adlnv- Sgd trdq ne onvdq KDC 2// rgntic l `jd rtdq sg`s hs hr trdc in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If ònnqr `qd bnudqdc vlsq rxmsgdsib l `sdq`k sgd qdk`siud humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short hmsdqqtosinmr `mc òtbtst `sinmr of the supply voltage acc. to IEC 61000-4-11	Voltage Dip: Dip to 0 % for 1 cycles @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100 % for 250/300 cycles	Voltage Dip: Dip to 0 % for 1 cycles @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the power LED 300 require continued operation during power mains interruptions, it is recommended that the power LED 300 be powered from an uninterruptible power supply or a battery.
L `fmdsib edic `s sgd r tooix frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the power LED 300 further from sources of onvdq eqdptdmbx l `fmdsib edicr nq sn hmrs`k l `fmdsib rgldcimf- Sgd onvdq eqdptdmbx l `fmdsib edic rgntic be measured in the intended installation location to `rrtdq sg`s hs hr r teebidmsix in v-
Radiated, radio-frequency, dldbsqn l `fmdsib edic l l t msx test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Oqnrwl l tsx QE edic sdrs ldudlr	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced ax q`cin,eqdptdmbx edicr acc. to IEC 61000-4-6	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	6 V _{rms} on 150 kHz to 80 MHz 1 kHz 80 % AM modulation 6 V _{rms} in ISM bands	

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Table 3 – Test levels for Radiated and conducted Immunity Tests
Guidance and manufacturer's declaration – electromagnetic immunity

Sgd onvdq KDC 2// hr hmsdmcdc enq trd hm sgd dldbsqn I `fmdslb dmulqnm I dms rodbiedc adlnv- Sgd trdq ne onvdq KDC 2// rgntlc I `jd rtqd sg`s hs hr trdc in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	6 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the power LED 300, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: f where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Eidic rsqdmfsg r eqn I ewdc QE sq`mr I lssdqr+ `r cdsdq I hmdc ax `m electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b . $f \quad 80 \text{ MHz to } 800 \text{ MHz}$ $h \quad 800 \text{ MHz to } 2.5 \text{ GHz}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	2 U. I	Interference may occur in the vicinity of equipment marked with the following symbol: 7

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Mnsd9 Sgd rd f tldclmdr I `x mns `oolx hm `lk rist `sinmr- Sgd oqno` f`sinm ne dldbsqn I `fmdslb pt`mslsidr hr `eedbsdc ax `arnqosinmr `mc qd0dbsinmr ne atlkclmfr+ objects, and persons.

a Eidic rsqdmfsg r eqn I ewdc sq`mr I lssdqr+ rtbg `r a`rd rs`sinmr enq q`cln `bdltk`q.bngcldr r(sdldognm dr `mc k`mc I naikd q`cln+ `I `sdtq q`cln+ @L `mc EL q`cln aqn`cb`rs `mc SU aqn`cb`rs b`mms ad oqdcibsd csgdnqdsib`lxx vlsq `bbtq`bx- Sn `rrdr rsgd dldbsqn I `fmdslb dmulqnm I dms ctd sn ewdc QE sq`mr I lssdqr+ `m dldbsqn I `fmdslb risd rtqudx rgntlc ad bnmrlcdqdc- He sgd I d`rtqdc edic rsqdmfsg hm sgd lnb`sinm hm vgibg sgd power LED 300 is used exceeds the applicable RF compliance level above, the power LED 300 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the power LED 300.

b Nudq sgd eqdptdmx q`mfd 04/ jGy sn 7/ LGy+ edic rsqdmfsg r rgntlc ad ldr sg`m 2 U. I -

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
Table 4 – Emission class and group		
Guidelines and manufacturer's declaration – electromagnetic emissions		
Sgd onvdq KDC 2// hr hmsdmcdc enq t rd hm sgd dldbsqn l ` fmdslb dmlqnm l dms rodblędc adknv- Sgd bt rsn l dq nq t rdq ne sgd onvdq KDC 2// rgntk dmr tqd that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The power LED 300 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class B	The power LED 300 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Unis` fd Ôtbst`sinmr.Ôibjdq D l hrrlnmr `r odq IEC 61000-3-3	Complies	

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
Table 5			
Recommended separation distances between portable and mobile RF communications equipment and the power LED 300			
The power LED 300 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the power LED 300 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the power LED 300 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz f	80 MHz to 800 MHz f	800 MHz to 2.5 GHz h
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies. Mnsd9 Sgd9d f t1cd9mdr l `x mns `oolx lm `kk rlst`sinmr- Sgd oqno` f`sinm ne dldbsqn l `fmdstb pt`mslsidr lr `aedbsdc ax `arngosinmr `mc qd0dbsinmr ne rsqtbstqdr+ objects and people.			

This equipment has been evaluated for compatibility with high-frequency surgical equipment. It has been tested according to IEC 60601-2-2 Annex BB.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlántic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada, Del. Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamveien1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Sörsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestučio st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-it-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B.V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N.V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S.A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G/1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S.A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: info-pl@karlstorz.com

KARL STORZ Endoskop Magyarország Kft.
Toberek utca 2. fsz. 17/b
HU-1112 Budapest, Hungary
Phone: +36 195 096 31, Fax: +36 195 096 31
E-Mail: info-hu@karlstorz.com

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: info-ro@karlstorz.com

KARL STORZ Endoskope Greece M.E.P.E.*
Patriarhou Grigoriou E' 34
54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
E-Mail: info-gr@karlstorz.com

*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9000, Fax: +90 216 442 9030

**Sales for Industrial Endoscopy

000 KARL STORZ Endoscopy – WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: Info-ru@karlstorz.com

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20
E-Mail: marketing@karlstorz.com.ua

KARL STORZ SE & Co. KG
Representation Office
Sabir Orudschow 1184, apt. 23
1025 Baku, Azerbaijan
Phone: +99 450 613 30 60
E-Mail: info-az@karlstorz.com

KARL STORZ ENDOSKOPE –
East Mediterranean and Gulf (Offshore) S.A.L.
Spark Tower 1st floor
Charles Helou St., Horch Tabet – Sin El Fil
Beirut, Lebanon
Phone: +961 1 501105, Fax: +961 1 501950
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
Saryarka, 6, BC "Arman", off. 910
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: infoasia@karlstorz.in

KARL STORZ SE & CO. KG
Interchange 21 Tower, Level 33
399 Sukhumvit Road,
North Klongtoey, Wattana,
10110 Bangkok, Thailand
Phone: +84 28 3823 8000
Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
200030, Shanghai, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

www.karlstorz.com

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com