
INSTRUCTION MANUAL

ULTRASONIC A/B SCANNER
AND PACHYMETER

UD-800

Read this manual thoroughly before using the instrument to ensure proper and safe operation.

Contact Tomey Corporation or your local distributor if you have any questions or you encounter any problems during operation.

Note

- *Always follow the operation procedures described in this manual.*
- *Keep this manual in a readily accessible location while operating the instrument.*
- *Contact your local distributor if you lose this manual.*

76AA9090-0L

i. Important safety information

- ***Do not install this instrument in a location where explosives or flammable substances are used or stored. Otherwise, fire or explosion may occur.***
- ***Do not remove the cover of the instrument. Otherwise electric shock may result.***
- ***Do not disassemble or modify the instrument. Otherwise electric shock may result.***

- ***Disconnect the power cord from all the connected devices before installing and/or servicing the instrument. Otherwise, electric shock may result.***
- ***Sterilize probes and the eye cup before use.***
- ***Do not use probes and the eye cup with its contact section damaged. Measurement error may occur, as well as the cornea or eyelid may be injured.***
- ***Do not commence service or maintenance work while the instrument is being used for a patient.***

- ***Do not place water or chemicals on the instrument. Any water or chemicals entering the instrument may cause electric shock or failure.***
- ***Do not allow water or chemicals to splash on probes, except for the tip cap. Any water or chemicals entering the instrument may cause electric shock or failure.***
- ***Use only the specified data terminal for connection of the instrument. Using another type of terminal may result in failure of the instrument.***
- ***This instrument is a diagnostic/measuring device specially designed for ophthalmology. Never use the instrument for other purposes.***
- ***The terminal for connecting the instrument to external devices is not isolated from the internal circuit. Inappropriate wiring may damage the internal circuit. Never touch any of these terminals and the patient at the same time. Be sure to contact Tomey or your local distributor before using the external output terminal.***
- ***When operating this instrument connected to other devices, use only devices that satisfy IEC60601-1 or equivalent safety requirements, or that conform to IEC60950 and whose power source is isolated with an isolation transformer, in order to satisfy the relevant safety requirements for medical equipment.***

This page is intentionally left blank.

ii. How to read this manual

Outline

This manual is structured as follows.

1. PRIOR TO USE

Describes safety precautions and important information to be understood before installing and using the instrument.

2. NAMES AND FUNCTIONS

Describes names and functions of each section of the instrument.

3. OPERATION PROCEDURES

Describes information required for installing and using the instrument.

4. TECHNICAL INFORMATION

Describes useful technical information about the instrument.

5. INSPECTION AND MAINTENANCE

Describes procedures for replacing consumable parts, etc. that the user of the instrument should normally conduct.

6. TROUBLESHOOTING

Describes how to solve problems.

7. CONSUMABLES AND OPTIONAL EQUIPMENT

Describes consumable parts and optional equipment.

8. SPECIFICATIONS

Describes the specifications of the instrument.

Symbols used in this manual

Sentences accompanied with the symbols below indicate the following:

- *This is a precaution that, if unheeded, will result in a hazardous situation where there is an imminent danger of serious injury or death.*

- *This is a precaution that, if unheeded, could result in a hazardous situation where there is a possibility of serious injury or death.*

- *This is a precaution that, if unheeded, may result in a situation where there is a possibility of minor or moderate injury or damage to property.*

- *This is an additional instruction which may contain a special precaution on company policy related, either directly or indirectly, to the safety of personnel or to the protection of property.*

iii. Contents

i. Important safety information.....	i-1
ii. How to read this manual.....	ii-1
Outline.....	ii-1
Symbols used in this manual.....	ii-2
iii. Contents	iii-1
1. PRIOR TO USE	1-1
1.1 Precautions for operation	1-1
1.2 Checking package contents	1-5
1.3 Glossary.....	1-6
1.4 Overview	1-9
2. NAMES AND FUNCTIONS.....	2-1
2.1 Front of the main unit.....	2-1
2.2 Back of the main unit	2-2
2.3 Sides of the main unit	2-3
2.4 Screen.....	2-4
2.4.1 Common items.....	2-4
2.4.2 B mode image diagnosis screen	2-5
2.4.3 Axial length measurement/view screen.....	2-7
2.4.4 Pachymetry screen.....	2-9
2.4.5 A-Diag measurement screen.....	2-10
2.5 Symbols used for marking	2-12
3. OPERATION PROCEDURES.....	3-1
3.1 Preparation before use	3-1
3.1.1 Connections.....	3-1
3.1.2 Turning the power on and adjustment after turning the power on	3-5
3.1.3 Switching modes.....	3-5
3.1.4 Entering the patient data	3-6
3.1.5 Deleting all measurement data (measurement preparation for other patient)	3-9
3.1.6 Selecting an eye	3-9
3.2 B mode image diagnosis	3-10
3.2.1 Settings related to measurement	3-10
3.2.2 Displaying probe angle	3-10
3.2.3 Gain adjustment.....	3-12
3.2.4 Patient preparation	3-13

3.2.5 Obtaining ultrasonic image	3-14
3.2.6 Various functions	3-15
3.3 Biometry function	3-22
3.3.1 Setting measurement conditions	3-22
3.3.2 Operation check.....	3-25
3.3.3 Preparation for measurement.....	3-25
3.3.4 Precautions for measurement	3-26
3.3.5 Measurement:.....	3-32
3.3.6 Checking waveforms after measurement.....	3-35
3.4 IOL power calculation	3-38
3.4.1 Calculation	3-38
3.4.2 Entering calculation parameters.....	3-38
3.4.3 Setting calculation formula	3-43
3.4.4 Entering data after surgery.....	3-44
3.5 Statistical processing	3-45
3.6 Pachymetry function	3-47
3.6.1 Setting the data type to be displayed	3-47
3.6.2 Setting measurement conditions	3-47
3.6.3 Displaying and setting measurement points	3-49
3.6.4 Operation check.....	3-50
3.6.5 Preparation for measurement.....	3-50
3.6.6 Measurement:.....	3-51
3.6.7 Checking measurement data	3-52
3.7 A mode diagnosis.....	3-56
3.7.1 Connecting the probe	3-56
3.7.2 Switching to A mode diagnosis function.....	3-56
3.7.3 Registering Ref. Data	3-56
3.7.4 Measurement:.....	3-57
3.7.5 Setting gain.....	3-57
3.7.6 Switching data display	3-58
3.7.7 Setting the beam direction.....	3-58
3.7.8 Gradation characteristics function	3-59
3.7.9 Analysis functions	3-59
3.8 Export, print, and save.....	3-61
3.8.1 Export.....	3-61
3.8.2 Print.....	3-63
3.8.3 Saved data management	3-70
3.9 System setup	3-73

3.9.1 General	3-74
3.9.2 Measurement.....	3-76
3.9.3 Application	3-81
3.9.4 Connecting & Printing.....	3-83
4. TECHNICAL INFORMATION	4-1
4.1 IOL power calculation formula	4-1
4.1.1 Haigis optimized / Haigis Standard	4-1
4.1.2 Hoffer® Q formula	4-3
4.1.3 Holladay 1	4-5
4.1.4 SRK-II formula	4-7
4.1.5 SRK/T formula	4-9
4.1.6 SRK SHOWA formula.....	4-11
4.1.7 Shammas-PL formula.....	4-12
4.1.8 SRK/T Double K	4-14
4.2 How to calculate axial length using the biometry function	4-16
4.2.1 Phakic eye	4-16
4.2.2 Dense cataract eye.....	4-16
4.2.3 Aphakic eye	4-16
4.2.4 IOL eye	4-16
4.3 Amplifier characteristics of the A-scan diagnosis	4-18
4.4 Ultrasonic wave output	4-18
5. INSPECTION AND MAINTENANCE	5-1
5.1 Warranty	5-1
5.2 Operation life	5-2
5.3 Inspection	5-2
5.4 Routine maintenance.....	5-2
5.4.1 Maintenance of measurement probes.....	5-2
5.4.2 Cleaning and Disinfection of measurement probes (for Europe).....	5-4
5.4.3 Maintenance of main unit	5-6
5.5 Replacing consumables	5-7
5.5.1 Fuses	5-7
5.5.2 Printer paper	5-7
5.6 Storing.....	5-8
5.7 Disposal	5-8
6. TROUBLESHOOTING	6-1
6.1 Common items.....	6-1
6.2 B mode image diagnosis	6-3

6.3 Biometry function	6-5
6.4 Pachymetry function	6-8
7. CONSUMABLES AND OPTIONAL EQUIPMENT	7-1
7.1 Spare parts	7-1
7.2 Optional parts	7-1
8. SPECIFICATIONS	8-1
8.1 Specifications	8-1
8.2 Quantity of ultrasonic energy	8-5
8.3 Noise	8-7
8.4 Operating environment	8-7
8.5 Classification of ME equipment	8-7
8.6 Declaration of conformity to EMC	8-9

1. PRIOR TO USE

- *Read this manual thoroughly before using this instrument to ensure proper and safe operation.*
- *Always follow the operation procedures described in this manual.*
- *Check that there are no devices that generate a strong magnetic field near the instrument. A strong magnetic field may cause noise and affect measurement.*
- *Do not place any other objects on this instrument.*

1.1 Precautions for operation

- *Do not use the instrument to the patients who has ophthalmic disease or systemic disease which will have aggravation by contact to cornea, conjunctiva or eyelid.*
- *Allow only qualified operators to use the instrument.*
- *When using the IOL calculation result to select intraocular lenses, also thoroughly examine cataract surgery methods and other inspections to determine the selection.*
- *The ultrasound images may contain artifacts. If the image is doubtful, consider other examinations to carefully determine whether to adopt the value.*
- *Precautions when installing the instrument and accessories*

- *Install the instrument in a location free of water or chemicals. Any water or chemicals entering the instrument or devices may cause electric shock or failure.*
- *Do not install the instrument in a location where chemicals are stored or gases may occur. Any chemicals or vapor entering the devices may result in fire.*
- *Check the frequency, voltage, and allowable current (or power consumption) of the power source. Otherwise, fire or electric shock may occur.*
- *Connect the power plug to a grounded 3-pin outlet. Otherwise, a short circuit due to failure of the instrument may result in electric shock.*
- *Do not place any heavy object on the power cord or squash the power cord. Fire or electric shock may occur.*

-
- *Fully insert the power plug into the outlet. Faulty contact, allowing any metal to contact the exposed terminal of the plug, or dust accumulated on the exposed terminal of the plug may result in fire or electric shock.*
 - *Conduct grounding work correctly. Otherwise, you may get electric shock.*
 - *Do not connect a device with data transmission specifications that are not compatible. Fire or electric shock may occur. Contact Tomey Corporation or your local distributor before using the instrument while connected to another device.*
 - *If there is a problem with the instrument, turn off the power switch and disconnect the power plug. Install the instrument in a place where this can be performed smoothly.*
 - *Install the instrument in a location not subject to direct sunlight, high temperature and humidity, or air containing dust, salt, and/or sulfur. Otherwise, failure or malfunction may occur.*
 - *Install the instrument in a leveled stable location free of vibration or mechanical impact. Otherwise, measurement cannot be conducted correctly. Also, the instrument may topple or fall, resulting in fire or a serious accident.*
 - *Install the instrument in a location with ample clearance from other devices to allow smooth inspection.*

■ **Precautions during operation**

- *Do not allow water or chemicals to splash onto the B mode probe, except for the tip cap. Any water or chemicals entering the instrument may cause electric shock or failure.*

■ **Precautions before operation**

- *Check the electrical contact of switches, polarity, dial setting, and meters, and that the instrument is working correctly.*
- *Check that all cables are connected correctly.*
- *Since simultaneous use of multiple devices can cause misdiagnosis or result in a hazardous situation, exercise caution when using this instrument.*
- *Check the sections that the patient will directly touch.*
- *Check that the instrument is correctly grounded.*
- *Check that the date set in the instrument conforms to the actual operation date and time.*

■ Precautions during operation

- Do not place any container with liquid in it on the instrument. Any liquid entering the instrument may cause electric shock or failure.
- Do not allow the patient to touch the terminals for connecting the instrument to external devices.

- Do not lean on the instrument or press the instrument from the top. The instrument may topple over, resulting in mechanical failure or injuries.
- Be extremely careful not to take too much time or too many measurements, as this can stress the patient.
- Constantly observe that both the instrument and patient are free of problems.
- If a problem with the instrument or the patient occurs, take appropriate action, such as stopping the device, to ensure the patient's safety.
- Do not allow the patient to touch the instrument.
- If any smoke, offensive odor, or abnormal sound occurs, turn off the instrument immediately, disconnect the power plug from the outlet, and contact your local distributor or Tomey Corporation.

■ Precautions after operation

- Do not place any container with liquid in it on the instrument. Any liquid entering the instrument may cause electric shock or failure.
- Do not use organic solvents such as thinner, benzene, or acetone to clean the instrument. Fire or electric shock may occur. (These solvents may also corrode the resin or coating on the cover of the instrument.)
- Follow the specified procedures to return the operation switch, dial, etc. to their original positions and turn the instrument off.
- Hold the plug when disconnecting the power plug from the outlet to avoid placing excessive force on the cord. Pulling the cord may damage the inner core wires, resulting in electric shock or fire.
- When disconnecting cords, do not apply too much force to them, for example, do not hold and pull the cord.
- Refer to "5.6 Storing" for instructions on the storage of the instrument.
- Clean the instrument at the end of operation in preparation for the next use.

-
- *Clean and neatly arrange the accessories and cables.*

■ ***If any failure occurs in the instrument, immediately stop operation, indicate the failure in the instrument, and contact your local distributor for repairs.***

- *Do not modify the instrument. Doing so may cause electric shock or failure of the instrument. There is a high-voltage section in the instrument. Touching this section will result in death or serious injuries.*

- *Disconnect the power cord from the outlet when replacing fuses. Otherwise, you may get electric shock, resulting in death or serious injuries.*

- *Use the power cord and fuses provided with the instrument or specified by Tomey to ensure safety. Also, do not use the accessories provided with the instrument for other equipment.*

- *Conduct regular inspections of the instrument and components.*
- *When any instrument failure occurs, indicate the failure, and contact your local distributor to request inspection and repair. Do not attempt to repair the instrument yourself.*
- *When the instrument has not been in use for 1 month or longer, check that the instrument is operating correctly and safely before starting operation. Refer to "5.3 Inspection" in this manual for the relevant procedures.*

1.2 Checking package contents

Open the package and check that the required quantity of the following items is included and they are not damaged. If any item is missing or damaged, contact your local distributor as soon as possible.

Note ■ *Keep the box and packing materials for use when moving or transporting the instrument.*

■ *When taking the instrument out of the box, pull the outer box upward and then remove the packing materials. Be careful not to lift the instrument up by directly holding the head, chin rest, forehead pad, joystick, or cables. The instrument may be damaged.*

- Main unit
- 10MHz B probe 1
- Power code 1
- Foot switch 1
- Printer paper 1
- Ultrasound gel 1
- Fuse 4
- (Two are included as accessories)
- Dust cover 1
- Instruction manual (this manual) 1
- DATA Transfer installation CD 1
- DATA Transfer startup guide 1

1.3 Glossary

[ACD]	: Anterior chamber measurement
[Acrylic]	: Acrylic lens
[Aphakic]	: Aphakic eye
[AVG]	: Average value
[Axial]	: Axial length measurement
[B-Mode]	: B mode image diagnosis
[Calibration]	: Calibration of B mode probes and pachymetry probes
[CCT]	: Corneal center thickness
[Contact]	: Method of taking measurements by directly applying the contact section of the biometry probe to the center of the cornea
[Clinical History Method]	: Method of calculating the corneal refractive power before refractive correction surgery. This is used together with the Double K Method.
[D]	: Diopter. The unit of refractive power representing the level of myopia, hyperopia, or astigmatism. Reciprocal of the focal length measured in meters.
[DATA Transfer]	: System to output the examination data created by TOMEY products as a file
[Export]	: To send measurement data to DATA Transfer or Tomey Link.
[FREEZE]	: Status where ultrasound waves are not transmitted
[Immersion]	: Measurement mode when applying cornea protective agent between the contact section of the axial length probe and the cornea.
[IOL]	: Intraocular lens
[K1]	: Radius of corneal curvature of minor meridian [mm or D]
[K2]	: Radius of corneal curvature of major meridian [mm or D]
[Obj.]	: Section to be measured
[Pachy]	: Measurement of corneal thickness
[Pre]	: Before refractive power correction surgery
[Phakic]	: Phakic eye

[PL]	: Abbreviation of POST LASIK
[PMMA]	: PMMA lens
[Post]	: After refractive power correction surgery
[Ref.]	: Reference section
[Silicone]	: Silicone lens
[SNR]	: S/N (signal to noise) ratio A larger value represents less noise, meaning higher quality signals can be obtained.
[TOMEY Link]	: Digital medical record system to manage data measured with Tomey products
[VD] (Vertex distance)	: Distance between corneal vertexes [mm] Represents the distance between the corneal vertex and the posterior surface of the lenses to be prescribed. When (VD=CL), vertex distance is calculated as 0 mm.
[Auto Power Off function]	: The LCD turns off automatically when no operation is conducted for a specified time. (Auto power off mode) Touch any button to return to Normal mode.
[Caliper]	: Manually move the measurement point.
[Gray scale]	: Indicates the range of gradation values to be measured in images represented in shades from white to black.
[Gain]	: Adjusts amplification, brightness of images, and amplitude of waveforms.
[Gate]	: Specifies the range of waveforms to be detected.
[Smoothing]	: This function smoothes the continuation of images and reduces noises, improving the image quality.
[Measurement point]	: Indicates the point where corneal thickness is measured, specified by diameter and angle.
[Touch panel]	: Allows you to make various settings, etc. by directly touching the monitor.
[Harmonic]	: Image composition method for ultrasound diagnosis which uses harmonic components of ultrasonic waves.
[Bias value]	: Value obtained by converting actual measurements using a specified percentage or correction value

[Vector A mode]	: Indicates the direction of ultrasonic wave transmission and displays the A mode waveform.
[Real time]	: Status where ultrasound waves are transmitted

1.4 Overview

This instrument is designed as an ophthalmology device with functions to retrieve ultrasound topographic images of living tissues using ultrasound waves generated from an ultrasound transducer built into the probe, to retrieve A mode waveforms and to measure axial length and corneal thickness.

B mode image diagnosis function

- From the B mode probe connected to the main unit, images and displays the tissue of the scanned eyeball cross section generating and scanning ultrasound waves into the eyeball and using ultrasonic echoes reflected from each tissue inside the eyeball.

Biometry function

- The measurement probe connected to this instrument transmits ultrasound waves to the inside of an eyeball and then receives ultrasound echoes reflected by the tissue boundary surface in the eyeball.
- The time required for the ultrasound wave to be transmitted from the front of the cornea to the measurement section is measured, and the length of the living tissue is calculated based on the preset converted acoustic velocity.

L: Length of living tissue

V: Converted sound speed

t: Measurement time

- With an IOL power calculation function, this instrument has various IOL power formulae applicable to ordinary cataract surgery and/or cataract surgery of eyes with corrected corneal refractive power, and provides data necessary to determine IOL power.

Pachymetry function

- Function to measure corneal thickness in the same manner as measurement of axial length


A mode diagnosis

- Function to identify areas with pathological changes by observing A mode waveforms

This page is intentionally left blank.


2. NAMES AND FUNCTIONS

2.1 Front of the main unit

- (1) LCD and touch panel
- (2) Foot switch
- (3) Power LED
- (4) Power switch
- (5) Biometry probe connector / Pachymetry probe connector
- (6) Probe
- (7) Probe holder

2.2 Back of the main unit

(1) Maintenance switch

Our service personnel use this switch for maintenance. Never touch this switch.

(2) LAN connector

(3) USB-D connector (PC) Connect a personal computer here.

(4) Foot switch connector

(5) Inlet


(6) Fuse holder

(7) Connector for external storage media such as USB flash memory, a video printer, and an external ID input device, etc.

(8) Connector for fixlight

Connect the external fixlight of the special chin rest here.


2.3 Sides of the main unit

(1) Built-in printer

2.4 Screen


2.4.1 Common items


- (1) Eye button (eye selection button)
- (2) Patient information button Moves to the patient information input screen.
- (3) Patient information field
- (4) Version Displays the version of the system.
- (5) "Setup" button Moves to the Setup screen.
- (6) Print / Export / Save buttons
- (7) Database button Moves to the Patient List screen.
- (8) Touch and hold this button to delete the examination data of both eyes and the patient's information.

2.4.2 B mode image diagnosis screen

(Real-time screen)

(FREEZE screen)


- (1) Image monitor Displays an image and the following information. Patient ID, inspection date/time, eye to be examined (R/L), image page No., converted acoustic velocity Total gain (TG) / Dynamic range (DR) / Near gain (NG) Settings for smoothing, harmonic, gradation characteristics, frequencies and scope of the image<
- (2) Buttons for playing movie
- (3) Vector- A
Select ON to display a vector-A line and vector waveforms.
- (4) Smoothing
Makes links between screens smooth.
- (5) Harmonic
Mode to extract harmonic components and create the image
- (6) Scope

Select the image display depth. Standard / Wide

- (7) Vector- A tool Operate the Vector-A line.
- (8) Image Quality Preset values of gains and dynamic ranges. Standard / High contrast / Wide dynamic range
- (9) Total Gain 20-80dB
- (10) Dynamic Range 30-70dB
- (11) Near Gain 1-60step
- (12) FREEZE / FREEZE cancel button
- (13) Probe Icon button Displays the probe set angle.
- (14) "Comment" button Allows you to enter a comment.
- (15) Zoom function Enlarges and displays the image being shown.
- (16) Save JPEG button Saves the image in Jpeg format.
- (17) FullSize button Zooms the image to full screen.


2.4.3 Axial length measurement/view screen

- (1) Gain display/adjustment field
- (2) Mode display – Contact/Immersion
- (3) Waveform display area
- (4) The waveform between these two cursors is the waveform of the lens.
- (5) The waveform on the right of this cursor position is measured as the waveform of the retina.
- (6) Gate cursor movement button Moves the gate cursor.
- (7) Measurement data will be taken when the waveform rises above this cursor/line position.
- (8) Measurement cursor/measurement line The distance at this cursor/line should be taken as the measurement data.
- (9) Measurement conditions display field
- (10) "View" button Displays the view screen.

-
- (11) "IOL" button
Opens the IOL power calculation screen.
 - (12) Deletes the measurement data currently displayed and measures the same eye again.
 - (13) Axial length measurement data display Displays measurement conditions and measured axial length.
 - (14) Measurement data display The axial length, anterior chamber depth, lens, and their average values are listed.
 - (15) "Delete/Restore" button Deletes and recovers the measurement data.
 - (16) The measurement data at the cursor position is selected to be used for calculating IOL power.
The selected measurement data is indicated by an "**".
 - (17) Selection cursor UP/DOWN buttons
 - (18) "Gate" button
 - (19) "Caliper" button
 - (20) "Measure" button


2.4.4 Pachymetry screen


- (1) Displays the actual measurements of the loaded measurement points. The latest loaded data are displayed during a measurement process.
- (2) Corneal thickness data display
- (3) Measurement data display
- (4) "Delete/Restore" button
- (5) "Caliper" button
- (6) "Meas. Point" button
- (7) "Change Meas. Point" button
- (8) Measurement waveform display field
- (9) Measurement conditions display field
- (10) "IOP" button
- (11) "Differential" button
- (12) "Calib." button

Calibrates the sensitivity of pachymetry probes.

2.4.5 A-Diag measurement screen

(FREEZE screen)

- (1) Gain
- (2) Displays tools for analysis of the line or point.
- (3) Hold this button briefly to delete "Ref." data and open the loading screen.
- (4) "Ref. Wave" button
- (5) Waveform display area Displays a waveform and the following information. Patient ID, inspection date/time, eye to be examined (R/L), image page No., converted acoustic velocity Point to which the probe is applied, direction of the ultrasound beams, probe type, beam direction, comment
- (6) "Comment" button
Allows you to enter a comment.
- (7) Gain adjustment 0 - 80dB

-
- (8) Scale selection Log / Linear / S
 - (9) Select the point to which the probe is applied and the direction of beams that are to be displayed in the waveform display field.
 - (10) FREEZE / FREEZE cancel button
 - (11) "Frame-by-frame" button


2.5 Symbols used for marking

Refer to instruction manual.

“ON” (power)

“OFF” (power)

Type B attached part

Grounding (earth)

3. OPERATION PROCEDURES

3.1 Preparation before use

3.1.1 Connections


- *The power connector completely isolates the instrument from the commercial power source. If there is a problem with the instrument, turn off the power switch and disconnect the power connector. Install the instrument in a place where this can be performed smoothly.*

- *Faulty connection of the instrument and relevant accessories may result in fire, electric shock, or fatal accident. Turn off all devices and disconnect the power plugs from the outlet before connection.*
- *The connector needs to be inserted in a specific direction. Check the direction and firmly insert the connector.*


a) Connection of the power cord


- 1) Insert the connector of the power cord into the power socket (1) in the correct orientation. Connect all three pins of the plug.

b) Connecting the probe cable

Fully insert it in the correct direction. Press the unlock button when disconnecting the probe cable.


c) Connecting the foot switch


- Insert the footswitch connector into the footswitch connector on the back of the instrument, aligning the notches of the connectors.
- Turn the fixture (3) of the connector clockwise until it stops to secure the connector.

d) Connecting the biometry probe / A-scan diagnosis probe

Insert the connector of the biometry probe into the biometry probe/A-scan diagnosis probe connector, indicated by “Axial” on the front of the instrument, in the correct orientation.

e) Connecting the pachymetry probe

Insert the connector of the pachymetry probe into the pachymetry probe connector, indicated by “Pachy” on the front of the instrument, in the correct orientation.

f) Connecting an external digital printer

Note

■ **Use only the specified printer and cables. Use a video signal cable or remote cable to connect this instrument to a device conforming to IEC60601-1 or JIS T 0601-1, or conforming to IEC60950-1 and with a power source isolated by an isolation transformer.**


• Video printer

Insert the cable plug A (1) of the USB cable into the USB connector on the side of the main unit in the correct orientation.

Connect the other cable plug B (2) of the USB cable to the video printer.


Follow the instruction manual of the video printer for details on how to connect the printer.

g) External ID input device

Plug the connector of the external ID input device (e.g. barcode reader, card reader, keypad, and keyboard) into the USB-H terminal (1) on the side of the main unit, checking the correct orientation.

h) Connecting the personal computer


- Note** ■ The network support system “TOMEY Link” (optional) or the examination data receiving software “DATA Transfer” (included in the package) is required for data communication with the instrument.
- Refer to the corresponding operation manual for installation, settings, and operation of TOMEY Link and DATA Transfer.
- Connection settings on the Setup screen must be completed in advance to connect with TOMEY Link and DATA Transfer. Refer to “3.9.4 Connecting & Printing” for the setting method.
- Be sure to make network settings under the consent of your network administrator.
- For LAN connection, be sure to connect the unit to a computer via a hub. The unit does not work correctly if directly connected to a computer.

• Connecting LAN cables


Prepare the following items.

- LAN cables (straight type, category 5 or higher) : 2
- A network hub (A 100MHz switching hub recommended) : 1
- A computer with TOMEY Link or DATA Transfer installed : 1

(Wiring example)

i) Connecting of the power plug of the Fixation Lamp for Chin Rest

- 1) The Fixation Lamp is interlocked with the ON/OFF power source for the instrument.
- 2) Insert the Power plug (1) for Chinrest fixation lamp into the terminal marked with "FIX LIGHT" provided at the rear side of the main unit.

3.1.2 Turning the power on and adjustment after turning the power on

a) Turning the power on

- 1) Turn on the power switch (1) on the front of the main unit.
- 2) The startup screen (Fig. 1) appears and then the measurement screen in the mode used last time appears.

(Fig. 1)

b) Adjustment after turning power on

The brightness of the monitor can be adjusted to suit the illumination in the examination room (see "3.9.1 General").

3.1.3 Switching modes

Note

■ *When the mode is changed, all data captured in the previous measurement mode is deleted.*

Touching the "Setup" button (1) on each screen displays the setup screen. Select any of the mode selection buttons (2) as desired to enter each mode.

There are following modes:


- B-Mode: B mode image diagnosis
- Axial : Axial length measurement
- Pachy : Measurement of corneal thickness
- A-Diag : A-scan diagnosis

3.1.4 Entering the patient data

Note


- When entering new patient information, carefully check that the same ID is not registered in the database or other electronic medical records.
- Only the first 14 digits of the ID number are displayed in the patient information field. Check that the ID number is correct on the patient information entry screen.
- The network support system “TOMEY Link” (optional) is required to use the patient information query function. “DATA Transfer” provided with the instrument is not available for query of the patient data.
- For details on the TOMEY Link Server settings, refer to the TOMEY Link instruction manual.
- Appropriate connection settings are required for connection with TOMEY Link. Refer to “3.9.4 c) PC setting” to make settings.
- Once examination data is output (printed, exported or saved), the ID cannot be changed.

- 1) Touch the patient information button (1) to open the Patient Information screen (Fig. 2).
 - When both the patient information and examination data exist, a warning message (Fig. 3) appears indicating you are about to overwrite data when the patient information button is touched. Touch the “Yes” button (7) to overwrite the data.

(Fig. 1)

(Fig. 2)

(Fig. 3) (7)

2) Enter necessary data.

– Entry method 1: Software keyboard

Patient's ID, name, sex, and date of birth can be entered using the software keyboard on the Patient Information screen (Fig. 2).

If the entered ID is already recorded in the memory, the corresponding name, sex, and date of birth is automatically recalled from the memory in the main unit when the ID is entered.

– Entry method 2: Database (USB memory stick)

The patient can be selected from the list saved in the database (USB memory stick).

Touch the "Patient Search" button (2) on the Patient Information screen (Fig. 2) to open the Patient List screen (Fig. 4). Touch the desired patient in the list (8) and then the "Select" button (9).


– Entry method 3: Inquiry to TOMEY Link

Enter the ID and touch the "Tomey Link Inquiry" button (3) on the Patient Information screen (Fig. 2). The patient information read from TOMEY Link in the personal computer appears. The inquiry reference can be changed to "DiscomWorkList" on the System Setup screen.

- Entry method 4: External ID input device

Data can be entered using a barcode reader, card reader, or keypad.

The barcode reader and card reader are also available on the measurement screen in addition to the Patient Information screen (Fig. 2)

(Fig. 4)

- 3) Touch the "OK" button on the Patient Information screen (Fig. 2), designate the patient information, and return to the previous screen.
 - Touch the "Cancel" button (5) to discard the edited information and return to the previous screen.


3.1.5 Deleting all measurement data (measurement preparation for other patient)

Note ■ *The deleted examination data cannot be restored. Carefully check the data before deleting it.*

■ *Be sure to touch the “New” button to delete the data for the previous patient before measuring a new patient. If a new measurement is started without deleting the previous data, the data for the previous patient may be included.*

1) Touch the “New” button for a while.

The patient information (ID, name, sex, and date of birth) and relevant examination data is deleted, and the measurement screen appears.

3.1.6 Selecting an eye

Select an eye using the eye selection button (1).

3.2 B mode image diagnosis

Refer to “3.1.3 Switching modes” for how to enter B mode image diagnosis.


When the probe is changed, all unsaved data captured by the previous probe is deleted. Be very careful when performing initialization.

The instrument automatically detects the connected probe.

3.2.1 Settings related to measurement

Touch the "Setup" button (1) to display the Setup (Ultrasound Pachy) screen (Fig. 2). Touch the "Exit" button (2) after settings are completed to apply the selected contents and return to the measurement screen (Fig. 1).

a) Acoustic velocity during measurement

Set the acoustic velocity used during measurement by the measurement function.

b) Settings for items displayed in image

Selects the items to be displayed in the image monitor.

Patient ID / Patient's eye / Measurement date/time / Comment / Focus marker (40 MHz only)

3.2.2 Displaying probe angle


Indicates the angle that the probe was set against the eye.


Confirm that the displayed probe angle matches the probe angle captured in the picture.

- 1) Touch the "Probe Icon" button (1) or the probe icon shown on the image to open the probe angle selection screen. The selected angle is shown in the lower right of the monitor. "N" means nasal (nose) and "T" means temporal (ear).


2) Touch the “OFF” button to hide the probe icon.

<Relationship between the displayed image and the probe angle>


The B mode probe has a mark (1) that indicates the scanning direction. The relationship between the scan direction mark and the image captured by the probe is as shown below:

<Understanding the probe angle mark>


The scan direction mark (1) on the B mode probe is indicated by the extending section (2) of the probe angle mark.

3.2.3 Gain adjustment

Note ■ Settings made here are only effective for the eye currently selected.

Examine the ultrasound image when you make any fine adjustment of the gain.

a) TOTAL Gain (1)

Sets the echo sensitivity of the overall image. (Adjustment range: 20 - 80 dB)

b) Dynamic Range (2)

Adjusts the dynamic range of the image. (Adjustment range: 30 - 70 dB)

A higher dynamic range makes the black and white contrast weaker. A lower dynamic range makes the black and white contrast stronger and creates a high contrast image.

c) Near Gain (3)

Sets the echo reflection sensitivity in the neighborhood of the anterior eye and vitreous.


3.2.4 Patient preparation

Note ■ *Cooperation by the patient is required to perform examination smoothly. Explain the examination method to the patient before starting the examination to help them relax.*

■ *Make sure that the air is not mixed in when applying the ultrasound diagnosis gel. When you observe unnecessary echoes, apply sufficient ultrasound diagnosis gel again.*

■ *When you use an eye cup, use one with the size appropriate for the patient. Leaked saline solution may cause an indirect contact infection.*

■ *Make sure that air bubbles do not mix in when you inject ultrasound media in the eye cup.*

(Fig. 1)

Follow the procedure below to use an eye cup.

- 1) Lay the patient on the back and apply topical anesthesia.
- 2) Apply the contacting part of the eye cup (1) to the eye (2).
- 3) Pour the saline solution (3) gently to the eye cup.


3.2.5 Obtaining ultrasonic image

- **Do not allow water or chemicals to splash on probes, except for the tip cap. Any water or chemicals entering the instrument may cause electric shock or failure.**

- **Be sure that the screw is not loose when you use the arm. Otherwise the probe attached to the arm may drop to be damaged or may drop to damage an eye ball or eye lid of the patient. Operating the arm directly with your hand may injure your finger pinched by the arm.**
- **Be sure to touch the “New” button to delete all the measurement data of the previous patient before measuring another patient. If new measurement is started without deleting the previous data, the measurement data of the previous patient may be included.**
- **Releasing a freeze by touching the FREEZE button or stepping on the foot switch deletes all images previously obtained. Be very careful when performing initialization.**

(1)

- 1) When the instrument is in FREEZE mode, touch the “FREEZE cancel” button (1) or step on the foot switch to release FREEZE mode.
- 2) Apply an appropriate amount of the ultrasound diagnosis gel to the contact section of the B mode probe.
- 3) Apply the B mode probe directly to the lid of the examined eye.

3.2.6 Various functions

Available functions vary depending on the measurement status. See the following table:

10MHz			
Playing Movie	FREEZE	Comment	FREEZE
Vector-A	FREEZE Real-time	Smoothing	FREEZE Real-time
Zoom	FREEZE	Measurement	FREEZE
Export	FREEZE	Print	FREEZE
Save in USB memory stick	FREEZE	Save in JPEG format	FREEZE
Measurement	FREEZE	Scope	Real-time
Full screen	FREEZE		

a) Playing Movie


Up to 200 images per eye for the patient are saved in the 10-MHz B probe. This function executes video playback operations and frame advance/reverse using the playback buttons. The “Stop” button is shown while a video is playing.

Playbacktool

b) Commenting


- 1) Touching the "Comment" button (1) displays the input comment window and the software keyboard.
- 2) Enter a comment and touch the “OK” button (2) to save the entered comment and close the comment input window.
- 3) Touching the “Cancel” button (3) will discard the changes made.

c) Vector mode-A

This function displays captured mode-A waveforms of the image.


Please note that the zoom and measurement functions cannot be used when vector mode-A is on.

- 1) Selecting "ON" for "Vector-A" (1) displays the cursor line and the mode-A waveform for the line position on the monitor.
- 2) The mode-A waveform of the cursor line position is shown on a real-time basis as you move the cursor line by touching the "Move Cursor Line up (or down)" button (2).
- 3) Touching a position on the screen moves the cursor to that position. Touching the "HRZN" button moves the cursor line to the center.


d) Smoothing

This function smoothes the continuation of images and reduces noises, improving the image quality. This function does not suit for viewing videos that include fast movements.

e) Zoom

This function magnifies the displayed image. Vector mode-A cannot be used when zoom is used.

The navigation monitor displays the whole image, and the area zoomed in the monitor is indicated by a frame. Moving the frame moves the displayed area accordingly.

Touching the “Zoom” button (1) alters the magnification between 100% and 200%.

f) Measurement

This function measures lengths, angles, and areas of the image displayed on the monitor.

■ ***When you use the result of this function for operations and medical treatments, make a comprehensive decision considering the results of other measurement devices and differences between left and right.***


■ ***The zoom function is not available during measurement. Apply the zoom function and then take measurement to measure on the enlarged image.***

■ ***The Vector-A display and smoothing function are not available during measurement.***

■ ***The measurement points and results are all deleted when the measurement tool is closed, and they cannot be recovered.***

- 1) Open the measuring tool by touching the “Open” button (1) in “Measurement” in the FREEZE menu.
- 2) Select a tab to display a measurement tool (length, angle, or area).
- 3) Touching the “Close” button (2) closes the measuring tool screen.

[Length measurement]

This function measures and displays the distance between two arbitrary points in an image.

Touching two points in the image automatically calculates the slant distance between the two points and displays the result in the image.

Touching the screen again changes the position of the active caliper (blue) and the angle is re-calculated.

The following tools are available.

SW-C button : Switches the active caliper (blue).

Cursor button : Moves the active caliper (blue).


Line button : Shows/hides the line between the caliper marks.

Perpendicular Line button : Shows/hides a line perpendicular to the above line.

Clear button : Deletes the measurement result.

[Angle measurement]

This function measures and displays the angle determined by three arbitrary points in an image.

Touching three points in the image automatically calculates an angle (acute angle) regarding the second point as the center and displays the resultant angle in the image.

Touching the screen again changes the position of the active caliper (blue) and the angle is re-calculated.

The following tools are available.

SW-C button : Switches the active caliper (blue).


Cursor button : Moves the active caliper (blue).

“Line” button : Shows/hides the line between the caliper marks.

“Clear” button : Deletes the measurement result.

[Area measurement]

This function measures a target area and displays the result. The target area can be defined freehand in the screen. The gradation value range can also be specified. Measurement is not performed if the specified area does not form a closed shape.

The following tools are available.

Gray level setting : Set the range of gradation.

Calculate button : Re-calculates the area.

Clear button : Deletes the measurement result.

g) Harmonic

This mode transmits 10MHz ultrasound (fundamental wave) and extracts twofold frequencies (harmonic ingredients) for creation of the image.

This mode creates images with better resolution in the azimuth direction (vertical direction in the image) and with less artifacts, providing a more faithful rendering of the tissue.

This mode, therefore, suits for obtaining clear structure images of an optic papilla, macular area, and tumor.

h) Scoping

This function selects the image display depth from “Standard” and “Wide.”

Choose one that suits the observation target area.

i) Preset values of gains and dynamic ranges

Preset values of each gain can be set in this instrument. Select a preset setting from Image Quality.

[Standard]

Each gain is set to the standard value.

[High Contrast]

This preset setting creates images of high contrast. However, its gradation granularity is low, not suited for observing subtle differences of tissues. This mode mainly suits observation of a retina.

[Wide Dynamic Range]


Its gradation granularity is high, so that it suits observation of subtle differences of tissues. This mode, however, tends to capture unnecessary echoes and artifacts (virtual images created by inflections and reflections of ultrasound beams). This mode mainly suits the observation of opacity and bleeding in a corpus vitreum.

j) Full screen display

This function zooms the measurement data to the maximum size on the screen.

Note ■ *The zoom function is not available while the full screen function is activated. Return to the measurement screen to edit the data.*

- 1) Touch the "FullSize" button (1) to open the Full Size screen.
- 2) Select a tab to display a measurement tool (length, angle, or area).
- 3) Touch the "Back" button (2) to return to the measurement screen.

• IOL eye (acryl)

Select this when there is an acrylic IOL implant.

• User Setting

Select this when registering optional materials for IOL eyes.

Register the settings on the converted acoustic velocity setting screen. (Refer to "3.3.1 b)

Setting converted acoustic velocity.")

"Anterior chamber depth" and "lens" are not measured for an aphakic eye and "lens" is not measured for an IOL eye.

The instrument may not recognize waveforms on the back of the crystal lens due to multiple echoes in the crystal lens in an eye with dense cataracts.

Inspected eye type	Axial length	ACD (anterior chamber depth)	Lens
Phakic eye (Average acoustic velocity, divisional acoustic velocity)	○	○	○
Dense cataract eye	○	○	△
Aphakic eye	○	*	*
IOL eye (All materials)	○	○	*

3.3 Biometry function

Refer to “3.1.3 Switching modes” for how to enter biometry mode.

3.3.1 Setting measurement conditions

Note ■ *Settings made here are only effective for the eye currently selected. Settings cannot be made for both eyes simultaneously. Complete necessary settings for each eye.*

Touch the "Setup" button (1) to display the Setup (Ultrasound Pachy) screen (Fig. 2). Set items related to operation conditions.

Touch the “Exit” button (2) after setting is completed to apply selected contents and return to the measurement screen (Fig. 1).

a) Selecting the eye type to be examined

Select the eye to be examined from the following 7 options according to the patient's eye.

- **Phakic eye**

Select this for a normal eye or when the crystalline lens nucleus is comparatively soft, such as an eye in the initial stages of cataracts.

- **Dense cataract eye**

Select this when the crystalline lens nucleus of the eye to be examined is comparatively hard, such as an eye with dense cataracts, and it is difficult to take measurements in the mode used for a phakic eye due to an echo reflected in the crystal lens.

- **Aphakic eye**

Select this in the case of aphakic eyes.

- **IOL eye (PMMA)**

Select this when a PMMA IOL is implanted.

- **IOL eye (silicon)**

Select this when a silicon IOL is implanted.

- **IOL eye (acryl)**

Select this when an acrylic IOL is implanted.

- **User Setting**

Select this when registering optional materials for IOL eyes.

Register the settings on the converted acoustic velocity setting screen.

(Refer to “3.3.1 b) Setting converted acoustic velocity.”)

“Anterior chamber depth” and “lens” are not measured for an aphakic eye and

“lens” is not measured for an IOL eye.

The instrument may not recognize waveforms on the back of the crystal lens due to multiechoes in the crystal lens in an eye with dense cataracts.

Inspected eye type	Axial length	ACD (anterior chamber depth)	Lens
Phakic eye (Average acoustic velocity, divisional acoustic velocity)	○	○	○
Dense cataract eye	○	○	△
Aphakic eye	○	*	*
IOL eye (All materials)	○	○	*


○ : Measurements are displayed.

△ : Measurements may not be displayed.

* : Measurements are not displayed

b) Setting converted acoustic velocity

Touch the “Change Setting” button (1) to open the setting change screen (Fig. 1).

(Fig. 1)

- 1) Entry boxes (2) for the acoustic velocities required for the eye currently selected appear. Enter the converted acoustic velocity using the software keyboard.
- 2) After entry, touch the "Enter" key to move to the next entry box. Also, by touching the "Incorporate in L" button (3) or the "Incorporate in R" button (4), the data entered will be applied to the other eye.
- 3) Touch the "Initial Setting" button (5) to enter the initial settings.
- 4) Touch the "OK" button (6) to apply the entry and return to the Setup (Ultrasound Axial) screen. Touch the "Cancel" button (7) to discard the entry and return to the Setup (Ultrasound Axial) screen with the previous sonic velocity still set.

(Initial setting and setting range of converted acoustic velocity)

- Phakic eye

Average acoustic velocity for axial length	:	1,550m/s, 1,500 - 1,600m/s
Crystal lens acoustic velocity	:	1,641m/s, 1,540 - 1,740m/s
Anterior chamber depth acoustic velocity	:	1,532 m/s, 1,430 - 1,630m/s
Vitreous acoustic velocity (divisional velocity only)	:	1,532 m/s, 800 - 2,000 m/s
- Dense cataract eye

Average acoustic velocity for axial length	:	1,548m/s, 1,500 - 1,600m/s
Crystal lens acoustic velocity	:	1,629m/s, 1,540 - 1,740m/s
Anterior chamber depth acoustic velocity	:	1,532 m/s, 1,430 - 1,630m/s
- Aphakic eye

Average acoustic velocity for axial length	:	1,532 m/s, 1,430 - 1,630m/s
--	---	-----------------------------
- IOL eye

Implanted lens acoustic velocity	(Acryl)	:	2,200m/s, 800 - 3,000m/s
	(Silicon)	:	1,049m/s, 800 - 3,000m/s
	(PMMA)	:	2,718m/s, 800 - 3,000m/s
Anterior chamber depth acoustic velocity	:	1,532 m/s, 1,460 - 1,630m/s	
Vitreous acoustic velocity	:	1,532 m/s, 800 - 2,000 m/s	
IOL thickness	(Acryl)	:	0.80 mm, 0.10 - 4.00 mm
	(Silicon)	:	1.00mm, 0.10 - 4.00 mm
	(PMMA)	:	0.80 mm, 0.10 - 4.00 mm

c) Measurement mode


Set the measurement mode.

- Auto
Select this for normal measurements.
- Auto Quick
Select this when measurement is difficult. The conditions for "Measurement Complete" will have greater variation than those of Auto.
- Manual
Select this when measurement is difficult in Auto and Auto quick modes.
A foot switch is required for measurement.

d) Setting contact/immersion mode

Complete the settings for contact mode and immersion mode referring to "3.9.2 Measurement."

3.3.2 Operation check

(Fig. 1)

Note *This biometry test piece cannot measure or calibrate the accuracy of the instrument.*

- 1) The biometry test piece provided with the instrument allows you to check the operation of the instrument.
- 2) Select "Aphakic eye" for the eye type to be measured.
- 3) Wet the surface of the test piece and set the biometry probe perpendicular to the surface to start measurement.


3.3.3 Preparation for measurement

Note *Cooperation by the patient is required to perform measurement smoothly. Explain the measurement method to the patient before starting measurement to help them relax.*

- 1) Check the measurement condition settings and operation as needed.
- 2) Apply topical anesthesia to the eye to be examined.
- 3) Have the patient sit down on a chair or lie on their back so that they are in a comfortable position while measurement is performed.


3.3.4 Precautions for measurement

a) Handling the biometry probe in contact/immersion modes

<Contact mode>

Directly apply the contact section of the biometry probe perpendicularly to the center of the cornea.


<Immersion mode>

When the immersion attachment is not used

Apply the axial length probe perpendicular to the corneal center where corneal protective media, etc. is applied.

Make sure that the clearance between the cornea and the contact section of the axial length probe is approximately 2.0 – 5.0 mm.


When the immersion attachment is used

■ **When putting the corneal protective media in the cup at the tip of the immersion attachment, pay extreme attention not to allow air to enter the cup.**


If air is in the cup, measurement cannot be performed or correct measurement cannot be obtained.

(Fig. 1)

(Fig. 2)


- 1) Lightly apply the corneal protective media etc. to the contact section (1) of the biometry probe and cover the probe with the immersion attachment (Fig. 1).
- 2) Put the corneal protective media (2) in the cup at the tip of the immersion attachment until it rises slightly.
- 3) Set the tip of the immersion attachment perpendicular to the center of the cornea, referring to the biometry probe axis and visual axis. The appropriate position is where the tip of the immersion attachment does not contact the cornea (Fig. 3).


(Fig. 3)

b) Attaching the biometry probe to the applanation tonometer

Use the attachment for applanation tonometer provided with the instrument when measuring axial length with the applanation tonometer.

(Fig. 1)

(Fig. 2)

- 1) Firmly hold the prism attachment section (1) of the applanation tonometer with your fingers so that it does not move.
- 2) Insert the applanation tonometer attachment (2) into the prism attachment section from the patient side (Fig.1).
- 3) Insert the biometry probe (3) into the attachment from the physician side, tighten the plastic screw (4) with a screwdriver and secure the biometry probe to the attachment (Fig. 2).
- 4) Secure the measurement knob of the applanation tonometer to the following positions to prevent pressure on the cornea.
 - 10 mmHg when the intraocular pressure of the eye is 20 mmHg or lower
 - 10 mmHg lower than the intraocular pressure when the intraocular pressure of the eye is 20 mmHg or higher

c) Selecting the retina waveform


Note

■ **Settings made here are only effective for the eye currently selected. Settings cannot be made for both eyes simultaneously. Complete necessary settings for each eye.**

Set the retina gate cursor to the left of the retina waveform when there is a waveform between the waveform on the back of the crystal lens and retina waveform and the instrument cannot recognize the original retina waveform. The instrument recognizes a waveform on the right of the retina gate cursor as the retina waveform and conducts measurement.

Select "Retina" and touch the gate cursor movement button to move the retina gate cursor to the left of the retina waveform.

Waveform generated
in the vitreous

Move the cursor to the left
of the retina waveform.


d) Selecting the waveform on the back of the lens

Note ■ *Settings made here are only effective for the eye currently selected. Settings cannot be made for both eyes simultaneously. Complete necessary settings for each eye.*

Set the gate cursor for the back of the lens to the right of the waveform on the back of the lens when there are multiple waveforms in the crystal lens and the instrument cannot recognize the original waveform on the back of the lens. The instrument recognizes a waveform on the left of the gate cursor for the back of the lens as the waveform on the back of the lens and conducts measurement. The gate cursor on the back of the lens can be set only when the eye type is set to “Phakic eye” or “Dense cataract.”


Select “LENS Rear” and touch the cursor movement button to move the lens rear gate cursor to the right of the lens rear waveform.

Waveform generated
in the crystal lens

Move the cursor to the right of
the waveform of the back of the lens.

e) What is a good waveform?

[Contact mode]

The instrument recognizes a waveform as a good one if conditions 1) to 3) listed below are satisfied and captures the measurement data in auto measurement mode.

- 1) Each of the following waveforms rises above the level cursor/line position.

- Eye with initial stage cataracts

Waveform on the front of crystal lens / Waveform on the back of crystal lens /

Retina waveform

- Dense cataract eye

Waveform on the front of crystal lens / Retina waveform

- Aphakic eye

Retina waveform

- IOL eye

Waveform on the front of crystal lens / Retina waveform

- 2) The retina waveform rises vertically. (Ultrasound wave reaches the retina perpendicularly.)

When the waveform slopes upward, the measurement tends to be long.


- 3) Within the measurement range

The following is not a requirement for capturing the measurement data, but a method to confirm that the ultrasound wave correctly captures the geometrical axis and ideal waveforms are taken.

- i. When the retina waveform spikes, the ultrasound wave reaches the retina perpendicularly.
- ii. When the waveforms on the front/back of the crystal lens spike, the ultrasound wave correctly captures the geometrical axis.
- iii. When the retina waveform can be distinguished from the sclera waveform, the ultrasound wave reaches the retina perpendicularly. Because the depression between two waveforms (choroid) cannot be confirmed when the gain is set high, etc., these waveforms may not need to be distinguished from each other.
- iv. When there is no tailing that follows the initial waveform (cornea waveform), the biometry probe directly contacts the retina. When there is a layer of tears or cornea protective agent between the probe and cornea, tailing occurs after the cornea waveform. In this case, the measurement data may be unstable or tends to be longer.

[Immersion mode]

The following conditions are added to conditions (1) to (3) in contact mode for automatic measurement in immersion mode.

-
- 4) The cornea waveform is formed within 2.0 – 5.0 mm from the initial waveform.

(Range between the dotted lines in the illustration)

The following is not a requirement for capturing the measurement data, but a method to confirm that ideal waveforms are taken. Check the following together with items i – iv in contact mode.

- v. There are no unnecessary waveforms between the initial waveform and cornea waveform. An unnecessary waveform appears when air bubbles are contained in the ultrasound media in the cup at the tip of the immersion attachment or the media placed between the biometry probe and cornea.

3.3.5 Measurement:

- *Sterilize the biometry probe before use.*
- *Do not use the biometry probe with its contact section damaged. Measurement error may occur, and the cornea may be injured.*
- *Immersion attachment is a disposable part. Do not reuse it. Otherwise, you may contract diseases.*


- *The converted acoustic velocity directly affects the measurement data. Check that the desired value is set before starting measurement.*
- *Check the settings in contact mode and immersion mode before starting measurement. In immersion mode, apply an ultrasound media such as cornea protective agent so that the clearance between the cornea and the contact section of the biometry probe is approximately 2.0 – 5.0 mm.*
- *Be very careful not to apply pressure to the cornea when taking measurement. Also, do not apply too much cornea protective agent if used. The measurement data may be affected.*
- *The automatic measurement function is an auxiliary function to take measurements more easily. This is not a function for making actual clinical judgments.*
- *Be sure to touch the “New” button to delete all the measurement data of the previous patient before measuring another patient. If new measurement is started without deleting the previous data, the measurement data of the previous patient may be included.*

a) Auto measurement

A waveform that satisfies the conditions is captured and up to 10 values are displayed. 20 or more data sets are captured and measurement is completed when variations are within a specific value.

b) Manual measurement

- 1) The instrument beeps when the capturing conditions are satisfied. Press the foot switch to capture the measurement data.
- 2) Press the foot switch again to proceed to the next measurement and take measurements in a similar manner. When 10 data sets are taken, a beep goes off and measurement is completed.

c) Setting gain

Note ■ *Settings made here are only effective for the eye currently selected. Settings cannot be made for both eyes simultaneously. Complete necessary settings for each eye.*


Adjust the waveform height according to the gain settings. The gain and the waveform height increases as the gain value increases.

Touch the gain adjustment buttons (1) to change the gain.

d) Setting gate cursor

Note ■ *Settings made here are only effective for the eye currently selected. Settings cannot be made for both eyes simultaneously. Complete necessary settings for each eye.*

-
- 1) Select a gate cursor to be adjusted by touching the gate selection buttons (1).
The selected and active cursor is displayed in red but the other cursors are displayed in white.
 - 2) Touch the gate cursor movement buttons (2) to set the active gate cursor position.

e) Settings for re-measurement of the same patient

Note ■ *The deleted data cannot be restored. Carefully check the data before deleting it.*

Retain the patient information and delete only the measurement data for each eye.
Delete all data referring to “3.1.5 Deleting all measurement data (measurement preparation for other patient)” when the patient is changed.

Touch the “Retake” button (1) for approximately 1 second until a beep sounds to delete the data. It retains the patient information and the measurement data of the eye other than the inspected eye being displayed.

3.3.6 Checking waveforms after measurement

a) Displaying waveform of arbitrary measurement data

After all necessary measurement data is captured or when the “View” button (1) is touched, the edit screen (Fig. 1) opens.

- 1) The waveform measured at the selection cursor is displayed. Touch the selection cursor movement buttons (2) to move the cursor up or down.

The following marks are assigned to the measurement data.


- * : Data used for calculating IOL power
- L : Longest axial length
- S : Shortest axial length
- C : Callipered data

b) Deleting and recovering measurement values

- 1) Touch the “Delete” button (1) to delete the data at the selection cursor. The data deleted here is excluded from the average.
- 2) To recover the deleted data, move the cursor to the data. Touch the “Recall” button (2) to cancel deletion. However, once the screen returns to the measurement screen, the deleted data cannot be restored by using the “Return” button even when the screen is switched to the edit screen again.

c) Gate change function


- 1) Touch the "Gate" button (1) to open the change gate screen (Fig. 2).
- 2) Touch the gate cursor selection buttons (2) to select the gate cursor to be changed.
- 3) Touch the "gate cursor movement" buttons (3) to change the position of active gate cursor. The modified measurement data is displayed in the edit data display field along with the movement of the gate cursor.
- 4) Touch the "Initial position" button (4) to discard changes.
- 5) Touch the "Apply" button (5) to define the changes and closes the change gate screen. If the "Cancel" button (6) is touched, the changes made are discarded and the change gate screen is closed.

d) Caliper function

Note ■ Values measured and displayed using the caliper function are rough estimates and may differ from the actual measurement result.


This function is used to measure the distance of the measurement waveform. Dotted caliper lines appear and the distance to the selected caliper line is displayed.

- 1) Touch the "Caliper" button (1) to open the caliper screen (Fig. 3).
- 2) Touch the "Switch" button (2) to select the caliper line to be changed.
- 3) Touch the caliper line movement buttons (3) to change the active caliper line position. The modified measurement data is displayed in the edit data display field (4) along with the movement of the caliper line.
- 4) Touch the "Initial position" button (4) to discard changes.
- 5) Touching the "Apply" button (5) will apply the changes made and close the caliper screen. Touching the "Cancel" button (6) will discard the changes made, and the caliper screen is closed.

e) Selecting specific measurement data to be used for calculating IOL power

The average value is normally adopted for the measurement data to be used in calculating the IOL power, but specific data can also be selected.

- 1) Move the selection cursor to the data which is to be used in calculating the IOL power. To use the average value in IOL power calculation, move the selection cursor to the average data at the bottom.
- 2) When the "Select" button (1) is touched, it is selected as the data to be used in calculating the IOL power, and a "*" appears on the left of the data.

3.4 IOL power calculation

Note ■ When using measurement results for calculation of the IOL power, the physician must examine the measurement result beforehand.

■ Refer to "3.8 Export, print, and save" for waveforms printed and exported from the IOL calculation screen.

Touch the "IOL" button (1) located on each screen to open the IOL power calculation screen.

3.4.1 Calculation

Note ■ When using measurement results for calculation of the IOL power, the physician must examine the measurement result beforehand.

■ Calculation by this instrument may cause some errors due to the number of significant digits in internal calculations.

■ Complex numbers may be generated for the SRK/T formula. In this case, the “√ section” is calculated as zero and an asterisk “*” is shown on the right of the calculation result.


This instrument automatically starts calculation and displays the result when all items required for IOL power calculation are set.

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When an axial length and anterior chamber depth measurement is completed, the measurement data is already entered automatically and cannot be entered manually.

When measurement has not yet been performed, enter the data as follows.

- 1) Touch the axial length input field (1) or the ACD input field (2) to activate it. Enter the data and touch the “Enter” key to apply the data.

[Input range]

Axial length : 13.00 - 45.00 mm

Anterior chamber depth : 0.00 - 10.00 mm

b) Corneal refractive power and radius of corneal curvature (K1/K2)

- 1) Touch the K1 input field (1) or K2 input field (2) to activate it. Enter the data and touch the “Enter” key to apply the data.

[Input range]

Corneal refractive power : 30.00 - 60.00D

Radius of corneal curvature : 5.00 - 11.00 mm

c) Target Ref.


- 1) When the converted sonic velocity entry field (1) is touched, it is activated. Enter the data and touch the “Enter” key to apply the data.

[Input range] Angle : -30.00 - +10.00

The entered value is stored in the main unit, and is not cleared even when the power is turned off.

d) Lens constants (A-constant/SF/ACD-constant/a0·a1·a2)

Enter various lens constants for IOL according to the formula. This instrument is able to calculate up to 4 constants for 1 formula simultaneously. The formula (1) corresponds to lenses (2) and (3), and the formula (4) corresponds to lenses (5) and (6). There are 2 ways to enter the values. One is to enter them directly using the keypad, and the other is by selecting them from the IOL data list.

• Entering using keypad

The lens constants a0, a1, and a2 in the Haigis optimized formula cannot be entered directly using the keypad. Refer to “Entering through IOL data list” below.


- 1) Touch an input field (1) to activate it. The keypad (2) appears. Enter the data and touch the “Enter” key to apply the data.

[Input range]

A-constant : 100.00 - 130.00
 ACD-constant : 0.00 - + 10.00
 SF : - 5.00 - + 10.00

■ Entering through IOL data list

When IOL data has already been registered in the “IOL data registration,” you can select data from the IOL data list.

- 1) Touch an input field (1) to activate it. The IOL data list (2) appears.
- 2) You can check the registration content by sliding the IOL data list being displayed with the scroll bar (3).
- 3) Select data using the “IOL data selection” buttons (4).
- 4) Touch the "Edit" button (5) to display the screen for editing the registration content.
- 5) Touch the “Apply” button (6) to apply the lens constants and to display them on the previous screen.

e) Entering parameters for Clinical History Method

Activated only when SRK/T Double K is selected as the formula.

The screenshot shows the TOMY IOL Master software interface. At the top, patient information is displayed: ID: 1234567890, Name: TOMY TARD, Phys.: AAA, Time: 2015.Feb.06 09:44, Ver.000. The main input area is divided into several sections. The 'Ultrasound' section includes 'Axial' (24.95) and 'ACD' (3.63). The 'K' section includes 'K1' (40.00 D), 'K2' (40.00 D), 'K1pre' (1.3375), 'K2pre' (S), and 'Kpost' (12.0). The 'Ref.' section includes 'Ref. pre' (S) and 'Ref. post' (C). A pull-down menu for 'V.D.' is set to '12.0'. The 'IOL Power Formula' section shows 'Double K SRK/T' selected. The 'Haigis optimized' section shows parameters for 'a0', 'a1', 'a2', 'M4', and 'E4'. A numeric keypad is visible at the bottom left. Four numbered callouts point to specific fields: (1) points to 'K1pre', (2) points to 'Ref. post', (3) points to 'Kpost', and (4) points to the 'V.D.' pull-down menu.

- 1) Enter the corneal refractive power or radius of corneal curvature before the surgery for correcting the refractive power in K1pre and K2pre fields.
- 2) Manually enter eye refractive power before refractive correction surgery in Ref.pre (S,C) (1) and eye refractive power after refractive correction surgery in Ref.post (S,C) (2).
- 3) When values are entered in K1Pre, K2Pre, Ref.pre (S,C), and Ref.post (S,C), the calculation starts automatically and the corneal refractive power or radius of corneal curvature after the refractive power correction surgery is displayed in Kpost (3). The data can be entered directly. However, Ref.pre (S,C) and Ref.post (S,C) cannot be entered.
- 4) Select the distance between the vertexes (V.D.) from the pull-down menu (4).
C.L./12.0/13.5/14.0/15.5/16.0

3.4.3 Setting calculation formula

Touch the calculation formula pull-down button (1) and select a formula.

The following 9 types of IOL power calculation formulae are provided with this instrument.


The IOL calculation formulae listed in the pull-down menu can be set as described in “3.9.3 b) Selecting IOL power formula.”

- Haigis standard formula
- Haigis optimized formula
- Hoffer® Q formula
- Holladay 1 formula
- SRK-II formula
- SRK/T formula
- SRK SHOWA formula
- Shammas-PL formula
- SRK/T Double K formula

3.4.4 Entering data after surgery

Select the name of the IOL type that was actually implanted in the surgery.

- 1) Select the type of implanted IOL. When the entry field (1) is touched, it is activated. Select the name of the IOL type that was actually implanted in the surgery from the IOL list.
- 2) Enter the implanted IOL power (2), date of surgery (3), and refractive power after surgery (4). The keypad or calendar appears when the input field is touched.
- 3) When the “Personal Constant” button (5) is touched, the Personal Value screen (Fig. 1) appears. Touch the “Print” button (6) to print the personal constants only.

(Fig. 1)


(6)

3.5 Statistical processing

Note ■ The lens constant registered on this screen will be used for the next IOL calculation. Carefully check the data before registration.

Statistical processing can be performed on this screen using the examination data accumulated in the facilities. Predictive errors when using the registered lens constant are calculated and the histogram is displayed. In addition, the optimized lens constant can be calculated and registered according to the statistical processing result.

Calculation results when lens constants are entered manually are excluded from statistical processing.

Set filter conditions (1) for targets subject to statistical processing.

- IOL : Select IOL. Touch the “Select IOL” button. The IOL data selection screen appears to allow you to select an IOL model.
- Term (data extraction term) : Set the term subject to statistical processing.
- Physician : The registered physician IDs are displayed in the pull-down menu. “All” includes all physicians.

-
- 1) After all the conditions are set, touch the “OK” button (2). Statistical processing is performed and the result appears on the right side of the screen.

The histogram (3) in the upper section shows registered lens constants and the statistical processing results of those lens constants. When lens constants are not registered, the result does not appear.

 - Registered lens constant
 - Histogram of predictive error in refractive power
 - Average of predictive error in refractive power (Avg)
 - Standard deviation of predictive error in refractive power (SD)
 - No. of data sets (n)

The histogram (4) in the lower section shows the average of lens constants optimized by statistical processing and predictive results of processing when the average is used.

 - Average (Avg) of optimized lens constant
 - Histogram of predictive error in refractive power
 - Average of predictive error in refractive power (Avg)
 - Standard deviation of predictive error in refractive power (SD)
 - Number of data (n)
 - Percentage that predictive errors in refractive power fall under ± 0.5 [D], ± 1.0 [D], or ± 1.5 [D]
 - 2) Select an IOL power formula (5). Available formulae are as follows.
 - Haigis standard formula
 - Haigis optimized formula
 - Hoffer[®] Q formula
 - Holladay 1 formula
 - SRK-II formula
 - SRK/T formula
 - SRK SHOWA formula
 - Shammas-PL formula
 - SRK/T Double K formula
 - 3) Click the “Update IOL Constant” button (6) to register the optimized lens constant for the target manufacturer, model, physician, and formula. When the button is touched, the IOL data change confirmation screen appears. Click the “OK” button to complete registration.
 - 4) Touch the “Setup” button (7) to re-set filter conditions.
 - 5) Touch the “Print” button (8) to print the displayed results.
 - 6) Touch the “Exit” button (9) to close the statistical processing screen.

3.6 Pachymetry function

Refer to “3.1.3 Switching modes” for how to enter pachymetry mode.

3.6.1 Setting the data type to be displayed

Select the measurement data type to be displayed from the following three options.


- Latest : Displays the last measurement data taken.
- Minimum : Displays the minimum value of measurement data.
- Average : Displays the average of measurement data.

Make settings, referring to the instrument settings “measurement settings.”


3.6.2 Setting measurement conditions

Note ■ *Settings made here are only effective for the eye currently selected. Settings cannot be made for both eyes simultaneously. Complete necessary settings for each eye.*

Touch the "Setup" button (1) to display the Setup (Ultrasound Pachy) screen (Fig. 2). Set items related to operation conditions. Touch the “Exit” button (2) after setting is completed to apply selected contents and return to the measurement screen (Fig. 1).

(Fig. 1)

(Fig. 2)

a) Measurement range

Select the measurement range from the following 3 options.

- 150 - 350μm
- 300 - 1000μm
- 900 - 1500μm

b) Selecting measurement values to be displayed

Select measurement values to be displayed from the following 2 options.


- Actual measurement
- Bias value

c) How to display converted acoustic velocity and bias values

Select how to display the bias value from the following two options.

- Percentage bias: Converts the actual measurement using the preset bias rate (percentage) and displays the result.
- Plus or minus bias: Adds subtracts the preset correction value to from the actual measurement, and displays the result.

- 1) Touch the “Change Setting” button (1) on the Setup (Ultrasound Pachy) screen to display the screen for setting converted sonic velocity and bias value (Fig. 3).

(Fig. 3)


- 2) When the converted sonic velocity entry field (2) is touched, it is activated. Enter the converted sonic velocity using the software keyboard. Input range of converted acoustic velocity: 1400 - 2000 m/s
- 3) Touch the "Display of Bias" button (3) to select the percentage bias (%) or plus/minus bias (μm).
Using the software keyboard, enter the bias rate when percentage bias is selected, or the correction (4) when plus/minus bias is selected.
Input range of bias rate: 60 - 130%
Input range of correction: -600 - +450 μm
- 4) Touch the “Initial Setting” button (5) to enter the initial settings.
- 5) Touch the “OK” button (6) to define the entered settings. Touch the “Cancel” button (7) to discard the entered settings and close the setting change screen.

d) Measurement mode

Set the measurement mode.

- Auto mode : Select this for measurement in Auto mode.
- Manual mode : Select this when measurement is difficult in Auto mode.

3.6.3 Displaying and setting measurement points


The measurement point indicates “diameter – angle – S/I.”


- Radius : Radius from the center (mm)
- Angle : The upper and lower angles (°) when regarding the horizontal axis as zero degrees
- S/I : Section above the horizontal axis (Superior), section below the horizontal axis (Inferior)

The center point is indicated as “CCT.”

- 1) Touch the "Meas Point" button (1) to display measurement points in the measurement data display field (2).

- 2) Touch the “Change Meas. Point” button (3) to open the measurement point change screen (Fig. 1) to change the measurement point.


- 3) Touch the input field to activate it. Set values using the keypad and the like. As for the setting of “S/I,” the indication switches between “S” and “I” each time the S/I button (4) is touched.
- 4) Touch the Preset button (5) to set preset values. (Refer to “3.9.2 d) Measurement setting Axial.”)
- 5) Touch the “Set” button (6) to define the settings. Touch the “Cancel” button (7) to discard the settings and close the measurement point change screen.

3.6.4 Operation check

■ **The corneal thickness test piece is for operation checks. This test piece cannot measure or calibrate the accuracy of the instrument.**

(Fig. 1)

The corneal thickness test piece allows you to check the operation of the instrument.

- 1) Select “Auto” for the measurement mode.
- 2) Place the test piece in water and set the pachymetry probe perpendicularly to start measurement.

3.6.5 Preparation for measurement

■ **Cooperation by the patient is required to perform measurement smoothly. Explain the purpose of the examination and measurement method to the patient before starting measurement to help them relax.**

- 1) Check the measurement condition settings and operation as needed.
- 2) Apply topical anesthesia to the eye and have the patient sit down on the chair or lie on their back so that they are in a comfortable position while measurement is performed.

3.6.6 Measurement:

■ *Sterilize the pachymetry probe before use.*

■ *Do not use the pachymetry probe with its contact section damaged. Measurement error may occur, and the cornea may be injured.*

■ *The converted acoustic velocity directly affects the measurement data. Check that the desired value is set before starting measurement.*

■ *The measurement data can be displayed as the actual measurement or a bias value. Check the setting of the data display mode.*

■ *When the surface of the cornea is dry, wet it with saline solution.*

■ *The automatic measurement function is an auxiliary function to take measurements more easily and is not a function used to actually make clinical judgment. The physician must examine the measurement result before using it.*

a) Auto measurement

- 1) Set the pachymetry probe perpendicular to the point to be measured.
- 2) The instrument automatically starts capturing measurement data when capturing conditions are satisfied.
- 3) A beep sounds when the measurement data is captured.
- 4) When 20 sets of measurement data are captured, a beep sounds twice to indicate completion of measurement.
- 5) When the pachymetry probe is removed from the cornea, the next memory slot enters standby mode. Measurement is completed and the next memory slot enters standby mode when the probe is removed from the cornea even before 20 data sets are captured.
- 6) Take measurements for all 10 memory allotments in the same manner.
- 7) When retaking measurements in a specific memory slot during or after measurement, touch the desired memory number to enter standby mode and perform measurement again.

b) Manual measurement

- 1) Set the pachymetry probe perpendicular to the point to be measured.
- 2) The instrument beeps when the capturing conditions are satisfied.

- 3) Step on the foot switch to capture measurement data. A beep sounds when the measurement data is captured.
- 4) Repeat measurements in the same manner. When 10 sets of measurement data are taken, a beep sounds to indicate completion of measurement.
- 5) To take measurement data in the next memory slot before 10 sets of measurement data have been taken, touch the next memory number to enter standby mode and perform measurement.
- 6) When the pachymetry probe is removed from the cornea, the next memory slot enters standby mode.

c) Re-measurement of the same patient

Note ■ *The deleted data cannot be restored. Carefully check the data before deleting it.*

Retain the patient information and delete only the measurement data for each eye. Delete all data referring to “3.1.5 Deleting all measurement data (measurement preparation for other patient)” when the patient is changed.


Touch the “Retake” button (1) for approximately 1 second until a beep sounds to delete the data. The patient information is retained and the measurement data of the eye other than the inspected eye being displayed.

3.6.7 Checking measurement data

a) Intraocular pressure correction

This instrument automatically starts calculation and displays the result when all items required for calculating the intraocular pressure correction are set.

- 1) Touch the "IOP" button (1) to open the intraocular pressure correction screen (Fig. 1). When a measurement has already been completed at this point, the measured value is displayed in the CCT field (2).

However, if measurement points are set and two or more "CCTs" are included in them, the average of the CCTs is displayed. When measurement points are not set, the average of all memory numbers is displayed.

- 2) The input fields of the parameter 1 (3) and parameter 2 (4) for the intraocular pressure correction formula become active when touched, and the keypad appears. Also, when a formula selection button (5) is touched, the parameters already set are displayed in the intraocular pressure correction formula.

[Input range]

- Parameter 1 : 0 - 1500
- Parameter 2 : 0.0000 - 1.0000


- 3) Enter the intraocular pressure data. Touch the input field (6). The keypad appears.

[Input range]

- Intraocular pressure data : 1 - 60.0 (mmHg)
1.33 - 79.99 (hPa)

- 4) Touch the "Back" button (7) to close the intraocular pressure correction screen (Fig. 1) and return to the previous screen.

b) Deleting part of measurement data


- 1) Touch a memory number button (1) to select the measurement data to be deleted.
- 2) Touch the "Delete" button (2) to delete the selected measurement data, and the average and the standard deviation of the corneal thickness are recalculated.
- 3) If you have deleted data mistakenly or want to cancel the deletion, select that data number and you will see that the "Delete" button (2) will change to the

“Return” button (3). Touch the “Return” button to cancel the deletion. When new measurement data is captured, however, the deleted data cannot be restored by pressing the “Return” button.

c) Caliper function

The caliper function is for temporary reference and cannot store the callipered data. The caliper screen cannot be opened while a bias value is selected.


(Fig. 2)

- 1) Touch the "Caliper" button (1) to display the caliper screen (Fig. 2).
- 2) Touch the “Switch” button (2) to select a caliper line to be adjusted.
- 3) Touch the “active caliper line movement” button (3) to change the position of the active caliper line. In conjunction with the movement of the caliper line, the distance between the caliper lines (4) is displayed.
- 4) Touch the “Close” button (5) to close the caliper screen (Fig. 2) and return to the previous screen.

d) Difference display function

This function displays the difference value between the 2 selected measurement data.


(Fig. 3)

- 1) Touch the “Differential” button (1) to open the difference display screen (Fig. 3).
- 3). When a measurement has already been completed at this point, the

measured values are displayed in the data 1 display field (2).

- 2) The measurement data being stored can be selected as Data 1 or Data 2. Touch the "Memory Data" button (3) to display a list of measurement data stored in a USB memory. (Fig. 4) If a patient ID has not been selected, however, the Preservation Data List screen (a list of patient IDs stored in a USB memory) opens. Refer to "3.8.3 Saved data management."

(Fig. 4)

- 3) Touch the "Data Display" button (4) to display the selected measurement data in the subtraction display screen.
- 4) When measurement data are selected in Data 1 and Data 2, the difference value between Data 1 and Data 2 is displayed in the difference display field (5). If measurement points have been selected, however, only the difference values of the memory numbers whose measurement points match in both data are displayed.
- 5) When the "Data 1 - Data 2" button (6) is touched, the indication on the button will change to "Data 2 - Data 1," and the difference value of Data 2 - Data 1 will be displayed.

3.7 A mode diagnosis


3.7.1 Connecting the probe

Connect the A-scan diagnosis probe to the instrument.

3.7.2 Switching to A mode diagnosis function

Refer to “3.1.3 Switching modes” for how to enter A mode diagnosis mode.

3.7.3 Registering Ref. Data



Load the reference data for line analysis. Loading of data and adjustment of gains will be described later. Adjust the line position with the “HRZN” (select line position) button (3) as needed and touch the “Set Ref.” button (1) to register the line position after freezing. The registered Ref. data will be retained even after turning power off. Hold the “Reset” button (2) briefly to delete the data.

After registration, touching the “Ref. Wave” button (4) displays the message “The ΔdB is calculated by using Ref. data” and reference waveforms.

3.7.4 Measurement:

Note ■ *CVCooperation by the patient is required to perform examination smoothly. Explain the examination method to the patient before starting the examination to help them relax.*


(1)

Apply an appropriate amount of cornea protective agent to the contact section of the A-scan diagnosis probe, and apply the probe to the eye to be examined. Adjust the total gain, and touch the “FREEZE” button (1) or step on the foot switch to capture waveforms.

3.7.5 Setting gain

Note ■ *Settings made are only effective for the eye currently selected.*


Adjust the waveform height according to the gain settings.

The gain value (1) shown in the upper left of the measurement screen (Fig. 1) can be adjusted with the gain adjustment buttons (2).

The applicable setting range is 0 - 80 dB.


3.7.6 Switching data display

This instrument stores waveforms for up to 100 data sets. Captured waveform data is displayed in order every time the frame advance/retract buttons (1) is touched.

3.7.7 Setting the beam direction

Use the beam direction tabs (1) to change the position to apply the probe or the direction of the ultrasound beam, or to display the position of a point of note (with pathological changes, etc.).

Probe set angle selection button: 9 directions


Beam direction 1: This helps indicate the angle of the meridian.

1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12


Beam direction 2: This helps to specify a notable point in an eyeball.

AX, P, PE, EP, E, EA, O, CB


(Probe set angle)

(Beam direction 1)


(Beam direction 2)

3.7.8 Gradation characteristics function


When executing line analysis, select “Log, Linear or S” (1) for the gain curve that converts ultrasound waves obtained by receiving ultrasound wave echoes generated from the ultrasound transducer.

3.7.9 Analysis functions

a) Line analysis

- 1) Select “Line Analysis” (1) for Analysis Method.

- 2) Select “LINE” for Analysis Method.


If the reference waveform is not yet registered, register it referring to the aforementioned “Registering Ref. Data.”


- 3) Take a measurement. Adjust the gain so that the peak of the target section is aligned with the light height and take measurements.
The gain when measuring the target section, reference gain, and the difference between them are shown (2).

b) Point analysis

- 1) Select "Point Analysis" (1) for Analysis Method.

- 2) After measurement is completed, the "X" cursor appears on the waveform monitor. Move the cursor by touching the screen or using the cursor movement button (2), and touch the "SetP1" button to set the gain value at the cursor position as P1.

- 3) Move the cursor to the next point position and touch the "SetP2" button (4) to set the gain as P2. Two gains and the difference between them are shown.
- 4) When the point position needs to be corrected, move the cursor and touch the "SetP1" or "SetP2" button again. The point position is overwritten and the difference is re-calculated.


3.8 Export, print, and save

Set desired functions to 2 buttons for Export, Print, and Save provided on the operation screen before executing the export, print, and save functions. (Refer to “3.9.4 a) Setting functions of save, export, and print buttons.”)

The print, export, and save buttons work according to the preset functions when touched after measurement.

Export

Print

Save

Export & Save

Export & Print

Save & Print

ALL (internal)

3.8.1 Export

Note ■ The network support system “TOMEY Link” (optional) or the examination data transfer software “DATA Transfer” (included in the package) is required for data communication.

■ Refer to the corresponding operation manual for the TOMEY Link and DATA Transfer settings.


■ Connection setting is required for connecting with TOMEY Link and DATA Transfer. Refer to “3.9.4 c) PC setting” to make settings.

■ **Data cannot be sent in some configurations if the ID number is not entered. Refer to “3.9.1 General.”**

■ **It is recommended to send data after entering the ID number. The data can be sent without ID numbers entered; however, the corresponding patient may not be traced from the inspection file output from DATA Transfer. When the data was sent without entering the ID number, check the content of the inspection file immediately after sending the data, and move and save the file in an appropriate location. In addition, when data is sent to TOMEY Link, all data without ID numbers is regarded as data of the same patient, so must be handled with great care to avoid mishandling.**

■ **Waveforms of the same form as those to be printed out are output. Refer to “3.8.2 Print” for details.**

- 1) Touch the “Export” button (1) to open the Transmission confirmation screen (Fig. 2).

(Fig. 2).

- 2) The transfer confirmation screen (Fig. 1) displays patient ID, patient name, and sex. Touch the "Export" button (2) to start sending data. Touching the "Cancel" button (3) returns you to the previous screen without sending data.

3.8.2 Print

The printer and printing type can be selected in the system setup. Refer to “3.9.4 b) Selecting and setting printer.”

[Printing by built-in printer]

Note ■ *Different waveforms are printed depending on the screen in biometry mode.*

Waveforms to be printed on each screen are as below.

Measurement screen:

The waveform which is the closest to the average is output.

View screen:

Waveform displayed on the screen is printed.

IOL screen:

The waveform the closest to the average value is printed when entering the IOL screen from the measurement screen. The waveform displayed on the view screen is printed if entering the IOL screen from the view screen.

a) Common items

(1) — YMD:2015/Jun/17 19:47:09
 (2) — UD-800 Ver.00
 (3) — Clinic ID/Name:
 (4) — TOMEY CLINIC
 Clinic Address:
 ABCDEFGHIJKLMNOPQRSTU
 (5) — ID : 123456789012
 (6) — Name : TOMEY TARO
 (7) — Sex : Female
 (8) — Phys : YAMADA

RIGHT

Auto

Phakic(M/V)

Velocity : Avg. 1550m/s

Velocity : ACD 1532m/s

Velocity : LENS 1641m/s

Gain : 20dB

Avg. AXIAL : 25.12mm

SD : 0.07mm

RANGE : 0.20mm

Avg. ACD : 3.38mm

Avg. LENS : 4.35mm

NO AXIAL ACD LENS

1 25.13 3.39 4.70

25.12 3.38 4.35

Avg 25.12 3.38 4.35

No.6

AXIAL 24.13

ACD 3.39

LENS 4.31

Print:2015/Jun/19 19:55

TOMEY CORP.

Probe:Axial 10MHz

Frequency : 10MHz

Measuring Fault : 0.1mm


Resolution : 0.6mm

Time Distance : 0.1s

Running Time : 0.13us

- 1) Measurement date and time
- 2) Product Name
- 3) Clinic ID/Name
- 4) Address of clinic
- 5) Patient ID
- 6) Patient's name
- 7) Sex
- 8) Physician's name
- 9) Printing date/time
- 10) Character string defined by user
- 11) Measurement information

b) Example printout of axial length

c) Example printout of corneal thickness

(1) —	RIGHT	1) Measurement eye (right eye/left eye)
(2) —	Auto	2) Measurement mode / Type of displayed data
(3) —	Minimum	
	Velocity : 1640m/s	3) Converted acoustic velocity
	----- Actual -----	
(4) {	1 477µm CCT	4) Actual measurement in each memory slot
	2 472µm 8.0mm 0°S	5) Measurement point of measurement data
	3 490µm 8.0mm 30°S	6) Average of actual measurements
	4 479µm 8.0mm 60°S	7) Standard deviation of actual measurements
	5 485µm 8.0mm 90°S	
	6 476µm 8.0mm 120°S	8) Bias value display method and bias rate or correction
	7 470µm 8.0mm 150°S	9) Bias Value (for each memory)
	8 484µm 8.0mm 170°S	10) Average of bias values
	9 475µm 10.0mm 0°S	11) Standard deviation of bias values
(6) —	10 476µm 10.0mm 170°S	12) Intraocular pressure correction parameter
	Avg. : 478µm SD : 5.9µm	13) Average corneal thickness
(7) }		14) Intraocular pressure
(8) —	----- 70% Biased -----	15) Corrected intraocular pressure
(9) {	1 333µm CCT	
	2 330µm 8.0mm 0°S	
	3 343µm 8.0mm 30°S	
	4 335µm 8.0mm 60°S	
	5 339µm 8.0mm 90°S	
	6 333µm 8.0mm 120°S	
	7 329µm 8.0mm 150°S	
	8 338µm 8.0mm 170°S	
	9 332µm 10.0mm 0°S	
(10) —	10 333µm 10.0mm 170°S	
	Avg. : 335µm SD : 4.1µm	
(11) }		
(12) {	IOP CALCULATION	
	PARAMETER1 : 554	
	PARAMETER2 : 0.0450	
(13) —	CCT(Avg.) : 478µm	
(14) —	Meas. IOP : 25.0mmHg	
(15) —	Cor. IOP : 28.4mmHg	

d) Printout of IOL power calculation

(1)	RIGHT
	Input Parameters (US)
(2)	AXIAL : 23.84mm
(3)	ACD : 2.47mm
(4)	K1 : 30.00D
	K2 : 35.00D
(5)	Keratometric Index : 1.3375
(6)	Target Ref. : -1.00D
	(Clinical History)
(7)	PreK1 : 32.00D
	PreK2 : 34.21D
	PostK : 31.00D
	Pre REF S : 5.00D
	Pre REF C : -1.00D
	Post REF S : 2.00D
	Post REF C : 0.00D
	V.D. : 12.0mm
(8)	SRK/T
(9)	IOL-01
	MAKER 1
(10)	Aconst IOL Ref.
	117.00 13.00 1.74
(11)	Power 13.50 1.34
	16.43 14.00 0.94
	14.50 0.54
	15.00 0.14
	15.50 -0.26
	16.00 -0.66
	16.50 -1.06
	17.00 -1.46
	17.50 -1.86
	18.00 -2.26
	18.50 -2.66
	19.00 -3.06
	19.50 -3.46
	20.00 -3.86
	SRK/T
	IOL-02
	MAKER 2
	Aconst IOL Ref.
	117.00 13.00 1.74
	Power 13.50 1.34
	16.43 14.00 0.94
	14.50 0.54
	15.00 0.14
	15.50 -0.26
	16.00 -0.66
	16.50 -1.06
	17.00 -1.46
	17.50 -1.86
	18.00 -2.26
	18.50 -2.66
	19.00 -3.06
	19.50 -3.46
	20.00 -3.86

(12)


- 1) Measurement eye (right eye/left eye)
- 2) Axial length
- 3) ACD (anterior chamber depth)
- 4) Corneal refractive power or radius of corneal curvature
- 5) Cornea equivalent refractive index
- 6) Target Ref. after surgery
- 7) Parameters printed when SRK/T Double K is selected
- 8) Formula name
- 9) IOL model name, manufacturer name
- 10) Lens constant
- 11) Calculation result
- 12) IOL standard (15 levels) and estimated refractive power after surgery

e) Printout of personal constants

(1)	Personal Value	
	RIGHT	
	Input Parameters	
(2)	Post-Op Ref.	: -0.50D
(3)	Imp. Diopter	: 17.00D
(4)	AXIAL	: 24.50mm
(5)	ACD	: 3.45mm
(6)	K1	: 42.00D
	K2	: 43.00D
(7)	Aconst.	: 112.00
	(Clinical History)	
(8)	PreK1	: 45.00D
	PreK2	: 46.00D
	PostK	: 41.47D
	PreREF S	: -5.00D
	PreREF C	: -1.00D
	PostREF S	: -1.00D
	PostREF C	: -0.50D
	V.D.	: 12.0mm
	Solution	
	Haigis std	
	Personal Aconst.	= 114.99
	Hoffer Q	
	Personal ACD	= 3.86
	Holladay 1	
	Personal SF	= -0.11
	SRK II	
	Personal Aconst.	= 116.38
	SRK/T	
	Personal Aconst.	= 115.99
	SRK SHOWA	
	Personal Aconst.	= 115.17
	Shammas-PL	
	Personal Aconst.	= 114.37
	SRK/T Double K	
	Personal Aconst.	= 113.48

- 1) Measurement eye (right eye/left eye)
- 2) Post-operation eye refractive power
- 3) Power of the implanted IOL
- 4) Axial length
- 5) ACD (anterior chamber depth)
- 6) Corneal refractive power or radius of corneal curvature
- 7) Lens constant
- 8) Parameters for SRK/T Double K
- 9) Calculation result

f) Example printout for A-diagnosis

- 1) Measurement eye (right eye/left eye)
- 2) Gradation characteristics (Log / Liner / S)
- 3) Analysis method (Line / Point)
- 4) Beam direction
- 5) Acoustic velocity
- 6) Total gain
- 7) Analysis result (Ref, Gain, Obj. Gain, ΔdB)
- 8) Ref. Data waveform
- 9) Obj. Data waveform


3.8.3 Saved data management

a) Saving

Note


- **Always enter the ID number when saving the measurement data.**
- **Data saved in a USB memory device may be automatically erased after a long period of use. Consider the USB memory as a temporary storage, and make a backup for important data.**
- **Do not remove the external storage medium during data saving. This may destroy not only the data being saved, but also other data in the external storage medium.**

Touch the “Save” button (1) on each screen or step on the foot switch to save the measurement data. On the B mode screen, touching the “Save JPEG” button (2) saves only an image in JPEG format.


b) Browsing saved data

- 1) Touch the “Database” button (1) on each screen to open the Database screen (Fig. 1).

(Fig. 1)

- 2) Touch the “Patient List” button (2) to open the Patient List screen (Fig. 2) saved in the USB memory stick.

- 3) When characters are entered in the search field (3), the patient ID and name are searched and the result appears. Touch the patient information button (4) to open the Patient Information screen and check details of the patient information. Touch the title of the line (5) to change the order of the list. Touching this button sorts the data in ascending order or descending order.
- 4) Touch the "Measurement List" button (6) to open the Preservation Data List (measurement data list of selected patient ID) screen (Fig. 3).
- 5) Touch the "View" button (7) to display the view screen (Fig. 4). The data can be viewed and edited on this screen.

- 6) Touch the "Meas. List" button (8) to return to the Preservation Data List screen (Fig. 3). Touch the "Patient List" button (9) to return to the Patient List screen (Fig. 2).

c) Deleting saved data

Note ■ When any patient information is deleted, all the measurement data saved with the patient ID is also deleted, regardless of the mode being displayed. Be sure to check the contents saved before deleting them.

[Delete all on Database screen]

Refer to “d) Database setting” in “3.9.4 Connecting & Printing.”


[Deleting data on Patient List screen]

Touch the “Delete” button (1) to delete the selected patient ID and measurement data.

[Deleting data on Preservation Data List screen]

Touch the “All Delete” button (1) to delete the selected patient ID and measurement data, and return to the Patient List screen. Touch the “Delete” button (2) to delete the selected measurement data only.

3.9 System setup

Touch the “System Setup” button (1) on the Setup screen (Fig. 1) in each mode to open the System Setup screen (Fig. 2).


- 1) The system setup consists of 4 major categories.
 - General ... Time, version information, etc.
 - Measurement ... Time, version information, etc.
 - Application ... Calculation, analysis, and correction
 - Connection & Print ... Connection and printing

Touch the menu tab (2) to open the corresponding setup screen. The setup screen has a tree configuration. The displayed location is shown in the address bar (3) and you can jump to the specified screen by touching the displayed characters.

- 2) Touching the “Exit” button (4) will reflect the changes made and open the measurement screen.

3.9.1 General

Set common items here. Select the General tab.

(1) System Ver.

Displays the version information.

(2) Date & Time

Set the display format of the date and the date and time. The key pad appears when the input field is touched.

(3) Sound

Set whether to emit sounds when operating the screen.

(4) Measurement Sound

Set whether to emit sounds and the sound volume when taking measurements.

(5) Power Saving

This function automatically turns off the LCD when the instrument is not operated for a specified time. Touch the monitor screen to return to normal status.

Sets the time to automatically turn off the LCD.

5 mins / 10 mins / OFF

(6) Brightness

Set the brightness of the monitor screen.

(7) Machine No. Setting

Set the number to identify each instrument when multiple instruments of the same model are installed.

(8) ID

Open the setup screen for patient IDs.

(9) ID Necessity SAVING/EXPORT

Set how the patient ID must be handled when saving and sending examination data to Tomey Link.

User Discretion : Select this when ID numbers are always required.
Data cannot be saved or sent when the ID is not entered.

Yes : Select this when the ID is usually entered, but data can also be saved and sent without entering the ID. Data is saved and sent with the ID as "No ID" when the ID is not entered.

No : Select this when entry of the ID is optional. The transfer confirmation screen appears, but the examination data is automatically saved and sent. Data is saved and sent with the ID as "No ID" when the ID is not entered.

(10) ID INPUT

Set whether to display the patient information input screen before measurement when starting measurement for a new patient.

Yes : Pre-measurement: Displays the patient information entry screen before starting measurement.

No : Displays the measurement screen.

(11) Input Mask

Set which part of the character string should be recognized as the ID number when entering the ID number from a barcode. If it starts from "0" or "1", the characters from the very beginning are identified as the ID number.

(12) Confirmation

Select whether to open the patient information confirmation screen when saving or sending the examination data to TOMEY Link.

(13) Keyboard

Open the keyboard setup screen for patient's name and physician's name. The data can be entered with the selected keyboard.

(14) Key arrangement setting

Select the software keyboard type.

(15) Language Setting


Select the display language for the instrument.

3.9.2 Measurement

Make settings related to measurement on this screen. Select the Measurement tab.

a) Measurement Selection: B-Mode 10MHz

Touch the “B-Mode 10MHz” button (1) to make settings for 10 MHz

- Starting Condition

Set the default condition of the instrument when power is turned on and when the measurement mode is switched.

Latest Setting : The system starts with values set when power was last turned off.

Favorite Setting : The system starts with setting values shown in the setting display.

Touch the “Apply Current Setting” button (2) to show setting values on the last real-time screen in the setting display.

- Foot Switch (Two Button)

Select functions of the two-button foot switch.

All functions selected by one operation of the foot switch are executed.

Export : ON/OFF


Print : ON/OFF

Save : ON/OFF

Save JPEG : ON/OFF

- Image Quality Setting

Open the Image Quality setting screen and set preset values of each gain and dynamic range. Select Standard, High Contrast or Wide Dynamic Range (1). The corresponding control tools appear on the right. Touch the “Apply Current Setting” button (2) to apply the setting values shown on the current screen for each box. Touch the “Default” button to return to the default settings.

b) Measurement Selection: B-Mode 40MHz


Touch the “B-Mode 40MHz” button (1) to make settings for 40MHz. Settings are the same as for “Measurement Setting B-Mode 10_10MH.”

c) Measurement Setting: Axial

Touch the “Axial” button (1) to make settings for “Axial.”

(1) Default Measurement Method

Set the measurement method normally used.

Contact / Immersion

(2) Fix light

Select whether the fix light in the biometry probe is turned on or off.

ON / OFF

(3) Setting Velocity


Mean Velocity (M/V) / Sectional Velocity (S/V)

(4) Initial Gain

Set the gain when axial length mode is started and when a new patient is examined.

d) Measurement Setting: Pachy

Touch the “Pachy” button (1) to make settings for “Pachy.”


(1) Data Selection Method

Select the measurement data type to be displayed.

Latest / Minimum / Average

(2) Measure Point

Opens the measurement setup screen (Fig. 1).

(Fig. 1)

(3) PreSet1 / PreSet2

Changes the selected preset measurement point.

(4) Eye selection

Changes the preset measurement point of the selected eye.

(5) Radius input field

Specify the radius from the center of the measurement point. Touch the input field. The keypad appears. Touch the “▼” on the right of the input field. The pull down menu appears. Select CCT or enter the desired value using the keypad.

Input range: 0.1 - 15.0 mm

(6) Angle input field

Enter the angle from the horizontal axis at the measurement point. The keypad appears when the input field is touched.

Input range: 0 - 179

(7) S/I selection button

The mode alternates between “S” and “I” every time you touch the button.

(8) “Apply Current Value” button

Enters the measurement points currently being used in each field on the measurement screen.

e) Measurement Setting: A-Diag

Touch the “A-Diag” button (1) to make settings for “A-Diag.”

- Foot Switch (Two Button)

Select functions of the two-button foot switch.

All functions selected by one operation of the foot switch are executed.

Export: ON/OFF

Print: ON/OFF

Save: ON/OFF

Save JPEG: ON/OFF


3.9.3 Application

Make settings related to measurement of calculation, correction, etc. Select the Application tab.

a) Registering IOL data

Touch the IOL Constant “Setting” button to open the System Setup – IOL Constant screen (Fig. 1).

(1) Data list No.

Set the registration number of the IOL data list.

(2) “Data List” button

The IOL data list (Fig. 2) is displayed.

(3) IOL data input field

Enter each item. The software keyboard appears when the entry field is touched.

(4) “Calculated from A-Const” button

Calculates the values corresponding to SF and ACD constants. The calculated values are displayed in the respective input boxes.

(5) Warning message

Displays warnings for determination of a_0 , a_1 , and a_2 .

When the “OK” button is touched, the warning message disappears and the input field appears.

b) Selecting IOL power formula

Select the IOL formula to be used on the IOL power calculation screen. The selected formulae will be displayed in the pull-down menu on the IOL power calculation screen. Touch the “Selection” button to open the IOL Formula Selection screen (Fig. 3).

(Fig. 3)

When Haigis Standard is selected, the a_1 and a_2 setting fields are displayed. When SRK/T Double K is selected, the setting field for cornea equivalent refractive index (Keratometric Index KI) is displayed.

c) Setting IOP formula

Register the intraocular pressure correction formula to be used for the intraocular pressure correction. The IOP correction formula set here can be selected on the IOL correction screen. Touch the “Selection” button to open the IOL Formula Selection screen (Fig. 4).

(Fig. 4)

(1) Intraocular pressure correction formulas

Register the intraocular pressure correction formula to be used for the intraocular pressure correction.


	Input range
IOP correction formula name	: Up to 10 characters
Parameter 1	: 0 - 1500
Parameter 2	: 0.0000 - 1.0000

(2) Indication Unit mmHg / hPa

3.9.4 Connecting & Printing

Make settings related to connection, saving, export, and printing on this screen.

Select the Connection & Print tab.

a) Setting functions of export, print, and save buttons

Assign functions to 2 buttons for saving, export, and print provided on each screen (Fig. 1).

Button 1

Button 2

Touch the "Setting" button for Function settings for button to open the System Setup screen (Fig. 2).


Select the function to be assigned to each of two buttons. The button icon changes according to the selections.

b) Selecting and setting printer

The built-in printer, video monochrome printer, video color printer, network printer or data transfer printer can be selected. Touch the “Setting” button to open the Setup screen for each printer.

[Built-in]

(Fig. 1)

(Fig. 2)

(Fig. 3)

(1) Single Eye / Both Eyes

Set whether to print the data of the displayed eye only or the data of both eyes.

Single Eye / Both Eyes

(2) Print Form

Select the display form of printing.

Standard / Simple

(3) Header Info.

Select items to always be printed.

Inspection date and time / Version / Patient's name / Sex /

Physician's name / Printing date and time / Patient ID / User-defined character

(4) IOL Window (IOL power calculation screen)

Make settings for printing from the IOL power calculation screen.

Touch the "Setting" button to open the System Setup screen (Fig. 3).

IOL data number : 3 / 9 / 15

Appending Measurement : Yes / No
data

Select Calculation Result at : Yes / No

Printing

(5) IOP Window (intraocular pressure correction screen)

Set whether to print the measurement data on the intraocular pressure correction screen.

Yes / No

(6) Pachymetry measurement

Make settings for printing on the pachymetry screen.

Touch the "Setting" button to open the System Setup screen (Fig. 2).

Measurement Point (for each : ON / OFF
memory)

Bias Value (for each memory) : ON / OFF

Average & SD of Bias (average : ON / OFF
and standard deviation of bias
value)

[Video]

(1) B/W Reverse (black/white reverse printing)

ON / OFF

(2) Print Size

Small / Large

Small

Large


(3) "Print" button

Set the area to be printed in B-Mode

Image / Screen

[Network printer setting]

(1) IP Address

Set the IP address of the printer to be connected.

Touch the input field to display the keyboard.

(2) Paper size

L / A4 / B5

(3) Print Area


Set the area to be printed in B-Mode.

Image / Screen

(4) IP Check button

Performs a connection test.

[Data transfer printer setting]


- (1) Connection permission and selection of cable type
- (2) Login to Tomey Link
- (3) Host IP address
- (4) “Connection Test” button
- (5) Port number
- (6) Paper size
L / A4 / B5
Set this according to the printer setting.
- (7) B-Mode printing area
Set the area to be printed in B mode.
Image / Screen

*Refer to “d) PC setting” for details on functions (1) to (5).

c) PC setting

Make connection settings and DICOM connection settings to export data to a personal computer with TOMEY Link or Data Transfer installed. Touch the “DICOM Setting” button to move to the DICOM setup screen (Fig. 2).

(Fig. 1)

(Fig. 2)

(1) Connection permission and selection of cable type

LAN : Use a LAN cable for connection.

When “LAN” is selected here, all items listed on this screen need to be set.

USB : Use a USB cable for connection.

Set the Login User ID and Password. Items other than the above do not need to be set.

Disable : Not connected with PC.

(2) Auto Inquiry

Set the auto inquiry reference for the patient's ID.

(3) Login Tomey Link

Set the user ID and password to log in to the TOMEY Link. Maximum number of characters to be entered for the user ID or password is 16.

(4) Host IP address

Set the IP address of the PC to be connected. Touch the “DNS” button to specify the PC name. A DNS server must be running in the LAN environment to enable connection using the DNS.

Touch the input field to display the keyboard.

(5) Port number

A port number can be set.

Input range: 0 - 65535

Default setting: 80

(6) "Connection Test" button

Performs a connection test.

(7) DICOM - Settings

Moves to the DICOM settings screen.

(8) SCU AE Name

Title of DICOM Application Entity.

Conforms to the DICOM standard, and is used for communicating connection between AEs for bilateral certification.

When registered on the server (SCP), unauthorized access is prevented, and registered users can identify this application on the network.

(9) SCP AE Name

Set the title of the DICOM Application Entity compatible with the DICOM server (SCP).

DICOM communication is generally conducted by combinations of IP address, port number and AE title. Be sure to set this SCP AE name.

Check the DICOM Conformance Statement before setting the AE title of the communication destination.

(10) SCP IP Address

Set the IP address of the DICOM compatible server (SCP).

(11) Port Number

Set the port number used for DICOM communication.

(12) "Connection Test" button

Sends a "C-ECHO-RQ" message to the DICOM compatible server (SCP).

During DICOM communication, C-ECHO service is conducted to check connection between the server and terminal ends.

d) Database setting

The external storage media is initialized.

(1) External memory

Initialization

[About initialization function]

- The initialization function of this instrument deletes all data in the external storage media. Deleted data cannot be recovered. Check the data carefully before deleting it.
- The initialization function of this instrument only deletes the data. It does not format the external storage media. Use a PC when formatting is necessary.
- The type of format available on the UD-800 is "FAT32."

File system: FAT32、


Allocation unit size: Standard

Format option: Uncheck "Quick format."

Note: "Quick format" is not recommended as defective sectors cannot be checked.

e) Extension Export

Sets whether to output information related to the clinic and measurement techniques.

(1) Outputting

ON / OFF

(2) Clinic ID / Name

Set the clinic ID and name.

(3) Clinic Address

Set the address of the clinic.

f) My IP

(1) IP Address

Selects either of DHCP (dynamic IP) or MANUAL (static IP) in the IP Setting Type. When “Manual” is selected, set My IP Address, Sub Net Mask, and Default Gateway.

(2) IP Check button

Performs a duplicate IP address check.

Setting Example

■ **When only this instrument and the Tomey Link Server or a computer running Data Transfer are connected to the network, the following settings can be used.**

1) Check the computer IP address.

Check and record the IP address and subnet mask of the computer with DATA Transfer installed. Refer to the DATA Transfer operation manual for details.

2) Settings on the Instrument

The settings shown in the table below describe an example when the IP address of the computer with DATA Transfer installed is "192.168.2.128" and the subnet mask is "255.255.255.0." The IP setting method in this case is "Manual."

	Computer settings				Instrument settings			
(Local) IP address	192	168	2	128	192	168	2	128
	Check these on the DATA Transfer screen.				Same value as computer setting			(*1)
Subnet mask	255	255	255	0	255	255	255	0
	Check these on the DATA Transfer screen.				Same value as computer setting			
Default gateway	-----				0	0	0	0
					All "0"			
Host IP address	-----				192	168	2	128
					IP address of computer			

*1: An arbitrary number from 1 to 255 excluding numbers used by the computer (128 in this example)

This page is intentionally left blank.

4. TECHNICAL INFORMATION

4.1 IOL power calculation formula

4.1.1 Haigis optimized / Haigis Standard

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$$P = \frac{1000na}{L - d} - \frac{na}{\frac{na}{z} - \frac{d}{1000}}$$

Where:

$$Z = DC + \frac{REF}{1 - \frac{REF \cdot V}{1000}}$$

$$d = a0 + a1 \cdot ACD + a2 \cdot L \quad (ACD \neq 0)$$

$$d = (a0 - 0.241 \cdot a1) + (a2 + 0.139 \cdot a1)L \quad (ACD = 0)$$

$$a0 = 0.62467A - 72.434 \dots^*$$

$$DC = \frac{1000(nc - 1)}{RC}$$

2. Estimated refractive power after surgery (D)

$$REF_{iol} = \frac{1000(1000Y - DC \cdot X)}{V(1000Y - DC \cdot X) + 1000X}$$

Where:

$$X = d^2 \cdot P + 1000L \cdot na - d \cdot L \cdot P$$

$$Y = na (1000 \cdot na - L \cdot P + d \cdot P)$$

$$na : \text{Aqueous humor and vitreous refractive index} \\ = 1.336$$

$$nc : \text{Corneal reflective index} \\ = 1.3315$$

$$A : \text{A-Constant}$$

$$RC : \text{Average radius of corneal curvature (mm)} \\ = (r1 + r2)/2$$

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a1 : 0.4...*

a2 : 0.1...*

* The HAIGIS optimized formula uses the registered a0, a1, and a2 for calculation.

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Where:

$$R = \frac{Rx}{1 - 0.012Rx}$$

2. Target Ref. after surgery (D) when wearing glasses

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r : Radius of corneal curvature (mm)
= $\{(K1 - 1) \times 1000\}/K$

K1 : Cornea equivalent refractive index

SF : Distance from the iris surface to the optical center of the
implanted IOL (mm)

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$$\begin{aligned} \text{Ofst} &: \text{Calculated distance from the iris surface to the optical} \\ &\text{surface of the implanted IOL (including corneal thickness)} \\ &\text{(mm)} \end{aligned}$$

$$\begin{aligned} &= ACD \text{ const} - 3.336 \\ &= (0.62467A - 68.747) - 3.336 \end{aligned}$$

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= {(K1 – 1) × 1000}/K

K1 : Cornea equivalent refractive index

REFc : Refractive power at corneal refractive surface

$$= \frac{1000}{\frac{1000}{REF} - VD}$$

REF : Target Ref. after surgery (D)

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Lc = L – 0.1(L – 23)

VD : Vertex distance (mm)

= 12

4. Personal A-Constant

$$A = \frac{1}{0.5835} \left[\frac{1336 \cdot a}{2} + \frac{Lc}{2} - \frac{1336}{2} \sqrt{\left(a + \frac{Lc}{1336}\right)^2 - \frac{4}{1336} \left\{ Lc \left(a + \frac{1}{P}\right) - \frac{1336 \cdot a}{P} \right\}} - 0.05 + 64.4 \right]$$

$$a = \frac{1.0125}{Kc + AREFc}$$

AREFc: Refractive power of the eye at corneal refractive surface

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C1 : Estimated anterior chamber depth after surgery (mm)

$$= H + Ofst.$$

Ofst: Calculated distance from the iris surface to the optical surface of the implanted IOL (including corneal thickness) (mm)

$$= ACD_{const} - 3.336$$

$$= (0.62467A - 68.747) - 3.336$$

H : Height of cornea dome (mm)

$$= r_{pre} - \sqrt{r_{pre}^2 - W^2 / 4}$$

W : Calculated corneal diameter (mm)

$$= -5.41 + 0.58412LC + 0.098K$$

LC : Corrected axial length (mm)

$$L \leq 24.2 \rightarrow LC = L$$

$$L > 24.2 \rightarrow LC = -3.446 + 1.716L - 0.0237L^2$$

A : A-Constant

K_{pre} : Average corneal refractive power (D) before refractive correction surgery
 $= (K1 + K2)/2$
 r_{pre} : Average radius of corneal curvature (mm) before refractive correction surgery
 $= \{(K1 - 1) \times 1000\} / K_{pre}$
 $K1$: Cornea equivalent refractive index
 K_{post} : Average corneal refractive power (D) after refractive correction surgery
 $= (K1post + K2post)/2$
 r_{post} : Average radius of corneal curvature (mm) after refractive correction surgery
 $= \{(K1 - 1) \times 1000\} / K_{post}$
 P : Implanted IOL power (D)
 na : Aqueous humor and vitreous refractive index
 $= 1.336$
 nc : Corneal reflective index
 $= 1.333$
 REF : Target Ref. after surgery (D)
 V : Vertex distance (mm)
 $= 12$

4. Personal A-Constant

$$A = \left(-b + \sqrt{b^2 - 4ac} \right) / 2a$$

Where:

a : $0.62467^2 \alpha$
 b : $0.62467 \{2\alpha (H - 72.083)\} + \beta$
 c : $\alpha (H - 72.083)^2 + \beta (H - 72.083) + \gamma$
 α : $P (1 - nc) + 0.001P \cdot AREF \{V (nc - 1) - r\}$
 β : $P [na \cdot r + L1 (nc - 1) + 0.001AREF \{L1 \cdot r + V \cdot L1 (1 - nc) - na \cdot V \cdot r\}]$
 γ : $na [1000X - P \cdot L1 \cdot r + AREF \{0.001P \cdot V \cdot L1 \cdot r - (V \cdot X + L1 \cdot r)\}]$

$AREF$: Refractive power of eye after surgery (D)

4.2 How to calculate axial length using the biometry function

4.2.1 Phakic eye

The axial length measurement is obtained using the following formulae.

<When using the average acoustic velocity>

$$L = \frac{V \cdot t}{2}$$

L : Axial length measurement
V : Average acoustic velocity for axial length
t : Measurement time

<When using the divisional acoustic velocity>

$$L = \frac{V_{acd} \cdot t_{acd}}{2} + \frac{V_{Ins} \cdot t_{Ins}}{2} + \frac{V_{vit} \cdot t_{vit}}{2}$$

L : Axial length measurement
 V_{acd} : Anterior chamber depth acoustic velocity
 V_{Ins} : Crystal lens acoustic velocity
 V_{vit} : Vitreous acoustic velocity
 t_{acd} : Time for measuring anterior chamber depth
 t_{Ins} : Time for measuring crystal lens
 t_{vit} : Time for measuring vitreous

4.2.2 Dense cataract eye

Take the axial length measurement in the same manner as when using the average acoustic velocity for a phakic eye.

4.2.3 Aphakic eye

Take the axial length measurement in the same manner as when using the average acoustic velocity for a phakic eye.

4.2.4 IOL eye

[Default setting for acoustic velocity at delivery (changeable)]

- IOL acoustic velocity (reference)

Acrylic	: 2,200m/s
Silicone	: 1,049m/s
PMMA	: 2,718m/s

* The acoustic velocity varies depending on materials, makers, and temperature.

- Vitreous acoustic velocity : 1,532 m/s

- Anterior chamber depth acoustic velocity : 1,532 m/s

The axial length measurement is obtained using the following formulae.

$$Lp = \frac{V_{acd} \cdot t_{acd}}{2} + Th + \frac{V_{vit} \cdot t_{vit}}{2}$$

$$t_{vit} = t_{axl} - t_{acd} - \frac{2Th}{V_{IOL}}$$

Lp : Axial length measurement

Th : Center thickness of IOL (entered by user)

V_{IOL} : IOL acoustic velocity

V_{acd} : Anterior chamber depth acoustic velocity

V_{vit} : Vitreous acoustic velocity


t_{axl} : Measurement time

t_{acd} : Time for measuring anterior chamber death

t_{vit} : Time for measuring a vitreous

4.3 Amplifier characteristics of the A-scan diagnosis

The I/O characteristics of the log amplifier, linear amplifier, and S-curve amplifier selected for the A-scan diagnosis are as shown below.

4.4 Ultrasonic wave output

- MI (mechanical index)

MI is a parameter that indicates the mechanical effects of ultrasound waves on biological tissue. Mechanical effects are, for example, movement (or flow) around air bubbles compressed when an acoustic pressure wave of the ultrasonic wave passes through biological tissue, or energy released when cavitation temporarily bursts the air bubbles.

- TIS (thermal index of soft tissue)

TIS is a parameter that indicates the thermal effects of ultrasound waves on biological tissue. Thermal effects mean an increase in the tissue temperature due to emission of ultrasonic waves.

5. INSPECTION AND MAINTENANCE

5.1 Warranty

One-Year Limited Warranty

The seller warrants this product to be free from defects in material and workmanship under the normal use of this product for one (1) year or other term complying with local regulations from the date of invoice issued by Seller to the original purchaser.

Lamps, paper and other consumable items shall not be covered by this warranty. This warranty also shall NOT apply if the product has not been installed, operated or maintained in accordance with the INSTRUCTION MANUAL of Tomey Corporation (here in after called "Tomey"). Neither seller not Tomey shall be liable for any damages caused by purchaser's failure to follow instruction for proper installation, use and maintenance of product.

This warranty is only applicable to the new product and DOES NOT cover any damage resulting from or caused by accident or negligence, abuse, misuse, mishandling, improper modification of this product, by persons other than personnel duly authorized by Tomey, not to a product whose serial number or batch number is removed, altered or effaced.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ANY AND ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED (INCLUDING SPECIFICALLY, WITHOUT LIMITING THE GENERALITY OF THE FOREGOING, ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE), AND ALL OTHER OBLIGATION AND LIABILITY ON THE PART OF SELLER AND TOMEY. NEITHER SELLER NOR TOMEY SHALL BE LIABLE FOR INCIDENTAL, CONSEQUENTIAL OR SPECIAL DAMAGES UNDER ANY CIRCUMSTANCES OR FOR MORE THAN REPAIR, REPLACEMENT OR REFUND OF THE PURCHASE PRICE OF DEFECTIVE GOODS.

5.2 Operation life

This instrument is designed to have an operation life of 8 years when operated under the appropriate environment and it is adequately inspected and serviced.

5.3 Inspection

- *When there is a problem, correct waveforms or images may not be obtained correctly. Contact Tomey or your local distributor immediately for repair or replacement.*
- *Be sure to hold the connector when removing the measurement probe from the main unit. Holding and pulling the cord may damage the inner core wires.*
- *Do not touch the connector of the measurement probe, otherwise the probe may be damaged.*

5.4 Routine maintenance

- *Hold the plug when disconnecting the power cord from the outlet to avoid applying excessive force on the cord. Pulling the cord may damage the inner core wires, resulting in electric shock or fire.*

- *Do not use organic solvents such as thinner, benzene, or acetone to clean the instrument and accessories. These solvents may damage the surface of the instrument.*
- *Disconnect the power plug from the outlet and place the dust cover over the main unit when not used for long periods.*
- *Place the dust cover over the main unit when not being used.*

5.4.1 Maintenance of measurement probes

- *Be sure to hold the connector when removing the probe from the main unit. Holding and pulling the cord may damage the inner core wires.*
- *Never touch the connector, otherwise the probe will be damaged.*
- *Attach the protective cap to the probe and place it in the probe holder when not used.*

a) Appearance check

- Check that the contact section of the probe is not damaged.
- Check that the connector is firmly inserted and the cable is not damaged.

b) Cleaning

Note

■ **Use absorbent cotton to clean the contact section of the measurement probe. Otherwise, the contact section may be damaged.**

■ **Holding and pulling the cord may damage the inner core wires.**

■ **Do not clean the connecting section of the connector.**

- When cleaning the main body of the measurement probe, wipe it with absorbent cotton dampened with antiseptic soap and water or alcohol, and then with cotton dampened with water.
- When cleaning the tip of the measurement probe, wipe it with absorbent cotton dampened with alcohol, and then with cotton dampened with water.
- Wipe cords with a well wrung cloth dampened with water.
- Completely remove water after cleaning.

c) Sterilization

Note


■ **Effective chemicals vary depending on the types of germs and viruses. Select appropriate chemicals and immersion time based on the guidelines, articles, literature, etc. provided by local administrative organizations and/or societies.**

■ **Do not use EOG sterilizer for the biometry probe and eye cup.**


■ **Never use a high-pressure steam sterilizer.**

■ **Follow the instruction manual for the EOG sterilizing conditions. Otherwise, safe sterilization cannot be assured.**

- The B mode probe can be sterilized by wiping or immersing the section up to 10 mm from the contact section using the following chemicals. As for the biometry and pachymetry probes, the section up to 15 mm from the contact section can be wiped and immersed in the following chemicals. Thoroughly rinse the probe with purified water after sterilization.

- Alcohol
- Sodium hypochloride (0.5% solution)

- Immerse the whole eye cup in alcohol or sodium hypochloride solution (0.5%), rinse them with purified water to completely remove chemicals, and thoroughly dry them.

When sterilizing the pachymetry probe with EOG, remove the protective cap from the tip of the probe and put the whole probe in a sterilizing pack.

<Sterilization conditions>

EOG : Mixture of ethylene oxide gas (20%) and carbonic acid gas (80%)

Sterilization temperature : $50 \pm 5^\circ\text{C}$

Sterilization pressure : $(9.8 \pm 0.98) \times 10^4 \text{ Pa}$

Sterilization time : $5 \pm 1 \text{ hours}$

Aeration time : 10 hours or longer

5.4.2 Cleaning and Disinfection of measurement probes (for Europe)

[Specific Section for European continent]

Note

■ Follow the legal requirements and regulation in the country.

■ Refer to the instruction manual concerning test-fluid for detail.

a) Cleaning

Note

■ Holding and pulling the cord may damage the inner core wires.

■ Do not clean the connecting section of the connector.

- When cleaning the main body of the measurement probe and, wipe it with a suitable fluid eg. "Cleanisept® Wipes" and then wipe with cotton dampened with water
- Wipe the cord with a well wrung cloth dampened with water and remove water after cleaning.

b) Disinfection


Note

■ Effective chemicals vary depending on the types of germs and viruses. Select appropriate chemicals and immersion time based on the guidelines, articles, literature, etc. provided by local administrative organizations and/or societies.

■ **A high-pressure steam sterilizer is acceptable only for the immersion attachment.**

- The Ultrasound probe (biometry probe, pachymetry probe, and A-scan diagnostic probe) can be disinfected by immersing the contact section of the probe tip in suitable cleaning fluids.

©Biometer / A-diagnostic probe/Pachy probe

One clinically tested cleaning product would be:

- PERFEKTAN®TB
(4% concentration and an exposure times of 30 min)

Antimicrobial Properties Perfektan® TB is:

- bactericidal (incl. MRSA=multi resistant Staphylococcus aureus) -
tuberculocidal - fungicidal (C. albicans), levurocidal - virus inactivating incl.
HBV/HIV/BVDV/HCV/Papova - Polyoma/Adeno virus - sporicidal Thoroughly
rinse the probe with purified water and thoroughly dry it after cleaning.

5.4.3 Maintenance of main unit

Note ■ *When disconnecting cords, do not apply excessive force on the cord; for example, do not hold and pull the cord to disconnect it.*

■ *Do not touch the connection terminal surface of the main unit or measurement probe.*

- Wipe off dirt on the main unit by lightly rubbing with a well wrung cloth dampened with water and then with a dry cloth. When it is very dirty, clean it by lightly rubbing with a well wrung cloth dampened with diluted neutral detergent. Then, wipe it with a cloth dampened with water and then with a dry cloth.

Note ■ *Do not leave any water, chemicals, etc. remaining on the touch panel as the surface of the touch panel is vulnerable to moisture. The surface of the touch panel may be damaged and unable to maintain the appropriate performance.*

- Wipe off dirt on the surface of the touch panel with a dry soft cloth. When it is very dirty, lightly wipe it with a soft cloth dampened with alcohol. Further, wipe it with a dry soft cloth to completely remove moisture.

5.5 Replacing consumables


5.5.1 Fuses

■ **Disconnect the power cord from the outlet when replacing fuses. Otherwise you may get electric shock, resulting in death or serious injuries.**

■ **Use fuses specifically designed for this instrument.**

■ **When the instrument does not work correctly after fuses have been replaced, there may be other causes of the problem. Turn off the instrument immediately and contact your local distributor.**


- 1) Turn power off and disconnect the power plug from the outlet.
- 2) Disconnect the power plug from the power terminal.
- 3) Insert a flat screwdriver into the groove of the fuse holder under the power terminal and remove the holder case.
- 4) Replace the blown fuse with a new one.
- 5) Install the fuse holder in the reverse order of removal.

5.5.2 Printer paper

■ **Always use genuine TOMEY paper for the printer. Using other types of paper may cause printer failure.**

■ **Do not pull paper forcibly. Trying to pull the paper out may cause printer failure. Touch the "PRINT" button without any data stored to feed blank paper.**

Replace the paper roll according to the procedure below when red lines appear on the edges of the printer paper.

- 1) Press the button (1) to open the printer cover.
- 2) Remove the old printer paper roll and install a new roll. If the paper roll is installed in the wrong direction, nothing will be printed.
- 3) Close the printer cover with the end of the paper protruding 2 mm from the outlet. Press the cover firmly until you hear a click.

5.6 Storing

- *Install the instrument in a location free of water or chemicals. Any water or chemicals entering the instrument may cause electric shock or failure.*

- *Do not store the instrument in a location where chemicals are stored or gases may occur. Spilt chemicals or vapor may enter the instrument and result in fire.*
- *Disconnect the power cord from the outlet to ensure safety when the instrument is not operated for 1 month or longer.*
- *Store the instrument in a location not subject to direct sunlight, high temperature and humidity, or air containing dust, salt and/or sulfur. Otherwise, failure or malfunction may occur.*
- *Store the instrument in a leveled stable place free of vibration or mechanical impact. Otherwise, measurement cannot be conducted correctly. The instrument may also topple or fall, resulting in fire or a fatal accident.*
- Clean and dry the measurement probe, attach the protective cap, and place it in the probe case for storage.
- When the instrument is not used for a while, check that the instrument is operating correctly before starting operation again.

Storage conditions

- Ambient temperature range $-20 - +50^{\circ}\text{C}$
- Relative humidity range 45 - 85%

5.7 Disposal

- *Keep the box and packing materials for use when moving or transporting the instrument.*
- *Keep the packing materials and box together.*
- *When disposing of the packing materials, sort them by type and abide by local government rules and regulations.*
- *A board with a lithium battery is used in the instrument. Since handling of lithium batteries varies depending on governing bodies, follow relevant laws and local rules.*

6. TROUBLESHOOTING

Check the following first when you encounter any problems.

If the problem is not solved even after checking the applicable item listed below, contact your local distributor to request inspection and/or repair.

■ ***Do not remove the cover of the instrument. Otherwise you may get electric shock.***

■ ***Do not take any actions other than those specified below.***

6.1 Common items

- The instrument does not start when the power switch is turned on.

?Cause 1 Problem with the power plug

Solution ➡ Check that the power plug is firmly connected to the outlet. Check that there are no flaws in the power cord, such as cracks or tears.

?Cause 2 Problem with the power outlet

Solution ➡ Check that power is supplied to the outlet to which the power cord is connected.

?Cause 3 Blown fuse

Solution ➡ Check that fuses are not blown. If blown, replace the fuse (refer to “5.5.1 Fuses”). When the new fuse is blown again, the instrument may be faulty. Contact your local distributor to request inspection and/or repair.

- Nothing appears on the monitor screen.

?Cause 1 The auto power off function, which automatically turns off the screen when the instrument is not operated for the specified time, has been activated.

Solution ➡ Touch the monitor screen.

? Cause 2 The maintenance switch on the side is in the P (up) position.

Solution ➡ Turn off the power, return the maintenance switch to the down position, and then turn on the power.

- The whole monitor screen is dark and not easy to see.

? Cause 1 The brightness of the monitor is low.

Solution ➡ Adjust the brightness of the screen using the screen brightness setting described in "3.9.1 General."

- The data cannot be printed by the built-in printer.

? Cause 1 Printer paper.

Solution ➡ Check the remaining printer paper.

? Cause 2 Setting the output printer

Solution ➡ Set "built-in printer" to the destination to which output is directed when "Print button" is pushed. Refer to "3.9.4 b) Selecting and setting printer" for setting details.

? Cause 3 The printer cover is open.

Solution ➡ Check that the printer cover is completely closed.

? Cause 4 The printer paper is not set in the correct direction.

Solution ➡ Check that the printer paper is set correctly. (Refer to "5.5.2 Printer paper")

- A different button to the one touched on the touch panel becomes active.

? Cause 1 The touch panel is not calibrated correctly.

Solution ➡ Contact your local distributor.

-
- Nothing is displayed on the screen of this instrument when inputting IDs in the external ID input device.

?Cause 1 USB cable.

Solution ➡ Verify that the USB cable of the external ID input instrument is correctly connected to the main unit. For correct connection, refer to "3.1.1 g) External ID input device."

?Cause 2 A screen which does not receive an ID input appears.

Solution ➡ Entry of ID from an external ID input instrument is accepted only on a screen where the patient information is displayed. Input ID after displaying such screens.

- Data cannot be saved.

?Cause 1 The ID is not entered.

Solution ➡ The instrument uses the ID as a file name.
Enter the ID and then save the data.

?Cause 2 There is not enough space available in the memory.
(The confirmation screen appears when attempting to save the data.)

Solution ➡ Delete unnecessary data in the memory and save the data again.

6.2 B mode image diagnosis

- B mode images are dark and not easy to see.

?Cause 1 The setting of the total gain is low.

Solution ➡ Increase the total gain value.

- Ultrasound images are not clear or exact.

?Cause 1 The B mode probe is not connected correctly.

Solution ➡ Check that the B mode probe is firmly connected.

? Cause 2 Gel for ultrasound diagnosis is not used or is insufficient.

Solution ➡ Apply gel for ultrasound diagnosis to the tip of the probe. If the problem is not solved even after the protective agent is applied, apply more.

? Cause 3 Noise is generated in the peripheral area

Solution ➡ If there is any source of noise (devices such as a motor, laser surgical equipment, etc.) near the instrument, move it away from the instrument.

- The vector A mode waveform is not displayed.

? Cause 1 Functions for which cannot be used concurrently are attempted to be effective.

While the following functions are effective, the vector A mode function cannot be enabled.

- Zoom function
- Measurement function

Solution ➡ Disable the above functions.

- A probe error is displayed when starting operation or switching to the B mode image diagnostic function.

? Cause 1 The B mode probe is not connected correctly.

Solution ➡ Check that the B mode probe is firmly connected.

- No comment is displayed on images.

? Cause 1 The display setting of comments is "hidden".

Solution ➡ Turn on the display of comments as described in "3.2.1 b) Settings for items displayed in image."

- Freezing cannot be canceled when stepping on the foot switch.

? Cause 1 The foot switch is not connected correctly.

Solution ➡ Verify that the foot switch is connected firmly and locked.

-
- The instrument enters FREEZE mode although the foot switch is not pressed.

? Cause 1 The auto freezing function is activated.

Solution ➡ When there is no entering from a user in the real time mode for a long time, freeze mode is automatically entered. When you want to return to real time mode, cancel freezing by touching the “FREEZE” button or stepping on the foot switch.

6.3 Biometry function

- The axial length cannot be measured.

? Cause 1 The biometry probe is not connected correctly.

Solution ➡ Correctly connect the probe until it is locked.

? Cause 2 The measurement conditions in immersion mode are not satisfied. (In the case of immersion mode)

Solution ➡ In immersion mode, data is captured when a corneal waveform is within 2.0 - 5.0 mm from the initial waveform. Place the cornea protective agent or immersion attachment between the cornea and the tip of the biometry probe to satisfy these conditions.

? Cause 3 The biometry probe is directly applied to the cornea without using a cornea protective agent or immersion attachment. (In the case of immersion mode)

Solution ➡ Use the immersion attachment or cornea protective agent for measurement in immersion mode.

- Measurement cannot be performed in Auto mode.

? Cause 1 The mode is set to Manual.

Solution ➡ Set to Auto mode.

? Cause 2 Noise is generated in the peripheral area.

Solution ➡ If there is any source of noise (devices such as a motor, laser surgical equipment, etc.) near the instrument, move it away from the instrument.

?Cause 3 Patient's sight is unstable.

Solution ➡ Use the fixlight in the biometry probe or on the chin rest to guide the sight of the patient.

?Cause 4 The contact section of the biometry probe is damaged.

Solution ➡ If damaged, immediately stop measurement and contact your local distributor.

● Measurement cannot be performed in Manual mode.

?Cause 1 The foot switch is not connected correctly.

Solution ➡ Connect the foot switch correctly until it is locked as described in “3.1.1 c) Connecting the foot switch.”

?Cause 2 The mode is set to Auto.

Solution ➡ Set to Manual mode.

● Measurements are unstable or inappropriate.

?Cause 1 The biometry probe is not connected correctly.

Solution ➡ Correctly connect the probe until it is locked.

?Cause 2 The retina gate cursor is not set in an appropriate position.

Solution ➡ Set the retina gate cursor on the immediate left of the actual retina waveform. Make sure that there are no unnecessary waveforms between the retina gate cursor and retina waveform.

?Cause 3 The biometry probe is not applied appropriately.

Solution ➡ Apply the probe to the center of the cornea perpendicularly so that the retina waveform spikes.

?Cause 4 Too much corneal protective agent is applied. (In the case of contact mode)

Solution ➡ When too much cornea protective agent is applied, a gap is formed between the biometry probe and cornea, and the measurement result may tend to be longer. Apply appropriate amount of cornea protective agent.

? Cause 5 The immersion attachment is used. (In the case of contact mode)

Solution ➡ Select immersion mode and take measurements.

? Cause 6 The biometry probe is directly applied to the cornea without using the cornea protective agent. (In the case of immersion mode)

Solution ➡ Use the immersion attachment or cornea protective agent for measurement in immersion mode.

? Cause 7 The converted acoustic velocity is not set appropriately.

Solution ➡ Check the setting of the converted acoustic velocity. A different converted acoustic velocity can be set for right and left eyes. Check the setting for both eyes.

? Cause 8 The contact section of the biometry probe is damaged.

Solution ➡ If damaged, immediately stop measurement and contact your local distributor.

● Monitoring sound does not stop.

Data is captured even while the probe does not contact the eyeball.

? Cause 1 The tip of the biometry probe is wet.

Solution ➡ Remove drops of water.

● Monitoring sound does not go off.

? Cause 1 The measurement sound is set to “OFF.”

Solution ➡ Set the volume to an option other than “OFF” as described in “3.9.1 General.”

-
- The probe error confirmation screen appears when the power is turned on.

? Cause 1 The biometry probe is not connected correctly.

Solution ➡ Correctly connect the probe until it is locked.

- The fixlight in the biometry probe is turned off.

? Cause 1 Lighting of the fixlight is set to “OFF.”

Solution ➡ Set to ON as described in “3.9.2 c) Measurement Setting : Axial.”

6.4 Pachymetry function

- The probe error confirmation screen appears when the power is turned on or the pachymetry probe is calibrated.

? Cause 1 The pachymetry probe is not connected correctly.

Solution ➡ Correctly connect the probe until it is locked.

- The error message “Out of range!” appears and measurement cannot be performed.

? Cause 1 The measurement range setting is not correct.

Solution ➡ Touch the “Meas Range” button to set the measurement range appropriately.

- Stable measurements cannot be obtained. Appropriate measurements cannot be obtained.

? Cause 1 The pachymetry probe is not connected correctly.

Solution ➡ Correctly connect the probe until it is locked.

? Cause 2 The pachymetry probe is not applied appropriately.

Solution ➡ Apply the pachymetry probe perpendicular to the cornea at the position to be measured.

? Cause 3 The converted acoustic velocity is not set appropriately.

Solution ➡ Check that the converted acoustic velocity is set correctly. (Refer to “3.6.2 c) How to display converted acoustic velocity and bias values.”)

?Cause 4 The contact section of the pachymetry probe is damaged.

Solution ➡ If damaged, immediately stop measurement and contact your local distributor.

?Cause 5 The pachymetry probe is not calibrated correctly.

Solution ➡ Thoroughly wipe the tip of the pachymetry probe and press the sensitivity calibration button to re-calibrate the probe.

● Measurement cannot be performed in Auto mode.

?Cause 1 The measurement mode is set to Manual.

Solution ➡ Switch the mode to auto measurement mode.

● The measureable point does not move immediately when the pachymetry probe is removed from the cornea.

?Cause 1 The tip of the pachymetry probe is wet.

Solution ➡ Thoroughly remove any drops of water.

● Monitoring sound does not stop. Data is captured even while the probe does not contact the eyeball.

?Cause 1 The tip of the pachymetry probe is wet.

Solution ➡ Thoroughly remove any drops of water.

This page is intentionally left blank.

7. CONSUMABLES AND OPTIONAL EQUIPMENT

7.1 Spare parts

Ultrasound gel

Fuse

Specify the fuse type as "Fuse for UD-800."

Built-in printer paper

Specify the paper type as "TOMEY common printer (58 mm)."

7.2 Optional parts

40MHz UBM probe

Probe cap membrane

Syringe

PROBE Arm

Eye cup

Biometry probe (with case and test piece)

Immersion attachment

Pachymetry probe (with case and test piece)

A-scan diagnosis probe (with case)

This page is intentionally left blank.

8. SPECIFICATIONS

8.1 Specifications

B mode image diagnosis 10MHz B-Probe

- Frame rate
 - Basic mode : 20 frame / sec
- Maximum number of pages in a movie
 - 200 pages x 2
- Image display range
 - Standard : 35.2mm / 52° (at ultrasound velocity=1550 m/sec)
 - Wide : 48.0mm / 52°(at ultrasound velocity=1550 m/sec)
- Resolution
 - Lateral accuracy : 0.6mm
 - Axial accuracy : 0.6mm
- Distance measurement
 - [Distance accuracy between cursors]**
 - Lateral accuracy : ±0.5mm
 - Axial accuracy : ±0.5mm
 - [Minimum distance unit between cursors]**
 - Normal: : 0.055mm step(velocity=1550 m/sec)
 - Wide : 0.075mm step(velocity=1550 m/sec)
- Acoustic lines
 - 131 lines(step by 0.4°)
- Color scale
 - 256 scale level
- Scan Type : Sector scanning
- Transducer type : 2-ring
- Transducer frequency : 10MHz

Biometry / IOL power calculations

●Measurement range

■ Axial length	: 15.00 - 45.00mm
■ ACD	: 1.80 - 7.00mm
■ Lens thickness	: 2.00 - 6.00mm

●Accuracy

■ Measurement Accuracy	: ± 0.1 mm
■ Measurement Resolution	: 0.01 mm

●Factory setting of the converted velocity

■ Average velocity of the axial length	
Phakic Eye	: 1,550m/s
Dense Cataract	: 1,548m/s
Aphakic	: 1,532 m/s
■ Lens velocity	
Phakic Eye	: 1,641m/s
Dense Cataract	: 1,629m/s
■ ACD velocity	
Phakic Eye	: 1,532 m/s
Dense Cataract	: 1,532 m/s
IOL Eye	: 1,532 m/s
■ Velocity of vitreous	
Phakic Eye	: 1,532 m/s
(In case of sectional velocity)	
IOL Eye	: 1,532 m/s
■ IOL velocity	
Acrylic	: 2,200m/s
Silicone	: 1,049m/s
PMMA	: 2,718m/s
User	: 2,200m/s
■ IOL thickness	
Acrylic	: 0.80mm
Silicone	: 1.00mm
PMMA	: 0.80mm
User	: 0.80mm

- IOL power calculation

- Haigis standard
- Haigis optimized
- Hoffer [®]Q
- Holladay 1
- SRK II
- SRK/T
- SRK SHOWA
- Shammas-PL
- SRK/T Double K

- Biometry probe

- Type : Solid state
- Fix light : Built in the probe, Red LED
- Transducer frequency : 10MHz
- Tip diameter : 5.3mmφ
- Dimension / Weight : 8mmφ × 100mm / 30g

Pachymetry

- Measurement range : 150 to 1,500um

- Accuracy

- Measurement Accuracy : ±5μm
- Measurement Resolution : 1μm

- Percent bias/plus minus bias setting range

- Percent bias : 60 - 130%
- Plus minus bias : -600 to +450um

- Factory setting and setting range of the converted velocity

: 1,640m/s / 1,400 to 2,000m/s

- Pachymetry probe

- Type : Solid state
- Transducer frequency : 20MHz
- Tip diameter : 1.5mm with an angle of 45°
- Dimension / Weight : 8.8mmφ × 90mm / 40g

A-scan Diagnosis

- Measurement range : 60mm
- A-Scan diagnostic probe
 - Type : Solid state
 - Transducer frequency : 10MHz
 - Tip diameter : 6.0mmφ
 - Dimension / Weight : 8mmφ × 100mm / 30g

Dimensions and weights

- Dimension : 310(W)×326(H)×214(D) mm
- Weight : 6.0kg

Display

- TFT LCD : 10.4 inches, color touch screen

Power source

- Input Voltage : 100-240VAC
- Frequency : 50/60 Hz
- Power Consumption : 120 VA

8.2 Quantity of ultrasonic energy

Ispta.3 : 17 mW/cm² or less
 Isppa.3 : 28 W/cm² or less
 MI : 0.23 or lower

Standards for ultrasonic wave output in an eyeball based on the 510(k) guidelines of FDA

Ispta.3 : Spatial peak temporal average strength

Isppa.3 : Spatial peak pulse average strength

MI : Mechanical index

The ultrasound wave output table based on IEC60601-2-37 is as shown below.

Ultrasound wave output table (IEC60601-2-37)

Index label			MI:	TIS Scan
			10M B probe	10M B probe
Maximum index value			0.207	0.06481
Associated acoustic parameters	pr,α[Mpa]		0.648	N/A
	P[mW]		N/A	0.06671
	Min.of [Pa(zs),lta,α(zs)]		N/A	N/A
	Zs		N/A	N/A
	Zbp		N/A	N/A
	Zb		N/A	N/A
	Z at max.lpi,α [cm]		1.85	N/A
	deq(zb)		N/A	N/A
	fawf[MHz]		9.83	9.83
	Dim of Aaprt [cm]	X	N/A	0.52
		Y	N/A	0.52
Other information	Td[us]		0.131	N/A
	Prr[Hz]		2620	N/A
	Pr at max.lpi [MPa]		1.2113	N/A
	deq at max.lpi		N/A	N/A
	lpa,α at max.MI [W/ cm ²]		12.227	N/A
Operating control conditions	Control		N/A	N/A

Index label		MI			TIS Non-scan($A_{aprt} \leq 1 \text{ cm}^2$)		
		Axial	Pachy	A-Diag	Axial	Pachy	A-Diag
Maximum index value		0.175	0.136	0.129	1.8×10^{-3}	0.449×10^{-3}	3.6×10^{-3}
Associated acoustic parameters	pr, α [Mpa]	0.557	0.629	0.386	N/A	N/A	N/A
	P[mW]	N/A	N/A	N/A	0.0373	4.42×10^{-3}	0.0542
	Min.of [P α (zs),lta, α (zs)]	N/A	N/A	N/A	N/A	N/A	N/A
	Zs	N/A	N/A	N/A	N/A	N/A	N/A
	Zbp	N/A	N/A	N/A	N/A	N/A	N/A
	Zb	N/A	N/A	N/A	N/A	N/A	N/A
	Z at max.lpi, α [cm]	1.36	0.05	1.36	N/A	N/A	N/A
	deq(zb)	N/A	N/A	N/A	N/A	N/A	N/A
	fawf[MHz]	10.147	21.362	8.91	10.147	21.362	8.91
	Dim of Aaprt [cm]	X	N/A	N/A	0.4	0.19	0.48
		Y	N/A	N/A	0.4	0.19	0.48
Other information	Td[us]	0.303	0.0815	0.171	N/A	N/A	N/A
	Prr[Hz]	1200	4000	640	N/A	N/A	N/A
	Pr at max.lpi [MPa]	0.897	0.653	0.586	N/A	N/A	N/A
	deq at max.lpi	N/A	N/A	N/A	N/A	N/A	N/A
	lpa, α at max.MI [W/ cm ²]	5.427	8.893	2.606	N/A	N/A	N/A
Operating control conditions	Control	N/A	N/A	N/A	N/A	N/A	N/A

8.3 Noise

The instrument generates machine noise while operating when:

- Turning power on
- Performing self diagnosis
- Printing out data
- Using the keys
- Inserting, detecting, storing data to, and removing a USB memory
- Generating errors

8.4 Operating environment

Operate the instrument under the following environmental conditions.

- Installation : Indoors, not in direct sunlight
- Temperature : +10°C - +35°C
- Humidity : 30 – 75%
- Atmospheric pressure : 800 -1060 hPa
- Power fluctuation : less than ±10% of normal voltage

Store and/or transport the instrument in the instrument's box under the following environmental conditions.

- Temperature : -20°C - +50°C
- Humidity : 45 - 85%

8.5 Classification of ME equipment

Protection against electric shock
: Class I ME equipment

Applied parts
: B applied parts (the tip of 10MHz B probe, the tip of the 40MHz UBM probe, the tip of the biometry probe, the tip of the pachymetry probe, the tip of the A-Scan diagnostic probe)

IP Code
: IP20 (Main unit) IP41 (10MHz B probe except connector, biometry

probe except connector, pachymetry probe except connector, A-Scan diagnostic probe except connector) IP47 (the tip of 10MHz B probe, the tip of the biometry probe, the tip of the pachymetry probe, the tip of the A-Scan diagnostic probe) IP21 (40MHz UBM probe except connector, Foot switch) IP27 (the tip of 40MHz UBM probe)

Mode of Operation

: Continuous operation

8.6 Declaration of Conformity to EMC

- *When using specified accessories or cables excluding internal parts sold as repair parts of TOMEY Corporation, electromagnetic emission of the UD-800 may increase or electromagnetic immunity may decrease.*
- *When using handheld or portable communication equipment, locate it more than 30cm(12inches) away from the UD-800. It may degrade the performance of the UD-800 if the equipment is located closer.*
- *Do not use the UD-800 adjacent to or stacked on other equipment. If it needs to be adjacent to or placed on other equipment, it must be confirmed that UD-800 operates correctly in such location.*

- *The UD-800 falls under Group 1, Class B according to EN 55011(CISPR 11). This means that the UD-800 does not generate RF energy intentionally in the form of electromagnetic irradiation, inductive coupling and/or capacitive coupling for processing of materials or inspection/analysis, and that it is suitable for use in facilities directly connected to a low-voltage power network that supplies power to general household facilities and buildings used for home use.*
- *The UD-800 can be used in a specialized health care facility such as hospitals and clinics (except anywhere near the source of a strong EMI) and a general commercial facility such as eye glasses stores.*
- *The UD-800 requires special attention for EMC and needs to be installed, serviced and used based on the following information.*
- *Do not use cables other than those provided or specified by TOMEY Corporation.*
- *The person who connects additional equipment to the signal I/O section as part of the medical system, shall take responsibility for ensuring the system conforms to IEC/EN 60601-1-2 requirements.*

Cable list

Name	Specification	Length	Note	Essential
AC Power Cord	Unshielded	2.5m	Accessories	
Foot Switch Cable	shielded	1.6m	Accessories	
USB cable	shielded	1.8m	Not accessories	
USB cable	shielded	1.2m	Not accessories	
USB cable	shielded	1.8m	Not accessories	
B Probe Cable	shielded	1.6m	Accessories	
Biometry Probe Cable	shielded	1.8m	Not accessories	
PACHY Probe Cable	shielded	1.8m	Not accessories	

performance :

Display B mode image of the retina, vitreous body and anterior chamber that enable the physician to diagnose the presence of diseases.

Display the correctly measured axial length and perform IOL calculation.

Display the correctly measured corneal thickness.

Display A mode waveforms of the retina and vitreous body that enable the physician to diagnose the presence of diseases.

<EMISSIONS>

Test	Standard	Compliance
		HOME HEALTHCARE ENVIRONMENT
Conducted Emissions	CISPR 11	Group 1 Class B
Radiated Emissions	CISPR 11	Group 1 Class B
Harmonic distortion	IEC 61000-3-2	Harmonic Class A
Voltage fluctuations and flicker	IEC 61000-3-3	Complies

<IMM
UNITY
>
Table
1

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Test	Standard	IMMUNITY TEST LEVELS
		HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Radiated RF EM fields	IEC 61000-4-3	10 V/m 80MHz – 2.7GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See Table 2
RATED power frequency magnetic fields	IEC 61000-4-8	30 A / m 50Hz or 60Hz

Table 2 【Test specifications for ENCLOSURE PORT IMMUNITY to
RF wireless communications equipment】

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V / m)
385	380 - 390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
450	430 - 470	GMRS 460, FRS 460	FM ±5kHz deviation 1kHz sine	2	0.3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800 / 900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
870						
930						
1 720	1 700 - 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	2	0.3	28
1 845						
1 970						
2 450	2 400 - 2 570	Bluetooth, WLAN, 802.11 b / g / n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28
5 240	5 100 - 5 800	WLAN 802.11 a / n	Pulse modulation 217Hz	0.2	0.3	9
5 500						
5 785						

Table 3 【Input a.c. power PORT】

Test	Standard	IMMUNITY TEST LEVELS
		HOME HEALTHCARE ENVIRONMENT
Electrical fast transients / bursts	IEC 61000-4-4	± 2 kV 100kHz repetition frequency
Surges Line-to-line	IEC 61000-4-5	± 0.5 kV, ± 1 kV
Surges Line-to-ground	IEC 61000-4-5	± 0.5 kV, ± 1 kV, ± 2 kV
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V 0.15MHz - 80MHz 6V in ISM and amateur radio bands between 0.15MHz and 80MHz 80% AM at 1kHz
Voltage dips	IEC 61000-4-11	0% U_T ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
		0% U_T ; 1 cycle and 70% U_T ; 25 / 30 cycles Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	0% U_T ; 250 / 300 cycle

Table 4 【Signal input / output parts PORT】

Test	Standard	IMMUNITY TEST LEVELS
		HOME HEALTHCARE ENVIRONMENT
Electrical fast transients / bursts	IEC 61000-4-4	± 1 kV 100 kHz repetition frequency
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0.15 MHz - 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz

Table 5 【PATIENT coupling PORT】

Test	Standard	IMMUNITY TEST LEVELS
		HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0.15 MHz - 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz

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

BIO & PACHY METER

AL-4000

Read this manual thoroughly before using the instrument to ensure proper and safe operation. Contact Tomey Corporation or our local distributor if you have any questions or you encounter any problems during operation.

Note

- Always follow the operation procedures described in this manual.
- Keep this manual in a readily accessible place while operating this instrument.
- Contact our local distributor if you lose this manual.

716A9090-1G

i. Important safety information

- ***Do not install this instrument in a location where explosives or inflammable substances are used or stored. Fires or explosion may occur.***
- ***Do not remove the cover of the instrument. You may be directly exposed to high voltage sections.***
- ***Do not disassemble or modify the instrument. You may be directly exposed to high voltage sections.***


- ***Disconnect the AC cord from all the connected devices before installing and/or servicing the instrument. Otherwise, you may get an electric shock.***
- ***Sterilize measurement probes before use.***
- ***Do not use a measurement probe with a damaged contact section. Measurement error may occur, as well as the corneal or eyelid may be injured.***
- ***Immersion attachment and disposable tip are a disposable part. Do not reuse these parts. Otherwise, you may contract diseases.***
- ***You are cautioned that changes or modifications not expressly approved by the party responsible for compliance could void your authority to operate the equipment.***
- ***Do not commence service or maintenance work while the instrument is being used for a patient.***

- ***Do not place water or chemicals on the instrument. Any water or chemicals entering the instrument may cause an electric shock or failure.***
- ***Do not allow water or chemicals to splash on measurement probes except for the tips. Any water or chemicals entering the probes may cause an electric shock or failure.***
- ***Only use the specified terminal for connection of the instrument. Using another type of terminal may result in failure of the instrument.***
- ***This instrument is a measuring device specially designed for ophthalmology. Never use the instrument for other purposes.***
- ***The terminal for connecting the instrument to external devices is not isolated from the internal circuit. Inappropriate wiring may damage the internal circuit. Never touch any of these terminals and the patient at the same time. Be sure to contact Tomey or our local distributor before using the external output terminal.***
- ***When operating this instrument connected to other devices, only use devices that satisfy IEC60601-1 or equivalent safety requirements, or that conform to IEC60950 and whose power source is insulated with an insulation transformer, in order to satisfy the relevant safety requirements for medical equipment.***

- **Never mark or damage the caution labels on the instrument. The caution labels are provided at two locations as shown below.**
 - **If a label is damaged or becomes illegible, please contact Tomey Corporation or our local distributor.**
-

- *This device complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. The antenna used for this transmitter must be installed to provide a separation distance of at least 20 cm from all persons and must not be co-located or operated in conjunction with any other antenna or transmitter.*
-

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- *This is a precaution that, if unheeded, will result in a hazardous situation where there is imminent danger of serious injury or death.*

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iii. Table of contents

i. Important safety information	1
ii. How to read this manual.....	1
Outline.....	1
Symbols used in this manual.....	2
iii. Table of contents	1
1. PRIOR TO USE.....	1-1
1.1 Precautions for operation.....	1-1
1.2 Glossary	1-6
1.3 Checking package contents.....	1-8
1.4 Overview	1-9
2. NAMES AND FUNCTIONS	2-1
2.1 Front of the main unit	2-1
2.2 Back of the main unit.....	2-3
2.3 Sides of the main unit.....	2-4
2.4 Top of the main unit.....	2-5
2.5 Bottom of the main unit	2-6
2.6 Screen	2-7
2.6.1 Basic structure and common items	2-7
2.6.2 Axial Length Measurement (Biometry)	2-9
2.6.3 Corneal thickness measurement.....	2-15
2.6.4 Diagnosis in mode A.....	2-17
2.7 Auto power off function	2-20
2.8 Auto power off function	2-21
3. OPERATION PROCEDURES	3-1
3.1 Safety precautions	3-1
3.1.1 Precautions for installing the instrument.....	3-2
3.1.2 Precautions for connecting the power	3-3
3.2 Preparation before use.....	3-4
3.2.1 Connections	3-4
3.2.2 Wireless connection with UD-8000 and IOL calculation unit.....	3-12
3.2.3 Turning the power on and adjustment after turning the power on.....	3-16
3.2.4 Switching modes.....	3-18
3.2.5 Setting measurement conditions	3-19

3.2.6	Entering the patient data.....	3-19
3.2.7	Clear all measurement data (preparation for measuring a new patient).....	3-22
3.2.8	Selecting the eye to be measured	3-23
3.3	Axial length measurement function.....	3-24
3.3.1	Connecting the biometry probe	3-24
3.3.2	Setting axial length calculation method and converted acoustic velocity for eye to be measured	3-24
3.3.3	Contact/immersion mode settings	3-29
3.3.4	Operation check.....	3-29
3.3.5	Preparation for measurement.....	3-29
3.3.6	Precautions for measurement	3-30
3.3.7	Measurement	3-40
3.3.8	Checking waveforms after measurement.....	3-45
3.4	IOL power calculation.....	3-53
3.4.1	Calculation	3-53
3.4.2	Setting calculation conditions	3-53
3.4.3	Entering calculation parameters	3-53
3.4.4	Setting calculation formula.....	3-59
3.4.5	Saving measurement data.....	3-60
3.5	Pachymetry function.....	3-61
3.5.1	Connecting the pachymetry probe.....	3-61
3.5.2	Setting the data type to be displayed	3-61
3.5.3	Setting how to display the converted acoustic velocity and bias value.....	3-61
3.5.4	Setting the measurement mode/range and data display.....	3-63
3.5.5	Operation check.....	3-65
3.5.6	Preparation for measurement.....	3-65
3.5.7	Measurement	3-66
3.5.8	Checking the measurement data.....	3-73
3.5.9	Saving the measurement data.....	3-76
3.6	A-scan diagnosis	3-77
3.6.1	Connecting the probe	3-77
3.6.2	Switching to the A-scan diagnosis	3-77
3.6.3	Gain setting.....	3-79
3.6.4	Setting analysis method.....	3-80
3.6.5	Measurement	3-81
3.6.6	Various functions in FREEZE mode	3-82
3.6.7	Various functions in real time mode.....	3-85

3.6.8	Analysis function	3-86
3.7	Utilities	3-90
3.7.1	Buttons and icons on the screen	3-91
3.7.2	Operations for each setting.....	3-95
4.	TECHNICAL INFORMATION	4-1
4.1	IOL power calculation formula.....	4-1
4.1.1	SRK-II formula	4-1
4.1.2	SRK/T formula	4-2
4.1.3	HOLLADAY formula	4-4
4.1.4	Hoffer Q formula	4-6
4.1.5	HAIGIS standard formula.....	4-7
4.2	How to calculate axial length using the axial length measurement function	4-9
4.2.1	Phakic eye	4-9
4.2.2	Dense cataract eye	4-9
4.2.3	Aphakic eye	4-9
4.2.4	IOL eye.....	4-10
4.3	Amplifier characteristics of the A-scan diagnosis	4-11
4.4	Checking the software version.....	4-11
4.5	Ultrasound wave output	4-12
4.5.1	MI (mechanical index).....	4-12
4.5.2	TIS (thermal index of soft tissue).....	4-12
5.	INSPECTION AND MAINTENANCE	5-1
5.1	Warranty	5-1
5.2	Durable years	5-2
5.3	Inspection	5-2
5.4	Routine maintenance	5-2
5.4.1	Maintenance of measurement probes	5-3
5.4.2	Cleaning and Disinfection of measurement probes (for Europe)	5-6
5.4.3	Maintenance of main unit.....	5-8
5.5	Replacing consumables	5-9
5.6	Storing	5-9
5.7	Disposal.....	5-10
6.	TROUBLESHOOTING	6-1
7.	SPARE PARTS AND OPTIONAL PARTS	7-1
7.1	Spare parts	7-1

7.2	Optional parts	7-1
8.	SPECIFICATIONS.....	8-1
8.1	Specifications	8-1
8.1.1	Biometry/IOL power calculations	8-1
8.1.2	Function of pachymetry	8-2
8.1.3	A-scan Diagnosis	8-3
8.1.4	Main unit.....	8-3
8.1.5	AC adapter	8-3
8.1.6	Battery.....	8-3
8.2	Energy and other consumption	8-4
8.2.1	Influences of ultrasound energy on the human body	8-4
8.2.2	Ultrasound energy.....	8-5
8.3	Noise	8-10
8.4	Environment conditions.....	8-10
8.5	Classification of ME Equipment	8-11
8.6	Declaration of Conformity to EMC	8-12
9.	INDEX.....	9-1

1. PRIOR TO USE

Note

- *Read this manual thoroughly before using this instrument to ensure proper and safe operation.*
- *Always follow the operation procedures described in this manual.*
- *Do not place any other objects on this instrument.*
- *Check that there are no devices that generate strong magnetic field near the instrument. A strong magnetic field may cause noise and affect measurement.*

1.1 Precautions for operation

- *When measuring the axial length, fully examine the measured data for waveforms and variations. If the measurement result is doubtful, perform measurement again or another inspection to review the inspection result. If incorrect measurement data is used to select intraocular lenses, further surgery might be required.*
- *When using the IOL calculation result to select intraocular lenses, thoroughly determine the selection by also examining cataract surgery methods and other inspections.*
- *When measuring corneal thickness, the measurement accuracy may exceed $\pm 5 \mu\text{m}$ depending on the conditions of measurement or patient's corneal. Carefully consider the history of corneal diseases and surgeries, and review the inspection result by performing measurement again or another inspection if the measurement result is doubtful.*
- *Only allow qualified operators to use the instrument.*
- *Artifacts may occur in the ultrasonic measurements. If measurement values are questionable, determine the values by overall judgment (using the wave forms or results from other equipment).*

■ Precautions when installing the instrument and accessories

- *Install the instrument in a location free of water or chemicals. Any water or chemicals entering the instrument or devices may cause an electric shock or failure.*
- *Do not install the instrument in a location where chemicals are stored or gases may occur. Any chemicals or vapor entering the devices may result in fire.*
- *Check that the frequency, voltage, and allowable current (or power consumption) of the power source are appropriate. Otherwise, fire or electric shock may occur.*
- *Do not place any heavy object on the AC cord or squash the power cord. Fire or electric shock may occur.*
- *Fully insert the power plug into the outlet. Faulty contact, allowing any metal to contact the exposed terminal of the plug, or dust accumulated on the exposed terminal of the plug may result in fire or electric shock.*
- *Do not connect a device with data transmission specifications that are not compatible. Fire or electric shock may occur. Contact Tomey Corporation or our local distributor before using the instrument while connected to another device.*

- *Install the instrument in a location not subject to direct sunlight, high temperature and humidity, or air containing dust, salt, and/or sulfur. Otherwise, failure or malfunction may occur.*
- *Install the instrument in a leveled stable location free of vibration or mechanical impact. Otherwise, measurement cannot be conducted correctly. Also, the instrument may topple over or fall down, resulting in fire or a serious accident.*
- *Install the instrument in a location with ample clearance from other devices to allow smooth inspection.*
- *Check that the frequency, voltage, and allowable current (or power consumption) of the power source are appropriate.*
- *Check the power source (discharge condition, polarity, etc.).*

-
- *A magnet is built into the instrument. Do not place any electronic device, magnetic card, and other objects that may be affected by magnetism near the instrument. These may become unusable.*

■ **Precautions before operation**

- *Check the electrical contact of switches, polarity, dial setting, and meters, and that the instrument is working correctly.*
- *Check that all cables are connected correctly.*
- *Do not allow water or chemicals to splash on measurement probes, except for the tips. Any water or chemicals entering the probes may cause an electric shock or failure.*
- *Check the sections that the patient will directly touch.*
- *Check that the date set in the instrument conforms to the actual operation date and time.*

■ **Precautions during operation**

- *Do not place any container with liquid in it on the instrument. Any liquid entering the instrument may cause electric shock or failure.*
- *Be sure to touch the “New” button on the main unit to delete the measurement data for the previous patient before measuring a new patient. If new measurement is started without deleting the previous data, the measurement data for the previous patient may be included.*
- *Do not allow the patient to touch the terminals for connecting the instrument to external devices.*

- *Do not press on the instrument from the top. Doing so may cause failure of the instrument.*
- *Do not allow water or chemicals to splash on probes, except for the tips. Any water or chemicals entering the probes may cause an electric shock or failure.*
- *Be extremely careful not to take too much time or too many measurements, as this can stress the patient.*

-
- *Constantly observe that both the instrument and patient are free of problems.*
 - *If a problem with the instrument or the patient occurs, take appropriate action, such as stopping the device, to ensure the patient's safety.*
 - *Do not allow the patient to touch the instrument.*
 - *If any smoke, offensive odor, or abnormal sound occurs, turn off the instrument immediately, disconnect the power plug from the outlet, and contact our local distributor or Tomey Corporation.*

■ **Precautions after operation**

- *Do not place any container with liquid in it on the instrument. Any liquid entering the instrument may cause electric shock or failure.*

- *Do not use organic solvents such as thinner, benzene, or acetone to clean the instrument. Fire or electric shock may occur. (These solvents may also corrode the resin or coating on the cover of the instrument.)*
- *Follow the specified procedures to return the operation switch, dial, etc. to their original positions and turn the instrument off.*
- *Hold the plug when disconnecting the power plug from the outlet to avoid placing excessive force on the cord. Pulling the cord may damage the inner core wires, resulting in electric shock or fire.*
- *When disconnecting cords, do not apply too much force to them, for example, do not hold and pull the cord.*
- *Refer to “5.6 Storing” for instructions on the storage of the instrument.*
- *Clean the instrument at the end of operation in preparation for the next use.*
- *Clean and neatly arrange the accessories and cables.*

-
- ***If any failure occurs in the instrument, immediately stop operation, identify the failure in the instrument, and contact our local distributor for repairs.***

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- *Use the AC cord provided with the instrument or one specified by Tomey to ensure safety. Also, do not use the accessories provided with the instrument for other equipment.*

- *Conduct regular inspections of the instrument and components. When the instrument is not used for 1 month or longer, check that the instrument is operating correctly and safely before starting operation. Refer to “5.3 Inspection” in this manual for the relevant procedures.*

1.2 Glossary

[ACD]	:	Anterior chamber depth
[Acrylic]	:	Acrylic lens
[Aphakic]	:	Aphakic eye
[AVG]	:	Average value
[Axial]	:	Axial length measurement (biometry)
[A-Diag]	:	A-scan diagnosis
[BD address]	:	Bluetooth® Device Address Number set to each device to identify the connection
[Bluetooth ®]	:	Registered mark of Bluetooth SIG, Inc.
[Calibration]	:	Calibration of the pachymetry probe
[CCT]	:	Corneal center thickness
[Contact]	:	Method of taking measurements by directly applying the contact section of the biometry probe to the center of the corneal
[D]	:	Diopter The unit of refractive power representing the level of myopia, hyperopia, or astigmatism. Reciprocal of the focal length measured in meters.
[FREEZE]	:	Status where ultrasound waves are not transmitted
[Immersion]	:	Method of taking measurements by applying corneal protective agent between the corneal and contact section of the biometry probe, or attaching a special tool to the eye and using water between the corneal and contact section.
[IOL]	:	Intraocular lens
[K1]	:	Radius of corneal curvature of minor meridian [mm]
[K2]	:	Radius of corneal curvature of major meridian [mm]
[Obj.]	:	Section to be measured
[Pachy]	:	Corneal thickness measurement
[Pre]	:	Before refractive power correction surgery
[Phakic]	:	Phakic eye
[PMMA]	:	PMMA lens
[Post]	:	After refractive power correction surgery
[Ref.]	:	Reference section
[Silicone]	:	Silicone lens
[US]	:	Ultrasound biometry

[Auto power off function]	:	The power turns off automatically when no operation is conducted for a specified time.
[Caliper]	:	Manually moves the measurement point.
[Gain]	:	Adjusts amplification and amplitude of waveforms.
[Gate]	:	Specifies the range of waveforms to be detected.
[Measurement point]	:	Indicates the point where corneal thickness is measured by diameters and angles.
[Touch panel]	:	Allows you to make various settings, etc. by directly touching the monitor.
[Bias value]	:	Value obtained by converting actual measurements using a specified percentage or correction value

1.3 Checking package contents

Open the package and check that the required quantity of the following items is included and they are not damaged. If any item is missing or damaged, contact our local distributor as soon as possible.

■ ***Keep the box and packing materials for use when moving or transporting the instrument.***

- Main unit BIO & PACHY METER AL-4000..... 1
- Biometry probe (with protective cap and the case)..... 1
- Footswitch..... 1
- AC adapter 1
- AC cord..... 1
- Biometry Test piece 1
- Immersion attachment (20pcs/box)..... 1
- USB cable..... 1
- Sleeve for eye axis probe 1
- Probe holder 1
- Touch pen 1
- Dust cover 1
- Instruction manual (this manual) 1
- Installation manual for TB-1000 1
- Installation USB flash drive for TB-1000 1

1.4 Overview

- This instrument is a device used in ophthalmology to measure the distance between living tissue, utilizing ultrasound waves generated from the ultrasound transducer built into the measurement probe. This device is designed to measure the axial length and corneal thickness at optical clinics. The A-scan diagnosis allows you to identify areas with pathological changes by observing mode-A waveforms.
- Transmission of ultrasound waves is stopped when waveforms freeze.
- The measurement probe connected to this instrument transmits ultrasound waves to the inside of the eyeball and then receives ultrasound echoes reflected by the tissue boundary surface in the eyeball.
- The time required for the ultrasound wave to be transmitted from the front of the corneal to the measurement section is measured, and the length of the living tissue is calculated based on the preset converted acoustic velocity.

$$L = \frac{V \cdot t}{2}$$


L: Length of living tissue
V: Converted acoustic velocity
t: Measurement time

- This instrument is equipped with an automatic measurement function and sound monitor function, as well as functions for measuring axial length and corneal thickness, for ease of operation by the physician.
- This instrument adopts the SRK/T II, SRK/T, HOLLADAY, Hoffer Q, and HAIGIS standard as the IOL power formula, and the IOL power can be calculated immediately after measuring axial length.
- The measurement data can be saved in the internal memory.

This page is intentionally blank

2. NAMES AND FUNCTIONS

2.1 Front of the main unit

(1) LCD and touch panel

Measurements, waveforms, and IOL calculation results are displayed on the 3.5" TFT color LCD (320 x 240 dots). Operations are performed by touching buttons on the touch panel attached over the LCD.

(2) Pachymetry probe holder

Holder for the pachymetry probe

(3) Probe holder (small)

Holder for the biometry/A-scan diagnosis probe

(4) "Mode" button

Switches the axial length measurement function and pachymetry function.

(5) Power LED

Turns on when power is supplied to the main unit.

Lit orange while charging.


(6) Power switch

Turns the instrument on and off.

(7) Gain adjustment knob

Knob to adjust the gain

2.2 Back of the main unit

(1) Battery cover

Open this cover when replacing the built-in battery.


(2) Stand

Stand to hold the main unit at an angle

(3) Touch pen

Used when operating the touch panel.

2.3 Sides of the main unit

(1) AC adapter/footswitch connector

Connect the AC adapter and footswitch here.


(2) Maintenance switch

Our service personnel use this switch for maintenance.

(3) Cable hook

Wind the probe cable that can be stored in the main unit here.

2.4 Top of the main unit


(1) USB connector

Connect the personal computer and IOL calculation unit here.

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(2) Biometry probe connector

Connect the biometry probe here.


(3) Corneal thickness/A-scan diagnosis probe connector

Connect the pachymetry probe or A-scan diagnosis probe here.

2.6 Screen

2.6.1 Basic structure and common items

The dark area varies depending on operation modes.

(1) Inspected eye

Touch the “R” or “L” button to display the measurement screen for the selected eye (right or left).

(2) Patient information field

Displays the patient ID entered on the patient information entry screen.

(3) “ID” button

Opens the patient information entry screen.

(4) “New” button

Deletes all measurement data and patient information for both eyes, and returns to the initial measurement screen.

(5) Communication status button

Displays the communication status as an icon.

Communicating via wireless access	:		
Communicating via USB	:		
Communication stopped	:		
Communication not possible	:		* Operations not possible.

Utility operations are required for connection.

(6) "Save" button

Saves the measurement data in the AL-4000 measurement unit.

(7) "Utility" button


Opens the utility screen.

(8) Battery indicator

Shows the battery status.

2.6.2 Axial Length Measurement (Biometry)

a) Measurement screen

(1) Gain indicator

Displays the gain.

(2) Mode display

Displays the mode set on the utility screen (contact or immersion).

(3) "Setting" button

Displays the measurement condition setting menu.

(4) Axial length data display

Displays the average, standard deviation, and data range of the axial length.

(5) "Retake" button

Deletes the measurement data currently displayed and opens the measurement screen for the same eye.

(6) Level cursor/level line

Measurement data will be taken when the waveform rises above this cursor/line position.

(7) Measurement cursor/measurement line

The distance at this cursor/line should be taken as the measurement data.

(8) Average acoustic velocity

Displays the average acoustic velocity when calculating the axial length.

(9) Cursor button

Switches the active gate cursor.

Corneal -> LensF (front of lens) -> LensR (back of lens) -> Retina

(10) Gate cursor movement button

Moves the active gate cursor.

(11) Gate cursor for the front of lens

The waveform on the right of this cursor position is measured as the waveform of the front of the lens.

(12) Gate cursor for the rear of lens

The waveform on the left of this cursor position is measured as the waveform of the rear of the lens.

(13) Retina gate cursor

The waveform on the right of this cursor position is measured as the waveform of the retina.

(14) "Meas (measurement)" button

Displays the measurement screen.

(15) "Edit" button

Opens the edit screen.

(16) "IOL" button

Opens the IOL power calculation screen.

b) Edit screen (data list)

(1) Measurement data display

Measurement data of the axial length, anterior chamber depth, and crystal lens thickness, and their averages are listed. The following marks are assigned to the measurement data.

“*”: Data used for calculating IOL

L : Longest axial length

S : Shortest axial length

C : Calipered data

“!” : Indicates only when wave form appears at more than 1mm away from Retina spike wave form.

(2) “Retake” button

Deletes only the displayed measurement data and takes another measurement.

(3) Eye type to be measured

Displays the current eye type to be measured.

(4) “Wave Form” button

Opens the edit screen (waveform).


(5) “Del/Ret” button

Deletes or retakes the data.

(6) "Sel (select)" button

Selects the data at the cursor position to be used for calculating the IOL. "*" is assigned to the selected data.

c) Edit screen (waveform)

(1) Memory number selection key

Touch these buttons to select the memory number of the measurement data.

(2) Waveform, memory number, mean velocity, and eye type of the selected measurement data

The waveform, memory number, mean velocity, and eye type of the selected measurement data is displayed.

(3) Level cursor/level line

Measurement data will be taken when the waveform is above this cursor position in auto measurement mode.

(4) Measurement cursor/measurement line

The distance at this measurement position should be the measurement data. The measurement point is indicated with “●.”

(5) “Gate” button

Opens the gate change screen.

(6) “Calip (caliper)” button

Opens the caliper screen.

(7) “Return” button

Returns to the edit screen (data list).

d) IOL calculation screen

The screenshot shows the IOL calculation screen with the following fields and values:

- ID:** 1234
- AXIAL:** 24.90
- Desired Ref.:** -1.00
- K1:** 43.00
- K2:** 44.00
- SRK / T:** (selected formula)
- A-Const (SRK/T):** 110.00
- IOL Power:** 11.66
- IOL Ref Table:**

IOL	Ref
10.50	0.15
11.00	-0.34
11.50	-0.83
12.00	-1.34
12.50	-1.85

Callouts (1) to (4) point to the following fields:

- (1) Parameter field (AXIAL, K1, K2, Desired Ref.)
- (2) A-Const (SRK/T)
- (3) IOL Power
- (4) IOL Ref Table

(1) Parameter field

Enter the following parameters used for calculating IOL power.

Axial length / anterior chamber depth / K1 (mm or D) / K2 (mm or D) / Desired refractive power.

The anterior chamber depth is displayed only when “HAIGIS standard formula” is selected.

(2) A-Const (lens constant)

Enter the appropriate lens constant for the formula.

(3) IOL Power


The IOL power that satisfies the “desired refractive power after surgery” is displayed.

(4) IOL standard and estimated refractive power after surgery

The IOL standards based on calculation results and estimated refractive power after surgery to implant the IOL are listed.

2.6.3 Corneal thickness measurement

a) Measurement screen

(1) Setting display

The bias value and types of data to be displayed are listed.

(2) Average and standard deviation of corneal thickness

The average and standard deviation of the corneal thickness data is displayed.

(3) Measurement data display

The memory number, measurement data, and measurement point are displayed.

(4) "Setting" button

Shows the measurement condition setting menu.

(5) "Retake" button

Deletes the measurement data currently displayed and opens the measurement screen for the same eye.

(6) "Wave Form" button

Displays the waveform of the selected memory number.

(7) "Del/Ret" button

Deletes or retakes the data.

(8) "Meas Point (measurement point)" button

Shows and hides measurement points.

(9) "Pre 1"/"Pre 2" buttons

These buttons appear when the "Meas Point" button is touched and allow you to show the contents preset for the selected measurement point in the measurement data field.

(10) "M/E" button

Displays the measurement/edit screen.

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
(12) "IOP" button

Displays the intraocular pressure correction screen.

2.6.4 Diagnosis in mode A

a) Measurement screen (real time)

(1) Gain display

The gain setting is displayed.

(2) Log (gradation characteristics)

The gradation characteristics are displayed.

(3) "Setting" button

Shows the measurement condition setting menu.

(4) Analysis tool

Analyzes waveforms of the measurement data.

(5) Converted acoustic velocity

The acoustic velocity for conversion is displayed.

(6) Probe position

The point to which the probe is applied and the direction of ultrasound beams are displayed.

(7) Measurement waveform

The mode-A waveform is displayed.

(8) "Save" button

Saves the measurement data in the AL-4000 measurement unit.

(9) "FREEZE" button

Changes the mode from real time to freeze.

(10) Probe display

The probe type currently applied is shown.


(11) Scale 1

The 20-mm point is indicated with a 1-dot-wide yellow line.

(12) Scale 2

The 40-mm point is indicated with a 2-dot-wide yellow line.

b) Measurement screen (FREEZE)

(1) Page button and page number display

Changes each data of the captured mode-A waveform. The page number currently displayed is shown.


(2) "Edit" button

Displays the edit menu.

(3) "Release" button

Changes the mode from freeze to real time.

2.7 Auto power off function

(Fig. 1)

This function automatically turns off the LCD when the instrument is not operated for a specified time. Touch the monitor to return to normal status. When the AC adapter is not connected, the power is automatically turned off when the specified time has passed after the auto power off function was activated. Turn on the power to resume measurement. Refer to “3.2.3 Turning the power on and adjustment after turning the power on.”

2.8 Auto power off function

Caution

Stand-by

TYPE B APPLIED PART

Foot switch

Operating instructions

This page is intentionally blank

3. OPERATION PROCEDURES

3.1 Safety precautions

- *Sterilize measurement probes before use.*
- *Do not use a measurement probe with its contact section damaged. Measurement error may occur, additionally the corneal or eyelid may be injured.*
- *The immersion attachment is a disposable part. Do not reuse it. Otherwise, you may contract diseases.*

- *This instrument is a measurement device specifically designed for ophthalmology. Never use the instrument for other purposes.*
- *Only use the specified terminal for connection of the instrument.*
- *Attach the protective cap to the measurement probe and place it in the probe holder when not in use.*

3.1.1 Precautions for installing the instrument

- *Install the instrument in a location free of water or chemicals. Any water or chemicals entering the instrument may cause an electric shock or failure.*

- *Do not install the instrument in a location where chemicals are stored or gases may occur. Spilt chemicals or vapor may enter the instrument and result in fire.*

- *Install the instrument in a location not subject to direct sunlight, high temperature and humidity, or air containing dust, salt, and/or sulfur. Otherwise, failure or malfunction may occur.*

- *Install the instrument in a leveled stable location free of vibration or mechanical impact. Otherwise, measurement cannot be conducted correctly. Also, the instrument may topple over or fall down, resulting in fire or a serious accident.*

- *Install the instrument so that the user can use the system under optimal conditions. Carefully connect the devices so that the wiring is not disconnected unintentionally during operation, and does not hinder operation of the instrument.*

3.1.2 Precautions for connecting the power

- *Check that the frequency, voltage, and allowable current (or power consumption) of the power source are appropriate. Otherwise, fire or electric shock may occur.*
- *Do not place any heavy object on the AC cord or squash the power cord. Fire or electric shock may occur.*
- *Fully insert the power plug into the outlet. Faulty contact, allowing any metal to contact the exposed terminal of the plug, or dust accumulated on the exposed terminal of the plug may result in fire or electric shock.*
- *The power plug completely isolates the instrument from the commercial power source. If there is a problem with the instrument, turn off the power switch and disconnect the power plug. Install the instrument in a place where this can be performed smoothly.*
- *Charge the instrument in a stable location. When the instrument is placed on an angle or in an unstable location, the instrument may fall, resulting in injuries or breakage.*


3.2 Preparation before use

3.2.1 Connections

- Note** ■ *Faulty connection of the instrument and relevant accessories may result in fire, electric shock, or fatal accident. Turn off all devices and disconnect the power plugs from the outlet before connection.*
- *The connector needs to be inserted in a specific direction. Check the direction and firmly insert the connector.*

a) Connecting the AC adapter


- Note** ■ *Always use the specified AC adapter.*



- 1) Insert the AC adapter connector (1) to the AC adapter/footswitch connector (2) on the side of the instrument.
- 2) Connect the plug (3) to the power outlet.

b) Replacing the probe holder

The probe holder (small) can be replaced with the probe holder (large) provided with the instrument. The probe holder (large) is for the biometry probe to which the probe sleeve is attached.


- 1) Remove the screws (2) that secure the probe holder (small) (1).

- 2) Secure the probe holder (large) (3) on the side of the instrument with screws so that the arrow (4) on the back of the holder points upward (toward the LCD).


c) Connecting the biometry probe

(Fig. 1)

Align the red marks and insert the connector of the biometry probe (1) to the biometry probe connector (2) at the bottom of the instrument until a click is heard.


d) Attaching the probe sleeve

(Fig. 1)

- 1) Hold the probe sleeve so that the slit (1) is aligned with the probe and slide it to cover the biometry probe (2) from the connector side.
- 2) Align the slit of the probe sleeve with the groove of the biometry probe (3) and fully insert the probe into the sleeve.


e) Connecting the footswitch

(Fig. 1)

- 1) Insert the footswitch connector (1) to the AC adapter/footswitch connector (2) on the side of the instrument.

[When using the AC adapter and footswitch together]

(Fig. 2)


- 2) Insert the connector of the AC adapter (1) to the AC adapter connector (2) of the footswitch.
- 3) Turn the connector counterclockwise and clamp the cable with the fixture (3) to fix it.

f) Connecting the corneal thickness/A-scan diagnosis probe

Note

■ *When using the pachymetry probe with a disposable tip, pay extreme attention not to allow air to enter the disposable tip. If air enters the tip, measurement cannot be performed or correct measurement cannot be obtained.*

■ *When using the pachymetry probe with a disposable tip, press the “Calib (calibration)” button and calibrate the sensitivity of the probe after the disposable tip is replaced. Otherwise, measurement cannot be performed or correct measurement cannot be obtained.*

(Fig. 1)

- 1) Align the red marks and insert the connector of the corneal thickness/A-scan diagnosis probe (1) to the corneal thickness/A-scan diagnosis probe connector (2) at the bottom of the instrument until a click is heard.

g) Connecting UD-8000, IOL calculation unit, and personal computer

Note ■ *Check that either or both of the connected devices is in new patient mode before starting to communicate the patient information and measurement data. Communication cannot be started while both devices contain the information and measurement data of the last patient.*


This instrument can be connected to the UD-8000, IOL calculation unit, and personal computer.

Install the AL-4000 PC KIT (optional) software on the personal computer when connecting the instrument to the personal computer.

USB connection (wired) and wireless connection are available.



Direct connection (built-in connection) is also available to connected the UD-8000.

[Wire connection]

1) Insert the USB cable plug A (3) into the USB connector of the mating unit in the correct direction.

2) Insert the USB cable plug B (4) into the USB connector (2) on the top of the AL-4000.


[Wireless connection]

Complete settings according to “3.7.2 c) Wireless Communication Setting.”

[Built-in connection]

Directly connect the instrument to UD-8000.

(Fig. 1)

Slowly push and slide the instrument to the back.

When the instrument is assembled in the UD-8000 main unit, power is supplied to the instrument via the UD-8000 power terminal (1).

Slowly pull out the instrument toward you to remove it.

The measurement data is communicated via wired connection or wireless connection. Refer to “Wired connection” and “Wireless connection” in this section for details.

3.2.2 Wireless connection with UD-8000 and IOL calculation unit

■ *This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:*

- Reorient or relocate the receiving antenna.*
- Increase the separation between the equipment and receiver.*
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.*
- Consult the dealer or an experienced radio/TV technician for help.*

a) Wireless communication setting

Complete settings referring to “3.7.2 c) Wireless Communication Setting.”

b) Connection for wireless communication

- Note** ■ *Check that both or either of the connected devices is in new patient mode before starting to communicate the patient information and measurement data. Communication cannot be started while both devices maintain the information and measurement data of the last patient.*
- *When connecting a device with a USB cable during wireless communicating, wireless connection is changed to wired connection.*

[Automatic connection upon startup]

Wireless communication automatically starts with the first device that establishes a connection when this instrument is turned on. Refer to “3.7.2 c) Wireless Communication Setting” for how to set “automatic connection upon startup.” Refer to the instruction manual of each device for settings of the device.


[No communication state]

Note ■ *The confirmation screen appears when the data between the partner device and this instrument does not match.*

- 1) The confirmation screen appears if measurement data remains in the partner device when establishing a connection in new patient mode. (Fig. 1)
Click the “OK” button (1) to delete the data in the partner device and establish connection.
Click the “Cancel” button (2) to stop communication. (Fig. 2)


(Fig. 1) (1) (2)

(Fig. 2)

[Resuming communication from the no communication state]

- 1) Click the “New” button of the partner device to clear the data in that device.


(Fig. 1)


- 2) Select the same measurement mode as that of the partner device in the “mode selection field.”
- 3) Click the Communication check” button (1) to resume communication.

3.2.3 Turning the power on and adjustment after turning the power on

a) Turning the power on


- 1) Turn on the power switch (1) on the front of the main unit.

(Fig. 2)

- 2) The startup screen (Fig. 2) appears and allows you to make adjustments. When the pachymetry probe or A-scan diagnosis probe are connected together with the biometry probe, first priority is given to adjustments for the mode used when the power was last turned off.

Wireless communication starts only when the use of wireless communication is permitted and the device set in the utility is found.

(Fig. 3)

- 3) The patient's eye selection screen (Fig. 3) appears when adjustments are completed correctly. When an eye is selected, the measurement screen for the current mode appears.

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
The brightness of the monitor can be adjusted according to the illumination in the examination room.

Also, the sound volume can be adjusted. (Refer to "3.7.2 a) Common Setting.")


3.2.4 Switching modes

Note ■ *When the mode is changed, all data captured in the previous measurement mode is deleted. Pay careful attention when changing the mode.*

(Fig. 1)

- 1) Press the MODE button (1) on the front of the main unit for two seconds to change modes. Modes are changed from the biometry mode to pachymetry mode.

(Fig. 2)

- 2) When the mode is changed, the patient's eye selection screen (Fig. 2) appears. Select an eye and start new measurements.

3.2.5 Setting measurement conditions


Refer to “3.7 Utility” and set the measurement conditions.

3.2.6 Entering the patient data

Note ■ *If the patient information and measurement data already exists, make sure that both the new patient information and existing measurement data belong to the same patient before replacing the existing patient information with the new data.*


■ *The patient's ID can be entered only when the patient information entry screen is displayed.*

■ *Also, entry of the patient's ID is accepted only when this instrument is used independently.*

(Fig. 1)

- 1) Touch the “ID” button in the upper left of the measurement screen in each mode to open the patient information entry screen (Fig. 2).


(Fig. 2)

(7)

(8)

- 2) Enter the ID (3) and name (4) using the software keyboard (2). Periods (.) and spaces are not acceptable for the ID number.
- 3) Touch "Male" button or "Female" button to set the gender (5).
- 4) Enter the patient's name (6) using the software keyboard (2).
- 5) After all data is entered, touch the "OK" button (7) to set the entered data and return to the previous screen (Fig. 1). Touch the "Cancel" button to ignore the entered data and return to the previous screen (Fig. 1).
- 6) If an ID number is entered while the patient information and measurement data already exist, the confirmation screen (Fig. 3) appears.

(Fig. 3)

(9)

(10)


-
- 7) Touch the “OK” button (9) to retain the measurement data, overwrite the existing patient information with the newly entered data, and return to the previous screen (Fig. 1).
Touch the “Cancel” button (10) to ignore the entered data and return to the previous screen (Fig. 1).

3.2.7 Clear all measurement data (preparation for measuring a new patient)

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■ *Be sure to touch the “New” button to delete all the measurement data for the previous patient before measuring another patient. If the data of the next patient is captured without touching the “New” button, the patient information will not match the captured data.*

■ *The “New” button is effective only when this instrument is used independently.*

(Fig. 1)

When the “New” button (1) in the upper right of the screen in each mode is held for approximately 1 second until a beep sounds, the patient information (ID, name, and sex) and measurement data is all deleted and the measurement screen for a new patient appears.

3.2.8 Selecting the eye to be measured

(Fig. 1)

Touch either the “R (right)” (1) or “L (left)” (2) eye display button.

3.3 Axial length measurement function

3.3.1 Connecting the biometry probe

Connect the biometry probe to the instrument. Refer to “3.2.1 c) Connecting the biometry probe.”

3.3.2 Setting axial length calculation method and converted acoustic velocity for eye to be measured

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Select the eye to be measured from the following 4 options according to the patient's eye.

Select the material of the implanted IOL in the case of an eye with IOL.

■ Phakic eye

Select this for a normal eye or when the crystalline lens nucleus is comparatively soft, such as an eye in the initial stage of cataracts.

■ Dense cataract eye

Select this when the crystalline lens nucleus is comparatively hard, such as an eye with dense cataracts, and when it is difficult to take measurements in the mode for an eye with initial stage cataracts using an echo reflected from the back of the crystal lens.

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IOL eye: Acrylic / Silicone / PMMA / user

Select this when an IOL is implanted in the patient's eye.

“Anterior chamber depth” and “lens” are not measured for an aphakic eye and “lens” is not measured for an IOL eye.

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
Eye to be measured	Axial length	Anterior chamber depth	Lens
Phakic eye (Average acoustic velocity, divisional acoustic velocity)	○	○	○
Dense cataract eye	○	○	△
Aphakic eye	○	*	*
IOL eye (All materials)	○	○	*

○: Measurements are displayed.


△: Measurements may not be displayed.

*: Measurements are not displayed

(Fig. 1)

1) Touch the "Setting" button (1) to open the setting screen (Fig. 2).


(Fig. 2)

- 2) Touch one of the “MODE” buttons (2) to select the measurement mode.
 - Auto
Select this for normal measurements.
 - Auto quick
Select this when measurement is difficult.
 - Manual
Select this when measurement is difficult in Auto and Auto quick modes.
- 3) Touch any of the “Eye Type” buttons (3) to set the eye to be measured.
- 4) Touch the “CLOSE” button (4) to return to the capturing screen (Fig. 1) when the converted acoustic velocity is not changed. Touch the “Change Setting” button (5) to open the change setting screen (Fig. 3) to change the converted sonic velocity.
- 5) Similarly, edit the setting in the change setting screen (Fig. 3) to change the IOL materials. However, only the “User” can be changed.

The screenshot shows a 'Setting change' dialog box with two panels. The left panel is for the right eye (R) and the right panel is for the left eye (L). Both panels have the same settings: Eye Type: Phakic, Avg.V: 1550, LENS.V: 1641, Vit.V: 1532. Below these settings are buttons 'Incorporate in L' and 'Incorporate in R'. At the bottom of the dialog are three buttons: 'Initial Setting', 'Registration', and 'Cancel'.


(Fig. 3)

- 6) Entry boxes for the acoustic velocities required for the measurement conditions currently set appear. Enter the converted acoustic velocity for the corresponding eye type using the keypad (6) that appears when the entry box is touched.

(Fig. 4)

- 7) Touch the "Ent (enter)" key to apply the setting and proceed to the next entry box.

(Fig. 5)

(9) (7) (10) (8) (11)

- 8) Touch the "Incorporate in L" button (7) or "Incorporate in R" button (8) to apply the data entered for one eye to the other eye. Touch the "Initial Setting" button (9) to enter the initial settings.

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■ Phakic eye

Average acoustic velocity for axial length	: 1,550 m/s, 1,500 - 1,600 m/s
Crystal lens acoustic velocity	: 1,641 m/s, 1,540 - 1,740 m/s
Anterior chamber depth acoustic velocity	: 1,532 m/s, 1,430 - 1,630 m/s
Vitreous acoustic velocity (divisional velocity only)	: 1,532 m/s, 800 - 2,000 m/s

■ Dense cataract eye

Average acoustic velocity for axial length	: 1,548 m/s, 1,500 - 1,600 m/s
Crystal lens acoustic velocity	: 1,629 m/s, 1,540 - 1,740 m/s
Anterior chamber depth acoustic velocity	: 1,532 m/s, 1,430 - 1,630 m/s

■ Aphakic eye

Average acoustic velocity for axial length	: 1,532 m/s, 1,430 - 1,630 m/s
--	--------------------------------

■ IOL eye

Implanted lens acoustic velocity	(Acrylic)	: 2,200 m/s, 800 - 3,000 m/s
	(Silicone)	: 1,049 m/s, 800 - 3,000 m/s
	(PMMA)	: 2,718 m/s, 800 - 3,000 m/s
Anterior chamber depth acoustic velocity		: 1,532 m/s, 1,430 - 1,630 m/s
Vitreous acoustic velocity		: 1,532 m/s, 800 - 2,000 m/s
IOL thickness	(Acrylic)	: 0.80 mm, 0.10 - 4.00 mm
	(Silicone)	: 1.00 mm, 0.10 - 4.00 mm
	(PMMA)	: 0.80 mm, 0.10 - 4.00 mm


9) Touch the “Registration” button (10) to apply the entered settings and close the setting change screen.

10) Touch the “Cancel” button (11) to ignore the entered settings and close the setting change screen.

3.3.3 Contact/immersion mode settings

Complete the settings for contact mode and immersion mode referring to “3.7.2 d) Eye Axis Setting.”

3.3.4 Operation check

(Fig. 1)

Note ■ *This biometry test piece cannot measure or calibrate the accuracy of the instrument.*

- 1) The biometry test piece provided with the instrument allows you to check the operation of the instrument.
- 2) Select “Aphakic eye” for the eye type to be measured.
- 3) Wet the surface of the test piece and set the biometry probe perpendicular to the surface to start measurement.


3.3.5 Preparation for measurement

Note ■ *Cooperation by the patient is required to perform measurement smoothly. Explain the measurement method to the patient before starting measurement to help them relax.*

- 1) Check the measurement condition settings and operation as needed.
- 2) Apply topical anesthesia for the eye and, in chin-rest measurement mode, adjust the height of the chair, optical table, and chin rest.
- 3) For measurement in manual mode, have the patient sit down on the chair or lie on their back so that they are in a comfortable position while measurement is performed.
- 4) When the other eye needs to be guided with the fix light in chin-rest measurement mode, place the lamp on the chin rest in front of the other eye.

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[Contact mode]

Directly apply the contact section of the biometry probe perpendicularly to the center of the corneal.


[Immersion mode]

The biometry probe can be used while its tip is immersed in ultrasound media (water, corneal protective agent, etc.).

■ When the immersion attachment is not used

Note ■ *The portion to be immersed in the ultrasound media should be within 10 mm from the probe tip. Exceeding this portion may cause electric shock or failure of the probe.*


- 1) Apply the biometry probe perpendicularly to the corneal center where ultrasound media is applied.
- 2) Apply an ultrasound media such as corneal protective agent between the contact section of the biometry probe and corneal so that the gap between the corneal and contact section of the probe is always approximately 2.0 – 5.0 mm.

- When the immersion attachment is used

Note ■ *When putting the ultrasound media in the cup at the tip of the immersion attachment, pay extreme attention not to allow air to enter the cup. If air is in the cup, measurement cannot be performed or correct measurement cannot be obtained.*


- *When using the immersion attachment, be very careful not to allow the tip to contact the corneal. The corneal may be damaged or the axial length is measured shorter due to pressure.*

(Fig. 3)


Lightly apply the ultrasound media to the contact section (1) of the biometry probe and cover the probe with the immersion attachment. (Fig. 3)

Put the ultrasound media in the cup at the tip of the immersion attachment until it rises slightly. (Fig. 4)


(Fig. 4)

Set the tip of the immersion attachment perpendicular to the center of the corneal, referring to the biometry probe axis and visual axis. The appropriate position is so that the tip of the immersion attachment does not come in contact with the corneal. (Fig. 5)

(Fig. 5)


b) Handling the slider in chin-rest measurement mode

(Fig. 1)

- 1) Attach the biometry probe. When the cord length from the biometry probe to the cord hook is insufficient, the corneal may experience pressure.
- 2) Firmly press the biometry probe so that the surface of the moving section (1) of the slider is almost flush with the surface of the fixed section (2) as shown in the illustration when performing measurement. When the probe is not pressed firmly, the axial length may be measured longer due to corneal protective agent or tears. However, if the probe is pressed too far beyond this position, the corneal may experience pressure. Do not press the probe excessively.

c) Attaching the biometry probe to the applanation tonometer


(Fig. 1)

- 1) Insert the biometry probe (1) to the prism attachment section (2) of the applanation tonometer from the patient's side.
- 2) Set the measurement knob of the applanation tonometer to the following positions to prevent pressure on the corneal.
 - 10 mmHg when the intraocular pressure of the eye is 20 mmHg or lower
 - 10 mmHg lower than the intraocular pressure when the intraocular pressure of the eye is 20 mmHg or higher

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
(Fig. 1)

- 1) Check that the retina gate cursor (1) is active. If not, touch the "Retina" button (2) to switch the active gate cursor.
- 2) Touch the gate cursor movement buttons (3) to move the retina gate cursor (1) to the left of the retina waveform. (Fig. 2)

Waveform generated in the vitreous

Move the cursor to the left of the retina waveform.


(Fig. 2)

e) Selecting the waveform on the back of the lens

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
Set the gate cursor for the back of the lens to the right of the waveform on the back of the lens when there are multiple waveforms in the crystal lens and the instrument cannot recognize the original waveform on the back of the lens.


The instrument recognizes a waveform on the left of the gate cursor for the back of the lens as the waveform on the back of the lens and conducts measurement. The gate cursor on the back of the lens can be set only when the eye type is set to “Phakic eye” or “Dense cataract.”

- 1) Check that the gate cursor for the back of the lens (1) is active. If not, touch the “Lens. R” button (2) to switch the active gate cursor.
- 2) Touch the gate cursor movement buttons (3) to move the gate cursor for the back of the lens (1) to the right of the waveform on the back of the lens. (Fig. 2)

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- 1) Each of the following waveforms rises above the level cursor/line position.
 - Eye with initial stage cataracts
Waveform on the front of crystal lens / Waveform on the back of crystal lens / Retina waveform
 - Dense cataract eye
Waveform on the front of crystal lens / Retina waveform
 - Aphakic eye
Retina waveform
 - IOL eye
Waveform on the front of crystal lens / Retina waveform
- 2) The retina waveform rises vertically. (Ultrasound wave reaches the retina perpendicularly.)

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
- i. When the retina waveform spikes, the ultrasound wave reaches the retina perpendicularly.
- ii. When the waveforms on the front/back of the crystal lens spike, the ultrasound wave correctly captures the geometrical axis.
- iii. When the retina waveform can be distinguished from the sclera waveform, the ultrasound wave reaches the retina perpendicularly. Because the depression between two waveforms (choroid) cannot be confirmed when the gain is set high, etc., these waveforms may not need to be distinguished from each other.
- iv. When there is no tailing that follows the initial waveform (corneal waveform), the biometry probe directly comes in contact with the retina. When there is a layer of tears or corneal protective agent between the probe and corneal, tailing occurs after the corneal waveform. In this case, the measurement data may be unstable or tends to be longer.

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(Fig. 2)

- 4) The corneal waveform is formed within 2.0 – 5.0 mm from the initial waveform.
(Range between the dotted lines in the illustration)

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- v. There are no unnecessary waveforms between the initial waveform and corneal waveform. An unnecessary waveform appears when air bubbles are contained in the ultrasound media in the cup at the tip of the immersion attachment or the media placed between the biometry probe and corneal.

3.3.7 Measurement

- *Sterilize the biometry probe before use.*
- *Do not use the biometry probe with its contact section damaged. Measurement error may occur, and the corneal may be injured.*
- *Immersion attachment is a disposable part. Do not reuse it. Otherwise, you may contract diseases.*

- *The converted acoustic velocity directly affects the measurement data. Check that the desired value is set before starting measurement.*
- *Check the settings in contact mode and immersion mode before starting measurement. In immersion mode, apply an ultrasound media such as corneal protective agent so that the clearance between the corneal and the contact section of the biometry probe is approximately 2.0 – 5.0 mm.*
- *Be very careful not to apply pressure to the corneal when taking measurement. In chin-rest measurement mode, be very careful not to apply pressure to the corneal by moving the slider beyond the allowable range. Also, do not apply too much corneal protective agent if used. The measurement data may be affected.*
- *The automatic measurement function is an auxiliary function to take measurements more easily. This is not a function for making actual clinical judgments.*
- *Be sure to touch the “New” button to delete all the measurement data of the previous patient before measuring another patient. If new measurement is started without deleting the previous data, the measurement data of the previous patient may be included.*

a) Automatic measurement (Auto / Auto quick)

[Auto]


A waveform that satisfies the conditions is captured and its values are displayed. Up to 10 values are displayed and measurement is completed when variations are within a specific value.

[Auto quick]

A waveform that satisfies the conditions is captured and its values are displayed. Up to 10 values are displayed and measurement is completed when variations are within a specific value. (The conditions are different to those for Auto mode.)

- 1) The instrument automatically starts capturing when a waveform satisfies the conditions. When the measurement data is captured, a beep sounds.
- 2) At least 20 data sets are captured and measurement is completed when variations are within a specific value. The waveform displayed at this time is the measurement data closest to the average value.

b) Manual measurement

(Fig. 1)

- 1) The instrument beeps when the capturing conditions are satisfied.
- 2) Press the footswitch (1) to capture the measurement data.
- 3) Press the footswitch (1) again to start the next measurement.
- 4) Repeat measurements in the same manner. When 10 sets of measurement data are taken, a beep sounds again and measurement is completed.

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- Note** ■ *The deleted data cannot be restored. Carefully check the data before deleting it.*
- *Be sure to touch the “New” button to delete all the measurement data of the previous patient before measuring another patient. If the data of the next patient is captured without touching the “New” button, the patient information does not match the captured data.*

Retain the patient information and delete only the measurement data for each eye. Delete all data referring to “3.2.7 Clear all measurement data (preparation for measuring a new patient)” when the patient is changed.


(Fig. 1)

The data of only one eye is deleted.
Hold the “Retake” button (1) for approximately 1 second until a beep sounds to delete the data.
The patient information is not deleted.

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
Adjust the waveform height according to the gain settings. The gain and the waveform height increase as the gain value increases.

(Fig. 1)


- 1) Turn the gain adjustment knob (1) to set the gain. The setting is shown in the gain display field (2) on the monitor.

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(Fig. 1)

- 1) Check that the gate cursor to be adjusted is active. The active gate cursor is displayed in red and the other cursors in white. When the cursor is not active, touch the “Retake” button (1) to switch the active gate cursor. The “Retake” button appears for 5 seconds when the gate cursor is touched.
- 2) Touch the gate cursor movement buttons (2) to set the gate cursor position.

3.3.8 Checking waveforms after measurement

a) Displaying the optional waveform data


An optional waveform can be recalled from the measurement data and checked.

- 1) Touch the “Edit” button (1) when all necessary measurement data is captured or during measurement to open the measurement data check screen. (Fig. 2)

- 2) Directly touch the measurement data to select it and touch the “Wave Form” button (2). The waveform of the selected measurement data appears. (Fig. 3)

(Fig. 3)

b) Selecting specific axial length data to be used for calculating IOL power

The average value is normally adopted for the measurement data to be saved or used for calculating the IOL power, but specific data can also be selected.

(Fig. 1)

- 1) Highlight the specific data on the edit screen.
- 2) Touch the "Sel (select)" button (1) to select the data used for calculating the IOL power. An asterisk "*" (2) moves to the left of the selected data.
- 3) To use the average value for IOL power calculation, highlight the average axial length data at the bottom and touch the "Sel" button (1).

c) Deleting part of measurement data

(Fig. 1)

- 1) Highlight the specific data on the edit screen.
- 2) When the "Del" button (1) is touched to delete the data, the button indication changes to "Ret." (Fig. 2)

(Fig. 2)

- 3) When data is deleted by mistake or to cancel deletion of data, touch the "Ret" button (2) to restore the data. However, once the screen returns to the measurement screen, the deleted data cannot be restored by the "Ret" button even when the screen is switched to the edit screen again.


d) Gate change function

(Fig. 1)

(1)

- 1) Touch the “Gate” button (1) on the edit (waveform display) screen to open the gate change screen. (Fig. 2)

(Fig. 2)

(3)

(2)

(6)

(7)

- 2) Check that the gate cursor to be adjusted is active. When the cursor is not active, touch the “Retina” button (2) to switch the active gate cursor.
- 3) Touch the gate cursor movement buttons (3) to set the gate cursor position.

-
- 4) The changed measurement data is displayed in the edit data display field (4) along with the movement of the gate cursor. When the gate cursor position is changed by mistake, touch the “Initial Position” button (5) to return the gate cursor to the original position.
 - 5) Touch the “OK” button (6) to apply the edited contents and return to the previous screen (Fig. 1). Touch the “Cancel” button (7) to ignore the edited contents and return to the previous screen (Fig. 1).

e) Caliper function

Note ■ *Values measured and displayed using the caliper function are rough estimates and may differ from the actual measurement result.*

This function is used to measure the distance at an arbitrary section of the measurement waveform. Dotted caliper lines appear and the distance between these lines is displayed. Three caliper lines are shown in contact mode and four lines in immersion mode. One caliper line is displayed in red (active caliper) and the other lines in yellow.

(Fig. 1)

(1)

- 1) Touch the "Calip" button (1) on the edit screen (waveform display) to open the caliper screen. (Fig. 2)


(Fig. 2)

(6)

(2)

- 2) An enlarged waveform centered around the active caliper line is displayed. Check that the caliper line to be adjusted is active. If not, touch the “Cornea” button (2) to switch the active caliper line.
- 3) Touch the caliper line movement buttons (3) to change the caliper line position. The changed measurement data is displayed in the edit data display field (4) along with the movement of the caliper line. When the line position is changed by mistake, touch the “Initial Position” button (5) to return the caliper line to the previous position.
- 4) Touch the “ALL” button (6) to display the waveforms in normal magnification. (Fig. 3)

- 5) The distance between the start line (7) and the retina caliper line (10) is shown in the axial length field, that between the start line (7) and the caliper line for the front of the crystal lens (8) is shown in the ACD (anterior chamber depth) field, and that between the caliper lines for the front of the crystal lens (8) and for the back of the crystal lens (9) is shown in the Lens field.
- 6) Touch the “OK” button (11) to apply the edited contents and return to the previous screen (Fig. 1). Touch the “Cancel” button (12) to ignore the edited contents and return to the previous screen (Fig. 1).

3.4 IOL power calculation

3.4.1 Calculation

- Note** ■ *When using the result of axial length for calculation of the IOL power, the physician must examine the measurement result.*
- *Calculation by this instrument may cause some errors due to the number of significant digits in internal calculations.*
- *Complex numbers may be generated for the SRK/T formula. In this case, the “√ section” is calculated as zero and an asterisk “*” is shown on the right of the calculation result.*

This instrument automatically calculates the IOL and displays the result when all necessary items described in “3.4.2 Setting calculation conditions” and “3.4.3 Entering calculation parameters” are set.

3.4.2 Setting calculation conditions


Refer to “3.2.8 Selecting the eye to be measured.”

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
When measurement has already been completed, the axial length and anterior chamber depth are automatically entered and cannot be entered manually. (Refer to “3.3.8 b Selecting specific axial length data to be used for calculating IOL power.”) When measurement has not yet been performed, enter the data as follows.


(Fig. 1)

- 1) Touch the AXIAL input field (1) to activate it. The keypad (2) appears. (Fig. 2)

(Fig. 2)

- 2) Enter the data using the keypad (2).
[Acceptable range]
Axial length: 13.00 - 45.00 mm
- 3) Touch the "Ent (enter)" key to apply the data


(Fig. 3)

- 4) When HAIGIS standard is selected, the ACD field (3) also appears. Enter the appropriate value in the same manner.

[Acceptable range]

ACD : 0.00 - 10.00 mm

b) Corneal refractive power and radius of corneal curvature (K1/K2)


(Fig. 1)

- 1) Touch the K1 input field (1) or K2 input field (2) to activate it. The keypad (3) appears.

- 2) Enter the data using the keypad (3).
[Acceptable range]
Corneal refractive power: 30.00 D - 60.00 D
Radius of corneal curvature: 5.00 - 11.00 mm
- 3) Touch the “Ent (enter)” key to apply the data.

c) Desired refractive power

(Fig. 1)


- 1) Touch the Desired Ref. input field (1) to activate it. The keypad (2) appears.

Enter the data using the keypad (2).
[Acceptable range]
-30.00 D - +10.00 D
- 2) Touch the “Ent (enter)” key to apply the data.

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d) Lens constants (A-constant/SF/ACD-constant/ACD constant)

Enter lens constants for IOL according to the formula. This instrument is able to internally save up to 5 constants for 1 formula for the right and left eyes as input history. There are 2 ways to enter values. One is to enter it directly using the keypad, and the other is by selecting it from the list.

(Fig. 1)

● Entering using keypad:

- 1) Touch the A-Const field (1) to activate it. The keypad (2) appears.
- 2) Enter the data using the keypad (2).

[Acceptable range]

A-constant : 100.00 - 130.00


ACD-constant : 0.00 - +10.00

SF :-5.00 - +10.00

- 3) Touch the “Ent (enter)” key to apply the data.

● Selecting from the list:

The last 5 constants entered are saved in the instrument and these can be selected from the list.


(Fig. 1)

- 1) Touch the arrow button (1) on the right of the A-Const input field. The list appears.
- 2) Touch the desired constant to select it. The list closes and the selected value is defined.

3.4.4 Setting calculation formula


Select the calculation formula from the following 5 options.

- SRK-II
- SRK/T
- HOLLADAY
- Hoffer Q
- HAIGIS standard

(Fig. 1)

- 1) Holding the formula selection button (1) for approximately 1 second displays the formula list (2).


(Fig. 2)

- 2) Touch the desired formula to select it. The list closes and the selected formula is defined.

3.4.5 Saving measurement data

Note ■ Always enter the ID number when saving the measurement data in the internal memory.

(1)

Touch the “Save” button (1) on each screen.

The “Save” button is inactive on the IOL power calculation screen and utility screen.

[Data to be saved]

- Data of 10 measurements
- Waveforms captured by 10 measurements

3.5 Pachymetry function

3.5.1 Connecting the pachymetry probe

Connect the pachymetry probe to this instrument. Refer to “3.2.1 f) Connecting the corneal thickness/A-scan diagnosis probe.”

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- Latest: Displays the last measurement data taken.
- Minimum: Displays the minimum value of measurement data.
- Average: Displays the average of measurement data.

Refer to “3.7.2 e) Corneal Thickness Setting” for setting method.

3.5.3 Setting how to display the converted acoustic velocity and bias value

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

(Fig. 1)

- 1) Touch the “Setting” button (1) to open the setting screen (Fig. 2).

(Fig. 2)

- 2) Touch the “Setting change” button (1) to open the setting change screen (Fig. 3).

(Fig. 3)

- 3) Touch the Converted acoustic velocity Value input field (3) to activate it. The keypad (4) appears.
- 4) Enter the converted acoustic velocity using the keypad (4) and touch the “Ent” key to apply the data.
[Acceptable range]
Conversion acoustic velocity 1400 - 2000 m/s

-
- 5) Touch the “Bias Value” button (5) to select the percentage bias (%) or plus/minus bias (μm).
 - 6) Touch the bias value input field (6) to activate it. The keypad (4) appears. Using the keypad (4), enter the bias rate when the percentage bias is selected, and the correction when the plus/minus bias is selected.
[Acceptable range]
Bias rate: 60 - 130%
Correction: -600 – +450 μm
 - 7) When a value is entered for one eye, the setting for the other eye is automatically changed.
 - 8) Touch the “Initial Setting” button (7) to enter the initial settings.
 - 9) Touch the “Cancel” button (9) to ignore the entered settings and close the setting change screen.
 - 10) Touch the “Registration” button (8) to apply the entered settings and close the setting change screen.

3.5.4 Setting the measurement mode/range and data display

Select the measurement mode from the following two options.

■ Auto

Select this for measurement in Auto mode.

■ Manual

Select this when measurement is difficult in Auto mode.


Set the measurement range according to the object to be measured.

Also, set whether to display the measurement data as an actual measurement or bias value.

(Fig. 1)

- 1) Touch the "Setting" button (1) to display the setting screen (Fig. 2).

(Fig. 2)

- 2) Touch one of the "Meas Mode" buttons (2) to select the measurement mode.
- 3) Select one of the "Meas Range" buttons (3) to set the measurement range.
- 4) Select either of the "Data Display" buttons (4) to set the data display mode.
- 5) Touch the "Close" button (5) to return to the measurement screen (Fig. 1).

3.5.5 Operation check

(Fig. 1)

Note ■ *The corneal thickness test piece is for operation checks. This test piece cannot measure or calibrate the accuracy of the instrument.*

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3.5.6 Preparation for measurement

Note ■ *Cooperation by the patient is required to perform measurement smoothly. Explain the purpose of the examination and measurement method to the patient before starting measurement to help them relax.*

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3.5.7 Measurement

- *Sterilize the pachymetry probe before use.*
- *Do not use the pachymetry probe with its contact section damaged. Measurement error may occur, and the corneal may be injured.*
- *The disposable tip is a disposable part. Do not reuse it. Otherwise, you may contract diseases.*

- *The converted acoustic velocity directly affects the measurement data. Check that the desired value is set before starting measurement.*
- *The measurement data can be displayed as the actual measurement or a bias value. Check the setting of the data display mode.*
- *When the surface of the corneal is dry, wet it with saline solution.*
- *The automatic measurement function is an auxiliary function to take measurements more easily and is not a function used to actually make clinical judgment. The physician must examine the measurement result before using it.*

a) Auto measurement

- 1) Set the pachymetry probe perpendicular to the point to be measured.
- 2) The instrument automatically starts capturing measurement data when the measurement conditions are satisfied.
- 3) A beep sounds when the measurement data is captured.
- 4) When 20 sets of measurement data are captured, a beep sounds twice to indicate completion of measurement.
- 5) When the pachymetry probe is removed from the corneal, the next memory enters standby mode. Measurement is completed and the next memory allotment enters standby mode when the probe is removed from the corneal even before 20 data sets are captured.
- 6) Take measurements for all 10 memory allotments in the same manner.
- 7) When retaking measurements in a specific memory allotment during or after measurement, touch the desired memory number to enter standby mode and perform measurement again.

b) Manual measurement


- 1) Set the pachymetry probe perpendicular to the point to be measured.
- 2) The instrument beeps when the capturing conditions are satisfied.
- 3) Step on the footswitch to capture measurement data. A beep sounds when the measurement data is captured.
- 4) Repeat measurements in the same manner. When 10 sets of measurement data are taken, a beep sounds to indicate completion of measurement.
- 5) To take measurement data in the next memory allotment before 10 sets of measurement data have been taken, touch the next memory number to enter standby mode and perform measurement.
- 6) When the pachymetry probe is removed from the corneal, the next memory allotment enters standby mode.

c) Waveform display


(Fig. 1)

- 1) Measurements and waveforms are displayed while data is retrieved.

(Fig. 2)

- 2) When data capture is completed for the selected memory allotment, the memory number and measurement data (1) is displayed.
- 3) When data capture is completed or the memory allotment is in standby mode, select the memory number (1) for which measurement is completed and then touch the "Wave Form" button (2) to display the measured waveform of the selected memory number. (Fig. 3)


(Fig. 3)

- 4) Waveforms cannot be displayed when the data selection mode is set to "Average."

d) Measurement point display

This function displays the measurement point.

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(Fig. 2)

-
- 1) Touch the “Meas Point” button (1) to display the measurement point. This button is not available after measurement is started.


Touch the “Pre 1” or “Pre 2” button (2) to apply the contents preset on the utility screen to the measurement point. The contents of the selected presetting are displayed in the measurement data field. (Refer to “3.7.2 e) Corneal Thickness Setting.”)

e) Re-measurement of the same patient

Note *The deleted data cannot be restored. Carefully check the data before deleting it.*

Be sure to touch the “New” button to delete all the measurement data for the previous patient before measuring another patient. If the data of the next patient is captured without touching the “New” button, the patient information will not match the captured data.

Retain the patient information and delete only the measurement data for each eye. Delete all data referring to “3.2.7 Clear all measurement data (preparation for measuring a new patient)” when the patient is changed.

(Fig. 1)

The data for one eye is deleted.


Hold the “Retake” button (1) for approximately 1 second until a beep sounds to delete the data.

The ID, patient’s name, and gender are retained.

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
(Fig. 1)


(1)

- 1) Touch the "IOP" button (1) to open the intraocular pressure correction screen. (Fig. 2)

(Fig. 2)

- 2) Each input field for parameter 1 (2), parameter 2 (3), and IOP data (4) becomes active when touched, and the keypad (7) appears. (Fig. 3)


(Fig. 3)

- 3) Enter each parameter using the keypad.
[Acceptable range]
 - Parameter 1: 0 - 1500
 - Parameter 2: 0.0000 - 1.0000
 - Intraocular pressure data: 1.0 - 60.0 (mmHg)
1.33 – 79.99 (hPa)
- 4) The measured value is displayed in the CCT field (5) when the intraocular pressure correction screen is opened after measurement. However, if measurement points are set and one or more “CCT” is included in them, the average CCT is displayed. When measurement points are not set, the average of all memory numbers is displayed.
- 5) When all parameters are entered, the instrument automatically starts calculation and displays the corrected intraocular pressure (6).

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
Note ■ The measurement data cannot be restored once deleted. Be very careful when deleting the data.

(Fig. 1)

- 1) Touch the desired memory number button (1) to select the measurement data to be deleted.
- 2) When the "Del" button (2) is touched to delete the data, the button changes to "Ret." (Fig. 2)


(Fig. 2)

- 3) When data is deleted by mistake or to cancel deletion of data, touch the "Ret" button (3) to restore the data. When new data is captured, the deleted data cannot be restored by pressing the "Ret" button.

3.5.9 Saving the measurement data

Note ■ Always enter the ID number when saving the measurement data in the internal memory.

(Fig. 1) (1)

Touch the "Save" button (1).

[Data to be saved]

- Data of 10 measurements
- Waveforms captured by 10 measurements
- Each measurement point (only when measurements are taken with the measurement points displayed)


3.6 A-scan diagnosis

Note ■ *Both the A-scan diagnosis probe and biometry probe can be used for A-scan diagnosis, but the A-scan diagnosis probe is selected for operation when both probes are connected.*

3.6.1 Connecting the probe


Connect the A-scan diagnosis probe or biometry probe to this instrument. Refer to “3.2.1 c) Connecting the biometry probe” or “3.2.1 f) Connecting the corneal thickness probe/A-scan diagnosis probe.”

3.6.2 Switching to the A-scan diagnosis

(Fig. 1)

- 1) Touch the “To A-Diag” button (1) on the utility screen of the axial length measurement function.
- 2) The initial screen to allow you to check the Ref. data when switching modes is displayed. (Fig. 1)


- 3) Hold the “Ref” button (2) and activate the “FREEZE” button (3) to start measurement using the registered Ref. data. Touch the “FREEZE” button (3) to start measurement..
- 4) Touch the “Reset” button (4) for two seconds to delete the Ref. data. The measurement screen appears after the Ref. data is deleted
- 5) The same initial screen appears when all measurement data is cleared (preparation for measuring a new patient).
(Refer to “3.2.7 Clear all measurement data (preparation for measuring a new patient).”)

3.6.3 Gain setting

Note *Settings made here are only effective for the eye currently selected. Settings cannot be made for both eyes simultaneously. Complete necessary settings for each eye.*

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
(Fig. 1)

- 1) Turn the gain adjustment knob (1) on the exterior of the instrument to set the gain. The setting is shown in the gain field (2) on the monitor.


(Fig. 2)

3.6.4 Setting analysis method

(Fig. 1)

- 1) Touch the "Setting" button (1) to display the measurement condition setting menu.


(Fig. 2)

- 2) Select either analysis method (2). When the analysis method is selected, the analysis tool changes according to the selection.

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
(1)


(Fig. 2)

- 1) When the instrument is in FREEZE mode, touch the "Release" button (1) or step on the FREEZE pedal (2) on the footswitch to release FREEZE mode.
- 2) Apply an appropriate amount of corneal protective agent to the contact section of the A-scan diagnosis probe.
- 3) Set the A-scan diagnosis probe to the patient's eye to be examined.
- 4) Adjust the total gain using the rotary encoder.
- 5) Touch the "FREEZE" button (1) or step on the footswitch (2) to enter FREEZE mode when capturing mode-A waveforms.

3.6.6 Various functions in FREEZE mode

a) Switching data display


(Fig. 1)

This instrument stores waveforms for up to 100 data sets. Captured waveform data is displayed in order every time the “Page” button (1) is touched.


b) Setting the beam direction

Use this function to change the position to apply the probe or the direction of the ultrasound beam, or to display the position of a point of note (with pathological changes, etc.).

(Fig. 1)

- 1) Touch the “Setting” button (1) to display the measurement condition setting menu.


(Fig. 2)

(2)


(3)

(4)


(Fig. 3)

- 2) One of the nine symbols is shown each time the probe position button (2) is touched. (Fig. 3)
- 3) Numbers from 1 to 12 appear in order every time the beam direction 1 button (3) is touched. This helps indicate the angle of the meridian. (Fig. 4)
- 4) The indicator cycles through AX → P → PE → EP → E → EA → O → CB in this order every time the beam direction 2 button (4) is touched. This helps to specify a notable point in an eyeball. (Fig. 5)

(Fig. 4)


(Fig. 5)

c) Saving data

This function saves captured mode-A waveforms.

Note ■ Always enter the ID number when saving the measurement data in the internal memory.

(Fig. 1)
(1)

- 1) Display the waveform data to be saved on the FREEZE screen.
- 2) Touch the "Save" button (1) to save the image displayed on the screen.

[Data to be saved]

■ Data of 1 measurement

- Waveform captured by 1 measurement (460 bytes/waveform)
- Analysis result (Line analysis or point analysis)

3.6.7 Various functions in real time mode

a) Switching gradation characteristics

(Fig. 1)

- 1) Touch the "Setting" button (1) to display the measurement condition setting menu.

(Fig. 2)


- 2) Select a gradation characteristic (2).
The gradation characteristic can be changed only when the Ref. data is not set.

b) Setting the beam direction


Refer to "3.6.6 b) Setting the beam direction."

3.6.8 Analysis function

a) Line analysis

- 1) Set the line (1) position in real time mode using the line position button (2).
- 2) Touch the “Def.” button (3) to return the line to the original position.
- 3) Adjust the gain to align the peak of the reference section (Ref. peak) with the line height.

- 4) Touch the “FREEZE” button (4) to perform measurement and then touch the “Set Ref” button (5). The value displayed in the gain field is set to the Ref. gain.

- 5) Release FREEZE mode again and align the peak of the object section (Obj. peak) with the line height.

(Fig. 3)

- 6) After measurement is completed by touching the “FREEZE” button (4), the gain of the waveform currently displayed automatically appears beside “Set Obj.” When changing the page by touching the page button, the gain shown beside “Set Obj” also changes to the gain of the waveform displayed on the current page.
- 7) The difference between the Ref. gain and Obj. gain is displayed (6).
- 8) When FREEZE mode is released after the result is displayed, the Obj. gain is cleared.
- 9) Touch the “Ref” button (7) to open the Ref. Data display screen to set the Ref. gain again.
Refer to “3.6.2 Switching to the A-scan diagnosis” for how to delete the Ref. Data.


b) Point analysis

(Fig. 1)

(1)

- 1) When measurement is completed by touching the “FREEZE” button (1), the cursor (2) appears in the waveform monitor screen. (Fig. 2)


(Fig. 2)

(2)


(3)

- 2) Move the cursor using the cursor movement button (3) to the position to be set in P1.
- 3) When the “Set P1” button is touched, the P1 mark (4) appears at the current cursor position and the gain at the cursor position is set as P1.(5)

(Fig. 2)


- 4) Move the cursor again to the position to be set in P2.
- 5) When the "Set P2" button is touched, the P2 mark (6) appears at the current cursor position and the gain at the cursor position is set as P2.(7)
- 6) The difference between P1 and P2 appears when P1 and P2 are set.(8)

(Fig. 3)

- 7) To correct the position of P1 or P2, move the cursor and touch the "Set P1" or "Set P2" button again. The setting is overwritten and the difference between the new P1 and P2 appears.
- 8) When FREEZE mode is released while analysis is performed or after analysis is completed, the analysis result will be deleted.


3.7 Utilities

(Fig. 1)

(1)

- 1) Touch the “Utility” button (1) on each screen to open the utility menu screen (Fig. 2).


(Fig. 2)

(3)

- 2) Touch any of the setting item buttons (2) on the screen to open the corresponding setting screen. Refer to 3.7.1 and the following sections for details.
- 3) Touch the target screen button (3) to apply the setting and close the utility screen.

3.7.1 Buttons and icons on the screen

a) Utility menu screen of axial length measurement function

(1) Common Setting

Set the sound volume, screen brightness, auto power off, and date/time.

Refer to “3.7.2 a) Common Setting” for setting method.

(2) Language Setting

Set the language.

Refer to “3.7.2 b) Language Setting” for setting method.

(3) Radio Comm. Setting

Set the wireless communication.

Refer to “3.7.2 c) Wireless Communication Setting” for setting method.

(4) Eye Axis Setting

Make the settings related to measurement of the axial length.

Refer to “3.7.2 d) Eye Axis Setting” for setting method.

(5) Internal Memory Management

Set the internal memory in this instrument.

Refer to “3.7.2 f) Internal Memory Management” for setting method.

(6) Version Information

Displays the version information.

Refer to “3.7.2 g) Version Information” for setting method.

(7) IOL Parameter Setup


Registers the parameters that are used in IOL calculation.

Refer to “3.7.2 h) Registering IOL parameters” for the setting method.

(8) To A-Diag

Hold this button for two seconds to go to the A-scan diagnosis.

b) Utility menu screen of corneal thickness measurement function

(1) Common Setting

Set the sound volume, screen brightness, auto power off, and date/time.

Refer to “3.7.2 a) Common Setting” for setting method.

(2) Language Setting

Set the language.

Refer to “3.7.2 b) Language Setting” for setting method.

(3) Radio Comm. Setting

Set the wireless communication.

Refer to “3.7.2 c) Wireless Communication Setting” for setting method.

(4) Cornea Thickness

Set the method for selecting the corneal thickness data.

Refer to “3.7.2 e) Corneal Thickness Setting” for setting method.

(5) Internal Memory Management

Set the internal memory in this instrument.


Refer to “3.7.2 f) Internal Memory Management” for setting method.

(6) Version Information

Displays the version information.

Refer to “3.7.2 g) Version Information” for setting method.

c) Utility screen of A-scan diagnosis

(1) Common Setting

Set the sound volume, screen brightness, auto power off, and date/time.
Refer to “3.7.2 a) Common Setting” for setting method.

(2) Language Setting

Set the language.
Refer to “3.7.2 b) Language Setting” for setting method.

(3) Radio Comm. Setting

Set the wireless communication.
Refer to “3.7.2 c) Wireless Communication Setting” for setting method.

(4) Internal Memory Management

Set the internal memory in this instrument.
Refer to “3.7.2 f) Internal Memory Management” for setting method.

(5) Version Information

Displays the version information.
Refer to “3.7.2 g) Version Information” for setting method.

(6) To Axial


Hold this button for two seconds to go to the axial length measurement function.

3.7.2 Operations for each setting

“Common Settings” are parameters common to all modes.

a) Common Setting

[Page 1/3: Volume and Screen Brightness]

(1) Volume

Set the sound volume. Touch the “Mute” button to set the volume to “0.”

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+: Bright

(3) Page No.

The page number currently displayed (page number of “Common Setting”) is shown.


(4) “Next” button

Calls the next page in “Common Setting.”

(5) “Menu” button

Opens the utility menu screen.

[Page 2/3: Auto Power Off]

(1) Auto Power Off


Set the time until Auto Power Off mode (to turn power off) is activated.

3 min / 6 min / 10 min

When measurement data is retained when the Auto Power Off function is activated, the measurement data is automatically saved.

Refer to “3.7.2 f) Internal Memory Management” for how to check the saved data.

[Page 3/3: Date & Time Setting]


(1) Date & Time Setting

Set the date and time. Set the parameters for date and time using the keypad (2).

Y/M/D

Hour: Minute: Second

b) Language Setting


(1) Language Setting

Touching the "Language Setting" button changes the language between "Japanese" and "English."

c) Wireless Communication Setting

[Page 1/2]

Note ■ While the utility screen is opened on a device during wireless communication, operation of the other device is disabled. Disable wireless communication to open the utility screen on both devices at the same time.

Settings for wireless communication with UD-8000, IOL calculation unit, and PC KIT are described in this section.

Settings for wireless communication need to be made on each device to be connected. Refer to the instruction manual of each device for these settings.

(1) Wireless Communication Button

Enables or disables wireless communication.

Permit: Wireless communication is enabled.

Prohibit: Wireless waves for wireless communication are stopped and wireless communication is disabled.


(2) BD Address Display of Bureau

The BD address of the device is displayed.

(3) Target Place

The target place name and BD address are displayed. These fields are blank when nothing is connected.

[Page 2/2]

(1) Connection place list

Registered devices to be connected and their BD addresses are listed.

(2) "Setting" button

Opens the connection setting screen.

[Connection setting screen]

The screenshot shows a 'Connection setting screen' with a light beige background. At the top, there are two input fields: 'Name' containing 'IOL_CAL' and 'BD Address' containing '00A0962F9010'. To the right of the 'Name' field is a label 'No. 1'. Below these fields are two buttons: 'Apply' and 'Cancel'. At the bottom of the screen is a virtual keyboard with keys for numbers 1-0, letters QWERTYUIOP, ASDFGHJKL, ZXCVBNM, and special keys: 'Char/Sign', 'Space', 'A / a', 'Enter', 'Delete', 'BS', and 'Clear'. Annotations with numbers in parentheses point to specific elements: (1) points to the 'No. 1' label, (2) points to the 'BD Address' field, (3) points to the 'Space' key on the keyboard, and (4) points to the 'A / a' key on the keyboard.

(1) Connection place name

Enter the name for the connection place.

The desired display name can be set for each connection place. This name is displayed in the connection list on the wireless communication setting screen or connection selection screen.

(2) BD Address

The BD address of the connected device is displayed.

<BD address setting examples>

The following describes the settings on each device when the BD address of each device to be connected is as listed in Table 1.

Table 1: BD address examples of each device

BD address example of this instrument	00A0962F952 0	BD address example of the device to be connected Example: IOL calculation unit	00A0962F9010
---------------------------------------	------------------	---	--------------

Table 2: Setting contents on “connection setting screen” of each device

Settings of this instrument		Settings of the IOL calculation unit (connection place)	
Connection place name	IOL_CAL	Connection place name	AL4000 MEAS2
Automatically obtain BD address	00A0962F9010	BD address	00A0962F9520

- 1) When starting connection with the settings as listed in Table 2 for the partner device, the BD address of that device is automatically set in this instrument.
 - Connection place name: Set the desired name.
 - BD address: The BD address of the connected device is automatically set.
- 2) Select this instrument as a “connection place” on the wireless communication setting screen of the partner device.
- 3) Select “Permit” for “Wireless Communication” on the wireless communication setting screen on both devices.


(3) “Apply” button

Applies the entered settings and returns to the top page of the wireless communication setting screen.

(4) “Cancel” button

Cancels the entered settings and returns to the top page of the wireless communication setting screen.

d) Eye Axis Setting

(1) Select contact mode / immersion mode

Select the mode for measuring the axial length.

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(3) Calculation Method

Select the method to calculate the axial length.

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The key pad appears when the entry field is touched.

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

Latest / Minimum / Average

(2) IOP

Select the units for intraocular pressure correction.

mmHg / hPa


[Page 2/2: Measurement Selection]

(1) Measurement point

Selects measurement points to be edited and opens the measurement point edit screen.

[Measurement point edit screen]

(2) Radius

Specify the distance from the center of the measurement point. Select CCT or enter the desired value using the keypad (8).

[Acceptable range]

Radius: 0.00 - 15.00 mm

(3) Angle

Enter the angle from the horizontal axis at the measurement point.

[Acceptable range]

Angle: 0 - 179

(4) S/I selection button

Touching this button changes the setting between “S” and “I.”

(5) UP/DOWN arrow buttons

Switches memory number groups (Nos. 1 – 5 and 6 – 10) to be listed.

(6) “Confirm” button

Applies the changes and returns to the previous screen.

(7) “Cancel” button


Ignores changes and returns to the previous screen.

f) Internal Memory Management

Note ■ *The measurement data cannot be restored once deleted. Be very careful when deleting the data.*

■ *When initialization is performed, all data in the internal memory is cleared and cannot be restored. Be very careful when performing initialization.*

The patient data saved in the internal memory is displayed and edited on this screen.


(1) Internal Memory Initialization

Hold this button for two seconds to initialize the internal memory.

(2) Data Arrange

Displays a list of data in this instrument.

(1) Saved data list

The data saved in the internal memory of the measurement unit is listed. Data for up to 5 measurements can be saved.

Patient ID / Patient Name / Sex / MODE / Measurement date and time / Physician

(2) Data saved by Auto Power Off

The data automatically saved when the Auto Power Off function is activated is displayed. "*" is assigned to the data saved by the Auto Power Off function. Only 1 measurement can be saved automatically.

(3) "Delete" button

Deletes the selected data from the internal memory in the measurement unit. The data for both eyes is deleted.

(4) Mode-specific operation buttons

Different buttons appear depending on the mode of the selected data. Refer to the explanation of the mode-specific buttons described below.


(5) "List" button

Displays the data list screen.

(6) "To Utility" button

Returns to the internal memory setting screen of the utility.

[When selecting axial length measurement data]


(1) "Edit" button

Displays the selected measurement data on the edit screen.

(2) "IOL" button

Displays the selected measurement data on the IOL power calculation screen.

[When selecting corneal thickness measurement data]

(1) "Edit" button

Displays the selected measurement data to the measurement/edit screen.

(2) "IOP" button

Displays the selected measurement data on the intraocular pressure correction screen.

[When selecting A-scan diagnosis data]

Patient ID	Patient Name	Sex	MODE
00001	Tomey Taro	Male	A-DIAG
00001	Tomey Taro	Male	PACHY
00002	Tomey Hanako	Female	AXIAL
00002	Tomey Hanako	Female	PACHY
00003	Tomey Jiro	Male	AXIAL
* 00003	Tomey Jiro	Male	PACHY


Delete Meas. List Utility

(1)

(1) "Meas." Button


Displays the selected measurement data.

g) Version Information

Displays the version information of this instrument.

h) Registering IOL parameters

Parameter setting for HAIGIS standard formula

Set parameters for the HAIGIS standard formula. Holding the “Edit” button (1) for approximately 1 second displays the keypad to allow you to enter a parameter. Holding the “Init” button (2) for approximately 1 second automatically enters the default parameter.

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$$= (K1 + K2) / 2$$

L : Axial length (mm)

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REF : Desired refractive power after surgery

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SF : Distance from the iris surface to the optical
center of the implanted IOL (mm)

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Where

$$R = \frac{Rx}{1 - 0.012Rx}$$

2. Desired refractive power after surgery (D) when wearing glasses

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$$= 1.3315$$

A : A-Constant

RC: Average radius of corneal curvature (mm)

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4.2 How to calculate axial length using the axial length measurement function

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<When using the average acoustic velocity>

$$L = \frac{V \cdot t}{2}$$

L: Axial length measurement
V: Average acoustic velocity for axial length
t: Measurement time

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
4.5 Ultrasound wave output

4.5.1 MI (mechanical index)

MI is a parameter that indicates the mechanical effects of ultrasound waves on biological tissue ultrasound.

Mechanical effects are, for example, movement (or flow) around air bubbles compressed when an acoustic pressure wave of the ultrasound wave passes through biological tissue, or energy released when cavitation temporarily bursts the air bubbles.

4.5.2 TIS (thermal index of soft tissue)

TIS is a parameter that indicates the thermal effects of ultrasound waves on biological tissue ultrasound.

Thermal effects mean an increase in the tissue temperature due to emission of ultrasound waves.

5. INSPECTION AND MAINTENANCE

5.1 Warranty

One-Year Limited Warranty

The seller guarantees this product to be free from defects in material and workmanship under the normal use of this product for one (1) year or other term complying with local regulations from the date of invoice issued by Seller to the original purchaser.

Lamps, paper and other consumable items shall not be covered by this warranty.

This warranty also shall NOT apply if the product has not been installed, operated or maintained in accordance with the INSTRUCTION MANUAL of Tomey Corporation (here in after called "Tomey"). Neither seller nor Tomey shall be liable for any damages caused by purchaser's failure to follow instruction for proper installation, use and maintenance of product.

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5.2 Durable years

This instrument is designed to be remain durable for 8 years when operated under the appropriate environment and adequately inspected and serviced.

5.3 Inspection

- *When there is a problem, waveforms may not be obtained correctly. Contact Tomey or our local distributor immediately for repair or replacement.*
- *Wipe the exterior of the instrument with a soft cloth dampened with diluted neutral detergent before transferring the instrument to Tomey Corporation for maintenance and inspection. Sterilize the measurement probes as described in 5.4.1 c) before returning them.*
- *Be sure to hold the connector when removing the measurement probe from the main unit. Holding and pulling the cord may damage the inner core wires.*

Check that the contact section of the measurement probe is not damaged. Check that the connector is firmly inserted and the cable is not damaged.

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- *Do not use organic solvents such as thinner, benzene, or acetone to clean the instrument and accessories. These solvents may damage the surface of the instrument.*

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Note

■ *Be sure to hold the connector when removing the probe from the main unit. Holding and pulling the cord may damage the inner core wires.*

■ *Attach the protective cap and place the probe in the probe holder when the biometry probe, pachy probe, and A-scan diagnosis probe are not used.*

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b) Cleaning

Note

■ *Use absorbent cotton to clean the contact section of the measurement probe. Otherwise, the contact section may be damaged.*

■ *Holding and pulling the cord may damage the inner core wires.*

■ *Do not clean the connecting section of the connector.*

- When cleaning the main body of measurement probe, wipe it with absorbent cotton dampened with antiseptic soap and water or alcohol, and then with cotton dampened with water.
- When cleaning the tip of the measurement probe, wipe it with absorbent cotton dampened with alcohol, and then with cotton dampened with water.
- Wipe cords with a well wrung cloth dampened with water.
- Completely remove water after cleaning.

c) Sterilization


Note ■ *Effective chemicals vary depending on the types of germs and viruses. Select appropriate chemicals and immersion time based on the guidelines, articles, literature, etc. provided by local administrative organizations and/or societies.*

■ *Follow the described EOG sterilizing conditions. Otherwise, safe sterilization cannot be assured.*

■ *Do not use EOG sterilizer for the biometry probe and A-scan diagnosis probe.*

■ *A high-pressure steam sterilizer is acceptable only for the immersion attachment.*

- The biometry probe, pachymetry probe, and A-scan diagnosis probe can be sterilized by wiping or immersing the section within 10 mm from the contact section using the following chemicals. Thoroughly rinse the probe with purified water and thoroughly dry it after sterilization.

- Immerse the whole immersion attachment and disposable tip in alcohol or sodium hypochloride solution (0.5%), rinse them with purified water to completely remove chemicals, and thoroughly dry them.

-
- When sterilizing the pachymetry probe with EOG, remove the protective cap from the tip of the probe and put the whole probe in a sterilizing pack.

<Sterilization conditions>

EOG: Mixture of ethylene oxide gas (20%) and carbonic acid gas (80%)

Sterilizing temperature: $50 \pm 5^{\circ}\text{C}$

Sterilizing pressure: $(9.8 \pm 0.98) \times 10^4 \text{ Pa}$

Sterilizing time: 5 ± 1 hours

Aeration time:

Corneal thickness probe: 10 hours or more (continuous)

Immersion attachment: 10 hours or more (continuous)

Disposable tip: 20 cycles of 30 minutes (supply air for ventilation after each cycle)

- * Use the probe after leaving for at least 7 hours, and the disposable tip for at least 6 days, after completing sterilizing with EOG.

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Note ■ *Follow the legal requirements and regulation in the country.*

■ *Refer to the instruction manual concerning test-fluid for detail.*

a) Cleaning

Note ■ *Holding and pulling the cord may damage the inner core wires.*

■ *Do not clean the connecting section of the connector.*

- When cleaning the main body of the measurement probe and, wipe it with a suitable fluid eg. “Cleanisept® Wipes” and then wipe with cotton dampened with water
- Wipe the cord with a well wrung cloth dampened with water and remove water after cleaning.


b) Disinfection

Note ■ *Effective chemicals vary depending on the types of germs and viruses. Select appropriate chemicals and immersion time based on the guidelines, articles, literature, etc. provided by local administrative organizations and/or societies.*

■ *A high-pressure steam sterilizer is acceptable only for the immersion attachment.*

-
- The Ultrasound probe (biometry probe, pachymetry probe, and A-scan diagnostic probe) can be disinfected by immersing the contact section of the probe tip in suitable cleaning fluids.

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- PERFEKTAN®TB

(4% concentration and an exposure times of 30 min)

Antimicrobial Properties Perfektan® TB is:

- bactericidal (incl. MRSA=multi resistant Staphylococcus aureus)
- tuberculocidal
- fungicidal (C. albicans), levurocidal
- virus inactivating incl. HBV/HIV/BVDV/HCV/Papova
- Polyoma/Adeno virus - sporicidal

Thoroughly rinse the probe with purified water and thoroughly dry it after cleaning.

5.4.3 Maintenance of main unit

- Note** ■ *Do not use organic solvents such as thinner to clean the main unit of the instrument. These solvents may damage the surface of the instrument.*
- *When disconnecting cords, do not apply excessive force on the cord; for example, do not hold and pull the cord to disconnect it.*
- *Do not touch the connection terminal surface of the main unit or measurement probe.*

- Wipe off dirt on the main unit by lightly rubbing with a well wrung cloth dampened with water and then with a dry cloth. When it is very dirty, clean it by lightly rubbing with a well wrung cloth dampened with diluted neutral detergent. Then, wipe it with a cloth dampened with water and then with a dry cloth.
- Wipe off dirt on the monitor screen with a dry cloth. When it is very dirty, clean it by lightly rubbing with a well wrung cloth dampened with diluted neutral detergent. Then, wipe it with a cloth dampened with water and then with a dry cloth.

- Note** ■ *Do not leave any water, chemicals, etc. remaining on the touch panel as the surface of the touch panel is vulnerable to moisture. The surface of the touch panel may be damaged and unable to maintain the appropriate performance.*

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5.5 Replacing consumables

- Contact Tomey Corporation or our local distributor for replacement of the battery.

5.6 Storing

- *Install the instrument in a location free of water or chemicals. Any water or chemicals entering the instrument may cause an electric shock or failure.*

- *Remove the battery from the main unit when not used for long periods.*

- *Do not store the instrument in a location where chemicals are stored or gases may occur. Spilt chemicals or vapor may enter the instrument and result in fire.*

- *Disconnect the power cord from the outlet to ensure safety when the instrument is not operation for 1 month or longer.*

- *Store the instrument in a location not subject to direct sunlight, high temperature and humidity, or air containing dust, salt and/or sulfur. Otherwise, failure or malfunction may occur.*

- *Store the instrument in a leveled stable place free of vibration or mechanical impact. Otherwise, measurement cannot be conducted correctly. The instrument may also topple over or fall down, resulting in fire or a fatal accident.*

- Clean and dry the measurement probe, attach the protective cap, and place it in the probe case for storage.

-
- When the instrument has not been used for a while, check that the instrument is operating correctly before starting operation again.
 - The service life of battery may be affected when the instrument is stored in a hot place. Store the instrument under the following conditions.

Storage conditions

- Ambient temperature range: -20 – +40°C
- Relative humidity range: 45 – 85%

5.7 Disposal

- Note** ■ *Keep the box and packing materials for use when moving or transporting the instrument.*
- *Keep the packing materials and box together.*
- *When disposing of the packing materials, sort them by type and abide by local government rules and regulations.*
- *A lithium battery is used in the instrument. Handling of the lithium battery varies depending on governing bodies. Follow relevant laws and local rules and regulations, or contact our local distributor or Tomey Corporation.*

6. TROUBLESHOOTING

Check the following first when you encounter any problems. If the problem is not solved even after checking the applicable item listed below, contact our local distributor to request inspection and/or repair.

■ ***Do not remove the cover of the instrument. You may be directly exposed to high voltage sections.***

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<Common items>

The power LED does not turn on when the power switch is turned on.

? Cause 1 Problem with the power plug

Solution ➡ Check that the power plug is firmly connected to the outlet. Check that the AC adapter is firmly connected to the main unit. Check that there is no flaws in the AC cord, such as cracks or tears.

? Cause 2 Problem with the power outlet

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Nothing appears on the monitor screen.

? Cause 1 The auto power off function, which automatically turns off the screen when the instrument is not operated for the specified time, has been activated.

Solution ➡ Touch the monitor screen.

? Cause 2 Built-in connection is established with UD-8000.

Solution ➡ Nothing appears on the monitor while connected to UD-8000. Disconnect AL-4000 from UD-8000.

A different button to the one touched on the touch panel becomes active.

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Solution ➡ Turn power on while pressing the MODE button on the front of the main unit to display the touch panel calibration screen. [+] appears on the screen. Touch the center. The point touched disappears and then the next [+] appears. Touch the center again. Repeat this step until [+] does not appear to complete calibration. Press the MODE button to return to the measurement screen. If this calibration of the touch panel problem occurs repeatedly during operation, contact our local distributor.

An error message “Connect the USB cable again” appears and measurement cannot be performed.

? Cause 1 The USB cable is disconnected.

Solution ➡ Connect the USB cable again. Communication is resumed and measurement can be started.

“Data remains on both devices and connection could not be established. Set either device to new patient mode and connect the devices again.” This error message appears and communication cannot be started.

? Cause 1 The patient information and measurement data is retained on both devices.

Solution ➡ Touch the “New” button on both or either device to delete the patient information and measurement data.

Wireless communication does not start even when power is turned on.

? Cause 1 The device to be connected is not turned on.

Solution ➡ Turn on the device to be connected.

? Cause 2 Settings for wireless communication are not set for this instrument.

Solution ➡ Set the wireless communication settings on the instrument, referring to “3.7.2 c) Wireless Communication Setting.”

? Cause 3 Settings for wireless communication are not set for the device to be connected.

Solution ➡ Set the wireless communication settings for the device to be connected.

The patient information cannot be entered or changed.

? Cause 1 The instrument is connected to another device.

Solution ➡ The patient information cannot be entered from the AL-4000 measurement unit while connected to another device. Cut the connection with the other device.

? Cause 2 The measurement result has already been saved.

Solution ➡ The patient information cannot be changed after the data is saved.

<Axial length measurement function>

The axial length cannot be measured.

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Solution ➡ Connect the probe correctly until it is locked as described in “3.2.1 b) Connecting the biometry probe.”

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Solution ➡ In immersion mode, data is captured when a corneal waveform is within 2.0 – 5.0 mm from the initial waveform. Place the corneal protective agent or immersion attachment between the corneal and the tip of the biometry probe to satisfy these conditions.

? Cause 3 The biometry probe is directly applied to the corneal without using a corneal protective agent or immersion attachment. (In the case of immersion mode)

Solution ➡ Use the immersion attachment or corneal protective agent for measurement in immersion mode as described in “3.3.5 a) Handling the biometry probe in contact/immersion mode.” Perform measurement in contact mode when directly applying the biometry probe to the corneal.

Measurement cannot be performed in Auto mode.

? Cause 1 The mode is set to Manual.

Solution ➡ Set to auto measurement (Auto or Auto quick). (Refer to 3.3.2 Setting the axial length calculation method and converted acoustic velocity for the eye to be measured.”)

? Cause 2 Noise generated in the surrounding area.

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? Cause 4 The contact section of the biometry probe is damaged.

Solution ➡ If damaged, immediately stop measurement and contact our local distributor.

Measurement cannot be performed in Manual mode.

? Cause 1 The footswitch is not connected correctly.

Solution ➡ Connect the footswitch correctly until it is locked as described in “3.2.1 c) Connecting the footswitch.”

? Cause 2 The mode is set to auto measurement (Auto or Auto quick).

Solution ➡ Set to Manual mode. (Refer to “3.3.2 Setting the axial length calculation method and converted acoustic velocity for the eye to be measured.”)

Measurements are unstable or inappropriate.

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Solution ➡ Connect the probe correctly until it is locked as described in “3.2.1 b) Connecting the biometry probe.”

? Cause 2 There is too much pressure on the corneal.

Solution ➡ When using the chin rest, ensure adequate length from the biometry probe to the cord hook so that the probe cord does not pull the slider of the chin rest. (Refer to “3.3.6 b) Handling the slider in chin-rest measurement mode.”)

? Cause 3 The retina gate cursor is not set in an appropriate position.

Solution ➡ Set the retina gate cursor on the immediate left of the actual retina waveform. Make sure that there are no

unnecessary waveforms between the retina gate cursor and retina waveform.

? Cause4 The biometry probe is not applied appropriately.

Solution ➡ Apply the probe to the center of the corneal perpendicularly so that the retina waveform spikes.

? Cause5 Too much corneal protective agent is applied. (In the case of contact mode)

Solution ➡ When too much corneal protective agent is applied, a gap is formed between the biometry probe and corneal, and the measurement result may tend to be longer. Apply appropriate amount of corneal protective agent.

? Cause6 The immersion attachment is used. Or, a large amount of corneal protective agent is applied. (In the case of contact mode)

Solution ➡ The measurement result tends to be longer only by the thickness of corneal protective agent. Complete settings in immersion mode and contact mode correctly, and take measurements as described in “3.3.5 a) Handling the biometry probe in contact/immersion mode.”

? Cause7 The biometry probe is directly applied to the corneal without using the corneal protective agent. (In the case of immersion mode)

Solution ➡ Use the immersion attachment or corneal protective agent for

measurement in immersion mode as described in “3.3.5 a) Handling the biometry probe in contact/immersion mode.” Perform measurement in contact mode when directly applying the biometry probe to the corneal.

? Cause 8 The converted acoustic velocity is not set appropriately.

Solution ➡ Check the setting of the converted acoustic velocity. A different converted acoustic velocity can be set for right and left eyes. Check the setting for both eyes. (Refer to “3.3.2 Setting the axial length calculation method and converted acoustic velocity for the eye to be measured.”)

? Cause 9 The contact section of the biometry probe is damaged.

Solution ➡ If damaged, immediately stop measurement and contact our local distributor.

Monitoring sound does not stop. Data is captured even while the probe does not contact the eyeball.

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Solution ➡ Remove drops of water.

Monitoring sound does not go off.

? Cause 1 The volume is set to level 0 or “Mute.”

Solution ➡ Set the volume to a level other than “0” or “Mute” as described in “3.7.2 a)

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Solution ➡ Connect the probe correctly until it is locked as described in “3.2.1 b) Connecting the biometry probe.”

? Cause 2 Noise is generated in the surrounding area.

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The probe error confirmation screen appears when the power is turned on.

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Solution ➡ Connect the probe correctly until it is locked as described in “3.2.1 b) Connecting the biometry probe.”

The fixlight in the biometry probe is turned off.

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Solution ➡ Set to ON as described in “3.7.2 d) Eye Axis Setting.”

<Pachymetry function>

The probe error confirmation screen appears when the power is turned on or the pachymetry probe is calibrated.

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Solution ➡ Connect the probe correctly until it is locked as described in “3.2.1 d) Connecting the corneal thickness/A-scan diagnosis probe.”
Touch the “Calib” button on the measurement screen or restart the instrument to calibrate the pachymetry probe after connecting the probe.

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Solution ➡ Connect the probe correctly until it is locked as described in “3.2.1 d) Connecting the corneal thickness/A-scan diagnosis probe.”

The error message “Out of range!” appears and measurement cannot be performed.

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Stable measurements cannot be obtained. Appropriate measurements cannot be obtained.

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Solution ➡ Connect the probe correctly until it is locked as described in “3.2.1 d) Connecting the corneal thickness/A-scan diagnosis probe.”

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Solution ➡ Apply the pachymetry probe perpendicular to the corneal at the position to be measured.

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Solution ➡ Check that the converted acoustic velocity is set correctly. (Refer to “3.5.3 Setting how to display the conversion sound velocity and bias value.”)

? Cause 4 The contact section of the pachymetry probe is damaged.

Solution ➡ If damaged, immediately stop measurement and contact our local distributor.

? Cause 5 The pachymetry probe is not calibrated correctly.

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Measurement cannot be performed in Auto mode.

? Cause 1 The pachymetry probe is not calibrated correctly.

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? Cause 2 Noise is generated in the surrounding area.

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? Cause 3 The contact section of the pachymetry probe is damaged.

Solution ➡ If damaged, immediately stop measurement and contact our local distributor.

? Cause 4 The measurement mode is set to Manual.

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The measureable point does not move immediately when the pachymetry probe is removed from the corneal.

? Cause 1 The tip of the pachymetry probe is wet.

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Monitoring sound does not stop. Data is captured even while the probe does not contact with the eyeball.

? Cause 1 The tip of the pachymetry probe is wet.

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<A-scan diagnosis>

The mode-A waveform is not displayed.

? Cause 1 The A-scan diagnosis probe is not connected correctly.

Solution ➡ Check that the A-scan diagnosis probe is connected until its connector is locked.

? Cause 2 There is no corneal protective agent applied or the amount is too small.

Solution ➡ Apply the appropriate amount of corneal protective agent to the tip of the A-scan diagnosis probe. If the problem is not solved even after the protective agent is applied, apply more.

The mode-A waveform is disarranged.

? Cause 1 Noise is generated in the surrounding area

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The function to display the saved data is not available.

? Cause 1 Data is not saved in the internal memory.

Solution ➡ After the waveform data is saved in the internal memory by the function to save data in the internal memory.

The probe error confirmation screen appears when switching to the A-scan diagnosis.

? Cause 1 The A-scan diagnosis probe is not connected correctly.

Solution ➡ Check that the A-scan diagnosis probe is firmly connected.

<Related to data saving>

Data cannot be saved in the internal memory.

? Cause 1 The ID is not entered.

Solution ➡ The instrument uses the ID as a file name. Enter the ID and then save the data.

? Cause 2 There is not enough space available in the internal memory. (The confirmation screen appears when attempting to save the data.)

Solution ➡ Delete unnecessary data in the internal memory and save the data again.

7. SPARE PARTS AND OPTIONAL PARTS

The following spare parts and accessories are available from our local distributor of this instrument.

Contact our local distributor to order them.

7.1 Spare parts

- Immersion attachment

Specify the part type as “AL-4000 immersion attachment.”

7.2 Optional parts

- Pachymetry probe (with case and test piece)

Specify the part type as “AL-4000 pachy probe.”

- A-scan diagnosis probe

Specify the part type as “AL-4000 A-scan diagnosis probe.”

- Chin rest AL-1100

Specify the part type as “AL-4000 chin rest.”

This page is intentionally blank

8. SPECIFICATIONS

8.1 Specifications

8.1.1 Biometry/IOL power calculations

- Measurement range
 - Axial length : 13.00 to 45.00mm
 - ACD : 1.50 to 7.00mm
 - Lens thickness : 2.00 to 6.00mm

- Accuracy
 - Measurement Accuracy : ± 0.1 mm
 - Measurement Resolution : 0.01mm

- Factory setting of the converted velocity
 - Average velocity of the axial length
 - Phakic Eye : 1,550m/s
 - Dense Cataract : 1,548m/s
 - Aphakic : 1,532m/s
 - Lens velocity
 - Phakic Eye : 1,641m/s
 - Dense Cataract : 1,629m/s
 - ACD velocity
 - Phakic Eye : 1,532m/s
 - Dense Cataract : 1,532m/s
 - IOL Eye : 1,532m/s
 - Velocity of vitreous
 - Phakic eye : 1,532m/s
 - (In case of sectional velocity)
 - IOL Eye : 1,532m/s
 - IOL velocity
 - Acrylic : 2,200m/s
 - Silicone : 1,049m/s
 - PMMA : 2,718m/s
 - User : 2,200m/s

- IOL thickness

Acrylic	: 0.80mm
Silicone	: 1.00mm
PMMA	: 0.80mm
User	: 0.80mm

- IOL power calculation

- SRK II
- SRK/T
- HOLLADAY
- Hoffer Q
- HAIGIS Standard

- Biometry probe

■ Type	: Solid state
■ Fixlight	: Built in the probe, Red LED
■ Transducer frequency	: 10MHz
■ Tip diameter	: 6.0mmφ(concave)
■ Dimension/Weight	: 14.0mmφ×98mm/40g

8.1.2 Function of pachymetry

- Measurement range : 150 to 1,500μm

- Accuracy

■ Measurement accuracy	: ±5μm
■ Measurement resolution	: 1μm

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■ Transducer frequency	: 20MHz
■ Tip diameter	: 1.5mmφ with an angle of 45°(flat)
■ Dimension/Weight	: 8.8mmφ×90mm/30g

8.1.3 A-scan Diagnosis

● Display Range	
■ Line analysis	: 64mm
■ Point analysis	: 45mm
● A-Scan diagnostic probe	
■ type	: Solid state
■ Transducer frequency	: 10MHz
■ Tip diameter	: 5mmφ(flat)
■ Dimension/Weight	: 8mmφ×97mm/30g

8.1.4 Main unit

● TFT LCD	: 3.5 inch color
● Dimension/Weight	: 109(W)×52(D)×166(H) 470g

With AC adapter

● Voltage	: 100-120/220-240 VAC
● Frequency	: 50/60Hz
● Power consumption	: 18/24 VA

Using built-in battery

● Voltage	: DC 3.7V
-----------	-----------

8.1.5 AC adapter

● Manufacturer	: SINPRO ELECTRONICS CO., LTD
● MODEL No.	: MPU16B-104-P25E001

8.1.6 Battery

● Battery	: Lithium ion battery
● Voltage	: DC 3.7V
● Manufacturer	: SANYO ELECTRIC CO., LTD.
● MODEL	: 1UF603443S

8.2 Energy and other consumption

8.2.1 Influences of ultrasound energy on the human body

***This unit is exclusively designed for ophthalmic use,
so do not use the instrument for any purposes
other than ophthalmic use.***

This unit is exclusively designed for ophthalmic biometry and pachymetry purposes.

The unit is set with a very weak ultrasound energy appropriate for ophthalmic use; however, minimize the duration of measurement and give minute caution not to cause the Patient to be exposed and with ultrasound energy while measurement.

8.2.2 Ultrasound energy

Acoustic Output Reporting Table(FDA)

Application(s):Biometry

Acoustic Output			MI	Ispta.3 (mW/cm2)	Isppa.3 (W/cm2)
Pre-Amendment Maximum Acoustic Output			0.23	17	28
Global Maximum Value			0.173	1.599	8.738
Associated Acoustic Parameter	pr.3(Mpa)		0.548		
	W0(mW)			0.0133	0.0133
	fc(MHz)		10.09	10.09	10.09
	Zsp(cm)		1.83	1.83	1.83
	Beam dimensions	x-6(cm)		0.161	0.161
		y-6(cm)		0.159	0.159
	PD(usec)		0.153		0.153
	PRF(Hz)		1200		1200
	EB	Az.(cm)		0.54	
	D	Ele.(cm)		0.54	

Application(s):Pachymetry

Acoustic Output			MI	Ispta.3 (mW/cm2)	Isppa.3 (W/cm2)
Pre-Amendment Maximum Acoustic Output			0.23	17	28
Global Maximum Value			0.148	11.185	12.710
Associated Acoustic Parameter	pr.3(Mpa)		0.666		
	W0(mW)			0.0208	0.0208
	fc(MHz)		20.14	20.14	20.14
	Zsp(cm)		0.05	0.05	0.05
	Beam dimensions	x-6(cm)		0.057	0.057
		y-6(cm)		0.055	0.055
	PD(usec)		0.110		0.110
	PRF(Hz)		8000		8000
	EB	Az.(cm)		0.19	
	D	Ele.(cm)		0.18	

Application(s):A-Scan

Acoustic Output			MI	Ispta.3 (mW/cm2)	Isppa.3 (W/cm2)
Pre-Amendment Maximum Acoustic Output			0.23	17	28
Global Maximum Value			0.0750	0.0127	1.747
Associated Acoustic Parameter	pr.3(Mpa)		0.242		
	W0(mW)			0.000517	0.000517
	fc(MHz)		10.44	10.44	10.44
	Zsp(cm)		2.149	2.149	2.149
	Beam dimensions	x-6(cm)		0.30	0.30
		y-6(cm)		0.24	0.24
	PD(usec)		0.2425		0.2425
	PRF(Hz)		30		30
	EB	Az.(cm)		0.4	
	D	Ele.(cm)		0.4	

MI:	the Mechanical index
Ispta.3:	the derated spatial-peak temporal-average intensity
Isppa.3:	the derated spatial-peak pulse-average intensity
Pr.3	the derated peak rarefactional pressure
W0:	the ultrasonic power
fc:	the center frequency
Zsp:	the axial distance at which the reported parameter is measured
x-6,y-6:	are respectively the in-plane(azimuthal) and out of plane(elevational) -6dB dimensions in the x-y plane where Zsp is found
PD:	the pulse duration
PRF:	the pulse repetition frequency
EBD:	the entrance beam dimensions for the azimuthal and elevational planes

Acoustic Output Report (IEC60601-2-37)

Index label		MI			TIS Non-Scan($A_{aprt} \leq 1 \text{ cm}^2$)		
		Bio	Pachy	A-Diag	Bio	Pachy	A-Diag
Maximum index value		0.173	0.148	0.0750	0.00228	0.00188	0.000121
Associated acoustic parameters	pr, α [Mpa]	0.548	0.666	0.242			
	P[mW]				0.0474	0.0223	0.00243
	Min.of [Pa(zs),Ita, α (z s)]						
	Zs						
	Zbp						
	Zb						
	Z at max.lpi, α [cm]	1.83	0.05	2.149			
	deq(zb)						
	fawf[MHz]	10.09	20.14	10.44	10.09	20.14	10.44
	Dim of Aaprt [cm]						
Other information	X				0.54	0.19	0.4
	Y				0.54	0.18	0.4
	Td[us]	0.153	0.110	0.2425			
	Prr[Hz]	1200	8000	30			
	Pr at max.lpi[MPa]	1.037	0.689	0.525			
Operating Control conditions	deq at max.lpi						
	lpa, α at max.MI [W/cm ²]	8.738	12.710	1.747			
Control							

α :	ACOUSTIC ATTENUATION COEFFICIENT
Aaprt:	-12dB OUTPUT BEAM AREA
Deq:	EQUIVALENT APERTURE DIAMETER
deq:	EQUIVALENT BEAM DIAMETER
fawf:	ACOUSTIC WORKING FREQUENCY
Ipa:	PULSE-AVERAGE INTENSITY
Ipa, α :	ATTENUATED PULSE-AVERAGE INTENSITY
Ipi:	PULSE-INTENSITY INTEGRAL
Ipi, α :	ATTENUATED PULSE-INTENSITY INTEGRAL
MI:	MECHANICAL INDEX
P:	OUTPUT POWER
P α :	ATTENUATED OUTPUT POWER
Pr:	PEAK-RAREFACTIONAL ACOUSTIC PRESSURE
Prr:	PULSE REPETITION RATE
TIS:	SOFT-TISSUE THERMAL INDEX
td:	PULSE DURATION
X,Y:	-12dB OUTPUT BEAM DIMENSIONS
z:	DISTANCE FROM THE SOURCE TO A SPECIFIED POINT
zb:	DEPTH FOR TIB
zbp:	BREAK-POINT DEPTH
zs:	DEPTH FOR TIS

8.3 Noise

The instrument makes monitoring sounds on these occasions.

- When turning the power on
- When self-diagnosing
- When using various keys
- When taking measurement(taking measurement data)
- When errors occur

8.4 Environment conditions

Use this unit under the environment conditions as followings

- Placeindoor without sunshine
- Temperature+10°C to +40°C
- Humidity30 to 75%
- Pressure700 to 1060hPa
- Power variationwithin ±10%

Store and transport this unit in the original packing box under the environment conditions as followings.

- Temperature-20°C to +50°C
- Humidity45 to 85%

8.5 Classification of ME Equipment

Protection against electrical shock:

- Class II ME equipment (When using the AC adaptor)
- Internally powered ME equipment (When using the built-in battery)

Applied parts

- B applied parts (the tip of the biometry probe, the tip of the pachymetry probe and the tip of the A-scan diagnosis probe)

IP Code:

- IP20 (Main unit)
- IP21 (Foot switch)
- IP47 (the tip of biometry probe, the tip of the pachymetry probe and the tip of the A-scan diagnosis probe)
- IP41 (parts of the biometry probe, the pachymetry probe and the A-scan diagnosis probe except each tip and each connector.)

Mode of Operation:

- Continuous operation

8.6 Declaration of Conformity to EMC

- *When using specified accessories or cables excluding internal parts sold as repair parts of TOMEY Corporation, electromagnetic emission of the AL-4000 may increase or electromagnetic immunity may decrease.*

- *When using handheld or portable communication equipment, locate it more than 30cm(12inches) away from the AL-4000. It may degrade the performance of the AL-4000 if the equipment is located closer.*

- *Do not use the AL-4000 adjacent to or stacked on other equipment. If it needs to be adjacent to or placed on other equipment, it must be confirmed that AL-4000 operates correctly in such location.*

- *The AL-4000 falls under Group 1, Class B according to EN 55011(CISPR 11). This means that the AL-4000 does not generate RF energy intentionally in the form of electromagnetic irradiation, inductive coupling and/or capacitive coupling for processing of materials or inspection/analysis, and that it is suitable for use in facilities directly connected to a low-voltage power network that supplies power to general household facilities and buildings used for home use.*

- *The AL-4000 can be used in a specialized health care facility such as hospitals and clinics (except anywhere near the source of a strong EMI) and a general commercial facility such as eye glasses stores.*

- *The AL-4000 requires special attention for EMC and needs to be installed, serviced and used based on the following information.*

- *Do not use cables other than those provided or specified by TOMEY Corporation.*

- *The person who connects additional equipment to the signal I/O section as part of the medical system, shall take responsibility for ensuring the system conforms to IEC/EN 60601-1-2 requirements.*

Cable List

Name	Specification	Length(m)	Note
AC power code	Unshielded	1.8m	Accessories
DC power code (AC adapter)	Unshielded	1.3m	Accessories
Axial Probe Cable	Shielded	1.6m	Accessories
Pachy Probe Cable	Shielded	1.8m	Accessories
USB cable	Shielded	1.8m	Accessories
Foot Switch Cable	Shielded	1.4m	Accessories

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Table 3 【Input a.c. power PORT】

Test	Standard	IMMUNITY TEST LEVELS
		HOME HEALTHCARE ENVIRONMENT
Electrical fast transients / bursts	IEC 61000-4-4	$\pm 2\text{ kV}$ 100kHz repetition frequency
Surges Line-to-line	IEC 61000-4-5	$\pm 0.5\text{ kV}$, $\pm 1\text{ kV}$
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V 0.15MHz - 80MHz 6V in ISM and amateur radio bands between 0.15MHz and 80MHz 80% AM at 1kHz
Voltage dips	IEC 61000-4-11	0% U_T , 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
		0% U_T , 1 cycle and 70% U_T , 25 / 30 cycles Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	0% U_T , 250 / 300 cycle

Table 4 【Signal input / output parts PORT】

Test	Standard	IMMUNITY TEST LEVELS
		HOME HEALTHCARE ENVIRONMENT
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0.15 MHz - 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz

Table 5 【PATIENT coupling PORT】

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A

AC adapter	8-4
Accuracy	8-1, 8-2
ACD	1-6, 8-1
ACD velocity	8-1
A-constant	3-57
Acrylic	1-6, 4-10, 8-1, 8-2
A-Diag	1-6
Amplifier characteristics of the A-scan diagnosis	4-11
Analysis function	3-86
Angle	3-104
Aphakic	1-6, 8-1
Aphakic eye	3-24, 3-28, 3-37, 4-9
Appearance check	5-3
A-scan diagnosis	3-77
A-scan Diagnosis	8-3
A-Scan diagnostic probe	8-3
Attaching the biometry probe to the applanation tonometer	3-33

Attaching the probe sleeve	3-7
Auto	3-26
Auto measurement	3-67
Auto Power Off	3-96
Auto power off function	1-7, 2-20
Auto quick	3-26
Automatic connection upon startup	3-13
Automatic measurement (Auto / Auto quick)	3-41
Average	3-61
Average velocity of the axial length	8-1
AVG	1-6
Axial	1-6
Axial length	3-53, 8-1
Axial Length Measurement (Biometry) ...	2-9
Axial length measurement function	3-24

B

Basic structure and common items	2-7
Battery	8-4
BD address	1-6
BD Address	3-99
BD Address of Bureau	3-98
BD address setting examples	3-100
Bias value	1-7
Biometry probe	8-2
Biometry/IOL power calculations	8-1
Buttons and icons on the screen	3-91

C

Calculation	3-53
Calculation Method	3-101
Calibration	1-6
Caliper	1-7
Caliper function	3-51
CCT	1-6

Checking the measurement data	3-73
Checking the software version	4-11
Checking waveforms after measurement..	3-45
Cleaning	5-3, 5-6
Clear all measurement data (preparation for measuring a new patient)	3-22
Common Setting	3-91, 3-93, 3-94, 3-95
Communicating via USB	2-8
Communicating via wireless access.....	2-8
Communication check/wireless communication selection button	2-8
Connecting the AC adapter	3-4
Connecting the A-scan diagnosis probe.....	3-77
Connecting the biometry probe	3-6, 3-24
Connecting the corneal thickness/A-scan diagnosis probe	3-9
Connecting the footswitch	3-8
Connecting the pachymetry probe	3-61
Connecting UD-8000, IOL calculation unit, and personal computer	3-10
Connection for wireless communication.....	3-13
Connection place list	3-98
Connection place name.....	3-99
Connection setting screen.....	3-99
Connections.....	3-4
Contact	1-6
Contact mode	3-30, 3-37
Contact/immersion mode settings	3-29
Corneal Thickness.....	3-93
Corneal Thickness Setting	3-102
Corneal refractive power and radius of corneal curvature (K1/K2).....	3-55
Corneal thickness measurement.....	2-15
Cut (communication available).....	2-8

D

D	1-6
---------	-----

Data Arrange	3-105
Data Selection Method.....	3-102
Date & Time Setting	3-96
Deleting part of measurement data ..	3-48, 3-75
Dense Cataract	8-1
Dense cataract eye	3-24, 3-28, 3-37, 4-9
Desired refractive power	3-56
Diagnosis in mode A	2-17
Displaying the optional waveform data ..	3-45

E

Edit screen (data list)	2-11
Edit screen (waveform).....	2-13
Entering calculation parameters	3-53
Entering the patient data	3-19
Eye Axis Setting	3-91, 3-101
Eye with initial stage cataracts.....	3-37

F

Factory setting and setting range of the converted velocity	8-2
Function of pachymetry.....	8-2

G

Gain.....	1-7
Gain setting	3-43, 3-79
Gate	1-7
Gate change function.....	3-49

H

Handling the biometry probe in contact/immersion modes	3-30
Handling the slider in chin-rest measurement mode	3-32
How to calculate axial length using the axial length measurement function	4-9

I

Immersion	1-6
Immersion mode	3-30, 3-39
Initial Gain	3-101
Inspection	5-2
INSPECTION AND MAINTENANCE	5-1
Internal Memory Initialization	3-105
Internal Memory Management	3-91, 3-93, 3-94, 3-105
Intraocular pressure correction	3-73
Intraocular pressure data	3-74
IOL	1-6
IOL calculation screen	2-14
IOL Ey	8-1
IOL eye	3-24, 3-28, 3-37, 4-10
IOL Eye	8-1
IOL power calculation	3-53, 8-2
IOL power calculation formula	4-1
IOL thickness	8-2
IOL velocity	8-1
IOP	3-102

K

K1	1-6
K2	1-6

L

Language Setting	3-91, 3-93, 3-94, 3-97
Latest	3-61
Lens thickness	8-1
Lens velocity	8-1

M

Main unit	8-3
Maintenance of main unit	5-8
Maintenance of measurement probes ..	5-3, 5-6
Manual	3-26

Manual measurement	3-41, 3-67
Measurement	3-40, 3-66, 3-81
Measurement accuracy	8-2
Measurement Accuracy	8-1
Measurement point	1-7, 3-103
Measurement point display	3-70
Measurement range	8-1
Measurement resolution	8-2
Measurement Resolution	8-1
Measurement screen	2-9, 2-15
Measurement screen (FREEZE)	2-19
Measurement screen (real time)	2-17
Measurement Selection	3-103
Measuring range	8-2
MI (mechanical index)	4-12
Minimum	3-61

N

NAMES AND FUNCTIONS	2-1
---------------------------	-----

O

Obj	1-6
Operation check	3-29, 3-65
OPERATION PROCEDURES	3-1
Operations for each setting	3-95

P

Pachy	1-6
Pachymetry function	3-61
Pachymetry probe	8-2
Page No.	3-95
Parameter 1	3-74
Parameter 2	3-74
Percent bias	8-2
Percent bias/plus minus bias setting range...	8-2
Percentage bias	3-61
Phakic	1-6

Phakic eye	3-24, 3-28, 4-9, 8-1
Phakic Eye.....	8-1
Plus minus bias	8-2
Plus/minus bias	3-61
PMMA.....	1-6, 4-10, 8-1, 8-2
Point analysis	3-88
Post	1-6
Precautions for measurement	3-30
Preparation before use.....	3-4
Preparation for measurement.....	3-29, 3-65
Probe Light	3-101

R

Radio Comm. Setting	3-91, 3-93, 3-94
Radius	3-103
Ref.	1-6
Re-measurement of the same patient	3-42, 3-72
Replacing consumables	5-9
Replacing the probe holder	3-5
Routine maintenance	5-2

S

S/I selection button.....	3-104
Saved data list ...3-106, 3-107, 3-108, 3-109	
Saving data.....	3-84
Saving measurement data	3-59
Saving the measurement data	3-76
Screen	2-7
Screen Brightness	3-95
Select contact mode / immersion mode ...	3-101
Selecting specific axial length data to be used for calculating IOL power	3-47
Selecting the eye to be measured.....	3-23
Selecting the retina waveform	3-34
Selecting the waveform on the back of the lens	3-35

Setting analysis method.....	3-80
Setting axial length calculation method and converted acoustic velocity for eye to be measured	3-24
Setting calculation conditions	3-53
Setting how to display the converted acoustic velocity and bias value.....	3-61
Setting measurement conditions	3-19
Setting the beam direction	3-82
Setting the data type to be displayed....	3-61
Setting the gate cursor.....	3-44
Setting the measurement mode/range and data display	3-63
Silicone	1-6, 4-10, 8-1, 8-2
SPARE PARTS AND OPTIONAL PARTS	7-1
SPECIFICATIONS	8-1, 9-1
SRK/T formula	4-1
Sterilization.....	5-4, 5-6
Sterilization conditions	5-5
Storing	5-9
Switching modes.....	3-18
Switching to the A-scan diagnosis	3-77

T

Target Place	3-98
TIS (thermal index of soft tissue)	4-12
Touch panel.....	1-7
TROUBLESHOOTING.....	6-1
Turning the power on and adjustment after turning the power on	3-16

U

Ultrasound wave output	4-12
UP/DOWN arrow buttons.....	3-104
US	1-6
User.....	8-1, 8-2
Utilities	3-90

Utility menu screen of axial length	
measurement function	3-91
Utility menu screen of corneal thickness	
measurement function	3-93
Utility screen of A-scan diagnosis.....	3-94
V	
Various functions in FREEZE mode	3-82
Various functions in real time mode	3-85
Velocity of vitreous	8-1
Version Information	
.....	3-91, 3-92, 3-93, 3-94, 3-110
Volume	3-95
W	
Waveform display.....	3-68
What is a good waveform?	3-37
Wireless Communication	3-98
Wireless communication setting	3-12
Wireless Communication Setting	3-97
Wireless connection with UD-8000 and	
IOL calculation unit.....	3-12

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The image shows a 'Setting change' window with two side-by-side panels for the Right (R) and Left (L) eyes. Both panels are currently set to 'Phakic' eye type. The parameters for both eyes are identical: Avg.V is 1550, LENS.V is 1641, and Vit.V is 1532. Each panel has a button to 'Incorporate' the settings into the opposite eye. At the bottom of the window are three buttons: 'Initial Setting', 'Registration', and 'Cancel'.

Parameter	Right Eye (R)	Left Eye (L)
Eye Type	Phakic	Phakic
Avg.V	1550	1550
LENS.V	1641	1641
Vit.V	1532	1532

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The screenshot shows a connection settings screen with a yellow background. At the top right, it says "No. 1". Below this, there are two input fields: "Name" with the text "IOL_GAL" and "BD Address" with the text "00A0962F9010". Below these fields are two buttons: "Apply" and "Cancel". A virtual keyboard is overlaid on the bottom half of the screen. The keyboard has keys for numbers 1-0, letters QWERTYUIOP, ASDFGHJKL, ZXCVBNM, and special keys like "Delete", "BS", "Enter", "Char/Sign", "Space", "A/a", and "Clear".

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Tomey GmbH