BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4324	Date: 06-DEC-2022	Proofread No.: 1
Designer: CDS	Page: 1 of 20	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4324	Date: 06-DEC-2022	Proofread No.: 1
Designer: CDS	Page: 2 of 20	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4324	Date: 06-DEC-2022	Proofread No.: 1
Designer: CDS	Page: 3 of 20	
Colour Reference:	BLACK	

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CLINICAL BENEFITS

- **BiCart** provides the required buffer in a dialysate proportioning system to remove waste products during hemodialysis when proportioned with an acid concentrate and ultrapure water as evidenced by standard minimal adequacy measurements of minimum of Kt/V of 1.2.
- **BiCart** is used as part of the dialysate system in hemodialysis to correct acid/base imbalances during hemodialysis therapy as evidenced by a target serum bicarbonate level of approximately 22 mEq/L.
- Essential component of a system designed for ultrapure dialysis fluid to avoid treatment-related inflammation.
- **BiCart** decreases the occurrence of microbial contamination using a cartridge directly connected to hemodialysis hardware and does not require the use of a loop system or containers in the mixing and deliver of buffers as evidenced by ISO and European Best Practices Group standards and for Dialysate.

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For Malaysia: The **BiCart** product is intended for preparation of hemodialysis solutions.

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4324		Date: 06-DEC-2022 Proofread No.: 1
Designer: CDS		Page: 4 of 20
Colour Reference: BLACK		

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en

es

DEFINITION OF SYMBOLS			
SYMBOL AND REFERENCE NUMBER	TITLE OF SYMBOL	DESCRIPTION OF SYMBOL	STANDARD
 5.1.1	Manufacturer	Indicates the medical device manufacturer.	ISO 15223-1
 5.1.8	Importer	Indicates the entity importing the medical device into the locale.	ISO 15223-1
 5.1.2	Authorized representative in the European Community/ European Union	Indicates the authorized representative in the European Community/European Union.	ISO 15223-1
 5.1.4	Use-by date	Indicates the date after which the medical device is not to be used.	ISO 15223-1
 5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1
 5.1.6	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1
 N/A	Content weight	Indicates the content weight in grams (g).	N/A
 5.2.8	Do not use if package is damaged and consult instructions for use	Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.	ISO 15223-1
 5.3.1	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully.	ISO 15223-1
 5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture.	ISO 15223-1
 5.3.6	Upper limit of temperature is +40°C (+104°F)	Indicates the upper limit of temperature to which the medical device can be safely exposed.	ISO 15223-1
 5.4.3	Consult instructions for use or consult electronic instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1
 5.4.4	Caution	Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.	ISO 15223-1
 5.7.7	Medical device	Indicates the item is a medical device.	ISO 15223-1
 5.7.10	Unique device identifier	Indicates a carrier that contains unique device identifier information.	ISO 15223-1
 N/A	Recyclable polypropylene	Indicates that the cartridge is made of recyclable polypropylene (PP).	N/A
 N/A	Prescription only	Substitute for statement "Caution: Federal law (USA) restricts this device to sale by or on the order of a [physician]."	21 CFR 801.109(b)(1)
 No. 4	Use no Hooks	To indicate that hooks shall not be used for handling the transport package.	ISO 780
 N/A	Do Not use Knife	To indicate that do not use knife while opening the package.	N/A
 5.1.11	Country of Manufacturer	To identify the country of manufacture of products.	ISO 15223-1
 No. 13	This way up	Indicated the correct, upright position of the package.	ISO 780

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4324	Date: 06-DEC-2022	Proofread No.: 1
Designer: CDS	Page: 20 of 20	
Colour Reference:	BLACK	

Printed on chlorine-free, recyclable paper
Baxter and BiCart are trademarks of Baxter International Inc., or its subsidiaries.

Baxter Healthcare SA,
Thurgauerstrasse 130,
8152 Glattpark (Opfikon),
SWITZERLAND

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MANUFACTURED BY
Gambro Dasco S.p.A.
Via Modenese 66
41036 MEDOLLA (MO)
ITALY

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4326	Date: 17-FEB-2023	Proofread No.: 1
Designer: CDS	Page: 2 of 48	
Colour Reference:	BLACK	

THIS PAGE IS INTENTIONALLY LEFT BLANK

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4326	Date: 17-FEB-2023	Proofread No.: 1
Designer: CDS	Page: 3 of 48	
Colour Reference:	BLACK	

GENERAL INFORMATION

The information in this manual must be carefully read and understood before using the product.

INTENDED PURPOSE

The **BiCart** product is intended for preparation of hemodialysis solutions. The **BiCart** cartridge is also intended for preparation of hemodiafiltration and hemofiltration substitution fluids. Preparation of substitution fluids for hemodiafiltration and hemofiltration has not been cleared in the USA.

The **BiCart** product must always be used together with a suitable acid concentrate.

INDICATIONS

To be used in a bicarbonate dialysis treatment for patients suffering from acute renal failure, chronic renal failure or acute intoxication with dialysable substances.

CONTRA-INDICATIONS

There are no absolute contraindications for bicarbonate dialysis.

DEVICE DESCRIPTION

The **BiCart** cartridge is a polypropylene (PP) cartridge containing dry sodium bicarbonate which enables on-line preparation of a liquid sodium bicarbonate solution used in a bicarbonate based dialysis fluid.

When the **BiCart** cartridge is attached to the cartridge holder, water is drawn by the dialysis machine through the cartridge, producing a saturated solution of sodium bicarbonate. The dialysis machine mixes the sodium bicarbonate solution with water and the acid concentrate to produce a bicarbonate-based dialysis fluid.

In some rare cases the bicarbonate powder may contain discoloured particles or powder being slightly discoloured when exposed to water. These discolourations originate from the raw materials used during the manufacturing process of the sodium bicarbonate powder. Both the powder and the liquid solution obtained from these **BiCart** cartridges meet the demands of the European and US Pharmacopoeias. The safety and efficacy of the **BiCart** cartridge are not affected by these discolourations. This product is non-sterile.

INTENDED USERS

The operation of the specially adapted dialysis machine is restricted to persons who are trained in dialysis and who are completely familiar with the operating instructions stated in the operator's manual of the corresponding dialysis machine and in this Instruction for Use.

INTENDED PATIENT POPULATION

BiCart acts as an essential component in the delivery of a prescribed dialysate composition for hemodialysis and replacement fluid for hemodiafiltration in patients with chronic kidney disease, acute kidney injury, or acute intoxication with dialyzable substances in accordance with hemodialysis parameters prescribed by a physician. It is the physician's decision as to whether the choice of hemodialysis should be limited to any patient population. While many factors must be considered when treating special populations, the normal range of plasma water electrolyte composition is not determined by pregnancy or age.

⚠ WARNINGS

- DO NOT use the products if any damage is found or suspected.
- DO NOT use if cap is missing.
- DO NOT use any other type of acid concentrate than the type specified under Specification.
- DO NOT initiate dialysis without having followed the instructions in the operator's manual of the corresponding dialysis machine regarding any recommended tests to verify the proper composition of the dialysis fluid.

- Proper composition of dialysis fluid must be confirmed prior to dialysis.
- The cartridges shall only be used on dialysis systems that are in compliance with IEC60601-2-16, or other equivalent local standards.
- The cartridges must be used only on the direction of a physician who has evaluated all the pertinent features of the product in relation to the individual patient.
- The Cartridges can only be used together with an appropriate dialysis system equipped with a compatible cartridge holder.
- Federal law (USA) restricts this device to sale by or on the order of a physician.
- DO NOT expose the cartridges to an overpressure higher than 100 kPa or a negative pressure lower than 78 kPa during treatment.
- DO NOT refill the cartridges.
- Between treatments, the cartridges shall be stored in room temperature with caps on the connectors. Aseptic handling is recommended.
- DO NOT perform any unauthorized modifications or alterations to the products or to the dialysis machine as this may result in a malfunction or have other serious consequences on the safe operation of the dialysis treatment.

⚠ CAUTIONS

- The concentrate in the cartridges must be used within 24 hours after being primed.

ADVERSE REACTIONS

Some undesirable effects related to dialysis can occur, for instance hypotension, hypertension, nausea, vomiting, muscle spasm/cramping, cardiac complications, weakness, altered level of consciousness, neurological symptoms, and headache.

SERIOUS ADVERSE EVENT REPORTING

For a user and/or patient if, during the use of this device, or as a result of its use, a serious incident has occurred, please report this incident to the manufacturer, and/or its authorized representative, and to the competent authority of the Member State in which the user and/or patient is established.

RESIDUAL RISKS AND SIDE-EFFECTS

The residual risks and side-effects are present as warnings, cautions and adverse reaction in the appropriate section of Instructions for Use. Apart from these, no additional risks or side effects are expected.

CLINICAL BENEFITS

- **BiCart** provides the required buffer in a dialysate proportioning system to remove waste products during hemodialysis when proportioned with an acid concentrate and ultrapure water as evidenced by standard minimal adequacy measurements of minimum of Kt/V of 1.2.
- **BiCart** is used as part of the dialysate system in hemodialysis to correct acid/base imbalances during hemodialysis therapy as evidenced by a target serum bicarbonate level of approximately 22 mEq/L.
- Essential component of a system designed for ultrapure dialysis fluid to avoid treatment-related inflammation.
- **BiCart** decreases the occurrence of microbial contamination using a cartridge directly connected to hemodialysis hardware and does not require the use of a loop system or containers in the mixing and deliver of buffers as evidenced by ISO and European Best

Practices Group standards and for Dialysate.

SPECIFICATION

Contents

The **BiCart** cartridge contains dry sodium bicarbonate powder in compliance with Ph. Eur. and USP.

Capacity

At a machine setting of 140 mmol/L sodium and 34 mmol/L bicarbonate and together with a suitable acid concentrate, as specified below, up to 32 litres of dialysis fluid can be produced per 100 grams of bicarbonate powder.

Requirements on the acid concentrate

The acid concentrate must be designed for use with a pure sodium bicarbonate concentrate. It must also be suitable for a dilution ratio as recommended by the manufacturer of the dialysis machine to obtain the final dialysis fluid composition.

As an example, when the SoftPac C204 is mixed with a dilution ratio of 1:45, this yields the following concentration in mmol/L:

Sodium	103.0	Potassium	2.0
Calcium	1.75	Magnesium	0.5
Chloride	109.5	Acetic acid	3.0

Together with the sodium bicarbonate concentrate from the **BiCart** cartridge, the final dialysis fluid will have a sodium concentration of 140 mmol/L and a bicarbonate concentration of 34 mmol/L, at a normal machine setting.

Requirements for water

Dialysis machine inlet water: Inlet water quality must comply with local regulations and if no such regulations are available follow ISO 23500-3.

STORAGE

Store below +40°C (+104°F).

SHELF-LIFE

The shelf-life of **BiCart** is 2 years. Refer to the expiry date stated on the unit and box labels.

DISPOSAL

Discard cartridge with cap after use as per local regulations.

INSTRUCTIONS FOR USE

This Instruction for Use must always be used together with the Operator's manual of the corresponding dialysis machine which describes the handling instructions in detail.

FOR BEST PERFORMANCE

Keep the **BiCart** cartridge at +20 to +25°C (+68 to +77°F) for 3-4 hours before use.

Gently shake the **BiCart** cartridge to release any caking of the powder so it is free flowing before use.

ELECTRONIC INSTRUCTIONS FOR USE

The electronic version of the Instructions for Use can be found at the following location: "http://edocs.baxter.com".

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4326	Date: 17-FEB-2023	Proofread No.: 1
Designer: CDS	Page: 42 of 48	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4326	Date: 17-FEB-2023	Proofread No.: 1
Designer: CDS	Page: 43 of 48	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4326	Date: 17-FEB-2023	Proofread No.: 1
Designer: CDS	Page: 44 of 48	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4326	Date: 17-FEB-2023	Proofread No.: 1
Designer: CDS	Page: 45 of 48	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4326	Date: 17-FEB-2023	Proofread No.: 1
Designer: CDS	Page: 46 of 48	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4326	Date: 17-FEB-2023	Proofread No.: 1
Designer: CDS	Page: 47 of 48	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 07-19-00-4326	Date: 17-FEB-2023	Proofread No.: 1
Designer: CDS	Page: 48 of 48	
Colour Reference:	BLACK	

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Baxter Deutschland GmbH