STORZ

KARL STORZ — ENDOSKOPE

GEBRAUCHSANWEISUNG
C-MAC® VS 10331 BX

INSTRUCTION MANUAL
C-MAC® VS 10331 BX

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It is recommended that the suitability of the products for the intended procedure be checked prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favour of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper use, cleaning, disinfection and, if necessary, sterilization of the intubation videoscope. All the necessary details and actions are clearly explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

WARNING: KARL STORZ instruments are not sterile when delivered, and must be cleaned, disinfected and, if necessary, sterilized prior to initial use and each subsequent use.

KARL STORZ is continuously working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures, and descriptions in this manual.

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1	Symbols employed 6
1.1	Packaging symbols.....6
1.2	Symbols for safety instructions 7
2	Intended use 8
2.1	Intended use.....8
2.2	Contraindications 8
2.3	Target groups 8
2.4	User qualification 8
2.5	Training in the operation and function of the device 9
2.6	User profile of physician and assistants 9
3	Safety instructions.....10
4	Handling 14
4.1	Unpacking 14
4.2	Basic equipment.....14
4.3	Possible combinations 14
4.4	Preparation before use 15
4.5	Inspecting the surfaces.....15
4.6	Holding the device correctly.....16
4.7	Assembling the adaptor for tube fixation 16
4.8	Initial operation 17
4.8.1	Connection to the C-MAC® monitor 8403 ZX.....17
4.8.2	Connection to the C-MAC® PM 8403 XD 19
4.8.3	Connection to the C-HUB® II.....20
4.8.4	Focus 20
4.9	Assembling the endotracheal tube.....21
4.10	Handling the deflection mechanism 22
4.11	Multifunction button 'BlueButton' on the C-MAC® VS.....23
4.12	Single image capture 24
4.13	Video recording 24
4.14	Disassembly 24
5	Reprocessing 25
5.1	Overview table 25
5.2	General warnings.....26
6	Precleaning directly after use.....28

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
7	Leakage test 28
7.1	Performing the leakage test – dry 28
7.2	Performing the leakage test – in liquid 29
8	Disassembly 30
9	Manual wipe-down disinfection 30
9.1	Manual reprocessing.....30
9.1.1	Manual cleaning.....30
9.1.2	Manual disinfection 31
9.1.3	Manual High Level Disinfection (HLD).....31
9.2	Drying 31
10	Machine reprocessing.....32
10.1	Connecting 32
10.2	Machine cleaning/chemical-thermal disinfection 32
11	Assembly, inspection and care 33
12	Packaging.....34
13	Sterilization 34
13.1	From SN PR 1191 onwards: Hydrogen peroxide (H ₂ O ₂) sterilization – ASP STERRAD® 35
13.2	From SN PR 1191 onwards: Hydrogen peroxide (H ₂ O ₂) sterilization – STERIS® V-PRO™ 35
13.3	Ethylene oxide sterilization (EO).....35
14	Limits of reprocessing 36
15	Storage 36
15.1	Storage.....37
16	Reprocessing the tube holder.....38
16.1	General warnings.....38
17	Precleaning directly after use.....39
18	Disassembly.....39
19	Manual reprocessing.....39
19.1	Manual cleaning.....39
19.1.1	Manual disinfection 40
19.1.2	Manual High Level Disinfection (HLD).....40
19.2	Drying 40

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20	Machine reprocessing.....	41
20.1	Connecting	41
20.1.1	Machine cleaning/ chemical-thermal disinfection.....	41
20.1.2	Machine cleaning/thermal disinfection	41
21	Assembly, inspection and care	42
22	Packaging.....	42
23	Sterilization	42
23.1	Steam sterilization using the fractionated prevacuum procedure.....	43
23.2	Hydrogen peroxide (H ₂ O ₂) sterilization – ASP STERRAD®	43
23.3	Hydrogen peroxide (H ₂ O ₂) sterilization – STERIS® V-PRO™	44
23.4	Ethylene oxide sterilization (EO).....	44
23.5	Limits of reprocessing.....	44
24	Technical description	45
24.1	System description	45
24.2	Technical data.....	46
24.2.1	C-MAC VS 10331 BX	46
25	Service, repair and disposal	47
25.1	Servicing and repair	47
25.2	Maintenance.....	47
25.3	Limitation of liability.....	47
25.4	Warranty	47
25.5	Electromagnetic Compatibility.....	48
25.6	Directive compliance.....	48
25.7	Important information on dispatch.....	48
25.8	Disposal.....	49
25.9	Troubleshooting	50
26	Appendix.....	52
26.1	Accessories	52
27	Subsidiaries.....	56

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



- ① Connection for leakage tester or pressure compensation cap (vent port)
- ② Handle
- ③ Socket for video connecting cable ②⑥
- ④ Multifunction button (BlueButton)
- ⑤ Flexible sheath for deflection
- ⑥ Distal tip
- ⑦ Deflection lever
- ⑧ Adaptor for tube fixation
- ⑨ O₂ insufflation connection
- ⑩ Leakage tester 13242 XL
- ⑪ Pressure compensation cap 11025 E
- ⑫ Monitor
- ⑬ Secondary connection for videoscope
- ⑭ Primary connection for videoscope
- ⑮ Mains connection
- ⑯ Reset button
- ⑰ Holder for SD memory card
- ⑱ Status bar
- ⑲ Pilot lamp
- ㉑ On/Off switch
- ㉑ Charging indicator lamp (orange)
- ㉒ Soft keys
- ㉓ Function symbols
- ㉔ HDMI port
- ㉕ USB port
- ㉖ Connecting cable
- ㉗ Connection socket on endoscope
- ㉘ Connection socket on monitor

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Fragile

Keep dry

The device is subject to the requirements of the WEEE Directive 2012/19/EU.

Type BF device

Keep away from sunlight

Manufacturer

Humidity limitation

Atmospheric pressure limitation

Number of products in the product packaging

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Scope of protection against moisture:
Protection against continuous immersion

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1.2 Symbols for safety instructions

WARNING: A Warning indicates that the personal safety of the patient or personnel may be endangered. Failure to observe a Warning could result in injury to the patient or personnel.

CAUTION: A caution indicates that particular procedures or precautions must be followed to avoid possible damage to the instrument set.

NOTE: A Note indicates special information regarding operation or clarifies important information.

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2 Intended use

2.1 Intended use

The intubation videoscope with deflecting tip 10331 BX is used in combination with the C-MAC® monitor 8403 ZX and the C-MAC® PM 8403 XD for endotracheal intubation and inspection of the oropharynx.

Videoscopes are applied parts, which must only be used with KARL STORZ supply units which are intended for this purpose.

2.2 Contraindications

No contraindications relating directly to the product are currently known. The intubation videoscope must not be used if, in the opinion of the responsible physician, the patient's health is endangered through its use, e.g., due to the patient's general condition, or if the endoscopic method as such is contraindicated.

2.3 Target groups

This instruction manual is intended for physicians, medical support staff and sterilization assistants.

2.4 User qualification

The intubation videoscope from KARL STORZ may only be used by persons with an appropriate medical qualification and who are acquainted with laryngoscopic techniques. The information given in these instructions only serves to instruct in the correct handling and preparation of the intubation videoscope. It is not intended as an introduction to the techniques of laryngoscopy.

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2.5 Training in the operation and function of the device

Your local representative or responsible KARL STORZ member of staff is available to provide training and information on further training alternatives.

2.6 User profile of physician and assistants

- Recognized medical skills in the field of application (specialist physician, qualified medical staff).
- Adequate powers of comprehension to rationally assess the surgical situation in hand.
- Adequate language skills in one of the languages used on the device and in the instruction manual.
- Be thoroughly trained in the operation of the instruments.
- Knowledge of the contents of the instruction manual.
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic).

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3 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

- Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

Combination of system components

- Combinations of medical devices are only assured to be technically safe if
 - they are listed as such in the respective instruction manuals or
 - the intended purpose and interface specifications of the devices used in combination permit this.The instruction manuals and interface specifications of the medical devices used in combination must be observed precisely.
- No high-frequency instruments may be used in combination with the intubation videoscope.
- Only KARL STORZ devices may be used. Disregarding this warning can endanger the patient and user.
- Only use accessories (tube holders, etc.) from KARL STORZ. Accessories from other manufacturers can damage the intubation videoscope.

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Risks of injury

- i** **NOTE:** All users must have been instructed in the operation and use of the C-MAC® VS, otherwise this endoscope may not be used.
- i** **NOTE:** The C-MAC® VS must be used in accordance with the airway algorithm in the respective country of application or the official specifications for airway management.
- Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.
- Optical radiation poses a risk of injury to eyes. Never look into the free end of either a connected light cable or endoscope. Always wear standardized protective equipment when using lasers.
- The high-intensity light emitted by the C-MAC® VS can result in increased temperatures at the light outlet. However, these always stay below the value which could cause burns.
- Risk of explosion if flammable anesthetics are used in the immediate proximity of the flexible intubation videoscope. The intubation videoscope must not be used in the presence of flammable or combustible media.
- If the flexible intubation videoscope malfunctions during use on the patient, discontinue use immediately. Place the distal tip of the flexible intubation videoscope in the neutral, non-deflected position and remove the intubation videoscope slowly and carefully from the patient.
- Danger of electric shock! Have servicing on these devices carried out by authorized personnel only. Any opening of the device by

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unqualified persons will result in the forfeiting of all warranty claims. Exceptions are those that are described in the instruction manual.

- Before performing any maintenance and cleaning work on the intubation videoscope, remove the connecting cable 8403 X and the tube holder.
- The intubation videoscope must not be used during the discharge of a cardiac defibrillator. Remove the intubation videoscope from the operation area before discharge takes place.
- Do not allow any live components or instruments to come into contact with the intubation videoscope.

Risks of infection

- The intubation videoscope is not sterile when delivered, and must therefore be properly cleaned and disinfected prior to initial use and each subsequent reuse.
- Incorrectly reprocessed instruments expose patients, users and third parties to a risk of infection as well as the danger of an instrument malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

Damage to the medical device

- Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.
- Intubation videoscopes must never be steam sterilized (autoclaved).
- The maximum permissible temperature for decontamination, sterilization, transport and storage is 65 °C.

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- Also do not place the intubation videoscope on other devices that become very hot. Do not expose the device to direct sunlight or excessive heat, to avoid damage to the device.
- Intubation videoscopes must never be exposed to X-rays, to avoid damage.
- Do not hit the distal tip of the flexible intubation videoscope against hard objects, and never use a sharp object to remove dirt.
- Never lay heavy objects on the intubation videoscope.
- To avoid damaging the monitor 8403 ZX, do not remove the SD card or connecting cable whilst video recording.
- Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided.
- The intubation videoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for intubation videoscopes that are repaired by unauthorized personnel.

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4 Handling

4.1 Unpacking

WARNING: The intubation videoscope is not sterile when delivered, and must therefore be properly cleaned/disinfected prior to initial use and each subsequent reuse.

Check for missing items and evidence of shipping damage.

File any complaints immediately with KARL STORZ or the supplier.

4.2 Basic equipment

The set 10331 BXK comprises:

- 1x C-MAC® VS 10331 BX
- 1x Instruction manual
- 1x Leakage tester 13242 XL
- 1x Pressure compensation cap 11025 E
- 1x Tube holder 10331 BAX

4.3 Possible combinations

The C-MAC® VS can be operated with the following monitors:

- C-MAC® motor 8403 ZX (software version: 704v300 or higher)
- C-MAC® PM 8403 XD
- Via the camera control unit C-HUB® II 20 2903 01 (software version: 721v105 or higher) with other monitors

NOTE: All users must have been instructed in the operation and use of the C-MAC® VS, otherwise this endoscope may not be used.

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4.4 Preparation before use

NOTE: Before the device is used for the first time, it must be subjected to a complete reprocessing cycle, excluding preliminary cleaning.

WARNING: Whenever the intubation videoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged intubation videoscopes or damaged accessories must not be used.

4.5 Inspecting the surfaces

WARNING: The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

CAUTION: Do not squeeze the distal tip of the flexible intubation videoscope too tightly.

Every time the intubation videoscope is used, you must first check that its surface, and above all the surface of the sheath and distal tip, is in perfect condition. Intubation videoscopes with sharp edges must not be used. They may cause serious injury to the patient or user.

Check the connecting cable 8403X for damage (such as fractures, cracks, twisted or crushed areas, etc.) and the connecting plug for loose connections.

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4.6 Holding the device correctly

CAUTION: The intubation videoscope must not be carried by the sheath. Do not pull, clamp or twist the sheath of the intubation videoscope.

Take the handle of the intubation videoscope in your right or left hand as shown in the image. With this kind of control grip, the intubation videoscope can be operated with one hand.

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4.7 Assembling the adaptor for tube fixation

If necessary, push the adaptor for tube fixation from distal onto the C-MAC® VS and fix in place by turning the knurled screw clockwise.

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
4.8 Initial operation

4.8.1 Connection to the C-MAC® monitor 8403 ZX

i **NOTE:** The instruction manual of the C-MAC® monitor must be observed! (Art. no. 96076008DER).

Plug cable into monitor

To switch on the device, insert the connecting cable of the C-MAC® VS into the corresponding connection socket ⑭ or ⑬ on the C-MAC® monitor 8403 ZX.

i **NOTE:** The monitor 8403 ZX makes it possible to connect two CMOS video units with Lemo connectors. All CMOS videoscopes from KARL STORZ can be connected to both connections. Video units from other manufacturers cannot be connected.

i **NOTE:** If two cameras are connected, the user must select the camera for which the images are to be displayed in the menu (refer to the instruction manual for the C-MAC® monitor 8403 ZX, Art. no. 96076008DER).

Plug cable into C-MAC® VS

Connect the plug of the connecting cable 8403 X to connection socket of the C-MAC® VS 10331 BX. Note the visible and tangible lug and push together both components until you feel resistance.

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Switch on monitor

Switch on the monitor by pressing the power switch ⑳. The green pilot lamp on the monitor lights up.

The C-MAC® VS is ready for immediate use. and the C-MAC® VS is immediately ready for use. Once the multifunction button on the C-MAC® VS lights up blue, videos and images can be recorded.

If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the sheath of the intubation videoscope out of the patient's body.

NOTE: The connecting cable can be removed during operation and connected to another videoscope.

NOTE: Make sure that the cable is fully inserted. There must not be a gap.

WARNING: It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other units at any time.

WARNING: Optical radiation poses a risk of injury to eyes. Never look directly into the light output at the distal tip.

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4.8.2 Connection to the C-MAC® PM 8403 XD

NOTE: The C-MAC® VS 10331 BX can only be connected to the C-MAC® PM 8403 XD. The 8401 XD model is not compatible.

NOTE: The instruction manual for the C-MAC® PM must be observed (Art. no. 96076020DER).

Connect the C-MAC® PM to the socket on the C-MAC® VS. The orientation pin on the plug and socket facilitates connection. Note the visible and tangible lug and push together both components until you feel resistance.

! WARNING: Optical radiation poses a risk of injury to eyes. Never look directly into the light output at the distal tip.

NOTE: Make sure that the plug is fully inserted. There must not be a gap.

NOTE: The monitor can also be opened up when connecting both components.

The C-MAC® VS is ready for immediate use. As soon as the multifunction button on the C-MAC® VS lights up blue, a blue information bar will appear on the C-MAC® PM monitor for a few seconds. It is now possible to record images and videos.

NOTE: It is not necessary to perform a white balance.

NOTE: The image focus is adjusted automatically. Manual adjustment is not necessary.

NOTE: The C-MAC® PM can be removed during operation and connected to another laryngoscope.

NOTE: After 10 minutes, the display of the C-MAC® PM turns off automatically and must be turned back on again by closing it and reopening it or by disconnecting and reconnecting it (switching it off and on again).

NOTE: It may be more difficult to read the display in bright light, e.g., in direct sunlight.

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4.8.3 Connection to the C-HUB® II

i NOTE: The C-MAC® VS can only be connected to the C-HUB® II 20 2903 01 (software version 721v105 or higher). Older models are not compatible.

i NOTE: The instruction manual of the C-HUB® II must be observed (Art. no. 96206234DER)!

Connect the plug of the connecting cable 8403 X to connection socket 3 of the C-MAC® VS 10331 BX. Note the visible and tangible lug and push together both components until you feel resistance.

i NOTE: Make sure that the cable is fully inserted. There must not be a gap.

4.8.4 Focus

The image focus is adjusted automatically. Manual adjustment is not necessary.

The C-MAC® VS is ready for immediate use. The multifunction button on the C-MAC® VS lights up blue at the same time.

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
4.10 Handling the deflection mechanism

NOTE: The neutral position of the distal tip is straight, without any deflection whatsoever.

WARNING: Risk of injury to users: Before each use, check the functionality of the deflection lever and the connected distal tip. Defective intubation videoscopes must not be used.

NOTE: Test deflection without the endotracheal tube threaded on.

WARNING: Risk of injury to the patient: When inserting the intubation videoscope into the patient, the distal tip must be straight and not deflected.

WARNING: Risk of injury to the patient: Control the deflection up to the stop and avoid rapid or jerky movement so as not to injure the patient or damage the instrument. The resistance that can be felt when deflecting the distal tip using the deflection lever increases as the angle of deflection increases.

WARNING: Before the threaded on endotracheal tube is positioned in the patient via the vocal folds, the distal tip must be brought into the neutral position as otherwise the instrument could be damaged. To do this, release the deflection lever.

The distal tip can be deflected.

- $\geq 90^\circ$ without tube
- $\geq 60^\circ$ with endotracheal tube threaded on

CAUTION: Damage to the instrument. When removing the endotracheal tube, the distal tip must be in the neutral position.

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4.11 Multifunction button 'BlueButton' on the C-MAC® VS

The intubation videoscope features a multifunction button ④.

Color coding meanings for the multifunction button 'BlueButton'

The multifunction button 'BlueButton' ④ on the C-MAC® VS can be used to activate single image capture and video recording. Documentation with the multifunction button on the C-MAC® VS can be performed on the C-MAC® PM 8403 XD and C-MAC® monitor 8403 ZX.

The colors of the multifunction button have the following meanings:

Blue: permanently illuminated. The recording function is enabled.

Green: temporarily illuminated when recording an image, or flashes during video recording.

ⓘ NOTE: The instruction manual for the C-MAC® PM must be observed (Art. no. 96076020DER).

ⓘ NOTE: The instruction manual of the C-MAC® monitor must be observed (Art. no. 96076008DER)!

Approx. 5 seconds after each time the C-MAC® PM is switched on, a blue bar is displayed at the top of the screen for a few seconds. The remaining memory capacity in minutes is shown on the left-hand side, while the number of the directory in which the saved data are stored is shown on the right and the current software version in the middle.

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4.12 Single image capture

Single image capture is activated when the multifunction button ④ is pressed down once briefly. The current image on the screen is saved on the SD card in the connected monitor 8403 ZX / C-MAC® PM 8403 XD. During single image capture, the multifunction button ④ briefly lights up green and then switches back to blue.

4.13 Video recording

Video recording is started by pressing and holding the multifunction button ④ (for approx. 2 seconds) and is stopped by briefly pressing the multifunction button once more. The video recording is saved on the SD card in the connected C-MAC® monitor 8403 ZX or on the internal storage of the C-MAC® PM 8403.

During video recording, the multifunction button ④ on the C-MAC® VS flashes green, switching back to blue when pressed again briefly at the end of recording.

NOTE: The video stream can only be viewed on a PC if an MPEG 4-Codec is installed on the PC.

NOTE: If memory capacity of the C-MAC® PM 8403 XD is less than 10 minutes, the bar is displayed in orange-red.

4.14 Disassembly

Disconnect the connecting cable 8403 X from the C-MAC® monitor or C-HUB II.

Disconnect or separate the connecting cable or C-MAC® PM from the C-MAC® VS.

If necessary, release the adaptor for tube fixation from the C-MAC® VS by turning the knurled screw counter-clockwise. Only guide the tube holder cautiously over the deflection if the knurled screw is undone.

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- This reprocessing step is part of the validated process.
- Only for serial numbers from PR 1191:
This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility, however must be validated by the operator on site in respect of effectiveness.
- x On the grounds of potential damage, this reprocessing step must not be performed.

NOTE: Empty fields are not part of the validated process.

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CAUTION: Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.

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NOTE: We recommend using a plastic basin for disinfection to avoid damaging the medical device. Furthermore, use of a plastic basin prevents electrolytic corrosion which can occur if there are different metals in a solution.

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NOTE: In principle, the C-MAC® VS and C-MAC® connecting cable 8403 X are resistant to manual and machine cleaning and disinfection. If the protective cap is not being used, then manual rinsing and drying of the electrical contacts may be necessary in order to ensure that functioning is not impaired. For this reason, the protective cap 8403 YZ is also offered. The user is responsible for deciding whether or not this cap needs to be used, depending on the individual process management (drying).

NOTE: For information on reprocessing of the C-MAC® monitor 8403ZX, please refer to the corresponding instruction manual 96076008DER.

NOTE: For information on reprocessing of the C-MAC® PM 8403 XD, please refer to the corresponding instruction manual 96076020DER.

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Accessories

Accessories required for carrying out reprocessing: pressure compensation cap, manual leakage tester

Place of application	Art. number brush/adaptor
Clean surface	27652
Ventilation connector	13242 XA
Contacts (optional)	8403 YZ

6 Precleaning directly after use

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the endoscope immediately after use.

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NOTE: Do not use disinfectants for precleaning; these cause fixation of proteins and organic residues.

7 Leakage test

The medical device should be tested to ensure it is leakproof after use and before reprocessing.

7.1 Performing the leakage test – dry

The connection tube of the leakage tester and the vent port on the medical device must be dry. Connect the leakage tester connector securely to the supply connector vent port and lock by turning the bayonet catch clockwise 25°.

NOTE: The rubber seal in the connection adaptor of the leakage tester should be treated regularly with instrument oil.

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Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal end may expand slightly. Using the deflection mechanism on the control body, move the deflection part in all possible directions as far as possible and observe the pointer on the manometer. The test pressure can drop by up to 10 mmHg until the medical device is completely filled with air.

After stabilizing the pressure on the manometer, wait at least 2 minutes to see whether the pressure remains constant or drops.

7.2 Performing the leakage test – in liquid

NOTE: To avoid damage from the penetration of liquid, a test in the dry state must be performed in advance.

Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in the chapter ‘Performing the leakage test – dry’ and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof.

Remove the endoscope from the liquid and press the vent button on the manometer. The test pressure must drop to 0 before the leakage tester can be removed.

NOTE: Never release the pressure in the liquid.

Cleaning can then begin.

If a pressure loss occurs, this means that the endoscope leaks and may no longer be used, immersed in liquid or reprocessed by machine.

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ATTENTION: In this case, fill out the accompanying repairs form and send the medical device, which has been cleaned by wiping to remove organic soiling/residue, for repair or, alternatively, contact KARL STORZ. The endoscope must be labeled as not disinfected. Line the storage/carrying case with the transport protection sheet 13990 SFN.

8 Disassembly

Remove the tube holder prior to cleaning and disinfection.

9 Manual wipe-down disinfection

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications concerning material compatibility and necessary exposure time for the disinfectant must be observed. At the end of the necessary exposure time, wipe the surface with a dry low lint cloth.

NOTE: Validation was performed with the Tristel Trio Wipes System from Tristel Solutions Limited.

NOTE: To avoid damaging the sheath, do not use pushing and pulling movements.

9.1 Manual reprocessing

9.1.1 Manual cleaning

The endoscope must be completely immersed in a cleaning solution. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush.

Finally, it must be rinsed with cold water to ensure neutralization.

ATTENTION: Never remove dirt with a needle or other sharp objects, since these can cause severe damage.

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ATTENTION: The pressure compensation cap must not be in place on the vent port for cleaning and disinfection.

9.1.2 Manual disinfection

The endoscope must be completely immersed in a disinfectant solution. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

9.1.3 Manual High Level Disinfection (HLD)

The endoscope must be completely immersed in a high level disinfectant solution. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

The endoscope can be disinfected with the following high level disinfection solutions which have been tested by way of example:

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9.2 Drying

Finally, the surfaces are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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10 Machine reprocessing

The following method for machine decontamination has been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

10.1 Connecting

In order to ensure effective machine cleaning and disinfection, the working channel of the endoscope must be connected up to the washer and disinfectant.

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10.2 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C and pressure of 100 kPa (1 bar/14.5 psi) must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

The selection of a suitable endoscope holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

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

11 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.

i NOTE: Clean the distal lens with a cotton tip applicator soaked in lens cleaner or 70% isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement.

- Damaged or corroded medical devices must be withdrawn from use.

i NOTE: Inspect plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.

- If there is still residual moisture on the medical device, it must be dried again manually. The surfaces (e.g., connections) must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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- Afterwards, a functional check must be carried out.

i NOTE: Check for a clear view and steerability.

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13 Sterilization

WARNING: Before sterilization, all parts of the medical device must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

WARNING: Risk of infection: Sterilize the medical device in a disassembled state.

ATTENTION: The medical device must never be steam sterilized (autoclaved).

ATTENTION: The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

NOTE: All deviations from the recommended sterilization parameters must be validated by the user.

NOTE: When using reprocessing equipment, please observe the equipment manufacturers' instructions.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

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13.1 From SN PR 1191 onwards: Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

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The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device **from serial number PR 1191** onwards:

- STERRAD® 100S Long Cycle
- STERRAD® NX® Standard Cycle
- STERRAD® 100NX® Standard Cycle

13.2 From SN PR 1191 onwards: Hydrogen peroxide (H₂O₂) sterilization – STERIS® V-PRO™

NOTE: The 'STERIS Low Temperature Reprocessing Device Compatibility Matrix' can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.

The following STERIS V-PRO® sterilization methods have been validated and approved by KARL STORZ for this medical device **from serial number PR 1191** onwards:

- V-PRO® 1 Plus non-lumen Cycle
- V-PRO® 1 lumen Cycle
- V-PRO® 1 Plus lumen Cycle
- V-PRO® maX lumen Cycle
- V-PRO® 60 Flexible Cycle

13.3 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100% ethylene oxide at 55°C and a hold time of 30-45 minutes.

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15 Storage

The storage method must be selected taking into account the respective applicable national requirements.

KARL STORZ only recommends dust-protected, suspended storage for decontaminated endoscopes.

Horizontal storage is possible for sterilized endoscopes.

ATTENTION: Do not expose the medical device to direct sunlight or excessive heat.

ATTENTION: The distal tip of the medical device must be in the neutral position during storage.

ATTENTION: Do not store the medical device in the case.

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10331 BXH

15.1 Storage

CAUTION: Do not expose intubation videoscopes to direct sunlight or excessive heat.

CAUTION: The distal tip of the intubation videoscope must be in the neutral position during storage.

CAUTION: The sheath of the intubation videoscope must never be kinked or coiled up too tightly. Intubation videoscopes should be packaged and stored in sterile cloths or suitable containers in accordance with EN 868, and must be absolutely dry before being put into storage. It is most suitable to store the intubation videoscope suspended vertically with the distal end pointing down. The room where it is stored should be at a normal temperature and humidity.

The intubation endoscope can be safely attached to the mobile stand 8401 YA using the holder 10331 BXH (see Figure). It is secured directly to the cart's equipment rail.

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16 Reprocessing the tube holder**16.1 General warnings**

WARNING: Risk of infection: Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.

WARNING: Risk of infection from non-sterile delivery. Reprocess the medical device before each use. Check the medical device for visible soiling before use.

WARNING: Before reusing the device, replace any defective parts or send the medical device for repair.

ATTENTION: Never use metal brushes to clean the medical device.

ATTENTION: Constant changing between various reprocessing methods is hard on the materials and should be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, is possible.

CAUTION: Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com

NOTE: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

NOTE: Observe the chemical manufacturer's specifications concerning concentration and exposure times. Bear in mind the microbiological range of action of the chemicals.

NOTE: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.

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Accessories

Accessories required for carrying out reprocessing:

Accessories: pressure compensation cap, manual leakage tester

Place of application	Art. number brush/adaptor
Clean surface	27652
Lumen	27650 A

17 Precleaning directly after use

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use.

Wipe-clean the exterior surfaces of the medical device with a damp cloth.

i NOTE: Do not use disinfectants for precleaning; these cause fixation of proteins and organic residues.

18 Disassembly

The tube holder must be removed from the endoscope prior to cleaning and disinfection.

19 Manual reprocessing

19.1 Manual cleaning

The medical device must be completely immersed in a cleaning solution. The lumina must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The lumina are brushed until the brush head is visibly clean. Never brush back and forth within the channel.

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19.1.1 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. The lumina must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the medical device must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

19.1.2 Manual High Level Disinfection (HLD)

The medical device must be completely immersed in a high level disinfectant solution. The lumina must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

The medical device can be disinfected with the following high level disinfection solutions which have been tested by way of example:

- 2.4% glutaraldehyde solution (e.g., CIDEX® Activated Dialdehyde Solution)
- 0.55% ortho-phthalaldehyde solution (e.g., CIDEX® OPA Solution)
- 2.0% accelerated hydrogen peroxide solution (e.g., Revital-Ox® Resert™ or Revital-Ox Resert XL HLD)

19.2 Drying

Finally, the surfaces and lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

ATTENTION: The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

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20 Machine reprocessing

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

i NOTE: Both machine decontamination methods (chemothermal and thermal decontamination) are possible.

20.1 Connecting

In order to ensure effective machine cleaning and disinfection, the medical instrument must be connected up to the washer and disinfectant.

The selection of a suitable device holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

i NOTE: If necessary, the instrument must be dried off manually afterwards.

20.1.1 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C and a maximum pressure of 100 kPa (1 bar/14.5 psi) must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

20.1.2 Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A₀ value must be taken into account when using this method.

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- If there is still residual moisture on or in the endoscope, it must be dried again manually. The surfaces and lumina must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.
- Afterwards, a functional check must be carried out.

22 Packaging

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2 – 10, EN ISO 11607 Parts 1 + 2, DIN 58953).

23 Sterilization

! WARNING: Before sterilization, all parts of the medical device must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

! WARNING: Risk of infection: Sterilize the medical device in a disassembled state.

i NOTE: All deviations from the recommended sterilization parameters must be validated by the user.

i NOTE: When using reprocessing equipment, please observe the equipment manufacturers' instructions.

i NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

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The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The procedure must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

The following sterilization procedures have been validated and approved by KARL STORZ for this medical device:

23.1 Steam sterilization using the fractionated prevacuum procedure

NOTE: Sterilize greased components in a disassembled state.

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132°C–137°C with a minimum exposure time of 4 minutes and maximum of 18 minutes.

NOTE: This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.

or

- 134°C–137°C with a minimum exposure time of 3 minutes and maximum of 18 minutes.

23.2 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

ATTENTION: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

NOTE: Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

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The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S Long Cycle
- STERRAD® NX® Standard Cycle
- STERRAD® 100NX® Standard Cycle

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The following STERIS V-PRO® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- V-PRO® 1 Plus non-lumen Cycle
- V-PRO® 1 lumen Cycle
- V-PRO® 1 Plus lumen Cycle
- V-PRO® maX lumen Cycle
- V-PRO® 60 Flexible Cycle

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23.5 Limits of reprocessing

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

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24 Technical description

24.1 System description

The intubation videoscope from KARL STORZ has a working length of 410 mm and a diameter of 5.5 mm, and is made of stainless steel and a flexible Pellethane sheath. The PEEK (polyetheretherketone) distal tip is 5.5 mm in diameter and houses the CMOS video module. The distal tip of the intubation videoscope can be deflected through $\geq 90^\circ$ without the tube. The handle houses the deflection lever and the function button.

The intubation videoscope from KARL STORZ is connected to the recommended C-MAC® monitor 8403 ZX. If an SD memory card is inserted in the C-MAC® monitor 8403 ZX, images and videos can be saved using the function button on the videoscope and the monitor which are assigned the functions of video recording and image capture.

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24.2 Technical data

24.2.1 C-MAC VS 10331 BX

Direction of view:	0°
Angle of view:	100°
Deflection:	≥ 90° without tube ≥ 60° with tube
Working length:	41 cm
Outer diameter:	5.5 mm
Illumination:	Integrated LED light source
Electrical power consumption:	1 W
Sensor:	CMOS technology
Weight	445 g
Drip water protection:	IPX8
Type of protection against electric shocks:	Protection Class II
Provides Type BF protection when used in conjunction with a KARL STORZ camera control unit with a BF symbol on the connection port.	

Operating/storage conditions

Operating conditions:	0°C to +40°C 30% to 70% rel. humidity
Storage conditions:	-20°C to +50°C 5% to 95% rel. humidity
Atmospheric pressure:	700 to 1080 hPa

CAUTION: The specifications in relation to atmospheric pressure can only be guaranteed if the pressure equalization cap is fitted to the vent port during transportation.

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25 Service, repair and disposal

25.1 Servicing and repair

Defective instruments must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts only.

25.2 Maintenance

Preventive maintenance is not essential. However, regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be requested from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

25.3 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this instrument only subject to the conditions that all assembly, system expansions, readjustments, modifications, or repairs to same have been performed by persons authorized by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

25.4 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see Chapter 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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25.5 Electromagnetic Compatibility

The applied part described in this instruction manual has been tested as a system with the following basic devices. For relevant electromagnetic compatibility (EMC) information, refer to 'Electromagnetic Compatibility (EMC) information' in the instruction manuals for these basic devices

Basic device	Instruction manual
C-MAC® monitor 8403 ZX	96076008 DER
C-MAC® PM 8403 XD	96076020 DER
C-HUB® II 20 2903 20	96206529 DER

The EMC warning statements, precautions, notes and Emission/Immunity limits specified in the instruction manuals for the basic devices also apply to the applied part described in this instruction manual.

25.6 Directive compliance

In accordance with the Medical Device Directive (MDD):
C-MAC® VS 10331 BX:

This medical device belongs to Class I

This medical device bears a CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

25.7 Important information on dispatch

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/devices to the sender. Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no warranty on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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25.8 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the device as electronic scrap. Please ask KARL STORZ SE & Co. KG, a KARL STORZ subsidiary, or your authorized dealer about the collection point in your area. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this unit.

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26 Appendix

26.1 Accessories

8401 YA: Stand,
height 120 cm, rollable, with five legs and anti-static castors, crossbar 25 cm x diameter 25 mm, for positioning the monitor, with tray for laryngoscopes, dimensions (w x d x h): 30 x 20 x 10 cm

8401 YAA: Crossbar,
for stands 8401 YA, 50 cm x diameter 25 mm, for positioning the C-MAC® monitor 8403 ZX

8401 YAB: Crossbar,
for stands 8401 YA, 70 cm x diameter 25 mm, for positioning the C-MAC® monitor 8401 ZX

8401 YB: Multifunctional clamp,
with VESA 75 standard mount. Suitable for mounting on tubes up to 24 mm in diameter and rails up to 25 mm in height

10331 BXH: Holder for C-MAC® VS,
endoscope holder for C-MAC® rigid intubation videoscope 10331 BX, made of Plexiglas, distal open end, for securing to a standard rail

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200200 85: Airway TROLL-E mobile stand,
on 4 antistatic double castors, 2 of which
are lockable, with VESA 75/100 holder for
mounting a monitor, integrated cable channel
in the central upright, max. load capacity of
monitor module 15 kg, dimensions of monitor
cart: 670 x 1660 x 670 mm (W x H x D), castor
diameter: 100 mm.
Mobile stand is not assembled on delivery.

200200 86: C-MAC® TROLL-E,
on 4 antistatic double castors, 2 of which are
lockable, with stainless steel tube,
dimensions of mobile stand: 670 x 1500 x
670 mm (w x h x d),
castor diameter: 100 mm
Mobile stand is not assembled on delivery.

29005AK: Storage basket,
for accessories, for mounting to equipment rail
29003 GS,
dimensions (w x h x d): 300 x 100 x 200 mm,
for use with equipment cart 29005

20020047: Console,
incl. installation accessories, max. capacity
12 kg, dimensions: 455 x 350 mm (w x d), for
use with subracks for mobile Stands 20020060
and 20020061

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopy México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorz@cuba@gmail.com

KARL STORZ Endoscopy México S.A. de C.V.
Lago Constanza No 326,
Col. Granada, Del. Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopy Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA – Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamveien 1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se
KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kęstućio st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-it-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & CO. KG
Interchange 21 Tower, Level 33
399 Sukhumvit Road,
North Klongtoey, Wattana,
10110 Bangkok, Thailand
Phone: +84 28 3823 8000
Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd.
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

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**It is recommended to check the suitability of
the products for the intended procedure prior
to use.**

Thank you for your expression of confidence in
the KARL STORZ brand name. Like all of our
other products, this product is the result of years
of experience and great care in manufacture. You
and your organization have decided in favor of
a modern, high-quality item of equipment from
KARL STORZ.

This instruction manual is intended to serve as an
aid in the proper use, cleaning and sterilization
of the C-MAC® video laryngoscope. All essential
details of the equipment and all actions required
on your part are clearly presented and explained.
We thus ask that you read this manual carefully
before proceeding to work with the instrument.
Keep this manual available for ready reference.

WARNING: KARL STORZ instruments are
provided unsterilized, and must be cleaned,
disinfected and, if necessary, sterilized
prior to the initial use and before each
subsequent use.

KARL STORZ is constantly working on the further
development of all products. Please appreciate
that changes to the scope of supply in form,
equipment and technology are possible for this
reason. Therefore, no claims may be deduced
from the information, figures and descriptions in
this manual.

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

Controls

- ① Laryngoscope blade with holder for electronic module
- ② Electronic module
- ③ Connecting cord for video monitor
- ④ Control key for video recording
- ⑤ Control key for still image capture
- ⑥ Video connection socket on laryngoscope blade
- ⑦ LED illuminator
- ⑧ Image sensor
- ⑨ Electronic module connection
- ⑩ Instrument coupling (C-CAM®)
- ⑪ Focus ring (C-CAM®)
- ⑫ C-MAC® pocket monitor
- ⑬ TFT monitor
- ⑭ Battery charge status display

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Avoidance of environmental pollution through the use of electronic equipment (China RoHS)

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Important: Observe the accompanying documentation

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Read accompanying documents before use

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Charging station for C-MAC® PM

WARNING: Do not touch the patient and the charger plug at the same time. For this reason, the battery should be recharged outside the patient environment.

NOTE: Use only the charging station Art. No. 8401 XDL for recharging. The procedure is described in the instruction manual of the charging station (Art. No. 96076005D).

NOTE: The battery will be approx. 80% charged after just 1 hour. It takes around 3 hours to charge the battery fully.

Insert the C-MAC® PM into the charging compartment.

The charge status display ⑭ of the C-MAC® PM lights up orange.

If the status display lights up green, remove the C-MAC® PM from the charging compartment. Perform a test for proper operation.

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Important note	II
Instrument illustrations	III
Controls	IV
Symbol description.....	VI
Charging station for C-MAC® PM	VII

Safety instructions	2
Normal use	2
User qualification	2
Safety precautions when using	3
Symbols and hazard information	4

Handling	6
Operating the C-MAC®	6
Operating the Electronic Module	8
C-CAM®	8
Function keys for the video laryngoscope.....	10
Function keys for the monitor	10
Inserting an oxygen or suction catheter	11
Disassembly	11
C-MAC® bag	11

Preparation	12
Cleaning and disinfection.....	12
Manual cleaning and disinfection	13
Machine cleaning and disinfection	14
Sterilization	15
Overview: Preparation	15

Maintenance	16
Servicing and repair.....	16
Limitation of liability.....	16
Warranty.....	16
Directive compliance	17
Standard compliance	17
Disposal	18
Repair program	18
Important information	18

Technical description	19
Technical data	19
Troubleshooting.....	21

Appendix	24
Accessories.....	24
EMC	25

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Normal use

The C-MAC® Video Laryngoscope 8401 AX/BX/HX/AXC/BXC/DXC/GXC/KXC/HXP is used in combination with Electronic Module 8401 X/8402 X and Monitor 8402ZX/8403 ZX or the C-MAC® PM for endotracheal intubation and the inspection of the oropharynx.

A separate instruction manual (96076008DER) is available for Monitor 8403 ZX.

Contraindications: No contraindications relating directly to the product are currently known. Use of the C-MAC® video laryngoscope is contraindicated if, in the opinion of the responsible physician, the health of the patient is endangered through its use, e.g., due to the general condition of the patient, or if laryngoscopy itself is contraindicated.

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User qualification

The C-MAC® video laryngoscope may only be used by persons with an appropriate medical qualification and who are acquainted with the laryngoscopic technique. The information given in these instructions only serves to instruct in the correct handling, cleaning, and sterilization of the C-MAC® video laryngoscope. It is not intended as an introduction to laryngoscopy.

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Safety precautions when using the C-MAC® video laryngoscope

The C-MAC® video laryngoscope must be used according to the medical rules and procedures of laryngoscopy recognized for this method.

NOTE: The C-MAC® video laryngoscope has been successfully tested for electro-magnetic compatibility as per DRF EMI Test REPORT (Rev. C) on the BK 117 B2, EC 135 and EC 145 helicopters.

NOTE: The C-MAC® video laryngoscope and the C-MAC® monitor 8402 ZX have been successfully tested for the RTCA/DO-160F (Section 21).

NOTE: The C-MAC® video laryngoscope and the C-MAC® monitor 8402 ZX have passed a fall and crash test successfully.

NOTE: The C-MAC® video laryngoscope and the C-MAC® monitor 8402 ZX have been successfully tested for the MIL-STD-461F, RE102 (Fixed Wing).

THE C-MAC® video laryngoscope and the C-MAC® monitor 8402 ZX have been successfully tested by the German technical center for aircraft "Wehrtechnische Dienststelle für Luftfahrzeuge" (WTD 61) on the Sea Lynx Mk. 88A and Sea King Mk.41 helicopters for use in SAR (Search and Rescue).

Safety precautions when using the C-CAM® and the C-MAC® PM

The C-CAM® with eyepiece telescopes and the C-MAC® PM with C-MAC® intubating video laryngoscope must be used according to the medical rules and procedures recognized for laryngoscopic methods.

NOTE: The C-CAM® and the C-MAC® PM have been successfully tested for the RTCA/DO-160F (Section 21).

NOTE: The C-CAM® and the C-MAC® PM have been successfully tested for the MIL-STD-461F, RE102 (Fixed Wing).

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Symbols and hazard information

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words Warning and Caution stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or personnel may be jeopardized. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular procedures or precautions must be followed to avoid possible damage to the instrument.

NOTE: A Note indicates special information about the operation of the instrument or clarifies important information.

WARNING: The C-MAC® video laryngoscope is delivered unsterilized and should be cleaned, disinfected and/or sterilized prior to initial use and subsequent reuse.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manual or
- the intended purpose and interface specifications of the devices used in combination permit this.

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WARNING: Check the C-MAC® video laryngoscope prior to each application to ensure that it functions correctly and has been properly cleaned, disinfected and/or sterilized.

WARNING: Before each use, the outer surface of the parts of the endoscope and every endoscopically used accessory which is intended to be inserted into the patient must be checked to ensure that there are no accidental uneven surfaces, sharp edges or protruding parts, which could cause harm to the patient.

WARNING: If damage is apparent, the C-MAC® video laryngoscope should no longer be used, if necessary it should be sent off to be repaired.

WARNING: Danger of glare! Never look directly into the light outlet at the tip of the blade.

WARNING: The high-intensity light emitted by the laryngoscope may lead to raised temperatures at the light outlet of the LED.

WARNING: C-MAC® blades contain a nickel compound and can trigger a nickel allergy in allergy-prone people.

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CAUTION: The C-MAC® video laryngoscope must not be cleaned in an ultrasound bath.

CAUTION: The C-MAC® video laryngoscope must not be sterilized with steam.

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Operating the C-MAC® PM

Push the C-MAC® PM ⑫ into the mount on the C-MAC® video laryngoscope ① (photo A).

Switching the OTI* display on

- Flip the OTI* display ⑬ upward (photo A, B).
- Flip the OTI* display to the left (photo C).
- Set the OTI* display so that you can see the image optimally (photo D).

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NOTE: Depending on the composition of the image, there may be a slight image reduction in the first 5-8 seconds.

The C-MAC® video laryngoscope with C-MAC® PM is ready for operation.

! WARNING: After 10 minutes, the display of the C-MAC® PM turns off automatically.

Switching the OTI* display off

- Put the OTI* display ⑬ into the vertical position (photo C).
 - Flip the OTI* display to the right (photo B).
 - Flip the OTI* display downward (photo A).
- The C-MAC® PM is now turned off.

NOTE: You can remove the C-MAC® PM from the C-MAC® video laryngoscope when in operation and insert it in another C-MAC® video laryngoscope.

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Energy management

-  **NOTE:** The battery charge status is shown on the display for the first 10 seconds of operation.
-  **NOTE:** The lithium ion battery can be used for approx. 1 hour when fully charged.
-  **WARNING:** After 10 minutes, the display of the C-MAC® PM turns off automatically and must be turned back on again by closing it and reopening it.
-  **WARNING:** If the charge status symbol flashes red, the charge of the C-MAC® PM is too low. In this state, it will only function for another 10 minutes.

Connect the C-MAC® PM to the corresponding charge station (8401 XDL) when the charge status symbol flashes red at the latest.

-  **NOTE:** A specially designed bag 8402 YE is available for the C-MAC® PM, which can be optionally ordered (see accessories, p. 24). In addition, the C-MAC® PM is part of a set which includes several different C-MAC® video laryngoscopes amongst other items. This concerns article no. 8400 B, "Intubation Set -C22-, ULM model".

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Operating the Electronic Module

Insert the electronic module ② into the receptacle on the video laryngoscope ①.

! NOTE: The electronic module can be removed from the laryngoscope and inserted into another laryngoscope during operation.

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C-CAM®

To connect an endoscope or fiberscope, rotate the outer ring of the coupling against the spring force (in the direction of the arrow), insert the DIN endoscope or fiberscope eyepiece, and release the outer ring. The spring force enables the instrument coupling to automatically click into place.

! NOTE: All KARL STORZ endoscopes and fiberscopes can be used.

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! CAUTION: No further tightening of the outer ring is required to position the eyepiece. This may cause damage to the instrument coupling!

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Connect the light cable to the endoscope (twist the knurled screw on the screw base through a quarter turn).

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Plug the connecting cable ③ of the C-MAC® electronics module into the corresponding connection socket of the monitor.

(Cf. further connection and application options of monitor in the instruction manual "C-MAC® Monitor 8403 ZX", art. no. 96076008DER).

The connecting cord ③ can be connected or disconnected during operation.

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Insert the cable from the power supply into the corresponding connection socket of the monitor and connect the power supply with the mains socket.

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Turn the monitor on using the power switch (Indicator light lights up).

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Focus

When using C-MAC® intubating video laryngoscopes and the C-MAC® S IMAGER with intubating video laryngoscopes for single use, the image is focused automatically. Manual focusing is not possible/necessary. When using a C-CAM®, the focus can be set by turning the silver focus ring ⑪ on the camera head. The system is now ready for use.

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Function keys for the video laryngoscope

Video recording: When the function key "video recording" ④ is pressed, the C-MAC® video laryngoscope saves a video stream on the SD memory card incorporated in the connected Monitor 8403 ZX. Pressing the key a second time stops the recording.

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NOTE: Video stream playback on a PC is only possible if an MPEG-4 codec is installed on the PC in question.

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Single image capture: When the function key "single image capture" ⑤ is pressed, the current image on the monitor is saved on the SD memory card incorporated in the connected Monitor 8403 ZX.

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Function keys for the monitor

The monitor is equipped with 4 softkeys that can be used to activate several other functions.

A more detailed description of the functions can be found in the instruction manual for the monitor "C-MAC® Monitor 8403 ZX" (art. no. 96076008DER).

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Inserting an oxygen or suction catheter (only with 8401 AX/BX/HX)

When using a blade with an oxygen or suction catheter guide, the probe or catheter is advanced so far beyond the guide until visibility in the image is guaranteed.

When using the blade 8401 AX, the 14 Fr (green) and the 16 Fr (orange) catheters can be used. When using the blade 8401 HX and 8401 BX, the 16 Fr (orange) and the 18 Fr (red) catheters can be used.

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Disassembly

Disconnect the connecting cord from the C-MAC® monitor port.

Push the C-MAC® electronics module 2 (Fig. A) or the C-MAC® PM ⑫ (Fig. B) out of the video laryngoscope.

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C-MAC® bag 8403 YD

The optionally available protective bag 8403 YD facilitates the mobile use of the C-MAC® system.

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Cleaning and disinfection

CAUTION: The instructions 'Cleaning, Disinfection, Sterilization and Care of KARL STORZ Instruments' must be followed. There, the cleaning, disinfection and sterilization procedures are explained in detail. The list of KARL STORZ approved materials can be found on the Internet at 'www.karlstorz.com'.

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
NOTE: As a general rule, both the connector and its associated sockets in the laryngoscope and in the electronics module are completely waterproof. All three connections are also resistant to the process chemicals of the validated preparation procedure used. During the validated preparation procedures, these connections can generally be left unprotected. No further accessories are required, particularly if the storage basket 39501 LC2 is used in a preparation machine.

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CAUTION: To avoid possible damages, these connections may only be assembled in dry conditions. Where a preparation is carried out which does not involve drying or during which no compressed air drying can be carried out, the reprocessing cap 8401 YZ must be used.

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The **C-MAC® laryngoscopes** and the **electronics modules 8401 X/8402 X** are suitable and validated for the following low-temperature preparation methods up to max. 65°C: Manual/machine cleaning and disinfection, sterilization with Steris® AMSCO V-PRO 1, Sterrad® (100/NX), and EO gas.

The **monitor 8402 ZX/8403 ZX** can be wipe-down disinfected.

The **C-MAC® PM 8401 XD** is suitable and validated for the following low-temperature preparation methods up to a maximum of 65 °C: Manual/machine cleaning and disinfection:

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NOTE: Preparation using STERRAD® can cause cosmetic changes to the surface of the blade.

CAUTION: The laryngoscope and the electronic module are fitted with gold contacts which prevent corrosion. Particularly with manual cleaning/disinfection, you must ensure that, in the last rinse, the contacts are rinsed thoroughly with microbiologically clean/sterile water. Then the contacts and instrument must be dried carefully, preferably using sterile compressed air.

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Manual cleaning and disinfection

Reprocessing of the instrument should begin as soon as possible after use, as incrustations may cause irreparable damage. Remove heavy soiling using an active cleaner and a disinfected sponge or soft brush.

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Immerse in a disinfectant solution, making sure that all surfaces are covered and that all air bubbles have been removed. Follow the chemical manufacturer's instructions on concentration and exposure time.

Rinse thoroughly with completely demineralized water and then dry.

NOTE: The electrical contacts in particular must be dried carefully.

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Machine cleaning/Disinfection

The choice of washing programs must be agreed with the washing machine manufacturer.

CAUTION: Clean and chemically disinfect the C-MAC® video laryngoscope in the washer-disinfector at a maximum temperature of 65 °C. The temperature of the drying program may not exceed 65 °C.

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In the case of the C-MAC® video laryngoscope, clean the optical surfaces of image chip ⑧ and LED light ⑦ with a cotton tip applicator soaked in 70% isopropyl alcohol.

NOTE: The electrical contacts in particular must be dried carefully.

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Sterilization

CAUTION: The C-MAC® video laryngoscope, the electronic module 8401 ZX/8402 ZX/8403 ZX and the C-MAC® PM 8401 XD must not be steam sterilized.

NOTE: In the case of gas or Steris System 1® sterilization, the country-specific laws and regulations must be observed.

CAUTION: Reprocessing cap 8401 YZ must be replaced before sterilization with the Steris System 1® and Steris System 1E®.

Gas or plasma sterilize the C-MAC® video laryngoscope and store disassembled.

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✓ Suitable procedures
---- Unsuitable procedures

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts.

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Limitation of liability

As the supplier of these medical devices we only regard ourselves responsible for the safety, reliability and performance of the products if assembly, system expansions, readjustments, modifications, or repairs have been performed by people authorized by KARL STORZ, and the device has been used in accordance with its operating instructions at all times.

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Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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Directive compliance

In accordance with the Medical Device Directive (MDD):

Medical device Class I
This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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In accordance with EN 60601-1, EN 60601-2-18, UL 2601, CSA 601.1:

- Type of protection against electric shock: Protection class II (not 8401 XD)
- Degree of protection against electric shock: Applied part of type BF

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DO-160F

In accordance with IEC 60601-1-2:

Please read the Electromagnetic Compatibility Information in the appendix.

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Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its service life, dispose of the device as electronic waste. Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information about your local collection point. Within the scope of application of this directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this device.

To dispose of C-MAC® laryngoscope blades for single use, no special measures are necessary. National laws and regulations must be observed.

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Repair program

To bridge the repair period, you will receive a unit on loan, which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Straße 34
78532 TUTTLINGEN

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

For hygiene reasons or to prevent infection, telescopes, instruments and units must be cleaned, disinfected/sterilized before they are sent for repair.

We reserve the right to return contaminated instruments/units to the sender.

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Technical data

Electronic module 8401 X/8402 X

Classification:	IEC Class II
Frequency:	DC
Voltage:	5V
Power consumption:	2 VA
Protection against submersion	IPX8

Blade 8401 AX/BX/HX/AXC/BXC/DXC/GXC/KXC/HXP

Protection against submersion	IPX8
Camera Technology	CMOS
Camera resolution	320 × 240
Illumination	LED, white, 1 W

Color settings Basic Setting

8403 ZX:

Brightness	40%
Contrast	35%
Saturation	40%
Hue	50%

8402 ZX:

Brightness	30%
Contrast	25%
Saturation	40%
Hue	100%

8401 ZX:

Brightness	56%
Contrast	50%
Saturation	20%
Hue	0%

Operating/storage conditions for C-MAC®

Operating conditions:	0 °C – +40 °C
Storage conditions:	-20 °C – +50 °C
	5% – 95% rel. humidity

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Charging station 8401XDL

Classification:	IEC Class II
Line frequency:	DC
Voltage:	5V
Power consumption:	2 VA
Weight:	approx. 230 g

C-MAC® PM 8401 XD

Frequency	DC
Voltage	3,7V
Power consumption	1.5 W
Protection against submersion	IPX8

C-CAM® 202901 31/ 202901 32

Camera technology	CMOS
Camera resolution	640 x 480
Protection against submersion	IPX8

Operating/ storage conditions for C-MAC® PM

Operating conditions:	0 °C – +50 °C 30% – 70% rel. humidity
Storage condiditons:	-20 °C – +50 °C 5% – 95% rel. humidity

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Error description:

Cloudy picture, stripes, streaks, or similar.

Possible causes:

- Telescope of the laryngoscope is soiled.
- Contacts of electronic module/C-MAC® PM are soiled.

Remedy:

- Clean with cotton swab and alcohol solution or special cleaning paste.
- Clean the contacts of the electronic module/C-MAC® PM.

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Color distortions.

Possible causes:

- White balance not carried out correctly.
- Electronic module defect.

Remedy:

- Carry out new white balance.
- Send supply unit to KARL STORZ for repair.

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- Camera connecting cable is defect.

Remedy:

- Have a new camera connecting cable fitted.

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Image cannot be saved.

Possible causes:

- No memory card inserted
- Memory card full
- Memory card is not recognized

Remedy:

- Insert memory card
- Change memory card
- Insert the memory card in the switched-off C-MAC® monitor and turn the monitor on again.
- Perform reset
- Format the memory card on the PC with FAT (standard).

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Abhilfe:

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Error description:

Video stream cannot be played back on the PC monitor.

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SD memory card cannot be written on.

Possible causes:

- SD memory card is formatted incorrectly

Remedy:

- Format the SD memory card on the PC with FAT (standard).

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Forceps designed by MAGILL, modified by BOEDEKER, length 25 cm, suitable for the endoscopic removal of foreign bodies, for use with video laryngoscopes of sizes 2 – 4

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
Forceps designed by MAGILL, for children, modified by BOEDEKER, length 20 cm, suitable for the endoscopic removal of foreign bodies, for use with video laryngoscopes of sizes 1 and 2

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Multifunctional cap for use with Steris System 1® to protect the contacts during reparation. Suitable for the C-MAC® video laryngoscope, the C-MAC® electronics module, and the C-MAC® PM. The cap is reusable.

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Bag for intubation set -C22-, ULM model, made from water resistant, hard-wearing material, washable, the bag comprises two compartments with several storage options for C-MAC® video laryngoscopes with C-MAC® pocket monitor and for conventional laryngoscopes, for use with C-MAC® pocket monitor 8401 XDK, C-MAC® video laryngoscopes and conventional laryngoscopes

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Instrument tray for the cleaning, sterilization and storage of two C-MAC® video-laryngoscope blades, includes an electronic module, with holder for securing and sealing the electronic connections.
Overall dimensions (w x d x h):
260 x 120 x 170 mm

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WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this Appendix during installation and commissioning.

The C-MAC® video laryngoscope 8401 AX/BX/HX/AXC/BXC/DXC/GXC/KXC/HXP in conjunction with the electronics module 8401 X /8402 X and the monitor 8401 ZX/8402 ZX/8403 ZX as well as in conjunction with the C-MAC® PM 8401 XD correspond to IEC 60601-1-2 : [CISPR 11 Class B] and therefore meet the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The aforementioned devices and combinations of devices are Group 1 devices (as per CISPR 11).

Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

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NOTE: The tables and guidelines that are included in this Appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

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WARNING: The use of accessories other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the C-MAC® video laryngoscope 8401 AX/BX/HX/AXC/BXC/DXC/GXC/KXC/HXP in conjunction with the electronics module 8401 X/8402 X and the monitor 8401 ZX/8402 ZX/8403 ZX as well as in conjunction with the C-MAC®-PM 8401 XD. The accessories listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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Table 200
Accessories and cables which have been shown to comply with EN/IEC 60601-1-2:

Cable type	Shielded	Length [m]	Ferrite	Used for
Charger 40 1500 31	No	2.0	No	Charging the internal battery of the C-MAC® monitor 8401ZX/8402ZX
Charger ET27-30-0004515	No	2.0	No	Charging the internal battery of the C-MAC® monitor 8403 ZX
Charging station 8401 XDL	not applicable	Approx. 1.8	No	Charging the internal battery of the C-MAC® PM

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Table 201 Guidance and manufacturer's declaration – electromagnetic emissions		
The C-MAC® video laryngoscope 8401 AX/BX/HX/AXC/BXC/DXC/GXC/KXC/HXP in conjunction with the electronics module 8401 X/8402 X /8403 ZX and the monitor 8401 ZX/8402 ZX as well as in conjunction with the C-MAC® PM 8401 XD is intended for use in the environment specified below. The customer or user of the system should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Complies Class 1	The aforementioned devices and combinations of devices only use HF energy for their internal functioning. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Complies Class B	The aforementioned devices and combinations of devices are suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A. Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 202
Guidance and manufacturer's declaration – electromagnetic immunity

The C-MAC® video laryngoscope 8401 AX/BX/HX/AXC/BXC/DXC/GXC/KXC/HXP in conjunction with the electronics module 8401 X/8402 X and the monitor 8401 ZX/8402 ZX/8403 ZX as well as in conjunction with the C-MAC® PM 8401 XD is intended for use in the electromagnetic environment specified below. The customer or user of the system should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. If image distortion occurs it may be necessary to further filter the power quality by using a AC line filter.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Complies <5% U_T * (>95% dip in U_T) for 0.5 cycle Complies 40% U_T (60% dip in U_T) for 5 cycles Complies 70% U_T (30% dip in U_T) for 25 cycles Complies <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

* Note: U_T is the a.c. mains voltage prior to application of the test level.

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Table 204 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
The C-MAC® video laryngoscope 8401 AX/BX/HX/AXC/BXC/DXC/GXC/KXC/HXP in conjunction with the electronics module 8401 X/8402 X and the monitor 8401 ZX/8402 ZX/8403 ZX as well as in conjunction with the C-MAC® PM 8401 XD is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			PoPortable and mobile RF communications equipment should be used no closer to any part of the aforementioned devices or combinations of devices, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	$d = [3.5/3]\sqrt{P}$ 150 kHz to 80 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = [3.5/3]\sqrt{P}$ 80 MHz to 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz
			where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.			

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Table 206

**Recommended separation distances between portable and mobile RF communications equipment and the
KARL STORZ C-MAC® video laryngoscope 8401 AX/BX/HX/AXC/BXC/DXC/GXC/KXC**

The C-MAC® video laryngoscope 8401 AX/BX/HX/AXC/BXC/DXC/GXC/KXC/HXP in conjunction with the electronics module 8401 X/8402 X and the monitor 8401 ZX/8402 ZX/8403 ZX as well as in conjunction with the C-MAC® PM 8401 XD is intended for use in an electromagnetic environment in which the HF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.1	0.1	0.2
0,1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@kse.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326 , Col. Chapultepec Morales
D.F.C.P. 11520, México, México
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14,
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4. kerros, 00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygang Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

KARL STORZ Endoscopy Kasachstan LLP
6, Saryarka str. BC „Arman“, office 430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) – Units 205 – 206
Dubai Health Care City, P.O.Box 118069
Dubai - United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 43743000-30,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Suite 2503-5, 25th Floor, AXA Tower, Landmark East,
100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016
People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisin Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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KARL STORZ Endoscopy Canada Ltd.

7171 Millcreek Drive, Mississauga, Ontario L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free (Canada only) Phone: 1-800-268-4880, Fax: 1-800-482-4198
(Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.

2151 East Grand Avenue, El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free (USA only) Phone: 800 421-0837, Fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.

1 South Los Carneros Road, Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopia Latino-America, Inc.

815 N. W. 57th Avenue, Suite 480, Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.

Edificio Atlantic, Oficina 3G, Calle D e/ 1ra y 3ra, 10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.

Av. Ejercito Nacional No. 453 Piso 2, Colonia Granada, Alcaldia Miguel Hidalgo, C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.

Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi, CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.

Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez, Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS

Stamveien 1, 1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB

Storsåtragränd 14, 127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY

Taivaltie 5, 01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG

Representation Office
Kęstučio st. 59 / Lenktoji st. 27, 08124 Vilnius, Lithuania
Phone: +370 5 272 0448, Mobile: +370 685 67 000
E-Mail: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S

Skovlytoften 33, 2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.

415 Perth Avenue, Slough, Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

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**Sales for Industrial Endoscopy

ООО KARL STORZ Endoscopy – WOSTOK

Derbenyevskaya nab. 7, building 4, 115114 Moscow, Russia

Phone: +7 495 983 02 40, Fax: +7 495 983 02 41

E-Mail: Info-ru@karlstorz.com

TOV KARL STORZ Ukraine

Obolonska naberezhna, 15 building 3, office 3, 04210 Kiev, Ukraine

Phone: +380 44 42668-14, -15, -19 Fax: +380 44 42668-20

E-Mail: marketing@karlstorz.com.ua

KARL STORZ SE & Co. KG Representation Office

Sabit Orudschow 1184, apt. 23, 1025 Baku, Azerbaijan

Phone: +99 450 613 30 60

E-Mail: info-az@karlstorz.com

KARL STORZ ENDOSKOPE – East Mediterranean and Gulf (Offshore) S.A.L.

Spark Tower 1st floor Charles Helou St., Horch Tabet – Sin El Fil, Beirut,
Lebanon

Phone: +961 1 501105, Fax: +961 1 501950

E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.

P.O. 6061, Roggebaai, 8012 Cape Town, South Africa

Phone: +27 21 417 2600, Fax: +27 21 421 5103

E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan

Saryarka, 6, BC "Arman", off. 910, 010000 Astana, Republic of Kazakhstan

Phone: +7 7172 552-549, 552-788, Fax: -444

E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf (branch)

Building West Side 7A – Unit 7WA – 3008, Dubai Airport Free Zone, P.O. Box

54983, Dubai - United Arab Emirates

Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282

Service Hotline: +971 (0)4 3415882

E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited

11th Floor, Dr. Gopal Das Bhawan, 28, Barakhamba Road, New Delhi

110001, India

Phone: +91 11 4374 3000, Fax: +91 11 4374 3010

E-Mail: corporate@karlstorz.in

KARL STORZ SE & CO. KG

Interchange 21 Tower, Level 33, 399 Sukhumvit Road, North Klongtoey,

Wattana, 10110 Bangkok, Thailand

Phone: +84 28 3823 8000 Fax: +84 28 3823 8039

E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG

Resident Representative Office

14th Floor, MPlaza Saigon, 39 Le Duan, District 1, Ho Chi Minh City, Vietnam

Phone: +84 28 3823 8000, Fax: +84 28 3823 8039

E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.

Room 2503-05, 25F AXA Tower, Landmark East, No. 100 How Ming Street,

Kwun Tong, Kowloon, Hong Kong, People's Republic of China

Phone: +852 28 65 2411, Fax: +852 28 65 4114

E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch

Room 1805-1807, Building B, 18F Beijing IFC, No. 8, Jianguomenwai Street,

Chaoyang District, 100022, Beijing, People's Republic of China

Phone: +86 10 5638188, Fax: +86 10 5638199

E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch

Room 701A Building 5 & Room 501 Building 7, No. 3000 Longdong Avenue,

Pilot Free Trade Zone, 201203, Shanghai, People's Republic of China

Phone: +86 21 60339888, Fax: +86 21 60339808

E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch

Room 803-805, 8F Jin Jiang International Building, No. 1 West Linjiang

Road, Wuhou District, 6100414, Chengdu, People's Republic of China

Phone: +86 28 86587977, Fax: +86 28 86587975

E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch

Room 2001-2005, 20F N-MEDIA International Center, No. 167 Youth Avenue,

Shenhe District, 110014, Shenyang, People's Republic of China

Phone: +86 24 23181118, Fax: +86 24 23181119

E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch

Room 02B & 03 & 04A, 35F Teem Tower, No. 208 Tianhe Road, Tianhe

District, 510620, Guangzhou, People's Republic of China

Phone: +86 20 87321281, Fax: +86 20 87321286

E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.

No. 8 Commonwealth Lane #03-02, Singapore 149555, Singapore

Phone: +65 69229150, Fax: +65 69229155

E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd

No. 8 Commonwealth Lane #03-02, Singapore 149555, Singapore

Phone: +65 69229150, Fax: +65 69229155

E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.

9F Hyowon-Building, 97, Jungdae-ro, Songpa-gu, 05719 Seoul, Korea

Phone: +82-70-4350-7474, Fax: +82-70-8277-3299

E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.

12F, No. 192, Sec. 2, Chung Hsin Rd., Sindian District, New Taipei City,

Taiwan

Phone: +886 933 014 160, Fax: +886 2 8672 6399

E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG Representative Office Philippines

1901 Picadilly Star Bldg., 4th Avenue, BGC, Taguig City 1636, Philippines

Phone: +63 2 317 45 00, Fax: +63 2 317 45 11

E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.

Stage Bldg. 8F, 2-7-2 Fujimi, Chiyoda-ku, Tokyo 102-0071, Japan

Phone: +81 3 6380-8622, Fax: +81 3 6380-8633

E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd .

68 Waterloo Road, Macquarie Park NSW 2113, P O Box 50 Lane Cove NSW

1595, Australia

Phone: +61 (0)2 9490 6700, Fax: +61 (0)2 9420 0695

Toll free: 1800 996 562 (Australia only)

E-Mail: info@karlstorz.au

www.karlstorz.com

KARL STORZ SE & Co. KG

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com

96076019FU • V4.3 • 12/2019

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Авторское право[©]