
NUGA
www.nuga.kr

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
T5

User Manual

Thank you for purchasing this product.
Before use...

Please read the cautions for safety in the
user manual to use this product safely.

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
NUGA MEDICAL CO., LTD.

Product Introduction

Product Components	8
Names and Functions of Product	9
Descriptions for Controller Buttons	10

How to Use Product

How to Use Product	12
How to Take Emergency Actions during Use	16

Product Storage and Care

How to Store and Clean after Use	18
Self-check Items	20
Before any A/S request	21
Check Indicator Window of Controller	22
Product Specifications	23
Product Warranty	24

Cautions for your safety

- These cautions include descriptions for the safety of users and prevention of any damage to properties among others.
- Read through the following details to use the product appropriately.
- This Manual classifies the risk level into following three degrees and the following graphical symbols remind users of cautions for matters and operations of concerns about any risk.

Interpretations on symbols

Danger

If there is high possibility of death due to non-compliance.

Warning

If there is a possibility of an injury or death due to non-compliance.

Caution

If there is a possibility of a minor injury or damage to goods due to non-compliance.

Disclaimer

- The manufacturer, importer and/or agent shall not be responsible for any problem and/or other damage to the product which has been caused from the user's fault or operation.
- Information in the user manual has been prepared on the basis of the current product specifications.
- Details in the user manual might be different from those on the product.
- The product manufacturer will be adding and supplementing new functions at any time and new technologies shall be applied to the product continuously in the future as well.
- All the standards and functions are subject to change without any notice to users to improve performance.


Danger

Do not fold the product in use.
(Otherwise, it might cause damage to the product.)


Avoid pulling out the power cord parts or touching the power plug with the wet hands.
(Otherwise, it might cause electrical shock or fire.)

Do not use it for other usages than the described purposes of uses.
(Otherwise, it might cause damage to the product.)

Anyone using a cardiac pacemaker should consult with a doctor before any use of this product.


Do not apply washing or dry-cleaning.
(Otherwise, it might cause damage to the product or any danger of electric shocks or fire.)

Anyone other than technicians qualified for our company's A/S should not disassemble, repair or modify the product for any reason.
(Otherwise, it might cause damage to the product or electric shock and/or fire.)


Warning

Do not allow children use the product alone.

Do not place any heavy object or other electric heater on this product.
(Otherwise, it might cause injuries or damage to the product.)

Do not place the controller covered with a blanket during use.
(Otherwise, it might cause damage to the product or a danger of fire.)

Avoid bending the power cord excessively or pressing it down with heavy objects.

Never use other temperature controllers with different specifications.
(Otherwise, it might cause damage to the product.)

Avoid placing inflammable substances like candlelight, cigarette light or liquors and/or vessels containing water, heavy objects or other objects on the product.
(Otherwise, it might cause fire.)

Caution


Please abide by defined time and methods.


Separate the power cord while the product is not being in use for a long time.

Keep any pin, wire or other metallic pieces from stabbing or inserting to the mat.
(Otherwise, it might cause damage to the product.)


Use the product in a proper range of temperature.
(Otherwise, it might cause burns.)

Do not expose the product to the direct sunlight.
(Otherwise, it might cause damage to the product.)

Check out the applied voltage before any use and apply the power plug only to a dedicated socket.

Product Introduction

Product Components 8

Names and Functions of Product 9

Descriptions for Controller Buttons 10


Product Components

i Introduction

Please check out whether all the components are unpacked in hand before installing the product.

Please inquire from of your supplier if there is any missing or damaged component.

Some specifications and/or designs are subject to change for improving the product performance.

Device


User manual


Device Bag

Names and Functions of Product


Name	Function
① Heater	The heater as the main component in the product improves the heating effects by combining the tourmanium mat.
② Power cord	It is connected to a socket to input the power supply to the product.
③ All-in-one controller	It is used to regulate the temperature of the heating unit.

Descriptions for Controller Buttons

Name	Function
① Power display lamp	A display that indicates whether the controller is connected to a power supply
② Power button	A button used to turn the power on or off for the device and controller
③ Time display lamp	An LED lamp that indicates the time set for the device operation
④ Time control button	A button used to set the time for the device operation (3, 6 or 9 hours)
⑤ Temperature display lamp	An LED lamp that indicates the temperature set for the device operation
⑥ Temperature control button	A button used to set the temperature for the device operation


How to Use Product

How to Use Product	12
How to Take Emergency Actions during Use	16

How to Use Product

i Before using product



- Do not use it for other usages than the described purposes of uses.
- Do not fold the product in use.
- Check out the voltage in use before any use.
- Anyone using a cardiac pacemaker should consult with a doctor before any use of this product.
- Do not place the controller covered with a blanket during use.
- Use the product in a proper range of temperature.

step. 1

When you connect the power cord to a power outlet, the power display lamp (LED) will turn red.



If the red LED lamp does not turn on, then check whether the power cord has been inserted into the power outlet properly.

step. 2

When you press the power button once, the power display lamp (LED) will turn green, and the device will enter the standby mode.

When the device is in standby mode, the power display lamp (LED) will be green.



step. 3

Press the time control button on the controller to set the desired device operation time.

A single press: 9 hours
Two presses: 6 hours
Three presses: 3 hours

The time display lamp (LED) will turn white.
The set time will be shown on the time display.


step. 4

Press the temperature control button on the controller to set the desired device operation temperature.

Temperature range
:40~70°C



The temperature display lamp (LED) will turn white.

The set temperature will be shown on the time display.

step. 5


After setting the desired temperature, place the device on the area of body that you wish to massage.

step. 6

To stop the device operation, press the power button and the device will enter the standby mode.

The power display lamp (LED) on the controller will turn red.


step. 7

Remove the power cord from the socket in any case after using the product.


i Notice

The power supply on the product is automatically turned off after 9 hours.

How to Take Emergency Actions during Use


Stop using the product immediately and remove the power plug if there is any lightning.


Stop using the product immediately and contact your supplier branch or our company if you smell any burn or detect any smoke from the product.


Stop using the product immediately and consult with a doctor if any abnormality occurs to the human body.

Turn off the power switch immediately, stop using the product and contact our company if you find out any failure or power failure from the product.

Turn off the power switch, remove the power plug immediately, and then contact our company if you have poured water into the product.

Product Storage and Care

How to Store and Clean after Use	18
Self-check Items	20
Before any A/S request	21
Check Indicator Window of Controller	22
Product Specifications	23
Product Warranty	24

How to Store and Clean after Use

i Introduction

1. Separate the power cord while the product is not in use for a long time.
2. Be aware of the following for the storage place.
 - ① Keep the product in a place with no water.
 - ② Keep the product in a place with no concern of any adverse effect such as the air pressure, temperature, humidity, ventilation, sunlight and/or air containing dust, salt, sulfur, etc.
 - ③ Be careful of safety conditions such as the slope, vibration, shock, etc. (including those during transportation).
 - ④ Do not keep the product in a place near any chemical or gas.
 - ⑤ Avoid any direct sunlight during the storage.
3. Store the product after cleaning the components.

Environmental conditions

- During operation: temperature (0~35 °C), humidity (20~80%)
- Transport and storage conditions: temperature (-25~70 °C), humidity (20~80%), atmospheric pressure (700hPa to 1060hPa)

Proper Disposal of the Device


- This system contains parts and components that must be disposed of as industrial wastes, and an inappropriate disposal thereof will lead to environmental pollution. Thus, do not discard this device together with general industrial wastes or household wastes.
- You must comply with the relevant laws and regulations when disposing of this system, in part or in whole. For inquiries regarding proper disposal, please contact the authorized retailer at a location near you.

How to clean


Caution

Ensure unplugging of the power plug before cleaning. Otherwise, it might cause electric shocks.


Proper Storage Method

step. 1


Unplug the device before storage.

step. 2


Plug the power cord properly into the hole on the back of the device.

Make sure that the power cord does not become twisted. There is a risk of wire disconnection.

step. 3

Roll up the mat.

step. 4

Place the rolled mat into the bag.

Make sure to roll the mat and place it into the bag to carry or store it.

Caution

There might be a danger of injuries or fire if you do not check out the following. Please make sure that you check them out.

Check before use

1. Is the power cord connected to an outlet? Yes ☒ No ☐
2. When you turn on the device, do the power, time and temperature display lamps turn on? Yes ☒ No ☐
3. Is the temperature increasing well in the product? Yes ☒ No ☐
4. Is the temperature in the product not too high? Yes ☐ No ☒
5. Did you check out conditions for any damage to the power cord or cable conditions? Yes ☒ No ☐

Check after use

1. Is the power supply automatically turned off after 9 hours? Yes ☒ No ☐
2. Is the power, time and temperature display lamps turn off? Yes ☒ No ☐
3. Is the power cord removed from the outlet? Yes ☒ No ☐
4. Is the place to keep the product exposed to the humidity dust and/or chemicals? Yes ☒ No ☐

Caution

Ensure checking out the followings before applying for A/S if any abnormality occurs in the product.
The product might fail to operate normally even without any failure if you use it without reading the manual.

Checklist for any abnormality in the power supply

If there is no power supply on the power lamp;

- ① Check out if the power plug on the controller is connected to an outlet.
- ② Check out if the power button in the controller is pressed to turn on the product.
- ③ Check the operating voltage.

Checklist for any abnormality in the heating

If the temperature is not increasing in the product;

- ① Check out if the power button in the controller is pressed to turn on the product.
- ② Check out if the temperature is set to too weak.
- ③ Do not disassemble the product but contact your supplier or our company.



Check Indicator Window of Controller

Caution

Do not disassemble the product arbitrarily if any error occurs, but remove the power plug, turn off the power switch, and then contact your supplier.

If you disassemble the product arbitrarily, it might cause serious injuries and/or damage to the product.

Display	Descriptions	Display locatio
 	Disconnection of temperature sensor	<p>Temperature Indicator Window</p> <p>Temperature display lamp (LED) blinking Buzzer sound will play when it turns on</p>
 	A short circuit in the temperature sensor Overheated (temperature exceeding 85 °C)	<p>Temperature Indicator Window</p> <p>Temperature display lamp (LED) blinking Buzzer sound will play when it turns on</p>

Product Specifications

Product name and model name	Thermal massager T5
Manufacturer	NUGA MEDICAL CO.,LTD. 185, Jiraeul-ro, Jijeong-myeon, Wonju-si, Gangwon-do, Korea
Date of manufacture / number of manufacture	Marked on the product
Product weight	2.2 kg
Product size	727 X 390 X 50 (mm)
Cautions during use	Refer to the user manual
Operating method	Refer to the user manual
Electric rating	AC 220-240 V~, 50/60 Hz
Power consumption	90 W 

Company	Product name	
	Model name	
	Warranty period	__ years from purchase
	Date of purchase	
	Place of purchase	
	Phone	
Customer	Name of customer	
	Customer number	
	Address	

- When damage is caused by repairs or modifications done by a person who is not a certified repair personnel of NUGA Medical's local office or service center.

<http://www.nuga.kr>