Diagnostic Guide Wires

INSTRUCTIONS FOR USE

DESCRIPTION

Merit Medical guide wires are fabricated from high quality stainless steel utilizing a sophisticated construction process and are available with or without a PTFE coating. Guide wires are supplied sterile, non-pyrogenic, and are intended for single use only.

Merit Medical guide wires are packaged in a plastic hoop, which is fitted with a luer hub. This packaging is provided to facilitate compliance with the manufacturer recommended guidelines that the wire be flushed with saline or heparinized saline prior to use (See directions for use - Note).

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CONTRAINDICATIONS

InQwire diagnostic guide wires, are contraindicated for use in the coronary and cerebral vasculature.

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PRECAUTIONS

Angiography should be undertaken only by an experienced angiographer.


FOR ONE TIME USE ONLY

Guide wires will collect blood and other foreign material in their lumens; neither autoclaving nor ultrasonic cleaning will completely remove foreign material, therefore guide wires are recommended for one-time use.

Inspect all guide wires prior to use. Do not use any unit if the package is open or damaged.

Employ an aseptic technique during removal from the package and during use.

All guide wires are secured in the hoop dispenser by the locking J-tip straightener. To avoid damaging the guide wire during removal from the flush hoop, grasp the J-tip straightener near the base and slide it forward approximately 5mm or until the J-tip straightener is no longer attached to the flush hoop adaptor. Holding both the guide wire and J-tip straightener, continue to dispense guide wire from the hoop.

Do not use excessive force to advance the moveable core while the guide wire is in a vessel. Advancement with excessive force may cause coil penetration and vessel damage. Never push, auger, or withdraw a guidewire which meets resistance as this could potentially affect other indwelling devices.

Avoid withdrawing PTFE coated guide wires back through a metal needle. The sharp edge of the needle may scrape the coating. It is suggested that a catheter or PTFE vessel dilator replace the needle as soon as the guide wire has reached the appropriate position.

During advancement of the catheter and guide wire within the aorta, it is recommended that the guide wire be removed at the appropriate level of the aorta.


Care should be taken when manipulating a catheter during placement and withdrawal to prevent possible intravascular tissue damage. If resistance is felt during advancement, manipulation or removal from the catheter, stop immediately and confirm the guide wire tip position under fluoroscopy. Note proximity of other potential indwelling devices within the patient's anatomy. Never push, auger, or withdraw a guidewire which meets resistance as this could potentially affect other indwelling devices. Resistance may be felt tactilely or noted by tip buckling during fluoroscopy.

A guide wire is a delicate instrument and remains the most fallible instrument used in a percutaneous procedure. Any time that a guide wire is used there is a possibility of thrombus formation/emboli, vessel wall damage, and plaque dislodgement, which could result in myocardial infarction, cardiac arrhythmia, stroke or death. The physician should be familiar with the use of angiography products and the literature concerning the complications of angiography.

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
FLUSH HOOP INSTRUCTION GUIDE

1. Attach flush filled syringe to flush port luer
2. Rotate syringe clockwise (as pictured)
3. Inject saline into hoop
4. Detach syringe from hoop

5. Dispense guide wire

Note: In order to reduce the potential of clot formation, it is recommended that the guide wire be flushed with saline or heparinized saline prior to use. Attach a filled syringe to the luer hub located at the end of the plastic hoop, flush several times. After flushing, remove guide wire from hoop and use as described above.


DIRECTIONS FOR USE

NOTE: Confirm guide wire and needle compatibility prior to use.


The following schematic shows a typical procedure for percutaneous entry utilizing the Seldinger technique. Variations in individual patient anatomy may preclude the utilization of this technique.

1. Vessel puncture with a two-part needle


VESSEL


2. When using a two-part needle remove stylet leaving cannula in place, insert flexible (distal) end of guide wire through cannula and into vessel


3. "J" guide wires are shipped with "J" straightener to aid in the insertion of the wire into the puncture needle. Advance the straightener until 2-3 mm of the tip extends from the tip. Insert wire into hub and through needle. Remove "J" straightener proximally and discard.


4. Remove needle cannula leaving the guide wire within the lumen of the vessel.

5. Pass dilator or catheter over the guide wire directly into the vessel.

6. Carefully remove the guide wire leaving the catheter in place.

NOTE: Device is sterilized as stated on the package label.

WARNING: After use, dispose of device in a manner consistent with standard protocols for biohazard waste disposal.

REUSE PRECAUTION STATEMENT: For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

NOTE: Sterile if package is unopened and undamaged.

NOTE: Device is non-pyrogenic.

HOW SUPPLIED: Individually packaged 5-10 per box, refer to catalog for ordering information.

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SYMBOL	DESIGNATION
	Single Use
	Do Not Re-sterilize
	Caution
	Sterilized Using Ethylene Oxide
	Do Not Use If Package is Damaged or Opened and Consult Instruction for Use
	Consult Instructions for Use For electronic copy scan QR code, or go to www.merit.com/ifu and enter IFU ID. For printed copy, call U.S.A. or EU Customer Service
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Non-Pyrogenic
	Medical Device
	Authorized Representative in European Community
	Manufacturer
	Catalog Number
	Date of Manufacture
	Lot Number
	Single Sterile Barrier System
	Use-By Date
	Unique Device Identifier

CE 2797

www.merit.com

Manufacturer:
Merit Medical Systems, Inc.
1600 West Merit Parkway,
South Jordan, Utah 84095 U.S.A.
1-801-253-1600
U.S.A Customer Service 1-800-356-3748

Authorized Representative:
Merit Medical Ireland Ltd,
Parkmore Business Park West,
Galway, Ireland
EU Customer Service +31 43 358 82 22

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INSTRUCTIONS FOR USE

INDICATION FOR USE:

The EN Snare® Endovascular Snare System is intended for use in the coronary and peripheral vascular system or hollow viscous to retrieve and manipulate foreign objects. Retrieval and manipulation procedures include, indwelling venous catheter repositioning, indwelling venous catheter fibrin sheath stripping, and central venous access venipuncture procedure assistance.

DESCRIPTION:

The EN Snare® system consists of three interlaced, cabled, super-elastic Nitinol, preformed loops. The super-elastic Nitinol construction enables the loops to be introduced through catheters without the risk of device deformation.

WARNINGS:

1. This device is not intended for the removal of foreign objects entrapped by tissue growth.
2. This device should not be used for fibrin sheath stripping in the presence of septal defects of Persistent Foramen Ovale.
3. This device is not intended for removal of implanted pacing leads.
4. Pull forces applied to catheters during fibrin sheath stripping may damage, stretch, or break indwelling catheters 6 French or smaller in diameter. Do not use excessive pull force when attempting fibrin sheath stripping of catheters 6 French or smaller in diameter.
5. Do not use excessive force when manipulating the catheter through an introducer. Excessive force may damage the catheter.
6. This device has been sterilized utilizing Ethylene Oxide and is considered sterile if the package is not opened or damaged. It is intended for Single Patient Use Only. Do not attempt to clean or re-sterilize the device. After use this device may be a potential biohazard. Handle in a manner that will prevent accidental contamination. Do not use a device that has been damaged or if the package is open or damaged.
7. For single patient use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
8. Nitinol is a nickel titanium alloy. Possible reaction may occur for those patients who exhibit sensitivity to nickel.

Rx ONLY

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician

Do not use if package is damaged and consult instructions for use.

Non-pyrogenic

PRECAUTION: Care should be observed when using this device for removal of a large fibrin sheath in order to minimize risk of pulmonary embolism.

POTENTIAL COMPLICATIONS:

1. Potential complications associated with foreign body retrieval devices in arterial vasculature include, but are not limited to:
 - Embolization
 - Stroke
 - Myocardial infarction (depending upon placement)
2. Potential complications associated with snare retrieval devices in venous vasculature include, but are not limited to:
 - Pulmonary embolism
3. Other potential complications associated with foreign body retrieval devices include, but are not limited to:
 - Vessel perforation
 - Device entrapment

Catheter damage can occur when attempting fibrin sheath stripping on small French size diameter catheters. (See WARNINGS) Incidence of pulmonary embolism after fibrin sheath stripping may occur. (See PRECAUTION).

Prepare the EN Snare® System:

Select the appropriate Snare diameter range for the site in which the foreign body is located. The Snare diameter range should approximate the size of the vessel in which it will be used.

1. Remove the Snare and Snare Catheter from their hoop holders and inspect for any damage.
2. Remove the Insertion Tool and Torque Device from the proximal end of the snare shaft.
3. Load the Snare into the Snare Catheter by inserting the proximal end of the snare into the distal (non-hubbed) end of the Snare Catheter, until the proximal end of the Snare shaft exits the hub and the loops can be retracted into the distal end of the Snare Catheter.
4. Test and inspect the device by extending and retracting the snare loops through the distal end of the snare catheter 2-3 times, while carefully examining the Snare Catheter and the device for any damage or defects.
5. When appropriate, the system (Snare and Snare Catheter) can be advanced to the desired site as a single unit assembled as described above.

Alternative Preparation of the EN Snare®:

If the Snare Catheter is already positioned within the vasculature, the provided Insertion Tool (located on the proximal end of the Snare and just distal to the Torque Device) may be used to

position the Snare in the indwelling Snare Catheter.

1. Remove the Snare from the protective holder and inspect for any damage.
2. Move the provided Insertion Tool distally until the loops of the Snare are enclosed within the tubing portion of the Insertion Tool.
3. Insert the distal end of the Insertion Tool into the hub of the indwelling Snare Catheter until resistance is felt. This will indicate the tip of the Insertion Tool is properly aligned with the inner lumen.
4. Hold the Insertion Tool as straight as possible, grasp the shaft of the Snare just proximal to the hub of the Insertion Tool and advance the Snare until it is secure within the lumen of the Snare Catheter. The insertion tool can be removed by grasping the blue tab and firmly peeling it away from the snare shaft.

Snare Assisted Retrieval and Manipulation Suggestions:

1. If present, remove the indwelling delivery catheter.
2. If a guidewire is in a patient at the location of a foreign body, advance a snare catheter over the guidewire to the desired location. Then remove the guidewire and advance the snare through the snare catheter. Alternatively, cinch one loop of the snare over the proximal end of the guidewire and advance the entire system (snare and snare catheter assembly) into a guide catheter or introducer sheath until the distal end of the snare catheter is positioned proximal to the foreign body.
3. If a guidewire is not present, pull the snare into the distal end of the snare catheter and advance through a guide catheter or introducer sheath until it is positioned proximal to the foreign body. Alternatively, collapse the snare loop by pulling the device into the distal end of the snare insertion tool. Place the tapered end of the snare insertion tool into the proximal (hub) end of the snare catheter, guide catheter or sheath and advance the snare forward maintaining constant contact between the insertion tool and snare catheter hub. NOTE: When attempting to utilize guide catheters or sheaths not specifically manufactured for use with the EN Snare® system, it is important to test product compatibility prior to use.
4. Gently push the Snare shaft forward to completely open the loops. The loops are then slowly advanced forward, and may be rotated if desired, around the proximal end of the foreign body. Alternatively, the Snare may be advanced beyond the target location and the loops brought back around the distal end of the foreign body.
5. By advancing the Snare Catheter, the loops of the device are closed to capture the foreign body. (Note that attempting to close the loops by pulling the Snare into the Snare Catheter will move the loops from their position around the foreign body.)
6. To manipulate a foreign body, maintain tension on the Snare Catheter to retain the hold on the foreign body, and move the Snare and Snare Catheter together to manipulate a foreign body to the desired position.
7. To retrieve a foreign body, maintain tension on the Snare Catheter and move the Snare and Snare Catheter assembly together proximally to, or into a guide catheter or sheath. The foreign body is then withdrawn through or together with the guiding catheter or introducer sheath. Withdrawal of large foreign bodies may require the insertion of larger sheaths, guiding catheters, or a cut-down at the peripheral site.

Snare Assisted Removal of Fibrin Sheaths from Indwelling Catheters:

1. Using standard technique, prepare a femoral vein approach, advance the selected Snare to the inferior vena cava or right atrium.
2. Advance a .035" guidewire through the end port (distal or venous port if more than one lumen) of the indwelling catheter and into the inferior vena cava or right atrium.
3. Position one of the Snare loops around the guidewire.
4. Advance the Snare over the distal end of the catheter to a position proximal to the fibrin sheath.
5. Close the Snare around the catheter and continue applying light traction while gently pulling the Snare down toward the distal end of the catheter over the end ports.
6. Repeat steps 4 & 5 until the catheter is free of fibrin sheath.

Snare Assisted Venous Canalization:

1. Introduce the Snare at a patent venous access site and position in the vasculature at the desired site.
2. Open the Snare loops to provide a target to guide an entry needle into the desired venous access site.
3. Introduce a guidewire through the needle and through the Snare loops. Remove the needle.
4. Close the Snare over the guidewire by advancing the Snare Catheter.
5. Pull the guidewire into the desired location.

EN Snare® FOREIGN BODY RETRIEVAL DEVICE SYSTEM*

Description	Snare Diameter Range	Snare Length	Catheter Size	Catheter Length
Mini EN Snare System	2-4 MM	175 CM	3.2 F	150 CM
Mini EN Snare System	4-8 MM	175 CM	3.2 F	150 CM

Description	Snare Diameter Range	Snare Length	Catheter Size	Catheter Length
EN Snare System	6-10 MM	120 CM	6 F	100 CM
EN Snare System	9-15 MM	120 CM	6 F	100 CM
EN Snare System	12-20 MM	120 CM	6 F	100 CM
EN Snare System	18-30 MM	120 CM	7 F	100 CM
EN Snare System	27-45 MM	120 CM	7 F	100 CM

* Every EN Snare System contains: (1) Snare, (1) Snare Catheter, (1) Insertion Tool and (1) Torque Device.

Manufacturer:

Merit Medical Systems, Inc.
1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A Customer Service 1-800-356-3748

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INSTRUCTIONS FOR USE
MODE D'EMPLOI
ISTRUZIONI PER L'USO
GEBRAUCHSANWEISUNG
INSTRUCCIONES DE USO
INSTRUÇÕES DE UTILIZAÇÃO
INSTRUÇÕES DE USO
GEBRUIKSAANWIJZING
BRUKSANVISNING
BRUKSANVISNING
BRUGSANVISNING
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
KULLANMA TALIMATLARI
INSTRUKCJA UŻYTKOWANIA
NÁVOD K POUŽITÍ
ИНСТРУКЦИИ ЗА УПОТРЕБА
HASZNÁLATI UTASÍTÁS
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
LIETOŠANAS PAM ĀCĪBA
KOMPRESIJOS PRIETAISAS
INSTRUCȚIUNI DE UTILIZARE
NÁVOD NA POUŽITIE
UPUTE ZA UPORABU
KASUTUSJUHEND
NAVODILA ZA UPORABO

Guide Wires

INSTRUCTIONS FOR USE

DESCRIPTION:

Merit Medical AmplatZ guide wires are fabricated from high quality stainless steel with PTFE coating. AmplatZ wires have increased stiffness which provides added strength and stability. Guide wires are packaged in a plastic hoop fitted with a luer hub. This packaging facilitates compliance with the manufacturer recommended guidelines that the guide wire be flushed with saline or heparinized saline prior to use. Guide wires are supplied sterile and non-pyrogenic.

INDICATIONS FOR USE:

Merit Medical guide wires are used to facilitate the placement of devices during diagnostic and interventional procedures.

CONTRAINDICATIONS:

InQwire diagnostic guide wires are contraindicated for use in the coronary and cerebral vasculature.

POTENTIAL COMPLICATIONS:

Potential complications which may result from the use of the device include but are not limited to: Air Embolism/Thromboembolism, Allergic Reaction, Cardiac Arrhythmia, Amputation, Arteriovenous (AV) Fistula, Breathing Difficulty, Death, Embolism, Hematoma, Hemorrhage, Hemoglobinuria, Infection or Sepsis/Infection, Myocardial Ischemia and/or Infarction, Pseudoaneurysm, Stroke (CVA)/Transient Ischemic Attacks (TIA), Thrombus, Vessel Occlusion, Vessel Perforation, Vessel Dissection, Vessel Trauma or Damage, Vessel Spasm, Wire Entrapment/Entanglement, Foreign body/Wire Fracture. Some of the stated potential adverse events may require additional surgical intervention.

INSPECTION PRIOR TO USE:

Product is sterile if package is unopened and undamaged. Prior to use, carefully examine all guide wires to verify that the sterile package or product has not been damaged in shipment. Prior to use and when possible during the procedure, inspect the guidewire carefully for coil separation, bends or kinks which may have occurred. Do not use a wire which has a damaged tip.

PREPARATION FOR USE:

Caution: A guide wire is a delicate instrument. Any time that a guide wire is used there is a possibility of thrombus formation/ emboli, vessel wall damage, and plaque dislodgement, which could result in adverse procedural complications and /or adverse patient outcomes. The physician should be familiar with the use of angiography products and the literature concerning the complications of angiography. Angiography should be undertaken only by an experienced angiographer.

Note: Distal tip of wire may be positioned inside the flush hoop to protect the fragile tip.

Note: Confirm guide wire and needle compatibility prior to use.

Caution: To avoid damaging the guide wire tip during removal from the flush hoop, slide proximal portion of guide wire body forward in the flush hoop allowing the distal wire tip to exit the flush hoop. Gently grasp guide wire tip and J straightener together as a unit and gently pull forward to withdraw the distal wire portion from the flush hoop.

1. Attach flush filled syringe to flush hoop luer.
2. Inspect and prepare the catheter or device to be used according to the manufacturer's instructions. This includes flushing the catheter to be used with saline solution.
3. Inject saline into hoop to completely fill loop until dripping out opposite end.

Note: In order to reduce the potential of clot formation, it is recommended that the guide wire be flushed with saline or heparinized saline prior to use.

4. Detach syringe from flush hoop luer.
5. Dispense guide wire into the port of the catheter.

Note: Employ an aseptic technique during removal from the package and during use.

Caution: Advancement with excessive force may cause coil penetration and vessel damage.

INSTRUCTIONS FOR USE:

Warning: Use extreme caution when withdrawing PTFE coated guide wires back through a metal needle. The sharp edge of the needle may scrape the coating. It is suggested that a catheter or PTFE vessel dilator replace the access needle as soon as the guide wire has reached the appropriate position.

1. Insert the guide wire J-straightener into the guide wire port of the intended catheter or device.
2. Carefully advance the distal guide wire tip through the J straightener and device lumen. Remove the J straightener by withdrawing it over the guidewire.

Warning: When reintroducing a guide wire into a catheter or device within a vessel, confirm that the catheter tip is free within the lumen (i.e. not against the vessel wall). Always advance or withdraw a wire slowly. Free movement of the guidewire within a catheter provides valuable tactile information. Never push, auger, or withdraw a guidewire which meets resistance as this could potentially affect other indwelling devices. Resistance may be felt tactfully or noted by tip buckling during fluoroscopy. Test all systems for resistance prior to use.

The guide wire and catheter should be removed as a unit when possible to prevent potential damage to the vessel wall.

3. Confirm guide wire tip placement under fluoroscopy to assure that the distal tip is intraluminal and in the intended vessel.
4. Hold the guide wire in position while manipulating the catheter over the guide wire to prevent unintended movement of the distal wire tip.

Warning: When reintroducing a guide wire into a catheter or device within a vessel, confirm that the catheter tip is free within the lumen (i.e. not against the vessel wall).

Warning: Always advance or withdraw a wire slowly. Free movement of the guidewire within a catheter provides valuable tactile information. Never push, auger, or withdraw a guidewire which meets resistance as this could potentially affect other indwelling devices. Resistance may be felt tactfully or noted by tip buckling during fluoroscopy. Test all systems for resistance prior to use.

REUSE PRECAUTION STATEMENT:

For Single Patient Use Only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DISPOSAL: After use, dispose of product and packaging in accordance with hospital protocol.

In the EU, any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the applicable Member State.

For a copy of this device's current European Summary of Safety and Clinical Performance (SSCP), please go to the European database on medical devices (Eudamed), where it is linked to the basic UDI-DI. <https://ec.europa.eu/tools/eudamed>.

Basic UDI-DI: 088445048761E4

SYMBOL	DESIGNATION
	Single Use
	Caution
	Sterilized Using Ethylene Oxide
	Do Not Use If Package is Damaged or Opened and Consult Instruction for Use
	Consult Instructions for Use For electronic copy scan QR code, or go to www.merit.com/ifu and enter IFU ID. For printed copy, call U.S.A. or EU Customer Service
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Non-Pyrogenic
	Medical Device
	Authorized Representative in European Community
	Manufacturer
	Catalog Number
	Date of Manufacture
	Lot Number
	Single Sterile Barrier System
	Use-By Date
	Unique Device Identifier

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