
CONTENT

Inhalt / Contenu / Inhoud / содержимое /
Contenuto / Contenido / İçerik / Conteúdo /
Treść / Зміст / / ЁЗЛW

Deutsch	4
English	11
Français	17
Nederlands	24
Русский	30
Italiano	37
Español	44
Türkçe	51
Português	58
Polski	65
Українська	73
.....	82
89	с/ОК

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Product Description

Visalys® Temp, developed by Kettenbach, is a highly esthetic, tooth-colored composite for fabricating temporary dental restorations that is also suitable for the anterior area. The two-component system is based on a multifunctional acrylic composite and is free of Bisphenol A and its derivatives. It has properties similar to filling materials and is available in the following shades: A1, A2, A3, A3.5, B1, BL. Temporary restorations fabricated with Visalys® Temp can also be supplemented and customized with low-viscosity, light-curing composites.

Indications / Techniques

- Fabrication of temporary crowns, partial crowns, bridges, inlays, onlays, and veneers
- Customization of preformed acrylic and metal crowns
- Fabrication of long-term temporary restorations

Warnings

Do not use if there is a known allergy to methacrylic acid esters. Allergic reactions to the product cannot be ruled out in sensitive persons. Sensitization as a result of skin contact is possible.

Visalys® Temp temporary restorations are intended to be used for temporary restorations and are not to be used as permanent restorations. Residual methacrylate-based material may change the setting behavior of silicone impression taking and bite registration materials. If a restoration has been fabricated beforehand, the preparation and surrounding tissue must be carefully cleaned with gauze soaked in alcohol and water spray before taking an impression. Eugenol-based temporary luting cements may inhibit the curing of the permanent composite cements that may subsequently be used for the final cementation.

Cautions

Do not use after expiration date. Do not leave any residue of the temporary crown and bridge material in the sulcus/on the tooth or oral cavity. Do not swallow material! If swallowed: Consult a medical doctor if problems arise or persist. Avoid contact with the eyes. If there is accidental contact with the eyes: Rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist. Do not inhale grinding dust. Use safety glasses and face mask! When taking impressions of areas with severe undercuts and wide interdental spaces, standard blocking-out techniques must be used.

Avoid contact with clothing, as the material cannot be removed by dry cleaning.

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Anatomical Impressions

To fabricate a temporary crown, bridge, inlay, or onlay, all standard impression materials can be used.

An A-silicone such as Silginat® is particularly well suited for anatomic impression taking because it remains dimensionally stable for several weeks and can be stored. When remaking a new temporary restoration, the impression taken with Silginat® can be reused.

Instructions for Use

Preparation

Take an anatomic impression of the dental arch beforehand. A laboratory fabricated, vacuum-formed splint or, in the case of single crowns, a prefabricated metal or plastic crown can be used instead of an impression for forming the shape. The impression can be cut out before use to eliminate any undercuts and to increase the stability of the temporary restoration. For bar connections of spaced abutment teeth, a channel can be cut into the impression. In the anterior area, denture teeth can be used as placeholders before taking the impression.

Before applying Visalys® Temp, prepared tooth surfaces of compomers, resin-modified glass ionomers or composites must be isolated, using petroleum jelly for example.

For optimal results we recommend wearing gloves made of vinyl or latex.

Dispensing and Mixing

The pastes are mixed and dispensed using Applyfix® 6 dispensing guns DS-50 4:1/10:1 and the corresponding mixing tips. Only use the enclosed Kettenbach mixing tips included.

The material requires a certain amount of time to flow out after the mixing tip has been attached. The material flow stops as soon as pressure is reduced on the dispenser handle. Any material already cured in the mixing tip must not be extruded with force because otherwise the cartridge and the mixing tip may be damaged.

Dispensing and Applying

Application

1. Select the desired shade, insert the cartridge in the dispensing gun, advance the plunger into the cartridge, remove the cap, and attach a new mixing tip (turn 90° clockwise until it clicks into place; Figure 1, page 95).

Before attaching a mixing tip onto a new cartridge for the first time, discharge a small amount of material until both materials extrude uniformly (the original cartridge cap can be replaced after use).

2. Before each use, extrude and discard a little material (Figure 2, page 95).

This point is defined as the start of mixing so start the stopwatch at the same time.

3. Fill the impression, vacuum-formed splint, or prefabricated crown with Visalys® Temp from the base upwards (Figure 3, page 95).

Keep the end of the mixing tip immersed in the material to avoid trapping any air bubbles.

The curing of the material is influenced by a number of different factors (e.g., mouth temperature, impression process, storage temperature, anatomic impression).

4. Insert and cement the filled mold, applying slight pressure, within 40 sec. of starting the mixing into the patient's mouth (Figure 4, page 95).

5. Remove the temporary restoration from the patient's mouth during the elastic phase between 1:30 min. and 2:30 min. after the start of mixing (Figure 5, page 95).

Visalys® Temp reaches an optimal elastic phase for easy removal on average 2 min. after the start of mixing. Use a probe to check the degree of hardness of excess material intraorally to determine whether the material can be removed from the mouth.

If this is not possible, we recommend applying a pea-sized quantity to the back of the hand directly after filling the impression. When the material has cured, the impression can be removed with the temporary restoration.

6. Finish the Visalys® Temp temporary restoration about 4 min. after the start of mixing (Figure 6, page 95).

After use, replace the sealing cap (alternatively, the used mixing tip can be used to seal the cartridge after carrying out hygiene measures). Before using the cartridge again, remove the cartridge cap (or the filled mixing tip), check the outlet of the cartridge to ensure that no cured material is blocking the openings (clear any blocked openings and discharge a small quantity of material until both materials extrude uniformly). Attach new mixing tips and proceed as described above.

Finishing

Completely remove the inhibition layer after the material has fully cured using alcohol, for example.

Once the material has fully cured (4 min. after the start of mixing), the temporary restoration can be finished using cross-cut burs and polishers, for example.

Cleaning / Disinfection

If necessary, Visalys® Temp can be disinfected with standard aqueous disinfectant.

Customization / Recontouring

Using a bonding agent is not necessary when using a filling composite in combination with Visalys® Temp for customization or recontouring. For other products, please follow the instructions for use.

Cementation

Remove any residual isolating material from the prepared tooth before temporary cementation.

Standard temporary luting materials can be used for cementation.

Residual eugenol-based temporary luting cement may negatively affect the curing of composite cements for the subsequent final cementation. If the use of composite cements

is planned, it is recommended to use a eugenol-free temporary luting cement.

After Cementation

Carefully check the sulcus and the surrounding areas and if necessary remove any residual temporary material or cement.

Repair

Roughen the fracture site and create mechanical retention on the adjacent areas.

Clean using Orthoskavident® C or alcohol and then air dry.

Join or add to newly fabricated restorations with Visalys® Temp or a flowable, light-curing filling composite.

Join or add to existing restorations with Visalys® Temp or a combination of a flowable, light-curing filling composite and a bonding agent.

Cure Visalys® Temp or the composite in accordance with the relevant instructions for use.

Finish as usual once the material has cured.

Cleaning and Disinfection of the Applyfix® 6

Remove the cartridge.

Remove any uncured paste using a cloth soaked in alcohol.

The Applyfix® 6 dispensing gun can be used several times if properly disinfected. Kettenbach recommends exchanging the gun if there are visible traces of use. Disinfection can be carried out using a commercially available dipping disinfectant. Kettenbach recommends the use of glutaraldehyde-based solutions.

Storage Conditions

Store in a dry place, protected against sunlight and at room temperature (18°C/64°F–25°C/77°F). Avoid exposure to extreme temperature fluctuations.

Do not store in a refrigerator.

Disposal

The empty cartridges/mixing tips and cured material can be disposed in the household rubbish.

Trademarks

Applyfix® 6, Visalys® Temp, Orthoskavident® C and Silginat® are registered trademarks of Kettenbach GmbH & Co. KG.

Note

Selected Kettenbach materials are available in certain markets only.

For professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

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Limitation of Liability

Except where prohibited by law, Kettenbach GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

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
Technical data

Working and setting times

00:00–00:40 min.	Pouring Visalys® Temp into the impression and placement in the mouth
	Curing in the mouth
01:30–02:30 min.	Removal from the mouth
	Full curing
From 04:00 min.	Finishing the temporary restoration

Note

Working times apply for products, which have been stored and processed at a temperature of 23°C / 74°F and a relative humidity of 50 %. Higher temperatures shorten the working time and lower temperatures extend the working time.

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VISALYS® TEMP

1


Select shade, insert cartridge, remove cap, attach new mixing tip.

2

Before each use, extrude and discard a little material, start a stopwatch.

3

Fill from the base upwards, keep the end of the mixing tip in the material. Tip: apply a pea-sized quantity of material to the back of the hand (see step 5).

4

Insert into the patient's mouth, applying slight pressure, within 40 sec. of starting the mixing.

5


Check the degree of hardness intraorally (or use the sample on your hand, see step 3). Remove the temporary restoration between 1:30 min. and 2:30 min. after the start of mixing.

6

Finish the Visalys® Temp about 4 min. after the start of mixing.

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Product Description

Visalys® CemCore is a dual-curing composite containing fluoride for adhesive cementation of indirect restorations and root posts. Moreover, Visalys® CemCore can be used to fabricate radiopaque core build-ups and build-up fillings.

Due to the Active-Connect-Technology (ACT) Visalys® CemCore includes an additional initiator system for the optimization of the adhesive bond.

Visalys® CemCore in the self-mixing 1:1 automix syringe is free of bisphenol A and its derivatives. BPA or BPA precursors are not used in the manufacturing process of this device.

Visalys® CemCore is easy to use and has high adhesive strength and good flowability with stable consistency.

Visalys® CemCore shall only be used in combination with self-etching Visalys® Tooth Primer. Visalys® Restorative Primer is recommended as the primer for the restoration surfaces.

Visalys® CemCore is offered in the shades Universal (A2 / A3), Dark (A4), Opaque, Bleach and Translucent.

Indications for Use

Permanent cementation of:

- crowns, bridges, inlays, onlays, veneers and adhesive bridges made of metal, ceramic, composite, hybrid and oxide ceramics
- root posts and indirect core build-ups

Composite for:

- adhesive core build-ups

Patient Target Group

Patients being treated in the context of a dental procedure.

Contraindications

Known allergies to methacrylic acid ester, benzotriazole compounds, benzoyl peroxide, camphorquinone, 4-methoxyphenol, or other constituents of Visalys® CemCore or Visalys® Tooth Primer (see Visalys® Tooth Primer instructions for use).

If there are no adequate means of ensuring drainage.

Warnings

Allergic reactions to the product cannot be ruled out in sensitive persons. Sensitization through contact with skin / mucous tissue is possible.

Materials containing eugenol can inhibit setting of composites. The use of hydrogen peroxide solutions can impair the adhesive strength. Retraction paste containing iron(III)sulfate can impair curing and adhesion of the material.

In the case of vital teeth, areas near the pulp may have to be treated with suitable pulp / dentin protection.

Keep out of reach for children.

Cautions

Do not use after the expiration date.

Do not swallow material. If the material is inadvertently swallowed: Consult a medical doctor if problems arise or persist.

Take the appropriate protective measures (e.g. wear gloves or safety glasses). Avoid contact of skin and mucous tissue with non-cured material.

Avoid contact with the eyes. If there is accidental contact with the eyes: Rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.

Avoid contact with clothing, as permanent damage to clothes cannot be excluded.

Safety data sheets are available at www.kettenbach.com.

General Notes

Given refrigerator storage, the product has to have reached room temperature prior to use.

In case of contact with Visalys® Tooth Primer and direct incident radiation of ambient light, curing of Visalys® CemCore accelerates.

For layers thicker than the specified polymerization depth, it is necessary to wait until the self-curing process is complete in all cases.

Working times apply for products, which are processed at a temperature of 23°C / 74°F and a relative humidity of 50 %. Higher temperatures reduce the working time and lower temperatures extend these times.

Application Instructions for Cementation

After removing the temporary appliance, carefully remove any residues of temporary cement from the core / cavity surface.

Checking the Fitting Accuracy and Occlusion of the Restoration

Thin and brittle restorations can break under excessive pressure during try-in. For this reason, do not apply high pressure or only check the occlusion after permanent cementation.

With the available Visalys® CemCore Try In Pastes, the shade effect of the selected composite shade can be checked on the slightly moist core prior to final cementation. Apply the paste in the selected shade on the inside surfaces of the restoration and check the shade effect on the restoration. The shade of the Visalys® CemCore Try In Paste matches the shade of the cured cementation composite.

Remove the Visalys® CemCore Try In Paste from the restoration and core surface with water.

After the try-in, the surfaces of the restoration that may be contaminated with saliva, blood or silicone residues, for example, have to be cleaned with water spray or in the ultrasonic bath with distilled water and then dried with water/oil-free air. All the surfaces must be free of contamination.

In the case of restorations made of non-precious metal or oxide ceramics, there may be no cleaning with phosphoric acid.

A rubber dam is recommended for a dry operating field.

Pretreatment of the Restoration, the Root Post or the Indirect Core Build-Up

Unless prescribed otherwise in the instructions for use for the restoration material, root post or material used for the indirect core build-up, the following procedure should be selected for pretreatment:

- 1.1 Adhesive surfaces of noble metal, metal and oxide ceramics (zirconium oxide and aluminum oxide ceramics)
 - Sand blasting of the surfaces with aluminum oxide powder ($\leq 50 \mu\text{m}$). Observe the restoration material manufacturer's blasting parameters.
 - Then rinse off thoroughly with water spray or clean in the ultrasonic bath.
DO NOT CLEAN WITH PHOSPHORIC ACID!
 - Dry with water / oil-free air.
- 1.2 Adhesive surfaces of silicate ceramic (feldspar ceramic, glass ceramics) and hybrid ceramics
 - Etch the surfaces with hydrofluoric acid according to the specifications in the instructions for use of the restoration material.
 - Then rinse off thoroughly with water spray or clean in the ultrasonic bath.
 - Dry with water / oil-free air.
- 1.3 Adhesive surfaces of composite
 - Sand blasting of the surfaces with aluminum oxide powder ($\leq 50 \mu\text{m}$). Observe the restoration material manufacturer's blasting parameters.

- Then clean thoroughly with alcohol.
 - Dry with water / oil-free air.
2. Drop Visalys® Restorative Primer into a preparation tray. Light protection is not necessary. Close primer bottle carefully after each use.
 3. Apply Visalys® Restorative Primer with a disposable applicator or brush onto the adhesive surfaces and allow it to work in for 60 sec. Then blow out and dry with a gentle jet of water / oil-free air. Visalys® Restorative Primer is non-light curing.

The surfaces of some types of posts can be attacked by sandblasting or hydrofluoric acid. The specifications of the root post manufacturer on surface treatment must be observed.

Pretreatment of the Core / Cavity Surface (Enamel, Dentin, Metal, Composite)

Visalys® CemCore may only be used together with Visalys® Tooth Primer. Visalys® Tooth Primer is a self-etching single-component primer and can therefore be used without additional phosphoric acid etching. Only with cementation of veneers and adhesive bridges and with uncut enamel is selective phosphoric acid etching necessary. Depending on the clinical conditions, in other cases too, selective phosphoric

acid etching may be optionally performed. Use the etching gel according to the instructions for use from the relevant manufacturer.

1. Following the try-in of the restoration, saliva, blood and any residues of the Visalys® CemCore Try In Paste must be removed from the adhesive surfaces of the dental hard substance and restorations surfaces using water spray and then keep the surfaces dry. The adhesive surfaces must be free of contamination. A rubber dam is recommended for a dry operating field.
2. Drop Visalys® Tooth Primer into a preparation tray. Light protection is not necessary. After applying the Visalys® Tooth Primer, start with the further processing immediately. Close primer bottle carefully after each use.
3. Apply Visalys® Tooth Primer to the entire adhesive surface of the tooth surface using a suitable applicator and rub in for 20 sec. Any excess can be removed with a cotton pellet and blown off with air. Then dry with a gentle jet of water / oil-free air. The prepared primer should be processed within 5 min., as constituents can evaporate. The Visalys® Tooth Primer is non light-curing.

For cementation on an implant abutment or on an indirect core build-up, the specifications of the respective material manufacturer must always be considered. Pretreat the entire

adhesive surfaces of the abutment/core build-up dependent on the material, as described under “Pretreatment of the restoration” in the Visalys® CemCore instructions for use, unless prescribed otherwise in the relevant material manufacturer’s instructions for use.

The emergence profiles and screw channels, as applicable, must be protected by suitable measures (e.g. wax, silicone covers).

Application of the Cementation Composite

1. Apply Visalys® CemCore on the inside surfaces of the restorations and on the core or cavity surfaces, as applicable. In the case of application in the mouth, the processing time is reduced and the restorations must be placed within 2 min.
2. Place the restoration on the core or in the cavity.
3. Removing excess
- 3.1. **Method 1:** Remove after initial curing through brief exposure to light or self-curing until the “gel phase”

The excess is either polymerized with a commercially available polymerization lamp (1200 mW/cm²) for 2–3 sec. oral and vestibular as well as mesial and distal – i.e. at four points – or wait for the gel phase in self-cur-

ing. The excess, which is still soft-elastic, can be easily removed with a scaler.

3.2. **Method 2:** Removing without initial curing

Remove the excess directly, e.g. with a cotton pellet, brush or foam pellet.

While removing the excess, the restoration has to be held in its final position. All excess, especially in difficult-to-reach approximal areas, has to be removed prior to final curing.

4. Final curing

After removal of excess, the cement joint can be covered with a glycerin gel or polyethylene glycol paste, e.g. Visalys® CemCore Try In Paste, to avoid the formation of an inhibition layer. This can then be removed with water.

Then, polymerize each surface / cement joint for 10 seconds with light.

For non-translucent restorations, the final curing must be allowed to complete beneath the restoration (see the table for technical data).

Translucent restorations can also be cured in several layers according to the specified polymerization times dependent on the size of the light aperture.

5. Finishing

Rough cement joints can be finished and polished with the appropriate instruments.

Special Procedure for Cementation of Root Posts and Indirect Core Build-Ups

1. Prepare the root canal, root post or indirect core build-up according to the instructions for use from the relevant manufacturer. If the manufacturer of the root post or indirect core build-up recommends silanization or pre-treatment with a primer for metal or oxide surfaces or a primer-adhesive / universal adhesive-activator combination, the Visalys® Restorative Primer can be used for this purpose.
2. Whetten the cleaned and dried surfaces of the root canal, the coronal dentin and the enamel parts, as applicable, completely with Visalys® Tooth Primer and rub in for 20 sec. Remove excess from the root canal with a paper tip and dry all surfaces with a gentle jet of water/oil-free air.
3. Apply Visalys® CemCore into the canal lumina and also on the root post. The end of the tip used should always remain in the material while applying. Do not use Lentulo spirals.

4. Immediately afterwards, insert the root post (the working time for Visalys® CemCore in the root canal is approx. 2 min.). It is recommended that any excess material is specifically distributed on the occlusal core surface.
5. To fix the post, the material can be light-cured with a suitable polymerization lamp for 10 sec.
6. Perform the build-up of the core as described under "Application instructions for core build-up".

Application Instructions for Core Build-Up

1. Clean and dry the cavity.
2. Areas close to the pulp must be protected with the appropriate means.
3. Ensure a clean and dry operating field (rubber dam) during the entire application.
4. Selective phosphoric acid etching of the enamel may be undertaken optionally.
5. Drop Visalys® Tooth Primer into a preparation tray. Light protection of the primer is not necessary.
6. Apply Visalys® Tooth Primer to the entire adhesive surface of the core or the cavity using a suitable applicator and rub in for 20 sec. Any excess can be removed with a cotton pellet and blown off with air. Then dry with a jet of water / oil-free air.

The prepared primer should be processed within 5 min., as constituents can evaporate.

7. Place a matrix and apply Visalys® CemCore directly into the cavity, as applicable. The use of a transparent matrix makes it easier for light to enter in the case of optional light curing. Alternatively, in the case of very badly damaged teeth, a core build-up mold can be used. The mold is filled with Visalys® CemCore and then placed on the residual tooth substance. Visalys® CemCore can still be modeled within the processing time.
8. Wait for final curing (see Technical Data). Optionally, the material can be light-cured for 20 sec. to accelerate curing down to a depth of approx. 2.5 mm. Light-curing should be performed with a commercially available polymerization lamp (1200 mW / cm²). In case of a lower light intensity, prolong the light curing duration accordingly. Place the lamp as close to the material as possible.

Deeper layers that are not accessible to light will cure in a purely chemical manner after approx. 5 min.
9. For layers thicker than the specified polymerization depth, it is necessary to wait until the self-curing process is complete in all cases.
10. After full curing, remove if applicable any matrices or core build-up formers.
11. Performing core preparation.

12. Isolation of the core prior to fabrication of a temporary restoration.

In the case of self-curing, lower temperatures, e.g. of the product or during work in the anterior region, lead to prolongation of the curing time.

Application

The paste is mixed and delivered directly from 9 g / 5 ml-auto-mix syringe and the appropriate mixing tips. Only use original mixing tips (ML Ø 2.5 mm, blue) and Intraoral / Endo tips!

The following points have to be observed:

1. Remove the syringe cap and prior to first attaching a mixing tip on the syringe, expel a small amount until both components are evenly delivered. Then discard the paste expelled.
2. Attaching the mixing tip. Turn the cap 90° clockwise to lock in place. If applicable, mount an Intraoral / Endo tip locking on the mixer.
3. The material requires a certain amount of time to flow out after the mixing tip has been attached. The material flow stops as soon as pressure is reduced on the syringe piston. Any material already cured in the mixing tip must not be extruded with force, because otherwise the syringe and the mixing tip may be damaged.

4. Dispose of used mixing tip, replace the syringe cap and check the outlet at the start of the next application.

Cleaning and Disinfection

The risk of contamination of the automix syringe can be reduced by covering the syringe with a disposable protective sleeve.

Clean the syringe with a disposable wipe and standard media.

Disinfection of the syringe in accordance with national guidelines.

Observe standard practice hygiene measures.

Storage Conditions

Store dry, protected from sunlight. The product has to be kept in the refrigerator at 2°C/36°F–8°C/46°F before first use.

Shelf-Life After First Use

- a) For continued storage in the refrigerator:
Until the end of the specified expiration date
- b) For continued storage at room temperature:
Use within six months

Avoid extreme temperatures. Given refrigerator storage, the product has to have reached room temperature prior to use.

Do not use after the expiration date.

Constituents

Main components of Visalys® CemCore are approx. 42 vol. % inorganic fillers in the size range 0.2–20 µm including ytterbium fluoride. The polymer base consists of aliphatic dimethacrylates.

Classification

ISO 4049, type 2, class 3.

Disposal

The empty automix syringes/mixing tips/applicator tips and cured material can be disposed in the household rubbish.

Trademarks

Visalys® CemCore, Visalys® CemCore Try In Paste, Visalys® Tooth Primer and Visalys® Restorative Primer are registered trademarks of Kettenbach GmbH & Co. KG.

Note

Not all Kettenbach products are available in all countries.

For dental professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

If serious incidents occur that are related to the product, please immediately contact the manufacturer and the competent authority.

Limitation of Liability

Except where prohibited by law, Kettenbach GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the legal basis asserted, including warranty, contract, negligence, or strict liability.

The information provided for Kettenbach products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use.

Updated: May 15, 2023

Color stability meets the requirements of ISO 4049.

* in a light wavelength range 400–500 nm

** Aluminum has a radiopacity equivalent to dentin. Therefore, 1 mm of a material that has a radiopacity equivalent to 1 mm of aluminum has a radiopacity equivalent to 1 mm of dentin, and 2 mm of aluminum is equivalent to the tooth enamel.

Technical Data

Working time (23°C/74°F)	approx. 4 minutes
Working time (intraoral)	approx. 2 minutes
Exposure time for initial curing (tack curing)	2–3 seconds per polymerization point
Time for purely chemical curing including working time (intraoral)	approx. 5 minutes
Light curing* (cementation)	10 seconds per surface / cement joint
Light curing* (core build-up)	20 seconds
Radiopacity**	approx. 2.5 mm Al
Depth of cure Translucent	approx. 2.5 mm
Depth of cure Universal (A2 / A3)	approx. 2.0 mm
Depth of cure Bleach, Dark (A4)	approx. 1.5 mm
Depth of cure, Opaque	approx. 0.5 mm
Light intensity	1200 mW/cm ²

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Product Description

Visalys® Core, developed by Kettenbach, is a dual-curing, fluoride-containing composite used for the fabrication of radiopaque core build-ups and core fillings. The automixing, two-component system is based on a multifunctional acryl composite. BPA or BPA precursors are not used in the manufacturing process of this device. Visalys® Core has a firm consistency with good flow behavior and is available in the following shades: white and dentine. Curing can be controlled as required at any time using optional light-curing.

Indications / Techniques

- Different types of core build-ups
- Luting of root posts

Warnings

Do not use in case of a known allergy to methacrylic acid ester or benzoyl peroxide.

Allergic reactions to the product cannot be ruled out with sensitive persons.

Sensitization is possible on contact with the skin / mucosa.

Materials containing eugenol could inhibit the curing of composites.

If necessary, pulp areas of vital teeth should be treated using suitable pulp / dentine protection.

Keep out of reach of children.

Cautions

Do not use after expiration date.

Do not swallow the material. If swallowed: Consult a medical doctor if problems arise or persist.

Avoid skin and mucosal contact with uncured material.

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Avoid contact with articles of clothing, as the material cannot be removed by dry cleaning.

Safety data sheets can also be found at www.kettenbach.com.

Note

The working times apply for products, which have been stored and processed at a temperature of 23° / 74°F and 50 % relative humidity. Higher temperatures shorten and lower temperatures extend these times.

Direct exposure to ambient light accelerates the curing of Visalys® Core.

Recommended Use

Core Build-Up

1. Remove existing fillings and caries.
2. Provide for a dry treatment site using appropriate measures (rubber dam).
3. If the bonding agent requires an additional etching stage, complete etching according to the manufacturer's instructions.
4. Apply the bonding agent according to the manufacturer's instructions.
5. Place a matrix around the tooth if the cavity wall is missing and apply Visalys® Core directly in the cavity. When using a core build-up former, fill it directly using Visalys® Core and place on the tooth. Visalys® Core can be sculpted during the working time.

6. Visalys® Core cures automatically after 5 minutes. The material can optionally be light cured for 20 seconds to accelerate curing. The curing depth with light-curing is 2 mm. The minimum light intensity should be 800 mW/cm². In case of lower light intensity, the light-curing time must be extended accordingly. Place the light as close as possible to the material.

When using AdheSE or AdheSE One F (Ivoclar Vivadent) and Adper Prompt L-Pop (3M ESPE) adhesive, apply the core build-up material Visalys® Core in incremental layers (2 mm), and light-cure.

During curing, lower temperatures, e.g. of the product or working in the anterior region, extend the curing time.

7. Remove the matrix or core build-up former, if necessary.
8. Prepare the core after fully curing.
9. Apply a separating medium on the core surface, e.g. using Vaseline, before fabricating a composite temporary restoration.

Luting Root Posts

For the luting of root posts it is recommended to use Visalys® Core with the 9 g / 5 ml-automix syringe.

1. Prepare the filled root canal for placing the post.
2. If the bonding agent requires an additional etching stage, complete etching according to the manufacturer's instructions.
3. Prepare the root canal, coronal dentine and, if applicable, the enamel section according to the bonding agent manufacturer's instructions. It is essential to remove any excess from the lumen using a paper tip to prevent blockage of the canal.
4. Condition the surface of the post according to the material manufacturer's instructions. Then avoid contamination of the post surface.
5. Apply Visalys® Core in the canal lumen and, if necessary, on the root post. The tip used should always remain in the material during application. Do not use a lentulo.
6. Immediately afterwards, insert the root post (the working time for Visalys® Core in the root canal is 45 seconds). It is recommended that any excess material is specifically distributed on the occlusal stump surface.
7. To fix the post, the material can be light-cured with a suitable polymerization lamp for 20 seconds.
8. Build up the core as described above.

The coronal section of the tooth can be prepared as a core build-up using Visalys® Core. Use of a matrix is recommended during the build-up due to the flowable consistency.

Application

The pastes should be mixed and dispensed directly from the 9 g/5 ml-automix syringe or using an Applyfix® 8 dispensing gun DS 24 1:1/2:1 and for the 45g/25 ml-cartridge corresponding mixing tips. Only use the original mixing tips MBT Ø 4.2 mm or MLT Ø 2.5 mm and Intraoral/Endo tips!

The following should be observed when using the 45g/25 ml-cartridge (when using the automix syringe, please start at Point 2):

1. Insert the cartridge in the Applyfix® 8 dispensing gun. Ensure that the notches of the cartridge base face downwards. If the cartridge is fitted incorrectly in the dispensing gun, the locking clip will not close.
2. Remove the cartridge/syringe cap and, before fitting a mixing tip for the first time, extrude a small amount of material until both materials are delivered uniformly. Then discard the extruded paste.
3. Fit the mixing tip. Turn the cap 90° in a clockwise direction until it engages. If applicable, mount an Intra-oral/Endo tip locking on the mixer.

4. Apply constant pressure while extruding material through the mixing tip. The material flow is interrupted as soon as the pressure is eased on the dispenser grip or syringe plunger. Do not use force to extrude material that has already cured in the mixing tips, as the cartridge/syringe and mixing tips could be damaged by creating excess back pressure.
5. Dispose of used mixing tip, replace the cartridge/syringe cap and check the outlet at the start of the next application.

Cleaning and Disinfection of Applyfix® 8

Remove the cartridge.

Remove uncured paste using a cloth soaked in alcohol.

The Applyfix® 8 dispensing gun can be used several times if properly disinfected. Kettenbach recommends exchanging the gun if there are visible traces of use. Disinfection can be carried out using a commercially available dipping disinfectant. Kettenbach recommends the use of glutaraldehyde-based solutions.

Storage Conditions

Store dry and protect from sunlight, storage in a refrigerator or at room temperature (2°C/36°F–24°C/75°F), avoid extreme temperatures. The product must be heated to room

temperature before use when it has been stored in a refrigerator.

Classification

ISO 4049, Type 2, Class 3.

Disposal

The empty automix syringes/cartridges/mixing tips/appliator tips and cured material can be disposed in the household rubbish.

Trademarks

Applyfix® 8 and Visalys® Core are registered trademarks of Kettenbach GmbH & Co. KG.

Note

Selected Kettenbach materials/products are available in certain markets only.

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Limitation of Liability

Except where prohibited by law, Kettenbach GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the legal theory asserted including warranty, contract, negligence or strict liability.

The information provided for Kettenbach products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply and consider all recommendations and information in connection with any use.

Updated: November 2, 2020

Technical Data

Working and curing times

Working time	90 seconds
Working time in the root canal	45 seconds
Light-curing (in a light wavelength range 400-500 nm)	20 seconds
Self-curing	5 minutes
Radiopacity*	2.5 mm Al
Flexural strength	> 100 MPa

Visalys® Core contains 42 vol% inorganic fillers with a size range of 0.2–20 µm.

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*Aluminum has a radiopacity equivalent to dentin. Therefore, 1 mm of a material that has a radiopacity equivalent to 1 mm of aluminum has a radiopacity equivalent to 1 mm of dentin, and 2 mm of aluminum is equivalent to the tooth enamel.

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