
Surgiwand™ II Auto Suture™

Suction and Irrigation Device



1301209





COVIDIEN

15 HAMPSHIRE STREET
MANSFIELD, MA 02048
U.S.A

[508]261-8000 [T]

2

3

4

5

6

7

EN**BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.****IMPORTANT!**

This booklet is designed to assist in using this product. It is not a reference to surgical techniques.

This device was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury. Reprocessing and/or resterilization of this device may create the risk of contamination and patient infection. Do not reuse, reprocess or resterilize this device.

Where indicated on the unit package, the Surgiwand™ II device contains DEHP. When used as indicated, transient exposure can be anticipated. In order to avoid undue risk of DEHP exposure in children and nursing or pregnant women, product should be used only as directed.

DESCRIPTION

The Surgiwand™ II single use suction and irrigation device without cautery consists of a 5 mm diameter shaft, 34 cm in length, with (9) orifices for pool suction and irrigation. This device also has a sealed portal to accommodate up to a 2.4 mm laser probe.

The Surgiwand™ II single use suction and irrigation device with unipolar cautery consists of a 5 mm diameter shaft with pool suction holes for enhanced suction and irrigation capabilities. The device is available with a nonremovable Spatula tip or L-Hook tip. The device can accommodate a gravity feed system, an electronic pump system and a nitrogen pump system.

INDICATIONS

The Surgiwand™ II device is designed for introduction and use through all appropriately sized Covidien™ trocar sleeves or larger sized trocar sleeves with the use of a converter. The device is intended for use in laparoscopic procedures as a multi-purpose operating tool with irrigation and suction capabilities.

CONTRAINDICATIONS

1. The Surgiwand™ II device with unipolar cautery is NOT intended for contraceptive coagulation of fallopian tissue, but may be used to achieve hemostasis following transection of the fallopian tube.
2. This device is not intended for use except as indicated.

WARNINGS AND PRECAUTIONS

1. This device is NOT intended to generate high pressure irrigant streams in hysteroscopic procedures or in encapsulated areas.
2. Endoscopic procedures should be performed only by physicians having adequate training and familiarity with endoscopic techniques. A thorough understanding of the operating principles, risks versus benefits, and the hazards involved in utilizing an endoscopic approach is necessary to avoid possible injury to the user and/or patient.
3. Verify mechanical and electrical compatibility of devices from different manufacturers prior to using them together in a procedure.
4. Always ensure that the outer sleeve is pushed fully forward over the tip of the device prior to insertion or removal of the device through the trocar sleeve. Failure to do so may cause damage to the device.
5. When not using suction, ensure that the trumpet valve providing suction does not remain depressed so as to minimize loss of pneumoperitoneum.
6. When using electrocautery, ensure that the outer sleeve is fully retracted to expose the tip of the device, as failure to do so may compromise discharge of energy from the electrode.
7. When using electrocautery, ensure that the electrode is not in contact with conductive irrigation fluid, as this may compromise discharge of energy from the electrode.
8. When using the device to suction or irrigate, ensure that the outer sleeve is pushed fully forward to cover the tip of the device.
9. This device is provided STERILE and is intended for use in a SINGLE procedure only. DISCARD AFTER USE. DO NOT RESTERILIZE.
10. Where indicated on the unit package, the Surgiwand™ II device contains DEHP. When used as indicated, transient exposure can be anticipated. In order to avoid undue risk of DEHP exposure in children and nursing or pregnant women, product should be used only as directed.
11. The following applies to the Surgiwand™ II device with unipolar cautery:

Regarding Patient Contact:

- The rated accessory voltage for these devices is 4500 volts (peak) or less. Do not use with output voltages greater than this rating. Excessive power levels may result in instrument malfunction and possible patient or user injury.
- The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth (for example operating table supports, etc.). The use of antistatic sheeting is recommended for this purpose.
- The cables to the surgical electrodes should be positioned in such a way that contact with the patient or other leads is avoided. Temporarily unused active electrodes should be stored in a location that is isolated from the patient.
- Skin-to-skin contact (for example between the arms and the body of the patient) should be avoided, for example by the insertion of dry gauze.
- For surgical procedures where the high-frequency current could flow through parts of the body having a relatively small cross sectional area, the use of bipolar techniques may be desirable in order to avoid unwanted tissue damage.

Regarding Fluids and Gases:

- The use of flammable anesthetics or oxidizing gases such as nitrous oxide (N₂O) and oxygen should be avoided if a surgical procedure is carried out in the region of the thorax or the head, unless these agents are sucked away.
- Over-insufflation of air or inert gas prior to high frequency surgery or laser assist gas could result in a gas embolism that creates a safety hazard.
- Explosive gas concentrations in the area of use of these devices could create a safety hazard and should be avoided.
- These devices are supplied sterile and intended for single use only. Should additional cleaning be or disinfection become necessary, use only non-flammable agents.
- Any flammable agents, solvents, or adhesives should be allowed to evaporate prior to high-frequency surgery. There is a risk of pooling of flammable solutions under the patient or in body depressions such as the umbilicus, and in body cavities such as the vagina. Any fluid pooled in these areas should be mopped up before these devices are used. Endogenous gases may ignite during high-frequency surgery, and some materials such as cotton, wool or gauze, when saturated with oxygen may be ignited by sparks produced during the normal use of these devices.

Regarding Equipment Compatibility:

- Allowing the energized tips of these devices to contact the uninsulated portion of other laparoscopic devices could result in patient or user injury.
- For patients with cardiac pacemakers or other active implants, a possible hazard exists because interference with the action of the pacemaker may occur, or the pacemaker may be damaged. In case of doubt, approved qualified advice should be obtained.
- When endoscopes are used with these devices, the patient leakage currents may be additive. This is particularly important if a Type CF Applied Part endoscope is used, in which case a Type CF Endoscopically-Used Accessory should be used in order to minimize total patient leakage current.
- When these devices and/or other endoscopic equipment are used with laser equipment, suitable eye protection or filtering should be employed in order to avoid potential eye damage to the operator.
- When these devices are used simultaneously with physiological monitoring equipment on the same patient, any monitoring electrodes should be placed as far as possible from the surgical electrodes. Needle monitoring electrodes are not recommended. In all cases, monitoring systems incorporating high frequency current limiting devices are recommended.

Regarding Use of the Equipment:

- The operator should regularly inspect these devices for possible damage, in particular, the electrode cables.
- To avoid shock and burn hazard to both the patient and operator(s) male cautery adaptor should remain covered when cautery not in use.
- Before each use, the outer surface of the portions of the endoscope and any endoscopically-used accessories which are intended to be inserted into a patient should be checked to ensure there are no unintended rough surfaces, sharp edges or protrusions which may cause a safety hazard.

1 SCHEMATIC VIEW

- A) SURGIWAND™ II DEVICE WITHOUT CAUTERY
- B) SHAFT OF DEVICE (34 cm in length)
- C) IRRIGATION VALVE
- D) SUCTION VALVE
- E) SURGIWAND™ II DEVICE WITH UNIPOLAR CAUTERY
- F) SHAFT/OUTER SLEEVE
- G) COLLAR
- H) CAUTERY ADAPTER
- I) LUER LOCK CONNECTOR
- J) SPATULA TIP
- K) L-HOOK TIP
- L) ELECTRONIC PUMP SYSTEM TUBING
- M) ELECTRONIC PUMP IRRIGATION TUBING
- N) OCCLUSION CLAMPS
- O) GRAVITY FEED IRRIGATION TUBING
- P) IRRIGATION BAG SPIKES
- Q) SUCTION TUBING

INSTRUCTIONS FOR USE

2 Attachment to the Suction Canister

When using the Surgiwand™ II device follow the instructions for attachment listed below.

1. Attach the suction tubing to the suction canister or wall unit by inserting the blue connector at the end of the tubing into the suction canister.

- Q) SUCTION TUBING
- R) BLUE CONNECTOR
- S) SUCTION CANISTER

3 Attachment to the Nitrogen Pump

When using the Surgiwand™ II device ONLY, follow the instructions for attachment listed below.

1. If using the nitrogen pump, attach the irrigation output from the pump to the luer lock connector on the instrument.

I) LUER LOCK CONNECTOR

2. Attach the suction tubing to the suction connection on the back of the instrument. If not using the Surgiwand™ II Device Suction Tubing, the suction tubing must have an inside diameter of 1/4". To attach the device to the suction canister, follow instruction #1 in Attachment to Suction Canister.

Q) SUCTION TUBING

4 Attachment to the Electronic Pump

When using the Surgiwand™ II device with Electronic Pump Tubing, follow the instructions for attachment listed below.

1. To hook up the irrigation tubing to the electronic pump, attach the irrigation tubes to the barb fitting on the electronic pump bottle.

T) IRRIGATION TUBING

U) BARB FITTING

V) ELECTRONIC PUMP BOTTLE

2. Attach the electronic pump tubing to the electronic pump bottle at the luer lock connector. The push-on connector attaches to the front of the electronic pump. If the second electronic pump bottle is not in use, clamp the irrigation tube which is not in use with the attached occlusion clamp, to eliminate the loss of irrigation fluid. To attach the device to the suction canister, follow instruction #1 in Attachment to Suction Canister.

W) ELECTRONIC PUMP

L) ELECTRONIC PUMP TUBING

I) LUER LOCK CONNECTOR

X) PUSH-ON CONNECTOR

N) OCCLUSION CLAMPS

Y) WALL SUCTION

5 Attachment to Gravity Feed System

When using the Surgiwand™ II device with gravity feed system (irrigation bags), follow the instructions for attachment listed below.

1. To hook up the irrigation tubing to the irrigation bag, insert the spikes into the irrigation bags. If the second irrigation bag is not in use, clamp the irrigation tube which is not in use with the attached occlusion clamp to eliminate the loss of irrigation fluid. To attach the device to the suction canister, follow instruction #1 in Attachment to Suction Canister.

Z) IRRIGATION DRIP BAGS

P) SPIKES

N) OCCLUSION CLAMPS

Y) WALL SUCTION

6 Laser Port Insertion

1. The laser port will accommodate a 2.4 mm or smaller diameter laser through the laser seal. The laser must have a minimum length of 45 cm to ensure that the laser tip will completely exit the shaft.

H) CAUTERY ADAPTER (ELECTROCAUTERY INSTRUMENT ONLY)

Aa) LASER PORT

Ab) SUCTION CONNECTION

Ac) IRRIGATION CONNECTION

NOTE: During initial insertion, additional force is required to puncture the seal.

The device can be held in two different ways:

2. Hold with thumb and forefinger on the side of the device like a pencil, or grasp the instrument in the palm of your hand.

3. Insert the device into a trocar sleeve.

4. For suction and irrigation, manually depress the desired suction or irrigation valve.

7 When using electrocautery with the Surgiwand™ II device with unipolar cautery

1. Ensure that the outer sleeve is fully retracted to expose the L-Hook tip or Spatula tip. To expose the tip, pull back on the white collar on the shaft of the device.

2. Attach the manufacturer's electrocautery cable to the cautery adapter provided on the back of the device.

3. To access the suction or irrigation capabilities with the device, ensure that the outer sleeve is pulled forward to cover the Spatula tip or the L-Hook tip. To uncover the tip of the device, pull the white collar back.

4. When removing or inserting the device from the trocar sleeve, ensure that the outer sleeve is pulled completely forward to cover the tip of the device. This will prevent any damage to the device or trocar sleeve.

STORE AT ROOM TEMPERATURE. AVOID PROLONGED EXPOSURE TO ELEVATED TEMPERATURES.

Blossom O'Bradovich

Blossom O'Bradovich
Regulatory Affairs – Contingent

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