sigvaris

Specialities

ARM SLEEVE ADVANCE

ARM SLEEVE TRADITIONAL

Gebruiksaanweisung	2	Инструкции за употреба	30
Instructions for use	3	Instrucțiuni de utilizare	31
Istruzioni per l'uso	4	Naudojimo instrukcija	33
Mode d'emploi	6	Lietošanas instrukcija	34
Brugsvejledning	7	Kasutusjuhend	35
Bruksanvisning	8	Інструкція з користування	37
Käyttöohjeet	10	Petunjuk penggunaan	38
Gebruiksaanwijzing	11	使用说明	39
Bruksanvisning	12	사용 설명서	40
Instrucciones de uso	14	使用方法	41
Instruções de uso	15	ຄຳມອນໃນການໃຊ້	42
Инструкция по применению	17	გამოყენების ინსტრუქცია	43
Instrukcja użytkowania	18	45	تعليمات الاستخدام
Návod k použití	20	46	دستورالعمل های استفاده
Upute za upotrebu	21	47	הנחיות שימוש
Návod na použitie	22		
Navodila za uporabo	23		
Uputstvo za upotrebu	25		
Οδηγίες χρήσης	26		
Kullanım talimatı	27		
	29		

Deutsch

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Specialities ARM SLEEVE ADVANCE /
Specialities ARM SLEEVE TRADITIONAL

English

Intended purpose

Specialities ARM SLEEVE ADVANCE / TRADITIONAL

Compression arm sleeves are aids for compression therapy. The products are intended for compression therapy for patients

with impaired venous and/or lymphatic outflow in the arm. The products are intended for repeat use and for application on the arms.

The effects of medical compression stockings

Compression therapy is used to improve quality of life in chronic venous diseases, for the prevention and therapy of venous edemas and venous skin changes, as well as for ulcer prevention, ulcer therapy, and thrombosis prevention.

Carefully controlled pressure on the tissue and veins improves the blood return in the extremities and increases lymphatic flow.

At the same time, it improves microcirculation and therefore the inflow to and outflow from the tissue.

Target group

The intended target group includes patients who have been prescribed an arm sleeve for individual compression therapy by the responsible physician or trained health professional. The target group also includes user groups from health care professions who apply the products to patients.

Indications & contraindications

Indications Specialities ARM SLEEVE TRADITIONAL:

For lymphedemas which regress easily and are soft; congestion and swelling following mammary amputation with and without postoperative radiation; post-traumatic and postoperative edemas following completion of basic treatment; polyarthritic tissue thickening; compression treatment for burns and for extensive scars; relapse prevention after successful decongestion with pronounced initial findings.

Indications Specialities ARM SLEEVE ADVANCE:

Prophylaxis of lymphedema following surgery, e.g. following amputation of the mamma with and without postoperative radiation; discrete lymphedema; relapse prevention of lymphedema following completed decongestive treatment; for mild and still soft lymphedemas.

Contraindications Specialities ARM SLEEVES ADVANCE / TRADITIONAL:

Medical arm sleeves must not be worn in the following cases:

- Advanced peripheral arterial occlusive disease PAD (radial arterial pressure < 60 mmHg, finger arterial pressure < 30 mmHg or TcPO₂ dorsum of the hand < 20 mmHg)
- Severe cardiac insufficiency (NYHA IV)
- Septic thrombophlebitis
- Compression of an existing arterial bypass (epifascial)

Relative contraindications, whereby benefits and risks must be weighed up in particular:

- Allergy against compression material
- Inflammatory processes (compression products to reduce inflammation, pain, or edema in leg erysipelas or cellulitis should only be used in combination with antibacterial treatment)

- Mild PAD (reduced radial arterial pressure in lateral comparison, but > 60 mmHg radial arterial pressure)
- Pronounced exuding dermatoses
- Blue or white coloration of the fingers, numbness of the hands
- Primary chronic polyarthritis
- Severe sensory disturbances of the extremity
- Advanced peripheral neuropathy (e.g., in diabetes mellitus)

Risks and side effects

Medical compression garments can cause skin necrosis and pressure damage to the peripheral nerves if handled improperly. Itching, scaling, and signs of inflammation may occur under the compression aids in cases of sensitive skin, so adequate skin care is advised.

A manifest lymphedema requires medically prescribed complex decongestive treatment. Regular medical check-ups are required to maintain successful therapy. If acute pain or skin irritation occurs while wearing, consult your physician immediately.

Special Instructions

Do not cut loose ends on the insides of the stocking. Protect the compression stockings from ointments, oils, and grease to preserve their service life. Compression stockings must not be worn over open wounds. When wearing compression stockings with grip tops, the skin must be regularly ventilated. Noncompliance with this information may limit safety and effectiveness of the product. All serious incidents that occur when using the product should be reported to the manufacturer and the competent supervisory authority.

Specialities ARM SLEEVE TRADITIONAL contains latex that does not come into direct contact with the skin.

Application

We recommend wearing special gloves to make it easier to put the stockings on and take them off. Sharp fingernails and rings can cause defects.

Donning

Make sure that your arm is dry.

1. Fold the upper half of the arm sleeve over the lower half.
2. Slip the sleeve with the overlap over the hand and over the arm until the overlap reaches the elbow. Do not pull, but spread the fabric evenly over the arm using gloves.
3. Pull the overlapped sleeve section over the elbow and push it gradually towards the shoulder.

Doffing

1. Take hold of the sleeve at the upper edge and pull downwards towards the elbow.
2. Work the fabric over your hand with the flat of the hand.

Donning and doffing aids

For easy donning and doffing of our arm sleeves, we recommend the donning aid DOFF N' DONNER.

Storage, wearing period, and service life

Unless otherwise prescribed by your physician, wear your compression stocking daily from morning to evening for optimal therapy results. Wearing and washing your compression stocking every day may cause the medical pressure gradient and elasticity of your compression stocking to diminish. The shelf life date is indicated on the packaging using the symbol of an hourglass. We guarantee maximum effect for the subsequent six months of wearing.

Specialities ARM SLEEVES ADVANCE / TRADITIONAL

Follow the instructions for care on the packaging and/or the textile labels. Wash inside a mesh laundry bag. The stocking must be stored at room temperature and in a dry place away from light. Keep away from direct heat.

Material abbreviations

PA	polyamide
EL	elastane
ED	elastodiene
CMD	modal (micro-modal)

Applies to the following products

Specialities ARM SLEEVE ADVANCE /
Specialities ARM SLEEVE TRADITIONAL

Italiano

Destinazione d'uso

Specialities ARM SLEEVE ADVANCE / TRADITIONAL

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SIGVARIS and SIGVARIS GROUP are registered trademarks of SIGVARIS AG,
© 2022 Copyright by SIGVARIS AG, CH-9014 St. Gallen/Switzerland
SIGVARIS GmbH, Dr.-Karl-Lenz-Strasse 35, 87700 Memmingen, **Germany**
SIGVARIS S.r.l., Via Galileo Galilei 2 / 2, 31057 Silea (TV), **Italy**
SIGVARIS S.A., ul. Powidzka 50, Gutowo Małe, 62-300, Września, **Poland**
SIGVARIS (Australia) Pty. Ltd., 8/53 Stanley Street, West Melbourne VIC 3003, **Australia**