DESTIKNEE™ Knee System

EN
(CR & PS)

Carefully read all instructions and be familiar with the surgical techniques prior to use.

STERILE EO

Description

The Meril Healthcares **Destiknee** Knee System is a system of components intended to replace the femoral and tibial articular surfaces of the knee joint. Components are available in many styles and sizes. The components are made from various types of metals and non-metallic materials. The component style, size, compatibility and specific component material is provided on the outside carton label. All implantable components are designed for single use only.

Product Selection Information

- Appropriate matching of the components will occur when the articular component is matched to the femoral component (by letter size designation and component style) and to the tibial base (by numerical size designation). For example a size B1+PS (Posterior Stabilized) articular component is appropriately matched to a size B PS femoral component and an extra size 1 or size 2 tibial base plate. Mismatching may result in poor surface contact and produce pain, decrease wear resistance, produce instability of the implant, or otherwise reduce implant life.
- Use only instruments and trials specifically designed for use with these devices to help ensure accurate surgical preparation, soft tissue balancing, and maximum motion function.
- Selections between the various sizes and options are matters of physician discretion.

Indications

- Severe knee joint pain, loss of mobility, and disability due to: rheumatoid arthritis, osteoarthritis, traumatic arthritis, polyarthritis.
- Correction of functional deformities.
- Post-traumatic loss of knee joint control, particularly when there is patellofemoral erosion, dysfunction, or prior patelloclay.
- Moderate valgus, varus, or flexion trauma.
- Knee fractures unreleasable by other methods.
- Misalignment of the knee joint and soft tissue integrity are present.

Contraindications

- History of infection in the affected joint that may affect the function of the implanted prosthesis.
- Less than optimal bone stock on femoral or tibial surfaces resulting from a history of disease, infection, or prior surgery.
- Compromised skeletal bone quality.
- Neuropathic disease that adversely affects the prosthetic joint.
- Osteoporosis or deficiency of musculature that compromises the affected limb.
- Pain free and stable arthrosis in an adequate functional position.
- Inactive knee joint secondary to negative collateral ligament integrity.
- Patients whose weight, age, or activity level may cause extreme loads and early failure of the system.

Precautions

- Before clinical use, the surgeon should thoroughly understand all aspects of the surgical procedure and the limitations of the product. Patients should also be instructed on the limitations of the product, including but not limited to the impact of patient weight and activity, and the weight to be placed on the knee joint.
- The prosthesis will not restore function to the level expected with normal healthy bone, and the patient should understand their expected level of function.
- As with all prosthetic implants, the durability of the components is affected by numerous biological, biomechanical, and other external factors which can limit their service life. Adherence to the indications, contraindications, precautions, and warnings for this product is essential for maximizing service life.

Adverse Effects

- Long term swelling or infection.
- No improvement in range of motion.
- Neuropathic disorders.
- Dislocations, bone fractures, and/or joint instability.
- For literature, there is a chance that wear of polyethylene components may result in bone resorption, loosening, and possible infection.
- Possibility for metal sensitivity reactions.
- Vascular thrombosis.
- Prolonged and excessive joint pain and/or inflammation.
- Aseptic loosening of implant.

Warnings

- This device is for single patient use only. Never reuse the implant, even though it may appear undamaged.
- Discard all damaged implants.
- Polished bearing areas must not come in contact with hard or abrasive surfaces.
- Bearing areas must be free of debris and clean prior to assembly.
- Contaminant or bending of the implant may reduce its fatigue strength and cause premature failure under load.
- All polyethylene tibial base plates should be limited to use in low demand (i.e. modest weight, activity level) patients with good bone quality.
- Return all packages with flaws in the sterile barrier to the supplier. Do not re-sterilize.

Materials

All of the materials used in the Meril Healthcares **Destiknee** Knee System meet ASTM standards for implants. The metal components comply with the F75 Standard Specification for Cobalt-28 Chromium-6 Molybdenum Alloy Castings and Casting Alloy for Surgical Implants (UNS R30075). The polyethylene components comply with the F663 Standard Specification for Ultra-High-Molecular-Weight Polyethylene Powder and Fabricated Form for Surgical Implants. The **Destiknee** Knee System has not been evaluated for safety and compatibility in the MR environment. The **Destiknee** Knee System has not been tested for imaging orientation in the MR environment.

Component	Material Grade	ASTM / ISO
Femoral	Cobalt Chromium-(F-75)	F-75 / ISO 5832-4
Tibial base Plate	Cobalt Chromium-(F-75)	F-75 / ISO 5832-4
UHMW-PE (GUR-1020)	F648 / ISO 5834-4	F648 / ISO 5834-4
Tibial Liner all Poly	UHMW-PE (GUR-1020)	F648 / ISO 5834-4
Patella	UHMW-PE (GUR-1020)	F648 / ISO 5834-4

Sterilization

- All implants are supplied sterile in protective packaging to a Sterility Assurance Level (SAL) of 10⁻⁶.
- The method of sterilization is provided on the outside carton label. Components sterilized through ETO gas.
- The packaging of all sterile products should be inspected for flaws in the sterile barrier before opening. In the event of such a flaw, the product must be returned non-sterile. Trial components should be used to avoid having to open any aspect of the sterile package prior to the component use.
- Care should be taken to prevent contamination of the component. In the event of contamination, this product should be discarded or returned to the supplier.
- If the package is opened, but the product is not used, the component should not be re-sterilized and should be discarded or returned to the supplier.

Patient Guidance

The probability of resulting complications and/or failure of knee prosthesis is increased in cases where the patient's physical and medical presentation (e.g. weight, or other diseases) is a detriment, the patient's functional outcome goals are above what is reasonably attainable, or the patient's expectations for knee joint function are unrealistic. Patients with unrealistic expectations for the outcome of their knee surgery are at a higher risk of post-operative injury can result in loosening, wear, and/or fracture of the knee implant. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail, or need to be replaced. The implant may not, and is not guaranteed to last the rest of the patient's life. Because prosthetic joints are not as strong, resilient, or durable as natural, healthy joints, all prosthetic knees may need to be replaced at some point.

WARNING: THIS DEVICE IS INTENDED FOR CEMENTED USE ONLY.

Meril Healthcare Pvt. Ltd.
Ground & First Floor, Survey No. 1734
and First Floor, H-145, Meril Park, Survey
No. 135/2/B & 1742, Mukundnagar Marg,
Chala, Vapi - 396191, Gujarat, India.
Customer Care Contact:
Tel. : +91 (260) 305 2100
E-mail: askinfo@merilite.com

Orthopedics

EC REP OBELIS S.A
Bd. General Vasile, 53
1030 Brussels, Belgium
T +32 2 732 59 54, F +32 2 732 60 03
E mail@obelis.net, W www.obelis.net

Rev. 02
Rev. 02 2019-12
MH/FUM/L004

Temperature limitation
Keep Dry
Keep Away from Sunlight
Do not Use If Package is damaged
Non Pyrogenic

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