SURE-Paste®

Dental Root Canal Filling Material / Calcium Hydroxide with Barium Sulfate

Description

Sure-Paste is a temporary root canal filling material for pulp capping and an infected root canal. It is also applied a good intracanal medicament. It contains Calcium Hydroxide and Barium Sulfate mainly, having good radiopacity and antibacterial characteristics. (pH≥12). Sure-Paste is very stable without any solidification or separation and is a flowable paste premixed in a convenient syringe. Also Sure-Paste has excellent accessibility to the apex and biocompatibility to the periapical tissue.

Composition 1. Calcium Hydroxide 2. Barium Sulfate 3. Thickening Agent 4. Glycol Derivatives 5. Water

Intended use To fill the root canal temporarily during endodontic treatment.

Indications

1. Treatment for an infected root canal
2. Pulp capping of exposed pulp
3. Disinfection after pulp extirpation treatment
4. Pulp capping of vital pulpotomy

Contra-Indications Do not apply as a permanent filling material or luting material.

Side effects When Sure-Paste is overfilled to beyond the apex, it may cause inflammation, tissue necrosis and nerve injury.

Recommended Procedure


1. After root canal preparation, rinse and dry the canal gently.
2. For your convenience, fix a free angle holder to the syringe.
3. Insert the tip of syringe into the root canal.
4. Apply paste by compressing the plunger gently and smoothly.
5. Withdraw the syringe gently when the paste fills the canal completely.
6. Remove excess paste on the surface.
7. Sure-Paste needs to be completely removed before final obturation.
8. Removal of Sure-Paste is accomplished using appropriate irrigation solution and instrument.
9. Remove the disposable tip of syringe, close the cap tightly and store in a cool place.

Precaution

1. For dental use only.
2. It should be used by a dentist and dental professional. Observe the procedure and the precautions.
3. Keep it away from a reach of non-dental users including children, senior citizens and others.
4. Do not press paste with excess force, or the patient may have some acute pain.
5. Do not overfill beyond the apex.
6. Avoid contact with skin, eye, and mucous membrane. Contact of calcium hydroxide with the skin, eyes and mucous membrane may cause local irritation and burns. If inadvertent contact occurs, rinse with copious amounts of water and seek medical advice.
7. Please keep the syringe cap tightly when not using Sure-Paste.
8. Do not reuse the disposable tip after use.
9. Do not use as an internal medicine.
10. Do not use after expiry date.
11. Do not use with the patients who have shown hypersensitivity or allergy to treatment materials.
12. Do not use if the paste has hardened, phase separation has occurred, or there are foreign substances.

Storage 1. Avoid exposure to direct sunlight. 2. Store in a cool place. (10°C-25°C, 50°F-77°F)

Symbols

Shelf Life : 3 years from the date of manufacture

Warranty

Sure Dent Corporation will replace a product that is proven to be defective. Sure Dent Corporation does not accept a liability for any damage or loss, direct or consequential, stemming from the use of or inability to use the product described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

SURE DENT CORPORATION

#809, 52 Sagimakgol-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, 13210, Republic of Korea
Tel) +82-31-737-4720 Fax) +82-31-737-4724 E-mail: suredent@suredent.com

Date of Issue: February 20, 2023, Instruction for Use of Sure-Paste - Rev.8 TFSOSP-IFU

Sure Dent Corporation

[Handwritten signature]
Hojin Bae / CEO

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Dental Root Canal Cleansers / EDTA Root Canal Preparation Cream

Description

Sure-Prep including Ethylenediamine Tetraacetic Acid(EDTA) is a chelating agent which allows efficient cleaning and easier preparation of the root canal. When Urea-peroxide in Sure-Prep reacts with sodium hypochlorite, pulp debris is lifted out the canal by the bubbling effect. Sure-Prep improves the tubular penetration of endodontic medicaments and filling materials as dentinal tubules are opened. Also, Sure-Prep improves the efficiency in injecting into the root canals by including a special lubricant.

Composition

Ethylenediamine Tetraacetic Acid (EDTA) 15%
Urea Peroxide 10%
PPG and PEG

Indications

1. Facile shaping and cleansing of the root canal
2. Easy and fast removal of debris
3. Encouraging decalcification of dentin
4. Simple handling of instruments in the canals due to lubricating effect

Recommended Procedure

1. Isolate the tooth with a rubber dam and open the pulp chamber.
2. Inject a proper amount of Sure-Prep into the pulp chamber or place Sure-Prep directly on the file.
3. File and ream canal in conventional manner.
4. Re-apply the above step if you needed.
5. With sodium hypochlorite, clean the root canal.


Storage

1. Keep Sure-Prep in a refrigerator to maintain the performance when not using.
2. Do not keep in a freezer not to freeze.
3. Avoid exposure to direct sunlight.

Precaution

1. For Dental Use Only
2. Avoid extrusion of the materials out of the apex.
3. Apply carefully not to have any hole in the root canal.
4. Always close the cap tightly after use.
5. Using protective gloves is recommended.
6. Do not use after expiry date.
7. If accidental contact with eyes or prolonged oral tissues occurs, flush immediately with large amounts of water.
8. It should be used by a dentist and dental professional. Observe the procedure and the precautions.
9. Keep it away from a reach of non-dental users including children, senior citizens and others.
10. Do not use this product in patients with known allergies to ingredients of the product.

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Date of Issue: February 20, 2023, Instruction for Use of Sure-Prep - Rev.7 TFSP-IFU

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Root Canal Sealing Material

Description

Sure-Seal Root is a convenient premixed ready-to use, injectable calcium aluminosilicate compound paste developed for permanent obturation of the root canal. **Sure-Seal Root** is based on a calcium aluminosilicate compound composition, which requires the presence of water to set and harden.

Sure-Seal Root does not shrink during setting and demonstrates excellent physical properties.

Setting time

Setting time is 25min. measured according to ISO 6876:2012 (100% humidity condition).
However, in normal root canals, the setting time can be more than 2.5 hours.

Composition Calcium Aluminosilicate Compound, Zirconium oxide, Filler and Thickening agent

Indications Permanent obturation of the root canal.

Recommended Procedure

1. Prepare and irrigate the root canal using standard endodontic techniques.
2. Remove the syringe cap from Sure-Seal Root syringe. Securely attach a disposable tip with a clockwise twist to the hub of the syringe. The disposable tip is flexible and can be bent to facilitate access to the root canal.
3. Insert the tip of the syringe into the root canal no deeper than the coronal one third(1/3).
Gently and smoothly dispense Sure-Seal Root into the apex of the root canal.
A radiograph can be taken at this stage to ensure an adequate apical seal. If the placement is not sufficient, thoroughly irrigate Sure-Seal Root out of the root canal and repeat the procedure.
4. Continue filling the root canal with Sure-Seal Root while withdrawing the disposable tip. Fill the root canal completely; avoid the formation of air bubbles and overfilling of the root canal.
5. Remove excess Sure-Seal Root from the top of the root canal with a moist cotton pellet.
6. Place the Gutta-Percha Points into the root canal using standard endodontic techniques.


Storage

1. Please keep the syringe cap tightly when not using.
2. Avoid exposure to direct sunlight.
3. Store in a cool place.

Precaution

1. For dental use only.
2. It should be used by a dentist and dental professional. Observe the procedure and the precautions.
3. Keep it away from a reach of non-dental users including children, senior citizens and others.
4. Do not press paste with excess force, because the patient may have some acute pain.
5. Do not overfill beyond the apex.
6. In accidental contact with eyes occurs, flush immediately with a large amount of water.
7. Do not reuse the disposable tip after use.
8. Do not use after expiry date.
9. Do not use as an internal medicine.

Symbols

Shelf Life : 2 years from the date of manufacture

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
Date of Issue: February 20, 2023, Instruction for Use of Sure-Seal Root - Rev.4 TFSSR-IFU

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