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
PERISO SA

Index

1	INTRODUCTION.....	4
1.1	Introduction to the CTU S wave device.....	4
1.2	Introduction to the manual	4
1.3	Safety Precautions and Warnings.....	5
1.4	Shipping damage	7
1.5	Packing contents.....	8
1.6	Warranty conditions	9
1.7	Jurisdiction	9
2	GLOSSARY OF SYMBOLS AND WRITING CONVENTIONS.....	9
2.1	Glossary of symbols.....	9
2.2	Writing conventions.....	12
3	INSTALLATION.....	13
3.1	Installation instructions.....	13
3.2	Operating handpieces	17
4	OPERATING INSTRUCTIONS.....	19
4.1	Skills required.....	19
4.2	Collateral effects.....	19
4.3	Operational section	20
4.4	Procedure for carrying out the treatment	25
5	CLINICAL INDICATIONS, CONTRAINDICATIONS AND PRECAUTIONS.....	25
5.1	Clinical indications	25
5.2	Contraindications	25
5.3	Precautions	26
6	MAINTENANCE and RESOLUTION OF THE PROBLEMS.....	28
6.1	Cleaning and disinfection	28
6.2	Guide to proper maintenance	28
6.3	Troubleshooting	29
7	TECHNICAL SPECIFICATIONS	31
8	INFORMATION PROVIDED BY THE MANUFACTURER.....	32

1 INTRODUCTION

1.1 Introduction to the CTU S wave device

CTU S wave is a physiotherapy and rehabilitation device that can be defined as a diamagnetic generator.

The diamagnetic generator exploits the principle of diamagnetism, that is a very weak effect of a quantum nature by virtue of which the materials involved, defined as diamagnetic, are characterized by magnetization with the opposite direction to the magnetic field, and therefore are weakly "rejected".

The effect of the diamagnetic repulsion, obtained with the use of high intensity and low frequency magnetic fields, determines the forward thrust of the lenses used for the therapy. This effect is combined with a particular focusing technique that uses Fresnel acoustic lenses, equipped with specific geometries and characteristics that allow an extremely precise focusing without aberrations of the generated waves. The CTU S wave diamagnetic generator allows to modify the temporal and spatial gradients thus obtaining a temporal focus and a spatial focus of the generated wave.

1.2 Introduction to the manual

This manual is an integral part of the equipment supplied with the CTU S wave equipment, constitutes documentation and support for use and is not intended for sale.

Reproduction, even partial, is prohibited without the explicit authorization of the manufacturer.

It is made by the technical department of the manufacturing company for the exclusive use with the CTU S wave equipment.

It is guaranteed that the manual delivered is the one relating to the equipment. Corresponding corrections will be included in new editions.

All the information contained may be subject, without notice, to changes by the Manufacturer for reasons of technical or commercial improvement.

The use of the CTU S wave implies the knowledge of this manual in all its parts. Use that does not comply with what is indicated in the manual automatically invalidates the guarantee and the responsibility of the manufacturer for what follows.

1.3 Safety Precautions and Warnings

- ❗ Before connecting the CTU S wave device to the power supply, it is necessary to make sure that it has the following characteristics:
 - Mains voltage 230 Volt AC (+/- 10%)
 - Network frequency 50-60 Hz
 - Power consumption 130vA (75W)
 - ❗ It is also necessary that the environmental conditions are the following:
 - Ambient temperature from 10 ° C to 30 ° C
 - Relative humidity from 30% to 75%
 - ❗ The environmental storage conditions must be as follows:
 - Ambient temperature from 0 ° C to 40 ° C
 - Relative humidity from 10% to 80%
 - ❗ The installation of the equipment should be avoided in the following cases:
 - Near heat sources
 - Exposed to rain or humidity
 - Exposed to direct sunlight
 - 🕒 The manufacturer declines any responsibility for damage or malfunction attributable to:
 - Power supply network not compliant with the above
 - Unsuitable environmental conditions
- Do not use the device:
- ❗ In the presence of combustible gases, flammable vapors, in oxygen chambers or in a detonating atmosphere.
 - ❗ In the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide.
 - ❗ Near equipment with strong electromagnetic emissions or equipment sensitive to electromagnetic emissions.
 - ❗ Near or on top of other appliances.
 - ❗ Keep the appliance away from water splashes, even accidental.
 - ❗ Avoid using the device in combination with monitoring devices.
 - ❗ Portable and mobile radio communication devices can influence the functioning of the electro-medical device.

Connection to the power supply

Connection of the appliance to the power supply network is simple and immediate and must be carried out using the power cable following the instructions below (chapter 3).

- ❗ To avoid the risk of electric shock, this appliance must only be connected to an earthed power supply.
- ❗ In the event of damage to the power supply cable, it must only be replaced by the manufacturer Periso or by personnel suitably trained by Periso.
- ❗ Any repair of the CTU S wave device must be carried out by the manufacturer Periso or by personnel suitably trained by Periso.
- ❗ The device is not permanently installed and uses the power plug as a means of isolation from the mains supply.
- ❗ Leave at least 90 cm of space between the machine and the power socket to be able to intervene in case of need to disconnect the device from the power socket.

Moving the device

The device is equipped with wheels to facilitate movement in various environments. To carry out any movement, it is advisable to grasp the device by the external frame making sure that the wheel brakes have been released, as shown in the figure.

- ❗ Make sure that the wheels are unlocked before moving the appliance in any way.

- ❗ Check the correct locking of the wheels every time you want to start a treatment.
- ❗ Before starting a treatment, check that the device is placed on a level floor.

1.4 Shipping damage

The equipment and all its components and accessories are subject to careful and accurate testing by the Quality Control Office, which certifies their integrity and perfect functioning. Any damage and / or failure or incorrect functioning of the same, must therefore be promptly communicated to the manufacturer and to be attributed exclusively to serious negligence on the part of the carrier. The original packaging must be kept in order for any future shipments to / from the manufacturer.

1.5 Packing contents

The CTU S wave device is carefully packaged and shipped with all the accessories and the manual, provided in the basic equipment. Below is a list of the items included in the shipping package.

Quantity	Description
1	User and Maintenance Manual
1	Service and Warranty Booklet
1	CE Declaration of Conformity
1	Method of Working
1	Trolley
1	USB cable
1	Single handpiece
1	Set of focal lenses (0-20-30-40-60-Mini13-Mini30)
1	Control module
1	Power module
1	Complex Gel Periso
1	Spare Parts – Screw kit

Optional

Quantity	Description
1	Double handpiece
2	Medium ranges
2	Large bands
1	Pair of double handpiece supports
1	Set of focal lenses (0-20-30-40-60-Mini13-Mini30)

1.6 Warranty conditions

Reference to the service and warranty booklet provided in the shipping wrapper. The expected service life of the CTU S wave device is 10 years old.


1.7 Jurisdiction

For any dispute concerning these conditions, the exclusive jurisdiction is that of Lugano.


2 GLOSSARY OF SYMBOLS AND WRITING CONVENTIONS

2.1 Glossary of symbols


Before starting, check for the check mark on "Server Online". Indicates the connection to the Server.


If checked, it indicates that the processor is functional.


If checked, it indicates that a patient username has been entered.


Patient username entry box.

Indicates that the patient's symptoms have been entered.


Patient symptoms, select based on the sensations and symptoms present.


It must be selected by the user, it indicates that the handpiece has been connected.

CHECK LIST incompleta


Error message that appears if any of the fields described above are not filled in or 'checked'.


Decreases or increases the intensity of the shock wave.


Decreases or increases the desired wave front.

Decrease or increase the hit rate.

Decreases or increases the intensity of the diamagnetic effect.

It decreases or increases the duration of therapy.

The "START" and "STOP" buttons start or stop the machine. After the "STOP" the therapy will continue for a few seconds.

"6 min" indicates the remaining therapy time.

Before starting a new therapy when the time comes goes to 0 you have to reset the duration again.

Label symbols

CE mark.

Temperature limits.

Indicates the name and address of the manufacturer of the medical device.

Humidity limits.

Consult instructions for use.

Atmospheric pressure limits.

Type BF applied component.

Keep away from sunlight and heat sources.

Non-ionizing emissions.

Keep away from humidity and rain.

It is mandatory to read the instructions for use.

It is forbidden to sit on the device.

Serial number.

Lot number.

It is forbidden to get on the device.

Waste electrical and electronic equipment.

It is forbidden to push the sideways device.

2.2 Writing conventions

Underline and / or bold

It highlights some parts of the document, of particular importance or in any case noteworthy.

Note

Notes highlight important information for specific use.

Caution

These messages appear before the description of operations which, if not observed, could cause damage to the equipment or its accessories, and to people.

Prohibition

This message is placed before operations that must never be performed.

3 INSTALLATION

The installation of the medical device must be carried out by Periso-trained technical personnel.

3.1 Installation instructions

The equipment must be extracted from the front of the packaging after having removed the safety locks. Then extract all the accessories and check their integrity.

The device is partially assembled, the only operations to be performed are listed below:

Switch
(Magnetothermic safety switch)

Pay attention to the yellow warning label.

PWR: tablet power cable.

USB or JACK: connection cable between the machine and the tablet.

NB: only one of the two is required, the cable must be a Periso product.

Back of the tablet

The connection cables with the machine are already connected.
Pay attention to the yellow warning label.


- ❗ The use of extension cords and / or multiple cables between the mains socket and the power cable is always not recommended as their possible malfunction could cause damage to the system. However, if it is necessary to use them, make sure that these components are fully functional and equipped with a grounding conductor.
- ❗ The device must be positioned in such a way that it is feasible to operate on the plug or switch if it is necessary to interrupt the connection to the electrical network.

Power on and off

To switch on the device, first use the switch located on the back of the device at the bottom of the power module.

To switch off the device, first make sure that the machine is not supplying power, then turn off the device using the switch on the back of the device.


- ❗ The appliance must always be switched on using the appropriate switch! The machine must never be turned off by unplugging the power cable to avoid software malfunctions or damage to the software.

Emergency button

If it is necessary to IMMEDIATELY stop the operation of the device, push the red button in the figure.

To unlock the device, turn the button.


- ⊙ Modifications of the device are not permitted.
- ❗ Never switch off the machine by unplugging the power cord, to avoid software malfunctions or possible damage.

Disconnection of the power and control module from the truck

If it is necessary to remove the power and control module from the truck proceed as follows:

First of all, disconnect the USB or JACK cable and the power cable, then unscrew the two black screws highlighted in the photo and then proceed to remove the power module from the trolley.

Remove the 4 screws shown in the photo (2 on the left and 2 on the right). Then remove the control module from the carriage.

These operations, or one of the two, may be necessary if the power module, or the control unit, must be sent to Technical Assistance.

- ❶ Such intervention may only be carried out by Periso-trained technicians.
- ❶ During these operations it is important to pay MAXIMUM ATTENTION and in case of doubts about the execution procedure, consult the Technical Assistance to request clarifications on the matter.

E-mail: assistancemed@periso.ch Telephone number: +41 919359400

3.2 Operating handpieces

The CTU S wave has two operating handpieces:

Single handpiece

Double handpiece (OPTIONAL)

The energy delivered by the device is divided between the two operating handpieces. The operator takes this into account in relation to the type of therapeutic treatment to be administered to the patient. The presence of the operator is always recommended during the development of the therapy.

The double handpiece is connected to a single connector. They can be placed in two different points on the body and work using a different focal lens for each point.

The operating handpieces are the only applied parts of the machine.


Use of operating handpieces

The therapy is administered in direct contact with the patient, applying the emitting head of the handpiece on the patient's skin. A lens for acoustic propagation is interposed between the head of the handpiece and the patient's skin and an adequate quantity of gel is interposed between the lens and the tissue.

Fresnel lenses must be seated and screwed on without tightening them too much. The choice of lenses depends on the depth of focus and therefore on the pathology to be treated.

The lenses are numbered as follows:

- 00 - defocused;
- 20 - focus with depth 2 cm;
- 30 - focus with depth 2.5 - 3.5 cm;
- 40 - focus with depth 3.5 - 4.5 cm;
- 60 - focus with depth 5.5 - 6.5 cm.

- ❗ The treatment must be performed always placing a lens between the handpiece and the skin and between the lens and the skin gel must necessarily be used.

4 OPERATING INSTRUCTIONS

4.1 Skills required

The CTU S wave is a device intended to be used by medical and / or physiotherapy personnel. It can be used in:

- Physiotherapy and hospital outpatient setting
- Rehabilitation in general

4.2 Collateral effects

Improper use of the device and / or excessive exposure of the patient to the treatment beyond the times and methods adequately adapted by the operator (doctor or paramedic and in any case qualified), based on the characteristics of the pathology and of the subject under examination , can cause cell damage as well as side effects of a nature not specifically known to date.

- ❗ The manufacturer is not responsible for the improper use of medicines, overdoses or inappropriate dosages.

4.3 Operational section

The CTU S wave device has an internet connection and, in order to perform treatment, can be used in online mode or offline mode. The internet connection allows the device to connect to an online database, on which each therapist/physician/clinic is registered with his/her own account. This account is associated exclusively with the purchased device, and the credentials provided (username and password) are used to register their patients on the database via codes, which are needed to initiate treatment in online mode.


The steps for starting treatment in both modes are outlined below.

ONLINE MODE


1. Turn on the device using the appropriate key. The following screen will appear:


- Before opening the "H4" application on the tablet, check that the internet connection is on. In the upper right corner, sliding your finger down will open the dropdown where you can click on the connection icon:


- Open the "H4" application. The following loading screen will appear:

Next, the application will open:


4. In order to begin treatment, all items on the checklist must be completed:


- Check from the checklist that the device is connected to the Internet: "Online server".
- Verify from the checklist that the power module is properly connected to the device: "Power module Online".

- In the space provided, enter the "patient code" you want to use for treatment:

When entered correctly, the "Select Patient" item appears complete on the checklist.


- Enter the symptoms:

Static pain: pain at rest
 Touching pain: pain on palpation
 Movement difficulties: motor difficulties
 Generic disfunction: general problems

subsequently appears on the complete checklist the item "Symptomatology".

- Click on "Probe Ready" when the handpiece is ready for treatment with the chosen lens properly inserted.
5. Select the desired action mode and treatment parameters, along with the minutes:


- Energy = shock wave power
- Rise = Rise time, speed at which the shock wave develops
- Repetition = pulses frequency
- DIA = diamagnetic effect intensity


6. Press start to start the treatment.
7. Refer to the document "MOD015 INSTRUCTION FOR PATIENT REGISTRATION IN PERISO DATABASE_FINAL USER" for the patient code registration procedure.

OFFLINE MODE

1. See section "online mode".
2. Before opening the "H4" application on the tablet, check that the internet connection is not active. In the top right corner, sliding your finger down will open the dropdown where you can click on the connection icon.
3. See section "online mode".
4. In order to start a treatment, even in offline mode, all items on the checklist must be completed:

Click on "Force OFFline". Automatically the checklist will appear as follows:

You have to complete the patient's symptoms anyway, so that you also complete "Symptomatology" on the checklist.

The red writing reminds you that you are using the device offline and that with this mode the number of treatments is limited. Once the therapies available in offline are finished, at least one therapy must be done in online mode to reload them.

- 5-6. See section "online mode".

4.4 Procedure for carrying out the treatment

The procedure for carrying out the treatment using the CTU S wave device can be divided into four steps:

- 1) Choose the area of the body to be treated. Prepare the patient for treatment in the most appropriate position.
- 2) Choose the lens, or lenses if you choose the two handpieces, to be mounted on the head of the handpiece. The lens, or lenses, will be chosen according to the area and the pathology to be treated.
- 3) Select the treatment parameters from the work screen.
- 4) Spread an adequate amount of gel on the area to be treated.
- 5) Place the handpiece on the area to be treated and start the treatment.
- 6) Adhere as closely as possible to the protocols provided in the Method of Working supplied for more effective treatment and to avoid possible overheating of the operating handpiece.

5 CLINICAL INDICATIONS, CONTRAINDICATIONS AND PRECAUTIONS

5.1 Clinical indications

The following are the therapeutic indications on the musculoskeletal system:

- **Bone tissues:** delayed consolidation, pseudo-arthritis, stress fractures.
- **Soft tissue:** chondrocalcinosis, joint stiffness, calcification / ossification, ossifying myositis, quadriceps tendinopathy, ligament fibromatosis, fascia fibromatosis.
- **Tendinopathies:** calcific tendinopathy of the shoulder, epicondylitis, trochanter tendinopathy, tendinopathy
- goose paw, knee tendinopathy, plantar fasciitis, patellar tendon, Achilles tendon.
- **Other pathologies:** skin ulcers, myofascial pain, scarring, spasticity.

5.2 Contraindications

IMPORTANT

Before carrying out one or more treatments with the CTU S wave device, the patient must be subjected to an anamnesis in order to detect:

- Unsuitable physiological states
- Presence of ferromagnetic material inside the areas of the body to be treated (this in the case in which the diamagnetic effect is used).

If the diamagnetic effect is used, the application is limited to the treated area, and it is therefore postponed to a clinical evaluation by the doctor, in the following subjects with:

- Fertile epiphyses
- Splinters, ferromagnetic metal fragments
- Clips on previous interventions
- Clips on aneurysms (blood vessels), aorta, brain
- Heart valves
- Spinal distractors
- Infusion pump for insulin or other drugs
- Heart pacemaker and other electronic devices
- Neurostimulators, electrodes, brain or subdural implants
- Spinal or ventricular lead
- Lens prosthesis

If the diamagnetic effect is used, the application is absolutely forbidden on subjects with:

- Severe heart disease
- Tumors
- Fertile epiphyses
- Pregnancies

If the diamagnetic effect is used, it is advisable for people with low blood pressure (hypotension) to check their blood pressure before and after the treatment and, if necessary, reduce the time of the session and allow more time to pass between sessions and the other one.

Therapy via the CTU S wave device cannot be performed in the following cases:

- Acute soft tissue / bone infection
- Pernicious primary diseases
- Epiphysiolysis in the focal point
- Blood clotting diseases (patients on anticoagulant therapy)
- Pregnancy
- Patients with pacemakers
- Lung tissue at the focal point (avoid the chest)
- Brain, spinal cord, large nerves in the focal point (neurocranium, spine, ribs).

5.3 Precautions

- Y Make sure you know the health status of the patient to be treated.
- Y If the equipment emits smoke, excessive heat or noise, stop the device immediately.
- Y If you notice an irregular thermal rise in the operating handpiece, stop the treatment immediately, as this indicates a malfunction.

- Y The device must only be operated by personnel trained in the use of the device and qualified to use it.
- Y In the absence of personnel trained to use the machine, switch off the device using the switch located on the back of the appliance.
- Y Treat the operating handpieces with care, taking care not to cause damage to the surface (scratches, falls, bumps, etc.).
- Y At each use, check the perfect integrity of the connection cables and their connectors, and of all moving or fixed parts.
- Y The use of accessories, transducers or cables, other than those supplied by the manufacturer, is prohibited, as they could cause an increase in emissions or a decrease in the immunity of the device.
- Y Make sure that the body of the equipment does not come into contact with liquid substances.
- Y In case of malfunctions, stop the treatment and notify the assistance.
- Y In case of damage or deterioration of the cables or the power supply, contact assistance.
- Y The device must be disposed of in accordance with the disposal of WEEE.

Products made by PERISO SA are marked with the following symbol:

On the basis of Art. 28 of Legislative Decree 14 March 2014, n. 49 all electrical and electronic equipment to be placed on the market after 13 August 2005 must bear a trademark, which must be visible, legible and indelible, which makes it possible to unequivocally identify the manufacturer of the EEE and the symbol indicating that the product must be disposed of separately.

Information to users on the disposal of WEEE waste

As expected in the Ordinance concerning the restitution, the resumption and the disposal of applianceselectrical and electronic, January 14, 1998 (ORSAE), His equipment cannot be disposed of as mixed urban waste but separate collection must be carried out.

Electrical and electronic appliances are those covered. The Ordinance obliges traders, manufacturers and importers to take back - at no additional cost - the appliances included in their assortment. In turn, consumers are obliged to return the appliances they intend to discard.

Used devices, in fact, must not be disposed of together with municipal waste or bulky waste but must be delivered, free of charge, to all official SENS and SWICO sales points or delivery centers for their disposal.

With regard to disposal in countries outside Switzerland, the user and the distributor must comply with the regulations in force in the country where the device is used.

6 MAINTENANCE and RESOLUTION OF THE PROBLEMS

6.1 Cleaning and disinfection

Cleaning the appliance body must be done with a damp cloth and using neutral and non-aggressive detergents.

Do not use alcohol to clean any part.

Focal lenses should be cleaned and disinfected with a damp cloth, lukewarm water and mild detergent after each use to prevent hard-to-remove gel deposits that can potentially damage accessories.

For a better cleaning of the lenses, use a soft-bristled toothbrush to remove the gel between the lens rings after each treatment to avoid caking.

After cleaning the lenses make sure to dry them completely.

Clean the single / double handpiece and focal lenses on first use.

- ❗ Make sure that the device is disconnected from the power supply during the cleaning phase.
- ❗ Use only the gel produced by Periso SA for treatments with the CTU S wave device.

6.2 Guide to proper maintenance

- ❗ The device must be reviewed annually to ensure that its characteristics and performance are not altered.
- ❗ The appliance must be installed and commissioned in accordance with the information contained in this manual.
- ❗ Do not leave the equipment exposed to humidity, rain, direct sunlight.
- ❗ Do not spill liquids on the equipment and do not place any weight on it.
- ❗ Do not leave live cables and / or in passageways, in order to avoid accidental breakage.
- ❗ Only connect the device to a system that complies with national standards.

Type of verification	Frequency	Connection to the power supply network
Power on	Daily	Yes
Cable integrity dispenser handpieces	Daily	No
Network cable integrity	Daily	No
Integrity of power-control module connection cable	Monthly	No
Integrity of dispensing handpieces	For each use	No
Focal lens integrity	For each use	No

6.3 Troubleshooting

Problem	Possible solution
Power Supply	
The device does not switch on.	Y Is the power cable connected to the 230V 50 / 60Hz mains?
	Y Is the magneto-thermal switch set to ON?
	Y Check that the emergency button is not pressed.
The display does not switch on.	Y Is the power cable connected to the 230V 50 / 60Hz mains?
	Y Is the magneto-thermal switch set to ON?
	Y Check that the emergency button is not pressed. Make sure the USB or JACK between the control and power module is correctly connected.
	Y Make sure the power cable (PWR) is properly connected to both the power module and the tablet.

Erogation

The device turns on but the handpiece does not deliver

- Y Check that a work program has been set, that START has been pressed on the display and the boxes are all ticked.
- Y Make sure the handpiece is properly connected to the power model. Check the integrity of the cable. Check the integrity of the connector.
- Y Make sure the USB cable between the control module and the power module is properly connected
- Y Make sure the handpiece is properly connected to the power model.
- Y Check the integrity of the cable.
- Y Check the integrity of the connector.

If the problem persists, contact PERISO SA support to send the device to the service center. It is mandatory to make a careful description of the problem by sending a written communication to the assistance center at the following e-mail address: assistancemed@periso.ch

7 TECHNICAL SPECIFICATIONS

Template	CTU S wave
Trading name	Diamagnetic Shock Wave
Power supply	50/60 Hz 230 Vac
Current consumption	0.57 A
Absorbed power	130 VA (75 W)
Conversion factor	0.56
Stored energy	70 J
Internal protection	1 fuse of 3.15 A + 1 fuse of 5 A 1 thermal magnetic breaker of 4.5 A
	<u>Not replaceable by the user</u>
Minimum inter-pulse time	10 sec
Visualization	DigitalTouch Screen display
Headquarters European authorized representative	Studio Murgia S.R.L. Via F. Benzi 1, 22100 Como (IT)
Producer	PERISO SA Via Senago 42d 6912 Pazzallo Switzerland CHE-101.514.101 VAT
Medical class	II A
Dimensions	48 x 45 x 121 cm
Weight	About 30 kg
Working temperature	From 10° C to 30 ° C
Storage temperature	0 ° C to 40 ° C
Working humidity	30% to 75% non-condensed
Storage humidity	10% to 80% non condensed
Atmospheric pressure	70 to 100 kPa

8 INFORMATION PROVIDED BY THE MANUFACTURER

8.1 Power module labels

UE labels A, values varies by country

PERISO^{SA}
technology that evolves

EC REP **Studio Murgia S.R.L.**
Via F.Benzi 1, 22100 Como (IT)

Manufacturer: PERISO sa - Via Senago 42/d, 6912 Pazzallo (CH)
Model: CTU S wave
General Power Supply: 220V 60Hz
Power Absorbed: 0.57A - 130VA (75W)
Power factor: 0.56
Storage Energy: 70J
Protection: Differential Thermal Magneto
Internal Protection: Fuses 1 x 3.15A, 1 x 5A
Minimum interpolation time: 10 sec

Ver. 3
Swiss Made
Software Ver. 1.0.0
Design Vers. 0
Weight: 30 kg

Warning: Annual Periodic Inspection

CE 279  10°C 30°C

A Medical Device II A Class

PERISO^{SA}
technology that evolves

EC REP **Studio Murgia S.R.L.**
Via F.Benzi 1, 22100 Como (IT)

Manufacturer: PERISO sa - Via Senago 42/d, 6912 Pazzallo (CH)
Model: CTU S wave
General Power Supply: 110V 60Hz
Power Absorbed: 0.57A - 130VA (75W)
Power factor: 0.56
Storage Energy: 70J
Protection: Differential Thermal Magneto
Internal Protection: Fuses 1 x 3.15A, 1 x 5A
Minimum interpolation time: 10 sec

Ver. 3
Swiss Made
Software Ver. 1.0.0
Design Vers. 0
Weight: 30 kg

Warning: Annual Periodic Inspection

CE 279  10°C 30°C

A Medical Device II A Class

Extra UE labels A, values varies by country

PERISO^{SA}
technology that evolves

PERISO sa
Via Senago 42 d , 6912 Pazzallo (CH)

Model: CTU S wave
General Power Supply: 110V 60Hz
Power Absorbed: 0.57A - 130VA (75W)
Power factor: 0.56
Storage Energy: 70J
Protection: Differential Thermal Magneto
Internal Protection: Fuses 1 x 3.15A, 1 x 5A
Minimum interpolation time: 10 sec

Ver. 2
Swiss Made
Software Ver. 1.0.0
Design Vers. 0
Weight: 30 kg

Warning: Annual Periodic Inspection

CE 279  10°C 30°C

A Medical Device II A Class

PERISO^{SA}
technology that evolves

PERISO sa
Via Senago 42 d , 6912 Pazzallo (CH)

Model: CTU S wave
General Power Supply: 220V 50/60Hz
Power Absorbed: 0.57A - 130VA (75W)
Power factor: 0.56
Storage Energy: 70J
Protection: Differential Thermal Magneto
Internal Protection: Fuses 1 x 3.15A, 1 x 5A
Minimum interpolation time: 10 sec

Ver. 2
Swiss Made
Software Ver. 1.0.0
Design Vers. 0
Weight: 30 kg

Warning: Annual Periodic Inspection

CE 279  10°C 30°C

A Medical Device II A Class

PERISO^{SA}
technology that evolves

PERISO sa
Via Senago 42 d , 6912 Pazzallo (CH)

Model: CTU S wave
General Power Supply: 230V 50/60Hz
Power Absorbed: 0.57A - 130VA (75W)
Power factor: 0.56
Storage Energy: 70J
Protection: Differential Thermal Magneto
Internal Protection: Fuses 1 x 3.15A, 1 x 5A
Minimum interpolation time: 10 sec

Ver. 3
Swiss Made
Software Ver. 1.0.0
Design Vers. 0
Weight: 30 kg

Warning: Annual Periodic Inspection

CE 279  10°C 30°C

A Medical Device II A Class

Annual Periodic Inspection

PERISO^{SA}
technology that evolves

Date of Production XX../.....20XX.....

Year Review Reference ID Signature

Year Review Reference ID Signature

Year Review Reference ID Signature

Year Review Reference ID Signature

Year Review Reference ID Signature

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LOT **09/2019**

SN **1B** **??????????** **D**

8.4 Operating handpiece connector label

8.5 Label on the cart


8.6 Applied parts label (handpiece)

8.7 Packaging label

Model: CTU S wave
 ID PERISO: 30.00.014
 HS Code: 9018.90
 Dimensions: 80x60xh135 cm
 Weight: 75.5 kg +/- 5%

PERISO^{SA}
 technology that evolves


EU Representative: Studio Murgia S.R.L.
 Via F. Benzi 1, 22100 Como (IT)
 Manufacturer: Periso SA
 Via Senago 42d, 6912 Pazzallo (CH)

PERISO CODE	DESCRIPTION
A388.050	○ 1 User and Maintenance Manual
A388.025	○ 1 Service and Warranty Booklet
A388.026	○ 1 CE Declaration of Conformity
A388.051	○ 1 Method of Working
30.04.012	○ 1 Trolley
A145.004	○ 1 USB cable
30.01.013	○ 1 Single handpiece
30.10.015	○ 1 Set of focal lenses (0-20-30-40-60-Mini13-Mini30)
30.03.029	○ 1 Control module
30.03.031	○ 1 Power module
A716.005	○ 1 Complex Gel Periso
30.10.033	○ 1 Spare Parts – Screw kit

**Environmental
conditions of
operation:**

**Environmental
conditions of
storage:**

**Avoid storage in
the following
cases:**

OPTIONAL

30.10.010	○ 1 Double handpiece
A1011.002	○ 2 Medium ranges
A1011.003	○ 2 Large bands
A495.595	○ 1 Pair of double handpiece supports
30.10.015	○ 1 Set of focal lenses (0-20-30-40-60-Mini13-Mini30)

Medical device class IIA

Rev.13 of 24/11/2023

Tables of electromagnetic emissions, electromagnetic immunity, recommended separation distances

TABLE 1		
GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS - FOR ALL EQUIPMENT AND SYSTEMS		
The CTU S wave are intended for use in the electromagnetic environment specified below. The customer or the CTU S wave should assure that it is used in such environment.		
Emissions Test	Compliance	Electromagnetic environment - Guidance
RF emissions CISPR 11	Group 1	The equipment CTU S wave uses RF energy only for its internal function, which is to connect to Periso's server via Wi-Fi of the Tablet. Frequency band of transmission: 2,4GHz band and 5GHz band. Radio module type: Wi-Fi 802.11 b/g/n/ac + Bluetooth 4.1 Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The equipment CTU S wave is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics emissions IEC 61000-3-2	Class A	
Voltage fluctuation / flicker emissions IEC 61000-3-3	Complies	

TABLE 2			
GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY - FOR ALL EQUIPMENT AND SYSTEMS			
The equipment CTU S wave is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment CTU S wave should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV in air-air	± 8 kV in air-air	The floors must be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	The mains voltage quality should be that of a typical commercial or hospital environment.
Over voltage Surge IEC 61000-4-5	± 1 kV line-line ± 2 kV line - earth	± 1 kV line-line ± 2 kV line - earth	Mains voltage quality should be that of a typical commercial or hospital environment. The mains voltage quality should be that of a typical commercial or hospital environment.

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T. (> 95% dip in U_T) for 0.5 cycle 40% U_T. (60% dip in U_V) for 5 cycles 70% U_T. (30% dip in U_T) for 25 cycles <5% U_T. (> 95% hole in dip in U_T) for 5 sec	<5% U_T. (> 95% dip in U_T) for 0.5 cycle 40% U_T. (60% dip in U_V) for 5 cycles 70% U_T. (30% dip in U_T) for 25 cycles <5% U_T. (> 95% dip in U_T) for 5 sec with $U_T = 230$ V ac 50 Hz	The mains voltage quality should be that of a typical commercial or hospital environment. If the user of the CTU S wave device requires continued operation also during the voltage mains interruptions, it is recommended to power the CTU S wave device with an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A / m	3 A / m	The mains frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note U_T is the ac mains voltage prior to application of the test level.			

TABLE 3			
GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY- FOR EQUIPMENT AND SYSTEMS THAT ARE NOT LIFE-SUPPORTING			
The equipment CTU S wave is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment CTU S wave should assure that it is used in such an environment.			
<i>Immunity Test</i>	<i>IEC 60601 test level</i>	<i>Compliance level</i>	<i>Electromagnetic environment - Guidance</i>
Conducted RF IEC 61000-4-6 Radiated RF Radiated RF IEC 61000-4-3	3 Veff_Vrms from 150 kHz to 80 MHz 150 kHz to 80 MHz 3 V / m 80 MHz to 2.5 GHz 80MHz to 2.5 GHz	3 V 3 V / m	Portable and mobile RF communications equipment should be used no closer to any part of the equipment CTUS wave, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d \geq \sqrt{1.17 P}$ $d \geq \sqrt{1.17 P}$ from 80 MHz to 800 MHz 80 MHz to 800 MHz $d \geq \sqrt{2.33 P}$ 800 MHz to 2.5 GHz 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

Notes:

- ! At 80 MHz and 800 MHz, the higher frequency range applies.
- ! These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment CTU S wave is used exceeds the applicable RF compliance level above, the equipment CTU S wave should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment CTU S wave
- b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V₁] V / m.

TABLE 4

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE EQUIPMENT CTU S wave - FOR EQUIPMENT AND SYSTEMS THAT ARE NOT LIFE-SUPPORTING

The equipment CTU S wave is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the equipment CTU S wave can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the equipment CTU S wave

as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter (m)		
	From 150 kHz to 80 MHz	From 80 MHz to 800 MHz	From 800 MHz to 2.5 GHz
	$d \geq 1.17 \sqrt{P}$	$d \geq 1.17 \sqrt{P}$	$d \geq 2.33 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.80	3.80	7.30
100	12.00	12.00	23.00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note

- At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

PERISO SA

Via Senago 42 d 6912 Pazzallo CH

Tel. +41 (0) 919359400

info@periso.ch - www.periso.ch

VAT number CHE-101.514.101 VAT