BeGraft

coronary

Coronary Stent Graft System

en	Coronary Stent Graft System
bg	Система с коронарна присадна протеза
cs	Systém koronárního stentgraftu
de	Koronares Stentgraftsystem
da	Koronart stentgraftsystem
es	Sistema de endoprótesis coronaria recubierta
et	Koronaarse stentsiiriku süsteem
fi	Sepelvaltimon stenttisiirrejärjestelmä
fr	Système d'endoprothèse coronaire couverte
el	Σύστημα στεφανιαίου μοσχεύματος στεντ
hr	Sustav koronarnog stent grafta
hu	Koronária sztentgraft rendszer
it	Sistema di endoprotesi coronarica
lt	Koronarinio stento implantavimo sistema
lv	Koronāro stentu implantātu sistēma
nl	Coronair stentgraftsystem
no	Koronart stentgraftsystem
pl	System stentgraftu wieńcowego
pt	Sistema endoprótese coronária coberta
ro	Sistem de endoproteză coronariană
ru	Стент-графт коронарный с системой доставки
sk	Systém koronárneho stentgraftu
sl	Sistem koronarnih žilnih opornic
sv	Koronarstentgraftsystem
tr	Koroner Stent Graft Sistemi

en Caution **bg** Внимание **cs** Pozor **de** Achtung **da** Forsigtig **es** Precaución **et** Ettevaatus! **fi** Huomio **fr** Mise en garde **el** Προσοχή **hr** Oprez **hu** Figyelem! **it** Attenzione **lt** Dėmesio! **lv** Uzmanību! **nl** Pas op **no** Forsiktig **pl** Przestroga **pt** Atenção **ro** Atenție **ru** Внимание! **sk** Pozor **sl** Pozor **sv** Försiktig **tr** Dikkat

en Consult instructions for use **bg** Направете справка с инструкциите за употреба **cs** Nahlédněte do návodu k použití **de** Gebrauchsanweisung beachten **da** Se brugsanvisningen **es** Consulte las instrucciones de uso **et** Vt kasutusjuhiseid **fi** Lue käyttöohjeet **fr** Consulter le mode d'emploi **el** Συμβουλευτείτε τις οδηγίες χρήσης **hr** Pogledajte upute za uporabu **hu** Lásd a használati útmutatót **it** Consultare le istruzioni per l'uso **lt** Skaitykite naudojimo instrukcijas **lv** Skatīt lietošanas instrukcijas **nl** Raadpleeg de gebruiksaanwijzingen **no** Se bruksanvisningen **pl** Patrz instrukcja użycia **pt** Consultar as instruções de utilização **ro** Consultati instrucțiunile de utilizare **ru** Ознакомьтесь с инструкцией по применению **sk** Pozri návod na použitie **sl** Glejte navodila za uporabo **sv** Läs bruksanvisningen **tr** Kullanım için talimatlara başvurun

en Medical Device **bg** Медицинско изделие **cs** Zdravotnický prostředek **de** Medizinprodukt **da** Medicinsk udstyr **es** Producto sanitario **et** Meditsiinisead **fi** Lääkinnällinen laite **fr** Dispositif médical **el** Ιατροτεχνολογικό προϊόν **hr** Medicinski proizvod **hu** Orvostechnikai eszköz **it** Dispositivo medico **lt** Medicinos priemonė **lv** Mediciniska ierīce **nl** Medisch hulpmiddel **no** Medicinsk utstyr **pl** Wyrób medyczny **pt** Dispositivo médico **ro** Dispozitiv medical **ru** Медицинское изделие **sk** Zdravotnícky pomôcka **sl** Medicinski pripomoček **sv** Medicinteknisk produkt **tr** Tıbbi cihaz

en Manufacturer **bg** Производител **cs** Výrobce **de** Hersteller **da** Producent **es** Fabricante **et** Tootja **fi** Valmistaja **fr** Fabricant **el** Κατασκευαστής **hr** Proizvođač **hu** Gyártó **it** Fabbricante **lt** Gamintojas **lv** Ražotājs **nl** Fabrikant **no** Tilvirker **pl** Producent **pt** Fabricante **ro** Producător **ru** Изготовитель **sk** Výrobca **sl** Proizvajalec **sv** Tillverkare **tr** Üretici

en Date of manufacture **bg** Дата на производство **cs** Datum výroby **de** Herstellungsdatum **da** Fremstillingsdato **es** Fecha de fabricación **et** Valmistamiskuupäev **fi** Valmistuspäivä **fr** Date de fabrication **el** Ημερομηνία κατασκευής **hr** Datum proizvodnje **hu** Gyártás dátuma **it** Data di fabbricazione **lt** Pagaminimo data **lv** Izgatavošanas datums **nl** Productiedatum **no** Fremstillingsdato **pl** Data produkcji **pt** Data de fabrico **ro** Data fabricației **ru** Дата изготовления **sk** Dátum výroby **sl** Datum proizvodnje **sv** Tillverkningsdatum **tr** Üretim tarihi

en Number of units **bg** Брой елементи **cs** Počet jednotek **de** Anzahl der Einheiten **da** Antal enheder **es** Cantidad de unidades **et** Ühikute arv **fi** Yksiköiden määrä **fr** Nombre d'unités **el** Αριθμός μονάδων **hr** Broj jedinica **hu** Az egységek darabszáma **it** Numero di unità **lt** Vienetų skaičius **lv** Vienību skaits **nl** Aantal eenheden **no** Antall enheter **pl** Liczba zestawów **pt** Número de unidades

ro Numărul de unități **ru** Количество **sk** Počet jednotiek **sl** Število enot **sv** Antal enheter **tr** Birim sayısı

en Keep away from sunlight **bg** Да не се излага на пряка слънчева светлина **cs** Chraňte před slunečním zářením **de** Vor Sonnenlicht schützen **da** Beskyttes mod sollys **es** Mantener alejado de la luz directa del sol **et** Hoida otsese päikesevalguse eest **fi** Pidettävä pois auringonvalosta **fr** Conserver à l'abri des rayons du soleil **el** Προστατέψτε από την άμεση έκθεση σε ηλιακό φως **hr** Ne izlažite sunčevu svjetlost **hu** Napfénytől távol tartandó **it** Tenere lontano dalla luce diretta del sole **lt** Saugokite nuo saulės šviesos **lv** Sargāt no saules gaismas **nl** Buiten direct zonlicht bewaren **no** Må ikke utsettes for sollys **pl** Chronić przed światłem słonecznym **pt** Manter afastado da luz solar **ro** A nu se expune direct la razele solare **ru** Защищать от воздействия прямых солнечных лучей **sk** Nevystavujte slnečným lúčom **sl** Hranite na temnem **sv** Skyddas mot solljus **tr** Güneş ışığından uzak tutun

en Keep dry **bg** Да се съхранява на сухо **cs** Uchovávejte v suchu **de** Trocken aufbewahren **da** Opbevares tørt **es** Mantener en un lugar seco **et** Hoida kuivaks **fi** Pidettävä kuivana **fr** Conserver au sec **el** Φυλάσσετε σε ξηρό χώρο **hr** Držati na suhom **hu** Szárazon tartandó **it** Tenere asciutto **lt** Laikykite sausi **lv** Uzglabāt sausā vietā **nl** Koel bewaren **no** Hold tørt **pl** Chronić przed wilgocią **pt** Conservar seco **ro** A se păstra într-un loc uscat **ru** Защищать от влаги **sk** Uchovávať v suchu **sl** Hranite na suhem **sv** Förvaras torrt **tr** Kuru yerde saklayın

en Use-by date **bg** Валидно до **cs** Datum spotřeby **de** Verwendbar bis **da** Anvendes inden **es** Fecha de caducidad **et** „Kõlblik kuni“ kuupäev **fi** Viimeinen käyttöpäivä **fr** Date limite d'utilisation **el** Ημερομηνία λήξης **hr** Upotrijebiti do **hu** Felhasználhatóság **it** Data di scadenza **lt** Naudokite iki **lv** Izlietot līdz **nl** Uiterste gebruiksdatum **no** Siste forbruksdato **pl** Termin ważności **pt** Prazo de validade **ro** A se folosi înainte de **ru** Годен до **sk** Dátum spotreby **sl** Rok uporabnosti **sv** Utgångsdatum **tr** Son kullanma tarihi

en Do not use if package is damaged **bg** Не използвайте, ако опаковката е повредена **cs** Nepoužívejte, je-li balení poškozeno **de** Bei beschädigter Verpackung nicht verwenden **da** Må ikke anvendes, hvis pakken er beskadiget **es** No utilizar si observa que el embalaje está dañado **et** Mitte kasutada juhul, kui pakend on vigastatud **fi** Ei saa käyttää, jos pakkaus on vahingoittunut **fr** Ne pas utiliser si l'emballage est endommagé **el** Μη χρησιμοποιείτε το προϊόν εάν η συσκευασία έχει υποστεί ζημιά **hr** Ne upotrebljavajte ako je pakiranje oštećeno **hu** Ne használja, ha a csomagolás sérült **it** Non usare se la confezione è danneggiata **lt** Nenaudokite, jei pakuotė yra pažeista **lv** Nelietot, ja iepakojums ir bojāts **nl** Niet gebruiken indien de verpakking is beschadigd **no** Må ikke brukes dersom pakningen er skadet **pl** Nie należy używać zestawu, jeżeli opakowanie jest uszkodzone **pt** Não utilizar se a embalagem estiver danificada **ro** A nu se folosi dacă ambalajul a fost deschis sau distrus **ru** Не использовать, если упаковка повреждена **sk** Nepoužívať, ak je balenie poškodené **sl** Ne uporabiti, če je embalaža poškodovana **sv** Använd inte produkten om förpackningen är skadad **tr** Ambalajı hasarlıysa kullanmayın

en Single use **bg** За еднократна употреба **cs** Pro jednorázové použití **de** Nur zum Einmalgebrauch **da** Må ikke genbruges **es** No reutilizar **et** Ühekordseks kasutamiseks **fi** Kertakäyttöinen **fr** Utilisation unique **el** Προϊόν μίας χρήσης **hr** Jednokratna uporaba **hu** Nem szabad újrafelhasználni **it** Monouso **lt** Naudokite vieną kartą **lv** Nelietot atkārtoti **nl** Voor eenmalig gebruik **no** Engangsbruk **pl** Wyrób jednorazowego użytku **pt** Utilização única **ro** De unică folosință **ru** Не использовать повторно **sk** Jednorazové použitie **sl** Za enkratno uporabo **sv** Engangsbruk **tr** Tek kullanımlık

en Do not resterilize **bg** Забранено повторно стерилизиране **cs** Nesterilizujte opakovaně **de** Nicht erneut sterilisieren **da** Må ikke reesteriliseres **es** No reesterilizar **et** Mitte reesteriliseerida **fi** Ei saa steriloida uudelleen **fr** Ne pas reesteriliser **el** Μην επαναποστεριώνετε **hr** Nemojte ponovno sterilizirati **hu** Nem szabad újraszterilizálni **it** Non risterrilizzare **lt** Nesterilizuokite pakartotinai **lv** Nesterilizēt atkārtoti **nl** Niet opnieuw steriliseren **no** Må ikke reesteriliseres **pl** Nie sterylizować ponownie **pt** Não reesterilizar **ro** Nu reesterilizați **ru** Не стерилизовать повторно **sk** Nesterilizujte opakovane **sl** Ne sterilizirajte ponovno **sv** Sterilisera inte om produkten **tr** Yeniden sterilize etmeyin

en Sterilized using ethylene oxide **bg** Стерилизирано спомоща на етиленов оксид **cs** Sterilizováno etylenoxidem **de** Sterilisiert mit Ethylenoxid **da** Steriliseret med ethylenoxid **es** Esterilizar mediante óxido de etileno **et** Steriliseeritud etüleenoksiidiga **fi** Steriloitu eteenioksidilla **fr** Stérilisé à l'oxyde d'éthylène **el** Αποστεριρωμένο με οξείδιο του αιθυλενίου **hr** Sterilizirano etilen oksidom **hu** Etilénoxiddal sterilizálva **it** Metodo di sterilizzazione: ossido di etilene **lt** Sterilizuota etileno oksidu **lv** Sterilizēts ar etilēna oksīdu **nl** Gesteriliseerd met gebruik van ethyleenoxide **no** Sterilisert med etylenoksid **pl** Sterylizowano tlénkiem etylenu **pt** Esterilizado em óxido de etileno **ro** Sterilizat cu oxid de etilen **ru** Стерилизовано окисью этилена **sk** Sterilizované etylénoxidom **sl** Sterilizirano s etilenoksidom **sv** Steriliserad med etylenoxid **tr** Etilen oksit kullanılarak sterilize edilmiştir

en Rapid Exchange **bg** Бърза смяна **cs** Rychlá výměna **de** „Rapid Exchange“ **da** Hurtig udskiftning **es** Intercambio rápido **et** „Rapid Exchange“ **fi** Pikavaihdettava **fr** Échange Rapide **el** Ταχεία εναλλαγή **hr** Brza izmjena **hu** Gyors csere **it** Scambio Rapido **lt** „Rapid Exchange“ **lv** Strauja maiņa **nl** Rapid Exchange **no** Hurtigutskifting **pl** Szybka wymiana **pt** Troca rápida **ro** Schimbare rapidă **ru** Катетер быстрой замены **sk** Rychlá výměna **sl** Hitra menjava **sv** Snabbt byte **tr** Hızlı Değiştirme

en Catalogue number **bg** Каталоген номер **cs** Katalogové číslo **de** Artikelnummer **da** Katalognummer **es** Número de catálogo **et** Kataloginumber **fi** Luettelonumero **fr** N° de référence **el** Αριθμός καταλόγου **hr** Kataloški broj **hu** Katalógusszám **it** Numero di catalogo **lt** Katalogo numeris **lv** Kataloga numurs **nl** Catalogusnummer **no** Katalognummer **pl** Numer katalogowy **pt** Número de catálogo **ro** Numărul de catalog **ru** Номер по каталогу **sk** Katalógové číslo **sl** Kataloška številka **sv** Katalognummer **tr** Katalog numarası

LOT en Lot number **bg** Номер на партида **cs** Číslo šarže **de** Fertigungslosnummer **da** Lotnummer **es** Número de lote **et** Partii number **fi** Eränumero **fr** Numéro de lot **el** Αριθμός παρτίδας **hr** Broj serije **hu** Tételszám **it** Numero di lotto **lt** Serijos numeris **lv** Partijas numurs **nl** Lotnummer **no** LOT-nummer **pl** Numer partii **pt** Número de lote **ro** Număr lot **ru** Номер партии **sk** Číslo šarže **sl** Številka serije **sv** Lotnummer **tr** Lot numarası

REC GW en Recommended Guide Wire **bg** Препоръчителен водач **cs** Doporučený vodič **de** Empfohlener Führungsdraht **da** Anbefalet guidewire **es** Alambre guía recomendado **et** Soovitavat juhttraat **fi** Suositeltava ohjauslanka **fr** Fil-guide recommandé **el** Συνιστώμενος οδηγός **hr** Preporučena žica vodilica **hu** Ajánlott vezetődrót **it** Filo guida consigliato **lt** Rekomenduojama pravedėjo viela **lv** Ieteicamā vadītājstīga **nl** Aanbevolen voerdraad **no** Anbefalt ledevaier **pl** Zalecany przewódnik **pt** Fio-guia recomendado **ro** Ghidajul recomandat **ru** Екомендуемый проводник **sk** Doporučený vodič **drót** **sl** Priporočena vodilna žica **sv** Rekommenderad ledare **tr** Önerilen Kılavuz Tel

MIN GC en Minimum Guiding Catheter Size **bg** Минимален размер на водещия катетър **cs** Minimální velikost vodičho katétru **de** Minimale Größe des Führungskatheters **da** Mindste guidekaterstørrelse **es** Tamaño mínimo de catéter guía **et** Minimaalne juhtekateetri suurus **fi** Ohjainkateetri minimikoko **fr** Taille minimum du cathéter-guide **el** Ελάχιστο μέγεθος οδηγού καθετήρα **hr** Minimalna veličina vodećeg katetera **hu** Maximális vezető-katéter mérete **it** Dimensione minima del catetere guida **lt** Mažiausias nukreipiamo kateterio dydis **lv** Virzošā katetra minimālais izmērs **nl** Minimumgrootte geleidekatheter **no** Minste ledetekaterstørrelse **pl** Minimalny rozmiar cewnika prowadzącego **pt** Tamanho mínimo do cateter-guia **ro** Dimensiune minimă cateter de ghidare **ru** Минимальный диаметр направляющего катетера **sk** Minimálna veľkosť vodiaceho katétra **sl** Minimalna velikost vodilnega katetra **sv** Införingskateterns minimistorlek **tr** Minimum Kılavuz Kateter Boyutu

P en Pressure **bg** Налягане **cs** Tlak **de** Druck **da** Tryk **es** Presión **et** Rõhk **fi** Paine **fr** Pression **el** Πίεση **hr** Tlak **hu** Nyomás **it** Pressione **lt** Slėgis **lv** Spiediens **nl** Druk **no** Trykk **pl** Ciśnienie **pt** Pressão **ro** Presiune **ru** Давление **sk** Tlak **sl** Tlak **sv** Tryck **tr** Basınç

NP en Nominal Pressure **bg** Номинално налягане **cs** Jmenovitý tlak **de** Nominaldruck **da** Nominelt tryk **es** Presión nominal **et** Nimirõhk **fi** Nimellispaine **fr** Pression nominale **el** Ονομαστική πίεση **hr** Nominalni tlak **hu** Névleges nyomás **it** Pressione nominale **lt** Nominalus spaudimas **lv** Nominālais spiediens **nl** Nominale druk **no** Nominelt trykk **pl** Ciśnienie nominalne **pt** Pressão nominal **ro** Resiunea nominală **ru** Номинальное давление **sk** Menovitý tlak **sl** Nominalni tlak **sv** Nominellt tryck **tr** Nominal basınç

RBP en Rated Burst Pressure **bg** Необходимо налягане за разширение **cs** Jmenovitý tlak roztržení **de** Nenn-Berstdruck **da** Nominelt sprængningstryk **es** Presión nominal de rotura **et** Hinnanguline lõhkemisrõhk **fi** Nimellinen puhkeamispaine **fr** Pression d'éclatement nominale **el** Ονομαστική πίεση ρήξης **hr** Maksimalni dopušteni tlak **hu** Névleges repedési nyomás **it** Pressione massima di rottura **lt** Nominalus balionėlio sprogimo spaudimas **lv** Aprēķinātais plīšanas spiediens **nl** Nominale barstdruk **no** Nominelt bristetrykk **pl** Nominalne ciśnienie rozerwania **pt** Pressão nominal de rebentamento **ro** Presiunea de spargere estimată **ru** Давление разрыва баллона **sk** Menovitý tlak roztrhnutia **sl** Nominalni razpočni tlak **sv** Nominellt bristningstryck **tr** Nominal patlama basıncı

 **STENT** en Expanded Stent Outer Diameter **bg** Външен диаметър на разгнат стент **cs** Vnější průměr prodlouženého stentu **de** Außendurchmesser des entfalteten Stents **da** Ydre diameter på udvidet stent **es** Diámetro exterior de la endoprótesis expandida **et** Laiendatud stendi välisläbimõõt **fi** Laajennetun stentin ulkohaikaisija **fr** Diamètre externe du stent déployé **el** Εξωτερική διάμετρος εκπτυγμένου στέντ **hr** Vanjski promjer proširenog stenta **hu** Kiterjesztett sztent külső átmérője **it** Diametro esterno dello stent espanso **lt** Išplėsto stento išorinis skersmuo **lv** Izplesta stenta ārējais diametrs **nl** Buitendiameter van geëxpandeerde stent **no** Ytre diameter på utvidet stent **pl** Zewnętrzna średnica rozprężonego stentu **pt** Diâmetro Exterior do Stent Expandido **ro** Diametrul exterior extins al stentului **ru** Наружный диаметр расширенного стента **sk** Vonkajší priemer rozťahnutého stentu **sl** Zunanji premer razširjene žilne opornice **sv** Ytterdiameter för utvidgad stent **tr** Genişletilmiş Stent Diş Çapı

 **STENT** en Labelled Stent Length **bg** Дължина на стента, означена върху етикета **cs** Délka stentu uvedená na štítku **de** Auf dem Etikett angegebene Stentlänge **da** Anført stentlængde **es** Longitud de la endoprótesis que se indica en la etiqueta **et** Stendi pikkus, mis on märgitud sildile **fi** Etikettiin merkitty stentin pituus **fr** Longueur du stent indiqué sur l'étiquette **el** Μήκος επισήμασμένου στέντ **hr** Označena dužina stenta **hu** A címkén feltüntetett sztenthosszúság **it** Lunghezza dello stent riportata in etichetta **lt** Užimė norādītais stenta garums **lv** Etikēte nurodytas stento ilgis **nl** Gelabelde stentlengte **no** Lengde på merket stent **pl** Długość stentu podana na etykiecie **pt** Comprimento do stent marcado **ro** Etiketėje nurodytas stento ilgis **ru** Длина стента, указанная на этикетке **sk** Dĺžka stentu na označení **sl** Oznaka za dolžino stenta **sv** Angiven stentlängd **tr** İşaretli Stent Uzunluğu

 **en** Non-pyrogen **bg** Негоримо **cs** Nepyrogenní **de** Nicht pyrogen **da** Ikke-pyrogen **es** Apirogén **et** Mittepyrogeenne **fi** Ei-pyrogeeninen **fr** Apyrogène **el** Μη πυρογόνο **hr** Nepirogeno **hu** Nem pirogén **it** Non pirogenico **lt** Nepirogeniškas **lv** Nepirogēns **nl** Niet pyrogeen **no** Ikke-pyrogen **pl** Produkt niepyrogeny **pt** Apirogénico **ro** Apirogen **ru** Апирогенно **sk** Nepyrogénne **sl** Nepirogeno **sv** Icke-pyrogen **tr** Apirojendir

 **MR** en MR Conditional **bg** Безопасно за ЯМР при определени условия **cs** Podmíněně použitelné v prostředí MR **de** Bedingt MRT-tauglich **da** MR-betinget **es** Compatibilidad condicional con la RM **et** MR-tingimuslik **fi** Ehdollisesti MR-yhteensopiva **fr** Compatibilité conditionnelle à la résonance magnétique (RM) **el** Υπό συνθήκες συμβατό με MR **hr** Sigurno za uporabu s MR pri točno propisanim uvjetima **hu** MR-kompatibilis **it** Compatibilità RMI condizionata **lt** Santykinai saugus naudoti MR aplinkoje **lv** Saderīgs ar MR **nl** MR-veilig **no** MR-betinget **pl** Produkt bezpieczny w określonych warunkach wykonywania badania MRI **pt** Condicional para RM **ro** Condiționat(ă) de RM **ru** Условно МР-совместимый **sk** Podmienenčne bezpečné v prostredí MR **sl** Slikanje MR dovoljeno pri določenih pogojih **sv** MR-villkorat **tr** MR koşullu

en	English	Instructions for use	5
bg	Български	Инструкции за употреба	7
cs	Česky	Návod na použití	11
de	Deutsch	Gebrauchsanweisung	14
da	Dansk	Brugsanvisning	17
es	Español	Instrucciones de uso	20
et	Eesti	Kasutusjuhised	24
fi	Suomi	Käyttöohjeet	26
fr	Français	Mode d'emploi	29
el	Ελληνικά	Οδηγίες χρήσης	33
hr	Hrvatski	Upute za uporabu	36
hu	Magyar	Használati útmutató	39
it	Italiano	Istruzioni per l'uso	42
lt	Lietuviškai	Naudojimo instrukcijos	46
lv	Latviešu	Lietošanas instrukcijas	48
nl	Nederlands	Gebruiksaanwijzingen	51
no	Norsk	Bruksanvisning	55
pl	Polski	Instrukcja użytkowania	58
pt	Português	Instruções de utilização	61
ro	Română	Instrucțiuni de utilizare	64
ru	Русский	Инструкция по применению	67
sk	Slovenčina	Návod na použitie	71
sl	Slovenščina	Navodila za uporabo	74
sv	Svenska	Bruksanvisning	77
tr	Türkçe	Kullanım talimatları	80

English

Note:

These instructions are intended for all personnel who will handle the „BeGraft Coronary Stent Graft System“ and should be read carefully before using this product!

BeGraft Coronary Stent Graft System

1 Description

The BeGraft Coronary Stent Graft System consists of:

- a balloon-expandable cobalt-chromium (L605) stent, covered with an ePTFE graft material, which is
- pre-mounted on the balloon of a rapid exchange (RX) stent delivery system.

The delivery system is compatible with 0.014" (0.36 mm) guide wires and is available in usable lengths of 143 cm (Ø 2.5-4.0 mm) and 153 cm (Ø 4.5/5.0 mm).

2 How Supplied

Content per Package:

- One "BeGraft Coronary Stent Graft System"
- One Compliance Card
- One Implant Information Card
- One Instruction for Use

Accessories as listed in chapter 8.1 are not included.

Storage: Store in a dry, dark and cool place. Do not use the system after the "use by" date (indicated on the product label).

Non Pyrogenic.

Sterile: Sterilized with ethylene oxide (EO). Do not use if the package is open or damaged.

Note: This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and/or chemical characteristics introduced under conditions of repeated use, cleaning, and/or re-sterilization may compromise the integrity of the design and/or materials, leading to contamination due to narrow gaps and/or spaces and diminished safety and/or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and/or user.

3 Indications

The BeGraft Coronary Stent Graft System is indicated for transluminal implantation in coronary arteries or aorto-coronary bypass grafts for the treatment of:

1. Acute perforation or rupture of coronary arteries.
 2. Aneurysm of coronary arteries or coronary bypass-vein graft.
- Patients being considered for stent graft implantation should be acceptable candidates for coronary artery bypass graft surgery.

4 Contraindications

- Patients with hemorrhagic diathesis or other disorders in which the application of anticoagulant or antiplatelet therapy is contraindicated and patients unable or unwilling to tolerate anticoagulant / antiplatelet therapy (e.g. septic ulceration, cerebrovascular complications).
- Patients with proximal atherosclerosis and/or significant vessel tortuosity in whom guide wire access and adequate guiding catheter support is prohibited.
- Patients with diffuse coronary artery disease resulting in low blood flow in the stented area.
- Patients with ostial lesions.
- Lesions with coronary thrombus formation.
- Patients with unprotected left main coronary vessel disease (i. e. > 50 % stenosis).
- Lesions that cannot be successfully predilated (i. e. < 20 % diameter stenosis reduction, and > 50 % residual diameter stenosis).
- Lesions with threatened or abrupt closure produced during attempted predilation prior to proposed stent implantation.
- Severe stenosis of left main coronary artery.
- Patients with totally obstructed coronary arteries.
- Diffuse lesions.
- Patients with known allergies to the used materials: stent material (L605) and/or to the graft material PTFE and/or delivery system (stainless steel) material.
- Severe allergic reaction to the contrast medium and/or the necessary pharmaceutical treatment (e.g. aspirin) or bleeding complications and patients unable or unwilling to tolerate anticoagulant / antiplatelet therapy.

5 Potential Adverse Effects

Potential adverse effects include, but may not be limited to the following:

- Dissection, perforation, injury, rupture, sudden and total occlusion or other damages to the coronary arteries
- Thrombosis of the coronary arteries
- Acute or subacute closure of stent graft
- Pseudoaneurysm formation
- Unstable angina
- Arterial spasm
- Acute myocardial infarction
- Cardiac dysrhythmia, including ventricular fibrillation
- Hypo / hypertension
- Restenosis of the treated artery
- Distal emboli (air, particle, thromboembolic emboli)
- Hemoptysis
- Hemorrhage or hematoma
- Hemorrhagic stroke
- Hematuria
- Bleeding from anticoagulation/antiplatelet medication
- Side effects caused by the concomitant medication
- Allergic reaction to contrast medium
- Allergic reaction to the implanted stent material (L605) and/or PTFE
- Infections, fever
- Arteriovenous fistula
- Emergency coronary vessel bypass graft surgery or surgical graft removal

- Deformation of stent symmetry on deployment or on re-crossing with another device (IVUS, balloon catheter, etc.)
- Stent graft migration
- Intimal tear
- Tissue reaction, necrosis
- Renal failure
- Death

6 Warnings and Precautions

The medical techniques and procedures described here do not entail all the procedures of PTCA acknowledged in medicine.

- It is explicitly advised that the product is not designed and tested to be used beyond the approved indications. The manufacturer shall not be liable if the product is used beyond the approved indications. In this case, the treatment is at the sole responsibility of the treating physician.
- The BeGraft Coronary Stent Graft System should only be used by physicians, who are well trained in interventional cardiology and who are experts in the techniques of coronary angioplasty and coronary stent implantation, as well as in management of possible complications.
- Stent graft implantation in patients who are not acceptable candidates for coronary artery bypass graft surgery require careful consideration, including possible hemodynamic support during stenting.
- The use of different stent types in the same vessel is not recommended. This may increase the risk of subacute thrombosis.
- During stent graft implantation, a cardiac surgery team should be on stand-by.
- Use aseptic techniques in order to avoid contamination of the system.
- Use the BeGraft Coronary Stent Graft System before the expiry date noted on the product label.
- Do not expose the system to organic solvents (e.g. alcohol).
- Do not use if the package is open or damaged.
- Ensure that labeled stent outer diameter and length of the stent graft are suitable for the intended procedure. The stent outer diameter must not exceed the diameter of the coronary artery proximally and distally to the aneurysm, rupture or perforation.
- If predilation of the lesion is required, do not overdilate the lesion site. The BeGraft Coronary Stent Graft System is not intended for use as primary dilatation device.
- Ensure that the stent graft is positioned between the RO-markers (don't touch) and undamaged after removal of the stent graft protector.
- Ensure that the shaft is not bent nor kinked. Do not attempt to straighten the shaft if it is bent or kinked, as it may break. Instead, prepare a new system.
- Administer appropriate anticoagulation and vasodilation therapy before insertion of the BeGraft Coronary Stent Graft System. Continue the anticoagulation treatment after the implantation according to national guide lines and international standards.
- Only use diluted contrast medium.

- Do not overtighten the hemostatic valve, as lumen constriction may occur, possibly impairing inflation/deflation of the balloon.
- Insert the BeGraft Coronary Stent Graft System into the vessel under fluoroscopic guidance only.
- Never attempt to advance the stent graft system against increased resistance.
- If increased resistance during advancement is detected, stop the procedure and only continue when the reason was found. If this is not possible, then carefully pull back the BeGraft Coronary Stent Graft System together with the guiding catheter as a unit. Do not attempt to retract the BeGraft Coronary Stent Graft System through the guiding catheter as this may cause the stent to dislodge and an embolization could occur.
- Only use an inflation device with a manometer to inflate the balloon and do not use air or any gaseous medium. Do not exceed the Rated Burst Pressure (RBP) mentioned on the label and compliance chart.
- A balloon burst during implantation may cause migration or deformation of the stent graft and problems to remove the catheter. In such a case, retract the catheter carefully, leaving the guide wire in place. If the balloon gets caught in the stent graft, the catheter should be removed immediately by cardiac surgery. A new PTCA catheter should be carefully inserted to reposition and completely dilate the stent graft.
- After stent graft delivery, manipulate the catheter only after the balloon has been fully deflated under vacuum.
- When passing an already implanted stent which is located proximally to the lesion, the BeGraft stent can be deformed or dislodged. In such a case, the correct position and expansion of the stent graft must be ensured.
- Recrossing a partially or fully deployed stent graft with adjunct device must be performed with extreme caution to avoid stent graft damage or deformation.
- The use of laser devices, mechanical atherectomy and hyperthermia are not recommended in the stented region.
- The long-term effectiveness of stent graft implantation has not been fully established.

7 MRI Safety and Compatibility

Non-clinical testing has demonstrated that the “BeGraft Coronary Stent Graft” is MR Conditional. A patient with this device can be safely scanned immediately after implantation in an MR system meeting the following conditions:

- static magnetic field of 1.5 Tesla and 3 Tesla, with
- maximum spatial field gradient of 6,500 G/cm (65 T/m)
- maximum force product of 117,000,000 G²/cm (117 T²/m)
- theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (normal operating mode)

Under the scan conditions defined above, the “BeGraft Coronary Stent Graft” is expected to produce a maximum temperature rise of less than

- **2.0°C** (2 W/kg, 1.5 Tesla) RF-related temperature increase with a background temperature increase of $\approx 1.6^\circ\text{C}$ (2 W/kg, 1.5 Tesla)

- **1.4°C** (2 W/kg, 3 Tesla) RF-related temperature increase with a background temperature increase of $\approx 1.0^\circ\text{C}$ (2 W/kg, 3 Tesla) after 15 minutes of continuous scanning.

The effect of heating in the MRI environment for overlapping stent grafts or stent grafts with fractured struts is not known.

In non-clinical testing, the image artifact caused by the device extends approximately 8.93 mm from the “BeGraft Coronary Stent Graft” when imaged with a spin echo pulse sequence and a 3 Tesla MR system.

8 Procedure

Note: The physician should consult literature on current medical practice of PTCA.

8.1 Accessories

- Femoral or brachial guiding catheter of appropriate size and configuration to match the coronary artery (5-9 F)
- Introducer sheath, compatible with the guiding catheter
- Hemostatic valve
- Sterile heparinized physiological saline
- 60 % contrast medium diluted 50:50 (per volume) with sterile heparinized physiological saline
- 20 ml Luer-Lock syringe (optional)
- Inflation device with manometer
- 3-way stopcock for a minimum pressure of 25 bar
- 0.014" (0.36 mm) guide wire (Attention: a hydrophilic coated guide wire can get stuck in the balloon catheter when the coating is swelling.)
- Guide wire introducer
- Torquer for the guide wire

8.2 Concomitant Medication

This information is a recommendation only. It is not intended to dictate medical practice. The conditions under which stent grafts have been used to date include:

- Administration of aspirin prior to stent graft implantation.
- Administration of heparin and nitro-glycerin during the procedure.
- Administration of antiplatelet therapy alone or in combination with other anticoagulants up to some months following stent graft implantation.

Most experience with PTCA procedures has been gained with anticoagulant/antiplatelet therapy up to now. Administration of antiplatelet therapy and anticoagulants are depending on the patients' condition and medical history. Prior anticoagulant/antiplatelet administration is an essential adjunct to stent graft implantation. Proper heparin treatment, before and after sheath removal, is also essential for the overall success of the procedure. Before, during and after stent implantation, specific blood values of the patient (e.g. prothrombin time PT, partial thromboplastin time PTT, platelet count etc.) must be determined in order to optimize the medication. The selection, duration and dosage of the medication are under the responsibility of the physician treating the patient, as well as the selection of long-term medication after stent graft implantation. Because of an extended endothelialization time it may be mandatory to administer a prolonged anticoagulant/antiplatelet therapy.

8.3 Prior to Use

Carefully inspect the entire product for any damage. For details please see “Warnings and Precautions” (section 6).

8.4 Preparation

1. Prepare the inflation device according to the manufacturer's instructions.
2. Prepare the sterile physiological saline and the diluted contrast medium.
3. Fill a 10 or 20 ml syringe with 4 ml of diluted contrast medium.
4. Carefully remove the stent graft system from the sterile packaging (aseptic technique!). Pull the stent graft system out of the transport coil.
5. In order to avoid loosening of the crimped stent graft from the balloon, hold the protector at its orange marker with one hand and the delivery system with the other hand, then pull the protector off.

Caution:

Ensure that the stent graft is not damaged (e.g. protruding struts) and that it is positioned between the RO-markers (don't touch the implant). Furthermore, the shaft has to be inspected on bends and kinks. For more details please see “Warnings and Precautions” (section 6). Return any damaged product to Bentley InnoMed GmbH.

6. Connect the 3-way stopcock (min. 25 bar) with the luer connector to the balloon catheter.
7. Connect the syringe to the stopcock. Position the syringe vertically with the piston pointing downwards.
8. Create a strong vacuum by pulling back the piston. Maintain the vacuum for 15 - 20 seconds and make sure that no more air bubbles can be observed in the diluted contrast medium.
9. Carefully release the vacuum, close the stopcock and remove the syringe.
10. Repeat this procedure until no more air bubbles can be observed.

Warning:

If it is not possible to create or to maintain a vacuum, check all connections. If you cannot find the cause, do not use this system under any circumstances and return the unit to Bentley InnoMed GmbH.

11. Create another vacuum in the system and close the 3-way stopcock.
12. Remove the air from the tubing and connector of the inflation device.
13. Connect the inflation device to the 3-way stopcock of the balloon catheter.
14. Open the stopcock carefully. The inflation device should be in a vertical position in order to avoid air entering the system.
15. Attach a needle to the syringe; insert it into the distal tip of the catheter and flush the guide wire lumen by applying gentle pressure to the syringe.
16. The BeGraft Coronary Stent Graft System is now ready for use.

8.5 Pre-Dilation

- If pre-dilation with a suitable PTCA dilatation catheter prior to stenting of the lesion is required, do not overdilate the lesion site.
- For the pre-dilatation of the vessel you must be aware of atherosclerotic plaques, which may inhibit advancement of the system.

8.6 Implantation of the Stent Graft

1. Insert a compatible introducer sheath in the femoral or brachial artery using conventional techniques.
2. Insert the guiding catheter into the vessel through the introducer sheath and fix the hemostatic valve to it.
3. Insert the guide wire into the guiding catheter and carefully push it forwards under fluoroscopic control until it is distally to the lesion/defect.
4. After the guide wire has passed the lesion/defect, close the hemostatic valve tightly to prevent displacement of the wire.
5. Advance the BeGraft Coronary Stent Graft System over the guide wire, ensuring that the guide wire exits through the side-port located approximately 25 cm from the catheter tip.
6. Carefully push the BeGraft Coronary Stent Graft System over the guide wire through the guiding catheter to reach the target lesion.
7. Continue the procedure under fluoroscopic control. The radiopaque markers facilitate the positioning of the implant in the area of the perforation, rupture or aneurysm.
8. After positioning, inflate the balloon to Nominal Pressure (NP). Keep this pressure for 10 to 90 seconds. Due to the restricted blood flow, the patient can feel heart pain. Adapt the inflation time to the patient's condition.
9. After deployment of the stent graft, deflate the balloon maintaining negative pressure, allowing adequate time for full deflation. After full deflation, open the stopcock to allow atmospheric pressure in the balloon.
10. The inflation can be repeated two to three times, but do not exceed the Rated Burst Pressure (RBP). Check the implant for full expansion as long as the balloon catheter has not been retracted (in order) to minimize the risk that part of the balloon got caught in the stent struts.

Caution:

- Do not retract the catheter prior to ensuring that the stent graft is implanted correctly.
Re-dilation with another balloon catheter is principally possible, but with increased risk.
- Do not overtighten the hemostatic valve, since this may restrict the flow of contrast medium into the balloon and thereby increase the inflation/deflation time.
- Do not overstretch the stent graft and do not exceed the labeled Rated Burst Pressure (RBP) of the balloon, as this may cause damage to the graft material and/or the vessel. Permanently check the correct expansion of the stent graft by monitoring with high resolution fluoroscopic equipment.
- Expansion of the stent graft should not be undertaken if the stent graft has not been properly positioned in the

area of perforation, rupture or aneurysm. If the position of the stent graft is not optimal, it should be repositioned before inflation.

Note:

- It is necessary that guide wire and/or catheter remain in the treated area until the procedure is completed.
- Different contrast media may affect the inflation/deflation time due to different viscosities.
- Multiple angiographic views and measurements should be taken in order to ensure the correct positioning of the implant.

8.7 Further Dilation of the Stent Graft

- If the primary angiographic result is not optimal, further expansion of the stent graft with the help of a PTCA catheter is possible.
- If the expanded stent graft needs to be recrossed with a guide wire, the distal tip of a flexible guide wire should be formed to a J-shape prior to reinsertion through the stent graft. This technique is recommended so the guidewire will more smoothly pass through the lumen of the stent graft minimizing the risk of getting caught in the stent struts.

8.8 Withdrawal of the Balloon Catheter

1. Open the hemostatic valve.
2. Hold the guide wire and the hemostatic valve firmly in one hand, while grasping the balloon catheter shaft with the other hand. Maintain the guide wire's position in the vessel by holding the wire stationary and begin pulling the balloon catheter carefully out of the guiding catheter.
3. Retract the balloon catheter until the opening for the guide wire lumen is reached. Then proceed in the same way as with an over-the-wire system until the balloon catheter has been totally removed.
4. Close the hemostatic valve.

Caution:

- Before retracting the balloon catheter, make sure that the balloon has been completely deflated.
- Observe the patient regularly during the first 30 minutes after stent placement and perform an angiographic examination for control. If stent graft placement is associated with the onset of thrombus in the stented region, intracoronary administration of a thrombolytic agent is recommended.

Note:

Once used, this product is potentially biohazardous. Handle and dispose of in accordance with current medical practices and applicable local and/or national laws and regulations.

9 Implant Information Card

An Implant Information Card is supplied with this product. This card must be carefully filled out and the peel-off label of the respective implant, indicating the Manufacturer / Product Name / REF# / Lot#, should be applied to the indicated area on the Implant Information Card (up to 4 peel-off labels can be applied). The patient should carry this Implant Information Card with him and show it to medical personnel in case of further treatments.

10 Absence of Defectiveness

Bentley InnoMed GmbH manufactured, sterilised and packaged this product with the due care described in the respective laws and regulations for Medical Devices. Each manufacturing lot is individually tested and inspected before it is released for distribution.

As a result of biological differences in individuals, no product is 100 % efficient under all circumstances. Therefore it is not possible to reliably predict product efficiency during use, nor is it possible to rule out unsatisfactory results. The product shall not be used if damaged during transport, storage or prior to use. Resterilisation or reuse of the product is prohibited. In case any difficulties arise during the use of the products of Bentley InnoMed GmbH, please send an e-mail to complaint@bentley.global. If the safety of users or patients is affected or if they have been harmed, the statutory notification obligations and deadlines shall be observed.

In the event of complaints or product failure, all components (including accessories required for the intervention such as guide wire, sheath, etc.) must be kept and stored. The manufacturer should be informed immediately after this and the relevant components returned to them.

For further information please contact:

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English

Note:

These instructions are intended for all personnel who will handle the “BeGraft Aortic Stent Graft System” and should be read carefully before using this product!

BeGraft Aortic Stent Graft System

1 Description

The “BeGraft Aortic Stent Graft System” consists of

- a balloon-expandable cobalt-chromium (L605) stent, covered with an ePTFE graft material, which is
- pre-mounted on the balloon of an over the wire (OTW) Stent Delivery System.

The “BeGraft Aortic Stent Graft” is available in 3 different stent designs (Ø 12 & 14 mm, Ø 16 & 18 mm and Ø 20 - 24 mm). The Delivery System is compatible with 0.035" (0.89 mm) guide wires and is available in usable lengths of 75 cm and / or 120 cm, depending on the Stent Graft diameter.

Two radiopaque markers located underneath the balloon fluoroscopically identify the position of the “BeGraft Aortic Stent Graft” on the balloon.

2 How Supplied

Contents per Package:

- One “BeGraft Aortic Stent Graft System”
- One yellow Applicator (for safe access of the “BeGraft Aortic Stent Graft System” into the introducer sheath)
- One Compliance Card
- One Implant Information Card
- One Instruction for Use

Accessories as listed in chapter 7.3 are not included.

Storage: Store in a dry, dark and cool place. Do not use the “BeGraft Aortic Stent Graft System” after the “use by” date (indicated on the product label).

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Note: This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and/or re-sterilization may compromise the integrity of the design and/or materials, leading to contamination due to narrow gaps and/or spaces and diminished safety and/or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and/or user.

3 Indications

The “BeGraft Aortic Stent Graft System” is indicated for

- the implantation in the native and/or recurrent coarctation of the aorta (CoA*) on adolescent or adult patients with the following clinical conditions:
 - Stenosis of the aorta resulting in significant anatomic narrowing as determined by angiography or non-invasive imaging, i.e. echocardiography, magnetic resonance imaging (MRI), CT scan;
 - Stenosis of the aorta resulting in hemodynamic alterations, resulting in systolic pressure gradient, systemic hypertension or altered left ventricular function;
 - Stenosis of the aorta where balloon angioplasty is ineffective or contraindicated;
 - Stenosis diameter < 20 % of the adjacent vessel diameter. Stenosis that would present increased risk of vascular damage of disruption.
- and furthermore for restoring and improving the patency of the iliac arteries**.

*) Sizes approved for CoA: Ø 12 - 16 mm in lengths of 19 - 49 mm, Ø 18 and 20 mm in lengths of 27 - 48 mm and Ø 22 and 24 mm in lengths of 37 and 48 mm.

**) Sizes approved for iliac arteries: Ø 12 - 16 mm in lengths of 29 - 59 mm.

4 Contraindications

Contraindications include, but may not be limited to:

1. When PTA is technically not possible (e. g. not feasible to access the lesion with guide wire or balloon catheter)
 2. Patients too small to allow safe delivery of the Stent Graft without compromise to the systemic artery used for delivery
 3. Clinical or biological signs of infection
 4. Active endocarditis
 5. Known allergy to aspirin, other antiplatelet agents, or heparin
 6. Pregnancy
 7. Aneurysms immediately adjacent to the site of Stent Graft implantation
 8. Stenosis distal to the site of Stent Graft implantation
 9. Lesions in or adjacent to essential collateral(s)
 10. Lesions in locations subject to external compression
 11. Heavily calcified lesions resistant to PTA
 12. Patients with diffuse distal disease resulting in poor Stent Graft outflow
 13. Patients with a history of coagulation disorders
 14. Patients with aspirin allergy or bleeding complications and patients unable or unwilling to tolerate anticoagulant / antiplatelet therapy and / or non-responders to anticoagulant / antiplatelet therapy
 15. Fresh thrombus formation
 16. Patients with known hypersensitivity to the stent material (cobalt-chromium (L605)) and / or ePTFE
- Allergic reaction to the implanted Stent Graft materials (cobalt-chromium (L605) and / or ePTFE)
 - Angina pectoris / abdominalis
 - Arterial spasm
 - Arteriovenous fistulae formation
 - Transitory arrhythmia
 - Bleeding or hematoma at the puncture site or from anticoagulation/antiplatelet medication (may require transfusion)
 - Blockage of major collateral artery or arterial branch
 - Cerebrovascular incident
 - Death
 - Distal emboli (air, particle, thromboembolic emboli)
 - Drug reactions with respect to medication
 - Edema
 - Emergency vascular surgery
 - Exacerbation of pre-existing conditions
 - Hematuria
 - Hemoptysis
 - Hemorrhage or hematoma requiring repair
 - Hemorrhagic stroke
 - Hypotension / hypertension
 - Infection, sepsis
 - Inflammation, fever
 - Injury, trauma, dissection, perforation, rupture, sudden and total occlusion of the vessel or other damages to the indicated arteries
 - Intimal flap, tear or trauma
 - Insufficient wall apposition of the Stent Graft
 - Ischemia / infarction of tissue / organ
 - Myocardial infarction
 - Nausea and vomiting
 - Pain
 - Pseudoaneurysm formation
 - Radiation injuries
 - Renal insufficiency / failure
 - Respiratory arrest
 - Restenosis of treated segment / edge of Stent Graft
 - Seizure
 - Sepsis
 - Shock
 - Short-term hemodynamic deterioration
 - Stroke / Transient Ischemic Attack (TIA)
 - Intraluminal thrombus / embolization / occlusion (Stent Graft or other)
 - Tissue reaction / cell necrosis at the site of implant

5 Potential Adverse Events

Potential adverse events include, but may not be limited to the following:

- Abscess
- Acute or subacute closure of Stent Graft
- Allergic reaction to contrast medium

6 Warnings and Precautions

The medical techniques and procedures described below do not entail all the PTA procedures acknowledged in medicine.

- It is explicitly advised that the product is not designed and tested to be used beyond the approved indications. The manufacturer shall not be liable if the product is used beyond the approved indications. In this case, the treatment is at the sole responsibility of the treating physician.
- The long-term effectiveness of the “BeGraft Aortic Stent Graft” has not been fully established.
- Do not use the system after the “use by” date (indicated on the product label).

- Do not expose this device to organic solvents, ionizing radiation or ultraviolet light.
- The “BeGraft Aortic Stent Graft System” should only be used by physicians trained in interventional techniques and percutaneous transluminal angioplasty.
- Assure that the labeled Stent Graft dimensions are suitable for the intended procedure.
- Do not use if the package is open or damaged.
- Employ aseptic techniques during removal of the product from the package and during use.
- Prior to use, please carefully inspect the product to verify that it has not been damaged during shipment (e.g. bends, kinks or other damage) and visually ensure that the crimped “BeGraft Aortic Stent Graft” is positioned between the balloon markers (do not touch the Stent Graft as this could decrease the retention force of the stent on the balloon).
- Administer appropriate anticoagulation / antiplatelet treatment before insertion of the “BeGraft Aortic Stent Graft System”. Continue the anticoagulation / antiplatelet treatment after the implantation according to current versions of national guide lines and international standards.
- For balloon inflation it is recommended to use a pressure monitoring device; optimal is the use of an inflation device with a calibrated manometer.
- Only use diluted contrast medium. Do not use air or any gaseous medium to inflate the balloon.
- The inflation pressure must not exceed the Rated Burst Pressure (RBP) indicated on the product label and compliance card.
- If a Tuohy-Borst valve is used for hemostasis, do not overtighten it after insertion of the “BeGraft Aortic Stent Graft System”, as lumen constriction may occur (possibly impairing inflation/deflation of the balloon).
- If an introducer sheath with a hemostatic valve is used, it is mandatory to use the yellow Applicator (for safe access of the “BeGraft Aortic Stent Graft System” into the introducer sheath).
- Use the recommended minimum introducer sheath size (refer to product label), as the integrity of the stent graft could be impaired (dislodgement of the stent graft and / or the ePTFE graft material), causing an embolism in the worst case.
- Do not attempt to straighten a catheter if the shaft is bent or kinked, as it may break. Instead, prepare a new catheter.
- Do not use the “BeGraft Aortic Delivery System” for post-dilatation of the “BeGraft Aortic Stent Graft”.
- The “BeGraft Aortic Stent Graft System” is intended to perform as a system. The Stent Graft should not be removed for use in conjunction with other dilatation catheters, nor should the Delivery System be used in conjunction with other stents / stent grafts.
- Only insert the “BeGraft Aortic Stent Graft System” under fluoroscopic guidance.
- After Stent Graft deployment, only move the Delivery System when the balloon is fully deflated.

- Prior to the expansion of the “BeGraft Aortic Stent Graft”, utilize fluoroscopy to verify that the Stent Graft has been properly positioned. Do not expand the “BeGraft Aortic Stent Graft” if it is not properly positioned.
- Once partly or fully deployed, the “BeGraft Aortic Stent Graft” cannot be repositioned or removed.
- Under-expansion of the “BeGraft Aortic Stent Graft” may result in migration.
- If the balloon of the Delivery System gets caught in the “BeGraft Aortic Stent Graft” during removal of the Delivery System, it should be removed immediately by surgery.
- A balloon burst during implantation may cause the migration or deformation of the “BeGraft Aortic Stent Graft”. In such a case, the Delivery System should be retracted carefully, leaving the guide wire in place. A new balloon catheter should be carefully inserted over the guide wire and the “BeGraft Aortic Stent Graft” can be expanded completely.
- Stenting across side-branches and the bifurcation will compromise blood flow with potential severe adverse events.
- A stent or stent graft previously implanted proximally to the lesion can be deformed or dislocated by contact, e.g. with the guide wire or the “BeGraft Aortic Delivery System”.
- Re-crossing a partially or fully deployed “BeGraft Aortic Stent Graft” with adjunct devices must be performed with extreme caution to avoid damage or deformation of the “BeGraft Aortic Stent Graft” (stent and cover material).
- The use of laser devices, mechanical atherectomy and hyperthermia are not recommended in the stented region.

6.1 MRI Safety and Compatibility

Non-clinical testing has demonstrated that the “**BeGraft Aortic Stent Graft**” is **MR Conditional**. A patient with this device can be safely scanned immediately after implantation in an MR system meeting the following conditions:

- static magnetic field of 1.5 Tesla and 3 Tesla, with
- maximum spatial field gradient of 3,400 G/cm (34 T/m)
- maximum force product of 62,000,000 G²/cm (62 T²/m)
- theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (normal operating mode)

Under the scan conditions defined above, the “BeGraft Aortic Stent Graft” is expected to produce a maximum temperature rise of less than

- **3.1°C** (2 W/kg, 1.5 Tesla) RF-related temperature increase with a background temperature increase of $\approx 1.3^\circ\text{C}$ (2 W/kg, 1.5 Tesla)
- **2.9°C** (2 W/kg, 3 Tesla) RF-related temperature increase with a background temperature increase of $\approx 0.8^\circ\text{C}$ (2 W/kg, 3 Tesla)

after 15 minutes of continuous scanning.

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In non-clinical testing, the image artifact caused by the device extends approximately 18.5 mm from the “BeGraft Aortic Stent Graft” when imaged with a spin echo pulse sequence and a 3 Tesla MR system.

7 Procedure

7.1 Inspection Procedure Prior To Use

Remove the “BeGraft Aortic Stent Graft System” from the package and carefully open the sterile pouch. Now take out the dispenser coil and slowly withdraw the catheter.

Caution: High resistance could be felt.

Note: A yellow Applicator for safe access of the “BeGraft Aortic Stent Graft System” into the introducer sheath is attached to the Compliance Card (which is attached to the dispenser coil).

Inspect the catheter shaft for bends, kinks, other damage and contamination. Do not use the system if any defects or contamination are noticed.

7.2 Stent Graft Protector Removal Procedure

In order to avoid loosening of the crimped Stent Graft from the balloon, special attention is warranted when the Stent Graft Protector (transparent with yellow distal end) is removed. For removal, pull the preformed flap (located proximally on the protector tube) distally. The protector tube will be torn open and therefore allow safe removal from the crimped Stent Graft.

Caution: Simply pulling on the protector tube may loosen the crimped Stent Graft from the balloon.

After removal of the Stent Graft Protector verify that the Stent Graft is still located between the radiopaque balloon markers. Do not use the system if this is not the case.

7.3 Accessories

- Introducer sheath in recommended minimum size (refer to product label).

Caution: Bench testing showed that at least the following introducer sheaths should NOT be used, as the integrity of the stent graft could be impaired (dislodgement of the stent graft and/or the ePTFE graft material), causing an embolism in the worst case: Medtronic “Sentrant™”, Medtronic “INPUT™ TS”, St. Jude Medical™ “Fast-Cath™”, Teleflex® “Super Arrow-Flex®”

- Sterile syringes
- 1000 u / 500 ml Heparinized Normal Saline (HepNS)
- One 0.035” (0.89 mm) guide wire of appropriate length
- 60 % contrast (e.g. with 300 mg Iodine / ml) diluted 1:1 with normal saline
- Inflation device with pressure monitoring capability (minimum volume: 60 ml)
- One three-way stopcock
- One torque device (if applicable)

7.4 Delivery System Preparation Procedure

1. Attach the syringe with HepNS to the straight guide wire port.
2. Flush until fluid exits at the distal tip.
3. Prepare the inflation device with diluted contrast medium.
4. Attach the inflation device to the stopcock; attach the stopcock to the inflation port of the Delivery System manifold (“Balloon”).
5. With the tip down, orient the Delivery System vertically.
6. Open the stopcock to the Delivery System; pull negative pressure for 30 seconds; release to neutral pressure for contrast fill.

- Caution:** Do not apply positive pressure to the balloon.
- Close the stopcock to the Delivery System; purge the inflation device of all air.
 - Repeat steps 3 through 5 until all air is expelled.
- Warning:** If it is not possible to create or to maintain a vacuum, check the connections. If you cannot find the cause, do not use this "BeGraft Aortic Stent Graft System" under any circumstances and return the unit to Bentley InnoMed GmbH.
- Open the stopcock to the Delivery System.
 - Leave on neutral pressure.

7.5 BeGraft Aortic Stent Graft Delivery Procedure

- Wipe the exposed guide wire with heparinized saline to remove residual blood or contrast medium. Ensure that the guide wire position across the target lesion is maintained during all further procedural steps.
- If a Tuohy-Borst valve is used for hemostasis, fully open the hemostatic valve and backload the "BeGraft Aortic Stent Graft System" onto the guide wire.
- If an introducer sheath with a hemostatic valve is used, it is **mandatory** to use the **yellow Applicator** (for safe access of the "BeGraft Aortic Stent Graft System" into the introducer sheath).
 - Advance the yellow Applicator with the flared end to the back over the guide wire and position it in front of the introducer sheath.
 - Advance the "BeGraft Aortic Stent Graft System" over the guide wire and place the Stent Graft completely inside the yellow Applicator.
 - Push the yellow Applicator with the "BeGraft Aortic Stent Graft System" gently through the hemostatic valve of the introducer sheath.
 - Further advance the "BeGraft Aortic Stent Graft System" into the introducer sheath.
 - Retract the yellow Applicator to achieve hemostasis.
- Maintain neutral pressure on the inflation device.
- Advance the Delivery System to the target lesion. Utilize the radiopaque balloon markers to position the "BeGraft Aortic Stent Graft" across the lesion. Perform angiography to confirm the position of the "BeGraft Aortic Stent Graft". If applicable tighten the hemostatic valve. Do not overtighten it as lumen constriction may occur.

Note: If you notice during the process of moving the Delivery System into position that the "BeGraft Aortic Stent Graft" has moved on the balloon, do not deploy the "BeGraft Aortic Stent Graft".

Caution: Do not attempt to pull an unexpanded Stent Graft back through the introducer sheath. If an unexpanded Stent Graft needs to be removed, the introducer sheath and the Delivery System should be removed as single unit (see "7.7 Delivery System Removal Procedure").

- The "BeGraft Aortic Stent Graft" is now ready to be deployed.

7.6 BeGraft Aortic Stent Graft Deployment Procedure

Caution: Refer to the product label for expanded outer diameter of the Stent Graft, nominal deployment pressure (NP) and rated burst pressure (RBP).

Inflate the Delivery System under fluoroscopic control to the specified nominal deployment pressure (NP) and confirm complete expansion and wall apposition. Do not exceed the RBP of the "BeGraft Aortic Delivery System".

7.7 Delivery System Removal Procedure

- With the guide wire remaining across the lesion, deflate the balloon by pulling negative pressure on the inflation device, remove the deflated balloon from the expanded Stent Graft and pull the Delivery System into the introducer sheath.
- Further withdraw the Delivery System through the introducer sheath.

Caution: Should unusual resistance be felt at any time during removal of the Delivery System post "BeGraft Aortic Stent Graft" implantation, the introducer sheath and Delivery System should be removed as a single unit. In this case, position the proximal balloon marker just distal to the tip of the introducer sheath, then advance the guide wire in the anatomy as far distally as safely possible and secure the Delivery System to the introducer sheath; then remove the introducer sheath and Delivery System as a single unit. Failure to follow these steps and / or applying excessive force to the Delivery System can potentially result in loss or damage of the Delivery System components.

Note: Once used, this product is potentially biohazardous. Handle and dispose of in accordance with current medical practices and applicable local and / or national laws and regulations.

7.8 BeGraft Aortic Post-Dilation

If the primary angiographic result is not optimal, the "BeGraft Aortic Stent Graft" can be post-dilated using a new, generic PTA balloon catheter.

Labelled diameter	Max. post-dilation diameter
Ø 12 and 14 mm	Ø 20 mm
Ø 16 and 18 mm	Ø 24 mm
Ø 20 - 24 mm	Ø 30 mm

Caution: Do not exceed the maximum post-dilation diameter of the "BeGraft Aortic Stent Graft".

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9 References

The physician should consult up-to-date literature on current medical practice of PTA.

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As a result of biological differences in individuals, no product is 100 % efficient under all circumstances. Therefore it

is not possible to reliably predict product efficiency during use, nor is it possible to rule out unsatisfactory results. The product shall not be used if damaged during transport, storage or prior to use. Resterilisation or reuse of the product is prohibited. In case any difficulties arise during the use of the products of Bentley InnoMed GmbH, please send an e-mail to complaint@bentley.global. If the safety of users or patients is affected or if they have been harmed, the statutory notification obligations and deadlines shall be observed. In the event of complaints or product failure, all components (including accessories required for the intervention such as guide wire, sheath, etc.) must be kept and stored. The manufacturer should be informed immediately after this and the relevant components returned to them.

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All rights reserved.
© 2019 Bentley InnoMed GmbH