APTIMAX™

Instrument Trays

- (fr) Plateaux à instruments
- (de) Instrumententrays
- (it) Vassoi per strumenti
- (es) Bandejas para instrumental
- (pt) Tabuleiros de Instrumentos
- (bg) Таби за инструменти
- (cs) Kontejnery na nástroje
- (da) Instrumentbakker
- (el) Δίσκοι εργαλείων
- (et) Instrumendikorvid
- (fi) Instrumenttialtaat
- (hr) Posude za instrumente
- (hu) Műszertálcák
- (id) Baki Instrumen
- (ja) 滅菌トレイ
- (lt) Instrumentų dėklai
- (lv) Instrumentu paplātes
- (mk) Кутии за инструменти
- (nl) Instrumenttrays
- (no) Instrumentbrett
- (pl) Tace na instrumenty
- (ro) Țăvi de instrumentar
- (ru) Лотки для инструментов
- (sk) Kontajnery na nástroje
- (sl) Pladnji za instrumente
- (sv) Instrumenttråg
- (tr) Alet Tepsileri
- (uk) Лотки для інструментів

ADVANCED STERILIZATION PRODUCTS, INC.

33 Technology Drive, Irvine, CA 92618, USA

+1(888)-STERRAD  www.asp.com

© ASP 2021

ASP, The Netherlands BV
BIC 1, 5657 BX, Eindhoven, The Netherlands

English	4
Français	6
Deutsch	8
Italiano	10
Español	12
Português	14
Български	16
Čeština	18
Dansk	20
Ελληνικά	22
Eesti	24
Suomi	26
Hrvatski	28
Magyar	30
Bahasa Indonesia	32
日本語	34
Lietuviškai	36
Latviešu	38
Македонски	40
Nederlands	42
Norsk	44
Polski	46
Română	48
Русский	50
Slovensky	52
Slovenščina	54
Svenska	56
Türkçe	58
Українська	60

APTIMAX™

Instrument Trays

en

INSTRUCTIONS FOR USE

INDICATIONS FOR USE

The APTIMAX Instrument Trays are designed to encase surgical instruments for sterilization in STERRAD™ Sterilization Systems:

- STERRAD™ 100S
- STERRAD NX™ (with or without ALLClear™ Technology)
- STERRAD™ 100NX (with or without ALLClear™ Technology)

GENERAL INFORMATION

1. Cleaning of the trays after the sterilization process is recommended after each use. A mild detergent is recommended.
2. Only STERRAD Instrument Mats and Holders have been validated by Advanced Sterilization Products for use with the APTIMAX Trays. (See product offering below.)

PRECAUTIONS

1. To enable maximum penetration/diffusion, do not stack the APTIMAX Trays in any STERRAD System.
2. Follow the load weight requirements found in your STERRAD System User's Guide
3. Discontinue use of the APTIMAX Trays when they become cracked or damaged in any way.
4. Only legally marketed (FDA cleared in US), polypropylene sterilization wraps should be used with the APTIMAX Trays in the STERRAD System.
5. APTIMAX Instrument Tray 13838 is not designed to fit in the STERRAD NX Sterilization System.

STERILIZATION INFORMATION

1. An appropriate sterile barrier will be required to maintain sterility post-processing. Consult the sterilization wrap manufacturer's instructions for use regarding compatible sterilization methods.
2. Use of the APTIMAX Instrument Trays must be in accordance with the operating instructions supplied with the sterilization equipment and in keeping with your hospital policy on sterilization validation. Several organizations (e.g., AAMI, AORN and ISO) publish guidance information on sterilization processing should you need additional information.
3. Complex instruments should be prepared and sterilized according to the device manufacturer's instructions.

DISPOSAL INSTRUCTIONS

1. Dispose of this product in accordance with the National, Federal, State, and local regulations, or per the healthcare facility's policies and procedures.

APTIMAX™

Instrument Trays

ACCESSORIES FOR APTIMAX INSTRUMENT TRAYS

STERRAD™ Instrument Tray Mats

REF Fits APTIMAX Instrument Tray #

99205	13832
99206	13833
99207	13834
99208	13835
99209	13829
99210	13830
99211	13831
99212	13836
99213	13837
No mat available	13838
15911	15910

APTIMAX™ Instrument Tray Holders

REF Hole Dimensions (millimeters)

99400	10 x 19 x 10
99401	19 x 10 x 19
99402	5 x 19 x 5
99403	19 x 5 x 19
99404	10 x 10 x 10
99405	15 x 15 x 15
99406	5 x 5 x 5
99407	Holding Bar

For a glossary of symbols used, go to eifu.asp.com

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