
C1® TECHNOLOGY GUIDE

READ THIS PRODUCT INSERT THOROUGHLY BEFORE USE

FIG. 1

- | | | |
|--------------|--------------|-------------------|
| ① NEEDLE | ④ DEVICE TIP | ⑦ BLUE LEVER |
| ② NEEDLE CAP | ⑤ SHAFT | ⑧ ROTATIONAL KNOB |
| ③ SUTURE | ⑥ HANDLE | ⑨ INDICATOR FIN |

C1® DEVICE - DESCRIPTION


Each sterile package contains one (1) SINGLE PATIENT USE C1® suturing device available with an angled (FIG. 1) shaft or with a straight (Fig. 2) shaft. The C1® DEVICE is used for the placement of surgical suture as supplied in an appropriate COR-SUTURE® or LS-5™ ePTFE QUICK LOAD® surgical suture. A short length of modified surgical stainless steel tubing, called a needle cap ②, is attached to each end of the suture ③. The needle cap is loaded into the needle cap compartment in the distal end of the device tip ④. Suture Placement and Suture Rearm (FIG. 6) are each achieved by sequentially squeezing and releasing the blue lever ⑦. During Suture Placement, the initial squeeze of the blue lever advances the retracted needle ① forward through the selected tissue placed in jaw of the device tip; the full squeeze advances the needle into the needle cap attached to its suture held in the device tip's distal end. Release of the blue lever retracts the needle, which pulls the now engaged needle cap and suture back through the tissue. Next, with the device tip oriented for Suture Rearm, a second squeeze of the blue lever advances the needle with its now engaged needle cap and suture forward through the empty jaw into the device tip's distal end, where a latch feature retains the needle cap and suture. Release of the blue lever returns the needle alone back to its retracted position in the distal shaft ready for repeat suture placement. The orientation of the device tip is set by rotating the rotational knob ⑧, which has six distinct positions, and an indicator fin ⑨ with an integrated suture management channel.

INDICATIONS

The C1® DEVICE is indicated for use in the approximation of soft tissue and prosthetic materials.

LS SOLUTIONS®

FIG. 2

The C1® DEVICE is available with a 5mm shaft in either an angled or straight configuration.

LOADING SUTURE


CAUTION: To avoid accidental needle exposure, **DO NOT** squeeze lever during suture loading.

FIG. 4

CAUTION: USE ONLY LSI SOLUTIONS COR-SUTURE® OR LS-5™ ePTFE QUICK LOAD® SURGICAL SUTURE WITH C1® DEVICE.

1 REMOVE Suture Tube From Tray

REMOVE suture tube from tray by grasping suture tube at suture release feature and pulling suture tube completely out of tray.

Straighten Tube

2 INSERT Suture Into Suture Track

INSERT one suture end into the suture track as shown; may require pulling needle cap and suture further out of suture tube.


3 PULL Needle Cap Into Compartment

PULL suture to seat needle cap into needle cap compartment in the distal end of the device tip. Make sure needle cap is fully seated in front of latch.

4 FIRE & REARM


Orient Suture As Shown, Squeeze And Release Lever To Ensure Suture Is Ready In Distal End Of Device Tip


FIRE & REARM to ensure suture is loaded properly. To avoid jamming the needle cap into the needle cap compartment, orient the suture directly away from the jaw as shown. Squeeze the lever to advance the needle through the jaw and into the newly loaded needle cap. Release the lever to pick up and retract the needle cap with attached suture back on the needle into the shaft. While continuing to orient the suture as shown, squeeze the lever to advance the needle, needle cap, and suture forward through the jaw to rearm the needle cap into its compartment. Release the lever again and retract back the now-empty needle, leaving the needle cap and suture ready for patient use.

OPTIONAL SUTURE MANAGEMENT CHANNEL

The opposite end of the suture tube can be passed through the fin's suture management channel and subsequently removed to have the suture captured in the fin.

Alternatively, the suture tube can be removed and the opposite end of the suture can be passed directly through the suture management channel.



COR-SUTURE® AND LS-5™ ePTFE QUICK LOAD® SURGICAL SUTURE

FIG. 3

COR-SUTURE® AND LS-5™ ePTFE QUICK LOAD® SURGICAL SUTURE

- ① SUTURE TRAY
- ② SUTURE RELEASE FEATURE
- ③ NEEDLE CAP
- ④ SUTURE
- ⑤ SUTURE TUBE

DESCRIPTION

Each LSI SOLUTIONS® COR-SUTURE® and LS-5™ ePTFE QUICK LOAD® surgical suture is held in a customized tray (1) with suture release feature (2) designed to enable the rapid, easy and reliable loading of suture into C1® devices. COR-SUTURE® QUICK LOAD® surgical suture is available as a non-absorbable monofilament and braided polyester surgical suture (FIG. 10). LS-5™ ePTFE QUICK LOAD® surgical suture is a non-absorbable monofilament expanded polytetrafluoroethylene (ePTFE) suture. A short length of modified surgical stainless steel tubing, called a “needle cap” (3), is attached to each end of the suture (4). The COR-SUTURE® and LS-5™ ePTFE QUICK LOAD® surgical suture also include a detachable clear suture tube (5) to keep the suture from tangling. Each sterile COR-SUTURE® and LS-5™ ePTFE QUICK LOAD® surgical suture is individually packaged for single patient use.

CAUTION:

To avoid accidental needle exposure, **DO NOT** squeeze lever during suture loading.


UNLOADING C1® DEVICE

UNLOADING NEEDLE CAP AND REMAINING SUTURE FROM C1® DEVICE

There is a simple and convenient method for removal of needle caps from the needles prior to reloading the C1® DEVICE. The AUTO-RELEASE Technique is illustrated below. This technique automatically removes the sutures and needle caps from the needle by simply pulling the lever fully forward.

FIG. 5

AUTO-RELEASE Technique *Retract Suture Back. Pull Lever Fully Forward. Remove Suture With Needle Cap.*


This unloading technique requires that the needle with its engaged needle cap and suture must first be retracted back into the distal end of the device shaft. If they are not, simply squeeze and release the blue lever to pick up and retract the needle cap, suture, and needle. Once the needle cap, suture, and needle are retracted into the shaft, pull the lever fully forward to automatically remove the needle cap and suture off of the needle. Inspect and discard the used needle cap and suture.


C1® DEVICE FUNCTION

FIG. 6

SUTURE PLACEMENT

SUTURE REARM

ACTIONS

To facilitate the placement of multiple stitches of the same suture (i.e., “running” the suture), the needle cap rearming mechanism of the *C1® DEVICE* enables the remote return and rearming of the needle cap back into the needle cap compartment. The operator selects the desired device tip orientation by turning the rotational knob and presents an appropriate tissue structure into the gap of the metal jaw in the device tip of the *C1® DEVICE*. During suture placement, the blue lever is squeezed to advance the retracted needle from the shaft of the device through the tissue in the jaw and into the needle cap. The distal end of the *C1® DEVICE* needle engages and captures the needle cap with its attached suture. Releasing the blue lever retracts the needle with attached needle cap and suture back through the tissue. For Suture Rearm, the device tip is then moved away from any tissue structures to clear the jaw for needle cap rearming. The blue lever is again fully squeezed to advance the needle, needle cap and suture forward through the empty jaw toward the needle cap compartment at the most distal end of the device tip. Care must be taken to ensure the suture crimped to the distal end of the needle cap is oriented to permit the suture to freely pass through the suture track opening in the needle cap compartment. Jamming misoriented suture into the needle cap compartment can damage suture and device. With the blue lever fully squeezed and the needle with its attached needle cap and suture completely forward, the now rotated needle will permit the needle cap rearm latch to engage the face of the needle cap and cause its release from the needle. With the first tissue suture placement complete, the needle cap rearmed into the needle cap compartment, and the blue lever back in its starting location, the *C1® DEVICE* is ready for another tissue bite. This sequence can be repeated for up to 24 bites (or 24 complete functional cycles) with up to 4 sutures per device.

CONTRAINDICATIONS

- Minimally invasive surgical procedures should only be performed by physicians having adequate training and familiarity with minimally invasive techniques. In addition, medical literature should be consulted relative to techniques, complications and hazards prior to the performance of minimally invasive procedures.
- The *C1® DEVICE* is not intended to be used with any suture other than a *COR-SUTURE®* or *LS-5™ QUICK LOAD®* surgical suture.
- Do not use this suture under conditions in which excessive suture tension can lead to tissue damage.

WARNINGS

- Federal (U.S.A.) law restricts this device to sale, distribution and use by, or on, the order of a physician.
- Do not resterilize.** The performance of the *C1® DEVICE* after cleaning or other reprocessing has not been verified and is not supported by LSI SOLUTIONS®.
- Discard open (unsealed), unused, expired or damaged devices or devices in damaged primary packaging.
- The *C1® DEVICE* will not fit through standard laparoscopic cannulas with seals.
- As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.
- Users should be familiar with surgical procedures and techniques involving suture before employing the *C1® DEVICE* for wound closure, as the risk of wound dehiscence may vary with the site of application.
- Acceptable surgical practice and good surgical judgement must be followed with respect to drainage and closure of infected or contaminated wounds.
- Redundant, cut-away suture remnants, used needle caps, and *C1® DEVICES*, along with packaging, must be inspected, handled and disposed of consistent with standard, accepted medical device disposal procedures.
- The *C1® DEVICE* is not for use in neurological procedures.
- Applications other than for soft tissue closure, or to anchor prosthetic materials, can result in failure to pick up suture or in damage to the device making it unsuitable for continued use.
- Never drive the needle into suture, bone, dense ligamentous tissue, severely calcific tissue or other instruments.
- Do not leave any foreign material (e.g. suture remnant, needle cap, etc.) unattached in areas potentially exposed to circulating blood.
- If securing the *COR-SUTURE®* or *LS-5™ QUICK LOAD®* surgical suture with any LSI SOLUTIONS® *COR-KNOT® DEVICE (COR-KNOT® MIS DEVICE, COR-KNOT MINI® DEVICE)*, or *MI-KNOT™ DEVICE* ensure any needle caps are removed from the suture ends being loaded, prior to loading the suture through the *COR-KNOT® DEVICE*.

PRECAUTIONS

- Check for hemostasis or leakage where appropriate.
- Always assure device tip is viewed under direct visualization before advancing the needle.
- Ensure obstructions do not interfere with the movement of the needle of *C1®*.
- Do not squeeze the blue lever of the *C1® DEVICE* while loading a surgical suture; squeezing the blue lever may expose the sharp needle and/or damage the needle.
- In handling the *C1® DEVICE*, care should be taken to avoid jamming the suture into the needle cap compartment and damage to the needle.
- Avoid damage to the needle, suture or needle caps due to direct application of surgical instruments, like forceps, needle holders, clamps, etc.
- Adequate knot security requires accurate completion of accepted surgical techniques for constructing surgically tied knots or the use of the *COR-KNOT® MIS DEVICE, COR-KNOT MINI® DEVICE*, and *COR-KNOT® QUICK LOAD® UNIT* or *INTRA-KNOT® QUICK LOAD® UNIT* or *MI-KNOT™ DEVICE* as warranted by surgical circumstances and the experience of the surgeon.
- Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure that mechanical function and electrical isolation and grounding are not compromised.
- Before loading the *C1® DEVICE* with another surgical suture, assure the remaining suture tail and needle cap from the previous load has been removed from the needle and needle cap compartment. Failure to appropriately remove used needle caps from the needle can result in damage to the device, including intracorporeal or extracorporeal fracturing off of the tip of the needle, making it unsuitable for continued use.
- After each loading and re-loading of a new suture into this device, squeeze the blue lever to drive the needle forward into the new needle cap loaded into the needle cap compartment. If the needle cap is picked-up by the needle, then squeeze the lever again to rearm the needle cap and suture. If the needle rotation is oriented to rearm the needle cap, then the needle will retract back without the needle cap and suture attached. This “cycling” (FIG. 4, Step 5) of the needle helps ensure that the previous needle cap was properly removed, the new needle cap is installed properly, and the operator receives the device with its needle oriented to pick-up the needle cap on its first needle advancement. If the previous needle cap was not properly removed from the needle prior to reloading the device, the needle will not fully advance into the new needle cap in the needle cap compartment. Failure to avoid driving a needle with a needle cap into another needle cap can lead to the breaking off of the tip of the needle.
- Do not use the *C1® DEVICE* to dissect or aggressively manipulate tissue structures.
- Verify that the needle cap is still retained within the needle cap compartment and the device has not been damaged or deformed before attempting to place a stitch.
- Do not manipulate the device at any time with the blue lever partially actuated. This may expose sharp surfaces that can cause trauma to the patient, the device operator or other staff, or damage the device.
- To avoid inadvertent suture damage, ensure the needle cap always enters the needle cap compartment with its suture oriented to freely pass through the needle cap compartment's suture track. Do not use damaged suture or expired suture.
- Ensure the advancing needle targets and enters the needle cap compartment. For example, during the suture placement, avoid using an extended needle to manipulate or lift tissue because such an action can cause the needle to deviate from its targeted course toward the needle cap compartment. A needle tip, not entering the needle cap compartment properly, can strike the distal tip of the device and lead to undesired outcomes, including needle tip fracture. For another example, during suture rearm, avoid applying tension to the suture from the needle cap on the needle. Tension on the suture can cause the needle to deviate off target and lead to the needle cap possibly striking the distal tip, which can cause needle tip fracture.
- Turning the rotational knob while squeezing the blue lever may damage tissue, and/or the *C1® DEVICE* and suture.

ADVERSE REACTIONS

Adverse effects associated with the use of suture include wound dehiscence, failure of adequate wound support in closure sites where expansion, stretching or distension occur, enhanced bacterial infectivity, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and pain, edema and erythema.

C1® TECHNIQUE PEARLS

If device fails to place a stitch, then under direct visualization and away from tissue, fire and rearm once or twice per Fig. 4, Step 4, then continue if appropriate. If suture placement is ineffective, try reloading with a new suture. If the device is still ineffective, replace device.

FIG. 10 C1® PRODUCT ORDERING

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The LSI logo, LSI SOLUTIONS, C1, Cor-Suture, Quick Load, LS-5, Cor-Knot, Cor-Knot Mini, Mi-Knot, Intra-Knot and Perfect Performance Policy are trademarks and registered trademarks of LSI SOLUTIONS, Inc.

Copyright © 2014, LSI SOLUTIONS®. All Rights Reserved.

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LSI SOLUTIONS®

7796 Victor-Mendon Road
Victor, New York 14564 U.S.A.
Phone: 585.869.6600
Customer Service: 866.575.3493
Technical Support: 866.428.9092
Fax: 585.742.8086
www.lsisolutions.com

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LSI SOLUTIONS®
7796 Victor-Mendon Road Victor,
New York 14564 U.S.A. Phone:
585.869.6600
Customer Service: 866.575.3493
Technical Support: 866.428.9092
Fax: 585.742.8086
www.lsisolutions.com


MADE IN THE USA

This Product Comes
with our LSI SOLUTIONS®
Perfect Performance Policy®
Call us at 866.575.3493 any time.

COR-KNOT MICRO® TECHNOLOGY GUIDE

 READ PRODUCT INSERT THOROUGHLY BEFORE USE

FIG. 1


COR-KNOT MICRO® DEVICE DESCRIPTION

Each sterile COR-KNOT MICRO® DEVICE package contains one single-patient-use COR-KNOT MICRO® DEVICE preloaded with a single COR-KNOT MICRO® FASTENER (1). Made from medical-grade titanium, the COR-KNOT MICRO® FASTENER is a hollow sleeve with a rounded base. A white target (2) (shown removed above) holds the loop shape of a wire snare (3). The wire snare passes through the COR-KNOT MICRO® FASTENER and is attached to a snare puller (4) knob. A suture slot (5) (not shown) in the device shaft (6) lies under the opening in the snare puller. The ends of a polypropylene suture (USP 6-0, 7-0, or 8-0) are passed through the wire snare and subsequently threaded into the titanium fastener. The snare puller is pulled up or retracted along the device shaft until it snaps onto the puller retainer (7) feature of the purple knob (8), which also has an integrated indicator fin (9). The suture slot and the indicator fin are located on the same side of the device shaft. The subsequently crimped fastener and remnant trimmed suture tails bend slightly in the direction away from or opposite the suture slot and indicator fin. By rotating the purple knob and the device's white handle (10), the surgeon can ergonomically orient the direction of the suture tails, if desired. A yellow lever stop (11) is located behind the purple lever (12) to restrict inadvertent squeezing of the lever during device handling before crimping. The lever stop is removed by pinching its sides together and pulling it out of the handle. By squeezing the purple lever, the COR-KNOT MICRO® DEVICE crimps the COR-KNOT MICRO® FASTENER to fasten together segments of suture and trims away excess suture.

INDICATIONS

The COR-KNOT MICRO® DEVICE with COR-KNOT MICRO® FASTENER when used in conjunction with USP 6-0, 7-0, or 8-0 polypropylene surgical suture, is indicated for use in the approximation of soft tissue.


FIG. 2


 **SOLUTIONS®**

PREPARATION

1 ENSURE


2 ROTATE


3 PUSH-OUT

Target Button

4 UNLOCK

5 INSPECT

PREPARING A COR-KNOT MICRO® DEVICE TO RECEIVE SUTURE

NOTE: The COR-KNOT MICRO® DEVICE does **NOT** require intraoperative loading of a titanium fastener. It is provided preloaded with a titanium fastener already positioned in the distal tip of the device.


Use proper operating room technique to pass the sterile COR-KNOT MICRO® DEVICE from its packaging. While maintaining appropriate sterile technique, follow the steps indicated in the illustrations. Apply appropriate hospital policy and practices to dispose of the components of this product.

1. **ENSURE** that both ends of the suture are approximately of equal length and attachments on the suture ends (e.g., needles, needle caps or ferrules, etc.) are cut away and/or removed, and no knots are in the sutures to be snared.
2. **ROTATE** the shaft's purple knob so that its indicator fin and suture slot are oriented approximately opposite the preferred direction of the trimmed suture tails.
3. **PUSH-OUT** and remove the white target from the wire snare by pressing against the target button.
4. **UNLOCK** the COR-KNOT MICRO® DEVICE by pinching the sides of the yellow lever stop and pulling it away from the handle.
5. **INSPECT** to ensure that the COR-KNOT MICRO® FASTENER is fully seated. If the fastener is not loaded properly, discard the COR-KNOT MICRO® DEVICE and obtain a new device.

FIG. 3

IMPORTANT : PULL KNOB ONLY AFTER SUTURE IS SNARED

DO NOT MOVE


NO

✗

Do not move snare puller knob before pushing-out the white target or passing suture.


NEVER ROTATE

NO


✗

Direct rotation of snare puller knob can break the wire snare.


NEAR SURGICAL SITE


1 LOWER


2 PASS

3 RETRACT/ SNAP

4 GRASP

5 INSPECT

THREADING SUTURE IN PRELOADED COR-KNOT MICRO® FASTENER

NOTE: Adjacent to the wound: Surgeons can use their dominant hand to hold the device's white handle and their nondominant hand to complete the suture threading technique.

1. **LOWER** the COR-KNOT MICRO® DEVICE tip near a sterile surface adjacent to the surgical site to avoid suture falling out of the wire snare.
2. **PASS** approximately 1in (2-3cm) of both ends of the suture through the open wire snare using fingers or forceps to draw suture through. Avoid tensioning suture at the closure site.
3. **RETRACT** the snare puller knob up the shaft while slightly lowering the device tip to thread the suture through the COR-KNOT MICRO® DEVICE and out of the suture slot. If desired, move snare puller knob slightly back down the shaft to free suture away from the shaft. **SNAP** the snare puller knob securely on the snare puller retainer. When fully engaged, the snare puller knob will not slide down the shaft.
4. **GRASP** both suture strands and assure both ends have exited the suture slot. If both suture ends are not threaded through the suture slot, do not use the COR-KNOT MICRO® DEVICE.
5. **INSPECT** to ensure that the COR-KNOT MICRO® FASTENER is fully seated in the distal device tip and retains 2 strands of suture segments. If fastener is not fully seated, discard device and replace.
6. **SLIDE** the COR-KNOT MICRO® DEVICE distal tip gently over lightly tensioned suture down to the targeted site. If desired, rotate the white handle to **ORIENT** the indicator fin and suture slot in a direction opposite subsequently trimmed suture tails.


NOTE: The orientation of the resultant angled crimped COR-KNOT MICRO® fastener is opposite the side of the indicator fin and the suture slot.

AT CLOSURE SITE

6 SLIDE & ORIENT

NOTE: Use one continuous motion when sliding the device to the targeted site. Failure to do so, or sliding with an in and out motion, could dislodge and unseat the fastener.

NOTE: The device tip and fastener should be stabilized and flush against the targeted tissue. Avoid over-compression of the tissue.

7 SQUEEZE (lever) & HOLD (1 second)

Fully squeeze **ONCE**.

8 TUG (suture) & RELEASE (lever)

Without moving the device tip, tug suture to trim. Provide full lever return.

9 PAUSE (1 second) & REMOVE (device)


Carefully and slowly withdraw device. If necessary, for easy release, gently rotate 45° in both directions about the shaft.

SQUEEZE & HOLD / TUG & RELEASE / PAUSE & REMOVE / INSPECT

With the distal tip on the targeted site, apply appropriate suture tension and:

7. **SQUEEZE** the purple lever once until it stops, maintain the device tip's position; **HOLD** lever for 1 second. **DO NOT RELEASE** or squeeze lever again. Fastener is now crimped. While the lever is being held, any additional tension will not be transferred to the tissue.
8. **TUG** the suture gently, without moving the device tip, to cut free both suture ends; fully **RELEASE** the purple lever.
9. **PAUSE** for 1 second; **REMOVE** device by slowly and carefully lifting up. If the crimped fastener does not readily release from the distal tip, ensure that purple lever is released, then gently rotate the white handle 45° in either direction about the shaft. If still necessary, turn the white handle back, then rotate 45° in the opposite direction. If the COR-KNOT MICRO® FASTENER still will not release, cut suture and ensure that fastener and any associated suture remnants are removed from patient.
10. **INSPECT** to ensure that the COR-KNOT MICRO® FASTENER and suture tails are complete and in an appropriate location.

10 INSPECT

IMPORTANT: “ONLY THE SURGEON TOUCHES THE PURPLE LEVER”**FIG. 4****DO NOT SQUEEZE BEFORE THE SNARE PULLER KNOB IS SNAPPED****SQUEEZE THE PURPLE LEVER ONLY ONCE!**

FIG. 4 Snapping the snare puller knob onto the snare puller retainer ensures that the wire snare and suture are completely through the titanium fastener and device tip. If the snare puller with its wire snare is not retracted fully up the device shaft and snapped on the snare puller retainer, the wire snare can remain inside of the titanium fastener and device tip, potentially leading to suture, titanium fastener, and cutting blade damage. When ready to use, squeeze the purple lever once only, hold during suture trimming, and then release the purple lever. After releasing the lever, pause for 1 second before removing the device.

FIG. 5

FIG. 5 Front and side enlarged views of a correctly crimped fastener


FIG. 6

FIG. 6 For comparison, front and side enlarged views of an incorrectly crimped fastener that was not fully seated in the device tip

CONTRAINDICATIONS

- The *COR-KNOT MICRO® DEVICE* is contraindicated for use in ophthalmic and neurologic surgery.
- The *COR-KNOT MICRO® DEVICE* is not intended to be used with any fastener other than the preloaded *COR-KNOT MICRO® FASTENER*.
- The *COR-KNOT MICRO® DEVICE* is not intended to be fired more than one (1) time.
- The *COR-KNOT MICRO® FASTENER* is not marketed for placement into circulating blood.

WARNINGS

- Federal law restricts this device to sale, distribution, and use by, or on, the order of a physician.
- Users should be familiar with standard procedures and techniques involving surgical suture and titanium usage before employing the *COR-KNOT MICRO® DEVICE* for fastening and trimming suture.
- Each *COR-KNOT MICRO® DEVICE* is intended to place a single *COR-KNOT MICRO® FASTENER*.
- If the yellow lever stop is not in place when removing the *COR-KNOT MICRO® DEVICE* from the package, discard the device.
- Do not resterilize. The *COR-KNOT MICRO® DEVICE* is designed and intended for single-patient use only. Do not reuse, reload, reprocess, or resterilize this product. The performance of the *COR-KNOT MICRO® DEVICE* after cleaning or other reprocessing has not been validated and is not supported by LSI SOLUTIONS®. Reuse, reprocessing, or resterilization may compromise the integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Discard any open (unsealed), unused, expired, or damaged *COR-KNOT MICRO® DEVICE*.
- Applications other than for soft tissue approximation can result in damage to the *COR-KNOT MICRO® DEVICE*, rendering it unsuitable for continued use.
- Do not leave any foreign material (e.g., suture fragment or fastener) unattached in areas potentially exposed to circulating blood.
- When securing suture with a *COR-KNOT MICRO® DEVICE*, ensure that attachments on the suture ends (e.g., needles, needle caps or ferrules, etc.) are removed prior to loading the suture through the *COR-KNOT MICRO® DEVICE*.
- Do not squeeze the lever of the *COR-KNOT MICRO® DEVICE* until the *COR-KNOT MICRO® FASTENER* has been appropriately positioned and the suture accurately tensioned at the closure site.
- Excessive suture tensioning can cause suture breakage or tissue deformation or necrosis.
- Direct contact between sensitive tissue structures and foreign materials can lead to tissue injury or damage, such as tissue erosion. Always orient *COR-KNOT MICRO® FASTENERS* and remnant suture tails to avoid direct contact with delicate tissue.
- As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.
- Ensure that the purple lever is squeezed until it stops when actuating the *COR-KNOT MICRO® DEVICE*. Failing to complete a full squeeze on a *COR-KNOT MICRO® DEVICE* can result in a partial crimp in the titanium fastener, which could lead to a reduced holding strength. Holding the lever squeezed for a second ensures full actuation.
- Adequate *COR-KNOT MICRO® FASTENER* security requires reasonable clinical judgment and appropriate surgical techniques as warranted by individual anatomy, surgical circumstances, and the experience of the surgeon.
- Each *COR-KNOT MICRO® DEVICE*, along with packaging, must be inspected, handled, and disposed of consistent with standard, accepted medical device disposal procedures.
- While the titanium of the *COR-KNOT MICRO® FASTENER* is physiologically very inert, routine surgical precautions must be employed whenever foreign materials are left in a patient.

PRECAUTIONS

- When handling the *COR-KNOT MICRO® DEVICE*, care should be taken to avoid damage.
- Surgical procedures should only be performed by physicians having adequate training and familiarity with such techniques. In addition, medical literature should be consulted relative to techniques, complications, and hazards prior to the performance of surgical procedures. Inappropriate actions or use can lead to patient harm or death.
- Avoid damage to the *COR-KNOT MICRO® DEVICE* due to inappropriate squeezing of the purple lever and/or due to application of surgical instruments like forceps, needle holders, clamps, etc.
- Do not use the *COR-KNOT MICRO® DEVICE* to aggressively manipulate tissue structures.
- Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure that electrical isolation or grounding are not compromised.
- Ensure that inappropriate obstructions do not interfere with the firing of the *COR-KNOT MICRO® DEVICE*; obstructions may cause damage or breakage.
- If the *COR-KNOT MICRO® FASTENER* falls out of the device tip or is not properly loaded or seated, discard the *COR-KNOT MICRO® DEVICE* and *COR-KNOT MICRO® FASTENER*.
- Irreparable damage to the *COR-KNOT MICRO® DEVICE* and fastener can occur if the purple lever is squeezed while the wire snare is in place at the tip of the device.
- Do not squeeze the purple lever on the same *COR-KNOT MICRO® DEVICE* more than once.
- If the purple lever of the *COR-KNOT MICRO® DEVICE* does not return completely forward on its own (i.e., without assistance), manually push the purple lever forward all the way until it stops to provide full lever release.
- Trim suture ends with scissors if the *COR-KNOT MICRO® DEVICE* makes a successful crimp but does not cut suture.
- If *COR-KNOT MICRO® DEVICE* makes an unsuccessful crimp, manually cut suture to remove fastener and suture.
- Check for hemostasis or leakage where appropriate.
- Avoid crushing or crimping damage to the *COR-KNOT MICRO® FASTENER* due to inappropriate squeezing of the *COR-KNOT MICRO® DEVICE* purple lever and/or due to application of surgical instruments like forceps, needle holders, clamps, etc.
- Inspect each *COR-KNOT MICRO® FASTENER* and its suture tails.

ACTIONS

When the *COR-KNOT MICRO® DEVICE* preloaded with a *COR-KNOT MICRO® FASTENER* is appropriately positioned at a suture closure site, squeezing the purple lever can instantly secure the fastener and trim the suture. LSI SOLUTIONS® has demonstrated that the *COR-KNOT MICRO® DEVICE* and *COR-KNOT MICRO® FASTENER* meet USP minimum values for tensile strength when used within the average range of USP diameters for 6-0, 7-0, and 8-0 polypropylene suture that complies with USP standards for surgical suture. The surgical titanium used in a *COR-KNOT MICRO® FASTENER* is not absorbed by the body; surgical titanium is generally not associated with significant inflammatory reactions.

ADVERSE REACTIONS

Adverse effects associated with the use of surgical suture and titanium can include, but are not limited to: wound dehiscence, thrombus formation, embolism, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infected wounds, minimal acute inflammatory tissue reaction, bleeding, and transitory local irritation. Surgical titanium is not absorbed by the body and is generally not associated with significant inflammatory reactions. Direct contact between sensitive tissue structures and foreign materials can lead to tissue injury or damage, such as tissue erosion.


MRI Testing

Based on MRI testing information, a titanium *COR-KNOT MICRO® FASTENER* will not present an additional hazard or risk to a patient undergoing an MRI procedure using a scanner operating with a static magnetic field of 3 Tesla or less and under the MRI-related heating conditions (MRI for 15 minutes at an MR system reported whole body averaged specific absorption rate (SAR) value of 2W/kg).

NOTE: The average weight of *COR-KNOT MICRO® FASTENER* is 0.002g; 1/6 the weight of a *COR-KNOT® FASTENER*.

MRI Safety Information – MR Conditional  Nonclinical testing demonstrated that the titanium *COR-KNOT MICRO® FASTENER* is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions: static magnetic field of 1.5 Tesla and 3 Tesla only; maximum spatial gradient magnetic field of 4,000 Gauss/cm (40T/m); and maximum MR system reported whole body averaged specific absorption rate (SAR) of 2W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode. Under the scan conditions defined, the titanium *COR-KNOT MICRO® FASTENER* is expected to produce a maximum temperature rise of 1.5°C after 15 minutes of continuous scanning (i.e., per pulse sequence). In nonclinical testing, the image artifact caused by the titanium *COR-KNOT MICRO® FASTENER* extends approximately 2mm from this implant when imaged using a gradient echo pulse sequence and a 3 Tesla MR system. Reference: ASTM 2503, *Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment*

COR-KNOT MICRO® FASTENER Knot Strength Relative to USP			
USP Size	USP Avg. Diameter Range for Polypropylene Suture	USP Knot Pull Tensile Strength (minimum)	<i>COR-KNOT MICRO® FASTENER</i> Knot Pull Tensile Strength Relative to USP Standard
6-0	0.070–0.099mm 0.0028–0.0039in	0.20kgf	>100%
7-0	0.050–0.069mm 0.0020–0.0027in	0.11kgf	>100%
8-0	0.040–0.049mm 0.0016–0.0019in	0.06kgf	>100%

COR-KNOT MICRO® PRODUCT ORDERING			SUPPLIED: STERILE
	REORDER	PRODUCT	DESCRIPTION
	REF 032500	COR-KNOT MICRO® DEVICE Pre-Loaded with (1) titanium COR-KNOT MICRO® FASTENER	Box of 6 Kits (1 Device per Kit)

The LSI logo, LSI SOLUTIONS, Cor-Knot, Cor-Knot Micro, and Perfect Performance Policy are trademarks and registered trademarks of LSI Solutions, Inc. Copyright © 2020, LSI SOLUTIONS®. All Rights Reserved.

Patents: www.lsisolutions.com/patents

Symbol Glossary: www.lsisolutions.com/symbols

LSI SOLUTIONS®
7796 Victor-Mendon Road
Victor, New York 14564 U.S.A.
Phone: 585.869.6600
Customer Service: 866.575.3493
Technical Support: 866.428.9092
Fax: 585.742.8086
www.lsisolutions.com

MADE IN THE USA

This Product Comes
with our LSI SOLUTIONS®
Perfect Performance Policy®
Call us at 866.575.3493 any time.

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
P/N 032477D 27 MAY 2022

miniRUMEL™ DEVICE TECHNOLOGY GUIDE

 READ THIS PRODUCT INSERT THOROUGHLY BEFORE USE

FIG. 1

miniRUMEL™ PRODUCT DESCRIPTION

Each shelf box provides 12 individually packaged SINGLE PATIENT USE sterile pouches with white, 15 cm miniRUMEL™ devices or 12 green, 25 cm miniRUMEL™ devices (FIG. 1).

miniRUMEL™ DEVICE DESCRIPTION

A miniRUMEL™ DEVICE has a white (15 cm) or green (25 cm) plastic target ① that is held by a wire snare ② connected to a twisted wire ③ that passes through the miniRUMEL™ tube ④. The tube and twisted wire run through the bottom base ⑤ and toggle top ⑥ that comprise the suture clamp ⑦. After suture has been placed through the wire snare, the snare knob ⑧ is pulled out away from the suture clamp to thread the suture in the wire snare through the miniRUMEL™ tube and suture clamp. After establishing the desired suture tension, the back close zone ⑨ can be pressed down to securely clamp the suture by engaging the features in the internal suture gripping area ⑩ (not shown here; see FIG. 2) as the toggle top rotates about the toggle joint ⑪, setting the miniRUMEL™ into its clamped, closed position. While in this position, the suture is held relative to the miniRUMEL™ tube and suture clamp. The forward release zone ⑫ is pressed down after pushing both twin tabs ⑬ forward to rotate the toggle top about the toggle joint, setting the miniRUMEL™ into its unclamped, open position. While in this position, the internal suture gripping area features are no longer holding the suture, so that it may be adjusted as necessary. The suture clamp can be repeatedly opened or closed as desired. The finger wall ⑭ on the bottom base helps stabilize the suture clamp during opening and closing.

INDICATIONS FOR USE


The miniRUMEL™ DEVICE is indicated for use to temporarily hold 2-0 polypropylene, polyester or silk suture, or 3-0 polypropylene suture during surgery.

LOADING


NOTE: Ensure that attachments on the suture ends (e.g., needles, suture caps or ferrules, etc.) are removed prior to threading suture through the *miniRUMEL*™ *DEVICE*.

NOTE: Remove device packaging and components from the surgical field when no longer in use and dispose per hospital policy.


1 **PUSH OUT** and remove the white or green plastic target from the wire snare of the open *miniRUMEL*™ *DEVICE*.


2 **LOWER** the *miniRUMEL*™ wire snare near a sterile surface adjacent to the surgical site to avoid suture falling out of the wire snare.

3 **PASS** about 1–2 in (2–5 cm) of both ends of the suture through the wire snare.


4 **PULL** the snare knob away from the suture clamp while slightly lowering the device to thread the suture through the *miniRUMEL*™ tube and the suture clamp. If both suture ends are not pulled through the suture clamp, do not use this *miniRUMEL*™ *DEVICE*.


CLOSING

5 SLIDE the miniRUMEL™ tube over the suture toward the tissue until appropriate tissue contact is established, then pull both suture ends to achieve the desired tension.


6 CLOSE the suture clamp by pressing down on the textured back close zone to snap closed the internal features in the suture gripping area and set the suture clamp into its clamped, closed position.


BACK DOWN  CLAMP DOWN

OPENING


7 OPEN the suture clamp by pressing down on the textured forward release zone while simultaneously pushing both white twin tabs forward to set the suture clamp into its unclamped, open position.


8 REMOVE the miniRUMEL™ DEVICE by sliding the miniRUMEL™ over the suture and away from the tissue.


BACK UP  OPEN UP

FIG. 2**OPEN**


Cut-out sections of suture clamp (shown above and below) are provided to illustrate the twin tab lock mechanism and the internal suture gripping area.

FIG. 3**CLOSED**

BACK UP  OPEN UP

BACK DOWN  CLAMP DOWN

If desired, the *miniRUMEL™ DEVICE* wire snare may be passed over the shaft of a suturing device before suturing to ensure that both ends of suture are snared. **NOTE:** Remember to remove any needles, ferrules or suture caps, etc., before threading suture through the *miniRUMEL™ DEVICE*.

FIG. 4



FIG. 5

CONTRAINDICATIONS


- None known.

WARNINGS

- Federal law restricts this device to sale, distribution and use by, or on, the order of a physician.
- Surgical procedures should only be performed by physicians having adequate training and familiarity with such techniques. In addition, medical literature should be consulted relative to techniques, complications and hazards prior to the performance of open or minimally invasive procedures.
- Read and become familiar with all instructions, warnings, and precautions before using this product.
- Do not exert excessive tension on the suture. Excessive suture tensioning can damage or tear tissue, cause bleeding, or adjacent structural injuries to tissue or cause suture breakage.
- Store at room temperature. Avoid prolonged exposure to elevated temperatures.
- Do not resterilize. The *miniRUMEL™ DEVICE* is designed and intended for single patient use only. Do not reuse, reload, reprocess, or resterilize this product. The performance of the *miniRUMEL™ DEVICE* after cleaning or other reprocessing has not been verified and is not supported by LSI SOLUTIONS®. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Do not use if package or product is damaged. Discard any open (unsealed), unused, expired or damaged product.
- Direct contact between sensitive tissue structures and foreign materials can lead to tissue injury or damage.
- Leaving tissue under tension for a prolonged period of time can cause tissue ischemia or damage. This device is for temporary use during a surgical procedure by operators trained and competent in related techniques.
- Ensure *miniRUMEL™* clamp is unclamped (i.e., fully released) and the suture slides freely through the *miniRUMEL™* tube prior to removing the *miniRUMEL™ DEVICE*.
- It is the responsibility of the user to dispose of the device and packaging in accordance with applicable regulations and hospital procedures.
- The *miniRUMEL™* suture lock must be fully closed to ensure tensioning and consistent engagement of the suture within the tube. Loss of suture tension can lead to loss of control at the tissue site and bleeding.

PRECAUTIONS

- Do not use the *miniRUMEL™ DEVICE* if it does not satisfactorily perform its intended function or if it has physical damage. Avoid mechanical impact or overstressing the devices.
- When handling the *miniRUMEL™ DEVICE*, care should be taken to avoid damage.
- Ensure proper surgical techniques, including adequate tissue bites, as needed for the handling of tissue and suture for use of *miniRUMEL™* tubes.
- Adequacy of the suture holding force should be confirmed by the user.
- Check for and ensure hemostasis and/or proper tissue control upon securing the suture and thereafter when appropriate.
- Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure that function is not compromised. Also, ensure that electrical isolation or grounding are not compromised.
- Use only with 2-0 Polypropylene, polyester or silk suture or 3-0 Polypropylene suture.
- For additional suture holding force, a surgical clamp can be applied onto the *miniRUMEL™* tube with its indwelling suture.
- A minimum suture length of 75cm (30 inches) is required for use with the 25 cm *miniRUMEL™ DEVICE* and a minimum suture length of 55cm (22 inches) is required for use with the 15 cm *miniRUMEL™ DEVICE*.
- Avoid holding excess tension on the suture, which can cause the suture to tear or cut through the tissue it is holding or, over time, cause ischemia or tissue strangulation by the tensioned suture.
- Inaccurately placed stay sutures can lead to tissue deformities, inadequate holding and/or tissue damage.
- Inadequate suture tensioning can lead to obstruction in viewing of the surgical view, improper setup of the surgical field, loss of tissue control, bleeding or leakage.

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The LSI logo, LSI SOLUTIONS, miniRUMEL, and Perfect Performance Policy are registered trademarks and trademarks of LSI Solutions, Inc. Copyright © 2023, LSI SOLUTIONS®. All Rights Reserved.

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
LSI SOLUTIONS®
 7796 Victor-Mendon Road
 Victor, New York 14564 U.S.A.
 Phone: 585.869.6600
 Customer Service: 866.575.3493
 Technical Support: 866.428.9092
 Fax: 585.742.8086
www.lsisolutions.com

MADE IN THE USA

This Product Comes
 with our LSI SOLUTIONS®
Perfect Performance Policy®
 Call us at 866.575.3493 any time.

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
MADE IN THE USA

This Product Comes with our LSI SOLUTIONS®
Perfect Performance Policy®
Call us at 866.575.3493 any time.

C1™ DEPTH LIMITER TECHNOLOGY GUIDE

FIG. 1 C1™ DEPTH LIMITER FEATURES

DESCRIPTION OF ACCESSORY

The C1™ Depth Limiter Accessory is a single-use metal component to be attached to the LSI 022000 C1™ Device. When attached to a C1™ Device, the C1™ Depth Limiter accessory limits the location of the target tissue within the C1 tip relative to the needle path thereby limiting the needle's bite depth. It is designed for applications, such as placing sutures, where a more controlled shallow tissue bite may be preferred. The accessory is supplied NON-STERILE by LSI and must be sterilized prior to use.

FIG. 2 SIDE VIEW OF DEPTH LIMITER IN PLACE ON C1™ DEVICE TIP



FIG. 3 DEPTH LIMITER INSTALLATION

The C1™ Depth Limiter is attached to the C1™ Device as shown in Fig. 3, Images A-D. Image A shows the Depth Limiter next to the C1™ Device tip. As shown in Images B - D, the Depth Limiter is inserted into the C1™ Device jaw by:

- Rotating the Depth Limiter 90° and placing the yoke, as defined in Fig. 1, into the opening in the C1™ Device jaw. (Fig. 3, Image B).
- Rotating the Depth Limiter back 90° and placing the yoke against the inside of the tip, above the Ferrule Compartment. (Fig. 3, Image C)
- Pushing the Flexures, as defined in Fig. 1, down so that they snap around the C1™ Device shaft (Fig. 3, Image D).
- Check that the C1™ Depth Limiter is securely attached to the C1™ Device shaft.
- Squeeze and release the C1™ Device lever to advance and retract the C1™ Needle. Check that the needle passes under the C1™ Depth Limiter Tissue Engagement Teeth, and that the C1™ Depth Limiter does not obstruct needle travel.

FIG. 4 DEPTH LIMITER REMOVAL

The C1™ Depth Limiter is removed from the C1™ Device as shown in Fig. 4, Images A-C.

- DO NOT advance the C1™ Needle while removing the C1™ Depth Limiter.
- Image A shows the Depth Limiter with the C1™ Device Jaw facing up. The C1™ Device should be placed on something firm, such as the surface of a sterile tray.
- Image B shows the next step, to lay a surgical instrument, such as grasper, alongside of the C1™ Device tip, with the instrument tip pushing against the Tab on the side of the Depth Limiter.
- Pushing down on the Tab with the grasper will dislodge the Depth Limiter from the C1™ Jaw and Shaft, as shown in Image C.
- The Depth Limiter can now be removed.

1. ACTIONS

The C1™ Depth Limiter restricts the amount of tissue that can fit into the C1™ Device jaw. The C1™ Depth Limiter is used with the C1™ Device to facilitate placing non-controlled shallower bites in tissue, for example, for non-transmural or non-full thickness tissue bites.

2. CONTRAINDICATIONS

- Minimally invasive surgical procedures should only be performed by physicians having adequate training and familiarity with minimally invasive techniques, including suturing techniques. In addition, medical literature should be consulted relative to techniques, complications and hazards prior to the performance of minimally invasive procedures.
- The C1™ Depth Limiter is not intended to be used with any device other than a C1™ DEVICE.

3. WARNINGS

- Do not resterilize or reuse. The performance of the C1™ Depth Limiter after cleaning or other reprocessing has not been verified and is not supported by LSI SOLUTIONS®.
- Discard open (unsealed), unused, or damaged devices or devices in damaged primary packaging.
- Users should be familiar with surgical procedures and techniques involving suture before employing the C1™ DEVICE for wound closure, as the risk of wound dehiscence may vary with the site of application.
- Acceptable surgical practice and good surgical judgment must be followed with respect to drainage and closure of infected or contaminated wounds.

4. PRECAUTIONS

- Check for hemostasis or leakage where appropriate.
- Always assure C1™ Depth Limiter and C1™ Device Tip are viewed under direct visualization before advancing the needle.
- Ensure obstructions do not interfere with C1™ Depth Limiter.
- Do not squeeze the blue lever of the C1™ DEVICE while attaching or detaching the C1™ Depth Limiter; squeezing the blue lever may expose the sharp needle and/or damage the needle.
- Avoid damage to or dislodgement of the C1™ Depth Limiter due to direct application of surgical instruments, like forceps, needle holders, clamps, etc.
- Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure that mechanical function and electrical isolation and grounding are not compromised.
- The surgeon's judgment regarding bite location remains imperative, for example, where avoidance of calcific tissue is usually warranted.

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6. Sterilization:

Make sure all exterior surfaces are free of dents, burrs, corrosion, or jagged edges. Do not use if device is bent or damaged.

C1™ Depth Limiter STERILIZATION INSTRUCTIONS:

The device must be properly cleaned and dried prior to sterilization.

STEAM AUTOCLAVE STERILIZATION:

- a. The C1™ Depth Limiter can be sterilized by steam autoclave with the following parameters (refer to AAMI TIR12): Double wrap the C1™ Depth Limiter in sterilization wrap.
- b. Perform sterilization cycle with the following parameters:
 - Pre-vacuum cycle: 132°C minimum, Full exposure time at 4 minutes
 - Drying time: 30 minutes

IMMEDIATE USE / “FLASH” STEAM STERILIZATION:

- a. Do NOT wrap C1™ Depth Limiter in sterilization wrap. Position devices directly into autoclave to allow direct steam exposure to all surfaces.
- b. Perform Immediate Use Sterilization with the following parameters:
 - Pre-vacuum cycle: 132°C minimum, Full exposure time at 4 minutes
 - Drying Time: N/A

Device(s) processed by Immediate Use Sterilization should be transferred immediately, using aseptic technique, from the sterilizer to the point of use.

Refer to ANSI/AAMI ST79 Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities.

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
FIG. 1

INTRA-KNOT® QUICK LOAD® UNIT DESCRIPTION

Each INTRA-KNOT® QUICK LOAD® UNIT consists of a blunt curved handle ①, one sterile COR-KNOT® FASTENER ② held in a customized loading unit consisting of a white target ③, a proximal wire snare ④, a connecting wire ⑤, a purple target ⑥, and a distal wire snare ⑦. Use only with LSI SOLUTIONS® specified 2-0 Polyester, 2-0 Polypropylene, or 3-0 Polypropylene suture.


INDICATIONS


The INTRA-KNOT® QUICK LOAD® UNIT, used in conjunction with the COR-KNOT® DEVICE, is indicated for use in the approximation of soft tissue and prosthetic materials. Refer to the Indications in the COR-KNOT® MINI® Technology Guide or the COR-KNOT® MIS Technology Guide for compatibility with INTRA-KNOT® QUICK LOAD® UNIT.

FIG. 2

INTRA-KNOT® PRODUCT ORDERING

SUPPLIED: STERILE

	REORDER	PRODUCT	DESCRIPTION
 x 12	REF 032040	INTRA-KNOT® QUICK LOAD® UNIT	Box of 12


ACTUAL LENGTH: 40cm

LSI SOLUTIONS®

LOADING AT SCRUB TABLE


Use proper operating room technique to pass the sterile *INTRA-KNOT*® QUICK LOAD® UNIT from its packaging. While maintaining appropriate sterile technique, follow the steps indicated in the illustrations.

1 INSERT


INSERT the blunt tip of the curved handle into the distal slot at the end of the *COR-KNOT*® *DEVICE* shaft. **ROTATE** the curved handle through the distal slot and out of the suture hole until the *COR-KNOT*® *FASTENER* occupies the shaft's distal slot. Fully **ENGAGE** the *COR-KNOT*® *FASTENER* within the tip of the *COR-KNOT*® *DEVICE* by pushing on the white target or by pulling on the curved handle.

2 PUSH-OUT


PUSH-OUT and remove the white target.

3 RELEASE

RELEASE the curved handle from the distal slot at the end of the *COR-KNOT*® *DEVICE* shaft.

4 INSPECT


INSPECT to ensure that the *COR-KNOT*® *FASTENER* is well loaded and fully seated.

INTRACORPOREAL MIS USE

For minimally invasive procedures, the distal snare of the *INTRA-KNOT*® QUICK LOAD® UNIT can be used to intracorporeally thread the suture. Follow steps 1 through 4 to load the *INTRA-KNOT*® QUICK LOAD® UNIT.

5 PUSH-OUT


6 FOLD

FOLD the distal wire snare back along the side of the *COR-KNOT*® shaft.


INSIDE OF PATIENT

7 INSERT


With the wire snare folded back against the shaft, **INSERT** the shaft through the cannula down to the targeted site.

8 PASS


PASS both ends of the suture through the unfolded distal wire snare.

9 THREAD


THREAD the suture through the *COR-KNOT*® FASTENER by continuing to pull the curved handle until the suture ends pass through the suture hole. Then **SLIDE** the *COR-KNOT*® DEVICE distal tip gently over the partially tensioned suture down to the targeted site.

10 SLIDE

11 SQUEEZE & HOLD 1 Second

To crimp *COR-KNOT*® FASTENER

SQUEEZE the purple lever to hard stop, maintain the device tip's position, **HOLD** it for one second.


12 TUG & RELEASE

Tug to trim suture

While still squeezing the purple lever, **TUG** all the sutures gently to cut them and **RELEASE** the purple lever fully to release the crimped *COR-KNOT*® FASTENER.


13 REMOVE DEVICE

Inspect *COR-KNOT*® FASTENER

REMOVE DEVICE Inspect to ensure *COR-KNOT*® FASTENER and suture tails are oriented away from delicate tissue and prosthetic structures.

14 PUSH & ROTATE


Rotate 1/4 turn

If the crimped *COR-KNOT*® FASTENER does not readily release from the distal tip, ensure purple lever is released, then gently **PUSH** inward and **ROTATE** the handle 90° about the shaft. If still necessary, **ROTATE** the handle back and then turn 90° in the opposite direction. **If *COR-KNOT*® FASTENER will still not release, cut suture.**


OUTSIDE OF PATIENT

5 PASS


PASS the free end of the suture through the proximal wire snare at the end of the shaft.

6 PULL


PULL the curved handle with its attached wire snare containing the suture end towards the purple lever to draw the snared bend in the suture into the COR-KNOT® FASTENER.

7 THREAD


Continue to pull the curved handle until the proximal wire snare and the end of the suture exit through the suture hole.

8 REMOVE


REMOVE the suture from the proximal wire snare. **GRASP** the suture with a clamp so that it doesn't become tangled.

9 PUSH-OUT


PUSH-OUT and remove the purple target.

10 PASS

PASS the end of the suture attached to the suturing device through the distal wire snare. If using a suturing device, fire device over the wire snare.


11 THREAD

Continue to pull the curved handle until the distal wire snare and suture exit through the suture hole. Be sure not to tangle the sutures with the **INTRA-KNOT®** wire.

INSIDE OF PATIENT

12 SLIDE


SLIDE the **COR-KNOT® DEVICE** distal tip gently over the partially tensioned suture down to the targeted site. **KEEP** suture next to shaft while pulling. With the distal tip on the tissue, maintain sufficient suture tension to hold tissue in appropriate apposition.

13 SQUEEZE & HOLD


1 Second to crimp
COR-KNOT® FASTENER

SQUEEZE the purple lever to hard stop, maintain the device tip's position, **HOLD** it for one second.


14 TUG & RELEASE

While still squeezing the purple lever, **TUG** all the sutures gently to cut them and **RELEASE** the purple lever fully to release the crimped **COR-KNOT® FASTENER**.


15 PUSH & ROTATE

If the crimped **COR-KNOT® FASTENER** does not readily release from the distal tip, ensure purple lever is released, then gently **PUSH** inward and **ROTATE** the handle 90° about the shaft. If still necessary, **ROTATE** the handle back and then turn 90° in the opposite direction. If **COR-KNOT® FASTENER** will still not release, cut suture.

16 RUN SUTURE

Continue running suture with the sewing device.

In the text below, the *COR-KNOT® DEVICE* refers to the LSI SOLUTIONS® *COR-KNOT®* family of suture fastening devices. Refer to *COR-KNOT MINI®* Technology Guide or the *COR-KNOT® MIS* Technology Guide Indications for compatibility with *INTRAKNOT® QUICK LOAD® UNIT*.

ACTIONS

When the *COR-KNOT® DEVICE* is loaded with a *COR-KNOT® FASTENER* using the *INTRAKNOT® QUICK LOAD® UNIT* and appropriately positioned at a suture closure site, squeezing the purple lever will instantly secure and trim the suture. The surgical titanium used in a *COR-KNOT® FASTENER* is not absorbed by the body and is generally not associated with significant inflammatory reactions. The *INTRAKNOT® QUICK LOAD® UNIT* can be used to secure 2, 3, or 4 strands of 3-0 suture or 2 strands of 2-0 suture.

CONTRAINDICATIONS

- Endoscopic procedures should only be performed by physicians having adequate training and familiarity with endoscopic techniques. Medical literature should be consulted relative to techniques, complications, and hazards prior to the performance of endoscopic procedures.
- The *INTRAKNOT® QUICK LOAD® UNIT* is not intended to be used with any device other than the *COR-KNOT MINI® DEVICE* or the *COR-KNOT® MIS DEVICE*.
- The *COR-KNOT MINI® DEVICE* and the *COR-KNOT® MIS DEVICE* are not intended to be loaded with anything other than a *COR-KNOT® QUICK LOAD® UNIT* or the *INTRAKNOT® QUICK LOAD® UNIT*.
- The *COR-KNOT® FASTENER* is NOT intended for placement into circulating blood unless used with compatible suture under conditions judged by the surgeon to be clinically appropriate.
- Use only with LSI SOLUTIONS® specified 2-0 Polyester, 2-0 Polypropylene, or 3-0 Polypropylene suture.
- Each *COR-KNOT® DEVICE* is not intended to be fired more than 12 times.

WARNINGS

- Federal (U.S.A.) law restricts this device to sale, distribution and use by, or on, the order of a physician.
- Users should be familiar with standard procedures and techniques involving surgical suture and titanium usage before employing the *COR-KNOT® DEVICE* with a *COR-KNOT® FASTENER* for fastening and trimming suture.
- Adequate *COR-KNOT® FASTENER* security requires reasonable clinical judgment and appropriate surgical techniques as warranted by surgical circumstances and the experience of the surgeon.
- When securing suture with an LSI *COR-KNOT® DEVICE*, ensure any ferrules or needles are removed from the suture ends to be loaded prior to loading the suture through the *COR-KNOT® DEVICE*.
- Excessive suture tensioning can cause suture breakage.
- Single patient use only. Do not reclean or resterilize. Adequate cleaning or removal of blood and other foreign materials from used *COR-KNOT®* products cannot be guaranteed. Validation of reesterilization is not established. Failure to eliminate inflammatory or infectious agents may cause patient harm. Product functional performance may be compromised in reprocessed devices or *COR-KNOT® FASTENERS*.
- Discard any open (unsealed), unused, expired or damaged *COR-KNOT®* product.
- *INTRAKNOT® QUICK LOAD® UNIT* components and each *COR-KNOT® DEVICE*, along with packaging, must be inspected, handled and disposed of consistent with standard, accepted medical device disposal procedures.
- Direct contact between sensitive tissue structures (e.g., pulsatile arteries, cardiac valve leaflets, valve chordae, etc.) and foreign materials can lead to tissue injury or damage, such as tissue erosion. Always orient *COR-KNOT® FASTENERS* and remnant suture tails to avoid direct contact between delicate tissue or prosthetic structures.
- As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.
- While the titanium of the *COR-KNOT® FASTENER* is physiologically very inert, routine surgical precautions must be employed whenever foreign materials are left in a patient.

PRECAUTIONS

- When handling the *INTRAKNOT® QUICK LOAD® UNIT* care should be taken to avoid damage.
- Do not squeeze the purple lever of the *COR-KNOT® DEVICE* while loading the *INTRAKNOT® QUICK LOAD® UNIT*.
- Irreparable damage to *COR-KNOT® DEVICE* suture cutting blade will occur if the purple lever is squeezed while the *INTRAKNOT® QUICK LOAD® UNIT* curved handle is in place at the tip of the instrument.
- Ensure that obstructions do not interfere with the firing of the *COR-KNOT® DEVICE*.
- Do not squeeze the purple lever of the loaded *COR-KNOT® DEVICE*, until the *COR-KNOT® FASTENER* has been appropriately positioned directly upon the tissue or prosthetic material and the suture accurately tensioned at the targeted site.
- Always squeeze and hold the purple lever and then fully release it before moving the *COR-KNOT® DEVICE* tip. Failure to appropriately release the purple lever can cause suture breakage. Inspect each *COR-KNOT® FASTENER* and its suture tails.
- Do not squeeze the purple lever on the same *COR-KNOT® FASTENER* more than once.
- Cut suture with scissors if the device fails to trim suture or release the *COR-KNOT® FASTENER*.
- Avoid crushing or crimping damage to the *COR-KNOT® FASTENER* due to inappropriate squeezing of the *COR-KNOT® DEVICE* purple lever and/or to the application of surgical instruments like forceps, needle holders, clamps, etc.
- If the *COR-KNOT® FASTENER* falls out of tip or is not properly loaded, retrieve loose fastener, reload with new fastener and start again.
- If the purple lever of the *COR-KNOT® DEVICE* does not return completely forward on its own (i.e., without assistance), manually push the lever forward all the way to release the *COR-KNOT® FASTENER*.
- Check for hemostasis or leakage where appropriate.
- Before endoscopic instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure that electrical isolation or grounding are not compromised.

ADVERSE REACTIONS

Adverse effects associated with the use of surgical suture and titanium can include, but are not limited to: wound dehiscence, thrombus formation, embolism, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infected wounds, minimal acute inflammatory tissue reaction, and transitory local irritation. Surgical titanium is not absorbed by the body and is generally not associated with inflammatory reactions.

MRI TESTING

Based on MRI testing information, titanium *COR-KNOT® FASTENERS* will not present an additional hazard or risk to a patient undergoing an MRI procedure using a scanner operating with a static magnetic field of 3-Tesla or less and under the MRI-related heating conditions (MRI for 15 min. at an MR system reported whole body averaged specific absorption rate, SAR, value of 3-W/kg).

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The LSI logo, LSI SOLUTIONS, Intra-Knot, Cor-Knot, Cor-Knot Mini, Quick Load, and Perfect Performance Policy are trademarks and registered trademarks of LSI Solutions, Inc. Copyright © 2014, LSI SOLUTIONS®. All Rights Reserved.

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P/N 032043F 8.1.19


COR-SUTURE® QUICK LOAD® SURGICAL SUTURE TECHNOLOGY GUIDE

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FIG. 1

COR-SUTURE® QUICK LOAD® SURGICAL SUTURE

- ① SUTURE TRAY
- ② SUTURE RELEASE FEATURE
- ③ NEEDLE CAPS
- ④ SUTURE
- ⑤ SUTURE TUBE

DESCRIPTION

Each LSI SOLUTIONS® COR-SUTURE® QUICK LOAD® surgical suture is held in a customized tray (①) with suture release feature (②) designed to enable the rapid, easy and reliable loading of suture into an appropriate LSI SOLUTIONS® suturing device. COR-SUTURE® QUICK LOAD® surgical suture is available as a non-absorbable monofilament polypropylene and braided polyester surgical suture (FIG. 2). A short length of modified surgical stainless steel tubing, called a needle cap (③), is attached to each end of the suture (④). The COR-SUTURE® QUICK LOAD® surgical suture also includes a detachable clear suture tube (⑤) to keep the suture from tangling. Each sterile COR-SUTURE® QUICK LOAD® surgical suture is individually packaged for single patient use.

INDICATIONS FOR USE

COR-SUTURE® QUICK LOAD® surgical suture is indicated for use in the approximation of soft tissue and prosthetic materials.

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The LSI logo, LSI SOLUTIONS, C1, Cor-Suture, Cor-Knot, Cor-Knot Mini, RD, Intra-Knot, Quick Load, RAM, Sew-Easy and Perfect Performance Policy are trademarks and registered trademarks of LSI SOLUTIONS, Inc. Copyright © 2014, LSI SOLUTIONS®. All Rights Reserved.

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ACTIONS

COR-SUTURE® QUICK LOAD® surgical suture typically elicits a minimal acute inflammatory reaction in tissues and gradual encapsulation of the suture by fibrous connective tissue. Both non-absorbable sutures, braided polyester and monofilament polypropylene suture, are not absorbed, nor is any significant change in tensile strength retention known to occur in vivo. The COR-SUTURE® QUICK LOAD® surgical suture is MR safe. The COR-SUTURE® QUICK LOAD® surgical suture must be used with a compatible LSI SOLUTIONS® sewing device. Compatible devices are called out as such in their respective technology guides.

CONTRAINDICATIONS

- Minimally invasive surgical procedures should only be performed by physicians having adequate training and familiarity with minimally invasive techniques. In addition, medical literature should be consulted relative to techniques, complications and hazards prior to the performance of minimally invasive procedures.
- Do not use this suture under conditions in which excessive suture tension can lead to tissue damage.

WARNINGS

- Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.
- Do not resterilize.** The performance of the COR-SUTURE® QUICK LOAD® surgical suture after cleaning or other reprocessing has not been verified and is not supported by LSI SOLUTIONS, Inc.
- Discard open (unsealed), unused, expired or damaged suture or suture in damaged primary packaging.
- As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.
- Users should be familiar with surgical procedures and techniques involving suture before employing the COR-SUTURE® QUICK LOAD® surgical suture for wound closure, as the risk of wound dehiscence may vary with the site of application.
- Acceptable surgical practice and good surgical judgement must be followed with respect to drainage and closure of infected or contaminated wounds.
- Redundant, cut-away suture remnants and used needle caps, along with packaging, must be inspected, handled and disposed of consistent with standard, accepted medical device disposal procedures.
- COR-SUTURE® QUICK LOAD® surgical suture is not for use in neurological procedures.
- Applications other than for soft tissue approximation, with or without appropriate prosthetic materials, can result in failure to pick up suture or in damage to the device rendering either unsuitable for continued use.
- Never drive the suturing device needle into suture, bone, dense ligamentous tissue, calcific tissue or other instruments.
- Do not leave any foreign material (e.g. suture remnant, needle cap, etc.) unattached in areas potentially exposed to circulating blood.
- If securing the COR-SUTURE® QUICK LOAD® surgical suture with any LSI SOLUTIONS® COR-KNOT® DEVICE (COR-KNOT® MIS DEVICE, COR-KNOT MINI® DEVICE), ensure any needle caps are removed from the suture ends being loaded, prior to loading the suture through the COR-KNOT® DEVICE.

PRECAUTIONS

- Check for hemostasis or leakage where appropriate.
- Always assure device tip is viewed under direct visualization before advancing the needle.
- Ensure obstructions do not interfere with the movement of the needle of suturing device.
- Do not squeeze the lever of the suturing device while loading a COR-SUTURE® QUICK LOAD® surgical suture; squeezing the lever may expose the sharp needle and/or damage the needle.
- In handling the suturing device, care should be taken to avoid jamming the suture into the needle cap compartment and damaging the needle.
- Avoid damage to the needle, suture or needle caps due to direct application of surgical instruments, like forceps, needle holders, clamps, etc.
- Adequate knot security requires accurate completion of accepted surgical techniques for constructing surgically tied knots or the use of the COR-KNOT® MIS DEVICE, COR-KNOT MINI® DEVICE, and COR-KNOT® QUICK LOAD® fastener or INTRA-KNOT® QUICK LOAD® fastener as warranted by surgical circumstances and the experience of the surgeon.
- Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure that mechanical function and electrical isolation and grounding are not compromised.
- Before loading the suturing device with another COR-SUTURE® QUICK LOAD® surgical suture, assure the remaining suture tail and needle cap from the previous load has been removed from the needle and needle cap compartment. Failure to appropriately remove used needle caps from the needle can result in damage to the device, including intracorporeal or extracorporeal fracturing off of the tip of the needle, making it unsuitable for continued use.
- Verify that the needle cap is retained within the needle cap compartment and the device has not been damaged or deformed before attempting to place a stitch.
- Do not use damaged suture or expired suture.

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Carefully read and understand the instructions for other Technology Guides for the related LSI SOLUTIONS® suture device products prior to using COR-SUTURE® QUICK LOAD® surgical suture.

FIG. 2 PRODUCT ORDERING			SUPPLIED: STERILE
	REORDER	PRODUCT	DESCRIPTION
 x 12	REF 021860	COR-SUTURE® QUICK LOAD® SURGICAL SUTURE 3-0 Polypropylene, Non-Absorbable, 53"	Box of 12 Sutures (1 Suture per Pouch)
 x 12	REF 021861	COR-SUTURE® QUICK LOAD® SURGICAL SUTURE 2-0 Polyester White, Non-Absorbable, 38"	Box of 12 Sutures (1 Suture per Pouch)
 x 12	REF 021862	COR-SUTURE® QUICK LOAD® SURGICAL SUTURE 2-0 Polyester Green, Non-Absorbable, 38"	Box of 12 Sutures (1 Suture per Pouch)

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P/N 022028E 07.14.2020



LS-5™ ePTFE QUICK LOAD® SUTURE TECHNOLOGY GUIDE

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FIG. 1

LS-5™ ePTFE QUICK LOAD® SUTURE

- ① SINGLE SUTURE TRAY
- ② SUTURE RELEASE FEATURE
- ③ FERRULE
- ④ SUTURE
- ⑤ SUTURE TUBE

DESCRIPTION

Each LSI SOLUTIONS® LS-5™ ePTFE QUICK LOAD® sterile surgical suture with 2 ferrules is held in a customized tray (①), with a suture release feature (②) and is intended for single patient use. The LS-5™ ePTFE QUICK LOAD® suture is designed to enable the rapid, easy and reliable loading of suture into compatible LSI SOLUTIONS® suturing devices. LS-5™ ePTFE QUICK LOAD® suture is a non-absorbable monofilament expanded polytetrafluoroethylene (ePTFE) suture. A short length of modified surgical stainless steel tubing, called a "ferrule" (③), is attached to each end of the suture (④) for interfacing with LSI suturing devices. The LS-5™ ePTFE QUICK LOAD® suture also includes a detachable clear suture tube (⑤) to keep the suture from tangling. The LS-5™ ePTFE QUICK LOAD® surgical suture is undyed and contains no additives.

INDICATIONS

LS-5™ ePTFE QUICK LOAD® surgical suture is indicated for use in all types of soft tissue approximation, including use in cardiovascular surgery.

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The LSI logo, LSI SOLUTIONS, LS-5, and Quick Load are registered trademarks and trademarks of LSI SOLUTIONS, Inc.
Copyright © 2016, LSI SOLUTIONS®. All Rights Reserved.

LSI SOLUTIONS®
INNOVATIVE SOLUTIONS FOR MINIMALLY INVASIVE SURGERY

Consult Operating
Instructions

STERILE EO

Sterilized Using Ethylene Oxide

Do Not Reuse

Do Not Resterilize

ACTIONS

ePTFE suture is a very inert suture material known to elicit minimal tissue reaction. LS-5™ ePTFE suture is not absorbed and is not known to degrade or lose strength over time. No significant loss in tensile strength retention occurs in vivo. The LS-5™ ePTFE QUICK LOAD® surgical suture is MR safe.

CONTRAINDICATIONS

- The LS-5™ ePTFE QUICK LOAD® surgical suture is contraindicated for use in ophthalmic surgery, microsurgery, and neural tissue.
- Do not use this suture under conditions in which excessive suture tension can lead to tissue damage.

WARNINGS

- Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.
- Do not resterilize.** The performance of the LS-5™ ePTFE QUICK LOAD® surgical suture after cleaning or other reprocessing has not been verified and is not supported by LSI SOLUTIONS, Inc.
- Discard open (unsealed), unused, or damaged suture or suture in damaged primary packaging.
- As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.
- Minimally invasive surgical procedures should only be performed by physicians having adequate training and familiarity with minimally invasive techniques. In addition, medical literature should be consulted relative to techniques, complications and hazards prior to the performance of minimally invasive procedures.
- Users should be familiar with surgical procedures and techniques involving suture before employing the LS-5™ ePTFE QUICK LOAD® surgical suture for wound closure, as the risk of wound dehiscence may vary with the site of application.
- Acceptable surgical practice and good surgical judgement must be followed with respect to drainage and closure of infected or contaminated wounds.
- Redundant, cut-away suture remnants and used ferrules, along with packaging, must be inspected, handled and disposed of consistent with standard, accepted medical device disposal procedures.
- Applications other than for soft tissue approximation, can result in failure to pick up suture or in damage to the suturing device rendering either unsuitable for continued use.
- Never drive the suturing needle into suture, bone, dense ligamentous tissue, severely calcific tissue or other instruments.
- Do not leave any foreign material (e.g. suture fragment, ferrule, etc.) unattached in areas potentially exposed to circulating blood.
- The LS-5™ ePTFE QUICK LOAD® surgical suture is intended for single patient use only.
- Tissue invasion of the LS-5™ ePTFE QUICK LOAD® surgical suture can result in attachment of the suture to the tissue it penetrates. Such attachment may make removal of the LS-5™ ePTFE QUICK LOAD® surgical suture difficult and, if attempted, may result in tissue damage.

PRECAUTIONS

- Check for hemostasis or leakage where appropriate.
- Always assure device tip is well visualized before advancing the needle.
- Ensure obstructions do not interfere with movement of the suturing device needle.
- Do not squeeze the lever of the suturing device while loading an LS-5™ ePTFE QUICK LOAD® surgical suture; squeezing the lever may expose the sharp needle and/or damage the needle.
- In handling the suture, care should be taken to avoid jamming the suture into the ferrule compartment and damaging the needle.
- Avoid damage to the needle, suture or ferrules due to direct application of surgical instruments, like forceps, needle holders, clamps, etc.
- Adequate knot security requires accurate completion of accepted surgical techniques for surgically constructed knots.
- Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure that mechanical function and electrical isolation and grounding are not compromised.
- Verify that the ferrule is retained within the ferrule compartment and the device has not been damaged or deformed before attempting to place a stitch.
- Do not use damaged suture.

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Carefully read and understand the Technology Guides for the related LSI SOLUTIONS® suturing device products prior to using LS-5™ ePTFE QUICK LOAD® surgical suture.

LS-5™ ePTFE sutures differ from and/or exceed USP requirements. See Fig. 2 and Fig. 3 below for LS-5™ suture characteristics.

FIG. 2 Comparison of LS-5™ ePTFE Suture to CV-5 GORE-TEX® Suture			* Per LSI documents and testing. ** Per GORE-TEX® IFU, SEPT. 2015
Property	LS-5™	CV-5	
Knot Pull Tensile Strength (kg)	> 0.60*	> 0.60**	
Average Diameter (mm)	0.305	0.246	

LSI LS-5™ ePTFE Suture is comparable to GORE-TEX® CV-5 suture, which has been marketed for decades as equivalent to USP 4-0 (FIG. 2). Both LS-5™ and GORE-TEX® CV-5 ePTFE suture are slightly oversized relative to the USP diameter range designation of 0.15 – 0.199 mm. Both sutures also exceed the minimum knot tensile pull apart parameters established by the USP for 4-0 of 0.60 kg. (FIG. 3).

FIG. 3 USP Requirements for Non-Absorbable Suture			
USP Size	Knot Pull Tensile Strength (kg)	Minimum Diameter (mm)	Maximum Diameter (mm)
2-0	1.44	0.30	0.399
3-0	0.96	0.20	0.249
4-0	0.60	0.15	0.199
5-0	0.40	0.10	0.149

FIG. 4 PRODUCT ORDERING			SUPPLIED: STERILE
	REORDER	PRODUCT	DESCRIPTION
 x 12	REF 022295	LS-5™ ePTFE QUICK LOAD® Suture, Monofilament, Non-Absorbable, 38", 2 Ferrules	Box of 12 Sutures (1 Suture per Pouch)

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Copyright © 2020, LSI SOLUTIONS®. Все права защищены. Патенты:
Патенты: www.lsisolutions.com/patents Логотип LSI, LSI SOLUTIONS, miniRUMEL и Perfect Performance Policy являются зарегистрированными торговыми марками и торговыми марками LSI Solutions, Inc. Copyright © 2023, LSI SOLUTIONS®. Все права защищены.
Логотип LSI, LSI SOLUTIONS, Intra-Knot, Cor-Knot, Cor-Knot Mini, Quick Load и Perfect Performance Policy являются товарными знаками и зарегистрированными товарными знаками компании LSI Solutions, Inc. Copyright © 2014, LSI SOLUTIONS *. Все права защищены.