The multiSTATION™ System: miniARM® triDOC™, monoDOC™, duoDOC™ L, duoDOC™ R RAIL CLAMP Technology Guide

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
FIG 1A: miniARM® triDOC™ RAIL CLAMP Front View

FIG 1B: miniARM® triDOC™ RAIL CLAMP Back View

SOLUTIONS®

miniARM® RAIL CLAMP DEVICE DESCRIPTION:

The miniARM® RAIL CLAMP device is a reusable, sterilizable, repositionable, multiple instrument holding metal base that clamps over the sterile drapes onto the underlying operating room table horizontal rail using a proprietary manually powered jaw mechanism for attachment of multiSTATION™ System accessories. FIG. 1 shows the front and back view of the triDOC™ RAIL CLAMP incorporating the crenellated central DOC 1, which is the receiving socket for the attachment of these accessories above a clamp knob 2 that rotates on the top of the fenestrated central stem 3. When viewed from the front, to the left of the central stem is the crenellated left DOC 4 and the left stem 5; to the right, the crenellated right DOC 6 and the right stem 7. The stems connect to the clamp body 8 at the front junction 9. A connector passing through the fenestrated central stem links the central DOC to a mobile upper jaw 10 positioned above two fixed lower jaws 11. Side DOCs also include removable DOC covers 12. The miniARM® RAIL CLAMP is available in four configurations based on the number of side DOCs and their locations: the miniARM® triDOC™ (FIG. 1, FIG. 2A), the monoDOC™ (FIG. 2B), the duoDOC™ L {left} (FIG. 2C) and the duoDOC™ R {right} (FIG. 2D).

{Continues on page 2}


FIG. 2A-D: Front View of the 4 different *miniARM*® RAIL CLAMP product offerings. 2A: *triDOC*™, 2B: *monoDOC*™, 2C: *duoDOC*™ L, 2D: *duoDOC*™ R

***miniARM*® RAIL CLAMP DEVICE**

DESCRIPTION (continued): FIG 3. shows a left cutaway of the front of a *triDOC*™ RAIL CLAMP with a purple colorized central DOC, connector link and mobile upper jaw in the “Up” position (FIG. 3A) compared to the upper jaw “Down” position (FIG. 3B) relative to a surgical table horizontal rail. With the fixed bottom jaws positioned (FIG. 3A) under the bottom surface of the targeted section of the draped operating room table rail, clockwise rotation of the clamp knob drives the mobile upper jaw down onto the upper surface of the table rail (FIG. 3B), thereby securely affixing the *miniARM*® RAIL CLAMP over the drape and onto the table. Each uncovered crenellated DOC is threaded and contoured to enable precise attachment and orientation of *multiSTATION*™ System accessories (e.g. *M45*™ CLICK RELEASE™ *miniARM*® HOLDER or *M55*™ RIGID *miniARM*® HOLDER). After removal of all *miniARM*® RAIL CLAMP accessories from each crenellated DOC, the *miniARM*® RAIL CLAMP can be repositioned or removed by counterclockwise rotation of the clamp mounting knob to pull the mobile upper jaw away (FIG. 3A) from its clamped position.


FIG. 3A: *triDOC*™ RAIL CLAMP in the “Up” Position

FIG. 3B: *triDOC*™ RAIL CLAMP in the “Down” Position

INDICATIONS FOR USE:


The *miniARM*® *monoDOC*™, the *duoDOC*™ L, the *duoDOC*™ R, and the *triDOC*™ RAIL CLAMPs are indicated for use by surgeons to hold instruments in a fixed position for a period of time.

INSTRUCTIONS FOR USE

GENERAL SETUP & USE:

The steps below show the *triDOC™ RAIL CLAMP*. The same attachment technique applies to the *monoDOC™*, *duoDOC™ L*, and *duoDOC™ R RAIL CLAMP* configurations. The *miniARM® RAIL CLAMP* is designed to attach to a surgical table's horizontal rail with dimensions of 9–10 mm in width x 25–30 mm in height (0.35–0.39 in. x 0.98–1.18 in.).


1 OPEN

ATTACHMENT TO A SURGICAL TABLE


1. **OPEN** the space between the mobile upper jaw and the fixed lower jaws of the clamp body to become wide enough to receive the selected portion of the operating room rail by turning the clamp knob counterclockwise to draw the mobile jaw away from the fixed lower jaws. DOC covers should be attached on side DOCs to reduce exposure of DOC crenellations.

2 POSITION


2. **POSITION** the jaws of the clamp body over the draped table rail with the front of the rail clamp facing away from the table and angled straight up. Avoid damaging the surgical drapes.

3 CLOSE

3. **CLOSE** the mobile upper jaw of the rail clamp down onto the fixed lower jaws by turning the clamp knob clockwise until it comes to a firm stop. Do not touch below level of sterility described in your hospital's policies. Care should be taken to not touch the crenellations within the DOCs, which could damage gloves or cause injury if inappropriately handled.


4 UNCAP


4. **UNCAP** crenellated side DOC(s) by turning the DOC cover(s) counterclockwise until fully unthreaded from the rail clamp. Keep the DOC cover(s) within the sterile field for reinstallation if accessory is detached and/or prior to rail clamp removal. Compatible accessories may now be attached per each component's instructions for use technology guide.

REMOVAL AND CLEANUP

Remove all previously attached accessories from the *miniARM*® *HOLDER* as well as each *miniARM*® *HOLDER* according to each component's instructions for use. Reattach sterile DOC cover(s) to side DOC(s).



1. **OPEN** the rail clamp jaws by turning the clamp knob counterclockwise until the rail clamp jaws fully release from the surgical table rail.


2. **REMOVE** rail clamp with attached DOC cover or covers by grasping the upper portion of a side stem and passing the rail clamp off from the surgical field.

miniARM® RAIL CLAMP REPROCESSING

- Disassemble and clean the device immediately after use. Do not allow a soiled device to dry.
- The device is not validated to be cleaned or sterilized with a holder and/or adapter attached including all DOC Covers.
- Cleaning agent used in validation: Steris Prolystica 2X (enzymatic, neutral pH).
- Perform the final rinse using only freshly prepared purified water/highly purified water.
- Never use metal brushes or steel wool for cleaning.
- Prepare and reprocess any holders or adapters according to each device's instruction for use.
- Avoid the risk of damaging gloves or injuring skin on the crenellations within the DOCs.

POINT OF USE	1. Disassemble and clean device immediately after use. 2. Do not allow soiled devices to dry.
PREPARATION	1. Turn the rail clamp knob counterclockwise until all threads (in the red circle) are visible and the rail clamp knob stops turning (FIG. 4) and the <i>miniARM</i>® RAIL CLAMP is as high as it goes (FIG. 5)   FIG. 4 FIG. 5

INSPECTION (CONTINUED)

3. Visually inspect the device for damage. If parts are damaged, discontinue use of the device in surgery, but complete reprocessing.
4. Perform a functional check for adequate tightening and release by turning the clamp knob clockwise and counterclockwise.

PACKAGING

1. Ensure the upper jaw is in the fully open reprocessing position (FIG. 9).
2. If using a sterilization pouch, place each instrument in its own individual pouch. If using a sterilization tray, read the sterilization tray instructions for use before proceeding. Ensure that the sterilization tray has been cleaned according to the sterilization tray instructions for use and load the tray bottom according to FIG. 9. Place detached DOC cover(s) in general instrument cleaning area.
3. Package the device according to TABLE 1. The barrier system for sterilized reusable instruments should meet the following requirements:
 - ISO 11607-1
 - Suitable for pre-vacuum steam sterilization
 - Appropriate for medical use
 - Grade appropriate for weight of loaded tray per sterilization tray instructions for use and facility procedures

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
TABLE 1: miniARM® RAIL CLAMP Sterile Packaging and Processing

Method	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Immediate Use Steam Sterilization according to ANSI/AAMI ST79
Container	<i>multiSTATION™ System Sterilization Tray</i> P/N 100034	No tray	No tray
Cycle	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)
Packaging	2-layer polypropylene wrap	Pouch	No packaging
Temperature	132°C (270°F) (minimum)	132°C (270°F) (minimum)	132°C (270°F) (minimum)
Exposure Time	4 minutes (minimum)	4 minutes (minimum)	4 minutes (minimum)
Dry Time	65 minutes (minimum)	25 minutes (minimum)	N/A

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CONTRAINDICATIONS

- Do not use with attachments other than accessories provided by LSI SOLUTIONS®.
- This device is not intended for use except as indicated.

WARNINGS


- Federal law restricts this device to sale, distribution and use by, or on, the order of a physician.
- Read and become familiar with all instructions, warnings, and cautions before using this product.
- The *miniARM*® *RAIL CLAMP* shall be used in accordance with these instructions for use.
- Improper use of the *miniARM*® *RAIL CLAMP* may cause serious injury. In addition, improper care and maintenance of the device may render the device non-sterile prior to patient use and may cause serious injury to the health care provider or the patient.
- When using the *miniARM*® *RAIL CLAMP*, patients must be immobilized or anesthetized.
- Discontinue use of the *miniARM*® *RAIL CLAMP* when moving patient or when patient is moving.
- Surgical or endoscopic procedures should be performed only by physicians having adequate training and familiarity with relevant techniques and anatomy. Medical literature relative to techniques, complications, and hazards should be consulted prior to use.
- Avoid touching with sterile gloved hands any part of the rail clamp that extends below the level considered sterile in your operating room. For example, many hospitals consider any item below a gowned person's waist to be no longer in the sterile field.
- Treat all breaches (e.g., holes, tears and perforations) in the protective sterile drape as a site that is no longer sterile. Any section of a breached drape or part of a rail clamp that contacts a rail under a breached drape should be treated as contaminated. Avoid contact by sterile gloves, instruments, cords, tubes, etc. with any potential contamination source.

PRECAUTIONS

- **The *miniARM*® *RAIL CLAMP* is packaged as non-sterile.** Clean and sterilize prior to use.
- Ensure the patient is well secured to the surgical table and never in direct contact with the rail clamp.
- If there are any variations between these instructions for use and either your facility's policies and/or your cleaning/sterilizing equipment manufacturer's instructions, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your device.
- Use of a *miniARM*® *RAIL CLAMP* for a task other than what it is intended for can result in a damaged or broken device.
- Prior to use, inspect the *miniARM*® *RAIL CLAMP* to ensure proper function and condition. Do not use the device if it does not satisfactorily perform its intended function or if it has physical damage.
- Surgical instruments vary between manufacturers. Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure electrical isolation or grounding are not compromised.
- Avoid mechanical shock or overstressing the *miniARM*® *RAIL CLAMP*.
- Avoid inappropriate contact with the crenellations within the DOCs, which could damage gloves or cause injury. Glove perforation can lead to contamination of the surgical field.
- Keep the DOC cover(s) within the sterile field for reinstallation if desired.
- Only the cleaning and sterilization processes which are defined within these instructions for use have been validated.
- Check stability of surgical table rails or rail adaptors before table mounting the *miniARM*® *RAIL CLAMP*. Only mount to secure surgical tables and well-affixed horizontal table rails; avoid attaching rail clamps to structures that are loose or not secure.
- The *miniARM*® *RAIL CLAMP* is compatible with USA, EU and Swiss operating room table rail sizes, 9–10 mm wide x 25–30 mm tall (0.35–0.39 in. wide x 0.98–1.18 in. tall).
- Store at room temperature.

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* The *multiSTATION*™ System Sterilization Tray is manufactured by Summit Medical, 815 Vikings Parkway, Suite 100, St. Paul, MN 55121 U.S.A.

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EC REP

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

Polystica is a registered trademark of Steris. Weiman is a registered trademark of PurposeBuilt Brands. Use of these names, trademarks and brands does not imply endorsement. All company, product and service names used in this technology guide are for identification purposes only.

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MADE IN THE USA

This Product Comes with our LSI SOLUTIONS®
Perfect Performance Policy®
Call us at 866.575.3493 any time.

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The LSI logo, LSI Solutions, Multistation, Miniarm, Monodoc, Duodoc, Tridoc, M45, Click Release, M55 and Perfect Performance Policy are registered trademarks and trademarks of LSI SOLUTIONS®.

Copyright © 2021, LSI SOLUTIONS®. All Rights Reserved.

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
FIG. 1 – *miniARM®* M45™ CLICK RELEASE™ *miniARM®* HOLDER

 SOLUTIONS®

M45™ CLICK RELEASE™ *miniARM®* HOLDER DEVICE DESCRIPTION:

The M45™ CLICK RELEASE™ *miniARM®* HOLDER (FIG. 1) is an accessory in the *multiSTATION™* System. It is a reusable, sterilizable, repositionable support member providing an exceptionally stable base, when used with a separately available *miniARM®* RAIL CLAMP, on which interchangeable adapters are connected. Adapters and devices can be attached at the quick connect post 1 as needed for surgery. The upper joint 2 allows the quick connect post to be positioned accurately relative to the upper arm 3. With the first manual squeeze of the upper arm's CLICK RELEASE™ lever 4, the integrated lever latch 5 (not shown) connects with the locking catch 6 and locks the holder in the selected orientation. The next squeeze of the CLICK RELEASE™ lever releases the holder for repositioning. A subsequent lever squeeze again locks the holder; this locking and unlocking pattern alternates with each lever squeeze. The middle joint 7 allows for the articulation of the upper arm relative to the lower arm 8. The lower joint 9 allows for the articulation of the lower arm relative to the angled connector 10, which terminates in a DOC ball 11 with a spherical articulating surface that engages the receiving socket of the *miniARM®* RAIL CLAMP.

The attachment nut 12 compresses the upper hemisphere of the DOC ball into the receiving socket when the attachment nut and its internal threads 13 (not shown) are rotated clockwise, threading into the receiving socket. As it is tightened into the socket, the DOC ball and the angled connector are compressed so that they can no longer move. The two lock wings 14 engraved with the lock symbol () can swivel approximately 90° away from their storage positions on the attachment nut to ease clockwise threading. Once the attachment nut is tightened, the lock wings can be swiveled back into their storage position. To loosen the attachment nut, both unlock wings 15 labeled with the unlock symbol () can swivel from their storage position to ease counterclockwise threading.

INDICATIONS FOR USE:


The M45™ CLICK RELEASE™ *miniARM®* HOLDER is indicated for use by surgeons to hold instruments in a fixed position for a period of time.

INSTRUCTIONS FOR USE

LOCKING, RELEASING, and PLACING the M45™ CLICK RELEASE™ miniARM® HOLDER


Prior to each use, test the M45™ CLICK RELEASE™ miniARM® HOLDER by confirming that acceptable holding strength is achieved. Discontinue use if holding strength is inadequate. These instructions for attaching an M45™ CLICK RELEASE™ miniARM® HOLDER to a receiving DOC apply to use with the receiving DOC of any separately available miniARM® RAIL CLAMP already attached to a draped table rail.

1 LOCK

1. **LOCK** the holder by squeezing the CLICK RELEASE™ lever down toward the lever catch. The CLICK RELEASE™ lever hinges at the middle joint until the CLICK RELEASE™ lever is in its locked position, parallel to the upper arm. An audible click indicates that the lever latch and locking catch are engaged.


2 RELEASE

2. **RELEASE** the holder by fully squeezing the CLICK RELEASE™ lever down toward the locking catch. An audible click indicates that the lever latch is now disengaged, allowing the CLICK RELEASE™ lever to hinge at the middle joint to reach its unlocked position. If necessary, manually raise the CLICK RELEASE™ lever away from the upper arm to the unlocked position.


NOTE: When pivoting the CLICK RELEASE™ lever from the locked to unlocked position, make sure to support the holder.

3 PLACE

3. **PLACE** the attachment nut and the DOC ball onto the receiving DOC of the miniARM® RAIL CLAMP to attach the holder to the DOC.


4 ROTATE

4. **ROTATE** the attachment nut onto the DOC using both opposing lock wings () to turn the attachment nut clockwise, threading it onto RAIL CLAMP DOC.

NOTE: If threading is difficult, turn the attachment nut counterclockwise to unscrew, then carefully realign the attachment nut with the DOC and repeat the clockwise rotation using the lock wings. Do not force the attachment nut to rotate, as doing so can cause the holder and RAIL CLAMP to bind.

5 POSITION


5. **POSITION** the holder while in its unlocked state by pushing, pulling, or rotating the quick connect post, upper arm, lower arm, and angled connector about the three joints and DOC ball for adjustability in all directions. When the holder has been positioned as desired, lock the system by squeezing the CLICK RELEASE™ lever onto the locking latch.

ATTACHING AN ADAPTER TO THE M45™ CLICK RELEASE™ miniARM® HOLDER


Steps 1–3 show one example of how *miniARM®* adapters (separately available) are attached and adjusted to the M45™ CLICK RELEASE™ *miniARM®* HOLDER. The adapter may be attached to the M45™ CLICK RELEASE™ *miniARM®* HOLDER either before or after attaching the holder to a *miniARM®* RAIL CLAMP.

1 SLIDE


1. **SLIDE** the desired adapter onto the quick connect post while the adapter locking lever is rotated to the unlocked position.

2 ROTATE

2. **ROTATE** the adapter locking lever to the locked position to secure the adapter onto the holder.

3 ADJUST


3. **ADJUST** the upper joint of the M45™ CLICK RELEASE™ *miniARM®* HOLDER manually for minor repositioning without unlocking the holder.

REMOVAL and CLEANUP


Removal of the *M45™ CLICK RELEASE™ miniARM® HOLDER* may be eased by keeping the holder in its locked position.

1 ROTATE


1. **ROTATE** the locking lever of the attached adapter to the unlocked position.

2 DETACH


2. **DETACH** the adapter by sliding it off the quick connect post.

3 ROTATE

3. **ROTATE** the attachment nut off the *miniARM® RAIL CLAMP DOC* using both opposing unlock wings (). Turn the attachment nut counterclockwise to unscrew from the DOC.

4 REMOVE


4. **REMOVE** the *M45™ CLICK RELEASE™ miniARM® HOLDER* from the RAIL CLAMP DOC.

M45™ CLICK RELEASE™ miniARM® HOLDER REPROCESSING

- Disassemble and clean the device immediately after use. Do not allow a soiled device to dry.
- The device is not validated to be cleaned or sterilized with an adapter and/or a rail clamp attached.
- Cleaning agent used in validation: Steris Prolystica 2X (enzymatic, neutral pH).
- Perform the final rinse using only freshly prepared purified water/highly purified water.
- Never use metal brushes or steel wool for cleaning.
- Prepare and reprocess any adapters or rail clamps according to each device's instruction for use.

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

PREPARATION

1. Remove the M45™ CLICK RELEASE™ miniARM® HOLDER from the rail clamp by rotating the unlock wings away from their storage position and then turning the attachment nut counterclockwise (FIG. 2).



2. Remove any adapters from the M45™ CLICK RELEASE™ miniARM® HOLDER by rotating the adapter locking lever and sliding the adapter off the quick connect post (FIG. 3).


NOTE: FIG. 3 shows the cannula adapter being removed; other adapters are removed in the same manner.


3. If the CLICK RELEASE™ lever is in the locked position (FIG. 4), release the lock by fully squeezing the lever (FIG. 5) until an audible click indicates that the lever latch is disengaged from the locking catch. If necessary, manually raise the CLICK RELEASE™ lever away from the upper arm.

4. Raise the CLICK RELEASE™ lever completely to ensure the M45™ CLICK RELEASE™ miniARM® HOLDER is fully open and in the reprocessing position (FIG. 6).


5. Rotate both sets of lock and unlock wings away from their storage positions to their open positions (FIG. 7).

MANUAL PRE-CLEANING

1. Soak in enzymatic, neutral pH cleaning solution for a minimum of 5 minutes. Refer to detergent manufacturer's instructions for temperature and concentration.
2. Use a plastic-bristled brush to thoroughly scrub the device exterior surfaces, paying special attention to the following:
 - a. Retaining ring and upper ball joint (FIG. 8)
 - b. CLICK RELEASE™ lever latch (FIG. 9) and locking catch (FIG. 10)
 - c. CLICK RELEASE™ lever pivot point (FIG. 11)
 - d. DOC Ball, attachment nut, and internal threads (FIG.12)
3. Rinse with warm water for a minimum of 2 minutes.
4. Proceed with one of the two required cleaning options: ultrasonic or automated.

ULTRASONIC CLEANING


1. Clean in ultrasonic bath with enzymatic, neutral pH cleaning solution for a minimum of 15 minutes. Refer to detergent manufacturer's instructions for temperature and concentration.
2. Rinse with warm water for a minimum of 4 minutes. To ensure a complete rinse, use a clean plastic-bristled brush to scrub the device.


OR

AUTOMATED CLEANING

1. A washer-disinfector with fundamentally approved efficiency (e.g., according to EN ISO 15883) is required and it must be properly installed, qualified, and regularly subjected to maintenance and testing.
2. Load the device into the washer-disinfector. Avoid contact between devices and arrange to allow for proper drainage.
3. Operate the washer-disinfector cycle with an additional rinse cycle.

NOTE: The sterilization tray is NOT designed for cleaning devices. It must be processed separately. The tray is only intended for sterilization, transport, and storage of reusable instruments. For more tray information, see the sterilization tray instructions for use.

↓ LUBRICATION	1. Apply instrument lubricant mixed to manufacturer's recommendations to prolong instrument life by submerging the entire device in the lubricant for a minimum of 30 seconds. If the hospital washer-disinfector has a lubrication cycle this can be used instead of manual lubrication. NOTE: LSI has validated the use of Weiman® Instrument Lubricant on this device. Other instrument lubricant brands have not been tested and performance and results cannot be guaranteed.
↓ INSPECTION	1. Carefully inspect the device to assure that all visible soil has been removed. Generally, unmagnified visual inspection under good light conditions is sufficient. Particular attention should be paid to the locations in FIGS. 8-12. Repeat cleaning process if soil is detected. 2. Visually inspect the device for mild or excessive corrosion. If corrosion is present, discontinue use of the device in surgery, but complete reprocessing. 3. Visually inspect the device for damaged latches or other parts. For example, see the broken lever latch outline of FIG. 13 and the bent lever latch of FIG. 14. If parts are damaged, discontinue use of the device in surgery, but complete reprocessing. 4. Perform a functional check for adequate device lock and release.  FIG. 13  FIG. 14
↓ PACKAGING	1. Ensure the <i>M45™ CLICK RELEASE™ miniARM® HOLDER</i> is in the fully open reprocessing position (refer to FIG. 6). 2. If using a sterilization pouch, place each instrument in its own individual pouch. If using a sterilization tray, read the sterilization tray instructions for use before proceeding. Ensure that the sterilization tray has been cleaned according to the sterilization tray instructions for use and load the tray bottom according to FIG. 15. 3. Package the device according to TABLE 1. The barrier system for sterilized re-usable instruments should meet the following requirements: • ISO 11607-1 • Suitable for pre-vacuum steam sterilization • Appropriate for medical use • Grade appropriate for weight of loaded tray per sterilization tray instructions for use and facility procedures  FIG. 15

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
TABLE 1: M45™ CLICK RELEASE™ miniARM® HOLDER Sterile Packaging and Processing			
Method	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Immediate use steam sterilization according to ANSI/AAMI ST79
Container	<i>multiSTATION™ System Sterilization Tray</i> P/N 100034	No tray	No tray
Cycle	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)
Packaging	2-layer polypropylene wrap	Pouch	No packaging
Temperature	132°C (270°F) (minimum)	132°C (270°F) (minimum)	132°C (270°F) (minimum)
Exposure Time	4 minutes (minimum)	4 minutes (minimum)	4 minutes (minimum)
Dry Time	65 minutes (minimum)	25 minutes (minimum)	N/A

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The M45™ CLICK RELEASE™ miniARM® HOLDER has been validated to 100 reprocessing cycles. The useful lifespan of a surgical instrument is largely dependent on the care and handling of the instrument. Careful inspection and functional testing of the instrument should be used to determine the end of its serviceable life.

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WARNINGS

- Federal law restricts this device to sale, distribution and use by, or on, the order of a physician.
- Read and become familiar with all instructions, warnings, and cautions before using this product.
- The M45™ CLICK RELEASE™ miniARM® HOLDER shall be used in accordance with these instructions for use.
- Improper use of the M45™ CLICK RELEASE™ miniARM® HOLDER may cause serious injury. In addition, improper care and maintenance of the device may render the device non-sterile prior to patient use and may cause serious injury to the health care provider or the patient.
- When using the M45™ CLICK RELEASE™ miniARM® HOLDER, patients must be immobilized or anesthetized.
- Discontinue use of the M45™ CLICK RELEASE™ miniARM® HOLDER when moving the patient or when the patient is moving.
- Surgical or endoscopic procedures should be performed only by physicians having adequate training and familiarity with relevant techniques and anatomy. Medical literature relating to techniques, complications, and hazards should be consulted prior to use.


PRECAUTIONS

- The **M45™ CLICK RELEASE™ miniARM® HOLDER** is packaged as **non-sterile**. Clean and sterilize prior to use.
- If there are any variations between these instructions for use and either your facility's policies and/or your cleaning/sterilizing equipment manufacturer's instructions, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your device.
- Use of the **M45™ CLICK RELEASE™ miniARM® HOLDER** for a task other than what it is intended for can result in a damaged or broken device.
- Prior to use, inspect the **M45™ CLICK RELEASE™ miniARM® HOLDER** to ensure proper function and condition. Do not use devices if they do not satisfactorily perform their intended function or if they have physical damage.
- Surgical instruments vary between manufacturers. Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure electrical isolation or grounding are not compromised.
- Avoid mechanical shock or overstressing the **M45™ CLICK RELEASE™ miniARM® HOLDER**.
- Only the cleaning and sterilization processes which are defined within these instructions for use have been validated.
- Check stability of surgical table accessory rails or rail adapters before table mounting the **M45™ CLICK RELEASE™ miniARM® HOLDER**. Only mount to secure surgical tables and well-affixed horizontal table rails; avoid attaching rail clamps to structures that are loose or not secure.
- When pivoting the **CLICK RELEASE™** lever from the locked to unlocked position, make sure to support the holder. An unlocked and unsupported holder can fall and may cause harm, damage to the device, damage to the affixed **miniARM®** adapter, or injury to the patient or user.
- Store at room temperature.

ADVERSE REACTIONS

- No documented adverse reactions.

ORDERING INFORMATION

TABLE 2: M45™ CLICK RELEASE™ miniARM® HOLDER PRODUCT ORDERING			
	REORDER	PRODUCT	DESCRIPTION
	REF 081010	M45™ CLICK RELEASE™ miniARM® HOLDER	1 Shelf Box
	REF 100034	multiSTATION™ System Sterilization Tray*	1 Shelf Box

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*The **multiSTATION™ System Sterilization Tray** is manufactured by Summit Medical, 815 Vikings Parkway, Suite 100, St. Paul, MN 55121 U.S.A.

MADE IN THE USA

This Product Comes with our LSI SOLUTIONS® **Perfect Performance Policy®**. Call us at 866.575.3493 any time.

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The LSI logo, LSI Solutions, Multistation, M45, Click Release, Miniarm, and Perfect Performance Policy are registered trademarks and trademarks of LSI SOLUTIONS®. Copyright © 2021, LSI SOLUTIONS®. All Rights Reserved.

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The LSI logo, LSI Solutions, Multistation, M45, Click Release, Miniarm, and Perfect Performance Policy are registered trademarks and trademarks of LSI SOLUTIONS®.
Copyright © 2021, LSI SOLUTIONS®. All Rights Reserved.

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
FIG. 1 – *M55™ RIGID miniARM® HOLDER*

 SOLUTIONS®

***M55™ RIGID miniARM® HOLDER* DEVICE DESCRIPTION:**

The *M55™ RIGID miniARM® HOLDER* (FIG. 1) is an accessory in the *multiSTATION™* System. It is a reusable, sterilizable, repositionable support member providing an exceptionally stable base, when used with a separately available *miniARM® RAIL CLAMP*, on which interchangeable adapters are connected. Adapters and devices can be attached at the quick connect post **1** as needed for surgery. The upper joint **2** allows the quick connect post to be positioned accurately relative to the upper arm **3**. Clockwise rotation of the adjustment knob **4** secures the holder in the selected orientation. Counterclockwise rotation of the adjustment knob releases the holder for repositioning. The middle joint **5** allows for both the articulation and secure positioning, through the locking teeth **6**, of the upper arm relative to the lower arm **7**. The lower joint **8** allows for the articulation and secure positioning, through the locking dimples (not shown) of the lower arm relative to the angled connector **9**. The angled connector terminates in a DOC coupler **10** that engages the receiving socket of the *miniARM® RAIL CLAMP*.

The attachment nut **11** locks the DOC coupler into the receiving socket when the attachment nut and its internal threads **12** (not shown) are rotated clockwise, threading into the receiving socket. As it is tightened into the socket, the DOC coupler and the angled connector are secured so that they can no longer move.

INDICATIONS FOR USE:


The *M55™ RIGID miniARM® HOLDER* is indicated for use by surgeons to hold instruments in a fixed position for a period of time.

INSTRUCTIONS FOR USE

LOCKING, RELEASING, and PLACING the M55™ RIGID miniARM® HOLDER

Prior to each use, test the M55™ RIGID miniARM® HOLDER by confirming that acceptable holding strength is achieved. Discontinue use if holding strength is inadequate. These instructions for attaching an M55™ RIGID miniARM® HOLDER to a receiving DOC apply to use with the receiving DOC of any separately available miniARM® RAIL CLAMP already attached to a draped table rail.


1 LOCK

1. **LOCK** the holder by rotating the adjustment knob clockwise until fully seated.

NOTE: If any undesirable movement of holder or attachment is observed, further tighten adjustment knob clockwise.



2 RELEASE

2. **RELEASE** the holder by rotating the adjustment knob counterclockwise.


NOTE: When rotating the adjustment knob, make sure to support the holder.

3 PLACE


3. **PLACE** the attachment nut and the DOC coupler onto the receiving DOC of the miniARM® RAIL CLAMP to attach the holder to the DOC.

4 ROTATE


4. **ROTATE** the attachment nut onto the DOC clockwise, threading it onto RAIL CLAMP DOC.

NOTE: If threading is difficult, turn the attachment nut counterclockwise to unscrew, then carefully realign the attachment nut with the DOC and repeat the clockwise rotation. Do not force the attachment nut to rotate, as doing so can cause the holder and RAIL CLAMP to bind.

5. ALIGN the locking teeth of the middle joint and ensure the lower arm is fully seated in a locking dimple of the lower joint.

NOTE: Full engagement of the lower arm into a locking dimple of the lower joint can be felt as the holder is positioned.


6. POSITION the holder while in its unlocked state by pushing, pulling, or rotating the quick connect post, upper arm, lower arm about the three joints for adjustability in all directions.


When the holder has been positioned as desired with the locking teeth aligned and the locking dimple seated, lock the system by rotating the adjustment knob clockwise until sufficient holding force has been achieved.

ATTACHING AN ADAPTER TO THE M55™ RIGID miniARM® HOLDER

Steps 1 and 2 show one example of how *miniARM*® adapters (separately available) are attached and adjusted to the M55™ RIGID *miniARM*® HOLDER. The adapter may be attached to the M55™ RIGID *miniARM*® HOLDER either before or after attaching the holder to a *miniARM*® RAIL CLAMP.

1. SLIDE the desired adapter onto the quick connect post while the adapter locking lever is rotated to the unlocked position.


2. ROTATE the adapter locking lever to the locked position to secure the adapter onto the holder.

REMOVAL and CLEANUP


Removal of the *M55™ RIGID miniARM® HOLDER* may be eased by keeping the holder in its locked position.

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

3 ROTATE

3. **ROTATE** the attachment nut off the *miniARM® RAIL CLAMP DOC*. Turn the attachment nut counterclockwise to unscrew from the DOC.

4 REMOVE

4. **REMOVE** the *M55™ RIGID miniARM® HOLDER* from the RAIL CLAMP DOC.

M55™ RIGID miniARM® HOLDER REPROCESSING

- Disassemble and clean the device immediately after use. Do not allow a soiled device to dry.
- The device is not validated to be cleaned or sterilized with an adapter and/or a rail clamp attached.
- Cleaning agent used in validation: Steris Prolystica 2X (enzymatic, neutral pH).
- Perform the final rinse using only freshly prepared purified water/highly purified water.
- Never use metal brushes or steel wool for cleaning.
- Prepare and reprocess any adapters or rail clamps according to each device's instructions for use.

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PREPARATION

1. Remove the M55™ RIGID miniARM® HOLDER from the rail clamp by turning the attachment nut counterclockwise (FIG. 2).
2. Remove any adapters from the M55™ RIGID miniARM® HOLDER by rotating the adapter locking lever and sliding the adapter off the quick connect post (FIG. 3).

NOTE: FIG. 3 shows the clamping adapter being removed; other adapters are removed in the same manner.

3. Rotate the adjustment knob counterclockwise until fully open such that the knob cannot be rotated anymore and every joint moves freely (FIG. 4).

4. If using the multiSTATION™ System Sterilization Tray (Part Number 100034), clean the tray separately according to the sterilization tray instructions for use.

FIG. 2


FIG. 3




FIG. 4


MANUAL PRE-CLEANING

1. Soak in enzymatic, neutral pH cleaning solution for a minimum of 5 minutes. Refer to detergent manufacturer's instructions for temperature and concentration.


2. Use a plastic-bristled brush to thoroughly scrub the device's exterior surfaces, paying special attention to the following:
 - a. Retaining ring, upper ball joint, and quick connect post (FIG. 5)
 - b. Locking teeth of the middle joint (FIG. 6)
 - c. Lower ball joint, locking dimples (FIG. 7)
 - d. DOC coupler, attachment nut, and internal threads (FIG. 8)

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
TABLE 1: M55™ RIGID miniARM® HOLDER Sterile Packaging and Processing			
Method	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Immediate use steam sterilization according to ANSI/AAMI ST79
Container	multiSTATION™ System Sterilization Tray P/N 100034	No tray	No tray
Cycle	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)
Packaging	2-layer polypropylene wrap	Pouch	No packaging
Temperature	132°C (270°F) (minimum)	132°C (270°F) (minimum)	132°C (270°F) (minimum)
Exposure Time	4 minutes (minimum)	4 minutes (minimum)	4 minutes (minimum)
Dry Time	65 minutes (minimum)	25 minutes (minimum)	N/A

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The M55™ RIGID miniARM® HOLDER has been validated to 100 reprocessing cycles. The useful lifespan of a surgical instrument is largely dependent on the care and handling of the instrument. Careful inspection and functional testing of the instrument should be used to determine the end of its serviceable life.

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WARNINGS

- Federal law restricts this device to sale, distribution and use by, or on, the order of a physician.
- Read and become familiar with all instructions, warnings, and cautions before using this product.
- The M55™ RIGID miniARM® HOLDER shall be used in accordance with these Instructions for Use.
- Improper use of the M55™ RIGID miniARM® HOLDER may cause serious injury. In addition, improper care and maintenance of the device may render the device non-sterile prior to patient use and may cause serious injury to the health care provider or the patient.
- When using the M55™ RIGID miniARM® HOLDER, patients must be immobilized or anesthetized.
- Discontinue use of the M55™ RIGID miniARM® HOLDER when moving the patient or when the patient is moving.
- Surgical or endoscopic procedures should be performed only by physicians having adequate training and familiarity with relevant techniques and anatomy. Medical literature relating to techniques, complications, and hazards should be consulted prior to use.


PRECAUTIONS

- The **M55™ RIGID miniARM® HOLDER** is packaged as non-sterile. Clean and sterilize prior to use.
- If there are any variations between these instructions for use and either your facility's policies and/or your cleaning/sterilizing equipment manufacturer's instructions, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your device.
- Use of the **M55™ RIGID miniARM® HOLDER** for a task other than what it is intended for can result in a damaged or broken device.
- Prior to use, inspect the **M55™ RIGID miniARM® HOLDER** to ensure proper function and condition. Do not use devices if they do not satisfactorily perform their intended function or if they have physical damage.
- Surgical instruments vary between manufacturers. Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure electrical isolation or grounding are not compromised.
- Avoid mechanical shock or overstressing the **M55™ RIGID miniARM® HOLDER**.
- Only the cleaning and sterilization processes which are defined within these instructions for use have been validated.
- Check stability of surgical table accessory rails or rail adapters before table mounting the **M55™ RIGID miniARM® HOLDER**. Only mount to secure surgical tables and well-affixed horizontal table rails; avoid attaching rail clamps to structures that are loose or not secure.
- When loosening the adjustment knob, make sure to support the holder. An unlocked and unsupported holder can fall and may cause harm, damage to the device, damage to the affixed **miniARM®** adapter, or injury to the patient or user.
- Store at room temperature.

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
The LSI logo, LSI Solutions, Multistation, M55, Miniarm, and Perfect Performance Policy are registered trademarks and trademarks of LSI SOLUTIONS®.
Copyright © 2021, LSI SOLUTIONS®. All Rights Reserved.

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
FIG. 1 – *miniARM®* ACCESS ADAPTER

 **SOLUTIONS®**

***miniARM®* ACCESS ADAPTER DEVICE DESCRIPTION:**

The *miniARM®* ACCESS ADAPTER (FIG. 1) is used with the *M55™ RIGID miniARM® HOLDER* or the *M45™ CLICK RELEASE™ miniARM® HOLDER*. The threaded receiver 1 provides a connection point onto which the mounting screw of a compatible LSI device is attached (refer to the Technology Guide of the selected device for more information before using the adapter). Rotating the locking lever 2 to the unlocked position allows the adapter to slide over the quick connect post of a *miniARM® HOLDER*. The adapter is locked into place by rotating the locking lever to the locked position. The adapter is shown attached to a holder with a compatible LSI device affixed to the threaded receiver in FIG. 2.

INDICATIONS FOR USE:

The *miniARM®* ACCESS ADAPTER is indicated for use by surgeons to hold instruments in a fixed position for a period of time.



FIG. 2 – *miniARM®* ACCESS ADAPTER 1 attached to a *miniARM®* HOLDER 2 with the *RAM®* Ring 3 attached to the threaded receiver.

FIG. 2


INSTRUCTIONS FOR USE

See the respective LSI SOLUTIONS® Technology Guides for full instructions for use and adjustment of LSI accessories in conjunction with the *miniARM® ACCESS ADAPTER*.


The *miniARM® ACCESS ADAPTER* may be attached to the *miniARM® HOLDER* either before or after attaching the holder to a *miniARM® RAIL CLAMP*.

1. SLIDE the adapter onto the quick connect post of the holder while the locking lever is rotated to the unlocked position.

2. ROTATE the locking lever to the locked position to secure the adapter onto the holder.


3. CONNECT the mounting screw of a compatible LSI accessory (RAM® Ring is shown here) to the adapter threaded receiver, using the 1/8" hex key (included with the RAM® Ring). Rotate the mounting screw clockwise to tighten..

REMOVAL AND CLEANUP

DISASSEMBLY


The *RAM® Ring* is a single-use device that should be disposed of in accordance with local regulations and hospital procedures. The *miniARM® ACCESS ADAPTER* can be reprocessed with the *multiSTATION™ System* using the instructions in this technology guide.


1. **DISCONNECT** the mounting screw of the LSI accessory (*RAM® Ring* is shown here) from the adapter threaded receiver, using the 1/8" hex key. Rotate the mounting screw counterclockwise.

REMOVING THE *miniARM® ACCESS ADAPTER* FROM THE *miniARM® HOLDER*

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
1. **ROTATE** the locking lever to the unlocked position.

2. **DETACH** the adapter by sliding it off the holder's quick connect post.

miniARM® ACCESS ADAPTER REPROCESSING

- Disassemble and clean the device immediately after use. Do not allow a soiled device to dry.
- The device is not validated to be cleaned or sterilized with a holder and/or rail clamp attached.
- Cleaning agent used in validation: Steris Prolystica 2X (enzymatic, neutral pH).
- Perform the final rinse using only freshly prepared purified water/highly purified water.
- Never use metal brushes or steel wool for cleaning.
- Prepare and reprocess any holders or rail clamps according to each device's instructions for use.

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LUBRICATION

1. Apply instrument lubricant mixed to manufacturer's recommendations to prolong instrument life by submerging the entire device in the lubricant for a minimum of 30 seconds. If the hospital washer-disinfector has a lubrication cycle this can be used instead of manual lubrication.

NOTE: LSI has validated the use of Weiman® Instrument Lubricant on this device. Other instrument lubricant brands have not been tested and performance and results cannot be guaranteed.

INSPECTION

1. Carefully inspect the device to assure that all visible soil has been removed. Generally, unmagnified visual inspection under good light conditions is sufficient. Particular attention should be paid to the locations in FIG. 4. Repeat cleaning process if soil is detected.

2. Visually inspect the device for mild or excessive corrosion. If corrosion is present, discontinue use of the device in surgery, but complete reprocessing.

3. Visually inspect the device for any damage. If parts are damaged, discontinue use of the device in surgery, but complete reprocessing.

PACKAGING

1. If using a sterilization pouch, place each instrument in its own individual pouch. If using a sterilization tray, read the sterilization tray instructions for use before proceeding. Ensure that the sterilization tray has been cleaned according to the sterilization tray instructions for use and load the tray bottom according to FIG. 5.

2. Package the device according to TABLE 1. The barrier system for sterilized re-usable instruments should meet the following requirements:

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
TABLE 1: <i>miniARM</i>® ACCESS ADAPTER Sterile Packaging and Processing			
Method	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Immediate use steam sterilization according to ANSI/AAMI ST79
Container	<i>multiSTATION</i> ™ System Sterilization Tray P/N 100034	No tray	No tray
Cycle	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)
Packaging	2-layer polypropylene wrap	Pouch	No packaging
Temperature	132°C (270°F) (minimum)	132°C (270°F) (minimum)	132°C (270°F) (minimum)
Exposure Time	4 minutes (minimum)	4 minutes (minimum)	4 minutes (minimum)
Dry Time	65 minutes (minimum)	25 minutes (minimum)	N/A

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The *miniARM*® ACCESS ADAPTER has been validated to 20 reprocessing cycles. The useful lifespan of a surgical instrument is largely dependent on the care and handling of the instrument. Careful inspection and functional testing of the instrument should be used to determine the end of its serviceable life.

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WARNINGS

- Federal law restricts this device to sale, distribution, and use by, or on, the order of a physician.
- Read and become familiar with all instructions, warnings, and cautions before using this product.
- The *miniARM*® ACCESS ADAPTER shall be used in accordance with these instructions for use.
- Improper use of the *miniARM*® ACCESS ADAPTER may cause serious injury. In addition, improper care and maintenance of the device may render the device non-sterile prior to patient use and may cause serious injury to the health care provider or the patient.
- When using the *miniARM*® ACCESS ADAPTER, patients must be immobilized or anesthetized.
- Discontinue use of the *miniARM*® ACCESS ADAPTER when patient is moving or being moved.
- Surgical or endoscopic procedures should be performed only by physicians having adequate training and familiarity with relevant techniques and anatomy. Medical literature relating to techniques, complications, and hazards should be consulted prior to use.

PRECAUTIONS

- **The *miniARM*® ACCESS ADAPTER is packaged as non-sterile.** Clean and sterilize prior to use.
- If there are any variations between these instructions for use and the policies of your facility and/or the instructions of your cleaning/sterilizing equipment manufacturer, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your device.

PRECAUTIONS (CONTINUED)


- Use of the *miniARM® ACCESS ADAPTER* for a task other than what it is intended for can result in a damaged or broken device.
- Prior to use, inspect the *miniARM® ACCESS ADAPTER* to ensure proper function and condition. Do not use devices if they do not satisfactorily perform their intended function or if they have physical damage.
- Surgical instruments vary between manufacturers. Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure electrical isolation or grounding are not compromised.
- Avoid mechanical shock or overstressing the *miniARM® ACCESS ADAPTER*.
- Only the cleaning and sterilization processes that are defined within these instructions for use have been validated.
- Store at room temperature.

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*The *multiSTATION™ System Sterilization Tray* is manufactured by Summit Medical, 815 Vikings Parkway, Suite 100, St. Paul, MN 55121 U.S.A.

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The LSI logo, LSI Solutions, Multistation, M45, M55, Click Release, Miniarm, Ram, and Perfect Performance Policy are registered trademarks and trademarks of LSI SOLUTIONS®. Copyright © 2021, LSI SOLUTIONS®. All Rights Reserved.

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
 SOLUTIONS®

DEVICE DESCRIPTION:

The *miniARM® CANNULA ADAPTER KIT* (FIG. 1) attaches to the *M55™ RIGID miniARM® HOLDER* or the *M45™ CLICK RELEASE™ miniARM® HOLDER*. The *miniARM® CANNULA ADAPTER KIT* has a miniARM® cannula adapter **1** that slides over the quick connect post of the miniARM® Holder and locks into place with the locking lever **4**. A miniARM® cannula **5**, for use with compatible endoscopes, fits into the miniARM® cannula adapter, which has a cannula latch **3** that secures the cannula. When depressed, the cannula latch releases the cannula. A miniARM® obturator **8** is placed into the proximal cannula opening to enable the cannula to be passed through tissue exposed by a skin incision. After the miniARM® obturator is removed, a scope may be placed into the cannula and held in place with the movable cannula ring **6**. A rotation dial **2** enables the cannula and scope to be rotated as desired.


INDICATIONS FOR USE:

The *miniARM® CANNULA ADAPTER KIT* is indicated for use by surgeons to hold instruments in a fixed position for a period of time.

miniARM® CANNULA ADAPTER ASSEMBLY AND POSITIONING


The miniARM® cannula adapter may be attached to the *miniARM® HOLDER* either before or after attaching the holder to a *miniARM® RAIL CLAMP*.

1 SLIDE

1. **SLIDE** the adapter onto the quick connect post of the holder while the locking lever is rotated to the unlocked position.

2 ROTATE


2. **ROTATE** the locking lever to the locked position to secure the adapter onto the holder.

3 POSITION


3. **POSITION** the adapter so that the word "SURGEON" is facing the user.

4 INSERT


4. **INSERT** the cannula from the "SURGEON" side until it is fully seated and secured by the cannula latch.

NOTE: cannula should be fully seated


If a gap is visible, cannula is NOT fully seated


miniARM® OBTURATOR USE:

1 ROTATE

1. **ROTATE** the cannula ring clockwise until the gap in the cannula ring aligns with the notch in the cannula.


2 INSERT

2. **INSERT** the obturator into the cannula with the obturator depth stop passing through the gap in the cannula ring into the notch in the cannula.

miniARM® OBTURATOR USE (continued):


3 **SECURE**

3. **SECURE** the obturator in the cannula by rotating the cannula ring counter-clockwise over the obturator depth stop.


ENDOSCOPE USE:

1 **ROTATE**


1. **ROTATE** the cannula ring clockwise until the gap in the cannula ring aligns with the notch in the cannula.

2 **INSERT**


2. **INSERT** the endoscope into the cannula with the light post passing through the gap in the cannula ring into the notch in the cannula.

3 **SECURE**


3. **SECURE** the endoscope in the cannula by rotating the cannula ring counter-clockwise over the light post.

4 **ROTATE**

4. **ROTATE** the endoscope with the cannula rotation dial to achieve desired viewing orientation.

5 **PULL**


5. **PULL** the rotation dial toward surgeon side to lock the rotation dial; push the rotation dial away from the surgeon side to unlock the rotation dial.

REMOVAL AND CLEANUP


This process may be made easier by having the *miniARM® HOLDER* locked while handling.

1 ROTATE


1. **ROTATE** the cannula ring clockwise until the gap in the cannula ring aligns with the notch in the cannula to release the obturator or endoscope.

2 REMOVE


2. **REMOVE** the obturator or endoscope.

3 REMOVE


3. **REMOVE** the cannula by pressing the cannula latch and pulling the cannula out of the adapter.

4 ROTATE

4. **ROTATE** the locking lever to the unlocked position.







5 DETACH


5. **DETACH** the adapter by sliding it off the holder's quick connect post..

miniARM® CANNULA ADAPTER KIT REPROCESSING

- Disassemble and clean the device immediately after use. Do not allow a soiled device to dry.
- The device is not validated to be cleaned or sterilized with the adapter, obturator, or cannula still assembled.
- Cleaning agent used in validation: Steris Prolystica® 2X (enzymatic, neutral pH).
- Perform the final rinse using only freshly prepared purified water/highly purified water.
- Never use metal brushes or steel wool for cleaning.
- Prepare and reprocess any holder or rail clamp according to each device's instructions for use.

POINT OF USE	1. Disassemble and clean device immediately after use. 2. Do not allow soiled devices to dry.
PREPARATION	1. Rotate the cannula ring clockwise until the gap in the cannula ring aligns with the notch in the cannula to release the obturator or endoscope. (FIG. 2) 2. Remove the obturator or endoscope (FIG. 3) 3. Remove the cannula by pressing the cannula latch and pulling the cannula out of the adapter. (FIG. 4) 4. Rotate the locking lever to the unlocked position and slide the adapter off the quick connect post. (FIG. 5).    
MANUAL PRE-CLEANING	1. Soak in enzymatic, neutral pH cleaning solution for a minimum of 5 minutes. Refer to detergent manufacturer's instructions for temperature and concentration. 2. Use a plastic-bristled brush to thoroughly scrub the device's exterior surfaces, paying special attention to the following: a. Cannula adapter attachment opening (FIG. 6) b. Length of cannula lumen (FIG. 7) 3. Rinse with warm water for a minimum of 2 minutes. 4. Proceed with one of the two required cleaning options: ultrasonic or automated.  

miniARM® CANNULA ADAPTER KIT REPROCESSING (continued)

miniARM® CANNULA ADAPTER KIT REPROCESSING (continued)

PACKAGING

1. Confirm that all components are present. If using a sterilization pouch, place one disassembled device (cannula adapter, cannula, and obturator) in its own individual pouch. If using a sterilization tray, read the sterilization tray instructions for use before proceeding. See FIG. 8 and FIG. 9 for placement of components in tray. Ensure that the sterilization tray has been cleaned according to the sterilization tray instructions for use.

2. Package the device according to TABLE 1 below. The barrier system for sterilized reusable instruments should meet the following requirements:

- ISO 11607-1
- Suitable for pre-vacuum steam sterilization
- Appropriate for medical use
- Grade appropriate for weight of loaded tray per sterilization tray instructions for use and facility procedures


FIG. 8

FIG. 9


STERILIZATION

1. The device must be properly cleaned prior to sterilization.

2. Perform sterilization cycle according to TABLE 1:

TABLE 1: miniARM® CANNULA ADAPTER KIT Sterile Packaging and Processing			
Method	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Immediate Use Steam Sterilization according to ANSI/AAMI ST79
Container	multiSTATION™ System Sterilization Tray P/N 100034	No tray	No tray
Cycle	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)
Packaging	2-layer polypropylene wrap	Pouch	No packaging
Temperature	132°C (270°F) (minimum)	132°C (270°F) (minimum)	132°C (270°F) (minimum)
Exposure Time	4 minutes (minimum)	4 minutes (minimum)	4 minutes (minimum)
Dry Time	65 minutes (minimum)	25 minutes (minimum)	N/A

miniARM® CANNULA ADAPTER KIT REPROCESSING (continued)

CONTRAINDICATIONS

- This device is not intended for use except as indicated.

WARNINGS


- Federal law restricts this device to sale, distribution and use by, or on, the order of a physician.
- Read and become familiar with all instructions, warnings, and cautions before using this product.
- The *miniARM® CANNULA ADAPTER KIT* shall be used in accordance with these instructions for use.
- Improper use of the device may cause serious injury. In addition, improper care and maintenance of the device may render the device non-sterile prior to patient use and may cause serious injury to the health care provider or the patient.
- When using the *miniARM® CANNULA ADAPTER KIT*, patients must be immobilized or anesthetized.
- Discontinue use of the *miniARM® CANNULA ADAPTER KIT* when the patient is moving or being moved.
- Surgical or endoscopic procedures should be performed only by physicians having adequate training and familiarity with relevant techniques and anatomy. Medical literature relative to techniques, complications, and hazards should be consulted prior to use.
- When using the *miniARM®* cannula/obturator, an insufficient skin incision may cause increased penetration force which may reduce the surgeon's control during entry.

PRECAUTIONS

- **The device is packaged as non-sterile.** Clean and sterilize prior to use.
- If there are any variations between these instructions for use and your facility's policies and/or your cleaning/sterilizing equipment manufacturer's instructions, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your device.
- Use of a *miniARM® CANNULA ADAPTER KIT* for a task other than what it is intended for can result in a damaged or broken device.
- Prior to use, inspect the *miniARM® CANNULA ADAPTER KIT* to ensure proper function and condition. Do not use the device if it does not satisfactorily perform its intended function or if it has physical damage.
- Surgical instruments vary between manufacturers. Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure electrical isolation or grounding are not compromised.
- Avoid mechanical shock or overstressing the *miniARM® CANNULA ADAPTER KIT*.
- Always use caution when inserting or removing devices through the cannula. Lateral pressure on a cannulated instrument during removal can damage the working tip and/or shaft of the device. Be sure the tips are closed and the device is pulled straight out until completely clear of the cannula.
- Only the cleaning and sterilization processes which are defined within these instructions for use have been validated.
- Store at room temperature.

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MADE IN THE USA

This Product Comes with our LSI SOLUTIONS® *Perfect Performance Policy**
Call us at 866.575.3493 any time.

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The LSI logo, LSI Solutions, Multistation, Miniarm, M45, Click Release, M55 and Perfect Performance Policy are registered trademarks and trademarks of LSI SOLUTIONS®.
Copyright © 2021, LSI SOLUTIONS®. All Rights Reserved.

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
FIG. 1 – *miniARM*® CLAMPING ADAPTER

 **SOLUTIONS**®

***miniARM*® CLAMPING ADAPTER DEVICE DESCRIPTION:**

The *miniARM*® CLAMPING ADAPTER (FIG. 1) is used with the *M55*™ RIGID *miniARM*® HOLDER or the *M45*™ CLICK RELEASE™ *miniARM*® HOLDER. The adjustable jaw 1 is opened or closed by rotating the actuation thumb screw 2 counterclockwise or clockwise, respectively. Opening and closing the adjustable jaw allows the user to secure or release a desired surgical device. Rotating the locking lever 3 to the unlocked position allows the adapter to slide over the quick connect post of the holder. The adapter is locked into place by rotating the locking lever to the locked position. The adapter is designed to hold surgical instruments with a shaft diameter of 4-14 mm.


INDICATIONS FOR USE:

The *miniARM*® CLAMPING ADAPTER is indicated for use by surgeons to hold instruments in a fixed position for a period of time.

INSTRUCTIONS FOR USE


The *miniARM*® CLAMPING ADAPTER may be attached to the *miniARM*® HOLDER either before or after attaching the holder to a *miniARM*® RAIL CLAMP.

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

3 OPEN

3. **OPEN** the adjustable jaw by rotating the actuation thumb screw counterclockwise.

4 SECURE


4. **SECURE** a surgical device in the adapter by rotating the actuation thumb screw clockwise, thus closing the adjustable jaw.


NOTE: Do not use excessive force when securing a surgical device in the adjustable jaw.

REMOVAL and CLEANUP

This process may be made easier by having the *miniARM*® HOLDER locked while handling.

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


PACKAGING

1. If using a sterilization pouch, place each instrument in its own individual pouch. If using a sterilization tray, read the sterilization tray instructions for use before proceeding. Ensure that the sterilization tray has been cleaned according to the sterilization tray instructions for use and load the tray bottom according to FIG. 4.

2. Package the device according to TABLE 1. The barrier system for sterilized reusable instruments should meet the following requirements:

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TABLE 1: <i>miniARM® CLAMPING ADAPTER</i> Sterile Packaging and Processing			
Method	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Immediate use steam sterilization according to ANSI/AAMI ST79
Container	<i>multiSTATION™ System Sterilization Tray</i> P/N 100034	No tray	No tray
Cycle	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)
Packaging	2-layer polypropylene wrap	Pouch	No packaging
Temperature	132°C (270°F) (minimum)	132°C (270°F) (minimum)	132°C (270°F) (minimum)
Exposure Time	4 minutes (minimum)	4 minutes (minimum)	4 minutes (minimum)
Dry Time	65 minutes (minimum)	25 minutes (minimum)	N/A

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The *miniARM® CLAMPING ADAPTER* has been validated to 20 reprocessing cycles. The useful lifespan of a surgical instrument is largely dependent on the care and handling of the instrument. Careful inspection and functional testing of the instrument should be used to determine the end of its serviceable life.

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CONTRAINDICATIONS

- These devices are not intended for use except as indicated.

WARNINGS

- Federal law restricts this device to sale, distribution, and use by, or on, the order of a physician.
- Read and become familiar with all instructions, warnings, and cautions before using this product.
- The *miniARM® CLAMPING ADAPTER* shall be used in accordance with these instructions for use.
- Improper use of the *miniARM® CLAMPING ADAPTER* may cause serious injury. In addition, improper care and maintenance of the device may render the device non-sterile prior to patient use and may cause serious injury to the health care provider or the patient.
- When using the *miniARM® CLAMPING ADAPTER*, patients must be immobilized or anesthetized.
- Discontinue use of the *miniARM® CLAMPING ADAPTER* when the patient is moving or being moved.
- Surgical or endoscopic procedures should be performed only by physicians having adequate training and familiarity with relevant techniques and anatomy. Medical literature relating to techniques, complications, and hazards should be consulted prior to use.


PRECAUTIONS

- **The *miniARM® CLAMPING ADAPTER* is packaged as non-sterile.** Clean and sterilize prior to use.
- If there are any variations between these instructions for use and the policies of your facility and/or the instructions of your cleaning/sterilizing equipment manufacturer, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your device.
- Use of the *miniARM® CLAMPING ADAPTER* for a task other than what it is intended for can result in a damaged or broken device.
- Prior to use, inspect the *miniARM® CLAMPING ADAPTER* to ensure proper function and condition. Do not use devices if they do not satisfactorily perform their intended function or if they have physical damage.
- Surgical instruments vary between manufacturers. Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure electrical isolation or grounding are not compromised.
- Avoid mechanical shock or overstressing the *miniARM® CLAMPING ADAPTER*.
- Only the cleaning and sterilization processes that are defined within these instructions for use have been validated.
- Be sure to adjust the *miniARM® CLAMPING ADAPTER* actuation thumb screw enough to hold the device securely without overtightening. Overtightening the *miniARM® CLAMPING ADAPTER* can result in damage to the held device while insufficient tightening could result in device slippage.
- Store at room temperature.

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
MADE IN THE USA

This Product Comes with our LSI SOLUTIONS® **Perfect Performance Policy®** Call us at 866.575.3493 any time.

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The LSI logo, LSI Solutions, Perfect Performance Policy, M45, M55, Click Release, Miniarm, and Multistation are registered trademarks and trademarks of LSI SOLUTIONS®. Copyright © 2021, LSI SOLUTIONS®. All Rights Reserved.

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
MADE IN THE USA

This Product Comes with our LSI SOLUTIONS® **Perfect Performance Policy®**
Call us at 866.575.3493 any time.

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 SOLUTIONS®

DEVICE DESCRIPTION:

The *miniARM® LOCKING COUPLER 1* is used with either the *miniARM® CANNULA ADAPTER KIT* or the *miniARM® CLAMPING ADAPTER* to securely hold the Olympus® ENDOEYE FLEX® 5mm HD surgical video camera (LTF-S190-5). The camera is released from the coupler by pressing the release tab **2**. The *miniARM® CLAMPING ADAPTER* secures the Coupler by gripping the locking coupler shaft **3**.

INDICATIONS FOR USE:


The *miniARM® LOCKING COUPLER* is indicated for use by surgeons to hold instruments in a fixed position for a period of time.

miniARM® LOCKING COUPLER ASSEMBLY AND POSITIONING

The *miniARM® LOCKING COUPLER* may be attached to the *miniARM® HOLDER* using either the *miniARM® CLAMPING ADAPTER* or the *miniARM® CANNULA ADAPTER KIT*.


USE WITH THE CLAMPING ADAPTER:

1 OPEN

1. **OPEN** the clamping adapter jaw by rotating the actuation thumb screw counterclockwise.


2 SECURE

2. **SECURE** the coupler in the clamping adapter by rotating the actuation thumb screw clockwise to close the jaw around the coupler shaft


USE WITH THE CANNULA ADAPTER:

1 ORIENT


1. **ORIENT** the cannula adapter so that the word "SURGEON" is facing the user.

2 INSERT


2. **INSERT** the coupler from the "SURGEON" side until it is fully seated and secured by the cannula latch.

NOTE: After insertion, the locking coupler can be rotated and secured as desired. See the *miniARM® CANNULA ADAPTER KIT* Instructions for Use.

Coupler should be fully seated


If a gap is visible, the coupler is NOT fully seated

INSERTING THE OLYMPUS® ENDOEYE FLEX® SURGICAL VIDEO CAMERA:


Note: To avoid damage to the camera's distal tip, ensure that the camera is not flexed before insertion.

1 ORIENT


1. **ORIENT** the coupler so that the word "TOP" is facing the user.

2 INSERT


2. **INSERT** the camera into the coupler, making sure that the flat portion of the camera aligns with the flat portion of the coupler labeled "TOP," as shown.

DISASSEMBLY:


NOTE: To avoid damage to the camera's distal tip, ensure that the camera is not flexed before removal.

1 **PULL**

1. **PULL** the release tab to disengage the camera from the coupler, then carefully remove the camera.


2a. **PRESS**

IF USING THE CANNULA ADAPTER:

2a. **PRESS** the cannula latch while pulling the locking coupler from the cannula adapter.

2b. **OPEN**





IF USING THE CLAMPING ADAPTER:

2b. **OPEN** the clamping adapter jaw by rotating the actuation thumb screw counterclockwise to release the coupler.

miniARM® LOCKING COUPLER REPROCESSING

- Disassemble and clean the device immediately after use. Do not allow a soiled device to dry.
- The device is not validated to be cleaned or sterilized with a camera attached or while still attached to a cannula adapter or clamping adapter.
- Cleaning agent used in validation: Steris Prolystica® 2X (enzymatic, neutral pH).
- Perform the final rinse using only freshly prepared purified water/highly purified water.
- Never use metal brushes or steel wool for cleaning.

POINT OF USE	1. Disassemble and clean device immediately after use. 2. Do not allow soiled devices to dry.
PREPARATION	NOTE: To avoid damage to the camera's distal tip, ensure that the camera is not flexed before removal. 1. Pull the release tab to disengage the camera from the coupler, then carefully remove the camera. (FIG. 2) If using the Cannula Adapter: 2a. PRESS the cannula latch while pulling the coupler from the cannula adapter. (FIG. 3) If using the Clamping Adapter: 2b. OPEN the clamping adapter jaw by rotating the actuation thumb screw counterclockwise to release the coupler. (FIG. 4)
MANUAL PRE-CLEANING	1. Soak in enzymatic, neutral pH cleaning solution for a minimum of 5 minutes. Refer to detergent manufacturer's instructions for temperature and concentration. 2. Use a plastic-bristled brush to thoroughly scrub the device's exterior surfaces, paying special attention to the following (FIG. 5): a. Release tab pin pocket, b. Length of coupler lumen. Use round brush 1 mm larger than lumen (e.g., use 7 mm brush for 6 mm lumen). 3. Rinse with warm water for a minimum of 2 minutes. 4. Proceed with one of the two required cleaning options: ultrasonic or automated.

miniARM® LOCKING COUPLER REPROCESSING (continued)



FIG. 6

miniARM® LOCKING COUPLER REPROCESSING (continued)

STERILIZATION

1. The device must be properly cleaned prior to sterilization.
2. Perform sterilization cycle according to TABLE 1.

TABLE 1: miniARM® LOCKING COUPLER Sterile Packaging and Processing

Method	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Immediate Use Steam Sterilization according to ANSI/AAMI ST79
Container	multiSTATION™ System Sterilization Tray P/N 100034	No tray	No tray
Cycle	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)
Packaging	2-layer polypropylene wrap	Pouch	No packaging
Temperature	132°C (270°F) (minimum)	132°C (270°F) (minimum)	132°C (270°F) (minimum)
Exposure Time	4 minutes (minimum)	4 minutes (minimum)	4 minutes (minimum)
Dry Time	65 minutes (minimum)	25 minutes (minimum)	N/A

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The miniARM® LOCKING COUPLER has been validated to 20 reprocessing cycles. The useful lifespan of a surgical instrument is largely dependent on the care and handling of the instrument. Careful inspection and functional testing of the instrument should be used to determine the end of its serviceable life.

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WARNINGS

- Federal law restricts this device to sale, distribution and use by, or on, the order of a physician.
- Read and become familiar with all instructions, warnings, and cautions before using this product.
- The miniARM® LOCKING COUPLER shall be used in accordance with these instructions for use.
- Improper use of the device may cause serious injury. In addition, improper care and maintenance of the device may render the device non-sterile prior to patient use and may cause serious injury to the health care provider or the patient.
- When using the miniARM® LOCKING COUPLER, patients must be immobilized or anesthetized.
- Discontinue use of the miniARM® LOCKING COUPLER when the patient is moving or being moved.
- Surgical or endoscopic procedures should be performed only by physicians having adequate training and familiarity with relevant techniques and anatomy. Medical literature relative to techniques, complications, and hazards should be consulted prior to use.
- Do not use the miniARM® LOCKING COUPLER as a cannula.


PRECAUTIONS

- **The device is packaged as non-sterile.** Clean and sterilize prior to use.
- If there are any variations between these instructions for use and either your facility's policies and/or your cleaning/sterilizing equipment manufacturer's instructions, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your device.
- Use of a *miniARM® LOCKING COUPLER* for a task other than what it is intended for can result in a damaged or broken device.
- Prior to use, inspect the *miniARM® LOCKING COUPLER* to ensure proper function and condition. Do not use the device if it does not satisfactorily perform its intended function or if it has physical damage.
- Surgical instruments vary between manufacturers. Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure electrical isolation or grounding are not compromised.
- Avoid mechanical shock or overstressing the *miniARM® LOCKING COUPLER*.
- To avoid damage to the camera's distal tip, ensure that the camera is not flexed before insertion or removal.
- Only the cleaning and sterilization processes which are defined within these instructions for use have been validated.
- Store at room temperature.

ADVERSE REACTIONS

No documented adverse reactions.

TABLE 2: miniARM® LOCKING COUPLER PRODUCT ORDERING

	REORDER	PRODUCT	DESCRIPTION
	REF 080555	<i>miniARM® LOCKING COUPLER*</i>	1 Shelf Box
	REF 100034	<i>multiSTATION™ System Sterilization Tray†</i>	1 Shelf Box

* For Use with the Olympus® ENDOEYE FLEX® 5mm HD surgical video camera (LTF-S190-5) (NOTE: the *miniARM® Locking Coupler* requires either the *miniARM® Cannula Adapter Kit* (#080415) OR the *miniARM® Clamping Adapter* (#080519) in order to be attached to the *miniARM® Instrument Holder*)

† The *multiSTATION™ System Sterilization Tray* is manufactured by: Summit Medical, 815 Vikings Parkway, Suite 100 St. Paul, MN 55121 U.S.A.

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The LSI logo, LSI Solutions, Multistation, Miniarm, and Perfect Performance Policy are registered trademarks and trademarks of LSI SOLUTIONS®.

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LSI SOLUTIONS®

LSI SOLUTIONS, Inc.
7796 Victor-Mendon Road
Victor, New York 14564 U.S.A.
Phone: +1 585.869.6600
Customer Service: +1 866.575.3493
Technical Support: +1 866.428.9092
Fax: +1 585.742.8086
www.lsisolutions.com

Patents: www.lsisolutions.com/patents
The LSI logo, LSI Solutions, Multistation, Miniarm, and Perfect Performance Policy are registered trademarks and trademarks of LSI SOLUTIONS®.
Copyright © 2021, LSI SOLUTIONS®. All Rights Reserved.

Symbol Glossary:
www.lsisolutions.com/symbols

OLYMPUS and the ENDOEYE FLEX are registered trademarks of the Olympus Corporation. Prolystica is a registered trademark of Steris. Weiman is a registered trademark of PurposeBuilt Brands. Use of these names, trademarks and brands does not imply endorsement. All company, product and service names used in this technology guide are for identification purposes only.

The *multiSTATION™* System: *SubX™ STERNAL ASCENDER™ SYSTEM* and *LIMA/RIMA RetroSterno™ Paddles* Technology Guide

 READ THIS PRODUCT INSERT THOROUGHLY BEFORE USE


FIG. 1 – SubX™ STERNAL ASCENDER™ SYSTEM 

SubX™ STERNAL ASCENDER™ SYSTEM DEVICE DESCRIPTION:

The non-sterile, reusable *SubX™ STERNAL ASCENDER™ SYSTEM* (FIG. 1) is a mechanical hoist mechanism that provides elevation of a selected accessory (sold separately) when configured with compatible *M55™ RIGID miniARM® HOLDERS* and *miniARM® RAIL CLAMPS*.

The selected accessory, *LIMA RetroSterno™ Paddle* 1 shown as an example above, is loaded into the accessory slot 2 and secured by pulling down. The hand grip 3 provides an ergonomic handle to insert and manipulate the accessory extracorporeally and create initial retraction. The paddle indicator 4 indicates the location of the tip of the *RetroSterno™ Paddle* during placement. The pivot button 5 is pressed to unlock the system and allow rotation of the rack 6 about the rack pivot 7. Once the rack is oriented as desired, the pivot button is pressed again to lock the angular position of the rack. The arch keystone 8 includes 2 attachment points into which 2 *M55™ Rigid miniARM® HOLDERS* (not shown) are inserted to establish the bridge configuration. Two locking levers 9 rotate to secure and release the *miniARM® HOLDER* quick connect posts to the arch keystone or remove them from the arch keystone. The swivel bar 10 is rotated counterclockwise to raise the system. Every half turn, the ratcheting latch 11 engages an anti-backdrive feature in the arch keystone. An audible click indicates that the locking feature has been activated and secure elevation has been achieved.

Upon completion of the procedure, the ratcheting latch is pressed and held to bypass the anti-backdrive mechanism. The swivel bar is rotated clockwise to release the load and lower the system.


INDICATIONS FOR USE:

The *miniARM® SubX™ STERNAL ASCENDER™ SYSTEM* is indicated for use by surgeons to hold instruments in a fixed position for a period of time.


INSTRUCTIONS FOR USE

ATTACHING, PIVOTING, RAISING, and LOWERING the SubX™ STERNAL ASCENDER™ SYSTEM


Prior to each use, test the SubX™ STERNAL ASCENDER™ SYSTEM by confirming that the arch keystone can travel the entire length of the rack. Discontinue use if the arch keystone binds at any point along the rack. Steps 1-4 show a LIMA RetroSterno™ Paddle as an example of a compatible accessory. Accessories available separately.

1. **ATTACH** the ascender to 2 holders by inserting the quick connect posts into the arch keystone and rotating the locking levers to the locked position.


2. **PIVOT** the hand grip about the rack while the pivot button is disengaged. Engage the pivot button to lock the angle of the hand grip.

3. **RAISE** the hoist mechanism by rotating the swivel bar counterclockwise.

NOTE: Ensure that the swivel bar is rotated until an audible click is heard, indicating that the arch keystone is locked in position.

4. **LOWER** the hoist mechanism by rotating the swivel bar clockwise while pressing the ratcheting latch.

NOTE: The swivel bar cannot rotate clockwise without the ratcheting latch disengaged.

The *multiSTATION™* System: *SubX™ STERNAL ASCENDER™ SYSTEM* and *LIMA/RIMA RetroSterno™ Paddles* Technology Guide

READ THIS PRODUCT INSERT THOROUGHLY BEFORE USE

1. LIMA RetroSterno™ Paddle
2. RIMA RetroSterno™ Paddle
3. Textured Zone
4. Paddle Post
5. Attachment Stud


FIG. 2 – LIMA/RIMA RetroSterno™ Paddles

SOLUTIONS®
INNOVATIVE SOLUTIONS FOR MINIMALLY INVASIVE SURGERY


LIMA/RIMA RetroSterno™ Paddle DEVICE DESCRIPTION:

The non-sterile, reusable *RIMA RetroSterno™ Paddle* and *LIMA RetroSterno™ Paddle* (FIG. 2) are accessories to the *SubX™ STERNAL ASCENDER™ SYSTEM*. The *LIMA RetroSterno™ paddle 1* provides retraction of the sternum and sub-xiphoid access to the left internal mammary artery. The *RIMA RetroSterno™ paddle 2* provides retraction of the sternum and sub-xiphoid access to the right internal mammary artery. The textured zone **3** provides grip on the sub-sternal tissue during retraction. The *RetroSterno™ Paddles* are offset laterally to maximize working space and concavely contoured to facilitate optimal visualization of the targeted vessel. The shape of the paddle post **4** interfaces with space lateral to the xiphoid process and opposite the targeted internal mammary artery to create a low-profile anatomical insertion. The attachment studs **5** allow for quick connection to the accessory slot of the *SubX™ STERNAL ASCENDER™ SYSTEM*.

INDICATIONS FOR USE:


The *miniARM® LIMA RetroSterno™ Paddle*/*miniARM® RIMA RetroSterno™ Paddle* is indicated for use for retraction of bone and soft tissue.

1 SELECT

1. **SELECT** the LIMA paddle for access to the left internal mammary artery or the RIMA paddle for access to the right internal mammary artery.


2 ATTACH

2. **ATTACH** a paddle by inserting the paddle post into the accessory slot of the hand grip and pulling down to secure.

INSTRUCTIONS FOR USE

SETTING UP, ATTACHING, and POSITIONING the multiSTATION™ Bridge

multiSTATION™ Bridge DESCRIPTION:

The *multiSTATION™ Bridge* (FIG. 3) comprises 2 M55™ RIGID miniARM® HOLDERS, 2 miniARM® RAIL CLAMPS, a SubX™ STERNAL ASCENDER™ SYSTEM, and a compatible accessory such as a LIMA RetroSterno™ Paddle.


The bridge is assembled by attaching the rail clamps to draped surgical table rails, connecting a holder to the rail clamp, and then connecting the ascender to the holders. This configuration provides a stable, rigid system with which compatible accessories may be securely raised and lowered.

The holders may be connected to either a side or center receiving DOC of the rail clamp to provide positional adjustability. The ascender, with a compatible accessory attached, is connected to the 2 holders and adjusted into position. Once the desired position of the bridge is reached, the system is locked into place. The swivel bar is used to raise the compatible accessory, providing a stable opening to the surgical site. The accessory paddle can be lowered by rotating the swivel bar clockwise while pressing the ratcheting latch. The bridge can then be repositioned, or a different accessory paddle can be installed.

INSTRUCTIONS FOR USE

SETTING UP, ATTACHING, and POSITIONING the multiSTATION™ Bridge


1 SET UP

1. **SET UP** 2 rail clamps and 2 holders on opposite surgical table rails. See each device's instructions for use for more information.

NOTE: Fully tighten the rail clamps to the surgical table rails. Do not fully lock the holders.


2 ATTACH

2. **ATTACH** a paddle to the ascender and connect the arch keystone to the holders.

NOTE: Ensure that the arch keystone is at the top of the rack to allow for maximum retraction.

3 POSITION


3. **POSITION** the paddle by adjusting the ascender and/or the holders. The paddle post should interface with the space lateral to the xiphoid process, opposite the targeted internal mammary artery.

NOTE: The paddle indicator on the top of the hand grip represents the location of the tip of the paddle. The paddle should always be kept parallel to the sternum.

INSTRUCTIONS FOR USE

LOCKING, LIFTING, AND RELEASING the SubX™ STERNAL ASCENDER™ SYSTEM


1 LOCK

1. **LOCK** the bridge in place by rotating the holder adjustment knobs clockwise until fully tightened and engaging the pivot button.

NOTE: Verify that the bridge's position is secure before proceeding.


2 LIFT

2. **LIFT** the sternum by rotating the swivel bar counterclockwise until sufficient visualization of the surgical site has been achieved.

NOTE: An audible click indicates that the locking feature is engaged.


3 RELEASE


3. **RELEASE** the system by rotating the swivel bar clockwise while pressing the ratcheting latch. With the system fully lowered, the paddle can be removed or swapped.

NOTE: Do not attempt to remove a paddle or reposition the bridge while the sternum is still under load.


REMOVAL and CLEANUP

1. **ROTATE** the locking levers of the ascender to the unlocked position.

2. **DETACH** the ascender by sliding it off both of the quick connect posts.

3. **REMOVE** the paddle by lifting the post up and out of the accessory slot of the hand grip.

4. **REMOVE** the holders and rail clamps according to each device's instructions for use.

NOTE: Fully support the holder while removing it from the rail clamps. Removal of the holder may be eased by keeping the holder in its locked position.

SubX™ STERNAL ASCENDER™ SYSTEM REPROCESSING




- Disassemble and clean the device immediately after use. Do not allow a soiled device to dry.
- The SubX™ STERNAL ASCENDER™ SYSTEM is not validated to be cleaned or sterilized with an accessory attached.
- Cleaning agent used in validation: Steris Prolystica 2X (enzymatic, neutral pH).
- Perform the final rinse using only freshly prepared purified water/highly purified water.
- Never use metal brushes or steel wool for cleaning.
- Prepare and reprocess any holder or rail clamp according to each device's instructions for use.

POINT OF USE

1. Disassemble and clean device immediately after use.
2. Do not allow soiled devices to dry.





PREPARATION

1. Remove the ascender from the 2 holders by rotating the locking levers and sliding the device off the quick connect posts (FIG. 4).
2. Push the pivot button to the unlocked position to allow the hand grip to freely rotate (FIG. 5).
3. Rotate the swivel bar to position the arch keystone approximately in the middle of travel along the rack. Press and hold the ratcheting latch to rotate the swivel bar clockwise. (FIG. 6).
4. If using the *multiSTATION™ System Sterilization Tray* (Part Number 100034), clean the tray separately according to the sterilization tray instructions for use.


MANUAL PRE-CLEANING

1. Soak in enzymatic, neutral pH cleaning solution for a minimum of 5 minutes. Refer to detergent manufacturer's instructions for temperature and concentration.
2. Use a plastic-bristled brush to thoroughly scrub the device's exterior surfaces, paying special attention to the following:
 - a. Accessory slot (FIG. 7)
 - b. Rack pivot and pivot button (FIG. 8)
 - c. Arch keystone, including swivel bar, ratcheting latch, and attachment opening (FIG. 9)
 - d. Front and back of the rack, rack teeth (FIG. 10)

3. Rinse with warm water for a minimum of 2 minutes.
4. Proceed with one of the two required cleaning options: ultrasonic or automated.



FIG. 11

STERILIZATION

1. The device must be properly cleaned prior to sterilization.
2. Perform sterilization cycle according to TABLE 1:


TABLE 1: SubX™ STERNAL ASCENDER™ SYSTEM, LIMA/RIMA RetroSterno™ Paddle Sterile Packaging and Processing

Method	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Immediate use steam sterilization according to ANSI/AAMI ST79
Container	<i>multiSTATION™ System Sterilization Tray P/N 100034</i>	No tray	No tray
Cycle	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)
Packaging	2-layer polypropylene wrap	Pouch	No packaging
Temperature	132°C (270°F) (minimum)	132°C (270°F) (minimum)	132°C (270°F) (minimum)
Exposure Time	4 minutes (minimum)	4 minutes (minimum)	4 minutes (minimum)
Dry Time	65 minutes (minimum)	25 minutes (minimum)	N/A

Device(s) processed by immediate use sterilization should be transferred immediately, using aseptic technique, from the sterilizer to the point of use.

Refer to ANSI/AAMI ST79, Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities.

The SubX™ STERNAL ASCENDER™ SYSTEM, LIMA RetroSterno™ Paddle, and RIMA RetroSterno™ Paddle have been validated to 20 reprocessing cycles. The useful lifespan of a surgical instrument is largely dependent on the care and handling of the instrument. Careful inspection and functional testing of the instrument should be used to determine the end of its serviceable life.

STORAGE

1. During storage, ensure the device remains in a sterile condition ready for reuse.
2. Shelf life is dependent on the sterile barrier employed, storage manner, and environmental and handling conditions.

CONTRAINDICATIONS

- Do not use with attachments other than accessories provided by LSI SOLUTIONS®.
- These devices are not intended for use except as indicated.

WARNINGS

- Federal law restricts this device to sale, distribution, and use by, or on, the order of a physician.
- Read and become familiar with all instructions, warnings, and cautions before using this product.
- The *SubX™ STERNAL ASCENDER™ SYSTEM*, *LIMA RetroSterno™ Paddle*, and *RIMA RetroSterno™ Paddle* shall be used in accordance with these instructions for use.
- Improper use of the *SubX™ STERNAL ASCENDER™ SYSTEM*, *LIMA RetroSterno™ Paddle*, and *RIMA RetroSterno™ Paddle* may cause serious injury or death. In addition, improper care and maintenance of the device may render the device non-sterile prior to patient use and may cause serious injury to the health care provider or the patient.
- When using the *SubX™ STERNAL ASCENDER™ SYSTEM*, *LIMA RetroSterno™ Paddle*, and *RIMA RetroSterno™ Paddle*, patients must be immobilized or anesthetized.
- Discontinue use of the *SubX™ STERNAL ASCENDER™ SYSTEM*, *LIMA RetroSterno™ Paddle*, and *RIMA RetroSterno™ Paddle* when patient is moving or being moved.
- Do not use this retraction system without adequate knowledge or experience regarding surgical retraction, including sternal and rib retraction. Avoid compressing, tensioning, displacing, or compromising tissue structures, such as the heart, that may be injured or lead to interoperative or postoperative loss of function.
- Surgical or endoscopic procedures should be performed only by physicians having adequate training and familiarity with relevant techniques and anatomy. Medical literature relating to techniques, complications, and hazards should be consulted prior to use.

PRECAUTIONS


- **The *SubX™ STERNAL ASCENDER™ SYSTEM*, *LIMA RetroSterno™ Paddle*, and *RIMA RetroSterno™ Paddle* are packaged as non-sterile.** Clean and sterilize prior to use.
- Proper cleaning, sterilization, packaging, storage and deployment are required to ensure sterility and safe use.
- If there are any variations between these instructions for use and the policies of your facility and/or the instructions of your cleaning/sterilizing equipment manufacturer, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your device.
- Use of the *SubX™ STERNAL ASCENDER™ SYSTEM*, *LIMA RetroSterno™ Paddle*, and *RIMA RetroSterno™ Paddle* for a task other than what it is intended for can result in a damaged or broken device, and/or injury or death.
- Prior to use, inspect the *SubX™ STERNAL ASCENDER™ SYSTEM*, *LIMA RetroSterno™ Paddle*, and *RIMA RetroSterno™ Paddle* to ensure proper function and condition. Do not use devices if they do not satisfactorily perform their intended function or if they have physical damage.
- Ensure operating room table rails are of adequate strength to support use of this system to lift tissue structures as indicated for the individual patient.
- Ensure system connections, fittings and engagement features are secured appropriately, fully and accurately positioned to confer the required function.
- Ensure placement of system components to not compromise patient access during surgery or injure patient or tissue structures.
- Ensure retractors are placed in contact with tissue structures under adequate visualization. Do not force system components against non-targeted tissue.
- Ensure all system components are appropriately removed from the patient and operating room table prior to patient movement.
- Surgical instruments vary between manufacturers. Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure electrical isolation or grounding are not compromised.
- Like all metal implements used in surgery, the components of this system conduct electricity and can be associated with electrical shock or arcing from cautery or other sources of current.
- Avoid mechanical shock or overstressing the *SubX™ STERNAL ASCENDER™ SYSTEM*, *LIMA RetroSterno™ Paddle*, and *RIMA RetroSterno™ Paddle*.
- Only the cleaning and sterilization processes that are defined within these instructions for use have been validated.
- Store at room temperature.

ADVERSE REACTIONS

- No documented adverse reactions.

ORDERING INFORMATION

TABLE 2: *SubX™ STERNAL ASCENDER™ SYSTEM* PRODUCT ORDERING

	REORDER	PRODUCT	DESCRIPTION
	REF 080900	<i>SubX™ STERNAL ASCENDER™ SYSTEM</i>	1 Shelf Box
	REF 081214	<i>LIMA RetroSterno™ Paddle</i>	1 Shelf Box
	REF 081241	<i>RIMA RetroSterno™ Paddle</i>	1 Shelf Box
	REF 081202	<i>RetroSterno™ Paddle Kit</i>	1 Shelf Box
	REF 100034	<i>multiSTATION™ System Sterilization Tray*</i>	1 Shelf Box

LSI SOLUTIONS®

LSI SOLUTIONS, Inc.
7796 Victor-Mendon Road
Victor, New York 14564 U.S.A.
Phone: +1 585.869.6600
Customer Service: +1 866.575.3493
Technical Support: +1 866.428.9092
Fax: +1 585.742.8086
www.lsisolutions.com

*The *multiSTATION™* System Sterilization Tray is manufactured by Summit Medical, 815 Vikings Parkway, Suite 100, St. Paul, MN 55121 U.S.A.

MADE IN THE USA

This Product Comes
with our LSI SOLUTIONS®
Perfect Performance Policy*
Call us at 866.575.3493 any time.

Patents: www.lsisolutions.com/patents

The LSI logo, LSI Solutions, Perfect Performance Policy, M45, M55, Click Release, Miniarm, and Multistation are registered trademarks and trademarks of LSI SOLUTIONS®.
Copyright © 2021, LSI SOLUTIONS®. All Rights Reserved.

Symbol Glossary:
www.lsisolutions.com/symbols

EC REP
Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

Prolystica is a registered trademark of Steris. Weiman is a registered trademark of PurposeBuilt Brands. Use of these names, trademarks and brands does not imply endorsement. All company, product and service names used in this technology guide are for identification purposes only.

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MADE IN THE USA

This Product Comes
with our LSI SOLUTIONS®
Perfect Performance Policy®
Call us at 866.575.3493 any time.

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
VALVEVAULT™ DEVICE TECHNOLOGY GUIDE

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VALVEVAULT™ DEVICE DESCRIPTION

SOLUTIONS®

The VALVEVAULT™ DEVICE (FIG. 1) is a reusable, sterilizable delivery device for transporting a prosthesis through a surgical incision. Each holder (1) has two differently sized chambers (2) designed to hold a prosthesis. A single letter "S", "M" or "L" for device size (3), found on both the holder and cover (5), and denotes Small, Medium and Large sized VALVEVAULT™ DEVICES, respectively. A chamber diameter (4) is marked within a heart shaped symbol on both the cover and bottom of the VALVEVAULT™ DEVICE Holder to denote the approximate diameter for each chamber. Alignment arrows (6) point towards the larger chamber and are used to orient and pair a corresponding holder and cover together while establishing enclosure of the prosthesis and corresponding suture. A shaft (7) contains a channel (8) connecting the two chambers of the holder. During use, annular facing suture from the surgical site is guided through the suture opening (10) of the selected chamber allowing the prosthesis sewing cuff to rest on or just below the inside shoulder (9). Free ends of the suture from the prosthesis are routed within the channel and over the vacant chamber to exit through the opposite suture opening. The cover alignment arrows are then matched to the holder alignment arrows to completely encase and protect the prosthesis and suture for delivery through the surgical incision.

INDICATIONS FOR USE:

The VALVEVAULT™ DEVICE is indicated for use to hold and deliver a prosthesis through a surgical incision.

CONTRAINDICATIONS

- The VALVEVAULT™ DEVICE is contraindicated for use with any balloon or self-expandable prosthesis.
- The VALVEVAULT™ DEVICE is contraindicated for prosthesis extraction.

CHOOSING THE SIZE OF THE VALVEVAULT™ DEVICE


NOTE: In addition to the size listing on the package, the photo below shows the letter “S” in red circles denoting the small size VALVEVAULT™ DEVICE. Medium and Large VALVEVAULT™ DEVICES are marked with “M” and “L” respectively; see below.

NOTE: For each chamber of the VALVEVAULT™ DEVICE, the approximate chamber diameter (measured in mm) is printed within a heart shape on the bottom of the holder and top of the cover. This numerical designation for chamber diameter is **not intended** to represent the prosthesis size.

TABLE 1. VALVEVAULT™ DEVICE SIZE GUIDE			
VALVEVAULT™ DEVICE Size	Reorder Number	Chamber Diameter (mm)	Chamber Depth (mm)
Small (S)	081025	19	11.5
		21	12.5
Medium (M)	081030	23	13.5
		25	14.5
Large (L)	081035	27	15.5
		29	16.2

NOTE: Surgeons should choose their preferred device size and chamber size based on surgical approach, prosthesis knowledge, and patient anatomy.

1. **SELECT** the appropriately sterilized VALVEVAULT™ DEVICE, (i.e., small, medium or large) based on the chamber most likely to readily fit the outer profile (i.e., height as well as width) of the sutured prosthesis.

2. **SIZE** the chamber if desired to estimate an appropriate fit for the prosthesis. Compare the selected valve sizer to the selected prosthesis. The valve sizer should fit easily into the selected chamber.

LOADING THE VALVEVAULT™ DEVICE

For the following loading steps, note that suture sections are named as shown:


FREE ENDS OF THE SUTURE (ABOVE THE PROSTHESIS)

SURGICAL SITE


NOTE: Avoid potential prosthesis damage by ensuring suture does not become entangled with prosthesis leaflets. Ensure all surgical needles are removed prior to VALVEVAULT™ DEVICE use to affirm that only the suture strands from the annulus and the free ends are associated with the prosthesis.

1 ORIENT / PLACE


1. **ORIENT** the selected VALVEVAULT™ DEVICE Chamber to be closest to the patient. With suture from the annulus through the prosthetic sewing cuff, gently **PLACE** the prosthesis into selected chamber with leaflets facing down, sewing cuff up. The prosthetic sewing cuff should rest on or under the inside shoulder. Ensure the suture from the annulus is gathered completely within the suture opening. Never force the prosthesis into the chamber or the cover onto the prosthesis.

2 INSPECT

2. **INSPECT** the prosthesis and suture to ensure leaflets are not compromised and suture lays correctly.


3 GUIDE

3. **GUIDE** the free ends of the suture from the prosthetic sewing cuff through the channel and unused chamber and out of the opposite suture opening.

4 CLOSE

4. **CLOSE** the VALVEVAULT™ DEVICE Holder by matching the alignment arrows on the cover to those on the holder. Ensure all suture is contained within the suture openings and the cover is flush with the holder's edges without compressing the prosthesis. **Do not** force the cover onto the prosthesis if it does not easily lay flush with the holder. If needed, select the next size larger chamber and try again.

DELIVERING THE PROSTHESIS WITH THE VALVEVAULT™ DEVICE

1. **DELIVER** the VALVEVAULT™ DEVICE through the surgical incision, preferably with the device cover facing the lower rib.

2. **OPEN** the VALVEVAULT™ DEVICE by removing the cover from the holder. **Do not** expose the prosthesis or compromise it by contact with other structures. Remove the cover from the surgical incision and away from the surgical site.


3. **TRANSFER** the suture and prosthesis from the VALVEVAULT™ DEVICE Holder to the surgical site by gently turning the shaft within the incision. If the prosthesis does not freely exit the chamber, the sutures can be used to gently pull the prosthesis out of the chamber. Care must be taken to avoid applying excess force to the prosthesis or sutures connected to the annulus.


4. **REMOVE** the VALVEVAULT™ DEVICE Holder from the surgical incision taking care not to contact/damage the prosthesis, especially avoid contact with the prosthetic leaflets, while ensuring all suture is away from both suture openings.


VALVEVAULT™ DEVICE REPROCESSING

- Disassemble and clean the device immediately after use. Do not allow a soiled device to dry.
- The VALVEVAULT™ DEVICE is not validated to be cleaned or sterilized in the fully assembled state. It should be disassembled into the components as shown in FIGS. 2 and 3.
- Cleaning agent used in validation: Steris Prolystica 2X (Enzymatic, Neutral pH)
- Perform the final rinse using only freshly prepared purified water/highly purified water.
- Never use metal brushes or steel wool for cleaning.

PACKAGING (CONTINUED)

NOTE: The tray lid only fits the tray bottom in one direction. If the side locks cannot both be securely latched, it may be necessary to reorient the lid to achieve a secure fit.

For more tray information see the sterilization tray instructions for use.

3. Package the VALVEVAULT™ DEVICE according to TABLE 2. The barrier system for sterilized re-usable instruments should meet the following requirements:

- ISO 11607-1
- Suitable for Pre-Vacuum steam sterilization
- Appropriate for medical use
- Grade appropriate for weight of loaded tray per sterilization tray IFU and facility procedures

STERILIZATION

1. The device must be properly cleaned before sterilization.
2. Perform sterilization cycle according to TABLE 2:

TABLE 2: VALVEVAULT™ DEVICE Packaging and Processing			
Method	Moist heat (steam) sterilization according to ANSI/AAMI ST 79	Moist heat (steam) sterilization according to ANSI/AAMI ST 79	Immediate Use Steam Sterilization according to ANSI/AAMI ST 79
Container	VALVEVAULT™ DEVICE Sterilization Tray Part Number 100035	No Tray	No Tray
Cycle	Pre-Vacuum (Pre-Vac)	Pre-Vacuum (Pre-Vac)	Pre-Vacuum (Pre-Vac)
Packaging	2-layer polypropylene wrap	Pouch	No Packaging
Temperature	132°C (270°F) (minimum)	132°C (270°F) (minimum)	132°C (270°F) (minimum)
Exposure Time	4 minutes (minimum)	4 minutes (minimum)	4 minutes (minimum)
Dry Time	30 minutes (minimum)	20 minutes (minimum)	N/A

Device(s) processed by Immediate Use Sterilization should be transferred immediately, using aseptic technique, from the sterilizer to the point of use.

Refer to ANSI/AAMI ST79 Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities.

The VALVEVAULT™ DEVICE has been validated to one hundred (100) reprocessing cycles. The useful lifespan of a surgical instrument is largely dependent on the care and handling of the instrument. Careful inspection and functional testing of the instrument should be used to determine the end of its serviceable life.

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WARNINGS

- Federal law restricts this device to sale, distribution and use by, or on, the order of a physician.
- Read and become familiar with all instructions, warnings and cautions before using this product.
- This device shall be used in accordance with these Instructions for Use.
- The chamber diameter designations of the VALVEVAULT™ DEVICE are not intended to represent valve size.
- Improper use of this device may cause serious injury. Improper care and maintenance of the device may render the device non-sterile prior to patient use and may cause serious injury to the health care provider or the patient.
- All actions beyond passing the chamber through the surgical access site (e.g., an intercostal space) should be smooth and easily performed. Avoid applying any excess compression or tensioning forces to the prosthesis and/or the suture. Avoid unnecessary stressors on VALVEVAULT™ components, the prosthesis, suture etc.
- Do not use a VALVEVAULT™ chamber that does not readily accommodate the prosthesis. If there is any doubt regarding the appropriate fit of the prosthesis in the selected chamber, use the next size up chamber.
- If passage of the VALVEVAULT™ through the surgical access site is not easily achieved, consider increasing the size of the access site or appropriate alternative routes.
- When using most surgical devices, harsh contact with tissue or forcing the device during its handling and use can damage or injure the prosthesis, suture or tissue structures like bone, blood vessels, nerves, soft tissue (e.g., cardiac valve tissue), etc. Inappropriate device handling and use can lead to bleeding, tissue structural damage, nerve impairment, infection, disability and even death. Special care should be taken to avoid damaging delicate tissue structures.
- Discontinue use of the VALVEVAULT™ DEVICE when moving patient or when patient is moving.
- Endoscopic procedures should be performed only by physicians having adequate training and familiarity with endoscopic techniques and relevant anatomy. Medical literature relative to techniques, complications and hazards should be consulted prior to use.
- When using the VALVEVAULT™ DEVICE, an insufficient surgical access site incision and dissection may cause increased penetration force which may reduce the surgeon's control during entry.
- Direct contact between sensitive tissue structures and foreign materials can lead to tissue injury or damage.
- When using the VALVEVAULT™ DEVICE, keep the prosthesis irrigated per Manufacturer's instructions.

PRECAUTIONS

- Devices are packaged as non-sterile.** Cleaning and sterilizing of this device must occur prior to use.
- If there are variations between these Instructions for Use and either your facility's policies and/or your cleaning-sterilizing equipment manufacturer's instructions, those variations should be brought to the attention of appropriate responsible hospital personnel for resolution before proceeding with device cleaning and sterilizing.
- Use of a device for a task other than what it is intended for will usually result in a damaged or broken device.
- Prior to use, inspect the cover and holder to ensure proper fit and condition. Do not use the cover or holder if they do not fit together or if they have physical damage.
- Surgical instruments vary between manufacturers. Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure electrical isolation or grounding are not compromised.
- Avoid mechanical shock or overstressing the devices.
- Only the cleaning and sterilization processes defined within these Instructions for Use have been validated.
- Care should be taken to ensure that all sutures are manually seated properly within the channel and the suture openings to avoid damage to the prosthesis or damage to the sutures.
- Store at room temperature.

ADVERSE REACTIONS

No documented adverse reactions.

ORDERING INFORMATION

TABLE 3. VALVEVAULT™ DEVICE		SUPPLIED: NONSTERILE, REUSABLE	
	REORDER	PRODUCT	DESCRIPTION
 x 1 S	REF 081025	VALVEVAULT™ DEVICE (S) 19mm / 21mm	Box of 1 Device
 x 1 M	REF 081030	VALVEVAULT™ DEVICE (M) 23mm / 25mm	Box of 1 Device
 x 1 L	REF 081035	VALVEVAULT™ DEVICE (L) 27mm / 29mm	Box of 1 Device
 x 1	REF 081040	VALVEVAULT™ DEVICE KIT	Combo Box of 3 Devices (1 of each size)
 x 1	REF 100035	VALVEVAULT™ DEVICE STERILIZATION TRAY	Box of 1 Device

LSI SOLUTIONS®
 LEADING SURGICAL INNOVATION
 LSI SOLUTIONS, Inc
 7796 Victor-Mendon Road
 Victor, New York 14564 U.S.A.
 Phone: +1 585.869.6600
 Customer Service: +1 866.575.3493
 Technical Support: +1 866.428.9092
 Fax: +1 585.742.8086
 www.lsi-solutions.com
 Регистрация МИ в Росздравнадзоре
 www.nevacert.ru | info@nevacert.ru

The VALVEVAULT™ Sterilization Tray is manufactured by Summit Medical, 815 Vikings Parkway, Suite 100, St. Paul, MN 55121 U.S.A.

MADE IN THE USA

This Product Comes with our LSI SOLUTIONS® Perfect Performance Policy® Call us at 866.575.3493 any time.

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The LSI logo, LSI SOLUTIONS, Perfect Performance Policy, and VALVEVAULT are registered trademarks and trademarks of LSI Solutions, Inc. Copyright © 2020, LSI SOLUTIONS®. All Rights Reserved.

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LSI SOLUTIONS®
LEADING SURGICAL INNOVATION

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
FIG. 1A 3D[®] mini ACCESS PORTAL

DESCRIPTION

- 1 – SOFT CRADLES
- 2 – EXPANSION ARMS
- 3 – HEIGHT ADJUSTMENT SCREWS
- 4 – EXPANSION INDICATOR
- 5 – ACCESSORY POST
- 6 – GEAR RECEIVER
- 7 – RELEASE LEVER

The 3D[®] mini ACCESS PORTAL is a sterile SINGLE PATIENT USE device that enables the controlled displacement of adjacent ribs and soft tissue in both transverse and anterior-posterior directions. To facilitate less traumatic surgical access to intrathoracic structures, 3D[®] mini ACCESS PORTAL soft cradles **1** are configured to engage and spread adjacent ribs and overlaying cutaneous tissue through a small skin incision. A removable control knob **8** sets the desired expansion width and height of the expansion arms **2**. When the control knob expansion drive gear **9** is coupled to the gear receiver **6**, rotation of the control knob handle spreads the expansion arms transversely apart; release lever **7** enables closing of expansion arms. When the control knob hexagonal key **10** is turned within either of the height adjustment screws **3**, the selected expansion arm moves to the desired height. Removal of the control knob optimizes visualization and accessibility of the surgical field. The 3D[®] mini ACCESS PORTAL also has an accessory post **5** that serves as an attachment point for optional accessories.

INDICATIONS

The 3D® mini ACCESS PORTAL is indicated for the retraction of bone and soft tissue.

3D® mini ACCESS PORTAL ADJUSTMENTS

EXPANSION

Manually


FIG. 2

SPREAD

the expansion arms manually by pulling the expansion arms apart. An internal latch maintains the expansion width.

OR

Using Control Knob

FIG. 3A

INSERT

control knob expansion drive gears fully into gear receiver.

FIG. 3B

ROTATE

the control knob handle counter-clockwise ("OPEN" Direction) to achieve desired arm expansion. Note: Expansion has reached its maximum limit when the expansion indicator marks line up.

EXPANSION RELEASE

PRESS

the release lever (in the direction shown) to defeat the locking mechanism, enabling the expansion arms to collapse.


FIG. 4A

HEIGHT ADJUSTMENT

FIG. 5A

INSERT

hexagonal key of the control knob into the desired height adjustment screw.


FIG. 5B

TURN

the control knob handle clockwise ("UP" direction) in order to raise the expansion arm in the upward direction.


FIG. 5C

TURN

the control knob counter-clockwise in order to lower the expansion arm in the downward direction.

REPEAT

adjusting the height adjustment screws as needed.

IN USE

1. Create desired skin incision above the selected intercostal space to accommodate soft cradles. Dissect underlying muscles and pleural structures per routine (FIG. 6A).


FIG. 6A


FIG. 6B

May require adjustment as soft and bony tissue settles

2. Position soft cradles through incision and around the ribs and/or soft tissue (FIG. 6B).
3. Adjust expansion arm width with control knob as desired (see FIGS. 2, 3A-3B).
4. Adjust expansion arm height with control knob as desired (see FIGS. 5A- 5C).
5. Remove control knob and attach accessories, such as the RAM® RING, to the accessory post using the hexagonal key (FIG. 6C).


FIG. 6C

6. When thoracic access is no longer needed, undo height adjustment with hexagonal key (see FIGS. 5A- 5C).
7. Collapse expansion arms with release lever (see FIGS. 4A- 4B).
8. Remove 3D® mini ACCESS PORTAL from the incision.

FIG. 7 RAM ACCESS PRODUCT ORDERING		
REORDER	PRODUCT	DESCRIPTION
110300	3D® mini ACCESS PORTAL	Box of 6 (1 Device per Box)
080215	RAM® RING Suture Management System	Box of 8 (1 Device per Box)
080228	RAM® RACK Suture Management System	12 Shelf Boxes (3 Devices per Shelf Box) (36 Devices per Ship Box)

ACTIONS

The 3D® mini ACCESS PORTAL facilitates surgical access to the organs and tissues within the thoracic cavity. With the device placed between ribs or tissue, the expansion arms and soft cradles can displace the ribs in the transverse and inward or outward directions. Accessories, such as suture management devices like the RAM® RING, can be attached to the 3D® mini ACCESS PORTAL device as desired.

CONTRAINDICATIONS

- Do not use this technology under conditions in which excessive tissue compression or tension can lead to tissue damage.
- Do not use with attachments or actuators not approved by LSI SOLUTIONS®.

WARNINGS

- Federal (U.S.A.) law restricts this device to sale, distribution and use by, or on, the order of a physician.
- Minimally invasive surgical procedures should only be performed by physicians having adequate training and familiarity with surgical and endoscopic techniques. In addition, medical literature should be consulted relative to techniques, complications, and hazards prior to the performance of minimally invasive surgical procedures.
- Read and become familiar with all instructions, warnings, and precautions before using this product.
- **Do not resterilize.** The 3D® mini ACCESS PORTAL is designed and intended for single patient use only. Do not reuse, reprocess, or resterilize this product. The performance of the 3D® mini ACCESS PORTAL after cleaning or other reprocessing has not been verified and is not supported by LSI SOLUTIONS®. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Do not use if package or product is damaged. Discard open (unsealed), unused, or damaged product.
- It is the responsibility of the user to dispose of the device in accordance with local regulations and hospital procedures.
- As with all retractors and surgical tissue manipulation devices, care should be taken to use only as much retraction as necessary to provide adequate access and visualization; rib fracture, vascular injury (bleeding), nerve damage, and/or other patient trauma or death may result.
- Excess pressure on the 3D® mini ACCESS PORTAL may lead to 3D® mini ACCESS PORTAL damage and/or patient injury.
- The 3D® mini ACCESS PORTAL is intended to provide surgical access for minimally invasive cardiothoracic procedures by retraction of soft and bony tissue. Applications other than these may damage the device making it unsuitable for continued use.

PRECAUTIONS

- Do not use devices if they do not satisfactorily perform their intended function or if they have physical damage. Avoid mechanical shock or overstressing the devices.
- Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure that function is not compromised.

ADVERSE REACTIONS

No documented adverse reactions.


LSI SOLUTIONS®
INNOVATIVE SOLUTIONS FOR MINIMALLY INVASIVE SURGERY

LSI Solutions, Inc.
7796 Victor-Mendon Road
Victor, New York 14564 U.S.A.
Phone: 585.869.6600
Customer Service: 866.575.3493
Technical Support: 866.428.9092
Fax: 585.742.8086
www.lsisolutions.com

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The LSI logo, LSI SOLUTIONS, Perfect Performance Policy, 3D and RAM are registered trademarks and trademarks of LSI Solutions, Inc. Copyright © 2017, LSI SOLUTIONS®. All Rights Reserved.

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

Symbol Glossary: www.lsisolutions.com/symbols

AVR NAVIGATOR™ TECHNOLOGY GUIDE

 READ PRODUCT INSERT THOROUGHLY BEFORE USE

FIG. 1

AVR NAVIGATOR™


FIG. 2

AORTIC ANNULUS

EXPOSE
enhanced
viewing

ALIGN
optimized
orientation

FIXATE
tissue
under tension



FIG. 3


AVR NAVIGATOR™ RETRACTOR WITH DELIVERY DEVICE OVERVIEW

Each sterile package contains one single-patient-use AVR NAVIGATOR™ (shown separated in FIG. 3). Each retractor incorporates three removable snare targets ①, which support wire snares ② (shown in closeup insert with target removed). These wire snares facilitate passage of stay sutures in the commissures between the right non-coronary (RN), left non-coronary (LN), and left-right (LR) sinuses of Valsalva through the collapsed retractor frame ③ and suspension tubes ④ to exit from their corresponding clamps ⑤. The snares are sequentially pulled out (RN, then LN, then LR) by withdrawing their attached suture knobs ⑥. After the retractor frame is appropriately deployed and the three commissural suture clamps are locked, pulling the cinch knob ⑦ tensions the green-and-white cinch suture within the suspension tube and the yellow cinch clamp ⑧ to secure the retractor frame in place. The delivery device is attached at the shaft tip ⑨ to the collapsed retractor frame by the integrated sliding pin ⑩ that traverses a distal jaw ⑪ and a proximal jaw ⑫. Four clamp holders ⑬ are incorporated onto the device's white handle ⑭. Clamps are removed from their holders during retractor deployment. After removal of the yellow lever stop ⑮, squeezing the purple lever ⑯ approximately halfway partially mobilizes the retractor within the aortic root. Fully squeezing the lever completely releases the retractor and enables delivery device removal from the surgical site.

INDICATIONS

The AVR NAVIGATOR™ is indicated for use in surgical procedures for the retraction and stabilization of cardiac tissue.

FIG. 4

AVR NAVIGATOR™ RETRACTOR FRAME DESCRIPTION

The self-retaining flexible contoured plastic retractor frame features three arches that gently fit into their corresponding sinuses of Valsalva. The cinch tube, with its integrated green-and-white suture, is used to hold the retractor frame open. The delivery device is attached to the retractor frame using the Attachment Site Hole located between the L-Arch and the LN tube.


COMMISSURAL STAY SUTURE PATHS

FIG. 5

RETRACTOR ENGAGEMENT

FIG. 6

3 STEPS: 1. **THREAD** 2. **INSERT** 3. **EXPAND**

Step 1: **THREAD** the **RN**, **LN**, and then **LR** commissural stay sutures

1.1 **PUSH-OUT** **RN** designated target

1.2 **PASS** **RN** suture tails through **RN** wire snare

1.3 **PULL KNOB** to thread **RN** stay suture through **RN** tube

*Repeat steps 1.1-1.3 above for **LN** and then **LR**. Do not tangle or intertwine sutures.*

Step 2: **INSERT** collapsed retractor frame

2.1 **GRASP** the ends of the three commissural stay sutures

2.2 **TENSION** the stay sutures and *take up slack* to guide retractor frame placement

2.3 **SLIDE** retractor frame gently over lightly tensioned suture toward aortic annulus

NOTE: Retractor frame's N-Arch typically enters the aortotomy first.

Step 3: **EXPAND** — lock **RN**, **LN**, then *mobilize* R-Arch, then lock **LR** and **CINCH**

3.1.A **ORIENT & VIEW** **RN** tube near **RN** commissure; N-Arch in non-coronary sinus

3.1.B **TENSION & LOCK** **RN** suture using suture clamp

3.1.C **DETACH** **RN** suture clamp

3.2.A **ORIENT & VIEW** **LN** tube near **LN** commissure; L-Arch in left coronary sinus

3.2.B **TENSION & LOCK** **LN** suture using suture clamp

3.2.C **DETACH** **LN** suture clamp

3.3.A **DETACH** **CINCH** clamp

3.3.B **DETACH** lever stop

3.3.C **HALF SQUEEZE** lever to mobilize the R-Arch by releasing delivery device from the Attachment Site Hole

3.3.D **RELEASE** lever after mobilizing R-Arch

3.4.A **ORIENT & VIEW** **LR** tube near **LR** commissure; R-Arch in right coronary sinus

3.4.B **TENSION & LOCK** **LR** suture using suture clamp

3.4.C **DETACH** **LR** suture clamp

3.5.A **FULLY SQUEEZE** lever


3.5.B **DISENGAGE** the retractor frame from the delivery device

3.5.C **REMOVE DELIVERY DEVICE** from surgical field


3.5.D **TENSION & LOCK** **CINCH** clamp to secure retractor frame in place

FIG. 7


STEP 1: *THREAD* commissural stay sutures **RN**, **LN**, and **LR**.
PUSH-OUT / PASS / PULL KNOB

PUSH-OUT white target.

PASS RN suture tails through **RN** wire snare.


PULL KNOB to thread suture in **RN** tube.


Repeat steps 1.1-1.3 for **LN** and for **LR** threading. Do not tangle or intertwine sutures.
NOTE: Discard targets and suture knobs with wire snares after removing them.

FIG. 8


STEP 2: *INSERT* collapsed retractor frame into prepared aorta.
GRASP / TENSION / SLIDE

GRASP the ends of the three commissural stay sutures.

TENSION the stay sutures and take up slack to guide retractor frame placement.

SLIDE retractor frame gently over lightly tensioned suture down toward aortic annulus.

NOTE: The retractor frame's N-Arch typically enters the aortotomy first.

STEP 3: **EXPAND** RETRACTOR FRAME

3.1 **LOCK** RN

3.2 **LOCK** LN


3.3 **MOBILIZE** R-Arch

3.4 **LOCK** LR


3.5 **LOCK** CINCH

FIG. 9


3.1 - **LOCK** RN

3.1.A ORIENT & VIEW RN tube near RN commissure (see cyan highlights); N-Arch in non-coronary sinus.

3.1.B TENSION & LOCK RN suture using suture clamp. **NOTE:** See FIG. 18 as a locking reference.


3.1.C DETACH RN suture clamp.


NOTE: The **DETACH** step may precede the **TENSION & LOCK** step if preferred.

FIG. 10


3.2 - **LOCK** LN

3.2.A ORIENT & VIEW LN tube near LN commissure; L-Arch in left coronary sinus.


3.2.B TENSION & LOCK LN suture using suture clamp.


3.2.C DETACH LN suture clamp.

FIG. 11


3.3 - MOBILIZE R-Arch

3.3.A **DETACH** CINCH clamp.

3.3.B **DETACH** lever stop.


3.3.C **HALF SQUEEZE** lever to mobilize the R-Arch by releasing delivery device from the Attachment Site Hole.


3.3.D **RELEASE** lever after mobilizing R-Arch.

FIG. 12


3.4 - LOCK LR

3.4.A **ORIENT & VIEW** LR tube near LR commissure; lift device tip right and up to place R-Arch at right coronary sinus.


3.4.B **TENSION & LOCK** LR suture using suture.

3.4.C **DETACH** LR suture clamp.

FIG. 13


3.5 - LOCK CINCH


3.5.A **FULLY SQUEEZE** lever.

3.5.B **DISENGAGE** the retractor frame from the delivery device.

3.5.C **REMOVE DELIVERY DEVICE** from the surgical field.


3.5.D **TENSION & LOCK CINCH** clamp to secure retractor frame in place.

FIG. 14

SUSPENSION

NOTE: RAM® RING suture management system not included.

Place suspension tubes under tension to optimize the orientation of the aortic annulus relative to the surgical access site.

FIG. 15

REMOVAL

15A **CUT**

15A. CUT all suspension tubes near their suture clamps.

15B **WITHDRAW**

15B. WITHDRAW the retractor frame over its commissural stay sutures.

15C **PULL OUT**


15C. PULL OUT the cinch suspension tube and the green and white cinch suture.

15D **FOLD OPEN**

15D. FOLD OPEN the retractor frame around the annular sutures. Fully remove retractor.

FIG. 16


Sound aortic annulus, then select corresponding size retractor frame.

1:1 Scale

FIG. 17

Identify and mark coronary ostia



Clockwise pattern for threading and locking commissural stay sutures

FIG. 18

SUTURE CLAMP: Holding mechanism is on the back

TO LOCK: Fingers on back, squeeze together back of clamp

TO UNLOCK: Fingers on front, while pushing both thumbnail latches forward, squeeze together front of clamp

FIG. 19

Purple Lever Utilization to Mobilize and Release Retractor Frame


Full Out
(Retractor fully retained)

Half Squeeze
(Opens distal jaw)

Release
(Readvances sliding pin)

Full Squeeze
(Fully disengage)

WARNINGS

- Federal (USA) law restricts this device to sale, distribution and use by, or on, the order of a physician.
- Surgical procedures should only be performed by physicians having adequate training and familiarity with surgical techniques. In addition, medical literature should be consulted relative to techniques, complications and hazards prior to the performance of surgical procedures.
- Read and become familiar with all instructions, warnings, and precautions before using this product.
- Do not resterilize. The AVR NAVIGATOR™ is designed and intended for single patient use only. Do not reuse, reprocess, or resterilize this product. The performance of the AVR NAVIGATOR™ after cleaning or other reprocessing has not been verified and is not supported by LSI SOLUTIONS®, Inc. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Do not use if package or product is damaged. Discard any open and unused (unsealed), expired or damaged devices or devices in damaged primary packaging.
- The AVR NAVIGATOR™ retractor frame must be carefully deployed, applied, and gently opened to avoid damaging the wall of the aorta.
- Do not leave any foreign material (e.g. suture remnant, targets, etc.) unattached in areas potentially exposed to circulating blood.
- Direct contact between sensitive tissue structures and foreign materials can lead to tissue injury or damage.
- Do not put excessive tension on the suture. Excessive suture tensioning can damage tissue or cause suture breakage.
- Leaving tissue under tension for a prolonged period of time can cause tissue ischemia or damage. This device is for temporary use during a surgical procedure by operators trained and competent in related techniques.
- It is the responsibility of the user to dispose of the device and packaging in accordance with applicable regulations and hospital procedures.

PRECAUTIONS

- Do not use devices if they do not satisfactorily perform their intended function or if they have physical damage. Avoid mechanical impact or overstressing the devices.
- Ensure proper surgical techniques as needed for the handling of tissue and suture for use of the AVR NAVIGATOR™ retractor frame.
- Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure that function is not compromised.
- Care should be taken to avoid tangling of sutures and suspension tubes for proper management.

CONTRAINDICATIONS

- Use only with LSI specified suture: 3-0 Polypropylene, 2-0 Polypropylene, 2-0 Polyester, or 2-0 Silk.
- This product is not intended for use outside of the aortic root.
- Do not use this product if patient has a very fragile, severely calcific, or porcelain aorta.

AVR NAVIGATOR™ PRODUCT ORDERING			SUPPLIED: STERILE
	REORDER	PRODUCT	DESCRIPTION
 x 6	 081446	AVR NAVIGATOR™ SMALL RETRACTOR WITH DELIVERY DEVICE	Box of 6
 x 6	 081452	AVR NAVIGATOR™ MEDIUM RETRACTOR WITH DELIVERY DEVICE	Box of 6
 x 6	 081458	AVR NAVIGATOR™ LARGE RETRACTOR WITH DELIVERY DEVICE	Box of 6

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The LSI logo, LSI SOLUTIONS, AVR Navigator, RAM and Perfect Performance Policy are registered trademarks and trademarks of LSI Solutions, Inc. Copyright © 2023, LSI SOLUTIONS®. All Rights Reserved.

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LSI SOLUTIONS®
7796 Victor-Mendon Road
Victor, New York 14564 U.S.A.
Phone: 585.869.6600
Customer Service: 866.575.3493
Technical Support: 866.428.9092
Fax: 585.742.8086
www.lsisolutions.com

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Solutions, Inc. Copyright © 2020, LSI SOLUTIC
Патенты: www.lsisolutions.com/patents Логотип LSI, LSI SOLUTIONS, Perfect Performance Policy, 3D и RAM являются зарегистрированными торговыми марками и торговыми марками LSI Solutions, Inc. Copyright © 2017, LSI SOLUTIONS®. Все права защищены.
Логотип LSI, LSI SOLUTIONS, AVR Navigator, RAM и Perfect Performance Policy являются зарегистрированными торговыми марками и торговыми марками LSI Solutions, Inc. Copyright © 2023, LSI SOLUTIONS®. Все права защищены.