Cobalt-Chromium Alloy Powder

【Product Name】Cobalt Chromium Alloy Powder

【Type Specification】DPR-M01

【Quantity】5kg、10kg

【Application Scope】

For the production of dental prostheses in the selective laser melting technology.

For the production of metal inner crowns and bridges for metal porcelain restorations.

For the Production of metal full crowns, bridges, inlays and posts for metal restorations.

For the Production of removable partial denture, clasps and full denture.

【Taboo Syndrome】It is prohibited for clinical use in patients with a history of allergy to one or more of the ingredients of this product.

【Chemical Composition】

1. The chemical composition (wt%) shall comply with the requirements of Table 1:

Table1 Chemical composition of alloy powder

Co	Cr	Mo	W	Si
64±2	25±2	5±1	5±1	1±1

2. The content of other residual elements should meet the requirements of Table 2.

Table 2 Other residual elements (wt%)

O	N	H	Fe
< 0.1	<0.2	< 0.1	<0.5

3.The content of harmful ingredients should meet the requirements of Table 3:

Table 3 harmful ingredients (wt%)

Ni	Cd	Be	Pb
< 0.1	<0.02	< 0.02	<0.02

【Performance Requirement】See product technical requirements for dental laser selective melting cobalt-chromium alloy powder for details.

【Print Device Parameters】

Table 4 3D-Print machine compatibility parameters

Parameters	Numerical Value
Laser Power	100~280W
Spot Diameter	40~110μm

Scan Spacing	0.075~0.11mm
Powder Thickness	20~40μm
Scanning Speed	600~1200mm/s

【Use Instructions / Process Sequence】

Dental processing steps 1: Printing, performed by a dental technician

Model design: Follow the design rules provided. Note that there are sufficient connectors between the crown and the bridge. If the patient has molar symptoms, it is recommended to increase the design thickness.

Repair the system as steps: Performed by trained personnel.

Document preparation: Use the software provided by the technician (such as Magics, Cambridge, etc.) or the processing software provided by the equipment manufacturer to prepare the restoration model, including the placement of parts, support design, slice allocation and printing parameters, and follow the user manual of the software.

Prosthesis printing:

- 1) check whether the scraper is intact.
- 2) set the voltage of the cycle power system according to the printing equipment operation manual.
- 3) clean all the surfaces and prisms in the equipment processing warehouse.
- 4) if the circulatory system gives an alarm, the filter element needs to be replaced.
- 5) put the selected dental powder into the processing cabin (use automatic screening system or manual screen to screen the powder).
- 6) the printing process is carried out according to the standard operation procedure in the operation manual provided by the equipment manufacturer, and the recommended material is used to print the parameters.
- 7) start the print program and start printing.

Note: if the restoration is discarded and reprinted if it is interrupted during the printing process (such as the interruption caused by lack of powder); if the restoration and support are broken during the printing process, stop printing and discard the restoration, design a stronger support structure and reprint. The use of contaminated powder is strictly prohibited.

Sandblasting before de-stress: the restoration on the substrate is sandblasted with spherical ceramic powder (diameter 0.125-0.250mm) at 2.5-3.5bar pressure until the surface is all bright.

De-stress annealing: The prosthesis is stress-relieved annealed in the annealing furnace, and the stress-relieved annealing process is as follows:

Note: Wear protective clothing and protective gloves when turning on the hot stove!

- 1) heating the furnace from room temperature to 450°C, 60min.
- 2) holding at 450°C, 45min.
- 3) heating the furnace to 750°C, 45min.

- 4) holding at 750°C, 60min.
- 5) auto-shut off heating devices.
- 6) when the furnace is cooled to 600 °C, open the furnace door.
- 7) when the furnace is cooled to 300 °C, take the restoration from the furnace.

Note: The use of higher de-stress annealing temperature or longer holding time may increase the brittleness of the material. Softening heat treatment (if necessary, such as printing denture): under the protection of vacuum or inert gas, the substrate is exothermically treated in the furnace.

Note: When turning on the hot stove, wear protective clothing and protective gloves!

- 1) heating the furnace from room temperature to 750°C, 90min.
- 2) holding at 750°C, 45min.
- 3) heat the furnace to 1190 °C, 60min.
- 4) keep it at 1190°C, 60min (the temperature in the atmosphere box is 1190°C ±10 °C, the holding time is 60min ±20min).
- 5) turn off the heating, cool the furnace to 500°C, open the furnace door.
- 6) remove the workpiece from the furnace and cool it to room temperature.

Remove the restoration from the substrate: After the heat treatment is completed and the substrate is cooled, the restoration are removed from the substrate using tools such as a saw, wire cutter, or pliers, and the supports on the restoration are removed using tools such as pliers. Remember: Do not remove the ID tag from the restoration before it is packaged or individually marked.

Dental processing steps 2: baked porcelain, performed by dental technologists.

Surface preparation of finishes: Trim the ceramic veneer surface with a cross drill to ensure that the ceramic surface is completely polished and that the remaining wall thickness of the polished restoration shall not be less than 0.3mm. Then, the restoration was sandblasted with alumina with a particle size of 125-250 μ m under 3-4bar pressure, and the surface of the restoration was cleaned with high pressure gas. After cleaning, use clamps or similar tools to hold the restoration, do not touch the restoration!

Finishing face : The facing porcelain material and process with similar thermal expansion coefficient to DPR-M01 Co-Cr alloy powder are used.

Operate in strict accordance with the manual provided by the porcelain powder manufacturer to ensure that the restoration is stably fixed on the porcelain plate during the whole process of porcelain faking.

Pre-oxidation: It is suggested that it should be carried out according to the steps and process parameters prescribed by porcelain manufacturers.

baked porcelain: it is recommended to follow the steps and process parameters prescribed by the porcelain manufacturer.

Warning: hydrofluoric acid and strong inorganic acid may corrode metal parts!

Finishing: it is recommended to follow the steps and process parameters prescribed by porcelain

manufacturers.

Soft soldering: If necessary, only materials with a matching coefficient of thermal expansion with DPR-M01 cobalt chromium alloy powder shall be used, and the appropriate process shall be adopted. Operate in strict accordance with the soldering material manufacturer's manual.

Laser welding: If necessary, only materials with a matching coefficient of thermal expansion with DPR-M01 cobalt chromium alloy powder shall be used, and the appropriate process shall be adopted. Operate in strict accordance with the laser welding material manufacturer's manual.

【Matters Need Attention】

1. This product is used with 3D printing equipment. Users should undergo necessary training before using it.

2. This product should not be exposed and stored in the air after opening. The lid should be tightened and kept sealed.

3. Dental materials such as porcelain powder and dental resin used in conjunction with this product should obtain a medical device product registration certificate.

4. Products should be stored and used away from food and beverages to avoid possible contamination.

5. Metal powder or dust is at risk of inhalation, sensitive skin contact, and eye contamination. When dealing with metal powder, there should be adequate ventilation, protective goggles, protective masks and other protective measures. After handling the powder, wash your hands.

6. The powder can be recycled for 30 times.

7. Please dispose of waste powder according to the local regulations.

【Storage Condition】 This product has good stability, storage environment temperature $\leq 35\text{ }^{\circ}\text{C}$, relative humidity is not more than 60%.

【Date of Manufacture】 See product label.

【Service life】 5 years

【Medical Device Registration Certificate Number】 NMPA Reg. No.20203170036

【Production License Number】 20200015

【 Registrant / Production Enterprise / Name of After-Sales Service Unit 】 DPR New Materials Technology CO., LTD.

【 Domicile/ Production Address 】 Room 301, Building 12-1, No. 26, Yongwang West Road, Daxing District, Beijing, China

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【Quantity】 2.5kg

【Application scope】

For the production of dental restorations in the selective laser melting technology

For the production of metal inner crowns and bridges for metal porcelain restorations

For the production of metal full crowns, bridges, inlays and posts for metal restorations

For the production of removable partial denture, clasps and full denture

【Contraindications】 It is prohibited for clinical use in patients with a history of allergies to one or more ingredients of this product.

【Chemical Composition】

1. The main chemical composition (wt.%) shall comply with the requirements of Table 1.

Table1 Chemical composition (wt.%)

Ti	Al	V
Bal.	5.5~6.75	3.5~4.5

2. The content of trace elements shall meet the requirements of Table 2.

Table 2 Trace elements (wt.%)

Fe	C	N	H	O	Other elements	
					single element	total
≤ 0.30	≤ 0.08	≤ 0.05	≤ 0.015	≤ 0.20	≤0.1	≤0.4

3. The content of hazardous elements shall meet the requirements of Table 3.

Table 3 Hazardous elements (wt.%)

Cd	Be	Pb
≤ 0.02	≤ 0.02	≤0.02

【Performance Requirement】 For more details, please see product technical requirements..

【3D Printing Machine Parameters】

Table 4 3D Printing machine parameters

Parameters	Range Value
Laser Power	90~340W
Spot Diameter	40~100μm

Scan Spacing	0.085~0.14mm
Layer Thickness	25~60μm
Scan Speed	800~1250mm/s

【Instructions for use】

Processing steps 1: 3D Printing, performed by a technician.

Model design: Follow the design rules of restorations. Please note that there should be sufficient connectors between the crown and the bridge. it is recommended to increase the design thickness if the patient has molar symptoms.

Preparation of the restorations: Performed by a technician.

3D printing file preparation: To prepare the restoration model by using the software provided by the lab (such as Magics, Cambridge, etc.) or the software provided by the 3D printing machine provider, including the placement, support design, assigning the slicing and 3D printing parameters, according to the manual of the software.

3D printing:

- 1) Check the operating status of the 3D printing machine.
- 2) Clean all the surfaces and lenses in the equipment processing cabin.
- 3) Check and replace the filter, if necessary.
- 4) Add the alloy powder into the processing chamber.
- 5) Assign the related 3D printing parameters recommended by the 3D printing machine manufacturer.
- 6) Start the machine.

Note: If the 3D printing process is interrupted (for example, lack of powder), need to stop the machine and prepare machine for reprinting. Do not use contaminated powder.

Sandblasting: Sandblast the restorations using spherical ceramic powder (diameter 0.125~0.250mm) at a pressure of 2.5~3.5bar.

Stress relief heat treatment: The restorations shall be stress-relieved using a suitable oven, and the annealing process is as follows:

Note: When opening the oven, please wear protective clothing and gloves!

- 1) Heating the oven from room temperature to 800°C, 60min.
- 2) Holding at 800°C, 60min.
- 3) Auto-shut off the heating oven.
- 4) When the temperature drops below 100 °C, open the oven door and take the restorations from the oven.

Remove the restoration from the substrate: After the stress relief heat treatment and cooling of the platform, remove the restorations using a band saw, wire cutting machine, rotary instruments or forceps, for example. Also remove the remaining supports using forceps.

Processing step 2: ceramic veneering, performed by a technician.

Surface preparation: Make sure that the surfaces to be veneered with porcelain are completely polished, and the remaining wall thickness of the polished restoration shall not be less than 0.3mm. Then, use alumina powder with a particle size of 125~250 µm to sandblast the restoration at a pressure of 3~4 bar, and use high-pressure gas to clean the surface of the restoration. After cleaning, use clamps or similar tools to hold the restoration without touching it!

Pre-oxidation: It is suggested that the peroxidation of the restorations should be carried out according to the process parameters prescribed by the porcelain manufacturers.

Ceramic veneer: Use veneering ceramics with a suitable CTE (ISO 9693-1), observe the instructions for use of the respective ceramics manufacturer.

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Soldering: Objects with implant interfaces may not be soldered! Parallel-walled soldering gap: max. 0.2 mm. Use suitable flux. Following the soldering, the flux residue and metal oxides must be acid cleaned and the surfaces should be cleaned with a steam blaster or by boiling in distilled water.

Laser welding: When possible, work with X-sutures and filler material. Please follow the instructions for use and safety instructions of the equipment manufacturer!

【Precautions】

- 1) Personnel who operating the 3D printing machine need to be trained.
- 2) Metal dust is harmful to health. Avoid dust formation! The opening of packages, filling of powders, grinding and blasting of dental restorations should be performed carefully and using an appropriate extraction system. Respiratory protection of type FFP3-EN 149, protective goggles with side protection (DIN EN 166), protective gloves and ESD-certified safety shoes are recommended.
- 3) In the event of contact with eyes, rinse with plenty of water. In the event of skin contact, wash with water and soap. If irritation persists, seek a physician's care.
- 4) Collect any spilled amounts mechanically with a moist rag (water or isopropanol) and dispose of in accordance with local or national statutory regulations.
- 5) Metal powders are combustible. Remove all sources of ignition.
- 6) Suitable extinguishing media: special powders against metal fires, sand.
- 7) Pay attention to safety data sheet!

【Storage Condition】 Store dry in tightly closed container, storage temperature should below 35 °C, relative humidity shall no more than 60%.

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【Address】 Room 301, Building 12-1, No. 26, Yongwang West Road, Daxing District Beijing, China

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