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This warranty is given by Abbott Medical. To make inquiries regarding this warranty, use the contact information for either the manufacturer or Australian Sponsor. Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

The benefits to you given by this warranty are in addition to your other rights and remedies under the Australian Consumer Law.

Abbott Medical warrants to the buyer that the goods, for a period equal to the validated shelf life of the goods (the "Warranty Period"), shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Unless otherwise obligated by law or mandatory consumer guarantees, under this warranty Abbott Medical will, at its option, replace or repair the goods at its factory, if the buyer, at the buyer's expense, returns the goods within the Warranty Period to Abbott Medical at the address provided with the product and after Abbott Medical confirms that the goods are defective.

To the maximum extent permitted by law, if a mandatory term is implied by law or a mandatory consumer guarantee applies to the goods, and the goods are not of a kind ordinarily acquired for personal, domestic or household use or consumption, Abbott Medical's liability for the breach of the term or guarantee is limited to, at Abbott Medical's option, either replacement or repair of the goods or payment of the costs of replacing or repairing the goods.

EXCEPT FOR MANDATORY TERMS IMPLIED BY LAW, MANDATORY CONSUMER GUARANTEES, OR AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED.

See the Terms and Conditions of Sale for further information.

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Proposition 65, a State of California voter initiative, requires the following notice:

**WARNING:** This product can expose you to chemicals including ethylene oxide, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third party trademark, which is property of its respective owner.

Pat. <http://www.abbott.com/patents>

© 2023 Abbott. All Rights Reserved.

# Response™ Electrophysiology Catheter

## Instructions for Use

### Description

The Response™ Electrophysiology Catheter is a sterile, single use, electrophysiological diagnostic catheter. The catheter has a fixed distal curve available in various styles, no lumen, and varying electrode numbers and spacing. The electrodes detect cardiac electrical signals that are transmitted to an electrocardiogram recorder for diagnosis or deliver the electrical signal from an external pacemaker to the heart for electrical pacing purposes. See label for outer diameter.

### Performance Characteristics

The Response™ Electrophysiology Catheter is designed to facilitate electrophysiology mapping of the heart chambers and transmit pacing stimuli to the catheter electrodes to induce cardiac stimulation. The Response™ Electrophysiology Catheter is compatible with introducers or sheaths with a diameter one size larger than its labeled French size.

### Intended Purpose

Electrophysiology diagnostic catheters are intended to be used for EP mapping, stimulation, and recording of a variety of cardiac arrhythmias from endocardial and intravascular sites.

### Indications for Use

Abbott electrophysiology catheters can be used in the evaluation of a variety of cardiac arrhythmias from endocardial and intravascular sites.

### Patient Target Group

The Response™ Electrophysiology Catheter is intended for patients with diagnosed or suspected cardiac arrhythmia requiring electrogram recording and/or cardiac stimulation during diagnostic EP studies. Patients must have a vascular anatomy compatible with introducers with a minimum inner diameter of 1F size larger than the catheter (5F or 6F).

**CAUTION: The safety and effectiveness of the device has not been established in pre-pubescent children. Careful consideration must therefore be given for the use of the device in pre-pubescent children.**

### Intended User

This device should be used by board-certified electrophysiologists, or EP fellows in training, in a fully-equipped and operational electrophysiology laboratory.

### Clinical Benefits

The intended clinical benefit of the Response™ Electrophysiology Catheter is to provide diagnostic information, measured by successful intracardiac electrogram recording and delivery of pacing stimuli, to aid the physician in the treatment of arrhythmias.

### Contraindications

- The catheter is contraindicated for patients with an active systemic infection, as this may increase the risk of cardiac infection.
- The catheter is contraindicated for patients with prosthetic valves, left atrial thrombus or myxoma, or interatrial baffle or patch via transseptal approach.

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### Warnings

- Misuse of this catheter and accessories may result in serious complications.
- This device should only be used with equipment that complies with international safety standards.
- Individual patient anatomy and physician technique may require procedural variations.
- This device should be used by board-certified electrophysiologists, or EP fellows in training, in a fully-equipped and operational electrophysiology laboratory.
- Any attempt to resterilize and reuse may compromise the integrity of this system. Adverse effects of using nonsterile components may include, but are not limited to:
  - Local and/or systemic infection

- Mechanical damage
- Inaccurate functionality
- Vascular perforation or dissection is an inherent risk of any electrode placement. Careful catheter manipulation must be performed in order to avoid device component damage, thromboembolism, cerebrovascular accident, cardiac damage, perforation, pericardial effusion, or tamponade.
- Catheter entrapment within the heart or blood vessels is a possible complication of electrophysiology procedures. If resistance is observed, do not force catheter. Withdraw catheter, correct difficulty, and reinsert.
- For both patients and laboratory staff, cardiac catheterization procedures present the potential for significant X-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects due to the X-ray beam intensity and duration of the fluoroscopic imaging. Carefully consider the use of this catheter in pregnant women.
- Maintain an activated clotting time (ACT) of greater than 300 sec at all times during the use of the catheter.
- The materials identified below may cause allergic or hypersensitive reactions. Persons allergic to these may suffer an allergic reaction. Before use on the patient, the patient should be counseled on the materials contained in the device and a thorough history of allergies must be discussed.
  - Polyether block amide
  - Polyurethane (epoxy) adhesive
  - Platinum iridium alloy

## Precautions

- Read directions prior to use.
- United States law restricts this device to sale by or on order of a physician.
- Do not alter this device.
- Store in a cool, dark, dry place.
- Use care to isolate any unused connector pins of the electrogram cable. This will reduce the chances of developing accidental current pathways to the heart.
- Inspect the package prior to use. Do not use if the package is open, damaged, or expired.
- Personnel handling the electrophysiology catheter should wear gloves.
- To maintain optimal patient safety and electrode catheter integrity, do not wipe this catheter with alcohol.
- Excessive bending or kinking of the catheter may cause damage to the catheter.
- Inspect the catheter for any damage or kinks. Do not use if the catheter is damaged, kinked, or does not hold its curve.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Catheter materials are not compatible with magnetic resonance imaging (MRI).
- For catheters equipped with a cable connector, use with the appropriate Abbott cable.
- Remove the catheter from its package. Inspect the electrodes and catheter carefully for integrity, continuity of leads, and overall condition.
- Catheter advancement should be performed under fluoroscopic guidance to minimize the risk of cardiac damage, perforation, or tamponade. Compatible navigation and real time visualization systems may not be considered.
- Not recommended for long term pacing. This device should not be used for life sustaining pacing.
- After use, this device, its accessories and packaging should be appropriately classified for disposal, e.g. biohazard, sharps, non-hazardous waste etc., and carefully disposed of in compliance with facility procedures and applicable laws and regulations.

## Potential Adverse Events

Risks associated with electrical stimulation may include, but are not limited to, the induction of arrhythmias, such as atrial fibrillation (AF), bradycardia, ventricular tachycardia (VT), ventricular fibrillation (VF), and the need for cardioversion.

The risks associated with using diagnostic electrophysiology catheters include:

| Adverse Event Category                                                           | Adverse Event Subcategory       |
|----------------------------------------------------------------------------------|---------------------------------|
| Hemorrhage/bleeding                                                              | N/A                             |
| Hematoma                                                                         | N/A                             |
| Ecchymosis                                                                       | N/A                             |
| Superficial tissue injury                                                        | N/A                             |
| Cardiac perforation                                                              | Cardiac tamponade               |
|                                                                                  | Pericardial effusion            |
| Cardiovascular injury                                                            | Valvular damage                 |
|                                                                                  | Great vessel perforation        |
| Cerebrovascular injury                                                           | Asymptomatic cerebral emboli    |
|                                                                                  | Stroke/cerebrovascular accident |
|                                                                                  | Transient ischemic attack       |
| Death                                                                            | N/A                             |
| Electric shock                                                                   | N/A                             |
| Embolism                                                                         | Air embolism                    |
|                                                                                  | Foreign body embolism           |
|                                                                                  | Pulmonary embolism              |
|                                                                                  | Thromboembolism                 |
| Immunological reaction                                                           | Anaphylaxis                     |
| Infections                                                                       | Pericarditis                    |
|                                                                                  | Endocarditis                    |
|                                                                                  | Pneumonia                       |
|                                                                                  | Sepsis                          |
| Myocardial ischemia                                                              | Myocardial infarction           |
|                                                                                  | Transient ST elevation          |
|                                                                                  | Angina                          |
|                                                                                  | Coronary artery spasm           |
| Peripheral vascular injury due to vascular access                                | Pseudoaneurysm                  |
|                                                                                  | Arteriovenous fistula           |
|                                                                                  | Vascular dissection             |
| Thrombosis/thrombus                                                              | N/A                             |
| Hypotension                                                                      | N/A                             |
| Radiation injury                                                                 | N/A                             |
| Dislodgement of implantable cardioverter defibrillator or permanent pacing leads | N/A                             |
| Pain                                                                             | N/A                             |

## Connection to Other Equipment

Pericardial Mapping System connected to a commercially available EP recording system and navigation and visualization system using the connection cable. All systems must be patient isolated. For instructions regarding the use of these systems with the catheter, refer to the system's instructions for use.

Table 1. Compatible Systems

| System Device       | Connect via (Model Number)                                        | Catheter Type                                                 |
|---------------------|-------------------------------------------------------------------|---------------------------------------------------------------|
| EP Recording System | Response™ Electrophysiology Extension Cable (401970) <sup>1</sup> | Bipolar and Quadripolar Response™ Electrophysiology Catheters |
|                     | Response™ Electrophysiology Extension Cable (401972) <sup>1</sup> | Hexapolar and Decapolar Response™ Electrophysiology Catheters |

Table 1. Compatible Systems

| System Device       | Connect via (Model Number)                                                | Catheter Type                                                 |
|---------------------|---------------------------------------------------------------------------|---------------------------------------------------------------|
|                     | Response™ Electrophysiology Extension Cable (401976) <sup>1</sup>         | Hexapolar and Decapolar Response™ Electrophysiology Catheters |
| EnSite™ X EP System | EnSite™ X EP System Catheter Connection Cable, 4 Pin (D-ENS-L/R-BQ-CBL)   | Bipolar and Quadripolar Response™ Electrophysiology Catheters |
|                     | EnSite™ X EP System Catheter Connection Cable, 10 Pin (D-ENS-L/R-HOD-CBL) | Hexapolar and Decapolar Response™ Electrophysiology Catheters |

<sup>1</sup> Alternate applicable models may be used as listed in Abbott's Electrophysiology Product Catalog.

Directions

- 1. Remove the catheter from its package. Inspect the electrodes and catheter carefully for integrity, continuity of leads, and overall condition.
- 2. Insert the catheter by using the appropriate percutaneous catheter introducer.
- 3. For catheters equipped with a cable connector, use with the appropriate Abbott Medical cable.
- 4. The catheter should be passed from a peripheral vessel to the desired intracardiac position using a navigation system, fluoroscopy, or intracardiac echocardiography (ICE). Manipulate the catheter to different parts of the anatomy to acquire desired intracardiac electrogram.

Packaging and Shelf-Life

The catheter packaging is designed to prevent crushing of the product, minimize product exposure to the atmosphere, and provide for aseptic product transfer. The catheter is sterilized with ethylene oxide and supplied sterile if the package is unopened and undamaged. Do not resterilize or reuse. It is recommended that the product remain in the unopened package until time of use. The expiration date is marked on the outside of the package. The instructions for use is recyclable.

Incident Reporting

If, in the course of use of this device, you have reason to believe that a serious incident occurred, please report it to the manufacturer. For customers in the European Union, report the serious incident to your national authority as well as to the manufacturer.

Clinical Data Summary

For European Union countries: A summary of the safety and clinical performance (SSCP) for these devices will be published at <https://ec.europa.eu/tools/eudamed> by the Notified Bodies in accordance with the requirements of EU MDR 2017/745 and associated timelines. Search for the device using the UDI-DI provided on the outer packaging of the device. The UDI-DI is the 14-digit number between the first two sets of parentheses. For more information about finding and using the UDI-DI, refer to the Frequently Asked Questions document found at [medical.abbott/manuals](https://medical.abbott/manuals).

Регистрация МИ в Росздравнадзоре  
[www.nevacert.ru](http://www.nevacert.ru) | [info@nevacert.ru](mailto:info@nevacert.ru)


Warranty


Abbott Medical warrants that its products shall be free from defects in materials and workmanship under normal use. This warranty does not exceed the “Expiration” date stated on any product labeling. The authorized uses and approved methods of use of each of our products is set forth in the related “Instructions for Use” that accompany each product. Abbott Medical disclaims any responsibility and liability for the use of its products in a manner that has not been authorized or approved. Abbott Medical’s liability under this warranty is limited to replacing its products. The foregoing warranty excludes and is in lieu of all other warranties whether expressed or implied including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Abbott Medical disclaims any liability for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this product, other than as expressly provided by specific law. Abbott Medical neither assumes nor authorizes any other person to assume for it any other or additional liability

for loss, damage, or expense in connection with this product. For more details please review complete Abbott Medical warranty policy available from Abbott Medical or on the back of an Abbott Medical invoice.

## Symbols

The symbols below and harmonized symbols may be found on the product or product label. For harmonized symbols, refer to the Universal Symbols Glossary at [medical.abbott/manuals](http://medical.abbott/manuals).

| Symbol                                                                                                                                                        | Definition                                          |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------|
|                                                                              | Manufacturer                                        |
|                                                                              | Date of manufacture                                 |
|                                                                              | Country of manufacture (Costa Rica)                 |
|                                                                              | Use-by date                                         |
|                                                                              | Authorized representative in the European Community |
|                                                                              | Batch code                                          |
|                                                                             | Catalog number                                      |
|                                                                            | Sterilized using ethylene oxide                     |
|                                                                            | Do not resterilize                                  |
|                                                                            | Do not use if package is damaged                    |
|                                                                            | Protect from heat and radioactive sources           |
|                                                                            | Keep dry; keep away from rain                       |
|                                                                            | Keep away from sunlight                             |
|                                                                            | Do not reuse                                        |
|                                                                            | Consult instructions for use                        |
|                                                                            | Quantity                                            |
| <div>Регистрация МИ в Росздравнадзоре<br/><a href="http://www.nevacert.ru">www.nevacert.ru</a>   <a href="mailto:info@nevacert.ru">info@nevacert.ru</a></div> | Length                                              |
|                                                                            | Electrophysiology Catheter                          |
|                                                                            | Double sterile barrier system                       |
|                                                                            | Fixed curve diagnostic catheter                     |
|                                                                            | Spacing                                             |
|                                                                            | Number of electrodes                                |

| Symbol                                                                                                        | Definition                                                                                                                                                                                                                                                        |
|---------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                              | Recommended cables                                                                                                                                                                                                                                                |
|                              | Medical Device                                                                                                                                                                                                                                                    |
|                              | Unique Device Identifier                                                                                                                                                                                                                                          |
|                              | CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a physician.                                                                                                                                                                      |
|                              | Importer                                                                                                                                                                                                                                                          |
|                              | Conformité Européenne (European Conformity). Affixed in accordance with European Council Regulation 2017/745. Hereby, Abbott Medical declares that this device is in compliance with the essential requirements and other relevant provisions of this regulation. |
|                             | United Kingdom Conformity Assessed                                                                                                                                                                                                                                |
|                            | Responsible person in the United Kingdom                                                                                                                                                                                                                          |
|                            | Ukrainian mark of conformity to the technical regulations, where 101 is the assessment body.                                                                                                                                                                      |
| <br>medical.abbott/manuals | Follow instructions for use on this website                                                                                                                                                                                                                       |

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Abbott Medical UK Limited