Epic™ Plus/Epic™ Plus Supra

Stented Porcine Tissue Valves

Instructions for Use

EN: English	1
BG: Български	15
CS: Čeština	30
DA: Dansk	43
DE: Deutsch	57
EL: Ελληνικά	71
ES: Español	86
ET: Eesti keel	100
FI: Suomi	113
FR : Français	127
HR: Hrvatski	141
HU: Magyar	154
IT: Italiano	167
LT: Lietuvių	181
LV: Latviski	194
NL: Nederlands	207
NO: Norsk	221
PL: Polski	235
PT: Português	249
RO: Română	263
RU: Русский	277
SK: Slovenčina	293
SL: Slovenščina	306
SR: Srpski	319
SV: Svenska	333
TR: Türkçe	347
UK: Українська	360

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

Pat. <http://www.abbott.com/patents>

© 2023 Abbott. All Rights Reserved.

EN: English

Epic™ Plus/Epic™ Plus Supra

Stented Porcine Tissue Valves

Instructions for Use

Device Description

The Epic™ Plus valve and the Epic™ Plus Supra valve are bioprosthetic heart valves manufactured from select porcine aortic valve cusps. The cusps are carefully matched for optimum leaflet coaptation and hemodynamics.

Figure 1. Epic™ Plus aortic and Epic™ Plus mitral valves



Figure 2. Epic™ Plus Supra aortic valve

Following tissue fixation, the tissue is mounted on the FlexFit™ polyester-covered flexible acetal copolymer stent. The stent is a low-profile design with a scalloped shape.

The stent is covered with knitted polyester fabric. The sewing cuff on Epic™ Plus standard valves is formed from a braided polyester filler material that is covered by the knitted polyester fabric cover. The sewing cuff on Epic™ Plus Supra valves is formed by enclosing a molded silicone elastomer within the same knitted polyester fabric as the Epic™ Plus standard valve. The sewing cuff also contains three suture markers to facilitate valve placement. A stainless-steel wire is located under the sewing cuff and at each post stent for radiopacity.

A bovine pericardial tissue strip is attached to the outflow edge of the valve. This strip protects the leaflets as they open and close. The pericardial strip and the porcine valve cusps are preserved and crosslinked in a glutaraldehyde solution. Glutaraldehyde, formaldehyde, and ethanol are used in the valve sterilization process. The valve is attached to a holder and packaged in a sealed jar containing a formaldehyde solution. The Epic™ Plus and Epic™ Plus Supra valves are supplied sterile and non-pyrogenic.

Non-clinical testing has demonstrated that the Epic™ Plus and Epic™ Plus Supra valves are MR Conditional. Refer to the MRI Safety Information section of this manual for further information.

The Epic™ Plus valves are available in the aortic and mitral sizes indicated in the following table. Epic™ Plus standard valves are designed to allow intra-annular placement of the inflow edge of the valve with supra-annular placement of the sewing cuff.

Epic™ Plus Supra valves are available in the aortic sizes indicated in the following table. Epic™ Plus Supra valves are designed for supra-annular implantation of both the valve and the sewing cuff. The sewing cuff is the only difference between the Epic™ Plus valve and the Epic™ Plus Supra valve.

The Epic™ Plus and Epic™ Plus Supra valves are processed using the Linx™ anticalcification treatment.

Table 1. Catalog Number Descriptions and Reference Dimensions

Catalog Number	Tissue Annulus Diameter (mm)	Aortic/Ventricular Protrusion (mm)	Total Height(mm)	Cuff Outer Diameter (mm)	Stent Inner Diameter (mm) ¹
Aortic Heart Valves, Epic™ Plus					
E200-21A	21	9	14	25	18.7
E200-23A	23	9	15	27	20.8
E200-25A	25	10	16	29	22.6
E200-27A	27	11	17	31	24.5
E200-29A	29	12	19	33	26.3
Aortic Heart Valves, Epic™ Plus Supra					
ESP200-19	19	11	14	25	18.7
ESP200-21	21	11	15	28	20.8
ESP200-23	23	13	16	29	22.6
ESP200-25	25	13	17	31	24.5
ESP200-27	27	14	19	33	26.3
ESP200-29	29	15	20	35	28.4
Mitral Heart Valves, Epic™ Plus					
E200-25M	25	9	16	33	22.6

¹ Stent Internal Diameter does not account for the tissue/fabric covering the stent which may vary in thickness between approximately 1 to 2 mm.

Table 1. Catalog Number Descriptions and Reference Dimensions

Catalog Number	Tissue Annulus Diameter (mm)	Aortic/Ventricular Protrusion (mm)	Total Height(mm)	Cuff Outer Diameter (mm)	Stent Inner Diameter (mm)
E200-27M	27	9	17	35	24.5
E200-29M	29	10	19	37	26.3
E200-31M	31	10	20	39	28.4
E200-33M	33	11	20	41	30.3

Intended Purpose

The Epic™ Plus Valve and Epic™ Plus Supra Valve are intended to replace native mitral or aortic heart valves or previously implanted mitral or aortic prosthetic valves.

Indications

The Epic™ Plus valve and Epic™ Plus Supra valve are indicated for use as a replacement valve in patients with a diseased, damaged, or malfunctioning mitral or aortic heart valve or as a replacement for a prosthetic heart valve.

Target Patient Population

Patients with a diseased, damaged, or malfunctioning aortic or mitral valve.

Intended User

Physicians trained in surgical mitral and/or aortic valve replacement.

Intended Clinical Benefits

Based on clinical evidence, the intended clinical benefits for the Epic™ Plus family of valves are:

- Restoration of hemodynamic valve performance resolving valvular stenosis and/or valvular regurgitation
- Improvement in NYHA Classification

Summary of Safety and Clinical Performance (SSCP)

A summary of the safety and clinical performance (SSCP) for this device is available at <https://ec.europa.eu/tools/eudamed>. Search for the device using the basic UDI-DI listed in the following table. This is the SSCP location after the launch of the European Database on Medical Devices/Eudamed.

Table 2. Epic™ Plus Tissue Heart Valves Basic UDI-DI

Basic UDI-DI	Basic UDI-DI Description
5415067SHV1001PD	Surgical Heart Valve Products, Tissue Heart Valves

Performance Characteristics

The Epic™ Plus Aortic, Supra and Mitral Valves consist of the following components:

- Stent
- Cuff
- Leaflets
- Pericardium shield
- Radiopaque elements

The Flexfit™ stent is flexible acetyl copolymer with radiopaque wires wrapped at the annulus of the stent and at each post. The stent contains three intra-annular leaflets made from non or left coronary porcine cusp segments which are assembled in composite fashion. The silicone-filled cuff allows for supra-annular implantation, and the flexible polyester-filled cuff mitigates paravalvular leak and easily fits patient anatomy for the aortic and mitral positions. The low overall valve height mitigates risk of coronary obstruction in the aortic/supra position and LVOT obstruction in the mitral position. The pericardial shield provides a tissue-to-tissue interface to help prevent the risk of abrasion. The Epic™ Plus valve design improves long-term hemodynamic function.

The clinical investigation of the first generation Epic valve supports the safety and effectiveness of the Epic™ Plus valve and the Epic™ Plus Supra valve.

Table 3. Effectiveness Outcomes at One Year Follow-Up Visit, Hemodynamic Results - All Epic™ Aortic Valves

All aortic subjects entered into study: N=586

Hemodynamic Parameter	21 mm	23 mm	25 mm	27 mm	29 mm*
Mean Gradient (mmHg)	n= 49	n=120	n=121	n= 36	n= 10
~ Mean ± SD	19.1 ± 8.2	13.9 ± 6.0	12.1 ± 5.1	11.4 ± 4.1	7.5 ± 3.3
~ Min, Max	3.1, 43.5	1.7, 35.0	3.7, 34.3	6.5, 26.3	2.7, 12.7
EOA (cm²)	n= 46	n=118	n=121	n= 35	n= 10
~ Mean ± SD	1.0 ± 0.3	1.4 ± 0.5	1.5 ± 0.5	1.6 ± 0.4	2.4 ± 1.1
~ Min, Max	0.5, 2.3	0.5, 3.5	0.2, 3.3	0.8, 2.7	1.2, 4.6
Regurgitation (n, %)	n= 56	n=130	n=128	n= 38	n= 10
~ None	47 (84%)	103 (79%)	92 (72%)	28 (74%)	9 (90%)
~ Trivial	7 (13%)	21 (16%)	29 (23%)	9 (24%)	1 (10%)
~ Mild	2 (4%)	6 (5%)	6 (5%)	1 (3%)	0 (0%)
~ Moderate	0 (0%)	0 (0%)	1 (1%)	0 (0%)	0 (0%)
~ Severe	0 (0%)	0 (0%)	0 (0%)	0 (0%)	0 (0%)

Table 3. Effectiveness Outcomes at One Year Follow-Up Visit, Hemodynamic Results - All Epic™ Aortic Valves

All aortic subjects entered into study: N=586

Hemodynamic Parameter	21 mm	23 mm	25 mm	27 mm	29 mm*
~ Unknown	0 (0%)	0 (0%)	0 (0%)	0 (0%)	0 (0%)

*Includes 2 subjects at greater than one year visit

n= number of subjects evaluated

Table 4. Effectiveness Outcomes at One Year Follow-Up Visit, Hemodynamic Results - All Epic™ Mitral Valves

All mitral subjects entered into study: N=205

Hemodynamic Parameter	27 mm*	29 mm	31 mm	33 mm*
Mean Gradient (mmHg)	n= 30	n= 41	n= 26	n= 24
~ Mean ± SD	6.1 ± 2.9	5.5 ± 1.7	4.8 ± 1.4	4.1 ± 1.6
~ Min, Max	2.7, 14.0	2.9, 10.0	2.6, 8.3	1.5, 7.9
EOA (cm²)	n= 16	n= 26	n= 15	n= 22
~ Mean ± SD	1.4 ± 0.7	1.5 ± 0.5	1.6 ± 0.3	1.5 ± 0.3
~ Min, Max	0.6, 3.1	0.6, 2.8	1.1, 2.4	1.1, 2.2
Regurgitation (n, %)	n= 30	n= 45	n= 28	n= 25
~ None	30 (100%)	41 (91%)	23 (82%)	23 (92%)
~ Trivial	0 (0%)	1 (2%)	0 (0%)	1 (4%)
~ Mild	0 (0%)	3 (7%)	5 (18%)	0 (0%)
~ Moderate	0 (0%)	0 (0%)	0 (0%)	1 (4%)
~ Severe	0 (0%)	0 (0%)	0 (0%)	0 (0%)
~ Unknown	0 (0%)	0 (0%)	0 (0%)	0 (0%)

*Includes 4 (27 mm) and 6 (33 mm) subjects at greater than one year visit.

n= number of subjects evaluated

Table 5. Biocor™/Epic™ Effectiveness Outcomes, 25mm Mitral Hemodynamic results (≥ 11 months)

All mitral 25 mm subjects: N= 15

Brazil Biocor™ Valve Study(n=12)

Epic™ Valve IDE (n=3)

Hemodynamic Parameter	25 mm
Mean Gradient (mmHg)	n= 15
~ Mean ± SD	6.8 ± 2.5
~ Min, Max	3.7, 13.7
EOA (cm²)	n=15
~ Mean ± SD	1.2 ± 0.3
~ Min, Max	0.5, 1.7
Regurgitation (n, %)	n= 15
~ None	11 (73.3%)
~ Trivial	0 (0.0%)
~ Mild	4 (26.7%)
~Moderate	0 (0.0%)
~Severe	0 (0.0%)
~ Unknown	0 (0.0)

n= number of subjects evaluated

Limitations

No known limitations.

Contraindications

None known.

Warnings and Precautions

Warnings

- Valve size selection is based on the size of the recipient annulus, and for supra-annular aortic placement, the anatomy of the sinotubular space. Implantation of an inappropriately large bioprosthesis may result in stent deformation, valvular incompetence, valve damage, diminished tissue durability, and/or damage to the surrounding tissues. The use of an inappropriately small bioprosthesis may result in suboptimal hemodynamics. Use only the E2000™ Epic™ Plus Heart Valve Sizer Set with the Epic™ Plus and Epic™ Plus Supra valves.
- Accelerated deterioration due to calcific degeneration of the Epic Plus and Epic Plus Supra valve may occur in:
 - children, adolescents, or young adults;
 - patients with altered calcium metabolism (e.g., patients with hyperparathyroidism or chronic renal failure); or
 - individuals requiring hemodialysis.

- For single use only. Do not reuse or resterilize. Attempts to resterilize the valve may result in valve malfunction, inadequate sterilization, or patient harm.
- Passage of a catheter or transvenous pacing lead through any bioprosthesis may damage the valve and is therefore not recommended.
- Persons with known history of allergies to any of the components of this device listed in the “Product Materials” section may suffer an allergic reaction to this device. Prior to its use on the patient, the patient should be counselled on the materials contained in the device and a thorough history of allergies must be discussed.

Do not use if:

- the valve has been dropped, damaged, or mishandled in any way, or if there is any sign of deterioration;
- the expiration date has elapsed;
- the tamper-evident container seal is damaged, broken, or missing, or if fluid is leaking from the packaging; or
- the storage solution does not completely cover the valve.

Precautions

- The safety and effectiveness of the Epic™ Plus and Epic™ Plus Supra valves has not been established for the following specific populations:
 - patients who are pregnant
 - nursing mothers
 - patients with chronic renal failure
 - patients with aneurysmal aortic degenerative conditions (e.g., cystic medial necrosis, Marfan’s syndrome)
 - patients with chronic endocarditis
 - patients requiring pulmonic or tricuspid valve replacement
 - children, adolescents, or young adults
- Sizers are supplied non-sterile, and must be cleaned and sterilized prior to each use. Do not use cracked, deformed, or damaged sizer set components.
- Do not pass the flanged portion of the valve replica sizing tool through the annulus.
- Do not place the non-sterile exterior of the valve container in the sterile field.
- Do not expose the valve to solutions other than the formaldehyde valve storage solution in which it was shipped, the sterile isotonic saline solution used during the rinsing procedure, or the sterile isotonic saline solution used to irrigate the valve.
- Do not add antibiotics to either the formaldehyde valve storage solution or the rinse solution.
- Do not apply antibiotics to the valve.
- Do not allow the valve tissue to dry. Place the valve in sterile isotonic saline rinse solution immediately upon removal from the valve storage solution. Once removed from this solution, the valve should be periodically irrigated during implantation.
- Do not use the valve if shipping temperature indicators on the product carton have turned red, or if the valve has been improperly stored in temperature conditions outside of the 5°C to 25°C (41°F to 77°F) range.
- Do not implant the valve without thoroughly rinsing as directed.
- Do not lacerate the valve tissue. If a valve is damaged, the valve must be explanted and replaced.
- Do not attempt to repair a valve. Damaged valves must not be used.
- Do not use cutting edge needles, unprotected forceps, or sharp instruments as they may cause structural damage to the valve.
- Never handle the leaflet tissue.
- Position the mitral valve in a manner to avoid commissure obstruction of the left ventricular outflow tract, and minimize any potential of commissure contact with the ventricular wall.
- Position the aortic valve so that the stent posts do not obstruct the coronary ostia.
- When implanting the Epic™ Plus heart valve, assess the suitability of the selected valve size and stent post position for a potential future valve-in-valve procedure and whether the transcatheter valve-in-valve procedure may result in left ventricular outflow tract or coronary ostia obstruction. For a future valve-in-valve procedure in an Epic™ Plus valve, refer to the instructions for use supplied with the transcatheter heart valve along with the reference dimensions in the Catalog Number Descriptions and Reference Dimensions Table (page 1) to determine compatibility. The safety and effectiveness of valve-in-valve procedures in an Epic™ Plus or an Epic™ Plus Supra valve have not been established.
- Avoid prolonged contact with the formaldehyde storage solution. Immediately after contact, thoroughly flush any skin exposed to the solution with water. In case of contact with eyes, flush with water and seek appropriate medical care.

Adverse Events

The clinical investigation of the first generation Epic™ valve supports the safety and effectiveness of the Epic™ Plus valve and the Epic™ Plus Supra valve. Between January 2003 and March 2006, seven-hundred and sixty-two (762) subjects were implanted with 791 Epic™ valve(s) at 19 investigational sites in the United States (U.S.), and three sites in Canada. Five-hundred and fifty-seven (557) subjects received isolated aortic replacement, 176 received isolated mitral replacement, and 29 received replacement of both the aortic and mitral valves. The cumulative follow-up for all subjects was 773.51 patient-years with a mean follow-up of 1.02 patient-years (s.d. = 0.71 patient-years, range 0 – 3.10 patient-years).

Potential Adverse Events

Adverse events potentially associated with the use of bioprosthetic heart valves (in alphabetical order) include:

- angina
- cardiac arrhythmias
- endocarditis
- heart failure
- hemolysis
- hemolytic anemia
- hemorrhage, anticoagulant/antiplatelet-related
- leak, transvalvular or paravalvular
- myocardial infarction
- nonstructural dysfunction (entrapment by pannus or suture, inappropriate sizing or positioning, or other)
- prosthesis regurgitation
- stroke
- structural deterioration (calcification, leaflet tear, or other)

- thromboembolism
- valve thrombosis

It is possible that these complications could lead to:

- reoperation
- explantation
- permanent disability
- death

Device Lifetime

The product is designed and tested to support a minimum lifetime of 5 years. In addition, long-term clinical data published in the literature indicates Abbott's tissue heart valves perform as intended for at least 10 years.

MRI Safety Information

Non-clinical testing has demonstrated that the Epic™ Plus heart valves are MR Conditional. A patient with one of these devices can be safely scanned in an MR system under the following conditions:

- Static magnetic field of 1.5 Tesla (1.5T) or 3.0 Tesla (3.0T)
- Maximum spatial gradient field of 19 T/m (1900 G/cm)
- Maximum MR system reported, whole-body averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode)

Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than or equal to 2°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends radially up to 0.6 cm from the device when imaged with a gradient echo pulse sequence in a 1.5T MR system.

Individualization of Treatment

Anticoagulant and/or Antiplatelet Therapy

Long-term, low-dose aspirin, unless contraindicated, is recommended for all patients with bioprosthetic valves. Long-term anticoagulant therapy, unless contraindicated, is recommended for all patients with bioprosthetic valves who have risk factors for thromboembolism.

Specific Patient Populations

Safety and effectiveness of Epic™ Plus and the Epic™ Plus Supra valves have not been established for the following specific populations:

- patients who are pregnant
- nursing mothers
- patients with chronic renal failure
- patients with aneurysmal aortic degenerative conditions (e.g., cystic medial necrosis, Marfan's syndrome)
- patients with chronic endocarditis
- patients requiring pulmonic or tricuspid valve replacement
- children, adolescents, or young adults.

Patient Counseling Information

Long-term, low-dose aspirin, unless contraindicated, is recommended for all patients with bioprosthetic valves. Long-term anticoagulant therapy, unless contraindicated, is recommended for all patients with bioprosthetic valves who have risk factors for thromboembolism.

Patients with bioprostheses who undergo dental or other procedures that are potentially bacteremic should receive endocarditis prophylactic antibiotic therapy.

Abbott Medical publishes a patient brochure. Copies of this booklet are available through your Abbott sales representative.

Packaging and Storage

As delivered, the valve is attached to a valve holder by three retaining sutures. A flexible plastic support surrounds the valve. The valve holder and support facilitate handling and manipulation of the valve during removal from the container, rinsing, and implantation.

The valve is packaged in a formaldehyde storage solution.

Store the valve in the upright position.

The valve should be stored in temperatures from 5°C to 25°C (41°F to 77°F). Do not store the valve where significant temperature fluctuations may occur.

CAUTION:

- Do not implant the valve without thoroughly rinsing as directed.
- Do not use the valve if shipping temperature indicators on the product carton have turned red, or if the valve has been improperly stored in temperature conditions outside the 5°C to 25°C (41°F to 77°F) range.

Directions for Use

Pre-Implant Handling

The Epic™ Plus and the Epic™ Plus Supra valves are supplied in a storage container with a screw-cap closure and tamper-evident seal. The contents of the container are sterile, and must be handled aseptically to prevent contamination.

WARNING:

- Do not use the valve if the expiration date has elapsed.
- Do not use the valve if fluid is leaking from the packaging.
- Do not resterilize the valve by any method.

Removing the Valve from the Outer Packaging

PRECAUTION:

- Do not place the non-sterile exterior of the valve container in the sterile field.
 - Do not expose the valve to solutions other than the formaldehyde valve storage solution in which it was shipped, the sterile isotonic saline solution used during the rinsing procedure, or the sterile isotonic saline solution used to irrigate the valve.
 - Do not add antibiotics to either the formaldehyde valve storage solution or the rinse solution.
 - Do not apply antibiotics to the valve.
1. After sizing, choose a valve of the appropriate size.
 2. Once the valve container has been removed from the outer packaging, examine the container for evidence of damage.


WARNING:

- The valve must not be implanted if the tamper-evident container seal is damaged, broken, or missing; or if fluid is leaking from the packaging.
 - The valve must not be implanted if the storage solution does not completely cover the valve.
3. Verify the valve size and expiration date on the label.
 4. To remove the valve from the container, break the seal and remove the screw-top closure.
CAUTION: Avoid prolonged contact with the formaldehyde storage solution. Immediately after contact, thoroughly flush any skin exposed to the solution with water. In case of contact with eyes, flush with water and seek appropriate medical care.
 5. Complete the medical device registration form and return it to Abbott Medical. Place one of the pull-off labels with the designated model and serial number in the patient's chart.

Removing the Valve from the Storage Container


1. Select the appropriate valve holder handle.
 - For aortic valves, select the aortic holder handle.
 - For mitral valves, select the mitral holder handle.
2. With the circulating nurse holding the container, press the valve holder handle into the valve holder, and remove the valve from the container.
 - For aortic valves, connect the aortic valve holder handle to the valve holder by pressing the aortic holder handle into the holder.

Figure 3. Press the aortic valve holder handle into the valve holder

- For mitral valves, connect the mitral valve holder handle to the valve holder by threading the mitral holder handle clockwise into the valve holder.

Figure 4. Thread the mitral valve holder handle clockwise into the valve holder

CAUTION: Do not use cutting edge needles, unprotected forceps, or sharp instruments as they may cause structural damage to the valve. Never handle the leaflet tissue.

1. Inspect the valve for damage.

WARNING: Do not implant the valve if it has been dropped, damaged, or mishandled in any way, or if there is any sign of deterioration.


Rinse Procedure

CAUTION: Do not implant the valve without thoroughly rinsing as directed.

1. Within the sterile field, prepare two sterile basins with a minimum of 500 mL of sterile isotonic saline solution in each basin.
2. Holding the valve by the handle, fully immerse the valve support, the valve, the valve holder, and the portion of the holder handle that was submerged in the valve storage solution, in the sterile isotonic saline solution in the first basin.
3. Continually rinse the valve for 10 seconds, using a gentle back-and-forth motion.
4. Repeat steps two and three in the second basin.
5. After rinsing, leave the valve immersed in the basin until required by the surgeon for implantation.
6. Prior to implantation, remove the valve support by depressing the three tabs below the level of the valve support ring.

CAUTION: Do not allow the valve tissue to dry. Place the valve in sterile isotonic saline rinse solution immediately upon removal from the valve storage solution.

Figure 5. Remove the valve support from the valve cuff

1. Support ring


Sizing Epic™ Plus Standard Aortic Valves

Epic Plus standard aortic valves are designed for intra-annular stent placement and supra-annular cuff placement.

To determine the correct standard aortic valve size, use the Model E2000 aortic sizers. The Model E2000 aortic sizer is a double-ended tool, with a supra valve replica end and an annular sizing end.

NOTE: Use only the annular sizing end of the E2000 sizer to size standard aortic valves.

Figure 6. Aortic sizer

Identify the sizer that fits snugly in the annulus and select the corresponding valve size.

WARNING: Standard aortic valve size selection is based on the size of the recipient annulus. Implantation of an inappropriately large bioprosthesis may result in stent deformation, valvular incompetence, and/or damage to the surrounding tissues. The use of an inappropriately small bioprosthesis may result in suboptimal hemodynamics. Use only the Epic™ Plus Heart Valve Sizer Model E2000 to size Epic Plus standard aortic valves.

CAUTION: Sizers are supplied non-sterile, and must be cleaned and sterilized prior to each use. Do not use cracked, deformed, or damaged sizer set components.

Sizing Epic™ Plus Supra Aortic Valves

Epic™ Plus Supra valves are designed for supra-annular aortic placement.


To determine the correct supra aortic valve size, use the Model E2000 aortic sizers. The aortic sizer is a double-ended tool, with a supra valve replica end and an annular sizing end. Use the annular sizing end to determine the size of the annulus. Insert the corresponding supra valve replica end in the supra-annular space to confirm placement and fit of the valve.

WARNING: Supra-annular aortic valve size selection is based on the size of the recipient annulus and the anatomy of the sinotubular space. Implantation of an inappropriately large bioprosthesis may result in stent deformation, valvular incompetence, and/or damage to the surrounding tissues. The use of an inappropriately small bioprosthesis may result in suboptimal hemodynamics. Use only the Epic™ Plus Heart Valve Sizer Model E2000 to size Epic™ Plus Supra aortic valves.

CAUTION:

- Sizers are supplied non-sterile, and must be cleaned and sterilized prior to each use. Do not use cracked, deformed, or damaged sizer set components.
- Do not pass the flanged portion of the supra valve replica sizing tool through the annulus.

Figure 7. Place the supra valve replica end in the supra-annular space

Sizing Epic™ Plus Mitral Valves


Epic™ Plus mitral valves are designed for intra-annular stent placement and supra-annular cuff placement.

To determine the correct standard mitral valve size, use the Model E2000 mitral sizers. Identify the sizer that fits snugly in the annulus and select the corresponding valve size.

WARNING: Mitral valve size selection is based on the size of the recipient annulus. Implantation of an inappropriately large bioprosthesis may result in stent deformation, valvular incompetence, and/or damage to the surrounding tissues. The use of an inappropriately small bioprosthesis may result in suboptimal hemodynamics. Use only the E2000 Epic™ Plus Heart Valve Sizer Set to size Epic Plus mitral valves.

CAUTION: Sizers are supplied non-sterile, and must be cleaned and sterilized prior to each use. Do not use cracked, deformed, or damaged sizer set components.

Figure 8. Mitral Sizer

Surgical Guidelines

The actual choice of surgical technique, modified in accordance with the instructions described herein, is left to the discretion of the individual surgeon. When implanting supra-annular valves, non-everting mattress sutures are recommended. Avoid any contact between the implantation sutures or knot tying technologies and the leaflets.

Precautions

- Do not allow the valve tissue to dry. Place the valve in isotonic sterile saline rinse solution immediately upon removal from the valve storage solution. Once removed from this solution, the valve should be periodically irrigated during implantation.
- Do not use cutting edge needles, unprotected forceps, or sharp instruments as they may cause structural damage to the valve.
- Never handle the leaflet tissue.
- Do not lacerate the valve tissue. If a valve is damaged, the valve must be explanted and replaced.
- Do not attempt to repair a valve. Damaged valves must not be used.

Aortic Valve Implantation


The Epic™ Plus aortic valve is designed to allow intra-annular placement of the inflow edge of the valve with supra-annular placement of the sewing cuff. The Epic™ Plus Supra valve is designed for supra-annular implantation of both the valve and the sewing cuff.

1. After sizing, choose a valve of the appropriate size.

CAUTION: Position the valve so that the stent posts do not obstruct the coronary ostia.

2. To facilitate implantation, the valve holder handle may be removed from the valve holder by depressing the release button on the valve holder.


Figure 9. Release the valve holder handle

1. Release button

3. To remove the valve holder from the valve, cut the three retaining sutures and pull the valve holder away from the valve. Examine the valve to ensure that there are no holder suture remnants.

Figure 10. Cut three sutures to remove the aortic valve

Mitral Valve Implantation


1. After sizing, choose a valve of the appropriate size.

CAUTION: Position the mitral valve in a manner to avoid commissure obstruction of the left ventricular outflow tract, and minimize any potential of commissure contact with the ventricular wall.

2. To facilitate insertion of the mitral valve into the annulus, the mitral valve stent posts may be temporarily deflected inward during implantation. To deflect the valve stent posts inward, rotate the valve holder handle in the clockwise direction while securing the valve holder. The valve stent posts will remain deflected until the valve retaining sutures are cut. To facilitate implantation, the valve holder handle may be removed from the valve holder by depressing the release button on the valve holder.


NOTE: Take care to avoid looping or entangling sutures around the commissural posts, as this may result in a compromise of leaflet function.

Figure 11. Rotate the valve holder handle to deflect the mitral stent posts inward

3. To remove the valve holder from the valve, make a single cut to the retaining sutures, and pull the valve holder away from the valve. Examine the valve to ensure that there are no holder suture remnants.

Figure 12. Cut the sutures to remove the mitral holder

Intra-Operative Assessment

The suggested method for assessing competence of the Epic™ Plus and the Epic™ Plus Supra valve is with intra-operative Doppler echocardiography.

Future Valve-in-Valve Considerations

A transcatheter heart valve implanted in a surgical aortic or mitral tissue valve (ViV) is emerging as a reasonable treatment option for severely symptomatic patients with a failed surgical valve who are judged by the heart team to be at high or prohibitive risk of surgical reoperation, and in whom improvement in hemodynamics is anticipated. The information in this section is intended to be helpful for future ViV procedures in Epic™ Plus valves. Reference the Catalog Number Descriptions and Reference Dimensions Table (page 1) for valve and stent properties.

CAUTION: When performing a ViV procedure, the instructions for use of the transcatheter heart valve must be followed. The safety and efficacy of ViV procedures in an Epic™ Plus valve has not been established.

Aortic and mitral Epic™ Plus valves contain radiopaque wire markers at the annulus and stent posts to aid the interventional cardiologist when placing the transcatheter prosthesis within the Epic Plus valve. Physicians may use the annulus and/or stent post markers as a reference point, in conjunction with known transcatheter stent foreshortening properties, to determine the proper location to place the transcatheter valve prior to delivery system actuation. Epic Plus markers can continue to be utilized while adjusting the transcatheter valve during the deployment process.

Surgical Valve Implantation in the Aortic Position

When implanting the Epic™ Plus or Epic™ Plus Supra valves in the aortic position, it is important to:

- Ensure that there is adequate clearance of the coronary artery openings
- Assess the suitability of the implanted surgical valve and anatomy for a future ViV intervention

The replica end of the sizer accurately represents the valve and should be used to assess the position of the valve relative to the encountered anatomy and thereby ensure proper coronary opening clearance and assess suitability of the valve for a future ViV intervention.

Surgical Valve Implantation in the Mitral Position

When implanting the Epic™ Plus valve in the mitral position, it is important to:

- Ensure that there is no obstruction of the left ventricular outflow tract and that there is no direct contact between the commissural posts and the ventricular wall
- Assess the suitability of the implanted surgical valve and anatomy for a future ViV intervention

Three-dimensional CT imaging is the recommended measurement method to assess a failed surgical valve dimensions when selecting the appropriate transcatheter heart valve size during a ViV procedure. The reference dimensions in the Catalog Number Descriptions and Reference Dimensions Table (page 1) do not account for valve tissue thickening, calcium burden/distribution, or other factors that vary by patient and are important when selecting a transcatheter heart valve size.

Patient Identification Card

A Patient Identification Card is present within the packaging of the valve. Hospital personnel should complete the front side of the Patient Identification Card by printing the patient's name, implant date, along with the physician's name and healthcare facility information and phone number in the designated areas. A peel-away device label is provided in the sterile packaging which includes device information (model number, serial number, UDI number, and device name). Apply the peel-away device label to the back of the card in the designated area. If a peel-away label is not available, write the information in the designated area on the back of the card. In procedures where multiple implantable devices are implanted, include information for all implants in the designated areas.

NOTE: Write the valve serial number in the SN/LOT field.

After the card is completed, separate the top portion from the Patient Identification Card by tearing along the perforated line. Provide the completed Patient Identification Card to the patient. Instruct the patient to carry the Patient Identification Card at all times and to present the implant card to hospital staff in emergency situations. To obtain a replacement card if a patient loses or damages their card, contact Abbott Medical Technical Support. Once the replacement card is received, handwrite the required information in the respective areas and provide to the patient.

Product Materials

The following product materials are intended to come into contact with tissue.

Table 6. Aortic Valves

Component	Material Description	Substance	Substance Concentration within each component (% w/w)(maximum)	Component Weight (g)
Cuff		Polyester	99.9 - 100	0.1544 - 0.2332

Table 6. Aortic Valves

Component	Material Description	Substance	Substance Concentration within each component (% w/w)(maximum)	Component Weight (g)
Stent	Teflon (PTFE) impregnated braided polyester (polyethylene terephthalate - PET)	Acetal Copolymer	99.9 - 100	0.3640 - 0.7120
Cuff Filler		Polyester	99.9 - 100	0.0524 - 0.0648
Suture				
Wire Wrap	Stainless Steel	Carbon	0.030	0.1052 - 0.1280
		Manganese	2.00	
		Phosphorous	0.025	
		Sulfur	0.010	
		Silicon	0.75	
		Chromium	17.00 - 19.00	
		Nickel	13.00 - 15.00	
		Molybdenum	2.25 - 3.00	
		Nitrogen	0.10	
		Copper	0.50	
		Cobalt	0.099	
		Iron	Balance	
Porcine Cusp	Porcine Cusp		99.9 - 100	1.7116 - 2.9054
Bovine Pericardium	Bovine Pericardium		99.9 - 100	0.6180 -1.0472
Chemical Residue on Valve After Physician's Rinse				
Formaldehyde				0.27555 (% w/w)

Table 7. Mitral Valves

Component	Material Description	Substance	Substance Concentration within each component (% w/w)(maximum)	Component Weight (g)
Cuff	Teflon (PTFE) impregnated braided polyester (polyethylene terephthalate - PET)	Polyester	99.9 - 100	0.2200 - 0.3228
Stent		Acetal Copolymer	99.9 - 100	0.5000 - 0.8748
Cuff Filler		Polyester	99.9 - 100	0.1368 - 0.1712
Suture			99.9 - 100	0.0812 - 0.0908
Wire Wrap	Stainless Steel	Carbon	0.030	0.1228 - 0.1456
		Manganese	2.00	
		Phosphorous	0.025	
		Sulfur	0.010	
		Silicon	0.75	
		Chromium	17.00 - 19.00	
		Nickel	13.00 - 15.00	
		Molybdenum	2.25 - 3.00	
		Nitrogen	0.10	
		Copper	0.50	
		Cobalt	0.099	
		Iron	Balance	
Porcine Cusp	Porcine Cusp		99.9 - 100	1.8076 - 3.4058
Bovine Pericardium	Bovine Pericardium		99.9 - 100	0.6180 - 0.9520
Chemical Residue on Valve After Physician's Rinse				
Formaldehyde				0.27555 (% w/w)

Table 8. Supra Valves

Component	Material Description	Substance	Substance Concentration within each component (% w/w)(maximum)	Component Weight (g)
Cuff	Polyester		99.9 - 100	0.1480 - 0.2508

Table 8. Supra Valves

Component	Material Description	Substance	Substance Concentration within each component (% w/w)(maximum)	Component Weight (g)
Stent		Acetal Copolymer	99.9 - 100	0.1348 - 0.2044
Cuff Filler		Silicone	99.9 - 100	0.3640 - 0.8456
Suture (green)		Polyester	99.9 - 100	0.0088 - 0.0096
Suture	Teflon (PTFE) impregnated braided polyester (polyethylene terephthalate - PET)		99.9 - 100	0.0816 - 0.0836
Wire Wrap	Stainless Steel	Carbon	0.030	0.1020 - 0.1444
		Manganese	2.00	
		Phosphorous	0.025	
		Sulfur	0.010	
		Silicon	0.75	
		Chromium	17.00 - 19.00	
		Nickel	13.00 - 15.00	
		Molybdenum	2.25 - 3.00	
		Nitrogen	0.10	
		Copper	0.50	
		Cobalt	0.099	
		Iron	Balance	
Porcine Cusp		Porcine Cusp	99.9 - 100	1.7238 - 2.7958
Bovine Pericardium		Bovine Pericardium	99.9 - 100	1.0472 - 0.9040
Chemical Residue on Valve After Physician's Rinse				
Formaldehyde				0.27555 (% w/w)

CAUTION: The Epic™ Plus valves contain Formaldehyde (Chemical Abstracts Service CAS No. 50-00-0; EC No. 200-001-8) which is defined as a CMR 1B (Category 1B carcinogen) in a concentration above 0.1% weight by weight.

CAUTION: One or more components of the Epic™ Plus valves contains the following substance defined as CMR 1B in a concentration above 0.1% weight by weight:

Cobalt; Chemical Abstracts Service CAS No.7440-48-4; EC No. 231-158-0

Current scientific evidence supports that medical devices manufactured from metal alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

Incident Reporting

If, in the course of use of this device, you have reason to believe that a serious incident occurred, please report it to the manufacturer. For customers in the European Union, report the serious incident to your national authority as well as to the manufacturer.

Disposal

This instructions for use is recyclable. Dispose of all packaging materials as appropriate. Dispose of valves and accessories per standard solid biohazard waste procedures. Dispose of the Formaldehyde storage solution following hospital standards/protocols for disposal of hazardous substances.

Limited Warranty


Abbott Medical warrants that reasonable care has been used in the manufacturing of this device. THIS WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES NOT EXPRESSLY SET FORTH HEREIN, WHETHER EXPRESS OR IMPLIED BY OPERATION OF LAW OR OTHERWISE INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, since handling, storage, cleaning, and sterilization of this device as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Abbott's control directly affect this device and the results obtained from its use. ABBOTT SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE directly or indirectly arising from the use of this device other than the replacement of all or part of it. Abbott neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

Some states in the United States do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you. This limited warranty gives you specific legal rights, and you may have other rights which vary from jurisdiction to jurisdiction.

Descriptions of specifications, appearing in Abbott literature, are meant solely to generally describe the device at the time of manufacture and do not constitute any express warranties.

Symbols

The symbols below and harmonized symbols may be found on the product or product label. For harmonized symbols, refer to the Universal Symbols Glossary at medical.abbott/manuals.

Symbol	Definition
	Stented Tissue Valve
	Rinse - 2 x 500 mL x 10 seconds
	Contents
	Aortic
	Mitral

	Temperature Indicator: If red, do not use
STERILE LC	Sterilized by liquid chemical sterilant
UDI	Unique Device Identifier
	Follow instructions for use on this website
medical.abbott/manuals	
	Contains hazardous substances (Formaldehyde CAS 50-00-0)
(50-00-0)	
	Contains hazardous substances (Cobalt CAS 7440-48-4)
(7440-48-4)	
	Date
	Healthcare center or physician
	Patient identification
	Physician telephone
	Ukrainian mark of conformity to the technical regulations, where 101 is the assessment body.
UA.TR.101	
	Contains biological material of animal origin
MD	Medical Device
	European conformity, affixed in accordance with the relevant provisions of European Council Regulation 2017/745 (NB 2797). Hereby, Abbott Medical declares that this device is in compliance with the relevant provisions of this regulation.
2797	
	Patient information website
medical.abbott/manuals	
	United Kingdom Conformity Assessed
0086	
UKRP	Responsible person in the United Kingdom
	Importer
ID-EFFECTIVE	Effective Orifice Diameter
ID-INFLOW	Inflow Orifice Diameter
OD-SEWINGRING	External sewing ring diameter
OD-HOUSING	Valve housing external diameter

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A technical line drawing of a mechanical assembly. It features a central vertical shaft passing through a multi-tiered base. The base has a circular, ribbed bottom flange. Four curved, ribbed components are arranged radially around the central shaft, each secured with a bolt. A callout line labeled '1' points to a small, rectangular component mounted on the right side of the base.

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