Constic

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Instructions for Use English

Product description

Constic is a self-etching, self-adhesive, radiopaque, flowable composite which is light-cured. The composite is immediately ready for use because the preparatory steps of etching, priming and bonding are contained in Constic.

Indications

- Small restorations of class I + III
- Restorations of class V
- Cavity liner for class I + II
- Fissure sealing
- Extended fissure sealing
- Repairs of composite restorations
- Modifications to temporaries and long-term temporaries
- Blocking out and filling of undercuts in cavities
- Small occlusal primary tooth cavities

Contra-indications

- Do not use Constic in the case of allergies to any of the components or in the event of contact allergies.
- Application of the material is contraindicated if dry isolation or the recommended application technique is not possible.
- Do not apply Constic to the exposed pulp.
- Do not use Constic on etched dentine.

Basic Safety Instructions

- For dental use only!
- Keep out of reach of children!
- Avoid contact with the skin! In the event of inadvertent skin contact, wash the affected area thoroughly with soap and water. Wear suitable gloves.
- Avoid contact with the eyes! In the event of inadvertent contact with the eyes, rinse thoroughly with plenty of water immediately and consult a physician if necessary.

Side effects

Isolated cases of hypersensitivity or contact allergy to components of the material cannot be excluded.

Interactions with other materials

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the point of contact. Avoid materials containing eugenol, moisture and air containing oil!
- The use of cationic mouthwashes, plaque indicators and chlorhexidine may lead to discoloration.

Instructions for use

- Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm². Observe the manufacturer's instructions for the light unit.
- Position the curing light as close to the material as possible.
- Do not remove the oxygen-inhibited layer because this is required for the bond with the next layer.
- When using Constic in conjunction with other materials, follow the relevant manufacturer's instructions.

Recommended Use

Minor restorations of class I, III, V and underfillings of class I and II

1. Determine tooth color and select the appropriate shade.
2. Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges if necessary. Protect pulp in deep cavities with appropriate preparations.

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3. Clean tooth, remove all cleaning residue with water spray.

Note: Cleaning with alcohol may cause overdrying of the dentine.

4. Dry the tooth using water- and oil-free air in order to avoid overdrying of the dentine. A moist layer must remain on the surface of the tooth.
5. Apply the material onto the cavity surface with the aid of the Luer-Lock-Tip by pressing the syringe and massage a thin layer (≈ 0.5 mm) into the entire surface of the cavity wall for 25 s using the brush.
6. Remove any excess, if necessary.
7. Light-cure for 20 s.
8. Layer the material as follows:
 - Small restorations of class I, III, V: Use the material to build up layers with a maximum 2 mm layer thickness and cure each layer for 20 s.
 - Underfillings of classes I and II: Use the material to build up 1 layer with a maximum 1 mm layer thickness and light-cure for 20 s.

Note: It is not necessary to repeat massaging of the individual layers.

9. Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru tour the filling surface and finish using suitable instruments.
10. Use polishing paste to polish the filling surface and neighboring surfaces. Use a polishing cup or polishing disk for polishing work.

Fissure sealing and extended fissure sealing

1. For extended fissure sealing, widen exceedingly narrow and ampullate fissures.
2. Clean enamel with fluoride-free paste, remove all cleaning residue with water spray and dry the tooth with water- and oil-free air. In doing so, leave a moist surface on the enamel.
3. Determine tooth color and select the appropriate shade.

Note: DMG recommends enamel etching in addition to cleaning. To do this, etch enamel for 15–30 s with 37% phosphoric acid gel, rinse for 15 s with water and dry with oil- and water-free air. In doing so, leave a moist surface on the enamel.

4. Apply the material onto the surface of the fissure or cavity wall with the aid of the Luer-Lock-Tip by pressing the syringe and massage a thin layer (≈ 0.5 mm) into the entire surface for 25 s using the brush.
5. Remove any excess, if necessary.
6. Light-cure for 20 s.
7. Where necessary, use the material to build up layers with a maximum 2 mm layer thickness and light-cure for 20 s.

Note: It is not necessary to repeat massaging of the individual layers.

8. Check occlusion and correct if necessary.

Disinfection

The brush holder is reusable after disinfection, although it should be replaced if visible signs of wear, e.g. scratches, appear. Disinfection can be carried out using a commercially available wiping or dipping disinfectant. DMG recommends using only disinfecting solutions listed by the RKI (Robert Koch Institute).

Composition

Barium glass in a Bis-GMA-based matrix from dental resins. Pigments, additives and catalyst. Filler content: 66 wt.-% = 43 vol.-%. The range of variation of the inorganic filler particles is between 0.02 and 2.3 μ m.

Storage and shelf life

- Store at room temperature (15 to 25 °C/59 to 77 °F)!
- Reseal syringes after use!
- Do not use after the expiration date!

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English

Product description

Ecosite Bulk Fill is a light-curing, radiopaque restorative material developed especially for use in the posterior region.

Indications

- Restorations in the posterior region (Classes I and II, including replacement of individual cusps)
- Restorations of class V (cervical caries, root erosions, wedge-shaped defects)
- Core build-up
- Extended fissure sealing on molars and premolars
- Restorations in deciduous teeth

Contraindications

- Application of the material is contraindicated if dry isolation or the recommended application technique is not possible.
- Do not use directly on exposed pulp or in cases of allergies to any of the ingredients.

Interactions

- Avoid materials containing eugenol, moisture and oily air, as these may inhibit polymerization at the contact area.
- The use of cationic oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

Side effects

There are no known side effects to date. Hypersensitivity to components of the material cannot be ruled out in individual cases.

Basic safety information

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.

Notes for use

- For optimum polymerization, position the light-curing unit as close to the restorative material as possible and observe the recommendations relating to the light-curing time and the light intensity (see table). Observe the manufacturer’s instructions for the light-curing unit used.
- Follow the manufacturer’s instructions for other materials that are used with Ecosite Bulk Fill.
- If Ecosite Bulk Fill is applied into the patient’s mouth straight from the Safetip, the Safetip must only be used on one patient for hygiene reasons.

Recommended use

- Clean the tooth. Remove all cleaning residue with water spray. It is advisable to use a rubber dam.
- Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Ensure that the areas close to the pulp are protected.
- Condition the dental enamel and apply the adhesion primer according to the manufacturer’s instructions for the material used. DMG advises using the universal adhesion primer Ecosite Bond.
- If necessary, apply an underfilling of low-viscosity composite or compomer. We recommend using a matrix.
- Apply Ecosite Bulk Fill in the cavity with a maximum layer thickness of 5 mm. Apply taller restorations one layer at a time. Light-cure each layer separately.

Light intensity	> 1,000 mW/cm²	> 600 mW/cm²
Curing time	20 s	40 s

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- Removal of excess material, contouring and fine processing are done using diamond finishing tools and flexible polishing discs of decreasing grain size, for example. Ecosite Bulk Fill can be given a high shine using polishing paste together with polishing cups or discs.

Disinfection of the Safetip-Dispenser

The Safetip-Dispenser is reusable after disinfection, although it should be replaced if visible signs of wear appear.

Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute). The Safetip-Dispenser is also autoclavable.

Composition

Barium glass in a Bis-GMA-based matrix of dental resins. Pigments, additives and catalyst. Filler content: 82% by weight = 65% by volume The range of variation of the inorganic filler particles is between 0.02 and 0.7 µm

Storage and shelf life

- Store at room temperature (2 to 25°C / 36 to 77°F).
- Close syringes / composite tip immediately after use so that the material does not cure in the ambient light.
- Do not use after the expiration date.

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Starter Kit:	
1 syringe containing 4 g of Light 2 syringes each containing 4 g of Universal 1 bottle @ 5 ml Ecosite Bond	REF 220750
16 Safetips each containing 0,25 g of Light 32 Safetips each containing 0,25 g of Universal 1 bottle @ 5 ml Ecosite Bond	REF 220751

Packaging: 1 syringe containing 4 g paste	
Light	REF 220752
Universal	REF 220753
Contrast	REF 220754
Packaging: 16 Safetips each containing 0.25 g paste	
Light	REF 220755
Universal	REF 220756
Contrast	REF 220757

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English

Product description

Ecosite Elements is a light-curing, esthetic, radiopaque nanohybrid composite system for the anterior and posterior region. The range of highly viscous composites includes 5 VITA shades (B1, A2, A3, A3.5 and A4) and 4 enamel shades (Enamel Bleach = EB, Enamel Light = EL, Enamel Medium = EM and Enamel Dark = ED) in different opacities. This, combined with their unique handling and excellent polishing properties contributes to the creation of highly-aesthetic restorations.

The shades OA2 (Opak A2) and INC (Incisal) are radiopaque composites with optimum flow composites. The shade OA2 is recommended for the veneering of slightly discolored areas and for creating ultimate tooth opacity. The shade INC is particularly suitable for an esthetic incisal edge design thanks to its high transparency. Both shades are suitable for fissure sealing.

The shades W (White) and B (Brown) are flowable composites for the customization of composite restorations. The shade W is used for the reproduction of white spots and as an opaque covering layer. The shade B is an ideal choice for fissures and the replication of brown spots.

Intended purpose

Ecosite Elements in the colors B1, A2, A3, A3.5, A4, EB, EL, EM, ED

- Restorations of all cavity classes of every dentition
- Direct inlays, onlays and veneers (vestibular and occlusal)
- Core build-up
- Splinting/blocking of loose teeth
- Repairs

Ecosite Elements in the colors OA2 and INC

- Restorations of all cavity classes III-V of every dentition
- Liner
- Fissure sealing
- Splinting/blocking of loose teeth
- Repairs

Ecosite Elements in the colors W and B:

- Customization of composite restorations

Limitations of use

Application of the material is contraindicated if dry isolation or the recommended application technique is not possible.

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Indications

- Loss of tooth substance due to caries, fractures, tooth abrasion or development disorders
- Diastema, tooth discoloration or slight misalignment of the anterior teeth
- Dental or periodontal diseases, e.g. trauma and periodontitis, which require the splinting for the affected teeth
- Primary and secondary prevention of caries

Contraindications

- Do not use directly on the exposed pulp.
- Do not use the material where there are known allergies to any of the components or contact allergies.

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Dentist

Notes for use

- For optimum polymerization, position the light-curing unit as close to the restorative material as possible and observe the recommendations relating to the light-curing time and the light intensity (see table). Observe the manufacturer's instructions for the light-curing unit used.
- Follow the manufacturer's instructions for other materials that are used with Ecosite Elements.
- If the material is applied into the patient's mouth with a single use applicator, the single use applicator must only be used on one patient for hygiene reasons.

Recommended use

(See the end of the instructions for use)

- Clean the tooth. Remove all cleaning residue with water spray. Determine tooth shade with the Ecosite Shade Guide Pure, the Ecosite Shade Guide Layer or the VITA tooth shade. Select suitable Ecosite Elements shade or shade combinations. Either choose a monochromatic Pure shade or combine it with a layer shade according to the outline. The specification corresponds to the use of a 0.5 mm-thick Ecosite Elements Layer:

Note: DMG recommends the use of a rubber dam.

- Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Ensure that the areas close to the pulp are protected.
- Condition the dental enamel and apply the adhesion primer according to the manufacturer's instructions for the material used. DMG advises using the universal adhesion primer Ecosite Bond.
- If necessary apply an underfilling of low-viscosity composite or compomer. We recommend using a matrix.
- Insert Ecosite Elements with a maximum layer thickness of 2 mm into the cavity. Apply taller restorations one layer at a time. Light-cure each layer separately.

Light intensity	> 1,000 mW/cm²	> 600 mW/cm²
Radiation time	20 s	40 s
Max. layer thickness	2 mm	2 mm
Wavelength	420 – 480 nm	420 – 480 nm

Note: Additional buccal or lingual light-curing is advisable. Remove metal matrices beforehand.

6. Removal of excess material, contouring and fine processing are done using diamond finishing tools and flexible polishing discs of decreasing grain size, for example. Ecosite Elements can be given a high shine using polishing paste together with polishing cups or discs.

Customization of composite restorations with W (White) and B (Brown)

1. Select the color required for the desired color shade.
2. Pour a sufficient amount of the selected shade onto a mixing pallet.
3. Apply the material to the composite with a fine brush in a thickness of ≤ 0.1 mm. Light-cure each layer separately.

Note: Always cover the customization layer with a composite layer.

Note: Mix the customization shades with one another (to obtain new tones), but never with other products.

Light intensity	> 1,000 mW/cm²	> 600 mW/cm²
Radiation time	20 s	40 s
Max. layer thickness	0.1 mm	0.1 mm
Wavelength	420 – 480 nm	420 – 480 nm

Disinfection of the Safetip-Dispenser Pro

The Safetip-Dispenser Pro is reusable after disinfection, although it should be replaced if visible signs of wear appear.
The Safetip-Dispenser Pro is suitable for steam sterilization in autoclaves. The procedure must be carried out at 134°C/3 bar for 5 minutes.

Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid this.
- The use of oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

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Warnings / precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Reuse: Standard surface disinfection as required. No specific maintenance required.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

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Ecosite Elements in the shades B1, A2, A3, A3.5, A4, EB, EL, EM, ED complies with DIN EN ISO 4049: 2009 (Type 1, Class 2, Group 1); polymerization depth: ≤ 2 mm with an light-curing time of 40 seconds with a light intensity of ≥ 600 mW/cm² and a light-curing time of 20 seconds with a light intensity of ≥ 1000 mW/cm²
Ecosite Elements in shades OA2 and INC correspond to DIN EN ISO 4049:2009 (Type 2, Class 2, Group 1); polymerization depth ≥ 2 mm with a light-curing time of 40 seconds with a light intensity of ≥ 600 mW/cm² and a light-curing time of 20 seconds with a light intensity of ≥ 1000 mW/cm²

Storage / disposal

- Store in a dry place (2-25°C/36-77°F) at room temperature.
- Dispose of in accordance with national regulations.

Composition

Ecosite Elements B1, A2, A3, A3.5, A4, EB, EL, EM, ED:
(Barium) dental glass (approx. 81 %) in a resin matrix based on EBPADMA and Bis-GMA (approx. 19%). Pigments, additives and catalysts.
Filler content: approx. 81 wt. % = 65 vol.%. The inorganic filler particle size varies between 0.02 and 0.7 µm.

Ecosite Elements OA2, INC, W and B:
(Barium) dental glass (approx. 65 %) in a resin matrix based on EBPADMA, Bis-GMA and TEDMA (about 34%)... Pigments, additives and catalysts.
Filler content: approx. 65. wt% = 38 vol.%. The inorganic filler particle size varies between 0.02 and 0.7 µm.

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Product description

Ionosit-Baseliner is a light-curing, radiopaque, compomer base liner. The controlled expansion of Ionosit-Baseliner counteracts the shrinkage of a composite filling.

Indications

- Base liner for composite restorations
- Blocking out of undercuts

Basic safety instructions

- For dental use only!
- Keep out of reach of children!

Side effects

There have been no systemic side effects to date. Isolated cases of hypersensitivity or contact allergy to components of the material cannot be excluded.

Instructions for use

- Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm². Place the light as close as possible to the material.
- Ionosit-Baseliner adheres to metal instruments. Therefore, clean instruments with a paper towel before the material sets.
- Ionosit-Baseliner is light-sensitive. Protect the material from premature activation by ambient light and close the syringe immediately after use.
- Do not remove the oxygen inhibited layer because it is required for bonding to the next layer.
- The use of a matrix is recommended. Follow the manufacturers’ instructions for other materials used.

Recommended use

1. Before preparation, clean the tooth and remove all cleaning residue with water spray. Use of a rubber dam is not necessary if moisture can be controlled in another way. Avoid contamination with saliva during the application.
2. Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Protect pulp in deep cavities with suitable materials.
3. Etch the cavity in accordance with the preferred technique (e.g. Total-Etch). Observe the manufacturer’s instructions for the material used.
4. Apply bonding agent (e.g. TECO) in accordance with the manufacturer’s instructions.
5. Ionosit-Baseliner can be used both as a conventional base liner (coating the floor of the cavity) and to coat the entire cavity (stress-breaking function). Apply Ionosit-Baseliner with a brush or ball-pointed instrument in a thickness of about 1 mm.
6. Cure for at least 20 seconds with a suitable light unit.
7. Place the composite filling material according to the manufacturer’s instructions. DMG recommends the use of EcuSphere.

Storage and shelf life

- Store in a dry place at room temperature (15 to 25 °C/59 to 77 °F)!
- Do not use after the expiration date.

Composition

Ionomer glass in a matrix of polymerizable oligo- and polycarbonic acids and other light-cure dental resins. The ionomer glass used contains fluoride and zinc ions. Filler content: 72 % by weight = 55 % by vol. (0.02–33 µm).

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Instructions for use English

Product description

LuxaBond-Total Etch is a dual-cure bonding agent for enamel and dentin, suitable for use with dual-cure and self-cure composites, such as LuxaCore Z-Dual. It is necessary to etch enamel and dentin surfaces prior to application. Due to the dual-cure nature of the material, LuxaBond-Total Etch is particularly useful when bonding posts in areas where light can't penetrate.

Indications

Dual-cure bonding agent used in conjunction with the total-etch technique for core build-ups, and adhesion of resin cements for posts, crowns, bridges, inlays and onlays.

Contraindications

Do not apply the material in cases of allergies to any of the components, in particular benzoyl peroxide, or in the event of contact allergies.

Basic safety instructions

- For dental use only!
- Keep out of reach of children!
- Avoid contact with eyes, skin and oral mucosa. In the event of inadvertent contact with eyes, immediately rinse the eyes thoroughly with plenty of water and consult a physician if necessary. In the event of inadvertent skin contact, wash the affected area thoroughly with soap and water.

Side effects

- Isolated cases of contact allergies to components of the material cannot be excluded.

Interactions

Products containing eugenol, for example temporary cements or temporary root canal filling materials, may inhibit the polymerization.

- Products containing eugenol and other substances that inhibit polymerization must not be used together with the material.

Timing

Working time of the mixture of bond A + bond B (storage in a dark place)	5 min
Reaction time of the etching gel to enamel	20 to 60 s
Reaction time of the etching gel to dentin	max. 15 s
Pre-bond work-in time	15 s
Mixing time of bond A + bond B	5 s
Bond A + bond B work-in time	20 s
LuxaBond-Total Etch exposure time	min. 10 s

Instructions for use

- Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm². Place the light as close as possible to the material.
- Always use etching gel on enamel and dentin surfaces.
- Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

 premature polymerization due to ambient light, the material must be added to the mixing palette immediately prior to application.
- Ensure that the etched or bonded surface does not become contaminated prior to each subsequent working step.
- Do not remove the oxygen inhibition layer caused by the application because this is important for the bonding process with the material applied on top.
- Observe the manufacturer's instructions for all materials that are used with LuxaBond-Total Etch.

Recommended use

Use in the root canal, for cementing root posts, for cementing inlays

1. Prepare the root canal or cavity in accordance with the general rules for adhesive procedures. We recommend the use of a rubber dam.
2. Etch dentin and enamel with 37 % phosphoric acid gel or similar. Start the etching process with the enamel and then expand this to include the dentin. The reaction time on the enamel should be 20 to 60 s. The reaction time on the dentin must not exceed 15 s.

Please note: When using LuxaBond-Total Etch it is essential to etch the affected enamel and dentin surfaces.

3. Rinse with water for a minimum of 15 seconds.
4. Remove general water excess. After the rinsing process, an even moist layer must remain on the dentin surface requiring treatment. If necessary, wet the dentin surface using a lightly moistened micro brush.
5. Use a brush to work 1 to 2 drops of pre-bond into the etched tooth structure for 15 s. Remove any excess, e.g. with a paper tip. Blow gently using oil and water-free air. The surface should appear evenly moist.
6. Mix 1 to 2 drops of bond A and bond B in a trough in the mixing palette at a ratio of 1:1 for approx. 5 s and work into the tooth structure for 20 s using a micro brush.
7. Carefully remove bond excess from the preparation, e.g. with a paper tip, and then gently blow with oil and water-free air until the surface appears evenly moist.
8. Apply the attaching cement in accordance with the manufacturer's instructions.

Please note: The working time of the attaching cement may reduce due to the contact with LuxaBond-Total Etch.

Please note: The use of LuxaCore Z-Dual is recommended for the cementation of root posts.

Please note: The use of Vitique is recommended for the cementation of inlays.

Core build-ups, cementation of crowns and bridges

Please note: Due to the fact that excess material may collect between the tooth and matrix (pooling effect), if possible use a matrix only after the bond has been applied.

1. Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. We recommend the use of a rubber dam. In deep cavitated areas located close to the pulp, provide protection in the form of a suitable liner, e.g. by using a small quantity of calcium hydroxide and a thin layer of glass ionomer cement.
2. Etch dentin and enamel with 37 % phosphoric acid gel. Start the etching process with the enamel and then expand this to include the dentin. The reaction time on the enamel should be 20 to 60 s. The reaction time on the dentin must not exceed 15 s.

Please note: When using LuxaBond-Total Etch, it is essential to etch the affected enamel and dentin surfaces.

3. Rinse with water for a minimum of 15 s.
4. Remove general water excess. After the rinsing process, an even moist layer must remain on the dentine surface requiring treatment. If necessary, wet the dentin surface using a lightly moistened micro brush.
5. Use a brush to work 1 to 2 drops of pre-bond into the etched tooth structure for 15 s. Remove any excess, e.g. with a paper tip. Blow gently using oil and water-free air. The surface should appear evenly moist.
6. Mix 1 to 2 drops of bond A and bond B in a trough in the mixing palette at a ratio of 1:1 for approx. 5 s and work into the tooth structure for 20 seconds using a micro brush.
7. Carefully remove any bond excess and then gently blow with oil and water-free air until the surface appears evenly moist.
8. If required, light-cure LuxaBond-Total Etch for a minimum of 10 s using a suitable light curing unit.
9. Apply core build-up material (e.g. LuxaCore Z-Dual) or resin cements (e.g. Vitique) in accordance with the manufacturer's instructions.

Please note: The working time of the resin cements or the core build-up material may reduce due to the contact with LuxaBond-Total Etch.

Please note: The use of Vitique is recommended for the cementation of onlays and crowns.

Composition

- Pre-bond: Ethanolic aryl sulfonate solution
- Bond A: Hydrophilic Bis-GMA based resin matrix, catalyst
- Bond B: Hydrophilic Bis-GMA based resin matrix, benzoyl peroxide

Storage

- Store LuxaBond-Total Etch in a refrigerator at 2–8 °C (35–46 °F).
- Close bottles immediately after use.
- Do not use the material past the stipulated shelf life or if the fluid assumes a gelatinous consistency.
- Ensure that the lids of bond bottles are not switched in order to avoid any cross contamination.

Available product sizes

Intro Kit	
1 bottle @ 5 ml Pre-Bond, 1 bottle @ 5 ml Bond A, 1 bottle @ 5 ml bond B, 1 syringe @ 2 ml DMG Etching Gel, 25 Endobrushes Red, 25 Endobrushes Black, 3 Luer-Lock-Tips, 1 mixing palette	REF 212384

Accessories	
Endobrushes Red	REF 212385
Endobrushes Black	REF 212386
Luer-Lock-Tips White	REF 109942
DMG Etching Gel	REF 212873
LuxaBond Pre-Bond Refill	REF 212389

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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English

Product description

LuxaBond Universal is a dual-curing 2-bottle bonding agent. The material can be used as an adhesive on enamel and dentine with the self-etch, total-etch or selective-etch technique in combination with dual-curing, self-curing or light-curing composites (e.g. PermaCem Universal, LuxaCore Z Dual, Ecosite Elements). Thanks to its dual-curing property, the material can be used in the root canal area where it is difficult for light to penetrate.

The material can also be used as a primer for indirect restorations made of metal (base metal, precious metal), oxide and silicate ceramics or composites.

Intended purpose

As an adhesive for:

- Light-curing fissure-sealing materials
- Direct restorations made of light-curing composites and compomers
- Light-curing, self-curing and dual-curing core build-up materials
- Cementing of indirect restorations (inlays, onlays, crowns and bridges) and root posts with light-curing, self-curing and dual-curing composite-based luting cements
- Cementing of veneers with light-curing luting composite cements in conjunction with the enamel etching technique (e.g. with Vitique Cement)

As a primer for indirect restorations made from:

- Metals (base metal, precious metal)
- Zircon/aluminum oxide
- Silicate ceramic
- Composite

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Do not use the material in a sufficiently dry working area or the recommended application technique is not possible.

Indications

- Loss of tooth substance/tooth due to caries, trauma, tooth abrasion or development disorders
- Diastema, tooth discoloration or slight misalignment of the anterior teeth
- Trauma or periodontitis that require composite splints
- Molars and premolars with pronounced deep fissures and increased caries risk

Contraindications

- Do not use directly on the exposed pulp.
- Do not use the material if pulpitis is present.
- Do not use the material in cases of known allergies to any of the components, in particular benzoyl peroxide, or in the event of contact allergies.

Patient target group

People treated in the course of a dental procedure.

Intended users

Dentist

Notes for use

- Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be at least 600 mW/cm². Place the light-curing unit as close as possible to the material. Observe the manufacturer's instructions.
- Bring the material to room temperature (15 – 25 °C/59 – 77 °F) before use.
- In order to prevent premature polymerization due to ambient light, the material should be placed onto the mixing pad immediately prior to the application.
- Ensure that the etched or bonded surface does not become contaminated prior to each subsequent working step.
- Do not remove the oxygen inhibition layer caused by the application because this is important for the bonding process with the material applied on top.
- In deep cavities, material excesses can form at the margin between tooth and matrix («pooling effect»). Therefore, place a matrix band, if possible, only after the bond is applied.
- It is advisable to use a rubber dam.
- If the material is applied into the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

Recommended use

Prepare tooth structure and etch if necessary

1. Carry out cavity or endodontic treatment and prepare the root canal for the root post insertion in accordance with the general rules of adhesive technology.
2. In deep cavities, protect areas close to pulp with a suitable underfilling.

Note: Clinically sufficient bonding values are achieved (self-etch) by applying LuxaBond Universal. Additional etching is no longer needed. The bonding values can be increased, if necessary, by etching the enamel and dentin.

3. Where necessary, etch as follows:

- Selective-etch technique: Etch affected enamel with 37% phosphoric acid gel. The contact time of the etching gel on the enamel should be 20 to 60 seconds. The etched enamel should

have a chalky white appearance. If this is not the case, etching must be repeated. Accidental etching of dentin has no effect on the bonding values, but the contact time of the etching gel on dentin should not exceed 15 seconds.

- Total-etch technique: Etch dentin and enamel with 37% phosphoric acid gel. Start the etching process with the enamel and then expand this to include the dentin. The contact time of the etching gel on the enamel should be 20 to 60 seconds. The etched enamel should have a chalky white appearance. If this is not the case, etching must be repeated. The contact time on the dentin must not exceed 15 seconds.
- 4. Rinse with water for at least 15 seconds.
- 5. After the rinsing process an even, moist layer must remain on the dentin surface to be treated. If necessary, wet the dentin surface using a lightly moistened micro brush.

Apply LuxaBond Universal as an adhesive to the tooth structure

1. Mix 1 to 2 drops each of bond A and bond B in a trough in the mixing pad in a 1:1 ratio for approx. 5 seconds.
2. Apply a layer of the bond mix onto the preparation using a micro brush and work into the tooth structure for 20 seconds.
3. Remove excess material carefully.
4. Direct a stream of oil-free and water-free air at the material for at least 10 seconds until the surface has an evenly moist and shiny appearance and is firm.
5. If the mixed material in the mixing palette is still liquid, use the liquid material. If the mixed material in the mixing pad has already hardened, remix 1 drop each of bond A and bond B in a trough in the mixing pad in a 1:1 ratio for approx. 5 seconds.
6. Apply another layer of the bond mix onto the preparation using a micro brush and work into the tooth structure for 20 seconds.
7. Remove excess material carefully.
8. Direct a stream of oil-free and water-free air at the material for at least 10 seconds until the surface has an evenly moist and shiny appearance and is firm.
9. Cure the material as follows:
 - For core build-up materials, light-curing filling materials: Light-cure LuxaBond Universal for at least 20 seconds with a suitable light-curing unit.
 - For use when cementing root posts, inlays, onlays, crowns, bridges and veneers: Allow LuxaBond Universal to set. After the working time is completed, the material is sufficiently cured. In addition, the material can be light-cured with a suitable light-curing unit for at least 20 seconds.

Note: LuxaCore Z Dual is recommended for core build-up and for root post cementation. PermaCem Universal is recommended for cementing indirect restorations.

Apply LuxaBond Universal as a primer for indirect restorations made of metals (non-noble metal, precious metal), oxide and silicate ceramics or composites

1. Prepare the clean and dry restoration in accordance with manufacturer's instructions. If not otherwise described, sandblast the surfaces to be bonded with aluminum oxide ($\leq 50 \mu\text{m}$).

Note: Phosphoric acid can hinder an optimal bond with zircon oxide aluminum oxide, and metal restorations.

2. Clean the restoration with alcohol and dry with air that is free of oil and water. Etch the silicate ceramic in advance with 5% hydrofluoric acid in accordance with the manufacturer's instructions and rinse thoroughly with water spray. Dry the restoration with oil-free and water-free air.
3. Mix 1 to 2 drops each of bond A and bond B in a trough in the mixing pad in a 1:1 ratio for approx. 5 seconds.
4. Apply the bond mix onto the restoration using a micro-brush and work in for 20 seconds.
5. Remove excess bond carefully from the restoration.
6. Direct a stream of oil-free and water-free air at the material for at least 10 seconds until the surface has an evenly moist and shiny appearance and is firm.
7. Additional light-curing is optional.

Note: The bonding values can be increased by the optional light-curing.

Timing

Optional: Contact time of the etching gel on enamel	20 – 60 s
Optional: Contact time of the etching gel on dentin	max. 15 s
Working time when keeping protected from light and incl. mixing time depending on the ambient temperature	max. 60 s
Mixing time	5 s
Work-in time as adhesive on tooth: 2 layers	Massage in for 2x 20 seconds + Blow with compressed air for 10 seconds
Working time as a primer on restorations: 1 layer	Massage in for 1x 20 seconds + Blow with compressed air for 10 seconds
Optional: Light-curing time	min. 20 seconds

Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.
- The working time of the luting cement or the core build-up material may be reduced if there is contact with LuxaBond Universal.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings / precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the

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in authorities.

Performance characteristics of the device

SBS bonding values on etched and non-etched tooth substance measured in accordance with ISO29022:2013 are $\geq 7 \text{ MPa}$

Storage / disposal

- Store in a dry place in the refrigerator at $2-8^\circ\text{C}/36-46^\circ\text{F}$.
- Close bottles immediately after use.
- Do not use the material past the stipulated shelf life or if the fluid assumes a gelatinous consistency.
- Ensure that the lids of bond bottles are not switched in order to avoid any cross contamination.

Composition

Bond A: Resin matrix (HEMA, bis-GMA, MDP) approx. 97 %, catalyst, stabilizer, additives.

Bond B: Ethanol approx. 57 %, water approx. 35 %, catalyst, additives.

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Luxatemp-Glaze & Bond is a light-curing, 1-component material based on multi-functional methacrylates for surface sealing and repair of temporary crowns and bridges, as well as for the pre-treatment of indirect plastic restorations.

Indications

- Surface sealer and repair material for temporary crowns and bridges
- As a primer/bonding agent for indirect restorations of PEEK, composite, and PMMA

Contraindications

- ▶ Not to be used if allergies are known to exist with components (especially methyl methacrylate).

Notes for use

- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. Place the light as close as possible to the material. The light intensity should be a minimum of 400 mW/cm².
- ▶ Luxatemp-Glaze & Bond binds to teeth and films made of polycarbonate or poly acetate.
- ▶ Do not use Luxatemp-Glaze & Bond as a dentine or enamel adhesive.
- ▶ Avoid contact with materials containing eugenol.
- ▶ Close the bottle immediately after use.

Recommended use

Surface coating

1. Apply a thin layer of the material without entrapped air bubbles to the required surfaces.
2. Allow the material to act for 20 s and then cure each element for 20 s with a halogen or LED lamp (light intensity > 1000 mW/cm²) . Alternatively cure the entire restoration for at least 90 s in the light furnace.

Note: To ensure that no oxygen inhibition layer is formed, the material should be light-cured within 2:00 min after being applied.

Adhesion for repair purposes

1. Roughen the surface to be repaired with a cutter.
2. Apply a thin layer of the material without entrapped air bubbles to the required surfaces.
3. Expose for 20 s with a halogen lamp or ≥10 s with an LED lamp (light intensity > 1000 mW/cm²) .
4. Apply freshly mixed temporary material, e.g. Luxatemp Star, or a light-curing composite, e.g. LuxaFlow Star, to the relevant surfaces. Observe the manufacturer's instructions for use for the material used.

Bonding agent for restorations

1. Prepare the restoration according to manufacturer's specifications.
2. Apply a thin layer of the material without entrapped air bubbles to the required surfaces. Expose for 20 s with a halogen lamp or ≥10 s with an LED lamp (light intensity >1000mW/cm²).
3. Lute restoration with a suitable adhesive dual- or self-curing composite luting cement according to the manufacturer's instructions.

Storage and shelf life

- ▶ Store at 2–25 °C (36–77 °F) in a dry place.
- ▶ If stored in the refrigerator (2–8 °C/ 36–46 °F), prior to applying bring the material to room temperature (15–25 °C / 59–77 °F).
- ▶ Do not use after the expiration date.

Composition

Multifunctional acrylates, methyl methacrylate, catalysts, stabilizers, additives

Packaging

Package	
1 bottle @ 5 ml liquid, 25 brushes	REF 211927
Регистрация МИ в Росздравнадзоре www.nevacert.ru info@nevacert.ru	
50 brushes	REF 110041

Keep out of the reach of children! For dental use only!

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Instructions for useEnglish

Indication

Silanization of siliceous ceramics for improved adhesion of resin based materials, Improved bond between the looting cement and the glass fiber enforced composite post, LuxaPost

Recommended application

1. Prepare ceramic restoration with hydrofluoric acid. Carefully follow the directions for using hydrofluoric acid. Clean root post with solvents such as acetone or ethanol.
2. Mix 1 to 2 drops each of Vitique Silane Adhesive and Vitique Silane Activator in a ratio of 1:1 for approx. 15 seconds. Use silane mixture immediately.
3. Apply a thin layer of the silane mixture with a brush onto the prepared ceramic surface and let set for approx. 10 seconds.
4. Blow off the solvent carefully with a gentle stream of air.
5. Now the tooth can be cemented.

Please note

- ▶ Do not apply the material to exposed dentine or enamel.
- ▶ Avoid contact with skin and mucous membranes. In case of contact, wash immediately with water and soap. In case of contact with eyes, rinse immediately with a lot of water and contact a physician.
- ▶ So far no side effects are known. In individual cases hypersensitivity to components of the material cannot be precluded.
- ▶ Hold the bottle vertically to draw out material.

Composition

Silane; activator: distilled water, ethanol, phosphoric acid and catalyst

Storage

Store at room temperature not exceeding 25°C/77°F. Recap bottles after use. Do not use after expiration date.

Available product sizes

1 x 3 ml/2.7 g bottle of Vitique Silane Adhesive	REF 213118
1 x 3 ml/2.7 g bottle of Vitique Silane Activator	
1 brush holder + 50 brushes	

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Keep away from children! For dental use only!

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English

Product description

Ecosite Bond is a light-curing single-component bonding agent. The material can be used as an adhesive on enamel and dentin with the self-etch, total-etch or selective-etch technique in combination with light-curing filling materials.

Intended purpose

Bonding agent for direct restorations of all cavity classes with light-curing composites, compomers, and ormocers and for adhesive luting of indirect restorations with light-curing veneer cements.

Limitations of use

Do not use the material if a sufficiently dry working area or the recommended application technique is not possible.

Indications

- Loss of tooth substance due to caries, trauma, tooth abrasion or development disorders: for direct restorations made of composites, compomers and ormocers or in adhesive luting of indirect restorations with light-curing veneer cements, e.g. with Vitique
- Tooth or periodontal injuries which require composite splints
- Esthetic requirements (e.g. Skyce tooth jewelry)

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Dentist

Notes for use

- Light-curing units should have an output of 450 nm and should be checked regularly.
- The light intensity should be a minimum of 400 mW/cm². Place the light-curing unit as close as possible to the material. Observe the manufacturer's instructions.
- In order to prevent premature polymerization due to ambient light, the material must be added to the mixing pad immediately prior to application.
- Ensure that the etched or bonded surface does not become contaminated prior to each subsequent working step.
- Do not remove the oxygen inhibition layer caused by the application because this is important for the bonding process with the material applied on top.
- Observe the manufacturer's instructions for all materials that are used with the product.
- Do not subject the dentin to excessive drying, particularly after etching.
- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

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1. Prepare the cavity in accordance with the general rules for adhesive technique. It is advisable to use a rubber dam. In deep cavities, protect areas close to pulp with a suitable underfilling.
2. Where necessary, etch as follows:
 - Selective-etch technique: Etch affected enamel with 37% phosphoric acid gel. The contact time of the etching gel on the enamel should be 20 to 60 seconds. The etched enamel should have a chalky white appearance. If this is not the case, etching must be repeated. Accidental etching of dentin has no effect on the bonding values, but the response time of the etching gel on dentin should not exceed 15 seconds.
 - Total-etch technique: Etch dentin and enamel with 37% phosphoric acid gel. Start the etching process with the enamel and then expand this to include the dentin. The response time on the enamel should be 20 to 60 seconds. The etched enamel should have a chalky white appearance. If this is not the case, etching must be repeated. The contact time on the dentin must not exceed 15 seconds.
3. Rinse with water for at least 15 seconds. Remove large areas of excess water.
4. After drying, an even, moist layer must remain on the etched dentin surface requiring treatment. If necessary, wet the dentin surface using a lightly moistened micro brush.

Apply adhesive to the tooth structure

1. Place 2–3 drops of Ecosite Bond on a mixing pad. Using a micro brush, rub the bond in for 10 seconds before leaving to work for 20 seconds.
2. Blow the materials for 5 seconds to remove the solvent.

Note: Ensure that enough material remains on the tooth.

3. Light-cure the bond with a suitable LED polymerization device or a halogen polymerization device for 10 seconds.
4. Apply light-curing composite material to the cured Ecosite Bond and light-cure in accordance with manufacturer's specifications.

Interactions

Moisture, materials containing eugenol and oily air may inhibit polymerization at the contact area.

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Warnings / precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.
- Avoid contamination of the bottle content from unclean disposable brushes.

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SBS bonding values on etched and non-etched tooth substance with light-curing composite materials measured in accordance with ISO 29022:2013 are ≥ 15 MPa

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English

Product description

Icon Proximal is an innovative product for the micro-invasive treatment of lesions in the proximal area.

Icon Proximal has been specially developed for the gentle treatment of incipient proximal caries.

Intended purpose (Fig. see »Lesion depth classification«)

Micro-invasive treatment of non-cavitated cariogenic lesions (lesion depth up to D1*)

* Radiographic lesion depth classification according to bitewing images

Limitations of use

The use of Icon is contraindicated if sufficient drying and protection of the gingiva is not possible.

Indications

Non-cavitated cariogenic lesions (lesion depth up to D1*)

* Radiographic lesion depth classification according to bitewing images

Contraindications

- Do not apply the material for deeper seated lesions (D2 – D3) (Fig. see »Lesion depth classification«) with cavitated enamel (enamel defects).
- Do not apply the material where there are known allergies or contact allergies to any of the components.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended users

Dentist

Notes for use

- Light-curing units should emit at 450 nm and should be checked regularly.
- The light intensity should be a minimum of 800 mW/cm².
- Place the light unit as close to the material as possible.
- The etched surface must not be touched or contaminated with saliva or moisture before the treatment resumes. If contamination occurs after drying, re-etch for approx. 10 seconds and repeat drying with Icon Dry.
- For effective treatment results a sufficiently dry working area is critical. Therefore, appropriate measures for complete isolation must be taken. Do not use rubber dams made from thermoplastic elastomers such as Flexi-Dam (Coltène Whaledent/Hygenic).
- Discard syringes after use.
- The syringes included in the treatment unit contain sufficient material for treating two proximal lesions of one patient.
- Additional proximal surfaces can be treated consecutively during one visit.

Recommended use (for image, see »Handling«)

Preparation

1. Before the start of treatment, clean the affected tooth and adjacent teeth. Remove all cleaning residue with water spray.
2. Ensure that the working area is sufficiently dry. For this purpose, appropriate measures such as a rubber dam should be taken and used according to the manufacturer's instructions.

Separation

1. Introduce one of the enclosed dental wedges into the interproximal space (Fig. 1).
 - For better access to the proximal region the wedge handle can be bent or detached with a twisting motion.
 - In order to obtain sufficient separation of the teeth the wedges must be seated at an adequate depth in the interproximal space.
 - A separation of approx. 50 µm is recommended to ensure a successful treatment. This can be performed gradually: Insert the wedge to the point of resistance. Wait 3 to 5 seconds and gradually move the wedge deeper until sufficient separation is achieved.
 - Leave the wedge in the proximal space during the entire treatment procedure.
 - For tighter proximal areas, the use of an orthodontic rubber ring or other separation method prior to treatment (up to several days prior) will ensure adequate separation.

Exposing the lesion body

For effective treatment results the hypermineralization must be removed.

1. Screw the Proximal Tip onto the Icon Vestibular Etch syringe, and introduce the application tip into the proximal area. Make sure that the green side of the Proximal Tip faces in the direction of the surface to be treated (Fig. 2).
 - The material only dispenses on the green side of the Proximal Tip.
2. Apply Icon Vestibular Etch with a slight excess to the lesion site. 1.5 to 2 turns of the syringe shaft correspond approximately to the amount to be applied. Allow Icon Vestibular Etch to act for 2 minutes. Move the applicator slightly so that the etching gel will be activated. Remove excess material (Fig. 3).
3. Remove the application tip from the interproximal space. Aspirate off Icon Vestibular Etch and rinse with water for at least 30 seconds. Then dry thoroughly with air that is free of oil and water (Fig. 4).

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For effective treatment results it is necessary to dry the lesion to its full depth.

1. Screw the application tip onto the Icon Dry syringe, apply an ample amount of material onto the lesion and allow to take effect for 30 seconds (Fig. 5). Dry carefully with oil-free and water-free air.

Infiltration

The lesion can be penetrated to its full depth by the Infiltrant after complete drying.

- Do not apply Icon Infiltrant under direct operating light as this may cause the material to set prematurely.
1. Screw a new Proximal Tip onto the Icon Infiltrant syringe, and introduce the application tip into the proximal area. Make sure that the green side of the Proximal Tip faces in the direction of the surface to be treated.
 - The material only dispenses on the green side of the Proximal Tip.
 2. Apply Icon Infiltrant with a slight excess to the lesion site. 1.5 to 2 turns of the syringe shaft correspond approximately to the amount to be applied.
 3. Allow Icon Infiltrant to penetrate for 3 minutes, occasionally activating the Infiltrant by gently moving the applicator and then top up as necessary (Fig. 6).
 4. Remove the application tip from the interproximal space. Remove excess material with dental floss.
 5. Light-cure Icon Infiltrant from all sides for at least 40 seconds.
 6. Screw a new Proximal Tip onto the Icon Infiltrant syringe, repeat the application, and let set for 1 minute. Remove the application tip from the interproximal space. Remove excess material with dental floss. Then light-cure from all sides for at least 40 seconds (Fig. 7).
 7. If a second lesion needs to be treated, repeat the recommended application described above for the affected site.
 8. Remove the wedge and rubber dam. Use polishing strips for the surface finish.

Interactions

Do not use together with substances which inhibit polymerization. This is especially valid for materials which contain eugenol, such as temporary cements.

Residual risks / side effects

- In some cases, contact allergies with similar composite products have been reported.
- Contact of Icon Vestibular Etch with the oral mucosa causes a white surface coloration. This will subside after a few days.

Warnings / precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Caution! Icon Vestibular Etch can cause chemical burns if not used properly.

- Avoid contact with soft oral tissue, eyes, and skin. In case of inadvertent contact rinse with copious amounts of water, and, if necessary, consult a physician.
- Limit the application to the surface to be treated. Avoid contact of the etching gel with gingiva, exposed dentin, root element and adjacent teeth.
- Use only the enclosed syringe and application tip to apply the etching gel in the patient's mouth.

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Penetration coefficient: approx. 147 [cm/s]

Storage / disposal

- Store at 2-25 °C (36–77 °F) in a dry place.
- Protect against exposure to direct sunlight!
- Do not use after the expiration date!

Patient card + sticker (Fig. see »Patient card«)

This material is not radiopaque. For the documentation of the infiltrated surfaces and lesion depths, the enclosed patient card can be filled out at the time of the treatment and follow-up visits and be given to the patient. In addition, the enclosed sticker facilitates the documentation in the patient’s record.

1. Front side: Enter the name of the patient and add the office stamp.
2. Tooth diagram: Mark the treated area.
3. Record the treated tooth [1], the treated surface [2] (d = distal, m = mesial, v = vestibular (buccal, cervical, labial), l = lingual or palatal), lesion depth*, and the treatment date [3].
4. Follow-up visit [4]: Record the lesion depth* and date of the follow-up appointment. Annual check-ups are recommended.

*see »Lesion depth classification«

Composition

- Icon Vestibular Etch: Hydrochloric acid, water, glycerol, pyrogenic silicic acid, surface-active substances
- Icon Dry: Ethanol
- Icon Infiltrant: Methacrylate-based resin matrix, initiators, additives

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English

Product description

The etching gel Icon Vestibular Etch forms part of infiltration treatment with Icon. For the entire infiltration treatment, the following instructions for use must be strictly observed: Icon Vestibular

Intended purpose

Removing the pseudointact surface layer in connection with infiltration treatment with Icon.

Limitations of use

Sufficient drying and protection of the gingiva not possible.

Indications

- Non-cavitated carious lesions
- Fluorosis
- Idiopathically or traumatically caused spots

Contraindications

Do not use this material:
in cases of allergies to any of the ingredients or in cases of contact allergies.

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Dentist

Notes for use

- If the material is applied into the patient's mouth with a single use applicator, the single use applicator must only be used on one patient for hygiene reasons.
- The etched surface must not be touched or contaminated with saliva or moisture until the treatment resumes. If contamination occurs after drying re-etch for approx. 10 seconds and repeat drying with Icon Dry.
- For effective treatment results a sufficiently dry working area is critical. Therefore, appropriate measures for a completely isolation must be taken. A liquid dam can be used as well. Do not use rubber dams made from thermoplastic elastomers such as Flexi-Dam (Coltène Whaledent/Hygenic).

Recommended use (for image, see »Handling«)

Application of Icon Vestibular Etch:

- Screw the Vestibular Tip onto the Icon Vestibular Etch syringe.
- Apply an ample amount of Icon Vestibular Etch onto the lesion site by turning the syringe shaft carefully, let set for 2 minutes and activate it by moving (Fig.1). Remove excess material with dental wadding.
- Suck off Icon Vestibular Etch and rinse with water for at least 30 s. Then dry using oil-free and water-free air (Fig. 2).

Interactions

Do not use rubber dams made from thermoplastic elastomers such as Flexi-Dam (Coltène Whaledent/Hygenic).

Residual risks/side effects

- There are no known side effects to date. The residual risk of hypersensitivity in relation to allergies to components of the material cannot be ruled out.
- Contact of Icon Vestibular Etch with the oral mucosa causes a white surface coloration. This will subside after a few days.

Warnings / precautions

- Keep out of the reach of children.
- Caution! Icon Vestibular Etch can cause chemical burns
- Avoid contact with soft oral tissue, eyes, and skin. In case of inadvertent contact rinse with copious amounts of water, and, if necessary, consult a physician.
- Limit the application to the surface to be treated. Avoid contact of the etching gel with the gingiva, exposed dentine, root element and adjacent teeth.
- Use only the enclosed syringe and application tip to apply the etching gel in the patient's mouth.
- Dispose of application tip after use. Reattach syringe closure.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

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pH value ≤ 0

Storage / disposal

- Store in a dry place at 2 to 25°C (36 to 77°F)!
- Protect against exposure to direct sunlight!
- Do not use after the expiration date!

Composition

Hydrochloric acid, pyrogenic silicic acid, surface-active substances

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