IT **CE** 0477

RS-OIL

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EN **CE** 0477


RS-OIL

WARNING: This instruction sheet contains important information and warnings on the correct use of the medical device. Carefully read it before starting any operation and closely follow the directions for use.

COMPOSITION AND PACKAGING

RS-OIL is a sterile and pyrogen-free, medical grade, pure silicone oil (100% polydimethylsiloxane) with a viscosity of 1,000, 1,300, 2,000, 5,000 or 5,700 cSt. Product is available in vials and syringes, packaged in individual pouches. Labels are included to allow product traceability.

Silicone oil	Refractive index (35 ± 2°C)
1,000	1,399
1,300	1,399
2,000	1,399
5,000	1,399
5,700	1,399

INDICATIONS FOR USE

Long-term intraocular tamponade in vitreoretinal surgery. Any other use must be considered improper. The amount of silicone oil to be injected is at the exclusive discretion of the surgeon; an underfill may result in an ineffective inferior tamponade and an overfill may result in corneal abnormalities and elevated intraocular pressure. The length of the indwelling period of the silicone oil in the patient's eye is at the exclusive discretion of the surgeon. The selection of patients, and the choice of surgical procedure are at the exclusive discretion of the surgeon.

DIRECTIONS FOR USE

Make sure that the package is intact.

Vials

Remove the vial from the pouch. Under aseptic conditions, transfer the content of the vial to a sterile device suitable for injection, taking care not to introduce air bubbles. Do not mix with any other substances prior to use.

Syringes

Remove the syringe from the pouch. **Manual injection:** turn the syringe adapter (1) counter-clockwise to remove it; screw the piston (2) in the rubber stopper (3) of the syringe. **Mechanical injection system:** discard the piston (2) and turn the syringe adapter (1) clockwise to firmly fix it to the syringe. Attach the syringe to the mechanical injection system oil by means of an appropriate connection. **BEFORE INJECTION, PURGE SYRINGE AND CANNULA BY PUSHING OUT APPROXIMATELY 0,5 ML OF SILICONE OIL, VERIFYING THE CORRECT SLIDING OF THE RUBBER STOPPER (3) INTO THE SYRINGE.**

IF PERFLUORO-N-OCTANE OR PERFLUORODECALIN ARE USED DURING THE SAME SURGERY, PERFORM FLUID-AIR EXCHANGES TO ENSURE COMPLETE REMOVAL OF PERFLUOROCARBONS, BEFORE INJECTING SILICONE OIL.

Inject silicone oil according to the accepted medical practice. Properly discard any remaining silicone oil. Once satisfactory outcomes are achieved, remove silicone oil from the patient's eye, and dispose of as sanitary waste.

CONTRAINDICATIONS

- Silicone oil is a relatively inert substance, nevertheless it is recommended that it is not used in patients with known hypersensitivity to the declared ingredients.
- Silicone oil is contraindicated in pseudophakic patients with silicone intraocular lenses or as a permanent vitreous substitute.

WARNINGS AND PRECAUTIONS

- Integrity is warranted only if package is intact. Do not use if package is open or damaged.
- The expiry date is shown on product label. Do not use after the expiry date.
- Do not shake and do not admix with any other substances.
- Single use only – Do not sterilize. Use or sterilization of single use products that have been already used may cause cross-contamination.
- BEFORE INJECTION, PURGE SYRINGE AND CANNULA BY PUSHING OUT APPROXIMATELY 0,5 ML OF SILICONE OIL, VERIFYING THE CORRECT SLIDING OF THE RUBBER STOPPER (3) INTO THE SYRINGE.**

- IN THE LITERATURE, THE PERMANENCE OF PERFLUOROCARBON RESIDUES HAS BEEN ASSOCIATED WITH SILICONE OIL STICKINESS AND FLOATING. PERFORM FLUID-AIR EXCHANGES TO ENSURE COMPLETE REMOVAL OF PERFLUOROCARBONS, BEFORE INJECTING SILICONE OIL.**
- Prior to inject silicone oil, make sure that there are no holes or retinal discontinuities, which can trap silicone oil. Avoid any contact between silicone oil and corneal endothelium.
- During injection pay attention that any residues are not being trapped within retinal discontinuities.
- Product shall be used by qualified medical personnel only, in surgery room. Any residues must be disposed of as sanitary waste. In case of accidental breakage or spillage of silicone oil before or during use, confine the spill and absorb in inert material, place in a plastic container and dispose of in compliance with procedures in force.
- Dispose of empty packages in compliance with regulations in force.

SIDE EFFECTS

- Temporary elevated intraocular pressure, hypotony, glaucoma and cataract can occur after surgery.
- The migration of silicone oil into the anterior chamber has been associated with corneal changes, which include edema, hazing, Descemet's folds, or corneal decompensation.

INTERACTIONS WITH OTHER SUBSTANCES

Silicone oil may interact with strong oxidizing agents.

STERILIZATION

Vials and syringes are sterilized by moist heat. Single use only - Do not sterilize.

STORAGE

Store at room temperature.

For any report on this product, quote the lot number on package. For any comments or suggestions, go to www.alchimiasrl.com

MANUFACTURER
AL.CH.I.M.I.A. S.r.l.
Viale Austria 14
35020 Ponte San Nicolò (PD) - Italy
e-mail: info@alchimiasrl.com

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ECT “easy connection tubing”

WARNING: This instruction sheet contains important information and warnings for the correct use of the medical device. Carefully read it before starting any operations and closely follow the directions for use.

COMPOSITION AND PACKAGING

ECT “easy connection tubing” are sterile polyethylene or PVC connection tubing with polythene, acetal or polycarbonate connectors, which are compatible with the following silicone oil mechanical injection systems:

ECT 002-00, accessory, compatible with Millennium™, Stellaris® PC;

ECT 003-00, accessory, compatible with Accurus®, Constellation®, Associate® 2500, EVA;

ECT 004-00, accessory, compatible with Orbit, OS3, OS4, FAROS™, R-Evolution.

ECT are packaged in individual sterile pouches and each box contains one instruction sheet.

INDICATIONS FOR USE

ECT “easy connection tubing” are accessories intended to connect **RS-OIL ECS**, a medical device made of fluid, sterile and apyrogenic silicone oil for long-term tissue tamponade during vitreoretinal surgery, to silicone oil mechanical injection systems. **ECT** must only be used together with **RS-OIL ECS** by qualified medical personnel in an operating theatre. Any other use must be considered improper.

DIRECTIONS FOR USE

1. Visually inspect the package for integrity.
2. Remove the **ECT** from the pouch within the sterile field.
3. Open the pouch of **RS-OIL ECS** and take the syringe out. Do not shake.
4. Turn the syringe adapter clockwise to fix it securely to the syringe.
5. Remove the caps on either end of the **ECT**, if present.
6. Turn the male connector (**A**) of the **ECT** clockwise to attach it to the syringe adapter.
7. Attach the other end of the **ECT (B)** to the mechanical injection system by exerting a slight pressure on the connector and a contemporaneous clockwise turn, if necessary.
8. Make sure that the **ECT** is correctly and firmly fixed to both the mechanical injection system and the syringe adapter.
9. See the instruction sheet of **RS-OIL ECS** for further instructions, if necessary.

WARNINGS

- Store at the temperatures indicated on label. Protect from light.
- Single use only – Do not resterilize. Use or sterilization of single use products that have been already used may cause cross-contamination.
- Do not use if packaging assuring sterility is damaged.
- Do not use after the expiry date.
- Dispose of in compliance with regulations in force – Do not litter the environment.

STERILIZATION

Medical devices are sterile by ethylene oxide. Medical devices are for single use only - Do not resterilize.

STORAGE

Store at room temperature.

For any report on this product, quote the lot number on package. For any comments or suggestions, go to www.alchimiasrl.com.

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AL.CHI.MI.A. S.r.l.