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WARNING

Read the Instructions for Use carefully before using the Pressio® Monitoring Kit, Parenchymal with Bolt (PSO-PB or PSO-PBT).

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

1. Figures

Figure 1. Pressio® ICP Monitoring Kit, Parenchymal With Bolt (Model PSO-PB)

A: Pressio® Catheter for Intracranial Pressure Measurement With Its Dongle B: Guide Stylet C: Bolt D: Allen Key E: Drill Bit With Adjustable Stop

Figure 2. Pressio® ICP and ICT Monitoring Kit, Parenchymal With Bolt (Model PSO-PBT)

A: Pressio® Catheter for Intracranial Pressure and Temperature Measurement With Its Dongle B: Guide Stylet C: Bolt D: Allen Key E: Drill Bit With Adjustable Stop

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Furthermore, the PSO-PBT model makes it possible to measure intracranial temperature (ICT).

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- Young children (< 3 years), because the bony plate of their skull is not thick enough to use the fixation with the bolt.

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- a Pressio ICP Monitor PSO-3000 or PSO-4000 (hereinafter referred to as Monitor),
- a Pressio Interface PSO-IN00 (hereinafter referred to as Interface).

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The Monitoring Kit, **Parenchymal with Bolt PSO-PB (Figure 1)** allows continuous monitoring of ICP.

The Monitoring Kit, **Parenchymal with Bolt PSO-PBT (Figure 2)** allows continuous monitoring of ICP and ICT.

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
4.2.2. Bolt

The bolt enables the introduction and fixation of the Catheter in the parenchyma.

- [4] Screw thread to be inserted into the burr hole.
- [5] Spacer ring, mounted on the bolt, to adjust the seating depth for the bolt, depending on the thickness of the skull bone of the patient.
- [6] Bolt shoulder.
- [7] Wing nut to help screw the bolt into the skull bone.
- [8] Tightening nut to fix the Catheter at the required depth in the parenchyma.

4: Screw thread 5: Spacer ring 6: Bolt shoulder 7: Wing nut 8: Tightening nut

4.2.3. Guide stylet

The guide stylet is inserted into the bolt (nut [5] unscrewed). It pierces the dura mater to facilitate the introduction of the Catheter.

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5. Fix the helical drill bit to a drill.

7.1. Preparing the Catheter Implantation Site

Observe aseptic neurosurgical techniques to prepare the Catheter implantation site and the insertion of the bolt.

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6. Perforate the internal and external skull plates avoiding any possibility of a parenchymal injury when passing through the internal plate.

WARNING

Make sure to make a completely round hole with the drill bit.

WARNING

Do not move the shaft of the drill during perforation, as this could affect the hold of the bolt in the burr hole.

7. Having crossed the internal plate, remove the drill bit.
8. Irrigate the burr hole with sterile saline solution and remove the debris.


7.1.3. Putting the Bolt in Place

Screw the bolt into the skull using the wing nut. The spacer ring may be used to reduce the implantation depth of the bolt, otherwise remove and discard it.

CAUTION

Make sure the bolt is not screwed too deeply. The bolt shoulder (8) and its spacer ring are a simple visual marker and not a stop. They will not prevent deeper perforation when screwing the bolt.

It is the surgeon's responsibility to judge the depth of the bolt seating depending on the thickness of the skull.

CAUTION

Do not unscrew and then retighten the bolt in the same burr hole. This could affect the hold of the bolt in the burr hole.

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
3. Observing aseptic technique, connect the Catheter extension cable (non-sterile) to the Catheter, aligning the blue arrow on the Catheter dongle [8] with the blue arrow on the Catheter extension cable connector [9].

The Catheter dongle must be completely inserted into the Catheter extension cable connector. A blue line on the Catheter dongle indicates the point of complete insertion.

When correctly inserted:


- if using a Monitor, a message displays on the screen,
- if using an Interface, the zeroing diode flashes quickly.

NOTE

The appearance of error codes, “E001”, “E002” or “E005” on the Monitor after connecting the Catheter to the Catheter extension cable may be a sign of an incorrect connection.

In this case, make sure that the Catheter dongle is pushed right up to the stop in the Catheter extension cable, and that the latter is correctly connected to the Monitor.

4. Within the sterile field, immerse the sensor (metal tip) of the Catheter [10] in the cup of sterile saline solution [11], without touching the cup.

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5. While the sensor is immersed, press the Zero button  on the Monitor touch screen, or on the Interface:
Zeroing the Catheter requires approximatively 4 seconds. Do not move the Catheter during the zeroing procedure.
 - If the zeroing procedure is successful:
 - the Monitor touch screen displays a message indicating that the Catheter is ready to be implanted,
 - the Interface zeroing diode remains lit.
 - If the zeroing procedure is not successful:
 - follow the instructions displayed on the Monitor touch screen,
 - the Interface emits a series of beeps and sends a pressure value of 360 mmHg to the patient bedside monitor. For more information, refer to the instructions for use of the Interface.

NOTE

Zeroing the Catheter can only be performed once. Upon completion, the zero calibration information is stored in the Catheter dogle. It allows the Catheter to be disconnected from a Monitoring System, and re-connected to any Monitoring System, without losing the zero calibration information.

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


1. Turn the tightening nut counter-clockwise to open the bolt.
2. Carefully introduce the guide stylet into the bolt to create a passage through the parenchyma, to facilitate the introduction of the Catheter.

3. Remove the guide stylet and irrigate the bolt with sterile saline solution.

4. Introduce the zeroed Catheter into the bolt.
5. Push the Catheter through the parenchyma. Position it at the desired depth using the centimeter marking.

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

NOTE

To check the implantation depth, the Catheter has centimeter markings (see figure below). A thicker marker located 4 cm from the Catheter end represents the length of the bolt (A). Once the Catheter is introduced into the bolt, this marker is no longer visible when the sensor at the end of the Catheter just emerges from the bolt (B).

Push the Catheter 1 to 3 centimeters further into the cerebral parenchyma in relation to the thicker marker.

6. Once the chosen introduction depth is reached, turn the tightening nut on the bolt clockwise to hold the Catheter in place.

NOTE

The tightening nut must be fully tightened to ensure optimal hold of the Catheter in the bolt.

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1 mmHg corresponds to 13.60 mmH₂O and to 133 Pa.

NOTE

The appearance of error codes “- - -”, “999”, “E001”, “E002”, “E005”, or the message “Connect sensor” on the Monitor after the Catheter implantation, when the Catheter is correctly connected may be a sign of damage to the dangle or to the sensor located on the tip of the Catheter.

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After the first 24 hours, the Monitor displays the duration of the Catheter implantation in the form of the message “**Implantation: X days**”, at the top left corner of the touch screen, for the entire duration of its implantation.

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Non-clinical testing has demonstrated the Monitoring Kit, Parenchymal with Bolt (PSO-PB) for parenchymal ICP monitoring and the Monitoring Kit, Parenchymal with Bolt (PSO-PBT) for parenchymal ICP and ICT monitoring are MR Conditional.

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Under the scan conditions defined above, the Monitoring Kit, Parenchymal with Bolt (PSO-PB) for parenchymal ICP monitoring or the Monitoring Kit, Parenchymal with Bolt (PSO-PBT) for parenchymal ICP and ICT monitoring is expected to produce a maximum temperature rise of less than 2.2 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the length of the implanted Catheter extends approximately 55 mm from the Monitoring Kit, Parenchymal with Bolt (PSO-PB) for parenchymal ICP monitoring or the Monitoring Kit, Parenchymal with Bolt (PSO-PBT) for parenchymal ICP and ICT monitoring when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

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The "MR Conditional" status only concerns Monitoring Kits for ICP or ICP and ICT. The Monitors, Interface, cables and the temperature module are considered to be *"MR Unsafe"* (not MRI-compatible) and must not be exposed to an MRI environment.

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- The non-implanted part of the Catheter is coiled around the MRI Support, which itself is centered in the MRI. For more information, refer to the indications mentioned in section *Positioning of the Monitoring Kits and the MRI Support* below.
- The MRI has a static magnetic field less than or equal to 3 Tesla.
- The spatial gradient of the static magnetic field is less than or equal to 19 T/m.
- The whole body averaged SAR (Specific Absorption Rate) is less than or equal to 2 W/kg (normal mode) for 15 minutes continuous exposure at 3 Tesla.
- The whole body averaged SAR (Specific Absorption Rate) is less than or equal to 2 W/kg (normal mode) for 15 minutes continuous exposure at 1.5 Tesla.
- No limitation in the magnetic field gradients.

- The size of the artifact on the images may reach about 5.5 cm depending on the implantation depth for MRI examinations performed in 3 Tesla echo sequence.
- Do not use an emission/reception head coil (RF emission coil for the head) or a simple RF emission head coil. Only use an RF emission/reception whole body coil or an RF emission whole body coil with a simple RF reception head coil.
- Do not scan a patient with a high body temperature.

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NOTE

Carefully remove the Catheter from the bolt.

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12. Symbols

	The Catheters are not protected against the effects of a cardiac defibrillator.
	Refer to the instructions for use
	Do not re-sterilize
	Do not re-use
	MR Conditional. It has been demonstrated that the use of ICP or ICP and ICT Monitoring Kits within an MRI setting does not pose any risks as long as the specific conditions of use are complied with.

	Serial number
	Lot number
	Catalog reference
	Sterilization Method using Ethylene Oxide
	Manufacturer
	Use until

	BF TYPE EQUIPMENT: Providing an appropriate degree of protection against electric shock, having a Type BF insulated applied section (floating).
	Temperature conditions for storage and transport
	Keep away from liquids

	Fragile, handle with care
	CE conformity marking

13. Technical Specifications

Table 1. General specifications

Item	Specification
Sensor capsule diameter	1.2 mm
Catheter diameter	2.1 Fr (0.7 mm) $D(mm) = Fr / 3$ $Fr = D(mm) \times 3$
Sensor capsule material	Titanium
Catheter length	1,000 mm (radiopaque on the entire length)
Catheter sheath material	Polymide
Depth markers	Marking every centimeter between 1 and 10 cm from the proximal end. Thicker marking at 4 cm. Additional marking at 15 cm.
Leakage current	Less than 10 μA to 120 V AC
Weight	11 g
Drill bit diameter	2.7 mm, adapted to the screw thread of the bolt
Note The Monitoring Kits are latex- and phthalate-free.	

Table 2. Pressure sensor specifications

Item	Specification
Type of sensor	Piezoresistive with strain gauge type on silicon
Reference pressure	Atmospheric pressure
Display range (complete system)	-40 mmHg to +100 mmHg (PSO-3000) -40 mmHg to +150 mmHg (PSO-4000)
Accuracy of the pressure measurement (linearity and hysteresis) (complete system)	$\pm 2\%$ reading in the range 0 to 100 mmHg
Bandwidth	> 100 Hz
Functional range of over-pressure without damage	-700 mmHg to +1,250 mmHg
Temperature coefficient	0.1 mmHg/°C max
Input resistance	667 Ohms
Output resistance	810 Ohms

Item	Specification
Excitation voltage	1 to 8 V AC or DC
Drift from zero	Less than 1 mmHg during the first day, at 37°C Less than 2 mmHg during the first week, at 37°C

Table 3. Temperature sensor specifications

Item	Specification
Type of sensor	Thermistor
Unit	Celsius (°C) or Fahrenheit (°F) $^{\circ}F = ^{\circ}C \times (9/5) + 32$ $^{\circ}C = (^{\circ}F - 32) \times (5/9)$ Equivalence in Kelvin (K): $K = ^{\circ}C + 273.15$ $K = ^{\circ}F \times (5/9) + 255.37$
Display range (complete system)	+20 °C to +45 °C (+68 °F to +113 °F)
Maximum tolerated error (accuracy) in the reference conditions (complete system)	± 0.2 °C max. from 25 °C to 45 °C (± 0.4 °F max. from 77 °F to 113 °F) ± 0.4 °C max. from 20 °C to 25 °C (± 0.7 °F max. from 68 °F to 77 °F)
Resolution	0.1 °C/0.1 °F

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Table 5. Catheter history specifications

Item	Specification
Data accessible	First 15 days of monitoring
ICP data	A value every 20 seconds. This value is an average of the previous 20 seconds.
ICT data	A value every 60 seconds. This value is an average of the previous 60 seconds.

Item	Specification
Note	
	Data can only be consulted with a Pressio 2 ICP Monitor.

14. References

Table 6. Pressio® Catheter Kits

PSO-PB	Pressio® ICP Monitoring Kit, Parenchymal with Bolt
PSO-PBT	Pressio® ICP and ICT Monitoring Kit, Parenchymal with Bolt
PSO-PT	Pressio® ICP Monitoring Kit, Parenchymal Tunneling
PSO-PTT	Pressio® ICP and ICT Monitoring Kit, Parenchymal Tunneling
PSO-VT	Pressio® ICP Monitoring Kit, Ventricular Tunneling with external CSF drainage function
PSO-VTT	Pressio® ICP and ICT Monitoring Kit, Ventricular Tunneling with external CSF drainage function

Table 7. Pressio® Monitoring System

PSO-3000	Pressio® ICP Monitor Power cable and Catheter extension cable (PSO-EC20) included
PSO-4000	Pressio® 2 ICP Monitor Power cable and Catheter extension cable (PSO-EC30) included
PSO-IN00	Pressio® ICP Monitoring Interface Catheter extension cable (PSO-EC20) included

Table 8. Pressio® Accessories

PSO-EC20	Catheter extension cable Use only with a Pressio® ICP Monitor (PSO-3000) or a Pressio® Interface (PSO-IN00)
PSO-EC30	Catheter extension cable Use only with a Pressio® 2 ICP Monitor (PSO-4000)

PSO-MCxx	Pressure cable Use only with a Pressio® 2 ICP Monitor (PSO-4000) <ul style="list-style-type: none"> – MC01: PHILIPS (AGILENT) - 12 pins – MC02: SIEMENS (SIRECUST) - 10 pins – MC03: SPACELABS & MINDRAY - 6 pins – MC04: GE DATEX-Ohmeda - 10 pins – MC05: GE Solar (MARQUETTE) - 11 pins – MC08: NIHON KOHDEN - 5 pins
PSO-MCT-Y	Temperature cable Use only with a Pressio® 2 ICP Monitor (PSO-4000) <ul style="list-style-type: none"> – MCT-A: PHILIPS (AGILENT) - 2 pins – MCT-B: SIEMENS - 7 pins – MCT-C: SPACELABS - 10 pins – MCT-E: GE Solar (MARQUETTE), GE DATEX-Ohmeda - 11 pins – MCT-F: HELIGE, DATEX-Ohmeda, NIHON KOHDEN, MINDRAY & DATASCOPE - JACK 6.35 mm
PSO-MT00	Intracranial temperature module Enables temperature values to be displayed on the Pressio® ICP Monitor (PSO-3000).
PSO-TX00	Standard transmitter Enables the pressure and temperature value to be transmitted to a computer. Use only with a Pressio® ICP Monitor (PSO-3000).
PSO-DR	Single use, disposable drill
PSO-MRI	Pressio® MRI Support For positioning Pressio® Catheter during MRI examination

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Table of Contents

Pressio® Monitoring kit, Parenchymal tunneling	7
1. Figures	7
2. Indications	7
3. Contraindications	7
4. Description	7
4.1. Catheters	7
4.2. Accessories	8
4.2.1. Drill bit	8
4.2.2. Tunneling needle	8
4.2.3. Fixation tab	8
4.2.4. Disposable drill	8
5. Sterilisation - Decontamination	9
6. Setting-up the Monitoring System	9
7. Preparing, Implanting and Monitoring	9
7.1. Preparing the Catheter Implantation Site	9
7.1.1. Choice of the Implantation Area	9
7.1.2. Intracranial Access	9
7.2. Zeroing the Catheter	10
7.2.1. Prerequisites	10
7.2.2. Procedure	10
7.3. Implanting the Catheter	10
7.3.1. Prerequisites	11
7.3.2. Procedure	11
7.3.3. Tunneling the Catheter	11
7.3.4. Implanting the Catheter	11
7.3.5. Fixing the Catheter	12
7.3.6. Specific case: Sub-dural Implantation	12
7.4. Monitoring, Caring and Nursing of the Patient	12
7.5. Explanting the Catheter	13
8. Behaviour During an MRI Examination	13
8.1. General Points	13
8.2. Positioning of the Monitoring Kits and the MRI Support	14
8.3. In Vitro Test Conclusions	14
9. Complications / Side Effects	14
9.1. Infection	15
9.2. Cerebral Hemorrhage	15
10. Destruction or return of product	15
10.1. Destruction After Use	15
10.2. Return of Products	15
11. Warranty	15
12. Symbols	15
13. Technical Specifications	16
14. References	17
Kit Pressio® pour Monitoring, Parenchymateux Tunnelisable	19
1. Figures	19
2. Indications	19
3. Contre-indications	19
4. Description	19
4.1. Cathéters	19
4.2. Accessoires	20
4.2.1. Trépan	20
4.2.2. Aiguille de tunnelisation	20

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WARNING

Read the Instructions for Use carefully before using the Pressio® Monitoring Kit, Parenchymal Tunneling (PSO-PT or PSO-PTT).

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

1. Figures

Figure 1. Pressio® ICP Monitoring Kit, Parenchymal Tunneling (Model PSO-PT)

A: Pressio® Catheter for Intracranial Pressure Measurement With Its Dongle B: Tunneling Needle C: Fixation Tab D: Allen Key E: Drill Bit With Adjustable Stop

Figure 2. Pressio® ICP and ICT Monitoring Kit, Parenchymal Tunneling (Model PSO-PTT)

A: Pressio® Catheter for Intracranial Pressure and Temperature Measurement With Its Dongle B: Tunneling Needle C: Fixation Tab D: Allen Key E: Drill Bit With Adjustable Stop

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Furthermore, the PSO-PTT model makes it possible to measure intracranial temperature (ICT).

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- Parenchymal tunneling is contraindicated in newborns (< 3.5 kg).

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The Monitoring Kit, **Parenchymal Tunneling PSO-PT (Figure 1)** allows continuous monitoring of ICP.

The Monitoring Kit, **Parenchymal Tunneling PSO-PTT (Figure 2)** allows continuous monitoring of ICP and ICT.

It is designed to be implanted in the cerebral parenchyma.

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In both kits, the Catheter has a marking every centimeter between 1 and 10 cm from the proximal end, and a marker at 15 cm. The marking at 4 cm is thicker.

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4.2.2. Tunneling needle

- [4] Tunneling needle to create a passage for the Catheter under the scalp.
- [5] Stylet to create a passage for the Catheter. It avoids coring during the insertion of the needle under the skin. The stylet is delivered already inserted into the needle.

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in decubitus without interfering with the intracranial pressure surveillance function.

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5. Fix the helical drill bit to a drill.
6. Perforate the internal and external skull plates avoiding any possibility of a parenchymal injury when passing through the internal plate.
7. Having crossed the internal plate, remove the drill bit.
8. Irrigate the burr hole with sterile saline solution and remove the debris.
9. Pierce the dura mater.

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7.1.1. Choice of the Implantation Area

As an example, the right and left prefrontal areas are the main implantation areas. These areas enable the patient to turn his/her head from one side to the other whilst remaining

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
3. Observing aseptic technique, connect the Catheter extension cable (non-sterile) to the Catheter, aligning the blue arrow on the Catheter dongle [7] with the blue arrow on the Catheter extension cable connector [8].

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
The appearance of error codes, "E001", "E002" or "E005" on the Monitor after connecting the Catheter to the Catheter extension cable may be a sign of an incorrect connection.

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5. While the sensor is immersed, press the Zero button  on the Monitor touch screen, or on the Interface:
Zeroing the Catheter requires approximately 4 seconds. Do not move the Catheter during the zeroing procedure.
 - If the zeroing procedure is successful:
 - the Monitor touch screen displays a message indicating that the Catheter is ready to be implanted,
 - the Interface zeroing diode remains lit.
 - If the zeroing procedure is not successful:
 - follow the instructions displayed on the Monitor touch screen,
 - the Interface emits a series of beeps and sends a pressure value of 360 mmHg to the patient bedside monitor. For more information, refer to the instructions for use of the Interface.

Zeroing the Catheter can only be performed once. Upon completion, the zero calibration information is stored in the Catheter dangle. It allows the Catheter to be disconnected from a Monitoring System, and re-connected to any Monitoring System, without losing the zero calibration information.

Before implanting the Catheter, check that the screen of the Monitor (or of the patient bedside monitor if the Interface is used) displays the value “0 mmHg”.

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The final implantation of the device must satisfy the conditions for optimal positioning of the sensor in the parenchyma.

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

7.3.3. Tunneling the Catheter


1. Make a small incision at the chosen site for the emergence of the Catheter (incision b).

2. Starting at the incision made for the burr hole (implantation site: incision a), insert the tunneling needle between the scalp and the skull in the direction of the emergence site (incision b).

3. Remove the stylet from the needle.


4. Delicately thread the Catheter into the needle from the emergence site (distal) to the implantation site (proximal).


CAUTION

The end of the tunneling needle is sharp, introduce the Catheter very carefully.

Make sure at least 15 cm of Catheter have been threaded through the needle to the implantation site.


5. Remove the needle from under the scalp, over the end of the Catheter.

**7.3.4. Implanting the Catheter****CAUTION**


To perform a parenchymal measurement of intracranial pressure, the implantation depth must be 1 to 3 cm in the cerebral parenchyma.

1. Introduce the zeroed Catheter into the parenchyma up to the desired implantation depth using the centimeter markings, taking the thickness of the skull into account.
2. Bend the Catheter slightly at the exit from the skull.

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
1 mmHg corresponds to 13.60 mmH₂O and to 133 Pa.

NOTE

The appearance of error codes “- - -”, “999”, “E001”, “E002”, “E005”, or the message “Connect sensor” on the Monitor after the Catheter implantation, when the Catheter is correctly connected may be a sign of damage to the dongle or to the sensor located on the tip of the Catheter.

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7.3.6. Specific case: Sub-dural Implantation

Although this is not normal use, it is also possible to implant the Catheter sub-durally.

In such cases, the implantation of the Catheter is performed under the responsibility of an anaesthetist or an experienced neurosurgeon. The clinical condition of the patient must, therefore, be subject to increased surveillance.

1. Fold the Catheter to make an elbow. The part between the elbow and the Catheter end will be implanted sub-durally. Make the elbow so that the sensor(s) is directed to face the cortex once the Catheter is implanted.
2. To implant the Catheter, insert its end through the pre-perforated dura mater and push until it reaches the interior of the sub-dural space.

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After the first 24 hours, the Monitor displays the duration of the Catheter implantation in the form of the message **"Implantation: X days"**, at the top left corner of the touch screen, for the entire duration of its implantation.

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When this duration is exceeded:

- on a Monitor, the message **"Implantation: X days"** flashes at the top left corner of the touch screen of the Monitor,
- on the Interface, the diode showing the implantation duration for the Catheter flashes.

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8. Behaviour During an MRI Examination**MRI Safety Information**

Non-clinical testing has demonstrated the Monitoring Kit, Parenchymal Tunneling (PSO-PT) for parenchymal ICP monitoring and the Monitoring Kit, Parenchymal Tunneling (PSO-PTT) for parenchymal ICP and ICT monitoring are MR Conditional.

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Under the scan conditions defined above, the Monitoring Kit, Parenchymal Tunneling (PSO-PT) for parenchymal ICP monitoring or the Monitoring Kit, Parenchymal Tunneling (PSO-PTT) for parenchymal ICP and ICT monitoring is expected to produce a maximum temperature rise of less than 2.2 °C after 15 minutes of continuous scanning.

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The "MR Conditional" status only concerns Monitoring Kits for ICP or ICP and ICT. The Monitors, Interface, cables and the temperature module are considered to be **"MR Unsafe"** (not MRI-compatible) and must not be exposed to an MRI environment.

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- The non-implanted part of the Catheter is coiled around the MRI Support, which itself is centered in the MRI. For more information, refer to the indications mentioned in section *Positioning of the Monitoring Kits and the MRI Support* below.
- The MRI has a static magnetic field less than or equal to 3 Tesla.
- The spatial gradient of the static magnetic field is less than or equal to 19 T/m.
- The whole body averaged SAR (Specific Absorption Rate) is less than or equal to 2 W/kg (normal mode) for 15 minutes continuous exposure at 1.5 Tesla.
- The whole body averaged SAR (Specific Absorption Rate) is less than or equal to 2 W/kg (normal mode) for 15 minutes continuous exposure at 3 Tesla.
- No limitation in the magnetic field gradients.
- The size of the artifact on the images may reach about 5.5 cm depending on the implantation depth for MRI examinations performed in 3 Tesla echo sequence.
- Do not use an emission/reception head coil (RF emission coil for the head) or a simple RF emission head coil. Only use an RF emission/reception whole body coil or an RF

emission whole body coil with a simple RF reception head coil.

- Do not scan a patient with a high body temperature.

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
The *in vitro* tests conducted in accordance with the ASTM F2213 standard have shown that the torque forces induced by a spatial gradient of the magnetic field less than or equal to 19 T/m do not present any risks for the patient.

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Sophysa warranties that the Monitoring Kits are free of any material and manufacturing defects. Apart from this warranty, Sophysa does not provide any other warranty, express

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12. Symbols

	The Catheters are not protected against the effects of a cardiac defibrillator.
	Refer to the instructions for use
	Do not re-sterilize
	Do not re-use
	MR Conditional. It has been demonstrated that the use of ICP or ICP and ICT Monitoring Kits within an MRI setting does not pose any risks as long as the specific conditions of use are complied with.
	Serial number
	Lot number

	Catalog reference
	Sterilization Method using Ethylene Oxide
	Manufacturer
	Use until
	BF TYPE EQUIPMENT: Providing an appropriate degree of protection against electric shock, having a Type BF insulated applied section (floating).
	Temperature conditions for storage and transport

	Keep away from liquids
	Fragile, handle with care

	CE conformity marking
---	-----------------------

13. Technical Specifications

Table 1. General specifications

Item	Specification
Sensor capsule diameter	1.2 mm
Catheter diameter	2.1 Fr (0.7 mm) $D(mm) = Fr / 3$ $Fr = D(mm) \times 3$
Sensor capsule material	Titanium
Catheter length	1,000 mm (radiopaque on the entire length)
Catheter sheath material	Polyamide
Depth markers	Marking every centimeter between 1 and 10 cm from the proximal end. Thicker marking at 4 cm. Additional marking at 15 cm.
Leakage current	Less than 10 μA to 120 V AC
Weight	11 g
Drill bit diameter	3.5 mm
Note	
The Monitoring Kits are latex- and phthalate-free.	

Table 2. Pressure sensor specifications

Item	Specification
Type of sensor	Piezoresistive with strain gauge type on silicon
Reference pressure	Atmospheric pressure
Display range (complete system)	-40 mmHg to +100 mmHg (PSO-3000) -40 mmHg to +150 mmHg (PSO-4000)
Accuracy of the pressure measurement (linearity and hysteresis) (complete system)	$\pm 2\%$ reading in the range 0 to 100 mmHg
Bandwidth	> 100 Hz
Functional range of over-pressure without damage	-700 mmHg to +1,250 mmHg
Temperature coefficient	0.1 mmHg/°C max
Input resistance	667 Ohms
Output resistance	810 Ohms
Excitation voltage	1 to 8 V AC or DC

Item	Specification
Drift from zero	Less than 1 mmHg during the first day, at 37°C Less than 2 mmHg during the first week, at 37°C

Table 3. Temperature sensor specifications

Item	Specification
Type of sensor	Thermistor
Unit	Celsius (°C) or Fahrenheit (°F) $^{\circ}F = ^{\circ}C \times (9/5) + 32$ $^{\circ}C = (^{\circ}F - 32) \times (5/9)$ Equivalence in Kelvin (K): $K = ^{\circ}C + 273.15$ $K = ^{\circ}F \times (5/9) + 255.37$
Display range (complete system)	+20 °C to +45 °C (+68 °F to +113 °F)
Maximum tolerated error (accuracy) in the reference conditions (complete system)	± 0.2 °C max. from 25 °C to 45 °C (± 0.4 °F max. from 77 °F to 113 °F) ± 0.4 °C max. from 20 °C to 25 °C (± 0.7 °F max. from 68 °F to 77 °F)
Resolution	0.1 °C/0.1 °F

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Item	Specification	
Temperature	Operation	+10 °C (50 °F) to +40 °C (104 °F)
	Storage, transport	0 °C (32 °F) to +50 °C (122 °F)
Humidity	Operation	30% to 75% relative humidity
	Storage, transport	20% to 95% relative humidity
Atmospheric pressure	Operation	500 hPa to 1,060 hPa
	Storage, transport	500 hPa to 1,060 hPa

Table 5. Catheter history specifications

Item	Specification
Data accessible	First 15 days of monitoring

Item	Specification
ICP data	A value every 20 seconds. This value is an average of the previous 20 seconds.
ICT data	A value every 60 seconds. This value is an average of the previous 60 seconds.

Item	Specification
Note	Data can only be consulted with a Pressio 2 ICP Monitor.

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Table 8. Pressio® Accessories

PSO-EC20	Catheter extension cable <i>Use only with a Pressio® ICP Monitor (PSO-3000) or a Pressio® Interface (PSO-IN00)</i>
PSO-EC30	Catheter extension cable <i>Use only with a Pressio® 2 ICP Monitor (PSO-4000)</i>

PSO-MCxx	Pressure cable <i>Use only with a Pressio® 2 ICP Monitor (PSO-4000)</i> <ul style="list-style-type: none"> – MC01: PHILIPS (AGILENT) - 12 pins – MC02: SIEMENS (SIRECUST) - 10 pins – MC03: SPACELABS & MINDRAY - 6 pins – MC04: GE DATEX-Ohmeda - 10 pins – MC05: GE Solar (MARQUETTE) - 11 pins – MC08: NIHON KOHDEN - 5 pins
PSO-MCT-Y	Temperature cable <i>Use only with a Pressio® 2 ICP Monitor (PSO-4000)</i> <ul style="list-style-type: none"> – MCT-A: PHILIPS (AGILENT) - 2 pins – MCT-B: SIEMENS - 7 pins – MCT-C: SPACELABS - 10 pins – MCT-E: GE Solar (MARQUETTE), GE DATEX-Ohmeda - 11 pins – MCT-F: HELLIGE, DATEX-Ohmeda, NIHON KOHDEN, MINDRAY & DATASCOPE - JACK 6.35 mm
PSO-MT00	Intracranial temperature module <i>Enables temperature values to be displayed on the Pressio® ICP Monitor (PSO-3000).</i>
PSO-TX00	Standard transmitter <i>Enables the pressure and temperature value to be transmitted to a computer. Use only with a Pressio® ICP Monitor (PSO-3000).</i>
PSO-DR	Single use, disposable drill
PSO-MRI	Pressio® MRI Support <i>For positioning Pressio® Catheter during MRI examination</i>

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contact@sophysa.com

Sophysa USA

503 E Summit Street, Suite 5
Crown Point, IN 46307
USA

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This page is intentionally left blank.

Table of Contents

Pressio® Monitoring kit, Ventricular tunneling	9
1. Figures	9
2. Indications	9
3. Contraindications	9
4. Description	9
4.1. Catheters	9
4.2. Accessories	10
4.2.1. Drill bit	10
4.2.2. Tunneling trocar	10
4.2.3. Fixation tab	10
4.2.4. Disposable drill	11
5. Sterilisation - Decontamination	11
6. Setting-up the Monitoring System	11
7. Preparing, Implanting and Monitoring	11
7.1. Preparing the Catheter Implantation Site	11
7.1.1. Choice of the Implantation Area	11
7.1.2. Intracranial Access	11
7.2. Zeroing the Catheter	12
7.2.1. Prerequisites	12
7.2.2. Procedure	12
7.3. Implanting the Catheter	13
7.3.1. Prerequisites	13
7.3.2. Procedure	13
7.3.3. Tunneling the Catheter	13
7.3.4. Implanting the Catheter	13
7.3.5. Connection to an External Drainage System	14
7.3.6. Fixing the Catheter	14
7.4. Monitoring, Caring and Nursing of the Patient	15
7.5. Explanting the Catheter	15
8. Behaviour During an MRI Examination	15
8.1. General Points	15
8.2. Positioning of the Monitoring Kits and the MRI Support	16
8.3. In Vitro Test Conclusions	16
9. Complications / Side Effects	17
9.1. Infection	17
9.2. Cerebral Hemorrhage	17
9.3. Drainage Catheter Obstruction	17
9.4. Overdrainage	17
10. Destruction or return of product	17
10.1. Destruction After Use	17
10.2. Return of Products	17
11. Warranty	17
12. Symbols	17
13. Technical Specifications	18
14. References	19
Kit Pressio® pour Monitoring, Ventriculaire Tunnelisable	21
1. Figures	21
2. Indications	21
3. Contre-indications	21
4. Description	21
4.1. Cathéters	22
4.2. Accessoires	22

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12. Símbolos 66

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14. Referências 67

This page is intentionally left blank.

WARNING

Read the Instructions for Use carefully before using the Pressio® Monitoring Kit, Ventricular Tunneling (PSO-VT or PSO-VTT).


Read the Instructions for Use before using the Pressio® ICP Monitor (PSO-3000), the Pressio® 2 ICP Monitor (PSO-4000), or the Pressio® Interface (PSO-IN00).

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.


1. Figures

Figure 1. Pressio® ICP Monitoring Kit, Ventricular Tunneling (Model PSO-VT)

A: Pressio® Catheter for Intracranial Pressure Measurement With Its Dongle B: Luer-Lock Connector C: Stylet D: Fixation Tab E: Allen Key F: Drill Bit With Adjustable Stop G: Tunneling Trocar With Sheath

Figure 2. Pressio® ICP and ICT Monitoring Kit, Ventricular Tunneling (Model PSO-VTT)

A: Pressio® Catheter for Intracranial Pressure and Temperature Measurement With Its Dongle B: Luer-Lock Connector C: Stylet D: Fixation Tab E: Allen Key F: Drill Bit With Adjustable Stop G: Tunneling Trocar With Sheath

2. Indications

The use of a Pressio® monitoring kit is indicated in patients requiring continuous monitoring of intracranial pressure (ICP).

It is also indicated for use on brain-injured patients with a risk of intracranial hypertension, who cannot be evaluated clinically.

Furthermore, the PSO-VTT model makes it possible to measure intracranial temperature (ICT).

NOTE

The intracranial temperature and pressure values indicated on the Pressio® monitoring systems do not pre-judge the condition of the patient.

3. Contraindications

Contraindications to the use of a Monitoring Kit are the following:

- Established or suspected infections in the tissues in direct contact with Monitoring System components (meningitis, ventriculitis, septicaemia or bacteraemia), or any infection present in any part whatsoever of the body.
- Patients on anticoagulant therapy or presenting with bleeding diathesis.

- Ventricular tunneling is contraindicated in children less than one year of age.

WARNING

Do not use the Monitoring Kit if there are no trained personnel available to provide continuous observation.

4. Description

The Monitoring Kits are made up of:

- a Pressio® Catheter (hereinafter referred to as Catheter),
- accessories allowing the Catheter to be implanted and kept in position.

4.1. Catheters

Catheters are single-use implantable devices intended to be used only with a Pressio Monitoring System (hereinafter referred to as Monitoring System):

- a Pressio ICP Monitor PSO-3000 or PSO-4000 (hereinafter referred to as Monitor),

- a Pressio Interface PSO-IN00 (hereinafter referred to as Interface).

The Catheters are of the “BF” type .

WARNING

Catheters are not protected against defibrillation and may be damaged as a result.

Before defibrillation, withdraw the Catheter. If this is not possible, for safety reasons, change the Catheter after defibrillation to continue monitoring.

The Monitoring Kit, **Ventricular Tunneling PSO-VT (Figure 1)** allows continuous monitoring of ICP.

The Monitoring Kit, **Ventricular Tunneling PSO-VTT (Figure 2)** allows continuous monitoring of ICP and ICT.

Both kits allow simultaneous drainage of the cerebrospinal fluid.

The Catheter, designed to be placed in a cerebral ventricle, has three distinct lumina:

- The first contains a Catheter equipped with a capsule at its end, enclosing a pressure sensor (and a temperature sensor for the PSO-VTT model).
- The second contains the pre-inserted introduction stylet.
- The third is for the drainage of the cerebrospinal fluid.

The two kits can be distinguished visually by their respective arrows on the Catheter dangle. The presence of a red arrow indicates that it is a Catheter for monitoring both ICP and ICT.

In both kits, the Catheter has a marking every centimeter between 5 and 10 cm from the proximal end, and a marker at 15 cm.

The external drainage Catheter has a Luer-Lock connector at its distal end, for connection to an external drainage system.

The end of the Catheter, opposite the sensor, has a connector to connect the Catheter to a Monitoring System via a Catheter extension cable:

- Catheter extension cable PSO-EC20 to connect a Catheter to a Pressio ICP Monitor or a Pressio Interface.
- Catheter extension cable PSO-EC30 to connect a Catheter to a Pressio 2 ICP Monitor.

The Catheter extension cable transmits the measurement signals from the sensor(s) to the Monitoring System, in the form of analog signals.

The Catheter also exchanges a digital signal with the Monitoring System. This signal contains data, such as sensor calibration, or value and date on which the Catheter was zeroed.

Data is memorised in the Catheter dangle and is thus independent of the Monitoring System used.

NOTE

If a Pressio 2 ICP Monitor is used, the ICP or ICP + ICT data for the first fifteen days of monitoring are recorded in the Catheter memory.

Beyond 15 days, no more new data is recorded in the Catheter memory. This data may be consulted only with a Pressio 2 ICP Monitor by using the “History” function.

Specific case for the PSO-VTT

To monitor both ICP and ICT, the Monitoring System must either be a Pressio ICP monitor, software version V2, or a Pressio 2 ICP Monitor.

NOTE

The software version can be identified when powering on the Monitor: the Sophysa logo and software version display on the Monitor touch screen.

NOTE


If you do not have a Pressio 2 ICP Monitor or a Pressio ICP monitor, software version V2, you can still use the Catheter on one of the other Monitoring Systems (Pressio ICP monitor, software version V1 or Pressio Interface). In this case, only the pressure will be displayed, the temperature will not be displayed. The performances of the pressure sensor will be identical.

4.2. Accessories

The different accessories described below enable the Catheters to be implanted.

4.2.1. Drill bit

- [1] Drill bit for intracranial access.
- [2] Adjustable stop and its locking screw to set a drilling depth.
- [3] Allen key to set the position of the adjustable stop.

1: Drill bit 2: Adjustable stop 3: Allen key

4.2.2. Tunneling trocar


- [4] Tunneling trocar to create a passage for the Catheter under the scalp
- [5] Sheath to help pass the Catheter under the scalp.

4: Tunneling trocar 5: Sheath

4.2.3. Fixation tab

The silicone fixation tab [6] enables the fixation of the Catheter to the scalp.

6: Fixation tab

4.2.4. Disposable drill

A sterile disposable drill (PSO-DR) to perforate the skull is also available to order. It can be used in combination with the drill bit supplied in each Monitoring Kit.

5. Sterilisation - Decontamination

The Monitoring Kits are packed individually in double peel-off, sterile, pyrogen-free packaging. They are sterilized with ethylene oxide.

WARNING

Do not use the products if the sterile packaging is open or damaged, or if the expiry date has passed.

WARNING

This product is **for single use only**. It is intended to be used once only for a single patient.

WARNING

Do not re-sterilize or re-use after opening the packaging and/or after explantation.

Resterilization may damage the product, which could cause lesions to the patient.

Reusing this device may change its mechanical or biological features and may provoke failure, or the development of allergic reactions or bacterial infections.

NOTE

Sophysa cannot be held responsible for the performance of any product that has been re-sterilized, or for any complications which might result from this.

6. Setting-up the Monitoring System

CAUTION

The Monitoring System should only be used by trained personnel.

- Power on the Monitoring System.
- Connect the Catheter extension cable to the Monitoring System.

7. Preparing, Implanting and Monitoring

WARNING

Read the instructions for use of the Monitoring System in use before using a Monitoring Kit.

CAUTION

The implantation of the Catheter should be performed **immediately** after the Catheter is zeroed. Therefore, prepare the Monitoring System and the implantation site **before** zeroing the Catheter.

CAUTION

Do not perform the implantation of a Catheter without having a replacement Monitoring Kit available in case it is required.

CAUTION

Only connect Catheters to calibrated Monitoring Systems.

7.1.1. Choice of the Implantation Area

As an example, the right and left prefrontal areas are the main implantation areas. These areas enable the patient to turn his/her head from one side to the other whilst remaining in decubitus without interfering with the intracranial pressure surveillance function.

In most cases, the incision is made behind the hairline, which is acceptable from an aesthetic point of view.

7.1.2. Intracranial Access

Once the implantation site has been chosen, the area is shaved and prepared aseptically. A local anesthetic is applied in the incision area, generally 2 to 3 centimeters in front of the coronal suture on the mid-pupillary line.

1. Make an incision of about 1 cm down to the bone.
2. Expose the bony plate and perform haemostasis on the wound edges.
3. Loosen the locking screw with the Allen key to set the position of the adjustable stop on the drill bit, depending upon the drilling depth chosen.
4. Once the adjustable stop is correctly positioned, retighten the locking screw.

WARNING

Incomplete tightening of the locking screw will prevent the adjustable stop from playing its role, with the risk of drilling too deeply.

5. Fix the helical drill bit to a drill.

7.1. Preparing the Catheter Implantation Site

Observe aseptic neurosurgical techniques to prepare the Catheter implantation site.

Surgeons will perform the most appropriate technique depending upon their experience and the clinical status of the patient.

Use the insertion accessories supplied by Sophysa in the Monitoring Kit.

6. Perforate the internal and external skull plates avoiding any possibility of a parenchymal injury when passing through the internal plate.
7. Having crossed the internal plate, remove the drill bit.
8. Irrigate the burr hole with sterile saline solution and remove the debris.
9. Pierce the dura mater.

7.2. Zeroing the Catheter

Zeroing the Catheter involves calibrating the Catheter in relation to atmospheric pressure.

WARNING

Perform the zeroing procedure on each new Catheter, **before** implanting it in the patient.

CAUTION

Observe aseptic technique throughout the entire zeroing procedure.

NOTICE

Do not allow the Catheter extension cable connector to come into contact with liquid.

NOTE

The ICT sensor is calibrated in the factory. Therefore, the temperature setting does not need prior zeroing.

7.2.1. Prerequisites

- Power on the Monitoring System.
- Connect the Catheter extension cable to the Monitoring System:
 - if using a Monitor, a message inviting you to connect the Catheter displays on the Monitor screen,
 - if using an Interface, the zeroing diode flashes slowly.


7.2.2. Procedure

1. Unpack the Catheter (sterile) within sterile field.
2. Prepare a shallow cup of sterile saline solution (less than 5 mm).

CAUTION

Using a deep receptacle could lead to hydrostatic pressure greater than atmospheric zero, which would then lead to an erroneous reference zero.

3. Observing aseptic technique, connect the Catheter extension cable (non-sterile) to the Catheter, aligning the blue arrow on the Catheter dangle [7] with the blue arrow on the Catheter extension cable connector [8].

The Catheter dangle must be completely inserted into the Catheter extension cable connector. A blue line on the Catheter dangle indicates the point of complete insertion.

When correctly inserted:


- if using a Monitor, a message displays on the screen,
- if using an Interface, the zeroing diode flashes quickly.

NOTE

The appearance of error codes, “E001”, “E002 or “E005” on the Monitor after connecting the Catheter to the Catheter extension cable may be a sign of an incorrect connection.

In this case, make sure that the Catheter dangle is pushed right up to the stop in the Catheter extension cable, and that the latter is correctly connected to the Monitor.

4. Within the sterile field, immerse the sensor (metal tip) of the Catheter [9] in the cup of sterile saline solution [10], without touching the cup.

CAUTION

Avoid any contact with the sensor during zeroing procedure, as this could adversely affect the calibration, resulting in an inaccurate zero reference.

Keep the sensor immersed in sterile saline solution during the zeroing procedure. Zeroing the sensor in conditions other than those recommended may cause inaccurate pressure readings.

5. While the sensor is immersed, press the Zero button  on the Monitor touch screen, or on the Interface: Zeroing the Catheter requires approximately 4 seconds. Do not move the Catheter during the zeroing procedure.
 - If the zeroing procedure is successful:
 - the Monitor touch screen displays a message indicating that the Catheter is ready to be implanted,
 - the Interface zeroing diode remains lit.
 - If the zeroing procedure is not successful:
 - follow the instructions displayed on the Monitor touch screen,
 - the Interface emits a series of beeps and sends a pressure value of 360 mmHg to the patient bedside monitor. For more information, refer to the instructions for use of the Interface.

NOTE

Zeroing the Catheter can only be performed once. Upon completion, the zero calibration information is stored in the Catheter dangle. It allows the Catheter to be disconnected from a Monitoring System, and reconnected to any Monitoring System, without losing the zero calibration information.

7.3. Implanting the Catheter

Before implanting the Catheter, check that the screen of the Monitor (or of the patient bedside monitor if the Interface is used) displays the value "0 mmHg".

WARNING

Limit repeated intracerebral implantation of Catheters. Frequent perforation of the brain to allow Catheter insertion may predispose the brain to oedema and intracerebral haemorrhage resulting in an increase in intracranial pressure.

CAUTION

Observe aseptic neurosurgical technique when implanting a Catheter.

Surgeons will perform the most appropriate technique depending upon their experience and the patient clinical status.

The final implantation of the device must satisfy the conditions for optimal positioning of the sensor in the ventricle.

7.3.1. Prerequisites

The Catheter must be zeroed.

7.3.2. Procedure

Use the insertion accessories supplied by Sophysa in the Monitoring Kit.

NOTE


The Catheter must be tunneled under the scalp to improve its fixation and reduce the risks of infection. The Catheter usually emerges 5 cm from the burr hole in a posterior position.

7.3.3. Tunneling the Catheter


CAUTION

The end of the tunneling trocar is sharp, introduce the Catheter very carefully.


1. Insert the proximal end of the Catheter in the tunneling trocar sheath.


2. Make a small incision at the site chosen for the emergence of the Catheter (incision a).

3. Starting at the emergence site, insert the tunneling trocar between the scalp and the skull towards the burr hole.


4. Pull the Catheter out of the tunnel making sure that a length of at least 30 cm comes out of the implantation site, and remove the tunneling trocar with the sheath.


7.3.4. Implanting the Catheter

1. Holding it by the stylet pre-inserted in the specifically designed lumen, implant the zeroed Catheter in the direction of the ventricle in accordance with standard techniques.

2. Remove the protective cap from the Luer-Lock connector at the end of the Catheter drainage tubing.
3. Check correct positioning in the ventricle by observing the return of the cerebrospinal fluid at the drainage connector opening.


- If a series of air bubbles and liquid segments appear in the Catheter, purge the Catheter letting possible air bubbles move through to the end of the drainage tubing.
- If the ventricles of the patient are swollen, move the Catheter forward by several millimeters beyond the point where the first sample of fluid was taken. This way, the Catheter end will remain in the ventricle during decompression.

NOTE

If the ventricular approach fails, it is still possible to monitor ICP by leaving the Catheter in place in the parenchyma it has crossed. The ICP values measured will be the intraparenchymal pressure values on the sensor.

In this case, close the Luer-Lock connector to limit the risks of infection.

4. While holding the Catheter in place at the implantation site, remove the stylet by gently bending the Catheter at the two-line marker level.

7.3.5. Connection to an External Drainage System


- The Catheter may be connected to different external CSF drainage systems. Connect the Catheter to the external drainage tubing using the Luer-Lock connector.
- If the Catheter is not connected to an external drainage system, close the Luer-Lock connector to limit the risk of infection.

NOTE

Refer to the external drainage system Instructions for Use.

7.3.6. Fixing the Catheter


1. Hold the Catheter in place at the implantation site and pull very gently on the end located at the side of the connector until it forms a right angle and rests flat against the skull.

2. Close the incision above the burr hole in compliance with standard hospital procedures.
3. Use the fixation tab supplied to fix the Catheter onto the scalp at the emergence site.
4. To keep the Catheter in place and reduce the tension, roll the Catheter and attach the loop thus formed. Make sure no traction is exerted on the fixation tab during these stages.

WARNING

Make sure the loop formed by the Catheter is not too tight (radius of more than 15 mm) or the drainage of the CSF may be slowed.

Once the Catheter is implanted, the Monitor (or patient bedside monitor) displays the mean ICP value in millimeters of mercury (mmHg).

1 mmHg corresponds to 13.60 mmH₂O and to 133 Pa.

NOTE

The appearance of error codes “- -”, “999”, “E001”, “E002”, “E005”, or the message “Connect sensor” on the Monitor after the Catheter implantation, when the Catheter is correctly connected may be a sign of damage to the dangle or to the sensor located on the tip of the Catheter.

Specific case for the PSO-VTT

The temperature is also displayed and will stabilise at its accurate value within a maximum of 150 seconds.

NOTE

A temperature measurement can only be accurate if the sensor is implanted. The sensor is not suitable for measuring the temperature in air.

When the temperature read by the sensor is less than 20 °C or higher than 45 °C, the Monitor displays: “- - -”. This means that the measurement is in the sensor measuring range but outside the display ranges of the Monitor.

7.4. Monitoring, Caring and Nursing of the Patient

Once the Catheter is correctly implanted and connected to the Monitoring System, take into account the precautions listed below to ensure optimum monitoring.

CAUTION

The patient should only be monitored by trained and qualified personnel.

CAUTION

Attach the Catheter extension cable fixation clips to the bed sheets, or to the patient's clothing to limit traction on the implanted Catheter, and reduce the risk of Catheter disconnection.

CAUTION

Move the patient with care in order to avoid disconnecting any cables, or any movement of the implanted Catheter.

After moving the patient, check the connection of the Catheter to the Catheter extension cable, and the connection of the Catheter extension cable to the Monitoring System.

CAUTION

Do not use a Monitoring System and the implanted Catheter at the same time as a high frequency electro-surgical instrument or a defibrillator. The Catheter and/or the Monitoring System could be damaged or have their operation disrupted.

CAUTION

Do not allow any of the connectors to come into contact with liquid, especially during care of the patient.

CAUTION

During monitoring, observe the Catheter implantation site in compliance with standard hospital procedures.

After the first 24 hours, the Monitor displays the duration of the Catheter implantation in the form of the message "**Implantation: X days**", at the top left corner of the touch screen, for the entire duration of its implantation.

CAUTION

Catheters are recommended to be implanted for up to 6 days (144 hours).

When this duration is exceeded:

- on a Monitor, the message "**Implantation: X days**" flashes at the top left corner of the touch screen of the Monitor,
- on the Interface, the diode showing the implantation duration for the Catheter flashes.

The accuracy of the displayed ICP value is no longer guaranteed.

NOTE

If a Pressio 2 ICP Monitor is used, it is possible to access the records for the last 15 days of monitoring by using the "History" function on the Monitor.

7.5. Explanting the Catheter

1. Stop the Monitor, or disconnect the Interface.
2. Disconnect the Catheter from the Catheter extension cable.
3. Carefully explant the Catheter.
4. Check the integrity of the explanted Catheter.

8. Behaviour During an MRI Examination

MRI Safety Information

Non-clinical testing has demonstrated the Monitoring Kit, Ventricular Tunneling (PSO-VT) for ventricular ICP monitoring and the Monitoring Kit, Ventricular Tunneling (PSO-VTT) for ventricular ICP and ICT monitoring are MR Conditional.

A patient with these devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T.
- Maximum spatial field gradient of 1,900 gauss/cm (19 T/m) (extrapolated).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode).
- The non-implanted part of the catheter must be wrapped around the Pressio® MRI Support (hereinafter referred to as MRI Support) and the MRI Support must be centered in the scanner bore. Refer to the information in section *Positioning of the Monitoring Kits for the MRI examination*.

Under the scan conditions defined above, the Monitoring Kit, Ventricular Tunneling (PSO-VT) for ventricular ICP monitoring or the Monitoring Kit, Ventricular Tunneling (PSO-VTT) for ventricular ICP and ICT monitoring is expected to produce

a maximum temperature rise of less than 2.2 °C after 15 minutes of continuous scanning.

WARNING

Disconnect the Catheter from the Monitor before any MRI examination.

The "MR Conditional" status only concerns Monitoring Kits for ICP or ICP and ICT. The Monitors, Interface, cables and the temperature module are considered to be "**MR Unsafe**" (not MRI-compatible) and must not be exposed to an MRI environment.

NOTE

The ICP or ICP and ICT Monitoring Kits are considered to be "MR Conditional" according to the definition of the ASTM F2503 standard.

8.1. General Points

A patient implanted with an ICP or ICP and ICT Monitoring Kit may undergo an MRI examination after the implantation of the ICP or ICP and ICT Monitoring Kit.

The results of *in-vitro* tests on ICP or ICP and ICT Monitoring Kit have shown that the ICP or ICP and ICT Monitoring Kits do not present any danger for patients in the following conditions:

- The non-implanted part of the Catheter is coiled around the MRI Support, which itself is centered in the MRI. For more information, refer to the indications mentioned in section *Positioning of the Monitoring Kits and the MRI Support* below.
- The MRI has a static magnetic field less than or equal to 3 Tesla.
- The spatial gradient of the static magnetic field is less than or equal to 19 T/m.
- The whole body averaged SAR (Specific Absorption Rate) is less than or equal to 2 W/kg (normal mode) for 15 minutes continuous exposure at 3 Tesla.
- The whole body averaged SAR (Specific Absorption Rate) is less than or equal to 2 W/kg (normal mode) for 15 minutes continuous exposure at 1.5 Tesla.
- No limitation in the magnetic field gradients.
- The size of the artifact on the images may reach about 5.5 cm depending on the implantation depth for MRI examinations performed in 3 Tesla echo sequence.
- Do not use an emission/reception head coil (RF emission coil for the head) or a simple RF emission head coil. Only use an RF emission/reception whole body coil or an RF emission whole body coil with a simple RF reception head coil.
- Do not scan a patient with a high body temperature.

8.2. Positioning of the Monitoring Kits and the MRI Support

WARNING

Read the Instructions for Use before using the MRI Support.

WARNING


The *MR conditional* status of the ICP or ICP and ICT Monitoring Kits is valid if, and only if, the MRI Support is used. Do not use another support for the coiling of the ICP or ICP and ICT Monitoring Kits.

WARNING


Do not perform an MRI leaving the Catheter uncoiled or in a straight line along the patient's body. This configuration could cause serious injuries to the patient.

Once the patient is installed on the MRI bed:

1. Place the MRI support on the bed, the base of the support in contact with the bed.
2. Place the support as close as possible to the patient's head.

3. Move the MRI Support so that it is in the **center** of the bed.


4. Pass the non-implanted part of the Catheter through one of the side notches on the body of the MRI Support, and coil it. Depending on the Catheter length, it may be coiled 4 or 5 times round the MRI Support.

CAUTION

Handle the Catheter with care when coiling it to avoid any excessive traction on the implanted Catheter.


CAUTION

When handling the Catheter, avoid any liquid spray on the connector as this may damage it.

5. Fix the Catheter dangle by inserting it, as shown below, between the flanges of the MRI Support spool.

click

6. Before starting the MRI examination, make sure that the MRI Support is correctly centered. If not centered, adjust its position.

8.3. In Vitro Test Conclusions

The *in vitro* tests conducted in accordance with the ASTM F2182 standard have shown that in the specific conditions mentioned previously, the increase in temperature resulting from exposure to an MRI examination at 1.5 and 3 Tesla is negligible and has no clinical consequences for the patient (maximum increase measured + 2°C +/- 0.3°C).

The *in vitro* tests conducted in accordance with the ASTM F2052 standard have shown that the displacement forces induced by a magnetic field of 3 Tesla or less do not present any risks for the patient.

The *in vitro* tests conducted in accordance with the ASTM F2213 standard have shown that the torque forces induced

by a spatial gradient of the magnetic field less than or equal to 19 T/m do not present any risks for the patient.

The size of the artifacts has been evaluated in accordance with the ASTM F2119 standard.

NOTE

The *in vitro* tests have shown that the operation of the ICP or ICP and ICT Monitoring Kits is not affected by exposure to a magnetic field of 1.5 to 3 Tesla (30,000 Gauss).

9. Complications / Side Effects

Complications, which may result from the implantation of a Catheter, include the inherent risks in any surgical procedure and the insertion of a foreign body.

These complications require the rapid intervention of a physician.

9.1. Infection

An infection is the major complication associated with this type of monitoring.

Risks of infection can be reduced by observing aseptic techniques when handling and implanting the Catheter, and by respecting the maximum duration for the Catheter implantation (6 days or 144 hours). If monitoring must be continued after this duration, implant a new Catheter on another site.

WARNING

In case of infection, remove the Catheter and start a specific treatment by a general or intrathecal route.

9.2. Cerebral Hemorrhage

Cerebral hemorrhage may also be observed with this type of monitoring.

Risks of hemorrhage may be reduced by limiting the number of cerebral incisions during the introduction procedure, and by ensuring that this procedure is only performed by trained and qualified personnel.

9.3. Drainage Catheter Obstruction

Obstruction is a complication specific to the use of an intra-ventricular ICP monitoring Catheter.

The implantation of the Catheter end close to the choroid plexus may cause obstruction of the CSF inlet holes. This phenomenon can also occur when the CSF has a high protein content, contains blood, or if there is intra-ventricular debris such as blood clots or tissues.

This obstruction may have an impact on CSF drainage but will not have any effect on the reliability on the ICP or ICT measurement.

9.4. Overdrainage

Overdrainage is a complication specific to the use of an intraventricular ICP monitoring Catheter.

Overdrainage can result in a collapse of the ventricles and the appearance of a subdural hematoma.

10. Destruction or return of product

10.1. Destruction After Use

An unpacked, used or explanted Monitoring Kit should be destroyed in accordance with the procedures in force in the medical establishment.

10.2. Return of Products

NOTE

For more effective analysis, avoid cleaning the device.

If an explanted Monitoring Kit needs to be returned to Sophysa for analysis, indicate whether cleaning has been performed.

In order to assess the returned product properly, it must be accompanied by an explanatory Return to Manufacturer Authorization form.

11. Warranty

The performance of the Monitoring Kits is only warranted with Monitoring Systems and accessories designed, tested and manufactured by Sophysa.

Sophysa warrants that the Monitoring Kits are free of any material and manufacturing defects. Apart from this warranty, Sophysa does not provide any other warranty, express

or implicit, including commercialization or adaptation for a particular use. Sophysa cannot be held responsible for any incident, complication, damage or prejudice occurring directly or indirectly from the use of this device. Sophysa does not authorize anyone whomsoever to take responsibility on its behalf for its products.

12. Symbols

The Catheters are not protected against the effects of a cardiac defibrillator.

Refer to the instructions for use

	Do not re-sterilize
	Do not re-use
	MR Conditional. It has been demonstrated that the use of ICP or ICP and ICT Monitoring Kits within an MRI setting does not pose any risks as long as the specific conditions of use are complied with.
	Serial number
	Lot number
	Catalog reference
	Sterilization Method using Ethylene Oxide
	Manufacturer

	Use until
	BF TYPE EQUIPMENT: Providing an appropriate degree of protection against electric shock, having a Type BF insulated applied section (floating).
	Temperature conditions for storage and transport
	Keep away from liquids
	Fragile, handle with care
	CE conformity marking

13. Technical Specifications

Table 1. General specifications

Item	Specification
Sensor capsule diameter	1.2 mm
Catheter diameter	2.1 Fr (0.7 mm) $D(mm) = Fr / 3$ $Fr = D(mm) \times 3$
Diameter of the catheter with 3 lumina	9.3 Fr (3.1 mm)
Sensor capsule material	Titanium
Catheter length	1,000 mm (radiopaque on the entire length)
Drainage catheter length	700 mm
Catheter sheath material	Polyamide
External drainage catheter sheath material	Silicone
Depth markers	Marking every centimeter between 5 and 10 cm from the proximal end. Additional marking at 15 cm.
Leakage current	Less than 10 μA to 120 V AC
Weight	11 g
Drill bit diameter	3.5 mm
Note	
The Monitoring Kits are latex- and phthalate-free.	

Table 2. Pressure sensor specifications

Item	Specification
Type of sensor	Piezoresistive with strain gauge type on silicon
Reference pressure	Atmospheric pressure
Display range (complete system)	-40 mmHg to +100 mmHg (PSO-3000) -40 mmHg to +150 mmHg (PSO-4000)
Accuracy of the pressure measurement (linearity and hysteresis) (complete system)	$\pm 2\%$ reading in the range 0 to 100 mmHg
Bandwidth	> 100 Hz
Functional range of over-pressure without damage	-700 mmHg to +1,250 mmHg
Temperature coefficient	0.1 mmHg/°C max
Input resistance	667 Ohms
Output resistance	810 Ohms
Excitation voltage	1 to 8 V AC or DC
Drift from zero	Less than 1 mmHg during the first day, at 37°C Less than 2 mmHg during the first week, at 37°C

Table 3. Temperature sensor specifications

Item	Specification
Type of sensor	Thermistor

Item	Specification
Unit	<p>Celsius (°C) or Fahrenheit (°F)</p> <p>°F = °C x (9/5) + 32</p> <p>°C = (°F - 32) x (5/9)</p> <p>Equivalence in Kelvin (K):</p> <p>K = °C + 273.15</p> <p>K = °F x (5/9) + 255.37</p>
Display range (complete system)	+20 °C to +45 °C (+68 °F to +113 °F)
Maximum tolerated error (accuracy) in the reference conditions (complete system)	<p>± 0.2 °C max. from 25 °C to 45 °C (± 0.4 °F max. from 77 °F to 113 °F)</p> <p>± 0.4 °C max. from 20 °C to 25 °C (± 0.7 °F max. from 68 °F to 77 °F)</p>
Resolution	0.1 °C/0.1 °F

Table 4. Environmental specifications

Item	Specification
Temperature	Operation +10 °C (50 °F) to +40 °C (104 °F)

14. References

Table 6. Pressio® Catheter Kits

PSO-PB	Pressio® ICP Monitoring Kit, Parenchymal with Bolt
PSO-PBT	Pressio® ICP and ICT Monitoring Kit, Parenchymal with Bolt
PSO-PT	Pressio® ICP Monitoring Kit, Parenchymal Tunneling
PSO-PTT	Pressio® ICP and ICT Monitoring Kit, Parenchymal Tunneling
PSO-VT	Pressio® ICP Monitoring Kit, Ventricular Tunneling <i>with external CSF drainage function</i>
PSO-VTT	Pressio® ICP and ICT Monitoring Kit, Ventricular Tunneling <i>with external CSF drainage function</i>

Table 7. Pressio® Monitoring System

PSO-3000	Pressio® ICP Monitor <i>Power cable and Catheter extension cable (PSO-EC20) included</i>
PSO-4000	Pressio® 2 ICP Monitor <i>Power cable and Catheter extension cable (PSO-EC30) included</i>
PSO-IN00	Pressio® ICP Monitoring Interface <i>Catheter extension cable (PSO-EC20) included</i>

Item	Specification
	<p>Storage, transport 0 °C (32 °F) to +50 °C (122 °F)</p>
Humidity	<p>Operation 30% to 75% relative humidity</p> <p>Storage, transport 20% to 95% relative humidity</p>
Atmospheric pressure	<p>Operation 500 hPa to 1,060 hPa</p> <p>Storage, transport 500 hPa to 1,060 hPa</p>

Table 5. Catheter history specifications

Item	Specification
Data accessible	First 15 days of monitoring
ICP data	A value every 20 seconds. This value is an average of the previous 20 seconds.
ICT data	A value every 60 seconds. This value is an average of the previous 60 seconds.
<p>Note</p> <p>Data can only be consulted with a Pressio 2 ICP Monitor.</p>	

Table 8. Pressio® Accessories

PSO-EC20	<p>Catheter extension cable</p> <p><i>Use only with a Pressio® ICP Monitor (PSO-3000) or a Pressio® Interface (PSO-IN00)</i></p>
PSO-EC30	<p>Catheter extension cable</p> <p><i>Use only with a Pressio® 2 ICP Monitor (PSO-4000)</i></p>
PSO-MCxx	<p>Pressure cable</p> <p><i>Use only with a Pressio® 2 ICP Monitor (PSO-4000)</i></p> <ul style="list-style-type: none"> - MC01: PHILIPS (AGILENT) - 12 pins - MC02: SIEMENS (SIRECUST) - 10 pins - MC03: SPACELABS & MINDRAY - 6 pins - MC04: GE DATEX-Ohmeda - 10 pins - MC05: GE Solar (MARQUETTE) - 11 pins - MC08: NIHON KOHDEN - 5 pins
PSO-MCT-Y	<p>Temperature cable</p> <p><i>Use only with a Pressio® 2 ICP Monitor (PSO-4000)</i></p> <ul style="list-style-type: none"> - MCT-A: PHILIPS (AGILENT) - 2 pins - MCT-B: SIEMENS - 7 pins - MCT-C: SPACELABS - 10 pins - MCT-E: GE Solar (MARQUETTE), GE DATEX-Ohmeda - 11 pins - MCT-F: HELBIGE, DATEX-Ohmeda, NIHON KOHDEN, MINDRAY & DATASCOPE - JACK 6.35 mm
PSO-MT00	<p>Intracranial temperature module</p> <p><i>Enables temperature values to be displayed on the Pressio® ICP Monitor (PSO-3000).</i></p>

PSO-TX00	Standard transmitter <i>Enables the pressure and temperature value to be transmitted to a computer. Use only with a Pressio® ICP Monitor (PSO-3000).</i>
PSO-DR	Single use, disposable drill
PSO-MRI	Pressio® MRI Support <i>For positioning Pressio® Catheter during MRI examination</i>

Technical specifications and list of product references may be modified without notice. Availability may vary according to country.

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©2024 Sophysa. All rights reserved.