Medtronic

Bio-Medicus™

Femoral Cannula and Introducer

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Kanyle og indfører til lårbenet

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
1 Product description

Femoral cannula and introducer. (All contents are sterile/single-use only.)

These instructions for use apply to both uncoated cannulae and cannulae coated with Cortiva™ BioActive Surface. Cannulae coated with Cortiva bioactive surface include a "CB" prefix in the model number. For more information regarding Cortiva bioactive surface, see the insert "Information for Products with Cortiva BioActive Surface."

2 Indications for use

These devices are to be used by a trained physician only. Cannulae are used to cannulate vessels, perfuse vessels or organs and/or connect with accessory extracorporeal equipment. The cannula introducer is intended to facilitate proper insertion and placement of the appropriate sized cannula within the vessel for cardiopulmonary bypass. These products are intended for use for up to 6 hours.

3 Contraindications

There are no known contraindications for the use of the cannula, other than those generally contraindicated for cardiopulmonary bypass. The cannula introducer is only to be used with the appropriate sized Bio-Medicus™ cannula. Use only as directed.

4 Warnings and precautions

- Read the instructions for use for all related extracorporeal devices prior to use.
- Each cannula and introducer is for single use only. Do not reuse or resterilize. Discard after use.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Do not autoclave.
- Do not use if the sterile package has been damaged or if the seal has been broken in any way.
- Do not use alcohol or alcohol-based fluids for lubrication.
- Attach and manipulate the cannula in such a manner as to prevent kinks or any restrictions that may alter flow such as false lumen cannulation.
- Use caution when attaching the cannula to other devices in the extracorporeal circuit so as not to damage the device being connected or the cannula.
- Clamping too close to the connector may damage the connector end.
- Do not clamp wire-reinforced sections. If wire-reinforced cannulae are used, clamping on reinforced sections may result in permanent wall distortion and/or lumen collapse.
- Because the cannula wall is opaque due to the presence of the spring, be sure to completely walk the air out of the cannula during the priming procedure.
- A cannula inserted into the left ventricle through the left atrium can perforate the left ventricular apex if undue pressure is exerted by the tip.
- Take care to avoid overinsertion of the cannula and possible impingement of the cannula tip against the vessel/ventricular/atrial walls.
- The use of cold solutions and cardioplegic agents, administered externally or via the coronary circulation, may increase the stiffness of urethane plastic and contribute to ventricular perforation under conditions of cannula tip pressure during manipulation of the heart.
- Forcing the introducer hub or hemostasis cap into the straight connector could result in stress cracking of the connector.
- Do not attempt to cross-clamp the cannula with the introducer in place.
- Carefully monitor for signs of occlusion on both the inlet and outlet cannulae.
- The user assumes responsibility for any variations in extracorporeal procedures that could compromise the intended use of this device.
- Thoroughly inspect the cannula prior to use to verify that the lumen and side holes (if applicable) are patent and that the cannula has not been damaged or kinked.
- As with all medical devices, this device is to be used by a trained physician only.
- Dispose of the product according to hospital policy.

5 Adverse effects

This device, as do all extracorporeal blood system devices, has possible adverse effects, which include, but are not limited to, embolic events and dislocation, infections, blood loss, and thrombus formation. Vessel damage and complications at the puncture site may occur if the instructions for use are not followed.

6 Sterilization

All Bio-Medicus cannulae are sterilized (EtO) prior to shipment. The lot number is designated on the outer package label. This device should not be used if the sterile package has been damaged or if the seal has been broken. Resterilization is not recommended under any circumstances as this device is intended for single use only.

7 Instructions for use

1. Sterile cannula should be placed into the vascular system according to standard surgical or percutaneous techniques. Various sizes and shapes of cannulae are available in order to provide the physician the opportunity of selecting the best cannula for each specific procedure. When it is determined that extracorporeal bypass is no longer required, the cannulae are carefully removed from the vessels and their placement sites are repaired according to standard surgical practice.
Note: The introducer has been designed to accommodate a 0.97 mm (0.038 in) diameter guidewire (not supplied with the device).
2. Remove the femoral cannula and introducer from the shipping carton and inspect the protective package for damage.
3. Using aseptic technique, open the protective package containing the cannula and introducer.
Note: Use the Bio-Medicus™ cannula insertion kit when performing percutaneous cannulation.
4. Insert the introducer into the lumen of the cannula through the straight connector.
5. Secure the introducer in place by firmly, but gently, pushing the hub of the Introducer into the straight connector.
Note: If a hemostasis cap is provided, firmly place it onto the straight connector of cannula. When the bold band closest to the introducer hub is in alignment with the edge of the straight connector, the introducer tip and cannula tip are at the same depth. (Refer to Figure 1)
6. When desired cannula placement has been achieved, remove the guidewire (if applicable) and the introducer.
Warning: Ensure the clamp is properly positioned on the cannula tubing between the tubing connector and wire-reinforced sections. Clamping too close to the connector may damage the connector end. Clamping on wire-reinforced sections may result in permanent wall distortion or lumen collapse.
Note: When the bold band closest to the vascular end of the introducer is visible at the edge of the straight connector, the cannula can be safely cross-clamped.
Warning: Do not attempt to cross-clamp the introducer.
7. Slowly fill the remainder of the cannula tube with sterile priming fluid, displacing any air out through the tube connector end before connecting fluid-filled tubing.
8. At the completion of the procedure, dispose of the product according to hospital policy.

8 Disclaimer of warranty (global)

ALTHOUGH THE PRODUCT HAS BEEN DESIGNED, MANUFACTURED, AND TESTED UNDER CAREFULLY CONTROLLED CONDITIONS PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION FOR A VARIETY OF REASONS. MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER THE CLAIM IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. EXCEPT AS OTHERWISE AGREED TO IN WRITING BY AUTHORIZED MEDTRONIC PERSONNEL, NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Authorized Representative in the European Community

Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
+31 45 566 8000

Europe/Middle East/Africa

Medtronic International Trading Sàrl
Route du Molliat 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland
+41 21 802 7000

Australia

Medtronic Australasia Pty Ltd
5 Alma Road
Macquarie Park, NSW 2113
Australia
1800 668 670

Canada

Medtronic of Canada Ltd
99 Hereford Street
Brampton, Ontario L6Y 0R3
Canada
+1 905 460 3800

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Bio-Medicus™

Pediatric Cannula and Introducers

1. Product Description

Pediatric cannula and introducers. Sterile, nonpyrogenic, disposable single use only.

2. Indications for Use

These devices are to be used by a trained physician only.

Cannulae are used to cannulate vessels, perfuse vessels or organs, and/or connect with accessory extracorporeal equipment. The cannula introducer is intended to facilitate proper insertion and placement of the appropriately-sized cannula within the vessel for cardiopulmonary bypass. These products are intended for use up to 6 hours.

3. Contraindications

Alone, the cannula and introducer are not medical treatment devices. The introducer is to be used only with the appropriately-sized Bio-Medicus™ cannula. These devices are not intended for use except as indicated above.

4. Warnings and Precautions

- Read all warnings, precautions, and Instructions for Use carefully prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- Verify Use By date.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Do not autoclave.
- Do not use if the sterile package has been damaged or if the seal has been broken.
- Do not use alcohol or alcohol-based fluids for lubrication.
- Thoroughly inspect the cannula prior to use to verify that the lumen and side holes (if applicable) are patent and that the cannula has not been damaged or kinked.
- Use caution when assembling or manipulating the introducer within the cannula, or moving the suture ring.
- Avoid false lumen cannulation.
- Never advance, position, or place the cannula without the introducer as this may kink and damage the cannula.
- Use care to protect the side hole areas from kinking as they are passed into the vessel.
- Manipulate cannula in such a manner as to prevent kinks or any restrictions that may alter flow.
- Manipulation of the heart with cannula in place may cause flow restrictions.
- Avoid kinking when attaching tubing and when anchoring or untethering the cannula body.
- Caution should be used when attaching the cannula to other devices in the extracorporeal circuit so as not to damage the cannula or the device being connected.
- Clamping too close to connector may damage connector end.
- Do not clamp wire-reinforced sections. Clamping on reinforced sections may result in permanent wall distortion and/or lumen collapse.
- Since the cannula wall is opaque due to the presence of the spring, be sure to completely remove air from the cannula during the priming procedure and before clamping.
- Care must be taken to avoid overinsertion of the cannula and possible impingement of the cannula tip against the vessel/ventricular/atrial walls.
- The use of cold solutions and cardioplegic agents, administered externally or via the coronary circulation, may increase the stiffness of urethane plastic and contribute to ventricular perforation under conditions of cannula tip pressure during manipulation of the heart.
- Do not attempt to cross-clamp the cannula with the introducer in place.
- Carefully monitor for signs of occlusion on both the inlet and outlet cannulae.
- User assumes responsibility for any variations in extracorporeal procedures that could compromise the intended use of this device.
- As with all medical devices, this device is to be used by a trained physician only.
- Initiate adequate anticoagulation therapy prior to patient cannulation.
- Only physicians trained and experienced in using percutaneous cannulation techniques (such as the Seldinger technique) and cardiopulmonary bypass should use this device.
- Ischemia may occur distal to the cannula placement site.
- Cannula tip placement should be confirmed by using fluoroscopy.

- Extreme care must be exercised to ensure that the arterial cannula and venous cannula are not confused while being connected to the perfusion circuit.
- In some cases, suture rings may not be detectable on x-ray due to the variances in anatomical densities and x-ray techniques. Modifying the x-ray angle or distance may aid in detection.
- When utilizing vacuum for augmented drainage, do not exceed -100 mmHg.
- Maximum pressure is 300 mmHg.
- Do not expose to MRI while cannula is inserted in patient.

5. Adverse Effects

This device, as do all extracorporeal blood system devices, has possible side effects which include, but are not limited to, embolic events and dislocation, infections, blood loss, and thrombus formation. Vessel damage and complications at the puncture site may occur if the Instructions for Use are not followed.

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7. Sterilization

All Bio-Medicus cannulae are sterilized (EtO) prior to shipment. The Lot Number is designated on the outer package label. This device should not be used if the sterile package has been damaged or if the seal has been broken. Resterilization is not recommended under any circumstances as this device is intended for single use only.

8. Instructions for Use

Note: Sterile cannula should be placed into the vascular system according to standard surgical or percutaneous technique. Various sizes of dilators and cannulae are available in order to select the best cannula for each specific procedure. When it is determined that extracorporeal bypass is no longer required, the cannulae are carefully removed from the vessels and their placement sites are repaired according to standard surgical practice.

Note: The introducer has been designed to accommodate a 0.025-in (0.635-mm) diameter guidewire (not supplied with the device).

1. Remove the cannula and introducers from the shipping carton and inspect the protective, sterile package for damage.
2. Using aseptic technique, open the protective, sterile package containing the cannula and the introducers.
3. Select either the standard solid core introducer for traditional placement method or the hollow core guidewire-capable introducer for Seldinger technique placement method.

Note: Use a pediatric-sized Bio-Medicus™ insertion kit when performing Seldinger or percutaneous cannulation.

4. Insert the introducer through the hemostasis cap and ensure the introducer hub is fully inserted into the hemostasis cap. Insert the introducer into the lumen of the cannula through the connector (Figure 1).
5. Firmly place hemostasis cap onto the connector of cannula.
6. Insert the cannula using a guidewire placement Seldinger technique. Determine the length of cannula to be introduced and note the insertion markings on the cannula.

Note: Some cannulae are labeled with length markings in 10-cm (3.94-in) increments for depth reference during cannula placement.

Note: Confirm the anatomical location of the guidewire before inserting the cannula. Anatomical location of guidewire is critical before cannula insertion.

Caution: Ensure that the vessel is of adequate size for cannula size selected to allow distal perfusion of the limb when the cannula is in place. Improper cannula size may be difficult to advance. The vessel must be large enough to ensure perfusion and venous return.

Caution: Do not kink guidewire during insertion, vessel dilation, and manipulation or it may be difficult to remove.

7. Position the adjustable radiopaque suture ring to the desired depth.
8. Pass the cannula/introducer assembly over the guidewire, advancing until the guidewire can be grasped beyond the introducer hub.
9. Insert tip of the cannula introducer into the vessel. Grasp the guidewire and hold stationary while advancing the cannula to the desired depth.

10. Blood leakage from the proximal drainage holes of the venous cannulae models may occur during catheter insertion until these holes are intravascular.
Note: Advancing the cannula and introducer with a rotating motion will make cannula placement easier.
Note: When the bold band closest to the introducer hub is in alignment with the edge of the hemostasis cap, the cannula tip and the introducer tip are at the same depth (Refer to Figure 1).
11. When the desired cannula placement has been achieved, remove the introducer (and guidewire if applicable).
Warning: To maintain hemostasis, do not completely remove the introducer or hemostasis cap until the clamp is applied.
Warning: Ensure clamp is properly positioned on cannula tubing between tubing connector and wire-reinforced sections. Clamping too close to the connector may damage the connector end. Clamping on wire-reinforced sections may result in permanent wall distortion and/or lumen collapse.
Warning: When the bold band closest to the vascular end of the introducer is visible at the edge of the hemostasis cap, the cannula can be safely cross-clamped. Do not attempt to cross-clamp the introducer.
Note: Perfusion tubing clamps must be used. Many clamps are too weak to ensure complete clamping.
12. Remove the hemostasis cap exposing the 1/4-in (0.635-cm) barbed connector.
Note: Arterial hemostasis caps are color-coded red. Venous caps are color-coded blue.
13. Slowly fill the remainder of the cannula tube with sterile priming fluid, displacing any air out through the tube connector end before connecting the fluid-filled tubing.
14. Attach the extracorporeal circuit, taking care to connect the venous cannula to the venous return line (inlet to pump) and the arterial cannula to the arterial return line (outlet from membrane oxygenator), ensuring air-free connection. The tubing may be further secured with cable ties at the barbed connector.
Note: Venous cannula can be identified by an arrow, in the clamp area next to the Bio-Medicus cannula logo, indicating flow out of the connector. See Figure 1.
15. Secure the device to the patient following standard medical practice. In order to avoid accidental displacement during extracorporeal circulation, suture may be wrapped around the suture ring to secure the cannula and avoid kinking the body of the device.
16. The retention clip may be used to secure cannula and/or tubing to the patient or drapes. The clip may be initially placed using the adhesive tape. To ensure that the cannula is secure during the procedure, clamp or suture the retention clip through the holes.
Warning: Do not rely on the adhesive alone to secure the clip.
17. When it is determined that extracorporeal bypass is no longer required, clamp the non wire-reinforced tube section, carefully remove the cannulae from the vessels and repair their placement sites according to standard surgical practice.
Warning: This product has been qualified for use up to 6 hours. Suitability for long-term use has not been determined.
18. After use, this product may be a potential biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, regional, and national laws and regulations.

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Disclaimer of Warranty

ALTHOUGH THE PEDIATRIC CANNULA AND INTRODUCERS, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

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Disclaimer of Warranty

ALTHOUGH THE PEDIATRIC CANNULA AND INTRODUCERS, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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www.medtronic.com

+1 763 514 4000

LifeLine Technical Support, 24-hour consultation
service:

1 877 526 7890

Medtronic Perfusion Systems

7611 Northland Drive
Minneapolis, MN 55428
USA

+1 763 391 9000

Customer service and product orders:

1 800 854 3570

www.perfusionsystems.com

Medtronic Mexico S. de R.L. de CV
Av. Paseo Cuapah 10510 El Lago
C.P. 22210 Tijuana, Baja California
Mexico

EC REP

Medtronic B.V.

Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
+31 45 566 8000

Europe/Africa/Middle East Headquarters

Medtronic International Trading Sàrl

Route du Molliat 31

Case Postale 84

CH - 1131 Tolochenaz

Switzerland

www.medtronic.co.uk

+41 21 802 7000

Australia

Medtronic Australasia Pty Ltd

5 Alma Road

Macquarie Park, NSW 2113

Australia

1800 668 670

Asia

Medtronic International Ltd.

49 Changi South Avenue 2

Singapore 486056

+65 6436 5000

Latin America

Medtronic Latin America

3750 NW 87th Avenue

Suite 700

Miami, FL 33178

USA

+1 305 500 9328

Canada

Medtronic of Canada Ltd

99 Hereford Street

Brampton, Ontario L6Y 0R3

Canada

1 800 268 5346

M950500A002

© 2015, 2017 Medtronic
M950500A002 Rev. 1A

master:20160719, alldocMod:20161025, alldocMod2:20161024, alldocStyle_graph:20160504, alldocStyle:20160504,
master_pageset:20160504, languageLookup:20161110, table:20170124, characterJoin:20110211, 3of9:20110211, 3of9_svg:
20110211

Medtronic

Bio-Medicus™

Adult Cannula Kit
Sada kanyly pro dospělé
Kanylesæt til voksne
Kanülenset für Erwachsene
Κιτ κάνουλας ενηλίκων
Kit de cânula para adulto
Aikuisten kanyylipakkaus
Kit de canule pour adulte
Kanüláló szett felnőttek számára
Kit per cannula per adulti
Canuleset voor volwassenen
Kanylesett for voksne
Zestaw kaniuli dla dorosłych
Conjunto de cânula para adulto
Conjunto de cânula para adulto
Набор канюли для взрослых
Kanylset för vuxna
Yetişkin Kanül Kiti

Instructions for Use • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης •
Instrucciones de uso • Käyttöohjeet • Mode d'emploi • Használati útmutató • Istruzioni per l'uso •
Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização •
Instruções de utilização • инструкция по применению • Bruksanvisning • Kullanım Talimatları

 **Caution:** Federal law (USA) restricts this device to sale by or on the order of
a physician.

master:20171221, alldocMod:20171221, alldocMod2:20171221, alldocStyle_graph:20180121, alldocStyle:20171221,
templates_global:20171025, master_pageset:20171011, pagespecs:20171024, body_extent:20171025, languageLookup:
20171025, table:20170615, characterJoin:20170404, 3of9:20110211, 3of9_svg:20110211

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Adult Cannula Kit

1. Product Description

Adult cannula with insertion set. Sterile, nonpyrogenic, disposable, single use only.

2. Indications for Use

These devices are to be used by a trained physician only. Cannulae are used to cannulate vessels, perfuse vessels or organs, and/or connect with accessory extracorporeal equipment. The cannula introducer is intended to facilitate proper insertion and placement of the appropriately-sized cannula within the vessel for cardiopulmonary bypass. The Bio-Medicus cannulae (18 cm [7.09 in] tip-length models) may be used in either the femoral position as an arterial delivery cannula or in the jugular position as a venous return cannula.

This product is intended for use up to 6 hours.

3. Contraindications

Alone, this cannula is not a medical treatment device. The cannula introducer is only to be used with the appropriately-sized Bio-Medicus cannula. Selection of patients as candidates for such procedures is the physician's responsibility. The outcome is dependent on many variables, including patient pathology, surgical procedure, and perfusion procedures. These devices are intended for use only as indicated in these instructions for use. Do not use if the patient has severe peripheral atherosclerosis or severe arterial dissection.

4. Warnings and Precautions

- Read all warnings, precautions, and instructions for use carefully prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- Verify the use-by date.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death. This product has been qualified for use up to 6 hours. Suitability for long-term use has not been determined.
- Do not autoclave.
- Do not use if the sterile package has been damaged or if the seal has been broken in any way.
- Do not use alcohol or alcohol-based fluids for lubrication.
- Thoroughly inspect the cannula prior to use to verify that the lumen and side holes (if applicable) are patent and that the cannula has not been damaged or kinked.
- Use caution when assembling or manipulating the introducer within the cannula, or when moving the suture ring.
- Avoid false lumen cannulation.
- Never advance, position, or place the cannula without the introducer as this may kink and damage the cannula.
- Carefully protect the side-hole areas from kinking as they are passed into the vessel.
- Manipulate cannula in such a manner as to prevent kinks or any restrictions that may alter flow.
- Manipulation of the heart with the cannula in place may cause flow restrictions.
- Avoid kinking when attaching tubing and when anchoring or untethering the cannula body.
- Use caution when attaching the cannula to other devices in the extracorporeal circuit so as not to damage the cannula or the device being connected.
- Clamping too close to connector may damage connector end.
- Do not clamp wire-reinforced sections. Clamping on reinforced sections may result in permanent wall distortion and/or lumen collapse.
- Since the cannula wall is opaque due to the presence of the spring, be sure to completely remove the air from the cannula during the priming procedure.
- Avoid overinsertion of the cannula and possible impingement of the cannula tip against the vessel walls, ventricular walls, or atrial walls.
- The use of cold solutions and cardioplegic agents, administered externally or via the coronary circulation, may increase the stiffness of urethane plastic and contribute to ventricular perforation under conditions of cannula tip pressure during manipulation of the heart.
- Do not attempt to cross-clamp the cannula with the introducer in place.
- Carefully monitor for signs of occlusion on both the inlet and outlet cannulae.
- User assumes responsibility for any variations in extracorporeal procedures that could compromise the intended use of this device.
- As with all medical devices, this device is to be used by a trained physician only.
- Initiate adequate anticoagulation therapy prior to patient cannulation.

- Only physicians trained and experienced in using percutaneous cannulation techniques (such as the Seldinger technique) and cardiopulmonary bypass should use this device.
- Ischemia may occur distal to the cannula placement site.
- Cannula tip placement should be confirmed by using fluoroscopy.
- Extreme care must be exercised to ensure that the arterial cannula and venous cannula are not confused while being connected to the perfusion circuit.
- In some cases, suture rings may not be detectable on x-ray due to the variances in anatomical densities and x-ray techniques. Modifying the x-ray angle or distance may aid in detection.
- When utilizing vacuum for augmented drainage, do not exceed ~100 mmHg.
- Maximum pressure is 300 mmHg.
- Do not expose to MRI while cannula is inserted in patient.
- Due to the presence of phthalates in the product, the clinician must weigh the medical benefits of product use against the drawbacks of phthalate exposure for male children and pregnant or nursing women.

5. Adverse Effects

The following known adverse effects are associated with the use of the product: abrasion, blood loss, death, dissection, exsanguination, hemolysis, hemolytic anemia, hypervolemia, hypotension, hypovolemia, infection, ischemia, neurological dysfunction, organ dysfunction, perforation, and tissue damage. Vessel damage and complications at the puncture site may occur if the instructions for use are not followed.

6. Complications

Possible complications include those normally associated with cardiopulmonary bypass, anticoagulation, and large bore venous and arterial catheterization.

Note: The benefits of cannulation for extracorporeal circulation must be weighed against the risk of systemic anticoagulation and subsequent propensity for hemorrhage.

7. Sterilization

All Medtronic Bio-Medicus cannulae are sterilized (EtO) prior to shipment. The lot number is designated on the outer package label. This device should not be used if the sterile package has been damaged or if the seal has been broken. Do not resterilize under any circumstances as this device is intended for single use only.

8. Instructions for Use

Note: The introducer has been designed to accommodate a 0.038 in (0.97 mm) diameter guidewire.

1. Select a cannula size that is appropriate for the patient based on criteria including desired flow rate, body surface area, vasculature size, and cannula flow charts. Various components are included in this kit to assist in the access of the vessel. Various sizes are available in order to select the appropriate cannula for each specific procedure.

Note: Flow curve data shown in Figure 2, Figure 3, and Figure 4 were obtained via bench testing using room temperature water. Other factors for a given patient could result in different outcomes than the data presented in the figures.

2. Remove the cannula kit from the carton and inspect the protective, sterile pouches for damage.
3. Using aseptic technique, carefully open the sterile pouches containing the cannula and components.
4. Using sterile technique, transfer the cannula and components into the sterile field.
5. Assemble the cannula/introducer. Insert introducer into the lumen of the cannula and advance it until the hemostasis cap is pushed firmly onto the barb connector.

Note: Tighten the luer fitting cap on the cannula connector (if present) to ensure that the luer fitting is completely sealed.

6. Insert the cannula using a guidewire placement Seldinger technique. Determine the length of cannula to be introduced and note the insertion markings on the cannula.

Note: Some cannulae are labeled with length markings in 10 cm (3.94 in) increments for depth reference during cannula placement.

Note: Confirm the anatomical location of the guidewire before cannula insertion. Anatomical location of the guidewire is critical before cannula insertion.

7. With the vessel identified, introduce the 18 ga (1.02 mm) introducer needle into the desired vessel and aspirate blood.

Note: If no free flow of blood is present, reposition the needle and aspirate until adequate blood flow is established.

Caution: Ensure that the cannula is of adequate size for the vessel to allow distal perfusion of the limb when the cannula is in place. Improper cannula size may be difficult to advance. The vessel must be large enough to ensure perfusion and venous return.

8. Insert the “J” tip of the guidewire through the needle into the vessel.
Note: The guidewire should advance freely without resistance.
Caution: Do not withdraw spring guidewire back into needle as shearing of the guidewire may result. Always remove the needle first.
Caution: Do not kink guidewire during insertion, vessel dilation, and manipulation or it may be difficult to remove.
9. Advance the guidewire into the vessel until it is positioned at a point beyond the desired cannula tip position. Confirm the position of the guidewire with transesophageal echocardiography, fluoroscopy, or other imaging technique.
10. Hold the guidewire in place and remove the introducer needle by sliding it away from the patient over the guidewire.
Note: If using percutaneous technique, it may be necessary to use the #11 scalpel blade to enlarge the skin puncture site.
11. Pass the smaller stepped vessel dilator over the guidewire. Be sure to hold the guidewire firmly in place.
12. Advance the dilator past the step approximately 3 cm (1.18 in) to 5 cm (1.97 in) into the vessel to enlarge the opening.
13. Remove the vessel dilator by sliding it away from the patient over the guidewire. Be sure to apply pressure at the vessel site to maintain hemostasis and to hold the guidewire firmly in place.
14. Repeat steps 10 through 12 with the larger vessel dilator(s).
Caution: Sequentially dilate vessels to match the cannula size. For larger sizes, the tapered introducer tip further dilates the vessel during cannula insertion. To avoid tissue tearing, do not omit a dilator size. Overdilation of the vessel may compromise hemostasis.
15. Position the adjustable radiopaque suture ring to the desired depth.
16. Pass the cannula/introducer assembly over the guidewire, advancing it until the guidewire can be grasped beyond the introducer hub.
17. Insert the tip of the cannula introducer into the vessel approximately 3 cm (1.18 in) to 5 cm (1.97 in). Grasp the guidewire and hold it stationary while advancing the cannula to the desired depth.
Note: Blood leakage from the proximal drainage holes may occur during catheter insertion until these holes are fully advanced into the vessel.
Note: Advancing the cannula and introducer with a rotating motion will make cannula placement easier.
Note: When the bold band closest to the introducer hub is in alignment with the edge of the hemostasis cap, the cannula tip and the introducer tip are at the same depth (Refer to Figure 1). For venous cannula, the alignment of the cannula and introducer ends is indicated by the tip of the triangle-shaped symbol.
18. When advancing the venous cannula into the superior vena cava (SVC), verify positioning with tactile feel, indirect visualization method, or echocardiographic imaging.
Note: For venous models, graduated marks on the hub end of the introducer allow users to approximate the introducer length beyond the tip of the cannula.
19. For proper intracardiac advancement, once the side holes are intravascular, withdraw the introducer to the base of the triangle-shaped symbol representing the tapered tip of the introducer.
Warning: Do not withdraw to the point of the triangle-shaped symbol or there is risk of SVC, innominate vein perforation, or vascular, cardiac trauma during intravascular advancement. Never force the cannula. It should advance easily.
20. When desired cannula placement has been achieved, remove the guidewire and retract the introducer clear of the clamp site.
Caution: To maintain hemostasis, do not completely remove the introducer or hemostasis cap until the clamp is applied.
Warning: Ensure clamp is properly positioned on cannula tubing between tubing connector and wire-reinforced sections. Clamping too close to the connector may damage the connector end. Clamping on wire-reinforced sections may result in permanent wall distortion and/or lumen collapse.
Warning: When the bold band closest to the vascular end of the introducer is visible at the edge of the hemostasis cap, the cannula can be safely cross-clamped. Do not attempt to cross-clamp the introducer.
Note: Perfusion tubing clamps must be used. Many clamps are too weak to ensure complete clamping.
21. Remove the hemostasis cap, exposing the 0.95 cm (0.375 in) barbed connector.
Note: Venous caps are color-coded blue. Since 18 cm (7.09 in) tip-length models may be used as either an arterial delivery cannula or as a venous return cannula, those hemostasis caps are clear, not color-coded.
Note: If needed, an adapter is included with 23 Fr (7.7 mm) and larger venous models for connection to 1.27 cm (0.5 in) tubing.

22. Slowly fill the remainder of the cannula tube with sterile priming fluid, displacing any air through the tube connector end before connecting the fluid-filled tubing.
23. Attach the extracorporeal circuit, taking care to connect the venous cannula to the venous return line (inlet to pump) and the arterial cannula to the arterial return line (outlet from oxygenator), ensuring air-free connection.
- Note:** Venous cannula can be identified by an arrow, in the clamp area next to the Bio-Medicus cannula logo, indicating flow out of the connector. See Figure 1.
24. The tubing may be further secured with cable ties at the barbed connector. Secure the device to the patient following standard medical practice.
25. In order to avoid accidental displacement during extracorporeal circulation, the suture may be wrapped around the suture ring to secure the cannula and avoid kinking the body of the device.
26. The retention clip may be used to secure cannula and/or tubing to the patient or drapes. The clip may be initially placed using the adhesive tape. To ensure that the cannula is secure during the procedure, clamp or suture the retention clip through the retention clip holes.
- Warning:** Do not rely on the adhesive alone to secure the clip.
27. When it is determined that extracorporeal bypass is no longer required, clamp the non wire-reinforced tube section, carefully remove the cannulae from the vessels, and repair their placement sites according to standard surgical procedure.
- Warning:** This product has been qualified for use up to 6 hours. Suitability for long-term use has not been determined.
28. After use, this product may be a potential biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, regional, and national laws and regulations.

9. The following disclaimer of warranty applies to United States customers only:

Disclaimer of Warranty

ALTHOUGH THE ADULT CANNULA KIT, HEREAFTER REFERRED TO AS “PRODUCT” HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

10. The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of Warranty

ALTHOUGH THE ADULT CANNULA KIT, HEREAFTER REFERRED TO AS “PRODUCT” HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be

construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Medtronic

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
+1 763 514 4000
LifeLine Technical Support, 24-hour consultation
service:
1 877 526 7890
Medtronic Perfusion Systems
7611 Northland Drive
Minneapolis, MN 55428
USA
+1 763 391 9000
Customer service and product orders:
1 800 854 3570
www.perfusionsystems.com

Medtronic Mexico S. de R.L. de CV
Av. Paseo Cuapah 10510 El Lago
C.P. 22210 Tijuana, Baja California
Mexico

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Insertion Kit

1. Product Description

This product contains the necessary components to achieve insertion of a Bio-Medicus™ cannula and introducer. The included items are: a Seldinger needle, a guidewire, a scalpel blade, stepped dilators, and a catheter tip syringe.

Stepped dilators by model

96551 and 96552 kits include 8 Fr/10 Fr, 12 Fr/14 Fr, and 16 Fr/18 Fr dilators.

96553 kit includes 8 Fr/10 Fr, 12 Fr/14 Fr dilators.

Guidewire size by model

96551 kit includes 0.038-in (0.965-mm) x 180-cm (70.87-in) guidewire.

96552 kit includes 0.038-in (0.965-mm) x 100-cm (39.37-in) guidewire.

96553 kit includes 0.025-in (0.635-mm) x 60-cm (23.62-in) guidewire.

Sterile, nonpyrogenic, disposable, single use only.

2. Indications for Use

This kit is intended for use by a trained physician only, to assist in vessel cannulation for cardiopulmonary bypass circulation. Standard surgical or percutaneous insertion techniques can be employed. This kit is intended for use for up to 6 hours.

3. Contraindications

Alone, this kit is not a medical treatment device. Selection of patients as candidates for such procedures is the physician's responsibility. The outcome is dependent on many variables, including patient pathology, surgical procedure, and perfusion procedures. These devices are not intended for use except as indicated above. Do not use the kit if the patient has severe peripheral atherosclerosis or severe arterial dissection.

4. Warnings and Precautions

- Read all warnings, precautions, and Instructions for Use carefully prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- Verify Use By date.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Do not autoclave.
- Do not use if the sterile package has been damaged or if the seal has been broken in any way.
- Thoroughly inspect the cannula prior to use to verify that the lumen and side holes (if applicable) are patent and that the cannula has not been damaged or kinked.
- Use caution when assembling or manipulating the introducer within the cannula, or when moving the suture ring.
- Avoid false lumen cannulation.
- Never advance, position, or place the cannula without the introducer as this may kink and damage the cannula.
- Use care to protect the side hole areas from kinking as they are passed into the vessel.
- Manipulate cannula in such a manner as to prevent kinks or any restrictions that may alter flow.
- Manipulation of the heart with cannula in place may cause flow restrictions.
- Avoid kinking when attaching tubing and when anchoring or untethering the cannula body.
- A guidewire inserted into the left ventricle through the left atrium can perforate the left ventricular apex if undue pressure is exerted by the tip.
- User assumes responsibility for any variations in extracorporeal procedures that could compromise the intended use of this device.
- Thoroughly inspect the components prior to use.
- As with all medical devices, this device is to be used by a trained physician only.
- Initiate adequate anticoagulation therapy prior to patient cannulation.
- Only physicians trained and experienced in using percutaneous cannulation techniques (such as the Seldinger technique) and cardiopulmonary bypass should use this device.

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7. Sterilization

All Bio-Medicus™ insertion kits are sterilized (EtO) prior to shipment. The Lot Number is designated on the outer package label. This device should not be used if the sterile package has been damaged or if the seal has been broken. Resterilization is not recommended under any circumstances as this device is intended for single use only.

8. Instructions for Use

Note: Only sterile cannula should be placed into the vascular system. Various components are included in this kit to assist in the access of the vessel. Various sizes of dilators and cannulae are also available in order to select the best cannula for each specific procedure. When it is determined that extracorporeal bypass is no longer required, the cannulae are carefully removed from the vessels, and their placement sites are repaired according to standard surgical practice.

- Using aseptic technique, carefully open the sterile insertion kit tray.
Note: Determine the length of cannula to be introduced and note the insertion markings on the cannula. Some cannulae are labeled with length markings in 10-cm (3.94-in) increments for depth reference during cannula placement.
- With the vessel identified, introduce the Seldinger needle into the desired vessel and aspirate blood.
Note: If no free flow of blood is present, reposition the needle and aspirate until adequate blood flow is established.
Caution: Ensure that the vessel is of adequate size for cannula size selected to allow distal perfusion of the limb when the cannula is in place.
- Insert the "J" tip of the guidewire through the needle into the vessel.
Note: Guidewire should advance freely without resistance.
Caution: Do not withdraw spring guidewire back into the needle as shearing of the guidewire may result. Always remove the needle first.
Caution: Do not kink the guidewire during insertion, vessel dilation, and manipulation or it may be difficult to remove.
- Advance the guidewire into the vessel until it is positioned at a point beyond the desired cannula tip position. Confirm the position of the guidewire with transesophageal echocardiography, fluoroscopy, or other imaging technique.
- Hold the guidewire in place and remove the introducer needle by sliding it away from the patient over the guidewire.
- When using percutaneous technique, it may be necessary to enlarge the skin puncture site with the #11 scalpel blade.
- Pass the smaller, stepped vessel dilator over the guidewire. Be sure to hold the guidewire firmly in place.
- Advance the dilator past the step approximately 3 cm (1.18 in) to 5 cm (1.97 in) into the vessel to enlarge the opening.
- Remove the vessel dilator by sliding it away from the patient over the guidewire. Be sure to apply pressure to the site to maintain hemostasis and hold the guidewire firmly in place.
- Repeat steps 7 through 9 with the larger vessel dilator(s).
Caution: Sequentially dilate the vessel to match cannula size. For larger sizes, the tapered introducer tip further dilates the vessel during cannula insertion. To avoid tissue tearing, do not omit a dilator size. Dilation too large may compromise hemostasis.
- Pass the cannula/introducer assembly over the guidewire, advancing until the guidewire can be grasped beyond the introducer hub.
- Insert the tip of the cannula introducer into the vessel approximately 3 cm (1.18 in) to 5 cm (1.97 in). Grasp the guidewire and hold it stationary while advancing the cannula to the desired depth.
Note: Advancing the cannula and introducer with a rotating motion will make cannula placement easier.

13. Remove the guidewire.
14. Remove the introducer so the cannula can be clamped.
15. For additional information, refer to cannula Instructions for Use.
16. After use, this product may be a potential biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, regional, and national laws and regulations.

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Disclaimer of Warranty

ALTHOUGH THE INSERTION KIT, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

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Disclaimer of Warranty

ALTHOUGH THE INSERTION KIT, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

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The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.
Trademarks may be registered and are the property of their respective owners.
©2018 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.