Instructions for Use

Yuranos™ AUI and AE Stent Graft System

Lifetech Scientific (Shenzhen) Co., Ltd

Contents

1	Device description.....	1
2	Intended Use.....	4
3	Clinical Benefits	4
4	Contraindication.....	4
5	Pre-implantation Preparation.....	5
6	Implantation Procedure	6
7	Management Suggestions	12
8	Potential Adverse Event.....	14
9	Summary of Clinical Data.....	15
10	Warnings and precautions	15
11	Follow-up	17
12	Shelf Life.....	17
13	Storage and Transportation Requirements.....	18
14	Package and Label	18

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1 Device description

The Yuranos™ AUI and AE Stent Graft System consists of a pre-loaded stent graft and a delivery system. The stent grafts include Aorto-Uni-Iliac Stent Graft and Aortic Extension Stent Graft, which are composed of a metal stent (nitinol alloy), a graft (PET), suture (PET), a stainless-steel sleeve at the stent connection (316L stainless steel), and markers (tantalum). There is a barb structure on the bare stent of Aorto-Uni-Iliac Stent Graft and Aortic Extension Stent Graft. The delivery system is mainly composed of a sheath, an inner core, grips, an immobile pole, and other parts, to realize the delivery, control and release of the stent grafts in human bodies, the delivery system has a proximal hydrophilic coated segment of 520mm.

1.1 Abdominal Aortic Stent Graft



Figure 1 Yuranos™ AUI and AE Stent Graft System

(1: “O”-shaped mark point; 2: “8”-shaped mark point)

The Yuranos™ Aorto-Uni-Iliac and Aortic Extension Stent Graft System adopts a split-type design, including Aorto-Uni-Iliac Stent Graft, and Aortic Extension Stent Graft. The Aorto-Uni-Iliac Stent Graft includes a proximal bare stent, a stent graft's main body and a single iliac branch. The Aortic Extension Stent Graft includes a proximal bare stent and a stent graft's main body. The main body of Aorto-Uni-Iliac Stent Graft can relate to the Aortic Extension Stent Graft to form a whole assembly. Each stent graft is separately loaded in its respective delivery system to be introduced into the human vascular system during a medical procedure.

The Stent Graft is composed of a superelastic nitinol stent and a PET fabric cover, which is fixed onto the metal stent with suture. After the stent graft is implanted into the human body, the superelastic metal stent can provide sufficient and stable expansion force to deploy the stent graft and anchor it in the corresponding position. The bare stent with a barb can provide enough anchoring force to prevent the stent graft from migration. The fabric graft can serve as an alternative conduit for blood flow within the patient's vasculature by excluding the aneurysmal sac from blood flow and pressure, thereby reducing the risk of aneurysm rupture.

1.2 The delivery system for the AUI and AE Stent Graft

- | | | | | | |
|----|-------------------|----|-----------------|----|------------------|
| 1 | Tapered Tip | 2 | AUI or AE | 3 | Markband |
| 4 | Fixation cap | 5 | Fixation anchor | 6 | Immobile Pole |
| 7 | Sheath | 8 | Front Grip | 9 | Sliding Grip |
| 10 | Trigger | 11 | Screw Gear | 12 | Rear Grip |
| 13 | Proximal releaser | 14 | Luer Connector | 15 | Hemostatis valve |
| 16 | Unlocking buckle | | | | |

Figure 2 The delivery system for AUI and AE stent graft

1.3 Device specifications


Both AUI and AE device are pre-loaded onto their respective delivery system and supplied sterile. The AUI System and AE System are available in the following lengths and diameters as shown in following table 1, 2 and 3.

Table1. Specifications Summary of Yuranos™ Aorto-Uni-Iliac Stent Graft (AUI)

Unit: mm

Specification*	Proximal Dia. of main body (D1)	Distal Dia. of long branch (D2)	Covered length of main body & long branch (L1)	Covered length of main body (L2)	Length of Bare stent (L3)	Delivery System
AUI-20XX-WW-YYY	20	10	100, 110, 120, 130, 140, 160	50, 70	16	18F
AUI-22XX-WW-YYY	22	10	100, 110, 120, 130, 140, 160	50, 70	16	18F
AUI-24XX-WW-YYY	24	12	100, 110, 120, 130, 140, 160	50, 70	16	18F
AUI-26XX-WW-YYY	26	12	100, 110, 120, 130, 140, 160	50, 70	16	20F
AUI-28XX-WW-YYY	28	12,14	100, 110, 120, 130, 140, 160	50, 70	16	20F

AUI-30XX-WW-YYY	30	12,14	100, 110, 120, 130, 140, 160	50, 70	16	20F
AUI-32XX-WW-YYY	32	12,14	100, 110, 120, 130, 140, 160	50, 70	16	20F
AUI-35XX-WW-YYY	35	12,14	100, 110, 120, 130, 140, 160	50, 70	16	20F


* **AUI-ZZXX-WW-YYY** is the specification of Yuranos™ AUI Stent Graft System which ZZ is Proximal Dia. of main body (D1), XX is Distal Dia. of long branch (D2), WW is Covered length of main body (L2), and YYY is Covered length of main body & long branch (L1).

Table2. Specifications Summary of Yuranos™ Aortic Extension Stent Graft ((AE)

Unit: mm

Specification*	Proximal Dia. (D1)	Distal Dia. (D2)	Covered length of stent graft (L1)	Length of bare stent graft(L2)	Delivery System
AE-20XX-YY	20	20	40,60	16	18F
AE-22XX-YY	22	22	40,60	16	18F
AE-24XX-YY	24	24	40,60	16	18F
AE-26XX-YY	26	26	40,60	16	20F
AE-28XX-YY	28	28	40,60	16	20F
AE-30XX-YY	30	30	40,60	16	20F
AE-32XX-YY	32	32	40,60	16	20F
AE-35XX-YY	35	35	40,60	16	20F

* **AE-ZZXX-YY** is the specification of Yuranos™ AE Stent Graft System which ZZ is Proximal Dia. (D1), XX is Distal Dia.

(D2), and YYY is Covered length of stent graft (L1).

Table3. Specifications Summary of Delivery System (in mm)

Name of delivery system for Abdominal Aortic Stent Graft	Specification	Outer diameter of sheath (OD)	Effective length (L0)
Delivery system (applies to AE/AUI)	18F	6.00	520
	20F	6.67	

Table4. Stent Graft Materials Composition

Stent graft component configurations	Materials composition
Yuranos™ Aorto-Uni-Iliac Stent Graft (AUI)	Nitinol (56%-59% w/w), Stainless steel (4%-7% w/w), Film covering (PET, 23%-26% w/w), Tantalum (5%-10%w/w) Suture (PP, 1%-9% w/w)
Yuranos™ Aortic Extension Stent Graft (AE)	Nitinol (51%-56% w/w) Stainless steel (8%-12% w/w) Film covering (PET, 24%-33% w/w) Tantalum (1%-4% w/w) Suture (PP, 5%-7%w/w)

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2 Intended Use

The Yuranos™ AUI and AE Stent Graft Systems is for the treatment of infrarenal abdominal aortic aneurysm.

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4 Contraindication

The contraindications of Yuranos™ AUI and AE Stent Graft Systems are:

- (1) Patients with acute systemic infection;
- (2) Patients with mesenteric blood flow mainly supplied by the inferior mesenteric artery;
- (3) Patients with accessory renal artery original from abdominal artery;
- (4) Patients who have allergic reaction to the device;
- (5) Patients who can not tolerate contrast agents due to renal insufficiency;
- (6) Patients who are allergic to contrast agents;
- (7) Aneurysms neck with thrombus;
- (8) Women who are pregnant or nursing.

5 Pre-implantation Preparation

Read the instruction for use carefully before using the device. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician rather than replacing the physician judgment.

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5.3 Patient preparation

- Measure proximal aortic neck length and vessel inner diameter;
- Measurement of the angle between aneurysm neck and aneurysm, i.e. aneurysm neck angle;
- Evaluation of the quality of aneurysm neck;
- Confirmation of the suitability of vascular morphology of patient for endovascular repair;
- All kinds of necessary routine tests and examinations, and monitoring of important vital signs;
- Vital signs related to anesthesia, anticoagulation, and angiography that need to be monitored based on routine practice;
- Adjustment of the patient's position so that the area from the aortic arch to the iliac artery bifurcation

can be fully exposed to the X-ray;

- Adequate iliac/femoral artery access and exposure of both femoral arteries with standard surgical techniques;
- Systematic anticoagulant preparation (if the patient is allergic to heparin, choose an appropriate alternative).

5.4 Products preparation

- The stent graft size should be selected according to the vessel inner diameter, and the proximal diameter and distal diameter of the selected stent graft should be 110% to 120% of the landing zone vessel inner diameter;
- The length of AUI stent graft should be selected based on the lengths of the infrarenal abdominal aorta and iliac artery. Generally, the distal anchoring area should be at least 15 mm;
- Check and make sure there is no package damage due to transportation. If the package is damaged or there is a potential risk of sterilization failure or damage to the device, please do not use it and return it to Lifetech Scientific (Shenzhen) Co., Ltd.;
- After the stent graft size is determined, unpack the outer packaging box of the corresponding size product, and take out the product together with the dialysis bag. Open the dialysis bag and take out the stent graft system under a sterile environment. Carefully check whether there is any crease on the surface of the sheath of the delivery system and whether the tapered tip of the front-end fits well with the sheath port. If the package is damaged or the product is damaged, please do not use it.

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- Connect a syringe to the hemostatis valve and Luer connector respectively, and rinse the sheath inside and guidewire cavity with sterile heparinized saline solution, as shown in Figure 3.

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- Wipe the outer surface of the sheath with sterile heparinized saline solution to activate the hydrophilic coating on the sheath, as shown in Figure 4.

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6.2 Direction for use for Aorto-Uni-Iliac Stent Graft System

- 1) Keep the hemostatis valve at the end of the delivery system toward the opposite side of the operator and insert the delivery system of the Aorto-Uni-Iliac Stent Graft along the guidewire (the single iliac branch and hemostasis valve are on the opposite side of the delivery system when the stent graft is loaded).

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- 2) Determine the directions of the single iliac branch of the stent graft. The directions of the branch of stent graft are adjusted by rotating the delivery system. The “O”-shaped and “8”-shaped marker points at the proximal end of the graft can be used as references for adjusting the positions of single iliac branch (the “O”-shaped marker point is located on the side of the single iliac branch), so that the “8”-shaped marker point at the proximal end of the stent graft is always on the side of the branch on the opposite side of the patient. The final positions of the single iliac branch should be confirmed by the “O”-shaped mark points at the branch (from the bifurcation position, there are two consecutive “O”-shaped mark points on the single iliac, as shown in Figure 1).

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- 3) Position the proximal end of the stent graft. Adjust the delivery system so that the “O”-shaped and “8”-shaped mark points at the beginning end of the stent graft are located on the lower edge of the lower renal artery opening.
- 4) Slowly release the stent graft. Under monitoring with DSA, hold the front grip firmly with the left hand, and rotate the sliding grip counterclockwise with the right hand to slowly release the bare stent graft and the front segment of the stent graft until at least one coil of the stent graft completely engages the vascular wall, as shown in Figure 6.

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Note 1: Before releasing the stent graft, keep the exposed sheath straight and ensure that the front

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- 5) Quickly release the stent graft. When a stent graft is to be released quickly, hold the front grip with the left hand to stabilize the delivery system, hold the trigger with the right thumb backward, remove the screw gear lock, and pull the sliding grip backward to release the stent graft quickly. At the same time, the process of stent graft release under DSA was observed until the single iliac stent was completely released. As shown in Figure 8.


Figure 8 The single iliac stent was fully released

- 6) Later release the proximal end of the stent graft. Fix the delivery system with the left hand and remove the unlocking buckle at the end of the delivery system with the right hand (the unlocking buckle is a blue translucent safety device, and its removal is as follows: gently press the end of the unlocking buckle with the thumb and index finger and pull it backward to remove it). Under monitoring with DSA, rotate the proximal releaser at the end of the delivery system clockwise, and

pay attention to the forward movement of the tapered tip until the proximal bare stent graft is fully deployed, as shown in Figure 9 (rotate in the direction of the arrow printed on the releaser, the tapered tip can be advanced to deploy the proximal bare stent graft).

Note 1: Please make **sure** that the relative position of the delivery system and the patient remain unchanged during the rotation of the releaser, to avoid stent graft migration.

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7) Reset the tapered tip of the delivery system for Bifurcation Stent Graft.

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a. After the proximal stent graft is released, the whole delivery system is advanced by more than 2 cm with the guidance of DSA, so that the tapered tip and fixation anchor are far away from the bare stent graft, and in a straight vascular cavity (to prevent the fixation cap on the tapered tip from scraping the stent graft or blood vessels when the tapered tip is withdrawn), as shown in Figure 10.

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b. Fix the delivery system and rotate the releaser counterclockwise (the opposite direction of the arrow printed on the knob) to withdraw the tapered tip slowly. Observe under DSA until the gap between the fixation cap and anchor of the tapered tip is completely closed when the resistance to the knob increases, and the knob is basically rotated to the farthest end, as shown in Figure 11.

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c. Hold the sliding grip with the right hand, pull the trigger backward to unlock the screw gear, hold the front grip with the left hand, and withdraw it until the tapered tip is reset to the sheath port (at this point, the gap between the front grip and the sliding grip is also basically closed), as shown in Figure 12.

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
Figure 10 Advance the delivery system as a whole


Figure 11 Close the gap between the tapered tip and the fixation anchor


Figure 12 Tapered tip reset

8) Place a plug in the other side of iliac artery, as shown in Figure 13


Figure 13 Plug placement

9) Perform angiography to assess the wall engagement for the proximal and distal stent graft and endoleak, and to confirm the positions of the stent components relative to the aneurysm and renal artery. Finally, withdraw the guidewire and catheter and close the femoral artery incisions with the standard surgical procedure.

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When Aortic Extension Stent Graft needs to be used, ensure that there is a minimum 1.5cm overlaps between the Aortic Extension Stent Graft and Aorto-Uni-Iliac Stent Graft. Completed the Aorto-Uni-Iliac Stent Graft deployment steps, except rotate the handle to open the Aortic Extension Stent Graft entirely before releasing the proximal end of the suprarenal stent of the Aortic Extension Stent Graft.

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- a. Keep the outer sheath of the conveyor delivery system still, press the back ends of the two shell parts of the front grip with fingers to enlarge the gap between the shells, insert the surgical scissors into the gap, pry the front grip open, and remove it;
- b. Use the surgical scissors to pry off the blue trigger on the sliding grip;
- c. With the same method as in step b, squeeze the tail of the sliding grip in dislocation to enlarge the gap between the two valve shells, pry the two-valve shells with the surgical scissors, and remove it;
- d. Pull back the blue ring trigger sleeve (the part before the metal spring) to compress the spring. After unlocking the screw gear, remove the screw gear sliding sleeve assembly from the screw gear. At this point, the screw gear is in a half-loosened state. Check whether there is foreign matter in the screw gear which may hinder the sheath movement;
- e. After removing the obstacles, fix the screw gear and pull back the sheath connector (a white translucent part) to withdraw the outer sheath and complete the further release of the stent graft;
- f. There is a small hole on the left and right sides of the hemostatis valve hose connection at the end of the screw gear. Use surgical scissors or hemostatic forceps to press the spring buckle in the small holes at the same time, and pull the screw gear forward to further remove the screw gear, facilitating further inspection and troubleshooting;
- g. After the screw gear is removed, restore the delivery system to a simple backward pull release system. Fix the transparent immobile pole joint (connected with the steel pipe) and withdraw the sheath joint (connected with the blue sheath) to release the stent graft.

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- a. Make sure that the releaser is rotated to the farthest end (toward the sheath) and kept there for period. If the stent graft still cannot be released, slightly advance or withdraw the delivery system (the operation needs to be carried out carefully under DSA monitoring to prevent stent graft migration), and observe the deployment of the tapered tip end;
- b. Under monitoring with DSA, carefully turn the releaser back to its original position in the opposite direction. Use a compliant or semi-compliant aorta balloon, expanding the balloon in the main segment of the Abdominal Aortic Bifurcation Stent Graft, anchor the stent graft, then try to release the proximal bare stent graft again according to Step 6) in Section 5.2;
- c. Introduce the capture device through the upper limb approach, try to capture the tapered tip of the delivery system, rotate the releaser of the delivery system to the farthest end (towards the sheath), fix the grip of the delivery system, and gently drag the tapered tip with the capture device to facilitate the release of the later release mechanism.

Note: The above suggestions are only used as emergency measures under abnormal circumstances and cannot replace the professional judgment of physicians. The above operations should not be attempted during normal stent graft release.

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After the stent graft is released, the balloons of corresponding compatible sizes can be delivered selectively to the proximal end of the main bifurcated stent graft, the junction of the branch of the bifurcated stent graft and the extension stent graft, and the distal end of the extension stent graft. The balloons can be expanded to the sizes of the blood vessels, to facilitate the stent graft to engage the vessel wall and stabilize the stent graft. A leak at the proximal or distal end of the stent or at the joint may be reconstructed to the vessel wall using a balloon, but if the leak is severe and cannot be corrected by balloon expansion, it may be managed by adding an extended stent to the already placed stent assembly.

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- Aneurysm enlargement
- Aneurysm rupture
- Thrombosis or vessel thrombosis
- Arterial injury
- Artery perforation
- Arteriovenous fistula
- Bleeding, hematoma, or coagulopathy
- Heart failure/myocardial infarction
- Turn to surgical operation
- Edema
- Sexual dysfunction
- Infections
- Nerve injuries, local or systemic injuries (e.g. stroke, paraplegia, paraparesis)
- Pulmonary Adverse Events
- Renal function injury or renal failure
- Wound infection and subsequent attendant problem
- Allergic reaction
- Anesthetic complications
- Arrhythmia
- Bowel complications
- Embolism
- Intercostal pain
- Shock
- Tissue necrosis
- Wound complications and subsequent attendant problems
- Death
- Device-related Adverse Events include but are not limited to:
 - ♦ Insertion and removal difficulties
 - ♦ Release difficulties
 - ♦ Failure to deploy
 - ♦ Inaccurate placement
 - ♦ Endoleak
 - ♦ Stent migration
 - ♦ Stent fracture
 - ♦ Rupture of graft material
 - ♦ Erosion
 - ♦ Stent thrombosis or embolization
 - ♦ Stent graft expansion
 - ♦ Stent graft twisting or kinking

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9 Summary of Clinical Data

“Clinical performance of the Yuranos™ AUI and AE Stent Graft Systems is supported by clinical data from the Yuranos™ AB Stent Graft System and by bench and simulated use testing; no independent clinical data are available for these components.”

10 Warnings and precautions

10.1 Summary of Clinical investigation in China

- Please read all instructions carefully. Failure to follow the instruction for use may cause serious injury to patients.
- This product (including the stent graft and the delivery system) is disposable and cannot be used again after sterilization or disinfection.
- Please check if the sterile packaging is damaged or unintentionally opened before use. If yes, DO NOT use damaged or open packages. And please follow your local or hospital's rules to handle such a situation.
- Do not use if the sterile pouch is compromised in any way.
- Please check the product carefully before use. If any damage is found or physicians think there is any problem affecting the safety of use, please do not use it. Please follow your local or hospital's rules to handle such a situation.
- Regarding safe disposal of the device or its package, please follow your local or hospital's rules.
- This kind of device is only used by physicians with endovascular interventional treatment skills (including those trained on using this device).
- The long-term efficacy and safety of endovascular graft exclusion have not been established. All patients with endovascular graft exclusion are required to receive long-term follow-up to evaluate their health and stent graft status.
- After use, discard the products and packaging according to the policies of the hospital, regulations, or local government.
- STERILIE- DO NOT RESTERILIZATION-SINGLE USE ONLY. Do not reuse, reprocess or resterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death.

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- Patient condition should be evaluated before operation, including heart failure/infarction, pulmonary Adverse Events, edema, and renal insufficiency.
- After stent graft implantation, all patients need to be regularly monitored for stent graft status (including stent graft endoleak, migration, shape and location changes, aneurysm enlargement, etc.) at least once a year, by X-ray or CT with or without contrast;
- There is a risk of continuous aneurysm dilation or other diseases at the stent graft placement site after stent graft implantation, which requires secondary interventional operation or open surgery. A persistent aneurysm dilation or leakage may cause aneurysm rupture.
- Patients with infection after device implantation need antibiotic treatment.
- During stent graft implantation or the second intervention, a standby vascular surgery team is needed to promptly perform open surgery on patients in case of intervention failure.
- The size and shape of the access vessels should conform to the general principles of interventional operation, and match the size of sheath of the stent graft delivery system. Excessive curving blood

vessels will increase the risk during stent graft delivery and release and the delivery system recovery, increase the difficulty of operation, and even lead to operation failure.

- Significant or circumferential calcification or mural thrombus within the treatment length may adversely impact device patency
- The Yuranos™ AUI and AE Stent Graft System products are not recommended for patients who weigh more than 150 kg because obesity may reduce the accuracy of DSA.
- Patients with systemic infection have an increased infection risk after intravascular stent implantation.
- Patients with coagulation disorders have an increased risk of endoleak and bleeding after stent implantation.

The safety and efficacy of AUI and AE Stent Graft have not been evaluated in the following patients:

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- During stent implantation, anticoagulant therapy should be administered according to hospital or physician advice. If patients are allergic to heparin, other appropriate anticoagulants should be considered;
- In order to reduce the risk of device contamination and its resulting infection, the time to stent implantation after removing the sterilization barrier should be shortened as far as possible;
- In the process of inserting the delivery system, it is necessary to keep the guidewire in position;
- Do not bend the sheath of the delivery system, otherwise, the delivery system/or stent graft may be damaged;
- The DSA equipment must be used to observe during the delivery system introduction and delivery, stent graft release, and delivery system withdrawal;
- During stent implantation, a contrast medium should be used, which will increase the risk of postoperative renal failure in patients with renal insufficiency;
- If balloon catheters are used, do not over-expand or expand them outside the graft. Follow all manufacturer's instructions for catheter operation;
- When the delivery system needs to be rotated to adjust the position of the stent graft, it should be rotated as a whole. It can only be rotated when the stent graft is not released completely. Once the stent graft starts to deploy, the delivery system cannot be rotated. When the sheath is in a twisted blood vessel, it may not be able to rotate along with the rotating grip. At this point, it is necessary to consider withdrawing the sheath to a relatively straight segment of the vessel and then try to rotate it;
- Incorrect stent graft release position may lead to endoleak or vascular blockage due to migration;
- Inadequate anchorage of stent graft increases the risk of stent graft displacement, and the incorrect release mode or stent migration may require surgical intervention;
- During the introduction or withdrawal of the guidewire or the delivery system, if abnormal resistance is encountered, please stop the operation and eliminate the cause of failure before proceeding to the next operation;

- The stent graft should not be used to block any vascular branches supplying blood to viscera or limbs unless it is medically proved to be applicable.

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Item	Magnetic field conditions	
Static magnetic field strength	1.5 Tesla	3.0Tesla
Spatial gradient magnetic field	2000G/m	1700G/m
Maximum whole-body average specific absorption rate (SAR)	2.0 W/kg with imaging for 15 minutes	

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11 Follow-up

Follow-up period

It is suggested that the clinical follow-up should be performed at 24 hours, 1 month, 3 months, 6 months, 12 months, and 24 months after implantation. The follow-up period can be adjusted as appropriate according to the circumstances of each individual patient.

Methods used during follow up

According to the patient's condition, blood pathology, X-ray fluoroscopy, B-US, CTA, MRI or DSA examinations can be performed.

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12 Shelf Life

Yuranos™ AUI and AE Stent Graft System is sterilized with ethylene oxide. The manufacture date and expiration date of the product are shown on the label. The product labels indicate the safe use period, and expired products shall not be used. Please use before the Expiration Date marked on the packaging label.

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Yuranos™ AUI and AE Stent Graft System products should be stored in a clean, room temperature in a cool place, keep away from sunlight, and keep dry.

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The Yuranos™ Abdominal Aortic Stent Graft and the delivery system should not contact with organic solvents (including alcohol).


14 Package and Label

Yuranos™ AUI and AE Stent Graft System are packed in PETG tray boxes and sealed with double-layer Tyvek pouches. The Aorto-Uni-Iliac Stent Graft System and Aortic Extension Stent Graft System products are respectively packed in outer box with an IFU, and customer information feedback sheet.

The label indicate the product specification, manufacture date (sterilization date), validity period, batch number, product barcode, sterilization mark, etc.

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	Date of manufacture
	Use-by date
	Serial number
	Batch code
	Sterilized using ethylene oxide
	Do not reuse
	Consult Instructions for Use
	Keep away from the sunlight
	Keep dry
	Non-pyrogenic
	MR conditional
	Do not use if package is damaged
	Do not re-sterilize
	Manufacturer

	Catalogue number
	MR Conditional, for implant only
	Importer
	Country of manufacture
	Double sterile barrier system
	Medical device
	Unique device identifier
	Temperature limit

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Lifetech Scientific (Shenzhen) Co., Ltd.

I Device description

The Thoracoabdominal Artery Stent Graft System is intended for endovascular repair of thoracoabdominal aortic aneurysms and aortic dissections. The system comprises three sub-systems: 1) the G-Branch™ Thoracoabdominal Aortic Stent Graft System, 2) the SilverFlow™ PV Peripheral Vascular Stent Graft System, and 3) the Aortic Extension Stent Graft System.

The G-Branch™ Thoracoabdominal Aortic Stent Graft System consists of the G-Branch™ Thoracoabdominal Aortic Stent Graft (TAAA stent graft) and its delivery system.

The SilverFlow™ PV Peripheral Vascular Stent Graft System consists of the SilverFlow™ PV Peripheral Vascular Stent Graft (PV stent graft) and its delivery system.


The Aortic Extension Stent Graft System consists of the G-Branch™ AE Main Body Extension Stent Graft (AE stent graft), the G-Branch™ AAA Bifurcated Stent Graft (AAA stent graft), the G-Branch™ IE Iliac Extension Stent Graft (IE stent graft) and their respective delivery systems.

All stent grafts are pre-loaded onto their respective delivery systems. During the procedure, each stent system is introduced via an appropriate vascular access route, navigated to the target lesion site under fluoroscopic guidance, and the corresponding stent graft is accurately deployed to achieve exclusion of the aneurysm or dissection segment.

1.1 Stent Graft Description

1.1.1 G-Branch™ Thoracoabdominal Aortic Stent Graft (TAAA Stent Graft)

The G-Branch™ Thoracoabdominal Aortic Stent Graft is a self-expanding stent graft constructed from PET polyester fabric and Nitinol (NiTi) stent frame. The Nitinol stent frame is securely fixed to the PET polyester fabric using polyester sutures. The proximal and distal stent graft junctions incorporate a conical transition. The TAAA stent graft includes two internal branches for the celiac trunk (CA) and superior mesenteric artery (SMA), and two external branches for the left and right renal arteries (LRA and RRA), as shown in Figure 1.

1	"O" Radiopaque Marker	2	"8" Radiopaque Marker
3	"G" Radiopaque Marker	4	"U" Radiopaque Marker
5	"O" Radiopaque Marker	6	Radiopaque Marker Band of Internal Branches
7	Left Internal Branch	8	Right Internal Branch
9	External Branches	10	Radiopaque Marker Band of External Branches

Figure 1: G-Branch™ Thoracoabdominal Aortic Stent Graft (TAAA Stent Graft)

1.1.2 SilverFlow™ PV Peripheral Vascular Stent Graft (PV Stent Graft)

The SilverFlow™ PV Peripheral Vascular Stent Graft is a self-expanding stent graft constructed from a NiTi alloy stent frame covered with PTFE membrane, as shown in Figure 2. The PV stent graft is used in conjunction with the TAAA stent graft to reconstruct the visceral branches of the thoracoabdominal aorta.



Figure 2: SilverFlow™ PV Peripheral Vascular Stent Graft (PV Stent Graft)

1.1.3 Aortic Extension Stent Graft (G-Branch™ AE, AAA, IE Stent Graft)

The G-Branch™ AE Main Body Extension Stent Graft, G-Branch™ AAA Bifurcated Stent Graft, and G-Branch™ IE Iliac Extension Stent Graft are all self-expanding stent grafts constructed from NiTi alloy stent frame covered with PTFE membrane. The AE stent graft is used in conjunction with the TAAA stent graft to extend its distal end, as shown in Figure 3(a). The AAA stent graft features long and short limbs configurations, which are designed for use with the corresponding IE stent graft, as shown in Figures 3(b) and 3(c). In combination, the AAA and IE stent grafts provide a continuous endovascular pathway extending from the distal end of the TAAA stent graft to the iliac arteries, thereby enabling effective treatment of thoracoabdominal aortic aneurysms or aortic dissections involving the aortic bifurcation and iliac segments.

(c) G-Branch™ IE Iliac Extension Stent Graft
(IE Stent Graft)

1	“O” Radiopaque Marker	2	“8” Radiopaque Marker	3	“V” Radiopaque Marker
4	“e” Radiopaque Marker	5	“√” Radiopaque Marker		

1.2 Delivery System Description

The delivery system for the TAAA stent graft primarily consists of a tapered tip, a hydrophilic coated sheath, preloaded guidewires and a handle assembly. The stent graft is preloaded onto the delivery sheath and positioned near the tapered tip. The 0.018-inch preloaded guidewires pass through the internal branches of the TAAA stent graft and are secured within the sheath by a guidewire locking mechanism, as shown in position ⑪ of Figure 4. The detailed components of the delivery system are listed in Table 1.


Figure 4: TAAA Delivery System

Table 1 Components of TAAA Delivery System

1	Tapered Tip	5	Trigger	9	Safety Buckle
2	Radiopaque Marker Band	6	Slider Grip	10	Haemostatic Valve
3	Sheath	7	Release Lock	11	Guidewire Locker
4	Front Grip	8	Proximal Releaser	12	Preloaded Guidewire

1.2.2 SilverFlow™ PV Delivery System

The delivery system for the PV stent graft primarily consists of a tapered tip, a hydrophilic coated sheath, and a handle assembly. The stent graft is preloaded within the sheath and positioned near the tapered tip. During the procedure, the stent graft is released by pressing the trigger and withdrawing the slider grip. The detailed components of the delivery system are listed in Table 2.


Figure 5: PV Delivery System

Table 2 Components of PV Delivery System

1	Tapered Tip	2	Radiopaque Marker Band	3	Sheath	4	Front Grip
5	Slider Grip	6	Luer Lock	7	Trigger		

1.2.3 G-Branch™ AE and G-Branch™ AAA Delivery System

The delivery system for both the AE and AAA stent grafts is the same. The delivery system primarily consists of a tapered tip, a hydrophilic coated sheath and a handle assembly. The stent graft is preloaded onto the delivery sheath and positioned near the tapered tip. The detailed components of the delivery system are listed in Table 3.


Figure 6: AE and AAA Delivery System

Table 3 Components of AE and AAA Delivery System

1	Tapered Tip	4	Front Grip	7	Proximal Releaser
2	Radiopaque Marker Band	5	Trigger	8	Safety Buckle
3	Sheath	6	Slider Grip	9	Haemostatic Valve

1.2.4 G-Branch™ IE Delivery System

The delivery system for the IE stent graft primarily consists of a tapered tip, a hydrophilic coated sheath and a handle assembly. The stent graft is preloaded onto the delivery sheath, positioned near the tapered tip. Deployment of the stent graft is achieved by withdrawing the sheath during the procedure. The detailed components of the delivery system are shown in Table 4.


Figure 7: IE Delivery System

Table 4 Components of IE Delivery System

1	Tapered Tip	4	Haemostatic Valve	7	Rear Grip
2	Radiopaque Marker Band	5	Front Grip	8	Rear Grip Locking Cap
3	Sheath	6	Front Grip Locking Cap	9	Luer Connector

1.3 Specifications and catalogue numbers

The G-Branch™ Thoracoabdominal Aortic Stent Graft System, the SilverFlow™ PV Peripheral Vascular Stent Graft System, and the Aortic Extension Stent Graft System which includes the G-Branch™ AE Main Body Extension Stent Graft System, the G-Branch™ AAA Bifurcated Stent Graft System, and the G-Branch™ IE Iliac Extension Stent Graft System, are available in the following lengths and diameters shown in Table 5 to Table 9.

Table 5 Dimensions of G-Branch™ Thoracoabdominal Aortic Stent Graft System

Unit: mm

Catalogue number ¹	Stent graft proximal diameter (D1)	Stent graft distal diameter (D2)	Internal branch diameter (d3)	External branch diameter (d4)	Stent graft length (L1)	Effective sheath length (L0)	Sheath outer diameter (D0)
TAAA24WWXX XeYiZZZZ	20	14, 16, 18, 20	6	6, 7, 8	110, 120, 130, 140, 150	680	22F
TAAA24WWXX XeYiZZZZ	22	14, 16, 18, 20	8				
TAAA24WWXX XeYiZZZZ	24	14, 16, 18, 20	8				
TAAA26WWXX XeYiZZZZ	26	14, 16, 18, 20	8, 10				
TAAA28WWXX XeYiZZZZ	28	14, 16, 18, 20	8, 10				
TAAA30WWXX XeYiZZZZ	30	14, 16, 18, 20, 22	8, 10				
TAAA32WWXX XeYiZZZZ	32	14, 16, 18, 20, 22	8, 10				
TAAA34WWXX XeYiZZZZ	34	14, 16, 18, 20, 22, 24, 26	8, 10				
TAAA36WWXX XeYiZZZZ	36	14, 16, 18, 20, 22, 24, 26	8, 10				
TAAA38WWXX XeYiZZZZ	38	14, 16, 18, 20, 22, 24, 26, 28, 30	8, 10				
TAAA40WWXX XeYiZZZZ	40	14, 16, 18, 20, 22, 24, 26, 28, 30	8, 10				

¹**TAAAVVWWXXXeYiZZZZ** is the G-Branch™ Thoracoabdominal Aortic Stent Graft System (“TAAA stent graft system”). “VV” is the proximal diameter of the stent graft (D1). “WW” is the distal diameter of the stent graft (D2). The “XXX” is the length of the stent graft (L1). “e” means the external branch. “Y” is external branch diameter (d4), “i” means the internal branch. “ZZZZ” is the internal branch diameter (d3).

Table 6 Dimensions of SilverFlow™ PV Peripheral Vascular Stent Graft System

Unit: mm

Catalogue number ²	Stent graft proximal diameter (D1)	Stent graft distal diameter (D2)	Stent graft length (L1)	Effective sheath length (L0)	Sheath outer diameter (D0)
PV-05WWYYYS/L/XL7	5	5	20, 30, 40, 50, 60, 70, 80, 90, 100, 110, 120	S/L/XL (800/1000/1200)	7F
PV-06WWYYYS/L/XL7	6	5, 6			
PV-07WWYYYS/L/XL7	7	5, 6, 7, 8			
PV-08WWYYYS/L/XL8	8	6, 7, 8, 9, 10, 12			8F
PV-09WWYYYS/L/XL8	9	6, 7, 8, 9, 10, 12			
PV-10WWYYYS/L/XL8	10	6, 7, 8, 9, 10, 12, 14			
PV-12WWYYYS/L/XL8	12	8, 9, 10, 12, 14			

* Stents with 20/30 mm length are only available with D1 = D2.

²**PV-VVWWYYYS/L/XLZ** is the SilverFlow™ PV Peripheral Vascular Stent Graft System ("PV stent graft system"). "VV" is the proximal diameter of the stent graft. "WW" is the distal diameter of stent graft. "YYY" is the covered length of stent graft. "S/L/XL" means the length of the delivery system, where S=800mm, L=1000mm, XL=1200mm. "Z" is the sheath outer diameter of the delivery system.

Table 7 Dimensions of G-Branch™ AE Main Body Extension Stent Graft System

Unit: mm

Catalogue number ³	Stent graft proximal diameter (D1)	Stent graft distal diameter (D2)	Stent graft length (L1)	Effective sheath length (L0)	Sheath outer diameter (D0)
AE16WWBYYYYZ	16	14, 16, 18, 20, 22, 24, 26, 28, 30	70, 80, 90, 100, 110, 120, 130, 140, 150, 160, 170, 180, 190, 200	760	18F
		32, 34, 36, 38, 40			20F/22F
AE18WWBYYYYZ	18	14, 16, 18, 20, 22, 24, 26, 28, 30	70, 80, 90, 100, 110, 120, 130, 140, 150, 160, 170, 180, 190, 200		18F
		32, 34, 36, 38, 40			20F/22F
AE20WWBYYYYZ	20	14, 16, 18, 20, 22, 24, 26, 28, 30	70, 80, 90, 100, 110, 120, 130, 140, 150, 160, 170, 180, 190, 200		18F
		32, 34, 36, 38, 40			20F/22F
		20			40, 60
AE22WWBYYYYZ	22	14, 16, 18, 20, 22, 24, 26, 28, 30	70, 80, 90, 100, 110, 120, 130, 140, 150, 160, 170, 180, 190, 200		18F
		32, 34, 36, 38, 40			20F/22F
		20, 22			40, 60

Catalogue number ³	Stent graft proximal diameter (D1)	Stent graft distal diameter (D2)	Stent graft length (L1)	Effective sheath length (L0)	Sheath outer diameter (D0)
AE24WWBYYYYZZ	24	14, 16, 18, 20, 22, 24, 26, 28, 30	70, 80, 90, 100, 110, 120, 130, 140, 150, 160, 170, 180, 190, 200		18F
		32, 34, 36, 38, 40			20F/22F
		24	40		18F
AE26WWBYYYYZZ	26	14, 16, 18, 20, 22, 24, 26, 28, 30	70, 80, 90, 100, 110, 120, 130, 140, 150, 160, 170, 180, 190, 200		18F
		32, 34, 36, 38, 40			20F/22F
		26	40		18F
AE28WWBYYYYZZ	28	14, 16, 18, 20, 22, 24, 26, 28, 30	70, 80, 90, 100, 110, 120, 130, 140, 150, 160, 170, 180, 190, 200		18F
		32, 34, 36, 38, 40			20F/22F
		28	40		18F
AE30WWBYYYYZZ	30	30	40		18F
AE32WWBYYYYZZ	32	32			20F/22F
AE34WWBYYYYZZ	34	34			
AE36WWBYYYYZZ	36	36			

³**AEVVWWBYYYYZZ** is the G-Branch™ AE Main Body Extension Stent Graft System (“AE stent graft system”). “VV” is the proximal diameter of the stent graft (D1). “WW” is the distal diameter of the stent graft (D2). “B” means the stent graft has a proximal bare stent. “YYY” is the covered length of the stent graft (L1). “ZZ” is the sheath outer diameter of the delivery system.

Table 8 Dimensions of G-Branch™ AAA Bifurcated Stent Graft System

Unit: mm

Catalogue number ⁴	Stent graft proximal diameter (D1)	Stent graft distal diameter (D2)	Stent graft short branch distal diameter (D3)	Stent graft length (L1)	Stent graft main body length (L2)	Stent short branch length (L3)	Effective sheath length (L0)	Sheath outer diameter (D0)
AAA18WWBYYYYZZ	18	09, 10, 12	09, 10, 12	120, 140, 160	50	30	520	18F
AAA20WWBYYYYZZ	20	10, 12	10, 12					
AAA22WWBYYYYZZ	22	10, 12	10, 12					
AAA24WWBYYYYZZ	24	12	12					
AAA26WWBYYYYZZ	26	12	12					
AAA28WWBYYYYZZ	28	12, 14	12, 14					
AAA30WWBYYYYZZ	30	12, 14	12, 14					
AAA32WWBYYYYZZ	32	12, 14	12, 14					
AAA34WWBYYYYZZ	34	12, 14	12, 14					20F/22F

Catalogue number ⁴	Stent graft proximal diameter (D1)	Stent graft distal diameter (D2)	Stent graft short branch distal diameter (D3)	Stent graft length (L1)	Stent graft main body length (L2)	Stent short branch length (L3)	Effective sheath length (L0)	Sheath outer diameter (D0)
AAA36WWBYYYYZ	36	12, 14	12, 14					

* AAA stent grafts are only available with D2= D3.

⁴**AAAVVWWBYYYYZ** is the G-Branch™ AAA Bifurcated Stent Graft System ("AAA stent graft system"). "VV" is the proximal diameter of the stent graft (D1). "WW" is the branch distal diameter of the stent graft (D2, D3). "B" means that the stent graft has a proximal bare stent. "YYY" is the covered length of the stent graft (L1). "ZZ" is the sheath outer diameter of the delivery system.

Table 9 Specifications of G-Branch™ IE Iliac Extension Stent Graft System

Unit: mm

Catalogue number ⁵	Stent graft proximal diameter (D1)	Stent graft distal diameter (D2)	Stent graft length (L1)	Effective sheath length (L0)	Sheath outer diameter (D0)
IE10WWCYYYYZ	10	08, 10, 12, 14	40, 60, 80, 90, 100, 120, 140, 160, 180	510	16F
IE10WWCYYYYZ	10	16, 18, 20	80, 90, 100, 120, 140, 160, 180		
IE12WWCYYYYZ	12	08, 10, 12, 14, 16	40, 60, 80, 90, 100, 120, 140, 160, 180		
IE12WWCYYYYZ	12	18, 20, 22, 24, 26, 28	80, 90, 100, 120, 140, 160, 180		
IE14WWCYYYYZ	14	10, 12, 14, 16, 18	40, 60, 80, 90, 100, 120, 140, 160, 180		
IE14WWCYYYYZ	14	20, 22, 24, 26, 28	80, 90, 100, 120, 140, 160, 180		
IE16WWCYYYYZ	16	12, 14, 16, 18, 20	40, 60, 80, 90, 100, 120, 140, 160, 180		
IE16WWCYYYYZ	16	22, 24, 26, 28	80, 90, 100, 120, 140, 160, 180		

⁵**IEVVWWCYYYYZ** is the G-Branch™ IE Iliac Extension Stent Graft System ("IE stent graft system"). "VV" represents the proximal diameter of the stent graft (D1). "WW" indicates the distal diameter of the stent graft (D2). "C" means that the stent graft does not have a proximal bare stent. "YYY" is the covered length of the stent graft (L1). "ZZ" is the outer diameter of the sheath of the delivery system.

II Intended use

The Thoracoabdominal Artery Stent Graft System is intended for the endovascular repair of thoracoabdominal aortic aneurysms and aortic dissections.

2.1 Indications

The Thoracoabdominal Artery Stent Graft System is indicated for the endovascular repair of thoracoabdominal aortic aneurysms and aortic dissections.

- The G-Branch™ TAAA Stent Graft System is used for endovascular repair of the thoracoabdominal aorta and aortic dissections in patients requiring:
 - Adequate femoral, iliac, and upper extremity vessel access compatible with the required delivery system.
 - The proximal landing zone of the thoracoabdominal aorta should have a diameter ranging from 17 mm to 36 mm, with a minimum landing length of 25 mm.
- The SilverFlow™ PV Stent Graft System is used for the endovascular repair of the celiac trunk, superior mesenteric artery, and renal arteries in patients requiring:
 - Adequate femoral, iliac, and upper extremity vessel access compatible with the required delivery system.
 - The landing zone of the visceral branch should have a diameter ranging from 6 mm to 13 mm, with a minimum landing length of 15 mm.
 - The landing zone of the renal artery should have a diameter ranging from 4.5 mm to 9 mm, with a minimum landing length of 15 mm.
- The G-Branch™ AAA and IE Stent Graft Systems are used for the endovascular repair of the abdominal aorta and iliac arteries in patients requiring:
 - Adequate femoral, iliac, and upper extremity vessel access compatible with the required delivery system.
 - The distal landing zone of the iliac artery should have a diameter ranging from 7 mm to 25 mm, with a minimum landing length of 15 mm.
- The G-Branch™ AE Stent Graft System is used for extending the proximal and distal landing zones in thoracoabdominal aortic aneurysms and aortic dissections in patients requiring:
 - Adequate femoral, iliac, and upper extremity vessel access compatible with the required delivery system.

- The distal landing zone of the abdominal aorta should have a diameter ranging from 12 mm to 36 mm, with a minimum landing length of 20 mm.

2.2 Intended users

The Thoracoabdominal Artery Stent Graft System shall be used only by physicians and clinical teams trained in vascular interventional techniques, including specific training in the use of the device.

2.3 Intended patient population

The Thoracoabdominal Artery Stent Graft System is intended for any patient requiring treatment of thoracoabdominal aortic aneurysms or aortic dissections with anatomic characteristics and vessel morphology compatible with the device.

2.4 Clinical benefit

The Thoracoabdominal Artery Stent Graft System, designed for endovascular treatment of thoracoabdominal aortic aneurysms (TAAs) and aortic dissections, offers significant clinical benefits by enabling effective aneurysm exclusion, restoring physiological blood flow, and preserving visceral perfusion through a minimally invasive approach. It also provides an endovascular option for patients with TAA or dissection who are at high risk of open surgical repair, which is associated with high morbidity and mortality rates.

III Contraindications

3.1 The Thoracoabdominal Artery Stent Graft System is contraindicated in:

- Patients whose anatomy is not compatible with available device dimensions.
- Patients with acute systemic infections.
- Patients with another implanted device that interferes with its use.
- Patients with contraindications to X-ray imaging.
- Patients with known hypersensitivity or allergy to stent graft materials.
- Patients who cannot tolerate contrast agents due to renal insufficiency.
- Patients allergic to the contrast medium.
- Patients with severe stenosis, calcification, or mural thrombosis at the planned landing zone that may hinder proper stent graft apposition.

3.2 The safety and effectiveness of the Thoracoabdominal Artery Stent Graft System have not been

evaluated in:

- Patients with acute rupture of thoracoabdominal aortic aneurysms.
- Patients with connective tissue disorders.
- Patients in a hypercoagulable state.
- Patients with concurrent mesenteric artery occlusive disease.
- Pregnant or nursing women.
- Patients <18 years.

IV Warnings

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- This device shall be used only by physicians and teams trained in vascular interventional techniques, including training in the use of the device.
- Potential device embolization may occur, and physicians need to be prepared to remove the embolized device.
- Do not continue advancing if resistance is felt during advancement of the guidewire and delivery system. Stop and assess the cause of resistance.
- Allergic reactions may occur in patients allergic to nitinol following the implantation of the stent graft.
- Do not use if the package is opened or damaged, or it is suspected that the sterility of the device has been compromised, as infection and related serious potential patient harms could occur including infection.
- For single use only. DO NOT REUSE, REPROCESS OR RESTERILIZE. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death.
- Do not use expired product, if any component of the device shows abnormalities such as damage, do not use the device.
- If a serious incident related to the device occurs, immediately report the incident to Lifetech and the applicable competent authority or regulatory body.

V Precautions

5.1 MRI Information:

Non-clinical testing demonstrated that the Thoracoabdominal Artery Stent Graft is MR Conditional. A patient implanted with this device may be safely scanned in an MR system meeting all of the following conditions:

- Static magnetic field of 1.5 T or 3.0 T
- Maximum spatial magnetic field gradient of 1,700 gauss/cm at static magnetic field of 3.0 T and 2,000 gauss/cm at static magnetic field of 1.5 T
- Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operation Mode) for maximum scan duration of 15 minutes.

Under the defined scan conditions, the device is expected to produce a maximum temperature rise of less than 2.63 °C.

In non-clinical testing, the largest image artifact measured was 6.5 mm (GRE sequence) at 1.5T and 7.5 mm (GRE sequence) at 3.0T. This image artifact does not obscure the device lumen.

Note: MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters to compensate for the presence of this implant.

5.2 Other matters

- This device should be stored in a non-corrosive gas environment, well-ventilated, and clean room.
- The device must not come into contact with organic solvents.
- During transportation of the device, crushing, wetting, or exposure to sunlight is strictly prohibited. Handle with care and gently place it down.
- If any questions arise before using this device, contact the distributor or manufacturer directly.
- Physicians should determine the use of anticoagulants or antiplatelet drugs before, during, and after the use of this device based on clinical conditions.
- Ensure that the size of the access vessel is adequate for the selected delivery sheath.
- The device should be used via standard interventional cardiovascular catheter techniques.

- The sheath's effective length surface is treated with a hydrophilic coating. Avoid contact with sharp objects to prevent coating detachment.
- To activate the hydrophilic coating on the sheath surface, the surface shall be wiped with sterile gauze pads moistened with saline solution.

5.3 Overall qualitative and quantitative information on the materials and substances to which patients can be exposed

- The materials and substances to which the patient can be exposed in the G-Branch™ Thoracoabdominal Aortic Stent Graft (in wt.%) are NiTi (Nickel titanium, 49-67%), stainless steel (7-17%), PET (Polyethylene terephthalate, 20-40%), TA (tantalum, <1%). During the intended use, the thoracic aorta, abdominal aorta, iliac artery, femoral artery as well the blood will come in contact with the device.
- The materials and substances to which the patient can be exposed in the SilverFlow™ PV Peripheral Vascular Stent Graft (in wt.%) are NiTi (Nickel titanium, 49-67%), PTFE (Polytetrafluoroethylene, 15-35%), TA (tantalum, <1%). During the intended use, the thoracic aorta, abdominal aorta, iliac artery, femoral artery as well the blood will come in contact with the device.
- The materials and substances to which the patient can be exposed in the Aortic Extension Stent Graft (in wt.%) are NiTi (Nickel titanium, 49-67%), stainless steel (7-17%), PTFE (Polytetrafluoroethylene, 20-40%), TA (tantalum, <1%). During the intended use, the thoracic aorta, abdominal aorta, iliac artery, femoral artery as well the blood will come in contact with the device.

VI Complications

Possible adverse events and complications that may occur with the use of the Thoracoabdominal Artery Stent Graft System or in any endovascular repair procedure, which may or may not require intervention include, but are not limited to:

- | | |
|--|--|
| ■ Access site Complication | ■ Hemothorax |
| ■ Anaphylaxis | ■ Hepatic Failure |
| ■ Anemia | ■ Hypertension/Hypotension |
| ■ Anesthetic Complications | ■ Infection |
| ■ Aneurysm Enlargement/Rupture | ■ Inflammation |
| ■ Aortic dissection or intramural hematoma | ■ Kidney atrophy/Infarction |
| ■ Arrhythmia | ■ Lymph Fistula |
| ■ Arteriovenous Fistula | ■ Myocardial Infarction |
| ■ Bleeding complication | ■ Neurologic Local or Systemic Complications (e.g., Stroke, Transient Ischemic Attack, |

- | | |
|--|--|
| <ul style="list-style-type: none"> ■ Cardiac Arrest ■ Cardiac Tamponade ■ Claudication (e.g., Buttock and Lower Limb) ■ Component Separation ■ Congestive Heart Failure | <ul style="list-style-type: none"> ■ Paraplegia, Paraparesis, Spinal Cord Shock, Paralysis) ■ Occlusion of branch vessel ■ Organ Failure ■ Pericardial Effusion ■ Pseudoaneurysm ■ Pulmonary Complications (e.g., Atelectasis, Pneumonia, Respiratory Failure) |
| <ul style="list-style-type: none"> ■ Death ■ Device Migration ■ Edema ■ Embolism ■ Erosion ■ Fever | <ul style="list-style-type: none"> ■ Renal artery stenosis ■ Renal Failure/Dysfunction ■ Reoperation/Reintervention ■ Retrograde dissection ■ Sepsis ■ Splenic Injury (e.g., Infarction, Ischemia) ■ Thrombosis |
| <ul style="list-style-type: none"> ■ Gastrointestinal Complications (e.g., Intestinal Ischemia/Infarction/Obstruction/Ischemia/Necrosis) ■ Hemorrhage | <ul style="list-style-type: none"> ■ Vascular Trauma/Injury (e.g., Vessel Dissection, Perforation, Rupture) |

WARNING: Even if the procedure is performed correctly, the aforementioned complications may still occur. Therefore, it is crucial to rigorously select patients, carefully assess their specific conditions, and determine the appropriate pre-operative planning measurements.

VII Operation guide

Note: This guideline refers to an example of endovascular repair of a type of Crawford III thoracoabdominal aortic aneurysm with the Thoracoabdominal Artery Stent Graft System.

7.1 Preparation of equipment and materials (Not included in this system, required separately)

- 0.035" Super Stiff Guidewire (260cm), at least two pieces, such as Cook TSMG-35-260-LES Lunderquist® Extra-Stiff Wire Guide.
- 0.035" Stiff Guidewire (260cm), at least two pieces, such as Cook THSCF-35-260-1.5-Rosen Wire Guide.
- 0.035" Hydrophilic Soft-tipped Guidewire (260cm), at least two pieces, such as Terumo RF*GA35263M.
- 8F long sheath, at least two sets, such as Cook KCFW -8.0 -35-70-RB-HFANL0-HC Flexor® Ansel

Guiding Sheath.

- Endovascular snare system, such as Lifetech LT-SG-20 SeQure™ Snare System.
- Steerable introducer, such as Lifetech SVA8F -700 Fustar™.
- Guiding catheter such as Cordis 6F JR4.
- 0.035" Balloon Dilatation Catheter, compatible with 0.035" guidewire.
- 0.018" Balloon Dilatation Catheter compatible with 0.018" guidewire.
- Standard surgical instruments for exposure and closure of iliac, femoral, and brachial arteries or vascular closure devices.
- Arterial puncture devices.
- Various angiographic catheters.
- Contrast media.
- Digital subtraction angiography capable fluoroscopy equipment.
- Power injector.
- Heparin saline solution.
- Blood pressure monitoring kit.
- Triple lumen stopcock.
- Any other instrument that the operating physician deems necessary for the surgery.

7.2 Patient preparation

- Measure the length and aortic diameter of the proximal landing zone.
- Measure the length and aortic diameter of the distal landing zone.
- Measure the diameter and length of the aneurysm.
- Measure the diameter of the celiac trunk, superior mesenteric artery, and left and right renal arteries.
- Measure the distance from each branch's opening to the landing zone.
- Measure the length from each branch's opening to the opening of the branch's accessory arteries.
- Measure the diameter of the opening of the brachial artery (axillary artery).
- Measure the angle of the proximal neck of the aneurysm.

- Determine that the patient's vascular anatomy is suitable for endovascular repair.
- Perform various necessary routine examinations and monitor vital signs related to anesthesia, anticoagulation, and contrast administration.
- Ensure adequate vascular access via the femoral arteries, brachial arteries, etc., using standard surgical exposure of both femoral arteries (iliac arteries), left brachial artery (left axillary artery), left common carotid artery, and left subclavian artery.
- Prepare for systemic anticoagulation (if the patient is allergic to heparin, choose an appropriate alternative medication).

7.3 Product selection and preparation

Based on the measurements, select the appropriate TAAA Stent Graft, the Aortic Extension Stent Grafts, the PV Stent Grafts, and the appropriate sheath. The stent length should be as long as possible to match the anatomy of the patient's thoracoabdominal aorta and its branches. It is recommended that the diameter of the TAAA stent graft and the Aortic Extension Stent Graft should be oversized by 10-20% of the corresponding aortic diameter, and that the diameter of the PV stent graft should be oversized by 5-15% of the corresponding artery diameter.

7.4 Operating Instructions

7.4.1 Preparation before implantation

- **Unpack the Devices:**

Under sterile conditions, remove the delivery system from its package. For the G-Branch™ Thoracoabdominal Aortic Stent Graft System, carefully remove the protective tube from the coiled portion of the preloaded guidewires at the proximal end of the delivery system, as illustrated in Figure 9(a).

CAUTION: Take care to prevent the protective tube from falling off when removing the delivery system from its packaging. After removing the protective tube, gently wipe the surface of the preloaded guidewires with sterile heparinized saline.


Figure 9(a): Uncoiling the Preloaded Guidewire and Removing the Protective Tube from the TAAA Delivery System

- **Prepare the G-Branch™ Thoracoabdominal Aortic Stent Graft System and the Aortic Extension Stent Graft System:**

For the G-Branch™ TAAA Delivery System and the Aortic Extension Delivery System, attach syringes to the haemostatic valve and luer lock, respectively. Thoroughly flush both the sheath lumen and the guidewire lumen with sterile heparinized saline, as illustrated in Figure 9(b) (TAAA delivery system as an example).


Figure 9(b): Flushing the sheath lumen and the guidewire lumen of the TAAA Delivery System

- **Prepare the SilverFlow™ PV Peripheral Vascular Stent Graft System:**

For the SilverFlow™ PV Peripheral Vascular Stent Graft System, after unpackaging, flush the delivery system from the luer lock at the distal end using a syringe, as shown in Figure 9(c).


Figure 9(c): Flushing the PV Delivery System

- Wipe the outer surface of the delivery system sheaths with sterile heparinized saline to lubricate the sheath surface and activate the hydrophilic coating on the sheaths of the TAAA delivery

system, the PV delivery system, the AE and AAA delivery system, and the IE delivery system, as shown in Figure 10.


Figure 10: Wiping the sheath with sterile heparinized saline

- Under DSA guidance, access the upper limb artery and advance two 8F long sheaths to the descending aorta segment. Remove the dilator for standby use. Introduce a 0.035" super stiff guidewire through the femoral artery opening to the segment of the descending aorta above the celiac trunk ostium, as shown in Figure 11. (For reference only)


Figure 11: Introducing a super stiff guidewire to the specified position

CAUTION:

Maintain cleanliness of the surgical field:

Before beginning the implantation of the stent graft, ensure that the surgical table is kept clean and free from any instruments or materials (such as surgical clamps, gauze, etc.) that could interfere with the introduction, deployment, and retrieval operations of the stent graft system.

Operation under DSA guidance:

Once the device enters the body, all procedures must be carried out cautiously under the guidance of DSA equipment. In case of encountering abnormal resistance, immediately stop the operation and only proceed after any issues have been resolved.

7.4.2 Implantation of the TAAA Stent Graft (partial deployment)

- **Introduce the TAAA Delivery System:**

Ensure that the haemostatic valve at the proximal end of the delivery system is facing away from the operator.

Advance the TAAA delivery system over the super stiff guidewire while gently pushing and adjusting as needed, and monitor the "G" radiopaque marker (Figure 1, component③) of the TAAA stent graft under DSA guidance.

Stop advancing when the "G" radiopaque marker aligns with the celiac trunk ostium, as illustrated in Figure 12.


Figure 12: Advancing the TAAA Delivery System to the required position

- **Partially release the proximal segment:**

Under DSA guidance, hold the front grip steady with the left hand while rotating the slider grip counterclockwise with the right hand to slowly release the bare stent and the proximal segment of

the TAAA stent graft.

Stop releasing once the proximal graft covered segment (including the "O" and "8" radiopaque markers shown in Figure 1, component①&②), is fully released from the sheath, as illustrated in Figure 13.


Figure 13: Releasing the proximal segment of the TAAA Stent Graft

- **Confirm orientation and adjust axial position:**

Use the "O" and "8" radiopaque markers at the proximal graft covered end to confirm correct orientation of the internal branch portals toward the celiac trunk (CA) ostium.

While holding the front grip of the delivery system, rotate or gently adjust the delivery system axially, ensuring that the "G" marker remains aligned with the celiac trunk (CA) ostium.

As illustrated in Figure 1, the "O" marker (left) and the "8" marker (right) indicate that the internal branch outlets are oriented anteriorly (i.e., facing upward under DSA view), corresponding to the direction of the visceral arteries.

The "8" and "O" radiopaque markers help verify that the distal portals of the internal branches are directed toward the visceral arteries. Verify both the axial position and the rotational orientation of the stent graft before proceeding to the next step.

- **Semi-release the proximal end of the TAAA Stent Graft:**

Hold the front grip of the delivery system with the left hand. Rotate the slider grip counterclockwise with the right hand to gradually release the stent graft under DSA guidance.

Stop releasing once the "G" marker at the distal graft covered segment becomes visible and the radiopaque marker band at the distal end of the external branches is also fully visible, as shown in

Figure 14.


Figure 14: Semi-released state of the TAAA Stent Graft

At this stage, the proximal segment of the TAAA stent graft is a semi-released state, with approximately 70–80% expansion of its proximal diameter. The constraining wire remains in place to allow for fine circumferential and axial adjustment.

Perform any necessary positional corrections by carefully manipulating the delivery system handle.

Once optimal placement is confirmed, proceed to establish the branch access pathway for celiac trunk and superior mesenteric artery.

7.4.3 Establishment of branch access pathway for the celiac trunk and superior mesenteric artery

- **Advance the preloaded guidewire:**

Hold the front grip of the delivery system. At the proximal end of delivery system, rotate one guidewire locker approximately half a turn counterclockwise to unlock the corresponding preloaded guidewire.

Advance the preloaded guidewire gently along the longitudinal axis of the sheath while monitoring its position within the aortic lumen under DSA guidance. Repeat the same procedure to advance both preloaded guidewires into the descending aorta, ensuring that their distal tips remain slightly offset longitudinally (i.e., not positioned at the same level), as shown in Figure 15.

Once both preloaded guidewires are properly positioned, rotate the guidewire lockers clockwise to securely clock them in place.


Figure 15: Advancement and positioning of the preloaded guidewires

CAUTION: Before advancing the preloaded guidewires, ensure that they have been properly unlocked, and the delivery system is firmly secured throughout the advancement to maintain stability.

Operational Tip: The numbers ① and ② in the figure indicate the sequence of operations.

- **Capture the preloaded guidewire:**

Introduce the snare system into the long sheath, and expose the snare from the sheath. Unlock the preloaded guidewire, and slowly advance it back and forth toward the snare. An assistant controls the snare to capture the soft tip of the preloaded guidewire. Once the soft tip enters the snare, as shown in Figure 16 (a), gently retract the snare to pull the soft tip fully into the snare catheter. At this point, a noticeable increase in resistance should be felt, indicating that the preloaded guidewire has been successfully captured.

Retrieve the preloaded guidewire from the body through the long sheath, and stop when approximately 30-50 cm of the guidewire remains outside the body near the guidewire locker. Then, rotate the locker clockwise to secure the preloaded guidewire.

Repeat the procedure with the second preloaded guidewire using another snare system, as shown in Figure 16 (b).


Figure 16: Capturing the preloaded guidewires

CAUTION: When pulling the preloaded guidewire outward, an assistant is required to push from the end near the guidewire locker. The assistant should continuously monitor the length of the guidewire near the locker in real time to prevent the guidewire from being completely pulled into the delivery system.

- **Advance the long sheath:**

Advance the two long sheaths respectively along the preloaded guidewires into the interior of the TAAA stent graft slowly. The long sheaths pass through the bare stent and internal branches. When the long sheaths are exposed from the outlet of the internal branches, advance them 5 mm to 10 mm and stop, as shown in Figure 17.


Figure 17: Advancing the long sheaths

Operational Tip: During advancement of the long sheath into the interior of the TAAA stent graft, it is recommended to use the dilator for additional support. If the tip of the dilator becomes engaged with or cannot pass the bare stent portion, the dilator may be replaced with a 0.018" balloon catheter (compatible with the preloaded 0.018" guidewire) to perform dilation. The balloon can facilitate advancement of the long sheath into the internal branch.

CAUTION: Obstruction by the bare stent may occur during the sheath advancement. Use a balloon to carefully traverse and dilate the bare stent, and advance the long sheath gradually during deflation.

- **Exchange the guidewires:**

Introduce the guiding catheter through the long sheath alongside the preloaded guidewire.

Insert a hydrophilic soft-tipped guidewire through the guiding catheter, adjusting the direction of the guiding catheter to cannulate celiac trunk (or SMA), as shown in Figure 18(a).

Slowly adjusting the direction of the soft tip of the guidewire to direct it into the branch artery. Once the guidewire is in place, advance the guiding catheter into the branch artery. Then, withdraw the hydrophilic soft-tipped guidewire, exchange it with a stiff guidewire.

Unlock the preloaded guidewire, pull it completely out of the body through the long sheath.

Advance the long sheath to cannulate the celiac trunk (or SMA), then remove the guiding catheter, as shown in Figure 18(b).


Figure 18: Establishment of branch access pathway (CA/SMA)

Operational Tip 1: An angiographic catheter can be used to perform angiography of branch arteries (e.g., the CA and SMA) before introducing the guiding catheter.

Operational Tip 2: If advancing the guidewire into the target branch through the guiding catheter is technically difficult, an appropriate steerable sheath may be used instead to cannulate the branch.

- **Release the proximal constraining wire of the TAAA Stent Graft:**

After completing the guidewire exchange, reconfirm the position of the TAAA stent graft. Once the position is confirmed, gently press the release lock at the proximal end of the TAAA delivery system with your thumb, and slowly withdraw the constraining wire from the proximal end of the delivery system to fully release and expand the proximal graft covered segment of the TAAA stent graft, as shown in Figure 19.


Figure 19: Releasing the proximal graft covered segment of the TAAA Stent Graft

CAUTION: Pull the constraining wire slowly and in a direction parallel to the axis of the delivery system.

7.4.4 Implantation of PV Stent Graft into the celiac trunk and superior mesenteric artery

- **Implant a PV Stent Graft in the celiac trunk:**

Introduce the PV delivery system along the stiff guidewire into the target branch artery, then slightly withdraw the long sheath until its distal end is retracted of the internal branch of the TAAA stent graft.

Adjust the position of the PV delivery system so that the two "■" radiopaque markers at the proximal end of the PV stent graft are aligned with the lower edge of the proximal radiopaque marker band of the internal branch of the TAAA stent graft.

Hold the front grip of the PV delivery system with the left hand to stabilize it. Hold the slider grip with the right hand, press the trigger, and slide the grip backward to deploy the PV stent graft.

The operation method and the fully released stent state are shown in Figure 20 (right) and Figure 20 (left), respectively.

Observe the deployment under DSA guidance until the stent is fully released.


Figure 20: Releasing the PV Stent Graft

CAUTION 1: Ensure that the long sheath is withdrawn to a position above the proximal opening of the internal branch before releasing the PV stent graft. Failure to do so may result in partial deployment of the proximal portion of the PV stent graft within the long sheath, inability to release the stent, or unintended stent migration.

CAUTION 2: Keep the sheath portion exposed outside of the body straight at the beginning of stent deployment.

- **Remove the PV Delivery System:**

Press the trigger and hold the slider grip with right hand. Hold the front grip with the left hand, and slowly retract the front grip toward the slider grip, thereby retracting the tapered tip assembly, until

the tapered tip is aligned with the distal end of the sheath.

Release the trigger to lock the delivery system. Keep the long sheath in position, and withdraw only the PV delivery system, leaving the stiff guidewire in place.

Post-dilation with an appropriately sized balloon may be performed over the stiff guidewire to ensure optimal lumen patency.

Angiography should be performed to confirm the morphology and blood flow of the deployed PV stent graft.

Operational Tip: The stiff guidewire may be left in place until any required balloon post-dilation is completed.

CAUTION 1: When withdrawing the tapered tip assembly, perform the action slowly and gently to prevent the tapered tip from catching on the stent, and causing migration.

CAUTION 2: If resistance is encountered when retracting the tapered tip—suggesting that the tip may have caught on the stent—one or more of the following approaches may be attempted:

- 1) Retract the guidewire until it clears the contact point between the tapered tip and the stent, then slowly retract the front grip to allow the tapered tip to pass smoothly through the stent lumen.
- 2) Slightly rotate the delivery system before retracting.
- 3) Alternatively, gently advance the entire delivery system before retracting.

Do not pull the delivery system forcefully, as this may cause stent migration or sheath damage.

- **Implant a PV Stent Graft in the superior mesenteric artery:**

Following the same procedure used for the celiac trunk revascularization, advance and deploy the PV stent graft in the superior mesenteric artery. The PV delivery system may be withdrawn as described above.

Post-dilation with an appropriate size balloon may be performed over the stiff guidewire to ensure optimal lumen patency.

Angiography should be performed to confirm the morphology and blood flow of the deployed PV stent graft.


Figure 21: Implantation of the PV Stent Graft in the CA and SMA

After completion of the PV stent graft implantation into the celiac trunk and the superior mesenteric artery, proceed to establish the access pathway for the renal arteries.

7.4.5 Establishment of branch access pathway for the renal arteries

- **Retract the long sheath:**

Slowly retract the two long sheaths away from the proximal end of the TAAA stent graft under DSA guidance until their distal ends are positioned above the tapered tip of the TAAA delivery system (see Figure 22).


Figure 22: Retracting the long sheath

- **Release the bare stent of the TAAA Stent Graft:**

Hold the TAAA delivery system with the left hand. Remove the safety buckle at the proximal end of the delivery handle using the right hand (press the safety buckle in the direction indicated by arrow ①).

Gently pull the proximal releaser backward through its full stroke (pinch the proximal releaser between the thumb and index finger and pull in the direction indicated by arrow ②) while monitoring the deployment of the proximal bare stent under fluoroscopy (Figure 23, left).

When deployment is complete, the stent configuration should match Figure 23 (right).


Figure 23: Releasing the bare stent of the TAAA Stent Graft

- **Advance the long sheath through the external branch channel:**

Introduce a guiding catheter through the long sheath and advance it to the proximal inlet of the target external branch. Orient the guiding catheter tip orientation so it aligns with the external branch outlet and secure the guiding catheter in position.

Advance a hydrophilic soft-tipped guidewire through the guiding catheter into the external branch until the guidewire exits the distal outlet of the branch. When the guidewire tip is approximately 5–10 mm above the renal artery ostium, secure the guidewire in position.

Advance the guiding catheter over the guidewire to the same level, then advance the long sheath through the external branch to a position 5–10 mm above the renal artery ostium. Confirm position under fluoroscopy (see Figure 24).

- **Exchange guidewires to establish the renal artery access:**

Position the guiding catheter tip flush with and oriented toward the renal artery ostium. Manipulate the hydrophilic soft-tipped guidewire to cannulate the renal artery and advance the guiding catheter into the renal artery.

Withdraw the hydrophilic soft-tipped guidewire and exchange it for a stiff guidewire positioned within the renal artery. Withdraw the guiding catheter, leaving the stiff guidewire in place as the established access for subsequent PV stent graft delivery (Figure 24).


Figure 24: Establishment of the PV Stent Graft access pathways (renal)

Operational tip: If difficulty is encountered when cannulating the renal artery, consider using an angiographic catheter for local angiography to clarify anatomy or a steerable sheath to facilitate access.

7.4.6 Implantation of PV Stent Graft into the first renal artery

Introduce the PV delivery system over the stiff guidewire into the renal artery.

Align the two proximal "•" radiopaque markers of the PV stent graft with the upper edge of the proximal radiopaque marker band of the corresponding TAAA stent graft external branch.

Hold the front grip steady with the left hand. With the right hand, press the trigger and pull the slider grip to deploy the PV stent graft. Observe the deployment under DSA guidance until the PV stent graft is fully released.

Retract the tapered tip into the sheath and remove the PV delivery system. Leave the long sheath and, if required, the stiff guidewire in place for post-dilation.

Operational tip: The stiff guidewire may be left in place until any required balloon post-dilation is completed.

Post-implantation check: Perform angiography to confirm stent morphology and blood flow. If indicated, perform balloon post-dilation over the stiff guidewire to optimize lumen patency.


Figure 25: Schematic diagram of the first renal artery revascularization

CAUTION: Ensure the long sheath is withdrawn above the internal branch proximal opening before releasing the PV stent graft. Failure to do so may result in partial deployment inside the sheath, inability to release, or unintended stent migration.

7.4.7 Fully release the TAAA stent graft (distal release)

After revascularization of the CA, SMA and one renal artery, position the PV delivery system for the contralateral renal artery as required. Then fully release the distal portion of the TAAA stent graft to restore full apposition.

Hold the front grip steady with the left hand. Operate the release mechanism with the right hand (pull the trigger and pull the slider) to rapidly release the distal graft covered segments, as illustrated in Figure 26. Confirm the deployed configuration under fluoroscopy.


Figure 26: Releasing the distal end of TAAA Stent Graft

7.4.8 Implantation of a PV Stent Graft into the other renal artery

Refer to Section 7.4.6 for PV stent graft implantation steps. Advance and deploy the PV stent graft into the contralateral renal artery using the same technique. Perform angiography and balloon post-dilation as required. See Figure 27 for the final bilateral renal revascularization configuration.


Figure 27: Schematic diagram of the both renal artery revascularization

Operational tip: If post-dilation is required, keep the stiff guidewire in place until dilation is complete; then withdraw the guidewire and sheath.

7.4.9 Removing the TAAA Delivery System

- **Retract the tapered tip into the delivery system sheath:**

Hold the slider grip with the right hand to stabilize the system, and the front grip with the left hand. Pull the trigger backward, then slowly retract the front grip toward the slider grip to withdraw the tapered tip into the distal end of the sheath (Figure 28).


Figure 28: Retracting the tapered tip of the TAAA Delivery System

CAUTION: If resistance is encountered and the tapered tip appears to catch on the stent, do not apply excessive force. Slightly rotate the delivery system and attempt gentle retraction. If the issue persists, partially retract guidewires to alter geometry and retry removal carefully.

- **Remove the TAAA Delivery System:**

Carefully withdraw the TAAA delivery system through the femoral access site. Perform angiography to confirm stent morphology and blood flow. Upon satisfactory confirmation, remove the angiographic catheter.

Operational Tip: If subsequent implantation of AE, AAA or IE stent grafts is planned, the super stiff guidewire may be left in place to maintain vascular access and minimize the need for repeated vascular access.

7.4.10 Implantation of the AE Stent Graft (when distal landing zone extension is required)

- **Introduce the AE Delivery System:**

Ensure the haemostatic valve faces away from the operator, then introduce the AE delivery system along the super stiff guidewire from the femoral access.

Slowly advance the tapered tip of the AE delivery system into the distal inner portion of the TAAA stent graft under DSA guidance (Figure 29).


Figure 29: Introducing the AE Delivery System

- **Release the proximal end of AE Stent Graft:**

Hold the front grip steady with the left hand and rotate the slider grip counterclockwise with the right hand to slowly release the bare stent and the proximal graft covered end.

Stop when the proximal graft covered segment (including the "O" and "8" radiopaque markers) is completely released from the sheath (Figure 30).


Figure 30: Releasing the proximal end of AE Stent Graft

Ensure that the proximal radiopaque marker ("O" or "8") of the AE stent graft axially overlaps the distal "O" radiopaque marker of the TAAA stent graft by at least 30 mm and that the overlapping length between the distal segment of the AE stent graft and abdominal aorta is at least 20 mm.

- **Release the fully graft covered segment rapidly:**

Hold the front grip steady. Pull the slider grip with the right hand while simultaneously pulling the trigger backward with the right thumb to rapidly release the distal graft covered segment (Figure 31).


Figure 31: AE distal segment rapid release

- **Release the proximal bare stent:**

Remove the safety buckle at the proximal handle, then gently pull the proximal releaser backward to its full stroke, while monitoring under fluoroscopy (Figure 32). When deployment is complete, confirm final state.


Figure 32: AE proximal bare stent release

- **Remove the AE Delivery System:**

Retract the tapered tip into the sheath and withdraw the AE delivery system (Figure 33). Perform angiography to confirm matching and, if required, post-dilate with an appropriate size balloon.


Figure 33: AE Delivery System withdrawal and final state

Note: AE and AAA delivery systems share the same operating mechanics; see Section 7.4.11 for AAA-specific short-branch steps.

7.4.11 Implantation of the AAA Stent Graft and the IE Stent Graft(s) (reconstruction of abdominal

aorta)

- **Introduce and position the AAA Delivery System:**

Introduce the AAA delivery system along the super stiff guidewire and advance the tapered tip into the inner distal portion of the TAAA stent graft.

- **Release the proximal graft covered segment and short branch:**

Slowly release the proximal graft covered segment (O and 8 markers). Ensure axial overlap between the AAA “8” marker and the distal TAAA “O” marker (≥ 30 mm). Adjust rotation slightly if necessary. Continue to release until the short branch is deployed (Figure 34).


Figure 34: Short branch deployment of AAA Stent Graft

CAUTION: Ensure the AAA stent graft and the tapered tip do not interfere with the existing TAAA stent graft to avoid migration during advancement and positioning.

- **Release the bare stent of the AAA Stent Graft:**

Remove the safety buckle and pull the proximal releaser to its full stroke to disengage the bare stent. If deployment is incomplete, perform gentle advance/retraction of the proximal releaser as needed (Figure 35).


Figure 35: AAA proximal bare stent release

- **Implant the IE Stent Graft into the short branch:**

From the contralateral femoral approach, introduce another super stiff guidewire through the short branch and advance the IE delivery system until the IE "O" marker aligns with the short-branch "e" marker (Figures 36). Confirm fluoroscopic alignment before deployment.


Figure 36: Introducing the IE Stent Graft to the short branch

Rotate the front grip locking cap counterclockwise approximately half a turn to release the rear grip lock while holding the front grip.

Secure the rear grip with the right hand and slowly retract the front grip toward the rear grip with the left hand until the IE stent graft is fully deployed under DSA guidance (Figure 37, left).

Hold the front grip with the left hand and retract the rear grip with the right hand until the tapered tip is fully within the sheath. Rotate the front grip locking cap clockwise to lock the delivery system, and remove the IE delivery system from the body (Figure 37, right).


Figure 37: Releasing the IE Stent Graft and removing the IE Delivery System

- **Fully release the AAA Stent Graft:**

Release the long branch to full deploy the AAA stent graft and remove the AAA delivery system from the body (Figures 38).


Figure 38: Releasing the long branch of the AAA Stent Graft and removing the AAA Delivery System

- **Implant the IE Stent Graft into the long branch:**

From the ipsilateral femoral artery approach, Introduce the IE delivery system into the long branch, and adjust the position to that the "O" radiopaque marker of the IE stent graft is aligned with the "O" radiopaque marker in the middle of the long branch (Figure 39).


Figure 39: Introducing IE Delivery System into the long branch

Repeat the same deployment procedure as for the contralateral IE stent graft, including unlocking, gradual retraction of the front grip for full deployment under fluoroscopy, and removal of the IE

delivery system. If apposition is suboptimal, post-dilate with an appropriate size balloon. The final state after implantation is shown in Figure 40.


Figure 40: Final state after implantation of AAA and IE Stent Grafts

7.4.12 Angiographic Confirmation and Post-procedure

Perform angiography of the entire Thoracoabdominal Artery Stent Grafts using an angiographic catheter to confirm patency and proper deployment. After angiography, withdraw the catheter and guidewire, and close all vascular access sites according to standard surgical procedures.

VIII Troubleshooting Techniques

8.1 Alternative Instruction for Delivery System Withdrawal

In the unlikely event of delivery system resistance or entrapment during withdrawal, the following alternative steps may be performed:

- 1) Verify that all guidewires and sheaths are free from tension.
- 2) Under fluoroscopic guidance, slightly advance the delivery system by 2–3 mm to relieve tension.
- 3) Gently rotate the front grip 10–20° to alter the tip alignment.
- 4) If the tapered tip remains caught, partially retract the guidewire within the lumen to reduce contact force, then reattempt withdrawal.
- 5) Do not apply excessive pulling force; forced withdrawal may result in vessel injury or stent

migration.

- 6) If withdrawal remains unsuccessful, fully retract the slider grip to its end position, then stabilize the delivery system and fix the front grip in place. While keeping the front grip stationary, gently move the slider grip toward the front grip (monitoring it and preventing the sheath tip from pressing against the stent and causing stent displacement). This motion advances the sheath forward relative to the tapered tip through the stent graft lumen, allowing the two tips to meet and form a smooth, continuous profile. Once a smooth profile is achieved, withdraw the entire delivery system from the patient under fluoroscopic guidance.

8.2 Handle Disassembly Technique for Partial Deployment

If the stent graft cannot be fully deployed due to partial graft cover severance or slider malfunction:

- 1) Pull back the trigger and fully retract the slider.
- 2) Stabilize the delivery system under DSA.
- 3) Insert a hemostat into the slider grip disassembly ports.
- 4) Manually separate the slider grip.
- 5) Gradually retract the sheath manually until the stent graft is completely released.
- 6) Withdraw the delivery system carefully from the patient.

8.3 Alternative Method for Releasing Proximal Bare Stent

In the unlikely event that the proximal bare stent fails to release due to proximal anchor malfunction:

- 1) Ensure the delivery system remains stable under fluoroscopy.
- 2) Remove the safety buckle and rotate counterclockwise the lock of the proximal releaser using forceps.
- 3) Use hemostats to grip the deployment shaft and gently retract it until the bare stent is released.
- 4) Confirm complete implantation and slowly remove the delivery system.

IX Patient follow-up

Aspirin and clopidogrel, or alternative anticoagulants, are recommended for 3 months after implantation. Continued aspirin is recommended.

9.1 Periods of follow up

Clinical follow-up is recommended at 24 hours, 1, 3-6, 12 and 24 months after the procedure. The follow-

up intervals may need to be adjusted based on the specific conditions of each patient.

9.2 Methods used during follow up

Depending on the patient, hematopathology, X-ray, CTA, MRI, or DSA may be performed.

9.3 Follow up items:

- Size of aneurysms or dissection.
- Embolization
- Changing of aneurysms pulsatility or dissection.
- Migration
- Endoleak
- Stent graft distortion.
- In addition, patients with following specific findings during following up, endovascular intervention or surgical treatment should be considered.
- Aneurysms enlarge 5mm and above.
- Changing of aneurysms or dissection pulsatility with or without aneurysms enlarging or endoleak or dissection tearing.
- Continual endoleak with or without aneurysms enlarging or dissection tearing.
- Stent graft migration results in endoleak.

X Storage and transportation

- Keep away from sunlight, keep dry.
- During transportation, the device must not be compressed, exposed to moisture, or subjected to direct sunlight.
- No exposure to organic solvent.

XI Shelf life

The Thoracoabdominal Artery Stent Graft System is sterilized by ethylene oxide and supplied sterile. The manufacture date and expiration date are printed on the label. Do not use the expired product.

XII Disposal

After use, treat all components of the Thoracoabdominal Artery Stent Graft System, including the delivery system and accessories, as potentially biohazardous medical waste. Dispose of them in accordance with



the healthcare facility's procedures, as well as all applicable local, regional, and national regulatory requirements for the medical waste disposal.

Exercise caution during disposal to prevent injury from sharp or rigid components and to mitigate infection or microbial hazards associated with blood-contacting devices.

XIII Packaging and labeling

- The Thoracoabdominal Artery Stent Graft System includes the G-Branch™ Thoracoabdominal Aortic Stent Graft System, the SilverFlow™ PV Peripheral Vascular Stent Graft System, the G-Branch™ AE Main Body Extension Stent Graft System, the G-Branch™ AAA Bifurcated Stent Graft System and the G-Branch™ IE Iliac Extension Stent Graft System. Each Stent Graft System are supplied sterile and packaged separately.
- Each stent graft is pre-loaded onto its delivery system, which is placed in a PETG tray and sealed with double Tyvek pouches bearing the primary label. The packaged device is then placed in an outer carton box together with the instructions for use (IFU), implant card, and other accompanying documents. An outer label is affixed to the carton box.
- The labels provide key information including device specification, sterilization date, expiration date, serial number, batch number and other relevant details as required for identification and traceability.
- Label Symbol Description:

	Date of manufacture
	Manufacturer
	Use-by date
	Serial number
	Catalogue number
	Batch code
	Consult instructions for use
	Sterilized using ethylene oxide

	Do not use if package is damaged and consult instructions for use
	Do not reuse
	Non-pyrogenic
	Do not resterilize
	Keep dry
	Keep away from sunlight
	MR conditional
	Made in China
	Double sterile barrier system
	Medical device
	Unique device identifier
	Importer
	Temperature limit

Manufacturer:

Lifetech Scientific (Shenzhen) Co., Ltd

8F, LifeTech Scientific Building, No.22, Keji 12th Road South,

High-tech Industrial Park, Yuehai Subdistrict, Nanshan District, Shenzhen 518063

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Lifetech Scientific (Shenzhen) Co., Ltd.

Contents


1	Device description.....	1
2	Intended Use.....	5
3	Clinical Benefit.....	5
4	Contraindication.....	5
5	Pre-implantation Preparation.....	6
6	Implantation Procedure	7
7	Management Suggestions	16
8	Potential Adverse Event.....	17
9	Summary of Clinical Data.....	18
10	Warnings and precautions	19
11	Follow-up	21
12	Shelf life	22
13	Storage and Transportation Requirements.....	22
14	Package and label.....	22

Please read this instruction for use carefully before using this product. Failure to follow the instruction for use may cause serious injury to patients.

1 Device description

The Yuranos™ AB, IE and AE Stent Graft System consist of a pre-loaded stent graft and a delivery system. The stent grafts include Abdominal Aortic Bifurcation Stent Graft (AB), Iliac Artery Extension Stent Graft (IE), and Aortic Extension Stent Graft (AE), which are composed of a metal stent (nitinol alloy), a fabric cover (PET), suture (PET), a steel jacket at the stent connection (316L stainless steel), and markers (tantalum). There is a barb structure on the bare stent of the Abdominal Aortic Bifurcation Stent Graft, and Aortic Extension Stent Graft and on the proximal end of the Iliac Artery Extension Stent Graft. The delivery system is mainly composed of a sheath, an inner core, grips, an immobile pole, and other parts, to realize the delivery, control and release of the stent grafts in human bodies. The delivery system has a proximal hydrophilic coated segment of 510mm/520mm.

1.1 Abdominal Aortic Stent Graft

(a) Abdominal Aortic Bifurcation Stent Graft (AB) (b) Iliac Artery Extension Stent (IE) (c) Aortic Extension Stent (AE)

Figure 1 Yuranos™ AB, IE and AE Stent Graft System

(1: "O"-shaped mark point; 2: "8"-shaped mark point)

The Yuranos™ Abdominal Aortic Stent Graft System adopts a split-type design, including Abdominal Aortic Bifurcation Stent Graft, Iliac Artery Extension Stent Graft, Aorto-Uni-Iliac Stent Graft, and Aortic Extension Stent Graft.

The Abdominal Aortic Bifurcation Stent Graft (Figure 1) is the primary component which is inserted into the patient's abdominal aorta. It is designed with bifurcation, including a proximal bare stent, a main body, a long branch and a short branch.


The Iliac Artery Extension Stent Graft (Figure 1) is designed in a straight cylinder structure. It is used to extend the length of implanted Abdominal Aortic Bifurcation Stent Graft or Aorto-Uni-Iliac Stent Graft as needed based on the patient's anatomy.

The Aortic Extension Stent Graft (Figure 1) includes a proximal bare stent graft and a stent graft's main body. It is used to extend the length of implanted Abdominal Aortic Bifurcation Stent Graft or Aorto-Uni-Iliac Stent Graft as needed based on the patient's anatomy.

Each stent graft is separately loaded in its respective delivery system to be introduced into the human vascular system during a medical procedure.


An Abdominal Aortic Stent Graft is composed of a superelastic nitinol stent and a PET fabric graft, which is fixed onto the metal stent with suture. After the stent is implanted into the human body, the superelastic metal stent can provide sufficient and stable expansion force to deploy the stent graft and anchor it in the corresponding position as a whole. The bare stent with a barb design can provide enough anchoring force to prevent the stent from migration. The fabric graft can serve as an alternative conduit for blood flow within the patient's vasculature by excluding the aneurysmal sac from blood flow and pressure, thereby reducing the risk of aneurysm rupture.

1.2 The delivery system for the Abdominal Aortic Bifurcation Stent Graft

1	Tapered Tip	2	AB or AE	3	Markband
4	Fixation cap	5	Fixation anchor	6	Immobile Pole
7	Sheath	8	Front Grip	9	Sliding Grip
10	Trigger	11	Screw Gear	12	Rear Grip
13	Proximal releaser	14	Luer Connector	15	Hemostatis valve
16	Unlocking buckle				

(a) Yuranos™ Delivery System (AB Stent Graft and AE Stent Graft)

1	Tapered Tip	2	IE	3	Sheath
4	Front Gripgrip	5	Immobile pole grip	6	Immobile pole
7	Inner core	8	Hemostatis valve	9	Luer Connector

(b) Yuranos™ Delivery System (IE Stent Graft)

Figure 2 Yuranos™ Delivery Systems


1.3 Device specifications

Both AB, IE and AE device are pre-loaded onto their respective delivery system and supplied sterile. The AB, IE System and AE System are available in the following lengths and diameters as shown in following table 1, 2 , 3 and 4.

Table 1 Specifications Summary of Abdominal Aortic Bifurcation Stent Graft (AB)

Unit: mm

Specification*	Proximal Dia.(D1)	Distal Dia.		Length of sent graft main body & long branch(L1)	Length of main body (L2)	Length of short branch (L3)	Length of bare stent(L4)	Delivery System
		long branch (D2)	short branch (D3)					
AB-20XX-YY-ZZZ-S	20	10	10	100,110, 120,130, 140,160	50,70	30	16	18F
AB-22XX-YY-ZZZ-S	22	10	10					18F
AB-24XX-YY-ZZZ-S	24	12	12					18F
AB-26XX-YY-ZZZ-S	26	12	12					20F
AB-28XX-YY-ZZZ-S	28	12	12					20F
AB-28XX-YY-ZZZ-S	28	14	14					20F
AB-30XX-YY-ZZZ-S	30	12	12					20F
AB-30XX-YY-ZZZ-S	30	14	14					20F
AB-32XX-YY-ZZZ-S	32	12	12					20F
AB-32XX-YY-ZZZ-S	32	14	14					20F
AB-35XX-YY-ZZZ-S	35	12	12					20F
AB-35XX-YY-ZZZ-S	35	14	14					20F


***AB-WWXX-YY-ZZZ-S** is the specification of Yuranos™ AB Stent Graft System which WW is Proximal Dia. (D1), XX is Proximal Distal Dia. (D2&D3), and YY is Length of main body (L2), ZZZ is the Length of sent graft main body & long branch (L1), S is the code of bare stent graft.

Table 2 Specifications Summary of Iliac Artery Extension Stent Graft (IE)

Unit: mm

Specification*	Proximal diameter (D4)	Distal diameter (D5)	Covered film length(L5)	Delivery System
IE-12XX-YYY	12	10, 12, 14, 16, 18	40,60,70,80,90, 100,	16F
		20, 22, 24, 26., 28	120, 140, 160, 180	18F
IE-14XX-YYY	14	10, 12, 14, 16, 18	40,60,70,80,90, 100,	16F

		20, 22, 24, 26, 28	120, 140, 160, 180	18F
IE-16XX-YYY	16	10, 12, 14, 16, 18	40,60,70,80,90, 100,	16F
		20, 22, 24, 26, 28	120, 140, 160, 180	18F


*IE-ZZXX-YYY is the specification of Yuranos™ IE Stent Graft System which ZZ is Proximal Dia. (D4), XX is Distal Dia. (D5), and YYY is Covered length of stent graft (L5).

Table 3 Specifications Summary of Aortic Extension Stent Graft (AE)

Unit: mm

Specification*	Proximal diameter (D1)	Distal diameter (D2)	Covered length of stent graft (L1)	Length of bare stent graft (L2)	Delivery System
AE-20XX-YY	20	20	40,60	16	18F
AE-22XX-YY	22	22	40,60	16	18F
AE-24XX-YY	24	24	40,60	16	18F
AE-26XX-YY	26	26	40,60	16	20F
AE-28XX-YY	28	28	40,60	16	20F
AE-30XX-YY	30	30	40,60	16	20F
AE-32XX-YY	32	32	40,60	16	20F
AE-35XX-YY	35	35	40,60	16	20F

*AE-ZZXX-YY is the specification of Yuranos™ AE Stent Graft System which ZZ is Proximal Dia. (D1), XX is Distal Dia. (D2), and YY is Covered length of stent graft (L1).

Table 4 Product specifications of Abdominal Aortic Stent Graft

Unit: mm

Name of delivery system for Abdominal Aortic Stent Graft	Specification	The outer diameter of sheath OD	Effective length L0
Delivery system for IE	16F	5.33	510

	18F	6.00	
Delivery system for AB/AE/AUI	18F	6.00	520
	20F	6.67	

Table 5 Stent Graft Materials Composition

Stent graft component configurations	Materials composition
Yuranos™ Abdominal Aortic Bifurcation Stent Graft (AB)	Nitinol (48%-53% w/w) Stainless steel (8%-13% w/w) Fabric graft (PET, 23%-32% w/w) Tantalum Marker (3%-7% w/w) Suture (PP, 5%-8% w/w)
Yuranos™ Iliac Artery Extension Stent Graft (IE)	Nitinol (41%-52% w/w) Stainless steel (7%-11% w/w) Fabric graft (PET, 16%-28% w/w) Tantalum Marker (2 %-11% w/w) Suture (PP, 9%-24% w/w)
Yuranos™ Aortic Extension Stent Graft (AE)	Nitinol (51%-56% w/w) Stainless steel (8%-12% w/w) Film covering (PET, 24%-33% w/w) Tantalum (1%-4% w/w) Suture (PP, 5%-7%w/w)

Physicians should strictly follow the indications and contraindications below, otherwise serious harm may be caused to patients.

2 Intended Use

The Yuranos™ AB, IE and AE Stent Graft System is for the treatment of infrarenal abdominal aortic aneurysm

The target population is as follows:

- (1) Adult patients, male or female;
- (2) The diameter of the proximal abdominal aortic aneurysm neck ranges from 17 mm to 32 mm;
- (3) The length of the proximal abdominal aortic aneurysm neck is 10 mm or more, and the angle of the proximal aneurysm neck is 75 degrees or less;
- (4) The diameter of the iliac artery anchorage area ranges from 8 mm to 22 mm;
- (5) The length of the iliac artery anchorage area is 15 mm or more;
- (6) Patients with adequate iliac or femoral access.

3 Clinical Benefit

Note: After the EUDAMED website is launched, the Summary of Safety and Clinical Performance (SSCP) can be found at <https://ec.europa.eu/tools/eudamed> using the Basic UDI-DI: 69383701004900VX.

4 Contraindication

The contraindications of Yuranos™ AB, IE and AE Stent Graft System are:

- (1) Patients with acute systemic infection;
- (2) Patients with mesenteric blood flow mainly supplied by the inferior mesenteric artery;

- (3) Patients with accessory renal artery original from abdominal artery;
- (4) Patients who have allergic reaction to the device;
- (5) Patients who can not tolerate contrast agents due to renal insufficiency;
- (6) Patients who are allergic to contrast agents;
- (7) Aneurysms neck with thrombus;
- (8) Women who are pregnant or nursing.

5 Pre-implantation Preparation

Before using the Abdominal Aortic Stent Graft, please read the instruction for use carefully. The following operation process is only for guidance, and appropriate changes are inevitable during the procedure. This guidance is only to assist physicians and cannot replace their judgment ability.

5.1 Materials Required

- 0.035" super-stiff guidewire (260 cm);
- Pigtail catheter
- Surgical instruments for conventional femoral artery incision;
- Puncture instruments;
- Angiographic catheters and labeling catheters;
- Various balloon catheters (for blood vessel dilation before stent implantation or stent dilation after stent implantation);
- Contrast agent;
- Fluoroscopy imaging instrument with digital subtraction angiography;
- High-pressure injection syringes;
- Heparin saline solution;
- Any devices that operational physician considers necessary for the operation.

5.2 A multidisciplinary team that has combined procedural experience with:

- Local anesthesia or general anesthesia;
- Femoral artery exposure, arteriotomy, and repair;
- Conventional guidewire and catheter technology;
- Fluoroscopic and angiographic image interpretation;
- Angioplasty;
- Endovascular stent placement;
- Appropriate use of radiographic contrast agent;
- Techniques to minimize radiation exposure;
- Expertise in necessary patient follow-up modalities;
- Femoral artery exposure, arteriotomy, and repair;
- Experience in traditional vascular surgery for any possibility of failure in EVAR.

5.3 Patient preparations

- Measure proximal aortic neck length and vessel inner diameter;
- Measurement of the angle between aneurysm neck and aneurysm, i.e. aneurysm neck angle;
- Evaluation of the quality of aneurysm neck;
- Confirmation of the suitability of vascular morphology of patient for endovascular repair;
- All kinds of necessary routine tests and examinations, and monitoring of important vital signs;
- Vital signs related to anesthesia, anticoagulation, and angiography that need to be monitored based on routine practice;
- Adjustment of the patient's position so that the area from the aortic arch to the iliac artery bifurcation can be fully exposed to the X-ray;
- Adequate iliac/femoral artery access and exposure of both femoral arteries with standard surgical

techniques;

- Systematic anticoagulant preparation (if the patient is allergic to heparin, choose an appropriate alternative).

5.4 Products preparation

- The stent graft size should be selected according to the vessel inner diameter, and the proximal diameter and distal diameter of the selected stent graft should be 110% to 120% of the landing zone vessel inner diameter;
- The diameter of the proximal end of the Iliac Artery Extension Stent Graft D4 should be at least 2 mm larger than the diameter of the short branch of the Abdominal Aortic Bifurcation Stent Graft D3. If multiple extension stent grafts are to be connected, the diameter of the proximal end of the second extension stent graft should be at least 2 mm larger than the diameter of the distal end of the first extension stent graft, and so on;
- The main covered bifurcation stent graft and Iliac Artery Extension Stent Graft of appropriate length should be selected based on the lengths of the infrarenal abdominal aorta and iliac artery. Generally, the distal anchoring area of the Iliac Artery Extension Stent Graft should be at least 15 mm;
- Check and make sure there is no package damage due to transportation. If the package is damaged or there is a potential risk of sterilization failure or damage to the device, please do not use it and return it to Lifetech Scientific (Shenzhen) Co., Ltd.;


- After the stent graft size is determined, unpack the outer packaging box of the corresponding size product, and take out the product together with the dialysis bag. Open the dialysis bag and take out the stent graft system under a sterile environment. Carefully check whether there is any crease on the surface of the sheath of the delivery system and whether the tapered tip of the front end fits well with the sheath port. If the package is damaged or the product is damaged, please do not use it.

Note: If there is a small gap between the tapered tip and sheath port (this phenomenon is caused by the natural shortening effect of sheath material), it can still be used normally after appropriate measures are taken. There are several solutions for the gap with the delivery system for Bifurcation Stent Graft: 1. Remove the unlocking buckle at the end of the delivery system, and rotate the releaser at the end counterclockwise (opposite to the direction of the arrow) to tighten the inner core; 2. If there is a gap between the front grip and the sliding grip, the sliding grip can also be rotated forward; 3. Alternatively, push the tapered tip gently into the sheath to close the gap between the tapered tip and the sheath, then it can be used normally. For Iliac Artery Extension Stent Graft, if there is a little gap between the tapered tip and sheath port, loosen the tapered tip cap at the end of the outer sheath grip, and push the tapered tip into the sheath slightly for normal use.


6 Implantation Procedure

6.1 Preparations before implantation

- Connect an injection syringe to the hemostasis valve and Luer connector respectively, and rinse the sheath inside and guidewire cavity with sterile heparinized saline solution, as shown in Figure 3.

(a) Rinsing method for Abdominal Aortic Bifurcation Stent Graft system

(b) Rinsing method for Iliac Artery Extension Stent Graft System

Figure 3 Rinsing the stent graft system with sterile heparinized saline solution

- Wipe the outer surface of the sheath with sterile heparinized saline solution to activate the hydrophilic coating on the sheath, as shown in Figure 4;


Figure 4 Wipe the stent graft system with sterile heparinized saline solution

- Expose the femoral artery using the standard surgical procedures. Insert the 0.035" super stiff guidewire into the femoral artery, and advance the guidewire under monitoring with DSA, until the front end of the guidewire reaches the abdominal aorta above the renal artery opening, as shown in Figure 5.


Figure 5 The super stiff guidewire reaches the desired position

Note 1: Before the implantation, the operating table should be kept clean and tidy to prevent other devices or articles (such as surgical forceps, gauze, etc.) on the operating table from interfering with the introduction, release, and recovery of the stent system.

Note 2: Once the device enters the human body, all operations should be carried out carefully under the guidance of the DSA instrument. When encountering abnormal resistance, stop the operation immediately, and proceed with the next step only after the abnormality is eliminated.

6.2 Direction for use for Abdominal Aortic Bifurcation Stent Graft system and Iliac Artery Extension Stent Graft system

- 1) Keep the hemostatis valve at the end of the delivery system toward the opposite side of the operator, and insert the delivery system of the Abdominal Aortic Bifurcation Stent Graft along the guidewire (the short branch and hemostatis valve are on the same side of the delivery system when the stent graft is loaded).


Figure 6 The location on the abdominal aortic aneurysm where the delivery system loaded with the stent graft enters

- 2) Determine the directions of the long and short branches of the stent graft. The directions of the long

and short branches of stent graft are adjusted by rotating the delivery system as a whole to make the short branch face the branch on the other side of the patient. The “O”-shaped and “8”-shaped mark points at the proximal end of the fabric cover can be used as references for adjusting the positions of the long and short branches (the “8”-shaped mark point is located on the side of the short branch), so that the “8”-shaped mark point at the proximal end of the stent graft is always on the side of the branch on the opposite side of the patient. The final positions of the long and short branches should be confirmed by the “O”-shaped mark points at the distal end of the long and short branches (from the bifurcation position, there are three consecutive “O”-shaped mark points on the short branch, as shown in Figure 1).

Note: When the tapered tip of the delivery system cannot rotate along, the whole delivery system can be partially withdrawn to locate the sheath of the delivery system in a relatively straight vascular segment, then you can try to rotate again to adjust the orientation of the stent graft.

- 3) Position the proximal end of the stent graft. Adjust the delivery system so that the “O”-shaped and “8”-shaped mark points at the beginning end of the stent graft are located on the lower edge of the lower renal artery opening.
- 4) Slowly release the stent graft. Under monitoring with DSA, hold the front grip firmly with the left hand, and rotate the sliding grip counterclockwise with the right hand to slowly release the bare stent graft and the front segment of the stent graft until at least one coil of the stent graft completely engages the vascular wall, as shown in Figure 7.

Before the stent graft engages the vessel wall, the proximal end of the stent graft can be adjusted slightly again. The adjustment should be carefully carried out under monitoring with DSA to prevent the stent graft or tip from scratching the endangium.

Note 1: Before releasing the stent graft, keep the exposed sheath straight and ensure that the front grip position does not move again.

Note 2: When the stent graft is released, the sheath will move toward the inside of the grip with the sliding grip. At this time, it is necessary to prevent foreign matters (such as gauze, surgical towel, etc.) from wrapping on the surface of the sheath and entering the gap between the sheath and the front grip to prevent release.


Figure 7 Slowly release the stent graft

- 5) Quickly release the stent graft. When a stent graft is to be released quickly, hold the front grip with the left hand to stabilize the delivery system, hold the trigger with the right thumb backward, remove the screw gear lock, and pull the sliding grip backward to release the stent graft quickly. At the same

time, observe the stent graft release process under DSA until the sheath is exposed from the short branch opening, and stop the release. As shown in Figure 8.


Figure 8 Quickly release the stent graft

- 6) Later release the proximal end of the stent graft. Fix the delivery system with the left hand and remove the unlocking buckle at the end of the delivery system with the right hand (the unlocking buckle is a blue translucent safety device, and its removal is as follows: gently press the end of the unlocking buckle with the thumb and index finger and pull it backward to remove it). Under monitoring with DSA, rotate the proximal releaser at the end of the delivery system clockwise, and pay attention to the forward movement of the tapered tip until the proximal bare stent graft is fully deployed, as shown in Figure 9 (rotate in the direction of the arrow printed on the releaser, the tapered tip can be advanced to deploy the proximal bare stent graft).

Note 1: Please make **sure** that the relative position of the delivery system and the patient remains unchanged during the rotation of the releaser, to avoid stent graft migration.

Note 2: Rotate the releaser to the farthest end (toward the sheath) to ensure that the proximal bare stent graft is fully deployed and detached from the anchoring mechanism.

Note 3: When the releaser is **rotated** to the farthest end, if part of the wave crest shows that the proximal bare stent graft is not detached from the anchoring mechanism, you can try to release the knob by reversing it slightly (in the opposite direction of the arrow on the knob), and then turn it forward again (in the direction of the arrow on the knob) to try to release the bare stent graft.


Figure 9 Later release the proximal end of stent graft

- 7) Introduce the contralateral Iliac Artery Extension Stent Graft. A guidewire is introduced from the contralateral femoral artery into the main stent graft through the short branch of the Abdominal Aortic Bifurcation Stent Graft. Next, introduce the Iliac Artery Extension Stent Graft System along the guidewire, and with the help of DSA, ensure it overlaps the short branch of the Abdominal Aortic Bifurcation Stent Graft by 2 to 3 coils (about 3 cm; the proximal “O”-shaped mark point on the Iliac Artery Extension Stent Graft exceeds the distal “O”-shaped mark point on the short branch by 2 to 3 coils and doesn’t exceed the mark point on the bifurcation). As shown in Figure 10.


Figure 10 Position the Iliac Artery Extension Stent Graft

- 8) Release the contralateral Iliac Artery Extension Stent Graft. After confirming the accurate positioning of the Iliac Artery Extension Stent Graft, rotate the end cap of the sheath grip about half a circle counterclockwise to unlock the immobile pole. Hold the immobile pole grip still with one hand, and withdraw the sheath grip with the other hand to pull back the outer sheath. At the same time, observe

the release of the iliac artery stent graft under DSA until the stent graft is completely deployed, and the distal stent graft engages the vessel wall well. As shown in Figure 11. After completing the implantation of the Iliac Artery Extension Stent Graft, withdraw the delivery system of the Iliac Artery Extension Stent Graft with the detailed operation as follows: Fix the sheath grip with one hand, and withdraw the immobile pole grip with the other hand. Observe through DSA until the gap between the tapered tip and the sheath end is closed, then rotate to tighten the end cap of the sheath grip to lock the system, and withdraw the delivery system for the whole extension stent graft out of the body.


Figure 11 Release the Iliac Artery Extension Stent Graft

- 9) Release the long branch extension stent graft of the ipsilateral Abdominal Aortic Bifurcation Stent Graft. After the implantation of the contralateral Iliac Artery Extension Stent Graft, release the long branch of the Abdominal Aortic Bifurcation Stent Graft according to the operation guidance in Steps 4) or 5), to deploy the whole main stent graft and make the distal long branch well engage the vascular wall, as shown in Figure 12.


Figure 12 Release the proximal bare stent graft of the main stent graft (later release)

- 10) Reset the tapered tip of the delivery system for Bifurcation Stent Graft.

Reset the tapered tip of the delivery system for Bifurcation Stent Graft by following the three steps below:

- a. After the proximal stent graft is released, the whole delivery system is advanced by more than 2 cm with the guidance of DSA, so that the tapered tip and fixation anchor are far away from the bare stent graft, and in a straight vascular cavity (to prevent the fixation cap on the tapered tip from scraping the stent graft or blood vessels when the tapered tip is withdrawn), as shown in Figure 13.

Note 1: Please keep the delivery system moving as a whole during the advancement of the delivery system.

Note 2: If any interference between the tapered tip or the later release mechanism and bare stent graft is observed under DSA, the operation should be stopped. The next operation can be continued after the interference is removed by slightly rotating, withdrawing, and advancing the delivery system.

- b. Fix the delivery system and rotate the releaser counterclockwise (the opposite direction of the arrow printed on the knob) to withdraw the tapered tip slowly. Observe under DSA until the gap between the fixation cap and anchor of the tapered tip is completely closed when the resistance to the knob increases, and the knob is basically rotated to the farthest end, as shown in Figure 14.

Note: Please make sure that the tapered tip of the delivery system is completely reset, that is, the fixation cap and anchor are fully engaged, and otherwise accidental injuries may be caused to patients during the withdrawal of the delivery system.

- c. Hold the sliding grip with the right hand, pull the trigger backward to unlock the screw gear, hold the front grip with the left hand, and withdraw it until the tapered tip is reset to the sheath port (at this point, the gap between the front grip and the sliding grip is also basically closed), as shown in Figure 15.

Note: If the tapered tip sticks to the stent graft, try the following solutions:

- Firstly, withdraw the guidewire until the end of the guidewire leaves the point where the tapered tip sticks to the stent graft, and then slowly withdraw the grip to make the tapered tip pass through the inner cavity of the stent graft;
- Alternatively, rotate the delivery system slightly, then withdraw it;
- Alternatively, advance the delivery system as a whole slightly, then withdraw it.

Do not drag or pull the delivery system forcibly, to avoid migration of the stent graft or damage to the sheath.


Figure 13 Advance the delivery system as a whole


Figure 14 Close the gap between the tapered tip and the fixation anchor


Figure 15 Tapered tip reset

- 11) Withdraw the delivery system of the Abdominal Aortic Bifurcation Stent Graft. According to clinical needs, another Iliac Artery Extension Stent Graft can be placed in the ipsilateral iliac artery. The procedure follows Steps 7) and 8). The result after completing the implantation of the stent graft is shown in Figure 16.


Figure 16 Result after stent graft implantation completion

- 12) Perform angiography to assess the wall engagement for the proximal and distal stent graft and endoleak, and to confirm the positions of the stent components relative to the aneurysm and renal artery. Finally, withdraw the guidewire and catheter and close the femoral artery incisions with the standard surgical procedure.

6.3 Direction for use for Aortic Extension Stent Graft System

When Aortic Extension Stent Graft needs to be used, ensure that there is a minimum 1.5cm overlaps between the Aortic Extension Stent Graft and Abdominal Aortic Bifurcation Stent Graft. Completed the Abdominal Aortic Bifurcation Stent Graft deployment steps, except rotate the handle to open the Aortic Extension Stent Graft entirely before releasing the proximal end of the suprarenal stent of the Aortic Extension Stent Graft.

7 Management Suggestions

7.1 Management suggestions for incomplete stent graft release

During the stent graft implantation process, if the rotation function and backward pull function of the delivery system of the main stent graft fail and the stent cannot be released completely, it is suggested to take the following steps to disassemble the grip of the delivery system, eliminate the fault, and restore to a simple backward pull release system:

- a. Keep the outer sheath of the conveyor delivery system still, press the back ends of the two shell parts of the front grip with fingers to enlarge the gap between the shells, insert the surgical scissors into the gap, pry the front grip open, and remove it;
- b. Use the surgical scissors to pry off the blue trigger on the sliding grip;
- c. With the same method as in step b, squeeze the tail of the sliding grip in dislocation to enlarge the gap between the two valve shells, pry the two-valve shells with the surgical scissors, and remove it;
- d. Pull back the blue ring trigger sleeve (the part before the metal spring) to compress the spring. After unlocking the screw gear, remove the screw gear sliding sleeve assembly from the screw gear. At this point, the screw gear is in a half-loosened state. Check whether there is foreign matter in the screw gear which may hinder the sheath movement;
- e. After removing the obstacles, fix the screw gear and pull back the sheath connector (a white translucent part) to withdraw the outer sheath and complete the further release of the stent graft;
- f. There is a small hole on the left and right sides of the hemostatic valve hose connection at the end of the screw gear. Use surgical scissors or hemostatic forceps to press the spring buckle in the small holes at the same time, and pull the screw gear forward to further remove the screw gear, facilitating further inspection and troubleshooting;
- g. After the screw gear is removed, restore the delivery system to a simple backward pull release system. Fix the transparent immobile pole joint (connected with the steel pipe) and withdraw the sheath joint (connected with the blue sheath) to release the stent graft.

7.2 Management suggestions for stent graft deployment failure

During the stent graft implantation process, if the releaser is rotated to the farthest end and the later release mechanism on the tapered tip end cannot be opened in unexpected cases, trying the following measures is recommended:

- a. Make sure that the releaser is rotated to the farthest end (toward the sheath) and kept there for a period of time. If the stent graft still cannot be released, slightly advance or withdraw the delivery system (the operation needs to be carried out carefully under DSA monitoring to prevent stent graft migration), and observe the deployment of the tapered tip end;
- b. Under monitoring with DSA, carefully turn the releaser back to its original position in the opposite direction. Use a compliant or semi-compliant aorta balloon, expanding the balloon in the main

segment of the Abdominal Aortic Bifurcation Stent Graft, anchor the stent graft, then try to release the proximal bare stent graft again according to Step 6) in Section 5.2;

- c. Introduce the capture device through the upper limb approach, try to capture the tapered tip of the delivery system, rotate the releaser of the delivery system to the farthest end (towards the sheath), fix the grip of the delivery system, and gently drag the tapered tip with the capture device to facilitate the release of the later release mechanism.

Note: The above suggestions are only used as emergency measures under abnormal circumstances, and cannot replace the professional judgment of physicians. The above operations should not be attempted during normal stent graft release.

7.3 Management suggestions for poor stent graft wall engagement or type I or III endoleak

After the stent graft is released, the balloons of corresponding compatible sizes can be delivered selectively to the proximal end of the main bifurcated stent graft, the junction of the branch of the bifurcated stent graft and the extension stent graft, and the distal end of the extension stent graft. The balloons can be expanded to the sizes of the blood vessels, to facilitate the stent graft to engage the vessel wall and stabilize the stent graft.

A leak at the proximal or distal end of the stent or at the joint may be reconstructed to the vessel wall using a balloon, but if the leak is severe and cannot be corrected by balloon expansion, it may be managed by adding an extended stent to the already placed stent assembly.

Note: Please refer to the device instruction manual for all operations.

8 Potential Adverse Event

Potential Adverse Events of the Yuranos™ AB, IE and AE Stent Graft System include but are not limited to the following:

- Aneurysm enlargement
- Aneurysm rupture
- Thrombosis or vessel thrombosis
- Arterial injury
- Artery perforation
- Arteriovenous fistula
- Bleeding, hematoma, or coagulopathy
- Heart failure/myocardial infarction
- Turn to surgical operation
- Edema
- Sexual dysfunction
- Infections
- Nerve injuries, local or systemic injuries (e.g. stroke, paraplegia, paraparesis)
- Pulmonary Adverse Events
- Renal function injury or renal failure
- Wound infection and subsequent attendant problem
- Allergic reaction
- Anesthetic complications
- Arrhythmia
- Bowel complications
- Embolism
- Intercostal pain
- Shock

- Tissue necrosis
- Wound complications and subsequent attendant problems
- Death
- Device-related Adverse Events include but are not limited to:
 - ♦ Insertion and removal difficulties
 - ♦ Release difficulties
 - ♦ Failure to deploy
 - ♦ Inaccurate placement
 - ♦ Endoleak
 - ♦ Stent migration
 - ♦ Stent fracture
 - ♦ Rupture of graft material
 - ♦ Erosion
 - ♦ Stent thrombosis or embolization
 - ♦ Stent graft expansion
 - ♦ Stent graft twisting or kinking

Note: For any serious incident, the user or patient should report it to Lifetech at quality@lifetechmed.com or the local sale representative of Lifetech and your local Competent Authority (CA).

9 Summary of Clinical Data

9.1 Summary of Clinical investigation in China

A clinical investigation to evaluate the safety and efficacy of Yuranos™ Stent Graft System in the endovascular treatment for infrarenal abdominal aortic aneurysms was conducted by Lifetech in China in 2018. The study adopted a prospective, multi-center, single-arm objective performance criteria (OPC) design.

The primary safety endpoint is:

- Incidence of major adverse events (MAE) at 30 days after surgery.
- MAE refers to all-cause death, myocardial infarction, renal failure, respiratory failure, paraplegia, stroke and intestinal ischemia.

And the primary efficacy endpoint is:

- Treatment success rate of abdominal aortic aneurysm (AAA) at 12 months post operation.
- Treatment success is defined as immediate technical success, increase of maximum diameter of AAA ≤ 5 mm at 12 months post operation compared with that before operation, migration of stent graft ≤ 10 mm at 12 months post operation compared with that before discharge, no type I or III endoleak requiring further intervention or rupture of AAA at 12 months post operation.

A total of 22 centers participated in the study. The study enrolled 156 subjects with a mean age of 68.91 ± 6.90 years, consisting of 87.82% males ($n = 137$) and 12.18% females ($n = 19$). Both the Full Analysis Set (FAS) and Safety Set (SS) comprised all 156 enrolled subjects. Among these, 12 subjects were missing primary efficacy endpoint data: 10 were lost to follow-up at the 12-month postoperative visit, and 2 completed the 12-month visit but failed to undergo computed tomographic angiography (CTA) as required by the protocol, causing the primary efficacy endpoint unevaluable. Therefore, the Per Protocol Set (PPS) included 144 cases.

In the FAS, one subject did not achieve immediate technical success due to the failure of the delivery system to reach the predetermined position. Consequently, the immediate technical success rates in the FAS and PPS were 99.36% (155/156) and 100% (144/144), respectively.

For primary efficacy endpoint (Treatment success rate of AAA at 12 months post operation), 6 subjects experienced type I/III endoleak requiring further intervention at 12 months post operation. No increase in

the maximum AAA diameter exceeding 5 mm was observed at 12 months post-operation, no stent migration exceeding 10 mm (compared with measurements taken before discharge) occurred, and no ruptured aneurysms were observed at 12 months post-operation. Therefore, 138 cases achieved treatment success at 12 months post-operation, resulting in a treatment success rate of 88.5% (138/156, 95% CI:83.5%~93.5%) in the FAS and 95.8% (138/144, 95% CI:91.2%~98.3%) in the PPS. Comparison of the treatment success rate at 12 months post operation with OPC of 80% showed a statistical significance ($Z=4.7500$, $P=0.0001$).

For primary safety endpoint (Incidence of MAEs within 30 days post-operation), 2 subjects had MAE within 30 days post-operation, specifically all-cause mortality due to coronary heart disease and intestinal obstruction, resulting an incidence rate of 1.28% (2/156, 95% CI: 0~3.05%) in the SS. The incidence of MAEs within 30 days post operation was compared with OPC of 20% and the result showed a statistical significance ($Z=-5.8447$, $P<0.0001$). The incidence of perioperative delivery-related complications and acute limb ischemia were both 0.00% (0/156).

Three subjects died due to coronary heart disease, intestinal obstruction, and bone metastases, resulting an all-cause mortality of 2.0% (3/154) at 6 months post-operation, with two cases lost to follow-up. One subject died from lung cancer between 180- and 360-days post-operation, resulting an all-cause mortality of 2.6% (4/153) at 12 months post-operation, with three cases lost to follow-up. The deaths of these four subjects were not related to the investigation device.

There was no aneurysm-related mortality in 180- and 360-days post-operation. No serious adverse events were reported as being related to the device under evaluation. No device-related complications, such as stent migration, conversion to open surgery, or secondary intervention, were observed at either 6- or 12-month post-operation.

This study demonstrated the efficacy and safety of Yuranos™ Stent Graft System in treating patients with infrarenal abdominal aortic aneurysms. The follow-up period for this study is five years to evaluate the long-term safety and performance of the device. Updated clinical results will be released upon completion of the follow-up period.

10 Warnings and precautions

10.1 General warnings

- Please read all instructions carefully. Failure to follow the instruction for use may cause serious injury to patients.
- This product (including the stent graft and the delivery system) is disposable and cannot be used again after sterilization or disinfection.
- Please check if the sterile packaging is damaged or unintentionally opened before use. If yes, DO NOT use damaged or open packages. And please follow your local or hospital's rules to handle such a situation.
- Do not use if the sterile pouch is compromised in any way.
- Please check the product carefully before use. If any damage is found or physicians think there is any problem affecting the safety of use, please do not use it. Please follow your local or hospital's rules to handle such a situation.
- Regarding safe disposal of the device or its package, please follow your local or hospital's rules.
- This kind of device is only used by physicians with endovascular interventional treatment skills (including those trained on using this device).
- The long-term efficacy and safety of endovascular graft exclusion have not been established. All patients with endovascular graft exclusion are required to receive long-term follow-up to evaluate their health and stent graft status.
- After use, discard the products and packaging according to the policies of the hospital, regulations,

or local government.

- **STERILIE- DO NOT RESTERILIZATION-SINGLE USE ONLY.** Do not reuse, reprocess or resterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death.

10.2 Patient selection, treatment and follow-up requirements

Clinical studies have shown that inappropriate patient selection can reduce the clinical treatment effect, even leading to endovascular treatment failure which requires surgical operation and causing an increased risk of Adverse Events.

- Patient condition should be evaluated before operation, including heart failure/infarction, pulmonary Adverse Events, edema, and renal insufficiency.
- After stent graft implantation, all patients need to be regularly monitored for stent graft status (including stent graft endoleak, migration, shape and location changes, aneurysm enlargement, etc.) at least once a year, by X-ray or CT with or without contrast.
- There is a risk of continuous aneurysm dilation or other diseases at the stent graft placement site after stent graft implantation, which requires secondary interventional operation or open surgery. A persistent aneurysm dilation or leakage may cause aneurysm rupture.
- Patients with infection after device implantation need antibiotic treatment.
- During stent graft implantation or the second intervention, a standby vascular surgery team is needed to promptly perform open surgery on patients in case of intervention failure.
- The size and shape of the access vessels should conform to the general principles of interventional operation and match the size of sheath of the stent graft delivery system. Excessive curving blood vessels will increase the risk during stent graft delivery and release and the delivery system recovery, increase the difficulty of operation, and even lead to operation failure.
- Significant or circumferential calcification or mural thrombus within the treatment length may adversely impact device patency
- The Yuranos™ AB, IE and AE Stent Graft System products are not recommended for patients who weigh more than 150 kg because obesity may reduce the accuracy of DSA.
- Patients with systemic infection have an increased infection risk after intravascular stent implantation.
- Patients with coagulation disorders have an increased risk of endoleak and bleeding after stent implantation.

The safety and efficacy of Abdominal Aortic Stent Graft have not been evaluated in the following patients:

- Patients with traumatic aortic disease
- Patients with coagulation disorders
- Patients with hereditary connective tissue disorders (e.g. Marfan syndrome, Ehlers-Danlos syndromes)
- Abdominal aortic aneurysm patients with thoracic aortic disease
- Patients with acute systemic infection
- Pregnant or lactating women
- Morbidly obese patients
- Patients under 18 years old
- Patients with life expectancy of less than one year

10.3 Precautions during stent graft implantation

- During stent implantation, anticoagulant therapy should be administered according to hospital or physician advice. If patients are allergic to heparin, other appropriate anticoagulants should be considered;

- In order to reduce the risk of device contamination and its resulting infection, the time to stent implantation after removing the sterilization barrier should be shortened as far as possible;
- In the process of inserting the delivery system, it is necessary to keep the guidewire in position;
- Do not bend the sheath of the delivery system, otherwise, the delivery system/or stent graft may be damaged;
- The DSA equipment must be used to observe during the delivery system introduction and delivery, stent graft release, and delivery system withdrawal;
- During stent implantation, a contrast medium should be used, which will increase the risk of postoperative renal failure in patients with renal insufficiency;
- If balloon catheters are used, do not over-expand or expand them outside the graft. Follow all manufacturer's instructions for catheter operation;
- When the delivery system needs to be rotated to adjust the position of the stent graft, it should be rotated as a whole. It can only be rotated when the stent graft is not released completely. Once the stent graft starts to deploy, the delivery system cannot be rotated. When the sheath is in a twisted blood vessel, it may not be able to rotate along with the rotating grip. At this point, it is necessary to consider withdrawing the sheath to a relatively straight segment of the vessel and then try to rotate it;
- Incorrect stent graft release position may lead to endoleak or vascular blockage due to migration;
- Inadequate anchorage of stent graft increases the risk of stent graft displacement, and the incorrect release mode or stent migration may require surgical intervention;
- During the introduction or withdrawal of the guidewire or the delivery system, if abnormal resistance is encountered, please stop the operation and eliminate the cause of failure before proceeding to the next operation;
- The stent graft should not be used to block any vascular branches supplying blood to viscera or limbs unless it is medically proved to be applicable.

10.4 Warnings on MRI safety and applicability

Non-clinical tests have shown that the Abdominal Aortic Stent Graft can be safely scanned on the following MRI equipment:

Item	Magnetic field conditions	
Static magnetic field strength	1.5 Tesla	3.0Tesla
Spatial gradient magnetic field	2000 G/cm	1700 G/cm
Maximum whole-body average specific absorption rate (SAR)	2.0 W/kg with imaging for 15 minutes	

Based on the non-clinical test results, with the 1.5T and 3.0 T MRI equipment, the maximum whole body average specific absorption rate (SAR) of patients is 2.0 W/kg, and the local temperature rise due to stent graft is no more than 2°C after 15 minutes of imaging. Under this condition, the one-sided artifact of the image obtained with spin-echo sequence is not more than 20 mm, and the one-sided artifact of the image obtained with gradient-echo sequence is not more than 20 mm.

11 Follow-up

Follow-up period

It is suggested that the clinical follow-up should be carried out at 24 hours, 1 month, 3 months, 6 months, 12 months and 24 months after implantation. The follow-up period can be adjusted as appropriate according to the patients' conditions.

Follow-up examination means

According to the patient's condition, blood pathology, X-ray fluoroscopy, B-US, CTA, MRI or DSA

examinations can be performed.

Follow-up examination indicators

- Aneurysm size
- Thrombosis
- Aneurysm pulsation
- Stent graft migration
- Endoleak
- Stent graft deformation or distortion
- A second interventional operation or surgical treatment can be considered for patients with the following cases during follow-up
 - Aneurysm enlargement by more than 5 mm
 - Serious type I and/or type III endoleak
 - Persistent endoleak with dilated or non-dilated aneurysm
 - Endoleak due to stent graft migration
 - Stent graft broken or damaged

12 Shelf life

The Yuranos™ AB, IE and AE Stent Graft System is sterilized with ethylene oxide. The manufacture date and expiration date of the product are shown on the label. The product labels indicate the safe use period, and expired products shall not be used. Please use before the Expiration Date marked on the packaging label.

13 Storage and Transportation Requirements

The Yuranos™ AB, IE and AE Stent Graft System products should be stored in a clean, room temperature in a cool place, keep away from sunlight, and keep dry.

During transportation, the products are strictly prohibited from being squeezed, wet or exposed to sunlight, and must be handled with care.

The Yuranos™ Abdominal Aortic Stent Graft and the delivery system should not contact organic solvents (Including alcohol).


14 Package and label

The Yuranos™ AB, IE and AE Stent Graft System products are packed in PETG tray boxes and sealed with double-layer pouches. Then the Abdominal Aortic Bifurcation Stent Graft System, Aortic Extension Stent Graft System and Iliac Artery Extension Stent Graft System products are respectively packed in outer box with an IFU, and customer information feedback sheet.

The labels indicate the product specification, manufacture date (sterilization date), validity period, batch number, product barcode, sterilization mark, etc.

Symbols explanations are shown below table:

	Date of manufacture
	Use-by date
	Serial number
	Batch code

	Sterilized using ethylene oxide
	Do not reuse
	Consult Instructions for Use
	Keep away from the sunlight
	Keep dry
	Non-pyrogenic
	MR conditional
	Do not use if package is damaged
	Do not re-sterilize
	Manufacturer
	Catalogue number
	MR Conditional, for implant only
	Importer
	Country of manufacture
	Double sterile barrier system
	Medical device
	Unique device identifier
	Temperature limit

**Manufacturer**

Lifetech Scientific (Shenzhen) Co., Ltd

8F, LifeTech Scientific Building, No.22, Keji 12th Road South, High-tech Industrial Park, Yuehai Subdistrict, Nanshan District, Shenzhen 518063, P. R. China.

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