iChemVELOCITY CalChek Reagent

REF 800-7703

Rx Only in the U.S.A.

INTENDED USE

The iChemVELOCITY CalChek Reagents are a set of buffer-based solutions intended for the *in vitro* monitoring of specific gravity, color, and clarity. The iChemVELOCITY CalChek Reagents are manufactured for use only with the iChemVELOCITY and DxU Urine Chemistry Systems and are intended for use exclusively by healthcare professionals.

For *in vitro* Diagnostic Use with the iChemVELOCITY and DxU Urine Chemistry Systems.

PRINCIPLE

Specific Gravity: Specific gravity is determined by refractometry in the Urine Chemistry System.

Color and Clarity: Color and Clarity are measured directly from transmitted and scattered light in the sample using proprietary algorithms embedded in the Urine Chemistry System.

WARNING AND PRECAUTIONS

- The iChemVELOCITY CalChek Reagents have been determined to be nonhazardous under the guidelines issued by OSHA in 29 CFR 1910. 1200(d).
- Appropriate precautions in the collection, handling, storage and disposal of specimens should be exercised accordingly.
- Wear appropriate personal protective equipment when handling.
- Avoid solution contact with eyes and skin. Refer to the Safety Data Sheet for first aid measures.
- Consult local or state authorities with regard to recommended method of disposal.

GHS HAZARD CLASSIFICATION

iChemVELOCITY CalChek
Reagent #1,#2.#3

EUH208

May produce an allergic reaction.

reaction mass of:
5-chloro-2-methyl-4-isothiazolin
-3-one [EC# 247-500-7] and
2-methyl-4-isothiazolin-3-one [EC#
220-239-6](3:1) < 0.05%

SDS Safety Data Sheet is available at beckmancoulter.com/techdocs

NON-REACTIVE INGREDIENTS

- Water
- Food coloring
- Salts
- Suspended particles
- Preservatives

STORAGE AND STABILITY

- Store at 18° C to 25° C (65° F to 77° F).
- CalChek reagents are stable for twelve months from the date of manufacture.
- Open expiration: Unused reagents may be used for a one-time repeat if necessary within eight hours of opening. Discard unused reagents after this period.
- Do not freeze.
- Protect reagents against light.
- Do not use reagents after the expiration date printed on the tube labels.

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- iChemVELOCITY or DxU Urine Chemistry Analyzer

- CalChek Racks
- Personal Protective Equipment

PROCEDURE

Mix by inverting CalChek tubes. Remove stopper. Place in appropriate CalChek racks. Follow the directions in the system's Instructions for Use. The system will automatically read each analyte at a specified time. iChemVELOCITY CalChek reagents come with a barcoded, convenient 10-vial package. The results are displayed on screen or printed as soon as they are available.

CALCHEK VALUES			
CalChek #	Type	Calibrated Values	Report Values
1	Specific Gravity	1.002 +/- 0.003	Pass/Fail
2	Specific Gravity	1.030 +/- 0.005	Pass/Fail
3	Specific Gravity	1.060 +/- 0.005	Pass/Fail
4	Color	Colorless	Pass/Fail
5	Color	Straw	Pass/Fail
6	Color	Normal Yellow	Pass/Fail
7	Color	Amber	Pass/Fail
8	Clarity	Slightly Cloudy	Pass/Fail
9	Clarity	Cloudy	Pass/Fail
10	Clarity	Turbid	Pass/Fail

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May be covered by one or more pat. - see www.beckmancoulter.com/patents

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Rx Only in the U.S.A.

INTENDED USE

iQ Calibrator is intended for *in vitro* diagnostic use with iQ 200 and DxU Microscopy Series Analyzers.

PRINCIPLE

iQ Calibrator is a suspension of fixed human red blood cells in a particulate-free solution. It is used to automatically calibrate the iQ200 and DxU Microscopy Series Analyzers. An accurate count is certified and the expected value (accurate to approximately 8.3%) is printed on each label. Barcodes are provided that contain Lot number, Expiration date and Expected value.

WARNING AND PRECAUTIONS

- The iQ Calibrator has been determined to be nonhazardous under the guidelines issued by OSHA in 29 CFR 1910.1200(d).
- Appropriate precautions in the collection, handling, storage and disposal of specimens should be exercised accordingly.
- Wear appropriate personal protective equipment when handling.
- Avoid skin and eye contact. Refer to the Safety Data Sheet for first aid measures.
- Consult local or state authorities with regard to recommended method of disposal.

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NON-REACTIVE INGREDIENTS

iQ Calibrator: Human erythrocytes, water, salts, dyes and preservatives. The concentration of stabilized human erythrocytes is <1% in an aqueous suspension. Non-reactive for HbsAg, HCV 3.0 Ab, HIV 1 / 2, HIV-1 p24 Ag, RPR and ALT <82 IU/mL. The stabilized human erythrocytes are chemically fixed, which thus minimizes any biohazardous risks.

STORAGE AND STABILITY

- Store at 2° C to 8° C (35.6° F to 46.4° F)
- iQ Calibrator is stable for eight months from date of manufacture
- Open expiration: 24 hours from bottle open time
- Do not freeze
- Do not store in direct sunlight
- Each bottle is expected to yield single calibration and should not be used again
- Do not mix the calibrator from different bottles
- Bring to room temperature before use
- Do not use the solution after the expiration date printed on the bottle label

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- Approved tubes listed in the iQ200 or DxU Microscopy Series Analyzer Instructions for Use. Do not use preservative tubes.
- iQ200 or DxU Microscopy Series Analyzer
- Personal Protective Equipment.

LIMITATION OF THE PROCEDURE

No interferences are reported.

PROCEDURE

Shake the bottle as directed in the iQ200 or DxU Microscopy Series Analyzer Instructions for Use before use and follow the directions for calibration. The included barcodes provide the lot number, expiration date, and expected value to the instrument. The passing results are displayed on screen as soon as they are available.

iQ CALIBRATOR NOMINAL ASSAY VALUES AND RANGES

CONTROL	NOMINAL ASSAY VALUE	RANGE
iQ CALIBRATOR	1200 (MEAN; particles/ μ L)	± 100

Important Note: The actual assay values for the iQ Calibrator may change from LOT number to LOT number.

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Rx Only in the U.S.A.

INTENDED USE

IRISpec CA/CB/CC is an assayed QC material for monitoring of urine chemistry analytes and devices as listed in this Instruction for Use.

For *in vitro* diagnostic use only.

WARNING AND PRECAUTIONS

- All biological materials used in this control set are of non-human origin.
- Appropriate precautions in the collection, handling storage and disposal of specimens should be exercised accordingly.
- Wear appropriate personal protective equipment when handling.
- Avoid skin and eye contact. Refer to the Safety Data Sheet for first aid measures.
- Consult local or state authorities with regard to recommended method of disposal.
- Do not pipette.

GHS HAZARD CLASSIFICATION

IRISpec CA Urine Chemistry DANGER
Control

H360

P201

P280

P308+P313

May damage fertility or the unborn child.
Obtain special instructions before use.
Wear protective gloves, protective clothing and eye/face protection.
IF exposed or concerned: Get medical advice/attention.
Sodium Borate Decahydrate 1 - 2%

IRISpec CB Urine Chemistry DANGER
Control

H360

EUH208
P201

P280

P308+P313

May damage fertility or the unborn child.
May produce an allergic reaction.
Obtain special instructions before use.
Wear protective gloves, protective clothing and eye/face protection.
IF exposed or concerned: Get medical advice/attention.
Boric Acid 0.1 - < 0.3%
Subtilisin < 0.1%

IRISpec CC Urine Chemistry DANGER
Control

H360

May damage fertility or the unborn child.

P201	Obtain special instructions before use.
P280	Wear protective gloves, protective clothing and eye/face protection.
P308+P313	IF exposed or concerned: Get medical advice/attention. Boric Acid 0.3 - 1%

Safety Data Sheet is available at beckmancoulter.com/techdocs

REACTIVE INGREDIENTS

- Glucose, 11 g/L
- Bilirubin (bovine), 80 mg/L
- Sodium nitrite, 6.7 mg/L
- Ascorbic acid, 100 mg/dL
- Hemoglobin (bovine) 14.3 mg/L
- 2,5-dimethylindole, 54 mg/L
- Albumin (bovine) 4.5 g/L
- Calcium chloride, 8.0 g/L
- Acetone, 16 mL/L
- Microbial esterase, 120 mg/L

NON-REACTIVE INGREDIENTS

- Stabilizing Agents
- Buffers
- Salts
- Preservatives
- Deionized Water

STORAGE AND STABILITY

- Store at 2°C to 8°C
- Refrigerate immediately upon receipt
- IRISpec CA, CB and CC are stable for six months from the date of manufacture
- Open expiration: 15 days at 2°C to 8°C
- Protect from light
- Do not use solutions after the expiration date printed on the bottle labels

MATERIALS PROVIDED

- 3 x 100 mL IRISpec CA Control
- 3 x 100 mL IRISpec CB Control
- 3 x 100 mL IRISpec CC Control

MATERIALS REQUIRED, BUT NOT PROVIDED

- Urine Chemistry Strips
- Sample Tubes
- Personal Protective Equipment

PROCEDURE

1. Remove CA/CB/CC set from the refrigerator; pour an aliquot from each bottle into a separate tube. Ensure sufficient volume to wet the strip if testing manually or follow instrument manufacturer's instructions. Return tightly capped control bottles to the refrigerator immediately.
2. Allow aliquots to warm to room temperature (20° to 28°C). Use control aliquots within one hour of pouring. Protect aliquots from light.
3. Treat each aliquot as a separate control material. If using visual read strips, do not dip more than one strip into an aliquot.
4. Compare results obtained for each control with the Assay Value – Table of Acceptable Result Ranges (shown below) for the appropriate strip or strip-instrument combination. Ranges given are intended only as guidelines. More stringent ranges may be validated by each laboratory. However, establishing more lenient ranges is not acceptable. Note: a value less than 40mg/dL for ascorbic acid should prompt the operator to troubleshoot the chemistry strip, the control or instrument (if applicable).
5. Discard aliquots after single use. Do not pour back into bottle.

ASSAY VALUES – TABLE OF ACCEPTABLE RESULT RANGES			
iChemVELOCITY and DxU Urine Chemistry Systems with iChemVELOCITY Urine Chemistry Strips			
Analyte	IRISpec CA	IRISpec CB	IRISpec CC
Bilirubin	0.5 – ≥ 8.0 mg/dL 8.5 – ≥ 140 µmol/L (+) – (+++)	NEG.	NEG.
Urobilinogen	NORM.	2 – ≥ 8 mg/dL 34 – ≥ 140 µmol/L (+) – (+++)	NORM.
Ketone	NEG.	10 – > 100 mg/dL 1 – > 10 mmol/L (+) – (+++)	NEG.
Ascorbic Acid	NEG.	NEG.	40 mg/dL 2.28 mmol/L ++
Glucose	200 – ≥ 1000 mg/dL 11 – ≥ 56 mmol/L (++) – (++++)	NA	NEG.
Protein	100 – ≥ 600 mg/dL 1 – ≥ 6 g/L (++) – (++++)	NEG.	NEG.
Blood	0.03 – 0.5 mg/dL .3 – 5 mg/L Trace – (++)	NEG.	NEG.
pH	8 – 9	5 – 6	5 – 6
Nitrite	NEG.	0.1 – 0.2 mg/dL (+) – (++)	NEG.
Leukocytes	NEG.	75 – 500 WBC/µL (+) – (+++)	NEG.
SG*	1.005 – 1.020	1.010 – 1.020	NA

LIMITATION OF THE PROCEDURE

Minimize oxygen exposure of product. Do not shake bottles. Improper storage and handling may reduce analyte concentration (specifically Bilirubin) and may prematurely result in a negative reading. A color change from orange to green for the CA control may indicate such a process has occurred. The CB control may become turbid with time; this does not impact the performance. Each laboratory should establish its own precision parameters.

TECHNICAL SUPPORT

For addition information, or if damaged product is received, call Beckman Coulter Technical Support at 800-854-3633 (USA or Canada) or contact your local Beckman Coulter Representative.

A Glossary of Symbols (PN C27648) is available at beckmancoulter.com/techdocs.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see www.beckmancoulter.com/patents.

REVISION HISTORY

REVISION	DATE	DESCRIPTION
AA	December 2018	Initial Release

REVISION	DATE	DESCRIPTION
AB	August 2021	Updated Proprietary Name
		Updated website address
		Updated Legal Manufacturer Phone Number
		Updated copyright year.

REVISION	DATE	DESCRIPTION
AC	December 2023	Updated PROCEDURE section
		Updated Brazilian Legal Manufacturer Address
		Updated copyright year.

REVISION	DATE	DESCRIPTION
AD	February 2025	Updated Assay Table

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Rx Only in the U.S.A.

INTENDED USE

iQ Control/Focus Set is intended for *in vitro* diagnostic use with iQ 200 and DxU Microscopy Series Analyzers.

PRINCIPLE

iQ Positive Control and iQ Focus are suspensions of fixed human red blood cells in a buffered, isotonic balanced solution. They are used to control and focus the iQ200 and DxU Microscopy Series Analyzers. iQ Positive Control is used as an abnormal, and iQ Negative Control is used as a normal control to verify accurate counting by the instrument. iQ Focus is used to check light level and focus the instrument.

Barcode labels that contain Lot number, Expiration date and Expected value are provided for each product. An accurate concentration count of these cells is certified and the expected counting reference value (accurate to $\pm 10\%$) is printed on each label.

WARNING AND PRECAUTIONS

- The iQ Control/Focus has been determined to be nonhazardous under the guidelines issued by OSHA in 29 CFR 1910. 1200(d).
- Appropriate precautions in the collection, handling, storage and disposal of specimens should be exercised accordingly.
- Wear appropriate personal protective equipment when handling.
- Avoid skin and eye contact. Refer to the Safety Data Sheet for first aid measures.
- Consult local or state authorities with regard to recommended method of disposal.

GHS HAZARD CLASSIFICATION

iQ Positive Control,
iQ Focus

EUH208

May produce an allergic reaction.
reaction mass of:
5-chloro-2-methyl-4-isothiazolin
-3-one [EC# 247-500-7] and
2-methyl-4-isothiazolin-3-one [EC#
220-239-6](3:1) < 0.05%

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NON-REACTIVE INGREDIENTS

iQ Negative Control: Water, salts and preservatives.

iQ Positive Control/iQ Focus: Human erythrocytes, water, salts, dyes and preservatives. The stabilized human erythrocytes are <1% in an aqueous suspension and are non-reactive for HbsAg, HCV 3.0 Ab, HIV 1 / 2, HIV-1 p24 Ag, RPR and ALT <82 IU/mL. Stabilized human erythrocytes are chemically fixed, which minimizes any hazardous characteristics.

STORAGE AND STABILITY

- Store iQ Focus and iQ Positive Control at 2° C to 8° C (35.6° F to 46.4° F)
- Store iQ Negative Control at 2° C to 28° C (35° F to 82.4° F)
- iQ Control/Focus is stable for eight months from the date of manufacture
- Open expiration: 30 days from bottle open date
- Do not freeze
- Do not store in direct sunlight
- Bring to room temperature before use
- Do not save leftover material or mix Control/Focus material from different bottles
- Do not use iQ Control/Focus after the expiration date printed on the bottle labels

MATERIALS PROVIDED

- 2 X 125 mL iQ Focus
- 1 X 125 mL iQ Positive Control
- 1 X 125 mL iQ Negative Control
- 1 X 40 iQ Positive Control Barcode Labels
- 1 X 40 iQ Negative Control Barcode Labels
- 1 X 50 iQ Focus Control Barcode Labels

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- Approved tubes listed in the iQ200 or DxU Microscopy Series Analyzer Instructions for Use. Do not use preservative tubes.
- iQ200 or DxU Microscopy Series Analyzer
- Personal Protective Equipment

PROCEDURE

Follow the directions for control rack, focus, and daily control as described in the iQ200 or DxU Microscopy Series Analyzer Instructions for Use. Shake the Positive Control and Focus prior to use. **DO NOT SHAKE THE NEGATIVE CONTROL.** The included barcodes inform the instrument of the lot number, expiration date, and expected value. If the results are within the allowable range, the date and time of the QC will appear on the “Last Micro QC” field of the Instrument screen. Daily control results and statistics can be printed from the iQ QC Review and QC Statistics sections.

iQ CONTROL/FOCUS SET NOMINAL ASSAY VALUES AND RANGES

CONTROL	NOMINAL ASSAY VALUE	RANGE
POS	1000 (particles/ μ L)	$\pm$ 100
FOCUS	1100 (particles/ μ L)	$\pm$ 100
NEG	0 (particles/ μ L)	+ 20

Important Note: The actual assay values for the iQ Control/Focus Set may change from LOT number to LOT number.

TECHNICAL SUPPORT

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-854-3633 (USA or Canada) or contact your local Beckman Coulter Representative.

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