COULTER 6C Cell Control

REF 628027, A59925

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INTENDED USE

6C Cell Control is a hematology quality control material used to monitor the performance of UniCel DxH Coulter Cellular Analysis Systems listed in the TABLE OF EXPECTED RESULTS, in conjunction with specific COULTER reagents.

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SUMMARY AND PRINCIPLE

6C Cell Control is a reference product prepared from stabilized human blood. By design, 6C Cell Control confirms and monitors instrument accuracy and precision performance by providing quantitative measurements for the following parameters: WBC, RBC, Hgb, HCT, MCV, MCH, MCHC, RDW, RDW-SD, PLT, MPV, NE, LY, MO, EO, BA, NRBC, NE#, LY#, MO#, EO#, BA#, NRBC#.

REAGENTS

6C Cell Control consists of treated, stabilized human erythrocytes in a buffered isotonic medium. 6C Cell Control also contains a stabilized, platelet-sized component, and fixed erythrocytes to simulate leukocytes and nucleated red blood cells.

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Avian < 0.3% w/w

NON-REACTIVE INGREDIENTS

Buffers, Stabilizers, and Preservatives

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STORAGE, STABILITY, AND DISPOSAL

6C Cell Control is shipped in a thermally insulated container designed to keep it cool. Reagent in unopened packs is stable for up to 95 days (until the expiration date stated on the label) when stored according to requirements. It is not recommended to use the reagent after shelf life expiration indicated on the label. When stored at 2 to 8°C, sealed/unopened tubes are stable until the expiration date shown on the TABLE OF EXPECTED RESULTS. Open vial stability is 16 days assuming that the Instructions for Use section of the package insert is performed a maximum of 18 times within 16 days.

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MATERIALS PROVIDED


- 6C Cell Control.
- Table of Expected Results containing Assigned Values and Expected Ranges.

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INSTRUCTIONS FOR USE

1. Remove 6C Cell Control tube(s) from refrigerator and warm at AMBIENT TEMPERATURE for 10 to 15 minutes.
2. Prepare your instrument: Follow the instructions for Setting Up Controls in System HELP or Instructions for Use when a new lot of controls is received.
3. After warming, mix by hand as follows:
 - A. Roll the tube slowly between the palms of the hands eight times in an upright position.
 - B. Invert the tube and slowly roll it between the palms eight times.
 - C. Gently invert the tube eight times.

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◆ DxH Reagents for use on DxH 800 / 600 / 900 / 690T 🌿 DxH ECO Reagents for use on DxH 900 / 690T

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REVISION HISTORY

REVISION	DATE	DESCRIPTION
AA	September 2024	Initial release to include documentation updates for compliance to European In Vitro Diagnostic Regulations (IVDR) (2017/746). No changes were made to the product.
AB	June 2025	Updated wording in EXPECTED RANGES

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INTENDED USE

6C Plus Cell Control is a hematology quality control material used to monitor the performance of UniCel DxH Coulter Cellular Analysis Systems listed in the TABLE OF EXPECTED RESULTS, in conjunction with specific COULTER reagents.

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SUMMARY AND PRINCIPLE

6C Plus Cell Control is a reference product prepared from stabilized human blood. By design, 6C Plus Cell Control confirms and monitors instrument accuracy and precision performance by providing quantitative measurements for the following parameters: WBC, RBC, HGB, HCT, MCV, MCH, MCHC, RDW, RDW-SD, PLT, MPV, NE, LY, MO, EO, BA, NRBC, NE#, LY#, MO#, EO#, BA#, NRBC#, MDW and LDW. MDW and LDW are only reported when the MDW parameter is enabled.

REAGENTS

6C Plus Cell Control consists of treated, stabilized human erythrocytes in a buffered isotonic medium. 6C Plus Cell Control also contains a stabilized, platelet-sized component, and fixed erythrocytes to simulate leukocytes and nucleated red blood cells.

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WARNING AND PRECAUTIONS

POTENTIAL BIOHAZARDOUS MATERIAL

Product contains biologically sourced materials: Human, Avian, Bovine, Reptile and Ungulate.

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STORAGE, STABILITY, AND DISPOSAL

6C Plus Cell Control is shipped in a thermally insulated container designed to keep it cool. Reagent in unopened packs is stable for up to 95 days (until the expiration date stated on the label) when stored according to requirements. It is not recommended to use the reagent after shelf life expiration indicated on the label. When stored at 2 to 8°C, sealed/unopened tubes are stable until the expiration date shown on the TABLE OF EXPECTED RESULTS. Open vial stability is 16 days assuming that the Instructions for Use section of the package insert is performed a maximum of 18 times within 16 days.

NOTE: The MCV, RDW and/or RDW-SD parameters may show trending through the product's shelf life. This is inherent to the product and should not be considered an indicator of product instability. 95% of recovered values for these parameters should remain within the stated EXPECTED RANGES.

Storage of product with cap down (inverted) may require additional mixing for complete resuspension of cellular components.

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EVIDENCE OF DETERIORATION

Inability to obtain expected values in the absence of known instrument problems or gross hemolysis (darkly colored supernatant) is indicative of product deterioration. However, a slight pink color to the supernatant is normal and should not be confused with deterioration of the product.

MATERIALS PROVIDED


- 6C Plus Cell Control
- Table of Expected Results containing Assigned Values and Expected Ranges.

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- Sufficient COULTER reagents to complete analysis. For the reagents used on the system refer to the specific UniCel DxH 900 and 690T Systems IFU.
- Online System HELP or Instructions for Use.

INSTRUCTIONS FOR USE

1. Remove 6C Plus Cell Control tube(s) from refrigerator and warm at AMBIENT TEMPERATURE for 10 to 15 minutes.
2. Prepare your instrument: Follow the instructions for Setting Up Controls in System HELP or Instructions for Use when a new lot of controls is received.
3. After warming, mix by hand as follows:
 - A. Roll the tube slowly between the palms of the hands eight times in an upright position.
 - B. Invert the tube and slowly roll it between the palms eight times.
 - C. Gently invert the tube eight times.

4. Repeat steps 3A to 3C.
NOTE: DO NOT use any type of mechanical mixer.
5. Run the control in Cassette Presentation on the instrument.
6. Return the tube(s) to the refrigerator **within 30 minutes**.
7. Compare instrument values to those given in the TABLE OF EXPECTED RESULTS.

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EXPECTED RANGES include variation between instruments and represent 95% confidence limits for well maintained instrument systems.

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- ≥95% of the recovered values fall within the EXPECTED RANGE of the Beckman Coulter Assigned Value.
- Recovered values do not trend OUTSIDE the EXPECTED RANGE.
- Control results with System messages and review flags, R, can be accepted except for the specific System Event messages that indicate a hardware parameter is out of limit for some item such as voltage, temperature or pressure. These include System Event: D (or HGB, N, PLT, R, RBC, TNC, WBC).

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- Instrument model and serial number.
- Expiration dates, lot numbers, and supporting data for all reagents, for control(s) in question and for other Cell Control lots currently in use.
- Previous Cell Control lot numbers and the data for these previous lots.
- Data history for COULTER LATRON CP-X control analyses.
- Data from current repeatability study (N=10) using a fresh whole blood specimen and performed according to System HELP and Instructions for Use.
- Data from the IQAP survey as well as data from your last instrument calibration, where appropriate.

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ADDITIONAL INFORMATION

Glossary of Symbols is available at beckmancoulter.com/techdocs by searching for part number C07163.

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Revision History

REVISION	DATE	DESCRIPTION
AA	September 2024	Initial release to include documentation updates for compliance to European In Vitro Diagnostic Regulations (IVDR) (2017/746). No changes were made to the product.
AB	May 2025	Updated SUMMARY AND PRINCIPLE section
		Updated EXPECTED RANGES

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INTENDED USE

Body Fluid Control is a hematology quality control material used to monitor the performance and verify the measuring range of the body fluid cycle of UniCel DxH Coulter Cellular Analysis Systems listed in the TABLE OF EXPECTED RESULTS, in conjunction with specific COULTER reagents. Refer to the System HELP or Instructions for Use.

The assigned values and expected ranges on the TABLE OF EXPECTED RESULTS can be used to monitor performance, and verify the measurable range. This product can also be used to establish your own laboratory mean.

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SUMMARY AND PRINCIPLE

The Body Fluid Control is a reference product prepared from stabilized human erythrocytes to confirm and monitor instrument accuracy and precision performance by providing quantitative measurements during counting.

The COULTER Body Fluid Control TNC and RBC concentrations span the instrument measurable range. Results for each concentration are compared to the established expected range, to assess the instrument performance and to verify the measurable range.

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REAGENTS

Body Fluid Control consists of treated, stabilized human erythrocytes in a buffered, isotonic medium. The control also contains fixed erythrocytes to simulate leukocytes.

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STORAGE, STABILITY, AND DISPOSAL

The Body Fluid Control is shipped in a thermally insulated container designed to keep it cool. Reagents in unopened packs are stable for up to 120 days (until the expiration date stated on the label) when stored according to requirements. It is not recommended to use the reagent after shelf life expiration indicated on the label. When stored at 2 to 8°C, sealed/unopened tubes are stable until the expiration date shown on the TABLE OF EXPECTED RESULTS. Opened vials are stable for 16 days assuming that the Instructions for Use section of the package insert is performed a maximum of 18 times within the 16 days. Storage of product with cap down (inverted) position may require additional mixing for complete resuspension of cellular components.

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EVIDENCE OF DETERIORATION

Inability to obtain expected values in the absence of known instrument problems or gross hemolysis (darkly colored supernatant) is indicative of product deterioration. A slight pink color to the supernatant is normal and should not be confused with deterioration of the product.

INSTRUCTIONS FOR USE: DAILY QUALITY CONTROL

Before analyzing Body Fluid Control, review the following to ensure accurate results are obtained:

MATERIALS PROVIDED

- Body Fluid Control.
- Table of EXPECTED RESULTS containing Assigned Values and EXPECTED RANGES.

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
- Sufficient COULTER reagents to complete analysis. For the reagents used on the system refer to the IFU specific for the UniCel DxH Coulter Cellular Analysis Systems in use.
- Online System HELP or Instructions for Use.

INSTRUMENT SET UP

Prepare your instrument: Follow the instructions for Setting Up Controls in System HELP or the Instructions for Use when a new lot of controls is received.

Analyzing Body Fluid Controls:

1. Remove Body Fluid Control tube(s) from refrigerator and warm at AMBIENT TEMPERATURE for 10 to 15 minutes.
2. After warming, mix by hand as follows:
 - A. Roll the tube slowly between the palms of the hands eight times in an upright position.
 - B. Invert the tube and slowly roll it between the palms eight times.
 - C. Gently invert the tube eight times.

3. Repeat steps 2A to 2C.

NOTE: DO NOT use any other type of mechanical mixer.

4. Run the control using the Single Tube Presentation on the instrument. Run Body Fluid Control tubes in consecutive order beginning with Level 1 and ending with Level 3.
5. Return the tube(s) to the refrigerator **within 30 minutes**.
6. Compare instrument values to those given in the TABLE OF EXPECTED RESULTS.

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◆ DxH Reagents for use on DxH 800 / 600 / 900 / 690T ♻ DxH ECO Reagents for use on DxH 900 / 690T

Whether you use the Beckman Coulter Assigned Value or your own laboratory mean, the instrument is considered well maintained and operating correctly (in control) if:

- ≥95% of the recovered values fall within the EXPECTED RANGE of the Beckman Coulter Assigned Value.
- Recovered values do not trend OUTSIDE the EXPECTED RANGE.

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IMPORTANT:

Flags, Codes and Messages

- DO NOT ignore System messages (R).
- Level 3 has a TNC high enough to trigger the TNC Carryover System Message if the TNC of the following sample is less than or equal to 200 cell/μL or the WBC Carryover System Message if the WBC of the following sample is less than or equal to 0.20×10^3 cell/μL. Samples giving carryover system messages should be repeated.

PERFORMANCE LIMITS

Individual laboratories can expect better repeatability than that shown in the EXPECTED RANGE column. Refer to your System HELP and Instructions for Use for precision performance specifications for your UniCel DxH Coulter Cellular Analysis Systems.

The report generated by IQAP/eIQAP will include an assessment of the measuring range.

TROUBLESHOOTING PROCEDURE

The System HELP and Instructions for Use contain troubleshooting procedures for resolving control recovery problems. If you need additional help, call Beckman Coulter Customer Service at 800-526-7694 (USA and Canada) or contact your local Beckman Coulter Representative.

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- Instrument model and serial number.
- Expiration dates, lot numbers, and supporting data for all reagents, for control(s) in question and for other Body Fluid Control and COULTER 6C Control lots currently in use.
- Previous Body Fluid Control lot numbers and the data for these previous lots.
- Data history for COULTER LATRON CP-X Control analyses.
- Data from the last calibration.
- Data from current repeatability study (N=10) using a fresh whole blood specimen and performed according to System HELP and Instructions for Use.

- Data from the IQAP survey as well as data from your last instrument calibration, where appropriate.

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For additional information or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

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NOTICE TO USER

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/746/EU on In vitro Diagnostic Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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INTENDED USE

COULTER LATRON CP-X Control is intended for use on UniCel DxH Coulter Cellular Analysis Systems that use VCSn technology, in conjunction with specific COULTER reagents. It is used to monitor the volume, conductivity and light scatter parameter measurements.

SUMMARY

Coulter Cellular Analysis Systems with VCSn technology require a stable quality control material of uniform size with appropriate light scattering characteristics. The control monitors the stability of the electrical processing and the fluidic flow rate systems used to measure the volume, conductivity and light scatter.

PRINCIPLE

The COULTER LATRON CP-X Control is supplied as a ready-to-use suspension of latex particles. These particles pass through the flow cell and produce characteristic electrical signals that are measured as volume, conductivity and various angles of light scatter and light loss, for quality control of Coulter Cellular Analysis Systems with VCSn technology.

REAGENTS

COULTER LATRON CP-X Control: A suspension of polystyrene particles in a buffered medium containing a surfactant.

WARNING AND PRECAUTIONS

1. Product contains <0.1% Sodium Azide. Sodium Azide preservative may form explosive compounds in metal drain lines. See National Institute for Occupational Safety and Health Bulletin: Explosive Azide Hazards (8/16/76).
2. May cause sensitization by skin contact.
3. Avoid contact with skin.
4. Wear suitable gloves.
5. After contact with skin, wash immediately with plenty of water.

GHS HAZARD CLASSIFICATION		
COULTER LATRON CP-X Control	WARNING 	May cause an allergic skin reaction.
		Wear protective gloves, protective clothing and eye/face protection.
		If skin irritation or rash occurs: Get medical advice/attention.
		Take off contaminated clothing and wash it before use.
		Reaction Mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) <0.05%

STORAGE, STABILITY, AND DISPOSAL

When stored at 2-30°C, unopened tubes of COULTER LATRON CP-X Control are stable until the expiration date indicated on the label. Opened tubes are stable for 30 days when stored at recommended temperatures.

Bring tube to AMBIENT TEMPERATURE before use.

Keep tubes tightly capped when not in use.

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EVIDENCE OF DETERIORATION

Inability to recover expected results may indicate product instability or deterioration due to improper storage, handling or contamination.

Product frozen and thawed or exposed to temperatures outside the recommended temperature range might not recover indicated values.

Tubes of product that become contaminated with debris might exhibit more than one population resulting in an increase in Coefficients of Variation.

MATERIALS PROVIDED

- LATRON CP-X Control - 8 x 4 mL
- Table of EXPECTED RESULTS containing Assigned Values and EXPECTED RANGES

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- Sufficient COULTER reagents to complete analysis
- System HELP or Instructions for Use

INSTRUMENT OPERATION

Refer to your Systems HELP or Instructions for Use for proper instrument operating instructions.

INSTRUCTIONS FOR USE

NOTE: Bring COULTER LATRON CP-X Control to AMBIENT TEMPERATURE before use.

1. Prepare your instrument: Follow the instructions for Setting Up Controls in System HELP or the Instructions for Use when a new lot of controls is received.
2. Perform Daily Checks according to System HELP or the Instructions for Use.
3. Gently mix a room temperature tube of COULTER LATRON CP-X Control by inversion five (5) to eight (8) times. Avoid foaming.
4. Run the control using Cassette Presentation on the instrument.
5. Aspirate the COULTER LATRON CP-X Control.
6. Examine the Mean/Mode Channel and Coefficient of Variation results for the volume, conductivity and various light scatter measurements. Compare the results to the TABLE OF EXPECTED RESULTS.

NOTE: If a result exceeds the assay limits, refer to your System HELP or Instructions for Use for suggested actions. Follow your laboratory's procedure for handling results that exceed limits.

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EXPECTED RANGES include variation between lots and between instruments, and represent 95% confidence limits for well maintained instrument systems.

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Revision AL, 09/2018

Changes were made to:

- INTENDED USE
- Removed European Hazard Classification
- Added new Patent Information Statement
- Added new Dutch, Traditional Chinese, and Ukrainian languages

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COULTER LIN-X Linearity Control

REF 628029

For *In Vitro* Diagnostic Use
Rx Only in the U.S.A.

INTENDED USE

COULTER LIN-X Linearity Controls are intended to assess calibration and verify the reportable range of UniCel DxH Coulter Cellular Analysis Systems listed in the TABLE OF EXPECTED RESULTS, in conjunction with specific COULTER reagents. Refer to the Instructions for Use or System HELP.

SUMMARY AND PRINCIPLE

CLIA requires that laboratories establish and verify the performance specification, accuracy, precision and reportable range, of new instrumentation upon installation, following significant preventive maintenance, unusual trends/shifts in control recovery or when recommended by the instrument manufacturer.

COULTER LIN-X Linearity Controls can be used on calibrated instruments to comply with these guidelines. The product's WBC, RBC, HGB and PLT concentrations span the instrument's reportable range. Results from repeated measurements for each concentration are compared to the established expected range, to assess the instrument's calibration and to verify the reportable range.

REAGENTS

LIN-X Linearity Controls consist of treated, stabilized human erythrocytes, a stabilized platelet-sized component, and fixed erythrocytes to simulate leukocytes.

WARNING AND PRECAUTIONS

POTENTIAL BIOHAZARDOUS MATERIAL

Product contains biologically sourced material: Human, Avian and Ungulate.

Each human donor unit used in preparation of this material was tested by an FDA approved method for the presence of the antibodies to Human Immunodeficiency Virus (HIV-1 and HIV-2) and Hepatitis C Virus (HCV) as well as for hepatitis B virus surface antigen and found to be negative (were not repeatedly reactive).

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GHS HAZARD CLASSIFICATION		
COULTER LIN-X Linearity Control		May produce an allergic reaction.
		Neomycin Sulfate <0.1%
		Penicillin-G, sodium salt <0.1%

STORAGE, STABILITY, AND DISPOSAL

IMPORTANT: Store LIN-X linearity control tubes in an UPRIGHT (cap up) position.

LIN-X Linearity Controls are shipped in thermally insulated containers designed to keep kits cool. Sealed/unopened tubes are stable until the expiration date shown on the TABLE OF EXPECTED RESULTS when stored between 2 and 8°C. For opened vial stability, refer to the TABLE OF EXPECTED RESULTS. Maximum product performance is achieved when tubes are stored in an upright position. Storage of the control product in a cap down (inverted) position might require additional mixing for complete resuspension of cellular components and may affect product performance.

Dispose of waste product, unused product, and contaminated packaging in compliance with federal, state and local regulations.

INDICATIONS OF INSTABILITY OR DETERIORATION

Inability to obtain expected values in the absence of known instrument problems or gross hemolysis (dark colored supernatant) is indicative of product deterioration. A slight pink color to the supernatant is normal and should not be confused with deterioration of the product.

MATERIALS PROVIDED

- LIN-X Linearity Controls. (Refer to the Product Availability section of this insert for quantities and fill volumes.)
- Table of Expected Results.
- Graphs.
- Worksheet.

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- Sufficient COULTER reagents to complete analysis.
- S-CAL Calibrator, if needed.
- 6C Cell Control.
- Instructions for Use or System HELP.

INSTRUCTIONS FOR USE

Before analyzing LIN-X Linearity Controls, review the following information to ensure accurate results are obtained:

IMPORTANT:

Flagging

When analyzing LIN-X Linearity Controls patient flags will occur. These flags should be ignored if a numerical result is obtained. When voteout (- - - -), incomplete computation (.....) or over range (+++++) flags occur, the sample should be repeated. Some levels may generate a non-blood (N) flag. DO NOT ignore System messages (R).

INSTRUMENT SET UP


Prior to running the LIN-X Linearity Controls on UniCel DxH Coulter Cellular Analysis Systems:

1. Ensure there is sufficient reagent to complete the analysis.
2. Perform quality control procedures/calibration according to the Instructions for Use or System HELP. If acceptable, proceed with running LIN-X.
3. Set the Test Mode to CBC.
4. Temporarily disable the NRBC module.

Refer to the Instructions for Use or System HELP for instructions to perform the functions noted above.

Analyzing LIN-X Linearity Controls:

1. Remove LIN-X Linearity Controls from the refrigerator and warm at AMBIENT TEMPERATURE for 15 minutes.
2. After warming, mix BY HAND as follows:
 - A. Roll tube slowly between the palms of the hands eight times in an upright position.
 - B. Invert the tube and slowly roll it between the palms eight times.
 - C. Gently invert the tube eight times.

3. Repeat steps 2A to 2C three times. Visually inspect the tube contents for complete resuspension. If necessary, repeat steps 2A to 2C until tube contents are uniformly suspended.
4. Run LIN-X linearity control tubes in consecutive order beginning with Level 0 and ending with Level 11.
NOTE: Run the controls in Cassette Presentation. You must first mix as directed above.
5. Run each level of linearity control six times. Repeat flagged samples with non-numeric values.
6. Record the results of Level 0 through Level 10 for applicable parameters on the Worksheet.
NOTE: Level 11 is used to clean the system. The values for this level should not be recorded on the worksheet. **ENABLE** the NRBC module.
7. Calculate the mean of runs 2 through 6 for each applicable parameter and record results on Worksheet. Verify the mean of the recovered values are within the Expected Ranges listed in the TABLE OF EXPECTED RESULTS.
NOTE: Run 1 for each level is used as a Prime. The results should not be used to calculate the mean. A minimum of three sample runs must be used to calculate the mean when repetitive results with non-numeric values occur. Levels that do not meet this minimum (N=3) should be excluded.
8. On the Worksheet, record the Low limit and High limit for the ranges provided on the TABLE OF EXPECTED RESULTS for Level 0 through Level 10. Plot mean recovered values and limits for Level 0 through Level 10 on the LIN-X graphs provided.
9. Beckman Coulter will prepare tabular summaries and graphic presentations of your data. Prior to submission, record values for all runs and complete all account and instrument information on Worksheet. Be sure to record the Lot Number and Expiration Date for each control level. Submit your data online at www.beckmancoulter.com/qap/index.jsp.

ASSIGNED VALUES AND EXPECTED RESULTS

Refer to TABLE OF EXPECTED RESULTS.

ASSIGNED VALUES are determined on validated systems using specific COULTER reagents. The values are confirmed by multiple analysis of the control.

EXPECTED RANGES include variation between lots and between instruments, and represent 99% confidence limits for well maintained instrument systems.

Mean values should recover within the ranges. The ranges established by Beckman Coulter should be used as a guideline. Each laboratory must establish its own criteria for acceptable results.

TROUBLESHOOTING PROCEDURE

Refer to your instrument Instructions for Use or System HELP for maintenance and troubleshooting procedures.

For additional information, call Beckman Coulter Customer Service at 800-526-7694 (USA and Canada) or contact your local Beckman Coulter Representative. To provide faster handling of your inquiry, please have the following information available:

- Instrument model and serial number.
- Lot numbers and expiration dates of reagents currently in use.
- Control recovery data for all levels of controls currently in use.
- Data from the most current calibration procedure.

PRODUCT AVAILABILITY

COULTER LIN-X Linearity Controls 12 x 3.3 mL kit contains 1 tube each of:

- | | |
|-----------|------------|
| – Level 0 | – Level 6 |
| – Level 1 | – Level 7 |
| – Level 2 | – Level 8 |
| – Level 3 | – Level 9 |
| – Level 4 | – Level 10 |
| – Level 5 | – Level 11 |

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ADDITIONAL INFORMATION

Glossary of Symbols is available at techdocs.beckmancoulter.com (PN C07163)

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Revision AP, 09/2018

- Added new Patent Information Statement
- Added new Dutch, Traditional Chinese, and Ukrainian languages

Revision AR, 05/2023

- Updated EC REP Address
- Updated copyright date

COULTER LIN-X Linearity Control Worksheet

ACCOUNT INFORMATION			
Name			
Address			
City / State		Zip Code	
Contact			
Telephone			

Level 0	LOT	WBC	RBC	HGB	PLT
Run 1 (Prime)					
Run 2					
Run 3					
Run 4					
Run 5					
Run 6					
Mean (2 thru 6)					
High Limit					
Low Limit					

Level 1	LOT	WBC	RBC	HGB	PLT
Run 1 (Prime)					
Run 2					
Run 3					
Run 4					
Run 5					
Run 6					
Mean (2 thru 6)					
High Limit					
Low Limit					

Level 2	LOT	WBC	RBC	HGB	PLT
Run 1 (Prime)					
Run 2					
Run 3					
Run 4					
Run 5					
Run 6					
Mean (2 thru 6)					
High Limit					
Low Limit					

Level 3	LOT	WBC	RBC	HGB	PLT
Run 1 (Prime)					
Run 2					
Run 3					
Run 4					
Run 5					
Run 6					
Mean (2 thru 6)					
High Limit					
Low Limit					

Level 4	LOT	WBC	RBC	HGB	PLT
Run 1 (Prime)					
Run 2					
Run 3					
Run 4					
Run 5					
Run 6					
Mean (2 thru 6)					
High Limit					
Low Limit					

Level 5	LOT	WBC	RBC	HGB	PLT
Run 1 (Prime)					
Run 2					
Run 3					
Run 4					
Run 5					
Run 6					
Mean (2 thru 6)					
High Limit					
Low Limit					

Instrument Model		
Serial Number		
IQAP Participant Number		
Test Date		OP ID

Level 6	LOT	RBC	HGB
Run 1 (Prime)			
Run 2			
Run 3			
Run 4			
Run 5			
Run 6			
Mean (2 thru 6)			
High Limit			
Low Limit			

Level 7	LOT	WBC	PLT
Run 1 (Prime)			
Run 2			
Run 3			
Run 4			
Run 5			
Run 6			
Mean (2 thru 6)			
High Limit			
Low Limit			

Level 8	LOT	WBC	PLT
Run 1 (Prime)			
Run 2			
Run 3			
Run 4			
Run 5			
Run 6			
Mean (2 thru 6)			
High Limit			
Low Limit			

Level 9	LOT	WBC	PLT
Run 1 (Prime)			
Run 2			
Run 3			
Run 4			
Run 5			
Run 6			
Mean (2 thru 6)			
High Limit			
Low Limit			

Level 10	LOT	PLT
Run 1 (Prime)		
Run 2		
Run 3		
Run 4		
Run 5		
Run 6		
Mean (2 thru 6)		
High Limit		
Low Limit		

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INTENDED USE

Retic-X Cell Control is a hematology quality control material used to monitor the performance of UniCel DxH Coulter Cellular Analysis Systems listed in the TABLE OF EXPECTED RESULTS, in conjunction with specific COULTER reagents.

The assigned values and expected ranges on the TABLE OF EXPECTED RESULTS can be used to monitor instrument performance. This product can also be used to establish your own laboratory mean.

INTENDED USER

This device is intended to be used by a laboratory professional.

SUMMARY AND PRINCIPLE

Retic-X Cell Control is a reference product prepared from stabilized erythrocytes. Retic-X Cell Control confirms and monitors instrument accuracy and precision performance by providing quantitative measurements for reticulocyte identification using VCSn technology.

REAGENTS

Retic-X Cell Control consists of treated, stabilized human erythrocytes in a buffered isotonic medium. Retic-X Cell Control also contains a stabilized, reticulocyte-like component.

REACTIVE INGREDIENTS

Stabilized and/or fixed erythrocytes from:

Human < 50% w/w

Avian < 5% w/w

NON-REACTIVE INGREDIENTS

Buffers, Stabilizers and Preservatives

WARNING AND PRECAUTIONS

POTENTIAL BIOHAZARDOUS MATERIAL

Product contains biologically sourced materials: Human, Avian, and Ungulate.

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Because no test method can offer complete assurance that hepatitis B virus, Human Immunodeficiency Virus (HIV-1 and HIV-2), or other infectious agents are absent, this specimen/reagent should be handled at Biosafety Level 2, as recommended for any potentially infectious human serum or blood specimen.¹

Product contains <0.1% Sodium Azide. Sodium Azide preservative may form explosive compounds in metal drain lines.²

STORAGE, STABILITY, AND DISPOSAL

Retic-X Cell Control is shipped in a thermally insulated container designed to keep it cool. When stored at 2 to 8°C, sealed/unopened tubes are stable until the expiration date shown on the TABLE OF EXPECTED RESULTS. Opened vials are stable for 16 days assuming that the Instructions for Use section of the package insert is performed a maximum of 18 times within the 16 days.

Reagent in unopened packs is stable for up to 120 days (until the expiration date stated on the label) when stored according to requirements. It is not recommended to use the reagent after shelf life expiration indicated on the label.

Storage of the control product in a cap down (inverted) position might require additional mixing for complete resuspension of cellular components. Over mixing causes hemolysis.

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EVIDENCE OF DETERIORATION

System performance should be considered acceptable when 95% of recovered values fall within the Expected Range. Slight hemolysis occurs over time through routine use and does not interfere with control performance. Failure to recover values within the Expected Range might indicate product deterioration in the absence of known instrument problems or gross hemolysis. A pink to light red color of the supernatant is normal and should not be considered deterioration of the product.

PERFORMANCE LIMITS

Refer to your System HELP and Instructions for Use for precision performance specifications for your UniCel DxH Coulter Cellular Analysis Systems.

MATERIALS PROVIDED

- Retic-X Cell Control. (Refer to the Product Availability section of this insert for quantities and fill volumes).
- Table of Expected Results containing Assigned Values and Expected Ranges.

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- Sufficient COULTER reagents to complete analysis. For the reagents used on the system refer to the IFU specific for the UniCel DxH Coulter Cellular Analysis Systems in use.

- Online System HELP or Instructions for Use.

INSTRUCTIONS FOR USE

1. Remove the COULTER Retic-X Cell Control tubes from the refrigerator and warm at AMBIENT TEMPERATURE for 15 minutes.
2. Prepare your instrument: Follow the instructions for Setting Up Controls in System HELP or Instructions for Use when a new lot of controls is received.
3. After warming, mix BY HAND as follows:
 - A. Roll tube slowly between the palms of the hands eight times in an upright position.
 - B. Invert the tube and slowly roll it between the palms eight times.
 - C. Gently invert the tube eight times.

4. Inspect the tube contents for complete resuspension. If necessary, repeat steps 3A and 3C until tube contents are uniformly suspended. **DO NOT OVERMIX. DO NOT** use any type of mechanical mixer.
5. Prepare and analyze COULTER Retic-X Cell Control in the same manner as a patient sample. Refer to the Instructions for Use or System HELP for sample preparation and instrument operation instructions.
6. Return the control tubes to the refrigerator **within 30 minutes**.
7. Compare instrument values to those given in the TABLE OF EXPECTED RESULTS.

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Refer to TABLE OF EXPECTED RESULTS.

The retic mode of the instrument DOES NOT report RBC Values. The RBC counts given here for each Level are not to be used for performance monitoring of the RBC parameter, but are to be entered into the instrument at the prompt so that RET # values are calculated. If your laboratory routinely reports the reticulocyte absolute number, enter the RBC count given for each Level.

NOTE: Beckman Coulter DOES NOT assign values for the MRV and IRF and makes no performance claims for these parameters for the Retic-X Cell Control.

ASSIGNED VALUES are determined on validated systems using specific COULTER reagents. The values are confirmed by multiple analyses of the control.

EXPECTED RANGES include variation between instruments, and represent 95% confidence limits for well maintained instrument systems.

◆ DxH Reagents for use on DxH 800 / 600 / 900 / 690T 🌿 DxH ECO Reagents for use on the DxH 900 / 690T

Whether you use the Beckman Coulter Assigned Value or your own laboratory mean, the instrument is considered well maintained and operating correctly if:

- ≥95% of the recovered values fall within the EXPECTED RANGE of the Beckman Coulter Assigned Value.
- Recovered values do not trend (five or more consecutive values drift in a single direction) OUTSIDE the EXPECTED RANGE.
- No more than three consecutive values fall outside the EXPECTED RANGE.

Before your current cell control lot(s) expire, perform the following on your new lot(s):

- Confirm that recovered values are within the TABLE OF EXPECTED RESULTS.
- OR
- Establish your own laboratory mean.

QUALITY CONTROL

This product is designed to monitor the reticulocyte capabilities of the operating system. Assuring the quality of the complete system, however, requires monitoring all reagents, controls and system indicators. Refer to the inserts for COULTER 6C Cell Control, COULTER LATRON CP-X Control, and Instructions for Use or System HELP. Inability to obtain any assigned value and/or electronic measurement could indicate poor system performance.

TROUBLESHOOTING PROCEDURE

The Instructions for Use or System HELP contain troubleshooting procedures for resolving control recovery problems. If you need additional help, call Beckman Coulter Customer Service at 800-526-7694 (US and Canada) or contact your local Beckman Coulter Representative.

To provide faster handling of your inquiry, please have the following information available when you call:

- Instrument model and serial number.
- Lot numbers, expiration dates and supporting data for all reagents, for control(s) in question and for other cell control lots currently in use.
- Previous cell control lot numbers and the data for these previous lots.
- Data history for COULTER LATRON CP-X Control.
- Data from your last IQAP Survey Report.

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ADDITIONAL INFORMATION

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (US or Canada) or contact your local Beckman Coulter Representative.

Glossary of Symbols is available at beckmancoulter.com/techdocs by searching for PN C07163.

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The Summary of Safety and Performance is available from the EUDAMED database: ec.europa.eu/tools/eudamed.

SDS Safety Data Sheet is available at beckmancoulter.com/techdocs

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REVISION HISTORY

REVISION	DATE	DESCRIPTION
AA	February 2025	Initial release to include documentation updates for compliance to European In Vitro Diagnostic Regulations (IVDR) (2017/746). No changes were made to the product.
AB	June 2025	Updated wording in EXPECTED RANGES

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COULTER S-CAL Calibrator Kit

REF 628026

**For *In Vitro* Diagnostic Use
Rx Only in the U.S.A.**

INTENDED USE

S-CAL Calibrator is intended for the determination of calibration factors for UniCel DxH Coulter Cellular Analysis Systems listed in the TABLE OF EXPECTED RESULTS, in conjunction with specific COULTER reagents. Refer to the Instructions for Use or System HELP.

INTENDED USER

The device is intended to be used by a laboratory professional.

SUMMARY AND PRINCIPLE

Coulter Cellular Analysis Systems, for which this kit was designed, require a calibrator to convert electronic measurements of each sample into accurate quantitative results expressed in clinical terms.

The parameters (WBC, RBC, HGB, MCV, PLT and MPV) are calibrated using S-CAL Calibrator, a stabilized human erythrocyte preparation. S-CAL Calibrator is an acceptable alternative to whole-blood calibration.

The calibration procedure uses replicate measurements of S-CAL Calibrator. The average result for each parameter is divided into the S-CAL Calibrator Assigned Value to give a calibration factor.

HCT, MCH, MCHC, RDW, NRBC, Retic and DIFF parameters do not require calibration.

REAGENTS

S-CAL Calibrator consists of treated, stabilized, human erythrocytes and platelet sized components in a buffered, isotonic medium. Fixed erythrocytes are added to simulate leukocytes.

REACTIVE INGREDIENTS

Stabilized Human Erythrocytes

Proportion: < 50% W/W

WARNING AND PRECAUTIONS

POTENTIAL BIOHAZARDOUS MATERIAL

Product contains biologically sourced materials: Human, Avian, Reptile and Ungulate.

Each human donor unit used in preparation of the material was tested by an FDA approved method for the presence of antibodies to Human Immunodeficiency Virus (HIV-1 and HIV-2) and Hepatitis C Virus (HCV) as well as for Hepatitis B virus surface antigen and found to be negative (were not repeatedly reactive.)

Because no test method can offer complete assurance that hepatitis B virus, Human Immunodeficiency Virus (HIV-1 and HIV-2), or other infectious agents are absent, this specimen/reagent should be handled at Biosafety Level 2, as recommended for any potentially infectious human serum or blood specimen in the Centers for Disease Control/National Institutes of Health manual "Biosafety in Microbiological and Biomedical Laboratories," 1988.

Product contains <0.1% Sodium Azide. Sodium Azide preservative may form explosive compounds in metal drain lines. See National Institute for Occupational Safety and Health Bulletin: Explosive Azide Hazards (8/16/76).

STORAGE, STABILITY, AND DISPOSAL

The S-CAL Calibrator kit is shipped in a thermally insulated container designed to keep it cool. When stored at 2-8°C, sealed/unopened tubes are stable until the expiration date shown on the TABLE OF EXPECTED RESULTS. When the tubes are opened, the product is stable for 1 hour.

Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations.

EVIDENCE OF DETERIORATION

Inability to obtain expected values in the absence of known instrument problems or gross hemolysis (darkly-colored supernatant) is indicative of product deterioration. However, a slight pink color to the supernatant is normal and should not be confused with deterioration of the product.

MATERIALS PROVIDED

- COULTER S-CAL Calibrator.
- Table of Expected Results.

MATERIALS REQUIRED, BUT NOT PROVIDED

- Sufficient COULTER reagents to complete calibration.
- Instructions for Use or System HELP.
- 6C Cell Control Kit.

INSTRUCTIONS FOR USE

NOTE: Keep S-CAL Calibrator in the refrigerator until ready to begin Section II of INSTRUCTIONS FOR USE.

1. Before calibrating the instrument, refer to the appropriate Instructions for Use or System HELP, then go on to Section II.
2. DO NOT PROCEED IF INSTRUMENT PERFORMANCE DOES NOT MEET SPECIFICATIONS. CALL BECKMAN COULTER CUSTOMER SERVICE OR CONTACT YOUR LOCAL BECKMAN COULTER REPRESENTATIVE.
3. Calibration errors can be caused by:
 - Inadequate mixing
 - Excessive handling
 - Exceeding time limits
 - Failure to analyze a sufficient number of samples

IMPORTANT: Perform the following procedures exactly as described. Complete the entire procedure within 1 hour.

4. Remove S-CAL calibrator kit from the refrigerator and warm at AMBIENT TEMPERATURE for 15 minutes.
5. While the calibrator is equilibrating to ambient temperature, prepare the analyzer as follows:
 - A. Select Calibration Setup from CBC Calibration Screen in QA Menu.
 - B. Scan 2D barcode from the ASSAY SHEET to upload the assigned values.
6. After warming the S-CAL Calibrator, mix by hand as follows:

NOTE: Do **NOT** use a mechanical mixer.

- A. Roll the tube slowly between the palms of the hands eight times in an upright position.
 - B. Invert the tube and slowly roll it between the palms eight times.
 - C. Gently invert the tube eight times.
- 
- D. Repeat steps 6A through 6C.
 - E. Inspect the tube contents to determine if all cells have been uniformly distributed.
 - F. Repeat the mixing procedure if tube contents have not been uniformly distributed.
7. Analyze S-CAL Calibrator 10 times in Cassette Presentation Mode, following Instructions for Use or System HELP.

CHECK CALIBRATION CRITERIA

The DxH systems check the following criteria automatically and flag the calibration results with a yellow background if you need to calibrate a parameter or with a red background if there may be a problem.

Check calibration results. Calibration is required if the Difference is greater than DIFFERENCE LOWER LIMIT MAX and the Factor % Diff is greater than the FACTOR % DIFFERENCE LOWER LIMIT MAX.

PARAMETER	CALIBRATE IF DIFFERENCE LOWER LIMIT MAX	CALIBRATE IF FACTOR % DIFF LOWER LIMIT MAX
WBC	>0.1	>1.25%
RBC	>0.03	>0.7%
Hgb	>0.1	>0.78%
MCV	>1.0	>1.18%
PLT	>6.0	>2.70%
MPV	>0.5	>5.0%

VERIFY CALIBRATION

1. Analyze 6C Cell Controls according to package insert instructions.
2. Your results for the cell control levels should be within the expected range as stated on the TABLE OF EXPECTED RESULTS in the 6C Cell Control kit.
3. If your cell control results are not within the expected range, repeat the sample. If the results of the second sample are not within the expected range, follow the Troubleshooting Procedure on the 6C Cell Control package insert and call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

ASSIGNED VALUES

The assigned values for WBC, RBC, HGB, MCV, PLT and MPV parameters are derived from sufficient samples to recover mean values with a high degree of statistical validity. Measurements are made on each of several COULTER instruments.

All instruments are traceable to whole-blood calibration and operated according to the appropriate Instructions for Use. Calibration stability is verified by comparison of instrument results and Reference Method results for fresh whole-blood specimens. System performance is monitored daily using 6C Cell Control at multiple levels.

Whole-blood calibration consists of K₃EDTA anticoagulated whole blood analyzed using the following reference methods:

WBC and RBC: A single-aperture impedance cell counter such as a COULTER Z Series cell counter and the manufacturer's recommended reagents. Macro dilutions are made using Class A glassware. Both WBC and RBC data are corrected for coincidence. The erythrocytes and leukocytes measurands (analytes) in this calibrator are traceable to Clin Lab Haemat, 1994, 16(2):131-138.

HGB: The hemoglobin measurand (analyte) in this calibrator is traceable to CLSI Standard H15-A3.¹ This hemoglobincyanide spectrophotometric procedure employs modified Drabkins (Zijlstra) Reagent and is referenced to NIST certified filters.

MCV: MCV is calculated: PCV/RBC x 10. Packed-cell volume (PCV) is measured by a hematocrit procedure according to CLSI Standard H7-A3.²

PLT: Phase contrast microscopy.

MPV: Latex particles (PN 7531820).

The assigned values were established using representative samples from this lot of calibrator and are specific to the assay methodologies of the applicable Beckman Coulter reagent systems indicated in the calibrator Table of Expected Results (Assay Sheet).

Values assigned by other methodologies may be different. Such differences, if present, may be caused by inter-method bias.

QUALITY CONTROL

After completing the calibration and verification procedures, the instrument is within the accuracy limit stated in your Product Manuals. Good laboratory practices recommend that each series of patient samples be preceded and followed by a quality-control check. Use COULTER cell controls to check the performance of the instrument.

PRODUCT AVAILABILITY

COULTER S-CAL Calibrator Kit

628026 - 1 x 3.3 mL

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

ADDITIONAL INFORMATION

Glossary of Symbols is available at beckmancoulter.com/techdocs (PN C07163)

May be covered by one or more pat. - see www.beckmancoulter.com/patents

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

Safety Data Sheet is available at beckmancoulter.com/techdocs

The Summary of Safety and Performance is available from the EUDAMED database: <https://ec.europa.eu/tools/eudamed>.

NOTICE TO USER

For a patient/user/third party in the European Union and in countries with identical regulatory regime (*Regulation 2017/746/EU on In vitro Diagnostic Medical Devices*); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Revision History

Revision AH, 02/2017

Changes were made to:

- ADDITIONAL INFORMATION
- Updated Brazilian address
- Updated Japanese address
- Added new Russian address

Revision AJ, 01/2018

Changes were made to:

- Added new Serbian language
- Added new Vietnamese language

Revision AK, 09/2018

Changes were made to:

- Added new Patent Information Statement

A59927-IVDR-AA

4 of 5

- Added new Dutch, Traditional Chinese, and Ukrainian languages

Revision IVDR, 03/2021

Changes were made to:

- Updated ADDITIONAL INFORMATION section
- Updated website address
- Added location for Summary of Safety and Performance
- Added Notified Body number to CE marking
- Updated EC REP Address
- Added REACTIVE INGREDIENTS section
- Updated SUMMARY AND PRINCIPLE section
- Updated REAGENTS section
- Updated Revision History section
- Updated copyright year
- Updated TRADEMARKS section
- Added INTENDED USER section
- Added NOTICE TO USER section

REFERENCES

1. Clinical and Laboratory Standards Institute. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood- Third Edition; Approved Standard. CLSI publication H15-A3. Wayne, PA, 2000.
2. Clinical and Laboratory Standards Institute. Procedure for Determining Packed Cell Volume by the Microhematocrit Method- Third Edition; Approved Standard. CLSI publication H7-A3. Wayne, PA, 2000.

Beckman Coulter Ireland Inc.,
Lismeehan, O'Callaghan's Mills
Co. Clare, Ireland +(353) (0) 65 683 1100

Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.
www.beckmancoulter.com

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