CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE.
OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS.
FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

1.0 Device Description

RevoEdge high pressure cutting balloon dilatation catheter with a balloon diameter of 1.5 to 4.0 mm and a balloon length of 10 to 20mm can be quickly exchanged with a standard length 0.014-inch guide wire. The balloon is made of a minimally compliance material with a rated burst pressure of 20 atmospheres. The proximal shaft of the balloon catheter consisted of a female luer connector bonded to a PTFE-coated stainless-steel tube. The proximal shaft of the catheter allows for superior push performance with a smooth transition to the distal assembly consisting of an outer tube and an inner tube with a balloon laser welded to both tubes at the distal tip. Two radiopaque marker bands are positioned at the shoulder of each end of the balloon. The inner tube accepts a standard 0.014-inch PTCA guide wire. The guide wire enters the tip of the balloon catheter and advances coaxially out the rapid exchange port of the distal assembly, allowing for coaxial guidance and rapid catheter exchange with the use of a standard-length guide wire. Three wedge-shaped nickel-titanium wires with an angle of 120° are distributed on the surface of the balloon, which are parallel to the axis of the balloon. Two marked sections, 3-5mm in length located on the proximal shaft, indicate the position of the catheter relative to the tip of the brachial or femoral artery guide catheter. The designed of this dilatation catheter does not incorporate a lumen for distal dye injections or distal pressure measurements.

• Performance Characteristics:

Balloon sizes: 1.5-4.0 x 10-20mm
Catheter working length: 140cm
Nominal Pressure (NP): 10 atm
Rated Burst Pressure (RBP): 20atm

2.0 Clinical benefit

The clinical benefit is to restore the patency of indicated vessel lumen of patients. The indicated vessels include coronary artery.

3.0 How supplied

Contents: One (1) Balloon Dilatation Catheter
Sterile: Sterilized with ethylene oxide gas. Non-pyrogenic
Storage: Keep away from sunlight, keep dry, and store within temperature limit.

4.0 Intended use

The balloon dilatation catheter is intended for dilatation of stenosis in the coronary vasculature.

5.0 Intended user

Intended users are the competent physicians who have the training of PTCA Balloon catheter management.

6.0 Intended patient population

The RevoEdge is intended to be used for adults and geriatrics patients needing PCI Procedure during treatment who do not meet contraindications.

7.0 Indications

• The balloon dilatation catheter is indicated for dilatation of coronary artery stenosis for the purpose of improving myocardial perfusion.

8.0 Contraindications

Contraindications for patients and conditions:

- Patients with unprotected left main coronary artery disease
- Patients with coronary spasm without organic stenosis
- Patients who are not suitable for coronary artery bypass graft surgery (CABG)
- Patients with severe allergic reactions or resistance to surgery-essential agents such as contrast media.
- Patients who are pregnant or suspected to be pregnant
- Patients with stent fracture lesions
- Distal lesion of stent after stent implantation or delivery through the side cell of a previously placed stent
- Other patients diagnosed by doctors as unsuitable for the use of the device in question.

9.0 Warnings

- For single patient, single procedure use only. Do NOT resterilize and/or reuse, as this can potentially result in compromised device performance and increase risk of inappropriate resterilization and cross contamination.
- Do NOT use the catheter if its package has been opened or damaged.
- To reduce the potential for vessel damage in the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.

- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.

- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do NOT advance or retract the catheter unless the balloon is fully deflated under vacuum as this can potentially result in damage to the vessel wall. If resistance is met during manipulation, determine the cause of the resistance before proceeding.

- Balloon pressure should not exceed the rated burst pressure indicated on the package label for each balloon. The rated burst pressure is based on the results of in vitro testing. Use of a pressure monitoring device is recommended to prevent over pressurization. PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potentially injurious or life-threatening complication.
- Use only the recommended balloon inflation medium. To prevent the possibility of an air embolus, never use air or any gaseous medium to inflate balloon.
- Use the catheter prior to the "Use by" date (Expiration Date) specified on the package.

10.0 Precautions

- Prior to angioplasty, the catheter should be examined to verify functionality and ensure that its size and shape are suitable for the procedure for which it is to be used.
- The catheter system should be used only by physicians trained in the performance of percutaneous transluminal coronary angioplasty.
- Once the balloon protective cover is removed, it cannot be installed again. (Otherwise, the product may be damaged and bent).
- Appropriate anticoagulation, antiplatelet and vasodilator therapy should be administered to the patient.
- When using two guide wires, care should be taken when introducing, torquing and removing one or both guide wires to avoid entanglement. It is recommended that one guide wire be completely withdrawn for the patient before removing any additional equipment.
- When feel resistance during the operation, please confirm the reason before the operation.
- Cobalt (CAS No. 7440-48-4) defined as a CMR substance of Category 1B may be present in the device.

11.0 Possible complications and adverse reactions

Complications and adverse reactions that may result from the use of this product include, but are not limited to, the following:

- Death
- Total occlusion of the coronary artery or bypass graft
- Acute coronary artery occlusion
- Coronary artery vessels tear, perforation, rupture or damage
- Dilated vessel stenosis
- Unstable angina pectoris
- Allergic reactions to contrast media
- Hypo/hypertension
- Infection
- Coronary artery spasm
- Arteriovenous fistula
- Thrombosis
- Coronary artery dissection (NHLBI grading)
- Arrhythmias that require medication or device treatment
- Sterile inflammation or granulomas at the access site
- Myocardial infarct
- Cerebral infarct or cerebral apoplexy
- Embolism
- Additional surgical intervention
- Tissue necrosis, transient or long-lasting
- Bleeding/hemorrhage or hematoma
- Cardiac tamponade/pericardial effusion
- Hemodynamic compromise including vasovagal reaction
- Radiation injury
- Renal failure/insufficiency
- Respiratory failure/insufficiency

Notice: any serious incident that has occurred *in relation to the device* should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

INSTRUCTIONS FOR USE

ENGLISH

12.0 Materials to be used in combination with a balloon catheter include:

- Suitable guide wire, refer to label claim
- 20cc syringe for balloon preparation
- Suitable guiding catheter, refer to label claim
- 10cc or smaller syringe for manual dye injections
- Appropriate inflation medium (e.g. 50:50 sterile mixture of a contrast medium and saline)
- Pressure-indicating inflation device
- Hemostasis valve

13.0 Preparation for Use

- Select an appropriate balloon catheter for the target vessel
- Remove the device from the sterile packaging
- Prior to use, examine all devices carefully for defects. Examine the dilatation catheter for bends, kinks, or any other damage. Do NOT use any defective device.
- Remove the protective balloon stylet and balloon protector
- Flush the guidewire lumen with saline
- Balloon Purging, purge air from the catheter using a 20cc syringe filled with 2 to 3ml of the inflation medium with the balloon catheter pointing downward. Attach an inflation device to the balloon inflation port. Ensure that a meniscus of contrast medium is evident in both the catheter luer connector and the inflation device. Apply negative pressure with the inflation device. Do NOT attempt Pre-Inflation technique to purge the balloon lumen.
- In order to activate coating, wet the balloon catheter with a sterile saline or wipe the balloon catheter with saline saturated gauze immediately prior its insertion in introducer sheath.

Caution: All air shall be removed from the balloon and displaced with contrast medium prior to inserting into the body. Otherwise complications may occur.

14.0 Instruction for Use

• Insertion Technique

- Place the guiding catheter, with a hemostasis valve attached, in the orifice of the target coronary artery
- Advance the guide wire through the guiding catheter to reach and cross the target lesion. Advance the distal tip of the balloon catheter over the proximal end of the guide wire. Ensure that the guide wire exits the balloon catheter through the guide wire exit location
- The hemostasis valve should be gradually tightened to control back flow. Excessive valve tightening may affect balloon inflation/deflation time as well as movement of the guide wire.
- Track the balloon catheter over the wire to cross the lesion using the radiopaque marker(s) to locate the balloon across the lesion

• Balloon Inflation

- Inflate the balloon to dilate the lesion using standard PTCA techniques.
- Begin inflation at 4 atm and hold at this pressure for 10 seconds.
- Increase the inflation pressure at a rate of 2 atm every 10 seconds until the indentation in the balloon disappears. Hold at this pressure for 10 seconds.
- After the indentation has disappeared, increase the inflation pressure at a rate of 2 atm every 10 seconds until the target balloon diameter is achieved by fluoroscopic or ultrasonic assessment. Hold at the final (maximum) pressure for 10 seconds.
- After each subsequent inflation, the distal blood flow should be assessed.
- If a significant stenosis persists, successive inflations may be required to resolve the stenosis. Do NOT exceed the rated burst pressure (see labeling).
- Confirm the results with fluoroscopy.

• Removing the Catheter

- Apply negative pressure of the inflation device and confirm that the balloon is fully deflated.
- Withdraw the balloon catheter into the guiding catheter while preserving guide wire position.
- After the deflated balloon dilatation catheter is withdrawn, it should be wiped clean with gauze soaked with sterile normal saline.
- Inspect the balloon catheter integrity.
- If reinserting the same balloon dilatation catheter, flush the guide wire lumen of the balloon dilatation catheter as described in the "Preparation for Use" section. Prior to reinsertion, the balloon dilatation catheter should be wiped clean with gauze soaked with sterile normal saline.

• Disposal

– After use, dispose and discard the product and packaging in accordance with hospital, administrative and/or local government policy.

15.0 Reference

- Physicians should consult recent literature on current medical practice on balloon dilatation, such as published by American College of Cardiology/ American Heart Association.
- The Summary of Safety and Clinical Performance (SSCP) for product is available in the Eudamed: <https://ec.europa.eu/tools/eudamed>. Search for the device using the Basic UDI-DI 695848140068MR.

(Before the Eudamed is fully functional, please contact the Manufacturer at sales@brosmed.com.)

16.0 Disclaimer of Warranty

ALTHOUGH THE CATHETER, HEREAFTER REFERRED TO AS "PRODUCT", HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, BROSMEDE MEDICAL CO., LTD AND ITS AFFILIATES HAVE NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. BROSMEDE MEDICAL CO., LTD AND ITS AFFILIATES, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. BROSMEDE MEDICAL CO., LTD AND ITS AFFILIATES SHALL NOT BE LIABLE TO ANY PERSONAL OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL

DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND BROSMEDE MEDICAL CO., LTD AND ITS AFFILIATES TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusion and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law.

If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court or competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

EU:  EU Authorized Representative
BrosMed Medical BV
Mgr. Buckestraat 8, 6134 AP Sittard, The Netherlands
+31 850 140 901
cs@brosmed.com

Manufacturer:
BrosMed Medical Co., Ltd
15th building, SMEs Venture Park
Songshan Lake Hi-Tech Industrial Development Zone
Dongguan, Guangdong, 523808, China
Office: +86 (769) 2289 2018
Email: sales@brosmed.com
www.brosmed.com

EXPLANATION OF SYMBOLS

Description	Symbol
Catalog Number	
Lot Number	
Single Sterile Barrier System With Protective Packaging Inside	
Sterilized Using Ethylene Oxide	
Use-by date	
Do not re-use	
Caution	
Consult Instructions For Use or consult electronic instruction for use on company website	
Do Not Resterilize	
Recommended Guiding Catheter	
Contents (numeral represents quantity of units inside)	
Do not use if package damaged	
CE mark	
Manufacturer	
Date of manufacture	
Authorized representative in the European Community	
Medical Device	
Unique device identifier	
Keep away from sunlight	
Keep dry	
Temperature limit	
Cobalt May Be Present in the Device	

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1.0 Device Description

The Mars NC Balloon Dilatation Catheter is designed to allow easy exchange of the catheter using a standard length guidewire. Balloon diameters range from 2.0mm to 6.0mm. The balloon material is made of a minimally compliant material with a rated burst pressure of 22 atmospheres for 02.0-4.0mm and 20 atmospheres for 04.5-6.0mm balloon respectively. The minimally compliant balloon material will allow high pressure dilatation while maintaining precise control of the balloon diameter and length. The proximal shaft of the catheter is composed of a female luer connector bonded to a stainless steel tube. The proximal shaft allows superior proximal pushability with a smooth transition to a distal shaft composed of an outer tube and an inner tube with a balloon welded to both tubes at the distal tip. Two radiopaque marker bands are positioned within the balloon shoulders. The inner tube accepts a standard 0.014 inch PTCA guidewire. The guidewire enters the catheter tip and advances coaxially out the distal Rx port, thereby allowing both coaxial guidance and rapid exchange of catheter with a single standard length guidewire. Two marked sections, 3-5mm in length located on the proximal shaft, indicate catheter position relative to the tip of either a brachial or femoral guiding catheter. The design of this dilatation catheter does not incorporate a lumen for distal dye injections or distal pressure measurements.

2.0 Clinical benefit

The clinical benefit is to restore the patency of indicated vessel lumen of patients. The indicated vessels include coronary artery, bypass graft and post stent dilatation.

3.0 How supplied

Contents: One (1) Balloon Dilatation Catheter
Sterile sterilized with ethylene oxide gas. Non-pyrogenic.
Storage Keep away from sunlight, keep dry, and store within temperature limit.

4.0 Intended use

The balloon dilatation catheter is intended for dilatation of stenosis and post-deployed stent in the coronary artery or bypass graft.

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6.0 Intended patient population

Patients need PTCA during treatment.

7.0 Indications

- The balloon dilatation catheter is indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion.
- Balloon dilatation of a stent after implantation (balloon models 2.0 mm - 6.0 mm only) **Note: Bench testing was conducted with the Mars NC Balloon Dilatation Catheter and marketed balloon expandable stents. Consideration should be taken when this device is used with different manufacturers' stents due to difference in stent design.**

8.0 Contraindications

The catheter is contraindicated for use in:

- Unprotected left main coronary artery
- Coronary artery spasm in the absence of significant stenosis

9.0 Warnings

- For single patient, single procedure use only. Do NOT resterilize and/or reuse it, as this can potentially result in compromised device performance and increase risk of inappropriate resterilization and cross contamination.
- Do NOT use the catheter if its package has been opened or damaged.
- Do not reduce the potential for vessel damage in the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do NOT advance or retract the catheter unless the balloon is fully deflated under vacuum as this can potentially result in damage to the vessel wall. If resistance is met during manipulation, determine the cause of the resistance before proceeding.

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10.0 Precautions

- Prior to angioplasty, the catheter should be examined to verify functionality and ensure that its size and shape are suitable for the procedure for which it is to be used.
- The catheter system should be used only by physicians trained in the performance of percutaneous transluminal coronary angioplasty.
- Appropriate anticoagulation, antiplatelet and vasodilator therapy should be administered to the patient.
- When using two guide wires, care should be taken when introducing, torquing and removing one or both guide wires to avoid entanglement. It is recommended that one guide wire be completely withdrawn for the patient before removing any additional equipment.

11.0 Adverse Events

Possible adverse effects include, but are not limited to the following

- Death
- Acute myocardial infarction
- Total occlusion of the coronary artery or bypass graft
- Coronary vessel dissection, perforation, rupture or injury
- Restenosis of the dilated vessel
- Hemorrhage or hematoma
- Unstable angina
- Arrhythmias, including ventricular fibrillation
- Drug reactions, allergic reaction to contrast medium
- Hypo/hypertension
- Infection
- Coronary artery spasm
- Arteriovenous fistula
- Embolism

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13.0 Preparation for Use

- Select an appropriate balloon catheter for the target vessel
- Remove the device from the sterile packaging
- Prior to use, examine all devices carefully for defects. Examine the dilatation catheter for bends, kinks, or any other damage. Do NOT use any defective device.
- Remove the protective balloon stylet and balloon protector
- Flush the guidewire lumen with saline
- Balloon Purging, purge air from the catheter using a 20cc syringe filled with 2 to 3ml of the inflation medium with the balloon catheter pointing downward. Attach an inflation device to the balloon inflation port. Ensure that a meniscus of contrast medium is evident in both the catheter luer connector and the inflation device. Apply negative pressure with the inflation device. Do NOT attempt Pre-Inflation technique to purge the balloon lumen.
- In order to active coating, wet the balloon catheter with a sterile saline or wipe the balloon catheter with saline saturated gauze immediately prior its insertion in guiding catheter.

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14.0 Instruction for Use

- Insertion Technique
 - Place the guiding catheter, with a hemostasis valve attached, in the orifice of the target coronary artery
 - Advance the guide wire through the guiding catheter to reach and cross the target lesion. Advance the distal tip of the balloon catheter over the proximal end of the guide wire. Ensure that the guide wire exits the balloon catheter through the guide wire exit location
 - The hemostasis valve should be gradually tightened to control back flow. Excessive valve tightening may affect balloon inflation/deflation time as well as movement of the guide wire.
 - Track the balloon catheter over the wire to cross the lesion using the radiopaque marker(s) to locate the balloon across the lesion
- Balloon Inflation
 - Inflate the balloon to dilate the lesion using standard PTCA techniques.
 - After each subsequent inflation, the distal blood flow should be assessed.
 - If a significant stenosis persists, successive inflations may be required to resolve the stenosis. Do NOT exceed the rated burst pressure (see labeling).
 - Confirm the results with fluoroscopy.
- Removing the Catheter
 - Apply negative pressure of the inflation device and confirm that the balloon is fully deflated.
 - Withdraw the balloon catheter into the guiding catheter while preserving guide wire position.
 - After the deflated balloon dilatation catheter is withdrawn, it should be wiped clean with gauze soaked with sterile normal saline.
 - Inspect the balloon catheter integrity.
 - If reinserting the same balloon dilatation catheter, flush the guide wire lumen of the balloon dilatation catheter as described in the "Preparation for Use" section. Prior to reinsertion, the balloon dilatation catheter should be wiped clean with gauze soaked with sterile normal saline.
- Disposal
 - After use, dispose and discard the product and packaging in accordance with hospital, administrative and/or local government policy.

15.0 Reference

- Physicians should consult recent literature on current medical practice on balloon dilatation, such as published by American College of Cardiology/ American Heart Association.
- The Summary of Safety and Clinical Performance (SSCP) for product is available in the Eudamed: <https://ec.europa.eu/tools/eudamed> (Before the Eudamed is fully functional, please contact the Manufacturer at sales@brosmed.com.)

16.0 Disclaimer of Warranty

ALTHOUGH THE CATHETER, HEREAFTER REFERRED TO AS "PRODUCT", HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, BROSMEDE MEDICAL CO., LTD AND ITS AFFILIATES HAVE NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. BROSMEDE MEDICAL CO., LTD AND ITS AFFILIATES, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. BROSMEDE MEDICAL CO., LTD AND ITS AFFILIATES SHALL NOT BE LIABLE TO ANY PERSONAL OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND BROSMEDE MEDICAL CO., LTD AND ITS AFFILIATES TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusion and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court or competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

Manufacturer:
BrosMed Medical Co., Ltd
15th building, SMEs Venture Park
Songshan Lake Hi-Tech Industrial Development Zone
Office: +86(769)2289 2018
Dongguan, Guangdong, 523808, China
Email: sales@brosmed.com
www.brosmed.com

EC REP EU Authorized Representative
BrosMed Medical BV,
Mgr. Buckestraat 8, 6134 AP Sittard, The Netherlands
+31 850 140 901
cs@brosmed.com

EXPLANATION OF SYMBOLS

Description	Symbol
Catalog Number	REF
Lot Number	LOT
Balloon Diameter	BALLOON
Balloon Length	+
Single Sterile Barrier System With Protective Packaging Inside	STERILE
Sterilized Using Ethylene Oxide	STERILE EO
Use-by date	EXP
Do Not Reuse	NO REUSE
Caution	CAUTION
Cobalt May Be Present in the Device	CAS No. 7440-48-4
Consult Instructions For Use or consult electronic instruction for use on company website	IFU
Do Not Resterilize	NO RSTERIL
Recommended Guiding Catheter	GC
Contents (numeral represents quantity of units inside)	CONT
Do not use if package damaged	NO DAM
CE mark	CE 2797
Manufacturer	MAN
Date of manufacture	DATE
Authorized representative in the European Community	EC REP
Medical Device	MD
Unique device identifier	UDI
Keep away from sunlight	SUN
Keep dry	DRY
Temperature limit	10 C - 30 C

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1.0 Device Description

The Mars (Rx type) device is a coronary dilatation catheter designed for easy guidewire exchange. Balloon diameters range from $\varnothing$ 1.0mm to $\varnothing$ 4.0mm. The balloon material is made of a semi-compliant material for diameter 1.0mm to 4.0mm with a rated burst pressure of 14 atmospheres. The proximal shaft of the catheter is composed of a female luer connector bonded to a stainless steel tube. The proximal shaft joins with a smooth transition to a distal shaft composed of an outer tube and an inner tube with a balloon welded to both tubes at the distal tip. Two radiopaque marker bands are positioned within the balloon shoulders. Within the balloon segment with the exception of balloon diameters less than 2.0mm which incorporate a centrally positioned single marker band. The inner tube accepts a standard 0.014 inch PTCA guide wire. The guide wire enters the catheter's tip and advances coaxially out the distal Rx port, thereby allowing both coaxial guidance and rapid exchange of catheter with a single standard length guide wire. Two marked sections of 5mm length each located on the proximal shaft indicate catheter position relative to the tip of either a brachial or femoral guiding catheter. The design of this dilatation catheter does not incorporate a lumen for distal dye injections or distal pressure measurements.

2.0 Clinical benefit

The clinical benefit is to restore the patency of indicated vessel lumen of patients.

The indicated vessels include coronary artery or bypass graft.

3.0 How supplied

Contents: One (1) Balloon Dilator Catheter
Sterile sterilized with ethylene oxide gas. Non-pyrogenic.
Storage Keep away from sunlight. Keep dry, and store within temperature limit.

4.0 Intended use

The balloon dilatation catheter is intended for dilatation of stenosis in the coronary artery or bypass graft.

5.0 Intended user

Intended users are the competent physicians who have the training of PTCA and Balloon catheter management.

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7.0 Indications

The balloon dilatation catheter is indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion.

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9.0 Warnings

- For single patient, single procedure use only. Do NOT resterilize and/or reuse it, as this can potentially result in compromised device performance and increase risk of inappropriate resterilization and cross contamination.
- Do NOT use the catheter if its package has been opened or damaged.
- To reduce the potential for vessel damage in the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do NOT advance or retract the catheter unless the balloon is fully deflated under vacuum as this can potentially result in damage to the vessel wall. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Balloon pressure should not exceed the rated burst pressure indicated on the package label for each balloon. The rated burst pressure is based on the results of in vitro testing. Use of a pressure monitoring device is recommended to prevent overpressurization. PTCA should only be performed at hospitals

where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potentially injurious or life-threatening complication.

- Use only the recommended balloon inflation medium. To prevent the possibility of an air embolus, never use air or any gaseous medium to inflate balloon.
- Use the catheter prior to the "Use by" date (Expiration Date) specified on the package.

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13.0 Preparation for Use

- Select an appropriate balloon catheter for the target vessel.
- Remove the device from the sterile packaging.
- Prior to use, examine all devices carefully for defects. Examine the dilatation catheter for bends, kinks, or any other damage. Do NOT use any defective device.
- Remove the protective balloon stylet and balloon protector.
- Flush the guidewire lumen with saline.
- Balloon Purging, purge air from the catheter using a 20cc syringe filled with 2 to 3ml of the inflation medium with the balloon catheter pointing downward. Attach an inflation device to the balloon inflation port. Ensure that a meniscus of contrast medium is evident in both the catheter luer connector and the inflation device. Apply negative pressure with the inflation device. Do NOT attempt Pre-Inflation technique to purge the balloon lumen.
- In order to active coating, wet the balloon catheter with a sterile saline or wipe the balloon catheter with saline saturated gauze immediately prior its insertion in guiding catheter.

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14.0 Instruction for Use

- Insertion Technique
 - Place the guiding catheter, with a hemostasis valve attached, in the orifice of the target coronary artery.

INSTRUCTIONS FOR USE

- Advance the guide wire through the guiding catheter to reach and cross the target lesion. Advance the distal tip of the balloon catheter over the proximal end of the guide wire. Ensure that the guide wire exits the balloon catheter through the guide wire exit location
- The hemostasis valve should be gradually tightened to control back flow. Excessive valve tightening may affect balloon inflation/deflation time as well as movement of the guide wire.
- Track the balloon catheter over the wire to cross the lesion using the radiopaque marker(s) to locate the balloon across the lesion
- Balloon Inflation
 - Inflate the balloon to dilate the lesion using standard PTCA techniques.
 - After each subsequent inflation, the distal blood flow should be assessed.
 - If a significant stenosis persists, successive inflations may be required to resolve the stenosis. Do NOT exceed the rated burst pressure (see labeling).
 - Confirm the results with fluoroscopy.
- Removing the Catheter
 - Apply negative pressure of the inflation device and confirm that the balloon is fully deflated.
 - Withdraw the balloon catheter into the guiding catheter while preserving guide wire position.
 - After the deflated balloon dilatation catheter is withdrawn, it should be wiped clean with gauze soaked with sterile normal saline.
 - Inspect the balloon catheter integrity.
 - If reinserting the same balloon dilatation catheter, flush the guide wire lumen of the balloon dilatation catheter as described in the "Preparation for Use" section. Prior to reinsertion, the balloon dilatation catheter should be wiped clean with gauze soaked with sterile normal saline.
- Disposal
 - After use, dispose and discard the product and packaging in accordance with hospital, administrative and/or local government policy.

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16.0 Disclaimer of Warranty

ALTHOUGH THE CATHETER, HEREAFTER REFERRED TO AS "PRODUCT", HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, BROSSED MEDICAL CO., LTD AND ITS AFFILIATES HAVE NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. BROSSED MEDICAL CO., LTD AND ITS AFFILIATES, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. BROSSED MEDICAL CO., LTD AND ITS AFFILIATES SHALL NOT BE LIABLE TO ANY PERSONAL OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND BROSSED MEDICAL CO., LTD AND ITS AFFILIATES TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

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ENGLISH

Manufacturer:

BrosMed Medical Co., Ltd
15th building, SMEs Venture Park
SongShan Lake Hi-Tech Industrial Development Zone
Dongguan, Guangdong, 523808, China
Office: +86 (769) 2289 2018
Email: sales@brosmed.com
www.brosmed.com

EU Authorized Representative

BrosMed Medical B.V.
Mgr. Buckestraat 8, 6134 AP Sittard, The Netherlands
+31 850 140 901
cs@brosmed.com

EXPLANATION OF SYMBOLS

Description	Symbol
Catalog Number	REF
Lot Number	LOT
Balloon Diameter	BALLOON
Balloon Length	BALLOON
Single Sterile Barrier System With Protective Packaging Inside	STERILE
Sterilized Using Ethylene Oxide	STERILE
Use-by date	
Do Not Reuse	
Caution	
Cobalt May Be Present in the Device	
Consult Instructions For Use or consult electronic instruction for use on company website	
Do Not Resterilize	
Recommended Guiding Catheter	
Contents (numerical represents quantity of units inside)	
Do not use if package damaged	
CE mark	
Manufacturer	
Date of manufacture	
Authorized representative in the European Community	
Medical Device	
Unique device identifier	
Keep away from sunlight	
Keep dry	
Temperature limit	

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