Medtronic MiniMed
13000 Devonshire Street,
Northridge, CA 91325, USA
Tel. +1 (866) 942-1223
www.medtronic.com

ATTESTATION STATEMENT

As a responsible member of Medtronic MiniMed, I hereby attest that I believe the Annex to IFU, issued by the Medtronic MiniMed, to be true and accurate copy.

Masae Tanemura
Senior Regulatory Specialist

Signature

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document

State of California


County of LOS ANGELES

On NOVEMBER 7, 17 before me, NAVIN K. GAUR
(date) (Insert name and title of the officer)


Personally appeared MASAE TANEMURA, who provide to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the instrument. I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

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Annex to IFUs № 6025871-021, Mp6025560-21A2_a, Mp6025713-21A1_B, Mp6026011-12B1_a
manufactured by Medtronic MiniMed, USA

Name of medical device	Guardian 2 Link Transmitter Kit (model MMT-7775)
Information about manufacturer of medical device (name, location address, mailing address, phone number, address of place of production)	Developer: 1. Medtronic MiniMed 18000 Devonshire Street, Northridge, CA 91325, USA Phone: +1-800-646-4633 www.medtronicdiabetes.com Manufacturer: Medtronic MiniMed 18000 Devonshire Street, Northridge, CA 91325, USA Phone: +1-800-646-4633 www.medtronicdiabetes.com Authorized Representative: ООО «Медтроник», 123317, Россия, г. Москва, Пресненская наб., д. 10, / «Medtronic» LLC, 123317, Moscow, Presnenskaya Naberezhnaya., 10 Phone.: (495) 580-73-77 Fax: (495) 580-73-78 www.medtronic.ru Addresses of manufacturing facilities of medical device: 1. Medtronic MiniMed 18000 Devonshire Street Northridge, CA 91325 USA Phone: +1-800-646-4633 www.medtronicdiabetes.com

2. IntriCor Corporation
1260 Red Fox Road
Arden Hills, MN 55112 USA
Phone.: 651-6049553

Classification of medical device

Class depending on the potential risk of a medical device, in accordance with the nomenclature classification

Health Authority/ Country	Class risk	Citation/ Device regulated as per	Document
European Commission, EU	Ila	e.g. 93/42/EEC as amended	Medical Device Directive
FDA, USA	N/A	N/A	N/A
Ministry of Health, Russia	2a	Russian MOH Order № 4n dated 06.06.2012	Russian MOH Order № 4n dated 06.06.2012

GMDN

44611

Class of electrical safety, laser safety, and other types of hazards (if any)

Transmitter

Type of protection against electric shock: Internally powered equipment
Degrees of Protection Provided by Enclosures per IEC 60529: IP48 - Protected against the effects of continuous immersion in water (2.4 meters (8 feet) immersion for 30 minutes), when connected to the Glucose Sensor
Degree of protection against electric shock: Type BF applied part
Software safety class per IEC 62304: Class C - Serious injury or death is possible

Indications

Transmitter

The transmitter is indicated for single-patient use as a component of select Medtronic continuous glucose sensing systems and MiniMed sensor-enabled pump systems.

Serter

The serter is used as an aid for inserting the Enlite sensor. It is indicated for single-patient use.

Tester

	The Watertight Tester is an accessory item included with select Medtronic Diabetes continuous glucose transmitters. The tester connects to a transmitter and is used to test the transmitter. The tester is also used to form a waterproof seal to protect the connector on the transmitter while the transmitter is being cleaned. Charger The charger is an accessory item included with select Medtronic continuous glucose recorder and transmitter devices. It is used to recharge the internal battery for these devices, as needed.
Contraindications	Transmitter Do not expose your transmitter to MRI equipment, diathermy devices, or other devices that generate strong magnetic fields. If your transmitter is inadvertently exposed to a strong magnetic field, discontinue use and contact the 24 Hour HelpLine or your local representative for further assistance. Serter Do not use the serter on products other than the Enlite sensor. Medtronic cannot guarantee the safety or efficacy of this product if used with other products. Tester, Charger None known.
Warnings and Precautions	Warnings This product contains small parts and may pose a choking hazard for young children. If the tester comes in contact with blood, the tester must be discarded. Dispose of the tester according to the local regulations for medical waste disposal. Bleeding may occur after inserting the sensor. Make sure that the site is not bleeding before connecting the transmitter to the sensor. If bleeding occurs, apply steady pressure with a sterile gauze or clean cloth at the insertion site until bleeding stops. After bleeding stops, connect the transmitter to the sensor. Contact the 24 Hour HelpLine or your local representative if you experience any adverse reactions associated with the transmitter or sensor. Magnetic fields Do not expose your transmitter to MRI equipment, diathermy devices, or other devices that generate strong magnetic fields. If your transmitter is inadvertently exposed to a strong magnetic field, discontinue use and contact the 24 Hour HelpLine or your local representative for further assistance. X-rays, MRIs, diathermy devices, and CT scans If you are going to have an x-ray, diathermy treatment, CT scan, MRI or other type of exposure to radiation, remove your sensor and transmitter before entering a room containing any of these equipment. Important information about airport security systems, and using your transmitter on an airplane, can be found on the

Emergency Card. Be sure to carry the Emergency Card provided with your device when you are traveling.

Precautions

Always use the tester when cleaning the transmitter. Do not use any other test plug with the transmitter.

Do not twist the tester or sensor while attached to the transmitter. This will damage the transmitter.

Do not allow water, or any other liquid, to come in contact with the tester when it is not connected to the transmitter.

A wet tester can cause damage to the transmitter.

Do not allow the transmitter to come in contact with any liquid when not connected to a sensor or to the tester.

Do not clean the o-rings on the tester, as this can damage the o-rings.

Notice

Caution: Any changes or modifications to the devices not expressly approved by Medtronic Diabetes could interfere with your ability to operate the equipment, cause injury, and void your warranty.

Warnings for Charger

This product contains small parts and may pose a choking hazard for young children.

Warnings for Tester

This product contains small parts and may pose a choking hazard for young children.

If the tester comes in contact with blood, the tester must be discarded. Dispose of the tester according to the local regulations for medical waste disposal.

General warnings for Serter

Read this entire user guide before attempting to insert the sensor. The serter does not work the same as other Medtronic insertion devices. Failure to follow directions may result in improper insertion, pain, or injury.

Contact the 24 Hour HelpLine or your local representative, if you have questions or concerns.

Do not make therapy decisions based on sensor glucose values because sensor glucose and blood glucose values may differ. If your sensor glucose reading is low or high, or if you feel symptoms of low or high glucose, do the following prior to making therapy decisions. Confirm your blood glucose with your BG meter using a fingerstick blood sample.

A retractable needle is attached to the sensor and minimal blood splatter may occur. If you are a healthcare professional or caregiver, wrap sterile gauze around the sensor to minimize contact with blood.

Keep as much distance as possible between you and the patient when removing the needle.

Keep the needle housing within sight at all times to avoid an accidental needlestick or puncture.

Always inspect the sensor packaging for damage before use. Sensors are sterile and non-pyrogenic, unless the package has been opened or damaged.

Do not use the sensor if the sterile package has been opened or damaged. Use of an unsterile sensor can cause site infection.

Watch for bleeding at the insertion site (under, around, or on top of the sensor).

If bleeding occurs, do the following:

1. Apply steady pressure using sterile gauze or a clean cloth, placed on top of the sensor for up to three minutes. The use of unsterile gauze can cause site infection.
2. If bleeding stops, connect the transmitter to the sensor.

If bleeding does not stop, do not connect the transmitter (or recorder) to the sensor because blood can get into the transmitter connector, and could damage the device.

If bleeding continues, do the following:

1. Remove the sensor and continue to apply steady pressure until the bleeding stops. Discard the sensor in a sharps container.
2. Check the site for redness, bleeding, irritation, pain, tenderness, or inflammation. Treat based on instructions from your healthcare professional.
3. Insert a new sensor in a different location.

Make sure the sensor is securely placed in theserter to avoid improper insertion, pain, or minor injury.

Keep theserter away from children. This product contains small parts and may pose a choking hazard.

General precautions for Sertser

Wash your hands with soap and water before inserting the sensor to help prevent site infection.

If you are a healthcare professional or caregiver, wear gloves when inserting the sensor into someone other than yourself to avoid contact with patient blood.

Do not insert the sensor through tape. Inserting the sensor through tape may cause improper sensor insertion and function.

Only use alcohol to prepare the insertion site, to ensure that residue is not left on the skin.


Rotate the sensor insertion sites so that they do not become overused.

Discard used sensors and needle housings in a sharps container after each use to avoid accidental needlestick or puncture.


Do not clean, resterilize, or try to extract the needle from the needle housing. An accidental needlestick or puncture may occur.

Do not reuse sensors. Reuse of a sensor may cause damage to the sensor surface and lead to inaccurate glucose values, site irritation, or infection.

Always make sure the sensor is firmly attached so that it does not come out of your body. Factors that may increase the likelihood of the sensor falling off include increased physical activity (particularly for younger patients) and improper taping technique.

Adverse effects	Adverse effects are unknown
Patient contact	Type of contact with the body for all device: Mode of use: Non-invasive (Contact intact skin only) Duration of use: Long term, re-usable device
Deployment method	Charger The transmitter contains a non-replaceable, rechargeable battery that you can recharge as needed with the charger. The charger has a green light that shows the charging status and a red light that communicates any problems during charging. If you see a red light, see the Troubleshooting section. The charger needs one AAA alkaline battery. <i>Note: If the battery is installed incorrectly or is low, the charger will not work. Repeat the battery installation steps using a new battery.</i> Installing a battery in the charger To install a battery in the charger: 1. Push the battery cover in and slide it off (as shown in the following illustration). 2. Insert a new alkaline AAA battery. Make sure the + and - symbols on the battery align with these same symbols shown on the charger. 3. Slide the cover back on the charger until it clicks into place.  Charging the transmitter Caution: Charge the transmitter after each sensor use. A fully charged transmitter works at least six days without recharging. A depleted transmitter can take up to one hour to recharge. To charge the transmitter: 1. Connect the transmitter to the charger by lining it up, flat side down, with the charger. Push the two components

- together fully.
2. Within 10 seconds after the transmitter is connected, a green light on the charger will flash for one to two seconds as the charger powers on. For the rest of the charging time, the charger's green light will continue to flash in a pattern of four flashes with a pause between the four flashes.
 3. When charging is complete, the green light on the charger will stay on, without flashing, for 15 to 20 seconds and then turn off.
 4. After the green charger light turns off, disconnect the transmitter from the charger. The green light on the transmitter will flash for about five seconds and then turn off.

Inserting the sensor

Where to insert the sensor

CAUTION: Make sure the sensor site is at least 2.5 cm (1 inch) away from the insulin pump infusion site or manual injection site. When replacing the sensor, select a new site that is at least 2.5 cm (1 inch) away from the previous site to allow the area to heal. Any scarring or hardening of the tissue can cause inaccurate insulin delivery or sensor performance.

CAUTION: Avoid the 5.0 cm (2 inch) area around the navel to help ensure a comfortable insertion site and to help with sensor adhesion.

Do not insert the sensor in muscle or areas constrained by clothing or accessories, areas with tough skin or scar tissue, sites subjected to rigorous movement during exercise, or in sites under a belt or on the waistline for best sensor performance and to avoid accidental sensor removal.

Inserting the sensor

WARNING: Wear gloves when inserting the sensor into someone other than yourself to avoid contact with patient blood. Minimal bleeding may occur. Contact with patient blood can cause infection.

1. Wash your hands.
 2. Choose an insertion site that has an adequate amount of fat, either on the abdomen, the lower back, or the upper buttock areas. Shaded areas show the best sites for sensor insertion.
 3. Clean the insertion site with alcohol. Let the area air dry.
 4. Open the sensor package.
 5. Remove the pedestal from the package. Place the pedestal on a flat surface.
 6. Make sure that the Enlite sensor's adhesive tab is tucked under the sensor connector and sensor snaps.
 7. Place your thumb on the thumbprint marking to hold theserter without touching the buttons.
 8. Grip theserter, placing your thumb on the thumbprint marking, without holding the buttons. Carefully push theserter down onto the pedestal until the base of theserter sits flat on the table and you hear a click.
 9. To detach theserter from the pedestal, using one hand, place your thumb on the thumbprint marking to grip theserter without holding the buttons.
- With your other hand, place two fingers on the pedestal arms and slowly pull theserter straight up. Do not detach the pedestal from theserter in midair, as this might damage the sensor.

Note: The arrow on the side of theserter aligns with the needle inside theserter.

WARNING: Never point a loadedserter toward any body part where insertion is not desired. An accidental button-push may cause the needle to inject the sensor in an undesired location, causing minor injury.

10. Hold theserter steady against your cleaned insertion site, without pushing theserter too deep into your skin. Press and release the bump on both buttons at the same time, while holding theserter flat against your body. Continue holding theserter flat against your body to allow the adhesive to stick to your skin. Slowly lift theserter away from your body, making sure that the buttons are not pressed.

Note: Failing to hold theserter securely flat against your body may allow theserter to spring back after pressing the buttons, and result in improper insertion of the sensor.

If you inserted the sensor into yourself, complete step 11a. If you are a healthcare professional or caregiver who inserted the sensor into a patient, complete step 11b.

Patient:

- 11a. Gently hold the base of the sensor against skin. Hold the needle housing at the top and slowly pull straight out, away from the sensor.

Healthcare professional or caregiver:

- 11b. Wrap sterile gauze around the sensor. Gently hold the base of the sensor against skin. Hold the needle housing at the top and slowly pull straight out, away from the sensor.

WARNING: Watch for bleeding at the insertion site. If bleeding occurs under, around, or on top of the sensor, apply

steady pressure using sterile gauze or a clean cloth placed on top of the sensor for up to three minutes. The use of unsterile gauze can cause an infection. If bleeding does not stop, remove the sensor and apply steady pressure until the bleeding stops.

Note: Medtronic adhesives are pressure-sensitive. Pressing the adhesive against the skin ensures that the sensor remains adhered to the skin throughout the wear period.

12a. Hold the sensor in place and gently remove the adhesive liner from under the adhesive pad. Do not remove the adhesive liner from the rectangular adhesive tab. This tab will be used to secure the transmitter in a later step.

12b. Firmly press the adhesive pad against the skin to make sure that the sensor remains adhered to the skin.

13a. Untuck the adhesive tab from under the connector.


13b. Straighten the sensor adhesive tab so that it lies flat against the skin.

Connecting the transmitter to the sensor

Before proceeding, have your pump user guide available.

To connect the transmitter to the sensor:

1. After the sensor is inserted, consult yourserter user guide for details on applying the required overtape.
2. Hold the rounded end of the inserted sensor to prevent it from moving during connection.
3. Hold the transmitter as shown. Line up the two notches on the transmitter with the side arms of the sensor. The flat side of the transmitter should face the skin.
4. Slide the transmitter onto the sensor until the sensor's flexible arms snap into the notches on the transmitter. If the transmitter is properly connected, and if the sensor has had enough time to become hydrated, the green light on the transmitter will flash within 10 seconds.

5. If the transmitter light does not flash, disconnect the transmitter from the sensor, wait for several seconds and then reconnect. If the transmitter light still does not flash, charge the transmitter.


6. When the transmitter light flashes green when connected to the sensor, use your pump to start the sensor. For more instructions, see your pump user guide.
7. After the transmitter successfully sends sensor data to the pump, attach the sensor's adhesive tab to the transmitter.
8. Follow the instructions that appear on the pump screen or follow the instructions in your pump user guide.

Disconnecting the transmitter from the sensor

Before proceeding, have your pump user guide available.

To disconnect the transmitter from the sensor:

1. Carefully remove any occlusive dressing from the transmitter and sensor.
2. For the Enlite sensor, remove the adhesive tab from the top of the transmitter.
3. Hold the transmitter as shown, and pinch the flexible side arms of the sensor between your thumb and forefinger.

4. Gently pull the transmitter away from the sensor.
5. Follow the instructions that appear on the pump or follow the instructions in your pump user guide.

Removing the sensor

Always refer to the sensor user guide for instructions on how to remove the sensor.

Bathing and swimming


After the transmitter and sensor are connected, they form a waterproof seal to a depth of 2.4 meters (eight feet) for up to 30 minutes. You can shower and swim without removing them. No occlusive dressing or overtape is needed.

Watertight Tester

The tester is used to test the transmitter to make sure it is working. It is also used as a required component for cleaning the transmitter. Properly connecting the tester to the transmitter will ensure that fluids do not come in contact with the


transmitter's connector pins. Fluids can cause connector pins to corrode and affect the transmitter's performance. Do not twist the tester while attached to the transmitter. This will damage the transmitter. The tester can be used for one year. If you continue to use the tester for more than one year, the transmitter's connector pins could be damaged, because the tester cannot continue to provide a waterproof seal.

Caution: Only use the tester with the transmitter. Do not use any other test plug.

Inspecting the transmitter connector pins

This image is an example of how the connector pins should look.

Look inside the transmitter's connector opening to make sure that the connector pins are not damaged or corroded. If the connector pins are damaged or corroded, the transmitter cannot communicate with the charger or pump. Contact the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.

Also look for moisture inside the connector opening. If you see any moisture, allow the transmitter to dry for at least one hour. Moisture inside the connector opening could cause the transmitter to not work properly, and could cause


corrosion and damage over time.

Connecting the tester for testing or cleaning

Before proceeding, have available your pump user guide.

To connect the tester:

- 1 Hold the transmitter and the tester as shown. Line up the flat side of the tester with the flat side of the transmitter.
- 2 Push the tester into the transmitter until the flexible side arms of the tester click into the notches on both sides of the transmitter.


- 3 Within five seconds, when properly connected, the green light on the transmitter flashes for about 10 seconds.
- 4 To test the transmitter, check the sensor icon on the pump to ensure that the transmitter is sending a signal (see your pump user guide).
- 5 After testing or cleaning, disconnect the tester from the transmitter.

Disconnecting the tester

To disconnect the tester:

- 1 Hold the transmitter body as shown and pinch the side arms of the tester.
- 2 With the tester arms pinched, gently pull the transmitter away from the tester.

Note: To save transmitter battery life, do NOT leave the tester connected after cleaning or testing.

Use conditions

The Guardian 2 Link Transmitter Kit is intended to be used by a patient or health professional. They are intended to be used in health care and home environments.

Principles of operation of medical device and their features


Transmitter

The Medtronic MiniMed development name for the Guardian 2 Link Transmitter is GST3B. The development name, GST3B, is referenced in supporting documentation. Guardian 2 Link Transmitter is indicated for single-patient use as a component of select Medtronic continuous glucose sensing systems and MiniMed sensor-enabled pump systems. The device is intended to process, store, and transmit glucose sensor values to data collection and display devices. Guardian 2 Link Transmitter is not intended to function as a “stand-alone” device.

Medtronic continuous glucose sensing systems or sensor-enabled pump systems consist of a glucose sensor, transmitter device, and receiving device (e.g. pump or monitor). The glucose sensor is inserted into the subcutaneous tissue and is affixed to the skin of the user. The transmitter is then connected to the glucose sensor and secured. The transmitter reads, processes the sensor current, applies the calibration algorithm, and transmits a glucose value via radio frequency (RF) communication to the receiving device which displays and stores the information. Data stored by the receiving device may be further uploaded to a therapy management software application residing on a personal computer for retrospective review and analysis.

Continuous glucose sensing systems and sensor-enabled pump systems that include the Guardian 2 Link Transmitter are intended to provide the user with the convenience of viewing real-time glucose values that can be analyzed to track patterns and improve overall diabetes management. Real-time glucose values are not intended to be used directly for making therapy adjustment but rather to provide an indication that confirmation finger stick with a home blood glucose monitor that may be required. All insulin therapy adjustments are based on the home blood glucose monitor, not on the glucose values displayed by the receiving device.

The Guardian Link Transmitter is a component of the Continuous Glucose Monitoring (CGM) system for the MiniMed 640G insulin pump. Figure below is a graphical representation of the Guardian 2 Link Transmitter and the MiniMed 640G insulin pump.

MD Models/ Versions

Guardian 2 Link Transmitter Kit (model MMT-7775), components:

- **Guardian 2 Link Transmitter (MMT-7731) – 1 pcs.;**
- **Charger with AAA Battery (MMT-7715) – 1 pcs.;**
- **Watertight Tester (MMT-7726) – 1 pcs.;**
- **One-Press Sertter Insertion Device (MMT-7512) – 1 pcs.;**
- **Product literature (User Guide Guardian 2 Link Transmitter, User Guide One-Press Sertter Insertion Device, Enlite Sensor Performance, warranty)**


Device components and list of parts and accessories


The Guardian 2 Link Transmitter Kit identified as model MMT-7775 includes the following device and accessory items listed in Table below:

Model/Item Number	Device/Accessory	Description
MMT-7731	Guardian 2 Link Transmitter	Guardian 2 Link Transmitter consists of an oval-shaped polycarbonate case, housing a battery, printed circuit board assembly (PCBA), and an antenna. The Transmitter incorporates a start-up circuit that detects when the transmitter is connected to an inserted sensor. After connection to the inserted sensor, the transmitter automatically detects the capacitance of the inserted sensor and then subsequently initializes the sensor. The Transmitter sends glucose data to the receiving device via RF signal after initialization has completed. The Guardian 2 Link Transmitter powers the glucose sensor and

measures the sensor signal current (I_{SIG}) from the glucose sensor. (The I_{SIG} is an electrical current level that is proportional to the glucose level in the subcutaneous interstitial fluid of the user.) The I_{SIG} value is converted to a digital signal and is filtered to reduce noise artifacts. The Transmitter software application then applies the sensor calibration algorithm to the I_{SIG} value to calculate sensor glucose values and sensor glucose rate of change values. Glucose alerts are also identified. The Transmitter software application converts 1-minute sensor current samples into 5-minute packets. The Transmitter communicates data packets which contain sensor values and alerts to the receiving device through a radio frequency (RF) link once every 5 minutes.

The Guardian 2 Link Transmitter receives calibration values (BG meter readings) via radiofrequency (RF) communication from the receiving device (pump). When a calibration value is received, the Guardian 2 Link Transmitter applies the calibration algorithm to calculate the updated Calibration Factor. The updated Calibration Factor is then applied to the I_{SIG} value to calculate sensor glucose values.


MMT-7715

Charger with AAA Battery

The Charger is an accessory device intended to charge the rechargeable lithium polymer battery contained within the Guardian 2 Link Transmitter. The Transmitter is charged automatically following connection to the Charger. The Charger automatically terminates charging upon full charge.


The Charger is powered by an internal, replaceable, commercially available AAA or LR-03 alkaline battery and is designed with sliding rear cover to facilitate changing the battery. The Charger uses visual indicators (green/red LEDs) to show the status of transmitter device charging. The indicator light also communicates any problem with the transmitter device battery, or the charger battery.

MMT-7726

Watertight Tester

The Watertight Tester is an accessory device intended to facilitate troubleshooting and the cleaning of Guardian 2 Link Transmitter. The

			Watertight Tester is assembled from a plastic housing and a printed circuit board (PCB). The Tester is designed with a molded handle and plug that mate with a transmitter device connector in the same manner that a sensor would be attached to the transmitter device. The Tester operates as a sensor simulator creating a signal current at a level that is within the range of an in-vivo sensor during normal operation. The Tester PCB contains a resistor. The Tester connects the device counter and reference connections. Initialization of operating conditions occurs in the Transmitter after being connected to the Tester. The Tester resistor is the connection between the reference and working electrodes through which a constant current value (I_{SIG}) of 53.5 nA (nanoamps) flows and is read by the Transmitter device. The known current value provides the ability to confirm function of the attached transmitter device. The Tester also provides a watertight plug for cleaning and/or disinfection of the transmitter device. The Tester is used as a handle to aid in holding the transmitter device during cleaning and disinfection. 
	MMT-7512	One-Press Sertter Insertion Device	The One-Press Sertter (Model MMT-7512) has been designed for aiding the insertion of the Enlite Glucose Sensor (Model MMT-7008). It is intended to be used by a patient or a clinician to introduce the Glucose Sensor into the subcutaneous tissue at a fixed depth, with minimal discomfort and technique dependency, and with minimal exposure of the sensor needle. The device consists of a plastic barrel containing a stainless steel spring

and sear springs on the other side. The Serter is loaded and locked by lowering it onto the sensor pedestal and pushing the striker into the frame until it locks in place. The serter is fired and the sensor is released by depressing both green buttons simultaneously, while pressing the bottom of the serter flat against the skin, over the insertion site.

The possibility and ways of integrating with other medical devices (device compatibility) (if applicable)

The Guardian 2 Link Transmitter is packaged as system and requires the use of the Charger (MMT-7715) and Watertight Tester (MMT-7726). The replacement transmitter is packaged as a stand-alone device. Below lists accessories and

compatible device for the Guardian 2 Link Transmitter:

Model Number	Device Name
MMT-1511, MMT-1512 MMT-1711, MMT-1712	MiniMed 640G Insulin Pump
MMT-7008	Glucose Sensor (Enhanced Enlite)
MMT-7715	Charger
MMT-7726	Watertight Tester

The One-Press Serter (MMT-7512) is sold independently and is intended for use with the following products:
Enlite Glucose Sensor (MMT-7008)

The Charger (MMT-7715) is intended to recharge the lithium polymer battery contain within the following devices:

Model Number	Device Name
MMT-7703	MiniLink Transmitter
MMT-7763	Guardian Link Transmitter
MMT-7731	Guardian 2 Link Transmitter

Other devices are in development at Medtronic MiniMed for use with the 7715 Charger.

The Tester (MMT-7726) is intended for use with the following devices:

Model Number	Device Name
MMT-7703	MiniLink Transmitter
MMT-7763	Guardian Link Transmitter
MMT-7731	Guardian 2 Link Transmitter

Other devices are in development at Medtronic MiniMed for use with the Tester.

The materials from which it is made or its medical device basic parts

The Guardian 2 Link Transmitter Hardware Materials

Component	Material, Material Specification	Manufacturer	Contact Type
HSG, Top	Polycarbonate Bayer Makrolon 2558 + colorant CCP03-124916SS-886975 LT GRAY	Clariant Corporation, USA RTP Company, USA	1) There is no contact with the patient (attached to the glucose sensor that is on the patient's skin 2) Short-term
HSG, Bottom	Polycarbonate Bayer	Clariant Corporation,	

	Makrolon 2558 + colorant CCP03-124916SS-886975 LT GRAY	USA RTP Company, USA	contact with undamaged skin of the hands (at the time of installation)
Adhesive	3M 467MP - 100%	3M, USA	

The One-Press Sertter Materials

Component	Material, Material Specification	Manufacturer	Contact Type
Housing (pressure construction)	Polycarbonate Bayer Makrolon 2558	Clariant Corporation, USA	Short-term contact with undamaged skin of hands
Mandrin (inside the product)	Polycarbonate Lubriloy D2000-GY3E269	SABIC Innovative Plastics, USA	
	+ white colorant Pantone White 1309	Stehlin Hostag Ink UK Ltd., UK	
Needle holder (inside the product)	Polycarbonate Lubriloy D2000-GY3E269	SABIC Innovative Plastics	
	+ white colorant Pantone White 1309	Stehlin Hostag Ink UK Ltd., UK	
Left, right button	Polycarbonate Bayer Makrolon 2558	Clariant Corporation, USA	
	+ green colorant 2A16632CMB GREEN 94	POLYONE CORPORATION, USA	
Trigger left, right (inside the product)	Polycarbonate Lubriloy D2000-GY3E269	SABIC Innovative Plastics, USA	
	+ white colorant Pantone White 1309	Stehlin Hostag Ink UK Ltd., UK	

The Charger Materials

Component	Material, Material Specification	Manufacturer	Contact Type
Top side of the case	ABS plastic, brand S-885547 RTP 600 ABS 348 DK GRAY + gray colorant 10013617 PANTONE GRAY 425C	RTP Company, USA POLYONE CORPORATION, USA	Short-term contact with undamaged skin of hands
Bottom of the case	ABS plastic, brand S-885547 RTP 600 ABS 348 DK GRAY + gray colorant 10013617 PANTONE GRAY 425C	RTP Company, USA POLYONE CORPORATION, USA	
Battery cover	ABS plastic, brand S-885547 RTP 600 ABS 348 DK GRAY + gray colorant 10013617 PANTONE GRAY 425C	RTP Company, USA POLYONE CORPORATION, USA	
The nameplate on the battery compartment cover	polyethylene terephthalate LaserLIGHTS + gray colorant 10013617 PANTONE GRAY 425C	Rowmark POLYONE CORPORATION, USA	
Light indicators	Polysulfone Udel P-1700 + red colorant APPLE RED + green colorant 2A16632CMB GREEN 94	Solvay Specialty Polymers, USA POLYONE CORPORATION, USA POLYONE CORPORATION, USA	
Adhesive	3M P1212 - 100%	3M, USA	

The Watertight Tester Materials

Component	Material, Material Specification	Manufacturer	Contact Type
Tester body	Polycarbonate Bayer Macrolon 2558 + gray colorant CCP03- 124916SS-886975 LT GRAY	Clariant Corporation, USA RTP Company, USA	Short-term contact with undamaged skin of hands
O-rings	каучук ELASTOSIL® E 41 + black colorant PX-V- 190 PANTONE BLACK NEW 180 DEGREE	Wacker Silicones, USA PolyOne, USA	

Labeling, Packaging and Sticker Information

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



Labelling of the Transmitter


Labelling of the Tester



The One-Press Sertter is packaged in a plastic cover, which does not have a label. The package labels are applied to the customer boxes, which include either the Guardian 2 Link Start-up kit or the One-Press Sertter Standalone (Replacement) Kit.

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Packaging


The Guardian 2 Link Transmitter and accessories are placed in a folded cardboard tray and packaged in a cardboard customer carton.

For customer convenience the Guardian 2 Link Transmitter and accessories are offered as a system kit and is identified as a separate model number (MMT-7775). The Guardian 2 Link Transmitter kit identified as model MMT-7775 includes the following device and accessory items listed below.

The package contains:

- Guardian 2 Link Transmitter (MMT-7731)
- Charger with AAA Battery (MMT-7715)
- Watertight Tester (MMT-7726)
- One-Press Sertter Insertion Device (MMT-7512)
- Product literature (User Guide Guardian 2 Link Transmitter, User Guide One-Press Sertter Insertion Device, Enlite Sensor Performance, warranty)*

* Appendix to the Manual (IFU) not in package content, printed and shipped separately on the territory of Russian Federation

Symbols

	Serial number
	Catalogue or model number
	One per container/package
	Date of manufacture
	Manufacturer
	Caution
	Follow instructions for use (appears blue on label)
	Consult instructions for use
	Storage temperature range

	CE 0459	Marking of Conformity: This symbol means the device fully complies with MDD 93/42/EEC (NB 0459).
		Signifies compliance with Australian EMC, EME and Radio communications requirements
		Radio communication
	CONF	Configuration or unique version identifier
		Degree of protection against electric shock: Type BF applied part
	IC	Signifies compliance with Industry Canada EMC and Radio communications requirements
	IP48	Transmitter: Protected against the effects of continuous immersion in water (2.4 meters (8 feet) immersion for 30 minutes).
		Storage humidity range
	EC REP	Authorized representative in the European community

LOT

Batch code

This information is provided by MDT Russia

Information about sticker:

The final packaging for the Russian market is accompanied by additional sticker:


1. name of medical device in Russian language;
2. name of the manufacturer;
3. country of origin;
4. code UPN (unique product number);
5. code CFN (customer facing number);
6. date and number of Registration Certificate;
7. sign of conformity in Russia ;
8. warning and informing symbols and inscriptions (if applicable).

Sample of the sticker (provided by MDT Russia)

	Product name	
	Name of the manufacturer	Code UPN
	Country of origin	Code CFN
	See IFU, Date and Number of Registration Certificate	

The main technical parameters and characteristics of medical device (rate-based specifications, lifetime, expiration date).

The Guardian 2 Link Transmitter Specifications




Dimensions	3.57 cm $\pm$ 0.20 cm (w) x 2.84 cm $\pm$ 0.15 cm (l) x 0.93 cm $\pm$ 0.01 cm (h) nominal
Weight	4.76 grams nominal ($\pm$ 0.3 g)
Data Storage	Transmitter shall store the data from at least the 10 most recent hours.
Accuracy	$\pm$ 6% between 2 nA and 10nA $\pm$ 3% from 10 nA to 200nA
Operating conditions	Temperature: -5°C to 45°C (23°F to 113°F) Caution: When operating the transmitter on a tester in air temperatures greater than 41°C (106°F), the temperature of the transmitter may exceed 43°C (109°F) Relative humidity: 5% to 95% with no condensation Pressure: 61.36 to 106.17 kPa (8.9 to 15.4 psi)
Storage conditions	Temperature: -25°C to 55°C (-13°F to 131°F) Relative humidity: 10% to 100% with condensation Pressure: 61.36 to 106.17 kPa (8.9 to 15.4 psi)
Minimum battery life (during normal use) before requiring recharge	Six days
Transmitter frequency	2.4 GHz, 2M65G1D modulation, less than 1mW ERP
Maximum output power (EIRP)	-0.63 dBm
Radio Frequency (RF) communications	Pump to transmitter frequency: 2.4 GHz; proprietary Medtronic protocol; range up to 1.8 meters (6 feet)

Utilizes the IEEE 802.15.4 protocol with proprietary data format.
 Operating frequency: 5 frequencies are used: 2420, 2435, 2450, 2465, and 2480 MHz
 Bandwidth: 5 MHz, which is allocated channel bandwidth per IEEE protocol

Transmitter expected service life The transmitter expected service life is 1 year depending on patient usage


The Charger Specifications

Dimensions	6.1 cm (± 0.3 cm) x 4.1 cm (± 0.2 cm) x 2.9 cm (± 0.1 cm) nominal
Weight	28.2 grams nominal (±1.5 g)
Operating conditions	Temperature: 10 °C to 40 °C Relative humidity: 30% to 75% with no condensation
Storage conditions	Temperature: -10 °C to +50 °C



	Relative humidity: 10% to 95% with no condensation
Shelf Life	There is no expiration date, as the device is not provided as sterile.
Operating Life	The Charger is intended to operate properly for one year excluding user caused damage.
Power Requirements	Standard 1.5V AAA or LR-03 battery.

The One-Press Sertter Specifications

Dimensions	Approximate dimensions: 7.85 cm (± 0.40 cm) x 6.91 cm (± 0.30 cm) x 5.77 cm (± 0.28 cm) nominal
Weight	Approximate weight: 112.5 grams (± 5.6 g)
Operating conditions	The Sertter shall operate within a normal room temperature of 10 °C to 40 °C
Storage conditions	Temperature: -20 °C ± 2 °C to 55 °C ± 2 °C Atmospheric Pressure: 49.64 kPa ± 3.5 kPa to 106.18 kPa ± 3.5 kPa (7.2 psia ± 0.5 psia to 15.4 psia ± 0.5 psia). Relative Humidity: 15% to 85% condensing
Shelf Life	There is no shelf life requirement associated with the Sertter.
Operating Life	The sertter has a minimum service life of one year.

The Watertight Tester Specifications

Dimensions	3.0 cm (± 0.2 cm) x 2.0 cm (± 0.1 cm) x 0.7 cm (± 0.1 cm) nominal
Weight	approx. 1.48 grams (± 0.1 g)
Operating conditions	Operating range, temperature: -5°C to $+45^{\circ}\text{C}$ Operating range, humidity: 5% to 95% RH
Storage conditions	Storage temperature: -30°C to $+70^{\circ}\text{C}$ Storage humidity: 10% to 100% RH with condensation
Shelf Life	There is no expiration date, as the device is not provided as sterile.
Operating Life	The Tester is intended to operate properly for 122 uses excluding user caused damage.

Software

The Guardian 2 Link Transmitter is intended to process, store, and transmit glucose sensor values to data collection and display devices. The device contains three software components which reside in separate microprocessors. Microprocessors STM32 and STM8L contain software modules which are built separately and run on separate microcontrollers but are designed as a single software application. Both modules work together to fulfill the main software application requirements. The third microprocessor (CC2430) contains software responsible for radiofrequency (RF) communication.

The STM8L is an 8-bit processor which acts as a co-processor to the STM32 in the Guardian 2 Link design. The STM8L is a low power processor that enables it to remain on permanently and provide functionality the STM32 cannot provide while maintaining a low power state. Responsibilities of the STM8L include: multiplexing external inputs for the STM32, counting the output from the I-to-F converter, controlling watchdogs for STM8L and STM32, and controlling the power domains of the transmitter. The STM8L software module is used as multiplex interrupts for the STM32 microprocessor and to perform long-running sensor data gathering tasks while the STM32 is in low power mode. The STM8L periodically wakes up the STM32 in order to perform certain tasks, including sending/receiving RF data to/from a compatible device (pump) and performing sensor calibration logic. When the STM8L wakes up the STM32, it communicates over an inter-integrated circuit (I²C) to exchange information about the reason for waking up, sensor data, etc.

The STM32 acts as the main microprocessor for the Guardian 2 Link Transmitter. Responsibilities of the STM32 include performing sensor calibration, communicating with the RF chip, and performing timekeeping to determine when periodic readings occur. The STM32 software module receives sensor data from the STM8L over an I²C bus. Sensor data received by the STM32 software module is then filtered and applied to the glucose calibration algorithm using BG meter values received from a compatible device (pump) via RF communication. Output of this filtering and calibration logic is sensor glucose value, which is transmitted every 5 minutes via RF to a compatible device. Sensor glucose values and alert flags are transmitted and used by a compatible device (pump) to trigger sensor glucose high/low/predictive display alerts.

The STM32 microprocessor includes the following two software components:

Main Application – microprocessor that executes sensor algorithm.

Bootloader - considered an off the shelf software or software of unknown province (SOUP) which is embedded in the internal boot ROM of the STM32 microcontroller. This software is used only during manufacturing to facilitate loading the STM32 application software. The software is not intended to perform any function with respect to the performance of the device.

The Tel-D (Telemetry for Diabetes) subsystem is based on a CC2430 system-on-a-chip processor which consists of a 2.4 GHz IEEE 802.15.4 compliant RF transceiver and a low power 8051 microcontroller core. The CC2430 communicates with the STM32 using a serial peripheral interface (SPI). The CC2430 contains the Tel-D software which controls the SPI communications with the host processors (STM32 in Guardian 2 Link Transmitter) and other Medtronic devices within the Tel-D subsystem over the 2.4GHz radio.

1.1 Overview of GST3B Algorithm

The GST3B calibration algorithm is illustrated in **Figure 3** below. Process steps on the bottom row improve reliability of sensor glucose values and trending information. From the left, sensor measurements are filtered to create a current value (I_{SIG}) every five minutes. After conversion of I_{SIG} data to five minute packets, the I_{SIG} signal is then checked to ensure it is in a reasonable range and then subsequently processed in the 2nd stage filter.

Final Calibration Factor calculations occur 10-15 minutes following a finger stick entry. In most cases, a temporary Cal Factor is calculated 0-5 minutes after the finger stick entry as well. Calculations use the most recent finger stick entries and F_{ISIG} values to determine the appropriate Cal Factor value. Once calculated, the same value is used until the next finger stick entry, unless it follows an abnormal calibration or is in first day calibration adjustment mode.


Figure 3— Overview of GST3B Calibration Algorithm

GST3B contains the following modifications / additions in relation to the MiniLink (MMT-7003) and the most recently CE Marked transmitter, GST3A2' (MMT-7763A) algorithm:

1. Tuning – accommodate higher I_{SIG} values generated by Enlite / Enlite Enhanced
 - Offset
 - Cal Ratio Limits / Calibration Error Thresholds
 - Calibration Factor Buffer Weighting Calculation
2. Calibration Error – improved error handling

- Calibration Error Timing
- 3. Calibration Updates – reduced delay in updates
 - Calibration Factor Update
 - Calibration Adjustment
 - First Day Calibration Adjustment
- 4. Artifact Detection / Response – faster response to rapid glucose excursions
 - Artifact Detection
 - I_{SIG} signal processing and filtering
 - Low I_{SIG} detection
 - Lag Correction
- 5. System Initialization Periods – shorter warm-up times
 - Initialization Period
- 6. Sensor Termination – identification of sensors that are no longer viable
 - Noise-based checks
 - Recovery Mode

1.2 GST3B Software Algorithm (Tuning Modifications)

Enlite / Enhanced Enlite sensors generate higher I_{SIG} currents than earlier generation glucose sensors, Sof Sensor (**MMT-7002 / MMT-7003**). Higher I_{SIG} currents result in lower Calibration Factors amplitudes. The equation below shows that an increase in I_{SIG} currents (denominator) results in a lower Calibration Ratio:

$$CR = (BG / I_{SIG})$$

Where,

CR = Calibration Ratio

BG = Blood glucose values

I_{SIG} = Current proportional to glucose values in patient's subcutaneous interstitial fluid.

Figure 4 below demonstrates the lower calibration factor amplitudes for Enlite / Enlite Enhanced in relation to Sof Sensor.


Figure 4 – Cal Factor Amplitudes (Enlite vs. Sof Sensor) – Low I_{SIG} Detection / Tighter Calibration Factor Limits Accommodate Larger I_{SIG} Currents with Enlite / Enlite Enhanced

Tuning modifications to accommodate the change in sensor characteristics between Enlite / Enlite Enhanced and Sof Sensor are listed in the **Table 3** below:

Table 3 – GST3B Tuning Modifications (Cal Factor Buffer Weighting Calculation)

Feature Description	MiniLink Transmitter	Guardian Link Prime	Guardian 2 Link

		Calibration Factor Buffer Weighting Calculation	Alpha = Exponentially Decaying Function Beta = Inverse Variance Look-Up Table	Alpha = Exponentially Decaying Function Beta = Inverse Variance Look-Up Table	Alpha = Constant Values Beta = Function of Blood Glucose and Expected Calibration Factor (5.0 Expected Slope). Modification is intended to take advantage of improved sensor stability (Enlite Enhanced) and provide more accurate data. The feature reduces errors due to inaccurate Blood Glucose (BGs) or temporary sensitivity changes.	
					Figure 5 below shows modifications to the weighting algorithm. Figure 6 demonstrates an improved glucose measurement profile with the new weighting logic.	

Offset	2nA subtracted from background I _{SIG} current	2nA subtracted from background I _{SIG} current	3nA subtracted from background I _{SIG} current Tuning modification matches sensor characteristics better, resulting in better detection of hypo- and hyperglycemic conditions.
Cal Ratio Limits / Calibration Error Thresholds	1.5 to 20 mg/dL/nA	1.5 to 20 mg/dL/nA	2.5 to 15 mg/dL/nA (2.5 to 17 mg/dL/nA for instant calibration error check) Thresholds were tightened to account for better sensor stability with Enlite / Enlite Enhanced. Prevents display of sensor glucose values and prevents calibrations when the measured sensitivity is high or low.

$$CF = 20\% \times CF_{expected} + 80\% \times \frac{\sum_{i=1}^n \alpha_i \times \beta_i \times (fisig_i + offset) \times BG_i}{\sum_{i=1}^n \alpha_i \times \beta_i \times (fisig_i + offset)^2}$$

New for GST3B
Modified for GST3B

Figure 5 – GST3B Calibration Factor Buffer Weighting Calculation Change


Figure 6 – GST3B CF Weighting Logic – Comparison of GST3B Algorithm with and without the Expected Cal Factor Weight)

1.3 **GST3B Software Algorithm (Artifact Detection)**

Artifact detection features prevent system noise related to rapid glucose excursions are listed in **Table 4** below:

Table 4 – GST3B Artifact Detection Features

		Feature Description	MiniLink Transmitter	Guardian Link	Guardian 2 Link
		Artifact Detection	No detection	Detection of sensor signals with larger amplitudes than system noise (artifacts) and are associated with non-physiologic glucose rates of change. This feature prevents calibration and Low Glucose Suspend (LGS) trigger on an artifact signal.	Detection of sensor signals with larger amplitudes than system noise (artifacts) and are associated with non-physiologic glucose rates of change. This feature prevents calibration and Low Glucose Suspend (LGS) trigger on an artifact signal.
		Low I _{SIG} Detection	I _{SIG} values less than 2nA are flagged as invalid to prevent use in calculating sensor glucose values.	I _{SIG} values less than 6nA are flagged as invalid to prevent use in calculating sensor glucose values.	I _{SIG} values less than 6nA are flagged as invalid to prevent use in calculating sensor glucose values. This feature prevents abnormally low I _{SIG} signals from being used due to transient low sensor sensitivity. A sensor error is triggered when low I _{SIG} is detected.

		Lag Correction	No lag correction	No lag correction	Lag correction
					Modified algorithm to account for delay and handling of small I_{SIG} dropouts (Refer to Figure 7).
		(I_{SIG}) Signal Processing and Filtering	Two (2) Stages of Filtering: Stage 1: eight (8) previous one-minute counts are converted into a five-minute I_{SIG}. Stage 2: eight (8) previous five-minute I_{SIG} are filtered ($f I_{SIG}$) using an 8 point Finite Impulse Response (FIR) Gaussian filter.	Two (2) Stages of Filtering: Stage 1: eight (8) previous one-minute counts are converted into a five-minute I_{SIG}. Stage 2: eight (8) previous five-minute I_{SIG} are filtered ($f I_{SIG}$) using an 8 point Finite Impulse Response (FIR) Gaussian filter.	Two (2) Stages of Filtering: Stage 1: eight (8) previous one-minute counts are converted into a five-minute I_{SIG}. Transmitter handles high noise conditions better, which reduces the need for an aggressive 1st stage filter and reduces delay. Stage 2: eight (8) previous five-minute I_{SIG} are filtered ($f I_{SIG}$) using an 8 point adaptive Kalman filter. The adaptive Kalman filter has improved processing of signal noise and rate of change better, which reduces the need for an aggressive 1st stage filter and a reduction in delay.


Figure 7 – GST3B (Lag Correction) Reduces Delay on Small I_{SG} Dropouts (Example of Signal Dropouts identified by comparing the I_{SG} to companion sensors)

1.4 GST3B Software Algorithm (Calibration Factor Update)

Updates related to the calibration factor for GST3B in relation to Enlite / Enlite Enhanced are listed in **Table 5** below:

Table 5 – GST3B Calibration Factor Update Features

Feature Description	MiniLink Transmitter	Guardian Link	Guardian 2 Link
Calibration Factor Update	Calibration Factor updated 10-20 minutes after BG entry	Calibration Factor updated 10-20 minutes after BG entry	Calibration updated 5 minutes after BG entry and again 10-15 minutes after BG entry This feature reduces the time required to update sensor glucose values following a calibration attempt. (Refer to Figure 8)
Calibration Adjustment	Calibration Factor unchanged until next calibration	Calibration Factor unchanged until next calibration	Calibration Factor adjustment between BGs following abnormal calibrations (Refer to Figure 9)

First
Calibration
Adjustment

No adjustments

No adjustments

First day calibration
adjustment – Selective
adjustment reduces
values in the first 10
hours due to unexpected
sensor sensitivity
following the
initialization period
(Refer to Figure 10).


Figure 8 – GST3B (Calibration Factor Update) – Comparison of Sensor Glucose Values (Early Calibration Update vs. No Early Calibration Update)


Figure 9 – GST3B (Calibration Adjustment) – Abnormal calibration at hour 6. GST3B Calibration Factor is gradually adjusted, resulting in a reduction in the error at hour 20. Black dots are BG measurements; Circles indicate a calibration on the point.


Figure 10 – GST3B (First Day Compensation – FDC) – Algorithm readjusts the CF on every sensor reading (top plot), while the CF for the A2 algorithm remains constant until the next calibration. After 6 hours, a BG measurement indicates the B algorithm reduced the error present in the A2 algorithm.

1.5 **GST3B Software Algorithm (Calibration Error Timing)**

Calibration error timing features in relation to the Enlite / Enhanced Enlite are listed in **Table 6** below:

Table 6 – GST3B Calibration Error Timing Features

Feature Description	MiniLink Transmitter	Guardian Link	Guardian 2 Link
Calibration Error Timing	Cal Errors sent 10-20mins after Blood Glucose (BG) entry. Change Sensor: 10-15 min after second BG entry	Cal Errors sent 10-20mins after Blood Glucose (BG) entry. Change Sensor: 10-15 min after second BG entry	Instant Calibration Error Alert: Errors sent immediately upon BG entry (Refer to Figure 11) Change Sensor: 10-15 min after second BG entry that would cause a “Change Sensor Alarm.” Modification to calibration error timing is intended to no notify user of error immediately, as opposed to waiting 10-20 minutes.


Figure 11 – GST3B (Instant Calibration Error Alert) – Immediately identifies a calibration error

1.6 GST3B Software Algorithm (Sensor Initialization Period)

Sensor initialization features in relation to the Enlite / Enhanced Enlite are listed in **Table 7** below:

Table 7 – GST3B Sensor Initialization Features

Feature Description	MiniLink Transmitter	Guardian Link	Guardian 2 Link
---------------------	----------------------	---------------	-----------------

Initialization Period	Fixed 2 hour	Fixed 2 hour	Variable Initialization Period: - Min = 35-40 minutes - Max = 2 hours Dependent on Isig current magnitude and change between consecutive packets. System looks for 3 consecutive high, stable signals to initialize at the minimum duration. Reduces delays to calibrate sensors and review sensor glucose values. CEP257 simulation demonstrates ~ 40% of sensors initialize prior to 2 hours.
-----------------------	--------------	--------------	--

1.7 GST3B Software Algorithm (Change Sensor Conditions)

Change Sensor conditions in relation to the Enlite / Enhanced Enlite are listed in **Table 8** below:

Table 8 – GST3B Change Sensor Features

Feature Description	MiniLink Transmitter	Guardian Link	Guardian 2 Link
---------------------	----------------------	---------------	-----------------

		Recovery Mode	No recovery mode	No recovery mode	Up to five (5) days from sensor start-up, if two (2) consecutive fingerstick readings result in calibration errors, sensor will enter recovery mode instead of initiating “Change Sensor Alarm.” Recovery mode results in “SG value not available” for up to two (2) hours. The intent of recovery mode is to allow the sensor time to return to a CF within threshold. Recovery mode may occur, at maximum, three (3) consecutive times. After the third consecutive recovery mode, the next calibration will result in a “Change Sensor Alarm.”	
--	--	---------------	------------------	------------------	--	--

	Sensor Conditions	result in a change sensor condition:	conditions result in a change sensor condition:	conditions result in a change sensor condition:
		1. Two (2) consecutive fingerstick entries, resulting in calibration errors (threshold or percentage errors) will initiate a "Change Sensor Alarm." 2. Six (6) days after sensor start-up	1. Two (2) consecutive fingerstick entries, resulting in calibration errors (threshold or percentage errors) will initiate a "Change Sensor Alarm." 2. Six (6) days after sensor start-up	1. After five (5) days from sensor start-up two (2) consecutive fingerstick entries, resulting in calibration errors (threshold or percentage errors) will initiate a "Change Sensor Alarm." 2. Six (6) days since sensor start-up 3. Algorithm analyzes signal magnitude and noise level. Prolonged periods of high noise and low signal will also result in "Change Sensor Alarm." Noise based sensor checks can terminate a sensor between calibrations (Refer to Figure 12). 4. After three (3) consecutive recovery modes, if the next

**Performance
methods and results.**

tests: Transmitter Radio Frequency (RF) communication

This device complies with the United States Federal Communications Commission (FCC) and international standards for electromagnetic compatibility.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

These standards are designed to provide reasonable protection against excessive radio frequency interference, and prevent undesirable operation of the devices from unwanted electromagnetic interference.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and the receiver.

This device can generate, use, and radiate radio frequency energy and, if installed and used in accordance with the instructions, may cause harmful interference to radio communications. If the device does cause interference to radio or television reception, you are encouraged to try to correct the interference by one or more of the following measures:

- Decrease the distance between the transmitter and the insulin pump to 1.8 meters (6 feet) or less.
- Increase the separation between the transmitter and the device that is receiving/emitting interference.

If other devices that employ radio frequencies are in use, such as cell phones, cordless phones, and wireless networks, they may prevent communication between the transmitter and the insulin pump. This interference does not cause any incorrect data to be sent and does not cause any harm to your devices. Moving away from, or turning off, these other devices may enable communication. If you continue to experience RF interference, please contact your local representative.

Caution: Changes or modifications to the internal RF transmitter or antenna not expressly approved by Medtronic could void the user's authority to operate this insulin delivery system.

**Guardian 2 Link wireless communication
Quality of service**

The Guardian 2 Link transmitter and the MiniMed 640G insulin pump are associated as part of an 802.15.4 network for which the pump functions as the coordinator and the transmitter as an end node. In an adverse RF environment the pump will assess channel changing needs based on "noise" levels detected during an energy scan. The pump will perform the energy scan if after 10 minutes no CGM transmitter signal has been received. If the channel change occurs the pump will send beacons on the new channel.

The Guardian 2 Link transmitter will initiate a channel search when beacon detection fails on the associated channel. The search will be conducted across all five channels. When the beacon is located the transmitter will rejoin on the identified channel. Upon re-association any missed packets (up to 10 hours) will be transmitted from the transmitter to the pump. In normal operation the transmitter will transmit a packet every 5 minutes and retransmit the packet if the data is corrupted or missed.

Data security

The MiniMed 640G insulin pump is designed to only accept radio frequency (RF) communications from recognized and linked devices (you must program your pump to accept information from a specific device).

The MiniMed 640G insulin pump and system components (meters and transmitters) ensure data security via proprietary means and ensures data integrity using error checking processes, such as cyclic redundancy checks.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The Guardian 2 Link is intended for use in the electromagnetic environment specified below. The customer or the user of the Guardian 2 Link should make sure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The Guardian 2 Link must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class B	The Guardian 2 Link is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guardian 2 and Manufacturer's Declaration - Electromagnetic Immunity

The Guardian 2 Link is intended for use in the electromagnetic environment specified below. The customer or the user of the Guardian 2 Link should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 2 kV, ± 4 kV, ± 8 kV Air ± 2 kV, ± 4 kV, ± 6 kV Indirect	± 8 kV Air ± 6 kV Indirect ± 22 kV Air, <5% RH	For use in a typical domestic, commercial, or hospital environment.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Requirement does not apply to this battery powered device.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable	Requirement does not apply to this battery powered device.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle	Not applicable	Requirement does not apply to this battery powered device.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The Guardian 2 Link is intended for use in the electromagnetic environment specified below. The customer or the user of the Guardian 2 Link should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	400 A/m 4000 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical domestic, commercial, or hospital environment.

Note: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The Guardian 2 Link is intended for use in the electromagnetic environment specified below. The customer or user of the Guardian 2 Link should assure that it is used in such an electromagnetic environment.

Immunity Test	IEC 60601 Level	Compliance Level	Electromagnetic Environment Guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Guardian 2 Link, including cables, than the recommended separation distance calculated from the equation applicable to the power of the transmitter. Refer to the recommended separation distance table for more information.
Conducted RF IEC 61000-4-6	3 V/m 150 kHz to 80 MHz	Not applicable	Not applicable

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The Guardian 2 Link is intended for use in the electromagnetic environment specified below. The customer or user of the Guardian 2 Link should assure that it is used in such an electromagnetic environment.

Immunity Test	IEC 60601 Level	Compliance Level	Electromagnetic Environment Guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 6 GHz	10 V/m 80 MHz to 6 GHz	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ 800 MHz to 6 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption, and reflection from structures, objects and people.

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcasts and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Guardian 2 Link is used exceeds the applicable RF compliance level above, the Guardian 2 Link should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Guardian 2 Link.

^bOver the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between the Guardian Link and common household radio transmitters

Household RF Transmitter	Frequency	Recommended Separation Distance (meter)	Recommended Separation Distance (inch)
Telephones			
Cordless Household	2.4 GHz	0.3	12
Cordless Household	5.8 GHz	0.3	12
TDMA-50 Hz (cell phone)	1.9 GHz	0.3	12
TDMA-50 Hz (cell phone)	800 MHz	0.3	12
PCS (cell phone)	1.9 MHz	0.3	12
DCS (cell phone)	1.8 MHz	0.3	12
GSM (cell phone)	900 MHz	0.3	12
GSM (cell phone)	850 MHz	0.3	12
CDMA (cell phone)	800 MHz	0.3	12
Analog (cell phone)	824 MHz	0.3	12
CDMA (cell phone)	1.9 MHz	0.3	12
WiFi Networks			
802.11b; 11Mbps maximum	2.4 GHz	1	39.5
802.11g; 54 Mbps maximum	2.4 GHz	1	39.5
802.11n; 11Mbps maximum	2.4 GHz	1	39.5
Bluetooth 500 kb/s	2.4 GHz	0.1	3.93
ZigBee 250 kb/s	2.4 GHz	0.1	3.93

Recommended separation distances between portable and mobile RF communications equipment and the Guardian 2 Link

The Guardian 2 Link is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the Guardian 2 Link users can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the Guardian 2 Link as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to the frequency of transmitter (m)		
	150 kHz to 80 MHz Not applicable	80MHz to 800MHz $d = 0.35 \sqrt{P}$	800MHz to 6.0GHz $d = 0.70 \sqrt{P}$
0.01	Not applicable	0.035	0.07
0.1	Not applicable	0.11	0.11
1	Not applicable	0.35	0.7
10	Not applicable	1.1	2.2
100	Not applicable	3.5	7

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Sensor/Charger Interface Connector	The transmitter connector shall show no distortions or damage following 244 charger insertion/removal cycles. The transmitter connector shall require an insertion force of less than 3 pounds when connecting to the sensor. The transmitter connector shall require an insertion force of less than 3 pounds when connecting to the charger. The transmitter connector shall have a charger retention force of more than 0.5 pounds and less than 3 pounds.
Drop	The transmitter withstand six physical drops from a height of one meter onto a 50mm thick (2 in) hardwood board (hardwood>700kg/m3) lying flat on a

	rigid base, one each from a starting altitude in each direction of the three mutually perpendicular axes (6 total drops) tested per procedure specified in EN 60601-1
Mechanical Shock	The transmitter withstand 3 repeated mechanical shocks in each direction of each axis of 15g with a 11 msec pulse duration and a half sine wave velocity change profile and tested per procedure IEC 60068- 2-27. Per IEC 60068-2-27 Section A.7, the additional 2 axis should be tested using the same protocol on additional devices if possible. The transmitter withstand 3 repeated mechanical shocks in each direction of each axis of 30g with a 6 msec pulse duration and a half sine wave velocity change profile and tested per procedure IEC 60068-2-27. Per IEC 60068-2-27 Section A.7, the additional 2 axis should be tested using the same protocol on additional devices if possible.
Mechanical Vibration	Wave Form Random Acceleration by Frequency Range 10 Hz - 100 Hz 1.0 (m/s²)²/Hz 100 Hz - 200 Hz -3dB/octave 200 Hz - 2000 Hz 0.5 (m/s²)²/Hz Duration 30 minutes for each axis
Push	The transmitter shall withstand a steady force of 250N (25.5Kg) over a 30 mm diameter contact surface for a period of 5 seconds per EN 60601-1 Section 15.3.2.
Molding Stress Relief	The transmitter exterior housing shall be so constructed that any shrinkage or distortion of the material due to release of internal stresses caused by the molding or forming operation does not result in an unacceptable risk when the housing is exposed to 70°C for 7 hours as defined in EN 60601-1 , Section 15.3.6.
Impact Test	The GST3 shall comply with the impact requirements per EN 60601-1 Section 15.3.3.
Electrical	The transmitter shall comply with the safety requirements for type BF electrical equipment as specified in EN 60601-1. The transmitter sensor connections shall survive a direct connection to the device's internal battery voltage, from a fully charged battery referenced to circuit ground, without degradation or damage to the circuitry.
Battery	Low battery warning shall occur at least 1 day before the occurrence of a depleted battery warning under normal use


Charger Test Tool


Test Description	Specification per Standard	Specification per ES9500 (Section)	Results
Storage Temperature	IEC 60601-1-11, Section 4.2.1	5.4.1	Pass
Storage Humidity	IEC 60601-1-11, Section 4.2.1	5.4.1	Pass
Operating Temp Shock	IEC 60601-1-11, Section 4.2.2	5.4.2	Pass
Operating Atmospheric Pressure	IEC 60601-1-11, Section 4.2.2	5.4.2	Pass
Operating Humidity	IEC 60601-1-11, Section 4.2.2	5.4.2	Pass
Mechanical Vibration	IEC 60601-1-11, Section 10.1.2b	5.6.2	Pass
Mechanical Drop	IEC 60601-1, Section 15.3.4.1	5.6.1	Pass
Mechanical Shock	IEC 60601-1-11, Section 10.1.2a	N/A	Pass
Connector Retention Force	N/A	5.3.7	Pass
Connector Cycling	N/A	5.3.5	Pass
Push Test	IEC 60601-1, Section 15.3.2	N/A	Pass
Molding Stress Relief	IEC 60601-1, Section 15.3.6	N/A	Pass


Manufacturing scheme or process description brief

The Guardian 2 Link Transmitter and accessories are manufactured in accordance with Good Manufacturing Practices as described in 21 CFR 820 of the Federal Code, and under guidelines laid out in ISO 13485. All changes to product and manufacturing procedures are performed in accordance with Medtronic MiniMed practices and qualified prior to incorporation into production.

Guardian 2 Link Transmitter Manufacturing Overview

Requirements for the protection of the environment in the application of medical devices. Procedure for utilization/liquidation of MD	Transmitter The device must be discarded according to local regulations for battery disposal (non-incineration). Do not discard the transmitter in a medical waste container or otherwise subject it to incineration. The transmitter contains a battery that may explode upon incineration. Discard the transmitter according to local regulations for battery disposal. Charger Dispose the charger according to the local regulations for battery disposal (non-incineration). Tester If the tester comes in contact with blood, the tester must be discarded. Dispose of the tester according to the local regulations for medical waste disposal.
Methods and apparatus for disinfection and sterilization cleaning (if applicable)	Cleaning the transmitter The transmitter is a single-patient use device and not intended for multi-patient use. Caution: Do not discard the transmitter in a medical waste container or otherwise subject it to incineration. The transmitter contains a battery that may explode upon incineration. <i>Note: The tester is a required component for cleaning the transmitter.</i> Always clean the transmitter after each use. To clean the transmitter, you will need the following materials: mild liquid soap, a softbristled toddler toothbrush, a container, 70% isopropyl alcohol, and a few clean, dry cloths. Warning: Cracking, flaking, or damage of the housing are signs of deterioration and the performance of the device may be compromised. This may affect the ability to properly clean and disinfect the transmitter. If these signs are noted, stop using the device and call the 24 Hour HelpLine or your local representative. The device must be discarded according to local regulations for battery disposal (non-incineration). To clean the transmitter: 1 Attach the tester to the transmitter. 2 If optional occlusive dressing or overtape was used and there is adhesive residue on the transmitter, see “Removing adhesive residue” below. 3 Rinse the transmitter under room temperature tap water for at least one minute, and until visibly clean. Make sure all hard-to-reach areas are rinsed completely. 4 Prepare a mild liquid soap solution using 5 milliliters (1 teaspoon) of liquid soap per 3.8 liters (1 gallon) of room temperature tap water. 5 With tester still attached, submerge the transmitter in the cleaning solution and soak for one minute. 6 Holding the tester, brush the entire surface of the transmitter using a soft-bristled toddler toothbrush. Make sure to brush all hard-to-reach areas until visibly clean.

- 7 Rinse the transmitter under running room temperature tap water for at least one minute, and until all visible liquid soap is gone.
- 8 Dry the transmitter and tester with a clean, dry cloth.
- 9 Place the transmitter and tester on a clean, dry cloth and air dry them completely.
- 10 Disconnect the tester from the transmitter.

Removing adhesive residue

You may need to perform this procedure only if you have used optional occlusive dressing, which may leave adhesive residue on the transmitter. If you visually inspect the transmitter and see adhesive residue on it, follow the instructions below.

To remove adhesive residue, you will need the following materials: Detachol medical adhesive remover and cotton swabs.

To remove adhesive residue:

1. Make sure the tester is attached to the transmitter.
2. Holding the tester, saturate a cotton swab in the Detachol solution and gently rub the adhesive residue on the transmitter until it is fully removed.
3. Continue with the cleaning procedure.

Cleaning the charger

This procedure is for general cleaning as required, based on physical appearance.

Caution: The charger is NOT waterproof. Do NOT immerse in water or any other cleaning agent.

Caution: Dispose the charger according to the local regulations for battery disposal (non-incineration).

To clean the charger:

1. Wash your hands thoroughly.
2. Use a damp cloth with mild cleaning solution, such as a dishwashing detergent, to clean any dirt or foreign material from the outside of the charger. Never use organic solvents, such as paint thinner or acetone, to clean the charger.
3. Place the charger on a clean, dry cloth and air dry for 2–3 minutes.

Cleaning the senter

To clean the senter, you will need the following materials: mild liquid soap, a soft-bristled toddler toothbrush, and a container.

	To clean the serter: 1. Prepare a mild liquid soap solution using 5 milliliters (1 teaspoon) of mild liquid soap per 3.8 liters (1 gallon, which is equal to 16 cups) of tap water. 2. Rinse the serter under running tap water at room temperature for at least one minute. Continue rinsing until the serter is visibly clean. Make sure all hard-to-reach areas are rinsed completely. 3. When rinsing the serter, press and release the mechanism on the underside of the serter to ensure that the entire device is completely rinsed with water. 4. Submerge the serter in the mild liquid soap solution and soak it for at least 10 minutes. 5. Brush the entire surface of the serter using a soft-bristled toddler toothbrush until visibly clean. Pay close attention to brushing the hard-to-reach areas. 6. Rinse the serter under running tap water at room temperature for at least one minute, and until all visible liquid soap is gone. 7. Shake off any excess water or moisture and place the serter upright on a clean, dry, cloth to air dry.
Methods and conditions of sterilization time maintaining sterility (for sterile medical devices produced)	Not applicable. The Guardian 2 Link Transmitter and accessory devices are not provided as sterile, nor are intended for sterilization.
Lifetime	Operation life: The transmitter expected service life is 1 year depending on patient usage. The Charger is intended to operate properly for one year excluding user caused damage. The Serter has a minimum service life of one year. The Tester is intended to operate properly for 122 uses excluding user caused damage.
Operating conditions:	Transmitter: Temperature: -5°C to 45°C (23°F to 113°F) Caution: When operating the transmitter on a tester in air temperatures greater than 41°C (106°F), the temperature of the transmitter may exceed 43°C (109°F) Relative humidity: 5% to 95% with no condensation Pressure: 61.36 to 106.17 kPa (8.9 to 15.4 psi) Charger: Temperature: 10 °C to 40 °C

	Relative humidity: 5% to 75% with no condensation		
	Serter: The Serter shall operate within a normal room temperature of 10 °C to 40 °C.		
	Tester: Operating range, temperature: -5°C to +45°C Operating range, humidity: 5% to 95% RH		
Assembly and installation requirements (if applicable);	Not applicable.		
Information about maintenance (if applicable);	Troubleshooting The following table contains troubleshooting information for the transmitter, charger, and tester. For more information about troubleshooting, see your pump user guide.		
	Problem	Likely Cause(s)	Resolution
	You connected the transmitter to the charger and no lights came on.	The transmitter connector pins are damaged or corroded. Your charger battery has no power.	1. Check the transmitter connector pins for damage or moisture. If the pins are damaged or corroded, contact the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter. 2. If there is no damage to the connector pins, replace the battery in the charger.
	During charging, the flashing green light on the charger turns off and you see a flashing red light on the charger.	Your charger battery is low on power.	Replace the battery in the charger.
	During charging, the flashing green light on the charger turns off and you see a series of quick flashing red lights on the charger.	Your transmitter is low on power.	1. Charge the transmitter continuously for one hour. If flashing does not stop, proceed to step 2. 2. Charge the transmitter continuously for eight hours. If flashing does not stop, call the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.
	During charging, a mix of quick	Your charger and your transmitter	1. Replace the battery in the charger.

	and long flashing red lights appear on the charger.	are low on power.	2. Charge the transmitter continuously for one hour. If flashing does not stop, proceed to step 3. 3. Charge the transmitter continuously for eight hours. If flashing does not stop, call the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.
	The green light on the transmitter does not flash when you connect it to the sensor.	Your transmitter is not fully connected. Your transmitter is low on power. Your sensor is not properly inserted into your body.	1. Disconnect the transmitter from the sensor. 2. Wait for five seconds and reconnect them. If the green light still does not flash, proceed to step 3. 3. Fully charge the transmitter. If the green light still does not flash, proceed to step 4. 4. The sensor may not be properly inserted into your body. Insert a new sensor.
	The green light on the transmitter does not flash when you connect it to the tester.	Your transmitter is low on power. Your transmitter is not fully connected.	1. Check the connection between the transmitter and the tester. If the green light still does not flash, proceed to step 2. 2. Fully charge the transmitter. 3. Test the transmitter with the tester again. If you still do not see the green light flash, call the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.
	Your transmitter battery does not last for six days.	Your transmitter is not fully charged when you connect it to the sensor. The transmitter and pump frequently lose wireless connection.	1. Fully charge the transmitter before connecting it to the sensor. If the transmitter battery still does not last for six days, proceed to step 2. 2. Move away from any device that can cause RF interference. For more information on RF interference, see Radio Frequency (RF) communication, on page 2. 3. Make sure your pump and your transmitter are located on the same side of

			your body to minimize any RF interference. If your fully charged transmitter battery continues to lose power before a full six days, call the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.
	Your transmitter has lost connection with your pump.	Your pump is out of range. There is RF interference from other devices.	1. Move away from any device that can cause RF interference. For more information on RF interference, see Radio Frequency (RF) communication, on page 2. If your transmitter is still not communicating with your pump, proceed to step 2. 2. Make sure your pump and your transmitter are located on the same side of your body to minimize any RF interference. If your transmitter is still not communicating with your pump, call the 24 Hour HelpLine or your local representative for assistance.
Conditions of transportation and storage	Storing the devices Store the transmitter, charger, and tester in a clean, dry location at room temperature. If the transmitter is not in use, you must charge the transmitter at least once every 60 days. Although not required, you may store the transmitter on the charger. If you are storing the transmitter on the charger, you must disconnect and reconnect the charger and the transmitter at least once every 60 days. Transportation the devices Transmitter: The packaged GST3 shall comply with the requirements of ASTM D4169-09 using hot, humid and extreme heat, moderate RH pre-conditioning. Charger: The charger device shall not be damaged or experience loss of function due to shipping, per ASTM D4169. Tester: ASTM D4169 The device is to be returned to within the operational temperature range before acceptance testing. Serter: After Thermal Conditioning, the Serter as packaged and packed in its shipping box shall withstand mechanical testing per ASTM D4169. Conditions of transportation and storage (Temperature, relative humidity, atmospheric pressure): Transmitter: Temperature: -25°C to 55°C (-13°F to 131°F)		

Relative humidity: 10% to 100% with condensation
Pressure: 61.36 to 106.17 kPa (8.9 to 15.4 psi)

Charger:

Temperature: -10 °C to +50 °C
Relative humidity: 10% to 95% with no condensation

Serter:

Temperature: -20 °C +0/-2 °C to 55 °C +2/-0 °C
Atmospheric Pressure: 49.64 kPa +0/-3.5 kPa to 106.18 kPa +3.5/-0 kPa (7.2 psia +0/-0.5 psia to 15.4 psia +0.5/-0 psia).
Relative Humidity: 15% to 85% condensing

Tester:

Storage temperature: -30 °C to +70 °C
Storage humidity: 10% to 100% RH with condensation

Warranty

Medtronic MiniMed, Inc. (or such other legal entity as may be referred to as manufacturer on the labeling of this device “Medtronic MiniMed”) warrants the Medtronic transmitter to the purchaser of the product against defects in material and workmanship for a period of one (1) year and the charger for up to one (1) year from the date of purchase.

During the warranty period, Medtronic MiniMed will replace or repair, at its discretion, any defective transmitter or charger, subject to the conditions and exclusions stated herein. This warranty applies only to new devices. In the event a transmitter or charger is replaced, the warranty period will not be extended past its original expiration date.

This warranty is valid only if the Medtronic transmitter or charger is used in accordance with the manufacturer’s instructions. Without limitation, this warranty will not apply:

- If damage results from changes or modifications made to the transmitter or charger by the user, or third persons, after the date of purchase;
- If damage results from service or repairs performed by any person or entity other than the manufacturer;
- If damage results from a *Force Majeure* or other event beyond the control of the manufacturer;
- If damage results from negligence or improper use, including but not limited to: improper storage, submersion in water, physical abuse (such as dropping);
- If damage results from use of the device in a manner other than according to the manufacturer’s product labeling, instructions for use, or regulatory notifications.

This warranty shall be personal to the original purchaser. Any sale, rental or other transfer or use of the product covered by this warranty to or by a user other than the original purchaser shall cause this warranty to immediately terminate. This warranty does not apply to Glucose Sensors and other accessories.

The remedies provided for in this warranty are the exclusive remedies available for any breach hereof. Neither Medtronic MiniMed nor its suppliers or distributors shall be liable for any incidental, consequential, or special damage of any nature or kind caused by or arising out of a defect in the product.

All other conditions and warranties, other than mandatory statutory warranties, expressed or implied, are excluded, including the warranties of merchantability and fitness for a particular purpose.

This warranty gives the purchaser specific legal rights, and the purchaser may also have other rights that vary under local law. This warranty does not affect the purchaser's statutory rights.

The list of national and international regulatory documents/standards, which matches medical device

Transmitter

Standards ID	Reference Document	Title (equivalent standard/comments)
EN 60601-1	BSI BS EN 60601-1 +A11: 2006.11.30	Medical Electrical Equipment - Part1: General Requirement for Basic Safety and Essential Performance
EN 60601-1-2	BSI BS EN 60601-1-2: 2007.09.28/C:2010	Medical Electrical Equipment – Part 1-2: General Requirements for Basic Safety and Essential Performance – Collateral Standard: Electromagnetic Compatibility – Requirements and Tests
EN 60601-1-6	BSI BS EN 60601-1-6: 2010.05.31	Medical electrical equipment. General requirements for basic safety and essential performance. Collateral standard. Usability
EN 60601-1-11	BSI BS EN 60601-1-11: 2010.07.31	Medical electrical equipment -- Part 1-11: General requirements for basic safety and essential performance -- Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
EN 62304	BSI BS EN 62304: 2006.11.30 /C:2011	Medical device software –Software life cycle processes
EN 62366	BSI BS EN 62366: 2008.04.30	Medical Devices – Application of Usability Engineering to Medical Devices
EN 980	BSI BS EN 980: 2008.06.30	Symbols for use in the Labeling of medical devices
EN 1041	BSI BS EN 1041: 2008.08.31	Information Supplied By The Manufacturers Of Medical Devices
EN ISO 10993-1	BSI BS EN ISO 10993-1: 2010.06.30/2C:2010	Biological evaluation of medical devices – Part 1: Evaluation and testing. (ISO 10993-1:2009) CORR: July 31, 2010; CORR: November 30, 2010

EN ISO 10993-5	BSI BS EN ISO 10993-5: 2009.07.31	Biological evaluation of medical devices Part 5: Tests for in vitro cytotoxicity
EN ISO 10993-12	BSI BS EN ISO 10993-12: 2012.10.31	Biological evaluation of medical devices- Part 12: Sample preparation and reference materials
EN ISO 10993-17	BSI BS EN ISO 10993-17: 2009.06.30	Biological evaluation of medical devices- Part 17: Establishment of allowable limits for leachable substances
EN ISO 10993-18	BSI BS EN ISO 10993-18: 2009.06.30	Biological evaluation of medical devices- Part 18: Chemical characterization of materials
EN ISO 13485	BSI BS EN ISO 13485: 2012.03.31	Medical Devices – Quality Management Systems Requirements for Regulatory Purposes
EN ISO 14971	BSI BS EN ISO 14971: 2012.07.31	Medical Devices – Application of Risk Management to Medical Devices

Serter:

Standards ID	Reference Document	Title (equivalent standard/comments)
EN 1041	BSI BS EN 1041:2008.08.31	Information supplied by the manufacturer with medical devices.
EN ISO 10993-1	BSI BS EN ISO 10993-1: 2009.07.31/R:2010	Biological evaluation of medical devices. Part 1: Evaluation and testing within a risk management process CORRIGENDUM July 31, 2010
EN ISO 10993-5	BSI BS EN ISO 10993-5: 2009.07.31	Biological evaluation of medical devices, Part 5: Tests for in vitro cytotoxicity
EN ISO 10993-10	BSI BS EN ISO 10993-10: 2014.2.28	Biological evaluation of medical devices Part 10: Tests for irritation and skin sensitization (ISO 10993-10:2010)
EN ISO 10993-12	BSI BS EN ISO 10993-12: 2012.10.31	Biological evaluation of medical devices Part 12: Sample preparation and reference materials
EN ISO 10993-17	BSI BS EN ISO 10993-17: 2009.06.30	Biological evaluation of medical devices Part 17: Establishment of allowable limits for leachable substances
EN ISO 10993-18	BSI BS EN ISO 10993-18: 2009.06.30	Biological evaluation of medical devices Part 18: Chemical characterization of materials
EN ISO 13485	BSI BS EN ISO 13485: 2012.03.31/C:2012	Medical Devices – Quality Management systems – Requirements for Regulatory Purposes
EN ISO 14971	BSI BS EN ISO 14971: 2012.07.31	Medical devices – Application of Risk Management to Medical Devices
EN ISO 15223-1	BSI BS EN ISO 15223-1: 2012.07.31	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements - CORR: July 31, 2012
EN 62366	BSI BS EN 62366:2008.04.30	Medical Devices - Application of usability engineering to medical devices.

Tester:		
Standards ID	Reference Document	Title (equivalent standard/comments)
EN 60601-1	BSI BS EN 60601-1:2006.11.30/C:2011	Medical Electrical Equipment - Part 1: General Requirement for Basic Safety and Essential Performance
EN 60601-1-11	BSI BS EN 60601-1-11:2010.07.31	Medical electrical equipment -- Part 1-11: General requirements for basic safety and essential performance -- Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
EN 60601-1-6	BSI BS EN 60601-1-6:2010.01.01	Medical Electrical Equipment -- Part 1-6: General Requirements for Safety -- Collateral Standard:
EN 980	BSI BS EN 980: 2008.06.30	Symbols for use in the Labeling of medical devices
EN 1041	BSI BS EN 1041: 2008.08.31	Information Supplied By The Manufacturers Of Medical Devices
EN ISO 13485	BSI BS EN ISO 13485:2012.03.31/C:2012	Medical Devices -- Quality Management Systems Requirements for Regulatory Purposes
EN ISO 14971	BSI BS EN ISO 14971:2012.07.31	Medical Devices -- Application of Risk Management to Medical Devices
EN 62366	BSI BS EN 62366: 2008.04.30	BSI BS EN 62366: 2008.04.30
Charger:		
Standards ID	Reference Document	Title (equivalent standard/comments)
EN 980	BSI BS EN 980: 2008.06.30	Symbols for use in the Labeling of medical devices
EN 1041	BSI BS EN 1041: 2008.08.31	Information Supplied By The Manufacturers Of Medical Devices
EN ISO 13485	BSI BS EN ISO 13485: 2012.03.31/C:2012	Medical Devices -- Quality Management Systems- Requirements for Regulatory Purposes
EN ISO 14971	BSI BS EN ISO 14971: 2012.07.31	Medical Devices -- Application of Risk Management to Medical Devices
EN 60601-1	BSI BS EN 60601-1: 2006.11.30/C:2011	Medical Electrical Equipment - Part 1: General Requirement for Basic Safety and Essential Performance
EN 60601-1-2:2007	BSI BS EN 60601-1-2: 2007.09.28/C:2010BSI BS EN 60601-1-2: 2007	Medical Electrical Equipment -- Part 1-2: General Requirements for Basic Safety and Essential Performance -- Collateral Standard: Electromagnetic Compatibility -- Requirements and Tests
EN 60601-1-63	BSI BS EN 60601-1-6: 2010.05.31	Medical Electrical Equipment -- Part 1-6: General Requirements for Basic Safety and Essential Performance -- Collateral standard: Usability

	EN 60601-1-11	BSI BS EN 60601-1-11: 2010.07.31	Medical electrical equipment -- Part 1-11: General requirements for basic safety and essential performance -- Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
	EN 62366	BSI BS EN 62366: 2008.04.30	Medical Devices – Application of Usability Engineering to Medical Devices
Customer Service	Contact for all questions, including complaints, repair and service of medical devices the Authorized Representative in Russia: Medtronic, LLC, 123317, Moscow, Presnenskaya Naberezhnaya., 10 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru		

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Prepare Battery

Verify that welds are uniform.

Perform peel strength: Black wire > 1.9 lbs., Red wire > 1.7 lbs.

Prepare, anneal, laser mark bottom housing

Visual inspection; Marking should be legible and housing free from damage

Antenna form rig

Visual inspection : antenna should be free of foreign material and scratches. Formed in proper location

Connector Build

Visual inspection for part alignment, parts should be free of foreign material

PCBA Build

Perform shear test and visual inspection

Assemble

Anneal top housing

Stake antenna to top housing

visual inspection for part alignment, parts should be free of foreign material

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
Manufacturing Steps

Inspection & Testing

Solder Connector to PCBA, apply epoxy to leads & seams

Visual inspection to confirm coverage of epoxy. Contacts should be epoxy covered

Place PCBA/Connector assembly into top housing

Visual inspection to confirm proper placement of tape and O-ring

Solder antenna & battery to PCBA

Visual inspection to confirm solder joint and flashing LEO

Apply UV silicone fill to PCBA, battery leads, tabs

Verify that battery leads, exposed pads, solder joints and tabs are completely covered

Glue bottom housing to top housing assembly

Verify that gluing seam is filled and that there are no voids in seam

Polish /clean housing

Perform automated leak test.

Acceptance criteria: Leak rate <0.34 mbar

Apply acrylic coating to bottom housing

Visual inspection to confirm coating coverage. Coating should be free of foreign materials

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
Manufacturing Site

Inspection & Testing

SW Programming and Final Quality Assurance

Subassembly received

Visual inspection' Print should be legible and housing free from damage

Charge device

Test Battery Capacity- Pass/Fail

Perform RFTest

Perform automated radiofrequency and serial number testing:

- Visually read and program serial number and model number from the device rear label.
- Confirm transmitter power, frequency and Error Vector Magnitude (EVM).
- Confirm TelD channels and perform an automated RF Link test on at least one of the five TelD channels

Perform Leak Test

Perform automated leak test: Leak rate < 0.40 mbar

Perform Accuracy Test

Perform automated calibration accuracy test and program applications software.

Acceptable ranges for applied current

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ATTESTATION STATEMENT

As a responsible member of Medtronic MiniMed, I hereby attest that I believe the One-press Sertter User Guide, issued by the Medtronic MiniMed with 6026011-12B1_a attached herein, to be a true and accurate copy.

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
Introduction

The one-pressserter (MMT-7512) for Enlite™ glucose sensors (MMT-7008A, MMT-7008B) is used as an aid for inserting the Enlite sensor.

Potential risks related to sserter use

General risks with sserter use may include skin infection around the area where the sserter is used.

Indications for use

The sserter is used as an aid for inserting the Enlite sensor. It is indicated for single-patient use.

Contraindications

Do not use the sserter on products other than the Enlite sensor. Medtronic cannot guarantee the safety or efficacy of this product if used with other products.

General warnings

Read this entire user guide before attempting to insert the sensor. The sserter does not work the same as other Medtronic insertion devices. Failure to follow directions may result in improper insertion, pain, or injury.

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A retractable needle is attached to the sensor and minimal blood splatter may occur. If you are a healthcare professional or caregiver, wrap sterile gauze around the sensor to minimize contact with blood. Keep as much distance as possible between you and the patient when removing the needle.

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Always inspect the sensor packaging for damage before use. Sensors are sterile and non-pyrogenic, unless the package has been opened or damaged. Do not use the sensor if the sterile package has been opened or damaged. Use of an unsterile sensor can cause site infection.

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1. Apply steady pressure using sterile gauze or a clean cloth, placed on top of the sensor for up to three minutes. The use of unsterile gauze can cause site infection.
2. If bleeding stops, connect the transmitter to the sensor.
If bleeding does not stop, do not connect the transmitter (or recorder) to the sensor because blood can get into the transmitter connector, and could damage the device.

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Make sure the sensor is securely placed in the sserter to avoid improper insertion, pain, or minor injury.

Keep the sserter away from children. This product contains small parts and may pose a choking hazard.

Refer to the *enlite user guide* for additional sensor related warnings.

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Choose an insertion site that has an adequate amount of subcutaneous fat. The Enlite sensor has been approved for use on the abdomen, the lower back, or the upper buttock areas. For a diagram of the area, see *Inserting the sensor, step 2*.

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Removing the sensor

Always refer to the *enlite user guide* for instructions on how to remove the sensor.

Storage and handling

Store the serter in the released position, at room temperature, to maintain optimum product performance and life.

Specifications

Approximate Dimensions	Approximate Weight
7.85 x 6.91 x 5.77 centimeters (3.09 x 2.72 x 2.27 inches)	112.5 grams (3.97 ounces)

Serter life of use

The serter has a minimum service life of one year.

Serter maintenance

The serter is intended for single-patient use. The serter needs to be cleaned after every use or whenever there is debris, blood, or other contaminants.

Cleaning


To clean the serter, you will need the following materials: mild liquid soap, a soft-bristled toddler toothbrush, and a container.

To clean the serter:


1. Prepare a mild liquid soap solution using 5 milliliters (1 teaspoon) of mild liquid soap per 3.8 liters (1 gallon, which is equal to 16 cups) of tap water.
2. Rinse the serter under running tap water at room temperature for at least one minute. Continue rinsing until the serter is visibly clean. Make sure all hard-to-reach areas are rinsed completely.
3. When rinsing the serter, press and release the mechanism on the underside of the serter to ensure that the entire device is completely rinsed with water.
4. Submerge the serter in the mild liquid soap solution and soak it for at least 10 minutes.
5. Brush the entire surface of the serter using a soft-bristled toddler toothbrush until visibly clean. Pay close attention to brushing the hard-to-reach areas.
6. Rinse the serter under running tap water at room temperature for at least one minute, and until all visible liquid soap is gone.
7. Shake off any excess water or moisture and place the serter upright on a clean, dry, cloth to air dry.

©2015, Medtronic MiniMed, Inc. All rights reserved. Enlife™ is a trademark of Medtronic MiniMed, Inc.

Components


One-pressserter
A. Bump on both buttons
B. Thumbprint marking

C. Pedestal
D. Needle housing
E. Sensor
F. Clear liner

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

1. Wash your hands.


2. Choose an insertion site that has an adequate amount of fat, either on the abdomen, the lower back, or the upper buttock areas. Shaded areas show the best sites for sensor insertion.


3. Clean the insertion site with alcohol. Let the area air dry.


4. Open the sensor package.


5. Remove the pedestal from the package. Place the pedestal on a flat surface.


6. Make sure that the Enlite sensor's adhesive tab is tucked under the sensor connector and sensor snaps.


7. Holding serter Correctly

Place your thumb on the thumbprint marking to hold the serter without touching the buttons.


Holding serter Incorrectly


Your fingers should not be touching the buttons.


8a–8b. Grip the serter, placing your thumb on the thumbprint marking, **without holding the buttons**. Carefully push the serter down onto the pedestal until the base of the serter sits flat on the table and you hear a click.

9a. To detach the serter from the pedestal, using one hand, place your thumb on the thumbprint marking to grip the serter **without holding the buttons**.


9b. With your other hand, place two fingers on the pedestal arms and slowly pull the serter straight up. Do not detach the pedestal from the serter in midair, as this might damage the sensor.

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10a. Hold theserter steady against your cleaned insertion site, without pushing theserter too deep into your skin.


10b-10c. Press and release the bump on both buttons at the same time, while holding theserter flat against your body.

10d. Continue holding theserter flat against your body to allow the adhesive to stick to your skin.

10e. Slowly lift theserter away from your body, making sure that the buttons are not pressed.

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
Patient:

11a. Gently hold the base of the sensor against skin. Hold the needle housing at the top and slowly pull straight out, away from the sensor.

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11b. Wrap sterile gauze around the sensor (as shown in image 11b). Gently hold the base of the sensor against skin. Hold the needle housing at the top and slowly pull straight out, away from the sensor.

WARNING: Watch for bleeding at the insertion site. If bleeding occurs under, around, or on top of the sensor, apply steady pressure using sterile gauze or a clean cloth placed on top of the sensor for up to three minutes. The use of unsterile gauze can cause an infection. If bleeding does not stop, remove the sensor and apply steady pressure until the bleeding stops.

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



Applying the overtape


This image shows an example of the Enlite overtape applied correctly. To apply the overtape, go to step 1.


Note: The sensor body and the skin must be covered partially as shown in the image. The sensor connector and sensor snaps must be uncovered and accessible.

1. Remove the adhesive liner marked "1" from the overtape. Do not remove the adhesive liner marked "2" from the sides of the overtape.

2. Apply the overtape to the rounded part of the sensor and the skin, making sure that the sensor connector and sensor snaps are sticking out through the hole in the overtape.


3. Stretch the overtape around the sensor connector and sensor snaps so that it lies flat against the sensor adhesive pad.


4. Remove the adhesive liner marked "2" from both sides of the overtape and smooth the overtape to the skin.


5. For details on how to enter your sensor settings into your pump or monitor, refer to the corresponding user guide.


6. Connect the transmitter to the sensor.


7. Wait for the green light on the transmitter to flash. If the green light on the transmitter does not flash, refer to the Troubleshooting section of your transmitter user guide.

8a. Remove the adhesive liner from the adhesive tab.


8b-8c. Cover the transmitter with the adhesive tab. Do not pull the tab too tightly.

9. To apply a second overtape, complete Option 9a or Option 9b.


If you have skin irritation due to moisture buildup, complete Option 9a. If the transmitter catches on your clothes, complete Option 9b.

Note: Use the Enlite overtape included in the sensor box to secure the transmitter. If you run out of Enlite overtape, you may use over the counter medical tape.

Option 9a

Position the overtape so that the end of the transmitter is exposed through the opening in the tape.

Option 9b

Position the opening of the overtape over the center of the transmitter. Make sure that the end of the transmitter is securely covered by the overtape.

OR

Icon table

	Consult instructions for use.		Caution
	Batch code		Humidity limitation
	One per container/package		Date of manufacture (YYYY-MM-DD)
	Temperature limit		Authorized representative in the European Community.
	Catalogue or model number		
	Manufacturer		

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ATTESTATION STATEMENT

As a responsible member of Medtronic MiniMed, I hereby attest that I believe the Guardian 2 Link Transmitter guide, issued by the Medtronic MiniMed with 6025871-11C1_a attached herein, to be a true and accurate copy.

Masae Tanamura
Senior Regulatory Specialist

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State of California


County of LOS ANGELES

On FEBRUARY 2, 2017 before me, NAVIN K. GAUR
(date) (insert name and title of the officer)

personally appeared MASAE TANEMURA, who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument. I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature Navin K. Gaur (Seal)

Medtronic MiniMed
18000 Devonshire Street
Northridge, CA 91325, USA

Contacts:

Africa: Medtronic Africa (Pty) Ltd.
Tel: +27 (0) 11 677 4800

Argentina: Corpomedica S.A.
Tel: +(11) 4 814 1333
Medtronic Directo 24/7: +0800 333 0752

Australia: Medtronic Australasia Pty. Ltd.
Tel: 1800 668 670 (product orders)
Tel: 1800 777 808 (customer help)

Azerbaijan: Isomed
Tel: +994 (12) 464 11 30

Bangladesh: Sonargaon Healthcare Pvt Ltd.
Mobile: (+91)-9903995417
or (+880)-1714217131

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Malaysia: Medtronic International Ltd.
Tel: +603 7946 9000

Middle East and North Africa: Regional Office
Tel: +961-1-370 670

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Sri Lanka: Swiss Biogenics Ltd.
Mobile: (+91)-9003077499
or (+94)-777256760

Suomi: Medtronic Finland Oy
Tel: +358 20 7281 200
Help line: +358 400 100 313

Sverige: Medtronic AB
Tel: +46 8 568 585 20
Fax: +46 8 568 585 11

Taiwan: Medtronic (Taiwan) Ltd.
Tel: 02-21836000
Toll free: +886-800-005285

Thailand: Medtronic (Thailand) Ltd.
Tel: +662 232 7400

Türkiye: Medtronic Medikal Teknoloji
Ticaret Ltd. Sirketi.
Tel: +90 216 4694330

USA: Medtronic Diabetes Global Headquarters
24 Hour HelpLine: +1-800-646-4633
To order supplies: +1-800-843-6687

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The Guardian™ 2 Link transmitter is a component of the continuous glucose monitoring system for the MiniMed™ 640G insulin pump. The transmitter collects data from the Enlite™ glucose sensor. The transmitter then wirelessly sends the data to the insulin pump.

Guardian 2 Link transmitter kit components

A complete Guardian 2 Link transmitter kit includes the following components:

- Guardian 2 Link transmitter (MMT-7731)
- Watertight Tester (MMT-7726)
- Charger (MMT-7715)
- One-press setter (MMT-7512)

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Contraindications

Do not expose your transmitter to MRI equipment, diathermy devices, or other devices that generate strong magnetic fields. If your transmitter is inadvertently exposed to a strong magnetic field, discontinue use and contact the 24 Hour HelpLine or your local representative for further assistance.

Warnings

This product contains small parts and may pose a choking hazard for young children. If the tester comes in contact with blood, the tester must be discarded. Dispose of the tester according to the local regulations for medical waste disposal.

Bleeding may occur after inserting the sensor. Make sure that the site is not bleeding before connecting the transmitter to the sensor. If bleeding occurs, apply steady pressure with a sterile gauze or clean cloth at the insertion site until bleeding stops. After bleeding stops, connect the transmitter to the sensor.

Contact the 24 Hour HelpLine or your local representative if you experience any adverse reactions associated with the transmitter or sensor.

Magnetic fields

Do not expose your transmitter to MRI equipment, diathermy devices, or other devices that generate strong magnetic fields. If your transmitter is inadvertently exposed to a strong magnetic field, discontinue use and contact the 24 Hour HelpLine or your local representative for further assistance.

X-rays, MRIs, diathermy devices, and CT scans

If you are going to have an x-ray, diathermy treatment, CT scan, MRI or other type of exposure to radiation, remove your sensor and transmitter before entering a room containing any of these equipment.

Important information about airport security systems, and using your transmitter on an airplane, can be found on the Emergency Card. Be sure to carry the Emergency Card provided with your device when you are traveling.

Precautions

Refer to the Enlite Sensor User Guide for all precautions, warnings, and instructions relating to the sensor.

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Do not allow water, or any other liquid, to come in contact with the tester when it is not connected to the transmitter. A wet tester can cause damage to the transmitter.

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- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and the receiver

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- Decrease the distance between the transmitter and the insulin pump to 1.8 meters (6 feet) or less.
- Increase the separation between the transmitter and the device that is receiving/ emitting interference.

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For Canada only

This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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Directive 1999/5/EC

Medtronic declares that this product is in conformity with the essential requirements of Directive 1999/5/EC on Radio and Telecommunications Terminal Equipment.

For additional information, contact Medtronic MiniMed at the address or phone number provided on the back cover

Assistance

Please contact your local representative for assistance. Refer to the Medtronic Diabetes International Contacts list at the beginning of this user guide for contact information.

Charger

The transmitter contains a non-replaceable, rechargeable battery that you can recharge as needed with the charger. The charger has a green light that shows the charging status and a red light that communicates any problems during charging. If you see a red light, see the Troubleshooting section. The charger needs one AAA alkaline battery.


Note: If the battery is installed incorrectly or is low, the charger will not work. Repeat the battery installation steps using a new battery.

Installing a battery in the charger

To install a battery in the charger:

- 1 Push the battery cover in and slide it off (as shown in the following illustration).
- 2 Insert a new alkaline AAA battery. Make sure the + and - symbols on the battery align with these same symbols shown on the charger.

- 3 Slide the cover back on the charger until it clicks into place.



Charging the transmitter

Caution: Charge the transmitter after each sensor use. A fully charged transmitter works at least six days without recharging. A depleted transmitter can take up to one hour to recharge.

To charge the transmitter:

- 1 Connect the transmitter to the charger by lining it up, flat side down, with the charger. Push the two components together fully.
- 2 Within 10 seconds after the transmitter is connected, a green light on the charger will flash for one to two seconds as the charger powers on. For the rest of the charging time, the charger's green light will continue to flash in a pattern of four flashes with a pause between the four flashes.
- 3 When charging is complete, the green light on the charger will stay on, without flashing, for 15 to 20 seconds and then turn off.
- 4 After the green charger light turns off, disconnect the transmitter from the charger. The green light on the transmitter will flash for about five seconds and then turn off.

Inserting the sensor


Always refer to the senter user guide for instructions on how to insert the sensor.

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- 1 After the sensor is inserted, consult your senter user guide for details on applying the required overtape.
- 2 Hold the rounded end of the inserted sensor to prevent it from moving during connection.
- 3 Hold the transmitter as shown. Line up the two notches on the transmitter with the side arms of the sensor. The flat side of the transmitter should face the skin.
- 4 Slide the transmitter onto the sensor until the sensor's flexible arms snap into the notches on the transmitter. If the transmitter is properly connected, and if the sensor has had enough time to become hydrated, the green light on the transmitter will flash within 10 seconds.
- 5 If the transmitter light does not flash, disconnect the transmitter from the sensor, wait for several seconds and then reconnect. If the transmitter light still does not flash, charge the transmitter.
- 6 When the transmitter light flashes green when connected to the sensor, use your pump to start the sensor. For more instructions, see your pump user guide.
- 7 After the transmitter successfully sends sensor data to the pump, attach the sensor's adhesive tab to the transmitter.
- 8 Follow the instructions that appear on the pump screen or follow the instructions in your pump user guide.

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
- 1 Carefully remove any occlusive dressing from the transmitter and sensor.
- 2 For the Enlite sensor, remove the adhesive tab from the top of the transmitter.
- 3 Hold the transmitter as shown, and pinch the flexible side arms of the sensor between your thumb and forefinger.
- 4 Gently pull the transmitter away from the sensor.
- 5 Follow the instructions that appear on the pump or follow the instructions in your pump user guide.

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Watertight Tester

The tester is used to test the transmitter to make sure it is working. It is also used as a required component for cleaning the transmitter. Properly connecting the tester to the transmitter will ensure that fluids do not come in contact with the transmitter's connector pins. Fluids can cause connector pins to corrode and affect the transmitter's performance.

Do not twist the tester while attached to the transmitter. This will damage the transmitter.


The tester can be used for one year. If you continue to use the tester for more than one year, the transmitter's connector pins could be damaged, because the tester cannot continue to provide a waterproof seal. For instructions on how to check the connector pins, see *Inspecting the transmitter connector pins*, on page 8.

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
Also look for moisture inside the connector opening. If you see any moisture, allow the transmitter to dry for at least one hour. Moisture inside the connector opening could cause the transmitter to not work properly, and could cause corrosion and damage over time.

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To connect the tester:

- 1 Hold the transmitter and the tester as shown. Line up the flat side of the tester with the flat side of the transmitter
- 2 Push the tester into the transmitter until the flexible side arms of the tester click into the notches on both sides of the transmitter
- 3 Within five seconds, when properly connected, the green light on the transmitter flashes for about 10 seconds.
- 4 To test the transmitter, check the sensor icon on the pump to ensure that the transmitter is sending a signal (see your pump user guide).
- 5 To clean the transmitter, see *Cleaning the transmitter*, on page 9.
- 6 After testing or cleaning, disconnect the tester from the transmitter

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
Cleaning the transmitter

The transmitter is a single-patient use device and not intended for multi-patient use.

Caution: Do not discard the transmitter in a medical waste container or otherwise subject it to incineration. The transmitter contains a battery that may explode upon incineration.

Note: The tester is a required component for cleaning the transmitter. For details, see *Watertight Tester*, on page 7.

Always clean the transmitter after each use.

To clean the transmitter, you will need the following materials: mild liquid soap, a soft-bristled toddler toothbrush, a container, 70% isopropyl alcohol, and a few clean, dry cloths.

Warning: Cracking, flaking, or damage of the housing are signs of deterioration and the performance of the device may be compromised. This may affect the ability to properly clean and disinfect the transmitter. If these signs are noted, stop using the device and call the 24 Hour HelpLine or your local representative. The device must be discarded according to local regulations for battery disposal (non-incineration).

To clean the transmitter:

- 1 Attach the tester to the transmitter.
- 2 If optional occlusive dressing or overtape was used and there is adhesive residue on the transmitter, see *Removing adhesive residue*, on page 10.
- 3 Rinse the transmitter under room temperature tap water for at least one minute, and until visibly clean. Make sure all hard-to-reach areas are rinsed completely.
- 4 Prepare a mild liquid soap solution using 5 milliliters (1 teaspoon) of liquid soap per 3.8 liters (1 gallon) of room temperature tap water.
- 5 With tester still attached, submerge the transmitter in the cleaning solution and soak for one minute.
- 6 Holding the tester, brush the entire surface of the transmitter using a soft-bristled toddler toothbrush. Make sure to brush all hard-to-reach areas until visibly clean.
- 7 Rinse the transmitter under running room temperature tap water for at least one minute, and until all visible liquid soap is gone.
- 8 Dry the transmitter and tester with a clean, dry cloth.
- 9 Place the transmitter and tester on a clean, dry cloth and air dry them completely.
- 10 Disconnect the tester from the transmitter.

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To remove adhesive residue, you will need the following materials: medical adhesive remover (such as Detachol® - which is a mineral spirit) and cotton swabs.

Note: During testing, Medtronic MiniMed used Detachol to remove the adhesive residue from the transmitter. Detachol is recommended for use but may not be available in all countries.

To remove adhesive residue:

- 1 Make sure the tester is attached to the transmitter.

- 2 Holding the tester, saturate a cotton swab in the Detachol solution and gently rub the adhesive residue on the transmitter until it is fully removed.
- 3 Continue with the cleaning procedure. See *Cleaning the transmitter*, on page 9 for details.

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To clean the charger:

- 1 Wash your hands thoroughly.
- 2 Use a damp cloth with mild cleaning solution, such as a dishwashing detergent, to clean any dirt or foreign material from the outside of the charger. Never use organic solvents, such as paint thinner or acetone, to clean the charger.
- 3 Place the charger on a clean, dry cloth and air dry for two to three minutes.

Troubleshooting

The following table contains troubleshooting information for the transmitter, charger, and tester. For more information about troubleshooting, see your pump user guide.

Problem	Likely Cause(s)	Resolution
You connected the transmitter to the charger and no lights came on.	The transmitter connector pins are damaged or corroded. Your charger battery has no power.	1 Check the transmitter connector pins for damage or corrosion. For more information about your connector pins, see <i>Inspecting the transmitter connector pins</i>, on page 8. If the pins are damaged or corroded, contact the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter. 2 If there is no damage to the connector pins, replace the battery in the charger. For instructions on replacing your charger battery, see <i>Installing a battery in the charger</i>, on page 4.

Problem	Likely Cause(s)	Resolution
During charging, the flashing green light on the charger turns off and you see a flashing red light on the charger.	Your charger battery is low on power.	Replace the battery in the charger. For instructions on replacing your charger battery, see <i>Installing a battery in the charger</i> , on page 4.
During charging, the flashing green light on the charger turns off and you see a series of quick flashing red lights on the charger	Your transmitter is low on power.	1 Charge the transmitter continuously for one hour. If flashing does not stop, proceed to step 2. 2 Charge the transmitter continuously for eight hours. If flashing does not stop, call the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.
During charging, a mix of quick and long flashing red lights appear on the charger.	Your charger and your transmitter are low on power.	1 Replace the battery in the charger. For instructions on replacing your charger battery, see <i>Installing a battery in the charger</i>, on page 4. 2 Charge the transmitter continuously for one hour. If flashing does not stop, proceed to step 3. 3 Charge the transmitter continuously for eight hours. If flashing does not stop, call the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.
The green light on the transmitter does not flash when you connect it to the sensor	Your transmitter is not fully connected. Your transmitter is low on power. Your sensor is not properly inserted into your body.	1 Disconnect the transmitter from the sensor 2 Wait for five seconds and reconnect them. If the green light still does not flash, proceed to step 3. 3 Fully charge the transmitter. If the green light still does not flash, proceed to step 4. 4 The sensor may not be properly inserted into your body. Insert a new sensor.

Problem	Likely Cause(s)	Resolution
The green light on the transmitter does not flash when you connect it to the tester	Your transmitter is low on power. Your transmitter is not fully connected.	1 Check the connection between the transmitter and the tester. If the green light still does not flash, proceed to step 2. 2 Fully charge the transmitter. 3 Test the transmitter with the tester again. If you still do not see the green light flash, call the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.
Your transmitter battery does not last for six days.	Your transmitter is not fully charged when you connect it to the sensor. The transmitter and pump frequently lose wireless connection.	1 Fully charge the transmitter before connecting it to the sensor. If the transmitter battery still does not last for six days, proceed to step 2. 2 Move away from any device that can cause RF interference. For more information on RF interference, see <i>Radio Frequency (RF) communication</i>, on page 2. 3 Make sure your pump and your transmitter are located on the same side of your body to minimize any RF interference. If your fully charged transmitter battery continues to lose power before a full six days, call the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.
Your transmitter has lost connection with your pump.	Your pump is out of range. There is RF interference from other devices.	1 Move away from any device that can cause RF interference. For more information on RF interference, see <i>Radio Frequency (RF) communication</i>, on page 2. If your transmitter is still not communicating with your pump, proceed to step 2. 2 Make sure your pump and your transmitter are located on the same side of your body to minimize any RF interference. If your transmitter is still not communicating with your pump, call the 24 Hour HelpLine or your local representative for assistance.

Storing the devices

Store the transmitter, charger, and tester in a clean, dry location at room temperature. If the transmitter is not in use, you must charge the transmitter at least once every 60 days. Although not required, you may store the transmitter on the charger. If you are storing the transmitter on the charger, you must disconnect and reconnect the charger and the transmitter at least once every 60 days.

Disposal

Discard the transmitter according to local regulations for battery disposal.

Specifications

Biocompatibility	Transmitter: Complies with EN ISO 10993-1
Applied parts	Transmitter Sensor
Operating conditions	Transmitter temperature: -5 °C to 45 °C (23 °F to 113 °F) Caution: When operating the transmitter on a tester in air temperatures greater than 41 °C (106 °F), the temperature of the transmitter may exceed 43 °C (109 °F) Transmitter relative humidity: 5% to 95% with no condensation Transmitter pressure: 61.36 kPa to 106.17 kPa (8.9 psi to 15.4 psi) Charger temperature: 10 °C to 40 °C (50 °F to 104 °F) Charger relative humidity: 30% to 75% with no condensation
Storage conditions	Transmitter temperature: -25 °C to 55 °C (-13 °F to 131 °F) Transmitter relative humidity: 10% to 100% with condensation Transmitter pressure: 61.36 kPa to 106.17 kPa (8.9 psi to 15.4 psi) Charger temperature: -10 °C to 50 °C (14 °F to 122 °F) Charger relative humidity: 10% to 95% with no condensation
Battery life	Uses one new AAA battery to charge the transmitter.
Transmitter frequency	2.4 GHz, 2M65G1D modulation, less than 1mW ERP
Maximum output power (EIRP)	-0.63 dBm

Radio Frequency (RF) communications	Pump to transmitter frequency: 2.4 Ghz; proprietary Medtronic protocol; range up to 1.8 meters (6 feet) Utilizes the IEEE 802.15.4 protocol with proprietary data format Operating frequency: 5 frequencies are used: 2420, 2435, 2450, 2465, and 2480 MHz Bandwidth: 5 MHz, which is allocated channel bandwidth per IEEE protocol
Transmitter expected service life	The transmitter expected service life is 1 year depending on patient usage

Guardian 2 Link wireless communication

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The Guardian 2 Link transmitter will initiate a channel search when beacon detection fails on the associated channel. The search will be conducted across all five channels. When the beacon is located the transmitter will rejoin on the identified channel. Upon re-association any missed packets (up to 10 hours) will be transmitted from the transmitter to the pump.

In normal operation the transmitter will transmit a packet every 5 minutes and retransmit the packet if the data is corrupted or missed.

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Guidance and manufacturer's declaration

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Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The Guardian 2 Link is intended for use in the electromagnetic environment specified below. The customer or the user of the Guardian 2 Link should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±2 kV, ±4 kV, ±8 kV Air ±2 kV, ±4 kV, ±6 kV Indirect	±8 kV Air ±6 kV Indirect ±22 kV Air, <5% RH	For use in a typical domestic, commercial, or hospital environment.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable	Requirement does not apply to this battery powered device.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Not applicable	Requirement does not apply to this battery powered device.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% U _T (>95% dip in U _T) for 0.5 cycle	Not applicable	Requirement does not apply to this battery powered device.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The Guardian 2 Link is intended for use in the electromagnetic environment specified below. The customer or the user of the Guardian 2 Link should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	400 A/m 4000 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical domestic, commercial, or hospital environment.
Note: U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The Guardian 2 Link is intended for use in the electromagnetic environment specified below. The customer or user of the Guardian 2 Link should assure that it is used in such an electromagnetic environment.			
Immunity Test	IEC 60601 Level	Compliance Level	Electromagnetic Environment Guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Guardian 2 Link, including cables, than the recommended separation distance calculated from the equation applicable to the power of the transmitter Refer to the recommended separation distance table for more information.
Conducted RF IEC 61000-4-6	3 V/m 150 kHz to 80 MHz	Not applicable	Not applicable

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The Guardian 2 Link is intended for use in the electromagnetic environment specified below. The customer or user of the Guardian 2 Link should assure that it is used in such an electromagnetic environment.			
Immunity Test	IEC 60601 Level	Compliance Level	Electromagnetic Environment Guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 6 GHz	10 V/m 80 MHz to 6 GHz	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ 800 MHz to 6 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Note: At 80 MHz and 800 MHz, the higher frequency range applies.			
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption, and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcasts and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Guardian 2 Link is used exceeds the applicable RF compliance level above, the Guardian 2 Link should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Guardian 2 Link. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between the Guardian 2 Link and common household radio transmitters			
Household RF Transmitter	Frequency	Recommended Separation Distance (meter)	Recommended Separation Distance (inch)
Telephones			
Cordless Household	2.4 GHz	0.3	12
Cordless Household	5.8 GHz	0.3	12
TDMA-50 Hz (cell phone)	1.9 GHz	0.3	12
TDMA-50 Hz (cell phone)	800 MHz	0.3	12
PCS (cell phone)	1.9 MHz	0.3	12
DCS (cell phone)	1.8 MHz	0.3	12
GSM (cell phone)	900 MHz	0.3	12
GSM (cell phone)	850 MHz	0.3	12
CDMA (cell phone)	800 MHz	0.3	12
Analog (cell phone)	824 MHz	0.3	12
CDMA (cell phone)	1.9 MHz	0.3	12
WiFi Networks			
802.11b; 11Mbps maximum	2.4 GHz	1	39.5
802.11g; 54 Mbps maximum	2.4 GHz	1	39.5
802.11n; 11Mbps maximum	2.4 GHz	1	39.5
Bluetooth 500 kb/s	2.4 GHz	0.1	3.93
ZigBee 250 kb/s	2.4 GHz	0.1	3.93

Recommended separation distances between portable and mobile RF communications equipment and the Guardian 2 Link			
The Guardian 2 Link is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the Guardian 2 Link users can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the Guardian 2 Link as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to the frequency of transmitter (m)		
	150 kHz to 80 MHz Not applicable	80MHz to 800MHz $d = 0.35 \sqrt{P}$	800MHz to 6.0GHz $d = 0.70 \sqrt{P}$
0.01	Not applicable	0.035	0.07
0.1	Not applicable	0.11	0.11
1	Not applicable	0.35	0.7
10	Not applicable	1.1	2.2
100	Not applicable	3.5	7
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer			
Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			

Icon Table

	Serial number
	Catalogue or model number
(1X)	One per container/package
	Date of manufacture
	Manufacturer
	Follow instructions for use (appears blue on label)
	Consult instructions for use
	Storage temperature range
	Radio communication
	Configuration or unique version identifier
	Degree of protection against electric shock: Type BF applied part
IP48	Transmitter: Protected against the effects of continuous immersion in water (2.4 meters (8 feet) immersion for 30 minutes).
	Storage humidity range
	RF device distributed in Australia
CE 0459	Marking of Conformity: This symbol means the device fully complies with MDD 93/42/EEC (NB 0459).
IC	RF device distributed in Canada
	Authorized representative in the European community

©2015 Medtronic MiniMed, Inc. All rights reserved.

MiniMed™, Enlite™, and Guardian™ are trademarks of Medtronic MiniMed, Inc.

Detachol® is a registered trademark of Ferndale Laboratories Inc.

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