Tyvek® Pouch with STERRAD® Chemical Indicator (Heat Seal (125XX) and Self Seal (123XX)) Tyvek® Roll with STERRAD® Chemical Indicator (124XX)

EN

INSTRUCTIONS FOR USE

INDICATIONS FOR USE:
Tyvek® Pouches and Rolls with STERRAD® Chemical Indicator are intended to be used to enclose medical devices that are to be terminally sterilized in the STERRAD® Sterilization Systems and Cycles, in a single or double pouch configuration, and to indicate, by color change, that the pouch has been exposed to sterilant. After completion of the sterilization process, the pouch/roll is intended to maintain sterility of the enclosed device until used.

The pouches and rolls are printed with a chemical indicator bar which is a process indicator (ISO 11140-1:2014 [Type I]) that changes from red to the color indicated on the comparator bar included on the shelf pack (or lighter) when exposed to hydrogen peroxide vapor during processing in the STERRAD® Sterilization Systems.

The pouches and rolls can be used in the STERRAD® Sterilization Systems and Cycles shown in Table 1.

Table 1: Indications for Use

STERRAD® / Cycle		Single-channel stainless steel lumens	Other lumens	Load requirements
STERRAD® 100S		≥1 mm ID x ≤125 mm long ≥2 mm ID x ≤250 mm long ≥3 mm ID x ≤400 mm long	Metal and nonmetal ≥6 mm ID x ≤310 mm long	≤10 lumens per load
STERRAD NX® with or without ALLClear™ Technology	STANDARD	≥1 mm ID x ≤150 mm long ≥2 mm ID x ≤400 mm long	N/A	≤10 lumens per load ≤10.7 lbs
	ADVANCED	≥1 mm ID x ≤500 mm long	N/A	≤10 lumens per load ≤10.7 lbs
STERRAD® 100NX with or without ALLClear™ Technology	STANDARD	≥0.7 mm ID x ≤500 mm long	N/A	≤10 lumens per load ≤21.4 lbs**
	EXPRESS	N/A	N/A	≤10.7 lbs on bottom shelf only
	DUO	N/A	N/A	Accessory devices (e.g., cameras and light cords) only 2 cameras maximum

**The validation testing for this lumen size was conducted using a maximum of 10 lumens in the USA; 40 lumens for markets outside the USA. Your loads should not exceed the validated maximum number of lumens.

DESCRIPTION:
A. The pouches and rolls are constructed from two films: a transparent, non-permeable film, and an opaque, permeable film of Tyvek® material.
B. Real-time testing supports maintenance of post-processing pouch/roll sterility for a period of 12 months.
C. The pouches and rolls are printed with a chemical indicator bar which is a process indicator (ISO 11140-1:2014 [Type I]) that changes from red to the color indicated on the comparator bar included on the shelf pack (or lighter) when exposed to hydrogen peroxide vapor during processing in the STERRAD® Sterilization Systems. Real-time testing has demonstrated stability of the post-processing chemical indicator color change for up to 6 months.
D. The Product Codes for all Tyvek® Pouches and Rolls with STERRAD® Chemical Indicator are listed in Table 2.

ADVANCED STERILIZATION PRODUCTS
Division of Ethicon, Inc.
a Johnson & Johnson company
33 Technology Drive, Irvine, CA 92618
EC REP
Johnson & Johnson MEDICAL GmbH
Robert-Koch-Str. 1, 22851 Norderstedt, Germany

US Irvine, CA 92618
FR 92787 Issy-les-Moulineaux
DE 22851 Norderstedt
CH 6300 Zug
NL 3800 AD Amersfoort
GR 15125 Maroussi, Athens
SE 19184 Sollentuna
BE 1831 Diegem
CZ Johnson & Johnson, s.r.o. Karla Engliše 3201/6 150 00 Praha 5
SK Karadžičova 12 821 08 Bratislava
IT Johnson & Johnson Medical S.p.A. Via del Mare, 56 - 00040 Pomezia, Roma
ES Johnson & Johnson Medical Paseo de las Doce Estrellas, 5-7 28042 Madrid
GB Wokingham, RG40 3EW
CA Johnson & Johnson Medical Products Markham, ON, L3R 0T5
JP ジョンソン・エンド・ジョンソン株式会社 〒101-0065 東京都千代田区西神田３丁目５番２号
AT 1190 Wien
BR Rodovia Presidente Dutra, Km 154 S.J. Campos-S.P 12240-908
PT 2745-555 Barcarena
PL Iłżecka 24, 02-135 Warszawa
HU Johnson & Johnson Kft. H-1125 Budapest, Nagyenyed u. 8-14

© Ethicon, Inc. 2004-2017
110085-04

Table 2: Product Codes for Tyvek® Pouches and Rolls with STERRAD® Chemical Indicator

Product Type	Product Codes
Tyvek® Self-Seal Pouches with STERRAD® Chemical Indicator (123XX)	12320, 12326, 12335, 12332, 12342, 12340, 12348, 12356
Tyvek® Rolls with STERRAD® Chemical Indicator (124XX)	12407, 12410, 12415, 12420, 12425, 12435, 12442, 12450
Tyvek® Heat-Seal Pouches with STERRAD® Chemical Indicator (125XX)	12521, 12526, 12532, 12541, 12548, 12543*, 12544*, 12557*, 12558*, 12559*, 12560*

**Indicates Product Codes exclusively for markets outside the USA.*

DIRECTIONS FOR USE:
A. Place medical devices(s) through an open end of pouch/roll. If using a self-seal pouch, peel off self-seal tape and seal the open end of the pouch. If using a heat-seal pouch, place in heat sealer following directions listed in the "NOTE" directly following. For roll, cut to desired length and seal open ends following heat sealer instructions.
NOTE: When heat sealing the Tyvek® pouch/roll for the first time, the operator should establish the seal temperature/dwell time range for their specific machine according to the peel strength required when opening the package. Typical conditions for heat sealer vary by machine, but are typically between 121°C-135°C [250°-275°F], and generally require 1 to 3 seconds of seal time. Seal temperature/time used should be recorded for future reference.
B. Do not reuse pouches/rolls if a sterilization cycle is aborted by the operator, or cancels.
C. Refer to the comparator bar on the shelf pack for the proper post-processing color change. The chemical indicator changes from red to the color indicated on the comparator bar (or lighter) when exposed to H_2O_2 vapor during processing in the STERRAD® Sterilization Systems.
D. Correct color change of the chemical indicator bar indicates exposure to hydrogen peroxide, but does not indicate sterility.

STORAGE:
A. Store at room temperature, 15°C-30°C (59°F-86°F), under dry conditions. Avoid exposing the chemical indicator side of pouch/roll to light.
B. Before and after processing, store away from STERRAD® Sterilizer or sources of hydrogen peroxide, acids, alkalis, or volatile antimicrobial agents such as glutaraldehyde or formaldehyde.

CAUTION:
A. After processing, the chemical indicator may not be suitable for use as permanent record due to its light and moisture sensitivity. If one should choose to store it, place pouch/roll in a plastic bag protected from light and moisture.
B. If the chemical indicator bar has not completely changed from red to color indicated by comparator bar, introduction of H_2O_2, an essential part of the STERRAD® Sterilization cycle, may not have taken place.
C. STERRAD® Chemical Indicators offer an additional way to verify processing in the sterilization cycle and should be used in addition to, not in place of, the biological indicator. STERRAD® Chemical Indicators do not signify sterilization; they solely confirm that the indicator has been exposed to the hydrogen peroxide.
D. Do not use after expiration date.
E. Single use only. Do not reuse as the pouches and rolls will not function properly if reused. Should a cycle cancellation occur, repackage the devices using new pouches and rolls when restarting the cycle.
F. Do not use if pouch or roll is damaged.
G. Users should perform their own determination for their heat sealer parameters.
H. Tyvek® pouches may not be appropriate for some delicate devices requiring protective packaging such as flexible endoscopes. Always refer to the device manufacturer's instructions for use for special packing instructions.

Tyvek® is a DuPont registered trademark.

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©Ethicon, Inc. 2004-2017 110085-04